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
Operations/Service Manual

 **smith&nephew**
LENS 4K
Camera Control Unit

Glossary of symbols

	Follow instructions for use
	CE Mark
	Catalog number
	On/Off
	Up and Down Buttons
	White balance
	Menu
	Select button
	Brightness
	Do not plug in wet
	This way up
	Consult instructions for use
	Electrical Equipment Safety System (EESS)

	Authorized representative in the European community
	Manufacturer
	TUV certification
	Video output
	Video input
	USB output
	Humidity limitation
	Power plug
	Equipotentiality
	Atmospheric pressure
	Wi-Fi
	Color Camera
	Keep dry
	MR Unsafe
	Serial number
	Type CF equipment

	Fuse Rating
	EU: Not for general waste
	Non-ionizing electromagnetic radiation
	Temperature limitation
	Ethernet
	Date of Manufacturer
	Caution
	General Warning Sign
	Medical Device
	Unique Device Identifier
	Prescription use only

Preface

Preface

This manual provides the information needed to operate and maintain the Smith & Nephew LENS 4K Camera Control Unit. Prior to using this equipment, read this manual thoroughly, paying particular attention to the operating instructions, warnings, and precautions. Smith & Nephew surgical equipment is designed for use only by medical professionals who are completely familiar with the appropriate surgical techniques and video procedures.

Table of Contents

Preface	4	Warnings	20	CHANGE PROCEDURE NAME	32
Device Description	5	Turn the Camera Control Unit On	21	BRIGHTNESS	32
Intended Users	5	Connect the iPad	21	ENHANCEMENT	32
Intended Use	5	Operation	21	CHROMA (SATI)	33
Indications for Use	5	Select a Camera setup	21	PHASE (HUE)	33
Contraindications	5	MAIN MENU	22	ZOOM	33
Target Patient Group	5	Brightness	22	CONFIGURE ICONS	33
Warnings	5	Zoom	22	MORE SETTINGS	33
If using the Wi-Fi enabled control unit	6	Enhancement	22	DISABLE BEEP ON CAPTURE	34
Precautions	7	White Balance	22	Disable Other Beeps	34
Adverse Event Reporting	7	Camera Setups	22	DISABLE PATIENT INFO	34
System Components	7	Color Bar	23	DEFAULT CHROMA MODE	34
Accessories	7	System Configuration	23	660 DUAL CAPTURE	34
System Controls	8	Buttons	23	660 DUAL VIDEO	34
Front Panel Overview	8	Exit	23	ELC Settings	34
Rear Panel Overview	9	Capture Images	23	CHROMA MODE VALUE	35
System Controls	10	To Capture a Still Image	23	Save and Exit	35
Camera Head Controls	10	To Capture Video Images	24	Set the Default Startup Camera setup	35
Button Function	10	Security Settings	24	Cleaning	35
General Screen Navigation	12	Resetting the system if the user password is lost	24	To Clean the LENS 4K Camera Control Unit	35
Installation	12	Alternative Power Up	24	To Clean and Sterilize the Camera Head	35
Unpacking and General Inspection	12	Menu Maps	25	Troubleshooting	36
LENS 4K Camera Control Unit	12	Postoperative	26	Expected Life	39
Setup	12	Customize the System	26	Service	39
Secure the Control Unit	13	Set/Change Password Protection	26	Maintenance	39
IEC 60601-1 Compliant System Configuration	13	Set/Change Admin Password	27	Recommended Annual Performance Checks	39
Warnings	13	Configure Peripheral Ports	28	To Replace Fuses	39
Setup Configurations Summary	13	Language	28	Warnings	39
Cable Connectors	13	Reset Factory Defaults	28	System Specifications	40
Recommended System Configuration for the Wi-Fi-Enabled System	14	Network	29	LENS 4K Camera Control Unit Specifications	40
Set-up	14	To change the network ID:	29	Environmental Conditions	40
Configure the peripheral ports	14	To change the network password:	29	Ordering Information	41
Recommended System Configuration, Non Wi-Fi	16	MAC address:	29	LENS 4K Camera Control Unit - Wi-Fi	41
Set-up	16	To change the network channel:	29	Compatible Camera Head	41
Configure the peripheral ports	16	IP Address:	29	Compatible Devices	41
Inspect the System Components	18	To activate FCC power:	30	LENS 4K Camera Control Unit - Non Wi-Fi	41
Electrical Connections	18	ADVANCED SETTINGS	30	Guidance and Manufacturer's Declaration	42
Inspect the Fiber Optic Light Cable	18	To configure HDSOI Outputs:	30	Electromagnetic Emissions	42
Service Life Checklist	18	To configure COMPRESSION:	30	Guidance and Manufacturer's Declaration - Guidance for Separation Distances	43
Apply Power to the Camera Control Unit	19	To configure AUTO START:	30	Electromagnetic Immunity	44
Perform the Network Setup	19	To configure the HEAD PLUGIN	30	Special Environment	45
Preoperative	20	To configure VIDEO INPUT:	31	For Further Information	46
Prepare the Camera Head	20	Versions	31	References	46
Connect the Camera Head	20	Customize Camera setup Settings	31		
Connect the Coupler/Videoendoscope to the Camera Head	20	Customize Button Settings	31		
Insert Light Guide	20	IMAGE CAPTURE OSD	32		
		DISPLAY IN STARTUP MENU?	32		
		PORTRAITS ENABLED?	32		

Device Description

The Smith & Nephew LENS 4K Camera Control Unit (CCU), Wi-Fi and non-Wi-Fi systems is a state-of-the-art video system designed for use in arthroscopic and endoscopic surgical procedures. The System enables control of camera features (enhancement features and zoom capabilities) and accessories (video printers, and image management systems) through either the Camera Head, Control Unit buttons, or the optional LENS Connected Tower Application (Wi-Fi enabled CCU version only). The System includes an integrated LED Light Source, which provides illumination during the procedure via designated light guides and endoscopes.

Note: This document will use the term Tablet Application when referencing the LENS Connected Tower Application.

The Wi-Fi enabled version of the System provides image management features in conjunction with the Tablet Application. These features include image capture, video capture, patient management, media object management, sharing, and archiving.

Refer to the Tablet Application Instructions for Use (REF 10601295) for information regarding the Tablet Application.

Refer to the LENS 4K Camera Head Instructions for Use (REF 10601349) or (REF 10601459 US Only) for information regarding the Camera Head.

Intended Users

The Smith & Nephew LENS 4K Camera Control Unit is intended for use by healthcare professionals in accordance with these instructions for use. The use environment is professional healthcare facility in accordance with guidance and manufacturer declaration EMC section.

Intended Use

The LENS 4K Camera Control Unit and Camera Head are intended to provide illumination, visualization and capture of still and motion pictures of surgical sites.

Indications for Use

The LENS 4K Camera Control Unit & Camera Head is used in diagnostic and operative endoscopic procedures such as arthroscopy to provide illumination, visualization and capture of still and motion pictures of surgical sites within articular cavities, hollow organs, body cavities and canals. Additionally, LENS 4K Camera Control Unit & Camera Head are indicated for use in endoscopic surgical procedures in the thoracic cavity when used with an appropriately indicated thoracoscope.

Contraindications

None known.

Target Patient Group

The targeted patient population for the LENS 4K Camera Control Unit is patients requiring endoscopic surgical procedures as determined necessary by the prescribing healthcare professionals.

Warnings

- It is the health care professional's responsibility to be familiar with the appropriate surgical techniques prior to use of this device.
- Read these instructions completely prior to use.
- **DANGER:** Risk of explosion if used in the presence of flammable anesthetics.
- Use aseptic technique in accordance with standard operating room procedures.
- To prevent electric shock, do not remove any covers or screws from the control unit. There are no user-serviceable components inside. Dismantling the equipment will void the warranty.
- To prevent electric shock, unplug the unit from the electrical outlet before attempting to replace the fuses.
- To avoid fire hazard, use only fuses of the correct type, voltage rating, and current rating.
- If this unit is configured as part of a system, the entire system should be tested for compliance with IEC 60601-1.
- If the leakage current of the configured system exceeds the limits of IEC 60601-1, install an appropriately rated IEC 60601-1 approved isolation transformer with a minimum of 1000 VA and retest the system.
- To avoid System power loss, use an uninterruptible power supply (UPS).
- The use of accessory equipment not complying with the equivalent safety requirements of this equipment may lead to a reduced level of safety. Consider the following when choosing accessory devices:
 - Use of the accessory in the patient vicinity.
 - Evidence that the safety certification of the accessory has been performed in accordance with IEC 60601-1.
- Use of accessories and cables with this equipment, other than those specified or provided by the Smith & Nephew, could result in increased electromagnetic emissions or decreased electromagnetic immunity and result in improper operation.

Warnings

- To prevent damage to internal optics only use compatible light cables.
- To avoid damaging the optics and light guide detector, do not use wet light cables.
- To avoid damage to the turret use only compatible light cables.
- When the light source is turned on, do not look directly at the LED light guide without protective goggles. The high- intensity of the light guide may cause burns to the skin, or permanent damage to eyes.
- Hazardous high voltage and energy are present at the output and in the internal circuitry of this unit.
- When inspecting the light cable, NEVER aim the light cable at, or peer directly into the light source. Retinal damage may occur.
- Connect only accessory items that have been specified as part of the System or have been specified as compatible. Use of the accessory, transducer, or cable with equipment and systems other than those specified may result in increased emissions or decreased immunity of the equipment or system.
- During operation, avoid prolonged contact of the scope tip to patient tissue or flammable materials. The scope tip may reach high temperatures due to high intensity light transmission.
- Do not leave the operating light cable on a patient or the drapes. Failure to observe this warning may result in burns to the patient and/or the surrounding drapes.
- Do not use this device in oxygen rich environments (>25%).
- To avoid the risk of electric shock, this equipment must only be connected to supply mains with protective earth.
- No modification of this equipment is allowed
- Electrical Interference. This equipment is designed to be used as a stand-alone, or in proximity to other medical equipment. Although this equipment is designed and tested to minimize interference with, or by, other electrical equipment; interference may occur. If interference occurs with other equipment it may be corrected by one or more of the following measures:
 - Reorient or relocate this equipment, the other equipment, or both.
 - Increase the separation between the pieces of equipment
 - Connect the pieces of equipment into different outlets or circuits.
 - Consult a biomedical engineer.
- Abnormal performance, due to Electromagnetic Disturbances, may be manifested by a failure to properly maintain illumination and visualization of the surgical field. If abnormal performance is observed, it may be corrected by one or more of the following measures:
 - Reorient or relocate this equipment, the other equipment, or both.
 - Increase the separation between the pieces of equipment.
 - Connect the pieces of equipment into different outlets or circuits.
 - Consult a biomedical engineer.
- High Frequency (HF) Interference: This equipment was designed and tested to be used with HF Electrical Surgical Unit (ESU) equipment. If interference occurs with other equipment it may be corrected by one or more of the following measures:
 - Reorient or relocate this equipment, the other equipment, or both.
 - Increase the separation between the pieces of equipment
 - Connect the pieces of equipment into different outlets or circuits.
 - Consult a biomedical engineer.

If using the Wi-Fi enabled control unit

- Portable and mobile RF communication equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the LENS 4K Camera Control Unit, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.
- Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.
- Use of accessories and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

Precautions

R U.S. Federal law restricts this device to sale by or on the order only of a physician.

- Prior to each use, inspect the device to ensure it is functioning properly and is not damaged. Do not use a damaged device.
- Do not spill liquids, or allow fluids to enter the device. If liquid enters the unit through the vents, permanent damage may result.
- Electrical safety testing should be performed by a biomedical engineer or other qualified person.
- All products used in conjunction with this control unit, such as video monitors, video printers, etc., are likely to have an inherent leakage current at their video connections. Connection to this control unit likely will increase the leakage current at the control unit's chassis.

When looping the video output to more than one piece of equipment, the termination switch on the last looping piece of equipment should be in the 75 Ω (ohm) or the ON position if the equipment provides that option.

- Do not plug the camera cable connector into the camera control unit if wet. Moisture on the camera cable connector, including the gold fingers, on the card edge of the camera head cable will damage the circuitry and void the warranty. Ensure that the camera cable connector is completely dry prior to plugging the camera head into the camera control unit.
- This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in an industrial installation. This equipment generates and can receive radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one of the following measures:
 - Reorient or relocate the receiving antenna.
 - Increase the separation between the equipment and receiver.
 - Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
 - Consult the dealer or an experienced radio/TV technician for help.
- This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.
- Maintain a minimum separation distance of 20 cm between the antenna and the patient's body.
- Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate this equipment.

- Environmental Protection: This equipment contains electronic printed circuit assemblies. At the end of the useful life of this equipment, it should be disposed of in accordance with any applicable national or institutional policy relating to obsolete electronic equipment.
- Products may wear over time as a result of normal use. Periodic retesting of the leakage current is recommended.
- Do not operate at line voltages other than those stipulated on the back of the unit.
- Handle the unit with care. If the unit is dropped or damaged in any way, it must be returned immediately for service.
- To avoid loss of or damage to the control unit, remove the antenna from the Wi-Fi enabled control unit prior to transport.
- Disconnect the power cord before cleaning the unit or light source.
- Ensure that the unit is installed in an environment in which the vents and fans are not obstructed. Insufficient ventilation can cause the unit to overheat and shut down.
- In the event that the unit overheats and shuts down but the cooling fans are still running, do not turn off or unplug the unit. Continue to run the fans to prevent further damage to the unit.
- Ensure that the light guide is inserted into the proper front panel port. Damage to the unit may occur if the light guide is not inserted into the port associated with the correct light guide manufacturer.
- The mains inlet is used as the means of isolation from mains and shouldn't be bypassed. Ensure the equipment allows access to the plug.
- This product contains Nickel (Ni) which may cause allergic reactions.

Adverse Event Reporting

If, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and to your national authority.

complaints@smith-nephew.com

System Components

The complete System includes the control unit, LENS 4K Camera Head (sold separately), and the optional Tablet Application. It is configured to be connected as a stand-alone system on an endoscopic cart or tower. The System can be used in conjunction with compatible HD/UHD flat panel displays, image capture devices, and print and video recording devices.

Accessories

Refer to the Ordering Information section in this manual for more information.

System Controls

System Controls

Front Panel Overview


Figure 1. LENS 4K Camera Control Unit Front Panel

	Front Panel Control/Connector	Function
1	Standby System Power button	Turns the control unit on and off. The rear power switch must be switched to the (I) position in order to provide mains power to the system.
2	Standby Illumination button	When the Standby Illumination button is not illuminated, the light source is in standby mode and the light source is off. When illuminated, the button turns standby off and turns the light source on.
3	Active Light Guide port indicator	The light guide port to be used must be aligned with this orange arrow.
4	Light Cable Selection Turret	Rotate and click into position to select the port for the Smith & Nephew/Wolf™, Olympus™, ACMI™, or Storz™ Light Guide.
5	Olympus Light Guide port	This self-closing light port is designed to accept Olympus Light Guides.
6	ACMI Light Guide port	This self-closing light port is designed to accept ACMI Light Guides.
7	Storz Light Guide port	This self-closing light port is designed to accept Storz Light Guides.
8	Smith & Nephew/Wolf Light Guide port	This self-closing light port is designed to accept the Smith & Nephew/Wolf Light Guide end-fittings.
9	Camera Cable Connector receptacle	Smith & Nephew 4K camera head is connected via the camera cable connector. Be sure the camera cable connector side with the triangle faces upward. Firmly push the camera cable connector into the keyed receptacle.
10	White Balance button	The White Balance button automatically adjusts the camera for optimal white balance when the button is depressed. The camera must be focused on a white object during the white balance process.
11	Up Arrow button	The Up Arrow button causes the cursor to scroll upward or left through the available menu selections after the OSD (On- Screen Display) menu has been accessed via the SELECT button. When the end of a menu is reached, the next Up Arrow button press will bring the user to the bottom of the menu.
12	SELECT/Menu Button	The SELECT button allows the user to either select a menu command or call up the MAIN MENU screen. When the MAIN MENU is present, the SELECT button selects the highlighted menu command.
13	Down Arrow button	The Down Arrow button causes the cursor to scroll down or right through the available menu selections after the OSD menu has been accessed via the SELECT button. When the end of a menu is reached, the next Down Arrow button press will bring the user to the top of the menu.
14	USB receptacle	The LENS 4K supports compatible USB storage devices with partition scheme as Master Boot Records (MBR) and keyboards.

Rear Panel Overview


Figure 2. LENS 4K Camera Control Unit Rear Panel (REF 72205059 - Wi-Fi enabled)

Note: The rear panel of the non-Wi-Fi (REF 72205185), not shown, is the same as the Wi-Fi; minus the antenna and Wi-Fi symbol.

Rear Panel Connections		Function
1	DP	DisplayPort Connector: Video output port supporting Native 4K 60fps output.
2	HD-SDI Output	Serial digital interface with component digital video output for HD-SDI-compatible video displays or video transmission equipment. The HD-SDI connector support 1080i and 1080p. To configure the connector, select the Advanced Settings icon from the SYSTEM CONFIGURATION MENU screen.
3	HD-SDI (A, B, C, D)	Quad 3G HD-SDI video output connectors supporting native 4K video resolutions.
4	Wi-Fi antenna connector	Accommodates the removable antenna for communication with supported Wi-Fi devices.
5	Top Ethernet Port	For service use only.
	Bottom Ethernet Port	For use with Integration Broker installation.
6	Peripheral Controls	Two standard 1/8" mini phono plug (3.5 mm) peripheral cable connectors allow control of peripheral devices such as printers and image capture devices from the Camera Head.
7	Exhaust port for power supply	Exhaust port for venting of air through the system. To avoid the risk of the System overheating, do not block.
8	Fuse Holder	The clip holds a dual fuse. Refer to the System Specifications section of this manual.
9	AC Mains Power	The switch applies power to the control unit. The switch can be set to off (O) or on (I).
10	Power Cord Connector	Accommodates the hospital grade power cord accessory. This receptacle is an integral part of the power input module.
11	Equipotential Compensator Terminal (case ground)	Brings other equipment to the same case potential as the device.
12	USB Receptacle	The LENS 4K supports compatible USB storage devices with partition scheme as Master Boot Records (MBR) and keyboards.
13	Composite Input	Standard composite video (CV) provides a video input for the supported composite output devices (C-arms, operating room cameras, etc.). This input supports NTSC or PAL, but it must be configured for the desired input. Refer to the Advanced Settings section of this manual for information about how to configure the composite input but must be configured from the Advanced Settings screen to match the input to be used.

System Controls

Camera Head Controls

The LENS 4K Camera Head (Figure 3) is a three button Camera Head. Each button can be customized to provide two separate functions as part of a customized Camera Setup. The two functions are activated by a short press (<1.0 seconds) and a long press (>1.0 seconds). These buttons can be customized based on the user's needs. Refer to the Button Settings section of this manual for instructions regarding how to customize the buttons.


Figure 3. Smith & Nephew LENS 4K Camera Head

Button Function

Table 1, Live Video Functions, summarizes the default button functions for the Camera Head in the Live Video mode. Refer to the Button Settings section of this Operations/Service Manual for more information regarding default and customized button settings.

Camera Head Control	Short Press (<1.0 seconds)	Long Press (>1.0 seconds)
Left Button	Frame Capture (Single from Primary)	Light Mode
Middle Button	Brightness (Up/Down)	Menu Access
Right Button	Zoom	White Balance

*Button mapping is context-sensitive. Button settings may change, depending on the menu screen that is displayed on the monitor.

Table 1. Default Live Video Button Functions*

Table 2, Optional Programmable Button Functions, provides an overview of the programmable custom button functions available on the system. Customization can be set with the front panel buttons, Camera Head buttons, or with the Tablet Application. Refer to the Tablet Application Instructions for Use (REF 10601295) for information about how to customize settings with the Tablet Application.

Function	Effect
BRIGHTNESS	Control the luminance level of the video output
ENHANCEMENT	Enhance the sharpness of the displayed video
ZOOM	Digital zoom of field of view
ALT PAUSE	Pause the video while recording to a device connected to a rear panel port
PAUSE VIDEO	Pause the video from Camera Head
ALT VIDEO	Start or stop video capture from the secondary input source and save to the Tablet Application or USB
ALT FRAME	Capture a still picture from the secondary input source and save to the Tablet Application
FRAME CAPTURE	Capture a still picture from the endoscopic field of view and save to the Tablet Application or USB
LIGHT MODE	Activates and deactivates the illumination device on the control unit
VIDEO CAPTURE	Start or stop video capture from the endoscopic field of view and save to the Tablet Application or USB
LAST IMAGE	Displays the last known frame captures to the Tablet Application for review
MENU ACCESS	Access the graphical menu system of the control unit (MAIN MENU)
WHITE BALANCE	Correct for ambient color temperature
Peripheral Port 1 (Top port)	Triggers the accessory connected to the top port
Peripheral Port 2 (Bottom port)	Triggers the accessory connected to the bottom port

* Note: Software updates to either the LENS 4K Camera Head or Camera Control Unit, may reset the button-mappable functions to the factory settings.

Table 2. Optional Programmable Button Functions*

For information regarding how to customize Camera Head button settings, refer to the Button Settings section of this manual.

Installation

General Screen Navigation

Figure 4 Indicates the general information areas of the on-screen display.


Figure 4. General Screen Navigation

1. Patient information box. This box appears when the Tablet Application and/or a USB storage device are connected to the System.
 - a. patient name
 - b. number of images and videos recorded
 - c. current camera setup name

Note: When both a Tablet Application and a USB are connected, the Tablet Application information will appear on the OSD. Images and videos will be saved to both devices.

2. Device status bar
3. Screen name
4. Active menu area
5. Name of highlighted function
6. Button navigation map (top = short press (<1 seconds), bottom = long press (>1 seconds)).

Installation

Unpacking and General Inspection

Prior to using the LENS 4K Camera Control Unit, it is essential that all system components be inspected for damage, which can negatively impact the system's performance. Inspection should include all equipment to be used in surgery, including cables and peripheral devices.

1. Remove the control unit and accessories from the shipping container. Be certain to open and inspect all boxes and sealed plastic bags to ensure that all components are accounted for.

If any parts are missing or damaged, contact an authorized Smith & Nephew representative. Save the carton and packing materials in the event that a component needs to be returned for repair.

2. Inspect all video components for damage during shipment. If damage is found, contact an authorized Smith & Nephew representative. For more information, consult the Instructions for Use provided with the component.

LENS 4K Camera Control Unit

The following parts should have been received:

REF	Description
72205059	Smith & Nephew LENS 4K Camera Control Unit (1 ea)
72205185	- Wi-Fi enabled
	- non Wi-Fi
(included with REF 72205059, 72205185)	USB Type A Male/Female cable 3 ft. (1 ea)
(72204938 included only with REF 72205059)	Wi-Fi antenna (1 ea)
Refer to the country kits for the Ref Number.	Hospital-grade power cord, country-specific (1 ea)
	Smith & Nephew LENS 4K Camera Control Unit Operations/Service Manual, country-specific version (1 ea)
72201420	HD-SDI coaxial cable, 10 ft. (5 ea)
7206083	1/8" (3.5 mm) mono mini-plug to 1/8" (3.5 mm) mono mini-plug accessory cable, 6 ft. (1.8 m) (2 ea)
72205276	DisplayPort, LENS 4K CCU cable 10 ft. (1 ea)
72205277	Cat 5e, LENS 4K CCU cable 30 ft. (1 ea)

Setup

Only a qualified technician should install the System.

Secure the Control Unit

If the control unit will be placed on a platform or cart shelf equipped with a screw mount, follow these instructions to secure the control unit to the platform or shelf.

1. Center the control unit on the shelf and ensure that the screw hole located on the base of the control unit aligns with the captured #10-32 UN screw in the center of the shelf.
2. Hand-tighten the screw to secure the control unit to the shelf.

IEC 60601-1 Compliant System Configuration

⚠ Warnings

- If this unit is configured as part of a system, the entire system should be tested for compliance with IEC 60601-1.
- If the leakage current of the configured system exceeds the limits of IEC 60601-1, install an appropriately rated IEC 60601-1 approved isolation transformer with a minimum of 1000 VA and retest the system.
- The use of accessory equipment not complying with the equivalent safety requirements of this equipment may lead to a reduced level of safety of the resulting system. Consider the following when choosing accessory devices:
 - Use of the accessory in the patient vicinity.
 - Evidence that the safety certification of the accessory has been performed in accordance with IEC 60601-1.

Setup Configurations Summary

The System can be used with a variety of peripheral components to broaden the capabilities of the video system. Each setup is dependent upon the peripherals included in the video system. Contact an authorized Smith & Nephew representative for assistance with assembling the system to optimize system function.

Note: For best image, use a Quad 3G HDSDI signal (control unit output to monitor input). A direct connection between the control unit and a 4K monitor will always produce the highest quality image.

Equipment Required

- LENS 4K Camera Control Unit
- 4K Flat panel display

Optional

- 660HD or 660HD-E Image Management System(s)
- iPad and Tablet Application
- USB storage device

Cable Connectors

The System uses the following types of cables: BNC/SDI and BNC/HDSI coaxial cables, S-VHS cables which have a 4-pin mini-DIN connector, HD-DVI, and HD-SDI connectors.

BNC/SDI (composite)

To connect, push in and turn.


BNC/HDSI

To connect, push in and turn.

Note: BNC/SDI and BNC/HDSI connectors appear identical. For the highest quality image, be sure to use the cable provided with the System.

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To connect, push in firmly.

DisplayPort

To connect, push in firmly.

Ethernet

To connect, push in firmly

Installation

Recommended System Configuration for the Wi-Fi-Enabled System

This section explains the recommended equipment connections for optimum performance. Perform this procedure each time that the System is installed or reconfigured (Figure 5).

Required Equipment:

- LENS 4K Camera Control Unit, Wi-Fi-enabled (REF 72205059)
- LENS 4K Camera Head (REF 72205058)
- Wi-Fi Antenna
- 4K Flat Panel Display

Optional Equipment:

- iPad and Tablet Application
- USB storage device. The USB can be connected to either the USB port on the front panel or the port on the rear panel of the control unit.

Set-up

- For best image, use a DisplayPort connector (control unit output to monitor input). A direct connection between the control unit and a high definition monitor will always produce the highest quality image.
 - It is strongly suggested that a minimum of two outputs be connected at all times to provide a back-up signal if necessary.
1. Attach the antenna.
 2. Connect the supplied hospital-grade power cord to the power receptacle on the control unit's rear panel.
 3. To record still or video images from the camera, connect the Tablet Application to the LENS 4K Camera Control Unit Wi-Fi network or a insert a USB storage device in one of the USB ports on the control unit. refer to the System Specification section listed in this manual for a list of the requirements and compatibility. Refer to the Tablet Application IFU (REF 10601295) for details.
 4. Connect the supplied hospital-grade power cord to the power receptacle on the flat panel display and plug it in.
 5. If the flat panel display has a termination switch, set it to the 75 Ω or ON position.
 6. After connecting all components, including the Camera Head, turn the video components on using the appropriate power switches.

Once the Wi-Fi-enabled System has been turned on, the iPad can be connected using the Tablet Application. If using the non-Wi-Fi version of the System, the iPad and Tablet Application cannot be used.

Configure the peripheral ports

The default settings for the peripheral ports on the rear of the control unit are FRAME CAPTURE (Port 1, or the top port) and VIDEO CAPTURE (Port 2, or the bottom port). Each of the peripheral control ports can be configured for FRAME CAPTURE, GENERAL PURPOSE, or VIDEO CAPTURE. For information on how to customize the peripheral ports, refer to the Configure Peripheral Ports section of this manual.

Note: A USB storage device, iPad are not required in order to use the System.


Figure 5. Recommended setup configuration, Wi-Fi enabled System

Installation

Recommended System Configuration, Non Wi-Fi

This section explains the recommended equipment connections for optimum performance. Perform this procedure each time that the System is installed or reconfigured (Figure 6).

Required Equipment:

- LENS 4K Camera Control Unit (REF 72205185)
- Flat Panel Display
- LENS 4K Camera Head (REF 72205058)

Optional Equipment:

- Smith & Nephew 660HD or 660HD-E Image Management System: required for connection to hospital system (PACS/EMR through the Integration Broker device)
- USB storage device. The USB can be connected to either the USB port on the front panel or the port on the rear panel of the control unit

Set-up

Notes:

- For best image, use a DisplayPort or Quad 3G HDSDI signal (control unit output to monitor input). A direct connection between the control unit and a high definition monitor will always produce the highest quality image.
 - It is strongly suggested that a minimum of two outputs be connected at all times to provide a back-up signal if necessary.
1. Connect the supplied hospital-grade power cord to the power receptacle on the control unit's rear panel.
 2. Connect the supplied hospital-grade power cord to the power receptacle on the 660HD or 660HD-E Image Management System and plug it in.
 3. Connect the HD-SDI output, refer to Figure 2, connection #2 from the control unit to the HD-SDI input of the 660HD or 660HD-E Image Management System. Then connect the 660HD or 660HD-E HD-SDI output to the flat panel display. The 660HD image management system is compatible only with 1080i output. Refer to the Advanced Settings section of this manual for information about how to configure input and output settings.
 4. To record still images from the camera, connect the 3.5 mm mini phono plug from one of the control unit's peripheral controls outputs (Peripheral 1 or 2, or top port or bottom port), to the image management system's Still Capture IN port.

Note: The Camera Head buttons must be mapped so that they will trigger the peripheral ports. Refer to the Button Settings section of this manual for information about how to map the Camera Head buttons.

5. To record live video images from the camera, connect the 3.5 mm mini phono plug from one of the control unit's peripheral controls outputs (Peripheral 1 or 2, or top port or bottom port) to the image management system's Video Capture IN port.
6. Connect the supplied hospital-grade power cord to the power receptacle on the flat panel display and plug it in.
7. If the flat panel display has a termination switch, set it to the 75 Ω or ON position.
8. After connecting all components, including the Camera Head, turn the video components on using the appropriate power switches.

Configure the peripheral ports

The default settings for the peripheral ports on the rear of the control unit are FRAME CAPTURE (Peripheral 1, or the top port) and VIDEO CAPTURE (Peripheral 2, or the bottom port). Each of the peripheral control ports can be configured for FRAME CAPTURE, GENERAL PURPOSE, or VIDEO CAPTURE. For information on how to customize the peripheral ports, refer to the Configure Peripheral Ports section of this manual.

Fluoroscopy,
In-room camera,
etc.

LENS 4K Camera Control Unit Non Wi-Fi


Figure 4. Recommended setup configuration, Non Wi-Fi enabled

Inspect the System Components

Inspect the System Components

Prior to using the LENS 4K Camera Control Unit, it is essential that all system components be inspected for damage which can negatively impact the system's performance. Inspection should include all equipment to be used in surgery, including cables and peripheral devices.

CAUTION: Prior to each use, inspect the device to ensure it is functioning properly and is not damaged. Do not use a damaged device.

Electrical Connections

Examine the electrical connections.

- Electrocautery and other electrical noise-inducing medical equipment can interfere with the performance of control units and monitors. To prevent interference, plug monitors and camera equipment into an outlet on a wall separate from noise-inducing equipment.
- Check that the electrical equipment is properly grounded (i.e., plugs contain a ground prong). The control unit must be plugged into a hospital-grade AC outlet.
- If the monitor has a termination switch, it needs to be set to 75 Ω . If two or more monitors are used, only the termination switch on the last monitor needs to be on. If there is no termination switch on the last monitor, the monitor is self-terminating.

Double-check the equipment setup diagrams to ensure that all connections are correct.

Inspect the Fiber Optic Light Cable

 **WARNING:** When inspecting the light cable, NEVER aim the light cable at, or peer directly into the light source. Retinal damage may occur.

- Check the sheath for damage. Cuts, abrasions, or tears in the cable's silicone sheath will reduce overall light transmission.
- Aim one end of the cable toward a bright light, such as a room light, and inspect the other end for damaged fibers, e.g., black dots or dark gray areas (Figure 7). A combination of broken fibers in the cable or endoscope will result in reduced light transmission. The combined percentage of dark spots viewed in the cable end should not exceed 20% of its total area.


Figure 7. Inspect the fiber optic light cable

- Perform a visual fiber bundle diameter comparison. For optimal light transmission, the diameter of the fiber bundles in the endoscope's fiber optic light post should match the fiber bundle diameter of the light cable. This will prevent unnecessary heat buildup at the scope/light guide connection

Service Life Checklist

Perform a performance check for the installation and before each use case begins.

Check the following items to ensure that they function properly and the system is ready to use.

- ☐ Plug in the light guide and ensure that the light illuminates.
- ☐ Confirm that the video works correctly and that visualization is correct.
- ☐ Check the OSD status bar to confirm connection to the Wi-Fi.
- ☐ Check front panel to ensure that all buttons and connectors are functioning properly.
- ☐ Check the OSD status bar to confirm that the light is turned on and detected.
- ☐ Check the OSD status bar to confirm that the USB device is detected.
- ☐ Check for any damage to the camera head receptacle (cracks, corrosion, buildup or damage to the card edge receptacle).
- ☐ Check that the labeling information on the Camera Control Unit is still readable.

Apply Power to the Camera Control Unit

Notes:

- Perform the following steps only upon the initial power up of the system, if the system is powered down, or if a loss of power occurs.
- If the power to the system is controlled by the power switch on a cart or tower, these steps will not need to be repeated each time the system is powered up.

1. Make sure the AC mains power switch on the back panel is set to off (O).
2. Plug the female end of the power cord into the power receptacle on control unit rear panel.
3. Plug the power cord into a hospital-grade AC power outlet. Verify that the third pin on the plug is properly grounded.

4. Turn on the monitor.

Set the power switch module on the rear panel to (I). The green LED on the front panel blinks green, and the monitor displays a progress bar as the system configures itself.

Notes:

- All front panel buttons are disabled until the system is configured. The green LED on the front of the control unit blinks while the system configures itself. The green LED illuminates a steady green when the buttons are enabled.
 - Refer to the Alternative Power Up section for ways to power up the system in the event of a power failure.
6. When the LED changes to a constant green and the progress bar disappears, the monitor will display the READY screen (Figure 8).


Figure 8. READY screen

The READY screen displays the status of main and secondary video inputs, Wi-Fi, USB, and the light source. When each feature is ready and confirmed to operate correctly, it will appear orange. If the feature icon is not orange, use the following key to identify the problem:

A. Primary (left) and secondary (right) inputs:

Color	Explanation
Orange	Ready
White	Input is not ready
Red	Recording or busy

For example, in Figure 8, the left side of the icon is orange, which indicates that the Camera Head port is available for connection. The right side of the icon is white, which indicates that the secondary video input (image management system, fluoroscopic device, etc.) is not operational or not connected and should be checked.

B. Wi-Fi:

Attempt to connect the Tablet Application to the control unit in order to ascertain the status of the Wi-Fi connection.

Color	Explanation
Orange	Connected, valid patient
Yellow	Tablet Application is connected, no patient selected
White	Input is not working, need to check
Gray	Not connected

Note: The Wi-Fi icon notifies the user that the Tablet Application is communicating or not communicating with the control unit. If the Wi-Fi connection between the Tablet Application and the control unit is lost, the Wi-Fi icon on both the OSD status bar and the iPad will change from orange to gray. Should the connection be lost or interrupted during surgery, all control unit control features can be accessed either from the front panel of the control unit or the Camera Head buttons.

C. USB storage device:

Insert a USB to check its status.

Color	Explanation
Orange	USB detected
Yellow	Reconfiguring/Not ready
White	No USB detected
Red	Recording or busy

D. Illumination:

Insert a light guide into the appropriate port and press the standby illumination button to ascertain illumination status.

Color	Explanation
Orange	On
Gray	Off

Perform the Network Setup

Refer to the Customize the System, or the Network sections of this manual for more information about network setup.

Preoperative

Preoperative

Prepare the Camera Head

Clean and sterilize the LENS 4K Camera Head as described in the cleaning and sterilization procedures in the Instructions for Use received with the Camera Head.

Notes:

- Allow the Camera Head to cool to room temperature after it has been autoclaved.
- The Camera Head may be used immediately after chemical sterilization.

Connect the Camera Head

CAUTION: Do not plug the camera cable connector into the camera control unit if wet. Moisture on the camera cable connector, including the gold fingers, on the card edge of the camera head cable will damage the circuitry and void the warranty. Ensure that the camera cable connector is completely dry prior to plugging the camera head into the camera control unit.

Note: The Camera Head is not specific to any single LENS 4K Camera Control Unit. It can be used with multiple LENS 4K Camera Control Units.

To connect the camera cable connector to the Control Unit, plug the card edge connector into the Camera Head receptacle on the front of the Control Unit and push in firmly (Figure 9).


Figure 9. Connect the camera head to the control unit

Endocoupler focal length determines the image size on the monitor. For optimal image size use the LENS Coupler, 19.5 (72200315).

Connect the Coupler/Videoscope to the Camera Head

The System is designed to be used with instruments such as a videolaparoscope, videoarthroscope, or the coupler.

WARNING: Use aseptic technique in accordance with standard operating room procedures.

1. Attach a videoscope or coupler with a sterilized direct-view endoscope to the Camera Head.
2. Attach the fiber optic light guide from the light source to the scope.

Insert Light Guide

Warnings

- During operation, avoid prolonged contact of the scope tip to patient tissue or flammable materials. The scope tip may reach high temperatures due to high intensity light transmission.
- Do not leave the operating light cable on a patient or the drapes. Failure to observe this warning may result in burns to the patient and/or the surrounding drapes.

Set the multiport light guide adaptor turret for the light guide that will be used. To set the turret to the appropriate light port, rotate the turret in either direction until the desired light port is aligned with the orange arrow to the left of the turret (Figure 10).

Insert the light cable into the appropriate light port by pushing in firmly. To remove the light cable, grasp the cable connector and pull straight out of the light port. Do not pull out by pulling on the cable.

CAUTION: Ensure that the light guide is inserted into the proper front panel port. Damage to the unit may occur if the light guide is not inserted into the port associated with the correct light guide manufacturer.


Figure 10. Multiport light guide adaptor turret

Turn the Camera Control Unit On

1. Press the **Standby Power** button on the front of the control unit. The button illuminates green. This Power button acts as a standby button for the System.

The buttons on the front panel of the control unit activate and the **STARTUP** menu (Figure 11) opens. The default camera setup is **Arthroscopy**.


Figure 11. Arthroscopy

The **STARTUP** menu remains on the screen until a camera setup or the **Menu Access** icon is selected.

Notes:

- When a camera setup is selected, the camera will automatically perform a white balance. If a Camera Head is connected, ensure that the Camera Head is pointed at a white object such as a piece of gauze or a piece of paper. At this point, the camera is in Live Video mode.
- Camera setups stored on the Tablet Application will not appear on the **STARTUP** menu. However, the name of the camera setup in use appears in the patient information box on the OSD.
- If a Camera Head is not plugged in, or if the Camera Head is unplugged while the control unit is on, the monitor will display a color bar pattern. When the Camera head is plugged back in, perform a White Balance.

Press the **Illumination** button on the front of the control unit to turn on the light source. The button illuminates orange.

Connect the iPad

If the control unit in use has a Wi-Fi feature, the Tablet Application can now be connected to the system. Refer to the Tablet Application Instructions for Use (REF 10601295) for information about how to download the Tablet Application and connect it to the LENS 4K Camera Control Unit.

Note: The Wi-Fi icon notifies the user that the iPad is communicating or not communicating with the control unit. If the Wi-Fi connection between the Tablet Application and the control unit is lost, the Wi-Fi icon on both the OSD status bar and the Tablet Application will change from orange to gray. The Tablet Application will automatically connect when the signal is available again. Should the connection be lost or interrupted during surgery, all control unit control features can be accessed either from the front panel of the control unit or the Camera Head buttons.

Operation

During a procedure, the user can access any of the functions available in the system either through the Camera Head button mappings or via the **Select/Menu Access** button on the front panel of the control unit. Any settings changed during the camera setup will be temporarily saved until the system is powered down. To customize camera setups and save settings permanently, refer to the **Customize the System** section of this manual. Refer to the menu maps for an overview of how to access System screens.

When setup is complete and service life checks have been performed, the System is ready for the surgeon to use. The list below summarizes the initial portion of the operative case:

1. Obtain a sterile Camera Head.
2. Plug the Camera Head into the control unit.
3. Connect the light guide to the appropriate light guide receptacle on the front of the control unit.
4. Turn on the illumination by pressing the **Illumination** button
5. White Balance from within the sterile field.

For detailed information about how to operate the System, read this **Operations** section.

Select a Camera setup

To select a camera setup, use the **Up** and **Down** arrow buttons or the left and right Camera Head buttons to highlight the desired camera setup. Press the **Select** button. The default camera setup is the **Arthroscopy** camera setup.

To change the default camera setup, refer to the **System Configuration** section of this manual. To customize a camera setup, refer to the **Customize a Camera setup** section of this manual.

To access the **MAIN MENU** screen, select the **Menu Access** icon on the **STARTUP** screen (Figure 12).


Figure 12. Select Camera setup Access

MAIN MENU

The MAIN MENU screen (Figure 13) allows the user to access and change settings for the various features and functions of the System. Use the Up and Down arrows on the control unit or the left and right Camera Head buttons to navigate to the desired feature or function.

Note: Settings changed directly from the MAIN MENU are temporary settings and will not be saved after the system is powered down or another camera setup is selected. To save settings permanently, refer to the Camera setups section of this manual.


Figure 13. MAIN MENU screen

Brightness

To adjust Brightness, navigate to the Brightness icon and press Select. Use the Up and Down arrows on the control unit or the left and right Camera Head buttons to adjust Brightness to the desired level.

Zoom

To adjust Zoom, navigate to the Zoom icon and press Select. Use the Up and Down arrows on the control unit or the left and right Camera Head buttons to adjust Zoom to the desired level.

Enhancement

To adjust Enhancement, navigate to the Enhancement icon and press Select. Use the Up and Down arrows on the control unit or the left and right Camera Head buttons to adjust Enhancement to the desired level.

White Balance

To perform a White Balance, a Camera Head and scope must be connected to the front of the control unit. Point the end of the scope at a white object (a 4" x 4" gauze pad or flat white piece of paper) and focus (Figure 14). Fill as much of the screen as possible with the white object, but do not hold the scope close enough to touch the white object. Press the White Balance button at the upper right corner of the Camera Head connector (Figure 11). When the camera completes the White Balance process, the message White Balance Complete appears briefly on the screen, and a long tone sounds to indicate that the White Balance is complete. If the White Balance fails, a chirping tone will sound.

Notes:

- Each time the White Balance button is pressed, a single tone will sound. Press Select to return to the MAIN MENU.

- When a camera setup is chosen from the STARTUP screen, the camera automatically performs a White Balance.


Figure 14. White Balance

Camera Setups

Camera setups consists of 1-3 screens:

- Setup Selection
- Customize Camera Setups
- Set Default Startup Camera setup

Setup Selection

To select or customize a camera setup, or to set the default camera setup, navigate to the Camera setups icon on the MAIN MENU and press Select. The Camera setups menu opens (Figure 15).

Figure 15. Setup Selection

Camera setup Selection

To select a camera setup from the Camera setup menu, navigate to the SETUP SELECTION icon and press Select. The CAMERA SETUPS menu (Figure 16) opens.


Figure 16. Camera Setup Selection

Because Arthroscopy is set as the default Camera Setup, it will blink. Including the Laparoscopy, ENT, MOIRE, ARTHROSCOPY, the system contains 10 customizable camera setups. The Moire camera setup defaults to the optimal settings for fibroscopes. For more information about Camera Setups, refer to the Camera setup section of this manual.

Use the Up and Down arrows on the control unit or the left and right Camera Head buttons to navigate to the desired camera setup and press **Select** to select that camera setup.

Customize a Camera setup

To customize a camera setup, highlight the **Customize Camera setups** icon on the CAMERA SETUPS menu (Figure 17) and press **Select**.

Figure 17. Customize Camera Setups

Press **Select** to open the SETUP SELECTION screen (Figure 18). Refer to the Customize Camera setup Settings section of this manual for information about how to customize or create a camera setup.

Figure 18. Select Setup Selection to Customize

To exit the SELECT Camera setup menu, highlight the **Menu Access** icon and press the **Select** button.

Set Default Startup Camera setup

Set Default CAMERA SETUP allows the user to set the camera setup that is highlighted on the STARTUP screen and so can be selected immediately. For information about how to customize a default camera setup, refer to the Customize the System section of this manual.

Color Bar

To calibrate the monitor, navigate to the **Color Bar** icon and press **Select**. The monitor screen shows the color bars. Adjust the color bars on the monitor according to that monitor's instructions for use. To exit, press any button or wait approximately 30 seconds.

System Configuration

The System is shipped with a default configuration. For information about how to change the configuration, refer to the Customize the System section of this manual.

Buttons

To view the current button configuration, navigate to the **Buttons** icon and press **Select**. The button map in the lower right corner of the monitor will indicate the current short (top setting) and long (bottom setting) press button configurations for the Camera Head. To change the button configurations, refer to the Button Settings section of this manual.

Exit

To exit the SYSTEM CONFIGURATION MENU screen, navigate to the **Exit** icon and press **Select**.

Capture Images

Ensure that the camera icon on the OSD status bar indicates **Ready** and that the button mappings are set for image capture. The System defaults to Live Video after a White Balance has been performed. Images can be captured when the System is in Live Video mode. The button map for the selected camera setup appears in the lower right corner of the OSD. Refer to Table 1 in the Camera Head section of the System Controls section of this manual for the System's default button settings for the Camera Head.

Note: When both an iPad and a USB or other image storage device are connected, the information displayed in the patient information box (Figure 4) is the information from the Tablet Application. The Tablet Application supersedes the USB information on the OSD, but information is recorded to both devices.

To Capture a Still Image

To record a still image, firmly press the button mapped to **Frame Capture** (or **ALT FRAME**). If using the Tablet Application or a USB for image storage, when the image is captured, the image counter in the upper left corner of the OSD increases by one. If the counter does not increase, the image has not been recorded. If the image counter does increase, the image is saved to the connected image storage devices. For information regarding how to capture still images from the Tablet Application, please refer to Tablet Application Instructions for Use (REF 10601295). For information about how to customize button settings, refer to the Button Settings section of this manual.

Note: If the Tablet Application and a USB are connected to the System, the captured image will be saved to both devices.

If using the 660HD or 660HD-E, images will be stored directly to that device. Buttons must be mapped to the appropriate peripheral port. The patient information box will not appear on the OSD.

Security Settings

To Capture Video Images

Tablet/USB capture only:

To record a video image, press the button mapped to **Video Capture** and release the button. To stop recording, press and release the **Video Capture** button again. The camera icon on the OSD remains red until the video has been transferred to the storage device(s) and confirmed as saved. The video counter in the upper left corner of the OSD increases by one. If the counter does not increase, the video has not been recorded.

For information regarding how to capture video images from the Tablet Application, please refer to Tablet Application Instructions for Use (REF 10601295). For information about how to customize button settings, refer to the Button Settings section of this manual.

If capturing video to the 660HD or the 660HD-E image management system, ensure that the ports are mapped to the Camera Head buttons (Acc1 or Acc2) correctly. Videos will be saved directly to the image management system.

Security Settings

The default factory settings do not require passwords and allow access to all system functions. Each site can choose to customize the security settings as part of its security processes.

The user password allows the user to prevent unwanted access to their customized camera setup settings. The administrator password allows the site to prevent any changes to the system settings.

For more information, refer to the Set/Change Password Protection and Set/Change Admin Password sections of Customize the System in this manual.

Resetting the system if the user password is lost

If the user password is lost, the system can be reset to the factory defaults from the **Reset Factory Defaults** icon on the **SYSTEM CONFIGURATION MENU** screen.

Note that if an administrator password has been set, the System will ask the user for the administrator password before allowing passwords to be reset.

Note: The system clears access either on power down or by accessing the power standby button on the front panel.

Alternative Power Up

The System can be configured to automatically restart under conditions set by the administrator. In the event of a mains power failure, it can be set to automatically restart, always, never, or only in the case of a mains power loss.


ALWAYS skips the READY screen and goes directly to the STARTUP screen.

NEVER skips the READY screen. The user needs to press the Standby button to turn on system power.

POWER LOSS only skips the READY screen upon power loss.

For information about how to configure **AUTO START**, refer to the To configure Automatic Start section of the Camera setups section of this manual.

Menu Maps

Postoperative

To avoid possible burns or fire, be sure to turn off the illumination. To turn off the illumination, either press the **Standby illumination** button or unplug the light guide.

Turn the System off to erase all temporary changes to the current camera setups and ensure security by ending administrative password access. To turn off the System, press the **Standby power** button or turn off the power switch on the tower.

To print:

- From the Tablet Application, refer to the Tablet Application's Instructions for Use (REF 10601295) for information.
- From the 660HD Image Management System, refer to that product's manual for information (REF 10600055).
- From the 660HD-E Image Management System, refer to that product's manual for information (REF 10601236)
- If a USB is used for data storage, all images are saved as MP4 files to a media directory. To print from the media directory, the user must have the directory structure and know how to navigate the directory (Figure 19).


Figure 19. Media Directory

If the Tablet Application is connected and a patient file was created, and a USB storage device is connected to the control unit, media files will be saved to that patient's MRN directory on the USB (Figure 20).


Figure 20. MRN Directory

Customize the System

The System is shipped with a default configuration. The System's default settings can be customized and saved for any procedure. Navigate to the **System Configuration** icon on the MAIN MENU and press **Select** to open the SELECT SYSTEM CONFIGURATION MENU screen.

1. From the MAIN MENU, select the **System Configuration** icon to open the SYSTEM CONFIGURATION MENU screen (Figure 21).

Figure 21. System Configuration Menu

- The SYSTEM CONFIGURATION MENU contains the following icons:
- Set/Change User Password Protection
 - Set/Change Admin Password.
 - Configure Accessory Ports
 - Language
 - Reset Factory Defaults
 - Network
 - Advanced Settings
 - Versions. Lists software versions of attached devices.
 - Exit
2. Navigate to the icon representing the feature to be customized and press **Select**.

Set/Change Password Protection

A user password enables an individual to change and protect the settings for all the custom camera setups. No one else will be able to access the camera setup's settings.

To set or change an individual's user password:

1. Select the **Set/Change User Password** icon on the SYSTEM CONFIGURATION MENU screen to open the CHANGE PASSWORD screen (Figure 22). **PASSWORD REQUIRED** is blinking. The default is set to **NO**.


Figure 22. Set/Change Password

Set/Change Password

2. Select **PASSWORD REQUIRED**. The current setting (YES or NO) will blink.
3. Use the Up and Down arrows on the control unit or the left and right Camera Head buttons to switch between YES or NO.
4. Press **Select** to select either Yes or NO. To save the setting, navigate to **SAVE AND EXIT** and press **Select**. To return to the System Configuration menu without saving any changes, select **CANCEL**.

If Yes is selected, a blinking **NEW PASSWORD =** option appears on the screen (Figure 23).


Figure 23. New Password

5. Press **Select** to select **NEW PASSWORD =**. An OSD keyboard appears below the **CHANGE PASSWORD** menu (Figure 24).


Figure 24. Change Password

6. Use the Up and Down arrows on the control unit or the left and right Camera Head buttons to navigate to the desired letters on the keyboard.
 - A short press moves the highlighted cursor left or right.
 - A long press moves the cursor up or down from line to line.
 - To select a letter, number, or symbol, press **Select**.
 - To exit the keyboard, navigate to the **Return** key and press **Select**.
7. To save the changes and exit the **CHANGE PASSWORD** screen, select **SAVE AND EXIT** to return to the **SYSTEM CONFIGURATION MENU** screen. To exit the screen without saving any changes, navigate to **CANCEL** and press **Select** to return to the **SYSTEM CONFIGURATION MENU** screen.

Set/Change Admin Password

The administrator password allows the site to prevent any changes to the system settings.

Note: A user will still be able to create a user password to allow controlled access to the **CAMERA SETUPS** menu.

To set or change the administrative password:

1. Select the **SET/CHANGE ADMIN PASSWORD** icon to open the **CHANGE PASSWORD** screen.
2. Select **PASSWORD REQUIRED**. The current setting (YES or NO) will blink.
3. Use the Up and Down arrows on the control unit or the left and right Camera Head buttons to switch between YES and NO.
4. Press **Select** to select either YES or NO. To save the setting, navigate to **SAVE AND EXIT** and press **Select**. To return to the System Configuration menu without saving any changes, select **CANCEL**.

If Yes is selected, a blinking **NEW PASSWORD =** option appears on the screen (Figure 25).

Figure 25. Set/Change Admin Password

5. Press **Select** to select **NEW PASSWORD =**. An OSD keyboard appears below the **CHANGE PASSWORD** menu (Figure 26).


Figure 26. New Admin Password

6. Use the Up and Down arrows on the control unit or the left and right Camera Head buttons to navigate to the desired letters on the keyboard.

Note: A short press moves the highlighted cursor left or right. A long press moves the cursor up or down from line to line. To select a letter, number, or symbol, press **Select**. To exit the keyboard, navigate to the **Return** key and press **Select**.

Customize the System

7. To save the changes and exit the CHANGE PASSWORD screen, select **SAVE AND EXIT** to return to the SYSTEM CONFIGURATION MENU screen. To exit the screen without saving any changes, navigate to **CANCEL** and press **Select** to return to the SYSTEM CONFIGURATION MENU screen.

Note: Do not lose the administrative password. This will prevent the recovery of passwords through the RESET FACTORY DEFAULTS function. If the administrator password is lost, please contact Smith & Nephew.

Configure Peripheral Ports

The CONFIGURE ACCESSORY PORTS MENU screen allows the user to configure the peripheral ports on the rear panel of the control unit to the desired function. To configure the peripheral ports:

1. Select the CONFIGURE ACCESSORY PORTS icon on the SYSTEM CONFIGURATION MENU screen to open the CONFIGURE ACCESSORY PORTS MENU screen (Figure 27). TOP PORT = will blink. The default setting is FRAME CAPTURE.


Figure 27. Configure Accessory Ports

2. To customize the TOP PORT, press **Select**. The current setting to the right of TOP PORT will begin to blink.
3. Use the Up and Down arrows on the control unit or the left and right Camera Head buttons to scroll through the available options
4. Navigate to the desired setting (Figure 28) and press **Select**:
 - Select **FRAME CAPTURE** the CCU will send a still frame of the digital image captured on the device.
 - Select **GENERAL PURPOSE** to send a still frame to the printer.
 - Select **VIDEO CAPTURE** to send video recordings to the digital image capture device. When VIDEO CAPTURE is in use, the first button press starts the recording, and the second button press stops the recording.

FRAME CAPTURE
GENERAL PURPOSE
VIDEO CAPTURE

Figure 28. Select Setting

5. To customize the settings for the BOTTOM PORT, navigate to **BOTTOM PORT** and press **Select**. Repeat Steps 2 and 3 to change the current settings for **BOTTOM PORT**.
6. To save changes and exit the screen, select **SAVE AND EXIT**. To exit the screen without saving any changes, return all settings to their original settings and select **Save and Exit**.

Language

To customize the language that appears on the OSD:

1. Select the LANGUAGE icon on the SYSTEM CONFIGURATION MENU screen to open the LANGUAGE MENU screen (Figure 29).

Figure 29. Language Menu

2. Use the Up and Down arrows on the control unit or the left and right Camera Head buttons to scroll through the available languages.
3. When the desired language blinks, press **Select** to select it. The System will return to the SYSTEM CONFIGURATION MENU screen. The OSD will change to the language that was selected.

Reset Factory Defaults

To reset the System to the original factory default settings:

1. Select the RESET FACTORY DEFAULTS icon on the SYSTEM CONFIGURATION MENU screen to open the RESET FACTORY DEFAULTS screen (Figure 30). RESET FACTORY DEFAULTS will blink.

Figure 30. Open Reset Factory Defaults

Select **EXIT** to exit the RESET FACTORY DEFAULTS screen and return to the SYSTEM CONFIGURATION MENU screen without resetting to the factory defaults

2. To reset the System to the original factory default settings, select **RESET FACTORY DEFAULTS**. A confirmation screen containing the question **ARE YOU SURE?** will open (Figure 31). **YES** will blink


Figure 31. Select Reset Factory Defaults

3. Select YES to restore the original factory default settings. The System will replace all previously customized settings for all camera setups and advanced settings with the original factory default settings and return to the SYSTEM CONFIGURATION MENU screen.

Select NO to retain previously saved settings and exit the RESET FACTORY DEFAULTS screen without making any changes. The System returns to the SYSTEM CONFIGURATION MENU screen.

Note: Ensure to power cycle the system, after re-setting the factory defaults.

Network

Note: This section applies only to the Wi-Fi enabled version of the control unit.

The NETWORK SETTINGS screen allows the user to change the settings that enable the System network to connect to the Tablet Application. The network ID and network password are required in order to connect the Tablet Application to the System network. To change the network settings:

1. Navigate to the NETWORK icon on the SYSTEM CONFIGURATION MENU screen. Press Select to open the NETWORK SETTINGS screen (Figure 32).


Figure 32. Network Settings

2. To reset the Wi-Fi or image management system, select RESET WIFI, whichever one appears on the screen.

To change the network ID:

1. Navigate to NETWORK ID = and press Select. An OSD keyboard appears below the NETWORK SETTINGS menu.
2. Use the Up and Down arrows on the control unit or the left and right Camera Head buttons to navigate to the desired letters on the keyboard.

Note: A short press moves the highlighted cursor left or right. A long press moves the cursor up or down from line to line. To select a letter, number, or symbol, press Select. To exit the keyboard, navigate to the Return key and press Select.

3. To save the changes and exit the NETWORK SETTINGS screen, select SAVE AND EXIT to return to the SYSTEM CONFIGURATION MENU screen. To exit the screen without saving any changes, navigate to CANCEL and press Select to return to the SYSTEM CONFIGURATION MENU screen.

To change the network password:

1. Navigate to NETWORK PW = and press Select. An OSD keyboard appears below the NETWORK SETTINGS menu.
2. Use the Up and Down arrows on the control unit or the left and right Camera Head buttons to navigate to the desired letters on the keyboard.

Note: A short press moves the highlighted cursor left or right. A long press moves the cursor up or down from line to line. To select a letter, number, or symbol, press Select. To exit the keyboard, navigate to the Return key and press Select.

3. To save the changes and exit the NETWORK SETTINGS screen, select SAVE AND EXIT to return to the SYSTEM CONFIGURATION MENU screen. To exit the screen without saving any changes, navigate to CANCEL and press Select to return to the SYSTEM CONFIGURATION MENU screen.

MAC address:

All LENS 4K System control unit's ships with a set MAC address and is not changeable.

The MAC address is the address that identifies a particular LENS 4K control unit on the network.

To change the network channel:

There are 12 available network channels available in the System: channels 0–11. If the System performance seems to be slowing or otherwise degrading, change to a different network channel to try to resolve the problem. For more information about how to troubleshoot a slowing system, refer to the Wi-Fi Troubleshooting section of this manual.

To change the network channel:

1. Navigate to NET CHANNEL =. When NET CHANNEL = blinks, press Select to select it. The current channel will blink.
2. Use the Up and Down arrows on the control unit or the left and right Camera Head buttons to navigate to the desired channel and press Select.
3. To save the changes and exit the NETWORK SETTINGS screen, select SAVE AND EXIT to return to the SYSTEM CONFIGURATION MENU screen. To exit the screen without saving any changes, navigate to CANCEL and press Select to return to the SYSTEM CONFIGURATION MENU screen.

IP Address:

This section is used to indicate what IP address has been assigned by the hospital network to the CCU if enabling the optional communication to Integration Broker. Refer to Network section listed in this document. If the system has not been given an IP address, this field will indicate NO IP.

Note: Refer to Configure the Settings for Integration Broker section in the Tablet Application document #10601295).

Customize the System

To activate FCC power:

FCC power is a regional setting that affects Wi-Fi signal strength. In the United States, FCC power is typically set to YES. In the European Union, FCC power is typically set to NO. If the System is located outside of the US and EU, contact a local sales representative for information regarding the optimal System setting.

To activate FCC power:

1. Navigate to FCC POWER =. When FCC POWER = blinks, press Select to select it. The current setting will blink.
2. Use the Up and Down arrows on the control unit or the left and right Camera Head buttons to switch between YES and NO. Press Select to select the desired setting.
3. To save the changes and exit the NETWORK SETTINGS screen, select SAVE AND EXIT to return to the SYSTEM CONFIGURATION MENU screen. To exit the screen without saving any changes, navigate to CANCEL and press Select to return to the SYSTEM CONFIGURATION MENU screen.

ADVANCED SETTINGS

The ADVANCED SETTINGS feature allows the user to configure the System's video signal settings. If an administrative password has been set up on the system, that password will be needed.

1. From the SYSTEM CONFIGURATION MENU screen (Figure 33), use the Up and Down arrows on the control unit or the left and right Camera Head buttons to navigate to the Advanced Settings icon. Press Select to open the ADVANCED SETTINGS screen (Figure 34). To exit the ADVANCED SETTINGS screen, navigate to SAVE AND EXIT and press Select.

Figure 33. System Configuration Menu


Figure 34. Advanced Settings

To configure HDSOI Outputs:

1. Press Select to select the blinking HDSOI OUTPUTS =. The current setting will blink.
2. Use the Up and Down arrows on the control unit or the left and right Camera Head buttons to switch to 1080i. Refer to Figure 2, connection #2. Press Select to activate the output.

Note: The System does not allow the selection of an incompatible output. If the monitor in use is not compatible with the selected output, the screen will flash and return to the CONFIGURE SETTINGS menu after approximately 10 seconds.

To configure COMPRESSION:

COMPRESSION provides the user with the ability to adjust the video compression bit rate.

Note: The higher the number, the less compression of the video, which will result in higher quality and longer transfer times.

1. Navigate to COMPRESSION. Press Select to select the blinking COMPRESSION.
2. Use the Up and Down arrows on the control unit or the left and right Camera Head buttons to select the video bit rate from 1 to 30. Press Select to activate the desired output.
3. Select No to return to the Advanced Settings screen without saving any changes. Select YES to save the change and return to the Advanced Settings screen.

To configure AUTO START:

The AUTO START option allows the user to configure how the system behaves on a mains power-up. To configure when the System should automatically restart:

1. Navigate to AUTO START. When AUTO START blinks, press Select to select it. The current setting will blink.
2. Use the Up and Down arrows on the control unit or the left and right Camera Head buttons to scroll through the available options (Figure 35). Press Select to activate the desired option. Press Select to activate the desired option.

Figure 35. AUTO START

To configure the HEAD PLUGIN

The HEAD PLUGIN = option allows the user to customize the System's behavior when the Camera Head plugged in.

To configure HEAD PLUGIN:

1. Navigate to HEAD PLUGIN =. When HEAD PLUGIN = blinks, press Select to select it. The current setting will blink.
2. Use the Up and Down arrows on the control unit or the left and right Camera Head buttons to switch between POWER and NOTHING..
3. Select POWER to automatically turn on the System upon connection of the Camera Head to the System.

Select NOTHING to not automatically turn the system on upon connection of the Camera Head to the System.

To configure VIDEO INPUT:

The VIDEO INPUT option allows the user to configure the video input format for either NTSC or PAL on the rear panel.

Note: NTSC inputs are used in North America and parts of South America. PAL inputs are typically used in the EU and other countries. Please contact an authorized Smith & Nephew representative with any questions regarding which input is correct for your location.

To configure VIDEO INPUT:

1. Navigate to VIDEO INPUT. When VIDEO INPUT blinks, press Select to select it. The current setting will blink.
2. Use the Up and Down arrows on the control unit or the left and right Camera Head buttons to switch between NTSC and PAL.
3. Select NTSC for an NTSC video input. Select PAL for a PAL video input.

To save the settings and exit the CONFIGURE SETTINGS screen, navigate to SAVE AND EXIT and press Select to return to the SYSTEM CONFIGURATION MENU screen.

Versions

Versions allows the user to view the software version currently in use on the system. To view the software versions currently in use:

1. Navigate to VERSIONS. When VERSIONS blinks, press Select to open the VERSIONS screen.
2. To exit the VERSIONS screen, press any button to return to the SYSTEM CONFIGURATION screen.

Customize Camera setup Settings

1. From the MAIN MENU, highlight the camera setups icon and select it to open the Camera Setups menu.
2. Use the Up and Down arrows on the control unit or the left and right Camera Head buttons to navigate to the CUSTOMIZE CAMERA SETUPS icon and press Select to open the SETUP SELECTION screen.
3. Navigate to the camera setup (Figure 36) to be customized and press Select to open that camera's setup customization screen. The CUSTOMIZE BUTTON SETTINGS option will blink (Figure 39).

Note: Changes to the camera setup are not saved until SAVE AND EXIT is selected.


Figure 36. CUSTOMIZE CAMERA Setups SETTINGS screen

The Camera setup screen contains a MORE SETTINGS menu; refer to the MORE SETTINGS section in this document.

Customize Button Settings

Camera Head button settings can be customized to the user's preferred settings for short and long button presses. To customize button settings:

1. Press Select to open that Customize button Settings screen (Figure 37).

Figure 37. Camera Head button settings

2. Use the Up and Down buttons on the control unit or the right and left buttons on the Camera Head to navigate to the button press to be customized (MIDDLE SHORT, RIGHT LONG, etc.).
3. Press Select to select the desired button press. To the right of the button press name, the current setting for that button press will begin to blink.
4. Use the Up and Down buttons on the control unit or the right and left buttons on the Camera Head to scroll through the list of options (Figure 38). Press Select to select the setting for that button press.

BRIGHTNESS
ENHANCEMENT
ZOOM
TOP PORT
BOTTOM PORT
PAUSE VIDEO
ALT PAUSE
ALT VIDEO
ALT FRAME
FRAME CAPTURE
LIGHT MODE
VIDEO CAPTURE
LAST IMAGE
MENU ACCESS
WHITE BALANCE

Figure 38. Up and Down buttons

Refer to Table 2 in the System Controls section of this manual for an explanation of each setting option.

5. When all button presses have been customized as desired, use the Up and Down buttons on the control unit or the right and left buttons on the Camera Head to navigate to the Exit button. Select Exit to save the settings, exit the button settings screen, and return to the camera setup's customization menu.

Customize Camera setup Settings

IMAGE CAPTURE OSD

The IMAGE CAPTURE OSD setting allows the user to determine if the on screen display graphics appear in the image and video captures. The default setting is OFF. To customize Image Capture OSD:

1. Navigate to IMAGE CAPTURE OSD (Figure 39). When IMAGE CAPTURE OSD blinks, press Select to select it. The current setting will blink.

Figure 39. Image Capture OSD

2. Use the Up and Down arrows on the control unit or the left and right Camera Head buttons to switch between OFF and ON. Press Select to select the desired setting.
3. To save the changes and exit the Camera setups Settings screen, select SAVE AND EXIT.

DISPLAY IN STARTUP MENU?

The DISPLAY IN STARTUP MENU? setting allows the user to choose whether or not to display the customized camera setup in the STARTUP menu. To customize DISPLAY IN STARTUP MENU?:

1. Navigate to DISPLAY IN STARTUP MENU?. When DISPLAY IN STARTUP MENU? blinks, press Select to select it. The current setting will blink.
2. Use the Up and Down arrows on the control unit or the left and right Camera Head buttons to switch between OFF and ON. Press Select to select the desired setting.
3. To save the changes and exit the CAMERA SETUPS SETTINGS screen, select SAVE AND EXIT.

Note: This setting cannot be changed for the ARTHROSCOPY or TOWER camera setup. The ARTHROSCOPY or TOWER camera setup will always appear on the STARTUP screen.

PORTRAITS ENABLED?

The PORTRAITS ENABLED? setting allows the user to choose whether or not to briefly display image portraits on the OSD during image capture.

1. Navigate to PORTRAITS ENABLED?. When PORTRAITS ENABLED? blinks, press Select to select it. The current setting will blink.
2. Use the Up and Down arrows on the control unit or the left and right Camera Head buttons to switch between YES and NO. Select YES to have image portraits appear in the upper right corner of the OSD. Select NO to prevent image portraits from appearing on the OSD.
3. To save the changes and exit the CAMERA SETUPS SETTINGS screen, select SAVE AND EXIT.

CHANGE PROCEDURE NAME

The CHANGE PROCEDURE NAME setting allows the user to rename a camera setup. To rename a camera setup:

1. Navigate to CHANGE PROCEDURE NAME. When CHANGE PROCEDURE NAME blinks, press Select to select it. The RENAME PROCEDURE NAME screen will open and the current camera setup name will blink.
2. Press Select to select the current camera setup name. An OSD keyboard will appear below the RENAME PROCEDURE NAME screen.
3. Use the Up and Down arrows on the control unit or the left and right Camera Head buttons to navigate to the desired letters on the keyboard.

Note: A short press moves the highlighted cursor left or right. A long press moves the cursor up or down from line to line. To select a letter, number, or symbol, press Select. To exit the keyboard, navigate to the Return key and press Select.

4. When the new name of the camera setup is spelled out, exit the keyboard. ENTER will blink.
5. Select ENTER to change the camera setup name and return to the CAMERA SETUPS SETTINGS screen. Select CANCEL to keep the original camera setup name and return to the CAMERA SETUPS SETTINGS screen.
6. To save the changes and exit the CAMERA SETUPS SETTINGS screen, select SAVE AND EXIT.

BRIGHTNESS

The BRIGHTNESS setting allows the user to set the starting luminance level of the video output during the camera setup. Increase BRIGHTNESS to better illuminate the parts of an image that are normally dark, such as a recess in a knee.

BRIGHTNESS settings range from 0 to 10 in whole number increments. The default value is +5. Once the starting BRIGHTNESS level has been set for a camera setup, BRIGHTNESS can be increased or decreased from that point within the available range. To customize the BRIGHTNESS:

1. Navigate to BRIGHTNESS. When BRIGHTNESS blinks, press Select to select it. The current setting will blink.
2. Use the Up and Down arrows on the control unit or the left and right Camera Head buttons to scroll through the BRIGHTNESS level settings. Press Select to select the desired level.
3. To save the changes and exit the CAMERA SETUPS SETTINGS screen, select SAVE AND EXIT.

ENHANCEMENT

The ENHANCEMENT setting allows the user to set the starting enhancement level of the video output during the camera setup. Increase ENHANCEMENT to sharpen the edges in the image and increase contrast. Decrease ENHANCEMENT to soften the edges in the image and decrease contrast.

ENHANCEMENT settings range from 0 to 6 in whole number increments. The default value is +3. Once the starting ENHANCEMENT level has been set for a camera setup, ENHANCEMENT can be increased or decreased from that point within the available range. To customize ENHANCEMENT:

1. Navigate to **ENHANCEMENT**. When **ENHANCEMENT** blinks, press **Select** to select it. The current setting will blink.
2. Use the Up and Down arrows on the control unit or the left and right Camera Head buttons to scroll through the **ENHANCEMENT** level settings. Press **Select** to select the desired level.
3. To save the changes and exit the **CAMERA SETUPS SETTINGS** screen, select **SAVE AND EXIT**.

CHROMA (SAT)

The **CHROMA (SAT)** setting allows the user to set the color saturation level of the video output during the camera setup. Increase saturation to increase the intensity of the colors. Decrease saturation to decrease the intensity of the color. A saturation level of zero produces a monochrome image.

CHROMA (SAT) settings range from 0 to 10 in whole number increments. The default value is +5. To customize the **CHROMA (SAT)** level:

1. Navigate to **CHROMA (SAT)**. When **CHROMA (SAT)** blinks, press **Select** to select it. The current setting will blink.
2. Use the Up and Down arrows on the control unit or the left and right Camera Head buttons to scroll through the **CHROMA (SAT)** level settings. Press **Select** to select the desired level.
3. To save the changes and exit the **CAMERA SETUPS SETTINGS** screen, select **SAVE AND EXIT**.

PHASE (HUE)

The **PHASE (HUE)** setting allows the user to set the color hue of the video output during the camera setup. **PHASE (HUE)** settings range from 0 to 10 in whole number increments. Higher settings (longer wavelengths) will contain more red. Mid-range settings (medium wavelengths) will contain more green, and lower range (shortest wavelengths) settings will contain more blue. The default value is +5. To customize the **PHASE (HUE)** level:

1. Navigate to **PHASE (HUE)**. When **PHASE (HUE)** blinks, press **Select** to select it. The current setting will blink.
2. Use the Up and Down arrows on the control unit or the left and right Camera Head buttons to scroll through the **PHASE (HUE)** level settings. Press **Select** to select the desired level.
3. To save the changes and exit the **CAMERA SETUPS SETTINGS** screen, select **SAVE AND EXIT**.

ZOOM

See trace for detail [560: ZOOM - increase or decrease from 1.0X to 2.5X. Each time the button is pressed, the zoom factor increases by approximately 10%, going through the following discrete zoom factors: 1.0X, 1.2X, 1.3X, 1.5X, 1.7X, 1.9X, 2.1X, 2.3X, 2.5X.]

CONFIGURE ICONS

The **CONFIGURE ICONS** setting allows the user to customize which icons appear in the status bar on the OSD. The default setting for all icons is **YES**. To customize the icons display in the status bar:

1. Navigate to the camera setup to be customized and press **Select**. The **CAMERA SETUPS** screen opens.

2. Navigate to **CONFIGURE ICONS** and press **Select** to select it. The **CONFIGURE ICONS** screen opens (Figure 40). Each icon listed can be configured to **YES** (display) or **NO** (do not display).


Figure 40. Configure Icons

3. Navigate to the desired icon. When the icon blinks, press **Select** to select it. The current setting will blink.
4. Use the Up and Down arrows on the control unit or the left and right Camera Head buttons to switch between **YES** and **NO**. Press **Select** to select the desired setting.

Note: The Icons Fluid Double, Robot and Other are not available at this time. Icons displayed can be enabled or disabled individually or as a group.

To **ENABLE** all icons in the status bar, navigate to **ENABLE ALL**. When **ENABLE ALL** blinks, press **Select** to select it. All icons listed on the screen will immediately be reconfigured to **YES** and will appear in the status bar on the OSD.

To **DISABLE** all icons in the status bar, navigate to **DISABLE ALL**. When **DISABLE ALL** blinks, press **Select** to select it. All icons listed on the screen will immediately be reconfigured to **NO** and will not appear in the status bar. This is not recommended.

To exit the **CONFIGURE ICONS** screen, highlight the **EXIT** option and press **Select** to return to the **CAMERA SETUP SETTINGS** screen.

MORE SETTINGS

The **MORE SETTINGS** menu (Figure 41) is an extension of the Camera setup Selection menu.


Figure 41. MORE SETTINGS Menu

1. From the **MAIN MENU**, highlight the camera setups icon and select it to open the Camera setup Selection menu.
2. Use the Up and Down arrows on the control unit or the left and right Camera Head buttons to navigate to the **Customize Camera Setups** icon and press **Select** to open the Setup Selection setup screen.

MORE SETTINGS

3. Navigate to the camera setup to be customized and press **Select** to open that camera's setup customization screen. The **CUSTOMIZE BUTTON SETTINGS** option will blink (Figure 37).

Note: Changes to the camera setup are not saved until **SAVE AND EXIT** is selected.

DISABLE BEEP ON CAPTURE

Disable beep is the post capture beeps that indicate whether it is successful or not. To customize **DISABLE BEEP ON CAPTURE**:

1. Navigate to **DISABLE BEEP ON CAPTURE**. When **DISABLE BEEP ON CAPTURE** blinks, press **Select** to select it. The current setting will blink.
2. Use the Up and Down arrows on the control unit or the left and right Camera Head buttons to switch between YES and NO. Press **Select** to select the desired setting.
3. To save the changes and exit the **CAMERA SETUPS SETTINGS** screen, select **SAVE AND EXIT**.

DISABLE OTHER BEEPS

The **DISABLE OTHER BEEPS** setting allows the user to choose whether or not to have any audible beeps occur. To customize **DISABLE OTHER BEEPS**:

1. Navigate to **DISABLE OTHER BEEPS**. When **DISABLE OTHER BEEPS** blinks, press **Select** to select it. The current setting will blink.
2. Use the Up and Down arrows on the control unit or the left and right Camera Head buttons to switch between YES and NO. Press **Select** to select the desired setting.
3. To save the changes and exit the **CAMERA SETUPS SETTINGS** screen, select **SAVE AND EXIT**.

DISABLE PATIENT INFO

The **DISABLE PATIENT INFO** allows the user to choose whether or not to display the upper left screen containing the patient information.

1. Navigate to **DISABLE PATIENT INFO**. When **DISABLE PATIENT INFO** blinks, press **Select** to select it. The current setting will blink.
2. Use the Up and Down arrows on the control unit or the left and right Camera Head buttons to switch between YES and NO. Press **Select** to select the desired setting.
3. To save the changes and exit the **CAMERA SETUPS SETTINGS** screen, select **SAVE AND EXIT**.

DEFAULT CHROMA MODE

Default **CHROMA Mode** allows the user to choose whether or not to enable or disable Chroma mode to access an alternate color rendition.

Refer to **CHROMA MODE VALUE** section in this document.

1. Navigate to **DEFAULT CHROMA MODE**. When **DEFAULT CHROMA MODE** blinks, press **Select** to select it. The current setting will blink.
2. Use the Up and Down arrows on the control unit or the left and right Camera Head buttons to switch between YES and NO. Press **Select** to select the desired setting.
3. To save the changes select **EXIT**.

660 DUAL CAPTURE

660 Dual Capture – when the 660 is connected and the Dual Capture is selected it will trigger (capture) the tablet to send the image to the 660 when triggered from the camera head or Tablet Application.

1. Navigate to **660 DUAL CAPTURE**. When **660 DUAL CAPTURE** blinks, press **Select** to select it. The current setting will blink.
2. Use the Up and Down arrows on the control unit or the left and right Camera Head buttons to switch between YES and NO. Press **Select** to select the desired setting.
3. To save the changes and exit the **CAMERA SETUPS SETTINGS** screen, select **SAVE AND EXIT**.

660 DUAL VIDEO

660 Dual Video - is the same as Capture and this will mirror the capture from the camera head buttons.

1. Navigate to **660 DUAL VIDEO**. When **660 DUAL VIDEO** blinks, press **Select** to select it. The current setting will blink.
2. Use the Up and Down arrows on the control unit or the left and right Camera Head buttons to switch between YES and NO. Press **Select** to select the desired setting.
3. To save the changes and exit the **CAMERA SETUPS SETTINGS** screen, select **SAVE AND EXIT**.

ELC Settings

ELC (Electric Light Control) Settings allows the user to customize the response of the ELC. The recommended setting is **AUTO**, which is also the default setting. It can also be set to sensitivity levels 1–16.

If flicker occurs while the ELC is set to **AUTO**, reset the ELC to 1 to remove the flicker, then increase the sensitivity level until the speed of the ELC response is sufficient. When sensitivity is set to 1, the response time of the ELC is slow, but more stable. At 16, the response time is quick, but less stable.

To customize **ELC Settings**:

1. Navigate to **ELC Settings** and press **Select** to select it. The **ELC Settings** screen (Figure 42) opens.


Figure 42. ELC Settings

2. Press **Select** to select **SENSITIVITY**. The current setting will blink.
3. Use the Up and Down arrows on the control unit or the left and right Camera Head buttons to scroll through the settings options. Press **Select** to select the desired sensitivity level or **AUTO**.
4. To exit the **ELC SETTINGS** screen, highlight the **EXIT** option and press **Select** to return to the **CAMERA SETUP SETTINGS** screen.

CHROMA MODE VALUE

Works in conjunction with the Default CHROMA Mode, when enabled, the user can adjust the CHROMA MODE VALUE level. Refer to DEFAULT CHROMA MODE section in this document.

1. To INCREASE or DECREASE the value, navigate to the CHROMA MODE VALUE. When the VALUE blinks, press Select to select.
2. Use the Up and Down arrows on the control unit or the left and right Camera Head buttons to navigate to the desired value on the keyboard.
3. To exit the CAMERA SETUPS SETTINGS screen, highlight the EXIT option and press Select.

Save and Exit

To save all changes to the camera setup, navigate to **SAVE AND EXIT**. Press Select to save changes and exit to OSD. To exit without saving changes, custom settings will need to be changed manually to their original settings prior to selecting **SAVE AND EXIT**.

Set the Default Startup Camera setup

Set Default Startup Camera setup allows the user to set the camera setup that is highlighted on the STARTUP screen and so can be selected immediately. To set the default startup camera setup:

1. Navigate to the Camera setup icon and press the Select button to select it.
2. Highlight the Set Default Startup Camera setup icon (Figure 43) and press the Select button to select it. A screen showing the available camera setup opens.
3. Use the Up and Down arrows on the control unit or the left and right Camera Head buttons to scroll to the camera setup to be set as the default camera setup. Press Select to select the camera setup and return to the OSD.


Figure 43. Set Default Startup Camera setup

Cleaning

To Clean the LENS 4K Camera Control Unit

CAUTION: Disconnect the power cord before cleaning the unit or light source.

- Prior to cleaning the equipment, turn off the control unit and disconnect the Camera Head.
- Wipe down the control unit with a clean, dry cloth after every procedure.
- Use a damp cloth or sponge to remove dirt or debris. Avoid getting liquid into the side vents, and never immerse the control unit in any solution.
- Wipe down the control unit with a neutral pH cleaner.
- After cleaning, store the equipment in a cool, dry place out of direct sunlight or excessive heat.

Note: It is important to periodically inspect and clean the fan grill on the rear of the unit and the exhaust holes on the side of the unit. The most effective method is to use a vacuum with a soft brush attached. A soft, damp cloth may also be used to remove any accumulations in these areas.

To Clean and Sterilize the Camera Head

Refer to the Instructions for Use provided with the LENS 4K Camera Head (REF 10601349).

Note: It is important to select the appropriate sterilization method for each type of equipment.

For sterilization of a light guide, refer to the Sterilization section of that light guide's Instructions for Use.

For sterilization of a coupler, refer to the Sterilization section of that coupler's Instructions for Use.

For sterilization of an endoscope, refer to the Sterilization section of that endoscope's Instructions for Use.

Troubleshooting

Troubleshooting

The Smith & Nephew LENS 4K Camera Control Unit is designed to provide the highest quality video image available. However, should problems occur, use the following guide to make troubleshooting easier. If the problem persists, contact an authorized Smith & Nephew representative.

Problem	Try This	Probable cause
Image too bright/excess glare	Decrease the Brightness setting. Use the BRIGHTNESS button on the camera head or access BRIGHTNESS from the MAIN MENU screen.	BRIGHTNESS was set too high
Image too grainy	Decrease the Enhancement setting. Use the ENHANCEMENT button on the camera head or access ENHANCEMENT from the MAIN MENU screen. or	ENHANCEMENT was set too high
	Increase the Brightness setting. Use the BRIGHTNESS button on the camera head or access BRIGHTNESS from the MAIN MENU screen.	BRIGHTNESS was set too low
Image is cloudy or hazy	Adjust the focus on the coupler or the video endoscope or	Camera was out of focus
	Clean both the scope and the coupler. Refer to the Instructions for Use for the respective scope and the coupler in use.	Dirty optics
	If the scope or coupler has fogged, use a suction tube to remove the moisture or dry the equipment using the manufacturer's recommended procedure.	Fogging results when the light in the scope warms or because the air space and trapped moisture has condensed on the glass window.
	Clean any anti-fogging agent off of the lenses. The use of an anti-fogging agent is not recommended	Anti-fogging agent had been applied to the lenses.
	Increase the Enhancement setting. Use the ENHANCEMENT button on the camera head, the ENHANCEMENT adjustment on the iPad, or access ENHANCEMENT from the MAIN MENU screen.	ENHANCEMENT was set too low.
	Ensure that the Moire procedures are not selected.	A fibroscope procedure type was selected for a non-fibroscope setup.
Intermittent picture	Disconnect the camera head from the control unit	If the intermittent behavior is resolved, then refer to the camera head's Instructions for Use to troubleshoot camera head issues. If buildup of contaminants on the proximal end of the camera connector appear, refer to the Cleaning and Sterilization section in the LENS 4K Camera Head for US Ref #10601459, or OUS Ref #106013491.
	Replace the cable on that video output and check the intermittent video output on another monitor.	Bad Cable
	Connect all video outputs to a compatible monitor directly from the control unit with known compatible cables, then switch between outputs at the monitor controls and monitor issue. Do not pass signals through to other devices, for example, the 660 & 660HD-E Image Management System.	
	Restore connections through other devices (for example, the 660 & 660HD-E Image Management System).	If intermittent video behavior persists, contact Smith & Nephew Customer Service.

Problem	Try This	Probable cause
Cannot establish a connection from the iPad to the control unit	1. Go to the VERSIONS screen on the control unit and determine the NETWORK ID. 2. Follow the instructions in the System's Operations/Service Manual to reconnect.	User attempted to connect to the wrong control unit.
	• Check uniqueness of SSID	Network conflict due to non-unique SSID (network name)
	• Change the Wi-Fi channel in the network menu in following the order 1, 4 or 11 • Restart the Camera Control Unit	Contested Wi-Fi channel, or missing Wi-Fi.
	• Ensure that the iPad is within visible range (within 150 ft) of the control unit. If this works but persists within the O.R., contact technical support.	The iPad is out of range of the control unit wireless network
	• Check Wi-Fi antenna is connected on back of system.	Damaged Wi-Fi antenna
Loss of Wi-Fi during a workflow pull	• Ensure the Integration Broker (REF 1061567) is set-up correctly and the configuration is the expected workflow on the Integration Broker server.	Corrupted Workflow list on the Integration Broker server
Cannot capture images/video to the App	• Ensure connection to the control unit Wi-Fi network.	iPad is not connected to the control unit
	• Check that Wi-Fi antenna is connected on to the back of control unit. – If yes, replace the Wi-Fi antenna. – If no, connect a Wi-Fi antenna.	The Wi-Fi antenna was damaged or missing
Cannot capture images/video to the USB storage device	• Refer to the compatible devices in this manual and attempt to reconnect.	Call tech support
	• Inspect for damage to the USB ports on the system (use the alternate port).	Physical damage to port
Camera color is not correct	• Perform a white balance per instructions in IFU.	White balance was not successfully completed
	• Restore factory defaults from system configuration menu (note this will erase all current settings and restore to factory).	Color phase or hue were incorrectly customized in the procedures
	• Unplug and plug in the camera head and perform white balance. If this works and is repeatable, call tech support	
System AC Mains LED is not lit	• Check AC mains switch on back of unit and ensure that the proper mains cable and outlet are used.	Mains switch should be enabled on (1)
	• Replace the fuse.	Fuse was blown
	If this problem persists contact local facilities, then call technical support.	

Troubleshooting

Problem	Try This	Probable cause
No video displayed on the monitor or 660HD or 660HD-E System	Check cables (correct cable, actual connection, integrity of the cables (cracks, breaks, etc.). Note: Ensure the correct functional cables are in use.	Bad cables.
	Bring up the color bar from MAIN MENU and check all video outputs on the monitor. Note: Try using another Camera Head and re-try.	Bad camera head.
	Check the video output format in Advanced Settings menu of the control unit and ensure that settings are compatible with monitor or 660HD or 660HD-E Image Management System. Refer to the Advanced Settings section of this Instructions for Use for more information. Note: The 660HD or 660HD-E System accepts only 1080i HD-SDI output.)	Mismatched video format setting
Arthroscope is too hot	Ensure that the proper light guide diameter is in use for the current operative procedure.	Incorrect light guide
	Ensure that the light guide to light post connection is secure.	Faulty connection
	Turn down the brightness of the control unit and contact Smith & Nephew Customer Service.	Not known

Expected Life

The Smith & Nephew LENS 4K Camera Control Unit is a reusable device. The functional life of these devices depends on various factors, including frequency of use, the method and duration of each use, as well as post-operative methods and handling. Proper maintenance of the equipment is essential to achieve optimal performance over time.

The LENS 4K Camera Control Unit should be inspected before each procedural use per the "Service Life Checklist", if any issues occur, perform "Troubleshooting". If the problem persists, contact Smith & Nephew customer service or an authorized Smith & Nephew representative.

Service

There are no user-serviceable components inside the control unit. Repairs and adjustments are to be performed only by Smith & Nephew authorized service centers.

If service becomes necessary, call an authorized Smith & Nephew Customer Service Representative prior to returning the device and request a Return Authorization (RA) number. A representative can also explain the available Service Replacement and Repair Programs.

Service items should be carefully disinfected, repackaged and returned post-paid to Smith & Nephew. A Smith & Nephew Customer Service Representative can provide additional instructions.

Note: Returned product that is found to have been serviced by an unauthorized third party repair facility and/or sterilized with a sterilization method other than one approved by Smith & Nephew will incur additional costs, regardless of warranty status.

It is not necessary to include accessory items (i.e., power cords, paper trays, footswitches, remote controls, etc.) when returning a device for service.

Do not remove any digital interface boards that may be installed in the unit.

Maintenance

Recommended Annual Performance Checks

Smith & Nephew recommends that dielectric strength, earth leakage current, and protective earth testing be performed annually to assure continued compliance with applicable safety requirements. These tests should be conducted in accordance with specifications IEC 60601-1:2005 + A1:2012 IEC 60601-2-18:2009, IEC 60601-1-2:2014

CAUTION: Electrical safety testing should be performed by a biomedical engineer or other qualified person.

To Replace Fuses

Warnings

- To prevent electric shock, unplug the unit from the electrical outlet before attempting to replace the fuses.
- To avoid fire hazard, use only fuses of the correct type, voltage rating, and current rating.

Note: Fuses that require frequent replacement may indicate an internal problem with the control unit. If fuses continue to blow, contact Smith & Nephew Customer Service for assistance.

To change the fuse carrier:

1. Unplug the power cord from the power outlet and from the control unit.
2. Use a screwdriver to open the fuse compartment door on the AC receptacle (Figure 44) and slide out the two fuse carriers. Refer to the Rear Panel Layout section to locate the AC receptacle on the rear panel.


Figure 44. AC receptacle with fuses

3. Replace fuses. Refer to the System Specifications section for replacement fuse types.
4. Reinsert the fuse carriers using the arrows on the inside of the fuse compartment door as a guide.
5. Snap the fuse compartment door closed.

System Specifications

System Specifications

LENS 4K Camera Control Unit Specifications

**Specifications subject to change without notice.*

Power Requirements	100–240VAC, 50/60 Hz; dual fuse
Equipment Classification	Class 1, type CF equipment.
Unit Dimensions	4.6 in. H x 14 in. W x 15.5 in. L / 11.684 cm. H x 35.56 cm. W x 39.37 cm L.
Unit Weight	13.7 pounds / 6.2kg
Resolution	3840 x 2160p (native) 4K
Shutter	1/60 second through 1/10,000 second shutter speed (60 Hz)
Wi-Fi	2.4Ghz 802.11 b/g/n WPA-2 RF Module
External Storage	8 GB minimum; 32 GB recommended. FAT32 (for drives smaller than 32 GB) or NTFS partitions. Note: Drives must be formatted to support MBR partition tables. The system does not support GPT formatted drives. The following devices have been validated for use with the System: • 8 GB USB flash drives (Transcend, Corsair™, Imation™, Kingston) • 16 GB USB flash drives (ADATA™, Mushkin, Silcon Power™) • 32 GB USB flash drives (Corsair, Patriot™, PNY™, Silcon Power) • 64 GB USB flash drives (Mushkin, & ADATA) • 1 TB external SSD hard drives (Buffalo & Western Digital)
USB Output	Two (2) Type A connectors; provides USB-compatible connection to flash or hard drives.
USB Type	USB 2.0 or greater
Video Input	NTSC or PAL composite input
Video Outputs	
HD Serial Digital Interface (HD-SDI)	One serial digital interfaces (HD-SDI) with 4:2:2 component digital video output on BNC connector
4K Serial Digital Interface (4K-SDI)	Four serial digital interfaces (HD-SDI) with 4:2:2 component digital video output on BNC connector
DisplayPort (DP)	Single stream transport (SST) version 1.2
Video Accessory Controls	Two mini-phono connectors (3.5 mm)
Compliance	IEC 60601-1:2005 + A1:2012, IEC 60601-2-18:2009, IEC 60601-1-2:2014 and EN 60601-1-2:2015
Fuses	F10AH250V; V and I rating: 250V, 10A
Modes of Operation	Continuous

Environmental Conditions

Condition	Transport and/or Storage	Use
Temperature	-20 +140° F (-29 +60° C)	+50 - 104° F (+10 -40° C)
Humidity	15–85% (no condensation allowed)	30–75%
Atmospheric Pressure	70–106 kPa	70–106 kPa
Altitude	3000 meters maximum	3000 meters maximum

Ordering Information

LENS 4K Camera Control Unit - Wi-Fi

Each Smith & Nephew LENS 4K Camera Control Unit includes a 10-foot HD-SDI cable, two mini-phono plug accessory cables, a DisplayPort cable, a Cat 5e Ethernet cable, an Operations/Service Manual, and a power cord. The International System (REF 72205222) includes a European Continental power cord.

REF	Description
72205059	LENS 4K Camera Control Unit, (Wi-Fi)
72205211	LENS 4K Camera Control Unit, U.S. (Wi-Fi)
72205212	LENS 4K Camera Control Unit, German (Wi-Fi)
72205213	LENS 4K Camera Control Unit, Spanish (Wi-Fi)
72205214	LENS 4K Camera Control Unit, French (Wi-Fi)
72205215	LENS 4K Camera Control Unit, Swedish (Wi-Fi)
72205216	LENS 4K Camera Control Unit, Portuguese (Wi-Fi)
72205217	LENS 4K Camera Control Unit, Danish (Wi-Fi)
72205218	LENS 4K Camera Control Unit, Norwegian (Wi-Fi)
72205219	LENS 4K Camera Control Unit, Dutch (Wi-Fi)
72205220	LENS 4K Camera Control Unit, Korean (Wi-Fi)
72205221	LENS 4K Camera Control Unit, Italian (Wi-Fi)
72205222	LENS 4K Camera Control Unit, International (Wi-Fi)
72205223	LENS 4K Camera Control Unit, Chinese (Wi-Fi)
72205224	LENS 4K Camera Control Unit, UK (Wi-Fi)

Compatible Camera Head

REF	Description
72205058	LENS 4K Camera Head
72205344	LENS 4K Camera Head (OUS)
72205343	LENS 4K Camera Head (US Only)

Compatible Devices

REF	Description
72204699	Tablet Application
72204938	LENS Wi-Fi Antenna
72200315	LENS Coupler, 19.5

LENS 4K Camera Control Unit - Non Wi-Fi

Each Smith & Nephew LENS 4K Camera Control Unit non Wi-Fi includes a 10-foot HD-SDI cable, two mini-phono plug accessory cables, a displayport cable, a Cat 5e Ethernet cable, an Operations/Service Manual, and a power cord. The International System (REF 72205224) includes a European Continental power cord.

REF	Description
72205185	LENS 4K Camera Control Unit, (Non-Wi-Fi)
72205239	LENS 4K Camera Control Unit, U.S. (Non-Wi-Fi)
72205240	LENS 4K Camera Control Unit, German (Non-Wi-Fi)
72205241	LENS 4K Camera Control Unit, Spanish (Non-Wi-Fi)
72205242	LENS 4K Camera Control Unit, French (Non-Wi-Fi)
72205243	LENS 4K Camera Control Unit, Italian (Non-Wi-Fi)
72205244	LENS 4K Camera Control Unit, Swedish (Non-Wi-Fi)
72205245	LENS 4K Camera Control Unit, Portuguese (Non-Wi-Fi)
72205246	LENS 4K Camera Control Unit, Danish (Non-Wi-Fi)
72205247	LENS 4K Camera Control Unit, Norwegian (Non-Wi-Fi)
72205248	LENS 4K Camera Control Unit, Dutch (Non-Wi-Fi)
72205249	LENS 4K Camera Control Unit, Korean (Non-Wi-Fi)
72205250	LENS 4K Camera Control Unit, Chinese (Non-Wi-Fi)
72205224	LENS 4K Camera Control Unit, International (Non-Wi-Fi)
72205252	LENS 4K Camera Control Unit, UK (Non-Wi-Fi)

Guidance and Manufacturer's Declaration

Guidance and Manufacturer's Declaration

Electromagnetic Emissions

The Smith & Nephew LENS° 4K Camera Control Unit is intended for use in the electromagnetic environment specified below. The customer or the user of the LENS° 4K Camera Control Unit should assure that it is used in such an environment.

The Smith & Nephew LENS 4K Camera Control Unit was tested to the following:

IEC/EN 60601-1-2: 2014/2015

Specification	IEC 60601-1-2 (2014) / EN 60601-1-2 (2015) MEDICAL ELECTRICAL EQUIPMENT - Part 1-2: General requirements for basic safety and essential performance - Collateral Standard: Electromagnetic disturbances - Requirements and tests	
Test Method	IEC CISPR 11 (2010) / EN 55011 (2010), Industrial, scientific and medical equipment - Radio-frequency disturbance characteristics - Limits and methods of measurement Radiated	
Emissions Test	Compliance	Electromagnetic Environment – Guidance
RF emissions CISPR 11	Group I Environment, Professional Healthcare Facility, and not compatible with a Magnetic Resonance (MR) environment.	The LENS 4K Camera Control Unit is suitable for use in a Professional Healthcare environment where utility mains are privately used by the Healthcare Facility. Note: The EMISSIONS characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A). Not to be used in a residential environment (for which CISPR 11 class B is normally required) as this equipment might not offer adequate protection to radio-frequency communication services.  WARNING: Use of accessories and cables with this equipment, other than those specified or provided by the Smith & Nephew, could result in increased electromagnetic emissions or decreased electromagnetic immunity and result in improper operation.
RF emissions CISPR 11	Class A	
Harmonic emissions EN 61000-3-2	N/A	
Voltage fluctuations/ flicker emissions EN 61000-3-3	N/A	

Guidance and Manufacturer's Declaration – Guidance for Separation Distances

Recommended Separation Distances Between Portable and Mobile RF Communications Equipment and the Smith & Nephew LENS 4K Camera Control Unit

The LENS Integrated System is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the LENS 4K Camera Control Unit can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the LENS Integrated System as recommended below, according to the maximum output power of the communications equipment.

Rated Maximum Output Power of Transmitter (W)	Separation Distance According to Frequency of Transmitter (m)		
	150 KHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

Guidance and Manufacturer's Declaration

Electromagnetic Immunity

The Smith & Nephew LENS 4K Camera Control Unit is intended for use in the electromagnetic environment specified below. The customer or the user of the LENS 4K Camera Control Unit should assure that it is used in such an environment.

Specification	IEC 60601-1-2 (2014) / EN 60601-1-2 (2015) MEDICAL ELECTRICAL EQUIPMENT - Part 1-2: General requirements for basic safety and essential performance - Collateral Standard: Electromagnetic disturbances - Requirements and tests.
Test Method	EN 61000-4-2 (2009) / IEC 61000-4-2 (2008), Electromagnetic compatibility (EMC) – Part 4-2: Testing and measurement techniques – Electrostatic discharge immunity test EN 61000-4-3 (2010) / IEC 61000-4-3 +A1+A2 (2010), Electromagnetic compatibility (EMC) – Part 4-3: Testing and measurement techniques – Radiated, radio-frequency, electromagnetic field immunity test EN 61000-4-4 (2012) / IEC 61000-4-4 (2012), Electromagnetic compatibility (EMC) – Part 4-4: Testing and measurement techniques – Electrical fast transient/burst immunity test EN 61000-4-5 (2006) / IEC 61000-4-5 (2005), Electromagnetic compatibility (EMC) – Part 4-5: Testing and measurement techniques – Surge immunity test EN 61000-4-6 (2014) / IEC 61000-4-6 (2013), Electromagnetic compatibility (EMC) – Part 4-6: Testing and measurement techniques – Immunity to conducted disturbances, induced by radio-frequency fields EN 61000-4-8 (2010) / IEC 61000-4-8 (2009), Electromagnetic compatibility (EMC) – Part 4-8: Testing and measurement techniques – Power frequency magnetic field immunity test EN 61000-4-11 (2004) / IEC 61000-4-11 (2004), Electromagnetic compatibility (EMC) – Part 4-11: Testing and measuring techniques – Voltage dips, short interruptions and voltage variations immunity tests

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment – Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	+/- 8 kV contact +/- 15 kV air	+/- 8 kV contact +/- 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	+/- 2 kV for power supply lines +/- 1 kV for input/output lines	+/- 2 kV for power supply lines +/- 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	+ 2 kV line to P.E. + 1 kV line to line	+ 2 kV line to P.E. + 1 kV line to line	
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% Unom, ½ cycle 40% Unom, 5 cycles 70% Unom, 25 cycles 0% Unom, 250 cycles	0% Unom, ½ cycle 40% Unom, 5 cycles 70% Unom, 25 cycles 0% Unom, 250 cycles	Mains power quality should be that of a typical commercial or hospital environment. If the user of the LENS 4K Camera Control Unit requires continued operation during power mains interruptions, it is recommended that the LENS 4K Camera Control Unit be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Conducted RF IEC 61000-4-6	3 ^μ RMS 150 kHz to 80 MHz	3 ^μ 3 ^μ /m	 WARNING: Portable and mobile RF communication equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the LENS 4K Camera Control Unit, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.
Radiated RF IEC 61000-4-3	3 ^μ /m 80 MHz to 2.7 GHz		

Guidance and Manufacturer's Declaration

		Field strengths from RF transmitters, such as base stations for radio, (cellular/cordless) telephones, land mobile radios, amateur radios, AM and FM radio broadcasts, and TV broadcasts cannot be predicted. The Professional Healthcare Facility must assess the electromagnetic environment due to RF transmitters, an electromagnetic site survey should be considered to ensure compatibility with the immunity levels contain in IEC/EN 60601-1-2. The LENS 4K Camera Control Unit normal operation is for general surgical purposes throughout the duration of orthopedic and arthroscopy applications.  WARNING: Abnormal performance, due to Electromagnetic Disturbances, may be manifested by a failure to properly maintain illumination and visualization of the surgical field. If abnormal performance is observed, it may be corrected by one or more of the following measures: — Reorient or relocate this equipment, the other equipment, or both. — Increase the separation between the pieces of equipment. — Connect the pieces of equipment into different outlets or circuits. — Consult a biomedical engineer.  WARNING: Use of accessories and cables with this equipment, other than those specified or provided by the Smith & Nephew, could result in increased electromagnetic emissions or decreased electromagnetic immunity and result in improper operation.
--	--	---

Note: Unom = Nominal AC input voltage (90, 240VAC/50Hz)

Special Environment

The LENS 4K Camera Control Unit is intended for use in a Special Environment using HF surgical equipment.

 **WARNING - High Frequency (HF) Interference:** This equipment was designed and tested to be used with HF Electrical Surgical Unit (ESU) equipment. If interference occurs with other equipment it may be corrected by one or more of the following measures:

- Reorient or relocate this equipment, the other equipment, or both.
- Increase the separation between the pieces of equipment.
- Connect the pieces of equipment into different outlets or circuits.
- Consult a biomedical engineer.

For Further Information

For Further Information

If further information, including warranty, on this product is needed, contact Smith & Nephew Customer Service at +1 800 343 5717 in the U.S., or an authorized representative.

References

Sterilization Standards

- ANSI/AAMI ST79. *Comprehensive guide to steam sterilization and sterility assurance in healthcare facilities.*
- ANSI/AAMI/ISO 17665-1. *Sterilization of healthcare products – Moist Heat Part 1 Requirements for development, validation and routine control of a sterilization process for medical devices.*
- ISO 14937. *Sterilization of Health Care Products – General Requirements for Characterization of a Sterilization Agent and the Development, Validation and Routine Control of a Sterilization Process for Medical Devices.*
- ISO 17664. *Sterilization of medical devices – information to be provided by the manufacturer for the processing of resterilizable medical devices.*
- AAMI TIR12. *Designing, Testing, and Labeling Reusable Medical Devices for Reprocessing in HealthCare Facilities: A Guide for Medical Device Manufacturers.*

Cleaning Standards

- Guidance for Industry and FDA Staff - *Reprocessing Medical Devices in Health Care Settings: Validation Methods and Labeling.* U.S. FDA March 17, 2015.
- American Society for Testing and Materials (ASTM). *Standard Guide for Selecting Test Soils for Validation of Cleaning Methods for Reusable Medical Devices, F3208-17.*
- American Society for Testing and Materials (ASTM). *Standard Guide for Application of Test Soils for the Validation of Cleaning Methods for Reusable Medical Devices, F3293-18.*
- AAMI TIR12. *Designing, testing, and labeling reusable medical devices for reprocessing in health care facilities: A guide for device manufacturers*
- AAMI TIR 30: *A compendium of processes, materials, test methods, and acceptance criteria for cleaning reusable medical devices.*

All trademarks acknowledged.

Covered by one or more of the following U.S. patent numbers: 5,510,070; 5,563,481; 5,602,449; 5,620,447; 5,630,826; 5,672,945; 5,712,543; 5,745,647; 5,749,885; 5,804,936; 5,833,692; 5,871,493; 5,913,867; 6,090,122; 6,328,752; 7,150,713; 7,510,563; Des. 381,425; other patents pending.

Endoscopy
Smith & Nephew, Inc.
150 Minutemen Road
Andover, MA 01810
USA


+1 978 749 1000 Telephone
+1 978 749 1108 Fax
+1 800 343 5717 U.S. Customer Service

www.smith-nephew.com
lu.smith-nephew.com

©2020 Smith & Nephew, Inc.
*Trademark of Smith & Nephew, registered
U.S. Patent & Trademark Office.
All rights reserved.
05/2020 10601350 Rev. D

Instructions for Use

LENS 4K Camera Head

EN English

CE

Endoscopy

Smith & Nephew, Inc.
150 Minuteman Road
Andover, MA 01810
USA

+1 978 749 1000 Telephone
+1 978 749 1108 Fax
+1 800 343 5717 U.S. Customer Service

www.smith-nephew.com
ifu.smith-nephew.com

©2020 Smith & Nephew, Inc. All rights reserved.
*Trademark of Smith & Nephew, registered U.S. Patent &
Trademark Office. All trademarks acknowledged.

05/2020 10601349 Rev. D

LENS 4K Camera Head

 **smith&nephew**
LENS 4K

Surgical Imaging System

Instructions for Use

Preface

The Smith & Nephew LENS 4K Camera Head is part of a family of products specifically designed for multimodal sterilization including autoclave, STERRAD® 100S Short Cycle, STERRAD® 100 NX All Clear DUO Cycle, STERRAD® 100 NX All Clear Express Cycle, STERRAD® 100 NX All Clear Standard Cycle, STERRAD® 100 NX DUO Cycle, STERRAD® 100 NX Express Cycle, STERRAD® 100 NX Standard Cycle, STERRAD NX® All Clear Standard Cycle, STERRAD NX® Standard Cycle Sterilization Systems. It may also be sterilized with V-PRO® 1 Cycle, V-PRO® 1 Plus Lumen Cycle, V-PRO® maX Lumen Cycle, V-PRO® 1 Plus Non Lumen Cycle, V-PRO® maX Non Lumen Cycle, V-PRO® maX Flexible Cycle Sterilization Systems.

Prior to using the equipment, read these instructions thoroughly. Pay particular attention to the Warnings and Precautions. These instructions will familiarize the user with the camera head and the steps required for proper setup and maintenance.

Device Description

The camera head is a state-of-the-art, 4K video camera head designed for use in arthroscopic and endoscopic surgical procedures. It is compatible with the LENS 4K Camera Control Unit REF 72205059 – LENS 4K Camera Control Unit, Wi-Fi; REF 72205185 – LENS 4K Camera Control Unit, non-Wi-Fi).

Compatible Endoscopes and Couplers (Must have an industry standard 1.00 – 32 UN 2B C-mount thread.)

Arthroscopes:	1.9, 2.7, and 4.0 mm sizes
Laparoscopes:	5.0, 5.5, and 10 mm size
Hysteroscopes:	4.0, 5.0, 5c, and 8.0 mm
Thoracoscopes:	2.8 mm

Intended Users

Smith & Nephew LENS 4K Camera Head is intended for use by healthcare professionals in accordance with these instructions for use. The use environment is a professional healthcare facility.

Intended Use

The LENS 4K Camera Control Unit and Camera Head are intended to provide illumination, visualization and capture of still and motion pictures of surgical sites.

Indications for Use

The LENS 4K Camera Head and Camera Control Unit are used in diagnostic and operative endoscopic procedures such as arthroscopy to provide illumination, visualization and capture of still and motion pictures of surgical sites within articular cavities, body cavities, hollow organs and canals. Additionally, LENS 4K Camera Head and Camera Control Unit are indicated for use in endoscopic surgical procedures in the thoracic cavity when used with an appropriately indicated thoracoscope.

Contraindications

None known.

Target Patient Group

The targeted patient population for the LENS 4K is patients requiring endoscopic surgical procedures as determined necessary by the prescribing healthcare professionals.

Warnings

- The camera head is shipped non-sterile. It must be cleaned and sterilized before the first use. It must be cleaned and sterilized before every subsequent use.
- It is the health care professional's responsibility to be familiar with the appropriate surgical techniques prior to use of this device.
- Read these instructions completely prior to use.
- Read the Operations/Service Manual for the applicable camera control unit and the Instructions for Use for other components used in the procedure.
- To prevent patient burns, do not place the camera head on the patient when it is not in use. The camera head contains heat-generating components that raise its surface temperature higher than that of the surrounding environment. Prolonged contact with the patient may cause burns.

- **DANGER:** Risk of explosion if used in the presence of flammable anesthetics.
- Do not use this device in oxygen rich environments (>25%).
- To prevent electric shock, do not disassemble the camera head. There are no user-serviceable parts. Dismantling the equipment will void the warranty.
- The camera head is autoclavable and compatible with only one chemical sterilant throughout its lifetime. Once a sterilant has been used, continue to use that same sterilant for the lifetime of that specific camera head. Do not apply multiple chemical sterilants during the lifetime of any single camera head. Use only chemical sterilants from the list of acceptable chemical sterilants.
- Do not modify this equipment. Attempts to modify this equipment will void the warranty.
- The LENS 4K Camera head is not compatible with any Liquid chemical sterilant processes.

Precautions

Rx U.S. Federal law restricts this device to sale by or on the order only of a physician.

- Prior to use, inspect the device to ensure it is not damaged. Do not use a damaged device.
- This product may be sterilized with chemical sterilants STERRAD® 100S Short Cycle, STERRAD® 100 NX All Clear DUO Cycle, STERRAD® 100 NX All Clear Express Cycle, STERRAD® 100 NX All Clear Standard Cycle, STERRAD® 100 NX DUO Cycle, STERRAD® 100 NX Express Cycle, STERRAD® 100 NX Standard Cycle, STERRAD NX® All Clear Standard Cycle, STERRAD NX® Standard Cycle Sterilization Systems. It may also be sterilized with V-PRO® 1 Cycle, V-PRO® 1 Plus Lumen Cycle, V-PRO® maX Lumen Cycle, V-PRO® 1 Plus Non Lumen Cycle, V-PRO® maX Non Lumen Cycle, V-PRO® maX Flexible Cycle Sterilization Systems. The product may be sterilized with any of the above specified chemical sterilization methods. Failure to comply with this caution may result in damage to the product and will void the warranty.
- The camera head cannot be sterilized with Steris 1E™ or any other wet sterilant.
- To avoid fogging during surgery, ensure that the VideoArthroscopes™ or coupler and the camera head are entirely free of moisture prior to connecting them.
- Use ENZOL™ to clean the camera head. If ENZOL is not available in your facility a neutral pH cleaner can be used. Non-neutral pH cleaners may cause residue buildup, which may cause interference or damage the optical seals. Such damage will result in internal fogging and image discoloration.
- **US Customers:** Users should only use sterilizers and accessories (such as sterilization wraps, sterilization pouches, chemical indicators, biological indicators, and sterilization containers) that have been cleared by the US FDA for the selected sterilization cycle specifications (time and temperature).
- Handle all endoscopic instruments with care. If the camera head is dropped or damaged in any way, it must be returned immediately for service.
- Electrical safety testing should be performed by a biomedical engineer or other qualified service person.
- Do not point the camera head toward direct sunlight, extremely bright lights, or laser beams for prolonged periods as damage may occur.
- Do not plug the camera cable connector into the camera control unit if wet. Moisture on the camera cable connector, including the gold fingers, on the card edge of the camera head cable will damage the circuitry and void the warranty. Ensure that the camera cable connector is completely dry prior to plugging the camera head into the camera control unit.
- Do not sterilize damaged equipment, as further damage may result. Clean and disinfect the camera head prior to returning it to Smith & Nephew.
- Ensure that there are no deposits on the optical surfaces prior to sterilization.
- Allow the camera head to cool to room temperature after it has been autoclaved.
- Extreme caution must be taken to dry out mating portions of both the scope and camera head; otherwise trapped moisture may condense on glass surfaces, blurring the image.
- Carry and store the cable in loose coils to prevent kinks and stress on the internal wires. Do not carry the camera head by the cable.
- Do not pinch or damage the camera head cable (i.e., rolling heavy equipment over it, securing with clamps, pinching in doors, etc.).
- Do not modify this equipment. Attempts to modify this equipment will void the warranty.
- To minimize any reciprocal interference from or to the System, use a dedicated power outlet during any RF treatment or investigation.
- **Environmental Protection:** This equipment contains electronic printed circuit assemblies. At the end of the useful life of the equipment it should be disposed of in accordance with any applicable national or institutional policy relating to obsolete electronic equipment.
- For information regarding electrical emissions and immunity testing, guidance separation, and precautions, refer to the camera control unit's Operations/Service Manual.
- After use, this device may be a potential biohazard and should be handled in accordance with accepted medical practice and applicable local and national requirements.
- This product contains Nickel (Ni) which may cause allergic reactions.

Unpacking and General Inspection

Carefully unpack and inspect all of the components shipped with the camera head. You should have received the following:

- 1 ea. LENS 4K Camera Head and Cable (REF 72205058)
- 1 ea. LENS 4K Camera Head Instructions for Use (REF 10601349)

If any parts are missing or damaged, contact an authorized Smith & Nephew representative. Save the carton and packing materials in the event that a component needs to be returned or repaired.

Button Function

Table 1, Live Video Functions and Table 2, OSD (On screen Display) Menu Functions summarize the default button functions for the camera head for both the Live Video and OSD Menu modes. Table 3 provides Optional Programmable Button Functions. Refer to the Operations/Service Manual for the LENS 4K camera control unit (REF 10601350) for more information regarding default and customized button settings and a complete list of the associated functions.

Camera Head Control	Short Press (<1.0 seconds)	Long Press (>1.0 seconds)
Left Button	Frame Capture (Single from Primary)	Light Mode
Middle Button	Brightness (Up/Down)	Menu Access
Right Button	Zoom	White Balance

Table 1. Live Video Functions

Camera Head Control	Default Button Function
Left Button	Move the cursor UP to highlight a menu choice
Middle Button	SELECT the highlighted menu choice
Right Button	Move the cursor DOWN to highlight a menu choice

Table 2. OSD Menu Functions

Function	Effect
BRIGHTNESS	Control the luminance level of the video output
ENHANCEMENT	Enhance the sharpness of the displayed video
ZOOM	Digital zoom of field of view
ALT PAUSE	Pause the video while recording to a device connected to a rear panel port
PAUSE VIDEO	Pause the video from camera head
ALT VIDEO	Start or stop video capture from the secondary input source and save to the Tablet Application or USB
ALT FRAME	Capture a still picture from the secondary input source and save to the Tablet Application
FRAME CAPTURE	Capture a still picture from the endoscopic field of view and save to the Tablet Application or USB
LIGHT MODE	Activates and deactivates the illumination device on the control unit
VIDEO CAPTURE	Start or stop video capture from the endoscopic field of view and save to the Tablet Application or USB
LAST IMAGE	Displays the last known frame captures to the Tablet Application for review
MENU ACCESS	Access the graphical menu system of the control unit [MAIN MENU]
WHITE BALANCE	Correct for ambient color temperature
Peripheral Port 1 (Top port)	Triggers the accessory connected to the top port
Peripheral Port 2 (Bottom port)	Triggers the accessory connected to the bottom port

Table 3. Optional Programmable Button Functions*

* **Note:** Software updates to either the LENS 4K Camera Head or Camera Control Unit, may reset the button-mappable functions to the factory settings.

Adverse Event Reporting

If, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and to your national authority.

complaints@smith-nephew.com

Instructions for Use

The camera head buttons perform a different set of functions depending upon whether the On Screen Display (OSD) Menu or Live Video is displayed on screen (Tables 1 and 2). When the OSD menu is displayed, use the camera head buttons to navigate through the menus.

When Live Video is displayed, each of the camera head buttons is assigned to default or customized functions for short (less than 1 second) or long (>1 second) presses. The default functions can be changed using the CUSTOMIZE PROCEDURE SETTINGS menu. Each registered press of a button is followed by a beep from the LENS 4K camera head. The beep is the same for both short and long presses.

Note: The camera head is not specific to any single LENS 4K camera control unit. It is interchangeable among multiple System control units.

Refer to the Operations/Service Manual for the applicable LENS 4K control unit for more information regarding the setup of the control unit and the LENS 4K camera head.


Figure 1. LENS 4K Camera Head

Connect the Camera Head

CAUTION: Do not plug the camera cable connector into the camera control unit if wet. Moisture on the camera cable connector, including the gold fingers, on the card edge of the camera head cable will damage the circuitry and void the warranty. Ensure that the camera cable connector is completely dry prior to plugging the camera head into the camera control unit.

To connect the camera cable connector to the control unit, plug the card edge connector into the appropriate receptacle on the front of the camera control unit and push in firmly (Figure 2).


Figure 2. Connect the camera head to the control unit

Endocoupler focal length determines the image size on the monitor. For optimal image size use the LENS Coupler, 19.5 (72200315).

Service Life Checklist

Before each case begins, check the following items to ensure that they function properly and the system is ready to use:

- ☐ Confirm that the video works correctly and that visualization is correct.
- ☐ Check for any damage to the camera head and receptacle (cracks, corrosion, buildup or damage to the card edge receptacle).
- ☐ Check that the labeling information on the camera head connector is still readable.
- ☐ Ensure there are no deposits, scratches, or fogging/humidity on the front window of the camera head.
- ☐ Check for any damage to the cable (cuts or ruptured).
- ☐ Check that all camera head buttons are functioning properly by pressing each button.

White Balance the Camera

When a procedure is selected from the LENS 4K System STARTUP menu, the system automatically starts the White Balance procedure to ensure proper color rendition. The scope and the fiber optic light guide should already be connected in order to perform the White Balance procedure.


Figure 3. White Balance the camera

1. Point the end of the scope at a white object (a 4" x 4" gauze pad or flat white piece of paper) and focus. Fill as much of the screen as possible with the white object, but do not hold the scope close enough to touch the white object.
2. Select the desired procedure from the STARTUP menu. The monitor will display the message **WHITE BALANCE IN PROGRESS**.
3. Hold the gauze pad or paper in front of the scope until the monitor displays the message **WHITE BALANCE COMPLETE** and a long tone sounds.

If an error occurs during White Balance, the message **WHITE BALANCE INCOMPLETE** will appear on the display and a chirping tone sounds. To White Balance without restarting system, select **WHITE BALANCE** from the MAIN MENU or press the **WHITE BALANCE** button on the front panel of the control unit.

Cleaning and Sterilization

Notes:

- Disinfection is not a substitution for sterilization. Always inspect for damage and debris prior to sterilization.
- The camera head requires sterilization because it is used in the sterile field. No protection is needed for the camera head connector during cleaning and sterilization. The camera head connector must be dry prior to plugging into the camera control unit.

CAUTION: US Customers: Users should only use sterilizers and accessories (such as sterilization wraps, sterilization pouches, chemical indicators, biological indicators, and sterilization containers) that have been cleared by the US FDA for the selected sterilization cycle specifications (time and temperature).

WARNING: The camera head is autoclavable and compatible with only one chemical sterilant throughout its lifetime. Use only chemical sterilants from the list of acceptable chemical sterilants. Once a sterilant has been used on the camera head, continue to use that same sterilant for the lifetime of that specific camera head. The use of multiple chemical sterilants during the lifetime of any one camera head will reduce the service lifetime of that camera head.

CAUTION: Prior to use, inspect the device to ensure it is not damaged. Do not use a damaged device.

CAUTION: Use ENZOL™ to clean the camera head. If ENZOL is not available in your facility a neutral pH cleaner can be used. Non-neutral pH cleaners may cause residue buildup, which may cause interference or damage the optical seals. Such damage will result in internal fogging and image discoloration.

CAUTION: To avoid fogging during surgery, ensure that the VideoArthroscope or coupler and the camera head are entirely free of moisture prior to connecting them.

CAUTION: Do not plug the camera cable connector into the camera control unit if wet. Moisture on the camera cable connector, including the gold fingers, on the card edge of the camera head cable will damage the circuitry and void the warranty. Ensure that the camera cable connector is completely dry prior to plugging the camera head into the camera control unit.

CAUTION: Do not sterilize damaged equipment, as further damage may result. Clean and disinfect the camera head prior to returning it to Smith & Nephew.

Cleaning and sterilization requires four steps in the following order:

- Cleaning
- Equipment Inspection
- Equipment preparation
- Sterilization

Preparations at the point of use prior to processing

Keep instruments moist after use to prevent soil from drying on them.

Follow Universal Precautions for handling and transporting contaminated instruments to the designated cleaning area. Contaminated instruments should be transported to the area for cleaning in a way that avoids contamination of personnel and hospital.

Prior to cleaning, remove visible debris from the surfaces of the device. Dried-on soil is difficult and sometimes impossible to remove with automatic washing.

Clean instruments as soon as possible after use to prevent blood from drying on the devices.

Approved Cleaners/Pre cleaners

Use ENZOL a neutral pH enzymatic cleaner to clean or pre-clean the camera head. If ENZOL is not available in your facility a neutral pH cleaner can be used.

Water Quality

Warm tap water is recommended for cleaning the camera head assembly. Distilled water is recommended for rinsing the camera head assembly.

Cleaning

Immediately after a procedure, manually clean the camera head to make it safe for further handling. The following steps constitute a general guideline for routine cleaning of the camera head and other equipment.

1. Turn the camera control unit off.
2. Disconnect the light guide from the endoscope and light source.
3. Unplug the camera head cable and unplug any remote controls that were used in the sterile field. Detach the endoscope or endocoupler from the camera head.
4. Contain contaminated instruments for transportation to the cleaning area and preparation for manual cleaning.
5. Do not use germicides or any other agents containing benzene or phenol derivatives. This will damage the camera head cable and will void the warranty.

● Clean the camera head

Manual cleaning

1. Rinse the camera head assembly in cold tap water for a minimum of 1 minute.
2. Soak in ENZOL™ a neutral pH solution for a minimum of 2.5 minutes to prevent coagulation and help remove protein materials and other resistant debris. If ENZOL is not available in your facility a neutral pH cleaner can be used.
3. Clean manually using soft cloths, soft bristled brushes or sponges for a minimum of 1 minute.

- To prevent the scratching of lenses and electrical components, limit the number of items in the sink at any one time.
- Never use scouring pads or abrasive powders.
- Pay particular attention to the button area. To avoid scratching lenses, wipe them in one direction; do not scrub.

Note: It is the responsibility of the hospital to ensure adequate cleaning of the brushes between each use.

4. Rinse in purified water for a minimum of 1 minute.
5. Use isopropyl alcohol on the lens for a minimum of 1 minute. Do not allow the exposed lenses and glass surfaces to air dry. Use a cotton applicator dipped in 70% isopropyl alcohol to remove stubborn spots, streaks, and residual debris.
6. Dry thoroughly with a clean applicator. Do not allow to air dry.
7. Repeat the cleaning process until the device is visually clean.

Automated cleaning

1. Place the camera directly on the washer rack being careful to arrange the cables so they have full exposure to the washing/rinsing action of the washer.
2. Run the automatic washer cycle. The minimum automatic washer parameters are:
 - Five (5) minute cold prewash (<43 °C)
 - Five (5) minute enzyme wash (neutral pH cleaner)
 - Two (2) minute rinse hot water
 - One (1) minute final rinse with purified water (90° C)

To clean the coupler

For information about how to clean and sterilize the coupler, refer to the Instructions for Use provided with that product.

To clean the card edge connector:

1. Use a soft-bristled brush to manually remove debris from the card edge connector and cable. Never use abrasive powders.
Note: It is the responsibility of the hospital to ensure adequate cleaning of the brushes between each use.
2. If residue is visible on the metal housing of the card edge connector, use a mildly abrasive plastic scrubbing pad to gently scrub both sides of the metal housing. Ensure that all visible debris is removed.

DO NOT scrub the gold fingers on the end of the card edge connector as this will damage the electronic connections and cause the camera head to fail. Figure 4 indicates which areas to scrub and which to avoid.


Figure 4. Card edge connector

3. Use distilled water to thoroughly rinse off any residue. Do not use saline as salt can increase pitting or corrosion of surfaces.
4. Use a soft, lint-free cloth to thoroughly dry the connector.
5. Wipe the camera head cable connector with isopropyl alcohol and air dry or dry thoroughly with a sterile, lint-free cloth. Ensure that the connector is completely dry. If moisture is present in the connector, use compressed air to blow out any fluid within the connector.
6. Visually inspect the device to ensure all visible tissue and debris have been removed. If not, repeat steps 1–5.
Note: Do not lubricate the camera head or cables.
7. Repeat the cleaning process until the device is visually clean.

Equipment Inspection

Prior to each use, inspect the camera head to ensure the following:

Camera head cable jacket: There are no burn holes, cracks, tears, swelling, etc.

Camera head cable connector: There is no corrosion, damaged card edge, etc.

Glass window: There are no scratches and no deposits.

If any of the listed conditions are present, contact Smith & Nephew as directed in the Warranty section of this Instructions for Use. The use of damaged equipment could compound the damage, harm interfacing equipment, or void the warranty.

Equipment Preparation

After every use, ensure that the camera head has been thoroughly cleaned, inspected, and sterilized according to the directions in this document.

Note: The camera head can be autoclaved individually.

Sterilization

READ ALL WARNINGS AND CAUTIONS BEFORE STERILIZING.

WARNING: The camera head is autoclavable and compatible with only one chemical sterilant throughout its lifetime. Once a sterilant has been used, continue to use that same sterilant for the lifetime of that specific camera head. Do not apply multiple chemical sterilants during the lifetime of any single camera head. Use only chemical sterilants from the list of acceptable chemical sterilants.

CAUTION: This product may be sterilized with chemical sterilants STERRAD® 100S Short Cycle, STERRAD® 100 NX All Clear DUO Cycle, STERRAD® 100 NX All Clear Express Cycle, STERRAD® 100 NX All Clear Standard Cycle, STERRAD® 100 NX DUO Cycle, STERRAD® 100 NX Express Cycle, STERRAD® 100 NX Standard Cycle, STERRAD NX® All Clear Standard Cycle, STERRAD NX® Standard Cycle Sterilization Systems. It may also be sterilized with V-PRO® 1 Cycle, V-PRO® 1 Plus Lumen Cycle, V-PRO® maX Lumen Cycle, V-PRO® 1 Plus Non Lumen Cycle, V-PRO® maX Non Lumen Cycle, V-PRO® maX Flexible Cycle Sterilization Systems. The product may be sterilized with any of the above-specified chemical sterilization methods. Failure to comply with this caution may result in damage to the product and will void the warranty.

CAUTION: Ensure that there are no deposits on the optical surfaces prior to sterilization.

CAUTION: Allow the camera head to cool at room temperature after it has been autoclaved.

CAUTION: Extreme caution must be taken to dry out mating portions of both the scope and camera head; otherwise trapped moisture may condense on glass surfaces, blurring the image.

Notes:

- Disinfection is not a substitute for sterilization.
- Always inspect for damage and debris prior to sterilization.

For sterilization instructions for light guides, couplers, and endoscopes, refer to the Sterilization section of the respective product's Instructions for Use.

Sterilization Methods

WARNING: The camera head is autoclavable and compatible with only one chemical sterilant throughout its lifetime. Once a sterilant has been used, continue to use that same sterilant for the lifetime of that specific camera head. Do not apply multiple chemical sterilants during the lifetime of any single camera head. Use only chemical sterilants from the list of acceptable chemical sterilants.

Use one of the following methods to sterilize the camera head:

WARNING: The LENS 4K Camera head is not compatible with any Liquid Chemical sterilant processes.

Autoclave (Steam Sterilization Parameters)

Method	Cycle	Exposure Temperature	Exposure Time (minutes)	Dry Time (minutes)
Steam (wrapped)	Prevacuum	270 °F (132 °C)	4	30
Steam (wrapped)	Gravity	270 °F (132 °C)	15	30

CAUTION: Allow the camera head to cool to room temperature after it has been autoclaved.

Non US Customer

Autoclave (Steam Sterilization Parameters)

Method	Cycle	Exposure Temperature	Exposure Time (minutes)	Dry Time (minutes)
UK Cycle (wrapped)	NA	134 °C	3	30
World Health Organization (WHO)	NA	134 °C	18	30

Note: The UK cycle is not an FDA cleared sterilization cycle and only for use Internationally.

Note: The WHO cycle is not an FDA cleared sterilization cycle and only for use Internationally. This sterilization cycle has not been validated for prion inactivation.

CAUTION: Allow the camera head to cool to room temperature after it has been autoclaved.

Chemical Sterilization

Sterilize the device using any one of the sterilization systems. The device load is the same for all sterilization systems listed.

Device Load: The LENS 4K camera should be wrapped in a FDA cleared sterilization wrap individually.

List of Sterilization System:

- STERRAD® 100S Short Cycle
- STERRAD® 100 NX All Clear DUO Cycle
- STERRAD® 100 NX All Clear Express Cycle
- STERRAD® 100 NX All Clear Standard Cycle
- STERRAD® 100 NX DUO Cycle
- STERRAD® 100 NX Express Cycle
- STERRAD® 100 NX Standard Cycle
- STERRAD NX® All Clear Standard Cycle
- STERRAD NX® Standard Cycle
- V-PRO® 1 Cycle
- V-PRO® 1 Plus Lumen Cycle
- V-PRO® maX Lumen Cycle
- V-PRO® 1 Plus Non Lumen Cycle
- V-PRO® maX Non Lumen Cycle
- V-PRO® maX Flexible Cycle

Follow all recommended directions and warnings included with the sterilization system.

CAUTION: The camera head cannot be sterilized with Steris System IETM or any other wet sterilant.

Notes:

- The camera head can be used immediately after chemical sterilization.
- These parameters have been validated to ensure sterility.
- The camera head is validated for sterilization assurance using the sterilant manufacturer's instructions.

Storage

To prevent damage, ensure that camera heads remain immobile while stored.

Expected Life

The Smith & Nephew LENS 4K Camera Head is a reusable instrument. The functional life of these devices depends on various factors, including frequency of use, the method and duration of each use, as well as post-operative methods and handling. Proper maintenance of the equipment is essential to achieve optimal performance over time.

The LENS 4K Camera Head should be inspected before each procedural use per the "Service Life Checklist". If any issues occur, contact Smith & Nephew customer service or an authorized Smith & Nephew representative.

Service

WARNING: Do not modify this equipment. Attempts to modify this equipment will void the warranty.

There are no user-serviceable components in the camera head. Repairs and adjustments are to be performed only by Smith & Nephew authorized service centers.

Camera Head Specifications*

Equipment Classification	Protection against electrical shock Class 1, type CF equipment
Video Sensor Technology	4K 3-Chip MOS
Video Data Transmission	High-speed, low-voltage differential serial (LVDS), digital transmission via 12-foot camera cable
Video Output Format	Full UHD 2160p, 100% compatible with SMPTE ST 425-5 when used with the LENS 4K Camera Control Unit
Video Field Rate	60 Hz
Transport Shock and Vibration Compliance	ISTA 3A
Camera Head Dimensions	1.7" W x 2.0" H x 3.5" L 43.18 mm W x 50.8 mm H x 88.9 mm L
Camera Head Weight	1.5 lbs (0.7 kg)
Camera Head Cable	12 ft (366 cm)
Lens Mount	1.00 - 32 UN 2B C-mount thread
Compatible Control Units	Smith & Nephew LENS 4K Camera Control Unit
Compliance	Complies with IEC 60601-1:2005 + A1:2012, IEC 60601-2-18:2009, IEC 60601-1-2:2014
Modes of Operation	Continuous

*Specifications subject to change without notice.

If service becomes necessary, call an authorized Smith & Nephew Customer Service Representative prior to returning the device and request a Return Authorization (RA) number. A representative can also explain the available Service Replacement and Repair programs.

Service items should be carefully cleaned and disinfected, repackaged, and returned postpaid to Smith & Nephew. Smith & Nephew Customer Service or an authorized representative can provide additional instructions.

Note: Product that has been returned that is found to have been serviced by an unauthorized third party repair facility and/or sterilized with a sterilization method other than the ones approved for it by Smith & Nephew will incur additional costs, regardless of warranty status.

Maintenance

Recommended Annual Performance Checks

Smith & Nephew recommends that Dielectric Strength Testing be performed annually to assure continued compliance with applicable safety requirements. These tests should be conducted in accordance with specifications IEC 60601-1:2005 + A1:2012, IEC 60601-2-18:2009, IEC 60601-1-2:2014.

CAUTION: Electrical safety testing should be performed by a biomedical engineer or other qualified service person.

CAUTION: For information regarding electrical emissions and immunity testing, guidance separation, and precautions, refer to the camera control unit's Operations/Service Manual.

Note: Service life is determined by adherence to general maintenance procedures.

Environmental Conditions

Transport and/or Storage Use		Use
Ambient Temperature	-20 – +140 °F -29 – +60 °C	+50 –77 °F +10 –25 °C
Relative Humidity	15–85%	30–75%
Atmospheric Pressure	70–106 kPa	70–106 kPa
Operating Altitude	3000 meters maximum	3000 meters maximum

Ordering Information

REF 72205344 Smith & Nephew LENS 4K Camera Head
 REF 10601349 Smith & Nephew LENS 4K Camera Head
 Instructions for Use (all countries)

For Further Information

If further information including warranty on this product is needed, contact Smith & Nephew Customer Service at +1 800 343 5717 in the U.S., or an authorized Smith & Nephew representative.

References




Sterilization Standards

- ANSI/AAMI ST79. *Comprehensive guide to steam sterilization and sterility assurance in healthcare facilities.*
- ANSI/AAMI/ISO 17665-1. *Sterilization of healthcare products – Moist Heat Part 1 Requirements for development, validation and routine control of a sterilization process for medical devices.*
- ISO 14937. *Sterilization of Health Care Products – General Requirements for Characterization of a Sterilization Agent and the Development, Validation and Routine Control of a Sterilization Process for Medical Devices.*
- ISO 17664. *Sterilization of medical devices – information to be provided by the manufacturer for the processing of resterilizable medical devices.*
- AAMI TIR12. *Designing, Testing, and Labeling Reusable Medical Devices for Reprocessing in HealthCare Facilities: A Guide for Medical Device Manufacturers.*

Cleaning Standards

- Guidance for Industry and FDA Staff - *Reprocessing Medical Devices in Health Care Settings: Validation Methods and Labeling.* U.S. FDA March 17, 2015.
- American Society for Testing and Materials (ASTM). *Standard Guide for Selecting Test Soils for Validation of Cleaning Methods for Reusable Medical Devices, F3208-17.*
- American Society for Testing and Materials (ASTM). *Standard Guide for Application of Test Soils for the Validation of Cleaning Methods for Reusable Medical Devices, F3293-18.*
- AAMI TIR12. *Designing, testing, and labeling reusable medical devices for reprocessing in health care facilities: A guide for device manufacturers*
- AAMI TIR 30: *A compendium of processes, materials, test methods, and acceptance criteria for cleaning reusable medical devices.*

Symbols Glossary

	Follow instructions for use		Catalog number		MR Unsafe
	Non-sterile		Serial number		Atmospheric Pressure
	CE mark		EU: Not for general waste		Medical Device
	Caution		Type CF Equipment		Consult instructions for use
	General Warning Sign		Temperature Limitation		Authorized representative in the European community
	Manufacturer		Humidity Limitation		Prescription Use Only
	Date of Manufacture		Do not plug in wet		Unique Device Identifier
	Batch Code		Keep Dry		

Endoscopy

Smith & Nephew, Inc.
150 Minuteman Road
Andover, MA 01810
USA

+1 978 749 1000 Telephone
+1 978 749 1108 Fax
+1 800 343 5717 U.S. Customer Service

www.smith-nephew.com
ifu.smith-nephew.com

©2020 Smith & Nephew, Inc.
All rights reserved.

*Trademark of Smith & Nephew, registered
U.S. Patent & Trademark Office.

All trademarks acknowledged.

05/2020 10601349 Rev. D

Instructions for Use

Couplers, Beamsplitters, and Video Adaptors

EN English

CE

Endoscopy
Smith & Nephew, Inc.
150 Minuteman Road
Andover, MA 01810
USA

+1 978 749 1000 Telephone
+1 978 749 1108 Fax
+1 800 343 5717 U.S. Customer Service

www.smith-nephew.com
ifu.smith-nephew.com

*Trademarks of Smith & Nephew.
©2002, 2020 Smith & Nephew, Inc. All rights reserved.

05/2020 1060329 Rev. M

Couplers, Beamsplitters, and Video Adaptors

Instructions for Use

Device Description

Camera couplers and beamsplitters are used to focus the image exiting the eyepiece of an endoscope onto the surface of a CCD camera head. Camera couplers and beamsplitters can each be used to provide a rotating connection between a direct view endoscope and a camera head. The endoscope attaches to the distal end and the camera head attaches to the proximal end.

Most camera couplers and beamsplitters have a focus knob to move the lens axially to produce the sharpest possible image. Beamsplitters differ from camera couplers in that they have an internal optical device that splits the single image coming from the endoscope into two complete images. This feature allows the image from the endoscope to be viewed straight through by the user's eye via an eyepiece while also providing a view ninety degrees from the axis of the endoscope onto the CCD camera.

For the purposes of these Instructions for Use Couplers, Beamsplitters, and Video Adaptors are referred to collectively as couplers.

Note: Multimode refers to sterilization methods. Multimode couplers may be autoclaved, STERRAD® 100S, STERRAD® NX, STERRAD® 100NX or STERIS® V-PRO 1 sterilized. These couplers are marked "Autoclavable" or "Multimode". Refer to the "Sterilization" section for parameters. Couplers not marked "Autoclavable" or "Multimode" cannot be autoclaved.

Note: Non-autoclave couplers may be STERRAD 100S, STERRAD NX, STERRAD® 100NX or STERIS® V-PRO 1 sterilized. They must not be autoclaved. Refer to the "Sterilization" section for parameters.

Intended Users

The Smith & Nephew couplers are intended for use by healthcare professionals in accordance with these instructions for use. The use environment is a professional healthcare facility.

Intended Use

The coupler is intended to be used as a video accessory to a camera head to allow endoscopic visualization in surgery procedures.

Indications for Use

Smith & Nephew couplers are intended to be used as video accessories to camera heads. They are indicated for use in endoscopic surgical procedures to allow visualization of the articular cavities, body cavities, hollow organs, and canals when used in conjunction with an appropriately indicated endoscope. Additionally, these couplers are indicated for use in endoscopic surgical procedures in the thoracic cavity when used with an appropriately indicated thoracoscope.

Contraindications

None known.

Target Patient Group

The targeted patient population for the Couplers is patients requiring endoscopic surgical procedures as determined to be necessary by the prescribing healthcare professionals.

Warnings

- These products are shipped non-sterile. They must be sterilized before the first use. They must be cleaned and sterilized before every subsequent use.
- It is the health care professional's responsibility to be familiar with the appropriate surgical techniques prior to use of this device.
- Read these instructions completely prior to use.
- Ensure that only medical equipment that complies with IEC 60601-1 and its relative standards is connected to, or used in conjunction with these devices.

- For steam autoclave method of sterilization applicable to couplers marked as "Autoclavable" or "Multimode" all individual parts (camera head, scope and coupler) **MUST** be detached prior to decontamination and sterilization. Failure to detach them from each other prior to decontamination and sterilization will compromise sterility and constitute a patient hazard.
- If coupler and camera head are sterilized as an assembled unit (applicable only to non-autoclave sterilization methods) be aware that the juncture between the coupler and the camera head is not sterile. In case failure of either component during use both the camera head and the coupler must be replaced.

Precautions

R U.S. Federal law restricts this device to sale by or on the order of a physician.

- Prior to each use, inspect the device to ensure it is functioning properly and not damaged. Do not use a damaged device.
- This product contains Nickel (Ni) which may cause allergic reactions.
- Only couplers marked "Autoclavable" or "Multimode" can be steam autoclaved. Steam autoclaving couplers not specifically designed for autoclaving would permanently damage the optical components and void the warranty.
- When autoclaving "Autoclavable" or "Multimode" couplers they must be separated from the camera for cleaning and sterilization to minimize the possibility of fogging.
- Before autoclaving, ensure that exposed optical surfaces of the coupler are thoroughly clean and free of residue. Failure to clean the coupler prior to steam autoclaving may result in permanent damage to the optical surfaces.
- To avoid fogging during surgery, both the coupler and the camera must be at room temperature and the optical surfaces entirely free of moisture before connecting the coupler to the camera.
- Do not immerse a hot coupler in cold water or liquid. Forced cooling can cause damage to the coupler and will void your warranty.
- Do not expose non-autoclavable couplers to temperatures over 140° F (60° C).
- Any mechanical manipulation, attempt to repair, or modification of the couplers may result in seal breakage. Do not attempt to disassemble, repair, or modify.
- Couplers are delicate instruments and must be handled carefully.

Adverse Event Reporting

If, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and to your national authority.

complaints@smith-nephew.com

Instructions for Use

Attachment to Camera (A)

Refer to Figure 1.

Coupler mounting threads are compatible with their respective cameras. Before attaching to the camera head, inspect the coupler to ensure the O-Ring is in its groove and is not nicked or damaged. Attachment is made to the camera head by aligning the threads and turning in a clockwise direction (tighten firmly by hand-turning only). The firm attachment allows rotation of the endoscope independent of the camera head.

CAUTION: To avoid fogging during surgery, both the coupler and the camera must be at room temperature and the optical surfaces entirely free of moisture before connecting the coupler to the camera.

Coupler Detachment

Couplers are removed from the camera head by turning in a counterclockwise direction.

Connection of Endoscope (B)

Refer to Figure 1.

Open coupling assembly by depressing latching mechanism. Insert endoscope eyepiece into coupler, seating the eyepiece against the back of the coupler. Ensure the endoscope eyepiece is centered within the coupler. Release the latching mechanism which holds the endoscope in place.

When connecting certain types of flexible endoscopes, you must use the appropriate intermediate adaptor between the scope and coupler. Connect the eyepiece-shaped end of the adaptor to the coupler. Then attach a flexible scope so that the radial pins on the flexible scope proximal end (if there are any) are fully engaged into respective cutouts on the adaptor.

Turn the scope clockwise and lock with the adaptor's thumbscrew.


Figure 1.

Endoscope Detachment

Disconnect the endoscope from the coupler by depressing the locking mechanism and removing the endoscope.

Zoom couplers, specifically designed for Smith & Nephew 3-CCD Camera Heads, have a zoom range of 22–45 mm. The silver scalloped knob is for focus and the black scalloped knob is zoom.

Compatible Devices

Smith & Nephew Direct View scopes

Smith & Nephew Camera Heads (for example, the LENS Camera Head, LENS 4K Camera Head, etc.)

Cleaning

Couplers have two or three optical surfaces that must be cleaned and checked routinely to ensure maximum transmission and image quality. These are the front and rear windows and the eyepiece on the beamsplitter.

1. Soak the device in a bath containing an enzymatic, neutral pH cleaner for 5 minutes.
2. Scrub all crevices using a cleaning brush to remove any visible debris from all crevices taking care to not scratch any of the optical surfaces. If a focus ring is present it should be turned fully in each direction and the exposed areas should be scrubbed.
3. Scrub each optical surface with a gauze pad soaked in an enzymatic, neutral pH cleaner.
4. Rinse and flush the device with warm tap water for 2 minutes.
5. Scrub each optical surface with a gauze pad soaked in isopropyl alcohol. Rinse thoroughly with distilled water.
6. Scrub each optical surface with a gauze pad soaked in acetone. Rinse thoroughly with distilled water.
7. Reflect light off all optical surfaces to detect foreign material remaining or scratches that may degrade image quality.
8. Thermal disinfection 196° F (91° C) for one minute minimum.

Sterilization

Disinfection is not a substitution for sterilization. Always inspect for damage and debris prior to sterilization.

Couplers should be placed in a sterilization tray and wrapped in a suitable packaging for the sterilization process and sterilized prior to use.

Wrapping the coupler (as a means of securing it in an unapproved container) is not acceptable, as contaminants from the wrap (e.g., a towel) may remain as soils on the coupler.

Non-Autoclavable Couplers

Read all Warnings and Precautions before sterilizing.

Autoclavable and Multimode Couplers

Read all Warnings and Precautions before sterilizing.

Note: The coupler and camera head are sterilized separately. Extreme caution must be taken to dry the mating portions of both coupler and camera; otherwise trapped moisture may fog during use and condense on glass surfaces, blurring the image.

Autoclave Wrapped

Follow standard hospital procedures

Pre-Vacuum methods	at "270° F" ("132° C") for 4 minutes
	30 minute dry time
or	
Gravity methods	at "270° F" ("132° C") for 15 minutes
	30 minute dry time

In addition to autoclave wrapped methods, Smith & Nephew Autoclavable and Multimode Couplers can also be sterilized by the following methods:

Non US Customers

Autoclave (Steam Sterilization Parameters)

Method	Cycle	Exposure Temperature	Exposure Time (minutes)	Dry Time (minutes)
UK Cycle (wrapped)	NA	134° C	3	30
World Health Organization (WHO)	NA	135° C	18	30

Note: The UK cycle is not an FDA cleared sterilization cycle and only for use Internationally.

Note: The WHO cycle is not an FDA cleared sterilization cycle and only for use Internationally. This sterilization cycle has not been validated for prion inactivation.

Note: The Coupler can be autoclaved individually.

Chemical Sterilization

Sterilize the device using any one of the following sterilization systems:

- STERRAD® 100NX Standard Cycle
- STERRAD 100S Short Cycle
- STERRAD NX Standard Cycle

Note: The device should be placed in a sterilization tray and wrapped with the STERRAD recommended non-woven polypropylene wrap.

Follow the STERRAD 100S, NX or 100NX System Operator's Manual Instructions for Use when processing these devices.

- STERIS V-PRO cycles:
 - STERIS V-PRO® 1 Standard cycle
 - STERIS V-PRO 1 Plus lumen, non lumen cycles
 - STERIS V-PRO MaX lumen, non lumen, flexible cycles

Note: The device should be placed in a sterilization tray and wrapped with the STERIS recommended non-woven polypropylene wrap

Note: Follow the STERIS V-PRO System Operator's Manual Instructions for Use for the V-PRO cycles when processing these devices

- Follow all recommended directions and warnings included with the sterilization system.

Note: The device can be used immediately after chemical sterilization.

Note: These parameters have been validated to ensure sterility.

Storage

Couplers should be stored either in their shipping box or in a sterilization tray. In either case, proper care should be taken to ensure that they remain immobile to prevent damage. When the couplers are sterilized after use, the devices should be stored in a designated, limited access area that is well ventilated and provides protection from dust, moisture, vermin, and temperature and humidity extremes. Storage is event related and not time related. Sterile packaged devices can be stored as long as sterile packaging is not breached, or until expiration date.

Note: to US customers: Only items sterilized and packaged in materials cleared by the FDA for maintenance of sterility can be stored.

Expected Life

Smith & Nephew couplers are reusable instruments. The functional life of these devices depends on various factors including frequency of use, the method and duration of each use, as well as post-operative methods and handling. Proper maintenance of your equipment is essential to achieve optimal performance over time.

The coupler should be inspected before procedural use.

If any of the below occur, contact Smith & Nephew customer service or an authorized Smith & Nephew representative:

- Device focus ring is locked, spins freely, or required excessive force to use
- Device has internal moisture
- Device lens is cracked
- Device is unable to attach to camera head or scope

Note: Clean & Sterilize this device prior to return or disposal

Environmental Conditions

Condition	Transport and/or Storage
Temperature	-20 – 140° F (-29 – 60° C)
Humidity	15-85% (no condensation allowed)

For Further Information

If further information, including warranty, on this product is needed, contact Smith & Nephew Customer Service at +1 800 343 5717 in the U.S., or an authorized representative.

Reference

Sterilization Standards

- ANSI/AAMI ST79. *Comprehensive guide to steam sterilization and sterility assurance in healthcare facilities*
- ANSI/AAMI/ISO 17665-1. *Sterilization of healthcare products- Moist Heat Part 1 Requirements for development, validation and routine control of a sterilization process for medical devices.*
- ISO 14937. *Sterilization of Health Care Products – General Requirements for Characterization of a Sterilization Agent and the Development, Validation and Routine Control of a Sterilization Process for Medical Devices*
- ISO 17664. *Sterilization of medical devices - information to be provided by the manufacturer for the processing of resterilizable medical devices*
- AAMI TIR12. *Designing, Testing, and Labeling Reusable Medical Devices for Reprocessing in HealthCare Facilities: A Guide for Medical Device Manufacturers*

Cleaning Standards

- Guidance for Industry and FDA Staff - *Reprocessing Medical Devices in Health Care Settings: Validation Methods and Labeling.* U.S. FDA March 17, 2015.
- American Society for Testing and Materials (ASTM). *Standard Guide for Selecting Test Soils for Validation of Cleaning Methods for Reusable Medical Devices, F3208-17.*
- American Society for Testing and Materials (ASTM). *Standard Guide for Application of Test Soils for the Validation of Cleaning Methods for Reusable Medical Devices, F3293-18*
- AAMI TIR12. *Designing, testing, and labeling reusable medical devices for reprocessing in health care facilities: A guide for device manufacturers*
- AAMI TIR 30. *A compendium of processes, materials, test methods, and acceptance criteria for cleaning reusable medical devices*

English

Symbols Glossary

	Medical Device
	EU: Not for general waste
	CE Mark
	Caution
	General Warning Sign
	Unique Device Identifier
	Non-Sterile
	Temperature Limitation
	Humidity limitation

	Authorized representative in the European community
	Device restricted to sale by or on the order of a physician.
	Manufacturer
	Date of manufacture
	Catalog Number
	Serial Number
	Batch code
	Consult instructions for use

Endoscopy
Smith & Nephew, Inc.
150 Minuteman Road
Andover, MA 01810
USA

+1 978 749 1000 Telephone
+1 978 749 1108 Fax
+1 800 343 5717 U.S. Customer Service

www.smith-nephew.com
ifu.smith-nephew.com

Trademark of Smith & Nephew.
©2002 2020 Smith & Nephew, Inc.
All rights reserved.

**ADDENDUM TO INSTRUCTION FOR USE ON THE
MEDICAL DEVICE
«ENDOSCOPIC VIDEOSYSTEM 4K LENS WITH
ACCESSORIES»**

When delivering a product to the Russian Federation, the information contained in this addendum to the instruction for use is a priority

MEDICAL DEVICE (MD) NAME
Endoscopic videosystem 4K LENS with accessories in variations: I. 4K LENS Endoscopic videosystem with Wi-Fi, in composition: 1. 4K LENS System Camera Control Unit with Wi-Fi - 1pc. 2. Manual - 1 pc. 3. Power cord - 1 pc. 4. Camera head - 1 pc. 5. Optical coupler 19.5 mm - 1 pc. 6. Optical coupler 26.5 mm - 1 pc. (optional). II. 4K LENS Endoscopic videosystem, in composition: 1. 4K LENS System Camera Control Unit - 1pc. 2. Manual - 1 pc. 3. Power cord - 1 pc. 4. Camera head - 1 pc. 5. Optical coupler 19.5 mm - 1 pc. 6. Optical coupler 26.5 mm - 1 pc. (optional). Accessories: 1. Kit of accessories, consisting of: 1.1. Cable, HD-SDI - 5 pcs. 1.2. Cable, 3.5 mm mini plug - 2 pcs. 1.3. USB cable - 1 pcs. 1.4. DisplayPort Cable - 1 pcs. 1.5. Cat 5e cable - 1 pcs.
INTENDED USE
Endoscopic videosystem 4K LENS is intended to provide illumination, visualization and capture of still and motion pictures of surgical sites.
INDICATIONS FOR USE
The LENS 4K Camera Control Unit & Camera head is used in diagnostic and operative endoscopic procedures such as arthroscopy to provide illumination, visualization and capture of still and motion pictures of surgical sites within articular cavities. Additionally, LENS 4K Camera Control Unit & Camera head is indicated for use in endoscopic surgical procedures in the thoracic cavity when used with an appropriately indicated thoracoscope. The LENS 4K Camera Head and Camera Control Unit are used in diagnostic and operative endoscopic procedures such as arthroscopy to provide illumination, visualization and capture of still and motion pictures of surgical sites within articular cavities, body cavities, hollow organs and canals. Additionally, LENS 4K Camera Head and Camera Control Unit are indicated for use in endoscopic surgical procedures in the thoracic cavity when used with an appropriately indicated thoracoscope. Smith & Nephew couplers are intended to be used as video accessories to camera heads. They are indicated for use in endoscopic surgical procedures to allow visualization of the articular cavities, body cavities, hollow organs, and canals when used in conjunction with an appropriately indicated endoscope. Additionally, these couplers are indicated for use in endoscopic surgical procedures in the thoracic cavity when used with an appropriately indicated thoracoscope.
POTENTIAL USERS OF MD
The Smith & Nephew LENS 4K Camera Control Unit is intended for use by healthcare professionals in accordance with these instructions for use. The use environment is a Professional healthcare facility in accordance with guidance and manufacturer declaration EMC section. The targeted patient population for the LENS 4K is patients requiring endoscopic surgical procedures as determined necessary by the prescribing healthcare professionals.
ADVERSE EVENT REPORTING
If, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and to your national authority. complaints@smith-nephew.com
THE POSSIBILITY AND WAYS OF INTEGRATING WITH OTHER MEDICAL DEVICES
This device is not compliant with MR.

Endoscopic videosystem LENS 4K

The entire System can be used in conjunction with compatible high-definition flat-panel displays, image capturing devices, as well as with printing devices and video recorders;

Camera head LENS 4K;

Video adapter;

LENS 4K Wi-Fi antenna.

Optical Lens Coupler

Compatible Endoscopes :

Arthroscopes: 1.9, 2.7, and 4.0 mm sizes

Laparoscopes: 5.0, 5.5, and 10 mm size

Hysteroscopes: 4.0, 5.0, 5c, and 8.0 mm

Thorascopes: 2.8 mm

Camera Head

It is compatible with the LENS 4K System Control Unit;

LENS 4K System Control Unit with Wi-Fi;

Compatible Endoscopes and Couplers (Must have an industry standard C-mount threading of 1-32-UN- 28)

Arthroscopes: 1.9, 2.7, and 4.0 mm sizes

Laparoscopes: 5.0, 5.5, and 10 mm size

Hysteroscopes: 4.0, 5.0, 5c, and 8.0 mm

TECHNICAL CHARACTERISTICS OF MD**4K LENS Camera Control Unit Specifications**

Dimensions (LxWxH) cm	39,37 x 35,56 x 11,68
Weight, kg	6,2
Display Output, Hz	60
Resolution, pixels	3840 x 2160 (native) 4K
Shutter	1/60 second through 1/10,000 second shutter speed (60 Hz)
AC Leakage Current, µA	420,1
Fuses	F10AH 250 V; 250 V, 10 A
Modes Operation	Continuous
Power requirements	100-240VAC, 50/60 Hz; dual fuse
Equipment Classification	Class 1, type CF.
Ingress Protection	IPX0
Wi-Fi	2.4Ghz 802.11 b/g/n WPA-2 RF Module
USB Storage	8 GB minimum; 32 GB recommended. FAT32 (for drives smaller than 32 GB) or NTFS partitions. Note: Drives must be formatted to support MBR partition tables. The system does not support GPT formatted drives. The following devices have been validated for use with the System: • 8 GB USB flash drives (Transcend, Corsair™, Imation™, Kingston) • 16 GB USB flash drives (ADATA™, Mushkin, Silicon Power™) • 32 GB USB flash drives (Corsair, Patriot™, PNY™, Silicon Power) • 64 GB USB flash drives (Mushkin, & ADATA) • 1 TB external SSD hard drives (Buffalo & Western Digital)
USB Output	Two (2) Type A connectors; provides USB-compatible connection to flash or hard drives.
USB Type	USB 2.0 or greater
Video Input	NTSC or PAL composite input
Video Outputs	
HD Serial Digital Interface (HD-SDI)	One serial digital interfaces (HD-SDI) with 4:2:2 component digital video output on BNC connector
4K Serial Digital Interface (4K-SDI)	Four serial digital interfaces (HD-SDI) with 4:2:2 component digital video output on BNC connector
DisplayPort (DP)	Single stream transport (SST) version 1.2
Video Accessory Controls	Two mini-phonon connectors (3.5 mm)

± 10% accuracy, unless otherwise specified

4K LENS Camera Head Specifications

Camera Head Dimensions (LxWxH), cm	1.7" x 2.0" x 3.5" (Width x Height x Length)
Weight, kg	0,7
Dimensions, cm, Camera Head cable	365,76 +/- 7,62
Equipment Classification	Protection against electrical shock Class 1, type CF equipment, applied part camera head bezel
Video Sensor Technology	4K 3-Chip MOS
Video Data Transmission	High-speed, low-voltage differential serial (LVDS), digital transmission via 12-foot camera cable
Video Output Format	Full UHD 2160p, 100% compatible with SMPTE ST 425-5 when used with the LENS 4K Camera Control Unit
Video Field Rate, Hz	60
Transport Shock and Vibration Compliance	ISTA 3A
Lens Mount	1.00 - 32 UN 2B C-mount thread
Modes of Operation	Continuous

± 10% accuracy, unless otherwise specified

Optical Lens Coupler 19.5 mm Specifications

Length, mm	50
Diameter, mm	50
Weight, g	110

± 10% accuracy, unless otherwise specified

Optical Lens Coupler 26.5 mm Specifications

Length, mm	50
Diameter, mm	50
Weight, g	110

± 10% accuracy, unless otherwise specified

The full list of technical specifications for optical lens is available per request.

ELECTROMAGNETIC EMISSIONS

Guidance and Manufacturer's Declaration – Electromagnetic Emissions.

R a n g e 1	Guidance and Manufacturer's Declaration – Electromagnetic Emissions		
	2 The Smith & Nephew 4K LENS System is intended for use in the electromagnetic environment specified below. The customer or the user of the 4K LENS System should assure that it is used in such an environment.		
3	Emissions Test	Compliance	Electromagnetic Environment – Guidance
4	RF emissions CISPR 11	Group 1	The 4K LENS System uses RF energy only for its internal functions. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
5	RF emissions CISPR 11	Class A	The 4K LENS endoscopic video system is suitable for use in any locations, except for residential houses and buildings directly connected to the public mains network that supplies energy to residential buildings. The system can be used in residential houses and buildings directly connected to the public mains network that supplies energy to residential houses provided that the following warning is in place:
6	Harmonic emissions 61000-3-2 EN	Class A	
7	Voltage fluctuations/	Compliance	Warning. This system is intended for use by healthcare professionals only. This system may cause poor radio reception and disrupt the operation of equipment in the vicinity. In this case, it may be necessary to take mitigation measures such as

flicker emissions EN 61000-3-3	reorienting, relocating the 4K LENS video endoscope system, or shielding the location area.
-----------------------------------	---

Guidance and Manufacturer's Declaration – Electromagnetic Immunity

Guidance and Manufacturer's Declaration – Electromagnetic Immunity

4K LENS System is intended for use in the electromagnetic environment specified below. The customer or user of the 4K LENS System should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment – Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	+/- 8 kV contact +/- 15 kV air	+/- 8 kV contact +/- 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	+/- 2 kV for power supply lines +/- 1 kV for input/output lines	+/- 2 kV for power supply lines +/- 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	+ 2 kV line to P.E. + 1 kV line to line	+ 2 kV line to P.E. + 1 kV line to line	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% Unom, ½ cycle 40% Unom, 5 cycles 70% Unom, 25 cycles 0% Unom, 250 cycles	0% Unom, ½ cycle 40% Unom, 5 cycles 70% Unom, 25 cycles 0% Unom, 250 cycles	Mains power quality should be that of a typical commercial or hospital environment. If the user of the LENS 4K Camera Control Unit requires continued operation during power mains interruptions, it is recommended that the LENS 4K Camera Control Unit be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Note: Unom = Nominal AC input voltage.

Guidance and manufacturer's declaration – Electromagnetic Immunity

Guidance and manufacturer's declaration – electromagnetic immunity

The Endoscopic 4K LENS Video System is designed for use in electromagnetic environments defined below. Buyer or user of the system should ensure its use in the specified electromagnetic environment.

Immunity test	IEC 60601 test level	Compliance Level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 V (root mean-square value) 150 kHz to 80 MHz	3 V (root mean-square value) 150 kHz to 80 MHz	Portable and mobile RF communications equipment should be used no closer to any part of the Controller, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance
	3 V/m	3 V/m	Recommended separation

Radiated RF IEC61000- 4-3	80 MHz to 2,7 GHz	80 MHz to 2,7 GHz	distance : $d = 1,2\sqrt{P}$ $d = 1,2\sqrt{P}$, (80 MHz to 800 MHz); $d = 2,3\sqrt{P}$, (800MHz to 2.7 GHz), Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strength from fixed RF transmitters, as determined by an electromagnetic site survey, a should be less than the compliance level in each frequency range. b Interference may occur in the vicinity of equipment marked with the following  symbol:
---------------------------------	-------------------------	-------------------------	---

a - Field strength from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Controller is used exceeds the applicable RF compliance level above, the Controller should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the Controller.

b - Over the frequency range 150 kHz to 80 MHz, field strength should be less than 3 V/m.

Notes:

1. At 80 MHz and 800 MHz, the higher frequency range applies.

2. These guidelines may not apply in all situations.

Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Recommended separation distance between portable and mobile RF communication equipment and the 4K LENS Video System

Recommended separation distance between portable and mobile RF communication equipment and the 4K LENS endoscopic video system

4K LENS endoscopic video system is intended for use in the electromagnetic environment in which radiated RF disturbances are controlled. The customer or the operator can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the System as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power P of transmitter, W	Separation distance d according to frequency of transmitter, m		
	$d = 1,2\sqrt{P}$, 150 kHz to 80 MHz	$d = 1,2\sqrt{P}$, 80 MHz to 800 MHz	$d = 2,3\sqrt{P}$, 800 MHz to 2.5 GHz
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23

Notes

1. At frequencies 80 and 800 MHz, a higher field strength value is applied.

2. The above expressions are not applicable in all cases. The spread of electromagnetic waves is affected by absorption or reflection from structures, objects and people.

3. When evaluation the recommended values of spatial diversity d for transmitters with rated maximum output power not specified in the table, the above expressions substitute the rated maximum output power P in watts specified in the transmitter manufacturer's documentation.

DEVICES INCORPORATING NON-VIABLE MATERIALS OF ANIMAL ORIGIN

Medical device do not contain human or animal-derived components.

DEVICES INCORPORATING A MEDICINAL SUBSTANCE

The subject Smith & Nephew products do not incorporate a medicinal substance, nor is it present in the manufacturing process or present in the final device, therefore this section does not apply.

PACKAGING AND LABELLING INFORMATION**Draft of labeling****Article:****Name of the medical device: (according to RU)****Size: (If applicable to this product).****Packs: (If applicable to this product).****Date of manufacture, Use before, lot/series number: See on the package.****Manufacturer's name and address:****Authorized representative on the territory of the Russian Federation:****Sterilization method: (If applicable to this product).****Registration certificate number:****Designation, method of application, contraindications: Refer to instructions****Terms and conditions of storage:****Precautions: (If applicable to this product).****Instructions for use:****Information on passing the mandatory certification (PCT Mark)****Information on compliance with the Technical Regulations of the Customs Union (EAC Mark);***** PCT/EAC Mark entered after passing the certification procedure (if there is a completed document confirming compliance).****Data of the packaging****4K LENS System Camera Control Unit with Wi-Fi or without Wi-Fi**

4K LENS System Camera Control Unit with a Wi-Fi module or without a Wi-Fi module is supplied in a corrugated cardboard shipping package that includes the 4K LENS video camera console with a power cable, Operation manual, Wi-Fi antenna (for an endovideocamera console with Wi-Fi module), foam separators located on the front and back of the product for safe transportation of the medical product, as well as transport packaging with accessories made of corrugated cardboard (if necessary). Transport packaging with accessories includes foam separators for safe transportation of accessories.

Lens optical coupler

The optical Lens coupler is delivered in a transport package made of corrugated cardboard, which includes an optical Lens coupler (with a protective cap), foam separators for safe transportation of the medical device, as well as instructions for use for the medical device.

Camera Head

The LENS 4K Camera Head is delivered in a corrugated cardboard transport package, which includes the LENS 4K system camera head, foam separators for the safety of transporting the medical device, as well as instructions for use for the medical device.

RISK INFORMATION

The performed risk analysis provides assurances that the relevant risks have been identified and either eliminated or reduced as far as possible. For the remaining risks, risk management measures are presented, for each risk appropriate, to minimize the risk of possible failures during the estimated expected life of the product.

The results of the risk analysis show that the use of the product for its intended purpose does not adversely affect the health or safety of the patient, user or other persons.

INFORMATION ABOUT SOFTWARE USED IN MD

The software used in a medical device is classified as hazard class B.

LIST OF MAIN APPLICABLE STANDARDS

Standard	Description
EN 60601-1	Medical Device Electrical equipment equipment-Part 1 General requirements for Basis Safety and Essential Performance
EN 60601-1-2	Medical Electrical Equipment Part 1-2: General Requirements for Basis Safety and Essential Performance – Collateral Standard: Electromagnetic disturbances
EN 60601-2-18	Medical Electrical Equipment: Particular Requirements for the Safety of Endoscopic Equipment
Directive 2015/863	RoHS-Restrictions of the Use of certain Hazardous
Directive 2011/65/EU	RoHS-Restrictions of the Use of certain Hazardous Substances
REACH 1907/2006 (33)	Registration, evaluation, authorization and restriction of chemicals
Directive 2014/53/EU	Radio Equipment Directive
EN ISO 14971	Application of risk management to medical devices
EN ISO 13485	Medical devices – Quality management system – Requirements for regulatory purposes

USE & OPERATION CONDITIONS OF MD

Endoscopic - 4K LENS videosystem

Environmental requirements during surgery

Temperature	From +10 °C to + 40 °C
Humidity	30 % - 75 %
Pressure	From 70 kPa to 106 kPa
Altitude	3000 meters maximum

The camera head and optical coupler lens

Environmental requirements during surgery

Temperature	From +10 °C to + 25 °C
Humidity	30 % - 75 %
Pressure	From 70 kPa to 106 kPa
Altitude	3000 meters maximum

SHELF LIFE

The Smith & Nephew LENS 4K Camera Control Unit is a reusable instrument. The functional life of these devices depends on various factors, including frequency of use, the method and duration of each use, as well as post-operative methods and handling. Proper maintenance of the equipment is essential to achieve optimal performance over time.

The LENS 4K Camera Control Unit should be inspected before each procedural use per the "Service Life Checklist", if any issues occur, perform "Troubleshooting". If the problem persists, contact Smith & Nephew customer service or an authorized Smith & Nephew representative.

The Smith & Nephew LENS 4K Camera Head is a reusable instrument. The functional life of these devices depends on various factors, including frequency of use, the method and duration of each use, as well as post-operative methods and handling. Proper maintenance of the equipment is essential to achieve optimal performance over time.

The LENS 4K Camera Head should be inspected before each procedural use per the "Service Life Checklist". If any issues occur, contact Smith & Nephew customer service or an authorized Smith & Nephew representative.

Smith & Nephew couplers are reusable instruments. The functional life of these devices depends on various factors including frequency of use, the method and duration of each use, as well as post-operative methods and handling. Proper maintenance of your equipment is essential to achieve optimal performance over time.

Medical equipment under normal operating conditions in compliance with the instruction manual is designed for an unlimited service life.

The expected life is 1 year.

HOW SUPPLIED

Contents :

I. 4K LENS Endoscopic videosystem with Wi-Fi, in composition:

1. 4K LENS System Camera Control Unit with Wi-Fi - 1pc.
2. Manual - 1 pc.
3. Power cord - 1 pc.
4. Camera head - 1 pc.
5. Optical coupler 19.5 mm - 1 pc.
6. Optical coupler 26.5 mm - 1 pc. (optional).

II. 4K LENS Endoscopic videosystem, in composition:

1. 4K LENS System Camera Control Unit - 1pc.
2. Manual - 1 pc.
3. Power cord - 1 pc.
4. Camera head - 1 pc.
5. Optical coupler 19.5 mm - 1 pc.
6. Optical coupler 26.5 mm - 1 pc. (optional).

Accessories:

1. Kit of accessories, consisting of:
 - 1.1. Cable, HD-SDI - 5 pcs.
 - 1.2. Cable, 3.5 mm mini plug - 2 pcs.
 - 1.3. USB cable - 1 pcs.

- 1.4. DisplayPort Cable - 1 pcs.
1.5. Cat 5e cable - 1 pcs.

DISPOSAL RULES AND CONDITIONS

The device is subject to the Directive 2002/96 / EC (European Community) for waste electrical and electronic equipment (WEEE). In accordance with the registration status of this directive, this device cannot be disposed of with ordinary electrical and electronic equipment.

After use, the product may be potentially hazardous. Dispose of the product in accordance with local laws and regulations. The manufacturer provides product safety information upon request.

Class of medical waste that had no contact with the internal environment of the human body - A.

WARRANTY AND CLAIM

Smith & Nephew LLC (S&N) warrants all of its products to the original purchaser (Customer) against defect in workmanship and materials for one hundred eighty (180) days from date of purchase from Smith & Nephew; except for NAVIO® products, one (1) year from date of installation (Warranty Period). For any products found to not be in conformance with this warranty during the Warranty Period, this warranty provides and is restricted to, as elected by S&N, either (i) repair or replacement of such products without charge and within a reasonable period of time or (ii) a refund or credit in the amount of the purchase price of such products. If S&N repairs or replaces product under this warranty and requests Customer to return such product, Customer must ship such product to S&N freight prepaid by Customer. Customer shall be invoiced for any replacement product if Customer does not return the requested replaced product within thirty (30) days after S&N's shipment of the replacement product. This warranty does not cover and is voided by any of the following: (i) product packaged or labeled by someone other than S&N or its authorized agents; (ii) product not used in compliance with the specifications, instructions or claims for use of the product; (iii) equipment product used in conjunction with disposables or accessories not specified for use with such equipment; (iv) equipment product used in conjunction with reprocessed disposables or accessories; (v) modification of product; (vi) product past its expiration date; (vii) normal wear and tear; (viii) damage due to misuse, reprocessing, alteration, unauthorized repair or negligent handling or damage due to lack of care by the owner, user or handler of the product including but not limited to storage, handling or cleaning; or (ix) any other damage inflicted to products by the owner, user or handler. A 90-day warranty on repairs applies to the defective component repaired. This warranty applies only to the original purchaser from S&N (or its authorized distributor) and is not transferable. EXCEPT TO THE EXTENT PROHIBITED OR OTHERWISE REQUIRED BY APPLICABLE LAW, THIS WARRANTY IS THE SOLE WARRANTY OF S&N, AND ALL OTHER WARRANTIES OF ANY KIND OR DESCRIPTION WHATSOEVER, INCLUDING WARRANTIES OF MERCHANTABILITY, SATISFACTORY QUALITY AND FITNESS FOR A PARTICULAR PURPOSE, EXPRESSED OR IMPLIED, ARE EXCLUDED.

For all questions relating to the quality of products, as well as repair and maintenance, please contact an **AUTHORIZED REPRESENTATIVE ON THE TERRITORY OF THE RUSSIAN FEDERATION:**

Smith & Nephew, LLC, Russia, 105120, Moscow, 2nd Syromyatnichesky lane., 1, 9 floor, flat 1.

Tel: +7 (495) 755-55-03

NAME AND LEGAL ADDRESS OF MANUFACTURER

DEVELOPER:

Smith & Nephew, Inc. Endoscopy Division
150 Minuteman Road, Andover, Massachusetts, 01810 - 1031, USA

MANUFACTURER:

Smith & Nephew, Inc. Endoscopy Division
150 Minuteman Road, Andover, Massachusetts, 01810 - 1031, USA

MANUFACTURING SITE:

Smith & Nephew, Inc. Endoscopy Division
76 South Meridian Ave, Oklahoma City, Oklahoma, 73107-6512, USA.

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