JAMES I. SEXTON

NOTARY PUBLIC

Bishopsgate,
Henry Street,
Limerick, P.O. Box 146,
Ireland.

Telephone + 353 61 445513
Facsimile + 353 61 316597
email james.sexton@homs.ie


NOTARIAL CERTIFICATE

TO ALL TO WHOM THESE PRESENTS SHALL COME: GREETING

I, **JAMES I. SEXTON**, Notary Public duly appointed by the Chief Justice of Ireland to practise for life as a Notary Public in the Counties of Limerick and Clare and in the adjoining counties, certify that on 8th July 2020, **Dermot Scanlon**, known to me as Chief Executive Officer of Serosep Limited of Annacotty Business Park, Annacotty, Co. Limerick, Ireland, personally in my presence acknowledged that he did duly sign on behalf of Serosep Limited the attached **Serosep Instruction Manual document** on four (4) pages, together with the accompanying translation into the Russian language on five (5) pages, by signing his name on the first page of each of the original version in English and of the translation version in Russian as the same now appears, and I certify that the signature subscribed thereto is the true and proper signature of the said **Dermot Scanlon** as Managing Director of Serosep Limited, and that I have signed and stamped each page for identification.

IN FAITH AND TESTIMONY WHEREOF I have hereunto subscribed my name and affixed my seal this 9th July 2020.

Signed:

James I. Sexton

JAMES I. SEXTON
Notary Public


I AFFIRM

Chief Executive Officer

Serosep Ltd.
Annacotty Business Park, Annacotty,
Co. Limerick. V94 FF83, Ireland
Tel: +353 61 358190
Email: info@serosep.com

Serosep Limited

 (Dermot Scanlon)

08 July 2020
COMPANY STAMP

Instruction manual for an in vitro diagnostic medical device
"10% Buffered Neutral Formalin / HistoPot"

Purpose

10% buffered neutral formalin / HistoPot is intended for use as a fixator when processing biological tissues or clinical samples.

Intended user

HistoPot pre-filled specimen containers are used by trained personnel, i.e. Medical Scientists & Clinical Theatre staff. It is intended to be used in vitro for the examination of specimens derived from the human body. The HistoPot pre-filled specimen containers are not to be used by laypersons or patients.

Composition of product

35% Formaldehyde, sodium dihydrogen phosphate monohydrate, sodium hydrophosphate heptahydrate, deionized water

Product Codes

S10-B-FOR-"XX" 10% Buffered Formalin

Note: "XX" indicates HistoPot pot size, e.g. 20ml for 20ml pot

Method of application

10% buffered neutral formalin / HistoPot is delivered in ready-to-use form complete with instructions that are put into group packaging.

When tissue is removed from the body, it undergoes a process of self-destruction or autolysis. In order to stop the changing of material, to maintain the specimen's structural features to enable a clear staining, the specimen must be fixed.

Tissue should be fixed before processing is initiated. Insufficient fixed tissues are not adequately protected against the effects of processing.

Formalin solutions are used to fix all kinds of histopathology specimens. After fixating the specimens, they can be processed and stained. Stained specimens can be examined by microscope.

Tissue samples should be placed in a formalin solution of 10% buffered neutral formalin / HistoPot immediately after being removed from the body.

Tissue containers should be labeled in accordance with standard laboratory procedures and delivered to the laboratory.

Warning and precautions

Use of 10% buffered neutral formalin / HistoPot should only be performed in well-ventilated areas by trained professionals wearing appropriate personal protective equipment (PPE): gloves, protective clothing and goggles to avoid contact with skin, eyes and mucous membranes.

Spillage of 10% buffered neutral formalin / HistoPot should also be avoided.

10% Buffered Neutral formalin / HistoPot contained in sample containers should not be removed from the container or allowed to come into contact with human skin.

In case of failure of the product or deviations in its functioning, it is necessary to contact the authorized representative of the manufacturer in the territory of the Russian Federation: JSC "AMTEO M," Russian Federation, 111539, Moscow, Veshnyakovskaya Str., house 11 A
 Phone/fax: (495) 921 35 78, e-mail: info@amteo.ru

The risks to human health from exposure to formalin solutions are outlined in the table below:

10% Formalin Solution	
Classification according to Regulation (EC) No 1272/2008	
	GHS08 health hazard
Muta. 2	H341 Suspected of causing genetic defects.
Carc. 1B	H350 May cause cancer.
	GHS07
Skin Sens. 1	H317 May cause an allergic skin reaction.

For full safety information please refer to the Material Safety Data Sheet (MSDS) (available on request).

About undesirable events that have signs of an adverse event (incident), it is necessary to send a message to the authorized representative of the manufacturer in the Russian Federation: JSC "AMTEO M", Russian Federation, 111539, Moscow, ul. Veshnyakovskaya, building 11 A
 Phone / Fax: (495) 921 35 78, e-mail: info@amteo.ru

Contraindications

Exposure to direct sunlight and freezing is not allowed. Do not use an expired product.

List of national standards applied by the manufacturer

- GOST 1625-2016 FORMALIN TECHNICAL. SPECIFICATIONS
- GOST R ISO 15223-1-2014 "Medical products. Symbols used for marking on medical products, labels and accompanying documentation. Part 1. Main requirements";
- GOST ISO 10993-1-2011 "Medical products. Evaluation of biological action of medical products. Part 1. Assessment and researches";
- GOST ISO 10993-9-2015 "Medical products. Evaluation of biological action of medical products. Part 9. Basic Principles for Identification and Quantification of Potential Destruction Products";
- GOST ISO 10993-10-2011 "Medical products. Evaluation of biological action of medical products. Part 10. Studies of Irritant and Sensitizing Effects";
- GOST ISO 10993-11-2011 "Medical products. Evaluation of biological action of medical products. Part 11. Researches of all-toxic action";
- GOST ISO 10993-12-2015 "Medical products. Evaluation of biological action of medical products. Part 12. Preparation of samples and control samples";
- GOST ISO 10993-18-2011 "Medical products. Evaluation of biological action of medical products. Part 18. Study of chemical properties of materials";
- GOST R ISO 10993-2-2009 "Medical products. Evaluation of biological action of medical products. Part 2. Requirements for the treatment of animals";
- GOST 31214-2003 "Medical products. Requirements for samples and documentation submitted for toxicological, sanitary and chemical tests, sterility and pyrogenicity tests";
- GOST R 52770-2007 "Medical products. Safety requirements. Methods of sanitary-chemical and toxicological tests."

GOST 31508-2012 "Medical products. Classification depending on potential application risk. General requirements".

Transportation conditions

Transportation of 10% buffered neutral formalin / HistoPot is performed by all modes of covered transport in accordance with the rules of cargo transportation in this mode of transport at temperatures of 2-30° C and relative air humidity not more than 80% at 25° C. Transportation is performed in a vertical position without excessive load from above. Direct sunlight as well as freezing should be avoided.

Storage conditions

10% buffered neutral formalin HistoPot is delivered ready to use and keeps the properties when stored up to three years at 2-30°C and has to be stored at relative air humidity no more than 80% at 25°C (in accordance with GOST P 15150-69) during shelf life. The covers of containers and canisters must be tightly closed. Direct sunlight as well as freezing is not allowed. Shelf life (including after opening): 3 years from the date of production. Items transported and stored in violation of temperature conditions shall not be used.

Procedure, conditions and precautions for disposal

After the expiration date, 10% buffered neutral formalin / HistoPot is subject to utilization. Disposal is carried out by special licensed organizations on specially equipped sites in accordance with the requirements stipulated by the existing Federal Laws, compliance with mandatory requirements for environmental protection. Used products are stored in waste collection tanks of class B and disposed of in accordance with SanPiN 2.1.7.2790.

It is not allowed to dispose of 10% of buffered neutral formalin / HistoPot into sewerage, surface or ground water, as well as disposal together with household garbage.

Symbols used in product labelling

- Name of manufacturer and its address: Serosep Limited, Annacotty Business Park, Annacotty, Limerick, V94FF83, Ireland Phone +353 (0) 61 358 190 www.serosep.com info@serosep.com	
- Batch number	
- Catalog number	
- Number of items in the package	pcs
- Name of product	10% buffered neutral formalin HistoPot
- Date of production	(date/month/year)
- Shelf life	
- Temperature limit	
- Do not reuse	

- "Refer to Instructions"	
- <i>in vitro</i> diagnostic product	
- Protect from damp	
- Fragile, handle carefully	
- Protect from sunlight	
- Limit on the number of tiers in the stack	
- Packaging of products in cardboard box	
- Possibility of recycling of polymer products	
- Do not use in case of damaged package	
- CE compliance	
- Compliance with technical regulations	
- Long-term health hazards	
- Sign of danger	

Authorized representative of the manufacturer in the territory of the Russian Federation:
 JSC AMTEO M,
 11 A Veshnyakovskaya str., Moscow, 111539, Russian Federation
 Phone/fax: (495) 921 35 78, e-mail: info@amteo.ru

Revision: February 2019

ZEP

*James I. Sexton*JAMES I. SEXTON
Notary Public

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