Medtronic

710 Medtronic Parkway
Minneapolis, MN 55432
USA
www.medtronic.com


tel: 001-763-514-4000
fax: 001-763-514-4879

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To Whom It May Concern

Medtronic Inc., a corporation organized under the laws of the United States, having its principal place of business at 710 Medtronic Parkway, Minneapolis, MN 55432, USA, as the legal manufacturer of 'Telescope Guide Extension Catheter' hereby confirm that the attached documentation is a true and original document.

Yours sincerely,

Signed on behalf of Medtronic Inc.

Mike Green


Senior Regulatory Affairs Manager

05 Dec 2020

Date:

As a registered Notary Public for the State of Minnesota,

I hereby certify that this is a true and accurate copy of the
original document. The original is filed at Medtronic, Inc.

 05 Dec 2020
 **JENNIFER KASPER SMITH**
NOTARY PUBLIC - MINNESOTA
MY COMMISSION EXPIRES 01/31/25

The information specified in the Appendix to the Instruction for use is a priority on the territory of the Russian Federation.

Appendix to the Instruction for use

Telescope Guide Extension Catheter

PRESCRIPTION OF MEDICAL DEVICE

Telescope Guide Extension Catheter is intended to be used in conjunction with guide catheters to access discrete regions of the coronary and/or peripheral vasculature, and to facilitate placement of interventional devices.

ADVERSE EVENTS (POTENTIAL COMPLICATIONS)

Potential issues that can be associated with the guide extension catheter include, but are not limited to, the following adverse effects:

- Death
- Embolism
- Slow-flow
- Foreign body in patient
- Stent dislodgement
- Thrombus
- Vessel dissection
- Vessel perforation
- Occlusion (CVA/MI/spasm)
- Aneurysm
- Avulsion
- Reaction¹
- Infection
- Blood loss/hemorrhage/hematoma (vascular access complication)
- Renal failure (contrast induced nephropathy)

The occurrence of the above listed complications may lead to the need for a surgical intervention.

PACKAGING CONTENTS/QUANTITY OF CHAMFERED SHELF CARTON:

1. Telescope Guide Extension Catheter – 1
2. Instructions for Use - 1

¹ An allergic reaction can result from a non-visible foreign matter, chemical, particulate, or drug during the procedure/use of Telescope. Irritation, itching, rash, local, reaction, skin inflammation, and skin irritation are all possible reactions to a realized hazard (eg, sterile barrier is compromised).

Available upon request

INFORMATION ABOUT STICKER

The final packaging for the Russian market is accompanied with additional sticker.

Devices supplied in packaging marked with the following possible information:

- Name of medical device in Russian language;
- Name of the manufacturer;
- Country of origin;
- Code UPN (unique product number);
- Code CFN (customer facing number);
- Date and number of Registration Certificate;
- Sign of conformity in Russia;
- Warning and informing symbols and inscriptions (if applicable)
- Non pyrogenic

Sample of the sticker:

	Name of medical device Telescope Guide Extension Catheter, variant	
	Manufacturer:	Code UPN
	Country of origin:	Code CFN
	Manufacturer's authorized representative::	
	Regestretion Certificate № _____ or _____	
	Date of manuf../Lot/UBD: See on pack	See IFU See packaging symbols

MAIN TECHNICAL DESCRIPTION

PRODUCT CHARACTERISTIC	6F (specifications with tolerance)	7F (specifications with tolerance)
Effective Length	150 cm +2 cm /- 5 cm	150 cm +2 cm /- 5 cm
Weight	3.5 g	3.5 g
Maximum Tip ID/ Inner Diameter (ID) distal section	0.056" (1.42 mm)	0.062" (1.57 mm)
Maximum Tip OD/ Outer Diameter (OD) distal section	0.070" (1.78 mm)	0.078" (1.98 mm)
Distal Catheter Length (Distal Tip to start of Entry Port)	251 mm + 9.6 mm /- 8.6 mm	251 mm + 9.6 mm /- 8.6 mm
Tip Length	2.2 mm ± 0.5 mm	2.2 mm ± 0.5 mm
Length of Entry Port	40 mm ± 5 mm	40 mm ± 5 mm

Available upon request

PRODUCT CHARACTERISTIC	6F (specifications with tolerance)	7F (specifications with tolerance)
Pushwire Main Body Profile (Fully tubular profile)	0.019" $\pm$ 0.0001" (0.48 mm $\pm$ 0.0025 mm)	0.019" $\pm$ 0.0001" (0.48 mm $\pm$ 0.0025 mm)
Guidewire Compatibility	0.014" (0.36 mm)	0.014" (0.36 mm)
Guide Catheter Compatibility	6 Fr (2.0 mm)	7 Fr (2.3 mm)
Marker Location (Femoral)	1000 mm +15.6/- 14.6 mm	1000 mm +15.6/- 14.6 mm
Marker Location (Brachial)	900 mm +15.6/-14.6 mm	900 mm +15.6/-14.6 mm
Length of Hydrophilic coating	213 mm $\pm$ 10.6 mm	213 mm $\pm$ 10.6 mm
Pushwire Integrity	$\geq$ 300 gf ($\geq$ 2.9 N)	$\geq$ 300 gf ($\geq$ 2.9 N)
Deliverability	$<$ 307 gf ($<$ 3.0 N)	$<$ 522 gf ($<$ 5.1 N)
Tip compression	$<$ 66 gf ($<$ 0.6 N)	$<$ 76 gf ($<$ 0.7 N)
Lumen Patency	$\geq$ 0.0540" ($\geq$ 1.4 mm)	$\geq$ 0.0600" ($\geq$ 1.5 mm)
Main to proximal jacket tensile	Min Peak Tensile: 10 N	Min Peak Tensile: 15 N
Proximal Jacket to Pushwire Bond Tensile	Min Peak Tensile: 10 N	Min Peak Tensile: 15 N
Kink resistance	GEC distal section will not kink around a bend radius of 15 mm	GEC distal section will not kink around a bend radius of 15 mm
Straight Line Torque	The Device must be capable of withstanding 0.5 rotations before failure occurs	The Device must be capable of withstanding 0.5 rotations before failure occurs
Torque to Kink	The Device must be capable of withstanding 0.5 rotations before failure occurs	The Device must be capable of withstanding 0.5 rotations before failure occurs
Radiographic	Visibility - Medical device is radiopaque at the tip and the entry port for identification under fluoroscopy	
Corrosion Resistance	The metal components of the device have corrosion resistance during use, transportation, storage in accordance with ISO 10555-1	

STERILIZATION

1. This device is single-use only.
2. Device is sterile.
3. Sterilized by EO.
4. Sterility Assurance Level (SAL) - 10^{-6}
5. Term preservation of sterility for 2 years.
6. Primary packaging is sterilized in accordance with ISO 11607.

CONDITIONS OF TRANSPORTATION AND STORAGE

The recommendation is to "Store in a cool, dry, and dark place."
Based on the temperatures the device was subjected to during testing, storage at between 15°C and 30°C is recommended for shelf life of the product. Temperature excursions down to -35°C for 2 days and up to 55°C for 75 days are supported by testing.

Available upon request

SHELF LIFE

2 years.

PROCEDURE AND CONDITIONS FOR THE DISPOSAL OR DESTRUCTION OF A MEDICAL DEVICE

After use, this product may be a potential biohazard (hazard class B in accordance with SanPiN 2.1.7.2790-10). Handle and dispose of all such device in accordance with accepted medical practice and applicable local, regional and national laws and regulations.

OPERATING CONDITIONS

The products can be used under the standard conditions of a hospital operating room environment. There are no additional assembly, installation or maintenance requirements.

The product is intended for use by trained medical personnel in the operating room: at room temperature (18°C) not higher than 25°C, relative humidity of 75%, a pressure of from 86.0 to 106.0 kPa.

THE LIST OF NATIONAL AND INTERNATIONAL REGULATORY DOCUMENTS/STANDARDS, WHICH MATCHES MEDICAL DEVICE

Standard #	Revision #	Title of Standard
EN ISO 14971	2012	Medical Devices – Application of Risk Management to Medical Devices (ISO 14971:2007, Corrected version 2007-10-01)
EN ISO 14644-1	2015	Cleanrooms and associated controlled environments – Part 1: Classification of air cleanliness by particle concentration (ISO 14644-1:2015)
EN ISO 14644-2	2015	Cleanrooms and associated controlled environments – Part 2: Monitoring to provide evidence of cleanroom performance related to air cleanliness by particle concentration (ISO 14644-2:2015)
EN ISO 13485	2016 & AC:2016 & AC:2018	Medical Devices - Quality Management Systems - Requirements for Regulatory Purposes (ISO 13485:2016)
EN ISO 10555-1	2013+A1:2017	Intravascular catheters - Sterile and single-use catheters – Part 1: General requirements (ISO 10555-1:2013/A1:2017)
EN ISO 11070	2014+A1:2018	Sterile single-use intravascular introducers, dilators and guidewires - Second Edition (ISO 11070:2014/A1:2018)
EN 62366-1	2015	Medical devices - Part 1: Application of usability engineering to medical devices (IEC 62366-1:2015)
EN ISO 11607-1	2017	Packaging for terminally sterilized medical devices -- Part 1: Requirements for materials, sterile barrier systems and packaging systems (ISO 11607-1:2017)

Available upon request

EN ISO 11607-2	2017	Packaging for terminally sterilized medical devices --Part 2: Validation requirements for forming, sealing and assembly processes (ISO 11607-2:2017)
EN ISO 15223-1	2016	Medical devices- Symbols to be used with medical device labels, labelling, and information to be supplied – Part 1: General requirements (ISO 15223-1:2016)
EN 1041	2008+A1:2013	Information supplied by the manufacturer of medical devices (EN 1041:2008+A1)
EN 556-1	2001 AC:2006	Sterilization of Medical Devices - Requirements for Medical Devices to be Designated "Sterile" – Part 1: Requirements for Terminally Sterilized Medical Devices
EN ISO 11135	2014	Sterilization of health-care products - Ethylene oxide - Requirements for the development, validation and routine control of a sterilization process for medical devices (ISO 11135:2014)
ANSI/AAMI ST72	2011/(R)2016	Bacterial endotoxins—Test methods, routine monitoring, and alternatives to batch testing
EN ISO 11737-1	2006/AC:2009	Sterilization of Medical Devices – Microbiological Methods – Part 1: Determination of a Population of Microorganisms on Products (ISO 11737-1: 2006/ Corr 1:2007)
EN ISO 11737-2	2009	Sterilization of Medical Devices – Microbiological Methods – Part 2: Tests of Sterility Performed in the Definition, Validation and Maintenance of a Sterilization Process (ISO 11737-2:2009)
ISO 14161	2009	Sterilization of Health Care Products – Biological Indicators – Guidance for the Selection, Use and Interpretation of Results (ISO 14161:2009)
EN ISO 10993-7	2008/AC:2009	Biological Evaluation of Medical Devices – Part 7: Ethylene Oxide Sterilization Residuals (ISO 10993-7:2008 (Corr 2009)
EN ISO 10993-1	2009/AC:2010	Biological Evaluation of Medical Devices – Part 1: Evaluation and Testing within a risk management process (ISO 10993-1:2009)
EN ISO 10993-4	2017	Biological evaluation of medical devices – Part 4: Selection of tests for interactions with blood (ISO 10993-4:2017)
EN ISO 10993-5	2009	Biological evaluation of medical devices – Part 5: Tests for in vitro cytotoxicity (ISO 10993-5:2009)
EN ISO 10993-10	2013	Biological evaluation of medical devices – Part 10: Tests for irritation and skin sensitization (ISO 10993-10:2010)
EN ISO 10993-11	2018	Biological evaluation of medical devices Part 11: Test for systemic toxicity (ISO 10993-11:2017)
EN ISO 10993-12	2012	Biological evaluation of medical devices – Part 12: Sample preparation and reference materials (ISO 10993-12:2012)
EN ISO 10993-17	2009	Biological evaluation of medical devices – Part 17: Establishment of allowable limits for leachable substances (ISO 10993-17:2002)
EN ISO 10993-18	2009	Biological evaluation of medical devices – Part 18: Chemical Characterization of materials (ISO 10993-18:2005)

Available upon request

DISCLAIMER OF WARRANTY

Shelf Life Per the Label - 2 Years

Disclaimer of Warranty Per the Instructions for Use:

ALTHOUGH THE GUIDE EXTENSION CATHETER, HEREAFTER REFERRED TO AS "PRODUCT," HAS BEEN MANUFACTURED UNDER CAREFULLY CONTROLLED CONDITIONS, MEDTRONIC AND THEIR RESPECTIVE AFFILIATES (COLLECTIVELY, "MEDTRONIC") HAVE NO CONTROL OVER CONDITIONS UNDER WHICH THIS PRODUCT IS USED. MEDTRONIC, THEREFORE, DISCLAIMS ALL WARRANTIES, BOTH EXPRESSED AND IMPLIED, WITH RESPECT TO THE PRODUCT, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.

MEDTRONIC SHALL NOT BE LIABLE TO ANY PERSON OR ENTITY FOR ANY MEDICAL EXPENSES OR ANY DIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, FAILURE OR MALFUNCTION OF THE PRODUCT, WHETHER A CLAIM FOR SUCH DAMAGES IS BASED UPON WARRANTY, CONTRACT, TORT OR OTHERWISE. NO PERSON HAS ANY AUTHORITY TO BIND MEDTRONIC TO ANY REPRESENTATION OR WARRANTY WITH RESPECT TO THE PRODUCT.

The exclusions and limitations set out above are not intended to and should not be construed so as to contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held to be illegal, unenforceable or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of this Disclaimer of Warranty shall not be affected.

CUSTOMER SERVICE

Contact information:

Contact for all questions, including complaints, repair and service of medical devices the Authorized Representative in Russia:

Medtronic, LLC,
123112, Moscow, Presnenskaya nab., 10, floor 9, office III, room 41
Tel.: (495) 580-73-77, Fax: (495) 580-73-78,

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