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
USER MANUAL

2013.12.03

Rm. 515.516, Ace High-End Tower,
235-2, Guro-Dong, Guro-Gu, Seoul, Korea

SEOIL PACIFIC CORP.

CONTENTS

1. DEVICE INTRODUCTION	2
1.1 DEVICE NAME	2
1.2 INTENDED PURPOSE	2
1.3 OVERVIEW	2
1.4 APPLICABLE DISEASE	2
2. IMPORTANT SAFEGUARDS	3
2.1 WARNINGS/CAUTIONS	3
3. GENERAL DESCRIPTION	5
3.1 FRONT VIEW	5
3.2 REAR VIEW	6
3.3 PATIENT PORT	6
4. ACCESSORIES	7
5. OPERATING INSTRUCTION	8
5.1 LCD SCREEN (AUTOMATIC MODE)	8
5.2 LCD SCREEN (MANUAL MODE)	9
5.3 CONTROL PANEL	10
6. HOW TO USE THE DEVICE	11
6.1 PREPARATION BEFORE USING	11
6.2 INITIAL SETTING	11
6.3 SETTING & USING THE AUTOMATIC MODE	12
6.4 SETTING & USING THE MANUAL MODE	13
6.5 SETTING THE CONTROL LOCK	14
6.6 RUNNING THE DEVICE	14
6.7 OPERATION VERIFICATION	14
7. CLEANING & DISINFECTION	15
7.1 EXTERNAL HOUSING	15
7.2 PATIENT CIRCUIT	15
7.2.1 INSTITUTIONAL (HOSPITAL) USE	15
7.2.2 HOME (INDIVIDUAL) USE	15
8. STORAGING	16
8.1 HOW TO STORE	16
8.2 HOW TO CLEAN	16
8.3 HOW TO MOVE	16
8.4 HOW TO DISPOSE OF	16
9. SPECIFICATIONS	17
10. SYMBOLS	20
11. TROUBLE SHOOTING & DEVICE A/S	21
12. LIMITED WARRANTY	22
13. QUALITY WARRANTY CERTIFICATIONS	23

1. DEVICE INTRODUCTION

 **Caution :** Federal law restricts this device is to sale by or on the order of a physician

1.1 DEVICE NAME

- Device Name : ComfortCough® (Mechanical In-Exsufflator)
- Model Name : SICC2001

1.2 INTENDED PURPOSE

For use on patients unable to cough or clear secretions effectively due to reduced peak cough expiratory flow, resulting from high spinal cord injuries, neuromuscular deficits or severe fatigue associated with intrinsic lung disease. It may be used either with a face mask, mouthpiece, or an adapter to patient's endotracheal tube or tracheostomy tube.

- Clinical Settings : For use in a hospital, institution setting or in a home, given adequate training and physician's prescription.
- Patient Population : For use on adults or pediatric patients.
- Contraindications : Any patient with a history of bullous emphysema, known susceptibility to pneumothorax or pneumomediastinum, or known to have had any recent barotrauma, should be carefully considered before use.

1.3 OVERVIEW

The ComfortCough® mechanical In-Exsufflator(MI-E) noninvasively helps to clear bronchial secretions when patients can't cough for themselves. The patients with pulmonary disorders associated with impaired secretion clearance function might benefit from the use of the ComfortCough®. The ComfortCough® helps patients to clear bronchial secretions by gradually applying the positive pressure to the airway followed by a rapid shift to negative pressure. These rapid changes in pressure during different phases help to make the cough stronger and more effective. This user manual includes information on how the device operates, as well as information on how to set up the device to optimize treatment. This user manual is supplied with the device. Please ensure the patient has this user manual.

1.4 APPLICABLE DISEASE

- | | |
|---------------------------------|---------------------------|
| • Muscular Dystrophy | • Spinal Muscular Atrophy |
| • Amyotrophic Lateral Sclerosis | • Myasthenia Gravis |
| • Guillain-Barre Syndrome | • Multiple Sclerosis |
| • Spinal Cord Injury | • Emphysema |
| • Cystic Fibrosis | • Bronchiectasis |

or any other neurological disorder with some paralysis of the respiratory muscles. This device can be used versatily as non-invasively with face mask or mouthpiece or invasively with endotracheal tube or tracheostomy tube.

 **Caution :** Please carefully read this User Manual before using to achieve the required result safely.

2. IMPORTANT SAFEGUARDS

Throughout this guide the following definitions apply.

2.1 WARNINGS/CAUTIONS

⚠ WARNINGS : A condition that could cause electrocution or injury to a user or operator if instructions are not followed.

- Always check time and pressure setting before each treatment.
- Always use a new circuit when using the device on a new patient.
- Patients known to have cardiac instability should be monitored for pulse and oxygen saturation very closely.
- The entire manual should be read before using the device.
- Soreness and/or pain in the chest from a pulled muscle may occur in patients using ComfortCough® for the first time if the positive pressure exceeds pressures, which the patient normally receives during Positive Pressure Therapy. Such patients should start at a lower positive pressure during treatment, and gradually (over several days, or as tolerated) increase the positive pressure.
Positive Pressure Therapy includes the use of a volume ventilator, nasal or mask ventilation or CPAP (Continuous Positive Airway Pressure), or IPPB (Intermittent Positive Pressure Breathing).
- Use the device only for its intended use as described in this manual.
- Do not use in the presence of flammable anesthetics.
- Connection should be made to a grounded outlet only.
- Do not place or store the device where it can fall or be pulled into a tub or sink.
- If the device comes into contact with water, unplug the unit.
- Never operate ComfortCough® if it has a damaged cord or plug, is not working properly, or has been dropped, damaged or immersed in water.
- The device is recommended to use the AC power cord supplied with the unit. If a replacement power cord is required, contact Seoil Pacific Service Centre.
- Any patient with a history of pneumothorax or mediastinal emphysema, or has had any recent barotrauma, should be carefully considered before using the device.
- Before connecting the patient to the device, check the whole of the tubing to eliminate any risk of leaks or defective parts.
- Do not use when there are any damages on the outside of the device or when the sign of malfunction is showing on the monitor.
- As a general safety precaution, the unit should be turned off and AC power should be disconnected before performing any cleaning procedures.
- The patient circuit is for single use and must not be re-used on other patients. This is to avoid the risk of cross-infection.

CAUTION

- U.S. Federal law restricts this device to sale by or on the order of a physician.
- Make sure that the Cooling Fan, Cooling Air Inlet, Patient Air Inlet are not to be blocked.
- Check if the device is properly equipped with the bacterial/viral filter before using.
- This device is designed for intermittent operation only and not for continuous use. The device should not be cycled continuously for more than 5 minutes. After such time, the unit should either be turned off or left idling with the blower on for at least 5 minutes.
- Do not use in the area with strong electromagnetic waves.
- Do not operate the device with wet hand.
- Turn the device off when not in use.
- Keep the cord away from heated surfaces.
- The device may be affected by electromagnetic fields greater than 10V/M.
- The device is not protected against water penetration (IPX0). During transportation, we recommend to take the device in water proof condition, to provide maximum protection against accidental splashes.
- The device must not be exposed to direct sunlight.
- The device should only be used by trained personnel.
- Repairing or checking-up should be carried only by authorized person or dealer who were educated by Seoil Pacific Corp.


3. GENERAL DESCRIPTION

3.1 FRONT VIEW

There are major components inside the device.

On the front side of the device, there is the Screen Panel, the Control Panel, the Patient Port.


The Control panel has Inhale/Exhale/Pause Pressure and Time Setting Functions and Manual Control Lever.

Part	Function
Screen Panel	shows all the set value and operating status.
Control Panel	sets in-exhale pressure, time duration, manual/automatic mode.
Patient Port	carries in-exhale pressure to the patients.

3.2 REAR VIEW


There is the AC Power Inlet, the Power Switch, the Cooling Fan on the back

Parts	Function
Air Vent	takes the air to cool the inside parts of the device.
AC Power Inlet	Operator can connect the main power supply to this inlet. Power on/off switch right beside this power inlet.
Handle	is to hold the device when operator moves the device.
Patient Air Inlet	takes the air and deliver it to the patient after filtering inside the device.

3.3 PATIENT PORT

Operator can use different size circuit as per the patients.







4. ACCESSORIES

ComfortCough® provides replacement accessories,

⚠ CAUTION : ComfortCough®, including the associated tubing, is for single-patient use and must not be re-used on other persons. This is to avoid the risk of cross-infection.

⚠ CAUTION : If individual uses, these parts must be completely dried before re-use.


⚠ CAUTION : Do not try to wash the Bacterial/Viral filter.

Parts Name	Part No.	Description		Reference
COMFORT COUGH™ Patient Circuit		Breathing Hose 1ea, Bacterial/Viral Filter 1ea, Face Mask(Adult) 1ea, Adapter 1ea		
Breathing Hose	38017	Flexible, corrugated tube with inside diameter 22mm and the total length of 1.2meters.		
Bacterial/Viral Filter	70015	Blocking the outside materials to avoid the risk of cross-infection.		
Face Mask	5353	Child(Optional)	O.D : 22mm	
	5354	Adult(Optional)		
	5355	Adult Large		
Adapter	74505	O.D : 22mm / I.D : 15mm		
Power Cable				

5. OPERATING INSTRUCTION

5.1 LCD SCREEN (AUTOMATIC MODE)


The LCD screen displays the various parameters and operation status. To assist you in adjusting the device, the LCD is equipped with a backlight. The LCD backlight comes on when the device is turned on and goes out when the device is not in use more than 3 minutes. The backlight is on again when you touch any key if the device is not turned off.

No.	Symbol or Text	Function	Function Description
1	AUTO/MANUAL	Current Mode	Indicates the current cycling mechanism.
2		Control Lock	Indicates the setting keys are unlocked.
3	PRESSURE BAR	Pressure Bar	Displays the level of generated inhale and exhale pressures.
4	XXMIN XXSEC	Operating Time	Displays the operating time.
5	IP	Inhale Pressure	Indicates the setting value of inhale pressure.
6	EP	Exhale Pressure	Indicates the setting value of exhale pressure.
7	HIGH/LOW	Inhale Flow	Indicates immediate Inhale Flow or gradual Inhale Flow.
8	TI	Inhale Time	Indicates the setting value of inhale time.
9	TE	Exhale Time	Indicates the setting value of exhale time.
10	TP	Pause Time	Indicates the setting value of pause time.


5.2 LCD SCREEN (MANUAL MODE)

The LCD screen displays the various parameters and operation status. To assist you in adjusting the device, the LCD is equipped with a backlight. The LCD backlight comes on when the device is turned on and goes out when the device is not in use more than 3 minutes. The backlight is on again when you touch any key if the device is not turned off.

No.	Symbol or Text	Function	Function Description
1	MANUAL	Current Mode	Indicates the current cycling mechanism.
2	🔒	Control Lock	Indicates the setting keys are locked.
3	PRESSURE BAR	Pressure Bar	Displays the level of generated inhale and exhale pressures.
4	STOP	Current Status	Indicates the current stopping device.
5	IP	Inhale Pressure	Indicates the setting value of inhale pressure.
6	EP	Exhale Pressure	Indicates the setting value of exhale pressure.
7	HIGH/LOW	Inhale Flow	Indicates immediate Inhale Flow or gradual Inhale Flow.

5.3 CONTROL PANEL

No.	Button	Function	Function Description
1		Start/Stop	Starts and stops the device. (In auto/manual mode)
2		Inhale pressure	-Lets you enter edit mode for adjusting inhalation pressure. -Options: 0 ~ 60 cmH ₂ O ; 1.0cmH ₂ O increment
3		Exhale pressure	-Lets you enter edit mode for adjusting exhalation pressure. -Options: 0 ~ 60cmH ₂ O ; 1.0cmH ₂ O increment
4		Inhale Time	-Lets you enter edit mode for adjusting time duration of inhale phase.(In automatic mode only) -Options: 0.5 ~ 5.0sec ; 0.1sec increment
5		Exhale Time	-Lets you enter edit mode for adjusting time duration of exhale phase.(In automatic mode only) -Options: 0.5 ~ 5.0sec ; 0.1sec increment
6		Pause Time	-Lets you enter edit mode for adjusting time duration of pause phase.(In automatic mode only) -Options: 0.0 ~ 5.0sec ; 0.1sec increment
7		Up/Down	-Adjusting numerical values of pressure and time.(Down to decrease up to increase) -‘Control Lock’ is activated when pushing ‘Up’ and ‘Down’ keys simultaneously for at least three seconds. (only in “Stop”status)
8		Auto/Manual	Changes the operation mode as Auto / Manual.
9		Inhale flow	-Sets the operation mode as immediate Inhale Flow(High) or gradual Inhale Flow(Low). -Options: High, Low
10		Manual control lever	-Manually changes the unit to inhale or exhale. (In manual mode only)

6. HOW TO USE THE DEVICE

- ⚠ **WARNING** : Do not use for more than five minutes without any prescription.
Patient injury like Hyperventilation could occur.

6.1 PREPARATION BEFORE USING

1. Please check the rated voltage of the power supply.
2. Please check if the power supply cord is correctly connected.
3. Combined use with other devices without prior approval of Seoil Pacific Corp. Can cause malfunction or unwanted results.
4. Do not operate this device if unskilled.
5. Please read this User Manual carefully before operating.
6. If more than 2 patients are using this device together, please change patient's circuit(tubing, bacterial/viral filter, face mask, connector) to avoid unexpected cross contamination.

6.2 INITIAL SETTING

1. Position the unit on a stable and hand surface within easy reach of the patient, or the operator of the unit.

- ⚠ **CAUTION** : Position the device where the air intake ports on the bottom of the unit are not blocked.

2. Connect the power cord connector to the socket on the rear of the device.
3. Plug the power cord into an AC grounded outlet of appropriate outlet.
4. Assemble the patient circuit (filter, breathing hose and patient interface) as follows ;
 - a. Attach the bacterial/viral filter to the patient port on the front panel.
 - b. Attach the 1.2meter x 22mm smooth bore breathing hose to the bacterial/viral filter.

- ⚠ **WARNING** : Air tube from Seoil Pacific should be used with the device. A different type of air tube may alter the pressure the patients actually receive and reduce the effectiveness of this device.

- ⚠ **WARNING** : Do not use conductive or anti-static hoses or tubes.

- c. Attach the appropriate patient interface to the breathing hose. Patient interface options include a face mask and adapter, mouthpiece, lip seal or tracheostomy tube adapter (A face mask and adapter are included in each unit)

5. Turn on the power switch. Power can be turned on/off by pressing the Power Switch. The 'I' Symbol on the power switch designates the 'ON' position. The 'O' symbol on the power switch designates the 'OFF' position.

- ⚠ **CAUTION** : Be careful not to place the device where it can be bumped or where someone is likely to trip over the power cord.


- ⚠ **CAUTION** : Make sure the area around the device is dry and clean. It should also be clear of bedding, clothes and other potential blockages and free of any objects that could block the dust filter of device's bottom.

- ⚠ **WARNING** : Make sure the power cord and plug are in good condition and the equipment is not damaged.

- ⚠ **WARNING** : The air filter cover protects the device in the event of accidental liquid spillage onto the device. Ensure that the air filter and air filter cover are fit at all times.

6. The device performs pre-test by itself before care-giver operates the device.
7. To set the pressure and time, press the 'START/STOP' button.
8. Block the end of patient port and wait until you have the desired pressure through pressure bar.

6.3 SETTING & USING THE AUTOMATIC MODE


1. Press 'AUTO/MANUAL' button to set the automatic mode.
Automatic mode consists of an inhalation phase, an exhalation phase and a pause phase, after which inhalation begins again. The time for each phase is set with the three buttons on the control panel. Normally, inhale time and exhale time are set to 1 ~ 3 seconds and the pause time can be set up to 5 seconds, or eliminated by setting the pause time button to 0 seconds, depending on the patient's preference.
2. Set the Pressure Value.
 - a. Press the 'IP' button to edit the 'IP' field.
Press the '▲ (UP)/ ▼ (DOWN)' button until the desired setting value appears.
 - b. Press the 'EP' button to edit the 'EP' field.
Press the '▲ (UP)/ ▼ (DOWN)' button until the desired setting value appears.
※Notice! Maximum value of Inhale Pressure cannot exceed set value of Exhale Pressure.
3. To change the inhalation flow, Press the 'HIGH/LOW' button to adjust the inhale flows.
 - a. HIGH : Full inhalation flow
 - b. LOW : Gradual inhalation flow (Usually half of Full Inhalation Flow)
※Notice! When using gradual inhalation flow (LOW) there might be a small reduction in inspiratory pressure.
4. Set the inhalation/exhalation/pause time.
 - a. Press the 'TI' button to edit the 'TI' field.
Press the '▲ (UP)/ ▼ (DOWN)' button until the desired setting value appears.
 - b. Press the 'TE' button to edit the 'TE' field.
Press the '▲ (UP)/ ▼ (DOWN)' button until the desired setting value appears.
 - c. Press the 'TP' button to edit the 'TP' field.
Press the '▲ (UP)/ ▼ (DOWN)' button until the desired setting value appears.
5. Press the 'START/STOP' button to start the device.
6. Block the end of patient port and wait until you have the desired pressure through pressure bar.
7. After 4 to 5 breathings, remove the patient interface from the patient and allow time for a normal breathing pattern to return (20 to 60 seconds), or place the patient back on the ventilator if currently in use. Avoid prolonged periods of connection to the device. During this resting period, clear secretions that may have become visible in the mouth, throat or tracheostomy tube.

⚠ **WARNING** : Soreness and/or pain in the chest from a pulled muscle may occur in patients using COMFORT COUGH™ for the first time if the positive pressure used exceeds pressures which the patient normally receives during Positive Pressure Therapy. Such patients should start at a lower positive pressure during treatment, and gradually (over several days, or as tolerated) increase the positive pressure used.

[Positive Pressure Therapy includes the use of a volume ventilator, nasal or mask ventilation or CPAP (Continuous Positive Airway Pressure), or IPPB (Intermittent Positive Pressure Breathing).]

6.4 SETTING & USING THE MANUAL MODE

1. Press 'AUTO/MANUAL' button to set the manual mode.
2. Set the Pressure Value.
 - ※Notice : Each individual patient may require special settings for the maximum positive(inhalation) and negative(exhalation) pressures. For a patient using this device for the first time, it is advisable to begin with lower pressures such as 10~15cmH₂O positive and negative, and to familiarize the patient with the feel of mechanical insufflation-exsufflation. During subsequent treatments, pressures can be increased if needed to achieve adequate secretion clearance. Note that at these lower pressures ComfortCough® may have limited effectiveness in clearing secretions. Increasing pressures should improve the effectiveness.
 - a. Press the 'IP' button to edit the 'IP' field.
Press the '▲ (UP)/ ▼ (DOWN)' button until the desired setting value appears.
 - b. Press the 'EP' button to edit the 'EP' field.
Press the '▲ (UP)/ ▼ (DOWN)' button until the desired setting value appears.
 - ※Notice : Maximum value of Inhale Pressure cannot exceed set value of Exhale Pressure.
3. To change the inhalation flow, Press the 'HIGH/LOW' button to adjust the inhale flows.
 - a. HIGH : Full inhalation flow
 - b. LOW : Gradual inhalation flow (Usually half of Full Inhalation Flow)
 - ※Notice : When using gradual inhalation flow (LOW) there might be a small reduction in inspiratory pressure.
4. Press the 'START/STOP' button to start the device.
5. Press the manual control lever to exhale position for 5 sec. (Check the bar graph.)
Press the manual control lever to inhale position for 5 sec. (Check the bar graph.)
and wait till you have the desired pressure through Pressure Bar.
 - ※Notice : If you need to adjust the pressure value, follow the above procedures again.
- ⚠ CAUTION : Begin with lower pressures like 10 cmH₂O so that patient can accommodate to the device.
6. Attach the appropriate patient interface to patient
7. Shift the manual control lever to the inhale position(to the right)and hold for 2 ~ 3 seconds.
Then observe the pressure bar to see the pressure build slowly.
8. Rapidly shift the manual control lever to the exhale position(to the left) to induce the cough, holding it there for 1~2 seconds.
9. The lever can then either be left in the neutral position for a few seconds, or can be shifted immediately to the positive pressure phase for another cough cycle, depending on the patient's preference.
10. After 4 to 5 breathings, remove the patient interface from the patient and allow time for a normal breathing pattern to return (20 to 60 seconds), or place the patient back on the ventilator if currently in use. Avoid prolonged periods of connection to the device. During this resting period, clear secretions that may have become visible in the mouth, throat or tracheostomy tube.

6.5 SETTING THE CONTROL LOCK

To prevent unintentional or accidental touch to the keys, you can lock the control keys by yourself.

1. Push the '▲ (UP)/ ▼ (DOWN)' keys simultaneously for at least three seconds.
2. The padlock symbol  is displayed on the screen.
3. To unlock the setting keys, push the '▲ (UP)/ ▼ (DOWN)' keys simultaneously for at least three seconds.
4. The open padlock symbol  is displayed on the screen.

※Notice! Power ON/OFF and 'START/STOP' button will not be locked at anytime.

6.6 RUNNING THE DEVICE

 **WARNING** : Always check time and pressure setting before each treatment.

Treatment usually consists of 4 or 5 coughings in succession. The patient is then allowed to rest for 20 ~ 60 seconds, which helps avoid hyperventilation. The cycles can then be repeated 4 ~ 6 times for a full treatment.

6.7 OPERATION VERIFICATION

It is recommended that ComfortCough® is to be periodically tested to ensure that the cycling valve returns to the neutral or pause position after either the inhale or exhale phase. To determine this, follow these steps.

1. Turn on the power switch.
2. Attach the patient circuit to the device and block the end of the patient port.
3. Press the 'AUTO/MANUAL' button to set the manual mode.
4. Set the maximum pressure value.
5. Press the 'START/STOP' key to start the device.
6. Move the manual control lever from exhale to inhale (repeat a couple of times) and observe the pressure bar to ensure that positive and negative pressure is being applied to the patient circuit.
7. Release the manual control lever from the inhale position and observe that the pressure immediately drops to 0 cmH₂O. Repeat this process from the exhale position also. In either case, if the pressure does not drop to 0 cmH₂O, the device should be returned for repair.

7. CLEANING & DISINFECTION

7.1 EXTERNAL HOUSING

The exterior of the device and housing may be washed with a mild detergent and water, or with a bactericidal cleaning solution such as 70% isopropyl alcohol.

⚠ CAUTION : Do not sterilize the device with ethylene oxide gas or steam.

7.2 PATIENT CIRCUIT

⚠ WARNING : The patient circuit should not be sterilized.

Always use a new circuit when using the device on a new patient.

※Notice : If you need additional accessories, please contact the distributor in your area.

7.2.1 INSTITUTIONAL (HOSPITAL) USE

a. Breathing Hose, Patient Interface and Adapters :

If the device is to be used by more than one patient, the circuit must be replaced.

b. Bacterial/Viral Filter :

If the device is to be used by more than one patient, the filter must be replaced to prevent cross contamination. Do not try to wash the filter.

7.2.2 HOME (INDIVIDUAL) USE

a. Breathing Hose, Patient Interface and Adapters :

After use, the breathing hose and patient interface should be washed thoroughly in soap and water. These parts must be completely dried before reuse.

b. Bacterial/Viral Filter :

The filter, which protects the device from entraining foreign material to the patient, can be left in place as long as it is not blocked by sputum or trapped moisture. Do not try to wash the filter.

8. STORAGING

8.1 HOW TO STORE

1. If there are any problems regarding the device, stop using the device immediately and contact Seoil Pacific Service Center.
2. Be sure to store the device away from water source.
3. Be sure to store the device to prevent dust from outside.
4. Remove the power line from the outlet after use.
5. Do not install the device on a slope or on an unstable ground.
(Be sure to install the device on horizontal surface).
6. Be sure to turn off the power when not in use.
7. Refrain from having direct ray of sunlight on the device.
8. Be sure to clean the device after use for storage.
9. Store it away from children's reach.
10. Do not remove the cover.
11. When using the device after a long storage period, make sure that the device functions properly before using the device.
12. Free of any objects that could block the dust filter on the bottom of device
13. Do not store the device nearby heating places.
14. Do not store the device in extremely high or low temperatures.
 - Operating temperature range : 10°C~35°C
 - Operating humidity range : 30%~75%, non-condensation
 - Operating pressure : 700-1060hPa
 - Storage temperature range : -20°C~50°C
 - Storage humidity range : 10%~90%, non-condensation

8.2 HOW TO CLEAN

1. Be sure to pull the plug out before cleaning the device.
2. Do not apply shock to the device.
3. Be sure to clean the device with dried towel.
4. Do not place drinks around the device.
5. Be sure to dry the device after cleaning

8.3 HOW TO MOVE

1. Be sure to pull the plug out before moving the device.
2. Pack the device if moving the device long distance.
3. Do not apply severe shock to the device during transportation. (It may lead to malfunction)

8.4 HOW TO DISPOSE OF

The device can be reused or recycled in accordance with WEEE/RoHs directive, and you should use appropriate collection, reuse and recycling systems available in your region. If you require information on collection and disposal of your device, please contact Seoil Pacific service center or local distributor at your area.

*WEEE(2002/96/EC) : European Directive that requires the proper disposal of electrical and electronic equipment.

9. SPECIFICATIONS

Device Name	ComfortCough® (Mechanical In-Exsufflator)
Model Name	SICC2001
Intended Purpose	Mechanical In-Exsufflator(MI-E) noninvasively helps to clear bronchial secretions when patients can't cough for themselves.
How to Use	Read the accompanying manuals.
Input Voltage	100-240VAC, 50/60Hz
Input Power	230VA
Electrical Protection	BF, Class I
Dimensions	288 x 330 x 212 (mm)
Weight	1set, 6.4Kg
Maximum positive pressure	60CmH ₂ O
Maximum negative pressure	-60CmH ₂ O
Maximum Inhalation Flow	3.3 liters/second with inhale flow set to minimum; if set to maximum inhalation, the flow is the same as the exhalation flow..
Maximum Exhalation Flow Capacity	10 liters/second; actual flow depends on maximum pressure and on airway resistance.
Mode of Operation	Automatic and Manual Timing.
Inhalation, Exhalation Times	0.5 to 5.0 seconds
Pause Times	0.0 to 5.0 seconds
Pressure Range	+5~+60cmH ₂ O(3cmH ₂ O or within 5%) -5~-60cmH ₂ O(3cmH ₂ O or within 5%)
Blower Type	Two-stage centrifugal blower with BLDC Motor
Enclosure Degree of Ingress Protection	IPX0
Mode of Operation	Intermittent operation only

Guidance and manufacturer's declaration – electromagnetic emissions

The Model SICC2001 is intended for use in the electromagnetic environment specified below. The customer or the user of the Model SICC2001 should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The model SICC2001 uses RF energy only for its internal functions. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The model SICC2001 is suitable for use in all establishments including domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.

Guidance and manufacturer's declaration – electromagnetic immunity

The model SICC2001 is intended for use in the electromagnetic environment specified below. The customer or the user of the model SICC2001 should assure that it is used in such an environment.

Immunity test	IEC60601 test level	Compliance level	Electromagnetic environment guidance
Electrostatic discharge (ESD) IEC 61000-4-2	6 kV contact 8 kV air	6 kV contact 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	2 kV for power supply lines 1 kV for input/output lines	2 kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	1 kV differential mode 2 kV common mode	1 kV differential mode 2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruption, and voltage variations on power supply input lines IEC 60601-4-11	< 5 % U_T (> 95 % dip in U_T) for 0.5 cycle 40 % U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) < 5 % U_T (> 95 % dip in U_T) for 5 sec	< 5 % U_T (> 95 % dip in U_T) for 0.5 cycle 40 % U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) < 5 % U_T (> 95 % dip in U_T) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the model SICC2001 requires continued operation during power mains interruptions, it is recommended that the model SICC2001 be powered from an uninterruptible power supply or battery.
Power frequency (50/60 Hz) IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Note : U is the A.C. mains voltage prior to application of the test level.

Guidance and manufacturer's declaration – electromagnetic immunity

The model SICC2001 is intended for use in the electromagnetic environment specified below.
The customer or the user of model SICC2001 should assure that it is used in such an environment.

Immunity test	IEC60601 test level	Compliance level	Electromagnetic environment guidance
Conducted RF IEC61000-4-6	3 Vrms 150 kHz to 80MHz	3 Vrms 150 kHz to 80MHz	Portable and mobile RF communications equipment should be used no closer to any part of the model SICC2001, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = [3.5/V1] P$ $d = [3.5/E1] P$ 80MHz to 800MHz $d = [7/E1] P$ 800MHz to 2.5GHz Where P is the maximum output power rating of the transmitter in watts(W) according to the transmitter manufacturer and d is the recommended separation distance in meters(m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC61000-4-3	3 V/m 80 MHz to 2.5GHz	3 V/m 80 MHz to 2.5GHz	

Note 1 At 80MHz and 800MHz, the higher frequency range applies.

Note 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the model ROOTOR is used exceeds the applicable RF compliance level above, the model ROOTOR should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the model ROOTOR.

^b Over the frequency range 150kHz to 80MHz, field strengths should be less than $[V1] V/m$.

10. SYMBOLS

The labels or labeling for this device may have the following symbols which have the meaning presented below.

Symbol	Description
	Mains On
	Mains Off
	Attention, consult instruction manual.
	Type BF Equipment
	Read the accompanying manuals
	IPX0
	AC voltage
	Manufacturer
	Date of Manufacture
	Authorized representative in the European community
	Notified body of CE/MDD(SGS)
	Locate next to the protective earth terminal inside the device, signifies that the device is grounded.
	Device subject to separate collection.
	Serial Number
	Caution : U.S.Federal law restricts this device to sale by or on the order of a physician.

11. TROUBLE SHOOTING & DEVICE A/S

Please check the following categories before applying for after service.

Status	Handling method
The power is not turned on even if I press the power switch.	-Please check the connection status of the power plug. -Check to see if the main switch is turned on. -Check to see if the fan is operating.
When the air is not coming out.	-Check the pressure setting.
When the air is weak.	-Check the pressure setting. -Check the filter and change if necessary.
When 'TEMP ERR' appears on the screen.	-Blower temperature is too high. -Device will stop automatically, and TEMP ERR will be disappeared when the temperature is in appropriate condition. If TEMP ERR is not shown, you can use this device again.
When 'E1' appears on the screen.	-The position of stepping motor which controls the valve is not correct. -Device will reset by itself and will start after checking automatically. -Call service center if still E1 shows on the monitor even you have restarted the device.
When 'E2' appears on the screen.	-The pressure exceeds the set value. -The device stop automatically. -Call service center if still E2 shows on the monitor even you have restarted the device.

- When you suspect any malfunctions, please stop using the device, turn off the main switch and pull the power plugs out. Then, contact our Seoil Pacific Service Center for repairs.
- If used in normal usage setting, our device is fully guaranteed for 1 year from the purchasing date.

12. LIMITED WARRANTY

Seoil Pacific Corp. warrants that your COMFORT COUGH™ shall be free from defects in material and workmanship for the period specified below. This warranty is not transferable.

Devices	Warranty Period
Seoil Pacific COMFORT COUGH™	1 Year
Accessories and mask systems (including mask frame, cushion, headgear and tubing). Excludes single-use devices.	90 Days

If the device fails under conditions of normal use, Seoil Pacific Corp. will repair or replace, at its option, the defective device or any of its components. This Limited Warranty does not cover:

- a) Any damage caused as a result of improper use, abuse, modification or alteration of the device;
- b) Repairs carried out by any service organization that has not been expressly authorized by Seoil Pacific Corp.;
- c) Any damage or contamination due to cigarette, pipe, cigar or other smoke;
- d) Any damage caused by water being spilled on or into the device.

Warranty is void on the device sold, or resold, outside the region of original purchase. Warranty claims on defective device must be made by the initial consumer at the point of purchase.

This warranty is in lieu of all other express or implied warranties, including any implied warranty of merchantability or fitness for a particular purpose. Some regions or states do not allow limitations on how long an implied warranty lasts, so the above limitation may not apply to you.

Seoil Pacific Corp. shall not be responsible for any incidental or consequential damages claimed to have occurred as a result of the sale, installation or use of any Seoil Pacific Corp. device. Some regions or states do not allow the exclusion or limitation of incidental or consequential damages, so the above limitation may not apply to you. This warranty gives you specific legal rights, and you may also have other rights which vary from region to region.

However, this warranty is not subject to cover the damages by accident, incorrect use, over use, transform, or defects on the parts which are not affecting the function of the device.

For further information on your warranty rights, contact your local Seoil Pacific Corp. dealers or authorized representatives.

13. QUALITY WARRANTY CERTIFICATIONS

Device number (S/N)	
Purchase date (dd/mm/yy)	
Free repair warranty period	1 year from purchase date
Warranty area	
Name	
Purchaser Address	□□□-□□□
Telephone	TEL : () - / H.P : - -
Branch name	
Seller Address	□□□-□□□
Telephone	

This certificate promises free repairs based on its terms.

When any malfunctions occur during the warranty period, please submit this warranty certificate with the device during warranty period.

-If the above blanks are left blank, this certificate loses its power. If they are not filled out, please contact your sales manager immediately.

-This certificate may not be reissued.

-Our device is manufactured under a strict inspection process.

< Summary of warranty >

- | | |
|--|---|
| <p>1. If any malfunction occurs during normal use condition mentioned in this manual and device label, customers are entitled to free repairs based on customer protection regulations.</p> <p>2. When receiving repair service during the warranty period, please submit this certificate with the device. Repair costs will be charged when the repair occurs after the warranty period is over.</p> <p>3. When moving from your house, please contact your sales manager or Seoil Pacific Service Centre.</p> <p>4. If the following incidents occur, customers will be charged with repair cost even during the warranty period.</p> <p>(1) When malfunction occurs due to user's misuse, unauthorized repairs and device adjustments.</p> <p>(2) When the device is damaged during transportation or by any means of shock.</p> | <p>(3) When the device is used in abnormal environment that is not covered within the manual.</p> <p>(4) When malfunction occurs due to foreign objects such as drinks, coffee and other objects.</p> <p>(5) When malfunction occurs due to combination use with other devices.</p> <p>(6) When the damaged parts are consumption parts.</p> <p>(7) When the device is damaged by force majeure such as fire accidents, earthquake, water accidents and lightning.</p> <p>(8) When this certificate is not presented.</p> <p>(9) When date of purchase, customer name and sales branch have been tampered.</p> <p>(10) When it has been damaged due to improper use (used in shops instead of at homes).</p> <p>5. Please keep this certificate for better A/S.</p> |
|--|---|

※This certificate promises free repairs based on mentioned period and terms. Therefore, this certificate does not limit customer's legal rights. If you have any questions about service and repairs after the warranty period, please contact Seoil Pacific Service Centre or your sales manager.

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