KOMMA

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Product no 211

Carton: C211EU-02

Label: L211EU-02

Sample Carton: SC211EU-02

GLOBAL ENGLISH MASTER TEXT

*****Artelac ® Night*****

Product number: 211

Format: single-dose units

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Legal Manufacturer: *as indicated on valid DoC*

Product # 211
Leaflet: I211EU-02

Carton: C211EU-02

Label: L211EU-02
20/08/2020

IFU Master Text

Please read these instructions for use carefully – they contain important information on how to use this product. If you have any questions, please ask your pharmacist, your doctor or contact us.

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(logotype) ***** Artelac ® Night EDO *****

0.24% hyaluronic acid (as sodium hyaluronate), carbomer, glycerol and a lipid component (medium chain triglycerides). Does not contain a medicinal product in its composition

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What is ***** Artelac ® Night EDO *****?

***** Artelac ® Night EDO ***** is an opaque ophthalmic solution that contains 0.24% hyaluronic acid (as sodium hyaluronate), carbomer, glycerol and a lipid component (medium chain triglycerides).

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What is the application of ***** Artelac ® Night EDO *****?

***** Artelac ® Night EDO ***** is an eye lubricant that supports all three layers of the tear film to improve the moistening of the ocular surface in the presence of dry eye sensation.

When you should use ***** Artelac ® Night EDO *****

***** Artelac ® Night EDO ***** is inspired by the design of a natural tear film. It contains the natural moisturizer hyaluronic acid (as sodium hyaluronate), which is found in the structures of a healthy human eye and the natural tear film and two additional moisturizers: Glycerol and Carbomer. Containing an additional lipid-component (medium chain triglycerides) ***** Artelac ® Night EDO ***** is supporting all three layers of the tear film: lipid, aqueous and mucin phase.

Product # 211

Leaflet: I211EU-02

Carton: C211EU-02

Label: L211EU-02

20/08/2020

***** Artelac ® Night EDO ***** is designed to remain on the surface of the eye providing moisturization for immediate, long-lasting symptom relief. It is an ideal product for everybody suffering from a range of symptoms of dry eyes such as tired, stressed, gritty, burning and/or tearing eyes. Dry eyes are caused by


- external factors, like extensive computer use, air conditioning, heating, cigarettes smoke, as well as windy outdoor environments as well as wearing contact lenses.
- Other factors, like certain drugs (such as the contraceptive pill, antihistamines, betablockers, antispasmodics and diuretics) or certain diseases that play a role in the development of ocular dryness.

How often and how long should you use * Artelac ® Night EDO ***?**

Instill one drop of ***** Artelac ® Night EDO ***** into the eye as the need arises. You can use it every day and as often as you want. There is no limit of duration of use.

Discard the single dose unit after use.


How to use * Artelac ® Night EDO ***?**

	1. Separate one single-dose unit from the strip carefully
	2. Twist (don't pull) the top cap off the single-dose unit.

Product # 211
Leaflet: I211EU-02

Carton: C211EU-02

Label: L211EU-02
20/08/2020

3. Lean your head back slightly and gently pull your lower eyelid down with a finger of your free hand. Hold the single-dose unit as vertically as possible above your eye with the nozzle pointing downwards. Make sure that the single-dose unit does not touch your eye to prevent any possible injury or damage. Then gently squeeze the single-dose unit to apply one drop.

Close your eye afterwards and slowly roll it in order to optimally spread the drop on the surface of your eye.

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Important information about using * Artelac ® Night EDO *****

- Local hypersensitivity reactions have been very rarely reported. Treatment should discontinue in such cases.
- This product may cause temporary blurred vision after instillation which can have negative impact on driving vehicles or operating machines.
- Do not touch the eye surface with the tip of the single dose to prevent any possible injury.
- Remove contact lenses before usage.
- Wait approximately 15 minutes before using another eye product or before replacing contact lenses.
- Do not re-use.
- Do not use after the use-by date.
- Do not use if the package or single-dose unit is damaged.

How to store * Artelac ® Night EDO ***?**

Store at a temperature of +1 - +25°C.

Keep away from sunlight.

Product # 211
Leaflet: I211EU-02

Carton: C211EU-02

Label: L211EU-02
20/08/2020

Presentation:


*** Artelac ® Night EDO *** is available in packs containing 2, 10, 30, 60 single dose units with 0.5 ml eye drops.

This product is also available in preservative-free multi-dose bottles.

Artelac ® Night EDO ophthalmic solution has the following temperature conditions and allowable durations.

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Symbols

<u>Symbol</u>	<u>Title</u>
	Sterilized using aseptic processing techniques
	Temperature limit
	Do not use if package or single dose unit is damaged
	Consult instructions for use
	Use-by date
	Batch code

Product # 211
Leaflet: I211EU-02

Carton: C211EU-02

Label: L211EU-02
20/08/2020

	Do not re-use
	Keep away from sunlight
	Manufacturer

Recycling

Artelac ® Night EDO ophthalmic solution is not toxic, do not create an excessive load on a person or the environment during its operation, storage or transportation. To prevent damage to the natural environment, it is necessary to comply with safety requirements, storage, transportation and operation conditions. Artelac ® Night EDO ophthalmic solution is sent for disposal in accordance with hygiene requirements of SanPin 2.1.7.2790-10 "Sanitary and epidemiological requirements for medical waste management", adopted in the Russian Federation.

Artelac ® Night EDO ophthalmic solution that is subject to disposal in connection with the expiration of the expiration date, as well as in the event of a breach of the package integrity, is classified as hazard class A. Medical devices, after using for the intended purposes, are classified as hazard class B. After the hardware decontamination methods using physical methods and changing the physical appearance of the waste product, which excludes the possibility of their reuse, waste of Class B can be accumulated, temporarily stored, transported, destroyed and disposed of together with Class A. Packaging of decontamination medical waste of Class B must be marked with evidence of medical wastes decontamination. Disposal must be carried out by accredited organizations.

Product # 211

Leaflet: I211EU-02

Carton: C211EU-02

Label: L211EU-02

20/08/2020

More information www.bausch.ru

Authorization representative in RU: Bausch Health LLC, Shabolovka St. 31, Bld. 5,
Moscow, 115162, Russia

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Address: Brunsbütteler Damm 165/173, 13581, Germany

Insert revision number and date of IFU in the format: Product name, revision
number of master text; national revision date

® / ™ are trademarks of Bausch & Lomb Incorporated or its affiliates.

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Product number: 999

Format: multi-dose bottle

Regulatory Status: Medical Device class IIa

Notified Body: CE 0483 (MDC) *(as indicated on valid DoC)*

Legal Manufacturer: Dr. Gerhard Mann Chem.-pharm. Fabrik GmbH, Germany,
Brunsbütteler Damm 165/173, 13581 Berlin, Germany

Product no 999

Leaflet: I999EU-03

Carton: C999EU-03

Label: L999EU-03

20/08/2020

IFU Master Text

Please read these instructions for use carefully – they contain important information on how to use this product. If you have any questions, please ask your pharmacist, doctor or contact us.

(logotype) *insert LOGO if applicable*

(logotype) *** Artelac®Night ***

0.24% sodium hyaluronate, carbomer, glycerol, lipid component (medium chain triglycerides), sodium hydroxide and water for injection. Does not contain a medicinal product in its composition

Intended use: medical devices for moisturizing the surface of the eyes with a feeling of dryness in the eyes

What is * Artelac®Night ***?**

*** Artelac®Night *** is an opaque ophthalmic solution that contains 0.24% hyaluronic acid (as sodium hyaluronate), carbomer, glycerol, a lipid component (medium chain triglycerides), sodium hydroxide and water for injection.

The advanced formulation possesses a very high tolerability as it is free of preservatives and free of buffers.

The ophthalmic solution is presented in an advanced and innovative preservative free multi-dose system (PFMD), which is easier to use and safe. Once opened it can be used for 6 months.

What is the application of * Artelac®Night ***?**

*** Artelac®Night *** is an eye lubricant that supports all three layers of the tear film to improve the moistening of the ocular surface in the presence of dry eye sensation.

When you should use * Artelac®Night ***?**

*** Artelac®Night *** is— inspired by the design of the natural tear film. It contains the natural moisturizer hyaluronic acid (as sodium hyaluronate), which is found in the structures of a healthy human eye and the natural tear film and two additional

Product no 999

Leaflet: I999EU-03

Carton: C999EU-03

Label: L999EU-03

20/08/2020

moisturizers: Glycerol and Carbomer. Containing an additional lipid-component (medium chain triglycerides) *** Artelac®Night *** is supporting all three layers of the tear film: lipid, aqueous and mucin phase.

*** Artelac®Night *** is designed to remain on the surface of the eye providing intensive moisturization for immediate, long-lasting symptom relief. It is an ideal product for everybody suffering from a range of symptoms of dry eyes such as tired, stressed, gritty, burning and/or tearing eyes.

Dry eyes are caused by

- external factors, like extensive computer use, air conditioning, heating, cigarettes smoke, as well as windy outdoor environments and/or wearing contact lenses.
- other factors, like certain drugs (such as the contraceptive pill, antihistamines, betablockers, antispasmodics and diuretics) or certain diseases that play a role in the development of ocular dryness.

How often and how long should you use * Artelac®Night ***?**

Instill one drop of *** Artelac®Night *** into the eye as the need arises. You can use it every day and as often as you want. There is no limit of duration of use.

How to use * Artelac®Night ***?**

Attention! Before first use and if the bottle has not been used for a while:

Remove the cap and point the bottle downwards, press repeatedly up to a maximum of 20 times until a drop is released. Now the bottle is ready for use.




Product no 999

Leaflet: I999EU-03

Carton: C999EU-03

Label: L999EU-03

20/08/2020

 1	1. Before each use, remove the translucent white protective cap from the tip of the bottle with a rotating movement.
 2	2. Attention! Before first use and if the bottle has not been used for a while: Point the bottle with the tip downwards and press repeatedly the pump mechanism until a drop is released.
 3	3. Hold the bottle as indicated: Place your thumb on the ergonomic thumb support and enclose the lower part with your fingers. Two fingers should be placed on the bottom of the bottle.
 4	4. Point the bottle with the tip downwards. Lean back slightly your head. Gently pull your lower eye-lid down with a finger of your free hand, hold the bottle in a vertical position over your eyes as. Press gently the pump mechanism together to apply one drop into the eye. Close your eye afterwards and blink in order to optimally spread the drop on the surface of your eye.

Product no 999

Leaflet: I999EU-03

Carton: C999EU-03

Label: L999EU-03

20/08/2020

The innovative pump mechanism of *** Artelac®Night *** allows the ophthalmic solution to be preservative free and allows a comfortable instillation. The declared filling volume can be Nightly removed from the bottle. Due to the design of the bottle a small remnant of the solution remains within the bottle.

When you should not use * Artelac®Night ***?**

Do not use this product if you are hypersensitive to any of the ingredients.

Important information about using * Artelac®Night ***?**

- Local hypersensitivity reactions have been very rarely reported. Treatment should discontinue in such cases.
- This product may cause temporary blurred vision after instillation which can have negative impact on driving vehicles or operating machines.
- Do not touch the eye surface with the tip of the bottle to prevent any possible injury of the eye surface or contamination of the solution.
- Shake-off any possible residual drops after application.
- Careless use of the product contrary to the instructions for use could result in slipping a product from hand and may cause eye injury
- Remove contact lenses before usage.
- Wait approximately 15 minutes before using another eye product or before replacing contact lenses.
- Do not use after the use-by date.
- Do not use if the package or bottle is damaged.
- Keep out of the reach of children

How to store * Artelac®Night ***?**

Store at a temperature of +1 - +25°C.

The design of the bottle allows you to use *** Artelac®Night *** for 6 months after opening.

Presentation:

Product no 999

Leaflet: I999EU-03

Carton: C999EU-03

Label: L999EU-03

20/08/2020


*** Artelac®Night *** is available in packs of 10 ml ophthalmic solution. ***
Artelac®Night EDO*** is available in convenient, preservative free single dose
units (SDU) as well.

Transportation conditions:

Artelac®Night ophthalmic solution has the following temperature conditions and
allowable durations.

Temperature conditions	Relative humidity	Additional requirements
1° C to 30° C	Not more than 80%	must be protected from heat and rain (moisture)

Symbols

<u>Symbol</u>	<u>Title</u>
	Sterilized using aseptic processing techniques
	Temperature limit
	Do not use if package is damaged
	Discard the medical device after 6 months once opened.

Product no 999

Leaflet: I999EU-03

Carton: C999EU-03

Label: L999EU-03

20/08/2020

	Consult instructions for use
 (YYYY-MM)	Used-by date
	Batch code
	Manufacturer

Recycling

Artelac®Night ophthalmic solution are not toxic, do not create an excessive load on a person or the environment during its operation, storage or transportation. To prevent damage to the natural environment, it is necessary to comply with safety requirements, storage, transportation and operation conditions. Artelac®Night ophthalmic solution is sent for disposal in accordance with hygiene requirements of SanPin 2.1.7.2790-10 "Sanitary and epidemiological requirements for medical waste management", adopted in the Russian Federation.

Artelac®Night ophthalmic solution that is subject to disposal in connection with the expiration of the expiration date, as well as in the event of a breach of the package integrity, is classified as hazard class A. Medical devices, after using for the intended purposes, are classified as hazard class B. After the hardware decontamination methods using physical methods and changing the physical appearance of the waste product,

Product no 999

Leaflet: I999EU-03

Carton: C999EU-03

Label: L999EU-03

20/08/2020

which excludes the possibility of their reuse, waste of Class B can be accumulated, temporarily stored, transported, destroyed and disposed of together with Class A. Packaging of decontamination medical waste of Class B must be marked with evidence of medical wastes decontamination. Disposal must be carried out by accredited organizations.

More information www.bausch.ru

Authorization representative in RU: Bausch Health LLC, Shabolovka St. 31, Bld. 5, Moscow, 115162, Russia

Manufacturer: Dr. Gerhard Mann Chem.-pharm. Fabrik GmbH, Germany

Address: Brunsbütteler Damm 165/173, 13581, Germany


insert revision number and date of IFU in the format: Product name, revision number of master text; national revision date

® / ™ are trademarks of Bausch & Lomb Incorporated or its affiliates.
© Bausch & Lomb Incorporated

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