
Eurotrol

Instruction for use

Eurotrol Plasma/Low Hb QCs for Low Hemoglobin measuring in variations

Manufactured by Eurotrol B.V., the Netherlands

Approved by

Eurotrol B.V.

Date: 2021-03-29

(Signature)

Daniel Philippens

QA/RA Director

Stamp

Eurotrol B.V.
PO Box 722, 6710 BS Ede
Keplerlaan 20, 6716 BS Ede
The Netherlands
www.eurotrol.com
T: +31 315 695777

APOSTILLE

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- 1 Country THE NETHERLANDS
This public document
- 2 has been signed by W.W. Blankestijn
- 3 acting in the capacity of official of the Chamber of
Commerce and Industry
- 4 bears the seal/stamp of aforesaid Chamber of
Commerce and Industry

Certified

- 5 in Arnhem
6. on 02-04-2021
- 7 by the registrar of the district court of Gelderland
- 8 no 21-882

9. Seal/stamp

10 Signature

C. Boers

Eurotrol

Eurotrol Plasma/Low Hb Control

Intended Use

Eurotrol Plasma/Low Hb Control is an assayed hematology control intended for professional use in the verification of the precision and accuracy of the HemoCue Plasma/Low Hb System.

IVD medical device

The Eurotrol Plasma/Low Hb Control complies with the IVD Medical Device Directive 98/79/EC and carries the CE mark. Eurotrol Plasma/Low Hb Control is for in vitro diagnostic use only.

Summary

Eurotrol Plasma/Low Hb Control is a quality control material with a known hemoglobin concentration. Eurotrol Plasma/Low Hb Control may be used in the periodic verification of the precision and accuracy of the HemoCue Plasma/Low Hb system. Please refer to local guidelines for recommended frequency of use.

Reagents

Eurotrol Plasma/Low Hb Control is produced at three physiologically relevant levels. Each vial containing 1.0 mL of purified bovine hemolysate.

Storage and Stability

Eurotrol Plasma/Low Hb Control should be stored unopened in the refrigerator at 2-8 °C (36-46 °F). Stored unopened at this temperature it is guaranteed stable as indicated until the expiration date on the outer box. After opening the vial, it is stable for 30 days when properly recapped and stored at 2-8 °C (36-46 °F) or at room temperature 15-30 °C (59-86 °F).

Procedures/instructions for use

1. Allow the vial to stand for 15 minutes at room temperature 15-30 °C (59-86 °F).
2. Gently mix the vial 8-10 times or until a homogeneous solution is attained, before sampling.
3. Do not fill the cuvette from the vial. Dispense a drop of the control material onto a hydrophobic surface, for example a plastic film. Fill the cuvette according to the manufacturer's instructions. Place the filled cuvette in the instrument for analysis within 60 seconds after filling.
4. Wipe any excess material from the vial and the cap with a clean tissue. Recap the vial tightly.
5. If the control does not perform as expected, review the instructions for use of the instrument to see if the test was performed correctly. Check the expiration date and storage conditions for the control and the cuvettes. Repeat the test. If the control still does not perform as expected, contact technical assistance.

Eurotrol

Precautions

1. For in vitro diagnostic use only.
2. Caution: Eurotrol Plasma/Low Hb Control is an Animal Blood Product. Bovine based materials do not carry biohazards for man, such as Hepatitis B Surface Antigen (HBsAg), Hepatitis C Virus (HCV) and Human Immunodeficiency Virus (HIV-1/HIV-2). This product is free from TSE.
3. It is good laboratory practice to follow the same precautions as with patient samples when handling or disposing vials.
4. This product should not be disposed of in general waste. Consult local environmental authorities for proper disposal.
5. Eurotrol Plasma/Low Hb Control is not to be used as a calibrator.

Expected results

The mean assay value provided on the enclosed assay sheet is derived from replicate analyses on factory calibrated HemoCue Plasma/Low Hb systems. The values are only valid on the HemoCue Plasma/Low Hb system. Results obtained should fall within the expected range specified for the control.

For further information or technical assistance, please contact Eurotrol B.V.:

Phone +31 318 695 777

Fax +31 318 695 770

E-mail: office@eurotrol.com

www.eurotrol.com

Eurotrol B.V., Keplerlaan 20, 6716 BS Ede, The Netherlands.

ASSAYED CONTROLS FOR IN VITRO DIAGNOSTIC USE

Eurotrol

Addendum to the operational documentation*
of medical device

Eurotrol Plasma/Low Hb QCs for Low Hemoglobin measuring in variations

Manufactured by Eurotrol B.V., the Netherlands

2021

* Available on request

Page 4 of 13

Eurotrol B.V. P.O. Box 720, 6710 BS Ede - Soesterlaan 20, 6716 BS Ede - The Netherlands
+31 318 695777 - office@eurotrol.com - www.eurotrol.com
BIC: ABNANL2A - IBAN: NL92ABNA0561838904 - Trade register Arnhem: 09067367 - VAT No. NL009139448901

Eurotrol

Table of Contents

Intended Use.....	6
Composition	6
Main consumer characteristics	7
Service maintenance and repair.....	7
The measurement method is based on the technology of two-wave spectrophotometry at wavelengths of 570 nm and 880 nm.....	8
Analytical and Diagnostic Performance.....	8
Precision performance.....	8
Metrological traceability	8
Measuring range of the assay	8
Clinical Performance and Evidence	8
Mode of application	9
Shelf life	10
The shelf life of Eurotrol Plasma/Low Hb QC's for Low Hemoglobin measuring in variations is 20 months. Expired quality controls are not to be used.	10
Storage and transportation.....	10
Manufacturer's warranty	10
Standards	11
Symbols used.....	13

Eurotrol

Intended Use

Eurotrol Plasma/Low Hb Control is an assayed hematology control intended for professional use in the verification of the precision and accuracy of the HemoCue Plasma/Low Hb System.

Indications for use

For the verification of the precision and accuracy of the HemoCue Plasma/Low Hb System. For in vitro diagnostic.

Application area

For in vitro laboratory diagnostic.

In Vitro Diagnostic Medical Device (IVD)

Eurotrol Plasma/Low Hb complies with the IVD Medical Device Directive 98/79/EC and carries the CE mark. Eurotrol Plasma/Low Hb is for in vitro diagnostic use only.

Description of the analyte (not part of the manufacturer's review)

Hemoglobin is an important component of red blood cells, which is responsible for the transport of oxygen and carbon dioxide mainly to body tissue. Measurement of low concentrations of hemoglobin in plasma, blood serum and irrigating fluids is used to assess the degree of hemolysis. Assessment of the degree of hemolysis is used for quality control in many areas, for example, for analyzing the quality of blood components prepared in banks, in preoperative and postoperative blood processing, as well as in many other areas.

Reagents

Eurotrol Plasma/Low Hb Control is produced at three physiologically relevant levels. Each vial containing 1.0 mL of purified bovine hemolysate.

Composition

Eurotrol Plasma/Low Hb QCs for Low Hemoglobin measuring in variations

Variations:

1. Level '1.0', 2 x 1.0 ml in a package, includes:

- Vial with a red cap – 2 pc.
- Instruction for use – 1 pc.
- Certificate of Analysis – 1 pc.

2. Level '5.0', 2 x 1.0 ml in a package, includes:

- Vial with a white cap – 2 pc.
- Instruction for use – 1 pc.
- Certificate of Analysis – 1 pc.

Eurotrol

3. Level '20.0', 2 x 1.0 ml in a package, includes:

- Vial with a blue cap – 2 pc.
- Instruction for use – 1 pc.
- Certificate of Analysis – 1 pc.

Values of quality controls and acceptable ranges are indicated in the certificate of analysis for each lot of the product.

Note:

1. Eurotrol Plasma/Low Hb QCs for Low Hemoglobin measuring, in variations, are a blood product of animals. Control solutions based on materials of animal origin (bovine blood products) are not characterized by the presence of biologically hazardous factors for humans, for example, the presence of hepatitis B virus surface antigen (HBsAg), hepatitis C virus (HCV) and human immunodeficiency virus (HIV-1 / HIV-2). This product does not contain the causative agents of transmissible spongiform encephalopathy (TSE).

2. The product does not contain materials that come into direct or indirect contact with the body of the patient and personnel using the product when meeting the requirements of operational documentation (instructions for use).

3. The product contains no medicines or pharmaceutical substances.

Main consumer characteristics

The following information is not part of the manufacturer's review and has not been proven by the manufacturer":

50 tests can be done using 1 vial of Eurotrol Plasma/Low Hb.

Sterilization methods of the device

The device does not require sterilization.

Software

N/A.

Service maintenance and repair

The device does not require maintenance and cannot be repaired.

Principal Method

The operational principal of Eurotrol Plasma/Low Hb QCs is based on Hemoglobin measuring, that contains in QCs, on HemoCue Plasma/Low Hb analyzer in the regime «Quality Control».

The obtained measurement results are compared with the certified values of the Eurotrol Plasma / Low Hb QCs for Low Hemoglobin measuring in variations, as a result of which the functional characteristics of the HemoCue Plasma / Low Hb analyzer are checked.

Eurotrol

The measurement method is based on the technology of two-wave spectrophotometry at wavelengths of 570 nm and 880 nm.

Analytical and Diagnostic Performance

Precision performance

The hemoglobin concentration of the three different levels of Eurotrol Plasma/Low Hb Control was measured frequently during a period of 20 (working) days on the HemoCue Plasma/Low Hb photometer. The hemoglobin concentration was measured in duplicate twice a day. For each day new dropper bottles were used from the same batch. Dropper bottles were stored according to the instructions for use of Eurotrol Plasma/Low Hb Control.

Based on the results of this evaluation, we conclude that the within-run precision and total precision is within the clinical acceptable performance criteria (Level '1.0': <10% (CV%), Level '5.0': <7% (CV%) and Level '20.0': <7% (CV%).).

Metrological traceability

At value assignment mean assay values are determined and are derived from replicated analyses on factory calibrated HemoCue systems. Per dropper bottle 20 measurements are performed on the HemoCue Plasma/Low Hb reference system. In addition, the HemoCue Plasma/Low Hb system is traceable to the ICSH method (International Committee for Standardization in Haematology).

Measuring range of the assay

To cover the low hemoglobin levels of (free) hemoglobin, Eurotrol Plasma/Low Hb Control is manufactured at three different levels covering the lower range of (free) hemoglobin.

Device name	Article number	Hemoglobin (g/L)
Eurotrol Plasma/Low Hb QC for Low Hemoglobin measuring Level '1.0'	053.001.002	~0,8
Eurotrol Plasma/Low Hb QC for Low Hemoglobin measuring Level '5.0'	053.002.002	~5,0
Eurotrol Plasma/Low Hb QC for Low Hemoglobin measuring Level '20.0'	053.003.002	~20,0

The measuring range of the HemoCue Plasma/Low Hb system is 0–30 g/L. The HemoCue Plasma/Low Hb system is based on a dual wavelength spectrophotometric technology. Measurements are performed at 570 nm and 880 nm.

Clinical Performance and Evidence

The clinical performance of a quality control such as Eurotrol Plasma/Low Hb Control cannot be evaluated based solely on the QC. Only the combination of the IVD analyzer, reagents and QC have a clinical performance as these need to be used, and can only be used together. Therefore the clinical performance of the IVD analyzer applies also to the QC and no separate data of the clinical performance of the QC Eurotrol Plasma/Low Hb is available.

Test samples

Sampling is not provided.

Eurotrol

Mode of application

Equipment and consumables (for manual testing)

- HemoCue Plasma/Low Hb analyzer, Registration Certificate №P3H 2016/3622 dated 03.02.2016;
- HemoCue Plasma/Low Hb microcuvettes, Registration Certificate №ΦC3 2011/10948 dated 29.01.2016;
- general laboratory equipment.

Procedures/instructions for use

1. Allow the vial to stand for 15 minutes at room temperature 15-30 °C (59-86 °F).
2. Gently mix the vial 8-10 times or until a homogeneous solution is attained, before sampling.
3. Do not fill the cuvette from the vial. Dispense a drop of the control material onto a hydrophobic surface, for example a plastic film. Fill the cuvette according to the manufacturer's instructions. Place the filled cuvette in the instrument for analysis within 60 seconds after filling.
4. Wipe any excess material from the vial and the cap with a clean tissue. Recap the vial tightly.
5. If the control does not perform as expected, review the instructions for use of the instrument to see if the test was performed correctly. Check the expiration date and storage conditions for the control and the cuvettes. Repeat the test. If the control still does not perform as expected, contact technical assistance.

Instructions for use and a list of the obtained values (depending on the batch), as well as the barcode, are contained in the document enclosed in the product box. Due to the use of different data sheets of obtained values depending on the lot and product level, the instructions for use (the same for all levels) are printed separately and attached to each level.

Contradictions

No.

Precautions

1. For in vitro diagnostic use only.
2. Caution: Eurotrol Plasma/Low Hb Control is an Animal Blood Product. Bovine based materials do not carry biohazards for man, such as Hepatitis B Surface Antigen (HBsAg), Hepatitis C Virus (HCV) and Human Immunodeficiency Virus (HIV-1/HIV-2). This product is free from TSE.
3. It is good laboratory practice to follow the same precautions as with patient samples when handling or disposing vials.
4. This product should not be disposed of in general waste. Consult local environmental authorities for proper disposal.
5. Eurotrol Plasma/Low Hb Control is not to be used as a calibrator.

Calculated test

Concentrations of test substances are calculated automatically by the HemoCue Plasma/Low Hb analyzer.

Eurotrol

Expected results

The average value of the analysis and shown in the attached sheet is based on a parallel analysis performed on factory calibrated HemoCue Plasma / Low Hb analyzers. Values are only valid for HemoCue Plasma / Low Hb analyzers. Results should be within the expected range for the control solution.

Shelf life

The shelf life of Eurotrol Plasma/Low Hb QCs for Low Hemoglobin measuring in variations is 20 months. Expired quality controls are not to be used.

Storage and transportation

Storage

Eurotrol Plasma/Low Hb QCs should be stored closed at 2–8 ° C (36–46 ° F) for the entire shelf life in refrigerated chambers or refrigerators providing a regulated temperature regime with daily temperature recording. When these storage conditions are met, the samples remain stable until the expiration date indicated on the outer packaging. After opening the vial, the material remains stable for 30 days, provided that the vial is properly closed and stored at 2–8 ° C (36–46 ° F) or at room temperature 15–30 ° C (59–86 ° F).

Transportation

Eurotrol Plasma/Low Hb QCs for Low Hemoglobin measuring in variations should be transported in the original packaging of the manufacturer, which protects the product from external influences, on all modes of transport in covered vehicles in accordance with the transport rules applicable to this type of transport at temperatures from 2 to 8 ° C for the entire shelf life.

When carrying out loading and unloading operations and transportation, the requirements of the handling signs printed on the packaging and transport packaging should be strictly followed.

Attention: *If the conditions of transportation and storage are not observed, the medical device cannot be used!*

Manufacturer's warranty

Claims about the quality of the medical device are sent to the manufacturer or an authorized representative in the Russian Federation.

Manufacturer:

Eurotrol B.V.

Phone: +31 318 695 777

Fax: +31 318 695 770

Email: office@eurotrol.com

www.eurotrol.com

Keplerlaan 20, 6716 BS, Ede, The Netherlands

Page 10 of 13

Eurotrol B.V. • P.O. Box 777, 6716 BS Ede • Keplerlaan 20, 6716 BS Ede • The Netherlands
+31 318 695 777 • office@eurotrol.com • www.eurotrol.com
BIC: ABNANL2A • IBAN: NL92ABNA0561838909 • Trade register Arnhem 09067367 • VAT No. NL009199448B03

Eurotrol

Authorized Representative in Russia:

LLC «RADIOMETER»

Phone: +7 495 228 92 06

Email: moscow@radiometer.ru

Website: www.radiometer.ru, www.radiometer.com

Address: Russian Federation, 109004, Moscow, Stanislavskogo street, 21, building 3, room I

Utilization

Eurotrol Plasma/Low Hb QCs for Low Hemoglobin measuring in variations must not be disposed of with normal waste. For proper disposal, consult your local environmental authority.

In order to avoid possible environmental pollution, the waste from the used Eurotrol Plasma/Low Hb QCs for Low Hemoglobin measuring in variations, as well as the quality controls with expired shelf life, must be disposed of in accordance with local state and/or hospital regulatory requirements.

The following information in this subsection is not part of the manufacturer's review:

All samples and materials used and generated during research, including quality controls with an expired shelf life, should be disposed of in accordance with SanPiN 2.1.7.2790-10 as potentially infected medical waste (class B).

When disposing of and destroying used Eurotrol Plasma/Low Hb QCs for Low Hemoglobin measuring in variations, as well as with expired shelf life, it is necessary to act in accordance with hygiene standards and follow official instructions for disposal and destruction of potentially infectious and chemical material accepted in the territory of Russian Federation and in this health care facility. Key issues in the waste management strategy are optimization of tactics and approaches to end products, as well as minimization of possible environmental pollution.

When using Eurotrol Plasma/Low Hb QCs for Low Hemoglobin measuring in variations, class B wastes are generated, which are classified and disposed of in accordance with SanPiN 2.1.7.2790-10 "Sanitary and epidemiological requirements for the treatment of medical waste".

Standards

The following is a list of standards according to which Eurotrol Plasma/Low Hb QCs for Low Hemoglobin measuring in variations are produced:

- CEN-EN-ISO 13485: 2016 Medical devices. Quality Management Systems. Regulatory requirements
- CEN-EN-ISO 14971: 2012 Medical devices. Application of risk management to medical devices
- CEN-EN-ISO 18113-1: 2011 Medical diagnostic products in vitro. Information supplied by the manufacturer (labeling). Part 1. Terms, definitions and general requirements

Eurotrol

- CEN-EN-ISO 18113-2: 2011 Medical diagnostic products in vitro. Information supplied by the manufacturer (labeling) - Part 2: In vitro diagnostic professional reagents
- CEN-EN-ISO 23640: 2015 In vitro diagnostic medical equipment
- CEN-EN 13612: 2002 Laboratory evaluation of the technical characteristics of medical devices
- CEN-EN-ISO 11737-2: 2009 Sterilization of medical devices. Microbiological methods. Part 2. Sterility tests carried out during validation of the sterilization process
- ISO 15223-1: 2016 Medical devices. Symbols used for labeling on medical devices, labels and accompanying documentation. Part 1. Basic requirements.
- ISTA-7D: 2007 Thermally controlled shipping packaging for systematic delivery of parcels of goods
- CLSI-EP25-A: 2009 Assessment of the stability of in vitro diagnostic reagents
- CLSI-EP05-A3: 2014 Accuracy Assessment for Quantitative Measurement Methods

The following list is not part of the manufacturer's essential principle review, the conclusions were drawn by Radiometer LLC:

Eurotrol Plasma/Low Hb QCs for Low Hemoglobin measuring in variations also comply with the following standards:

- GOST R 52905-2007 "Medical laboratories. Safety requirements".
- SanPiN 2.1.7.2790-10 "Sanitary and epidemiological requirements for the treatment of medical waste"
- MU 287-113 on disinfection, pre-sterilization cleaning and sterilization of medical devices.
- GOST R 51088-2013 Medical products for in vitro diagnostics. Reagents, reagent kits, test systems, control materials, culture media. Product Requirements and Supporting Documentation
- GOST R ISO 15223-1-2014 Medical devices. Symbols used for labeling on medical devices, labels and accompanying documentation. Part 1. Basic requirements
- GOST 15150-69. Machines, devices and other technical products.
- Versions for different climatic regions. Categories, operating conditions, storage and transportation regarding the impact of climatic environmental factors
- GOST R 51352-2013 Medical products for in vitro diagnostics. Test methods.
- GOST ISO 18153-2011 Medical devices for in vitro diagnostics. Measurement of values in biological samples. Metrological traceability of the values of the catalytic concentration of enzymes attributed to calibrators and control materials.
- GOST R ISO 23640-2015 Investigation of the stability of in vitro diagnostic reagents

Eurotrol

- GOST R ISO 18113-2-2015 Medical products for in vitro diagnostics. Information provided by the manufacturer (labeling). Part 2. In vitro diagnostic reagents for professional use
- Order of the Ministry of Health of Russia dated 06.06.2012 No. 4n (as amended according to the Order of the Ministry of Health of the Russian Federation of September 25, 2014 No. 557n).
- Order of the Ministry of Health of Russia dated 01.09.2014, No. 2n.
- Order of the Ministry of Health of Russia dated January 19, 2017 No. 11n.

Symbols used

Caution!

Use by

In Vitro Diagnostic Medical Device

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