COMMONWEALTH OF MASSACHUSETTS


NOTARY ACKNOWLEDGEMENT

On this 11th day of January, 20 21, before me, the undersigned Notary Public, personally appeared Keith Edwards, as SVP & General Manager of Hamilton Thorne, Inc.,

proved to me through satisfactory evidence of identification, which were PA Driver's License, to be the person whose name is signed on the preceding or attached document, and acknowledged to me that (he) (she) signed it voluntarily for its stated purpose.

Nicole Marie Thibodeau

Nicole Thibodeau
Notary Public
Essex County


CAUTION: Federal (US) law restricts this device to sale by, or on the order of a physician.

LYKOS® Clinical Laser

Instructions for Use

Software Version 6.3

OFFICIAL DOCUMENT
Keith Edwards
KEITH EDWARDS
HAMILTON THORNE, INC.

100 Cummings Center, Suite 465E, Beverly, MA 01915
978.921.2050, Fax: 978.921.0250, sales@hamiltonthorne.com

LYKOS Laser with Dynamic Targeting System Instructions for Use

The content of this manual is the sole property of Hamilton Thorne, Inc., and along with the products to which it applies, is provided for the exclusive use of Hamilton Thorne, Inc. clients. Prior to use of all Hamilton Thorne devices, read all manuals and follow all instructions for use.

In no event does Hamilton Thorne, Inc. assume the liability for any technical or editorial errors or omissions, nor is Hamilton Thorne, Inc. liable for direct, indirect, incidental, or consequential damages arising out of the use or inability to use this manual.

This manual is protected by copyright, all rights reserved, and no part hereof may be photocopied or reproduced in any form without the prior written consent of Hamilton Thorne, Inc.

Document #: M-75-02-32

Revision #: Software Version 6.3, Rev A

Effective Date: May 21, 2019

Federal law in the United States restricts this device to sale by, or on the order of a physician.

Hamilton Thorne, Inc.
100 Cummings Center, Suite 465E
Beverly, MA 01915
978-921-2050
www.hamiltonthorne.com

© Copyright 2019. Hamilton Thorne, Inc. All rights reserved. LYKOS® and RED-i® are registered trademarks of Hamilton Thorne, Inc. All other brand names are trademarks or registered trademarks of their respective companies.

Advena Ltd
Pure Offices, Plato Close
Warwick, CV34 6WE UK

CE Mark in accordance with Directive 93/42/EEC.

WARRANTY, SERVICE INFORMATION, RETURNS AND LIMITATION OF LIABILITY

WARRANTY

HAMILTON THORNE INC. (HT) warrants that its products will be free from defects in workmanship or materials and perform in accordance with HT published specifications (or the specifications agreed to, in writing, by the Buyer and HT and made a part of the sales contract) for a period of one year from the date of installation. This warranty does not cover lost parts or security keys, and shall not apply to damage to the equipment resulting from abuse, negligence, accident or loss due to fire, flood, theft, power fluctuations or power failures, lightning strikes, temperatures or humidity outside of HT published Operating Environments, storage or use in a corrosive environment, off-label use, user-directed system changes such as incompatible computer systems, computer virus or malware induced system changes, use of non-HT approved software on HT provided computer systems, damage in transit or any other damages covered under the Distributor's or Final User's Insurance Policy. The warranty may be voided should the Buyer attempt any repairs, alterations or additions, including installation of third party software, without prior written permission of HT. This warranty is not valid unless a completed Installation Checklist for the system is returned to HT within 30 days from the installation date.

SERVICE INFORMATION

During this warranty period, HT will, at no cost, repair or replace any defective equipment returned to HT. Transportation charges to return the equipment to HT will be prepaid by the sender. The shipping method and packaging are critical to the repair process. Consult HT before shipping. When the Buyer requests expedited shipping or special handling, the Buyer shall pay any associated charges.

HT provides Software Maintenance which includes updates to software such as patches and reliability enhancements during the warranty at no charge. The warranty does not include major software upgrades.

Backup all files before returning the equipment for repair or replacement. HT recommends that you have an external back-up system at all times to reconstruct lost or altered files, data or programs. HT IS NOT RESPONSIBLE FOR ANY LOSS OF DATA.

RETURNS

A Return Authorization Number must be obtained before returning any product to HT. Please call 1-800-323-0503 in the U.S., 1-978-921-2050 outside of the U.S., your local distributor or email support@hamiltonthorne.com for this Return Authorization Number. When calling or contacting HT, please have the serial number of your system available.

LIMITATION OF LIABILITY

HT makes no other warranty, expressed or implied, and HT disclaims any implied warranty of merchantability or fitness for a particular purpose.

The Buyer and HT agree that the sole and exclusive remedies for breach of any warranty concerning the goods shall be repair or replacement of defective parts upon the terms above described or, at HT's discretion, refund of the purchase price. HT shall not be liable for contingent or consequential damages to persons or property, and its sole liability is as above set forth. Any action by the Buyer for any alleged breach of the warranty set forth herein shall be brought to the attention of HT by the Buyer within the warranty period, but not later than 30 days after the alleged breach.

This statement of warranty and limitations of liability is a complete and exclusive statement of all warranty and liability representations of HT. It may not be varied, supplemented, qualified or interpreted by any prior dealings between the parties or by any usage of the trade or upon the face or reverse of any form to which this is attached or a part of, nor may it be modified by any agent, employee or representative of HT unless such modification or representation is made in writing and signed by a duly authorized officer of HT.

Repairs and/or replacement under the terms of this warranty shall not extend the warranty life of the original equipment supplied. All repairs and service must be performed by HT service engineers or by an authorized representative.

For information about purchasing an additional Service Contract, please contact sales@hamiltonthorne.com

QF-03-01-1 Rev. C
Effective Date: May 22, 2012

Securing and Protecting Your Computerized Hamilton Thorne Systems

Hamilton Thorne Inc. (HT) products are designed to operate with Microsoft Windows operating systems. Therefore, as the owner of a Hamilton Thorne system, you have a great deal of freedom in determining how and whether to connect your system to a network, load other programs on the computer, and manage and protect your data. There may also be regulatory standards which pertain to your environment such as risk management for IT networks.

General Use:

First and foremost, this is medical/research equipment and needs to be treated as such. Overloading the system with unnecessary programs, exposing unprotected systems to the Internet, inserting virus-infected memory sticks and other similar actions which routinely cause problems on personal computers are also risks to your HT system. It is recommended that additions or modifications be limited to those required for proper operation of your HT system. It is strongly recommended that you backup your system periodically and prior to making any major changes.

Data Backup:

It is recommended that important data be externally backed up frequently to a location such as an external hard drive or a cloud backup service. Microsoft Windows includes backup utilities and there are also third-party backup and restore solutions available so that data you create as well as the "state" of the system can be backed up and then used in emergency situations to restore the system to a previous operational state.

Networking:

Hamilton Thorne ships most of its systems "network capable". Hamilton Thorne is not responsible for networking the HT system. To connect to a network, contact the system administrator for that network. There are some basic guidelines which are recommended.

- a. First, the HT system operation should be validated prior to network setup.
- b. The system should be backed up/stored at this state. Only then should the networking steps be pursued.
- c. Be aware that although networking is the easiest way to share data, joining a domain can alter performance especially if the HT system is trying to boot from or run programs from a domain. HT programs typically require administrator access to the local hard drive.
- d. A preferred method of network operation is to map a shared network drive that the HT programs can access.

Virus and Malware Protection:

Unless an HT system is secured and isolated both physically and electronically from sources of computer virus or malware (sandboxed), the use of an anti-virus anti-malware program is strongly recommended. Windows Security Essentials and Windows Defender are compatible with HT systems and are free to download from Microsoft or included with Windows. Hamilton Thorne may opt to pre-install anti-virus software on HT systems.

HT Software Program Note:

Many HT programs must install or run with Administrator privileges and may require custom Windows User Account Control settings.

Third Party Hardware and Software:

1. Ensure that third-party software, such as anti-virus programs, do not use excessive CPU resources during times which may interfere with HT system use (e.g. running virus scans or live updates during the day).
2. Installing third-party programs such as spreadsheet, database or PDF capabilities is possible on HT systems (e.g. Open Office, Microsoft Office). This is generally acceptable but it is recommended to minimize the installation and use of additional programs to only what is required to run with the HT system.

Software Release Announcements

When Hamilton Thorne releases new software, a Software Release Announcement (SRA) is generated. The SRA contains a summary of changes to the software including new features added and changes to existing features, as well as any known issues. These announcements are typically included on the software update CD and may also be downloaded from the Distributor and Customer sections of the Hamilton Thorne website. Please read the appropriate SRA before updating to new versions of Hamilton Thorne software.

Table of Contents

Table of Contents	a
List of Figures	e
Chapter 1. Indications, Precautions & Warnings	1
Intended Use	1
Indications for Use	1
Contraindications	1
Conventions used in this Guide	2
Special Notations	2
Laser Safety	2
Installation Precautions	3
Computer Precautions	3
Application Precautions	4
Chapter 2. Laser System Components	5
Components List	5
Reference Images	6
Preparing for Installation	9
Chapter 3. Hardware Installation	11
LYKOS Objective Turret Attachment	11
LYKOS Laser Controller Connections	12
Digital Camera Installation	13
Orienting the Camera to the Stage Movement	14
Computer Connections	15
Desktop Computer	15
Laptop Computer	16
Chapter 4. Start Up and Software Overview	17
System Start Up	17
Software Layout	17
Video Window	17
Control Panel	17
Image & Video Control Panel	18
Main Menu	18
Chapter 5. Setup Screens	21
Laser	21
Auto Labels	21
Image Capture	23
Laser Settings	25
Overlay Settings	29
Time Lapse	30
Video Capture	31
Optics	32
Camera Settings	32
Objectives	33
Chapter 6. Target Alignment and Magnification Calibration	35

Target Alignment	35
RED-i Target Alignment	35
Magnification Calibration	37
Chapter 7. Patient & Embryo Information	39
Patient Info	39
Embryo Panel	39
Info	39
Measurements	40
Data Fields	42
Chapter 8. Working with Images & Videos	43
Thumbnail Gallery	43
Working with Images	43
Opening, Saving and Deleting	43
Capturing Images	44
Adding Images to Report	44
Viewing and Closing Images	44
Image Toolbox	45
Working with Videos	47
Opening, Saving and Deleting	47
Capturing Real Time Videos	47
Capturing Time Lapse Video	47
Playing Videos	48
Chapter 9. Laser Modes and System Operation	49
Working in Full Screen View	49
Clinical Laser Mode	50
Activating the Laser Energy Level	50
Firing the Laser	50
Validation Mode	51
Configuring Validation Mode	51
Validation Control Panel	51
Multipulse Mode	52
Multipulse Laser Settings	52
Multipulse Laser Panel	52
Firing the Laser	53
Built-in Maximum Limits for Multipulse Mode	54
Chapter 10. Dynamic Targeting System	55
Introduction	55
Software Overview	55
Path Type	56
Initialization and Verification	56
Objective Calibration	56
RED-i Target Alignment to Laser	56
Using DTS Option	58
Clinical Mode	58
Validation Mode	60
Working with Path Types	61
IVQC Report	62
IVQC Report Screen Layout	63

Report Panel	63
IVQC Trend Report Toolbar	64
Save and Print Toolbar	64
Chapter 11. Reporting & ASCII	65
Working with Reports	65
Report Template	65
Viewing Reports	66
Designing Reports	66
Report Property Descriptions	69
General	69
Working with ASCII Data	74
Chapter 12. Laser Applications	77
Getting Ready	77
Check Target Alignment	77
Focusing on the Embryo	77
RED-i	77
Laser-assisted Hatching (LAH)	77
Laser Positioning	77
Laser Zona Drilling Options	77
Clinical Laser Energy Level	78
Summarized Laser Assisted Hatching Procedure	79
Troubleshooting Failed Drilling Attempts	80
Laser-assisted Blastomere Biopsy	80
Laser-assisted Trophectoderm Biopsy	81
Using Manual Pulses to Remove Trophectoderm Cells	81
Using Multipulse Mode to Remove Trophectoderm Cells	82
Chapter 13. Help Menu	85
About	85
Database	85
Backup	85
Restore	85
Optimize	86
Save Support Package	86
Operating Environment	87
Rated Input Parameters	87
Laser Specifications	87
Dimensions	87
Electrical	88
Laser	88
Laser Safety	88
RED-i Target Locator	89
Objective	89
Cleaning & Disinfection	89
Index	I

Table of Contents

List of Figures

Figure 2-1. LYKOS Laser Objective	6
Figure 2-2. Laser Controller and Accessories	6
Figure 2-3. RMS Thread Adapters and Thread Spacers	6
Figure 2-4. 1.5mm Hex Head Drivers	7
Figure 2-5. USB Security Key	7
Figure 2-6. C-Mount Adapters (representative samples).....	7
Figure 2-7. USB 3.0 Digital Camera and Cable (if ordered)	7
Figure 2-8. Foot Switch	8
Figure 2-9. Stage Micrometer	8
Figure 2-10. HT Laser Target Slide	8
Figure 2-11. Optional Compensator Cap	8
Figure 3-1. LYKOS laser installation using Leica M25 RMS adapter	12
Figure 3-2. Laser Cable Attachment Reference Image	12
Figure 3-3. Laser Controller Rear Panel	13
Figure 3-4. Laser Controller LED Status Bar	13
Figure 3-5. JAI Digital Camera Installation using Leica OEM C-Mount Adapter (example).....	14
Figure 4-1. Software Layout	18
Figure 5-1. Laser Settings Main Screen	21
Figure 5-2. Auto Labels	22
Figure 5-3. Programmable Field for Label Text	23
Figure 5-4. Image Capture	23
Figure 5-5. Auto: Use Patient ID and Embryo ID	24
Figure 5-6. Save Mode	24
Figure 5-7. Laser Settings	25
Figure 5-8. Circle and Arrow Target Types	26
Figure 5-9. Target Settings	26
Figure 5-10. Isotherm Rings Target Type	27
Figure 5-11. Overlay Settings	29
Figure 5-12. Crosshairs Inner & Outer Radius	29
Figure 5-13. Time Lapse Settings	30
Figure 5-14. Video Settings	31
Figure 5-15. Camera Settings	32
Figure 5-16. Objectives Settings Panel	33
Figure 6-1. HT Laser Target Slide	35
Figure 6-2. RED-i Adjustment	36
Figure 6-3. Thread Spacer Rings for RED-i Alignment	36
Figure 6-4. X Calibration Screen	37
Figure 7-1. Patient Info	39
Figure 7-2. Embryo Info Panel	40
Figure 7-3. Embryo Measurements	41
Figure 7-4. Images Report Buttons	42
Figure 7-5. Change Data Field Label and Persist	42
Figure 8-1. Thumbnail Gallery	43
Figure 8-2. Capture Options	44
Figure 8-3. Report Images	44
Figure 8-4. Image Toolbox	45
Figure 8-5. Full Screen Mode: Image Toolbox	46

List of Figures

Figure 8-6. Video Playback Controls	48
Figure 9-1. Laser Mode Menu.	49
Figure 9-2. Full Screen View: Clinical Mode.	49
Figure 9-3. Laser Energy Level set to Off (left) and Medium (right).....	50
Figure 9-4. Validation Control Panel	51
Figure 9-5. Multipulse Laser Settings.....	52
Figure 9-6. Multipulse Control Panel	53
Figure 10-1. Clinical DTS and Validation DTS Panels	55
Figure 10-2. Dynamic Laser Info Panel	57
Figure 10-3. Verify DTS Screen.....	57
Figure 10-4. Verification Results	58
Figure 10-5. Clinical Pulse Settings and Hole Size	58
Figure 10-6. Clinical Mode with DTS Enabled	59
Figure 10-7. DTS Path Types	59
Figure 10-8. Validation Control Panel with DTS Enabled	60
Figure 10-9. Draw Line.....	61
Figure 10-10. Manipulation of Laser Path	62
Figure 10-11. IVQC Report Screen.....	62
Figure 10-12. Laser and RED-i Offset Data	63
Figure 10-13. IVQC Trend Chart	64
Figure 10-14. IVQC Trend Report Toolbar.....	64
Figure 11-1. Design View and Report View.....	65
Figure 11-2. Reports Template Panel	66
Figure 11-3. Report Query	66
Figure 11-4. Blank Report Template.....	67
Figure 11-5. Report Design Menu.....	67
Figure 11-6. Alignment Guide	68
Figure 11-7. Quick Properties Dialog.....	69
Figure 11-8. Report Properties	69
Figure 11-9. Header, Body, Footer Properties	70
Figure 11-10. Image Properties.....	71
Figure 11-11. Text Properties	72
Figure 11-12. ASCII Filter Panel	74
Figure 12-1. LAH: Full Zona Breach (left); Zona Thinning (right)	78
Figure 12-2. Examples of Zona Openings for Blastomere Biopsy	81
Figure 13-1. Save Support Package	86

Chapter 1. Indications, Precautions & Warnings

1.1 Intended Use

The microsurgical laser system ("Laser") is intended for in vitro fertilization (IVF) laboratory use to drill or thin the zona pellucida of an IVF gamete or embryo during assisted reproduction hatching procedures, and to biopsy cells for purposes of genetic diagnosis or screening.

Indications for Use

- Whenever assisted hatching of IVF embryos is indicated, including infertile mothers with poor prognosis for successful pregnancy outcomes due to advanced maternal age, prior unsuccessful IVF procedures, hardened or cryopreserved IVF embryos or abnormal zona pellucida morphology.
- IVF embryos subject to biopsy for preimplantation genetic diagnosis (PGD).
- IVF embryos subject to biopsy for preimplantation genetic screening (PGS).

Caution Federal (US) law restricts this device to sale by or on the order of a physician or licensed healthcare practitioner.

1.2 Contraindications

There are no known contraindications for the use of the Laser. There are no known circumstances where it interacts with a drug. No contraindications have been imposed by the manufacturer or regulatory agencies and no reports from the use of the device have resulted in evidence requiring contraindications.

1.3 Conventions used in this Guide

Some of the text in this manual uses special formatting to help indicate emphasis or keystrokes. The text conventions are as follows:

Table 1-1: Text Conventions

Convention	Example	Meaning
Bold	From the dropdown menu, select Calibrate this image	Bold text will indicate a menu item, control, function or setting
Bold with Arrows	Select Patient > New and enter patient information in the form fields.	Indicates menu selection and menu choice.

Special Notations

NOTE: *Indicates an important point to keep in mind.*

Precaution: *Indicates a condition that must be followed to ensure proper operation and performance.*

Radiation Warning: *Indicates a possible hazardous laser radiation exposure.*

ADMIN: *Indicates the users must have administrator privileges to perform a function.*

1.4 Laser Safety

The Laser system contains an infrared laser. As with any laser-based system, appropriate care must be taken to ensure safe and proper use.

The Laser system is a Class I laser product (U.S. Code of Federal Regulations 1040.10, IEC 60825-1, Am.2 2001). When installed and used as instructed, the user is not exposed to potentially harmful laser radiation during normal operation.

Complementary features and operating instructions have been designed to avoid exposure to potentially harmful laser radiation.

Radiation Warning: *Use of controls or adjustments or performance of procedures other than those specified in this manual may result in hazardous laser radiation exposure.*

DO NOT disassemble or modify the laser objective. Doing so may result in exposure to hazardous levels of infrared radiation and voids the system warranty and/or service agreement.

To reduce risk of injury:

DO NOT use a magnifying glass or eye loupe above the objective lens at any time.

DO NOT direct the laser beam at unprotected personnel.

DO NOT operate the laser unless properly mounted to a functioning microscope.

1.5 Installation Precautions

The Laser installation process requires the assistance of an HT certified or equivalently qualified service technician. Proper system installation and compliance with all precautions ensure there is minimal risk of exposure to hazardous levels of laser radiation during operation. However, incorrect installation could result in instrument failure and/or increased risk of exposure to hazardous laser radiation.

The Laser Controller should only be used with the power supply provided. DO NOT attempt to power or control the laser by any equipment other than that provided with the Laser system. Doing so may cause a hazard and voids the warranty and/or service agreement.

DO NOT exchange the Laser Controller received with the laser system with that from any other laser system. The serial number of the Laser and the serial number of the Laser Controller must match exactly.

To ensure proper operation and prevent damage to the system, the Laser Controller must be unplugged from the power supply and detached from the computer before connecting or disconnecting the laser cable to the Laser objective.

1.6 Computer Precautions

The Laser system is a dedicated medical device which is based on PC computer technology.

To ensure proper operation, DO NOT install any program which is not directly required for the operation of the device.

Some basic system precautions are as follows:

1. Backup the system hard disk periodically per your laboratory's protocol before making any significant changes.
2. Anti-virus software, anti-spyware, and anti-malware software are recommended anytime the system will be connected to other computers on a network or the Internet.
3. Avoid installing and running other programs or hardware which heavily tax the CPU concurrently with the systems application.

4. It is recommended that a surge protector is used between this instrument and the power socket.

In the event that other connections, programs, or hardware are required to meet lab needs, all precautions must be taken to verify proper operation and ensure that the device remains fully functional for its intended applications.

1.7 Application Precautions

No therapeutic or clinical laser procedures should be performed under fluorescent illumination.

The Time Lapse mode should not be used for therapeutic or diagnostic procedures.

To eliminate any significant effects from thermal lensing, keep the pulse duration below 667 μ sec.

Additional procedural precautions are provided along with the specific applications of laser-assisted hatching, laser-assisted blastomere biopsy and laser-assisted trophectoderm biopsy (see Chapter 12).

Chapter 2. Laser System Components

2.1 Components List

The components included with the laser system will vary depending on the camera type (digital or analog) and computer format (desktop or laptop).

When unpacking your laser system, please make sure all components are available before installing the laser. Please refer to "Reference Images" on page 6 to help identify the laser system components.


Table 2-1: Components List (n/s = not shown)

Ref #	Component		LYKOS	
			Desktop	Laptop
2-1	LYKOS Laser Objective		Yes	Yes
2-2a	Laser Controller and Accessories	Laser Controller Box	Yes	Yes
2-2b		HDMI cable	Yes	Yes
2-2c		Power Supply	Yes	Yes
2-2d		USB Cable	Yes	Yes
2-3	RMS Thread Adapters (3) and thread spacer rings (3)		Yes	Yes
2-4	1.5mm Hex Head Drivers (2) used for RED-i adjustment (1 long, 1 short)		Yes	Yes
2-5	Green USB Security Key		Yes	Yes
2-6	C-mount Adapter (type differs based on microscope model)		Yes	Yes
2-7	Digital camera with USB cable		Yes	Yes
2-8	Foot switch		Yes	Yes
2-9	Stage micrometer		Yes	Yes
2-10	HT Laser Target Slide		Yes	Yes
2-11	Compensator Cap **Optional**		No	No
n/s	Wireless Keyboard (additional to wired)		Yes	No
n/s	Wireless Mouse (additional to wired)		Yes	Yes
n/s	Software CDs (HT Laser, Windows OS)		Yes	Yes
n/s	Operator's Manual		Yes	Yes
n/s	Desktop Computer System including • PC mini tower & power cord • Keyboard & mouse (wired) • Speakers • Monitor, monitor cable, and power cord		Yes	No
			Desktop	Laptop
n/s	Laptop Computer System including • USB hub • Mouse (wired) • External USB DVD Drive • Power supply		No	Yes

2.2 Reference Images

Figure 2-1: LYKOS Laser Objective

- a. Laser Controller Box
- b. HDMI Cable
- c. Power Supply
- d. USB Cable

Figure 2-2: Laser Controller and Accessories


Figure 2-3: RMS Thread Adapters and Thread Spacers

Figure 2-4: 1.5mm Hex Head Drivers

Figure 2-5: USB Security Key

Figure 2-6: C-Mount Adapters (representative samples).


Figure 2-7: USB 3.0 Digital Camera and Cable (if ordered)

Figure 2-8: Foot Switch


Figure 2-9: Stage Micrometer

Figure 2-10: HT Laser Target Slide

Figure 2-11: Optional Compensator Cap

2.3 Preparing for Installation

All Laser system components are overpacked and shipped in a single carton. Unpack all boxes in the shipping carton and verify components against the packing list.

The laser objective, laser controller, green security key, adapters, spacers, shims, and hex key or hex drivers are packaged together inside hard-sided, black plastic case.

The computer and monitor are packed in their original boxes. All accessories, including camera and C-mount, are packed in a single box.

Chapter 3. Hardware Installation

3.1 LYKOS Objective Turret Attachment

Caution: Electronic components are sensitive to static discharge. To reduce the effects of static electric charges, touch a non-conductive, large, metal object prior to handling electronic components.

1. If needed, attach an RMS thread adapter to an open port on the microscope turret. Three types of RMS thread adapters are included (see Table 3-1) and use depends on the make and model of the inverted microscope.

Table 3-1: RMS Thread Adapters

Type	Description	Examples of Compatible Microscopes
	LYKOS RMS Thread	Nikon Diaphot models ZEISS Olympus Leitz Fluorovert All older RMS microscopes
	NIKON CF160 Large Chrome 25 mm threads	(Nikon CF160 optics) Nikon Eclipse TE 200/300 Nikon Eclipse TE 2000 Nikon Ti-2000 series Nikon Ti2
	LEICA M25 Medium Brass 25 mm threads	Leica DMIL Leica DMIRB (HC optics) Leica DMI 3000/4000/6000
	ZEISS M27 Medium Brass 27 mm threads	ZEISS Axio Observer 3/5 ZEISS Axio Vert.A1 ZEISS 200 Axiovert

2. Carefully remove LYKOS laser module from its plastic case.

3. Attach the LYKOS objective to the RMS thread adapter (if required) or directly to the turret.


Figure 3-1: LYKOS laser installation using Leica M25 RMS adapter

NOTE:

The LYKOS objective should be situated so that the mini-HDMI connector is facing inward (towards the center of the turret) and the RED-i adjustment ports towards the outside for easy access (see Figure on page 35). If necessary, Thread Spacer Rings may be inserted to obtain the proper orientation (see Figure 6-3 on page 36).

NOTE:

Ensure that both the RMS thread adapter and the laser objective are secured tightly to the turret. Any movement of the laser module relative to the turret may affect the target alignment.

3.2 LYKOS Laser Controller Connections

Caution:

To reduce the risk of shock and prevent damage to the system, the Laser Controller must be unplugged from the power supply and detached from the computer before connecting or disconnecting the HDMI laser cable to the Laser objective.

1. Make sure that the laser controller is not attached to the power supply or the computer.
2. Attach the standard connector of the HDMI cable to the controller box.
3. Attach the mini connector of the HDMI cable to the LYKOS objective installed on the microscope turret.

Figure 3-2: Laser Cable Attachment Reference Image

4. Connect the USB A/B cable between the laser controller and an available USB port on the computer.


Figure 3-3: Laser Controller Rear Panel

5. Attach the external power supply to the laser controller and plug into a surge protected power source.

Figure 3-4: Laser Controller LED Status Bar

3.3 Digital Camera Installation

1. Locate the camera port of the inverted microscope. Its location varies depending on microscope model.
2. Attach the C-mount adapter to the camera port.
3. Attach the digital camera to the C-mount adapter (see Figure 3-5).
4. Tighten the C-mount adapter set screws.
5. Attach the mini USB connector of camera cable to the camera. Tighten the pin screws to secure.
6. Attach the camera cable directly to the USB 3.0 port labeled **Camera** on the computer. To ensure proper operation, do not connect the camera cable to any other USB port or use a USB hub.
7. Once the system is on and software is running, refer to “Orienting the Camera to the Stage Movement” on page 14 to complete the camera installation process.

Figure 3-5: JAI Digital Camera Installation using Leica OEM C-Mount Adapter (example).

NOTE: Cameras will only operate correctly when connected to the USB 3.0 port labeled "Camera."

NOTE: The USB 3.0 connection powers the digital camera so there is no need for an external power supply.

Orienting the Camera to the Stage Movement

1. Once the system is completely installed, open the HT Laser 6 software so that the live camera image appears on the screen.
2. Loosen the bottom camera set screw and rotate the camera as needed to provide proper orientation with the video screen. Make sure that the video image is horizontal and moves in the same direction as the stage movement (left to right / up and down).
3. When the correct orientation is obtained, securely tighten the camera set screw so that the camera cannot rotate.

Caution: Ensure that the camera and C-Mount are adequately secured. Any camera movement may affect the target position.

3.4 Computer Connections

Desktop Computer

Digital Camera


- Attach the mini USB connector of the camera cable to the digital camera and the standard USB connector to the USB 3.0 port labeled Camera on the computer.

NOTE: To ensure proper operation, the Digital Camera must be connected directly to a USB 3.0 port on the computer.

Monitor

- Attach the 15-pin monitor cable between the monitor and the computer.
- Attach the power cord to the monitor and plug into a surge protected power source.
- If you wish to use the USB slots on the monitor, connect the USB A/B cable (shown right) between the computer (USB A-type) and the monitor (USB B-type).

Security Key

- Insert the green security key into a USB port on the computer or monitor.

NOTE: The USB Security key represents the license for the Hamilton Thorne laser software and is required to run the software. Loss of this key will incur a high replacement cost. If you move the system, be careful not to lose this key.

Keyboard and Mouse

- For wired components, connect the keyboard and mouse to an available USB port on the front or back of the computer, or on the monitor.
- If using wireless components, insert the included wireless USB receivers into available USB ports on the front or back of the computer.
 - Ensure batteries are inserted in the mouse and keyboard and power switches are set to On.
 - Ensure wireless USB receivers are connected (or download and run the Unifying software available on the Logitech website [<http://support.logitech.com>] to use just one receiver).

Speakers

- Attach the green mini stereo plug of the speaker cable to the Audio Out jack on the back or front of the computer.
- The speakers are powered through the monitor. Connect the black power plug to the power receptacle on the back of the monitor.

Foot Switch

- The foot switch may be connected to a USB port on the front or back of the computer, or the monitor.

Computer Power

- Connect the computer power cord to the computer plug into a surge protected power source.

Laptop Computer

Digital Camera

- Connect the mini USB plug of the camera cable to the digital camera and the standard plug into the USB 3.0 port labeled Camera on the computer.

To ensure proper operation, the Digital Camera must be connected directly to a USB 3 port on the computer. DO NOT CONNECT THROUGH A USB HUB.

Security Key

- Insert the green security key into a USB port on the computer or USB hub.

NOTE:

The USB Security key represents the license for the Hamilton Thorne laser software and is required to run the software. Loss of this key will incur a high replacement cost. If you move the system, be careful not to lose this key.

Mouse

- If using a wired mouse, connect it to an available USB port on the front or back of the laptop computer, or the USB hub.
- If using the wireless mouse, connect the wireless receiver to available USB ports on the front or back of the computer, or USB hub.
 - Make sure that the batteries are inserted in the mouse.
 - Make sure the power switch on the mouse is set to On.

Foot Switch

- The foot switch may be attached to a USB connection on the front or back of the laptop computer, or the USB hub.

Computer Power

- Attach the power plug between the laptop computer and a surge protected power source.

Chapter 4. Start Up and Software Overview

4.1 System Start Up

WARNING: Before turning on the system, make sure that all components with manual voltage switches are set to the proper voltage. Failure to do so may cause irreparable damage to the system. Hamilton Thorne is not responsible for damage incurred if an improper voltage setting is used.

1. Before powering up, make sure all components are connected properly and plugged into a grounded surge protected power supply. For the highest level of protection, use a UPS (Uninterruptible Power Source) Battery Backup.
2. Turn on any components with an external ON/OFF control.
3. When the Windows desktop appears, double-click the HT Laser 6 software icon.
4. Check the status bar at the bottom left of the screen to confirm that the laser controller is properly connected.
 - a. If the message "ZTKMK3004 on COMx" appears, the laser is ready for use (x will be replaced by the COM port number).
 - b. If the message "Disconnected" appears, check that the laser serial or USB cable is connected to both the computer and laser controller (see Chapter 3) and verify that it is receiving good power (see Figure 3-4).

NOTE: As a safety precaution, laser operation will automatically disconnect when any of the setup, reporting or ASCII screens are open.

4.2 Software Layout

The software screen is divided into four distinct areas: the Video Window, the Control Panel, the Image and Video Control Panel, and the Menu Bar.

Video Window

- The Video Window displays either the live image from the camera or a saved image or video file.
- When viewing the live image, a Target Overlay is used to indicate the position of the laser beam in the field of view and, if activated, a set of Crosshairs indicates the center of the image area.

Control Panel

The Control Panel, located along the right hand side of the screen, provides access to all controls needed for operation and data entry. It is divided into five sections:

1. The Laser - Clinical panel provides controls to fire the laser. This panel will change based on the Laser Mode selected - Clinical, Validation or Multipulse (see "Laser Modes and System Operation" on page 49).
2. RED-i Enable control to turn RED-i target on and off and the RED-i Intensity Slider to adjust the brightness of the RED-i target beam.

3. The **Patient Info** panel allows entry of patient details or selection of the patient from the internal database.
4. The **Embryo** panel allows input of embryo identification details, provides access to tools to measure structures directly on the captured image and permits addition of user defined **Data Fields**.
5. The **Images** panel provides controls to save up to four images to the report.


Figure 4-1: Software Layout

Image & Video Control Panel

The **Image and Video Control Panel**, located below the **Video Window**, provides access to controls for capturing still images and real time or time lapse videos. The **Thumbnail Gallery** displays the captured or opened images and videos.

Main Menu

The **Main Menu** located along the top of the screen provides access to the following options: **Laser Mode**, **Objective**, **Reports**, **ASCII**, **Settings** and **Help**.

Laser Mode

The three laser modes available are **Clinical**, **Validation** and **Multipulse**. Please see “Laser Settings” on page 25 for more information.

Objective

The options under **Objective** include access to **Calibrate** (see “Magnification Calibration” on page 39) objectives, to select previously configured objectives (see “Objectives” on page 33) and to **Align Laser Target** (“Target Alignment” on page 35).

Reporting

The **Reporting** control opens the **Report Panel** which provides access to design, view and save patient reports. Please see "Reporting & ASCII" on page 65 for more information.

Ascii

The **Ascii** control opens the **ASCII Panel** which provides access to filter, view and export saved data. Please see "Reporting & ASCII" on page 65 for more information.

Settings

The **Settings** control opens the **Settings Panel** which provides access to laser, software and optics settings. Please see "Setup Screens" on page 21 for more information.

Help

The **Help** control provides access to the **License Agreement**, database tools, and support package creation. Please see "Help Menu" on page 85 for more details.

Chapter 5. Setup Screens

Select **Settings** from the **Menu Bar** to open the **Settings** dialog box. Selecting the different options from the setup tree panel on the left changes the controls available in the main section of the dialog box.

- To activate the new settings entered, the **OK** or **Apply** button must be selected.
- To exit the Setup screen without saving, press **Cancel**.


Figure 5-1: Laser Settings Main Screen

5.1 Laser

Auto Labels

The **Auto Labels** field provides the ability to customize labels that may be overlaid automatically on every image captured with editable options like the position, name, font, color, etc. The label may also contain variable field names. Multiple labels may be created and saved in the system.


Figure 5-2: Auto Labels

Auto Label Name

- To create a new label, select **New** and enter the label name.
- To edit an existing label name, select the name from the list and then **Edit**.
- To delete an existing label, select the name from the list and then **Delete**.

Font Properties

- **Font:** Select **Change** to update the font, size, and style of text to be used for the selected label.
- **Color:** Open the color palette to change the text color of the currently selected label.

Label

- **Programmable Fields:** Labels may be designed with a combination of manually entered text and **Programmable Field** variables (see Figure 5-3). To add a variable, select it from the list and press **Add**.
- **Edit Label Location:** Each label can be positioned as needed on the image area. Using the mouse, click and drag the label to the desired location. Label location editing requires a live video feed.
- **Enable:** Check to activate the label.


Figure 5-3: Programmable Field for Label Text

Image Capture


Figure 5-4: Image Capture

Save Location

The **Save Location** is the path/folder where files will be saved.

Image Type

Image files may be stored in JPG or TIF format. When saving images, the objective calibration factor (see “Magnification Calibration” on page 39) of the currently selected objective is automatically included in the image header information. This allows for accurate measurements to be made on recalled images.

Save Mode

When saving images, files may be named using one of three modes:

- **Manual: Ask for File Name:** If this option is used, a standard Save As dialog box will open and the file name is input from the keyboard.
- **Auto: Use Patient ID and Embryo ID:** If this option is selected, a file name will be automatically generated using the Patient ID, Retrieval Date, Current Date, Embryo ID and Count value. These values must be added to the Patient and Embryo Info prior to capturing images.

Patient ID	Retrieval Date	Current Date	Embryo ID	Counter
PID1234	2017-12-08	2017-12-13	EID1234	16.jpg

Figure 5-5: Auto: Use Patient ID and Embryo ID

- **Auto: Use Root Name:** If this option is selected, a file name will be automatically generated using the name entered in the Root Name text field and Count value.

Figure 5-6: Save Mode

NOTE:

The Count number is NOT DERIVED AUTOMATICALLY. It is generated from the manually entered Count value (see Figure 5-6), which increments as each image is saved. Changing the Count value will change the starting value appended to the image file name.

Autosave on Capture

For Individual File Storage, enabling Autosave on Capture causes the program to immediately save the captured image to the Save Location. It is advised that one of the Auto methods of saving files is used to avoid having to enter the file names manually.

Image File Name Label

To embed an identification label consisting of either the Full Path + File Name or File Name only onto the saved images, enable Image File Name Label. The color of the embedded label may be selected using the color chart palette.

Watermark

To save the Hamilton Thorne watermark on all saved images, check the Enabled box.

Laser Settings


Figure 5-7: Laser Settings

Laser

The **Laser Enable** toggle control enables/disables the laser. The laser may not be fired if the box is unchecked.

Remote Fire

The **Remote Fire Enable** toggle control enables/disables the firing of the laser via the foot switch. **Remote Fire** is grayed out if **Laser Enable** is not selected.

Auto Detect

The **Auto Detect** toggle control enables/disables automatic location of the USB port to which the laser is connected.

Target Type

The **Target Type** sets the graphic overlay used to designate the laser spot on the live image. Three options are available: **Circle**, **Arrow** or **Isotherms**.

Circle

- When using the **Circle** target, it should be aligned so that the laser fires in the middle of the **Circle** (see "Target Alignment" on page 35).
- The **Circle** size is NOT representative of the hole size.
- The color of the **Circle** is red and cannot be changed.

Figure 5-8: Circle and Arrow Target Types

Arrow

- When using the **Arrow** target, it should be aligned so that the laser fires at the tip of the **Arrow** (see “Target Alignment” on page 35).
- The color of the **Arrow** is red and cannot be changed.

Isotherms

Isotherms appear on the screen as a series of up to six concentric circles, with colors corresponding to temperature (see Table 5-1). Each circle indicates the peak temperature reached at that radius of the circle when the laser is fired at its current power. The size of each circle is derived by the software as a function of the laser power and laser pulse duration and will scale accordingly. A color key reference appears at the upper left corner of the live video window when **Isotherms** are enabled.

Table 5-1: Isotherm Color and Peak Temperature Key

Red	Laser Beam Diameter
Orange	140° C (approximate hole size in zona)
Yellow	100° C
Green	80° C
Blue	60° C
Violet	50° C

**NOTE:**

When operating in Clinical Mode, the Isotherm Rings only appear after the High, Medium or Low is selected on the Clinical Control Panel.


Figure 5-9: Target Settings

The **Isotherms** which appear on the screen may be activated or deactivated using the check boxes. In addition, one custom ring may be generated by entering the temperature value in the **Custom Temp** box. A custom color may also be selected for that ring. The central red ring that designates the laser beam is always enabled.

Figure 5-10 shows three examples of **Isotherm** activation. The image on the left shows all **Isotherms** activated. The center image shows only the 50° C **Isotherm** circle, while the image on the right shows only the 100° C **Isotherm** circle.


Figure 5-10: Isotherm Rings Target Type

Blanking

The **Black**, **Hide** and **None** options indicate what happens to the **Target** during the **Duration** period or the built-in safety delay (see NOTE under **Duration** below).

- To still see the **Target** on the image but see a visual indicator that the laser has been fired, set the **Target** to change to **Black** after laser firing.
- To turn off the **Target** after laser firing, select **Hide**.
- To make no change to the **Target** after laser firing, select **None**. If set to **None**, the laser still cannot be fired for the **Blanking Duration** or built-in safety delay, whichever is longer.

NOTE:

If you are using the Capture Overlay feature for image capture, Blanking should be set to None.

Duration

The **Duration** is the time the **Target** turns black or is hidden after the laser is fired. The **Duration** may be set from 0 to 10 seconds. This feature provides visual confirmation that the laser has been fired. The laser may not be fired during this **Duration** period.

NOTE:

Even if the Duration is set to 0 seconds, there is a built-in safety delay for Class 1 lasers in which the laser may not be fired. The longer the pulse duration, the longer the safety delay. However, the Duration, if set for a longer period than the safety delay, will take precedence.

Fire Animation

Upon firing the laser, animation of the target may be set to **None**, **Pulse** or **Ping**. This is an alternative method to **Blanking** to visually confirm that the laser has been fired.

- **None:** No target animation will occur.
- **Pulse:** The target itself will pulse.
- **Ping:** Circles will radiate out from the target.

NOTE: *If Target Animation is set to Pulse or Ping, it is recommended that Blanking be set to None.*

Pulse Resolution

The **Pulse Resolution** controls the interval the **Pulse Slider** value increases or decreases on the **Validation Control Panel**.

Multipulse

When in **Multipulse** mode, **Enable Duration (sec)** and **Mean Power (mw)** options are used when firing the laser. Please see "Multipulse Laser Settings" on page 52 for more information.

Clinical Pulse Settings

When in **Clinical** mode, the laser operates only at three pulse settings of **Low**, **Medium** or **High** with the laser power fixed at 100%. In addition, the estimated hole size (in μm) based on the selected pulse setting is displayed.

The pulse durations for each level may be set as follows:

- **High:** from 600 μs to 3000 μs
- **Medium:** from 200 μs to 600 μs
- **Low:** from 0 μs to 200 μs

Values may be adjusted by using the scroll arrows or entering values from the keyboard. If a value outside the limits is entered, the system will automatically reset it to the closest acceptable value.

NOTE: *Hamilton Thorne recommends pulse settings of 600 μs , 400 μs , or 200 μs .*

For each entered value, the built-in safety delay (the time that must pass between firing the laser again) is adjusted accordingly. The safety delays for the default values are:

- **600 μs :** 18 milliseconds
- **400 μs :** 12 milliseconds
- **200 μs :** 6 milliseconds

Audio Feedback

When **Audio Feedback** is enabled, a sound file (.wav) may be selected to annunciate upon laser fire.

NOTE: *External desktop speakers may be required on older systems for audio.*

Overlay Settings


Figure 5-11: Overlay Settings

Crosshairs

The **Crosshairs** are useful for positioning the embryo in the field of view. When the embryo is placed beneath the **Crosshairs** at low magnification, it can more easily be located when switching to a higher power objective.

Crosshair options that may be configured include:

- Enabling or disabling
- Dimensions of the **Inner** and **Outer Radius** in pixels
- Thickness of the **Crosshairs**
- Color of the **Crosshairs**
- Appearance of a dot in the center of the **Crosshairs**
- Alpha (transparency) of the **Crosshairs**


Figure 5-12: Crosshairs Inner & Outer Radius

Video Capture


Figure 5-14: Video Settings

Save Location

The Save Location is the path/folder where files will be saved.

Save Mode

When saving images from the Thumbnail Gallery (see “Working in Full Screen View” on page 46), users have the option to save individual video files or save all captured video files. For either method, files may be named using one of three methods:

- **Manual: Ask for file name:** If this option is used, a standard Save As dialog box (see “Working in Full Screen View” on page 46) will open and the file name can be input from the keyboard.
- **Auto: Patient ID and Embryo ID:** If this option is selected, a file name will be automatically generated using the Patient ID, Date of Oocyte Retrieval, Current Date, Embryo ID and Count value. These values must be added to the report prior to capturing images.
- **Auto: Use Root name:** If this option is selected, a file name will be automatically generated using the name entered in the Root Name text field and Count value.

NOTE:

The Count number is NOT DERIVED AUTOMATICALLY. It is generated from the manually entered Count value (see Figure 5-14), which increments as each image is saved. Changing the Count value will change the starting value appended to the image file name.

5.2 Optics

Camera Settings


Figure 5-15: Camera Settings

The following **Camera Settings** are available:

- **Device:** Select the installed camera from the Device drop-down list.
- **FPS Set:** Enter the desired number of frames per second. It is recommended to use FPS of 30 or 60.
- **FPS Actual:** Based on the FPS Set value, this is the actual FPS.
- **Resize:** Enabling resize will resize the image to fit within the Live Video Window
- **Flip Horizontal / Flip Vertical:** Enable to flip the orientation of the displayed image from the camera.
- **Settings:** Provides access to set the **Black Level** and **Gain Raw** camera values. The exposure time cannot be changed as auto exposure is used.
- **White Balance:** Set the **Blue**, **Green** and **Red Gain** values or use **Auto White Balance**.

Objectives

Select **Objectives** from the navigation panel to open the **Objectives** setup panel.


Figure 5-16: Objectives Settings Panel

Working with Objective Setups

- **New:** Adds a new **Objective** to the list.
- **Copy:** Creates a clone (copy) of the currently selected **Objective**.
- **Delete:** Deletes the currently selected **Objective**.

Objective Info

- **Name:** Enter a short name for the **Objective**. This name will also appear under **Objectives** on the **Menu Bar**.
- **Description:** If desired, enter a more detailed description for the **Objective**.
- **Type:** Select either **Standard** or **Laser** from the drop-down menu.

Magnification Calibration

If the magnification has already been calibrated, the X and Y values will appear. If a new **Objective** is being added, position the objective and then proceed with the calibration procedure outlined in "Magnification Calibration" on page 39 to set the **X Micron/Pixel** and **Y Micron/Pixel** values.

NOTE:

Proper setting of the magnification is required for accurate display of Isotherms and accurate measurements. Calibration information is stored with captured images for future recall and measuring.

Laser Objectives

If the Type selected is Laser, two additional settings for Laser Target and Com Port will appear.

Laser Target

The Laser Target X and Y values represent the position of the laser target in pixels from the upper left hand corner of the live video window. This target position is linked specifically to the selected laser objective.

Com Port

The Com Port to which the laser controller is attached must be selected for the system to function properly. This may be done by selecting the Com Port from the drop down list or by pressing the Auto Detect button.

Chapter 6. Target Alignment and Magnification Calibration

6.1 Target Alignment

In order to perform the target alignment, basic knowledge of how to activate and fire the laser is required. Please see "Validation Mode" on page 51 for more information. The target alignment must be checked before the laser is used in Clinical or Multipulse mode.

To reduce risk of laser misalignment and damage to the cell, the laser's target alignment must be verified prior to use.

To maintain target alignment, the laser objective must be securely attached to the turret and the camera to the microscope. If the objective and/or camera are loosened and re-tightened, the target alignment should be verified prior to use.

RED-i Target Alignment

1. Gather the following tools:
 - HT Laser Target alignment slide: A Laser Target alignment slide (HT part #500997) is included with all LYKOS systems and available to purchase from Hamilton Thorne. The Laser Target is a specially-coated, gridded slide for verifying target alignment.
 - Plastic petri dish


Figure 6-1: HT Laser Target Slide

2. Place the Laser Target slide, face up, in a dry petri dish on the microscope stage.
3. Focus the image so the grid lines on the Laser Target slide are sharp. Position the slide to an unused black area.
4. From the Objective menu on the Menu Bar, select Align Target. The Align Target screen appears and Validation Mode is automatically enabled.
5. On the Validation panel, *deselect* Enable DTS.
6. Enable RED-i and using the RED-i Intensity Slider, set the intensity to a lower level for more precise positioning.

7. Set **Pulse Duration** to 10 - 100 microseconds.
8. Fire the laser. The laser spot shows as a clear (white) mark resulting from the localized vaporization of the special coating. If needed, adjust the illumination on the microscope to increase the size of the laser spot.
9. To align the **RED-i**, insert a 1.5 mm hex head driver into the upper **RED-i** adjustment screw. Rotate the screw to move the **RED-i** spot on the screen. The lower **RED-i** adjustment screw will move the spot in a perpendicular direction. Using both screws, center the **RED-i** spot in the center of the vaporized laser spot.
10. Confirm alignment: Reposition **Laser Target** slide to an unused area and fire laser. Laser spot should appear superimposed on RED-i spot.
11. If using the eyepieces, ensure that the eyepiece view is aligned by firing on the **Laser Target** again. While looking through the microscope eyepieces, the **RED-i** dot appears as a red disk. This red dot should be centered on the laser burn mark.

Figure 6-2: RED-i Adjustment.

NOTE:

IMPORTANT: Calibration and the resultant measurements are functions of the magnification of the objective lens through to the camera-imaging device. If the setup of the microscope changes with respect to the laser-objective or other magnification multiplier between the laser module and the camera [e.g. 1.0x – 1.5x], then the magnification must be re-calibrated.

NOTE:

If the RED-i adjustment screws are not facing to the outside while the objective is screwed onto the turret, the included thread spacer rings (see Figure 6-3) may be used between the turret and the objective. Each spacer ring equals one quarter (1/4) turn of the LYKOS. These may stay in place during normal operation.


Figure 6-3: Thread Spacer Rings for RED-i Alignment

6.2 Magnification Calibration

Proper setting of the Magnification Calibration is required for accurate measurements using the measurement toolbox feature and for accurate display of the Isotherm Rings. The computer calculates measurements by translating pixels into micrometers (μm). To ensure that this translation is accurate, the computer must be calibrated against a standard 100 micrometer scale.

Caution: *IMPORTANT: Calibration and the resultant measurements are functions of the magnification of the objective lens through to the camera-imaging device. If the setup of the microscope changes with respect to the laser-objective or other magnification multiplier between the laser module and the camera [e.g. 1.0x – 1.5x], then the magnification must be re-calibrated.*

Calibration Procedure

1. Position the objective to be calibrated and select it from the Objectives menu, or create a new objective under **Settings > Objective**.
2. Place a 100 micrometer scale on the microscope stage. Focus on the scale to obtain an extremely sharp image.
3. From the **Main Bar**, select **Objectives > Calibrate**. The Calibration command is also accessible from **Objective** settings panel. The **X Magnification Calibration** screen appears.


Figure 6-4: X Calibration Screen

4. Using the mouse, adjust the **yellow vertical calibration lines** so they are 100 microns apart. Make sure that the yellow lines are placed on the same side of the micrometer scale lines. Press **Accept** to save the X magnification value.
5. The **Y Magnification Calibration** screen appears. Repeat the procedure outline in Step 4 for the Y axis. Press **Accept** to save the Y magnification value.

6. The **Magnification Calibration** is now complete.
7. For future reference, note the camera type, objective lens identification and the calibration in the table below.

Camera Type	Objective Lens ID	Current Calibration (microns/pixel)

Using Measurement Tools to Confirm Calibration

1. After the calibration is complete, capture an image of the 100 micron scale.
2. Select the image from the thumbnail gallery.
3. Select the **Zona** measurement tool and press Start.
4. Using the mouse, measure the distance between the calibration lines used to set the calibration. The value should be 100 μm .

Chapter 7. Patient & Embryo Information

7.1 Patient Info

The following Patient Info may be entered and saved to the internal database:

- Patient ID
- First Name
- Last Name
- Date of Birth
- Physician

Patient Info	
Patient ID	
98333083	
First Name	Last Name
Leslie	Smylee
Date Of Birth	Physician
11/ 3/1981	Dr. Jack Smothers

Figure 7-1: Patient Info

After the Patient Info has been entered or selected, the Embryo information may then be entered.

7.2 Embryo Panel

Info

NOTE:

A Patient must be selected or entered prior to entering Embryo Info. Entering or selecting a patient after embryo details are entered will erase the Embryo Info.

The following Embryo Info may be entered and saved to the internal database:

- Embryo ID
- Retrieval Date
- Fertilization Method
- Cycle Number
- Embryos in Cycle
- Note

Figure 7-2: Embryo Info Panel

To save the entered embryo details to the selected patient record, press the Save button at the top of the panel.

A previously saved Embryo ID for the currently selected patient may be chosen from the drop-down list. This is a useful tool for recalling and measuring images at a later time.

Measurements

The Measurements panel provides access to six assigned measurement Tool Types including: Embryo, Zona, Drill, PN 1 and PN 2.

NOTE:

Measurements may only be performed on captured or saved images. The calibration information is included in the image header for future recall and measurement.

Accuracy of measurements is reliant on use of a calibrated objective. Please see "Magnification Calibration" on page 39 for more information.

1. Click on a Thumbnail image to open in Video Window.
2. Select a tool from the drop-down list.
3. Click Start to initiate the first measurement.
4. Position the cursor at the starting point of the measurement, and click and release the mouse button.
5. Position the cursor at the ending point of the measurement, and click and release the mouse button.
6. The measurement line is drawn, the measured value appears on the image and the value is inserted into the appropriate field in the Measurement panel.
7. Repeat Steps 4 & 5 until the required values are acquired.
8. To undo the previous measurement, click Undo. Both the line drawn and the value of the most recent measurement is removed.
9. To delete all measurements and start fresh, press Clear (this deletes all current measurements).
10. To save the measurements to the report, press Save.

	um
1	108.7
2	101.5
3	102.1
4	101.4
5	101.6
Mean	101.5
SD	0.5

Number of Blastomeres: 0

Figure 7-3: Embryo Measurements

**NOTE:**

Not all positions need to be filled with measurements. Empty positions are not included mean and standard deviation (SD) calculations.

**NOTE:**

Once Save is selected, the measurements on the image can only be cleared by closing the image and reopening it.

Tool Type**Embryo**

The Embryo measurement tool is used to measure the diameter of the embryo. Up to five measurements may be recorded. Also, the Number of Blastomeres present may be recorded. The mean and the standard deviation of the diameter measurements are also calculated.

Zona

The Zona measurement tool is used to record up to five measurements of the zona pellucida's thickness. To minimize variability and reduce bias, it is recommended that measurements be made in at least four locations: such as 12, 3, 6 and 9 o'clock. Up to five measurements may be recorded. The mean and the standard deviation of the entered measurements are also calculated.

Drill

The Drill measurement tool is used to measure the size of the hole made by the laser. Up to five measurements may be recorded. The mean and the standard deviation of the entered measurements are also calculated.

PN (Pronucleus) 1 and 2

The PN 1 and 2 measurement tools are used to measure the diameter of the two pronuclei. Two measurements may be recorded: D1 and D2. The mean and the standard deviation of the two diameter measurements are also calculated.

Images

To temporarily save a copy of an image with its measurement overlays into the **Thumbnail Gallery**, press the **Save** button in the **Image Toolbox**. In order to save to disk, the image must be saved using the **Thumbnail Gallery** menu options (see “Opening, Saving and Deleting” on page 43).

NOTE: *Images saved with overlays are added as new images in the Thumbnail Gallery; the original image files are not overwritten.*

To save an image to the report, select the image from the **Thumbnail Gallery** and then select one of the four **Images Report Buttons** on the **Embryo** panel. A copy of the selected image appears on the button to which it was attached.

Figure 7-4: Images Report Buttons

- To change an image saved to the report, select another thumbnail image and either press the button to replace the current image or right-click on the button and choose **Select** from the menu.
- To add overlays to an image, right-click on the button and select **Edit**.
- To remove an image from the report, right-click on the button and select **Delete**.

Data Fields

Located under **Data Fields** are 25 data input fields with customizable labels that may be added to the report. To change the label text, click on the label text, enter a new descriptor in the text box, and press **OK**.


Figure 7-5: Change Data Field Label and Persist

Checking the box next to **Persist** causes the input data to remain when the analysis data is cleared.

Chapter 8. Working with Images & Videos

8.1 Thumbnail Gallery

The **Thumbnail Gallery** at the bottom of the screen provides access to image and video capture controls.


Figure 8-1: Thumbnail Gallery

Thumbnail representations of captured images and videos appear in the **Thumbnail Gallery**, under the corresponding tab. The number of currently captured images and videos is also indicated.

8.2 Working with Images


NOTE: Prior to capturing and saving images, the **Image Capture** options under **Settings** should be configured (see "Image Capture" on page 23).

Select the **Images** button to open the **Thumbnail Gallery** of saved images.

Opening, Saving and Deleting

- **Open:** To open a saved image file, press the **Open** button and select the file. Only .jpg or .tif images may be opened. The image will be added to the **Thumbnail Gallery**. Double click on the thumbnail image to view the full image in the **Video Window**.
- **Save:** If **Autosave on Capture** is not enabled, images must be manually saved.
 - **Save Selected Image:** Select one image and press **Save > Save Selected** to store to disk.
 - **To Save All Images:** First, make sure that **Images** is the active gallery. Select **Save > Save All** to store all images in the **Thumbnail Gallery** to disk
- **Delete:**
 - **Delete Selected Image:** Select one image and press **Delete > Delete Selected** to remove it from the **Thumbnail Gallery**.
 - **To Delete All Images:** First, make sure that **Images** is the active gallery. Select **Delete > Delete All** to delete all images in the **Thumbnail Gallery**.

NOTE: Ensure that **Images** is the active Thumbnail Gallery before **Save All** or **Delete All** is selected.

NOTE: If an image has already been saved to disk, deleting it from the Thumbnail Gallery does not delete it from the disk.

Capturing Images

To manually capture an image, press the **Capture** button. A new thumbnail image will appear in the gallery.

Two additional image capture options are accessible from the **Thumbnail Gallery** menu from the drop down menu on the **Capture** button:

- **With Video Overlays:** When enabled, the **Target**, **Crosshairs** and **Scale Bar**, will be captured in addition to the live video image.
- **Capture on Laser Fire:** When enabled, images will be captured each time the laser is fired in **Clinical** or **Validation** mode.


Figure 8-2: Capture Options

To enable either option, click the item in the menu. A check mark indicates the item is enabled.

Adding Images to Report

To save an image to the report, select the image from the **Thumbnail Gallery** and then select one of the four **Report Images** on the **Embryo** panel. A thumbnail copy of the image appears on the panel and the full-size image is saved to the report database.


Figure 8-3: Report Images

- To change an image saved to the report database, select another image from the **Thumbnail Gallery** and select the **Report Image** to be replaced or right-click on the **Report Image** and choose **Select** from the menu.
- To add overlays to an image, right-click on the **Report Image** and select **Edit**. This opens the **Image Toolbox** (see “Image Toolbox” on page 45). Please note that when the image is saved with the overlays, it does not automatically replace the Report Image; this will need to be done as directed above.
- To delete a Report Image from the report database, right-click on the **Report Image** and select **Edit**. This only removes the image from the report; it does not delete it from the **Thumbnail Gallery**.

Viewing and Closing Images

- To view a full size image in the **Video Window**, double-click the image in **Thumbnail Gallery**.
- To close the full image, click the red **Close Window** icon in the upper right hand corner of the **Video Window**.

Image Toolbox

When an image file is opened in the Video Window, the Image Toolbox is available to add a variety of overlays to the image.


Figure 8-4: Image Toolbox

The Image Toolbox provides tools to add text, shapes, and lines to an image as well as a ruler to make measurements.

Add Text

1. Use the color palette to select the text color.
2. Press the Text icon.
3. Click on the image where the text should start and move the mouse to where the text should end and click again (do not hold the mouse button while drawing the text box).
4. In the Image Label box that appears, enter the text.
5. Select the font to use.
6. Select OK to add the text to the image.

Add Arrow, Line, Ellipse, or Rectangle

1. Use the color palette to select the outline color of the graphic.
2. Select the Thickness of the outline.
3. Press either the Arrow, Line, Rectangle, or Ellipse.
4. On the image area, click the mouse to start drawing.
5. Move the mouse to draw the graphic and click again to complete the shape (do not hold the button down while drawing the shape).

Make a Measurement

NOTE:

The measurements made using the Ruler tool are added only to the image and not incorporated into the report data. For measurements to be added to the report, the Measurement feature (see "Measurements" on page 40) must be used.

1. Use the color palette to select the outline color of the ruler line.
2. Select the Thickness of the line.
3. Press the Ruler icon.
4. On the image area, click the mouse to start drawing.

5. Move the mouse to draw the line and click again to complete the line (do not hold the button down while drawing the shape).
6. The calculated distance is displayed at the end of the ruler line.

Undo

Pressing the **Undo** button will erase the last action performed. To erase additional actions, press the **Undo** button again.

Redo

Pressing the **Redo** button will recreate the action or actions previously undone.

Add Image with Graphic Overlay to Thumbnail Gallery


Figure 8-5: Full Screen Mode: Image Toolbox

To add an image with graphic overlays to the **Thumbnail Gallery**, click **Save**. This adds the new image to the **Thumbnail Gallery** only. It does not overwrite the original image, nor does it save the image to disk. See “Opening, Saving and Deleting” on page 43 for details on saving images to disk.

Working in Full Screen View

- With an image file open, click the **Full Screen Icon** at the bottom right corner of the **Video Window**. The image will open in full screen view, and the **Image Toolbox** will appear at the upper left-hand corner. All other controls are hidden.
- To return to standard view, click the **Reduce Screen Icon** at the bottom right corner of the **Video Window**.

8.3 Working with Videos

Select the **Videos** button to open the Thumbnail Gallery of saved videos.

Opening, Saving and Deleting

- **Open:** To open a saved video file, press the **Open** button and select the file. Only .avi or .mp4 videos may be opened. The video will be added to the Thumbnail Gallery.
- **Save:**
 - **Save Selected Video:** Select one video and press **Save > Save Selected** to store to disk.
 - **To Save All Videos:** First, make sure that **Videos** is the active gallery. Select **Save > Save All** to store all videos in the Thumbnail Gallery to disk.
- **Delete:**
 - **Delete Selected Video:** Select one image and press **Delete > Delete Selected** to remove it from the Thumbnail Gallery.
 - **To Delete All Videos:** First, make sure that **Videos** is the active gallery. Select **Delete > Delete All** to delete all videos in the Thumbnail Gallery.

NOTE: Ensure that **Videos** is the active Thumbnail Gallery before **Save All** or **Delete All** is selected.

NOTE: If a video has already been saved to disk, deleting it from the Thumbnail Gallery does not delete it from the disk.

Capturing Real Time Videos

NOTE: Prior to capturing and saving real time or Time Lapse videos, the **Video Capture** options under **Settings** should be configured (see "Video Capture" on page 31).

- **To Start Recording:** Press the **Record** button to start recording the live video stream. A counter showing the **Time Elapsed** appears.
- **To Stop Recording:** Press the **Stop** button.

Capturing Time Lapse Video

NOTE: Prior to capturing and saving Time Lapse videos, the **Time Lapse** options under **Settings** should be configured (see "Time Lapse" on page 30).

- **To Start Recording:** Press the **Time Lapse** button to start recording the live video stream. A counter showing the **Time Elapsed** and **Completed %** appears.
- **To Stop Recording:** Press the **Stop** button.

Playing Videos

Standard video controls are available to play the videos.

Figure 8-6: Video Playback Controls

- **To Play a Video:** Double-click on the video thumbnail. The video will open in full screen mode and start playing.
- **To Close a Video:** Either press Stop or close the window.
- **To Pause a Video:** Press the Pause button.
- **To Rewind a Video:** Press the Back button to jump back 10 seconds.
- **To Fast Forward a Video:** Press the Forward button to jump 10 seconds forward.

Chapter 9. Laser Modes and System Operation

The software features three modes of laser operation: **Clinical**, **Validation** and **Multipulse**. Depending on which mode is activated, the available controls will differ.

The mode is activated by selecting it from the **Laser Mode Menu** list from the **Menu Bar**.


Figure 9-1: Laser Mode Menu

9.1 Working in Full Screen View


Figure 9-2: Full Screen View: Clinical Mode

- Click the **Full Screen Icon** at the bottom right corner of the **Video Window**. The live image will open in full screen view, and the activated **Laser Energy Panel** will appear at the upper left hand corner. All other controls are hidden.
- To return to standard view, click the **Reduce Screen Icon** at the bottom right corner of the **Video Window**.

9.2 Clinical Laser Mode

To activate **Clinical** mode, select **Clinical** from the **Menu Bar** drop-down list.

The **Clinical Laser** panel houses the controls to select one of three preset laser energy levels (**Low**, **Medium** or **High**) and fire the laser. Upon entry in the program, the laser should automatically be powered **Off** (**Off** button highlighted) and **Laser Fire** button grayed out.


Figure 9-3: Laser Energy Level set to Off (left) and Medium (right)

Activating the Laser Energy Level

The laser energy level (**Low**, **Medium** or **High**) must be manually activated before firing the laser. These energy levels correspond to the **Low**, **Medium** or **High** settings entered on the **Settings > Laser Settings** screen (see “Clinical Pulse Settings” on page 28).

To activate the laser, select the **Low**, **Medium** or **High** button with the mouse. Once a level is selected, the **Fire** button becomes active and a 15 second countdown timer appears.

NOTE:

The laser target selected under Laser settings (see “Target Type” on page 25) will only appear after the Low, Medium or High button is pressed.

The laser may be manually deactivated by selecting the **Off** button.


Firing the Laser

Caution:

Before firing the laser, make sure that the target has been properly aligned to the laser spot as described in “Target Alignment” on page 35.

After the laser energy level is selected, the **Laser Fire** button becomes active and the laser may be fired either by pressing the **Laser Fire** button or using the foot switch.

The laser must be fired within 15 seconds of being activated or it will automatically power **Off**. In addition, the number of times the laser may be fired within the 15 seconds is determined by the level selected:

- **Low:** laser may be fired up to 15 times within the 15 seconds.
- **Medium:** laser may be fired up to 10 times within the 15 seconds.
- **High:** laser may only be fired once within the 15 seconds.

9.3 Validation Mode

Validation mode is not for use on human embryos.

To activate Validation mode, select Validation from the Menu Bar drop-down list.

Configuring Validation Mode

On the Settings> Laser Settings screen, select the Pulse Resolution (see "Pulse Resolution" on page 28). The Pulse Slider Resolution is set to 1 μ s, 10 μ s or 100 μ s and relates to the adjustment increments of the laser pulse on the Validation Control Panel.

Validation Control Panel

When Validation Mode is activated, the Pulse Duration (μ s) and Power (%) of the laser may be adjusted manually from the control panel.


Figure 9-4: Validation Control Panel

Setting the Power

The Power may be set from 1% to 100% by using the adjustment arrows or by typing a new value in the field.

Setting the Pulse

The Pulse may be set from 1 μ s to 3000 μ s by using the slider, the adjustment arrows or by typing a new value in the field.

NOTE: The Pulse cannot be increased above 3000 μ s on the Validation Mode Control Panel. If a higher value is entered in the Pulse field, the Pulse value automatically reverts to 3000 μ s. The built-in delay between pulses is directly related to the pulse value, in order to maintain Class I standards.

Firing the Laser

Before firing the laser, make sure that the target has been properly aligned to the laser spot as described in "Target Alignment" on page 35.

The laser may be fired either by pressing the Laser Fire button on the Validation Control Panel or using the foot switch.

9.4 Multipulse Mode

The Multipulse Mode allows the use of a series of laser pulses. To activate Multipulse mode, select Multipulse from the Menu Bar drop-down list.

Multipulse Mode is for use with trophectoderm biopsy procedures only.

To prevent excessive heating, always target the laser away from the inner cell mass. The laser should not be focused closer than 60 μm to the inner cell mass. To reduce risk of targeting the inner cell mass, ensure the manipulators are correctly oriented so only the trophectoderm cells are exposed to the laser pulses during the biopsy procedure.

Multipulse Laser Settings

From Settings > Laser set the Enable Duration and Mean Power.

- **Enable Duration (Sec):** The Enable Duration is the length of time during which, once enabled, the Laser Fire Button will remain active allowing the laser to be fired. The maximum allowable time is 600 seconds (10 minutes). This allows operators time to position the embryo for laser application, without the need to relinquish control of the micromanipulators to reactivate the laser.
- **Power Limit (mW):** The Power Limit is the maximum mean power of the pulse train during multipulse application. The maximum allowable setting is 10 mW. The Power Limit is used to determine the interdependent values used during multipulse laser application (see "Built-in Maximum Limits for Multipulse Mode" on page 54).

Figure 9-5: Multipulse Laser Settings

Multipulse Laser Panel

The Multipulse panel provides controls to customize Power, Pulse/Sec, Total Duration and Pulse Duration.


Figure 9-6: Multipulse Control Panel

- **Power:** Percent maximum laser power to be used (maximum 100%)
- **Pulses/Sec:** The number of laser pulses to be fired per second (maximum 1000)
- **Total Duration:** The time length in seconds that the pulses will fire in multipulse mode (maximum 10 seconds)
- **Pulse Duration:** The time length in microseconds of each laser pulse (maximum 3000 microseconds)
- **Enable Button:** The Multipulse mode is enabled only after the Enable button has been clicked, which starts the timer based on the Enable Duration (see "Multipulse" on page 28). The Enable button will change to a Disable button and show the remaining seconds on the timer. Once the timer expires, the Enable button must be pressed again.
- **Halt Button:** The Halt button is an emergency override that immediately stops the multipulse mode train of laser pulses, even if the foot switch is depressed. It is not intended to be used for normal operation.

**NOTE:**

The settings of Power, Pulses/Sec, and Pulse Duration are interdependent and settings are restricted to be less than or equal to the Mean Power. Please see "Built-in Maximum Limits for Multipulse Mode" on page 54 for a detailed explanation.

Firing the Laser

When in Multipulse mode, the laser can only be fired by using the foot switch.

When depressed, the foot switch starts the Multipulse laser fire mode. When the foot switch is released, the Multipulse mode is halted.

The laser pulsing will automatically halt at the set Total Duration value. As long as the timer is still active, the Multipulse mode can be restarted (after a blanking duration equivalent to previous total pulse train duration) by releasing and depressing the foot switch.

To deactivate the laser before the Enable Duration (see "Multipulse Laser Settings" on page 52) is reached, press the Disable button on the Multipulse panel.

Built-in Maximum Limits for Multipulse Mode

There are maximum limits on each of the four Multipulse settings, as indicated in the table below. In addition, the **Power Limit** of the pulse train is limited to 10 mW.

Table 9-1: Max Limits on Multipulse Settings

Settings	Maximum Limit
Power	100 (300 mW)
Pulse Duration	3000 microseconds
Pulses per Second (PPS)	1000 PPS
Total Duration	10 seconds

The three values of **Power**, **Pulse Duration** and **Pulses per Second** are interdependent. The product of these three settings must be equal to or less than the **Maximum Mean Power**.

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The built-in algorithm permits the user to set any two of the three inputs to any value, provided it is lower than the maximum limit (see Table 9-1). However, the third value will be restricted to a value so that the **Power Limit** value is not exceeded.

For example, if the **Power Limit** is set to 10 mW, and the **Power** is set to 100% and the **Pulse per Second** is set to 200, then the **Pulse Duration** will be limited to 167 μ s or less.

Chapter 10. Dynamic Targeting System

10.1 Introduction

The Dynamic Targeting System (DTS) is an option for the LYKOS laser system which enables a computer-controlled targeting system for precise targeting of the laser. In addition to single-shot positioning, it allows selection of a continuous path (line, curve, freehand or rectangle) or a multi-point arrangement along which the laser will fire.

10.2 Software Overview

The DTS option may be enabled in both **Clinical** and **Validation Modes** (not for use with human embryos).


Figure 10-1: Clinical DTS and Validation DTS Panels



- In **Clinical Mode**, only Medium and Low laser levels may be used in DTS mode.
- In **Validation Mode**, the pulse/hole size, laser power, and spacing of the holes may be set (Validation Mode is not for use with human embryos).

When the **DTS** option is disabled, the LYKOS operates as a fixed laser. Also, the **DTS** option is not available under **Multipulse Mode**.

Path Type

When **DTS** is enabled, a variety of paths may be selected and drawn on the screen. The mouse is used to draw the path or select the points. Please see "Working with Path Types" on page 61 for more details.

Table 10-1: Laser Path Types

 Single Shot	Single Shot: The user may click on any area of the field of view to position the laser target and initiate a single laser shot upon firing.
 Line / Curve	Line/Curve: A line of any length may be drawn in the field of view. This line may then be selected and "dragged" to form an arc.
 Freehand	Freehand: A freeform path may be drawn on the field of view, and then edited as needed by dragging points to reposition
 Multipoint	Multipoint: Individual target spots in the field of view may be designated.
 Rectangle	Rectangle: A rectangle may be drawn and edited as needed.

10.3 Initialization and Verification

Objective Calibration

The LYKOS objective must be calibrated prior to using the laser with **DTS** mode. Please refer to "Magnification Calibration" on page 39 for more detail.

RED-i Target Alignment to Laser

The **RED-i** target must be aligned to the laser spot prior to using the **DTS** function. Please refer to see "RED-i Target Alignment" on page 35. After the **RED-i** is aligned, the **Initialization** and **Verification** routines (see below) must be performed.

Initialize and Verify

Caution: To ensure proper operation, Initialization and Verification should be performed every 24 hours.

When in DTS mode, the RED-i target and laser are programmed to match the target path drawn on the screen. During Initialization, the system learns and maps the position of the RED-i target in relation to the camera's field of view. During Verification, the system then verifies that both the RED-i target and the laser precisely target 15 fixed coordinates in the mapped field on the display.

Initialize DTS

Caution

To ensure proper operation, the microscope light **MUST** be turned OFF before Initialize DTS button is selected.

1. Turn OFF microscope illumination.
2. Select Validation mode.
3. Ensure that the LYKOS objective is in place.
4. On the Dynamic Laser Info panel, select Initialize DTS.

Figure 10-2: Dynamic Laser Info Panel

5. A reminder to turn off illumination on microscope appears. Select OK to continue or Cancel to exit. The initialization screen appears, and the system begins the mapping process.
6. After 2 to 3 minutes the initialization process completes and Verify DTS can be performed.

Verify DTS

1. Turn ON microscope illumination.
2. Place the HT Laser Target slide (HT part # 500997) face-up in a dry Petri dish on the microscope stage.
3. Position the stage to the HT Laser Target slide and focus on a grid line.
4. Move the stage so that a fully black area of the HT Laser Target slide is in view.
5. Select Verify DTS.
6. The Verify DTS screen consisting of 15 target points will appear. The RED-i illuminates at each location and the laser automatically fires at each target.

Figure 10-3: Verify DTS Screen

7. The RED-i and laser spot should align within the limits of accuracy for the system.

- a. For each target where the RED-I and laser spot are within agreement, a green X is shown.
 - b. If either the RED-I or laser spot are outside the acceptable range, a red X appears.
8. Once the verification process completes, the Verification Results are shown:

Figure 10-4: Verification Results

9. If less than 15/15 Passed, the initialization process must be performed again.

To review the Verification results in report format see "IVQC Report" on page 62.

10.4 Using DTS Option

Clinical Mode

1. Select Clinical from the Main Menu.
2. Select Settings > Laser > Laser Settings.
3. Under Clinical Pulse Settings, set the desired Medium (200 μ s to 600 μ s) and Low (0 μ s to 200 μ s) pulse values. For reference, the estimated Hole Size, in μ m, for the given Pulse is shown. High Pulse is not available in DTS mode.

Clinical Pulse Settings		
	Pulse (μ s)	Hole Size (μ m)
High	600	14
Medium	400	11
Low	200	8

Figure 10-5: Clinical Pulse Settings and Hole Size

4. Make sure that Blanking Duration is set to None.
5. Save settings and exit.
6. Select DTS Enabled on the Clinical Control Panel.


Figure 10-6: Clinical Mode with DTS Enabled

NOTE: When DTS is enabled in Clinical Mode, only the Medium and Low laser power may be used.

7. Select the touch-key under **DTS Enabled** to open the drop-down options and select the target path type.


Figure 10-7: DTS Path Types

8. Select pulse setting: **Medium** or **Low** (**High** is disabled in DTS mode).
9. Using the mouse, draw the intended target path on the screen (see “Working with Path Types” on page 61).
10. When the path is acceptable and the embryo is positioned, press **Enable Laser**. A 15 second count-down timer begins, signifying the time interval the laser is enabled and may be fired. A **STOP** button appears.
11. To fire the laser, press **Fire** within the 15-second interval.
12. To terminate the laser at any point, press **STOP**.

Caution: The **STOP** function terminates the laser and programmed path. Rearming and firing the laser will reset the target path to the beginning. The target path will remain enabled until cleared.

13. To remove the path from the screen, select **Clear Path**.

NOTE:

The programmed target path may be used for subsequent procedures and remains active until cleared.

Validation Mode

1. Select **Validation** from the Main Menu.
2. Select **Settings > Laser > Laser Settings**. Under **Pulse Resolution**, select either 1, 10, or 100. This controls the interval that the **Pulse** value increases or decreases when the **Pulse Slider** is moved.
3. Make sure the **Blanking Duration** is set to **None**.
4. Save settings and exit.
5. Select **DTS Enabled** on the Validation Control Panel.


Figure 10-8: Validation Control Panel with DTS Enabled

6. Select the touch-key under **DTS Enabled** to open the drop-down options and select the target path type.
7. Set the **Pulse** and **Power** (it is recommended that 100% Power is always used).
8. Using the mouse, draw the path on the screen and adjust the spacing as needed (see "Working with Path Types" on page 61).
9. When the path is acceptable and the embryo is positioned, press **Fire**.
10. To terminate the laser at any point, press **STOP**.

The STOP function terminates the laser and programmed path. Firing the laser again will reset the target path to the beginning. The target path will remain enabled until cleared.

11. To remove the path from the screen, select **Clear Path**.

NOTE:

The programmed target path may be used for subsequent procedures and remains active until cleared.

Working with Path Types

Single Shot

To use **Single Shot**, click the mouse on the screen to position the target. This is where the laser will fire when activated.

NOTE:

Under Single Shot, if the target type is set to Circle under Laser Settings, this will not be representative of the hole size. In order to see the Hole Size, Isotherm Rings must be the target type and the 50° circle activated (see "Isotherms" on page 26).

Multipoint

Click the mouse on the image area to position a target. Continue adding up to 10 target points as needed.

1. Each point may be selected and the **Pulse** and **Power** for that point may be adjusted. Thus it is possible to have 10 targets reflecting varying settings. To select individual target points, click on the target using the left mouse button. The target will turn green and the settings may be adjusted. To change position, click and hold the left mouse button and drag the target to a new location.
2. To change all targets, click any target using the right mouse button. All targets will turn green. Adjust settings from the control panel. All targets will simultaneously reflect any changes in settings. To reposition, click and hold the right mouse button and drag the group of targets to a new location.

Line/Curve, Freehand and Rectangle

If drawing a **Line/Curve** or **Rectangle**, the form is restricted. Using **Freehand** mode, any type of shape may be drawn.

To draw a path, position the mouse at the start of the path, then press and hold the left or right mouse button and drag the mouse. When finished drawing release the mouse button. The line turns into a series of red circles, with the circle size determined by the current laser settings.


Figure 10-9: Draw Line

After the path is drawn, it may be adjusted in several ways:

1. Left-click the path to select it; a selected path will turn green (see Figure 10-10a).
2. In **Clinical Mode**, only the selected pulse (**Low** or **Medium**) may be changed.
3. In **Validation Mode** (not for use with human embryos), both the **Pulse** and **Power** may be adjusted. The circles will increase or decrease in size and number as changes are made. In addition, the **Spacing** may be adjusted to achieve the desired overlap (see Figure 10-10b & c).

4. Position the mouse over the path so it turns into a hand icon, click and hold the left mouse button and drag to adjust the shape/size. Either the middle or ends of the line may be repositioned (see Figure 10-10d).
5. Position the mouse over the path so it turns into a hand icon, click and hold the right mouse button to reposition the entire path. When done, release the mouse button.

Once the path is acceptable, click the mouse anywhere in the image area. The path will turn red and the laser may then be fired.


Figure 10-10: Manipulation of Laser Path

10.5 IVQC Report

To view recent and historical results from the **Initialize DTS** and **Verify DTS** routines in report format, select the **IVQC Report** button. The report screen shown in Figure 10-11 opens.


Figure 10-11: IVQC Report Screen

IVQC Report Screen Layout

The IVQC Report screen consists of three main areas:

1. Report Panel
2. IVQC Trend Report Toolbar
3. Save and Print Toolbar

Report Panel

When opened from the Dynamic Laser Info panel, the IVQC Report panel presents the following quality control data (based on the current date and time):

- Last Initialization Date
- Most Recent Verification
- Laser Offset and RED-i Offset data for the Most Recent Verification
- 4. IVQC Trend Chart for the past month of Verification time points

Laser and RED-i Offset Data

The Laser and RED-i Offset Data is presented in both numerical and graphical formats.


Figure 10-12: Laser and RED-i Offset Data

The data displayed in both charts always represents the data obtained for the Most Recent Verification, based on the End Date of the active date range ("IVQC Trend Report Toolbar" on page 64).

During the Verification process, the Laser Offset and RED-i Offset from the target are measured across 15 points (see "Verify DTS Screen" on page 57). In the numerical results, the Mean and Standard Deviation are presented along with the Max Offset detected across all points.

In addition, the Pass/Fail Status is given. A Max Offset exceeding 4µm results in a Fail.

The **Verification** graph plots the measured **Laser Offset** and the **RED-i Offset** across all 15 points. The **Laser Offset** is plotted in blue and the **RED-i Offset** is plotted in red.

IVQC Trend Chart

The **IVQC Trend Chart** shows the **Mean Offset** for the **Laser** and **RED-i** for each **Verification** performed over the selected time frame ("IVQC Trend Report Toolbar" on page 64).


Figure 10-13: IVQC Trend Chart

IVQC Trend Report Toolbar

The **IVQC Trend Report Toolbar** allows the user select a date range to view a **Trend Chart** of stored **Verifications**. By using the **Date Range**, **Start Date/Time** and **End Date/Time**, the **Verification Data** and **IVQC Trend Chart** may be viewed for different time frames.


Figure 10-14: IVQC Trend Report Toolbar.

NOTE:

The Last Initialization and Most Recent Verification are based on the selected End Date/Time.

The **IVQC Trend Report Toolbar** also provides access to designing and editing the layout of the report. The process of editing and designing reports is discussed in "Working with Reports" on page 65.

Save and Print Toolbar

The displayed **IVQC Report** may be saved as a PDF or TIFF file or printed. Using the **Fit Width** toggle, the report may also be viewed so it fills the **Report Panel**.

Chapter 11. Reporting & ASCII

11.1 Working with Reports

Select **Reports** on the **Main Menu** to open the **Reports** window. The **Reports** function allows the creation, copying, editing, viewing, printing or deleting of reports.

Figure 11-1 below illustrates the capabilities of the report by showing the same report in both the **Design View** (showing field placeholders) and the **Report View** (showing actual data).

New Valley Fertility Clinic
123 Elm St. Beverly, MA 01915
978-555-1212

Embryology Report
Tuesday, December 26, 2017 2:52 PM

Patient Information

ID: Patient Birthdate:
 First Name: Last Name:
 Patient Physician:
 Embryologist:

Embryo Information

ID: Retrieval date:
 Cycle Number: Embryos in Cycle:
 Fertilization Method: Note:

Embryo Measurements

Date:
 Number Of Blastocysts:

	Embryo	Zone	Drill	PN 1	PN2
1					
2					
3					
4					
5					
Mean					
SD					

Footer Page #

New Valley Fertility Clinic
123 Elm St. Beverly, MA 01915
978-555-1212

Embryology Report
Tuesday, December 26, 2017 2:42 PM

Patient Information

ID: 96773993 Patient Birthdate: 11/13/1981
 First Name: Leslie Last Name: Smyth
 Patient Physician: Dr. Jack Smothers
 Embryologist: J.B. Hubbs

Embryo Information

ID: 4.11172017 Retrieval date: 11/13/2017
 Cycle Number: 1 Embryos in Cycle: 6
 Fertilization Method: KCSI Note:

Embryo Measurements

Date: Friday, November 17, 2017 1:20 PM
 Number Of Blastocysts: 0

	Embryo	Zone	Drill	PN 1	PN2
1	154.3	12.3	9.7	11.7	13.2
2	100.5	12.8	9.8	12.5	13.3
3	101.6	11.8	9.7		
4	106.4	12.3	10.5		
5	99.2	13.1	10.6		
Mean	102.4	12.5	10.1	12.1	13.2
SD	2.9	0.5	0.3		

1

Figure 11-1: Design View and Report View

Report Template

Report Templates provide access to the design and customization of new or existing reports and the viewing of reports. Available functions are:

- **New Template:** Opens the **New Report** dialog box. Enter the name of the **New Report** and press **OK**.
- **New Folder:** Creates a new folder within **Report Templates** panel.
- **Copy:** Opens a **New Report** dialog box. Enter a **New Report** name and select the report to be copied.
- **Delete:** Select the template to be deleted from the **Report Template** panel and press **Delete**.
- **Design:** Select the template to be edited from the **Report Template** panel and press **Design**.
- **View:** Select the report to view from the **Report Template** panel and then press the **View** button.


Figure 11-2: Reports Template Panel

Viewing Reports

1. From the **Report Templates** panel, select the report to be viewed and press **View**, or double-click the **Report Name**.
2. From the header section of the **Report Filter** panel, initiate a query based on one or more of the following criteria:
 - a. All Patients
 - b. Patient ID
 - c. First and Last Names
 - d. Start and End Dates
3. Press **Refresh** to view a list of records meeting the search criteria.
4. Select a record within the resulting query table to view the report.


Figure 11-3: Report Query

Designing Reports

The **Design** option permits customization of the pre-designed forms or creation of entirely new forms. The user-friendly, “drag and drop” designer gives you complete control over the look and content of the report. Select **New Template** to open a blank report design window.


Figure 11-4: Blank Report Template

Design Area

The design area is divided into three parts, **Header**, **Body** and **Footer**. Items included in these sections are printed out as a report. Any input or output data may be added to the report. Free-form fields also allow the inclusion of non-analysis data such as contact information or company logo.

The menu items above the design area provides access to basic commands used in designing reports.


Figure 11-5: Report Design Menu

Adding Items to the Report

There are two methods for adding data fields to the report design area:

- Select the header, body or footer section of the report and then double click on the **Standard** or **Data Item** to be added. The field appears at the upper left corner of the selected section.
- Click on the item, hold the left mouse button and drag the field to the design area (Click Hold + Drag).

Please note that adding an item only brings the placeholder for the actual variable. Labels for the variables must be created separately using the **Toolbox/Text Label** tool (see "Standard Items" on page 68.)

Report Layout Aids

As fields are being dragged within the design area, blue alignment guidelines appear to help align and position the fields. Also, if you select items and right click, additional layout controls are available (i.e., center align, send to back, bring to front).


Figure 11-6: Alignment Guide

Toolbox

The Toolbox is divided into Standard Items and Data Items, which are used to add data results to the report.

See “Report Property Descriptions” on page 69 for information on setting various properties for the Standard and Data Items.

Standard Items

- **Pointer:** Standard mouse pointer
- **Image:** Insert an image into the report
- **Text Label:** Add a text label or text field
- **Report Date & Time:** Insert the date and time the report was generated
- **Page Number:** Insert a page number for multi-page reports

Data Items

The Data Items toolbox area consists of all data from the system’s CASA database. These items can be click-dragged into the desired position on the report. In the actual report, the corresponding value of that variable will appear.

Properties

Quick Properties Tool

When a Standard or Data Item is inserted into the report, the Properties may be set from the Properties Toolbox or the Quick Properties Tool located on each field box frame. This provides a shortcut to enter parameters that will design the item (e.g., select image, formatting, alignment, etc.).


Figure 11-7: Quick Properties Dialog

Properties Panel

The Properties toolbox area lets you format the report layout, **Standard Items** and **Data Items**. The variables listed under Properties will change depending on what is selected in the design window. For detailed descriptions of the various property values, please see “Report Property Descriptions” on page 69.

11.2 Report Property Descriptions

General

Report Layout

To view the **Report Properties**, click on the “light gray” area between the header and body or body and footer in the design area.


Figure 11-8: Report Properties

- **Appearance**
 - **Font:** Controls the typeface, size, and style of the default font.
- **Design**
 - **(Name):** The tag applied to the report. This can be left as is or changed.
- **Layout**
 - **Size**
 - *Width:* Report width based on the set margins.
 - *Height:* Report width based on the set margins.
- **Misc.**
 - *Landscape:* Determines paper orientation. Set to false to print standard portrait orientation.
 - *Margins:* Set desired margins.
 - *PaperType:* Select desired paper type (e.g., letter, legal, A3, A4, etc.).

Header, Body and Footer Properties

From within the design window, select either the header, body or footer section of the report. The default font and layout properties appear.


Figure 11-9: Header, Body, Footer Properties

- **Appearance**
 - **Font:** Controls the typeface, size, style of the default font.
- **Design**
 - **(Name):** The tag applied to the selected item. This can be left as is or changed.
- **Layout**
 - **Size**
 - *Width:* Based on the report width.
 - *Height:* Maximum header and footer height is 125. These may be adjusted by entering new values or by dragging the header or footer border. Body height is determined by report height minus the header and footer heights.

- **Misc** (header and footer properties only)
 - *ShowOnFirstPage*: Select true to show on first page.
 - *ShowOnLastPage*: Select true to show on last page.

Image Properties

Two types of images may be added to the report:

- Under **Standard Items**, the **Image** control allows the import of any compatible image.
- Under **Data Items**, inserting the **embryo_image_path_X** (1 through 4) fields pulls the saved **Report Images** (see “Images” on page 42) into the final report.

When the image box is selected, the options under Properties reflect Image settings.


Figure 11-10: Image Properties

- **Appearance**
 - *Image*: Browse to locate the image file. Compatible image files are: .bmp, .jpg, .png, .gif, .ico, .emf and .wmf. Image files are embedded in the report template.
- **Behavior**
 - *SizeMode*
 - *Normal*: inserts image into box with no resizing or repositioning.
 - *StretchImage*: stretches image to fit the current image box (non-proportional).
 - *AutoSize*: adjusts the size of the image box to fit the image.
 - *CenterImage*: positions image to the center of the box.
 - *Zoom*: proportionally resizes image to fit the image box.
- **Design**
 - *(Name)*: The tag applied to the selected item. This can be left as is or changed.

- **Layout**

- *Location*: The X,Y location of the upper left corner of the image box. This may be changed by adjusting the numerical value or by dragging the box to a new location.
- *Margin*: The "snap-to" margin between the image box and other control boxes.
- *Size*: The size of the image box. This may be changed by adjusting the numerical value or by dragging the box borders.

- **Misc.**

- *DBPath*: This field is automatically populated when the image is added to the report.

Text Properties

The Text Properties apply to the following:

- Text Labels
- Report Date & Time
- Page #
- Data Items

Common text box properties are shown below:

Properties	
Appearance	
BackColor	Window
Font	Arial, 9.75pt
ForeColor	WindowText
TextAlign	Left
Barcode	
Print As Barcode	False
Behavior	
Multiline	False
Design	
(Name)	ReportDataBox89
Layout	
Location	164, 151
Margin	3, 3, 3, 3
Size	100, 22
Misc	
DBPath	Report Query 1: Human Clinical Embryo Re
Precision	2

Figure 11-11: Text Properties

- **Appearance**

- *BackColor*: The background color of the text box.
- *Font*: Controls the typeface, size, style of the default font.
- *ForeColor*: The color of the text.
- *Text*: The text that appears in the text box in the design layout. For Text Labels, this will also appear on the final report. For Data Items, the text will be replaced with the correct value upon viewing report.
- *Text Align*: Left, right, center.

- **Barcode**
 - *Print as Barcode*: To print the text as a barcode, set to true. To use standard font, set to false. When using barcodes, make sure the box is long enough to accept the complete code. This may take some trial and error.
- **Behavior**
 - *Multiline*: To extend printing to additional lines within the selected text box, set to true.
- **Design**
 - *(Name)*: The tag applied to the selected item. This can be left as is or changed.
- **Layout**
 - *Location*: The X,Y location of the upper left corner of the text box. This may be changed by adjusting the numerical value or by dragging the box to a new location.
 - *Margin*: The "snap-to" margin between the text box and other control boxes.
 - *Size*: The size of the text box. This may be changed by adjusting the numerical value or by dragging the box borders.

The Report Date & Time Standard Items and the analysis_date_time Data Item also have the following Misc options:

- **DateFormat**
 - Year_Month_Day: 2013-05-06
 - Month_Day_Year: 05-06-2013
 - Day_Month_Year: 06-05-2013
 - Long: Monday, May 6, 2013
 - Short: 5/6/2013
- **DisplayFormat**
 - Date: Shows date only
 - Time: Shows time only
 - Date_Time: Shows date then time
 - Time_Date: Shows time then date
- **TimeFormat (the outputs may vary based on geographic locale)**
 - TwelveHour: 2:38:29 PM
 - Twentyfour Hour: 14:38:29
 - Short: 2:38 PM
 - Long: 2:38:29 PM

NOTE: The Short and Long date and time formats are set through the Windows operating system and may differ from the above.

- **Misc.**
 - *DBPath*: The database path for the selected Data Item. This is for reference purposes only. This field does not apply to Standard Items.
 - *Precision*: Number of decimal places to display for floating point data. This setting is only available for Data Items with numerical output.

11.3 Working with ASCII Data

The **ASCII** section provides ability to manually select data fields for ASCII file output. All data stored is available for inclusion. This allows you to filter data selectively and create multiple ASCII files based on the same data set.


Figure 11-12: ASCII Filter Panel

- From the header section of the **ASCII Filter** panel, initiate a query based on the following identifiers:
 - Patient ID
 - First Name
 - Last Name
 - Start and End Dates
- Press **Refresh** to view a list of records meeting the search criteria.
- Under the **Fields** section are the list of variables available which may be included in the ASCII output.
 - **Select All:** Will check all the variables that will go in the report.
 - **Clear Selection:** Will clear all the checked variables.
 - Typing the desired variable in the **Fields** box will filter variables from the field.
 - Click the check box of variables to be included in the output.
 - As these items are checked, a table will appear under the **Sample Data** section listing their values under their output names.
- Under the **Sample Data:**
 - Click the drop-down list and select the type of delimiter to be used: **Comma**, **Semicolon**, **Space** or **Tab**.
 - Set the data qualifier to **None**, **Double Quote** or **Single Quote**.

- Click **Refresh** to reload the current changes made to the sample data. To have the data refresh automatically, check the **Auto Refresh** box.

Click **Save ASCII** to save the ASCII table in the command line. The **Save As** dialog box appears, prompting you to provide a name and folder location to save the file. The file is saved in **.tab** format.

Chapter 12. Laser Applications

12.1 Getting Ready

Check Target Alignment

Using an inked test slide, check the alignment of the target and RED-i spot to the laser beam. Please see "Target Alignment" on page 35 for detailed instructions.

Focusing on the Embryo

Place a culture dish containing the embryo on the microscope slide. Using a low magnification objective (e.g., 4x) locate the embryo in the media droplet and then switch to laser objective. Focus the image at the largest diameter of the embryo.

TIP: *To determine precise focus, observe the periphery of the embryo. Using bright field illumination, clear or dark halos may be visible surrounding the inner periphery of the zona during focusing. The clear or dark halos are seen when the focus is adjusted for one hemisphere or the other, above or below the equator of the embryo. When the proper focus is obtained, neither halo can be seen.*

RED-i

Enable RED-i on the Control Panel and adjust the brightness of the RED-i spot using the Intensity Slider.

12.2 Laser-assisted Hatching (LAH)

Multiple openings around the zona pellucida may lead to embryo loss during hatching and/or abnormal embryo development. Follow laboratory Hatching Protocol and best clinical practices.

Laser Positioning

Using either the on-screen target or RED-i spot, position the embryo using either the microscope stage controls or holding pipette so the laser beam is tangential to the zona pellucida.

To minimize the possible heating of adjacent blastomeres during treatment, the laser beam should be directed toward a section of the zona pellucida where the adjacent perivitelline space is widest or next to an area of fragmentation, when visible.

Laser Zona Drilling Options

Full Zona Breaching

The laser beam should be delivered in such a way that the axis of the drilled trench intersects the zona at the site of the desired perforation. Only a single opening should be formed; however, several laser shots may be used to produce or enlarge a single opening.

For LAH applications, the hole should be wide enough to breach the zona pellucida (as seen at the focal plane, when the embryo is observed at its largest diameter). The laser treatment may or may not lead to the perforation of the perivitelline layer.

NOTE:

The optimal diameter of a hole is determined by the thickness and/or hardness of the zona. Holes of wide diameter may be necessary to accommodate a thick zona. Holes of smaller diameter may be preferred for a thin zona.

Zona Thinning

Some users prefer to ablate only a *fraction* of the zona thickness, yet across a wider area. This process is commonly referred to as zona thinning. Zona thinning does not open the perivitelline space.


Figure 12-1: LAH: Full Zona Breach (left); Zona Thinning (right)

Clinical Laser Energy Level

The three preset laser energy levels (**Low**, **Medium**, and **High**) should permit the operator to perform zona drilling on a wide range of embryos.

- **Low:** use to perforate a thin zona ($<10\ \mu\text{m}$) or as a first attempt for drilling. Low pulse duration minimizes exposure and risk to embryos.
- **Medium:** commonly used for drilling the zona of most embryos.
- **High:** recommended only for drilling particularly thick ($>17\ \mu\text{m}$) or hard zona, for example as occurs with frozen/thawed embryos.

Recommendations for Embryos

Control over the diameter of the trenches created is achieved using the three laser pulse durations indicated above. The intensity of the laser is fixed at 100% power in order to minimize pulse duration required and therefore minimize heat conduction into neighboring cells. A large trench can be obtained in the zona with one single laser pulse using the **High** laser pulse. However, the invisible thermal effect of these laser parameters could potentially result in damage to the blastomeres adjacent to the hole. It is recommended, for safe

operation of the laser with most embryos, to activate the laser several times using the Low or Medium laser pulse and enlarge the hole gradually, thereby reducing embryo risk.

**NOTE:**

The thermal effect may be previewed by activating Isotherms under the Target Type. The Isotherms show the maximum temperature reached during pulse as function of position around the laser beam, for the low, medium, and high pulse durations.

Multiple Manual Pulses

Full zona breaching may not be accomplished after a single laser pulse. The thickness or hardness of the zona of most embryos may require the manual application of two or more laser pulses to the same area. When this is required, the conditions for the additional laser pulses should be as follows:

- The additional laser pulses are of the same or shorter duration as the first one and applied to a contiguous location to complete zona breaching (resulting in the formation of an oval shaped opening).
- The additional laser pulse is of a higher duration and applied to the same location to ablate more zona material.

Very small openings may lead to embryo constriction and abnormal development.

Summarized Laser Assisted Hatching Procedure**NOTE:**

The following procedure assumes that the system is properly set up and the target or RED-I aligned to the laser spot.

1. Position the laser objective module.
2. Adjust microscope and focus on the largest diameter of the zona.
3. Position specimen so a portion of the zona is in the path of the laser beam (target).
4. Select the laser intensity (Low, Med, High) on the Clinical Laser Energy panel.
5. Fire the laser using either the mouse or foot switch.
6. For single hole full zona breaching, if the hole in the zona must be enlarged, apply an additional laser pulse according to the recommendations under "Multiple Manual Pulses" above.
7. For zona thinning, reposition the embryo and fire again. Repeat until the desired portion of the zona is thinned.

Troubleshooting Failed Drilling Attempts

If the first attempt to drill the zona fails, it may be related to several factors. A number of probable causes and corrective actions are listed in the table below.

Table 12-1: Troubleshooting Failed Drilling Attempts

Effect	Cause	Corrective Action
Limited or no drilling effect	Embryo is raised away from bottom of the culture dish	Adjust micro-manipulators to lower embryo so it is sitting on the bottom of the culture dish
Limited drilling effect	Insufficient laser energy	• Select higher laser energy setting (longer pulse duration) • Fine tune focus
Laser focused on the culture medium	Focus not properly adjusted to largest diameter of embryo midsection	Fine tune focus
Laser does not hit target zona	Laser target not properly aligned to laser spot	Align laser target marker (see "Target Alignment" on page 35)

12.3 Laser-assisted Blastomere Biopsy

Caution: To remove the risk of damage to or excessive heating of blastomeres, administer laser pulses at the lowest power levels necessary to achieve prescribed zona drilling or thinning effects.

For laser-assisted blastomere biopsy the drilling of the zona pellucida is performed in the same way as with LAH. In this case the drilled hole should penetrate the entire thickness of the zona pellucida.

Blastocyst biopsy is usually performed on a Day 3 embryo. The embryo is supported on a holding pipette and held as close as possible to the floor of the Petri dish. The hole may be drilled adjacent to the cell to be extracted. Various users report that the diameter required for a biopsy drilled hole is between 10 μm and 25-30 μm . The required drilled diameter is best obtained by repeated adjacent pulses if necessary (most users report using between 2 and 4 pulses at the **High** setting [500-600 μs pulse duration]). It is less thermally disruptive to nearby cells if multiple short pulses are used rather than a single long pulse to drill out a required diameter.


Figure 12-2: Examples of Zona Openings for Blastomere Biopsy

12.4 Laser-assisted Trophectoderm Biopsy

Several protocols may be followed in preparation for trophectoderm biopsy.

1. Some users employ the laser to create a small opening in the zona pellucida on Day 3 of development. The embryo is then returned to the incubator to develop further. This allows the embryo to begin the hatching process early on Day 5 of development.
2. Another method is to laser drill the zona of the blastocyst early on Day 5, on the opposite side from the inner cell mass, and return the embryo to the incubator for several hours to promote hatching.
3. A third option is to use the laser to create an opening in the zona of an intact Day 5 embryo (on the opposite side of the inner cell mass) and using microtools, the trophectoderm cells are gently pulled out through the hole.

In all cases, the goal is for the trophectoderm cells to herniate through the opening in zona. Using microtools to hold the embryo and herniated trophectoderm cells, the trophectoderm cells are dissected using the laser. This may be done either with multiple manual pulses or using the **Multipulse** mode.

Using Manual Pulses to Remove Trophectoderm Cells

NOTE:

The following procedure assumes that the system is properly set up and the target or RED-i aligned to the laser spot.

NOTE:

The number of trophectoderm cells removed is based on your laboratory protocol and/or the number of cells requested by the testing facility.

1. Position the laser objective module.
2. The embryo and herniated trophectoderm cells are held on the bottom of the culture dish by the holding pipette and biopsy pipette, respectively.
3. Adjust microscope and focus.
4. Using the microtools, gently pull the trophectoderm cells into the biopsy pipette.
5. Position specimen so the cell junction between 2 trophectoderm cells is in the path of the laser beam (target).

6. Select the laser intensity (Low, Med, High) on the Clinical Laser Energy panel.
7. Fire the laser using either the mouse or foot switch while continuing to gently pull the trophectoderm cells into the biopsy pipette.
8. Repeat steps 5 through 7 until the trophectoderm cells separate from the embryo.

Using Multipulse Mode to Remove Trophectoderm Cells

Multipulse mode should only be used when performing trophectoderm biopsy procedures.

To prevent excessive heating of the inner cell mass when using Multipulse mode, laser application should not be closer than 60 μ m to the inner cell mass.

When using Multipulse mode, the laser can only be fired using the foot switch.

NOTE:

The following procedure assumes that the system is properly set up and the target or RED-i aligned to the laser spot.

1. Select the Multipulse mode and set the Max Mean Power to 1mW under Setup > Laser.
2. Press Save and Exit.
3. On the Multipulse Laser panel set the following parameters:
 - **Total Duration:** 1 second.
 - **Power:** 100%
 - **Pulse/sec:** 6
 - **Pulse Duration:** 500 μ s

NOTE:

The software will impose dynamic limits on the Power, Pulse/sec and Pulse Duration settings not to exceed a mean power of 1mW.

4. Position the laser objective module.
5. The embryo and herniated trophectoderm cells are held on the bottom of the culture dish by the holding pipette and biopsy pipette, respectively.
6. Adjust microscope and focus.
7. Using the microtools, gently pull the trophectoderm cells into the biopsy pipette.

NOTE:

The number of trophectoderm cells removed is based on your laboratory protocol and/or the number of cells requested by the testing facility.

8. Position specimen so the "waist" between the cells to be removed and the embryo is in the path of the laser beam (target).
9. Using the micromanipulator controls, practice moving the embryo across the target ring or RED-i prior to enabling the laser. This will help the user to confirm the right direction of manipulator in order to perform trophectoderm biopsy safely.
10. Press the Enable button on the Multipulse Laser panel to activate the laser and the timer (the trophectoderm biopsy must be started and completed within that time).

11. Press the foot switch to fire the pulsed laser and move the trophoctoderm (optimally at the cell junctions) across the pulsed laser beam while gently pulling the trophoctoderm cells into the biopsy pipette. The trophoctoderm cells will separate from the embryo.
12. Release the foot switch to stop the laser pulses.
13. Expel the biopsied cells from the pipette.

**NOTE:**

The user may perform another Multipulse only after a blanking duration equal to the total duration of the previous pulse train.

Additional Notes on Multipulse Mode

Each pulse has virtually the same effect as if it were a single pulse, since the cooling between pulses is so rapid and complete that the pulses are essentially independent. The result is that in the case of an embryo with particularly strongly attached cells, the herniated biopsy may still be separated, without damaging the biopsied cells or harming the rest of the embryo. In the Multipulse mode the average power of the multipulse laser beam is limited to 10 mW. The laser wavelength used is 1460 nm, which is infrared and is too low in energy to cause any actinic effects. Consequently, the only effect of the laser on the embryo's material is to heat it. Keeping the mean power of the laser to 10 mW maximum prevents the temperature at the embryo from rising significantly during the pulse train.

- In certain cases the separation of the biopsy specimen is harder to effect. For example, the presence of actin filaments connecting the biopsy specimen to the rest of the embryo can make the dissection of the biopsy much more difficult. Actin is more strongly bonded and repeated pulses may be required to weaken the actin bonds. The Max Mean Power setting of 1 mW can then be increased to a maximum of 10 mW on the Setup > Laser screen.
- The total pulse train duration can be increased to a maximum of 10 seconds if required.

Chapter 13. Help Menu

The **Help** menu provides access to information on the HT Laser 6 software and database maintenance utilities.

13.1 About

- The **About** dialog box includes Version of the software.
- Click the **License Agreement** button to review and agree to Hamilton Thorne's standard end user license agreement (EULA).

NOTE:

If you do not agree to the terms and the license agreement, you will not be able to use the program.

- Click the **Advanced** button to see more information on **Plugins** and **Options** of the software.

13.2 Database

Database provides access to the following **Backup**, **Restore** and **Optimize** database utilities.

Backup

The backup utility creates a copy of the database.

1. Select **Database > Backup** from the **Help** menu.
2. Select **Backup** from the **DB Maintenance** dialog box that appears.
3. Select a location, enter a file name for the backup file, and click **Save**. The file is saved with the extension **.backup**.
4. When the **Backup Complete** message appears, select **Exit**.

Restore

NOTE:

IMPORTANT: Any new data added to the HT Laser 6 database since the backup file was created will be lost if the database is restored from a backup file.

1. Select **Database > Restore** from the **Help** menu.
2. A notice stating that the program must be exited appears. Select **Yes** to continue with the restore process or **No** to return to the program.
3. Select **Restore** from the **DB Maintenance** dialog.
4. Locate the **.backup** file from which to restore and click **Open**.
5. When the restore process is complete, the HT Laser 6 program will automatically restart.

Optimize

Database optimization occurs automatically each time the HT Laser 6 program is initiated; however, if the system is left running for an extended period of time, select **Database > Optimize** to make the database most efficient.

13.3 Save Support Package

In the event of an error, a detailed set of data may be saved and sent to Hamilton Thorne to help determine the cause of the problem. In addition to mandatory storage of the **Event Log** and **Application Settings**, **Report Templates** and **Database** may also be saved. All data is compressed into a single zip file.

Figure 13-1: Save Support Package

Appendix A. Laser System Specifications

A.1 Operating Environment

Factor	Metric	USA
Temperature		
Operating	10°C - 37°C	50°F - 98°F
Non-operating	-10°C - 60°C	14°F - 140°F
Relative Humidity (non-condensing)		
Operating	10% - 90%	10% - 90%
Non-operating	5% - 95%	5% - 95%
Maximum Altitude (unpressurized)		
Operating	4,000 m	13,000 ft
Non-operating	9,000 m	30,000 ft

A.2 Rated Input Parameters

Inputs	Rating
Operating Voltage (power supply)	100 - 240 VAC RMS
Operating Frequency Range (power supply)	50/60 Hz AC
Operating Current (desktop system)	3.9 A RMS (max)
Operating Current (laptop system)	3.3 A RMS (max)
Input Voltage (laser controller)	12 V DC

A.3 Laser Specifications

Dimensions

Dimensions	H in. (mm)	W in. (mm)	D in. (mm)	lb. (kg)
LYKOS Objective	1.65 (42)	1.26 (32)	1.26 (32)	0.28 (0.13)
Laser Controller	3.8 (95)	2.5 (63.5)	1.8 (46)	1.0 (0.45)

Electrical

Electrical	Laser Controller	Desktop	Monitor	Laptop
Input Voltage	12v	100 - 240 VAC	110 - 240 VAC	110 - 240 VAC
Power	1 W	300 W	30 W	90 W
Line Frequency	50/60 Hz	50/60 Hz	50/60 Hz	50/60 Hz

Laser

Laser	Specification
Type	1460 ± 20nm, Infrared Solid State Diode
Maximum Power	@ Focus = 300 mW (Class 1)
Laser Modes	Clinical, Validation, Multipulse
Laser Modes	• Clinical: Low, Medium, High • Validation: Adjustable 0.1 to 3.0 ms in increments of 0.1 ms • Multipulse: - Normal default mean power: 1 mW - Normal pulse train duration: 1 second - Maximum mean power: 10 mW - Maximum Pulse train duration: 10 seconds
Firing	Mouse or foot switch
Target Marker	Circle, Arrow, Isotherm Rings. Isotherm Rings showing peak temperature.
Crosshairs	Activation, size and color set by user (used for positioning embryo)
Target Alignment	Performed interactively on screen

Laser Safety

Laser	Specification
Class 1 Laser	Laser pulse and energy cannot exceed Class 1 limits
Timer	Laser must be discharged within 15 seconds of activation while in Clinical mode
Energy	Once energy level is selected in Clinical mode, other energy levels are deactivated

RED-i Target Locator

RED-i	Specification
Source	Red LED built into laser module
Alignment	X-Y mechanical adjustment to center of laser beam using Isotherm Rings as a guide
Intensity	On screen intensity adjustment
Wavelength	665 nm and power 3.96 mW

Objective

Objective	Specification
Standard	40x, N.A. = 0.60, High transmission in visible, IR and UV, long working distance
Scale Calibration	Performed interactively on screen
Thread	LYKOS: RMS
Adapters	LYKOS: Thread adapters available for M25, M25 Nikon, M27
Working Distance	In Air: 2.92 mm To Petri (1mm and Heated plate (0.5mm)): 1.73 mm

A.4 Cleaning & Disinfection

The system may be cleaned and disinfected as often as needed. Routine maintenance will ensure the laser functions with optimal performance and extends the laser's long-term use.

Cleaning: The objective lens should be cleaned occasionally with a dry lens cloth only.

Caution: *The laser objective top lens should be cleaned only when needed and with extreme care not to damage, scratch, or remove the anti-reflective coating. At no time should any liquid be allowed to penetrate any part of the system.*

Disinfection: All exposed parts may be disinfected by cleaning with a mild, nonabrasive cleaner such as methanol or ethanol. Verify that the cleaning agent has evaporated from the outer surfaces and completely dried before powering the system on.

Hamilton Thorne encourages routine preventive maintenance service visits to help maximize the utilization of the product. We recommend having a certified Hamilton Thorne technician examine the laser's performance, professionally clean the system, and replace common wear and tear parts at least once a year. Annual preventive maintenance visits increase both system lifespan and reduce the risk of costly repairs. Contact your Hamilton Thorne representative to explore our multiple service plan options. For additional guidance on caring for your system, please refer to the customers only section on www.hamiltonthorne.com.

Index

A

- activating the laser, 50
 - 15 second timer, 50
- add graphic elements to an image, 45
- add image to gallery
 - with graphic overlay, 46
- ADMIN notation, 2
- Alpha, 29
- altitude
 - operating environment, 87
- applications, 77
 - laser-assisted blastomere biopsy, 80
 - laser-assisted hatching (LAH), 77
 - laser-assisted trophectoderm biopsy, 81
- arrow target/Target Type
 - arrow, 26
- arrow tool, 45
- ASCII Export, 74
- ASCII Filter, 74
- Audio Feedback, 28
- auto file storage
 - use report file name, 24, 31
 - use root name, 24, 31
- Auto Labels, 21
- auto labels
 - label properties, 22
 - name, 22
- Auto White Balance, 32
- Autosave on Capture, 24

B

- blastocyst biopsy, 80
- Body Properties, 70

C

- calibration, 37
- calibration procedure, 37
- camera
 - digital, 7, 13, 15, 16
 - orienting to stage, 14
- camera image
 - flip, 32
- camera port, 13

- Camera Settings, 32
- capture images, 44
- capture images
 - Capture on Laser Fire, 44
- capture images
 - with video overlay, 44
- Capture on Laser Fire, 44
- caution, 2, 11, 12
- circle target, 25
- cleaning, 89
- Clear Path, 60
- Clinical Laser panel, 50
- Clinical Mode
 - DTS, 58
- Clinical Pulse Settings, 28
- c-mount adapter, 5, 7, 13
- Com port, 34
- compensator cap, 5, 8
- components, 5
- computer connections
 - desktop, 15
 - laptop, 16
- control panel, 17
- controller, 5
- Copy report template, 65
- Copy setup, 33
- Crosshairs, 29
- Crosshairs Setup
 - color, 29
 - inner radius, 29
 - outer radius, 29
- Cycle Number, 39

D

- Data Fields, 42
- Data Items, 68
- Database, 85
 - Backup, 85
 - Optimize, 86
 - Restore, 85
- Date of Birth, 39
- DBPath, 72
- delete images, 43, 47
- Delete report template, 65
- Delete setup, 33

- deleting videos, 47
- delimiter, 74
- Design report template, 65
- desktop computer connections, 15
 - digital camera, 15
 - foot switch, 15
 - keyboard, 15
 - monitor, 15
 - mouse, 15
 - power, 16
 - security key, 15
 - speakers, 15
- digital camera, 7, 15, 16
 - desktop computer connection, 15
 - installation, 13
 - laptop computer connection, 16
- disinfection, 89
- drill tool type, 41
- dry erase marker, 8
- DTS

- Clinical Mode, 58
- Enabled, 59, 60
- Freehand, 56, 61
- Hole Size, 58
- Initialization, 56
- Initialize, 57
- Line/Curve, 56, 61
- Multipoint, 56, 61
- Path Type, 56
- Rectangle, 56, 61
- Single Shot, 56, 61
- Validation Mode, 60
- Verify, 57
- Verification, 56

- Dynamic Targeting System, 55

E

- Edit Label Location, 22
- electrical specifications
 - input voltage, 88
 - line frequency, 88
 - power, 88
- ellipse tool, 45
- Embryo ID, 39
- Embryo Info, 39
- Embryo panel, 18
- embryo tool type, 41
- Embryos in Cycle, 39
- enable button
 - Multipulse Mode, 53

- Enable Duration, 52
- Enable label, 22

F

- failed drilling, 80
- Fertilization Method, 39
- Fire Animation, 27
- Fire button, 50
- firing the laser, 50
 - 15 second timer, 50
 - Multipulse Mode, 53
- First Name, 39
- Flip Horizontal, 32
- Flip Vertical, 32
- focus, 77
- foot switch, 5, 8
 - desktop computer connection, 15
- Footer Properties, 70
- FPS, 32
- FPS Actual, 32
- frames per second, 32
- Freehand, 56, 61
- Full Screen Mode, 49
 - Image Toolbox, 46
- full zona breaching, 77

H

- Halt button
 - Multipulse Mode, 53
- Header Properties, 70
- Help, 85
- hex head drivers, 7
- high laser energy level, 50
 - laser-assisted hatching, 78
- high pulse duration, 28
- high pulse setting, 28
- Hole Size
 - DTS, 58

I

- Image, 68
- Image and Video Control Panel, 18
- Image Capture
 - Autosave on Capture, 24
 - Image File Name Label, 24
 - Image Type, 23
 - Save Location, 23, 31

Save Mode, 24, 31
 Watermark, 24
 Image Capture Settings, 23
 Image Capture Setup, 23
 auto, 24, 31
 manual, 24, 31
 Image File Name Label, 24
 Image Properties, 71
 Image Toolbox, 45
 arrow tool, 45
 ellipse tool, 45
 line tool, 45
 measurement tool, 45
 rectangle tool, 45
 text tool, 45
 undo button, 46
 Image Type, 23
 images
 adding to report, 44
 capturing, 44
 closing
 closing images, 44
 deleting, 43, 47
 Image Toolbox, 45
 opening, 43, 47
 save with overlays, 42
 saving, 43, 47
 viewing, 44
 Images panel, 18
 Info
 Embryo, 39
 Initialize DTS, 57
 Initialize DTS, 57
 inner radius, 29
 input voltage, 88
 installation
 preparing for, 9
 Isotherms
 activating rings, 27
 colors, 26
 peak temperature, 26
 Isotherms target, 26
 IVQC Report, 58

K

keyboard
 desktop computer, 15

L

label properties, 22
 LAH, 77
 laptop computer connections, 16
 digital camera, 16
 mouse, 16
 power, 16
 security key, 16
 Laser
 Enable, 25
 Laser - Clinical panel, 17
 laser applications, 77
 laser controller, 5, 6
 LED status bar, 13
 LYKOS installation, 11, 12
 serial cable connection, 12
 serial to USB connection, 12
 USB cable connection, 12
 Laser Fire button, 52
 Laser Mode Menu, 49
 laser objective, 6
 Laser Objectives, 34
 laser panel
 Clinical, 50
 Multipulse, 52
 Validation, 51
 laser safety, 88
 Laser Settings, 25
 Audio Feedback, 28
 Clinical Pulse Settings, 28
 Fire Animation, 27
 Laser, 25
 Pulse Resolution, 28
 Remote Fire, 25
 Scale, 30
 Target Type, 25
 Laser Setup
 max mean power, 52
 Multipulse Mode, 52
 target type, 25
 laser specifications, 88
 Laser Target Slide, 8
 RED-i target alignment, 35
 Verify DTS, 57
 Laser Target X and Y values, 34
 laser-assisted blastomere biopsy, 80
 laser-assisted hatching, 77
 full zona breaching, 77
 laser positioning, 77
 procedure summary, 79

- zona thinning, 78
- laser-assisted trophectoderm biopsy, 81
 - manual pulses, 81
 - multipulse mode, 82
- Last Name, 39
- line frequency, 88
- line tool, 45
- Line/Curve, 56, 61
- low laser energy level, 50
 - laser-assisted hatching, 78
- low pulse duration, 28
- low pulse setting, 28
- LYKOS installation
 - laser controller, 11, 12
 - RMS thread adapters, 11

M

- Magnification Calibration, 33, 37
 - X Micron/Pixel, 33
 - Y Micron/Pixel, 33
- Main Menu, 18
- Manual file storage
 - ask for file name, 24, 31
- measurement tool, 45
- medium laser energy level, 50
 - laser-assisted hatching, 78
- medium pulse duration, 28
- medium pulse setting, 28
- micrometer, 5, 8
- microscope
 - camera port, 13
- monitor
 - desktop computer connection, 15
- mouse
 - desktop computer, 15
 - laptop computer connection, 16
- multiple manual pulses, 79
- Multipoint, 56, 61
- Multipulse Laser panel, 52
 - enable, 52
 - halt button, 52
 - power, 52
 - pulse duration, 52
 - pulses/sec, 52
 - total duration, 52
- Multipulse Mode, 52
 - firing the laser, 53
 - notes, 83

N

- New Folder, 65
- New setup, 33
- New Template, 65
 - designing new report, 66
- Note, 39
 - note, 2
- Number, 41

O

- Objective Info, 33
 - Description, 33
 - Name, 33
 - Type, 33
- objective specifications, 89
- Objectives Settings, 33
- Objectives Setup
 - Com port, 34
- open images, 43, 47
- opening videos, 47
- operating current, 87
- operating environment, 87
 - altitude, 87
 - relative humidity, 87
 - temperature, 87
- operating frequency, 87
- operating voltage, 87
- Optics Settings, 32
 - Camera, 32
 - Objective, 33
- outer radius, 29
- Overlay Settings, 29
 - Crosshairs, 29

P

- Page Number, 68
- Path Type, 56
- Patient ID, 39
- Patient Info, 39
- Patient Info panel, 18
- Persist, 42
- Physician, 39
- Playback Speed, 30
- playing videos, 48
- PN 1
 - PN 2, 41
- power, 16

- electrical specifications, 88
- Multipulse Mode, 53
- Validation Mode, 51
- Power Limit, 52, 54
- Programmable Fields, 22
- pronucleus tool type, 41
- Properties Panel, 69
- pulse duration
 - Multipulse Mode, 53
 - Validation Mode, 51
- Pulse Resolution, 28
- Pulse Settings
 - Clinical Pulse Settings, 28
 - Pulse Resolution, 28
- pulses/sec
 - Multipulse Mode, 53

Q

- Quick Properties Tool, 68

R

- real time video capture, 47
- Rectangle, 56, 61
- rectangle tool, 45
- RED-i
 - brightness, 17
 - enable, 17
 - Intensity Slider, 17, 35
- RED-i specifications, 89
- RED-i target alignment
 - Laser Target slide, 35
- relative humidity
 - operating environment, 87
- Remote Fire, 25
- Report Date & Time, 68
- report images, 44
- Report Layout, 69
- Report Properties, 69
- Report Template
 - Body, 70
 - Copy, 65
 - Delete, 65
 - Design, 65
 - design area, 67
 - Footer, 70
 - Header, 70
 - Image Properties, 71
 - New Folder, 65

- New Template, 65
- Properties, 68
- Properties Panel, 69
- Text Properties, 72
- Toolbox, 68
- View, 65
- Report template
 - adding items, 67
- Report Templates
 - layout aids, 68
- Resize, 32
- Retrieval Date, 39
- RMS thread adapters, 5, 6, 11

S

- safety delay, 27
- save images, 43, 47
- Save Location, 23, 31
- Save Mode, 24, 31
- Save Support Package, 86
- saving videos, 47
- Scale, 30
- security key, 5, 7, 15
 - desktop computer connection, 15
 - laptop computer connection, 16
- setting the laser power
 - Validation Mode, 51
- setting the pulse
 - Validation Mode, 51
- Settings
 - Auto Labels, 21
 - Image Capture, 23
 - Laser, 25
 - Optics, 32
 - Overlay, 29
 - Time Lapse, 30
- Setup
 - Delete, 33
- Setups
 - Copy, 33
 - New, 33
- Single Shot, 56, 61
- software layout, 17
- speakers, 15
- special notations
 - ADMIN, 2
 - caution, 2, 11, 12
 - note, 2
- specifications, 87
- stage micrometer, 5, 8

Standard Items, 68
system
 magnification calibration, 37

T

target alignment, 35, 77
Target Type, 25
 circle, 25
 Isotherms, 26
target type, 25
temperature
 operating environment, 87
Text Label, 68
Text Properties, 72
text tool, 45
thread adapters
 RMS, 5, 6, 11
Thumbnail Gallery, 18, 43
Time Lapse
 Playback Speed, 30
Time Lapse Duration, 30
Time Lapse Interval, 30
Time Lapse Settings, 30
Time Lapse videos, 47
tool type
 drill, 41
 embryo, 41
 pronucleus, 41
 zona, 41
total duration
 Multipulse Mode, 53
trophectoderm biopsy, 81
 manual pulses, 81
 multipulse mode, 82
troubleshooting
 failed drilling, 80

U

undo button, 46
USB security key, 5, 7, 15

V

Validation Laser panel, 51
Validation Mode, 51
 DTS, 60
Verify DTS, 57
video window, 17

videos
 capturing real time videos, 47
 capturing time lapse videos, 47
 deleting, 47
 opening, 47
 saving, 47
video
 playing, 48
View report template, 65
viewing images, 44

W

Watermark, 24
White Balance, 32
With Video Overlays, 44

X

X Magnification Calibration, 37
X Micron/Pixel, 33

Y

Y Magnification Calibration screen, 37
Y Micron/Pixel, 33

Z

zona thinning, 78
zona tool type, 41

Application No.1 to the User Manual of LYKOS Clinical Laser with the Dynamic Targeting System option (DTS) for In Vitro Fertilization (IVF) Laboratories in the following assembly

Information for users in the Russian Federation

Name and supply set of the medical device

LYKOS Clinical Laser with the Dynamic Targeting System option (DTS) for In Vitro Fertilization (IVF) Laboratories in the following assembly:

- 1. LYKOS laser objective.**
- 2. Laser controller in the following assembly:**
 - a. Laser controller unit.**
 - b. Interface cable.**
 - c. Source of supply for laser controller unit.**
 - d. Power cord supply source.**
 - e. USB cable.**
- 3. Threaded adapters equipped with a setting tool and spacers, 3 pcs. (if any).**
- 4. Screwdriver Adjusting Tool Set:**
 - a. Long hex driver 1.5 mm to adjust RED-i.**
 - b. Short hex driver 1.5 mm to adjust RED-i.**
- 5. USB access key.**
- 6. C-mount objective adapter (if any).**
- 7. Digital colored camera (if any).**
- 8. Cable to connect the digital camera (if any).**
- 9. Foot switch (pedal) with connecting cable.**
- 10. Micro-graduated glass.**
- 11. Laser targeting adjusting slide.**
- 12. DVD with HT Laser software.**
- 13. User manual.**
- 14. Wired mouse.**
- 15. Wireless mouse (if any).**
- 16. USB concentrator (if any).**
- 17. DVD-driving unit (if necessary).**
- 18. Plastic case for device components.**
- 19. Laptop with HT LASER predetermined software (if any).**
- 20. Desktop system unit (if any).**

Intended purpose of the medical device set by the manufacturer

The device is used as a portable laser assembly energized with the alternating current where the incoming energy excites the diode to form the low level laser beam in order to ablate the zona pellucida at IVF-based laser-assisted hatching or embryo biopsy. It is intended for use with a compatible inverted light microscope. It is in assembly with a combined microscope objective, laser controller unit, camera and additional devices for auxiliary use (adjusting devices, software, foot switch).

Scope of application

The device is intended for use in the area of assisted reproductive technologies (ART) at IVF.

Potential consumers

LYKOS Clinical Laser is intended for use in medical and health facilities that are certified to provide help in the field of assisted reproductive technologies (ART), such as clinics, laboratories, and IVF centers.

The device is to be handled by qualified embryologists with the necessary specialization and relevant professional knowledge and skills, as well as experience in the field of assisted reproductive technologies (ART).

Special training on how to operate LYKOS Clinical Laser is recommended as the method has many technical aspects essential for optimal IVF results.

Information on availability of a medicinal product for human use, any materials of animal and (or) human origin in the medical device

LYKOS Clinical Laser with the Dynamic Targeting System option (DTS) for In Vitro Fertilization (IVF) Laboratories does not contain any medicinal products or any materials of animal and (or) human origin.

Information on a sterile condition of the medical device, method of sterilization thereof and on order of actions in case if the sterile package is damaged (if the medical device is supplied in a sterile condition) or, if the medical device is not supplied in a sterile condition, indication on the necessity to sterilize it before use

LYKOS Clinical Laser with the Dynamic Targeting System option (DTS) for In Vitro Fertilization (IVF) Laboratories and its components/accessories are not sterile.

Disposal

Grade in accordance with Sanitary Regulations and Norms 2.1.7.2790-10 "Sanitary and epidemiological requirements for handling medical waste — Grade A".

It is prohibited to utilize the device with household waste. The user must contact the authorized representative or distributor of the company.

Health professionals must follow national/ministerial guidance to treat medical waste and general safety precautions. Protective clothing should be worn. It is prohibited to dispose of the device as household waste. An authorized representative or a distributor should be contacted.

Health professionals must follow national/ministerial guidance to treat medical waste and general safety precautions. Protective clothing should be worn.

Lifetime of the device

The lifetime mentioned by the manufacturer is 7 years. The issue of prolonging the lifetime can be decided based on the results of respective tests and technical maintenance after the nominal lifecycle is over.

The guaranteed service life is presented in the operational manual and specification. The expected service life is taken into account when analyzing the associated risks in accordance with EN ISO 14971:2012.

Storage conditions

Temperature: -10 to +70 °C

Relative humidity: 10 to 90 %

Pressure: 710 to 1,060 hPa

Transport conditions

Temperature: -10 to +70 °C

Relative humidity: 10 to 90 %

Pressure: 710 to 1,060 hPa

Service and support

The device must be installed, maintained and repaired by the service provider authorized by the manufacturer.

Hamilton Thorne encourages regular preventive maintenance by its own personnel. Then the product will be used to the maximum extent. We recommend that a certified technician of Hamilton Thorne check the laser working capacity, clean the system professionally and replace the details subject to deterioration at least once a year. Annual preventive visits extend the service life and decrease the risk of costly repair. Contact your representative of Hamilton Thorne to examine various types of maintenance plan. To obtain further maintenance instructions see Customer Only section at www.hamiltonthorne.com.

Warranty liability

The devices are covered by a standard 12-month warranty (warrantied shelf life and operation period). The warranty covers the devices that are stored and used in accordance with the user manual.

Authorized representative of the manufacturer of the medical device

INTERGEN OOO (ООО «ИНТЕРГЕН»),
legal and postal address: 4A, ul. Nagornaya, 117186 Moscow (117186, Москва,
ул. Нагорная, д. 4А)
INN 7727507861/KPP 772701001, OGRN 1047796214828,
tel.: (495) 724-38-62, tel./fax: (499) 969-83-38

Information on original release or the latest review of the operational documents
This application was drawn up in September, 2020.

Change history of the operational documents:		
Date	Change	Content/Comment:

OFFICIAL DOCUMENT
Keith Edwards
KEITH EDWARDS
HAMILTON THORNE INC

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New Valley Fertility Clinic
123 Elm St. Beverly, MA 01915
978-555-1212

Embryology Report
Tuesday, December 26, 2017 2:53 PM

Patient Information

ID: patient_id Patient Birthdate: patient_birthdate
First Name: patient_first_name Last Name: patient_last_name
Patient Physician: patient_physician
Embryologist: data_first_data_01

Embryo Information

ID: embryo_id Retrieval date: embryo_date_of_retrieval
Cycle Number: embryo_cycle_number Embryos in Cycle: embryo_number_of_embryos_in_cycle
Fertilization Method: embryo_fertilization_method Note: embryo_note

Embryo Measurements

Date: embryo_measurements_date
Number Of Blastomeres: embryo_measure

Embryo	Zone	Dish	PN 1	PN2
1	embryo_measure	embryo_measure	embryo_measure	embryo_measure
2	embryo_measure	embryo_measure	embryo_measure	embryo_measure
3	embryo_measure	embryo_measure	embryo_measure	embryo_measure
4	embryo_measure	embryo_measure	embryo_measure	embryo_measure
5	embryo_measure	embryo_measure	embryo_measure	embryo_measure
Mean	embryo_measure	embryo_measure	embryo_measure	embryo_measure
SD	embryo_measure	embryo_measure	embryo_measure	embryo_measure

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Patient Information

ID: 963330K3 Patient Birthdate: 11/3/1981
First Name: Leslie Last Name: Smythe
Patient Physician: Dr. Jack Smothers
Embryologist: J.B. Hubbs

Embryo Information

ID: 4-11172017 Retrieval date: 11/13/2017
Cycle Number: 1 Embryos in Cycle: 6
Fertilization Method: ICSI Note: None

Embryo Measurements

Date: Friday, November 17, 2017 1:20 PM
Number Of Blastomeres: 0

Embryo	Zone	Dish	PN 1	PN2
1	104.3	12.3	9.7	11.7
2	100.5	12.6	9.8	12.5
3	101.6	11.8	9.7	
4	106.4	12.3	10.5	
5	99.2	13.1	10.6	
Mean	102.4	12.5	10.1	12.1
SD	2.9	0.5	0.5	

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© Copyright 2019. Hamilton Thorne, Inc. Все права защищены. LYKOS® и RED-i® являются зарегистрированными товарными знаками компании Hamilton Thorne, Inc. Все другие торговые марки являются товарными знаками или зарегистрированными товарными знаками соответствующих компаний.