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Approved by
Technical Article Manager
Ken Olsson

Essity Hygiene and Health AB

Ken Olsson

Essity Hygiene and Health AB
405 03 GÖTEBORG
SWEDEN

Instruction for use

Urological pads TENA®
TENA Lady Slim Mini, TENA Lady Slim Normal, TENA Lady Normal,
TENA Lady Normal Night, TENA Lady Super, TENA Lady Maxi Night

Essity Hygiene and Health AB, SE-405 03 Goteborg, Sweden

I, the undersigned, Adrienne Bonde, Deputy Notary Public
in the City of Stockholm, Sweden, certify that

Ken Birger Olsson---

has issued and signed the foregoing document.

Fee
Crowns

Stockholm, 2020-12-01
Ex officio:

Adrienne Bonde

2020

1. Name of MD.

Urological pads TENA®:

TENA Lady Slim Mini, TENA Lady Slim Normal, TENA Lady Normal, TENA Lady Normal Night, TENA Lady Super, TENA Lady Maxi Night

2.Manufacturer: Essity Hygiene and Health AB, SE-405 03 Goteborg, Sweden
Production site: Essity Operations Hoogezand B.V., Abramskade 6, 9601 KM, Hoogezand, The Netherlands.

3.Intended use of MD

Urological pads TENA® protect against drop, light and moderate urine leakage.

They are developed so that women can stay confident and secure every day, in all situations. The TENA Lady products provide Triple Protection against leaks, odour and moisture, helping women stay dry, secure and confident every day.

TENA Lady Slim Mini protect against drop and light incontinence.

TENA Lady Slim Normal and TENA Lady Normal protect against light incontinence.

TENA Lady Super и TENA Lady Maxi Night protect against light and moderate incontinence. TENA Lady Normal Night purpose made products for night, designed for light and moderate incontinence; The pads are unsterile, 1st risk class.

4.INDICATIONS

For women suffering from drop, light and moderate incontinence.

5.CONTRAINDICATIONS

No contraindications found.

6.Side effects

No side effects found

7. Precautions

N/A

8. Information about potential customers

This medical device is intended for women living an active lifestyle, suffering from drop, light and moderate incontinence.

9. Functional characteristics and weight of MD

Model	Weight (g)± 15%	Average value of absorption capacity (g) Rothwell
TENA Lady Slim Mini	7,14	169

TENA Lady Slim Normal	13,64	349
TENA Lady Normal	15,43	350
TENA Lady Super	38	860
TENA Lady Normal Night	22,4	530
TENA Lady Maxi Night	35,4	914

Rewet under pressure:

A certain amount of liquid (0,9% NaCl) is added to the product and after a certain resting time the product is subjected to a pressure and rewet is measured

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Lady Slim Normal V= 50 ml

Lady Normal V= 50 ml

Lady Normal Night V= 70 ml

Lady Super V= 150 ml

Lady Maxi Night V= 150 ml

Acquisition:

Liquid added in two increments with a resting time of 10 minutes between.

Lady Slim Mini: First dose V= 10 ml, second dose V= 10 ml

Lady Slim Normal: First dose V= 20 ml, second dose V= 20 ml

Lady Normal: First dose V= 20 ml, second dose V= 20 ml

Lady Normal Night: First dose V=50 ml, second dose V= 50 ml

Lady Super: First dose V= 80 ml, second dose V= 80 ml

Lady Maxi Night: First dose V= 80 ml, second dose V= 80 ml

Centrifuge absorption capacity:

This parameter is very dependent on the type of the centrifuge used for the test and is the reason for different results in different laboratories. The important parameter to be known is the g-force and the time the product is subjected to during the centrifugation. Therefore, similar values between laboratories are difficult to reach.

<i>Name of product</i>	<i>Rewet</i>	<i>Inlet time</i>	<i>Sorption capacity after centrifugation</i>	<i>Distribution of the SAP in the structure of the absorption layer (digital centrifugation)</i>	<i>pH value of the product extract</i>
TENA® Lady Slim Mini	One dose of 20 ml; average 2,7 g; stdev 0,8	One dose of 10 ml: first inlet average 1,6 s stdev: 0,2 Second dose of 10 ml: average 1,5 s stdev: 0,2	N/A	moisture-retention property of symmetrical sections of the samples on the left and right is the same	6 stdv 0,05
TENA® Lady Slim Normal	One dose of 50 ml; average 7 g; stdev 1,7	One dose of 50 ml: first inlet average 2,3 s stdev: 0,2 Second dose of 50 ml: average 2,8 s stdev: 0,2	N/A	moisture-retention property of symmetrical sections of the samples on the left and right is the same	6 stdv 0,05

TENA® Lady Normal	One dose of 50 ml; average 13 g; stdev 1,5	One dose of 50 ml: first inlet average 1,9 s stdev: 0,13 Second dose of 50 ml: average 3,2 s stdev: 0,4	N/A	moisture-retention property of symmetrical sections of the samples on the left and right is the same	6 stdv 0,05
TENA® Lady Normal Night	One dose of 70 ml; average 19 g; stdev 3,1	One dose of 50 ml: first inlet average 5,2 s stdev: 0,5 Second dose of 50 ml: average 8,5 s stdev: 1,4	N/A	moisture-retention property of symmetrical sections of the samples on the left and right is the same	6 stdv 0,05
TENA® Lady Super	One dose of 150 ml; average 48 g; stdev 4,7	One dose of 80 ml: first inlet average 7,9 s stdev: 0,5 Second dose of 80 ml: average 12,6 s stdev: 0,7	N/A	moisture-retention property of symmetrical sections of the samples on the left and right is the same	6 stdv 0,05
TENA® Lady Maxi Night	One dose of 150 ml; average 33 g; stdev 3,2	One dose of 80 ml: first inlet average 9,4 s stdev: 1,0 Second dose of 80 ml: average 15,7 s stdev: 1,4	N/A	moisture-retention property of symmetrical sections of the samples on the left and right is the same	6 stdv 0,05

10. Risks of using a medical product, contraindications, expected and predictable side effects associated with the use of a medical product for the intended purpose.

None

11. How to apply MD

Open individual pack (if available) and take out a pad; Remove release paper
Fix a pad by adhesive tape to the underwear, Change pad when appropriate

12. Handling conditions

To handle product under normal conditions


Requirements for environmental conditions during handling	
Air temperature, °C	+15 C- +35C
Humidity (relative), without condensation	not more 40%

13. Technical characteristics of MD TENA Lady Slim Mini


№	Measurement point	mm± 15%
1	Length of pad (A)	219
2	Width of pad (B)	96
3	Thickness of pad (C)	3,95
4	Length of transfer layer (D)	219
5	Width of transfer layer (E)	50
6	Length of absorption core (F)	196
7	Width of absorption core (G)	76
8	Length of back sheet (H)	219
9	Width of back sheet (I)	96


Scheme of MD

TENA Lady Slim Normal


№	Measurement point	mm ± 15%
1	Length of pad (A)	277
2	Width of pad (B)	106
3	Thickness of pad (C)	7,3
4	Length of transfer layer (D)	160
5	Width of transfer layer (E)	50
6	Length of absorption core (F)	248
7	Width of absorption core (G)	85
8	Length of back sheet (H)	277
9	Width of back sheet (I)	106


TENA Lady Normal

№	Measurement point	mm ± 15%
1	Length of pad (A)	277
2	Width of pad (B)	106
3	Thickness of pad (C)	9,5
4	Length of transfer layer (D)	170
5	Width of transfer layer (E)	50
6	Length of absorption core (F)	250
7	Width of absorption core (G)	85
8	Length of back sheet (H)	277
9	Width of back sheet (I)	106

TENA Lady Normal Night


№	Measurement point	mm ± 15%
1	Length of pad (A)	347
2	Width of pad (B)	121
3	Thickness of pad (C)	10
4	Length of transfer layer (D)	170
5	Width of transfer layer (E)	50
6	Length of absorption core (F)	315
7	Width of absorption core (G)	92
8	Length of back sheet (H)	347
9	Width of back sheet (I)	121

TENA Lady Super


№	Measurement point	mm ± 15%
1	Length of pad (A)	430
2	Width of pad (B)	196
3	Thickness of pad (C)	13
4	Length of transfer layer (D)	170
5	Width of transfer layer (E)	70
6	Length of absorption core (F)	394
7	Width of absorption core (G)	edges 161/middle 85
8	Length of back sheet (H)	430
9	Width of back sheet (I)	196


TENA Lady Maxi Night

№	Measurement point	mm $\pm$ 15%
1	Length of pad (A)	394
2	Width of pad (B)	182 back 146 front
3	Thickness of pad (C)	13,5
4	Length of transfer layer (D)	180
5	Width of transfer layer (E)	50
6	Length of absorption core (F)	360
7	Width of absorption core (G)	155 back 119 front
8	Length of back sheet (H)	394
9	Width of back sheet (I)	182

12. How to apply MD.

Open individual pack (if available) and take out a pad; Remove release paper
Fix a pad by adhesive tape to the underwear, Change pad when appropriate.

13. Materials used

Components	Function	Material	Brand	Supplier	MSDS (date or №)
Absorbent core	Absorb and store urine	Fluff pulp	FR Supersoft Plus 8%	International Paper;	Safety data Sheet 25/09/2017
		Super absorbent polymer	SXM 9465	Evonik Industries	Safety data sheet 29/09/2016
Topsheet	User comfort/dry skin	Non woven; printed, silk	NW-SB, philic 111018HID 2593/U	Berry France	Safety data sheet 25/04/2018
Transfer layer	Transport urine	Acquisition layer (high loft)	AB50E334TPS	Texus	Safety data sheet MSAB0059 10/09/2009
Glue	Joining	Blend of polymers (Hot melt)	Kmelt PERFECT CARE W510	Colquimica S.A.	Safety data sheet 13/01/2016
Backsheet	Prevent leakage	Polyethylene film (Mini & Super)	10203 Hyfol PP Film Extra BTXB, 25 gsm FPL M25L	RKW	Safety data sheet 1/1/2014
		Polyethylene/Non-Woven		Aschersleben GmbH	Safety data sheet Oct 2018
Fixation glue	Fixates the product in the underwear	Pressure Sensitive Hot Melt	NW 1208	HB Fuller	Safety data sheet 09/10/2017
Release paper	Covers fixation glue	Kraft paper coated on one side with silicone	37 gsm	Mondi Herlen	Safety data sheet 01/01/2012

14. Labelling of consumer pack and MD

Back sheet of TENA Lady Slim Mini, TENA Lady Slim Normal, TENA Lady Normal, TENA Lady Super, TENA Lady Normal Night, TENA Lady Maxi Night contains LOT_Country_YYMMDDHH:MM.

Consumer pack contains:

- name of product, a trademark;

- country of manufacture;
- name of the legal manufacturer, address.
- number of pads in the package;
- date (YYMMDD HH:MM) of production
- shelf life; (BBD) YYMMDD
- storage conditions
- disposal instructions: "Do not throw into drains";
- absorbency of pads
- bar code;
- Lot number
- article;
- number of registration certificate and date of issue

Graphic symbols on consumer pack of MD

Symbol	Explanation
	The symbol of recycling of pulp products
	The possibility of processing low density polyethylene (LDPE).
	Conformity symbol when declaring conformity
	Throw in the ballot box
	Prohibition of re-use
	Do not dispose into drains
	Absorbency index
	CE Marking

	Standard icons for EU
	How to apply pads
	Date of production

15. Information about transportation of MD

Pads are transported only in a dry, clean transport that does not have extraneous odors
Temperature -30°C to +35°C; Relative humidity of air no more than 80%

16. Storage and Shelf life

Environmental conditions during storage

Requirements for environmental conditions during storage	
Air temperature, °C	+ 15°C - +35°C
Air relative humidity	not more 40%

Shelf life of the pads is 3 years from the production date indicated on the package.

17. UTILIZATION

The product must be disposed of as medical waste of hazard class "B"
Do not throw into drains.

18. Manufacturer's warranty


The manufacturer guarantees the quality of pads under the conditions of operation, storage and transportation within 3 years from the date of manufacture date indicated on the package.

19. Authorized representative in Russia

During the warranty period, claims are made to the manufacturer at the address: LLC Essity Russia, 117218 Moscow, ul.Krzhizhanovskogo, 14, building 3, Tel: +7 (495) 967-33-67.

20. List of applicable national standards

Harmonized and other International Standards that are relevant to prove conformity to Medical Devices Regulation (EU) 2017/745 by Essity Hygiene and Health AB applied in full or in part.

ISO 13485:2016	Medical devices - Quality management systems – Systems requirements for regulatory purposes This international standard specifies requirements for a quality management system that can be used by an organization for the design and development, production, installation and servicing of medical devices, and the design, development, and provision of related services.
ISO 9001:2015	Quality management systems This international standard specifies requirements for a quality management system where an organization needs to demonstrate its ability to consistently provide product that meets customer and applicable regulatory requirements, and aims to enhance customer satisfaction through the effective application of the system, including processes for continual improvement of the system and the assurance of conformity to customer and applicable regulatory requirements.
ISO 14971:2012	Risk management This international standard is a framework for effective management by the manufacturer of the risks associated with the use of medical devices. It includes a procedure to investigate, using available information, the safety of a medical device, including accessories, by identifying hazards and estimating the risks associated with the device.
EN ISO 10993:2009	Biological evaluation a) the fundamental principles governing the biological evaluation of medical devices the definition of categories of devices based on the nature and duration of contact with the body; the selection of appropriate tests.
EN ISO 17025:2005	General Requirements for the competence of testing and calibrating laboratories This international standard specifies the general requirements for the competence to carry out tests and/or calibrations, including sampling. It covers testing and calibration performed using standard methods, non-standard methods and laboratory-developed methods.
EN 1041:2008	Information supplied by the manufacturer of medical devices This standard specifies the information to be supplied by a manufacturer for different categories of medical devices, as required by the relevant EEC Directives. It does not specify the language to be used for such information.
ISO 15223:2016	Medical devices – Symbols to be used with medical device labels, labeling and information supplied. This international standard identifies symbols conventionally used to convey information essential for proper use to the user and others for safe and effective use of medical devices.
ISO 7000:2012	Graphical symbols for use on equipment – Index and synopsis. This international standard provides a synopsis of graphical symbols, which are placed on equipment or parts of equipment of any kind in order to instruct the persons handling the equipment as to its use and operation.

EN ISO 14155:2011

Clinical investigation of medical devices for human subjects

This international standard defines procedures for the conduct and performance of clinical investigations of medical devices. It specifies general requirements intended to --protect human subjects, --ensure the scientific conduct of the clinical investigation, --assist sponsors, monitors, investigators, ethics committees, regulatory authorities and bodies involved in the conformity assessment of medical devices.

Part 1 of ISO 14155 a) specifies requirements for the conduct of a clinical investigation such that it establishes the performance of the medical device during the clinical investigation intended to mimic normal clinical use, reveals adverse events under normal conditions of use, and permits assessment of the acceptable risks having regard to the intended performance of the medical device, b) specifies requirements for the organization, conduct, monitoring, data collection and documentation of the clinical investigation of a medical device, c) is applicable to all clinical investigation(s) of medical devices whose clinical performance and safety is being assessed in human subjects.

Part 2 of this standard provides requirements for the preparation of a Clinical Investigation Plan (CIP) for the clinical investigation of medical devices. The compilation of a CIP in accordance with the requirements of this standard and adherence to it will help in optimising the scientific validity and reproducibility of the results of a clinical investigation.

MEDDEV 2.4/1 rev. 9	Guidelines to the Classification of Medical Devices
MEDDEV 2.7.1 rev. 4	Guidelines on Medical Devices, Clinical Evaluation
MEDDEV 2.12/1 rev.8	Guidelines on a Medical Device Vigilance System
Commission Decision 2000/532/EC	European List of Waste

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