Approved by

Frank Ruseler

Date: August 11, 2020

(Signature)

Name of Responsible Person

Stamp

Instruction for use:

«Disposable Ophthalmic Surgical Forceps»

Manufactured by:

Vitrex B.V., The Netherlands

Vitrex B.V.
Seggelant-Noord 2
3237 MG Vierpolders
The Netherlands
T: +31 181 745 090
www.vitrex.com

INSTRUCTION FOR USE

Disposable Ophthalmic Surgical Forceps IFU for Russian Submission

2020

Approved by

Frank Ruseler

Date: _____

(Signature)

Name of Responsible Person

Stamp

Instruction for use:

«Disposable Ophthalmic Surgical Forceps»

Manufactured by:

Vitrex B.V., The Netherlands

INSTRUCTION FOR USE

Disposable Ophthalmic Surgical Forceps IFU for Russian Submission

2020

Name of the medical device

Disposable Ophthalmic Surgical Forceps, versions (hereinafter may be referred to as the forceps, the medical device):

1. Disposable Ophthalmic Surgical 20G Forceps, Serrated Gripping
2. Disposable Ophthalmic Surgical 25G Forceps, Serrated Gripping
3. Disposable Ophthalmic Surgical 23G Forceps, Serrated Gripping
4. Disposable Ophthalmic Surgical 20G Forceps, ILM Endgripping
5. Disposable Ophthalmic Surgical 23G Forceps, ILM Endgripping
6. Disposable Ophthalmic Surgical 25G Forceps, ILM Endgripping
7. Disposable Ophthalmic Surgical 27G Forceps, ILM Endgripping
8. Disposable Ophthalmic Surgical 20G Forceps, Asymmetrical Endgripping
9. Disposable Ophthalmic Surgical 23G Forceps, Asymmetrical Endgripping
10. Disposable Ophthalmic Surgical 25G Forceps, Asymmetrical Endgripping
11. Disposable Ophthalmic Surgical 20G Forceps, Eckardt Endgripping
12. Disposable Ophthalmic Surgical 23G Forceps, Eckardt Endgripping
13. Disposable Ophthalmic Surgical 25G Forceps, Eckardt Endgripping
14. Disposable Ophthalmic Surgical 27G Forceps, Eckardt Endgripping
15. Disposable Ophthalmic Surgical 23G Forceps, Kuhn 45° Power Gripping
16. Disposable Ophthalmic Surgical 25G Forceps, Kuhn 45° Power Gripping
17. Disposable Ophthalmic Surgical 23G Forceps, Shah Xtra Grip
18. Disposable Ophthalmic Surgical 25G Forceps, Shah Xtra Grip
19. Disposable Ophthalmic Surgical 27G Forceps, Shah Xtra Grip

Legal manufacturer: Vitreq B.V. / Barpex E.B., Seggelant-Noord 2, 3237 MG Vierpolders, The Netherlands.

Site of manufacture: Vitreq B.V. / Barpex E.B., Seggelant-Noord 2, 3237 MG Vierpolders, The Netherlands.

Intended use of the medical device

A hand-held manual instrument designed to grasp and manipulate intraocular tissue during ophthalmic surgery.

Indications:

Damaged, detached, dead or foreign tissue in the retina area, which need to be repaired, removed or repositioned

Contraindications:

- Resterilization and/or reuse of the device by the buyer. Resterilization and/or reuse of the device may result in cross-contamination and/or significant deterioration of the product performance;
- Use of the device if the packaging integrity has been compromised;
- Use of the expired device.

Information on potential users of the medical device

These instruments are intended for use by an ophthalmic surgeon during ophthalmic surgery. The devices are intended for professional use only.

Description of the medical device

Forceps consist of a body, a rod and a tip. The body is the same for all instruments. When closing the body shanks, the rod moves forward and covers the tip. The instrument tip determines its function. Forceps are usually used to grip the mucous membranes/internal limiting membrane/retina/blood clots/foreign bodies in the posterior segment of the eye. Instruments are available in a range of gauges: 20, 23, 25, 27 G.

3

2

1


Figure 1. Picture of disposable forceps: 1 – body, 2 – rod, 3 – tip

The tip design is similar to that used in commercially available products, apart from the Eckardt Endgripping tip for smaller models. The 20G and 23G models are equivalent to similar devices from other manufacturers, but for the 25G and 27G gauges, the tip width at the distal end corresponds to the 23G tip width, so this tip has a larger gripping area.


Tip styles:

Forceps, Shah Xtra Grip	
Forceps, Kuhn 45° Power Gripping	
Forceps, Eckardt Endgripping	

Forceps, Asymmetrical Endgripping	
Forceps, ILM Endgripping	
Forceps, Serrated Gripping	

Schematic images of the medical device:

**The tolerance is $\pm 5\%$, unless otherwise specified.*

A	Handle length
B	Length of the working part
C1	Width in the activated position
C2	Width in open position
D	Diameter of the working part
H	Activation force (Newtons)
I	Weight [grams]

Table. Overall dimensions of the medical device.

	A	B	C1	C2	D	H	I
20G	110 mm	36 mm	8.9 mm	14.6 mm	0.9 mm	≤ 5 N	4.6 g
23G	109.8 mm	35 mm	8.9 mm	14.6 mm	0.6 mm	≤ 5 N	4.4 g
25G	109.7 mm	33 mm	8.9 mm	14.6 mm	0.5 mm	≤ 5 N	4.8 g
27G	109.9 mm	35 mm	8.9 mm	14.6 mm	0.4 mm	≤ 5 N	4.5 g

Table. Technical specifications.

Roughness of metal parts, μm , max		Rockwell hardness
Forceps, Serrated Gripping		
20G	0.16 - for external surfaces	Min 40 HRC
23G		
25G		
Forceps, ILM Endgripping		
20G	0.16 - for external surfaces	Min 40 HRC
23G		
25G		
27G		
Forceps, Asymmetrical Endgripping		
20G	0.16 - for external surfaces	Min 40 HRC
23G		
25G		
Forceps, Eckardt Endgripping		
20G	0.16 - for external surfaces	Min 40 HRC
23G		
25G		
27G		
Forceps, Kuhn 45° Power Gripping		
23G	0.16 - for external surfaces	Min 40 HRC
25G		
Forceps, Shah Xtra Grip		
23G	0.16 - for external surfaces	Min 40 HRC
25G		
27G		

Package contents

Scope of delivery:

- Disposable ophthalmic surgical forceps, 5 pcs.
- Instructions for Use.

List and description of medical device materials

Component	Material	Brand	Manufacturer
Guide	Polypropylene	Eltex MED 100-MG03 (PP)	Eltex
Handle	Dyed ABS, white	Novodur HD M 203 FC (ABS) 6% Colorant MEVOPUR-WHITE SB0M176080-ZN	Styrolution
Screw	Stainless steel	AISI304	Fox Wire Limited
Base	Dyed ABS 20G - 2% - yellow	Novadur HD M 203 FC (ABS) Yellow MEVOPUR-YELLOW SB1M176037-ZN	Styrolution

	23G - 2% - blue 25G - 2% - orange 27G - 2% - gray	Blue MEVOPUR-BLUE SB5M176202-ZN Orange MEVOPUR-ORANGE SB2M176069-ZN Gray MEVOPUR-GREY SB7M176140-ZN	
Front side	Dyed ABS 20G - 2% - yellow 23G - 2% - blue 25G - 2% - orange 27G - 2% - gray	Novadur HD M 203 FC (ABS) Yellow MEVOPUR-YELLOW SB1M176037-ZN Blue MEVOPUR-BLUE SB5M176202-ZN Orange MEVOPUR-ORANGE SB2M176069-ZN Gray MEVOPUR-GREY SB7M176140-ZN	Styrolution
Capillary tube	Cobalt-nickel alloy or stainless steel	MP35N ASTM F562 / ISO 5832-6 AISI 304	Fox Wire Limited
Forceps tip	Spring stainless steel	AISI 301	Fox Wire Limited
Spring	Spring stainless steel	AISI 301	Fox Wire Limited

The adhesive used is Loctite 4601, made by Henkel Nederland

Medical device labelling and packaging data

- The device is placed in a PETG tray; the tray protects the device from impacts, damages and other effects.
- The tray containing the device is placed in a Steriking paper pouch, which consists of an inner blister, an outer blister and a Tyvek protective film: the paper pouch ensures the device sterility.
- 5 devices placed in the tray and the pouch are packaged together in a cardboard box with an IFU: the cardboard box protects the paper bags from damages, contains information on the product name and the main specifications of the packaged devices, and is used for commercial purposes.

Packaging specifications*:

**The tolerance is $\pm 10\%$, unless otherwise specified.*

Package	Length, mm	Width, mm	Height, mm	Material/manufacturer	Manufacturer
Primary	11	60	165	PETG	Wipak
Sterile	-	10	270	Medical paper, polyethylene/polypropylene	Inpac
Multiple-unit	75	100	270	Cardboard	Oto9
Transport	The transport package specifications depend on the scope of delivery and the quantity of the devices transported.				Corrugated fiberboard

The multiple-unit package has a label in Russian, which may contain the following information:

- Device name in Russian;
- Manufacturer;
- Country of origin;
- Catalogue No.;
- Batch number;
- Registration certificate date and No.;
- Compliance mark, if declared compliant;
- Warning and informing symbols and inscriptions.

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Description of symbols

Symbol	Meaning
	Manufacturer
REF	Catalogue No.
LOT	Lot number
	Use by...
STERILE EO	Sterilized using ethylene oxide
	Do not reuse
Rx ONLY	Caution Law restricts this device to sale by or on the order of a physician
	Do not resterilize
	Do not use if package is damaged
	Consult instructions for use
	Warning Consult instructions for use

Sterilization

This medical device is sterilized using ethylene oxide and remains sterile until the expiration date indicated on the product label.

Procedure for recycling or destruction of the medical device

After use, the instrument shall be disposed of in accordance with the regulations of the medical facility regarding biologically hazardous or acute items.

Class of medical waste in accordance with SanPiN 2.1.7.2790-10:

Class A (Unused devices with expired shelf life, unused products with the sterile packaging integrity compromised);

Class B (After contact with the patient's body)

Application method and operation principles of the medical device

Forceps consist of a body, a rod and a tip. The body is the same for all instruments. When closing the body shanks, the rod moves forward and covers the tip. The instrument tip determines its function. Forceps are usually used to grip the mucous membranes/internal limiting membrane/retina/blood clots/foreign bodies in the posterior segment of the eye. Instruments are available in a range of gauges: 20, 23, 25, 27 G.

Special applications for different types of forceps:


- Forceps, Serrated Gripping: for universal (non-specific) gripping of the mucous membrane.
- Forceps, ILM Endgripping: for precision gripping and extraction of the ILM with improved visualization of thin membranes in the macular region.
- Forceps, Asymmetrical Endgripping: for precision gripping and improved visualization of the ILM, thin epiretinal membranes and membranes in patients with diabetic retinopathy and proliferative vitreoretinopathy.
- Forceps, Eckardt Endgripping: for precision gripping and improved visualization of thin membranes in the macular region.
- Forceps, Kuhn 45° Power Gripping: for removal of proliferative intravitreal, epiretinal and subretinal membranes, lens membrane, iris, small foreign bodies in the eye, fibrinous membranes and other elements. The forceps are also used to remove the ILM and epiretinal membrane and attached posterior hyaloid membranes.
- Forceps, Shah X-tra Grip: thin tips allow the ILM to be lifted, while the toothed platform grips and holds the tissue to prevent slipping. The forceps are efficient for all tissue manipulations, including the ILM tissues, epiretinal membrane, in diabetic retinopathy and proliferative vitreoretinopathy.

Preparing for Use

- Please visually inspect the packaging of the sterile disposable forceps for possible damage and expiry date. Devices from previously opened or damaged packing must be deemed unsterile.
- Unpack the disposable forceps for introduction into a sterile environment.

Directions for Use

Insert the device in the posterior segment of the eye through the entry port or incision. Insert into or retract from the eye in the closed position as visualized in the picture below.

Shelf life

5 years from the date of sterilization

Operating conditions

Requirements for environmental conditions during operation	Value
Air temperature, °C	From +15 to +30
Air humidity (relative), without condensation	30-75%

Requirements for environmental conditions during surgeries	Value
Air temperature, °C	From +32 to +42
Air humidity (relative), without condensation	100%

Transportation and storage conditions

Requirements for environmental conditions during operation
Air temperature, °C
Air humidity (relative), without condensation

Store at room temperature. Avoid prolonged exposure to temperatures higher than room temperature

Requirements for environmental conditions during transportation	Value
Air temperature, °C	From -35 to +40
Air humidity (relative), without condensation	30-95%

Warranty

VitreQ makes no warranties, expressed or implied, relating to sales of forceps and disclaims warranty liability with regard to merchantability or suitability for use. VitreQ will not be liable for injury or damage resulting from selection and use of any disposable forceps for any particular patient or any patient's condition.

Guaranteed shelf life: 5 years

Customer service

For any questions, including complaints about products, repairs and maintenance, please contact the Authorized Representative in the Russian Federation:

Company name: OpusMed Limited Liability Company

Address: 79, Ostrovskogo str., Kazan, 420107, Republic of Tatarstan

APOSTILLE

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1. Country: THE NETHERLANDS
This public document
2. has been signed by **mr. E.D. Smit**
3. acting in the capacity of notary at 's-Gravenhage
4. bears the seal/stamp of aforesaid notary

Certified

5. in Den Haag
6. on 17-08-2020
7. by the registrar of the district court of Den Haag
8. no. 2020/6384

Seal/stamp:

10. Signature:

S. Petronilia

I, EELKO DREWES SMIT, civil law notary, residing at The Hague, The Netherlands, do hereby certify that the attached document is issued for a true and correct copy of the original as referred to in article 49 paragraph 3 of the "Wet op het Notarisambt" (law regulating the notary's office) by me, which document I have seen.

Signed at The Hague, The Netherlands,
on the 14th day of August, 2020.

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