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"Operating documents" for submission in the Russian Federation:

Monitor INVOS

Manufactured by:

Covidien llc., USA. 15 Hampshire Street, Mansfield MA 02048, USA

COUNTY/CITY OF Boulder
COMMONWEALTH STATE OF Colorado
THE FOREGOING INSTRUMENTS WAS ACKNOWLEDGED
BEFORE ME THIS 20 DAY OF November 2019

Greeshma Kayala
(NAME OF PERSON SEEKING ACKNOWLEDGMENT)
Darlene Joann Johnson 29 March 2023
NOTARY PUBLIC MY COMMISSION EXPIRES

DARLENE JOANN JOHNSON
NOTARY PUBLIC
STATE OF COLORADO
NOTARY ID 20194012270
MY COMMISSION EXPIRES MARCH 29, 2023

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1. Cho H, Nemoto EM, Yonas H, Balzer J. Cerebral monitoring by oximetry And Somatosensory Evoked Potentials (SSEP) during Carotid Endarterectomy. *J Neurosurg*. 89:533-38, 1998.
2. Roberts KW, Crnkovic AP, Linneman LJ. Near infrared spectroscopy detects critical cerebral hypoxia during carotid endarterectomy in awake patients. *Anesthesiology* 1998;89(3A):A934.
3. Samra S, Dy E, Welch K, Dorje P, Zelenock G, Stanley J. Evaluation of a cerebral oximeter as a monitor of cerebral ischemia during carotid endarterectomy. *Anesthesiology* 2000;93:964-70.
4. Murkin JM, Adams SJ, Novick RJ, Quantz M, Bainbridge D, Iglesias I, Cleland A, Schaefer B, Irwin B, Fox S. Monitoring brain oxygen saturation during coronary bypass surgery: a randomized, prospective study. *Anesth Analg*. 2007 Jan;104(1):51-8.
5. Casati A, Fanelli G, Pietropaoli P, Proietti R, Tufano R, Montanini S. Monitoring cerebral oxygen saturation in elderly patients undergoing general abdominal surgery: a prospective cohort study. *Eur J Anaesthesiol*. 2006 Jul 7;1-7 [электронная публикация до печати].
6. Casati A, Fanelli G, Pietropaoli P, Proietti R, Tufano R, Danelli G, Fierro G, De Cosmo G, Servillo G. Continuous monitoring of cerebral oxygen saturation in elderly patients undergoing major abdominal surgery minimizes brain exposure to potential hypoxia. *Anesth Analg*. 2005 Sep;101(3):740-7.
7. Edmonds HL Jr, Ganzel BL, Austin EH 3rd. Cerebral oximetry for cardiac and vascular surgery. *Semin Cardiothorac Vasc Anesth*. 2004;8:147-66.
8. Goldman S, Sutter F, Ferdinand F, Trace C. Optimizing intraoperative cerebral oxygen delivery using noninvasive cerebral oximetry decreases the incidence of stroke for cardiac surgical patients. *Heart Surgery Forum* 2004;7(5):#2004-1062.
9. Murkin JM, Iglesias I, Bainbridge D, Adams S, Schaefer B, Irwin B, Fox S. Monitoring cerebral oxygen saturation significantly decreases major organ morbidity in CABG patients: A randomized blinded study. *The Heart Surgery Forum* 2004;7(6):515.
10. Yao FSF, Tseng CCA, Ho CYA, Levin SK, Illner P. Cerebral oxygen desaturation is associated with early postoperative neuropsychological dysfunction in patients undergoing cardiac surgery. *J Cardiothorac Vasc Anesth* 2004;18(5):552-558.

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1. Kim MB, Ward DS, Cartwright CR, Kolano J, Chlebowski S, Henson LC. Estimation of jugular venous O₂ saturation from cerebral oximetry or arterial O₂ saturation during isocapnic hypoxia. *J Clin Monit Comput.* 2000;16(3):191-9.
2. Kolb JC, Ainslie PN, Ide K, Poulin MJ. Effects of five consecutive nocturnal hypoxic exposures on the cerebrovascular responses to acute hypoxia and hypercapnia in humans. *J Appl Physiol.* 2004 May;96(5):1745-54. Epub 2004 Jan 16.

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 Covidien Inc.

15 Hampshire Street, Mansfield, MA 02048 USA.

 Covidien Ireland Limited, IDA Business & Technology Park, Tallamore.

www.Medtronic.com





Medtronic

REF PMPAMP71

INVOS™

Preamplifier

For use with PM7100

PT00017912 Rev C 2018-06

Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic.

Other brands are trademarks of a Medtronic company.

© 2016 Medtronic. All rights reserved.

Covidien Inc, 15 Hampshire Street, Mansfield, MA 02048 USA.
Covidien Ireland Limited, IDA Business & Technology Park,
Tullamore, Ireland.

www.medtronic.com

Rx
ONLY


Medtronic

REF PMAC71RSC

INVOS™

Reusable Sensor Cable for PM7100

Instructions for Use

Part No. 10147644 Rev E 2018-06

Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic.

Other brands are trademarks of a Medtronic company.

© 2016 Medtronic. All rights reserved.
Covidien Inc, 15 Hampshire Street, Mansfield, MA 02048 USA.
Covidien Ireland Limited, IDA Business & Technology Park,
Tullamore, Ireland.

www.medtronic.com

Rx
ONLY


Medtronic







INVOS™

REF PMSENS71-A

Adult rSO₂ Sensor

en	Instructions For Use
en-GB	Instructions For Use
fr	Mode d'emploi
de	Gebrauchsanleitung
nl	Gebruiksaanwijzing
it	Istruzioni per l'uso
es	Instrucciones de uso
sv	Bruksanvisning
da	Brugsanvisning
no	Bruksanvisning
fi	Käyttöohjeet
pt	Instruções de utilização
ru	Инструкция по применению
pl	Instrukcja użycia
cs	Návod k použití
sk	Návod na použitie
hu	Használati utasítás
el	Οδηγίες χρήσης
tr	Kullanma Talimatı
ro	Instrucțiuni de utilizare

INVOS™ Adult rSO₂ Sensor

The INVOS™ adult rSO₂ sensor has been designed for use with the INVOS™ 7100 and 5100C cerebral/somatic monitoring systems for monitoring of site-specific regional oxygen saturation (rSO₂) in adult patients weighing > 40 kilograms.

Indications for Use

The INVOS™ Adult rSO₂ Sensor is indicated for single patient use when cerebral/somatic monitoring of site-specific regional oxygen saturation (rSO₂) is required in patients weighing >40 kilograms. This sensor is only intended to be used with INVOS™ Near Infrared Spectroscopy (NIRS) technology including monitoring systems and devices integrated with INVOS™ NIRS technology. For additional information regarding setup and use of the INVOS™ System including indications for use, contraindications, warnings and cautions, consult the *Monitoring System Operator's Manual*.

Contraindications

The INVOS™ Adult rSO₂ sensor is contraindicated for use on patients who exhibit allergic reactions to the adhesive tape.

WARNINGS

- The sensor is designed for external use only as described in the instructions. Do not use the sensors internally or over compromised skin for any reason.
- Do not use the INVOS™ System in the presence of flammable anesthetics or in other flammable environments.
- The INVOS™ System should not be used as the sole basis for diagnosis or therapy. It is intended only as an adjunct in patient assessment. It must be used in conjunction with clinical signs and symptoms.
- Disconnect the sensor from the patient during magnetic resonance imaging (MRI) scanning. Induced current may cause burns.
- Environments with excessive ambient light such as bright sunlight or strong operating room lighting may require loosely covering the area of the sensor with an opaque drape.
- Select sensor site over tissue area of interest (site selection will determine which body region is monitored). Avoid placing the sensor over thick fatty deposits, hair, or bony protuberances. Do not place the sensor over nevi, hematomas, or broken skin, as this may cause readings that are not reflective of tissue or no readings at all. Sensor location is at the clinician's discretion, provided it adheres to the criteria noted in this Instruction For Use.

CAUTIONS

- Make sure all connectors are fully engaged and free from moisture. Moisture intrusion may cause poor performance or no readings at all.
- Use of an electrosurgical/electrocautery instrument in the vicinity of the INVOS™ System may interfere with the signal and result in no readings.
- Use care when placing or removing sensor. Do not place on broken or undeveloped skin.
- To avoid pressure sores do not apply pressure (e.g. headbands, wraps, tape, or any additional adhesives) to sensor.
- Do not immerse the INVOS™ System or sensor in any liquids as they may cause electric shock hazard or damage the device.
- If two sensors are placed in close proximity to each other on a patient, the same preamplifier should be connected to both sensors to avoid poor performance.

Instructions for Use


Pre-use Check

Remove the sensor from the package and examine for visual signs of damage. If any signs of damage are observed, select another sensor.

Cerebral Site Selection

Select sensor site on the right and left side of forehead. Placement of the sensor in other cerebral locations, or over hair, may cause poor performance or no readings at all. Do not place the sensor over nevi, sinus cavities, the superior sagittal sinus, subdural or epidural hematomas, or other anomalies such as arteriovenous malformations, as this may cause readings that are not reflective of brain tissue or no readings at all.

EEG and level of consciousness sensors can be used simultaneously with the INVOS™ adult rSO₂ sensors, overlapping sensor pads if necessary, as long as the three optical windows on the INVOS™ sensor are not blocked.

If space is limited, the INVOS™ adult rSO₂ sensor pads may overlap as long as the three optical windows on each sensor pad remain unobstructed.

Somatic Site Selection

Select sensor site over tissue area of interest (site selection will determine which body region is monitored). Avoid placing the sensor over thick fatty deposits, hair, or bony protuberances. Do not place the sensor over nevi, hematomas, or broken skin, as this may cause readings that are not reflective of tissue or no readings at all. Sensor location is at the clinician's discretion, provided it adheres to the criteria noted in this *Instruction For Use*.

Placements may include:

- Over the spine (upper and lower) **4**
- On the lateral or posterior calf **5**
- On the upper arm or leg
- On the forearm
- On the chest

Patient Preparation

Remove any moisture or perspiration from the patient's skin with a dry gauze pad. Then, degrease the skin with isopropyl alcohol pad. Ensure patient's skin is completely dry and remove degreaser residue, if any, with a dry gauze pad **1**.

Sensor Placement

Slowly peel the protective backing **2** label from the adhesive side of the sensor, and apply the sensor to the skin **3**. Continue applying the sensor by smoothing it to the skin from the center outward. Ensure edges of the sensor are sealed to prevent light from entering.

Monitoring

Plug the sensor into the cable connector **6**. Secure the sensor cable to a fixed object to avoid strain on the sensor-to-skin interface using strain-relief clips. Ensure the reusable cable is properly inserted into the preamplifier. Calibration is automatic and monitoring will begin momentarily.

Status messages on the INVOS™ system display will appear if monitoring conditions are compromised. Check the sensor periodically to assess skin integrity (irritation or injury) and the security of sensor placement. Always ensure sensor is properly sealed to skin to reduce light interference. **For extended monitoring, Medtronic recommends using a new sensor every 24 hours or if adhesive is inadequate to seal the sensor to the skin.**

If present, the following may cause poor performance:

- "Cardiogreen," "Indigo Carmine," "Methylene Blue" or other intravascular dyes
- Carboxyhemoglobin or other dyshemoglobins
- Hemoglobinopathies
- Conjugated hyperbilirubinemia (direct)
- Myoglobin (Mb), hemoglobin from muscular tissue, in the blood.
- Dark skin pigment
- Externally applied coloring agents (dye, pigmented cream)

Cleaning/Disinfection/Sterilization

The sensor is intended for single use only and cannot be adequately cleaned, disinfected, or sterilized for safe reuse. Attempts to clean, disinfect, or sterilize the sensor for reuse may result in product failure and/or safety risk to the patient.

Disposal

Dispose of the sensor after a single use. Follow local governing ordinances regarding disposal of biologically hazardous waste.

Environmental Specification	
Operational temperature range	10°C to 35°C
Storage temperature range	10°C to 30°C
Transport temperature range	-10°C to 50°C up to 72 hours
Operational relative humidity range	15% to 95%
Storage relative humidity range	10% to 95%
Atmospheric pressure range	61.6 kPa to 107.5 kPa
Storage relative pressure range	46.5 kPa to 107.5 kPa

INVOS™

Adult rSO₂ Sensor

The INVOS™ adult rSO₂ sensor has been designed for use with the INVOS™ 7100 and 5100C cerebral/somatic monitoring systems for monitoring of site-specific regional oxygen saturation (rSO₂) in adult patients weighing > 40 kilograms.

Indications for Use

The INVOS™ Adult rSO₂ Sensor is indicated for single patient use when cerebral/somatic monitoring of site-specific regional oxygen saturation (rSO₂) is required in patients weighing > 40 kilograms. This sensor is only intended to be used with INVOS™ Near Infrared Spectroscopy (NIRS) technology including monitoring systems and devices integrated with INVOS™ NIRS technology. For additional information regarding set-up and use of the INVOS™ System including indications for use, contraindications, warnings and cautions, consult the *Monitoring System Operator's Manual*.

Contraindications

The INVOS™ Adult rSO₂ sensor is contraindicated for use on patients who exhibit allergic reactions to the adhesive tape.

WARNINGS

- The sensor is designed for external use only as described in the instructions. Do not use the sensors internally or over compromised skin for any reason.
- Do not use the INVOS™ System in the presence of flammable anaesthetics or in other flammable environments.
- The INVOS™ System should not be used as the sole basis for diagnosis or therapy. It is intended only as an adjunct in patient assessment. It must be used in conjunction with clinical signs and symptoms.
- Disconnect the sensor from the patient during magnetic resonance imaging (MRI) scanning. Induced current may cause burns.
- Environments with excessive ambient light such as bright sunlight or strong operating room lighting may require loosely covering the area of the sensor with an opaque drape.
- Select sensor site over tissue area of interest (site selection will determine which body region is monitored). Avoid placing the sensor over thick fatty deposits, hair or bony protuberances. Do not place the sensor over nevi, hematomas or broken skin, as this may cause readings that are not reflective of tissue or no readings at all. Sensor location is at the clinician's discretion, provided it adheres to the criteria noted in these Instructions For Use.

CAUTIONS

- Make sure all connectors are fully engaged and free from moisture. Moisture intrusion may cause poor performance or no readings at all.
- Use of an electrosurgical/electrocautery instrument in the vicinity of the INVOS™ System may interfere with the signal and result in no readings.
- Use care when placing or removing the sensor. Do not place on broken or undeveloped skin.
- To avoid pressure sores, do not apply pressure (e.g. headbands, wraps, tape or any additional adhesives) to the sensor.
- Do not immerse the INVOS™ System or sensor in any liquids as they may cause electric shock hazard or damage the device.
- If two sensors are placed in close proximity to each other on a patient, the same preamplifier should be connected to both sensors to avoid poor performance.


Instructions for Use**Pre-use Check**

Remove the sensor from the package and examine for visual signs of damage. If any signs of damage are observed, select another sensor.

Cerebral Site Selection

Select sensor site on the right and left sides of the forehead. Placement of the sensor in other cerebral locations or over hair may cause poor performance or no readings at all. Do not place the sensor over nevi, sinus cavities, the superior sagittal sinus, subdural or epidural hematomas or other anomalies such as arteriovenous malformations, as this may cause readings that are not reflective of brain tissue or no readings at all.

EEG and level of consciousness sensors can be used simultaneously with the INVOS™ adult rSO₂ sensors, overlapping sensor pads if necessary, as long as the three optical windows on the INVOS™ sensor are not blocked.

If space is limited, the INVOS™ adult rSO₂ sensor pads may overlap as long as the three optical windows on each sensor pad remain unobstructed.

Somatic Site Selection

Select sensor site over tissue area of interest (site selection will determine which body region is monitored). Avoid placing the sensor over thick fatty deposits, hair or bony protuberances. Do not place the sensor over nevi, hematomas or broken skin, as this may cause readings that are not reflective of tissue or no readings at all. Sensor location is at the clinician's discretion, provided it adheres to the criteria noted in these *Instructions For Use*.

Placements may include:

- Over the spine (upper and lower) **4**
- On the lateral or posterior calf **5**
- On the upper arm or leg
- On the forearm
- On the chest

Patient Preparation

Remove any moisture or perspiration from the patient's skin with a dry gauze pad. Then, degrease the skin with an isopropyl alcohol pad. Ensure that the patient's skin is completely dry and remove degreaser residue, if any, with a dry gauze pad **1**.

Sensor Placement

Slowly peel the protective backing **2** label from the adhesive side of the sensor and apply the sensor to the skin **3**. Continue applying the sensor by smoothing it to the skin from the centre outward. Ensure that the edges of the sensor are sealed to prevent light from entering.

Monitoring

Plug the sensor into the cable connector **6**. Secure the sensor cable to a fixed object to avoid strain on the sensor-to-skin interface using strain-relief clips. Ensure that the reusable cable is properly inserted into the preamplifier. Calibration is automatic and monitoring will begin shortly.

Status messages on the INVOS™ system display will appear if monitoring conditions are compromised. Check the sensor periodically to assess skin integrity (irritation or injury) and the security of sensor placement. Always ensure that the sensor is properly sealed to the skin to reduce light interference. **For extended monitoring, Medtronic recommends using a new sensor every 24 hours or if adhesive is inadequate to seal the sensor to the skin.**

If present, the following may cause poor performance:

- "Cardiogreen", "Indigo Carmine", "Methylene Blue" or other intravascular dyes
- Carboxyhaemoglobin or other dysaemoglobins
- Haemoglobinopathies
- Conjugated hyperbilirubinemia (direct)
- Myoglobin (Mb), haemoglobin from muscular tissue, in the blood.
- Dark skin pigment
- Externally applied colouring agents (dye, pigmented cream)

Cleaning/Disinfection/Sterilisation

The sensor is intended for single use only and cannot be adequately cleaned, disinfected or sterilised for safe reuse. Attempts to clean, disinfect or sterilise the sensor for reuse may result in product failure and/or a safety risk to the patient.

Disposal

Dispose of the sensor after a single use. Follow local governing ordinances regarding disposal of biologically hazardous waste.

Environmental Specification	
Operational temperature range	10°C to 35°C
Storage temperature range	10°C to 30°C
Transport temperature range	-10°C to 50°C up to 72 hours
Operational relative humidity range	15% to 95%
Storage relative humidity range	10% to 95%
Atmospheric pressure range	61.6 kPa to 107.5 kPa
Storage relative pressure range	46.5 kPa to 107.5 kPa

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