PEDRO-ÁNGEL CASADO MARTÍN
PUBLIC NOTARY
Registered with number 43/2021 in
my Registered Book of Certifications.

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CAUTIONARY MEASURES.

- Perform appropriate treatment of the patient using anticoagulants, vasodilators, etc., according to the procedure for using intravascular catheters.
- Before use, check the size and suitability of the catheter for its intended use.
- Proceed with extreme care so as not to damage the catheter while advancing it.
- The balloon must be inflated with a mixture of saline solution and contrast medium (preferably 1/1, though the proportion of contrast medium can be reduced in large balloons).
- Do not exceed the calculated burst pressure as the balloon could burst.
- If any resistance is noticed during advance, stop and determine its cause before continuing.
- Use with standard 0.014" guidewires and standard introducers. If the introducer length goes above 25cm or the braided type is used, it may be necessary to increase its size to avoid catheter friction.
- If any resistance to withdrawing the catheter is noted, it is recommended that the balloon catheter, the guidewire and introducer be removed as a single unit.
- After use, this product may represent a biohazard. Handle and dispose it in accordance with accepted medical practices and pertinent local, state or federal laws and regulations.
- Store in a cool, dry place away from direct sunlight.

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Catheter balloon (diameter x length in mm)	Number of radiopaque markers/marks	Distance from the distal end of the balloon to the distal mark (mm)	Distance from the balloon distal end to the proximal mark (mm)
1,25 x 10; 1,25 x 15; 1,25 x 20; 1,5 x 10; 1,5 x 15; 1,5 x 20; 2,0 x 10; 2,0 x 15; 2,0 x 20	One mark in the center	One mark in the center	One mark in the center
1,5 x 40	Two marks	1,42	41,42
1,5 x 60	Two marks	1,42	61,42
1,5 x 80	Two marks	1,42	81,42
2,0 x 40	Two marks	2,21	42,21
2,0 x 60	Two marks	2,21	62,21
2,0 x 80	Two marks	2,21	82,21
2,0 x 120	Two marks	2,21	122,21
2,0 x 150	Two marks	2,21	152,21
2,0 x 200	Two marks	2,21	202,21
2,5 x 40	Two marks	2,95	42,95
2,5 x 60	Two marks	2,95	62,95
2,5 x 80	Two marks	2,95	82,95
2,5 x 120	Two marks	2,95	122,95
2,5 x 150	Two marks	2,95	152,95
2,5 x 200	Two marks	2,95	202,95
3,0 x 40	Two marks	3,73	43,73
3,0 x 60	Two marks	3,73	63,73
3,0 x 80	Two marks	3,73	83,73
3,0 x 120	Two marks	3,73	123,73
3,0 x 150	Two marks	3,73	153,73
3,0 x 200	Two marks	3,73	203,73
3,5 x 40	Two marks	4,52	44,52
3,5 x 60	Two marks	4,52	64,52

3,5 x 80	Two marks	4,52	84,52
3,5 x 120	Two marks	4,52	124,52
3,5 x 150	Two marks	4,52	154,52
3,5 x 200	Two marks	4,52	204,52
4,0 x 40	Two marks	5,28	45,28
4,0 x 60	Two marks	5,28	65,28
4,0 x 80	Two marks	5,28	85,28
4,0 x 120	Two marks	5,28	125,28

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Part name of the device	Material of manufacture and manufacturer
1) Luer lock port	Polycarbonate (color: clear) Mark on the port surface: (laser application)
2) Connecting tube/strain relief	Pebax, color: white;
3) Catheter shaft: Front end Internal lumen and lumen for guidewire insertion Tip	Nylon 12; color: blue; Nylon 12/Plexar/HDPE; color: clear; Nylon 12; color: blue;
4) Balloon	Nylon 12; color: clear;
5) Protective balloon tube	Nylon 12; color: clear; Plexar/HDPE
6) Balloon radiopaque marker	Alloy platinum and iridium Tungsten + Pebax 63333
7) Reinforcing shaft	Stainless steel
8) Catheter coating: the hydrophilic coating of the catheter	HYDRAX PLUS
PACKAGING:	
9) Catheter protective case	High-density polyethylene HDPE

10) Catheter protective case clip	High-density polyethylene HDPE
11) Blister	LLDPE/Tyvek
12) Consumer packaging	Coated cardboard

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plunger on the inflation device back to deflate the balloon. Keep up the negative pressure for between 20 seconds depending on the size of the balloon. Make sure the balloon is fully deflated (by means of fluoroscopy) removing the catheter.

negative pressure in the inflation device and with the guidewire in position, carefully retract the catheter until through the introducer sheath.

Inject contrast liquid through the introducer to check the treated lesion before removing the guidewire.

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17.1. Transportation.

PTA balloon dilation catheter iVascular OCEANUS 14 PRO is transported by all modes of transport in covered vehicles in accordance with the rules of cargo transportation applicable to each type of transport.

Transport temperature from 0°C to +40 °C, at a relative humidity no more than 75%.

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THE MANUFACTURER GUARANTEES: sterility of the product, as well as its compliance with the requirements of international and national standards during the entire shelf life of the product, subject to the rules of transportation,

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 iVascular	- brand logo and manufacturer's name	0.00 x 00mm	The size of the balloon catheter format: balloon diameter (mm) x balloon length in (mm)
		LOT	The digital code of the lot.

WILSON

TITLES

Information about the number of radioopaque markers of the balloon	Hydrophilic coating	Number of devices per package	OTW describes the method of delivery of the catheter through the guide-wire.
Barcode or QR code	Information about the manufacturer's authorized representative in the Russian Federation	Registration certificate number and date	
THE FOLLOWING INFORMATION IS GIVEN IN TABLE FORM			
The extensibility curve, in the form of values of the balloon diameter, achieved with each step-by-step increase in the created pressure.			

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Martín, on 20 of
November of the year
2013, with number of
official record 2281. --

I, the Public
Notary, attest that I
understand the English
language.-----

Likewise I, the
Notary, certify the
identity of the
signatory party and by
its virtue I state that
the previous signature
belongs to Mrs.
M^aEUGENIA VILLANUEVA
PÉREZ.-----

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**PTA balloon dilatation catheter iVascular
OCEANUS 35****Manufacturing company:**

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Date of entry into force:

Version: -28-/July/2020

Action	FULL NAME	Position	Date	Signature
Prepared by	Aurea Roca	Regulatory Affairs Director	28/07/2020	
Approved by	M. Eugenia Villanueva	Regulatory Compliance Manager	28/07/2020	M ^a Eugenia Villanueva

LS

LIFE VASCULAR DEVICES BIOTECH S.L.

2020

1. PURPOSE OF THE DEVICE..

PTA balloon dilation catheter iVascular OCEANUS 35 is intended for percutaneous transluminal angioplasty of peripheral arteries.

2. RISK CLASS OF THE MEDICAL DEVICE.

The potential risk class for the use of a medical device is 3, in accordance with Directive 93/42/EEC. Type of contact with the human body is a product that briefly (for less than 24 hours) comes into contact with blood.

3. LIST OF INTERNATIONAL REGULATIONS/STANDARDS THAT THE MEDICAL DEVICE MEETS.

STANDARD CODE DESIGNATION	STANDARD NAME
ISO 13485	Medical devices. Quality management system. Requirements for regulatory purposes.
EN ISO 14971	Medical devices. The application of risk management to medical devices
UNE-EN 1041+A1	Information prepared by the manufacturer accompanying medical devices
EN 556-1	Sterilization of medical devices. Requirements for medical devices of the "sterile" category. Part 1. Requirements for medical devices subject to final sterilization.
UNE-EN ISO 11135	Sterilization of medical devices. Ethylene oxide. Part 1: Requirements for the development, validation, and ongoing monitoring of the sterilization process for medical devices.
UNE-EN ISO 11737-2	Sterilization of medical devices. Microbiological method. Part 2: Sterility tests performed at the stages of technical development, validation and continuous maintenance of the sterilization process.
UNE-EN ISO 14644-1	Clean rooms and their associated controlled areas. Part 1: Classification of air cleanliness.
UNE-EN ISO 14644-2	Clean rooms and their associated controlled areas. Part 2: Specifications for testing and monitoring to ensure continuous compliance with ISO 14644-1.
UNE-EN ISO 14698-2	Clean rooms and their associated controlled areas. Control of bio-contamination. Part 2: Analysis and interpretation of biosafety data.
UNE-EN ISO 14644-3	Clean rooms and their associated controlled areas. Part 3: Methods of control.
UNE-EN ISO 14698-1	Clean rooms and their associated controlled areas. Control of bio-contamination. Part 1: General principles and methods.
MDD 93/42/EEC	Directive 93/42/EEC Medical devices, devices and equipment.
UNE EN ISO 10555-1	Single-use sterile intravascular catheters. Part 1: General requirements.
UNE EN ISO 10555-4	Single-use sterile intravascular catheters. Part 4: Balloon dilation catheters.
UNE EN 20594-1	Conical fittings with a 6% Luer lock cone for syringes, needles and some other medical equipment. Part 1: General requirements.
UNE EN 11607-1	Packaging for medical devices that have undergone final sterilization. Part 1: requirements for materials, sterility maintenance systems, and packaging systems.
UNE-EN 62366	Medical devices: Application of serviceability design to medical devices.

UNE EN ISO 10993-1	Assessment of the biological effect of medical devices. Part 1. Analysis and testing in the process of implementing risk management.
UNE EN ISO 10993-4	Assessment of the biological effect of medical devices. Part 4. Selection of tests for interaction with blood.
UNE EN ISO 10993-5	Assessment of the biological effect of medical devices. Part 5. In vitro cytotoxicity tests.
UNE EN ISO 10993-7	Assessment of the biological effect of medical devices. Part 7. Residual amounts of ethylene oxide after the sterilization process.
UNE EN ISO 10993-9	Assessment of the biological effect of medical devices. Part 9. Standards for authenticating and quantifying potential degradation products.
UNE EN ISO 10993-10	Assessment of the biological effect of medical devices. Part 10. Skin irritation and sensitization tests.
UNE EN ISO 10993-11	Assessment of the biological effect of medical devices. Part 11. Tests for general toxic effects.
UNE EN ISO 10993-12	Assessment of the biological effect of medical devices. Part 12. Preparation of test samples and standard samples.
UNE EN ISO 10993-13	Assessment of the biological effect of medical devices. Part 13. Authenticity and quantification of degradation products when using polymer medical devices.
UNE EN ISO 10993-18	Assessment of the biological effect of medical devices. Part 18. Chemical characteristics of materials.
UNE EN ISO 11737-1	Sterilization of medical devices. Microbiological method. Part 1. Determination of microbial populations on products.

4. INFORMATION ON THE MANUFACTURER OF THE MEDICAL DEVICE

Manufacturing company:

LIFE VASCULAR DEVICES BIOTECH, S.L.
(LIFE VASCULAR DIVISION BIOTECH, Spain)
Camí de Can Ubach, 11; Pol. Ind. Les Fallulles
08620 Sant Vicenç Dels Horts- Barcelona, Spain.

Manufacturing organizations:

LIFE VASCULAR DEVICES BIOTECH, S.L.
(LIFE VASCULAR DIVISION BIOTECH, Spain)
Camí de Can Ubach, 11; Pol. Ind. Les Fallulles
08620 Sant Vicenç Dels Horts- Barcelona, Spain.
Tel/fax +34936724711,
E-mail: info@ivascular.es

5. DELIVERY SET.

- 1) Balloon catheter - 1 pc.
- 2) Instructions for use - 1 pc.
- 3) Tracking label - 6 pcs

6. INFORMATION ABOUT POTENTIAL CONSUMERS (USERS) OF A MEDICAL DEVICE

Qualified personnel of vascular surgery departments of medical treatment and preventive care institutions.

7. INFORMATION ABOUT THE VARIANTS OF THE MEDICAL DEVICE.

I.PTA balloon dilatation catheter iVascular OCEANUS 35:

1. PTA balloon dilatation catheter iVascular OCEANUS 35, with the guidewire of 0,035 inch, catheter length of 80 cm, balloon dimensions: 3,0 x 20mm; 3,0 x 40mm; 3,0 x 60mm; 3,0 x 80mm; 3,0 x 120mm; 3,0 x 150mm; 3,0 x 200mm; 4,0 x 20mm; 4,0 x 40mm; 4,0 x 60mm; 4,0 x 80mm; 4,0 x 120mm; 4,0 x 150mm; 4,0 x 200mm; 5,0 x

20mm; 5,0 x 40mm; 5,0 x 60mm; 5,0 x 80mm; 5,0 x 120mm; 5,0 x 150mm; 5,0 x 200mm; 6,0 x 20mm; 6,0 x 40mm; 6,0 x 60mm; 6,0 x 80mm; 6,0 x 120mm; 6,0 x 150mm; 6,0 x 200mm; 7,0 x 20mm; 7,0 x 40mm; 7,0 x 60mm; 7,0 x 80mm; 7,0 x 120mm; 8,0 x 20mm; 8,0 x 40mm; 8,0 x 60mm; 8,0 x 80mm; 9,0 x 20mm; 9,0 x 40mm; 9,0 x 60mm; 9,0 x 80mm; 10,0 x 20mm; 10,0 x 40mm; 10,0 x 60mm; 12,0 x 20mm; 12,0 x 40mm; 12,0 x 60mm.

2. PTA balloon dilatation catheter iVascular OCEANUS 35, with the guidewire of 0,035 inch, length of 140 cm, balloon dimensions: 3,0 x 20mm; 3,0 x 40mm; 3,0 x 60mm; 3,0 x 80mm; 3,0 x 120mm; 3,0 x 150mm; 3,0 x 200mm; 4,0 x 20mm; 4,0 x 40mm; 4,0 x 60mm; 4,0 x 80mm; 4,0 x 120mm; 4,0 x 150mm; 4,0 x 200mm; 5,0 x 20mm; 5,0 x 40mm; 5,0 x 60mm; 5,0 x 80mm; 5,0 x 120mm; 5,0 x 150mm; 5,0 x 200mm; 6,0 x 20mm; 6,0 x 40mm; 6,0 x 60mm; 6,0 x 80mm; 6,0 x 120mm; 6,0 x 150mm; 6,0 x 200mm; 7,0 x 20mm; 7,0 x 40mm; 7,0 x 60mm; 7,0 x 80mm; 7,0 x 120mm; 8,0 x 20mm; 8,0 x 40mm; 8,0 x 60mm; 8,0 x 80mm; 9,0 x 20mm; 9,0 x 40mm; 9,0 x 60mm; 9,0 x 80mm; 10,0 x 20mm; 10,0 x 40mm; 10,0 x 60mm; 12,0 x 20mm; 12,0 x 40mm; 12,0 x 60mm.
3. PTA balloon dilatation catheter iVascular OCEANUS 35 with the guidewire of 0,035 inches, catheter length of 200 cm, balloon dimensions: 6,0 x 20mm; 7,0 x 20mm; 8,0 x 20mm; 9,0 x 20mm; 10,0 x 20mm; 12,0 x 20mm.

8. DESCRIPTION OF THE MEDICAL DEVICE AND ITS PHOTOGRAPHIC IMAGES.

PTA balloon dilation catheter iVascular OCEANUS 35 is used for dilation and is used to inflate the balloon for stenosis or occlusion of peripheral arteries in order to recanalize the artery.

This procedure is performed using minimally invasive techniques, also called percutaneous transluminal angiography (PTA) of the peripheral arteries in which the puncture of the outer skin allows access to the blood flow and to the affected area.

The catheter has a coaxial design and is a dual-lumen catheter (from the luer connector to the tip), also called an OTW catheter, i.e. a catheter with a stent delivery system along the guidewire.


Fig. 1. Schematic representation of a catheter.

The tip is at the distal end of the catheter, it is atraumatic in shape to avoid damage to the arteries during its insertion. In the distal part of the catheter, in front of its tip, there is a balloon, or an inflatable segment that allows you to pump the pressure to various values. The use of a balloon leads to the expansion of the artery by pumping a contrast medium into it.

One or two radiopaque markers are placed: one marker - in the center of the balloon, two markers: one on each of the tapered tips of the balloon, in order to mark its length and help the user determine the location of the catheter when it is inside the patient (when displayed on the x-ray screen).

Between the Luer lock connector and the catheter body is a connecting tube (strain relief), which is necessary to prevent the catheter from bending at this point.

The connector has a branched (Y-shaped) shape and is provided with two ports directly connected to the catheter shaft:

- the side port (inflating port) allows for transmission contrast fluid to the catheter's inflation and deflating lumen (the lumen for inflating);
- the direct port (guidewire insertion port) is intended for insertion of the catheter guidewire. A guidewire is inserted

the straight port of the connector and exits the catheter tip in order to reach the affected area.
OCEANUS 35, the maximum guidewire diameter should not exceed 0.89 mm (0.035 in).


Fig. 2. Enlarged image of the proximal part of the PTA balloon dilation catheter iVascular OCEANUS 35.

Fig. 3. Enlarged image of the middle part of the PTA balloon dilation catheter iVascular OCEANUS 35.


Fig. 4. Enlarged image of the part at the junction of the tube with the balloon of the PTA balloon dilation catheter iVascular OCEANUS 35.


Fig. 5. Enlarged image of the distal part of the PTA balloon dilation catheter iVascular OCEANUS 14 PRO


Fig. 6. Image of a general view of the medical device: PTA balloon dilation catheter iVascular OCEANUS 35 in a protective case.

9. TERMS OF USE.

Only qualified specialists of the surgical departments of the treatment and prevention institutions.

10. SCOPE OF APPLICATION.

Vascular surgery.

11. INDICATIONS FOR USE. CONTRAINDICATIONS. PRECAUTIONARY MEASURES.

INDICATIONS.

PTA balloon dilation catheter iVascular OCEANUS 35, with hydrophilic coating is used for dilation of stenoses in the femoral, popliteal arteries, in the end branches of the popliteal arteries, in the renal arteries, as well as for the treatment of obstructive lesions of native or synthetic arteriovenous dialysis fistulae, whose diameters from 3,0 to 12,0 mm, and lengths from 20 to 200 mm.

CONTRAINDICATIONS.

Inability to pass the conductor through the stenosed area.

POSSIBLE SIDE EFFECTS/COMPLICATIONS.

Possible side effects and/or complications that may occur before, during, or after the catheterization procedure are:

- Death

- Stroke/embolism/thrombosis
- Haemodynamic deterioration
- Sharp elastic rebound
- Arrhythmias
- Aneurysm or pseudo-aneurysm
- Infections
- Complete arterial occlusion
- Reocclusion of the treated area; restenosis
- Spasm
- Perforation or dissection of the treated area
- Arteriovenous fistula
- Local hemorrhage with hematoma at the access point
- Hypo/hypertension
- Allergic reactions to materials
- Increase / decrease in blood pressure

PRECAUTIONARY MEASURES.

- Perform appropriate treatment of the patient using anticoagulants, vasodilators, etc., according to the procedure for inserting intravascular catheters.
- Before use, check the size and suitability of the catheter for its intended use.
- Proceed with extreme care so as not to damage the catheter while advancing it.
- The balloon must be inflated with a mixture of saline solution and contrast medium (preferably 1/1, though the proportion of contrast medium can be reduced in large balloons).
- Do not exceed the calculated burst pressure as the balloon could burst.
- If any resistance is noticed during advance, stop and determine its cause before continuing.
- Use with standard 0,035" guidewires and standard introducers. If the introducer length goes above 25cm or the braided type is used, it may be necessary to increase its size to avoid catheter friction.
- If any resistance to withdrawing the catheter is noted, it is recommended that the balloon catheter, the guidewire and introducer be removed as a single unit.
- After use, this product may represent a biohazard. Handle and dispose it in accordance with accepted medical practices and pertinent local, state or federal laws and regulations.
- Store in a cool, dry place away from direct sunlight.

WARNINGS.

- The device should only be used by experienced doctors familiar in the percutaneous transluminal angioplasty (PTA) technique.
- The percutaneous procedure must be carried out in hospitals equipped with emergency facilities to conduct open heart surgery or hospitals with easy access to other hospitals where this type of operation can be performed.
- The product is for single use. It must not be re-sterilised or re-used once the procedure has been completed. Re-using the product in another patient may lead to cross contamination, infection or transmission of infectious diseases from one patient to another. Re-use of the product may cause alterations to it and affect its effectiveness.
- The product is supplied sterile. Check the expiry date and do not use products which are past this date.
- Inspect the packaging before opening it. If there are any defects or the packaging is damaged, do not use the product.
- Use aseptic techniques when the product is taken out of the primary packaging.
- Select the balloon diameter according to the diameter of the artery to be treated. Do not insert a balloon with a diameter larger than the artery to be treated.
- Do not dry with gauze.
- Do not expose the product to organic solvents or to contrast media not suitable for intravascular use.
- Do not use air or gaseous media to inflate the balloon.
- Advance the product on the guidewire using fluoroscopy. Do not allow the product to advance without the guidewire inside it.
- Do not manipulate, advance or retract either the catheter or the wire guide while the balloon is inflated.

Table 1. The principal dimensions of the catheter.

Effective catheter length (mm)	Balloon diameter (mm)	Effective catheter length (mm)
800 1400	3,0	20 ⁺
	4,0	40
	5,0	60
	6,0	80
		120
		150
800 1400	7,0	200
		20
		40
		60
		80
800 1400	8,0	120
		20
		40
		60
800 1400	9,0	80
		20
		40
800 1400	10,0	60
		20
		40
		60
		20
2000	6,0	20
	7,0	
	8,0	
	9,0	
	10,0	
	12,0	

Table 2. Data on external and internal catheter diameters, depending on the parameters of the balloon.

The cylinders of the catheter (diameter x length mm)	External proximal catheter diameter (mm) [Fr]	External distal catheter diameter (mm) [Fr]	Internal catheter diameter (mm) [Inches]
3,0 x 20; 3,0 x 40; 3,0 x 60; 3,0 x 80; 3,0 x 120; 3,0 x 150; 3,0 x 200; 4,0 x 20; 4,0 x 40; 4,0 x 60; 4,0 x 80; 4,0 x 120; 4,0 x 150; 4,0 x 200; 5,0 x 20; 5,0 x 40; 5,0 x 60; 5,0 x 80; 5,0 x 120; 5,0 x 150; 5,0 x 200; 6,0 x 20; 6,0 x 40; 6,0 x 60; 6,0 x 80; 6,0 x 120; 6,0 x 150; 6,0 x 200; 7,0 x 20; 7,0 x 40; 7,0 x 60; 7,0 x 80;	1,64 [5]	1,64 [5]	0,95 [0,0374]
7,0 x 120; 8,0 x 20; 8,0 x 40; 8,0 x 60; 8,0 x 80; 9,0 x 20; 9,0 x 40; 9,0 x 60; 9,0 x 80; 10,0 x 20; 10,0 x 40; 10,0 x 60; 12,0 x 20; 12,0 x 40; 12,0 x 60;	1,75 [5,3]	1,75 [5,3]	0,95 [0,0374]

Table 3. Data of the balloon cross-profile depending on the diameter of the balloon.

Balloon diameter (mm)	Inches	mm	Fr
3,0	0,057	1,45	4,3
4,0	0,059	1,50	4,5

5,0	0,061	1,54	4,6
6,0	0,065	1,66	4,9
7,0	0,070	1,79	5,4
8,0	0,072	1,82	5,5
9,0	0,073	1,85	5,5
10,0	0,077	1,95	5,8
12,0	0,083	2,10	6,3

Table 4. Data on the location of radiopaque markers/marks of catheter balloons.

Catheter balloon (diameter x length in mm)	Number of radiopaque markers/marks	Distance from the distal end of the balloon to the distal mark (mm)	Distance from the balloon distal end to the proximal mark (mm)
3,0 x 20;	Two marks	2,08	22,08
3,0 x 40;	Two marks	2,08	42,08
3,0 x 60;	Two marks	2,08	62,08
3,0 x 80;	Two marks	2,08	82,08
3,0 x 120;	Two marks	2,08	122,08
3,0 x 150;	Two marks	2,08	152,08
3,0 x 200;	Two marks	2,08	202,08
4,0 x 20;	Two marks	3,34	23,34
4,0 x 40;	Two marks	3,34	43,34
4,0 x 60;	Two marks	3,34	63,34
4,0 x 80;	Two marks	3,34	83,34
4,0 x 120;	Two marks	3,34	123,34
4,0 x 150;	Two marks	3,34	153,34
4,0 x 200;	Two marks	3,34	203,34
5,0 x 20;	Two marks	4,83	24,83
5,0 x 40;	Two marks	4,83	44,83
5,0 x 60;	Two marks	4,83	64,83
5,0 x 80;	Two marks	4,83	84,83
5,0 x 120;	Two marks	4,83	124,83
5,0 x 150;	Two marks	4,83	154,83
5,0 x 200;	Two marks	4,83	204,83
6,0 x 20;	Two marks	6,15	26,15
6,0 x 40;	Two marks	6,15	46,15
6,0 x 60;	Two marks	6,15	66,15
6,0 x 80;	Two marks	6,15	86,15
6,0 x 120;	Two marks	6,15	126,15
6,0 x 150;	Two marks	6,15	156,15
6,0 x 200;	Two marks	6,15	206,15
7,0 x 20;	Two marks	7,58	27,58
7,0 x 40;	Two marks	7,58	47,58
7,0 x 60;	Two marks	7,58	67,58
7,0 x 80;	Two marks	7,58	87,58
7,0 x 120;	Two marks	7,58	127,58
8,0 x 20;	Two marks	8,85	28,85
8,0 x 40;	Two marks	8,85	48,85
8,0 x 60;	Two marks	8,85	68,85
8,0 x 80;	Two marks	8,85	88,85
9,0 x 20;	Two marks	10,21	30,21
9,0 x 40;	Two marks	10,21	50,21
9,0 x 60;	Two marks	10,21	70,21
9,0 x 80;	Two marks	10,21	90,21
10,0 x 20;	Two marks	11,55	31,55
10,0 x 40;	Two marks	11,55	51,55
10,0 x 60;	Two marks	11,55	71,55
12,0 x 20;	Two marks	14,19	34,19
12,0 x 40;	Two marks	14,19	54,19
12,0 x 60;	Two marks	14,19	74,19

Table 5. Data on the balloon pressure taking into account the length and diameter of the balloon.

Balloon length (mm)	Balloon diameter (mm)	Rated pressure (Atm) [kPa]	Calculated burst pressure (Atm) [kPa]
20 ÷ 60	3,0	6 [607,95]; 7 [709,275]	20 [2026,5]
80 ÷ 200 20 ÷ 200 20 ÷ 200 20 ÷ 200 20 ÷ 120 20 ÷ 80 20 ÷ 80 20 ÷ 60 20 ÷ 60	3,0 4,0 5,0 6,0 7,0 8,0 9,0 10,0 12,0	6 [607,95]; 7 [709,275]	13[1317,22] ÷ 16[1621,2]

Note: nominal pressure corresponds to the pressure at which the balloon reaches the diameter indicated in the labeling; calculated burst pressure corresponds to the maximum pressure to which it is possible to inflate the balloon.

Table 6. Nominal pressure and calculated burst pressure of the balloon (depending on the diameter).

Balloon diameter (mm)	3,0 (L≤60mm)	3,0 (L>60mm)	4,0	5,0	6,0	7,0	8,0	9,0	10,0	12,0
Nominal pressure (Atm)	7	7	7	7	6	7	7	6	6	6
Calculated burst pressure (Atm)	20	16	16	16	16	16	15	14	13	13

Secondary parameters:

Compatible Introducer (mm) [Fr]: for diameters 3.0 ÷ 6.0 = 1.64 [5]; for other diameters: 2.0[6] ÷ 2.3 [7]

Permissible conductor diameter (mm) [Inch]: 0,89 [0,035];

Maximum deflation time of the balloon for diameters 3,0 ÷ 9,0: no more than 8s; for others: no more than 10s.

Load on the gap in the connecting points of the catheter (N): not less than: 10N.

Mass of catheter (g): no more than 6,00.

Table 6. Medical device manufacturing materials

Part name of the device	Material of manufacture and manufacturer
1) Luer lock port	Polycarbonate (color: clear) Paint on the port surface: (laser application)
2) Connecting tube/strain relief	Pebax, color: white;
3) Catheter shaft: Front end Internal lumen and lumen for guidewire insertion Tip	Nylon 12; color: blue; Nylon 12/Plexar/HDPE; color: clear; Pebax 7033; color: blue;

4) Balloon	Nylon 12; color: clear;
5) Protective balloon tube	Nylon 12; color: clear; Plexar/HDPE
6) Balloon radiopaque marker	Alloy platinum and iridium Tungsten + Pebax 63333
7) Catheter coating:	Dispersion based on pharmaceutical grade silicone for medical use
PACKAGING:	
8) Catheter protective case	High-density polyethylene HDPE
9) Catheter protective case clip	High-density polyethylene HDPE
10) Blister	LLDPE/Tyvek
11) Consumer packaging	Coated cardboard

13. HOW TO WORK WITH A MEDICAL DEVICE.

Attention! Please refer to the instructions for use before use!

Warning. Do not use the device after the expiration date. If there is obvious damage to the device, it is forbidden to use it!

13.1. Preparing the catheter

Remove the device from its packaging. Check that the balloon procedure and the selected additional equipment the catheter is indicated on the label attached to the protective case and the stickers placed on individual and consumer packages.

Gently remove the protective sheath and stylet from the balloon.

Connect a 10ml syringe containing sterile saline solution to the port for the guidewire (straight part of the Y-connector) and irrigate the lumen until the liquid comes out of the catheter's tip.

Purge the whole air system:

- attach a three-way stopcock to the contrast liquid port (the side part of the connector). An extension line could be optionally placed between the balloon connector and the three-way stopcock;
- close the air passage through the balloon;
- attach a 10-20ml syringe with a one-third part of saline solution to the three-way stopcock;
- open the passage in the three-way stopcock between the syringe and the catheter;
- with the syringe in a vertical position, withdraw the plunger upwards allowing the air bubbles to exit to the liquid;
- when bubbles stop coming in, close the three-way stopcock at the catheter end and remove the syringe;

ATTENTION: If bubbles do not stop coming from the catheter into the syringe after 3 minutes of negative pressure, this may be a clear indication that the balloon catheter has leaks, is broken or the connections between the syringe and the three-way stopcock are not properly closed. If bubbles continue to be seen after checking the connections, do not use the device. Return it to the manufacturer or distributor for inspection.

13.2 Technique for Insertion / Treatment

- Position the 0.035" guidewire across the lesion using PTA techniques, using fluoroscopy to determine its position at all times.
- Insert the proximal end of the guidewire into the distal tip of the balloon catheter.

- Carefully advance the balloon catheter on the guidewire, through the introducer and with the aid of fluoroscopy until it reaches the area of stenosis. In the event that the balloon catheter selected cannot be placed inside the lesion, take a smaller size to predilate the lesion.
- Connect the inflation device (1/3 of its capacity of contrast liquid mixture: saline solution) to the three-way stopcock and open the passage between the balloon and the inflation device.
- Apply pressure to the inflation device so that the balloon inflates. Do not exceed the calculated burst pressure on the label and the compliance curve.
- Apply pressure for 30-60 seconds depending on the lesion.
- Draw the plunger on the inflation device back to deflate the balloon. Keep up the negative pressure for between 20 and 60 seconds depending on the size of the balloon. Make sure the balloon is fully deflated (by means of fluoroscopy) before moving the catheter.
- With negative pressure in the inflation device and with the guidewire in position, carefully retract the catheter until through the introducer sheath.
- Inject contrast liquid through the introducer to check the treated lesion before removing the guidewire.

14. INFORMATION ABOUT ADDITIONAL EQUIPMENT REQUIRED BUT NOT SUPPLIED.

- Heparinized saline solution;
- Contrast agents (use contrast agents intended for intravascular use);
- One guidewire with a diameter of 0.035", do not use conductors with other dimensions;
- One Introducer with a hemostatic valve, the size of which is indicated on the label. Do not use smaller introducer that might damage the catheter;
- One three-way tap;
- One device for inflating the balloon (hand pump with built-in pressure gauge).

15. INFORMATION ABOUT MAINTENANCE AND REPAIR REQUIREMENTS.

This device is not applicable because this medical device is intended for single use.

16. OPERATING CONDITIONS

The device is resistant to climatic factors during operation for temperature conditions from +15 to +35 °C, at a relative humidity of 45-80%, at atmospheric pressure from 84.0 to 106.7 kPa.

17. TRANSPORTATION AND STORAGE CONDITIONS

17.1. Transportation.

PTA balloon dilation catheter iVascular OCEANUS 35 is transported by all modes of transport in covered vehicles in accordance with the rules of cargo transportation applicable to each type of transport. Transport temperature from 0°C to +40 °C, at a relative humidity no more than 75%.

17.2. Storage.

The device must be stored in the warehouses of the supplier and consumer in a cool dry place (0°-40°C), away from direct sunlight, at a relative humidity no more than 75%.

18. INFORMATION ABOUT STERILIZATION.

The device is sterilized with ethylene oxide.

19. THE EXPIRATION DATE.

Expiration date of a medical device:

Storage conditions are met, the manufacturer sets the following expiration date of the medical device in its original packaging: 3 years from the date of production.

20. THE MANUFACTURER'S WARRANTY.

This product and its components have been designed, manufactured, tested, and packed with all appropriate precautions. LIFE VASCULAR DEVICES BIOTECH, S. L. guarantees the product until its expiration date, provided that its packaging is not damaged or torn.

THE MANUFACTURER GUARANTEES: sterility of the product, as well as its compliance with the requirements of international and national standards during the entire shelf life of the product, subject to the rules of transportation, storage and operating conditions established by the manufacturer and when used in accordance with the instructions for use.

21. REQUIREMENTS FOR ENVIRONMENTAL PROTECTION.

The materials used to make this product do not have a harmful effect on the environment, provided that it is disposed of in accordance with local biosecurity standards.

22. DISPOSAL OF MEDICAL DEVICES.

After use, the product presents a potential biological hazard. The product must be handled and disposed of in accordance with medical practice, as well as the relevant local, state, and Federal laws and regulations.

Hazard class of medical waste-B. (Sanitary regulations and standards 2.1.7.2790-10)

23. CLAIMS AND COMPLAINTS.

All claims and complaints related to a medical device should be sent to the authorized representative of the manufacturer. The address is listed below. In the letter of claim, you must specify your full name, first name and patronymic, as well as your contact phone number, the name of a legal entity or individual entrepreneur.

24. PACKAGING AND LABELING.

24.1. Packaging.

The medical device is delivered one piece per package. Each primary package contains a balloon catheter for dilation.

The PTA balloon dilation catheter iVascular OCEANUS 35 is placed inside the protective tube, which supports its folded state (the balloon dilation catheter iVascular OCEANUS 35 also has a reinforcing rod of the conductor lumen); the device is then placed in a special protective case in the form of a hollow HDPE tube, fixed to special clips to avoid possible damage.

Inside the primary package (blister) is also a card with the extensibility curve, which is a description of the characteristic elasticity function in the form of values of the cylinder diameter, achieved with each step-by-step increase in the created pressure. For product identification purposes, a sticker with a label is placed on the primary packaging (blister). Sterility indicators are placed on the reverse side of the blister.

The primary packaging (the blister that encloses this device, which is necessary to preserve the sterility of the device) is placed in an individual consumer cardboard package (together with a booklet with instructions for use) to avoid damage during transport and storage.



Further, in the amount of 42 units, the medical device in individual consumer cardboard packages is placed in a transport box made of corrugated cardboard, which is sealed with transparent adhesive tape and sent to the warehouse for shipment.

24.2. Labeling.

The following labels and symbols are applied to the port of a medical device with a Luer type lock:

	- brand logo and manufacturer's name	0.00 x 00mm	The size of the balloon catheter format: balloon diameter (mm) x balloon length in (mm)
Vascular		LOT	The digital code of the lot.

The following labels and symbols are applied to the individual blister pack:

SYMBOLS (including ISO 15223-1)			
REF	Device reference number		Lot code
	Do not re-use		Do not use with damaged packaging
	Use by...		Avoid direct sunlight
	Keep away from moisture		Apyrogenic
	Attention! View the attachment.		Manufacturer
	Do not re-sterilize		Sterilization using ethylene oxide
CE	Complies with European Directive 93/42/EEC		Do not dispose of as household waste
 Diagram showing the catheter with dimensions: Catheter length L, and diameter Ø GW 0.035" 0.89 mm.			Schematic image of the catheter

	Manufacturing date		Storage temperature limit
	Read the instructions for use		
TITLES			
Name of the medical device	Catheter size	Compatible introducer size	Compatible guide-wire parameters
Name of the manufacturer			
Information about the number of radioopaque markers of the balloon	Hydrophilic coating	Number of devices per package	OTW describes the method of delivery of the catheter through the guide-wire.
Barcode or QR code	Information about the manufacturer's authorized representative in the Russian Federation	Registration certificate number and date	
THE FOLLOWING INFORMATION IS GIVEN IN TABLE FORM			
The extensibility curve, in the form of values of the balloon diameter, achieved with each step-by-step increase in the created pressure.			

The transport package of a medical device contains the following labels and symbols:

SYMBOLS (including ISO 15223-1)			
REF	Device reference number		Lot code
	Use by...		Avoid direct sunlight
	Keep away from moisture		Manufacturer
CE	Complies with European Directive 93/42/EEC		Temperature limit during transportation
	Do not dispose of as household waste		Careful handling
	Top		

<i>TITLES</i>
Number of devices per package.
Information about the manufacturer's authorized representative in the Russian Federation

25. INFORMATION ABOUT THE MANUFACTURER'S AUTHORIZED REPRESENTATIVE ON THE TERRITORY OF THE RUSSIAN FEDERATION.

RAYMED TG PVT LTD

8, 1st Vokzalnaya Str., Mytishchi, Mytishchi district, Moscow region, 141001

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Date of entry into force:

Version: 28-/July /2020

Action	FULL NAME	Position	Date	Signature
Prepared by	Aurea Roca	Regulatory Affairs Director	28/07/2020	
Approved by	M. Eugenia Villanueva	Regulatory Com- pliance Manager	28/07/2020	

LS

2020

1. PURPOSE OF THE DEVICE.

PTA balloon dilation catheter iVascular OCEANUS 14 PRO is intended for percutaneous transluminal angioplasty of peripheral arteries.

2. RISK CLASS OF THE MEDICAL DEVICE.

The potential risk class for the use of a medical device is 3, in accordance with Directive 93/42/EEC. Type of contact with the human body is a product that briefly (for less than 24 hours) comes into contact with blood.

3. LIST OF INTERNATIONAL REGULATIONS/STANDARDS THAT THE MEDICAL DEVICE MEETS.

STANDARD CODE DESIGNATION	STANDARD NAME
ISO 13485	Medical devices. Quality management system. Requirements for regulatory purposes.
EN ISO 14971	Medical devices. The application of risk management to medical devices
UNE-EN 1041+A1	Information prepared by the manufacturer accompanying medical devices
EN 556-1	Sterilization of medical devices. Requirements for medical devices of the "sterile" category. Part 1. Requirements for medical devices subject to final sterilization.
UNE-EN ISO 11135	Sterilization of medical devices. Ethylene oxide. Part 1: Requirements for the development, validation, and ongoing monitoring of the sterilization process for medical devices.
UNE-EN ISO 11737-2	Sterilization of medical devices. Microbiological method. Part 2: Sterility tests performed at the stages of technical development, validation and continuous maintenance of the sterilization process.
UNE-EN ISO 14644-1	Clean rooms and their associated controlled areas. Part 1: Classification of air cleanliness.
UNE-EN ISO 14644-2	Clean rooms and their associated controlled areas. Part 2: Specifications for testing and monitoring to ensure continuous compliance with ISO 14644-1.
UNE-EN ISO 14698-2	Clean rooms and their associated controlled areas. Control of bio-contamination. Part 2: Analysis and interpretation of biosafety data.
UNE-EN ISO 14644-3	Clean rooms and their associated controlled areas. Part 3: Methods of control.
UNE-EN ISO 14698-1	Clean rooms and their associated controlled areas. Control of bio-contamination. Part 1: General principles and methods.
MDD 93/42/EEC	Directive 93/42/EEC Medical devices, devices and equipment.
UNE EN ISO 10555-1	Single-use sterile intravascular catheters. Part 1: General requirements.
UNE EN ISO 10555-4	Single-use sterile intravascular catheters. Part 4: Balloon dilation catheters.
UNE EN 20594-1	Conical fittings with a 6% Luer lock cone for syringes, needles and some other medical equipment. Part 1: General requirements.
UNE EN 11607-1	Packaging for medical devices that have undergone final sterilization. Part 1: requirements for materials, sterility maintenance systems, and packaging systems.
UNE-EN 62366	Medical devices: Application of serviceability design to medical devices.

UNE EN ISO 10993-1	Assessment of the biological effect of medical devices. Part 1. Analysis and testing in the process of implementing risk management.
UNE EN ISO 10993-4	Assessment of the biological effect of medical devices. Part 4. Selection of tests for interaction with blood.
UNE EN ISO 10993-5	Assessment of the biological effect of medical devices. Part 5. In vitro cytotoxicity tests.
UNE EN ISO 10993-7	Assessment of the biological effect of medical devices. Part 7. Residual amounts of ethylene oxide after the sterilization process.
UNE EN ISO 10993-9	Assessment of the biological effect of medical devices. Part 9. Standards for authenticating and quantifying potential degradation products.
UNE EN ISO 10993-10	Assessment of the biological effect of medical devices. Part 10. Skin irritation and sensitization tests.
UNE EN ISO 10993-11	Assessment of the biological effect of medical devices. Part 11. Tests for general toxic effects.
UNE EN ISO 10993-12	Assessment of the biological effect of medical devices. Part 12. Preparation of test samples and standard samples.
UNE EN ISO 10993-13	Assessment of the biological effect of medical devices. Part 13. Authenticity and quantification of degradation products when using polymer medical devices.
UNE EN ISO 10993-18	Assessment of the biological effect of medical devices. Part 18. Chemical characteristics of materials.
UNE EN ISO 11737-1	Sterilization of medical devices. Microbiological method. Part 1. Determination of microbial populations on products.

4. INFORMATION ON THE MANUFACTURER OF THE MEDICAL DEVICE

Manufacturing company:

LIFE VASCULAR DEVICES BIOTECH, S.L.
(LIFE VASCULAR DIVISION BIOTECH, Spain)
Camí de Can Ubach, 11; Pol. Ind. Les Fallulles
08620 Sant Vicenç Dels Horts- Barcelona, Spain.

Manufacturing organizations:

LIFE VASCULAR DEVICES BIOTECH, S.L.
(LIFE VASCULAR DIVISION BIOTECH, Spain)
Camí de Can Ubach, 11; Pol. Ind. Les Fallulles
08620 Sant Vicenç Dels Horts- Barcelona, Spain.
Tel/fax +34936724711,
E-mail: info@ivascular.es

5. DELIVERY SET.

- 1) Balloon catheter - 1 pc.
- 2) Instructions for use - 1 pc.
- 3) Tracking label - 6 pcs

6. INFORMATION ABOUT POTENTIAL CONSUMERS (USERS) OF A MEDICAL DEVICE

Qualified personnel of vascular surgery departments of medical treatment and preventive care institutions.

7. INFORMATION ABOUT THE VARIANTS OF THE MEDICAL DEVICE.

I.PTA balloon dilatation catheter iVascular OCEANUS 14 PRO:

1. PTA balloon dilatation catheter iVascular OCEANUS 14 PRO, with guidewire of 0,014 inch, hydrophilic coating, catheter length of 100 cm, balloon dimensions: 1,25 x 10mm; 1,25 x 15mm; 1,25 x 20mm; 1,5 x 10mm; 1,5 x

15mm; 1,5 x 20mm; 1,5 x 40mm; 1,5 x 60mm; 1,5 x 80mm; 2,0 x 10mm; 2,0 x 15mm; 2,0 x 20mm; 2,0 x 40mm; 2,0 x 60mm; 2,0 x 80mm; 2,0 x 120mm; 2,0 x 150mm; 2,0 x 200mm; 2,5 x 40mm; 2,5 x 60mm; 2,5 x 80mm; 2,5 x 120mm; 2,5 x 150mm; 2,5 x 200mm; 3,0 x 40mm; 3,0 x 60mm; 3,0 x 80mm; 3,0 x 120mm; 3,0 x 150mm; 3,0 x 200mm; 3,5 x 40mm; 3,5 x 60mm; 3,5 x 80mm; 3,5 x 120mm; 3,5 x 150mm; 3,5 x 200mm; 4,0 x 40mm; 4,0 x 60mm; 4,0 x 80mm; 4,0 x 120mm.

2. PTA balloon dilatation catheter iVascular OCEANUS 14 PRO, with guidewire of 0,014 inch, with hydrophilic coating, catheter length of 150 cm, balloon dimensions: 1,25 x 10mm; 1,25 x 15mm; 1,25 x 20mm; 1,5 x 10mm; 1,5 x 15mm; 1,5 x 20mm; 1,5 x 40mm; 1,5 x 60mm; 1,5 x 80mm; 2,0 x 10mm; 2,0 x 15mm; 2,0 x 20mm; 2,0 x 40mm; 2,0 x 60mm; 2,0 x 80mm; 2,0 x 120mm; 2,0 x 150mm; 2,0 x 200mm; 2,5 x 40mm; 2,5 x 60mm; 2,5 x 80mm; 2,5 x 120mm; 2,5 x 150mm; 2,5 x 200mm; 3,0 x 40mm; 3,0 x 60mm; 3,0 x 80mm; 3,0 x 120mm; 3,0 x 150mm; 3,0 x 200mm; 3,5 x 40mm; 3,5 x 60mm; 3,5 x 80mm; 3,5 x 120mm; 3,5 x 150mm; 3,5 x 200mm; 4,0 x 40mm; 4,0 x 60mm; 4,0 x 80mm; 4,0 x 120mm

8. DESCRIPTION OF THE MEDICAL DEVICE AND ITS PHOTOGRAPHIC IMAGES.

PTA balloon dilation catheter iVascular OCEANUS 14 PRO is used for dilation and is used to inflate the balloon for stenosis or occlusion of peripheral arteries in order to recanalize the artery.

This procedure is performed using minimally invasive techniques, also called percutaneous transluminal angiography (PTA) of the peripheral arteries in which the puncture of the outer skin allows access to the blood flow and to the affected area.

The catheter has a coaxial design and is a dual-lumen catheter (from the luer connector to the tip), also called an OTW catheter, i.e. a catheter with a stent delivery system along the guidewire.


Fig. 1. Schematic representation of a catheter.

The tip is at the distal end of the catheter, it is atraumatic in shape to avoid damage to the arteries during its insertion. In the distal part of the catheter, in front of its tip, there is a balloon, or an inflatable segment that allows you to pump the pressure to various values. The use of a balloon leads to the expansion of the artery by pumping a contrast medium into it.

One or two radiopaque markers are placed: one marker - in the center of the balloon, two markers: one on each of the tapered tips of the balloon, in order to labeling its length and help the user determine the location of the catheter when it is inside the patient (when displayed on the x-ray screen).

Between the Luer lock connector and the catheter body is a connecting tube (strain relief), which is necessary to prevent the catheter from bending at this point.

The connector has a branched (Y-shaped) shape and is provided with two ports directly connected to the catheter shaft:

- the side port (inflating port) allows for transmission contrast fluid to the catheter's inflation and deflating lumen (the lumen for inflating);
- the direct port (guidewire insertion port) is intended for insertion of the catheter guidewire. A guidewire is inserted into the straight port of the connector and exits the catheter tip in order to reach the affected area.

For OCEANUS 14 PRO, the maximum guidewire diameter must not exceed 0.36 mm (0.014 in).


Fig. 2. Enlarged image of the proximal part of the PTA balloon dilation catheter iVascular OCEANUS 14 PRO.


Fig. 3. Enlarged image of the middle part of the PTA balloon dilation catheter iVascular OCEANUS 14 PRO.


Fig. 4. Enlarged image of the part at the junction of the tube with the balloon of the PTA balloon dilation catheter iVascular OCEANUS 14 PRO.


Fig. 5. Enlarged image of the distal part of the PTA balloon dilation catheter iVascular OCEANUS 14 PRO


Fig. 6. Image of a general view of the medical device: PTA balloon dilation catheter iVascular OCEANUS 14 PRO in a protective case.

9. TERMS OF USE.

Only qualified specialists of the surgical departments of the treatment and prevention institutions.

10. SCOPE OF APPLICATION.

Vascular surgery.

11. INDICATIONS FOR USE. CONTRAINDICATIONS. PRECAUTIONARY MEASURES.

INDICATIONS.

PTA balloon dilation catheter iVascular OCEANUS 14 PRO, with hydrophilic coating is used for dilation of stenoses in the femoral, popliteal arteries, in the end branches of the popliteal arteries, in the renal arteries, as well as for the treatment of obstructive lesions of native or synthetic arteriovenous dialysis fistulae, whose diameters from 1.25 to 4.0 mm, and lengths from 10 to 200 mm.

CONTRAINDICATIONS.

Inability to pass the conductor through the stenosed area.

POSSIBLE SIDE EFFECTS/COMPLICATIONS.

Possible side effects and/or complications that may occur before, during, or after the catheterization procedure are:

- Death
- Stroke/embolism/thrombosis
- Haemodynamic deterioration
- Sharp elastic rebound
- Arrhythmias
- Aneurysm or pseudo-aneurysm
- Infections
- Complete arterial occlusion
- Reocclusion of the treated area: restenosis
- Spasm
- Perforation or dissection of the treated area
- Arteriovenous fistula
- Local hemorrhage with hematoma at the access point
- Hypo/hypertension
- Allergic reactions to materials

CAUTIONARY MEASURES.

Perform appropriate treatment of the patient using anticoagulants, vasodilators, etc., according to the procedure for inserting intravascular catheters.

- Before use, check the size and suitability of the catheter for its intended use.
- Proceed with extreme care so as not to damage the catheter while advancing it.
- The balloon must be inflated with a mixture of saline solution and contrast medium (preferably 1/1, though the proportion of contrast medium can be reduced in large balloons).
- Do not exceed the calculated burst pressure as the balloon could burst.
- If any resistance is noticed during advance, stop and determine its cause before continuing.
- Use with standard 0.014" guidewires and standard introducers. If the introducer length goes above 25cm or the braided type is used, it may be necessary to increase its size to avoid catheter friction.
- If any resistance to withdrawing the catheter is noted, it is recommended that the balloon catheter, the guidewire and introducer be removed as a single unit.
- After use, this product may represent a biohazard. Handle and dispose it in accordance with accepted medical practices and pertinent local, state or federal laws and regulations.
- Store in a cool, dry place away from direct sunlight.

WARNINGS.

- The device should only be used by experienced doctors familiar in the percutaneous transluminal angioplasty (PTA) technique.
- The percutaneous procedure must be carried out in hospitals equipped with emergency facilities to conduct open heart surgery or hospitals with easy access to other hospitals where this type of operation can be performed.
- The product is for single use. It must not be re-sterilised or re-used once the procedure has been completed. Re-using the product in another patient may lead to cross contamination, infection or transmission of infectious diseases from one patient to another. Re-use of the product may cause alterations to it and affect its effectiveness.
- The product is supplied sterile. Check the expiry date and do not use products which are past this date.
- Inspect the packaging before opening it. If there are any defects or the packaging is damaged, do not use the product.
- Use aseptic techniques when the product is taken out of the primary packaging.
- Select the balloon diameter according to the diameter of the artery to be treated. Do not insert a balloon with a diameter larger than the artery to be treated.
- Do not dry with gauze.
- Do not expose the product to organic solvents or to contrast media not suitable for intravascular use.
- Do not use air or gaseous media to inflate the balloon.
- Advance the product on the guidewire using fluoroscopy. Do not allow the product to advance without the guidewire inside it.
- Do not manipulate, advance or retract either the catheter or the wire guide while the balloon is inflated.

12. TECHNICAL CHARACTERISTICS AND MEDICAL DEVICE MANUFACTURING MATERIALS.

Table 1. The principal dimensions of the catheter.

Effective catheter length (mm)	Balloon diameter (mm)	Effective balloon-length (mm)	Reinforcing bar length (mm)
1000 1500	1,25	10	283
	1,5	15	
	2,0	20	
	2,5	40	
	3,0	60	
	3,5	80	
	4,0	120	
		150	
		200	

Table 2. Data on external and internal catheter diameters, depending on the parameters of the balloon.

The cylinders of the catheter (di- ameter x length mm)	External proximal catheter diameter (mm) [Fr]	External distal catheter diame- ter (mm) [Fr]	Internal catheter diameter (mm) [Inches]
1,25 x 10; 1,25 x 15; 1,25 x 20; 1,5 x 10; 1,5 x 15; 1,5 x 20; 2,0 x 10; 2,0 x 15; 2,0 x 20	1,16 [3,5]	0,75 [2,3]	0,43 [0,017]
Other balloon sizes	1,31 [3,9]	0,98 [2,9]	

Table 3. Data of the balloon cross-profile depending on the diameter of the balloon.

Balloon di- ameter (mm)	Inches	mm	Fr
1.25	0,021	0,53	1,6F
1.5	0,023	0,60	1,8F
2.0	0,025	0,65	1,9F
2.5	0,027	0,70	2,1F
3.0	0,029	0,75	2,2F
3.5	0,031	0,80	2,4F
4.0	0,033	0,85	2,5F

Table 4. Data on the location of radiopaque markers/marks of catheter balloons.

Catheter balloon (diameter x length in mm)	Number of radiopaque markers/marks	Distance from the distal end of the balloon to the distal mark (mm)	Distance from the bal- loon distal end to the proximal mark (mm)
1,25 x 10; 1,25 x 15; 1,25 x 20; 1,5 x 10; 1,5 x 15; 1,5 x 20; 2,0 x 10; 2,0 x 15; 2,0 x 20	One mark in the center	One mark in the center	One mark in the center
1,5 x 40	Two marks	1,42	41,42
1,5 x 60	Two marks	1,42	61,42
1,5 x 80	Two marks	1,42	81,42
2,0 x 40	Two marks	2,21	42,21
2,0 x 60	Two marks	2,21	62,21
2,0 x 80	Two marks	2,21	82,21
2,0 x 120	Two marks	2,21	122,21
2,0 x 150	Two marks	2,21	152,21
2,0 x 200	Two marks	2,21	202,21
2,5 x 40	Two marks	2,95	42,95
2,5 x 60	Two marks	2,95	62,95
2,5 x 80	Two marks	2,95	82,95
2,5 x 120	Two marks	2,95	122,95
2,5 x 150	Two marks	2,95	152,95
2,5 x 200	Two marks	2,95	202,95
3,0 x 40	Two marks	3,73	43,73
3,0 x 60	Two marks	3,73	63,73
3,0 x 80	Two marks	3,73	83,73
3,0 x 120	Two marks	3,73	123,73
3,0 x 150	Two marks	3,73	153,73
3,0 x 200	Two marks	3,73	203,73
3,5 x 40	Two marks	4,52	44,52
3,5 x 60	Two marks	4,52	64,52

3,5 x 80	Two marks	4,52	84,52
3,5 x 120	Two marks	4,52	124,52
3,5 x 150	Two marks	4,52	154,52
3,5 x 200	Two marks	4,52	204,52
4,0 x 40	Two marks	5,28	45,28
4,0 x 60	Two marks	5,28	65,28
4,0 x 80	Two marks	5,28	85,28
4,0 x 120	Two marks	5,28	125,28

Table 5. Data on the balloon pressure taking into account the length and diameter of the balloon.

Balloon length (mm)	Balloon diameter (mm)	Rated pressure (Atm) [kPa]	Calculated burst pressure (Atm) [kPa]
10 ÷ 200	1,25 ÷ 4,0	7 [709,275]	16 [1621,2]

Note: nominal pressure corresponds to the pressure at which the balloon reaches the diameter indicated in the labeling; calculated burst pressure corresponds to the maximum pressure to which it is possible to inflate the balloon.

Secondary parameters:

Compatible Introducer (mm) [Fr]: 1,3 [4];

Permissible guidewire diameter (mm) [Inch]: 0,36 [0,014];

Maximum deflation time of the balloon for all diameters and lengths: no more than 10s.

Load on the gap in the connecting points of the catheter (N): not less than: 10N.

Mass of catheter (g): no more than 6,00.

Table 6. Medical device manufacturing materials

Part name of the device	Material of manufacture and manufacturer
1) Luer lock port	Polycarbonate (color: clear) Mark on the port surface: (laser application)
2) Connecting tube/strain relief	Pebax, color: white;
3) Catheter shaft:	
Front end	Nylon 12; color: blue;
Internal lumen and lumen for guidewire insertion	Nylon 12/Plexar/HDPE; color: clear;
Tip	Nylon12; color: blue;
4) Balloon	Nylon 12; color: clear;
5) Protective balloon tube	Nylon 12; color: clear; Plexar/HDPE
6) Balloon radiopaque marker	Alloy platinum and iridium Tungsten + Pebax 63333
7) Reinforcing shaft	Stainless steel
8) Catheter coating: the hydrophilic coating of the catheter	HYDRAX PLUS
PACKAGING:	
9) Catheter protective case	High-density polyethylene HDPE

10) Catheter protective case clip	High-density polyethylene HDPE
11) Blister	LLDPE/Tyvek
12) Consumer packaging	Coated cardboard

13. HOW TO WORK WITH A MEDICAL DEVICE.

Attention! Please refer to the instructions for use before use!

Warning. Do not use the device after the expiration date. If there is obvious damage to the device, it is forbidden to use it!

13.1. Preparing the catheter

Remove the device from its packaging. Check that the balloon procedure and the selected additional equipment the catheter is indicated on the label attached to the protective case and the stickers placed on individual and consumer packages.

Gently remove the protective sheath from the balloon and the balloon's reinforcing shaft.

Connect a 10ml syringe containing sterile saline solution to the port for the guidewire (straight part of the Y-connector) and irrigate the lumen until the liquid comes out of the catheter's tip.

Pump the air system:

- attach a three-way stopcock to the contrast liquid port (the side part of the connector). An extension line could be optionally placed between the balloon connector and the three-way stopcock;
- close the air passage through the balloon;
- attach a 10-20ml syringe with a one-third part of saline solution to the three-way stopcock;
- open the passage in the three-way stopcock between the syringe and the catheter;
- with the syringe in a vertical position, withdraw the plunger upwards allowing the air bubbles to exit to the liquid;
- when bubbles stop coming in, close the three-way stopcock at the catheter end and remove the syringe.

ATTENTION: if bubbles do not stop coming from the catheter into the syringe after 3 minutes of negative pressure, this may be a clear indication that the balloon catheter has leaks, is broken or the connections between the syringe and the three-way stopcock are not properly closed. If bubbles continue to be seen after checking the connections, do not use the device. Return it to the manufacturer or distributor for inspection.

In order to activate the hydrophilic coating, it is recommended to wet the catheter with a sterile saline solution immediately prior to its insertion through the introducer.

13.2 Technique for Insertion / Treatment

- Position the 0.014" guidewire across the lesion using PTA techniques, using fluoroscopy to determine its position at all times.
- Insert the proximal end of the guidewire into the distal tip of the balloon catheter.
- Carefully advance the balloon catheter on the guidewire, through the introducer and with the aid of fluoroscopy, until it reaches the area of stenosis. In the event that the balloon catheter selected cannot be placed inside the lesion, take a smaller size to predilate the lesion.
- Connect the inflation device (1/3 of its capacity of contrast liquid mixture: saline solution) to the three-way stopcock and open the passage between the balloon and the inflation device.
- Apply pressure to the inflation device so that the balloon inflates. Do not exceed the rated burst pressure on the label and the compliance curve.
- Apply pressure for 30-60 seconds depending on the lesion.

Draw the plunger on the inflation device back to deflate the balloon. Keep up the negative pressure for between 20 and 60 seconds depending on the size of the balloon. Make sure the balloon is fully deflated (by means of fluoroscopy) before moving the catheter.

- With negative pressure in the inflation device and with the guidewire in position, carefully retract the catheter until through the introducer sheath.
- Inject contrast liquid through the introducer to check the treated lesion before removing the guidewire.

14. INFORMATION ABOUT ADDITIONAL EQUIPMENT REQUIRED BUT NOT SUPPLIED.

- Heparinized saline solution;
- Contrast agents (use contrast agents intended for intravascular use);
- One guidewire with a diameter of 0.014", do not use conductors with other dimensions;
- One Introducer with a hemostatic valve, the size of which is indicated on the label. Do not use smaller introducer that might damage the catheter;
- One three-way tap;
- One device for inflating the balloon (hand pump with built-in pressure gauge).

15. INFORMATION ABOUT MAINTENANCE AND REPAIR REQUIREMENTS.

This device is not applicable because this medical device is intended for single use.

16. OPERATING CONDITIONS

The device is resistant to climatic factors during operation for temperature conditions from +15 to +35 °C, at a relative humidity of 45-80%, at atmospheric pressure from 84.0 to 106.7 kPa.

17. TRANSPORTATION AND STORAGE CONDITIONS

17.1. Transportation.

PTA balloon dilation catheter iVascular OCEANUS 14 PRO is transported by all modes of transport in covered vehicles in accordance with the rules of cargo transportation applicable to each type of transport. Transport temperature from 0°C to +40 °C, at a relative humidity no more than 75%.

17.2. Storage.

The device must be stored in the warehouses of the supplier and consumer in a cool dry place (0°-40°C), away from direct sunlight, at a relative humidity no more than 75%.

18. INFORMATION ABOUT STERILIZATION.

The device is sterilized with ethylene oxide.

19. THE EXPIRATION DATE.

Expiration date of a medical device:

If the storage conditions are met, the manufacturer sets the following expiration date of the medical device in its original packaging: 3 years from the date of production.

20. THE MANUFACTURER'S WARRANTY.

This product and its components have been designed, manufactured, tested, and packed with all appropriate precautions. LIFE VASCULAR DEVICES BIOTECH, S. L. guarantees the product until its expiration date, provided that its packaging is not damaged or torn.

THE MANUFACTURER GUARANTEES: sterility of the product, as well as its compliance with the requirements of international and national standards during the entire shelf life of the product, subject to the rules of transportation.

storage and operating conditions established by the manufacturer and when used in accordance with the instruction use.

21. REQUIREMENTS FOR ENVIRONMENTAL PROTECTION.

The materials used to make this product do not have a harmful effect on the environment, provided that it is disposed of in accordance with local biosecurity standards.

22. DISPOSAL OF MEDICAL DEVICES.

After use, the product presents a potential biological hazard. The product must be handled and disposed of in accordance with medical practice, as well as the relevant local, state, and Federal laws and regulations.

Hazard class of medical waste-B. (Sanitary regulations and standards 2.1.7.2790-10)

23. CLAIMS AND COMPLAINTS.

All claims and complaints related to a medical device should be sent to the authorized representative of the manufacturer. The address is listed below. In the letter of claim, you must specify your full name, first name and patronymic, as well as your contact phone number, the name of a legal entity or individual entrepreneur.

24. PACKAGING AND LABELING.

24.1. Packaging.

The medical device is delivered one piece per package. Each primary package contains a balloon catheter for dilation.

The PTA balloon dilation catheter iVascular OCEANUS 14 PRO is placed inside the protective tube, which supports its folded state (the balloon dilation catheter iVascular OCEANUS 14 PRO also has a reinforcing rod of the conductor lumen); the device is then placed in a special protective case in the form of a hollow HDPE tube, fixed to special clips to avoid possible damage.


Inside the primary package (blister) is also a card with the extensibility curve, which is a description of the characteristic elasticity function in the form of values of the cylinder diameter, achieved with each step-by-step increase in the created pressure. For product identification purposes, a sticker with a label is placed on the primary packaging (blister). Sterility indicators are placed on the reverse side of the blister.

The primary packaging (the blister that encloses this device, which is necessary to preserve the sterility of the device) is placed in an individual consumer cardboard package (together with a booklet with instructions for use) to avoid damage during transport and storage.

Further, in the amount of 42 units, the medical device in individual consumer cardboard packages is placed in a transport box made of corrugated cardboard, which is sealed with transparent adhesive tape and sent to the warehouse for shipment.

24.2. Labeling.

The following labels and symbols are applied to the port of a medical device with a Luer type lock:

	- brand logo and manufacturer's name	0.00 x 00mm	The size of the balloon catheter format: balloon diameter (mm) x balloon length in (mm)
		LOT	The digital code of the lot.

The following labels and symbols are applied to the individual blister pack:


SYMBOLS (including ISO 15223-1)			
REF	Device reference number	LOT	Lot code
	Do not re-use		Do not use with damaged packaging
	Use by...		Avoid direct sunlight
	Keep away from moisture		Apyrogenic
	Attention! View the attachment.		Manufacturer
	Do not re-sterilize	STERILE EO	Sterilization using ethylene oxide
CE	Complies with European Directive 93/42/EEC		Do not dispose of as household waste
			Schematic image of the catheter
	Manufacturing date		Storage temperature limit
	Read the instructions for use		
TITLES			
Name of the medical device	Catheter size	Compatible introducer size	Compatible guide-wire parameters
Name of the manufacturer			

Information about the number of radioopaque markers of the balloon	Hydrophilic coating	Number of devices per package	OTW describes the method of delivery of the catheter through the guide-wire.
Barcode or QR code	Information about the manufacturer's authorized representative in the Russian Federation	Registration certificate number and date	

THE FOLLOWING INFORMATION IS GIVEN IN TABLE FORM

The extensibility curve, in the form of values of the balloon diameter, achieved with each step-by-step increase in the created pressure.

The transport package of a medical device contains the following labels and symbols:

SYMBOLS (including ISO 15223-1)			
REF	Device reference number		Lot code
	Use by...		Avoid direct sunlight
	Keep away from moisture		Manufacturer
CE	Complies with European Directive 93/42/EEC		Temperature limit during transportation
	Do not dispose of as household waste		Careful handling
	Top		
TITLES			
Number of devices per package.			
Information about the manufacturer's authorized representative in the Russian Federation			

25. INFORMATION ABOUT THE MANUFACTURER'S AUTHORIZED REPRESENTATIVE ON THE TERRITORY OF THE RUSSIAN FEDERATION.

RAYMED TG PVT LTD
8, 1st Vokzalnaya Str., Mytishchi, Mytishchi district, Moscow region, 141001
Phone/Fax: (495) 956-83-16
E-mail: raymedtgpvt@mail.ru

03/2020

I, **PEDRO-ANGEL CASADO MARTÍN**, Notary of the Catalanian Association of Notaries, living in the city of Barcelona (SPAIN),
CERTIFY:-----

That the previous signature belongs to **M. M^a EUGENIA VILLANUEVA PÉREZ**, of legal age, of Spanish nationality and with DNI/NIF. (National Identity Card) number 37323062L.-----

M. M^a EUGENIA VILLANUEVA PÉREZ acts in name and on behalf of the company **"LIFE VASCULAR DEVICES BIOTECH, S.L."**, a company validly existing and incorporated under the laws of Spain; with registered address at Sant Vicenç dels Horts, industrial state Les Fullulles, Camí de Can Ubach, number 11, C.P. 08620; with Tax Identification Number ("CIF") B-65405169; and registered in the Commercial Registry of Barcelona Volume 42171, folio 171 and sheet number B-402656.-----

She acts as attorneys by virtue of the power of attorney granted by the notary public of Barcelona, Mr. Pedro Angel Casado

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