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
Zumax Medical Co., Ltd.
Legal address and address of manufacture (plant):
5 Zhiying Street, Suzhou New District, Suzhou
Jiangsu, China
Tel: +86 512 6665 0512 ext.502/504
Fax: +86 512 6690 9655
<http://www.zumaxmedical.com>
E-mail: sales@zumaxmedical.com

To whom It May Concern

To be submitted to the Competent Health Authorities in Russian Federation

I herewith confirm, that the content of the «Binocular medical loupe, with accessories»
Operational documentation is valid and true.

Zumax Medical CO., LTD, 5 Zhiying Street Suzhou New District, Suzhou Jiangsu, China
is original manufacturer of the device mentioned in Operational documentation.


Zumax Medical CO., LTD
CEO

Xiangdong Li

Date:

2021.2.19

Medical Binocular Loupes
SLE, SLF, SLH
User's Guide

SLE

SLF

SLH

The CE mark on this product indicates it has been tested to and conforms with the provisions noted within the 93/42/EEC Medical Device Directive.

Zumax Medical Co., Ltd.

Thank you for purchasing this product. To prevent damage to your product or injury to yourself or to others, read the following safety precaution in their entirety before using this equipment. Keep these safety instructions where all those who use the product will read them.

Symbols:

Operating instructions

Warning

Attention. Read operating manual for cautions and instructions for use.

Warnings and Cautions

Warning: Do not expose your product to liquid or moisture or to humidity.

Warning: Do not touch the surface of lens with hand or hard objects.

Warning: Never clean the surface of lens with any hard and dirty materials.

Warning: Do not look into the sun or into a powerful light source with this loupe.

Warning: Federal law restricts this device to sale or order of a physician.

1. Purpose and features of the product

1.1 Purpose

The binocular loupes and headlights are used in Medicine.

Prescription

Designed to increase the image space and / or illuminate the field of view of the dentist during a patient examination or surgery

1.2 Features:

- Multi-coated, high grade glass lenses produce crisp, clear images.
- The Loupes which employ eyepiece of very large view field diameter. The depth of focus is deep and object can be seen stereoscopically with actual distance feeling.
- The angle of the loupes can be adjusted between 0° ~ 40° with ease and fatigue is minimized even after prolonged use.
- 6mm up and down displacement. Pupil distance can be adjusted individually which could be suited for all kind of faces.
- Adapts to any Head-Band and Headlight.
- Smart design, convenient for mobile diagnosis.

2. Operating instruction

2.1 SLE and SLF operating,

First step, please fasten the strap to the frame of glasses.

Put on the loupe and check that the pupil distance is symmetrically set.

Now look at a well it three-dimensional object. Adjust one eyepiece after the other until you see a round image with no shadows and a focused object in each. Then screw in both set screws. The position of the optics is adjusted over a wide range by means of the hinge mount. The ease of movement of the hinges can be adjusted by means of the hinge screws.

Adjusting headband

Pull the headline button to the end and then put on the glasses to adjust to the most comfortable position. Tighten the head cord button to ensure that the surgical magnifying glass is stable on the head. Figure 1


Figure 1

Adjust the interpupillary distance:

Loosen the Wrench to make the magnifying glass active. The user rotates the surgical magnifying lens barrel along the middle hinge according to his or her own interpupillary distance, and then puts on the surgical magnifying glass and ensures that the two lens barrels are distributed on both sides of the eye according to the interpupillary distance, so that the two eyepieces are Observe the field of view coincidence. Figure 2


Figure 2

Adjust position :

After adjusting the height and pitch angle of the surgical magnifying glass to the appropriate position, the Wrench is fixed and the Wrench is removed, and the user can clearly observe the object in the field of view. Figure 3


Figure 3

2.2 SLH operating

First step, please adjust the knob of headband, wearing the binocular loupes. Then check that the pupil distance is symmetrically set.

Now look at a well it three-dimensional object. Adjust one eyepiece after the other until you see a round image with no shadows and a focused object in each. Then screw in both set screws. The position of the optics is adjusted over a wide range by means of the hinge mount. The ease of movement of the hinges can be adjusted by means of the hinge screws.

Adjusting headband

Adjust the middle headband to the comfortable position in the middle, and adjust the two knobs to fix the headband to ensure that the surgical zoom lens is stable on the head.


Figure 4

Adjust the interpupillary distance:

Loosen the lever and the knob to make the magnifying glass active. The user rotates the surgical magnifying lens barrel along the middle hinge according to his or her own distance, and then puts on the surgical magnifying tube and ensures that the two barrels are distributed on both sides of the eye according to the interpupillary distance, so that two The observation fields of the eyepieces coincide. Figure 5


Figure 5

Adjust the height:

After adjusting the height of the surgical magnifying glass to a suitable height, tighten the knob to fix the height of the surgical magnifying glass. Figure 6


Figure 6

Adjust position :

After adjusting the height and pitch angle of the surgical magnifying glass to the appropriate position, tighten the lever to fix the position of the surgical magnifying glass, and the user can clearly observe the objects in the field of view. Figure 7


Figure 7

2.3 Optional accessory installation method

Secure the magnetic housing to the surgical magnifying glass with screws; Figer 4


Figure 8

Magnetic HL8200 can be directly attracted to the magnetic base;


Figure 9

3. Maintenance of product

3.1 Maintenance of lens

- The instrument not in use should be covered and put into the packing box.
- Do not disassemble the optical lens by unskilled or unqualified technicians.
- If the product does not work properly, contact an Authorized Representative of the manufacturer. Do not try to repair the product yourself.

3.2 Cleaning

- The outer surface of the instrument may be cleaned with wet cloth, the remaining stains can be cleaned off with the mixture of 50% alcohol and 50% distilled water. Do not wipe with any corrosive detergent lest that the surface should be damaged.
- Cleaning the lenses : Mirror-cleaning paper or a drop of liquid solvent (50% alcohol and 50% aether) may be used to clean the lens, blow it carefully. If there is dust stained on the lens, blow them with a blowball or brush them with a dust pen.

Caution : Be careful to prevent the solvent into the instrument.

4. Specifications

Optical specifications, SLE		
Magnification	Work distance (mm)	Field of view (mm)
Galilean Loupes		
2.5 x	340 / 420 / 460 / 500	80 / 90 / 100 / 110
3.0 x	340 / 420 / 460 / 500	60 / 65 / 68 / 70
3.5 x	340 / 420 / 460 / 500	50 / 60 / 70 / 80
Prismatic Loupes		
4.0 x	340 / 420 / 460 / 500	55 / 65 / 75 / 85
5.0 x	340 / 420 / 460 / 500	45 / 55 / 60 / 65
6.0 x	340 / 420 / 460 / 500	40 / 47 / 50 / 58




Optical specifications, SLF		
Magnification	Work distance (mm)	Field of view (mm)
Galilean Loupes		
2.5 x	340 / 420 / 460 / 500	80 / 90 / 100 / 110
3.0 x	340 / 420 / 460 / 500	60 / 65 / 68 / 70
3.5 x	340 / 420 / 460 / 500	50 / 60 / 70 / 80

Optical specifications, SLH		
Magnification	Work distance (mm)	Field of view (mm)
Galilean Loupes		
2.5 x	340 / 420 / 460 / 500	80 / 90 / 100 / 110
3.0 x	340 / 420 / 460 / 500	60 / 65 / 68 / 70
3.5 x	340 / 420 / 460 / 500	50 / 60 / 70 / 80
Prismatic Loupes		
4.0 x	340 / 420 / 460 / 500	55 / 65 / 75 / 85
5.0 x	340 / 420 / 460 / 500	45 / 55 / 60 / 65
6.0 x	340 / 420 / 460 / 500	40 / 47 / 50 / 58

Environment requirements

Operation	Temperature	+10 °C – +40 °C
	Relative humidity	30% – 75%
	Atmospheric pressure	700 hPa – 1060 hPa
Shipping	Temperature	-40 °C – +55 °C
	Relative humidity	10% – 90%
	Atmospheric pressure	500 hPa – 1060 hPa
Storage	Temperature	-40 °C – +55 °C
	Relative humidity	10% – 90%
	Atmospheric pressure	500 hPa – 1060 hPa

Applicable symbols

	Manufacturer
	Authorized Representative in the EU
	Manufacturing date
	Serial number
	Storage Temperature
	Storage Relative humidity
	Storage Atmospheric pressure
	Keep away from moisture
	Fragile
	CE mark
	Caution
	Consult instructions for use
	Non Sterile

Production date and warranty period

Production date on the body label.

Uselife of the device is 5 years.

The actual warranty period is 1 year.

Limited Warranty

«Zumax Medical CO., LTD» (Zumax) warrants that reasonable care has been used in the design and manufacture of this instrument. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose. Handling, storage and cleaning of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond Zumax's control directly affect the instrument and the results obtained from its use. Zumax's obligation under this warranty is limited to the repair or replacement of this instrument and Zumax shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this instrument. Zumax neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument.

Dispose in accordance with accepted medical practice and applicable local, state, and federal laws and regulations.

Class of hazard for medical wastes by reference to specification of morphological composition (in accordance with Russian Sanitary Regulations and Standards 2.1.7.2790-10) Class A

Zumax Medical Co., Ltd.


5 Zhiying Street, Suzhou New District, Suzhou, Jiangsu, China 215129

Phone : +86 512 66650502

Fax : +86 512 66909655

E-mail : sales@zumaxmedical.com

[Http://www.zumaxmedical.com](http://www.zumaxmedical.com)

MedNet EC-REP GmbH Borkstrasse 10, 48163 Münster, Germany

Authorized representative of the manufacturer in Russia

LLC "N. CELLA"

111033, Moscow, Volochaevskaya st., 12A, bldg. 1


+7 (495) 771-75-39

manager@nsella.ru

Medical Binocular Loupes SLT

User's Guide

CE The CE mark on this product indicates it has been tested to and conforms with the provisions noted within the 93/42/EEC Medical Device Directive.

Zumax Medical Co., Ltd.

Thank you for purchasing this product. To prevent damage to your product or injury to yourself or to others, read the following safety precaution in their entirety before using this equipment. Keep these safety instructions where all those who use the product will read them.

Symbols:

Operating instructions

Warning

Attention. Read operating manual for cautions and instructions for use.

Warnings and Cautions

Warning: Do not expose your product to liquid or moisture or to humidity.

Warning: Do not touch the surface of lens with hand or hard objects.

Warning: Never clean the surface of lens with any hard and dirty materials.

Warning: Do not look into the sun or into a powerful light source with this loupe.

Warning: Federal law restricts this device to sale or order of a physician.

1. Purpose and features of the product

1.1 Purpose

The binocular loupes and headlights are used in Medicine.

Prescription

Designed to increase the image space and / or illuminate the field of view of the dentist during a patient examination or surgery

1.2 Features:

- Multi-coated, high grade glass lenses produce crisp, clear images.
- The Loupes which employ eyepiece of very large view field diameter. The depth of focus is deep and object can be seen stereoscopically with actual distance feeling.
- Smart design, convenient for mobile diagnosis.

2. Instructions for use

SLT loupe is personalized to meet the requirements of each surgeon. The following parameters are preset according to the data on the customer sheet.

- Inter-ocular distance
- Working distance
- Lens correction
- Inclination or declination drilling position depending on the inter ocular distance and position of the eyes
- Position of the nose rest

The device requires fine adjustment before using the first time.

Adjust the sides if necessary and then try it on. SLT should be used with the headband fastened to keep it in position. The side guards clip onto the frame.

3. Maintenance of product

3.1 Maintenance of lens

- The instrument not in use should be covered and put into the packing box.
- Do not disassemble the optical lens by unskilled or unqualified technicians.
- If the product does not work properly, contact an Authorized Representative of the manufacturer. Do not try to repair the product yourself.

3.2 Cleaning

- The outer surface of the instrument may be cleaned with wet cloth, the remaining stains can be cleaned off with the mixture of 50% alcohol and 50% distilled water. Do not wipe with any corrosive detergent lest that the surface should be damaged.
- Cleaning the lenses : Mirror-cleaning paper or a drop of liquid solvent (50 % alcohol and 50 % aether) may be used to clean the lens, blow it carefully. If there is dust stained on the lens, blow them with a blowball or brush them with a dust pen.

Caution : Be careful to prevent the solvent into the instrument.

4. Specifications

Optical specifications, SLT		
Magnification	Work distance (mm)	Field of view (mm)
TTL Loupes		
2.0 x	340/ 420/ 460/ 500	100 / 110 / 130 / 140
2.5 x	340/ 420/ 460/ 500/550	70 /80 / 90 /100 / 120
3.0 x	340/ 420/ 460/ 500/550	60 /70 / 80 / 90 /100
PD		
PD, mm	55; 55,5; 56; 56,5; 57; 57,5; 58; 58,5; 59; 59,5; 60; 60,5; 61; 61,5; 62; 62,5; 63; 63,5; 64; 64,5; 65; 65,5; 66; 66,5; 67; 67,5; 68; 68,5;	

Environment requirements

Operation	Temperature	+10 °C – +40 °C
	Relative humidity	30% – 75%
	Atmospheric pressure	700 hPa – 1060 hPa
Shipping	Temperature	-40 °C – +55 °C
	Relative humidity	10% – 90%
	Atmospheric pressure	500 hPa – 1060 hPa
Storage	Temperature	-40 °C – +55 °C
	Relative humidity	10% – 90%
	Atmospheric pressure	500 hPa – 1060 hPa

5. Product type

The instrument is divided into a frame Loupes and a sports Loupes.

Titanium Frame (figure 1) Sport Frame (figure 2)



Figure 1

Figure 2

Applicable symbols

	Manufacturer
	Authorized Representative in the EU
	Manufacturing date
	Serial number
	Storage Temperature
	Storage Relative humidity
	Storage Atmospheric pressure
	Keep away from moisture
	Fragile
	CE mark
	Caution
	Consult instructions for use
	Non Sterile

Production date and warranty period

Production date on the body label.

Uselife of the device is 5 years.

The actual warranty period is 1 year.

Limited Warranty

«Zumax Medical CO., LTD» (Zumax) warrants that reasonable care has been used in the design and manufacture of this instrument. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose. Handling, storage and cleaning of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond Zumax's control directly affect the instrument and the results obtained from its use. Zumax's obligation under this warranty is limited to the repair or replacement of this instrument and Zumax shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this instrument. Zumax neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument.

Dispose in accordance with accepted medical practice and applicable local, state, and federal laws and regulations.

Class of hazard for medical wastes by reference to specification of morphological composition (in accordance with Russian Sanitary Regulations and Standards 2.1.7.2790-10) Class A

Zumax Medical Co., Ltd.

5 Zhiying Street, Suzhou New District, Suzhou, Jiangsu, China 215129

Phone : +86 512 66650502

Fax : +86 512 66909655

E-mail : sales@zumaxmedical.com

[Http://www.zumaxmedical.com](http://www.zumaxmedical.com)

MedNet EC-REP GmbH Borkstrasse 10, 48163 Münster, Germany


Authorized representative of the manufacturer in Russia

LLC "N. CELLA"

111033, Moscow, Volochaevskaya st., 12A, bldg. 1


+7 (495) 771-75-39

manager@nsella.ru

**Headlight HL8000
to Medical Binocular Loupes**

User's Guide

The CE mark on this product indicates it has been tested to and conforms with the provisions noted within the 93/42/EEC Medical Device Directive.

Zumax Medical Co., Ltd.

Thank you for purchasing this product. To prevent damage to your product or injury to yourself or to others, read the following safety precaution in their entirety before using this equipment. Keep these safety instructions where all those who use the product will read them.

Symbols:

Consult Article

Internal electrical power source

Illumination adjustment mark

Class II

Cautions

1. Read this article for Cautions and instructions for use.
2. Do not use this instrument in damp rooms. Do not expose the instrument to water splashes, dripping water or sprayed water.
3. This product must not be used in the presence of flammable gases.
4. Do not shine the light directly into the eyes.
5. Ensure that the headband power cord is located in the clips provided. The headlight body should not be left covered with a cloth or similar.
6. Please use special certificated adapter.
7. Never scratch the glass cover with fingers or any other hard materials.
8. Federal law restricts this device to sale or order of a physician.

1. Purpose and Features of the product

1.1 Purpose

The binocular loupes and headlights are used in Medicine.

Prescription

Designed to increase the image space and / or illuminate the field of view of the dentist during a patient examination or surgery

1.2 Features

- Light Source: LED source with color temperature of 5000 K .White, shadow-free light allowing the doctor to see tissue characteristics without distortion.
- Enjoy the flexibility of operating the headlight using the rechargeable Lithium battery, can use without the direct power source. It is one of the most comfortable lightweight headlights and suitable for all kinds of surgery.
- The headlight could be used when it is charging.
- Battery indicator lights.
- Adapts to any Binocular Loupes of Zumax.
- Smart design, convenient for mobile diagnosis.

2. Operating Instructions

Attach the power cord by means of the clip provided. Place the device on your head and adjust for width and height to that it fits comfortably.

Connect the Headlight to the power adapter. Use the headband right or left brightness control knob to select the required brightness.

Loosen the fixing screw to adjust the light head up and down and backwards or forwards. After adjustment, retighten the screw. The illuminated field of view can be selected by means of the slide control.

3. Maintenance

Caution: Qualified trained professional persons are needed to maintain the product. Please contact with authorized dealer or manufacture.

Caution: If the product does not work properly, contact an Authorized Representative of the manufacturer. Do not try to repair the product yourself.

3.1 Battery charging and replacing

3.1.1 Battery case BT552 is charged


Battery box BT552 (hereinafter referred to as the battery) yellow light, prompted to charge or replace as soon as possible to recharge the battery backup. If the power indicator only yellow light and continue to use, the battery will start self-protection, cut off the power to prevent the battery over-discharge.

When charging, connect one end of the power adapter output line to the round hole on the side of the battery, and the power adapter is plugged into the network power socket. At this point the charging indicator will flash from yellow to green followed by charging. Charging lights are all off, indicating that the battery is full. When charging, the headset inspection lamp can continue to use.

3.1.2 Battery case BT556 is charged

The headlamp can be plugged directly into the charger to charge when the battery power is depleted.

When charging, you should put the headlight in the charger and connect adapter one end connected to the charger and the other end connected to power supply.

Caution: Please use the appropriate battery. The battery should never be placed in municipal waste. Please check local regulations for disposal of batteries.

BT576

BT556

3.2 Battery use and care

The capacity of Lithium battery will be changed by it's using condition. For maximum battery capacity, use the battery in room temperature(10 °C~40 °C). Do not expose the battery to extreme temperatures, if the battery is used in high or low temperatures, the battery capacity will be reduced.

3.3 Cleaning Recommendations

- The outer surface of the instrument may be cleaned with wet cloth, the remaining stains can be cleaned off with the mixture of 50% alcohol and 50% distilled water. Do not wipe with any corrosive detergent lest that the surface should be damaged.
- Cleaning the lenses: Mirror-cleaning paper or a drop of liquid solvent (50% alcohol and 50 % aether) may be used to clean the lens, blow it carefully. If there is dust stained on the lens, blow them with a blowball or brush them with a dust pen.


Caution : Be careful to prevent the solvent into the device.

4. Installation

4.1 Installing the cradle

According to the battery box mounting hole in the wall, fighting two expansion screw holes and then typing the expansion screw sets and then installing the device as follow the picture.


4.2 Installing the battery into the headband

Installing as shown in below picture, the battery cartridge into the adjusting knob cover, then rotate along time until the level stops.

Removing, counter-clockwise rotation to 90 degree position and then pull it out.

There are two types of battery installation, respectively, up or back.

5. Technical Specification

Optical specifications	
illumination light spot at 420mm	20~80 mm
Illumination in 200mm working distance	≥ 75000 lx
Adjustable range	
Fluctuation distance	≥ 12.5 mm
Angle range	Vertical $\pm 25^\circ$
Weight	
Headlight (excluding head band)	69 g
Option parts	Red color filter
Electrical specifications	
Illumination Source	5W LED, life 20000 hours continuous burn, do not need to change the LED lamp.
Electrical Safety	EN60601-1, EN60601-1-2 Class II and internal electrical power source
Accessories	
Battery BT552	7.4 V Recharged Lithium battery (4400mAh)
Battery BT556	7.4 V Recharged Lithium battery (1800mAh)
Power adapter input voltage	100VAC ~ 240VAC, 50Hz or 60Hz.

Power adapter Output voltage and power	Output voltage 12 V, maximum input power ≤ 30 VA
Battery Run-time	≥ 7 h (Maximum luminance) BT552 ≥ 3 h (Maximum luminance) BT556
Battery Charging-time	≤ 5 h BT552 ≤ 2.5 h BT556
Protection	Short circuit protection, power protection and over current protection.
Longer spring cable (if necessary)	Length ≤ 200 cm
Cable (if necessary)	Length ≤ 19 cm

Electromagnetic compatibility

The use of non-original cables and batteries may result in increased emissions and decreased electromagnetic immunity of the medical device.

The medical electrical device should not be located or used near other devices. If the device is located or used near another device, monitor its operation.

Portable and mobile RF communications equipment can interfere with the Headlights. To avoid exposure, follow the guidelines in the tables below.

Guidance and Manufacturer's Declaration - Electromagnetic Emissions		
The device is intended for use in the electromagnetic environment specified below. The user of the product should ensure that it is used in the specified electromagnetic environment.		
Initial measurements	Requirements fulfilled	Electromagnetic environment - guidelines
RF emissions according to CISPR 11	Group 1	This product only uses very low levels of RF energy and therefore does not cause interference to nearby electronic equipment.
RF emissions according to CISPR 11	Class A	The device intended for use in all premises not used for domestic purposes and not connected to low-voltage distribution networks.
Harmonic emissions according to IEC 61000-3-2	Class A	
Voltage fluctuations / flickering emissions IEC 61000-3-3	Compliant	
NOTE. The emission characteristics of this device are suitable for use in industrial areas and hospitals (standard CISPR 11, class A). If this device is used in a residential environment (which typically requires CISPR 11 Class B), it may not provide adequate protection for radio frequency communications. The user may need to take mitigation measures such as moving or reorienting the device.		

Guidance and Manufacturer's Declaration - Immunity			
The product is intended for use in the electromagnetic environment specified below. The user of the product should ensure that it is used in the specified electromagnetic environment.			
Stability test	Test level according to IEC	Compliance level	Electromagnetic environment - guidelines

	60601		
Electrostatic discharge (ESD) according to IEC 61000-4-2	± 6 kV - contact discharge ± 8 kV - air discharge	± 6 kV - contact discharge ± 8 kV - air discharge	The floor in the room is made of wood, concrete or ceramic tiles. With floors covered with synthetic material, the relative humidity is at least 30%
Nanosecond impulse noise per IEC 61000-4-4	± 2 kV - for power lines ± 1 kV - for input / output lines	± 2 kV - for power lines ± 1 kV - for input / output lines	On-grid power quality based on typical commercial or hospital conditions
High energy microsecond impulse noise according to IEC 61000-4-5	± 1 kV with wire-to-wire interference ± 2 kV with wire-to-ground interference	± 1 kV with wire-to-wire interference ± 2 kV with wire-to-ground interference	The quality of electrical power in the electrical network should be ensured in accordance with a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on input power lines according to IEC 61000-4-11	$<5\%$ (voltage dip $>95\%$) for 0.5 period 40% (voltage dip 60%) for 5 periods 70% (voltage dip 30%) for 25 periods $<5\%$ (voltage dip $>95\%$) for 5 seconds	5% (voltage dip $>95\%$) for 0.5 period 40% (voltage dip 60%) for 5 periods 70% (voltage dip 30%) for 25 periods $<5\%$ (voltage dip $>95\%$) for 5 seconds	On-grid power quality - in line with typical commercial or hospital conditions. If the user needs to ensure continuous operation in conditions of possible power interruptions, it is recommended that the product be powered from an uninterruptible power supply or battery
Industrial magnetic field frequency (50/60 Hz) IEC 61000-4-8	3 A / m	3 A / m	Power frequency magnetic field levels should be provided in accordance with typical commercial or hospital conditions

Guidance and manufacturer's declaration - electromagnetic immunity

The product is intended for use in the electromagnetic environment specified below. The user of the product should ensure that it is used in the specified electromagnetic environment.

Stability test	Test level according to IEC 60601	Compliance level	Electromagnetic environment - guidelines
Conducted disturbance by radio frequency electromagnetic fields according to IEC 61000-4-6	3 V 150 kHz to 80 MHz	3 V	The distance between the mobile radiotelephone communication systems in use and any part of the product, including cables, must be at least the recommended separation separation, which is calculated in accordance with the expression below for the frequency

			of the transmitter.
RF electromagnetic field IEC 61000-4-3	3 V / m 80MHz to 2.5GHz	3 V/m	Recommended distance: $d=3,5/3 * \text{SQRT}(P)$ $d=3,5/3 * \text{SQRT}(P)$ 80MHz to 800MHz $d=7/3 * \text{SQRT}(P)$ 80MHz to 2,5 80GHz Where P is the maximum rated output power of the transmitter in watts (W) according to the manufacturer's data, and d is the recommended distance in meters (m). The field strength during propagation of radio waves from stationary radio transmitters, based on the results of observations of the electromagnetic environment, should be lower than the compliance level in each frequency band. Interference can occur in the vicinity of marked equipment. 
Notes 1. At 80 MHz and 800 MHz, the higher frequency range applies. 2. These guidelines may not apply in some situations. Electromagnetic propagation is affected by absorption or reflection from structures, objects and people.			

Recommended separation distances between portable and mobile RF communications equipment and the product			
The product is intended for use in an electromagnetic environment in which radiated disturbances are controlled. The user of the product can avoid the effects of electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the product as recommended below, taking into account the maximum output power of the communications equipment.			
Rated maximum output power of the transmitter, (W)	Recommended distance (m), depending on the frequency of the transmitter		
	150 kHz to 80 MHz $d=3,5/3 * \text{SQRT}(P)$	80 to 800 MHz $d=3,5/3 * \text{SQRT}(P)$	800 MHz to 2.5 GHz $d=7/3 * \text{SQRT}(P)$
0,01	0,1	0,1	0,2
0,1	0,4	0,4	0,7

1	1,2	1,2	2,3
10	3,7	3,7	7,4
100	11,7	11,7	23,3

When determining the recommended distance d for transmitters with a rated maximum output power not listed in the table, the expressions given are substituted for the rated maximum output power, P , in watts, as stated in the transmitter manufacturer's documentation.

Notes

1. At 80 MHz and 800 MHz, the higher frequency range applies.



2. These guidelines may not apply in some situations.

Electromagnetic propagation is affected by absorption or reflection from structures, objects and people.

Environment requirements

Operation	Temperature	+10 °C – +40 °C
	Relative humidity	30% – 75%
	Atmospheric pressure	700 hPa – 1060 hPa
Shipping	Temperature	-40 °C – +55 °C
	Relative humidity	10% – 90%
	Atmospheric pressure	500 hPa – 1060 hPa
Storage	Temperature	-40 °C – +55 °C
	Relative humidity	10% – 90%
	Atmospheric pressure	500 hPa – 1060 hPa

Applicable symbols

	Manufacturer
	Authorized Representative in the EU
	Manufacturing date
	Serial number
	Storage Temperature
	Storage Relative humidity
	Storage Atmospheric pressure
	Keep away from moisture
	Fragile

	CE mark
	Internal electrical power source
	Illumination adjustment mark
	Class II
	Caution
	Consult instructions for use
	Non Sterile

Production date and warranty period

Production date on the body label.

Useful life of the device is 5 years.

The actual warranty period is 1 year.

Limited Warranty

«Zumax Medical CO., LTD» (Zumax) warrants that reasonable care has been used in the design and manufacture of this instrument. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose. Handling, storage and cleaning of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond Zumax's control directly affect the instrument and the results obtained from its use. Zumax's obligation under this warranty is limited to the repair or replacement of this instrument and Zumax shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this instrument. Zumax neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument.

Dispose in accordance with accepted medical practice and applicable local, state, and federal laws and regulations.

Class of hazard for medical wastes by reference to specification of morphological composition (in accordance with Russian Sanitary Regulations and Standards 2.1.7.2790-10) Class A

Zumax Medical Co., Ltd.


5 Zhiying Street, Suzhou New District, Suzhou, Jiangsu, China 215129

Phone : +86 512 66650502

Fax : +86 512 66909655

E-mail : sales@zumaxmedical.com

[Http://www.zumaxmedical.com](http://www.zumaxmedical.com)

MedNet EC-REP GmbH Borkstrasse 10, 48163 Münster, Germany

Authorized representative of the manufacturer in Russia

LLC "N. CELLA"

111033, Moscow, Volochaevskaya st., 12A, bldg. 1


+7 (495) 771-75-39

manager@nsella.ru

zumax

Headlight HL8200 to Medical Binocular Loupes

User's Guide

The CE mark on this product indicates it has been tested to and conforms with the provisions noted within the 93/42/EEC Medical Device Directive.


Zumax Medical Co., Ltd.

Thank you for purchasing this product. To prevent damage to your product or injury to yourself or to others, read the following safety precaution in their entirety before using this equipment. Keep these safety instructions where all those who use the product will read them.

Symbols:

Consult Article

Internal electrical power source

Illumination adjustment mark

Class II

Cautions

1. Read this article for Cautions and instructions for use.
2. Do not use this instrument in damp rooms. Do not expose the instrument to water splashes, dripping water or sprayed water.
3. This product must not be used in the presence of flammable gases.
4. Do not shine the light directly into the eyes.
5. Ensure that the headband power cord is located in the clips provided. The headlight body should not be left covered with a cloth or similar.
6. Please use special certificated adapter.
7. Never scratch the glass cover with fingers or any other hard materials.
8. Federal law restricts this device to sale or order of a physician.

1. Purpose and Features of the product

1.1 Purpose

The binocular loupes and headlights are used in Medicine.

Prescription

Designed to increase the image space and / or illuminate the field of view of the dentist during a patient examination or surgery

1.2 Features

- Light Source: LED source with color temperature of 5000 K .White, shadow-free light allowing the doctor to see tissue characteristics without distortion.
- Enjoy the flexibility of operating the headlight using the rechargeable Lithium battery, can use without the direct power source. It is one of the most comfortable lightweight headlights and suitable for all kinds of surgery.
- The headlight could be used when it is charging.
- Battery indicator lights.
- Adapts to any Binocular Loupes of Zumax.
- Smart design, convenient for mobile diagnosis.

2. Operating Instructions

Attach the power cord by means of the clip provided. Place the device on your head and adjust for width and height to that it fits comfortably.

Connect the Headlight to the power adapter. Use the headband right or left brightness control knob to select the required brightness.

Loosen the fixing screw to adjust the light head up and down and backwards or forwards. After adjustment, retighten the screw. The illuminated field of view can be selected by means of the slide control.

3. Maintenance

Caution: Qualified trained professional persons are needed to maintenance the product. Please contact with authorized dealer or manufacture.

Caution: If the product does not work properly, contact an Authorized Representative of the manufacturer. Do not try to repair the product yourself.

3.1 Battery charging and replacing

3.1.1 Battery case BT576 is charged


Battery box BT576 (hereinafter referred to as the battery) yellow light, prompted to charge or replace as soon as possible to recharge the battery backup. If the power indicator only yellow light and continue to use, the battery will start self-protection, cut off the power to prevent the battery over-discharge.

When charging, connect one end of the power adapter output line to the round hole on the side of the battery, and the power adapter is plugged into the network power socket. At this point the charging indicator will flash from yellow to green followed by charging. Charging lights are all off, indicating that the battery is full. When charging, the headset inspection lamp can continue to use.

3.1.2 Battery case BT556 is charged

The headlamp can be plugged directly into the charger to charge when the battery power is depleted.


When charging, you should put the headlight in the charger and connect adapter one end connected to the charger and the other end connected to power supply.

Caution: Please use the appropriate battery. The battery should never be placed in municipal waste. Please check local regulations for disposal of batteries.

BT576

BT556

3.2 Battery use and care

The capacity of Lithium battery will be changed by it's using condition. For maximum battery capacity, use the battery in room temperature(10 °C~40 °C). Do not expose the battery to extreme

temperatures, if the battery is used in high or low temperatures, the battery capacity will be reduced.

3.3 Cleaning Recommendations

- The outer surface of the instrument may be cleaned with wet cloth, the remaining stains can be cleaned off with the mixture of 50% alcohol and 50% distilled water. Do not wipe with any corrosive detergent lest that the surface should be damaged.
- Cleaning the lenses: Mirror-cleaning paper or a drop of liquid solvent (50% alcohol and 50% aether) may be used to clean the lens, blow it carefully. If there is dust stained on the lens, blow them with a blowball or brush them with a dust pen.


Caution : Be careful to prevent the solvent into the device.

4. Installation

4.1 Installing the cradle

According to the battery box mounting hole in the wall, fighting two expansion screw holes and then typing the expansion screw sets and then installing the device as follow the picture.


4.2 Installing the battery into the headband

Installing as shown in below picture, the battery cartridge into the adjusting knob cover, then rotate along time until the level stops.


Removing, counter-clockwise rotation to 90 degree position and then pull it out.

There are two types of battery installation, respectively, up or back.


4.3 Headlight HL8200 installation instructions


Headlight HL8200 can be easily mounted on R type headbands, SLE loupes, STL loupes.


SLE loupes + HL8200

SLF loupes + HL8200

R type headbands + HL8200

SLT loupes + HL8200

5. Technical Specification

Optical specifications	
Illumination light spot at 420mm	120 mm
Illumination in 200mm working distance	≥ 35000 lx
Adjustable range	
Fluctuation distance	≥ 32 mm
Angle range	Vertical $\pm 45^\circ$
Weight	
Headlight (excluding head band)	16 g
Option parts	Orange color filter
Electrical specifications	
Illumination Source	3W LED, life 20000 hours continuous burn, do not need to change the LED lamp.
Electrical Safety	EN60601-1, EN60601-1-2 Class II and internal electrical power source
Accessories	
Battery BT576	7.4 V Recharged Lithium battery (2400mAh)
Battery BT556	7.4 V Recharged Lithium battery (1800mAh)
Power adapter input voltage	100VAC ~ 240VAC, 50Hz or 60Hz.
Power adapter Output voltage and power	Output voltage 12 V, maximum input power ≤ 30 VA
Battery Run-time	≥ 7 h (Maximum luminance) BT576 ≥ 3 h (Maximum luminance) BT556
Battery Charging-time	≤ 3.5 h BT576 ≤ 3 h BT556
Protection	Short circuit protection, power protection and over current protection
Longer spring cable (if necessary)	Length ≤ 200 cm
Cable (if necessary)	Length ≤ 19 cm

Electromagnetic compatibility

The use of non-original cables and batteries may result in increased emissions and decreased electromagnetic immunity of the medical device.

The medical electrical device should not be located or used near other devices. If the device is located or used near another device, monitor its operation.

Portable and mobile RF communications equipment can interfere with the Headlights. To avoid exposure, follow the guidelines in the tables below.

Guidance and Manufacturer's Declaration - Electromagnetic Emissions		
The device is intended for use in the electromagnetic environment specified below. The user of the product should ensure that it is used in the specified electromagnetic environment.		
Initial measurements	Requirements fulfilled	Electromagnetic environment - guidelines
RF emissions according to CISPR 11	Group 1	This product only uses very low levels of RF energy and therefore does not cause interference to nearby electronic equipment.

RF emissions according to CISPR 11	Class B	The device is suitable for use in all establishments, including home use, as well as in direct connection with public low-voltage power supply networks of buildings.
Harmonic emissions according to IEC 61000-3-2	Class A	
Voltage fluctuations / flickering emissions IEC 61000-3-3	Compliant	

Guidance and Manufacturer's Declaration - Immunity

The product is intended for use in the electromagnetic environment specified below. The user of the product should ensure that it is used in the specified electromagnetic environment.

Stability test	Test level according to IEC 60601	Compliance level	Electromagnetic environment - guidelines
Electrostatic discharge (ESD) according to IEC 61000-4-2	± 6 kV - contact discharge ± 8 kV - air discharge	± 6 kV - contact discharge ± 8 kV - air discharge	The floor in the room is made of wood, concrete or ceramic tiles. With floors covered with synthetic material, the relative humidity is at least 30%
Nanosecond impulse noise per IEC 61000-4-4	± 2 kV - for power lines ± 1 kV - for input / output lines	± 2 kV - for power lines ± 1 kV - for input / output lines	On-grid power quality based on typical commercial or hospital conditions
High energy microsecond impulse noise according to IEC 61000-4-5	± 1 kV with wire-to-wire interference ± 2 kV with wire-to-ground interference	± 1 kV with wire-to-wire interference ± 2 kV with wire-to-ground interference	The quality of electrical power in the electrical network should be ensured in accordance with a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on input power lines according to IEC 61000-4-11	<5% (voltage dip > 95%) for 0.5 period 40% (voltage dip 60%) for 5 periods 70% (voltage dip 30%) for 25 periods <5% (voltage dip > 95%) for 5 seconds	5% (voltage dip > 95%) for 0.5 period 40% (voltage dip 60%) for 5 periods 70% (voltage dip 30%) for 25 periods <5% (voltage dip > 95%) for 5 seconds	On-grid power quality - in line with typical commercial or hospital conditions. If the user needs to ensure continuous operation in conditions of possible power interruptions, it is recommended that the product be powered from an uninterruptible power supply or battery
Industrial magnetic field frequency (50/60 Hz) IEC 61000-4-8	3 A / m	3 A / m	Power frequency magnetic field levels should be provided in accordance with typical commercial or hospital conditions

Guidance and manufacturer's declaration - electromagnetic immunity

The product is intended for use in the electromagnetic environment specified below. The user of the product should ensure that it is used in the specified electromagnetic environment.

Stability test	Test level according to IEC 60601	Compliance level	Electromagnetic environment - guidelines
Conducted disturbance by radio frequency electromagnetic fields according to IEC 61000-4-6	3 V 150 kHz to 80 MHz	3 V	The distance between the mobile radiotelephone communication systems in use and any part of the product, including cables, must be at least the recommended separation separation, which is calculated in accordance with the expression below for the frequency of the transmitter.
RF electromagnetic field IEC 61000-4-3	3 V / m 80MHz to 2.5GHz	3 V/m	Recommended distance: $d=3,5/3 * \text{SQRT}(P)$ $d=3,5/3 * \text{SQRT}(P)$ 80MHz to 80080MHz $d=7/3 * \text{SQRT}(P)$ 80MHz to 2,5 80GHz Where P is the maximum rated output power of the transmitter in watts (W) according to the manufacturer's data, and d is the recommended distance in meters (m). The field strength during propagation of radio waves from stationary radio transmitters, based on the results of observations of the electromagnetic environment, should be lower than the compliance level in each frequency band. Interference can occur in the vicinity of marked equipment. 

Notes

1. At 80 MHz and 800 MHz, the higher frequency range applies.

2. These guidelines may not apply in some situations.

Electromagnetic propagation is affected by absorption or reflection from structures, objects and people.




Recommended separation distances between portable and mobile RF communications

equipment and the product			
The product is intended for use in an electromagnetic environment in which radiated disturbances are controlled. The user of the product can avoid the effects of electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the product as recommended below, taking into account the maximum output power of the communications equipment.			
Rated maximum output power of the transmitter, (W)	Recommended distance (m), depending on the frequency of the transmitter		
	150 kHz to 80 MHz $d=3,5/3 * \text{SQRT}(P)$	80 to 800 MHz $d=3,5/3 * \text{SQRT}(P)$	800 MHz to 2.5 GHz $d=7/3 * \text{SQRT}(P)$
0,01	0,1	0,1	0,2
0,1	0,4	0,4	0,7
1	1,2	1,2	2,3
10	3,7	3,7	7,4
100	11,7	11,7	23,3
When determining the recommended distance d for transmitters with a rated maximum output power not listed in the table, the expressions given are substituted for the rated maximum output power, P, in watts, as stated in the transmitter manufacturer's documentation.			
Notes 1. At 80 MHz and 800 MHz, the higher frequency range applies. 2. These guidelines may not apply in some situations. Electromagnetic propagation is affected by absorption or reflection from structures, objects and people.			

Environment requirements

Operation	Temperature	+10 °C – +40 °C
	Relative humidity	30% – 75%
	Atmospheric pressure	700 hPa – 1060 hPa
Shipping	Temperature	-40 °C – +55 °C
	Relative humidity	10% – 90%
	Atmospheric pressure	500 hPa – 1060 hPa
Storage	Temperature	-40 °C – +55 °C
	Relative humidity	10% – 90%
	Atmospheric pressure	500 hPa – 1060 hPa

Applicable symbols

	Manufacturer
	Authorized Representative in the EU
	Manufacturing date
	Serial number
	Storage Temperature
	Storage Relative humidity
	Storage Atmospheric pressure
	Keep away from moisture
	Fragile
	CE mark
	Internal electrical power source
	Illumination adjustment mark
	Class II
	Caution
	Consult instructions for use
	Non Sterile

Production date and warranty period

Production date on the body label.

Uselife of the device is 5 years.

The actual warranty period is 1 year.

Limited Warranty

«Zumax Medical CO., LTD» (Zumax) warrants that reasonable care has been used in the design and manufacture of this instrument. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose. Handling, storage and cleaning of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond Zumax's control directly affect the instrument and the results obtained from its use. Zumax's obligation under this warranty is limited to the repair or replacement of this instrument and Zumax shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this instrument. Zumax neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument.

Dispose in accordance with accepted medical practice and applicable local, state, and federal laws and regulations.

Class of hazard for medical wastes by reference to specification of morphological composition (in accordance with Russian Sanitary Regulations and Standards 2.1.7.2790-10) Class A

Zumax Medical Co., Ltd.


5 Zhiying Street, Suzhou New District, Suzhou, Jiangsu, China 215129

Phone : +86 512 66650502

Fax : +86 512 66909655

E-mail : sales@zumaxmedical.com

[Http://www.zumaxmedical.com](http://www.zumaxmedical.com)

MedNet EC-REP GmbH Borkstrasse 10, 48163 Münster, Germany


Authorized representative of the manufacturer in Russia

LLC "N. CELLA"

111033, Moscow, Volochaevskaya st., 12A, bldg. 1


+7 (495) 771-75-39

manager@nsella.ru

**Headlight HL8300
to Medical Binocular Loupes**

User's Guide

CE The CE mark on this product indicates it has been tested to and conforms with the provisions noted within the 93/42/EEC Medical Device Directive.


Zumax Medical Co., Ltd.

Thank you for purchasing this product. To prevent damage to your product or injury to yourself or to others, read the following safety precaution in their entirety before using this equipment. Keep these safety instructions where all those who use the product will read them.

Symbols:

Consult Article

Internal electrical power source

Illumination adjustment mark

Class II

Cautions

1. Read this article for Cautions and instructions for use.
2. Do not use this instrument in damp rooms. Do not expose the instrument to water splashes, dripping water or sprayed water.
3. This product must not be used in the presence of flammable gases.
4. Do not shine the light directly into the eyes.
5. Ensure that the headband power cord is located in the clips provided. The headlight body should not be left covered with a cloth or similar.
6. Please use special certificated adapter.
7. Never scratch the glass cover with fingers or any other hard materials.
8. Federal law restricts this device to sale or order of a physician.

1. Purpose and Features of the product

1.1 Purpose

The binocular loupes and headlights are used in Medicine.

Prescription

Designed to increase the image space and / or illuminate the field of view of the dentist during a patient examination or surgery

1.2 Features

- Light Source: LED source with color temperature of 5000 K .White, shadow-free light allowing the doctor to see tissue characteristics without distortion.
- Enjoy the flexibility of operating the headlight using the rechargeable Lithium battery, can use without the direct power source. It is one of the most comfortable lightweight headlights and suitable for all kinds of surgery.
- The headlight could be used when it is charging.
- Battery indicator lights.
- Adapts to any Binocular Loupes of Zimax.
- Smart design, convenient for mobile diagnosis.

2. Operating Instructions

Attach the power cord by means of the clip provided. Place the device on your head and adjust for width and height to that it fits comfortably.

Connect the Headlight to the power adapter. Use the headband right or left brightness control knob to select the required brightness.

Loosen the fixing screw to adjust the light head up and down and backwards or forwards. After adjustment, retighten the screw. The illuminated field of view can be selected by means of the slide control.

3. Maintenance

Caution: Qualified trained professional persons are needed to maintenance the product. Please contact with authorized dealer or manufacture.

Caution: If the product does not work properly, contact an Authorized Representative of the manufacturer. Do not try to repair the product yourself.

3.1 Battery charging and replacing

3.1.1 Battery case BT576 is charged

Battery box BT576 (hereinafter referred to as the battery) yellow light, prompted to charge or replace as soon as possible to recharge the battery backup. If the power indicator only yellow light and continue to use, the battery will start self-protection, cut off the power to prevent the battery over-discharge.

When charging, connect one end of the power adapter output line to the round hole on the side of the battery, and the power adapter is plugged into the network power socket. At this point the charging indicator will flash from yellow to green followed by charging. Charging lights are all off, indicating that the battery is full. When charging, the headset inspection lamp can continue to use.


3.1.2 Battery case BT556 is charged

The headlamp can be plugged directly into the charger to charge when the battery power is depleted.

When charging, you should put the headlight in the charger and connect adapter one end connected to the charger and the other end connected to power supply.

Caution: Please use the appropriate battery. The battery should never be placed in municipal waste. Please check local regulations for disposal of batteries.

BT576

BT556

3.2 Battery use and care

The capacity of Lithium battery will be changed by it's using condition. For maximum battery capacity, use the battery in room temperature(10 °C~40 °C). Do not expose the battery to extreme

temperatures, if the battery is used in high or low temperatures, the battery capacity will be reduced.

3.3 Cleaning Recommendations

- The outer surface of the instrument may be cleaned with wet cloth, the remaining stains can be cleaned off with the mixture of 50% alcohol and 50% distilled water. Do not wipe with any corrosive detergent lest that the surface should be damaged.
- Cleaning the lenses: Mirror-cleaning paper or a drop of liquid solvent (50% alcohol and 50% ether) may be used to clean the lens, blow it carefully. If there is dust stained on the lens, blow them with a blowball or brush them with a dust pen.


Caution : Be careful to prevent the solvent into the device.

4. Installation

4.1 Installing the cradle

According to the battery box mounting hole in the wall, finding two expansion screw holes and then typing the expansion screw sets and then installing the device as follow the picture.





Installing the battery into the headband

Installing as shown in below picture, the battery cartridge into the adjusting knob cover, then rotate along time until the level stops.


Removing, counter-clockwise rotation to 90 degree position and then pull it out.

There are two types of battery installation, respectively, up or back.


4.3 Headlight HL8300 installation instructions


Headlights can be easily mounted on R type headbands, SLE loupes, TTL


SLE loupes + HL8300

SLF loupes + HL8300

R type headbands + HL8300

SLT loupes + HL8200

5. Technical Specification

Optical specifications	
illumination light spot at 420mm	80 mm
Illumination in 200mm working distance	≥ 35000 lx
Adjustable range	
Fluctuation distance	≥ 32 mm
Angle range	Vertical $\pm 45^\circ$
Weight	
Headlight (excluding head band)	27 ± 3 g
Option parts	Orange color filter
Electrical specifications	
Illumination Source	3W LED, life 20000 hours continuous burn, do not need to change the LED lamp.
Electrical Safety	EN60601-1, EN60601-1-2 Class II and internal electrical power source
Accessories	
Battery BT576	7.4 V Recharged Lithium battery (2400mAh)
Battery BT556	7.4 V Recharged Lithium battery (1800mAh)
Power adapter input voltage	100VAC ~ 240VAC, 50Hz or 60Hz.
Power adapter Output voltage and power	Output voltage 12 V, maximum input power ≤ 30 VA
Battery Run-time	≥ 7 h (Maximum luminance) BT576 ≥ 3 h (Maximum luminance) BT556
Battery Charging-time	≤ 3.5 h BT576 ≤ 3 h BT556
Protection	Short circuit protection, power protection and over current protection
Longer spring cable (if necessary)	Length ≤ 200 cm
Cable (if necessary)	Length ≤ 19 cm

Electromagnetic compatibility

The use of non-original cables and batteries may result in increased emissions and decreased electromagnetic immunity of the medical device.

The medical electrical device should not be located or used near other devices. If the device is located or used near another device, monitor its operation.

Portable and mobile RF communications equipment can interfere with the Headlights. To avoid exposure, follow the guidelines in the tables below.

Guidance and Manufacturer's Declaration - Electromagnetic Emissions		
The device is intended for use in the electromagnetic environment specified below. The user of the product should ensure that it is used in the specified electromagnetic environment.		
Initial measurements	Requirements fulfilled	Electromagnetic environment - guidelines
RF emissions according to CISPR 11	Group 1	This product only uses very low levels of RF energy and therefore does not cause interference to nearby electronic equipment.

RF emissions according to CISPR 11	Class B	The device is suitable for use in all establishments, including home use, as well as in direct connection with public low-voltage power supply networks of buildings.
Harmonic emissions according to IEC 61000-3-2	Class A	
Voltage fluctuations / flickering emissions IEC 61000-3-3	Compliant	

Guidance and Manufacturer's Declaration - Immunity			
The product is intended for use in the electromagnetic environment specified below. The user of the product should ensure that it is used in the specified electromagnetic environment.			
Stability test	Test level according to IEC 60601	Compliance level	Electromagnetic environment - guidelines
Electrostatic discharge (ESD) according to IEC 61000-4-2	± 6 kV - contact discharge ± 8 kV - air discharge	± 4 kV - contact discharge ± 8 kV - air discharge	The floor in the room is made of wood, concrete or ceramic tiles. With floors covered with synthetic material, the relative humidity is at least 30%. It is recommended to use antistatic materials or to operate the device at high relative humidity.
Nanosecond impulse noise per IEC 61000-4-4	± 2 kV - for power lines ± 1 kV - for input / output lines	± 2 kV - for power lines ± 1 kV - for input / output lines	On-grid power quality based on typical commercial or hospital conditions
High energy microsecond impulse noise according to IEC 61000-4-5	± 1 kV with wire-to-wire interference ± 2 kV with wire-to-ground interference	± 1 kV with wire-to-wire interference ± 2 kV with wire-to-ground interference	The quality of electrical power in the electrical network should be ensured in accordance with a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on input power lines according to IEC 61000-4-11	<5% (voltage dip > 95%) for 0.5 period 40% (voltage dip 60%) for 5 periods 70% (voltage dip 30%) for 25 periods <5% (voltage dip > 95%) for 5 seconds	5% (voltage dip > 95%) for 0.5 period 40% (voltage dip 60%) for 5 periods 70% (voltage dip 30%) for 25 periods <5% (voltage dip > 95%) for 5 seconds	On-grid power quality - in line with typical commercial or hospital conditions. If the user needs to ensure continuous operation in conditions of possible power interruptions, it is recommended that the product be powered from an uninterruptible power supply or battery
Industrial magnetic field	3 A / m	3 A / m	Power frequency magnetic field levels should be

frequency (50/60 Hz) IEC 61000-4-8			provided in accordance with typical commercial or hospital conditions
---------------------------------------	--	--	---

Guidance and manufacturer's declaration - electromagnetic immunity			
The product is intended for use in the electromagnetic environment specified below. The user of the product should ensure that it is used in the specified electromagnetic environment.			
Stability test	Test level according to IEC 60601	Compliance level	Electromagnetic environment - guidelines
Conducted disturbance by radio frequency electromagnetic fields according to IEC 61000-4-6	3 V 150 kHz to 80 MHz	3 V	The distance between the mobile radiotelephone communication systems in use and any part of the product, including cables, must be at least the recommended separation separation, which is calculated in accordance with the expression below for the frequency of the transmitter.
RF electromagnetic field IEC 61000-4-3	3 V / m 80MHz to 2.5GHz	3 V/m	Recommended distance: $d = 3,5/3 * \sqrt{P}$ $d = 3,5/3 * \sqrt{P}$ 80MHz to 80080MHz $d = 7/3 * \sqrt{P}$ 80MHz to 2,5 80GHz Where P is the maximum rated output power of the transmitter in watts (W) according to the manufacturer's data, and d is the recommended distance in meters (m). The field strength during propagation of radio waves from stationary radio transmitters, based on the results of observations of the electromagnetic environment, should be lower than the compliance level in each frequency band. Interference can occur in the vicinity of marked equipment. 
Notes 1. At 80 MHz and 800 MHz, the higher frequency range applies. 2. These guidelines may not apply in some situations.			

Electromagnetic propagation is affected by absorption or reflection from structures, objects and people.

Recommended separation distances between portable and mobile RF communications equipment and the product

The product is intended for use in an electromagnetic environment in which radiated disturbances are controlled. The user of the product can avoid the effects of electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the product as recommended below, taking into account the maximum output power of the communications equipment.

Rated maximum output power of the transmitter, (W)	Recommended distance (m), depending on the frequency of the transmitter		
	150 kHz to 80 MHz $d=3,5/3 * \text{SQRT}(P)$	80 to 800 MHz $d=3,5/3 * \text{SQRT}(P)$	800 MHz to 2.5 GHz $d=7/3 * \text{SQRT}(P)$
0,01	0,1	0,1	0,2
0,1	0,4	0,4	0,7
1	1,2	1,2	2,3
10	3,7	3,7	7,4
100	11,7	11,7	23,3

When determining the recommended distance d for transmitters with a rated maximum output power not listed in the table, the expressions given are substituted for the rated maximum output power, P , in watts, as stated in the transmitter manufacturer's documentation.

Notes



1. At 80 MHz and 800 MHz, the higher frequency range applies.
2. These guidelines may not apply in some situations.

Electromagnetic propagation is affected by absorption or reflection from structures, objects and people.

Environment requirements

Operation	Temperature	+10 °C – +40 °C
	Relative humidity	30% – 75%
	Atmospheric pressure	700 hPa – 1060 hPa
Shipping	Temperature	-40 °C – +55 °C
	Relative humidity	10% – 90%
	Atmospheric pressure	500 hPa – 1060 hPa
Storage	Temperature	-40 °C – +55 °C
	Relative humidity	10% – 90%
	Atmospheric pressure	500 hPa – 1060 hPa

Applicable symbols

	Manufacturer
	Authorized Representative in the EU
	Manufacturing date
	Serial number
	Storage Temperature
	Storage Relative humidity
	Storage Atmospheric pressure
	Keep away from moisture
	Fragile
	CE mark
	Internal electrical power source
	Illumination adjustment mark
	Class II
	Caution
	Consult instructions for use
	Non Sterile

Production date and warranty period

Production date on the body label.

Useful life of the device is 5 years.

The actual warranty period is 1 year.

Limited Warranty

«Zumax Medical CO., LTD» (Zumax) warrants that reasonable care has been used in the design and manufacture of this instrument. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose. Handling, storage and cleaning of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond Zumax's control directly affect the instrument and the results obtained from its use. Zumax's obligation under this warranty is limited to the repair or replacement of this instrument and Zumax shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this instrument. Zumax neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument.

Dispose in accordance with accepted medical practice and applicable local, state, and federal laws and regulations.

Class of hazard for medical wastes by reference to specification of morphological composition (in accordance with Russian Sanitary Regulations and Standards 2.1.7.2790-10) Class B

Zumax Medical Co., Ltd.


5 Zhiying Street, Suzhou New District, Suzhou, Jiangsu, China 215129

Phone : +86 512 66650502

Fax : +86 512 66909655

E-mail : sales@zumaxmedical.com

[Http://www.zumaxmedical.com](http://www.zumaxmedical.com)

MedNet EC-REP GmbH Borkstrasse 10, 48163 Münster, Germany

Authorized representative of the manufacturer in Russia

LLC "N. CELLA"

111033, Moscow, Volochaevskaya st., 12A, bldg. 1

+7 (495) 771-75-39

manager@nsella.ru

Versions and composition of medical device

Binocular medical loupe, with accessories

1. Binocular loupe, versions: SLE, SLT, SLF, composed of:
 1. Binocular loupe;
 2. Screwdriver (if necessary);
 3. Band;
 4. Cloth;
 5. Case;
 6. User manual.
2. Binocular loupe SLH on headband composed of:
 1. Binocular loupe;
 2. Removable pads (3 pcs. pack);
 3. Screwdriver (if necessary)
 4. Cloth;
 5. Case;
 6. Case key (if necessary);
 7. User manual.
3. Binocular loupe, versions: SLE, SLT, SLF with headlight on headband, composed of:
 1. Binocular loupe;
 2. Screwdriver (if necessary);
 3. Band;
 4. Cloth;
 5. Headlight HL8000 on headband;
 6. Charger;
 7. Battery;
 8. Longer spring cable;
 9. Safety cap (3 pcs.);
 10. Removable pads (3 pcs. pack);
 11. Binocular loupe fasteners (if necessary);
 12. Case;
 13. Case key;
 14. User manual.
4. Binocular loupe, versions: SLE, SLT, SLF with cordless headlight on headband, composed of:
 1. Binocular loupe;
 2. Screwdriver (if necessary);
 3. Band;
 4. Cloth;
 5. Cordless headlight HL8000 on headband;
 6. Charger;
 7. Battery;
 8. Cable;
 9. Safety cap (3 pcs.);

10. Removable pads (3 pcs. pack);
 11. Binocular loupe fasteners (if necessary);
 12. Case;
 13. Case key;
 14. Wall mount;
 15. User manual.
5. Binocular loupe SLH with headlight on headband, composed of:
1. Binocular loupe;
 2. Removable pads (3 pcs. pack);
 3. Cloth;
 4. Headlight HL8000 on headband;
 5. Charger;
 6. Battery;
 7. Longer spring cable;
 8. Safety cap (3 pcs.);
 9. Case;
 10. Case key;
 11. User manual.
6. Binocular loupe SLH with cordless headlight on headband, composed of:
1. Binocular loupe;
 2. Removable pads (3 pcs. pack);
 3. Cloth;
 4. Cordless headlight HL8000 on headband;
 5. Charger;
 6. Battery;
 7. Longer spring cable;
 8. Safety cap (3 pcs.);
 9. Binocular loupe fasteners (if necessary)
 10. Case;
 11. Case key;
 12. Wall mount;
 13. User manual.
7. Binocular loupe, versions: SLE, SLT, SLF with headlight, composed of:
1. Binocular loupe;
 2. Screwdriver (if necessary);
 3. Band;
 4. Cloth;
 5. Headlight, versions: HL8200, HL8300;
 6. Battery;
 7. Charger;
 8. Longer spring cable;
 9. Orange filter (if necessary);
 10. Binocular loupe fasteners (if necessary);
 11. Screwdriver (if necessary);
 12. Case;
 13. User manual.

8. Binocular loupe, versions: SLE, SLT, SLF with cordless headlight, composed of:

1. Binocular loupe;
2. Screwdriver (if necessary);
3. Band;
4. Cloth;
5. Cordless headlight, versions: HL8200, HL8300;
6. Battery;
7. Charger;
8. Cable;
9. Orange filter (if necessary);
10. Wall mount;
11. Case;
12. User manual.

9. Binocular loupe SLH on headband with headlight, composed of:

1. Binocular loupe;
2. Removable pads (3 pcs. pack);
3. Cloth;
4. Headlight, versions: HL8200, HL8300;
5. Battery;
6. Charger;
7. Longer spring cable;
8. Orange filter (if necessary);
9. Binocular loupe fasteners (if necessary);
10. Screwdriver (if necessary);
11. Case;
12. Case key;
13. User manual.

10. Binocular loupe SLH on headband with cordless headlight, composed of:

1. Binocular loupe;
2. Removable pads (3 pcs. pack);
3. Cloth;
4. Cordless headlight, versions: HL8200, HL8300;
5. Battery;
6. Charger;
7. Cable;
8. Orange filter (if necessary);
9. Case;
10. Case key;
11. Wall mount;
12. User manual.

11. Headlight HL8000 on headband, composed of:

1. Headlight;
2. Charger;
3. Battery;

4. Longer spring cable;
5. Safety cap (3 pcs.);
6. Removable pads (3 pcs. pack);
7. Binocular loupe fasteners (if necessary);
8. Case;
9. Case key;
10. User manual.

12. Cordless headlight HL8000 on headband, composed of:

1. Cordless headlight;
2. Charger;
3. Battery;
4. Cable;
5. Safety cap (3 pcs.);
6. Removable pads (3 pcs. pack);
7. Binocular loupe fasteners (if necessary);
8. Case;
9. Case key;
10. Wall mount;
11. User manual.

13. Headlight, versions: HL8200, HL8300, composed of:

1. Headlight;
2. Battery;
3. Charger;
4. Longer spring cable;
5. Orange filter (if necessary);
6. Binocular loupe fasteners (if necessary);
7. Screwdriver (if necessary);
8. Case;
9. User manual.

14. Cordless headlight, versions: HL8200, HL8300, composed of:

1. Cordless headlight;
2. Battery;
3. Charger;
4. Cable;
5. Orange filter (if necessary);
6. Wall mount;
7. Case;
8. Case key
9. User manual.

Accessories:

1. Removable pads (3 pcs. pack)
2. Screwdriver
3. Case key

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