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20007601-2 E D C

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CERTIFICATE

中国国际贸易促进委员会暨中国国际商会
China Council for the Promotion of International Trade is China Chamber of International Commerce

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China Council for the Promotion of International Trade
China Chamber of International Commerce

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China Council for the Promotion
of International Trade

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证书查询网址 Website for verifying the certificate: <http://www.rzccpit.com/validate.html>

«I CERTIFY»

Director
Wang Kefeng

LINYI 2H HEALTHCARE PRODUCTS CO.,
LTD

Date «01» June 2020 r.

Signature Wang Kefeng

Stamp

EXTRACT FROM CONSOLIDATED TECHNICAL DOCUMENTATION
of the medical device:

Medical disposable absorbent devices

2020

1

The following extract was drawn up to provide information for the purposes of state registration of the medical device in the Russian Federation.

The information is provided to the extent required in accordance with the Rules of state registration of medical devices (RF Government Decree No. 1416 dated December 27, 2012).

1. Name of the medical device and design variants

Medical disposable absorbent devices in the following versions:

1. Dr.Skipp pull-up diapers for adults, size M (Medium) (type 233900);
2. Dr.Skipp pull-up diapers for adults, size L (Large) (type 233900);
3. Dr.Skipp diapers for adults, size M (Medium) (type 280360);
4. Dr.Skipp diapers for adults, size L (Large) (type 280360).

Hereafter referred to as medical devices.

2. Information about MD manufacturer, MD developer, addresses of MD manufacturing sites, MD manufacturer authorized representative in the Russian Federation

2.1 Information about the manufacturer of the medical device

Name: LINYI 2H HEALTHCARE PRODUCTS CO., LTD.);

Place of business address: NO.58 LUOQI ROAD, SHENGZHUANG STREET, LUOZHUANG DISTRICT, LINYI CITY, SHANDONG PROVINCE, CHINA;

Phone: +86 (0539) 5635398;

e-mail: sales@ak2h.com.

2.2 Information about the developer of the medical device

Name: LINYI 2H HEALTHCARE PRODUCTS CO., LTD.

Place of business address: NO.58 LUOQI ROAD, SHENGZHUANG STREET, LUOZHUANG DISTRICT, LINYI CITY, SHANDONG PROVINCE, CHINA;

Phone: +86 (0539) 5635398;

e-mail: sales@ak2h.com.

2.3 Address of the manufacturing site of the medical device

LINYI 2H HEALTHCARE PRODUCTS CO., LTD, NO.58 LUOQI ROAD, SHENGZHUANG STREET, LUOZHUANG DISTRICT, LINYI CITY, SHANDONG PROVINCE, CHINA.

2.4 Authorized representative of the manufacturer in Russia

Name: RM-plus, Limited Liability Company (RM-plus, LLC);

Place of business address: 19, Nauchny proezd, 8/7 floor/office, Moscow, 117246

Phone: +7 (495) 642-82-60;

e-mail: msk@rmplus.ru.

3. Functional purpose of the medical device

3.1 Purpose

Medical disposable absorbent devices are designed to absorb and retain urine when caring for adult patients suffering from moderate to severe incontinence, bedridden patients.

3.2 Indications for use

Indications for prescription of products are persistent pronounced or significantly pronounced excretion function disorders resulting from:

- diseases, consequences of injuries, malformations of central, peripheral nervous system, including conditions resulting from cerebrovascular disorder;
- diseases, consequences of trauma, chronic infections and malformations of genitourinary system;
- malignant neoplasms (clinical group IV);
- mental disorders with severe or profound mental deficiency, dementia;
- other diseases accompanied by moderate to severe urinary incontinence.

3.3 Contraindications for use

An absolute medical contraindication for providing incontinent adults with medical devices is an allergic reaction to the material from which the products are made.

It is not recommended to use the products in case of diseases of the external genital organs or skin lesions in the places of body contact with the product.

3.4 Adverse effects

Skin irritation in the area of contact with the product is possible if the patient has an allergic reaction to the material from which the products are made.

3.5 Scope of application

Care for patients suffering from urinary incontinence, care for bedridden patients.

3.6 Frequency of application

Non-sterile and not subject to sterilization medical devices intended for single use. The product is disposable and is intended for single patient use. The recommended time for use of the product is no more than 8 hours (with polyuria syndrome — no more than 5 hours).


3.7 Instruction for use

The method of application of the medical device is described in the instruction for use of the medical device.

The rules for use of products (represented as pictures) are shown on the consumer packaging.


Algorithm for putting diapers on:

1. Remove the diaper from the consumer package, straighten all its parts, smooth out folds, stretch slightly;
2. Put the patient on the side, bending the knees;
3. Fold the product and pass it between slightly open legs, closing the perineal region - putting on should start from the front, pushing the diaper backwards and straightening the back;
4. Turn the patient onto the back and align the diaper on the abdomen leaving minimum folds;
5. The edges of the product should fit snugly to the body; this will help avoid leakage. Making sure of sufficient fit, fix the diaper with the Velcro - first with the lower ones to wrap the legs more tightly, then with the upper ones.

Algorithm for putting the pull-up diapers on:

1. Remove the pull-up diapers from the consumer package;
2. Put on as usual underwear, straighten out the protective barriers.
3. Take off the used pull-up diapers as normal underwear, or tear the side seams.

3.8 Potential users

Adult patients suffering from urinary incontinence.

3.9 Precautions, warnings

- 3.9.1 Products should not be stored in highly humid environment, near heating devices and in direct sunlight.
- 3.9.2 Products should not be stored with laundry detergents, machine oil and other strong smelling substances, as the products may be saturated with foreign odors, even in a sealed package.
- 3.9.3 Products should not be stored near open flames, as it may deform and discolor the product.
- 3.9.4 If the skin condition in the area of contact with the product deteriorates and irritation occurs, stop using the product and refer to a doctor.
- 3.9.5 The products are for external use only.
- 3.9.6 It is forbidden to use expired products. These products should be disposed of.

3.10 Conditions of use of the medical device

The products are intended for use in healthcare institutions and at home, at environmental temperature of +5 °C to +40 °C and relative air humidity up to 98 %.

4. Type of contact with the human body

Short-term contact with skin;

The products are resistant to biological fluids and body tissue secretions.

5. Completeness of delivery

Dr.Skipp pull-up diapers for adults, size M (Medium), are supplied in consumer packaging, 20 pcs per package.

Dr.Skipp pull-up diapers for adults, size L (Large), are supplied in consumer packaging, 20 pcs per package.

Dr.Skipp diapers for adults, size M (Medium), are supplied in consumer packaging, 30 pcs per package.

Dr.Skipp diapers for adults, size L (Large), are supplied in consumer packaging, 30 pcs per package.

instruction for use, 1 pc., is enclosed to each consumer package.

6. Medical device classification

6.1 Medical device risk class

The medical device belongs to class I according to Appendix IX to the current version of MDD/93/42/EEC.

6.3 Use of medicinal drugs

The products do not contain any medicinal drugs as an integral part required for their operation.

6.4 Use of products of animal origin

The products do not contain any derivatives of animal origin as an integral part required for their operation.

6.5 Use of blood and its components

The products do not contain any derivatives of human blood or its components as an integral part required for their operation.

6.6 Use with other medical devices

Medical disposable absorbent devices are used independently of other medical devices.

6.7 Risk analysis results

Risk analysis was performed in accordance with the EN ISO 14971 standard.

7. Description and operation principle of the medical device

The products are multi-layer absorbent disposable products with an absorbent layer made of fibrous semi-finished products containing moisture-absorbing substances (superabsorbents). The products allow incontinent patients to carry out independent active life (from creating an opportunity to perform family and household duties to carrying out professional and production activities), despite the pathological changes in the body.

Products provide dry skin, prevent reproduction of pathogenic flora, i.e. used as a means of rehabilitation and care for prevention of skin damage, creating physical and psychological comfort.

Medical disposable absorbent devices are manufactured:

- 1) in the form of ready-made briefs (pull-up diapers)
- 2) in the form of cutting of briefs with fixing elements (diapers).

Depending on the amount of absorbed fluid, the products can be used for people with moderate to severe incontinence.

The shape of the diapers and pull-up diapers should correspond to the development of a part of the torso of the person with an additional increase in the area for the wrapover of the side parts.

The absorbent layer of the products is shaped to suit both men and women. Diapers and pull-up diapers are a multi-layer product.

The design of the products includes (starting from the layer in contact with human skin):

- **top cover layer:** the layer that directly contacts the human skin and allows the penetration of liquid into the product.
- **absorbent layer:** the inner base layer of the product that absorbs and retains the absorbable liquid inside the product.
- **protective layer:** the layer that is located directly behind the absorbent layer of the product and prevents the passage of liquid outside the device.
- **barrier elements:** barriers and side ruffles in the form of an arc through the groin with elastic bands, preventing the penetration of liquid onto the human skin, as well as ensuring the anatomical shape of the product and the best fit.
- **fixing elements (only for diapers):** parts of the diaper on its front and back parts to fix it in the desired position and fasten it to the person's waist (reusable Velcro fasteners in the amount of 4 pieces, allowing to perfectly fit the diaper to the anatomical features of any patient). Diapers, in the version "Dr.Skipp diapers for adults, size L (Large)" have an elastic crimped tie on the lower back, which allows you to adapt the diaper to any shape.
- **diaper wetness Indicator:** signals in the form of colored stripes (for diapers – blue stripes, for pull-up diapers – red stripes) on the protective layer, fixing the full filling of the product with the absorbed liquid.

The protective layer and fixing elements (for diapers) of products are made without finishing.

External defects are not allowed in the products - mechanical damages (rupture of edges, cuts, damage to fixing elements, etc.), stains of different origin, foreign inclusions visible to the naked eye.

The layers of the products are bonded by means of heat treatment, which ensures the strength of the gluing of the layers (seams) of the device. The seams should be continuous.

Diapers are made in the following sizes depending on the waist/hip circumference:

- 1) Size M (Medium): Waist/hips 68-115 cm;
- 2) Size L (Large): Waist/hips 84-150 cm;

Pull-up diapers are made in the following sizes depending on the Waist/hip circumference:

- 1) Size M (Medium): Waist/hips 80-120 cm;
- 2) Size L (Large): Waist/hips 100-140 cm;

Photographic images of products are shown in Figure 1 and Figure 2.

Figure 1 – Dr Skipp pull-up diapers for adults, size M (Medium)

Figure 2 – Dr Skipp diapers for adults, size L (Large)

8. Technical characteristics

8.1 The parameters of the products should correspond to the parameters specified in Table 1.

Table 1.

Name of parameter	Design variant			
	Dr. Skipp pull-up diapers for adults, size M (Medium)	Dr. Skipp pull-up diapers for adults, size L (Large)	Dr. Skipp diapers for adults, size M (Medium)	Dr. Skipp diapers for adults, size L (Large)
Length, mm	320±10	350±10	750±10	870±10
Middle width, mm	180±10	200±10	290±10	290±10
Width on edges, mm	260±10	300±10	750±10	750±10
Width on edges in stretched state, mm	620±10	700±10	Not applicable	Not applicable
Waist/hip circumference, mm	800-1200	1000-1400	680-1150	840-1500
Absorbent layer length, mm	420±10	450±10	520±10	720±10
Absorbent layer width in middle, mm	80±10	80±10	150±10	170±10
Absorbent layer width on edges, mm	130±10	130±10	235±10	280±10
Quantity of tightening elastic bands on side ruffles in the form of an arc through groin, pcs.	2	2	2	2
Quantity of tightening elastic bands on barriers, pcs.	3	3	3	3
Quantity of tightening elastic bands on band, pcs.	18	18	Not applicable	Not applicable
Fixing element length, mm	Not applicable	Not applicable	80±5	80±5
Fixing element width, mm	Not applicable	Not applicable	26±5	26±5
Crimped tie length, mm	Not applicable	Not applicable	Not applicable	125±5
Crimped tie width, mm	Not applicable	Not applicable	Not applicable	45±5
Crimped tie length, in stretched state, mm	Not applicable	Not applicable	Not applicable	190±5
Crimped tie width, in stretched state, mm	Not applicable	Not applicable	Not applicable	40±5
Crimped tie elongation, %, no more	Not applicable	Not applicable	Not applicable	152
Band elongation, %, no more than	238	233	Not applicable	Not applicable
Weight, g	30±5	40±5	45±5	60±5

8.2 Adhesive layer of the fixing elements should ensure stable adhesion to the backing layer when fixing the diapers at least 5 times. The product is disposable, however, refixation may be required to fix it in position and secure to the waist.

8.3 The force to be applied in order to separate the fixing elements and the protective layer should be at least 6.5 N.

8.4 Parameters that ensure the functional purpose of the products should comply with the requirements specified in Table 2.

Table 2.

Name of parameter	Design variant			
	Dr. Skipp pull-up diapers for adults, size M (Medium)	Dr. Skipp pull-up diapers for adults, size L (Large)	Dr. Skipp diapers for adults, size M (Medium)	Dr. Skipp diapers for adults, size L (Large)
Complete moisture absorption, g, no less than	1800	2000	2000	2400
Sorption capacity after centrifugation, g, no less than	440	500	480	540
Resorption, g, no more than	4,4	4,4	4,4	4,4
Absorption rate, cm^3/c , no less than	2,3	2,3	2,3	2,3
Breaking load of cover layer in a dry state in the longitudinal direction, kgf, no less than	4,0	4,0	4,0	4,0
Breaking load of cover layer in a wet state in the longitudinal direction, kgf, no less than	3,0	3,0	3,0	3,0

9. Materials used in the manufacture of products

Detailed information on the composition of the materials from which the devices are manufactured is shown in Table 3.

Table 3.

Product part	Material
Top cover layer	Nonwoven fabric (polypropylene of TPP F79FB brand), manufactured by PetroChina Group, China
Absorbent layer	1. Fluff cellulose of Treated brand, manufactured by GP Cellulose, USA. 2. Superabsorbent (sodium polyacrylate) of GS-301 brand, manufactured by LG Chem, Republic of Korea.
Protective layer	Polyethylene of 7000 F brand, manufactured by PetroChina Group, China
Barrier elements	Nonwoven fabric (polypropylene of TPP F79FB brand), manufactured

	by PetroChina Group, China
Fixing elements (Velcro fasteners)	Single-layer hydrophilic polyurethane film Avery MED 5590U with acrylic adhesive 1959, manufactured by Avery Zweckform, Germany

10. Packaging

Products in the quantity determined by the manufacturer are packed in plastic film bags (consumer packaging) which ensures safety of products during transportation and storage.

Dr.Skipp pull-up diapers for adults, size M (Medium), are packed in plastic film bags, 20 pcs per bag

Dr.Skipp pull-up diapers for adults, size L (Large) are packed in plastic film bags, 20 pcs per bag

Dr.Skipp diapers for adults, size M (Medium), are packed in plastic film bags, 30 pcs per package.

Dr.Skipp diapers for adults, size L (Large) are packed in plastic film bags, 30 pcs per bag.

instruction for use, 1 pc., is enclosed to each consumer package.

Seams in plastic film bags should be welded.

Products of the same group, type, variant of size range, design, execution, made from the same materials, with the same quality indicators, the same manufacture date (month, year) are packed in one bag.

No mechanical damage to packaging that opens access to the product surface is allowed.

Products packed in consumer packaging are packed in a box manufactured of liner cardboard, which contains 5 individual (consumer) packages.

Products of the same type, variant of size range, design, execution, made from the same materials, with the same manufacture date (month, year) are packed in one box.

An example of the product consumer package is shown in Figure 3 and Figure 4.


Figure 3 – Consumer package of the product as illustrated by Dr.Skipp diapers for adults, size L (Large)


Figure 4 – Consumer package of the product as illustrated by Dr.Skipp pull-ip diapers for adults, size L (Large)

11. Labelling


11.1 Labelling should be reliable, verifiable and legible. Labelling is applied to package or label (tag) attached to the package. Labelling is applied in any way (by printing, pressing, stamping) that ensures its clarity, clearness and readability. When using the printed method of labelling, setting off is prohibited.

11.2 Labelling on the consumer package of the products should contain:

- name of the manufacturing company and its legal address;
- address of the product actual manufacture;
- number of registration certificate and date of issue;
- full name of the manufacturer's authorized representative in Russia;
- legal address of the manufacturer's authorized representative in Russia;
- trademark;
- name of the medical device;
- name of the product execution;
- type of the product depending on the purpose (degree of urinary incontinence);
- product sizes by waist/hip circumference;
- recommendations for the product use (in the form of figures or text);
- disposal instructions: the words "Do not throw into the sewer";
- quantity of products per pack, pcs.;

- date of manufacture (month, year);
- manufacturer's expiry date;
- symbols according to EN ISO 15223-1-2012;
- designation according to GOST R 55082-2012.

The corresponding warning symbols (see below) are printed on the consumer package or on the label attached to the consumer package in accordance with the normative document EN ISO 15223-1-2012.

Symbol according to EN ISO 15223-1-2012	Symbol explanation
	- Batch code
	- Date of manufacture
	- Best before
	- Re-use is prohibited
	- Protect from sunlight
	- Protect from humidity
	- Storage temperature range
	- See instructions for use
	- Information about the manufacturer.
	- Information about the manufacturer's representative.

The labelling applied to the transportation packaging should contain information according to item 11.2 with the addition:

- quantity of consumer packages;
- gross weight, kg;
- net weight, kg;
- dimensions of the transportation packaging (length x width x height);

11.3 Label layout for consumer package is provided below.

Dr. Skipp

MEDICAL DISPOSABLE ABSORBENT DEVICES

Design variant: XXXXXXXXX

Quantity of products per pack: XX pcs.

Product sizes by waist/hip circumference: XX - XX cm.

The products can be used both for people suffering from moderate to severe incontinence.

Manufacturer: LINYI 2H HEALTHCARE PRODUCTS CO., LTD

Address: NO.58 LUOQI ROAD, SHENGZHUANG STREET, LUOZHUANG DISTRICT, LINYI CITY, SHANDONG PROVINCE, CHINA

Manufacturing site: LINYI 2H HEALTHCARE PRODUCTS CO., LTD., NO.58 LUOQI ROAD, SHENGZHUANG STREET, LUOZHUANG DISTRICT, LINYI CITY, SHANDONG PROVINCE, CHINA

REP

Authorized representative of the manufacturer in Russia:

RM-plus, Limited Liability Company (RM-plus, LLC)

Address: 19, Nauchny proezd, 8/7 floor/office, Moscow, 117246

Tel.: +7 (495) 642-82-60

e-mail: msk@rmpplus.ru

LOT

XXXXXX

XX.XXXX

XX.XXXX


Method of use:

For diapers:

For pull-up diapers:

Disposal method: It is prohibited to use a medical device that has a damaged individual (consumer) packaging, or expired shelf life; it should be disposed of as class A medical waste (epidemiologically safe waste, compositionally similar to solid household waste) in accordance with SanPiN 2.1.7.2790-10. Used devices should be disposed of as class B medical waste (epidemiologically hazardous waste) in accordance with SanPiN 2.1.7.2790-10.

"Do not throw into the sewer!"

Attention! Keep dry! GOST R 55082-2012

Registration certificate No.

dated

12. Biological compatibility

The biocompatibility testing was performed for the finished products in accordance with ISO 10993-1-2009 – Biological Evaluation of Medical Devices. Part 1: Evaluation and Testing as Part of the Risk Management Process. In general, the biocompatibility review supports the conclusion that the products have an acceptable biocompatibility profile.

The products do not have any cytotoxic, irritating, sensitizing, toxic, or hemolytic effects and are compliant with the requirements applicable to products that have a short-term contact with intact skin. The products are also resistant to biological fluids they are or may be exposed to during their use. Based on the tests, the product does not cause any irritation and complies with the requirements of ISO 10993-1-2009.

13. Test and control methods

The company carries out control in the manufacture of products at various stages, including the control of raw materials for subsequent production of products, control of the production process itself, packaging, as well as control of products after production according to specifications.

1. There are four types of testing methods:

- A. Incoming control.
- B. Production inspection.
- C. Final inspection.
- D. Functional testing, if applicable.

2. All four types of test methods are based on written work instructions and on a sampling plan according to ANSI ASQZ 1.4 standard.

14. Risk analysis

14.1 Introduction

Product risk analysis was carried out by LINYI 2H HEALTHCARE PRODUCTS CO., LTD in accordance with EN ISO 14971: 2012 - Application of Risk Management System to Medical Devices. This analysis was developed to identify potential hazards associated with the application of the device and to assess risks related to any potential hazards, as well as to identify hazard mitigation measures included in the project in order to verify that the risks have been eliminated or mitigated to acceptable levels.

14.2 Risk management principles

The basic principles of risk management applied at LINYI 2H HEALTHCARE PRODUCTS CO., LTD.;

1. The risk management process is an integral part of the Company's management system.
2. All risks that may arise at any stage of the life cycle of a medical device manufactured by the Company are subject to identification, assessment, development of risk mitigation measures, and documentation.
3. The objectives of the risk management process are minimizing the consequences of the realized risks.
4. The risk management process in the Company implies application of a uniform and standardized approach to risk identification, assessment, and management.
5. Managers at all levels are responsible for timely identification of risks, their assessment, development of risk management measures, and for informing all the Company stakeholders about the risks affecting the achievement of the goals set.
6. The Company's staff incentive scheme takes the efficacy of risk management into account.
7. In the process of risk management, the Company observes a reasonable balance between the costs of risk management and the size of possible damage due to the occurrence of a risk event. The Company may not take any additional measures to mitigate a risk if its level is acceptable and the costs of risk management exceed the possible effect.

14.3. Risk management process

LINYI 2H HEALTHCARE PRODUCTS CO., LTD documents the risk management processes, which include:

- risk analysis;

- risk assessment;
- risk management;
- post-production information.

14.4 Objectives of risk management

To implement the risk management strategy, the Company employees at all levels must fulfill and maintain the following risk management objectives:

- All identified risks shall be kept record of to ensure their subsequent assessment.
- All risks shall be analyzed in order to determine the level of their acceptability.
- Appropriate measures shall be taken for the risks that exceed the acceptable levels (high and medium levels).
- Reported risks shall be monitored until their elimination, and re-assessed.
- The following staff categories shall be familiarized with the risk management strategy and objectives (to ensure they understand the strategy and objectives correctly):
 - current Company staff (the strategy shall be explained by the heads of the corresponding units);
 - new employees (during admission to work / induction).

14.5 Risk management group (qualified personnel)

Personnel involved in the risk management process shall have knowledge in the following areas:

- design of diapers/pull-up diapers;
- operating principle of diapers/pull-up diapers;
- procedure for manufacture of diapers/pull-up diapers;
- clinical use of diapers/pull-up diapers.

14.6 Management responsibilities

The Company's management shall ensure its engagement in the risk management process by providing sufficient resources and appointing qualified employees to the relevant positions.

The Company's management shall develop and document a policy establishing risk acceptability criteria. This policy shall ensure that the established criteria are based on applicable national and regional regulatory documents and relevant international standards and shall factor in available information such as the current level of research and technology development and the needs of the stakeholders.

The Company's management shall review the suitability of the risk management process at scheduled intervals to ensure the continued efficacy of the process and shall document all the decisions and actions taken. If a manufacturer has a quality management system in place, this review can be a part of this quality management system analysis.

14.7. Analysis of risk management activities

The analysis of risk management activities is performed as part of quarterly operational meetings. These meetings are held to plan the development sequence. They are used to prioritize development projects for the next quarter, to formulate tasks for global development projects, and to assign performers to specific tasks.

These meetings are attended by the employees authorized and responsible for the analysis of the risk management process or their delegates.

The meetings address all issues related to the risk management process that are submitted for consideration at the date of the meeting. The agenda is prepared as part of the approved general procedure for considering issues to be included in the agenda.

14.8 Risk acceptability criteria

The risk acceptability criteria are based on the Company's risk acceptability policy, including cases where the likelihood of harm cannot be determined.

14.8.1 Determination of harm likelihood

The following classification of likelihood levels by qualitative parameters is used:

Ranking	General characteristics	Likelihood range
4	Rare	Up to 100 cases per 100,000
3	Remote	Up to 50 cases per 100,000
2	Unlikely	Up to 10 cases per 100,000
1	Implausible	Up to 1 case per 100,000

14.8.2 Determination of risk severity

The following classification of risk severity by qualitative parameters is used:

Class	General characteristics	Description
4	Absolute	Lethal outcome
3	Critical	Causes traumas or injuries requiring professional medical attention (such as an electric shock)
2	Maximum permissible	Causes temporary traumas or injuries that do not require professional medical attention
1	Can be disregarded	Inconvenience, temporary discomfort, slowdown of work, distortion of information

14.8.4 Severity determination standard


The following matrix is used to determine the severity:

Likelihood	Rare [4]	4 (ACC)	8 (ALARP)	12 (N/ACC)	16 (N/ACC)
	Remote [3]	3 (ACC)	6 (ALARP)	9 (ALARP)	12 (N/ACC)
	Unlikely [2]	2 (ACC)	4 (ACC)	6 (ALARP)	8 (ALARP)
	Implausible [1]	1 (ACC)	2 (ACC)	3 (ACC)	4 (ACC)
Total		Can be disregarded (1)	Maximum permissible (2)	Critical (3)	Absolute (4)
		Severity			

1. Risks marked as ACC do not require further risk control.
2. Risks marked as ALARP (As Low as Reasonably Practicable) are undesirable and can be mitigated by implementation of risk control measures. Risks in this area shall be mitigated according to the ASAP principle without regard to the economic costs.
3. Risks marked as N/ACC require urgent elimination.

14.9. Risk analysis

14.9.1 Risk management procedure

14.9.2 Intended use and determination of characteristics associated with the safety of the medical device (Table 4)

Table 4.

Item No.	Questions according to GOST ISO 14971:2012	A or N/A (applicable / not applicable)	Comments
1.	What is the intended use and how is the medical device to be used?		
1.1.	Purpose of the medical device: For diagnosis, prevention, monitoring, treatment, or alleviation of diseases; For compensation of an injury or a physical disability; For replacement or modification of body parts or control of conception.	A	Medical products are designed to absorb and retain urine when caring for adult patients suffering from moderate to severe incontinence, bedridden patients.
1.2.	Indications for use		
1.2.1.	Anticipated groups of users	A	Adult patients suffering from urinary incontinence.
1.3.	Is the medical device intended for preserving or maintaining life?	N/A	
1.4.	Is specialist involvement necessary in case of the medical device failure?	N/A	
2.	Is the medical device intended for implantation?	N/A	
3.	Is the medical device intended to be in contact with the patient or other persons? The factors to be considered include the following: the character of the possible contact (surface-to-surface contact, invasive contact, or implantation) and in each case, the duration and frequency of the contact	A	Short-term contact with intact skin
4.	What materials or components are included in the medical device, are used with or come into contact with it? The factors to be considered include the following: - Compatibility with the substances concerned; - Compatibility with tissues and biological fluids; - Safety-related characteristics; - Any materials of animal origin included in the medical device.	A	The materials used in the manufacture of the medical device are specified in item 9. The medical device does not contain any derivatives of animal origin. The medical device does not contain any derivatives of human blood or its components.
5.	Can the medical device pass energy to the patient or cause the patient to generate	N/A	

	energy? The factors to be considered include the following: - Type of transmitted energy - Energy management, quality, quantity, intensity, and duration of exposure - Is the level of energy higher than that in the currently used devices?		
6.	Are any substances administered to the patient or removed from him/her? The factors to be considered include the following: - Information on substance administration or removal - Is it one substance or several substances? - Information on the maximum and minimum substance administration / removal rate and on the management of this process.	N/A	
7.	Is the medical device used to process biological substances for their subsequent use, transfusion, or transplantation? The factors to be considered include the following: type of processing and processed substance(s) (e.g. autotransfusion, dialysis, treatment of blood components, or cell-based therapy).	N/A	
8.	Is the supplied medical device sterile, or is it intended to be sterilized by the user, or is another microbiological treatment required? The factors to be considered include the following: - if the medical device is intended for single or repeated use; - storage conditions; - limitation of the number of repeated uses; - methods of sterilization; - exposure to sterilization methods not stipulated by the manufacturer.	N/A	The medical device is a non-sterile and not subject to sterilization medical device intended for single use.
9.	Is the medical device intended for routine cleaning and disinfection by the user? The factors to be considered include the following: types of cleaning and disinfecting agents used, and any limitations of the number of cleaning cycles. The design of the medical device can affect the efficacy of routine cleaning and disinfection. Moreover, effect of cleaning	N/A	Diapers and pull-up diapers are single-use medical devices. Not subject to cleaning and disinfection.

	agents and disinfectants on the safety or performance of the medical device should be considered.		
10.	Does the medical device affect the patient's environment? The factors to be considered include the following: - temperature; - humidity; - composition of atmospheric gases; - pressure; - lighting level.	N/A	
11.	Is the medical device intended for measurement? The factors to be considered include the following: measured variables, measurement correctness and accuracy.	N/A	
12.	Is the medical device interpretive? The factors to be considered include the following: ability of the medical device to give conclusion based on the input or application-derived data, used algorithms, and confidence intervals. Special attention should be paid to unintended use of data or algorithms.	N/A	
13.	Is the medical device intended for use with other medical devices, medicines, or other medical products? The factors to be considered include the following: identification of any other medical devices, medicines, or medical products used and any possible issues associated with their interaction; patient's reaction to treatment.	N/A	The medical device does not contain any medicinal drugs. Diapers and pull-up diapers are used independently of other medical devices.
14.	Is there an unwanted release of energy or substances? Energy-related factors to be considered include the following: noise and vibration, heat, radiation (including ionizing, non-ionizing, and ultraviolet / visible / infrared radiation), contact surface temperature, leakage currents, electric or magnetic fields. Substance-related factors to be considered include the following: substances used in the manufacture, cleaning, and testing that can have an adverse physiological effect if they remain on the device. Other substance-related factors to be considered are excretion of chemical	N/A	

	substances, metabolic by-products, and biological fluids.		
15	Is the medical device sensitive to environmental exposure? The factors to be considered include the following: manufacturing environment, means of transport, and storage conditions. These include: lighting level, temperature, vibrations, leaks, sensitivity to changes in power supply and cooling, electromagnetic interference.	A	The device can be transported by any means of transport, such as covered railroad wagons, containers, trucks, ship cargo holds, heated sealed aircraft modules, in accordance with the rules applicable to these modes of transport and in compliance with the transportation conditions without subjecting the packaging to the risk of damage. When transporting the products, protection against dust and atmospheric precipitation should be provided. Transportation should be carried out at an ambient temperature from -50°C to +50°C and relative humidity up to 100%. Protect from humidity!
16	Does the medical device affect the environment? The factors to be considered include the following: - impact on power supply and cooling facilities; - release of toxic substances; - generation of electromagnetic interferences.	N/A	
17	Are there any consumables or accessories associated with the medical device? The factors to be considered include the following: technical requirements for consumables or accessories associated with the medical device, and any restrictions for users concerning the choice of these materials or accessories.	N/A	
18	Does the medical device require maintenance or calibration? The factors to be considered include the following information: - whether maintenance or calibration should be performed by the operator, the user, or a specialist; - whether proper maintenance or calibration require special substances or equipment.	N/A	Diapers and pull-up diapers are single-use medical devices. Diapers and pull-up diapers have a completed design. Not subject to maintenance and repair.
19	Does the medical device include a software? The factors to be considered include the	N/A	

	following: whether a software needs to be installed, verified, modified, or replaced by the operator, the user, or a specialist.		
20	Does the medical device have a limited shelf life? The factors to be considered include the following: relevant labelling and indicators, and information on the disposal of the medical device after its shelf life expires.	A	The product shelf life is 3 years from the date of manufacture. Disposal information is indicated on the label
21	Are there any delayed or long-term consequences of using the medical device? The factors to be considered include the following: ergonomic and cumulative effects. These include corrosive impact on saline infusion pumps, mechanical fatigue, loosening of assemblies or fasteners, exposure to vibration, labels abrasion or peeling, degradation of material properties due to the long-term use of the medical device.	N/A	
22	What mechanical forces can affect the medical device? The factors to be considered include the following: if the mechanical forces affecting the medical device are managed only by the user or by the user and other persons.	N/A	
23	What factors determine the service life of the medical device? The factors to be considered include the following: device ageing and battery depletion.	A	Natural ageing
24	Is the medical device intended for single use? The factors to be considered include the following: if the medical device can self-destruct after use, if there are any indications that the device has already been used.	A	The medical device is a non-sterile and not subject to sterilization medical device intended for single use. The product is disposable and is intended for single patient use. The recommended time for using the product is no more than 8 hours.
25	Does the medical device require safe dismantling or disposal? The factors to be considered include the following: waste resulting from the disposal of the medical device. For example, it should be noted whether the device contains any toxic, hazardous, or recyclable substances.	A	The devices are disposed of in accordance with the rules for disposal of medical devices in force in the country of use at the date of disposal. It is prohibited to use the medical device that has a damaged individual (consumer) packaging, or the expired shelf life; it should be disposed of as class A waste - epidemiologically safe waste, compositionally similar to solid household waste.

			After use, the device disposal shall correspond to disposal of disposable medical devices classified as class B hazardous medical waste.
26	Does the installation or use of the medical device require special training or special skills? The factors to be considered include the following: the novelty of the medical device; skills and knowledge required for installation of the device.	N/A	
27	What information on the safe use will be provided to the user? The factors to be considered include the following: - whether the information will be provided to the end user directly by the manufacturer or other persons such as assemblers / installers, healthcare providers, healthcare professionals, or pharmacists will be involved and how this will affect training; - information on the pre-commissioning (commissioning) and transfer of the device to the end user, and whether it is reasonable / possible to perform installation by persons without special skills; - whether it is necessary to re-train or re-qualify the operators or the maintenance staff depending on the expected service life of the device.	A	The instruction for use contains information on the safe use of the device.
28	Is the development and implementation of new manufacturing processes necessary? The factors to be considered include the following: new technologies or new scale of production.	N/A	
29	Does successful use of the medical device depend on a human factor such as the user interface?	N/A	
29.1	Could user interface design features lead to application errors? The factors to be considered include the following: user interface design features that could lead to application errors. These include controls or indicators, symbols used, ergonomics, product design and layout of controls, hierarchy of actions, menus for programmed devices, clear warnings,	N/A	

	audible alarms, standard color codes.		
29.2	Is the medical device used in an environment where distractions could lead to misuse? The factors to be considered include the following: - misuse consequences; - whether the distractions are familiar to the user; - the likelihood of rare distractions.	N/A	
29.3	Does the medical device have fittings or accessories? The factors to be considered include the following: whether the fittings could be connected incorrectly; whether there is similarity to fittings of other devices; how strong the connection is; whether there is any feedback in case of incorrect, too strong, or not strong enough connection.	N/A	
29.4	Does the medical device have a control interface? The factors to be considered include the following: layout of controls, their coding, grouping, layout scheme, feedback modes, gross errors, minor slips, control differentiation, image quality, direction of activation or change, discreteness or continuity of control element operation, reversibility of settings or actions.	N/A	
29.5	Does the medical device have an information display function? The factors to be considered include the following: image quality under different conditions; orientation, visual ability of the user; filling extent of the monitor screen and perspective image; the clarity of the information presented; information display units; color codes; access to the most important information.	N/A	
29.6	Does the medical device have a control menu? The factors to be considered include the following: the complexity and the number of menu levels, the user awareness of the menu state, location of the settings, transfer mode between objects, number of steps per action, issues related to clear sequence of actions and memorization, the relation between the importance and availability of a	N/A	

	control function, and the impact of deviations from the established operating procedures.		
29.7	Will the medical device be used by people with special needs? The factors to be considered include the following: intended user, their mental and physical abilities, skills and training, ergonomics, environment of application, installation / assembly requirements, patient's ability to control or influence the use of the medical device. Particular attention should be paid to users with specific needs, such as persons with disabilities, the elderly, and children. Specific needs can also include the need for assistance from another person while using the medical device. It should be considered whether the medical device can be used by persons with different levels of skills and different cultural backgrounds.	A	Adult patients suffering from urinary incontinence. There is often a need for assistance from another person while using the medical device. The necessary information on the correct method of use is shown in the instructions and on the label, including by means of symbols.
29.8	Can the user interface be used to initiate user actions? The factors to be considered include the following: whether it is possible to initiate intended actions of the user to enter the controlled operation mode, which increases risks for the patient; whether the user is aware of what to do under such conditions.	N/A	
30	Does the medical device use an alarm system? The factors to be considered include the following: the risk of false alarms, malfunctioning of the alarm system, disconnection of the alarm system, unreliable performance of remote alarm systems, and whether the medical personnel is capable of understanding how the alarm system works.	N/A	
31	What are the possible misuses of the medical device? The factors to be considered include the following: improper use of the switches; failure of the alarm system or of the safety elements; non-compliance with the manufacturer's recommendations for maintenance.	N/A	
32	Does the medical device store data important to patient management? The factors to be considered include the	N/A	

	consequences of data change or distortion.		
33	Is the medical device mobile or portable? The factors to be considered include the following: necessary clamps, handles, wheels, brakes, mechanical stability and strength of the medical device.	N/A	
34	Does the use of the medical device depend on its critical performance characteristics? The factors to be considered include the following: output characteristics of life support devices; triggering of the alarm system.	N/A	

14.9.3 Hazard identification and risk assessment

Source of hazard / parameter	Severity	Likelihood	Risk number
Incomplete instruction for use	1	3	RM-01
Inadequate description of performance characteristics	1	2	RM-02
Inadequate description of intended use	1	2	RM-03
Inadequate information about limitations	1	2	RM-04
Violation of the production process leading to a non-conformance of the products to the stated characteristics	2	2	RM-05
Error in the configuration	1	3	RM-06
Packaging and labelling error	1	3	RM-07
Non-observance of the storage / transportation conditions	1	4	RM-08
Repeated use	2	2	RM-09

14.9.4 Risk assessment results:

Likelihood	Rare [4]	4 (ACC):1	8 (ALARP)	12 (N/ACC)	16 (N/ACC)
	Remote [3]	3 (ACC):3	6 (ALARP)	9 (ALARP)	12 (N/ACC)
	Unlikely [2]	2 (ACC): 3	4 (ACC):2	6 (ALARP)	8 (ALARP)
	Implausible [1]	1 (ACC)	2 (ACC)	3 (ACC)	4 (ACC)
Total = 8		Can be disregarded (1)	Maximum permissible (2)	Critical (3)	Absolute (4)
		Severity			

Area 1 (ACC): 9 cases out of 9.

Objectives regarding risk control measures: 8 cases in total.

Acceptable	ALARP (as low as reasonably practicable)	N/ACC (not acceptable)
9	0	0

Thus, all the risks (for all the identified hazardous situations) are acceptable (ACC). The residual risk can be disregarded. There is no need for risk mitigation according to the criteria established in the Risk Management Plan.

14.10 Post-production support information

Prior to the commercial release of the masks, LINYI 2H HEALTHCARE PRODUCTS CO., LTD., had followed the internal quality control procedures such as customer complaint management, customer satisfaction procedures, and procedures for error correction and prevention at the production and post-production stages. The information is assessed for potential importance for safety, with particular attention being paid to the following:

- if there are any unknown sources of hazard or any hazardous situations unaccounted for.
- if the certain level of risk associated with a specific hazardous situation is acceptable.

If any of the above safety issues becomes urgent, its impact on the previously implemented risk management measures shall be assessed. The corresponding changes shall be also made to the risk management documentation.

14.11. Activities related to information collection and analysis

The methods of collecting post-production information are part of the established quality management procedures. The information collection procedures include collecting data from the following sources:

- Notifications from the employees involved in the production process, in any convenient format (by mail, face to face communication, etc.).
- Notifications from the customers received via the corporate phone number of the manufacturing company or by e-mail.
- Notifications recorded in the unified notification recording system by the manufacturing company employees who are in charge of agenda preparation for the operational meetings.

Based on the risk analysis, a documented decision (protocol) is made concerning the monitoring method suitable for collecting post-production information about the device. The protocol includes the full name of the designated employee (an expert or a person with the approval authority), the monitoring method (whether feedback or active monitoring is necessary).

14.12. Verification activities

The Risk Management Report specifies that the cumulative risk for all the identified hazardous situations is acceptable. Therefore, it is planned to perform one procedure to verify the efficacy of the risk mitigation measures.

The procedure is aimed to ensure that the final project has taken the appropriate risk management measures. The verification procedure is documented as a protocol and includes the following: the full name of the designated employee (an expert or a person with the approval authority) and the conclusion on the efficacy of the implemented measures.

14.13. Risk management report

Prior to introduction of the medical device, the manufacturing company LINYI 2H HEALTHCARE PRODUCTS CO., LTD., performed an analysis of the risk management process. This analysis demonstrates that:

1. The cumulative risk for all the identified hazardous situations is acceptable.
2. Appropriate methods are used to collect the necessary production and post-production information. The Company has created, documented, and maintains a system for collection and analysis of information regarding the product at all the production and post-production stages.


14.14. Conclusion

Planned risk management of the medical device was carried out: **Medical disposable absorbent devices**. Based on the analysis results, the cumulative residual risk is acceptable.

15. List of applicable standards, applicable directives and harmonized standards

Number of standard	Name of standard
93/42/EEC (2007/47/EC as amended)	Directive on medical devices
EN ISO13485:2016	Medical devices. Quality management systems. Requirements for regulatory purposes.
ISO 11948-1-1996	Diapers for adults. Part 1. Testing the entire product.
EN ISO 14971: 2012	Medical devices. Application of risk management to medical devices
ISO10993-1:2009	Biological evaluation of medical devices. Part 1. Evaluation and testing in the process of risk management.
ISO 15223-1:2012	Medical devices. Symbols used on medical equipment labels during labelling. Part 1. General requirements
EN 1041:2008	Information provided by the manufacturer of medical devices

16. Manufacture process

17. After-sales supervision procedure

The company carries out control in manufacture of products at various stages, including the control of raw materials for subsequent manufacture of products, control of the production process itself, packaging, as well as control of products after production according to specifications. Complaints and indicators of possible unfavorable incidents are studied in detail for all company products, which is an important tool to assess the application of the product after release. The frequency of ongoing inspections will be determined based on the risk management principles. Standard Operation Procedure - product risk assessment procedure - describes the process of conducting periodic product risk reviews to confirm the constant safety and efficiency of the product on the market.

18. Clinical evaluation

The Clinical evaluation document contains the evaluation and analysis of clinical data related to the medical disposable absorbent devices in order to confirm the clinical safety and technical characteristics of the device under the conditions of use provided for by LINYI 2H HEALTHCARE PRODUCTS CO., LTD.

19. Terms of transportation

The device can be transported by any means of transport, such as covered railroad wagons, containers, trucks, ship cargo holds, heated sealed aircraft modules, in accordance with the rules applicable to these modes of transport and in compliance with the transportation conditions without subjecting the packaging to the risk of damage.

When transporting the products, protection against dust and atmospheric precipitation should be provided.

Transportation should be carried out at an ambient temperature from - 50°C to + 50°C and relative humidity up to 98%.

Protect from direct humidity!

20. Conditions of storage

The products should be stored in the shipping packaging in closed, clean, and dry warehouses protected from moisture and direct sunlight, at a distance of at least 1 m from heating devices and placed on racks at least 5 cm away from the floor under conditions that exclude exposure to corrosive media (acidic, alkaline, etc.), highly-flammable and combustible fluids.

Storage should be carried out in the shipping packaging in the manufacturer's warehouses at an ambient temperature from + 5°C to + 40°C and relative humidity up to 98%.

Protect from direct humidity!

21. Shelf life of the product

The product shelf life is 3 years from the date of manufacture.

22. Requirements for environmental protection

The company LINYI 2H HEALTHCARE PRODUCTS CO., LTD. carries out official risk assessments and analyses (in accordance with ISO 14971 standards) of the products at the design and development stage and before starting production. This procedure allows to recognize any possible risks for the user and the environment. The products are also closely monitored for any harmful effects after they are sold. The products do not have any harmful impact on the environment if the conditions of

transportation, storage and use are observed. During the manufacture of the products, environmental hazard by production wastes is excluded. Wastes generated during the manufacture of the products is recyclable. Wastes that are not suitable for recycling should be disposed of according to the established procedure and sanitary regulations.

23. Requirements for disposal

The devices are disposed of in accordance with the rules for disposal of medical devices in force in the country of use at the date of disposal.

It is prohibited to use a medical device that has a damaged individual (consumer) packaging, or expired shelf life; it should be disposed of as class A medical waste (epidemiologically safe waste, compositionally similar to solid household waste) in accordance with SanPiN 2.1.7.2790-10.

Used devices should be disposed of as class B medical waste (epidemiologically hazardous waste) in accordance with SanPiN 2.1.7.2790-10.

24. Cleaning and disinfection

The products are single-use medical devices. Not subject to cleaning and disinfection.

25. Maintenance and repair

The products are single-use medical devices. The products have a completed design. Not subject to maintenance and repair.

26. Warranty obligations

The manufacturer guarantees the quality, safety and characteristics of the disposable medical absorbent devices if they are used in accordance with the Instruction for use before the expiration date and are stored according to the conditions of transportation and storage established by this regulatory document. Mechanical damage to the products cannot constitute the basis for any claims to the manufacturer.

The product shelf life is 3 years from the date of manufacture.

If you have any questions concerning the quality of the medical device in Russia, please contact the authorized representative of the manufacturer:

Authorized representative of the manufacturer in Russia:

Name: RM-plus, Limited Liability Company (RM-plus, LLC);

Place of business address: 19, Nauchny proezd, 8/7 floor/office, Moscow, 117246

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China Council for the Promotion of International Trade is China Chamber of International Commerce

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查询网址 Website for verifying the certificate: <http://www.rzccpit.com/validate.html>

«I CERTIFY»

Director
Wang Kefeng

LINYI 2H HEALTHCARE PRODUCTS CO.,
LTD

Date «01» June 2020 г.

Signature Wang Kefeng

INSTRUCTION FOR USE
for the medical device
Medical disposable absorbent devices

2020

1. PRODUCT NAME:

Medical disposable absorbent devices.

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Table 1.

Name of parameter	Design variant			
	Dr.Skipp pull-up diapers for adults, size M (Medium)	Dr.Skipp pull-up diapers for adults, size L (Large)	Dr.Skipp diapers for adults, size M (Medium)	Dr.Skipp diapers for adults, size L (Large)
Length, mm	320±10	350±10	750±10	870±10
Middle width, mm	180±10	200±10	290±10	290±10
Width on edges, mm	260±10	300±10	750±10	750±10
Width on edges in stretched state, mm	620±10	700±10	Not applicable	Not applicable
Waist/hip circumference, mm	800-1200	1000-1400	680-1150	840-1500
Absorbent layer length, mm	420±10	450±10	520±10	720±10
Absorbent layer width in middle, mm	80±10	80±10	150±10	170±10
Absorbent layer width on edges, mm	130±10	130±10	235±10	280±10
Quantity of tightening elastic bands on side ruffles in the form of an arc through groin, pcs.	2	2	2	2
Quantity of tightening elastic bands on barriers, pcs.	3	3	3	3
Quantity of tightening elastic bands on band, pcs.	18	18	Not applicable	Not applicable
Fixing element length, mm	Not applicable	Not applicable	80±5	80±5
Fixing element width, mm	Not applicable	Not applicable	26±5	26±5
Crimped tie length, mm	Not applicable	Not applicable	Not applicable	125±5
Crimped tie width, mm	Not applicable	Not applicable	Not applicable	45±5
Crimped tie length, in stretched state, mm	Not applicable	Not applicable	Not applicable	190±5
Crimped tie width, in stretched state, mm	Not applicable	Not applicable	Not applicable	40±5
Crimped tie elongation, %, no more	Not applicable	Not applicable	Not applicable	152
Band elongation, %, no more than	238	233	Not applicable	Not applicable
Weight, g	30±5	40±5	45±5	60±5

10.2 Adhesive layer of the fixing elements should ensure stable adhesion to the backing layer when fixing the diapers at least 5 times. The product is disposable, however, refixation may be required to fix it in position and secure to the waist.

10.3 The force to be applied in order to separate the fixing elements and the protective layer should be at least 6.5 N.

10.4 Parameters that ensure the functional purpose of the products should comply with the requirements specified in Table 2.

Table 2. Name of parameter	Design variant			
	Dr. Skipp pull-up diapers for adults, size M (Medium)	Dr. Skipp pull-up diapers for adults, size L (Large)	Dr. Skipp diapers for adults, size M (Medium)	Dr. Skipp diapers for adults, size L (Large)
Complete moisture absorption, g, no less than	1800	2000	2000	2400
Sorption capacity after centrifugation, g, no less than	440	500	480	540
Resorption, g, no more than	4,4	4,4	4,4	4,4
Absorption rate, cm^3/c , no less than	2,3	2,3	2,3	2,3
Breaking load of cover layer in a dry state in the longitudinal direction, kgf, no less than	4,0	4,0	4,0	4,0
Breaking load of cover layer in a wet state in the longitudinal direction, kgf, no less than	3,0	3,0	3,0	3,0

11. SUPPLY PACKAGE

Dr. Skipp pull-up diapers for adults, size M (Medium), are supplied in consumer packaging, 20 pcs per package.

Dr. Skipp pull-up diapers for adults, size L (Large), are supplied in consumer packaging, 20 pcs per package.

Dr. Skipp diapers for adults, size M (Medium), are supplied in consumer packaging, 30 pcs per package.

Dr. Skipp diapers for adults, size L (Large), are supplied in consumer packaging, 30 pcs per package.

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12. LABELLING


12.1 Labelling should be reliable, verifiable and legible. Labelling is applied to package or label (tag) attached to the package. Labelling is applied in any way (by printing, pressing, stamping) that ensures its clarity, clearness and readability. When using the printed method of labelling, setting off is prohibited.

12.2 Labelling on the consumer package of the products should contain:

- name of the manufacturing company and its legal address;
- address of the product actual manufacture;
- number of registration certificate and date of issue;
- full name of the manufacturer's authorized representative in Russia;
- legal address of the manufacturer's authorized representative in Russia;
- trademark;
- name of the medical device;
- name of the product execution;
- type of the product depending on the purpose (degree of urinary incontinence);
- product sizes by waist/hip circumference;
- recommendations for the product use (in the form of figures or text);
- disposal instructions: the words "Do not throw into the sewer";
- quantity of products per pack, pcs.;
- date of manufacture (month, year);
- manufacturer's expiry date;
- symbols according to EN ISO 15223-1-2012;

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The labelling applied to the transportation packaging should contain information according to item 12.2 with the addition:

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2.3 Label layout for consumer package is provided below.

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Product sizes by waist/hip circumference: XX – XX cm.

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REP Authorized representative of the manufacturer in Russia:

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
Dr. Skipp pull-up diapers for adults, size M (Medium), are packed in plastic film bags, 20 pcs per bag.

Dr. Skipp pull-up diapers for adults, size L (Large) are packed in plastic film bags, 20 pcs per bag.

Dr. Skipp diapers for adults, size M (Medium) are packed in plastic film bags, 30 pcs per package.

Dr. Skipp diapers for adults, size L (Large) are packed in plastic film bags, 30 pcs per bag.

Instruction for use, 1 pc., is enclosed to each consumer package.

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An example of the product consumer package is shown in Figure 1.



Figure 1 - Consumer packaging of the product:
a) Dr.Skipp diapers for adults, size L (Large)
b) Dr.Skipp pull-up diapers for adults, size L (Large)

14. METHOD OF USE


The rules for use of products (represented as pictures) are shown on the consumer packaging. Algorithm for putting diapers on:

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15. INSTRUCTIONS FOR USE

The use of products must be carried out in accordance with this Instructions for use.

Users who have studied the Instructions for Use are allowed to use the products.

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IT IS FORBIDDEN TO USE THE PRODUCT IN THE FOLLOWING CASES:

External defects of the product were found: foreign inclusions, damages of the structure of the material used, tears, folds, delamination and other mechanical damage.

PRECAUTIONS, WARNINGS:

Products should not be stored in highly humid environment, near heating devices and in direct sunlight.

Products should not be stored with laundry detergents, machine oil and other strong smelling substances, as the products may be saturated with foreign odors, even in a sealed package.

Products should not be stored near open flames, as it may deform and discolor the product.

If the skin condition in the area of contact with the product deteriorates and irritation occurs, stop using the product and refer to a doctor.

The products are for external use only.

It is forbidden to use the expired products. These products should be disposed of.

CONDITIONS OF USE: The products are intended for use in healthcare institutions and at home, at environmental temperature of +5 °C to +40 °C and relative air humidity up to 98 %.

16. CLEANING AND DISINFECTION REQUIREMENTS

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18. TERMS OF TRANSPORTATION

The product can be transported by any means of transport, such as covered railroad wagons, containers, trucks, ship cargo holds, heated sealed aircraft modules, in accordance with the rules applicable to these modes of transport and in compliance with the transportation conditions without subjecting the packaging to the risk of damage.

When transporting the products, protection against dust and atmospheric precipitation should be provided. Transportation should be carried out at an ambient temperature from - 50°C to + 50°C and relative humidity up to 98%.

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22. MANUFACTURER INFORMATION

Name: LINYI 2H HEALTHCARE PRODUCTS CO.;

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Phone: +86 (0539) 5635398;

e-mail: sales@ak2h.com.

Address of the manufacturing site of the medical device: LINYI 2H HEALTHCARE PRODUCTS CO., LTD, NO.58 LUOQI ROAD, SHENGZHUANG STREET, LUOZHUANG DISTRICT, LINYI CITY, SHANDONG PROVINCE, CHINA.

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Phone: +7 (495) 642-82-60;

e-mail: msk@rmplus.ru.

LIST OF APPLICABLE BY THE MANUFACTURER STANDARDS, APPLICABLE DIRECTIVES AND HARMONIZED STANDARDS

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