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6743 Southpoint Drive North
Jacksonville, Florida 32216-0980
Tel: 904.296.9600

TO WHOM IT MAY CONCERN

It is hereby certified that the dossier of papers to which this present Certificate is attached, and each page of which bears the stamp of Medtronic Xomed, Inc. is the NIM EMG Endotracheal Tube Manual.

Amina Abdal

Amina Abdal
Regulatory Affairs Specialist
Medtronic Xomed, Inc.

24 Jul 2014

Date

Stamp:


United States of America
State of Florida
County of Duval

Subscribed and sworn to before me this 24th day of July, 2014

Personally known ✓

Or Produced Identification _____

Type of Identification Produced _____


Michelle Hughes

Medtronic Xomed
6743 Southpoint Drive North
Jacksonville, Florida 32216-0980
www.medtronic.com
Tel: 904-296-9600
Fax: 904-296-2386

Medtronic

NIM EMG Endotracheal Tube

Product Information and Instructions

Informations et instructions concernant le produit / Informazioni sul prodotto e istruzioni per l'uso / Produktinformationen und Anweisungen / Instrucciones e información sobre el producto / Productinformatie en instructies / Produktinformation og -vejledning / Tuotetiedot ja käyttöohjeet / Produktinformation och instruktioner / Informações do produto e instruções / Πληροφορίες προϊόντος και οδηγίες / Informacja o produkcie i instrukcje / Informace o výrobku a návod k použití / Termékinformáció és utasítások / Ürün ile ilgili Bilgi ve Talimatlar / Produktinformasjon og instruksjoner / Информация о продукте и инструкция по его эксплуатации / Informacije o proizvodu i uputstvo / Podaci o proizvodu i upute za uporabu

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English

DEVICE DESCRIPTION

The NIM EMG Endotracheal tube is a standard-size PVC (DEHP-free polyvinyl chloride) endotracheal tube with inflatable high-volume, low-pressure cuffs. Each tube is imprinted with four (two pairs) conductive silver electrodes. The conductive silver is printed on the PVC of the main shaft of the endotracheal tube and exposed only for a short distance, approximately 40 mm, slightly superior to the cuff, for contacting the vocal cords. The electrodes are designed to make contact with the patient's vocal cords to facilitate electromyographic (EMG) monitoring of the laryngeal musculature during surgery when connected to a multi-channel EMG neuromonitoring device.

INTENDED USE

The EMG Endotracheal tube is intended for use as a means of providing both an open airway for patient ventilation and for intraoperative monitoring of EMG activity of the laryngeal musculature when connected to an appropriate EMG monitor.

INDICATIONS FOR USE

The EMG Endotracheal tube is indicated for use where continuous monitoring of the nerves supplying the laryngeal musculature is required during surgical procedures.

CONTRAINDICATIONS

- The EMG Endotracheal tube is not intended for postoperative use. The EMG Endotracheal tube should be removed and replaced with a standard endotracheal tube if ventilation is required beyond the surgical procedure.
- The use of paralyzing agents will significantly reduce, if not completely eliminate EMG responses to direct or passive neural stimulation.
- The EMG Endotracheal tube is not intended for use in infant and neonatal patients.
- Do not use the EMG Endotracheal tube in surgical procedures with flammable gases or high temperature generating devices, such as surgical lasers or cauteries inside the airway, to avoid combustion.

WARNINGS

- Do not use local anesthetic gels or creams to lubricate the tube or apply topical anesthetic sprays on the vocal cords. The use of neuromuscular blocking agents will significantly reduce, if not completely eliminate, EMG responses to direct or passive neural stimulation. Whenever nerve paralysis is suspected, consult the anesthesiologist.
- Do not attempt to use an EMG Endotracheal tube without first performing a cuff inflation test. Inflate the cuff with air and visually inspect the cuff for any abnormality. Replace with another EMG Endotracheal tube if any adverse abnormality is found.
- Failure to use a NIM stimulating device may result in not identifying the location or detecting injury to the nerve. During EMG monitoring only, (passive monitoring), neuromuscular tissue may not elicit a response from sharp dissection or all types of manipulation. Use the stimulating device to confirm location and integrity of the nerve.
- Avoid false negative responses (failure to evoke a nerve response), which may be caused by:
 - Neuromuscular fatigue from prolonged or repeated exposure to electrical stimuli.
 - Deep anesthesia, which may suppress neuroelectrical activity of the recurrent laryngeal nerve.
 - Insufficient contact between the electrodes and vocal cords due to misplacement or using a tube size that is too small.
 - Displacement during the procedure when rotating, flexing, or extending the patient's head and neck.
- Do not use an EMG Endotracheal tube with a high electrode impedance. This may be due to inadequate electrode continuity. Replace with another EMG Endotracheal tube.
- Cuff deflation or diffusion of a gas or gas mixture (nitrous oxide, oxygen, air) might result in cuff pressure and volume variation. Monitor the cuff pressure during the procedure.
- Do not over inflate the cuff. Over inflation may cause cuff deformation, deflation or rupture resulting in tube blockage and/or tracheal damage.
- Avoid damage to the trachea, larynx or vocal cords during insertion by inspecting the tube, cuff, and monitoring connections for damage which may occur if the cuff is over inflated during set-up testing.
- Do not attempt to manipulate an EMG Endotracheal tube with an inflated cuff after insertion. Manipulating a tube with an inflated cuff may cause partial airway blockage at the tip and/or Murphy eye, cuff herniation or tip deflection. Ensure that the cuff is fully deflated before any manipulation and confirm that the airway is free of any potential occlusion after repositioning.
- Do not attempt to remove an inflated cuff from the trachea. This may result in cuff deformation and/or injury of the larynx or vocal cords. If resistance is encountered during removal, ensure the cuff is fully deflated.
- Avoid insertion of a suction tube or stylet in a tube that has been distorted by biting or other forces. This may further damage the tube and block the airway.
- Do not attempt to straighten bent needles, which may weaken the metal. Improperly placed, weak or bent needles increase the risk of the needle breaking off in the patient.
- Do not excessively bend or flex the electrodes or the electrical integrity may be compromised. Check electrode integrity after insertion.
- Avoid inadequate oxygenation and/or tracheobronchial injury, which may be caused by:
 - Excessively bending the tube causing it to kink or retracting the tube which may cause it to crush.
 - Placement in the right main bronchus (deep intubation) by confirming the correct positioning after intubation.
 - A collapsed or blocked tube after performing a test inflation, by inspecting the inner diameter of tube for patency after test inflation.
- Failure to use a bite block may result in patient injury or damage to the tube. The use of a bite block is essential during transcranial Motor Evoked Potentials (TcMEP) procedures.
- To avoid disrupted air supply, confirm, monitor, and maintain anesthesia gas delivery systems and connections at all times.

NIM EMG Endotracheal Tube

Medtronic Xomed
 6743 Southpoint Drive North
 Jacksonville, Florida 32216-0001
 www.medtronic.com
 Tel: 904-296-9600
 Fax: 904-296-2386

English

- Do not perform Magnetic Resonance Imaging (MRI) on a patient with an EMG Endotracheal tube in place.
- Do not activate the electrosurgical instruments while the stimulator is in contact with tissue or in the surgical field to avoid patient burns.
- Do not connect electrode leads to any other types of equipment or connections than EMG recording/monitoring equipment or electrical shock to the patient or user may occur.
- Avoid delivery of current to unintended tissues by using an EMG monitor equipped with active audio and/or visual feedback systems for stimulating.
- Do not use if sterile package has been opened or is damaged.
- Do not attempt to re-sterilize or reuse single use disposable devices. Re-sterilization or reuse may result in device damage and patient injury.
- Avoid injury by disposing of the devices in an appropriate sharps biohazard container.
- Avoid injury by using extreme care when handling and cleaning instruments with sharp points or edges.

PRECAUTIONS

- It is strongly recommended that the clinician(s) have a thorough understanding of and experience with intraoperative EMG Monitoring prior to using the EMG Endotracheal tube in a surgical procedure. This device is not intended to replace the surgeon's medical judgment or knowledge of neural anatomy and physiology. It is intended to provide the surgeon with an additional tool with which to make better-informed decisions regarding the surgical procedure.
- It is strongly recommended that the surgeon consult with the attending licensed medical practitioner who will be administering anesthesia prior to the use of EMG monitoring to review EMG monitoring techniques, goals and the effects of the administration of anesthesia on neuromuscular activity.
- A patient with an implanted electronic device should not be subjected to electrical stimulation unless specialist medical opinion has first been evaluated.
- Proper sizing, intubation and extubation should be in accordance with accepted medical techniques and expert clinical judgment. A tube that is one size larger than standard selection is recommended whenever possible to improve electrode contact with vocal cords. The proper size tube for the patient should be determined prior to intubation by the anesthesia provider and/or surgeon.
- For safe and accurate EMG monitoring, proper handling, insertion and placement of electrodes and probes is critical. Review the appropriate user's guide for the monitor being used for the procedure.

STERILITY

- The needle electrodes are provided sterile and are intended for single patient use only.
- Do not resterilize this product. Medtronic assumes no liability for products which have been resterilized.

INSTRUCTIONS FOR USE

Please read and understand all Warnings and Precautions before using this device.

The NIM nerve integrity monitor is the recommended EMG monitor for use with the EMG Endotracheal tube. All references to connections and technical specifications for EMG monitors contained in these instructions are made with references to the NIM. Medtronic recommends that the user consult the NIM User's Guide for determining proper settings of the NIM and instructions for intraoperative EMG monitoring and motor nerve location and stimulation.

Preparation

1. Determine the proper size tube for the individual patient prior to intubation by the anesthesia provider and/or surgeon. A tube that is one size larger than standard selection is recommended to improve electrode contact with vocal cords.
2. Test the cuff for leakage with 15cc to 20cc of air. Thoroughly evacuate all air before intubation.
3. Lubricate cuff with an aqueous lubricant for intubation. Avoid local anesthetic gas and solutions which may impair monitoring.
4. Prior to insertion of the EMG Endotracheal tube, visually inspect the larynx to ensure that the tube can be inserted to allow adequate electrode contact with the patient's vocal cords.

EMG Endotracheal Tube Placement

1. Intubate the patient per the procedure for non-reinforced endotracheal tubes. Confirm that the tube lumen is free of any potential occlusion. Place the tube so that the marking is anterior, with red wire on the right and blue wire on the left.
2. Visualize electrode contact with the true vocal cords. A blue cross-hair is printed, anteriorly, on the tube to aid placement. The tube is optimally placed with the vocal folds in contact with the wide blue cross-hair section.

Note: The optimum exposure can be obtained at the midpoint, with the midpoint of the cross-hair approximately being 80-106mm above the caudal tip of the tube.

3. Confirm the depth of intubation per the procedure for non-reinforced endotracheal tubes. The EMG Endotracheal tube has depth markings on the surface of the tube. Electrode depth and location should be checked against preoperative inspection using these markings.
4. Use the EMG monitor to measure electrode impedance and imbalance. Suggested impedance values for the EMG Endotracheal tube are less than 5 kOhms. Impedance imbalance values of less than 2 kOhms are recommended between the positive and negative electrodes of a channel. Reposition if necessary.

2

NIM EMG Endotracheal Tube

Medtronic Xomed
6743 Southpoint Drive North
Jacksonville, Florida 32216-0001
www.medtronic.com
Tel: 904-296-9600
Fax: 904-296-2006

Note: Although monitoring can proceed, one is likely to experience problems with electrical and bioelectric artifact and excessive background noise with large recording electrode imbalance levels. It is better to have a pair of recording electrodes that have high impedance values, but are relatively balanced, than to have one electrode with very low impedance and the other with very high impedance.

5. Secure the tube once positioned properly in the trachea.

Note: A bite block is recommended for use with the EMG Endotracheal tube to prevent damage to the tube.

6. Inflate the endotracheal tube cuff with a minimum volume of air or saline injected with a syringe through the cuff inflation valve, to create an airtight seal and ensure against slippage of the tube in situ.
7. Monitor the cuff volume routinely to ensure an adequate seal is maintained.

Connect EMG Tubes to Monitor

1. Connect the electrodes of the EMG Endotracheal tube to the NIM in accordance with the NIM User's Guide. The proper connection of the EMG Endotracheal tube to the NIM system is essential for accurately recording the EMG activity of the vocal cords.
2. Insert the green recording ground electrode into the sternum, and then insert the white electrode with red plug into the sternum a few centimeters below the green electrode. The white electrode is the monopolar stimulator return (anode, +) electrode.
3. Plug the connectors into the corresponding locations in the patient interface.
4. Press **ELECTRODE CHECK** on the NIM monitor to verify proper impedance.

Patient Monitoring

- After you confirm the position of the EMG Endotracheal tube, note the baseline activity of monitoring activity. Baseline electrical activity may vary from patient to patient and during the course of an individual procedure.

Medtronic recommends that the user consult the NIM User's Guide to determine proper settings and instructions for intraoperative EMG monitoring and motor nerve location and stimulation.

EMG Endotracheal Tube Removal

Fully deflate the cuff prior to extubation. Exercise caution while extubating the patient. Extubate by pulling on the tube, not on the harness.

CUSTOMER SERVICE

For further information regarding the use of this product or to report any problems, please contact Medtronic, Inc. using the appropriate information provided on the blue and white contact information card packaged with each device; or contact your local distributor.

LIMITED WARRANTY

- A. This LIMITED WARRANTY provides assurance for the customer who purchases a Medtronic Xomed Product (hereinafter the "Product") that should the Product fail to function to Medtronic Xomed's published specifications during the term of this LIMITED WARRANTY (one year from the date of shipment for new Product, 90 days from date of shipment for refurbished or used Product), Medtronic Xomed will either replace, repair, or issue a credit (adjusted to reflect the age of the Product) for the Product or any portion thereof. This LIMITED WARRANTY is extended only to the buyer purchasing the Product directly from Medtronic Xomed or from its affiliate or its authorized distributor or representative.
- B. To qualify for this LIMITED WARRANTY, the following conditions must be met:
 - (1) The Product must be used on or before its "Use By" or "Use Before" date, if applicable.
 - (2) The Product must be used in accordance with its labeling and may not be altered or subjected to misuse, abuse, accident or improper handling.
 - (3) Medtronic Xomed must be notified in writing within thirty (30) days following discovery of a defect.
 - (4) The Product must be returned to Medtronic Xomed within thirty (30) days of Medtronic Xomed receiving notice as provided for in (3) above.
 - (5) Upon examination of the Product by Medtronic Xomed, Medtronic Xomed shall have determined that: (i) the Product was not repaired or altered by anyone other than Medtronic Xomed or its authorized representative, (ii) the Product was not operated under conditions other than normal use, and (iii) the prescribed periodic maintenance and services have been performed on the Product.
- C. This LIMITED WARRANTY is limited to its express terms. THIS LIMITED WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED WHETHER STATUTORY OR OTHERWISE, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. In no event shall Medtronic Xomed be liable for any consequential, incidental, prospective or other similar damage resulting from a defect, failure, or malfunction of the Product, whether a claim for such damage is based upon the warranty, contract, negligence or otherwise.
- D. The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. Users may benefit from statutory warranty rights under legislation governing the sale of consumer goods. If any part or term of this LIMITED WARRANTY is held by any court of competent jurisdiction to be illegal, unenforceable, or in conflict with applicable law, the validity of the remaining portion of the LIMITED WARRANTY shall not be affected, and all rights and obligations shall be construed and enforced as if this LIMITED WARRANTY did not contain the particular part or term held to be invalid.

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
REF	Catalog number	STERILEO	Sterilized by ethylene oxide
LOT	Lot number	STERILER	Sterilized by radiation
!	Caution	STERILEM	Medicated using steam or dry heat
⌚	Use by date	ECTREP	Authorized Representative in the European Community
🏭	Date of Manufacture		Temperature Limitation
🏭	Manufacturer		Consult instructions for use
SN	Serial Number		

ECTREP
EU/UKA Authorized Representative: Medtronic Xomed, Inc.
Earl Bakken Inc. 4122 Franklin Ave. Minneapolis, MN 55425
Phone: +1 612 556 5500 Fax: +1 612 556 5500

Germany
phone: +49 211 429 4290
fax: +49 211 520 5100
Italy
phone: +39 22 266 66 66
fax: +39 22 266 66 66
Japan
phone: +81 3 5640 7110
fax: +81 3 5640 7110
Netherlands
phone: +31 43 300 3000
fax: +31 43 300 3000
United Kingdom
phone: +44 192 521 2215
fax: +44 192 521 2215

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Precaution

If the phthalates symbol is on the device label then this device contains phthalates. The European Union has categorized these substances as those that "may" impair fertility or cause harm to the unborn child. Therefore minimum treatment time should be used for children, pregnant or nursing women.

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© Medtronic 2011