
INSTRUCTIONS FOR USE


Spectacle lenses made of polymer material: colorless, tinted, photochromic, afocal, bifocal, transfocal, single-vision, polarizing with antireflex or reinforcing coating or without coating with accessories, execution options


Production: «BBGR», 22 rue de Montmorency. 75003, Paris, France.

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signature de M. 

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1. PURPOSE AND DESIGN OF THE MEDICAL DEVICE.

Spectacle lenses made of polymer material: colorless, colored, photochromic, afocal, bifocal, transfocal, single-vision, polarizing with antireflex or reinforcing coating or without coating with accessories (hereinafter - lenses, product) - an optical product intended for the manufacture of individual corrective glasses intended for the correction of visual defects, incl. various degrees of eye ametropia, accommodation disorders and correction of the spectral composition of light passing through the lens.

All types of lenses in a pair can be combined with afocal lenses (supplied in a pair of lens + afocal lens) for **amlyopia**.

Scope - optics, optometric and ophthalmological offices of hospitals, clinics, health institutions.

Indications

Depending on the lens design:

Afocal lenses are designed for:

- improving the contrast of the image,
- correction of violation / decrease in the function of accommodation,
- correction of violations of the binocular system of the eyes,
- correction of the spectral composition of the light flux entering the eye,
- correction of anisometropia, etc.

Single vision lenses are designed for:

- correcting and / or slowing down the progression of myopia,
- correction of hyperopia, astigmatism, presbyopia,
- correction of violations of the binocular system of the eyes.

Bifocal lenses are designed for:

- correcting and / or slowing down the progression of myopia,
- correction of hyperopia, astigmatism in case of violation / decrease in the function of accommodation,
- correction of presbyopia in the absence of the need for clarity of vision at medium distances.

Trifocal lenses are designed for:

- correction of myopia, hyperopia, astigmatism in case of violation / decrease in the function of accommodation,
- correction of presbyopia with special requirements for clear vision distances.

Transfocal (progressive) lenses are designed for:

- correcting and / or slowing down the progression of myopia,
- correction of hyperopia, astigmatism in case of violation / decrease in the function of accommodation,
- correction of presbyopia.

Contraindications

Absent

Side Effects : The only potential physical hazard is risk

laceration from broken lenses. However, lenses have excellent impact resistance, but they can break with significant physical force. Broken lenses should be disposed of properly.

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Pic.3 Napkin

1.2. How supply

Spectacle lenses made of polymer material: colorless, tinted, photochromic, afocal, bifocal, transfocal, single-vision, polarizing with antireflex or reinforcing coating or without coating with accessories, execution options:

Spectacle lenses with refraction from -40 to +30 diopters, diameter 50-80 mm:

- 1 Spectacle lenses made of polymer material, colorless, single vision
- 2 Spectacle lenses made of polymer material, colorless bifocal
- 3 Colorless transfocal spectacle lenses made of polymer material
- 4 Spectacle lenses made of polymer material colored single vision
- 5 Spectacle lenses made of polymer material, colored bifocal
- 6 Spectacle lenses made of polymer material, colored transfocal
- 7 Spectacle lenses made of polymer material, photochromic single vision
- 8 Spectacle lenses made of polymer material, photochromic bifocal
- 9 Spectacle lenses made of polymer material, photochromic transfocal
- 10 Spectacle lenses made of polymer material, single-vision polarizing
- 11 Spectacle lenses made of polymer material polarizing bifocal
- 12 Spectacle lenses made of polymer material polarizing zoom
- 13 Spectacle lenses made of polymer material, polarizing photochromic single-vision
- 14 Spectacle lenses made of polymer material polarizing photochromic transfocal
- 15 Spectacle lenses made of polymer material polarizing photochromic bifocal
- 16 Spectacle lenses made of polymer material polarized colored single vision
- 17 Spectacle lenses made of polymer material polarized colored transfocal
- 18 Spectacle lenses made of polymer material polarized colored bifocal

Accessories:

Napkin (if necessary)

Spectacle lenses with a refraction of 0 diopters, diameter 50-80 mm:

- 1 Spectacle lenses made of polymer material, colorless afocal
- 2 Spectacle lenses made of polymer material, painted
- 3 Spectacle lenses made of polymer material photochromic afocal
- 4 Spectacle lenses made of polymer material polarizing afocal
- 5 Spectacle lenses made of polymer material polarizing photochromic afocal
- 6 Spectacle lenses made of polymer material polarized colored afocal

Accessories:

Napkin (if necessary)

2. TECHNICAL CHARACTERISTICS

Spectacle lens characteristic name	Spectacle lens type and value of characteristics			
	Afocal	Single vision	Bifocal	Transfocal (progressive)
spherical (S) refraction value, diopters	0	-40 to +30, step 0.25	-40 to +30, step 0.25	-40 to +30, step 0.25
value of cylindrical (C) refraction, diopters	0	0 to 8.0 in 0.25 steps	0 to 8.0 in 0.25 steps	0 to 8.0 in 0.25 steps
position of the axis of the cylinder / base of the prism (°)	0	from 0 to 360, step 1.0	from 0 to 360, step 1.0	from 0 to 360, step 1.0

value of prismatic action (prism), pr. diopters	from 0 to 5, step 0.25	from 0 to 20, step 0.25	from 0 to 20, step 0.25	from 0 to 20, step 0.25
nominal size, diameter, mm	50-80, step 1mm	50-80, step 1mm	50-80, step 1mm	50-80, step 1mm
additional refraction for near, diopters	-	-	from 0.25 to 4.0, step 0.25	-
prismatic action of the lens (prism), pdpt.	from 0 to 5, step 0.25	from 0.25 to 20.0, step 0.25	from 0.25 to 20.0, step 0.25	from 0.25 to 20.0, step 0.25
prismatic action of progressive spectacle lenses, diopters	-	-	-	from 0.25 to 10.0, step 0.25
additional refraction for near progressive spectacle lenses, diopters	-	-	-	from 0.75 to 3.5, step 0.25
information about the thickness of the lenses in the center, mm	from 1 to 5, step 0.1	from 1 to 20, step 0.1	from 1 to 20, step 0.1	from 1 to 20, step 0.1
information about the thickness of the lenses along the edge, mm	from 0.5 to 7, step 0.1	from 0.5 to 18, step 0.1	from 0.5 to 18, step 0.1	from 0.5 to 18, step 0.1

The maximum permissible deviations of the refraction of single-vision lenses and zones for the distance of multifocal lenses from the nominal values should correspond to those given in Table 1.

Table 1

In diopters

Surface refraction at the <i>second</i> major meridian	Maximum deviation on the <i>first</i> major meridian	Maximum deviation of the absolute value of the astigmatic difference (cylinder)			
		from 0.00 to 0.75 incl.	from 0.75 to 4.00 incl.	from 4.00 to 6.00 incl.	over 6.00
0.00 to 3.00 incl.	$\pm 0,09$	$\pm 0,09$	$\pm 0,12$	$\pm 0,18$	-
3.00 to 6.00 incl.	$\pm 0,12$	$\pm 0,12$	$\pm 0,12$	$\pm 0,18$	$\pm 0,25$
6.00 to 9.00 incl.	$\pm 0,12$	$\pm 0,12$	$\pm 0,18$	$\pm 0,18$	$\pm 0,25$
From 9.00 to 12.00 incl.	$\pm 0,18$	$\pm 0,12$	$\pm 0,18$	$\pm 0,25$	$\pm 0,25$
12.00 to 20.00 incl.	$\pm 0,25$	$\pm 0,18$	$\pm 0,25$	$\pm 0,25$	$\pm 0,25$
More than 20.00	$\pm 0,37$	$\pm 0,25$	$\pm 0,25$	$\pm 0,37$	$\pm 0,37$

The maximum permissible deviations of the refraction zone for given progressive spectacle lenses *from the nominal values* should correspond to those given in table 2.

Table 2

In diopters

Surface refraction at the <i>second</i> major meridian	Maximum deviation on	Maximum deviation of the absolute value of the astigmatic difference (cylinder)
--	----------------------	---

	the <i>first</i> major meridian				
		from 0.00 to 0.75 incl.	from 0.75 to 4.00 incl.	from 4.00 to 6.00 incl.	over 6.00
0.00 to 6.00 incl.	±0,12	±0,12	±0,18	±0,18	±0,25
6.00 to 9.00 incl.	±0,18	±0,18	±0,18	±0,18	±0,25
9.00 to 12.00 incl.	±0,18	±0,18	±0,18	±0,25	±0,25
12.00 to 20.00 incl.	±0,25	±0,18	±0,25	±0,25	±0,25
More than 20.00	±0,37	±0,25	±0,25	±0,37	±0,37

The maximum permissible deviations of the position of the cylinder axis from the nominal values (refer to multifocal, progressive and single-vision spectacle lenses with a predetermined orientation, for example, the position of the base of the prism) should correspond to the values shown in Table 3.

Table 3

The absolute value of the astigmatic difference (cylinder), diopters	Less than 0.50	0.50 to 0.75 incl.	0.75 to 1.50 incl.	More than 1.50
Axis deviation limit	±7°	±5°	±3°	±2°

The maximum permissible deviations of the additional refraction of the zone for the vicinity of multifocal and progressive spectacle lenses from the nominal values should correspond to those given in Table 4.

Table 4

In diopters

The value of the additional refraction zone for near	Less than 4.00	More than 4.00
Maximum permissible deviation	±0,12	±0,18

At the base point for given, the maximum permissible deviations of the resulting prismatic action of the prescribed prism and, if present, a thinning prism, should not exceed the values given in Table 5. (For spectacle lenses without a specified prismatic action, the data in Table 5 represent the maximum permissible values of the undesirable prismatic action caused by the deviation of the position of the constructive base point from the calculated one.)

Table 5

In prism diopters

Prismatic action	Lenses		
	Single vision	Multifocal and progressive	
		Horizontally	Vertically
0.00 to 2.00	$\pm(0,25+0,1 F' V_{\max})$	$\pm(0,25+0,1 F' V_{\max})$	$\pm(0,25+0,05 F' V_{\max})$
2.00 to 10.00	$\pm(0,37+0,1 F' V_{\max})$	$\pm(0,37+0,1 F' V_{\max})$	$\pm(0,37+0,05 F' V_{\max})$
More than 10.00	$\pm(0,50+0,1 F' V_{\max})$	$\pm(0,50+0,1 F' V_{\max})$	$\pm(0,50+0,05 F' V_{\max})$

Note - $F' V_{\max}$ - the highest absolute value of refraction at the main meridians.

Note - An example of the application of the zone tolerances indicated in the table for the distance of a multifocal lens with a prescription: +0.50; -2.50; 20 ° with a prismatic action of no more than 2.00 pdptr.

With regard to spectacle lenses without a given prismatic action, the data in Table 5 represent the maximum permissible values of undesirable prismatic action caused by the deviation of the position of the constructive base point from the calculated one.

For this recipe, the refractions on the main meridians are +0.50 diopters and -2.00 diopters.
 With the highest absolute value of refraction equal to 2,00 diopters *deflection of prismatic action* horizontally equals $\pm (0,25+0,1 \times 2,00) = \pm 0,45$ pdpt, *deflection of prismatic action* vertically equals $\pm (0,25+0,05 \times 2,00) = \pm 0,35$ pdpt.

Dimensional tolerances:

dimensional tolerances (for diameter) are calculated by the formula:

effective size: from nominal size -1 to usable size +2 mm;

useful size: more useful size -2mm.

Thickness tolerance (at the reference point of the front surface for distance normal to this surface) : The thickness should not differ from the nominal value by more than ± 0.3 mm.

Segment tolerances for multifocal lenses:

Any of the segment dimensions (width, depth and depth of the intermediate zone) should not differ from the nominal value by more than ± 0.5 mm.

When a matched pair of spectacle lenses is sold, any of the segment dimensions (width, depth and depth of the intermediate zone) must not differ from each other by more than ± 0.7 mm.

The quality of processing of the refractive surfaces of the lens must meet the following requirements:

- The deviation of the shape of the lens surface (waves), distorting the image of the object in question, is not allowed.
- The width of the scratches (b) and their total length (l) must not exceed the values given in table 5:

Table 5

Lens area	b, mm	l, mm
Central, diameter 30mm	Over 0.006 to 0.01 incl.	10
Regional	Over 0.006 to 0.02 incl.	10

Scratches up to 0.006mm and dots with a diameter of 0.05mm are allowed if their area in a limited area with a diameter of 5mm does not exceed 0.1 mm^2 .

- The number (N) of bubbles, dots and other foreign inclusions, permissible only when the distance between them is more than 5 mm, and their diameter (d) should not exceed the values indicated in table 6:

Table 6

Lens area	d, mm	N, no more
Central, diameter 30mm	0.05 to 0.1 incl.	1
Regional	0.1 to 0.2 incl.	2

- In the annular zone 2 mm wide along the edge of the lens, defects, except cracks, are not standardized. Cracks are not allowed.

Lenses must withstand the impact of a steel ball with a diameter of 22 mm, to which a force of $(100 \pm 2) \text{ N}$ is applied

UVB-100% .

Lenses must be made from the following materials:

Material	Index	Abbe ne	Specific gravity	UV cut	Chemical components
1.5 (CR39)**	1.502	58	1.31 g/cm3	355nm	diethylene glycol bis (allyl carbonate)
1.5 Photochromic (CR607)**	1.499	60	1.27 g/cm3	383nm	
1.5 Photochromic (CR39)	1.502	58	1.31 g/cm3	393nm	
1.5 Polarized (CR39)	1.502	58	1.31 g/cm3	380nm	
1.6 (MR8)	1.601	42	1.30 g/cm3	398nm	polythiourethanes produced by a reaction of a polythiol with di-isocyanates
1.6 Photochromic (MR8)	1.601	42	1.30 g/cm3	400nm	
1.6 Polarized (MR8)	1.601	42	1.30 g/cm3	400nm	
1.67 (MR7)	1.665	32	1.36 g/cm3	395nm	
1.67 Photochromic (MR7)	1.665	32	1.36 g/cm3	400nm	episulphides
1.74 (N19)	1.736	31	1.46 g/cm3	400nm	
1.74 Photochromic (N19)	1.736	31	1.46 g/cm3	400nm	polycarbonate
Polycarbonate (PC)	1.591	31	1.20 g/cm3	380nm	
Polycarbonate Photochromic (PC)	1.591	31	1.20 g/cm4	390nm	
Polycarbonate Polarized (PC)	1.591	31	1.20 g/cm5	390nm	

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
The chromaticity of photochromic lenses is from 5% to 85%.
Outer lens coatings consist of:

1.5 UV CLEAR	1,5	59	2.48 g/cm ³	silicon dioxide (SiO ₂), sodium oxide (Na ₂ O) from sodium carbonate (Na ₂ CO ₃), calcium oxide (CaO)
1.5 UNICROWN	1,5	58	2.61 g/cm ³	
1.6/41 White	1,6	42	2.63 g/cm ³	
1.7/35 White	1,7	34	3.21 g/cm ³	
1.7/42 White	1,7	42	3.21 g/cm ³	
1.8/35 White	1,8	34	3.65 g/cm ³	
1.9/31 White	1,9	30	3.99 g/cm ³	
UV DG 37	1,5	60	2.48 g/cm ³	
PHOTOGRAY EXTRA	1,5	57	2.41 g/cm ³	
PHOTOBROWN EXTRA	1,5	46	2.41 g/cm ³	
FUSIBLE PHOTOGRAY 16	1,6	42	2.70 g/cm ³	
XDF DARK GRAY	1,5	56	2.41 g/cm ³	
XDF LIGHT GRAY	1,5	57	2.41 g/cm ³	

Grades of dyes and coatings produced by "BBGR".

Napkin material: - Polyester 70% Texwiel + Polyamide 30% Nylon 6. The napkin is packed in a polyethylene envelope. Polyethylene (PE) brand LUTENE FB3000.

Overall dimensions of the napkin (LxW): 195 × 146 mm ± 2mm

Napkin weight: no more than 5 g

Overall dimensions of a polyethylene envelope (LxW): 97 × 52 mm ± 2mm


Case weight: no more than 5 g

3. CONDITIONS OF USE . OPERATING RULES. CARE RULES.


3.1. Lenses must be suitable for use at an air temperature of -45 to $+40^{\circ}\text{C}$, with a relative humidity of up to 80%.

3.2. Operating rules:

From the lens in an electronic facetting machine (bevelling method), an eyeglass lens is machined into the shape of a spectacle frame.

The method of decorative processing of the edge of the lens by cutting the edge for a given frame.

After that, the frame is straightened in accordance with STB GOST, as well as the patient's anthropological data.

The basis for the manufacture of glasses is a prescription from an ophthalmologist.

An automatic dioptrimeter controls the glasses being manufactured, and prints a check based on the result of the check, which confirms that the glasses are made in accordance with the prescription.

Glasses touch the skin, we wear glasses in smoky rooms, we cook food in them and look into a boiling saucepan, we leave them without a case and dust settles on them. If you look closely, all these traces of wearing can be seen on the nose pads, where the lens touches the frame.

Therefore, every 2 years you need to change the frame or do ultrasonic cleaning of glasses and lenses.

Ultrasonic cleaning is a high-tech procedure that allows you to clean the lenses and frames of glasses from various types of dirt. This procedure will not do any damage to your lenses. 3.3. Правила ухода.

Reusable Microfiber Lens Care Cloths

They are composed of ultra-fine fibers. Thanks to their small diameter - only 0.06 mm - microfibers can penetrate the micropores of the lens, removing any liquid or particulate matter faster than any other tissue. Microfiber fibers literally suck dust and dirt into their matrix and hold them securely. This cloth can absorb and retain dirt and moisture up to 7 times its own weight. The main rule of effective use of these wipes is to remember to keep them clean, as a dirty wipe can scratch the lens coating.

Disposable wet wipes

They are packed in sealed bags, which is very convenient for those who save time. These wipes contain cleaner and highly evaporating alcohol, and are also composed of a thin non-woven material to prevent scratching of the lens coating.

Cleaning sprays

The cleaning ability of sprays is due to the content of surfactants - special chemical compounds. The molecules of these substances readily come into contact with pollutants, for example, fats and oils that are difficult to wash, and remove them without difficulty. After spray application, the lens should be wiped with a reusable microfiber cloth. It is important to note that lens cleaning sprays do not contain alkaline additives that can damage the surface of plastic, polycarbonate, and highly refractive organic materials. In addition to cleaning, lenses require daily care. There are a few simple rules to follow:

1. Store lenses only in their original packaging. Never put in any packaging next to other things, as this can scratch the lenses.
2. Do not put lenses down. This will keep them from being scratched.
3. Make sure that loose lenses are not exposed to direct sunlight. Heating the lens in direct sunlight can cause a fire.
4. Do not throw your lenses.
5. It is not allowed to lay lenses with their convex side on a hard surface;
6. It is not allowed to apply significant mechanical stress to the lenses;
7. It is not allowed to heat the lenses to temperatures above plus 60 ° C, including using them in a bath, sauna.

4. MAINTENANCE INFORMATION

The product is warranted for a period of 18 months from the date of purchase, it is guaranteed that the product is free from defects in workmanship and quality of workmanship.

Make sure that any deviations when working with the lenses and when using them do not occur due to non-compliance with the operating rules. In case of defects in compliance with the operating requirements in accordance with the instructions for use, the manufacturer replaces the product. For all questions, please contact: LLC "Essilor-LUYS-Optics", Legal address / actual address: 127273, Moscow, st. Otradnaya, 2B, bldg. 1, pom. No. 3.11, office. No. 51, fl. 3, INN 7743654045, OGRN 1077758508365 tel / fax +7 (495) 137-48-48, info@luis-optica.ru

5. ENVIRONMENTAL PROTECTION REQUIREMENTS.

The main type of possible hazardous impact on the environment is the pollution of atmospheric air in populated areas, soil and water as a result of:

- emergency leaks (placers) of production materials;
- unorganized waste disposal on the territory of the manufacturer or outside of it;
- their arbitrary dump in places not intended for this purpose.

The materials used in the manufacture of the product should not pose a danger to life, health of people and the environment both during operation and after the end of its life.

6. MARCHING AND PACKAGING

Products are supplied unsterile in 1 and 2 pcs. packaged.

6.1. Each package must be printed with the following information *:

- optical effect of the lens: the value of spherical (S) and cylindrical (C) refraction, diopters
- nominal size (D), mm;
- colour;
- identification of the coating;
- trade name and location of the manufacturer;
- refractive index;
- standard GOST R 53950-2010;

- number and date of the registration certificate.

For multifocal spectacle lenses, the following additional data must be specified:

- additional refraction for near, diopters;
- type and width of the segment, mm, for example, D;
- prismatic action of the lens (prism), PDPT;
- method of measuring additional refraction for near, if it is not measured from the side of the segment;
- the position of the base point for distance in aspherical multifocal lenses.

For progressive spectacle lenses, the following additional data must be specified:

- additional refraction for near, diopters;
- indication of the place of installation of the lens in the frame (right or left);
- the method of measuring the additional refraction for near, if it is not measured from the progressive side.

6.2. Requirements for marking a progressive spectacle lens:

Permanent (permanent) marking

The spectacle lens must be permanently marked as follows:

- a) two alignment points located at a distance of 34 mm from each other, equidistant from the vertical plane passing through the center of the alignment cross or the base point of the prism;
- b) the value of additional refraction for near, diopters;
- c) manufacturer's or supplier's logo or trade name or trade mark.

Additional (removable) marking

The following removable markings are recommended:

- a) designation of adjustment points;
- b) designation of the base point for distance;
- c) designation of a constructive base point for the near;
- d) designation of the installation point (installation crosshair);
- e) designation of the base point of the prism.

* The list of information on the primary packaging does not exclude the possibility of adding additional information by the manufacturer at the request of the supplier.

The product must be placed in the appropriate slot of individual packaging, and fixed in it during transportation (packaging material is cardboard of the T-23 brand, the material of the retainer (fixation during transportation) is cardboard of the T-23 brand: Each package with the product, together with operational documentation, must be packed in a box or wrapped in paper. At the customer's request, operational documentation can be provided in electronic form. The box or bag must be tied with twine or pasted over with paper tape, glue on a paper base or polyethylene tape with an adhesive layer so that the package cannot be opened without integrity violations.

For transportation, canned products in packages with operational documentation must be packed in boxes made of sheet wood materials. Board boxes should be lined with waterproof material inside.

Transport markings depending on the types of vehicles used. Marking is applied directly to containers or labels: plywood, metal, paper or fiberboard. In this case, the manipulation signs applied to the container must correspond to the values: "Fragile. Caution", "Protect from moisture".

Symbols that may appear on the packaging, including handling signs:

Symbol	Description
	Batch code
	Manufacturer

	date of manufacture
	Non-sterile
	Fragile, handle with care
	Keep away from moisture

7. TRANSPORTATION AND STORAGE

The product is transported by all types of covered vehicles in accordance with the requirements and rules for the carriage of goods in force for each type of transport.

Transportation of cargo by sea transport must be carried out in accordance with the "Safety Rules for the Carriage of General Cargo by Sea". Type of shipment - in containers and small shipment.

Transportation temperature from -50 to + 50 ° C

Relative humidity up to 100%, non-condensing.

Lenses are stored in warehouses, excluding the effects of atmospheric precipitation and aggressive environments, at a distance of at least 1 m from heating and heating devices for heated warehouses.

Storage temperature from +5 to + 40 ° C

Relative humidity up to 80%, non-condensing.

The lenses have a shelf life of 3 years.

8. MANUFACTURER'S WARRANTIES .

8.1. The manufacturer guarantees the compliance of the above lenses with the requirements of this technical document, subject to the conditions of transportation, storage and operation.

8.2. The warranty period is 18 months from the date of receipt of the lenses by the customer.

9. REPAIR INFORMATION. COMPLAINT INFORMATION

It cannot be repaired, in the case of a factory barque it must be replaced.

Complaints should be sent to:

LLC "Essilor-LUYS-Optics", Legal address / actual address: 127273, Moscow, st. Otradnaya, 2B, bldg. 1, pom. No. 3.11, office. No. 51, fl. 3, INN 7743654045, OGRN 1077758508365 tel / fax +7 (495) 137-48-48, info@luis-optica.ru

10. RECYCLING INFORMATION.

When disposing of the product, follow the requirements of SanPiN 2.1.3684-21, Class A. Also, regarding disposal, please contact the official representative of the manufacturer: OOO Essilor-LUYS-Optics, Legal address / actual address: 127273, Moscow, st. Otradnaya, 2B, bldg. 1, pom. No. 3.11, office. No. 51, fl. 3, INN 7743654045, OGRN 1077758508365 tel / fax +7 (495) 137-48-48, info@luis-optica.ru

11. PRODUCT MEETS STANDARDS:

- ISO 8980-1 "Ophthalmic optics - Uncut finished spectacle lenses - Part 1: Specifications for single-vision and multifocal lenses";
- ISO 8980-2 "Ophthalmic optics - Uncut finished spectacle lenses - Part 2: Specifications for progressive power lenses";
- ISO 14889 "Ophthalmic optics - Spectacle lenses - Fundamental requirements for uncut finished lenses"/

EXAMPLE OF A WARRANTY CARD:

Name _____

Model _____

Date of purchase _____ l.m.

Malfunction _____

Warranty center mark

model	malfunction	decision	result

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