Abbott

Abbott GmbH & Co. KG
Max-Planck-Ring 2
65205 Wiesbaden
Postfach 1303, 65011 Wiesbaden

Tel. +49 61 22 58 0
Fax +49 61 22 58 1244

I, Michael Fritz, Manager Regulatory Affairs EMEA and Pakistan, Manager International Regulatory Affairs & Division Labeling, Abbott GmbH & Co. KG, hereby certify that the attached documents that keep the number of pages indicated below are the Insert. These documents are in a full accordance with actual manuals of manufacturer.

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Date: _____

2018 - 01-25

Michael Fritz,
Manager Regulatory Manager Regulatory Affairs
EMA and Pakistan, International Regulatory
Affairs & Division Labeling

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HCV Ag
6L47

F5-Y406-2/R09
B6L4H0

Read Highlighted Changes: Revised February 2016.

Package insert instructions must be carefully followed. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this package insert.

NAME

ARCHITECT HCV Ag

INTENDED USE

The ARCHITECT HCV Ag assay is a chemiluminescent microparticle immunoassay (CMIA) for the quantitative determination of Hepatitis C virus core antigen in human serum and plasma.

SUMMARY AND EXPLANATION OF THE TEST

ARCHITECT HCV Ag is a chemiluminescent microparticle immunoassay (CMIA) using microparticles coated with monoclonal anti-HCV for the detection of HCV Ag.

HCV Ag assays are used as an aid in the diagnosis of suspected Hepatitis C viral (HCV) infection and to monitor the status of infected individuals, i.e., whether the patient's infection has resolved or the patient has become a chronic carrier of the virus. An HCV Ag assay can detect acute HCV infection in newly infected individuals who are seronegative for antibodies to HCV due to the delayed response of HCV specific antibodies. These may include patients with an elevated risk of HCV infection for example intravenous drug users or in patients with impaired immune function such as patients undergoing hemodialysis or suffering from HIV-HCV coinfections, where HCV Ag may be the only serological marker to detect HCV infection.^{1, 2} Recent studies suggest a testing algorithm using the HCV Ag test to confirm active viral replication in anti-HCV positive individuals.¹⁹⁻²¹ Using an algorithm such as those proposed in the Mederacke, et al.¹⁹, Ottiger, et al.²⁰ and Cloherty, et al.²² publications could be helpful to accelerate access to new, improved anti-viral therapies. For the diagnosis of acute or chronic hepatitis, HCV Ag reactivity should be correlated with patient history and the presence of other Hepatitis C serological markers.^{3, 4}

BIOLOGICAL PRINCIPLES OF THE PROCEDURE

The ARCHITECT HCV Ag assay is a two-step immunoassay for the quantitative determination of core antigen of Hepatitis C virus using CMIA technology with flexible assay protocols, referred to as Chemiflex.

1. Sample, Pre-Treatment Reagent 1 and Pre-Treatment Reagent 2 are combined. An aliquot of the pre-treated sample is aspirated and dispensed into a new reaction vessel. The pre-treated sample, Assay Specific Diluent and anti-HCV coated microparticles are combined. The HCV Ag present in the pre-treated sample binds to the anti-HCV coated microparticles.
2. After washing, acridinium-labeled anti-HCV conjugate is added to create a reaction mixture.
3. Following another wash cycle, Pre-Trigger and Trigger Solutions are added to the reaction mixture.
4. The resulting chemiluminescent reaction is measured as relative light units (RLUs). There is a direct relationship between the amount of HCV Ag in the sample and the RLUs detected by the ARCHITECT iSystem optics.

The concentration of Hepatitis C core antigen in the specimen is determined using a previously generated ARCHITECT HCV Ag calibration curve. If the concentration of the specimen is greater than or equal to 3.00 fmol/L, the specimen is considered reactive for HCV Ag.

For additional information on system and assay technology, refer to the ARCHITECT System Operations Manual, Section 3.

REAGENTS

Kit Contents

ARCHITECT HCV Ag 6L47

REF	6L47-29
	100
MICROPARTICLES	1 x 6.7 mL
CONJUGATE	1 x 6.1 mL
ASSAY SPECIFIC DILUENT	1 x 30.0 mL
PRE-TREATMENT REAGENT 1	1 x 14.5 mL
PRE-TREATMENT REAGENT 2	1 x 11.0 mL
SPECIMEN DILUENT	1 x 5.9 mL
MICROPARTICLES Murine anti-HCV antibody coated microparticles in 400 mM Bicine, 50 mM TRIS buffer with protein (bovine) stabilizer. Minimum concentration: 0.025% solids. Preservatives: sodium azide and antimicrobial agents.	
CONJUGATE Murine anti-HCV antibody acridinium-labeled conjugate in 80 mM BIS-TRIS with protein (bovine) stabilizer. Minimum concentration: 0.3 µg/mL. Preservatives: sodium azide and antimicrobial agents.	
ASSAY SPECIFIC DILUENT HCV Ag Assay Specific Diluent containing 1.46 N NaOH.	
PRE-TREATMENT REAGENT 1 HCV Ag Pre-Treatment Reagent 1 containing 0.83 N HCl.	
PRE-TREATMENT REAGENT 2 HCV Ag Pre-Treatment Reagent 2 containing 0.83 N HCl.	
SPECIMEN DILUENT HCV Ag Specimen Diluent containing phosphate buffer with protein (horse serum) stabilizer. Preservatives: sodium azide and antimicrobial agents.	

Other Reagents

PRE-TRIGGER SOLUTION ARCHITECT Pre-Trigger Solution containing 1.32% (w/v) hydrogen peroxide.

TRIGGER SOLUTION ARCHITECT Trigger Solution containing 0.35 N sodium hydroxide.

WASH BUFFER ARCHITECT Wash Buffer containing phosphate buffered saline solution. Preservatives: antimicrobial agents.

Warnings and Precautions

- IVD
- For In Vitro Diagnostic Use

Safety Precautions

CAUTION: This product requires the handling of human specimens. It is recommended that all human-sourced materials be considered potentially infectious and handled in accordance with the OSHA Standard on Bloodborne Pathogens. Biosafety Level 2 or other appropriate biosafety practices should be used for materials that contain or are suspected of containing infectious agents. 5-#

The following warnings and precautions apply to: **MICROPARTICLES / SPECIMEN DILUENT**

Contains sodium azide.

EUH032	Contact with acids liberates very toxic gas.
P501	Dispose of contents / container in accordance with local regulations.

The following warnings and precautions apply to: **CONJUGATE**

WARNING:	Contains sodium azide, sodium fluoride, and cetyltrimethylammonium bromide.
EUH032	Contact with acids liberates very toxic gas.
H410	Very toxic to aquatic life with long lasting effects.

Prevention

P273	Avoid release to the environment.
------	-----------------------------------

Response

P391	Collect spillage.
------	-------------------

Disposal

P501	Dispose of contents / container in accordance with local regulations.
------	---

The following warnings and precautions apply to:

ASSAY SPECIFIC DILUENT

DANGER:	Contains sodium hydroxide.
H314	Causes severe skin burns and eye damage.
H290	May be corrosive to metals.

Prevention

P234	Keep only in original container.
P260	Do not breathe mist / vapors / spray.
P264	Wash hands thoroughly after handling.
P280	Wear protective gloves / protective clothing / eye protection.

Response

P301+P330+P331	IF SWALLOWED: Rinse mouth. Do NOT induce vomiting.
P305+P351+P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P303+P361+P353	IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water / shower.
P310	Immediately call a POISON CENTER or doctor / physician.
P390	Absorb spillage to prevent material damage.

Disposal

P501	Dispose of contents / container in accordance with local regulations.
------	---

The following warnings and precautions apply to:

PRE-TREATMENT REAGENT 1 / **PRE-TREATMENT REAGENT 2**

DANGER:	Contains hydrochloric acid, dodecyltrimethylammonium bromide, and polyethylene glycol octylphenyl ether.
H314	Causes severe skin burns and eye damage.
H290	May be corrosive to metals.
H412	Harmful to aquatic life with long lasting effects.

Prevention

P234	Keep only in original container.
P260	Do not breathe mist / vapors / spray.
P264	Wash hands thoroughly after handling.
P280	Wear protective gloves / protective clothing / eye protection.
P273	Avoid release to the environment.

Response

P301+P330+P331	IF SWALLOWED: Rinse mouth. Do NOT induce vomiting.
P305+P351+P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P303+P361+P353	IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water / shower.
P310	Immediately call a POISON CENTER or doctor / physician.
P391	Collect spillage.

Disposal

P501	Dispose of contents / container in accordance with local regulations.
------	---

Safety Data Sheets are available at www.abbottdiagnostics.com or contact your local representative.

For a detailed discussion of safety precautions during system operation, refer to the ARCHITECT System Operations Manual, Section 8.

Reagent Handling

- Do not use reagent kits beyond the expiration date.
- Do not open the ASSAY SPECIFIC DILUENT plastic bag until ready to use.
- If leaking is observed with the Assay Specific Diluent bottle, the reagent kit cannot be used due to a lack of homogeneity which may impact results. The Assay Specific Diluent contains sodium hydroxide, and can cause severe skin and eye burns. Leaking bottles should be handled with appropriate safety precautions.
- Do not pool reagents within a kit or between kits.
- Before loading the reagent kit on the system for the first time, the reagents shipped on dry ice must be completely thawed and mixed thoroughly. For mixing instructions, refer to the PROCEDURE, Assay Procedure section of this package insert.

- **Septums MUST be used to prevent reagent evaporation and contamination and to ensure reagent integrity. Reliability of assay results cannot be guaranteed if septums are not used according to the instructions in this package insert.**
 - To avoid contamination, wear clean gloves when placing a septum on an uncapped reagent bottle.
 - Once a septum has been placed on an open reagent bottle, **do not invert the bottle** as this will result in reagent leakage and may compromise assay results.
 - Over time, residual liquids may dry on the septum surface. These are typically dried salts and have no effect on assay efficacy.

For a detailed discussion of handling precautions during system operation, refer to the ARCHITECT System Operations Manual, Section 7.

Reagent Storage

When stored and handled as directed, reagents are stable until the expiration date.

	Storage Temperature	Maximum Storage Time	Additional Storage Instructions
Unopened/Opened*	2-8°C	Until expiration date	Shipped on dry ice, must be stored at 2-8°C after receipt Store in upright position.
On board	System temperature	30 days	Discard after 30 days. For information on tracking onboard time, refer to the ARCHITECT System Operations Manual, Section 5.

* Reagents may be stored on or off the ARCHITECT iSystem. If reagents are removed from the system, store them at 2-8°C (with septums and replacement caps) in an upright position. For reagents stored off the system, it is recommended that they be stored in their original trays and boxes to ensure they remain upright. **If the microparticle bottle does not remain upright (with a septum installed) while in refrigerated storage off the system, the reagent kit must be discarded.** For information on unloading reagents, refer to the ARCHITECT System Operations Manual, Section 5.

Indications of Reagent Deterioration

When a control value is out of the specified range, it may indicate deterioration of the reagents or errors in technique. Associated test results are invalid, and samples must be retested. Assay recalibration may be necessary. For troubleshooting information, refer to the ARCHITECT System Operations Manual, Section 10.

INSTRUMENT PROCEDURE

The ARCHITECT HCV Ag assay file (801_004.iae or higher) must be installed on the ARCHITECT iSystem from an ARCHITECT iSystem Assay CD-ROM prior to performing the assay.

For detailed information on assay file installation and viewing and editing assay parameters, refer to the ARCHITECT System Operations Manual, Section 2.

If the ARCHITECT Anti-HCV assay and the ARCHITECT HCV Ag assay are run on the same ARCHITECT iSystem, the Anti-HCV assay file (161_020.iae or higher) must be installed from ARCHITECT iSystem e-Assay CD-ROM found on www.abbottdiagnostics.com, prior to the installation of the HCV Ag assay file.

For information on printing assay parameters, refer to the ARCHITECT System Operations Manual, Section 5.

For a detailed description of system procedures, refer to the ARCHITECT System Operations Manual.

Alternate Result Units

Edit assay parameter "Result concentration units" to select an alternate unit.

(Concentration in fmol/L) x (0.02) = (Concentration in pg/mL)

Default result unit	Conversion factor	Alternate result unit
fmol/L	0.02	pg/mL

SPECIMEN COLLECTION AND PREPARATION FOR ANALYSIS

Specimen Types

Verified specimen types to be used with this assay:

Specimen Types	Collection Tubes
Human serum	Serum Serum separator tubes
Human plasma	Sodium EDTA Potassium EDTA Lithium Heparin Sodium Heparin Sodium Citrate CPD

- Other specimen collection tube types have not been tested with this assay.
- Liquid anticoagulants may have a dilution effect resulting in lower concentrations for individual patient specimens.
- The instrument does not provide the capability to verify specimen type. It is the responsibility of the operator to verify that the correct specimen types are used in the assay.

Specimen Conditions

- Do not use specimens with the following conditions:
 - heat-inactivated
 - pooled
 - grossly hemolyzed
 - obvious microbial contamination
- Performance has not been established for the use of cadaveric specimens or the use of body fluids other than human serum and plasma.
- For optimal results, serum and plasma specimens should be free of fibrin, red blood cells, or other particulate matter. Such specimens may give inconsistent results and must be transferred to a centrifuge tube and centrifuged at 3000 x g for 10 minutes. Serum specimens from patients receiving anticoagulant or thrombolytic therapy may contain fibrin due to incomplete clot formation.
- When using plasma primary tubes, samples must be priority loaded immediately after centrifugation to minimize disturbance of the red blood cells.
- To prevent cross contamination, use of disposable pipettes or pipette tips is recommended.

Preparation for Analysis

- Follow the tube manufacturer's processing instructions for collection tubes. Gravity separation is not sufficient for specimen preparation.
- Mix thawed specimens thoroughly by low speed vortexing or by inverting 10 times. Visually inspect the specimens. If layering or stratification is observed, continue mixing until specimens are visibly homogeneous.
- To ensure consistency in results, specimens must be transferred to a centrifuge tube and **centrifuged at 3000 x g for 10 minutes** before testing if
 - they contain fibrin, red blood cells, or other particulate matter,
 - they require repeat testing, or
 - they were frozen and thawed.

- Transfer clarified specimen to a sample cup or secondary tube for testing. For centrifuged specimens with a lipid layer, transfer only the clarified specimen and not the lipemic material.
- Inspect all specimens for bubbles. Remove bubbles with an applicator stick before analysis. Use a new applicator stick for each specimen to prevent cross contamination.

Specimen Storage

Specimen Type	Storage Temperature	Maximum Storage Time
Serum or plasma	2-8°C	≤ 5 days
	- 20°C or colder	--

Specimen may be stored on or off the clot, red blood cells, or separator gel for up to 5 days refrigerated at 2-8°C.

If testing will be delayed more than 5 days, remove serum or plasma from the clot, red blood cells, or separator gel and store frozen at -20°C or colder.

Avoid more than two freeze/thaw cycles.

Specimen Shipping

- Package and label specimens in compliance with applicable state, federal, and international regulations covering the transport of clinical specimens and infectious substances.
- Do not exceed the storage limitations listed above.

PROCEDURE

Materials Provided

6L47 ARCHITECT HCV Ag Reagent Kit

Materials Required but not Provided

- ARCHITECT HCV Ag Assay file obtained from the ARCHITECT iSystem e-Assay CD-ROM found on www.abbottiagnostics.com.
- 6L47-02 ARCHITECT HCV Ag Calibrators
- 6L47-11 ARCHITECT HCV Ag Controls
- ARCHITECT Pre-Trigger Solution
- ARCHITECT Trigger Solution
- ARCHITECT Wash Buffer
- ARCHITECT Reaction Vessels
- ARCHITECT Sample Cups
- ARCHITECT Septum
- ARCHITECT Replacement Caps
- Pipettes or pipette tips (optional) to deliver the volumes specified on the patient or control order screen.

For information on materials required for maintenance procedures, refer to the ARCHITECT System Operations Manual, Section 9.

Assay Procedure

- Before loading the reagent kit on the system for the first time, **the reagents shipped on dry ice must be completely thawed and mixed thoroughly**. After the first time the reagents have been loaded, no further mixing is required.
 - Tear open the ASSAY SPECIFIC DILUENT plastic bag.
 - Invert **all the reagent bottles** 30 times.
 - Visually inspect the bottle to ensure microparticles are resuspended. If microparticles are still adhered to the bottle, continue to invert the bottle until the microparticles have been completely resuspended.
 - **If the microparticles do not resuspend, DO NOT USE. Contact your local Abbott representative.**
 - Once the microparticles have been resuspended, place a septum on the bottle. For instructions on placing septums on bottles, refer to the **Reagent Handling** section of this package insert.
- Load the reagent kit on the ARCHITECT iSystem.
 - Verify that all necessary reagents are present.
 - Ensure that septums are present on all reagent bottles.

- Order calibration, if necessary.
 - For information on ordering calibrations, refer to the ARCHITECT System Operations Manual, Section 6.
- Order tests.
 - For information on ordering patient specimens and controls and for general operating procedures, refer to the ARCHITECT System Operations Manual, Section 5.
- Minimum sample cup volume is calculated by the system and printed on the Orderlist report. To minimize the effects of evaporation, verify adequate sample cup volume is present prior to running the test.

Maximum number of replicates sampled from the same sample cup: 10

 - Priority:
 - Sample volume for first test: 158 µL
 - Sample volume for each additional test from same sample cup: 108 µL
 - ≤ 3 hours on board:
 - Sample volume for first test: 158 µL
 - Sample volume for each additional test from same sample cup: 108 µL
 - > 3 hours on board: replace with a fresh sample (patient specimens, controls, and calibrators).
 - If using primary or aliquot tubes, use the sample gauge to ensure sufficient patient specimen is present.
- Prepare ARCHITECT HCV Ag Calibrators and Controls.
 - Mix calibrator(s) and controls by gentle inversion before use.
 - Hold bottles **vertically** and dispense recommended volumes into each respective sample cup.
 - Recommended volumes:
 - for each calibrator: 12 drops
 - for each control: 7 drops
- Load samples.
 - For information on loading samples, refer to the ARCHITECT System Operations Manual, Section 5.
- Press RUN.
- For additional information on principles of operation, refer to the ARCHITECT System Operations Manual, Section 3.
- For optimal performance, it is important to perform routine maintenance as described in the ARCHITECT System Operations Manual, Section 9. Perform maintenance more frequently when required by laboratory procedures.

Specimen Dilution Procedures

Specimens with an HCV Ag concentration of > 20,000 fmol/L will be flagged as ">20,000" and may be diluted with either the Automated Dilution Protocol or the Manual Dilution Procedure.

Automated Dilution Protocol

The system performs a 1:9 dilution of the specimen and automatically calculates the concentration of the specimen before dilution and reports the result.

Manual Dilution Procedure

Suggested dilution: 1:20

1. Add 20 µL of the patient specimen to 380 µL of ARCHITECT HCV Ag Negative Control.
2. The operator must enter the dilution factor in the Patient or Control order screen. The system will use this dilution factor to automatically calculate the concentration of the sample before dilution and report the result. The result (before the dilution factor is applied) should be greater than 3.00 fmol/L.

For detailed information on ordering dilutions, refer to the ARCHITECT System Operations Manual, Section 5.

Calibration

- Test Calibrators A-F in duplicate and a single sample of each control level as below,
 1. Order the calibration. **DO NOT START THE RUN.**
 2. Order 1 rep of CONTROL+2, 1 rep of CONTROL+1 and 1 rep of CONTROL- in this order.
 3. Then, start the run.

The calibrators should be priority loaded. Ensure that assay control values are within the concentration ranges specified in the control package insert.
- Calibration Range: 0 – 20,000 fmol/L.
- Once an ARCHITECT HCV Ag calibration is accepted and stored, all subsequent samples may be tested without further calibration unless:
 - A reagent kit with a new lot number is used or
 - Controls are out of range.
- For detailed information on how to perform an assay calibration, refer to the ARCHITECT System Operations Manual, Section 6.

Quality Control Procedures

- The recommended control requirement for the ARCHITECT HCV Ag assay is that a single sample of each control level be tested once every 24 hours each day of use. If the quality control procedures in your laboratory require more frequent use of controls to verify test results, follow your laboratory-specific procedures.
- The ARCHITECT HCV Ag Control values must be within the acceptable ranges specified in the control package insert. If a control is out of its specified range, the associated test results are invalid and samples must be retested. Recalibration may be indicated.

Verification of Assay Claims

For protocols to verify package insert claims, refer to the ARCHITECT System Operations Manual, Appendix B.
The ARCHITECT HCV Ag assay belongs to method group 5.

RESULTS

The ARCHITECT HCV Ag assay utilizes a 4 Parameter Logistic Curve fit data reduction method (4PLC, Y-weighted) to generate a calibration curve.

Flags

Some results may contain information in the Flags field. For a description of the flags that may appear in this field, refer to the ARCHITECT System Operations Manual, Section 5.

Interpretation of Results

Concentration Value	Interpretation
< 3.00 fmol/L	Nonreactive for HCV Ag
≥ 3.00 fmol/L	Reactive for HCV Ag

Specimens with concentration values ≥ 3.00 fmol/L to < 10.00 fmol/L should be retested in duplicate.

Duplicate Retest Results	
Instrument Interpretation	Specimen Classification
Both retest values nonreactive	Specimen considered nonreactive for HCV Ag.
One or both of the duplicates is (are) ≥ 3.00 fmol/L	Specimen considered repeatedly reactive for HCV Ag, and the initial value is used as the final reported value.

For details on configuring the ARCHITECT iSystem to use grayzone interpretations, refer to the ARCHITECT System Operations Manual, Section 2.

LIMITATIONS OF THE PROCEDURE

- For diagnostic purposes, results should be used in conjunction with other data; e.g., symptoms, results of other tests, clinical impressions, etc.
- Specimens from patients who have received preparations of mouse monoclonal antibodies for diagnosis or therapy may contain human anti-mouse antibodies (HAMA). Such specimens may show either falsely elevated or depressed values when tested with assay kits such as ARCHITECT HCV Ag that employ mouse monoclonal antibodies. Additional information may be required for diagnosis.^{9, 10} The ARCHITECT HCV Ag reagents contain a component that reduces the effect of HAMA reactive specimens. Additional clinical or diagnostic information may be required to determine patient status.
- If the ARCHITECT HCV Ag results are inconsistent with clinical evidence, additional testing is suggested to confirm the result.
- Samples containing particulate matter or red blood cells must be centrifuged prior to running the assay. Insufficiently separated plasma specimens from clots or red blood cells must not be used.
- Specimens from patients with high levels of IgM, e.g., specimens from patients with multiple myeloma, may generate “3350 Unable to process test, aspiration error for (Sample Pipettor) at (RV 24)”.
- **The ARCHITECT Anti-HCV assay (LN 6C37) may cause falsely elevated HCV Ag results in the ARCHITECT HCV Ag assay when both assays are run on the same module. As a preventive measure to protect the integrity of test results, it is important that one of the following actions be taken:**
 - a. Multiple I2000/I2000SR module customers should run the ARCHITECT HCV Ag and ARCHITECT Anti-HCV assays on separate modules or,
 - b. If you are unable to run these assays on separate modules:
 - i. Segregate the ARCHITECT HCV Ag samples and run them immediately following the ARCHITECT 6041 Daily Maintenance procedure.
 - ii. After running the ARCHITECT HCV Ag samples, it is recommended to immediately perform the ARCHITECT 6041 Daily Maintenance procedure again, as an extra precaution.
 - iii. It is recommended that duplicate retest be performed on samples between the range of 3 - 15 fmol/L instead of the range of 3 - 10 fmol/L in the RESULTS section.
- Due to the naturally occurring diversity of hepatitis C viruses, including amino acid variations within the core gene, rare specimens from HCV infected patients may produce unexpectedly low or negative results in the ARCHITECT HCV Ag assay.¹²⁻¹⁴ Two studies noted reduced assay sensitivity among a subpopulation of individuals infected with HCV genotype 3a.^{19, 14} Additional studies have reported varying detection rates for genotype 3 samples on the ARCHITECT HCV Ag test when compared to HCV RNA testing.¹³⁻¹⁸
 - In patients that are anti-HCV positive but HCV Ag negative, it is advisable to perform additional clinical or diagnostic testing (e.g. molecular testing) to determine patient status.

SPECIFIC PERFORMANCE CHARACTERISTICS

Precision

The ARCHITECT HCV Ag assay precision is < 10% total CV. A study was performed as described in the National Committee for Clinical Laboratory Standards (NCCLS) Protocol EP5-A2.¹¹ Five samples consisting of two buffer protein based HCV Ag positive controls and three serum based panels were assayed in replicates of two at two separate times per day for twenty days (n=80 for each sample), using three lots of reagents. Data from this study are summarized in the following table.*

Table I: ARCHITECT HCV Ag Precision

Sample	Instrument	Reagent		Mean Value	Within Run		Total	
		Lot	n		SD	%CV	SD	%CV
Level 1	1	1	80	52.41	4.82	9.2	4.96	9.5
		2	80	50.08	3.76	7.5	4.07	8.1
		3	80	52.92	3.80	7.2	4.10	7.8
	2	1	80	53.57	4.40	8.2	4.73	8.8
		2	80	52.31	3.81	7.3	3.81	7.3
		3	80	51.98	3.21	6.2	3.28	6.3
Level 2	1	1	80	303.18	20.60	6.8	23.01	7.6
		2	80	288.91	17.27	6.0	18.91	6.5
		3	80	303.95	22.83	7.5	24.03	7.9
	2	1	80	320.62	20.53	6.4	20.53	6.4
		2	80	311.64	14.79	4.7	15.64	5.0
		3	80	306.22	11.98	3.9	13.32	4.3
Panel 1	1	1	80	67.54	4.74	7.0	4.90	7.3
		2	80	65.09	4.35	6.7	5.13	7.9
		3	80	67.14	3.54	5.3	4.64	6.9
	2	1	80	65.57	4.61	7.0	5.00	7.6
		2	80	63.70	4.12	6.5	4.47	7.0
		3	80	63.09	4.64	7.4	4.64	7.4
Panel 2	1	1	80	327.17	22.40	6.8	24.81	7.6
		2	80	310.52	16.11	5.2	26.28	8.5
		3	80	324.36	15.22	4.7	17.29	5.3
	2	1	80	337.09	21.37	6.3	22.06	6.5
		2	80	322.28	19.84	6.2	21.83	6.8
		3	80	310.09	14.46	4.7	17.65	5.7
Panel 3	1	1	80	7629.26	406.99	5.3	549.08	7.2
		2	80	7765.96	545.13	7.0	599.72	7.7
		3	80	7360.23	318.91	4.3	363.50	4.9
	2	1	80	8231.22	475.04	5.8	508.93	6.2
		2	80	7553.65	361.03	4.8	420.65	5.6
		3	80	7516.93	254.33	3.4	359.53	4.8

* Representative data; results in individual laboratories may vary from these data.

Specificity

The ARCHITECT HCV Ag assay demonstrated a specificity of $\geq 99.5\%$ in a study where specimens from a blood donor population, hospitalized patients and specimens containing potentially interfering substances were tested. This study includes the specimens from individuals with medical conditions unrelated to HCV infection (Table II). A total of 5027 serum and plasma specimens from blood donors were evaluated. The initial and repeat reactive rates were 0.24% (12/5027) and 0.02% (1/5027), respectively. Four of 250 specimens obtained from hospital patients were repeatedly reactive and confirmed positive for HCV infection. In 126 specimens from individuals with medical conditions unrelated to HCV infection and specimens containing potentially interfering substances, five specimens were repeatedly reactive and confirmed positive for HCV infection.

Table II: Reactivity of the ARCHITECT HCV Ag Assay in Specimens from Blood Donors, Hospital Patients, Individuals with Medical Conditions Unrelated to HCV Infection, and in Specimens Containing Potentially Interfering Substances

Category	Number Tested	IR (% of Total)	RR (% of Total)	Number of Positive by Supplemental Testing ^a (% of RR)
Blood Donors				
Serum	2256	7 (0.31)	1 (0.04)	0 (0.00)
EDTA Plasma	411	0 (0.00)	0 (0.00)	0 (0.00)
Na Citrate Plasma	1180	2 (0.17)	0 (0.00)	0 (0.00)
Heparin Plasma	1180	3 (0.25)	0 (0.00)	0 (0.00)
Total Donors	5027	12 (0.24)	1 (0.02)	0 (0.00)
Hospital Patients	250	4 (1.60)	4 (1.60)	4 (100.0)
Individuals with Medical Conditions Unrelated to HCV Infection and Specimens Containing Potentially Interfering Substances ^b	126	5 (3.97)	5 (3.97)	5 (100.0)

IR = Initially Reactive; RR = Repeatedly Reactive

^a A positive result was defined as reactive to Abbott AxSYM Anti-HCV (1D30), Roche Amplicor HCV, and/or Abbott RealTime HCV (4J86).

^b Category includes the following: anti-CMV positive (5), anti-EBV positive (5), anti-HAV positive (5), HBsAg positive (5), anti-HIV-1 positive (9), anti-HIV-2 positive (1), HIV Ag positive (5), anti-HTLV-1 positive (5), anti-HBc positive (5), Syphilis positive (5), rheumatoid factor (5), alcoholic liver disease (5), human anti-mouse antibody positive (5), *E. coli* (5), anti-HSV-1 positive (5), multiparous female (5), hemodialysis patients (5), anti-nuclear antibody positive (5), elevated IgG (5), elevated IgM (5), multiple myeloma (3), pregnant female (5), influenza vaccine recipients (5), Toxoplasma positive (3), West Nile virus positive (5), multiple transfusion recipients (5).

Sensitivity

The ARCHITECT HCV Ag assay has a sensitivity of ≤ 3.00 fmol/L.

A total of 452 serum and plasma specimens known to be positive for HCV RNA including genotypes 1a, 1b, 2a, 2b, 3a, 3k, 4a, 5a, 6a, and 6i, were tested. Of the 452 specimens, 97.8% (442/452) were reactive.

Seroconversion

Sensitivity of the ARCHITECT HCV Ag assay was evaluated utilizing 10 commercially available panels of sequential specimens from patients who seroconverted for the detection of anti-HCV antibodies. In each panel, a positive result in the ARCHITECT HCV Ag assay was obtained prior to detection of anti-HCV antibody, resulting in an average reduction between the times of infection and detection of 35.8 days.

Interference

The ARCHITECT HCV Ag assay is designed to have a mean interference of $\leq 10\%$ difference in concentration for patient samples with triglycerides (3000 mg/dL), bilirubin (20 mg/dL), hemoglobin (500 mg/dL), and protein (9.2 g/dL).

In a representative study, the interference from hemoglobin, bilirubin, triglycerides, and protein was evaluated in the ARCHITECT HCV Ag assay. The following interferences were obtained:

- Hemoglobin < 10% at 500 mg/dL
- Bilirubin < 10% at 20 mg/dL
- Triglycerides < 10% at 3000 mg/dL
- Protein < 10% at 9.2 g/dL

BIBLIOGRAPHY

1. Miedouge M, Saune K, Kamar N, Rieu M, Rostaing L, Izopet J. Analytical evaluation of HCV core antigen and interest for HCV screening in haemodialysis patients. *J Clin Virol* 2010;48:18–21.
2. Thomson EC, Nastouli E, Main J, Karayiannis P, Eliahoo J, Muir D et al. Delayed anti-HCV antibody response in HIV-positive men acutely infected with HCV. *AIDS* 2009;23:89–93.
3. Aoyagi K, Ohue C, Iida K, et al. Development of a simple and highly sensitive enzyme immunoassay for Hepatitis C virus core antigen. *J Clin Microbiol.* 1999;37:1802-1808.
4. Fabrizi F, Lunghi G, Aucella F, et al. Novel assay using total hepatitis C virus (HCV) core antigen quantification for diagnosis of HCV infection in dialysis patients. *J Clin Microbiol.* 2005;43:414–420.
5. US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
6. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
7. World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
8. Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.
9. Primus FJ, Kelley EA, Hansen HJ, et al. "Sandwich"-type immunoassay of carcinoembryonic antigen in patients receiving murine monoclonal antibodies for diagnosis and therapy. *Clin Chem* 1988;34(2):261-264.
10. Schroff RW, Foon KA, Beatty SM, et al. Human anti-murine immunoglobulin responses in patients receiving monoclonal antibody therapy. *Cancer Res* 1985;45(2):879-885.
11. National Committee for Clinical Laboratory Standards (NCCLS). *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline—Second Edition*. NCCLS Document EP5-A2. Wayne, PA: NCCLS; 2004.
12. Murayama et al., Japanese Reference Panel of Blood Specimens for Evaluation of Hepatitis C Virus RNA and Core Antigen Quantitative Assays. *J. Clin. Microbiol.* 50:1943. 2012.
13. Ottiger, C.; Gygli, N.; Huber, A.R., Detection limit of architect hepatitis C core antigen assay in correlation with HCV RNA, and renewed confirmation algorithm for reactive anti-HCV samples. *Journal of Clinical Virology* 58:535-540. 2013.
14. Garbuglia et al., HCV core antigen and HCV-RNA in HIV/HCV co-infected patients with different HCV genotypes. *BMC Infectious Diseases* 2014;14:222.
15. Morota et al., A new sensitive and automated chemiluminescent microparticle immunoassay for quantitative determination of hepatitis C virus core antigen. *J Virol Methods* 157:8. 2009.
16. Heidrich et al., Hepatitis C virus core antigen testing in liver and kidney transplant recipients. *J Viral Hepat* 21:769. 2014.
17. Mederacke et al., Performance and clinical utility of a novel fully automated quantitative HCV-core antigen assay. *Journal of Clinical Virology* 46:210-215. 2009.
18. Nguyen et al., Hepatitis C Virus Core Mutations Associated with False-Negative Serological Results for Genotype 3a Core Antigen. *Journal of clinical microbiology* 53 (8) p.2697-2700.
19. Mederacke I, Wedermeyer H, Ciesek S, et al. Performance and Clinical Utility of a Novel Fully Automated HCV-Core Antigen Assay. *Journal of Clinical Virology* 2009;46:210-215.
20. Ottiger, C.; Gygli, N.; Huber, A.R., Detection limit of architect hepatitis C core antigen assay in correlation with HCV RNA, and renewed confirmation algorithm for reactive anti-HCV samples. *Journal of Clinical Virology* 58:535-540. 2013.
21. Cloherty G, Cheng K, Chevaliez S et al., Clinical utility of the Abbott ARCHITECT Hepatitis C core antigen testing in the staging and monitoring of treatment naïve non-cirrhotic patients receiving an all-oral, interferon-free regimen. *Journal of Hepatology* Volume 62, Supplement 2, Page S617.
22. Cloherty G.; Talal A.; Collier K. et al., Role of serologic and molecular diagnostic assays in identification and management of hepatitis C virus infection. *J Clin Microbiol* 54:265–273. 2016.

Key to Symbols

	Consult instructions for use
	Manufacturer
	Sufficient for
	Temperature limitation
	Use by/Expiration date
ASSAY SPECIFIC DILUENT	Assay Specific Diluent
CAL A	Calibrator (A,B,C,D,E or F)
CONJUGATE	Conjugate
CONTAINS: AZIDE	Contains Sodium Azide. Contact with acids liberates very toxic gas.
CONTROL -	Negative Control
CONTROL + 1	Positive Control 1
CONTROL + 2	Positive Control 2
CONTROL NO.	Control Number
ECO HAZARD	Ecological hazard
IVD	In Vitro Diagnostic Medical Device
LOT	Lot Number
MICROPARTICLES	Microparticles
PRE-TREATMENT REAGENT 1	Pre-Treatment Reagent 1
PRE-TREATMENT REAGENT 2	Pre-Treatment Reagent 2
PRE-TRIGGER SOLUTION	Pre-Trigger Solution
PRODUCED FOR ABBOTT BY	Produced for Abbott by
PRODUCT OF JAPAN	Product of Japan
REACTION VESSELS	Reaction Vessels
REAGENT LOT	Reagent Lot
REF	List Number
REPLACEMENT CAPS	Replacement Caps
SAMPLE CUPS	Sample Cups
SEPTUM	Septum
SN	Serial number
SPECIMEN DILUENT	Specimen Diluent
TRIGGER SOLUTION	Trigger Solution
WASH BUFFER	Wash Buffer

ARCHITECT, Chemiflex and AxSYM are trademarks of Abbott Laboratories in various jurisdictions. All other trademarks are property of their respective owners.

Anti-HCV monoclonal antibodies are prepared under US license by Chiron Corporation under a shared manufacturing agreement. The ARCHITECT HCV Ag assay is manufactured under contract agreement from Ortho Diagnostic Systems and Chiron Corporation.

Abbott GmbH & Co. KG
Max-Planck-Ring 2
65205 Wiesbaden
Germany
+49-6122-580

CE
0088

PRODUCED FOR ABBOTT BY

DENKA SEIKEN CO., LTD. Tokyo, Japan

**Customer Service: Contact your local representative
or find country-specific contact information on
www.abbottdiagnostics.com**

Revised February 2016.

©2009, 2016 Abbott Laboratories

ADDENDUM

to the Instruction for Use

ARCHITECT HCV Ag

Reagent for Russia

SHIPPING

Recommended shipping conditions for device correspond to storage conditions provided in the Instructions for Use (2-8°C). Device is stable until the expiration date indicated on the package. Device is no longer suitable for use if shipping temperature limits have been violated (GOST R 51088-2013). Ship in an enclosed carrier vehicle of any kind in accordance with shipping rules, applicable to the respective carrier.

WARRANTY

Manufacturer guarantees that device remains stable until the expiration date when stored and handled as directed. Manufacturer set the warranty shelf life of device which corresponds to period between the date of manufacture and expiration date that are indicated on the label and/or in accompanying documentation. Device must be discarded if stored with deviations in manufacturer's requirements (GOST R 51088-2013).

MAINTENANCE

Not required.

REPAIRS

Not required.

POTENTIAL USER

Professional user

INFORMATION ABOUT PRESENCE OF DRUGS IN COMPOSITION OF MEDICAL DEVICE

Does not contain any drugs

STERILIZATION

Medical device is not required to be sterile.

INFORMATION ABOUT POSSIBILITY OF MULTIPLE APPLICATION

Medical device cannot be re-applied

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1. Miedouge M, Saune K, Kamar N, Rieu M, Rostaing L, Izopet J. Analytical evaluation of HCV core antigen and interest for HCV screening in haemodialysis patients. *J Clin Virol* 2010;48:18–21.
2. Thomson EC, Nastouli E, Main J, Karayiannis P, Eliahoo J, Muir D et al. Delayed anti-HCV antibody response in HIV-positive men acutely infected with HCV. *AIDS* 2009;23:89–93.
3. Aoyagi K, Ohue C, Iida K, et al. Development of a simple and highly sensitive enzyme immunoassay for Hepatitis C virus core antigen. *J Clin Microbiol.* 1999;37:1802–1808.
4. Fabrizi F, Lunghi G, Aucella F, et al. Novel assay using total hepatitis C virus (HCV) core antigen quantification for diagnosis of HCV infection in dialysis patients. *J Clin Microbiol.* 2005;43:414–420.
5. US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
6. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
7. World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
8. Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.
9. Primus FJ, Kelley EA, Hansen HJ, et al. "Sandwich"-type immunoassay of carcinoembryonic antigen in patients receiving murine monoclonal antibodies for diagnosis and therapy. *Clin Chem* 1988;34(2):261–264.
10. Schreff RW, Foon KA, Beatty SM, et al. Human anti-murine immunoglobulin responses in patients receiving monoclonal antibody therapy. *Cancer Res* 1985;45(2):879–885.
11. National Committee for Clinical Laboratory Standards (NCCLS). *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline—Second Edition*. NCCLS Document EP5-A2. Wayne, PA: NCCLS; 2004.
12. Murayama et al., Japanese Reference Panel of Blood Specimens for Evaluation of Hepatitis C Virus RNA and Core Antigen Quantitative Assays. *J. Clin. Microbiol.* 50:1943. 2012.
13. Ottiger, C.; Gygli, N.; Huber, A.R., Detection limit of architect hepatitis C core antigen assay in correlation with HCV RNA, and renewed confirmation algorithm for reactive anti-HCV samples. *Journal of Clinical Virology* 58:535–540. 2013.
14. Garbuglia et al., HCV core antigen and HCV-RNA in HIV/HCV co-infected patients with different HCV genotypes. *BMC Infectious Diseases* 2014;14:222.
15. Morota et al., A new sensitive and automated chemiluminescent microparticle immunoassay for quantitative determination of hepatitis C virus core antigen. *J Virol Methods* 157:8. 2009.
16. Heidrich et al., Hepatitis C virus core antigen testing in liver and kidney transplant recipients. *J Viral Hepat* 21:769. 2014.
17. Mederacke et al., Performance and clinical utility of a novel fully automated quantitative HCV-core antigen assay. *Journal of Clinical Virology* 46:210–215. 2009.
18. Nguyen et al., Hepatitis C Virus Core Mutations Associated with False-Negative Serological Results for Genotype 3a Core Antigen. *Journal of clinical microbiology* 53 (8) p.2697–2700.
19. Mederacke I, Wedermeyer H, Ciesek S, et al. Performance and Clinical Utility of a Novel Fully Automated HCV-Core Antigen Assay. *Journal of Clinical Virology* 2009;46:210–215.
20. Ottiger, C.; Gygli, N.; Huber, A.R., Detection limit of architect hepatitis C core antigen assay in correlation with HCV RNA, and renewed confirmation algorithm for reactive anti-HCV samples. *Journal of Clinical Virology* 58:535–540. 2013.
21. Cloherty G, Cheng K, Chevaliez S et al., Clinical utility of the Abbott ARCHITECT Hepatitis C core antigen testing in the staging and monitoring of treatment naïve non-cirrhotic patients receiving an all-oral, interferon-free regimen. *Journal of Hepatology* Volume 62, Supplement 2, Page S617.
22. Cloherty G.; Talal A.; Collier K. et al., Role of serologic and molecular diagnostic assays in identification and management of hepatitis C virus infection. *J Clin Microbiol* 54:265–273. 2016.

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INTENDED USE

The ARCHITECT HCV Ag Calibrators are for the calibration of the ARCHITECT iSystem when used for the quantitative determination of core antigen to Hepatitis C virus in human serum and plasma. Refer to the ARCHITECT HCV Ag reagent package insert and the ARCHITECT System Operations Manual for additional information.

CONTENTS

6 Bottles (4 mL each) of ARCHITECT HCV Ag Calibrators. Calibrator A contains citrate buffer. Calibrators B through F contain recombinant core antigen to Hepatitis C virus prepared in acetate buffer.

The calibrators yield the following concentrations:

Calibrator	Concentration	
	(fmol/L)	(pg/mL)
CAL A	0.00	0.00
CAL B	10.00	0.20
CAL C	80.00	1.60
CAL D	600.00	12.00
CAL E	4,500.00	90.00
CAL F	20,000.00	400.00

STANDARDIZATION

The ARCHITECT HCV Ag Calibrators are referenced to the manufacturer's Internal Reference standard.

PRECAUTIONS

- IVD
- For In Vitro Diagnostic Use

The following warnings and precautions apply to: CAL B - CAL F

DANGER:	Contains myristyltrimethylammonium bromide, dodecyltrimethylammonium bromide, polyethylene glycol octylphenyl ether, and cetyltrimethylammonium bromide.
H318	Causes serious eye damage.
H315	Causes skin irritation.
H410	Very toxic to aquatic life with long lasting effects.
Prevention	
P264	Wash hands thoroughly after handling.
P280	Wear protective gloves / protective clothing / eye protection.
P273	Avoid release to the environment.

Response

P305+P351+P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P310	Immediately call a POISON CENTER or doctor / physician.
P302+P352	IF ON SKIN: Wash with plenty of water.
P332+P313	If skin irritation occurs: Get medical advice / attention.
P362+P364	Take off contaminated clothing and wash it before reuse.
P391	Collect spillage.
Disposal	
P501	Dispose of contents / container in accordance with local regulations.

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STORAGE

- Calibrators are stable until the expiration date when stored and handled as directed.
- Do not use past expiration date.

PREPARATION FOR USE

Calibrators may be used immediately after removal from 2-8°C storage.

Prior to use, mix by gentle inversion (5-10 times).

After each use, tightly close the caps and return the calibrators to 2-8°C storage.

Key to Symbols

	Consult instructions for use
	Manufacturer
	Temperature limitation
	Use by/Expiration date
CAL A	Calibrator (A,B,C,D,E or F)
IVD	In Vitro Diagnostic Medical Device
LOT	Lot Number
PRODUCED FOR ABBOTT BY	Produced for Abbott by
PRODUCT OF JAPAN	Product of Japan
REF	List Number

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INTENDED USE

The ARCHITECT HCV Ag Controls are for the verification of the accuracy and precision of the ARCHITECT iSystem when used for the quantitative determination of core antigen to Hepatitis C virus in human serum and plasma. Refer to the ARCHITECT HCV Ag reagent package insert and the ARCHITECT System Operations Manual for additional information.

CONTENTS

3 Bottles (8 mL each) of ARCHITECT HCV Ag Controls. The Negative Control (1 Bottle) contains recalcified human plasma. Preservatives: sodium azide, ProClin 300 and ProClin 950. The Positive Controls 1 and 2 (2 Bottles) contain recombinant core antigen to Hepatitis C virus prepared in acetate buffer.

The controls are at the following concentrations:

Control	Color	Concentration		Range	
		(fmol/L)	(pg/mL)	(fmol/L)	(pg/mL)
CONTROL -	Natural	0.00	0.00	0.00 - 2.99	0.00 - 0.05
CONTROL +1	Blue*	50.00	1.00	35.00 - 65.00	0.70 - 1.30
CONTROL +2	Red**	300.00	6.00	210.00 - 390.00	4.20 - 7.80

* Dye: Acid Blue No. 9

** Dye: Red D&C No. 33

STANDARDIZATION

The ARCHITECT HCV Ag Controls are referenced to the manufacturer's Internal Reference standard.

PRECAUTIONS

- IVD
- For In Vitro Diagnostic Use

-  **CAUTION:** This product contains human-sourced and/or potentially infectious components. Refer to the CONTENTS section of this package insert. No known test method can offer complete assurance that products derived from human sources or inactivated microorganisms will not transmit infection. Therefore, all human-sourced materials should be considered potentially infectious. It is recommended that these reagents and human specimens be handled in accordance with the OSHA Standard on Bloodborne Pathogens. Biosafety Level 2 or other appropriate biosafety practices should be used for materials that contain or are suspected of containing infectious agents.¹⁻⁴
- The human plasma used in the Negative Control is nonreactive for HBsAg, HIV-1 Ag or HIV-1 RNA, anti-HCV, and anti-HIV-1/HIV-2.

The following warnings and precautions apply to: CONTROL -	
	
WARNING:	Contains methylisothiazolones and sodium azide.
H317	May cause an allergic skin reaction.
EUH032	Contact with acids liberates very toxic gas.

Prevention	
P261	Avoid breathing mist / vapors / spray.
P272	Contaminated work clothing should not be allowed out of the workplace.
P280	Wear protective gloves / protective clothing / eye protection.
Response	
P302+P352	IF ON SKIN: Wash with plenty of water.
P333+P313	If skin irritation or rash occurs: Get medical advice / attention.
P362+P364	Take off contaminated clothing and wash it before reuse.
Disposal	
P501	Dispose of contents / container in accordance with local regulations.

The following warnings and precautions apply to: CONTROL +1 / CONTROL +2	
	
DANGER:	Contains myristyltrimethylammonium bromide, dodecyltrimethylammonium bromide, polyethylene glycol octylphenyl ether and cetyltrimethylammonium bromide.
H318	Causes serious eye damage.
H315	Causes skin irritation.
H410	Very toxic to aquatic life with long lasting effects.
Prevention	
P264	Wash hands thoroughly after handling.
P280	Wear protective gloves / protective clothing / eye protection.
P273	Avoid release to the environment.
Response	
P305+P351+P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P310	Immediately call a POISON CENTER or doctor / physician.
P302+P352	IF ON SKIN: Wash with plenty of water.
P332+P313	If skin irritation occurs: Get medical advice / attention.
P362+P364	Take off contaminated clothing and wash it before reuse.
P391	Collect spillage.
Disposal	
P501	Dispose of contents / container in accordance with local regulations.

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STORAGE

- Controls are stable until the expiration date when stored and handled as directed.
- Do not use past expiration date.

PREPARATION FOR ANALYSIS

Controls may be used immediately after removal from 2-8°C storage. Prior to use, mix by gentle inversion (5-10 times).

After each use, tightly close the caps and return the controls to 2-8°C storage.

BIBLIOGRAPHY

- US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
- US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
- World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
- Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline—Third Edition*. CLSI Document M29-A3. Wayne, PA: CLSI; 2005.

Key to Symbols

	Caution
	Consult instructions for use
	Manufacturer
	Temperature limitation
	Use by/Expiration date
CONC	Concentration
CONTAINS: AZIDE	Contains Sodium Azide. Contact with acids liberates very toxic gas.
CONTROL -	Negative Control
CONTROL + 1	Positive Control 1
CONTROL + 2	Positive Control 2
IVD	In Vitro Diagnostic Medical Device
LOT	Lot Number
PRODUCED FOR ABBOTT BY	Produced for Abbott by
PRODUCT OF JAPAN	Product of Japan
RANGE	Range
REF	List Number
WARNING: SENSITIZER	Warning: May cause an allergic reaction.

ARCHITECT is a trademark of Abbott Laboratories in various jurisdictions.

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PRODUCED FOR ABBOTT BY

DENKA SEIKEN CO., LTD. Tokyo, Japan

Customer Service: Contact your local representative or find country-specific contact information on www.abbottdiagnostics.com

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Signed: _____

Date: _____

Michael Fritz,
Manager Regulatory Manager Regulatory Affairs
EMEA and Pakistan, International Regulatory
Affairs & Division Labeling

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Page 1 of 1

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Read Highlighted Changes: Revised November 2017.

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ARCHITECT rHTLV-I/II

INTENDED USE

The ARCHITECT rHTLV-I/II assay is a chemiluminescent microparticle immunoassay (CMIA) for the qualitative detection of antibodies to HTLV-I and HTLV-II in human serum and plasma including specimens collected post-mortem (non-heart-beating).

The ARCHITECT rHTLV-I/II assay is intended to be used as an aid in the diagnosis of HTLV-I and HTLV-II infection and as a screening test to prevent transmission of HTLV-I and HTLV-II to recipients of blood, blood components, cells, tissue and organs.

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HTLV-I and HTLV-II are closely related human type C retroviruses.¹ HTLV-I has been etiologically associated with neoplastic conditions and a variety of demyelinating neurologic disorders including: adult T-cell leukemia (ATL),³ tropical spastic paraparesis (TSP),⁴ and/or HTLV-I associated myelopathy (HAM)⁶ and more recently HTLV-I associated polymyositis, arthritis, and infective dermatitis.⁷⁻⁹ Association of HTLV-II with leukemia pathogenesis is not established; however, there is some evidence of an association with a neuro-degenerative disease similar to HAM/TSP¹⁰ and occasionally also with lymphoproliferative disease.¹¹

HTLV-I infection is endemic in south Japan¹², the Caribbean¹³, in some regions of Africa¹⁴, Central and South America¹⁵ and also found in Melanesia¹⁶ and central and northern Australia¹⁷, while HTLV-II is endemic to a number of indigenous American Indian populations.¹⁵ Both HTLV-I and HTLV-II are distributed worldwide among populations at high risk of infection, such as intravenous drug abusers, sex workers and patients attending sexually transmitted disease clinics.^{18, 19} Transmission of both HTLV-I and HTLV-II is by sexual contact, exposure to infected cellular blood components by transfusion or intravenous drug use or perinatally by breast feeding. Detection of antibodies against HTLV-I and HTLV-II serves to aid in the diagnosis of HTLV infection and to protect the safety of blood supply.

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The ARCHITECT rHTLV-I/II assay is a two-step immunoassay for the qualitative detection of antibodies to HTLV-I and HTLV-II in human serum and plasma using CMIA technology with flexible assay protocols, referred to as Chemiflex.

1. Sample and assay diluent are combined. An aliquot of the pretreated sample and HTLV-I/HTLV-II coated paramagnetic microparticles are combined in a new reaction vessel. Anti-HTLV-I/HTLV-II antibodies present in the sample bind to the HTLV-I/HTLV-II synthetic peptide and HTLV-II recombinant antigen coated microparticles.
2. After washing, acridinium-labeled conjugates (HTLV-I/HTLV-II synthetic peptide and HTLV-I recombinant antigen) bind to the anti-HTLV-I/HTLV-II antibodies.

3. Following another wash cycle, Pre-Trigger and Trigger Solutions are added to the reaction mixture.
4. The resulting chemiluminescent reaction is measured as relative light units (RLUs). There is a direct relationship between the amount of anti-HTLV-I/HTLV-II antibodies in the sample and the RLUs detected by the ARCHITECT iSystem optics.

The presence or absence of anti-HTLV-I/HTLV-II in the specimen is determined by comparing the chemiluminescent signal in the reaction to the cutoff signal determined from an active calibration. If the chemiluminescent signal in the reaction is greater than or equal to the cutoff signal, the specimen is considered reactive for anti-HTLV-I and/or anti-HTLV-II.

The immunodominant region of gp21 is 100% identical between HTLV-I and HTLV-II viral genome, and the overall homology between the two gp21 proteins is 86%. Additionally the homology between the two gp46 peptides is 65%, therefore the selected rare reagents (synthetic peptides (gp46) and recombinant antigens (gp21)) in this assay are capable to detect antibodies for HTLV-I and HTLV-II simultaneously.

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Kit Contents

ARCHITECT rHTLV-I/II 6L61

NOTE: Some kit sizes are not available in all countries or for use on all ARCHITECT iSystems. Please contact your local distributor.

REF	6L61-25	6L61-35	6L61-30
	100	500	2000
MICROPARTICLES	1 x 6.6 mL	1 x 27.0 mL	4 x 27.0 mL
CONJUGATE	1 x 5.9 mL	1 x 26.3 mL	4 x 26.3 mL
ASSAY DILUENT	1 x 7.7 mL	1 x 39.2 mL	4 x 39.2 mL

MICROPARTICLES HTLV-I/HTLV-II synthetic peptide and HTLV-II recombinant antigen coated microparticles in TRIS buffered saline. Minimum concentration: 0.08% solids. Preservatives: ProClin 950 and sodium azide.

CONJUGATE Acridinium-labeled HTLV-I/HTLV-II synthetic peptide (maximum concentration: 100 ng/mL each) and acridinium-labeled HTLV-I recombinant antigen (maximum concentration: 10 ng/mL) in HEPES buffer with protein stabilizers and detergent. Preservative: ProClin 950.

ASSAY DILUENT rHTLV-I/II assay diluent containing TRIS buffer and detergent. Preservatives: ProClin 950 and sodium azide.

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- **IVD**
- For In Vitro Diagnostic Use

Safety Precautions

CAUTION: This product requires the handling of human specimens. It is recommended that all human-sourced materials be considered potentially infectious and handled in accordance with the OSHA Standard on Bloodborne Pathogens. Biosafety Level 2 or other appropriate biosafety practices should be used for materials that contain or are suspected of containing infectious agents.²⁰⁻²³

The following warnings and precautions apply to: **MICROPARTICLES**

WARNING	Contain methylisothiazolone and sodium azide.
H317	May cause an allergic skin reaction.
EUH032	Contact with acids liberates very toxic gas.

Prevention

P261	Avoid breathing mist / vapors / spray.
P272	Contaminated work clothing should not be allowed out of the workplace.
P280	Wear protective gloves / protective clothing / eye protection.

Response

P302+P352	IF ON SKIN: Wash with plenty of water.
P333+P313	If skin irritation or rash occurs: Get medical advice / attention.
P362+P364	Take off contaminated clothing and wash it before reuse.

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WARNING	Contains methylisothiazolone.
H317	May cause an allergic skin reaction.

Prevention

P261	Avoid breathing mist / vapors / spray.
P272	Contaminated work clothing should not be allowed out of the workplace.
P280	Wear protective gloves / protective clothing / eye protection.

Response

P302+P352	IF ON SKIN: Wash with plenty of water.
P333+P313	If skin irritation or rash occurs: Get medical advice / attention.
P362+P364	Take off contaminated clothing and wash it before reuse.

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The following warnings and precautions apply to: **ASSAY DILUENT**

WARNING	Contains polyethylene glycol octylphenyl ether (Triton X-100), methylisothiazolone and sodium azide.
----------------	--

H317	May cause an allergic skin reaction.
------	--------------------------------------

H319	Causes serious eye irritation.
------	--------------------------------

EUH032	Contact with acids liberates very toxic gas.
--------	--

Prevention

P261	Avoid breathing mist / vapors / spray.
P264	Wash hands thoroughly after handling.
P272	Contaminated work clothing should not be allowed out of the workplace.
P280	Wear protective gloves / protective clothing / eye protection.

Response

P305+P351+P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
----------------	--

P337+P313	If eye irritation persists: Get medical advice / attention.
-----------	---

P302+P352	IF ON SKIN: Wash with plenty of water.
-----------	--

P333+P313	If skin irritation or rash occurs: Get medical advice / attention.
-----------	--

P362+P364	Take off contaminated clothing and wash it before reuse.
-----------	--

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Reagent Handling

- Do not use reagent kits beyond the expiration date.
- **Do not pool reagents within a kit or between kits.**
- Before loading the reagent kit on the system for the first time, the microparticle bottle requires mixing to resuspend microparticles that may have settled during shipment. For microparticle mixing instructions, refer to the **PROCEDURE, Assay Procedure** section of this package insert.
- **Septums MUST be used to prevent reagent evaporation and contamination and to ensure reagent integrity. Reliability of assay results cannot be guaranteed if septums are not used according to the instructions in this package insert.**
 - To avoid contamination, wear clean gloves when placing a septum on an uncapped reagent bottle.
 - Once a septum has been placed on an open reagent bottle, **do not invert the bottle** as this will result in reagent leakage and may compromise assay results.
 - Over time, residual liquids may dry on the septum surface. These are typically dried salts and have no effect on assay efficacy.

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	Storage Temperature	Maximum Storage Time	Additional Storage Instructions
Unopened/Opened*	2-8°C	Until expiration date	May be used immediately after removal from 2-8°C storage. Store in upright position.
On board	System temperature	30 days	Discard after 30 days. For information on tracking onboard time, refer to the ARCHITECT System Operations Manual, Section 5.

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INSTRUMENT PROCEDURE

The ARCHITECT rHTLV-1/II assay file must be installed on the ARCHITECT iSystem from an ARCHITECT iSystem Assay CD-ROM prior to performing the assay.

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- Other specimen collection tube types have not been tested with this assay.
- Liquid anticoagulants may have a dilution effect resulting in lower concentrations for individual patient specimens.

- The instrument does not provide the capability to verify specimen type. It is the responsibility of the operator to verify that the correct specimen types are used in the assay.
- Performance has been established for the use of cadaveric blood specimens (specimens collected post-mortem, non-heart-beating), for details refer to section **Testing of Cadaveric Blood Specimens**.

Specimen Conditions

- Do not use specimens with the following conditions:
 - heat-inactivated
 - pooled
 - grossly hemolyzed (> 500 mg/dL hemoglobin)
 - obvious microbial contamination
 - body fluids other than human serum and plasma
- For accurate results, serum and plasma specimens should be free of fibrin, red blood cells, and other particulate matter. Serum specimens from patients receiving anticoagulant or thrombolytic therapy may contain fibrin due to incomplete clot formation.
- To prevent cross contamination, use of disposable pipettes or pipette tips is recommended.
- No qualitative performance differences were observed between experimental controls and a minimum of 10 nonreactive and 10 spiked reactive specimens tested with elevated levels of bilirubin (20 mg/dL), triglycerides (3000 mg/dL), protein (12 g/dL), red blood cells (0.4% v/v), or hemoglobin (500 mg/dL).

Preparation for Analysis

- Follow the tube manufacturer's processing instructions for collection tubes. Gravity separation is not sufficient for specimen preparation.
- Mix thawed specimens thoroughly by low speed vortexing or by inverting 10 times. Visually inspect the specimens. If layering or stratification is observed, continue mixing until specimens are visibly homogeneous.
- To ensure consistency in results, specimens must be transferred to a centrifuge tube and centrifuged at $\geq 10,000$ RCF (Relative Centrifugal Force) for 10 minutes before testing if
 - they contain fibrin, red blood cells, or other particulate matter,
 - they require repeat testing, or
 - they were frozen and thawed.Transfer clarified specimen to a sample cup or secondary tube for testing.
- Centrifuged specimens with a lipid layer on the top must be transferred to a sample cup or secondary tube. Care must be taken to transfer only the clarified specimen without the lipemic material.
- Inspect all specimens for bubbles. Remove bubbles with an applicator stick before analysis. Use a new applicator stick for each specimen to prevent cross contamination.

Specimen Storage

Specimen Type	Storage Temperature	Maximum Storage Time
Serum/Plasma	15-30°C	≤ 3 days
	2-8°C	≤ 14 days
	-20°C or colder	∞

Specimens may be stored on or off the clot, red blood cells, or separator gel.

If testing will be delayed more than the recommended storage time, remove serum or plasma from the clot or red blood cells, and store frozen (-20°C or colder).

No qualitative performance differences were observed between experimental controls and a minimum of 10 nonreactive or spiked reactive specimens subjected to 6 freeze/thaw cycles; however, multiple freeze/thaw cycles should be avoided.

Specimen Shipping

- Package and label specimens in compliance with applicable state, federal, and international regulations covering the transport of clinical specimens and infectious substances.
- It is recommended that specimens be removed from the clot, red blood cells, or separator gel.
- Specimens may be shipped ambient.
- Do not exceed the storage time limitations listed above.

Testing of Cadaveric Blood Specimens

- Performance has been established for the use of cadaveric blood specimens (specimens collected post-mortem, non-heart-beating) that have been collected up to 18.5 hours after death. Performance was established using 50 spiked and 50 non-spiked cadaveric blood specimens. 
- Testing of cadaveric blood specimens from patients with plasma dilution due to transfusions of > 2000 mL of blood or colloids within 48 hours, or > 2000 mL of crystalloids within 1 hour (or any combination thereof) prior to collection of the specimens have not been validated.
- Follow general standards and/or regulations for collection, storage and handling.
- Follow the tube manufacturer's processing instructions for serum or plasma collection tubes. After initial centrifugation, transfer the supernatant to a centrifuge tube and centrifuge at 10,000 RCF (Relative Centrifugal Force) for 10 minutes. If specimens are not processed directly after initial centrifugation, it is recommended to remove the supernatant from the clot, red blood cells or separator gel until further processing.
- Cadaveric blood specimens can be stored for up to 7 days at 2-8°C or up to 3 days at 15-30°C following collection.
- No qualitative differences were observed for cadaveric blood specimens (nonreactive or spiked reactive) when subjected to up to 3 freeze/thaw cycles. However, multiple freeze/thaw cycles should be avoided.

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- ARCHITECT rHTLV-I/II Assay file obtained from the ARCHITECT iSystem e-Assay CD-ROM found on www.abbottdiagnostics.com.
- 6L61-01 ARCHITECT rHTLV-I/II Calibrator
- 6L61-10 ARCHITECT rHTLV-I/II Controls
- ARCHITECT Pre-Trigger Solution
- ARCHITECT Trigger Solution
- ARCHITECT Wash Buffer
- ARCHITECT Reaction Vessels
- ARCHITECT Sample Cups
- ARCHITECT Septum
- ARCHITECT Replacement Caps
- Pipettes or pipette tips (optional) to deliver the volumes specified on the patient or control order screen.

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Assay Procedure

- Before loading the reagent kit on the system for the first time, the microparticle bottle requires mixing to resuspend microparticles that may have settled during shipment. After the first time the microparticles have been loaded, no further mixing is required.
 - **Invert the microparticle bottle 30 times.**
 - Visually inspect the bottle to ensure microparticles are resuspended. If microparticles are still adhered to the bottle, continue to invert the bottle until the microparticles have been completely resuspended.
 - **If the microparticles do not resuspend, DO NOT USE. Contact your local Abbott representative.**
 - Once the microparticles have been resuspended, discard the cap and place a septum on the bottle. For instructions about placing septums on bottles, refer to the **Reagent Handling** section of this package insert.
- Load the reagent kit on the ARCHITECT iSystem.
 - Verify that all necessary reagents are present.
 - Ensure that septums are present on all reagent bottles.
- Order calibration, if necessary.
 - For information on ordering calibrations, refer to the ARCHITECT System Operations Manual, Section 6.
- Order tests.
 - For information on ordering patient specimens and controls and for general operating procedures, refer to the ARCHITECT System Operations Manual, Section 5.
- Minimum sample cup volume is calculated by the system and printed on the Orderlist report. To minimize the effects of evaporation, verify adequate sample cup volume is present prior to running the test.

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 - Priority:

Sample volume for first test: 183 µL

Sample volume for each additional test from same sample cup: 133 µL
 - ≤ 3 hours on board:

Sample volume for first test: 183 µL

Sample volume for each additional test from same sample cup: 133 µL
 - > 3 hours on board: Additional sample volume required. For information on sample evaporation and volumes, refer to the ARCHITECT System Operations Manual, Section 5.
 - If using primary or aliquot tubes, use the sample gauge to ensure sufficient patient specimen is present.
- Prepare ARCHITECT rHTLV-I/II Calibrator and Controls.
 - Mix calibrator(s) and controls by gentle inversion before use.
 - Hold bottles **vertically** and dispense recommended volumes into each respective sample cup.
 - Recommended volumes:

for each calibrator: 20 drops

for each control: 10 drops
- Load samples.
 - For information on loading samples, refer to the ARCHITECT System Operations Manual, Section 5.
- Press RUN.
- For additional information on principles of operation, refer to the ARCHITECT System Operations Manual, Section 3.
- For optimal performance, it is important to perform routine maintenance as described in the ARCHITECT System Operations Manual, Section 9. Perform maintenance more frequently when required by laboratory procedures.

Specimen Dilution Procedures

Specimens cannot be diluted for the ARCHITECT rHTLV-I/II assay.

Calibration

- Test Calibrator 1 in replicates of three. The calibrator should be priority loaded.
A single sample of each control level must be tested to evaluate the assay calibration. Ensure that assay control values are within the ranges specified in the respective control package insert.
- Once an ARCHITECT rHTLV-I/II calibration is accepted and stored, all subsequent samples may be tested without further calibration unless:
 - A reagent kit with a new lot number is used or
 - Controls are out of range.
- For detailed information on how to perform an assay calibration, refer to the ARCHITECT System Operations Manual, Section 6.

Quality Control Procedures

The recommended control requirement for the ARCHITECT rHTLV-I/II assay is that a single sample of each control be tested once every 24 hours each day of use. If laboratory quality control procedures require more frequent use of controls to verify test results, follow those procedures.

Each laboratory should establish control ranges to monitor the acceptable performance of the assay. The ARCHITECT rHTLV-I/II Control values must be within the acceptable ranges specified in the control package insert. If a control is out of its specified range, the associated test results are invalid and must be retested. Recalibration may be indicated.

Verification of Assay Claims

For protocols to verify package insert claims, refer to the ARCHITECT System Operations Manual, Appendix B.

The ARCHITECT rHTLV-I/II assay belongs to method group 5, except functional sensitivity.

RESULTS

The ARCHITECT iSystem calculates the cutoff (CO) using the mean chemiluminescent signal (RLU) from three replicates of the Calibrator 1 multiplied by 0.25 and stores the result.

Calculation

The ARCHITECT rHTLV-I/II assay calculates results based on a cutoff determined by the following calculation.

- Cutoff RLU = Calibrator 1 mean RLU value x 0.25
- S/CO = Sample RLU/Cutoff RLU
- The cutoff RLU is stored for each reagent lot calibration.
- The ARCHITECT iSystem calculates a result based on the ratio of the sample RLU(s) to the cutoff RLU for each specimen and control.

Interpretation of Results

- Specimens with S/CO values < 1.00 are considered nonreactive (NR).
- Specimens with S/CO values ≥ 1.00 are considered reactive (R).

NOTE: All specimens that are initially reactive must be centrifuged and retested in duplicate. Refer to the **SPECIMEN COLLECTION AND PREPARATION FOR ANALYSIS** section in this package insert.

ARCHITECT rHTLV-I/II Final Interpretation

Initial Result (S/CO)	Retest Results (S/CO)	Final Result
NR (< 1.00)	No retest required (NA)	NR
R (≥ 1.00)	Both tests are NR (< 1.00)	NR
R (≥ 1.00)	One or both tests are R (≥ 1.00)	R

For details on configuring the ARCHITECT iSystem to use grayzone interpretations, refer to the ARCHITECT System Operations Manual, Section 2. The grayzone interpretation from the ARCHITECT interpretations screen is not used by the ARCHITECT System unless a grayzone is configured.

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LIMITATIONS OF THE PROCEDURE

- If the assay results are inconsistent with clinical evidence, additional testing is suggested to confirm the result.
- Heterophilic antibodies in human serum can react with reagent immunoglobulins, interfering with *in vitro* immunoassays. Patients routinely exposed to animals or to animal serum products can be prone to this interference, and anomalous values may be observed. Additional information may be required for diagnosis.²⁵

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Precision

The ARCHITECT rHTLV-I/II assay is designed to have an imprecision of ≤ 10% total** CV for specimens within the range of 1.00 to 6.00 S/CO. The study was performed at one internal and two external evaluation sites each using one instrument. A panel consisting of three different control lots and two human plasma specimens was tested in replicates of four across three reagent lots and one calibrator lot per site. A total of three different calibrator lots were used. Each combination of instruments, panel members, and reagent lots was tested in four runs. Data from this study are summarized in Table 1*.

Table 1: ARCHITECT rHTLV-I/II Precision

Panel Member	n	Mean			Within Run		Total **	
		S/CO	SD	%CV	SD	%CV		
Negative Control	432	0.10	0.02	23.61	0.03	25.44		
Positive Control	432	4.05	0.15	3.80	0.17	4.31		
Human Plasma HTLV I	144	1.77	0.06	3.31	0.07	3.72		
Human Plasma HTLV II	144	1.82	0.07	3.66	0.07	3.98		

* Representative data; results in individual laboratories may vary from these data. Higher imprecision may be observed in high reactive samples without changes in qualitative results.

** Total is an accumulation of within run, between run and between day.

Specificity

The ARCHITECT rHTLV-I/II assay is designed to have a specificity of ≥ 99.5% on a blood donor population and ≥ 99.0% on a hospitalized/diagnostic population. A study was performed at two external sites and one internal site on a total of 5646 serum and plasma specimens collected from six blood-donation centers and 692 hospitalized/diagnostic specimens. Data from this study are summarized in Table 2*.

Table 2: ARCHITECT rHTLV-I/II Specificity

Category	n	IR [%]	RR [%]	Resolved Specificity	95% Confidence Interval
Plasma	3583	2 [0.06]	2 [0.06]	99.94% (3581/3583)	99.80% - 99.99%
Serum	2063	2 [0.10]	1 [0.05]	99.95% (2062/2063)	99.73% - 100%
Overall Blood Donors	5646	4 [0.07]	3 [0.05]	99.95% (5643/5646)	99.84% - 99.99%
HP/Diagnostics	692	1 [0.14]	1 [0.14]	99.86% (691/692)	99.20% - 100%

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Sensitivity

The ARCHITECT rHTLV-I/II assay is designed to have a sensitivity equivalent or greater than a commercially available diagnostic kit. A total of 301 anti-HTLV-I and 105 anti-HTLV-II reactive specimens were tested, resulting in a sensitivity of 100% (406/406), 95% confidence interval: 99.10% - 100%. (Representative data; results in individual laboratories may vary from these data). These values were determined in a study where pedigreed specimens from individuals with HTLV infection were tested.

Interference

Additional studies were performed to evaluate other potential interfering disease states on the ARCHITECT rHTLV-I/II assay. A total of 260 specimens were tested from following categories: Viral infection (HBV, HSV, CMV, HCV, EBV, HIV-1, HIV-2); fungal/yeast/protozoal/bacterial infection (*T. pallidum*, *T. gondii*, *E. coli*, *C. trachomatis*, *N. gonorrhoea*); autoimmune (rheumatoid factor [RF], antinuclear antibodies [ANA]), other conditions (pregnant females all trimesters, multiparous females, elevated IgG, elevated IgM, monoclonal and polyclonal gammopathy, influenza vaccine recipients, hemodialysis patients, hemophiliacs, multiple transfusion recipients). Data from this study are summarized in Table 3^a.

Table 3: ARCHITECT rHTLV-I/II Potential Interfering Substances

Interfering Substances	n	IR n	RR n	Apparent Specificity	95% Confidence Interval
Total	260	1 ^a	1 ^a	100% (259/259)	98.59% - 100%

^a One multiple transfusion sample initially and repeatedly reactive on ARCHITECT rHTLV-I/II was confirmed positive with an alternative enzyme immunoassay and a Western Blot.

• Representative data; results in individual laboratories may vary from these data.

BIBLIOGRAPHY

- Poiesz BJ, Ruscetti FW, Gazdar AF, et al. Detection and Isolation of Type C Retrovirus Particles from Fresh and Cultured Lymphocytes of a Patient with Cutaneous T-cell Lymphoma. *Proc Natl Acad Sci USA* 1980; 77:7415.
- Poiesz BJ, Ruscetti FW, Reitz MS, et al. Isolation of a New Type C Retrovirus (HTLV) in Primary Uncultured Cells of a Patient with Sezary T-cell Leukemia. *Nature (London)* 1981; 294:268.
- Blattner WA, Takatsuki K, Gallo RC. Human T-Cell Leukemia-Lymphoma Virus and Adult T-Cell Leukemia. *JAMA* 1983; 250:1074.
- Gessain A, Barin F, Vermant JC, et al. Antibodies to Human T-Lymphotropic Virus Type-I in Patients with Tropical Spastic Paraparesis. *Lancet* 1985; ii:407.
- Rodgers-Johnson P, Gajdusek DC, Morgan OSC, et al. HTLV-I and HTLV-II Antibodies and Tropical Spastic Paraparesis. *Lancet* 1985; ii:1247.
- Osame M, Usuku K, Izumo S, et al. HTLV-I Associated Myelopathy. A New Clinical Entity. *Lancet* 1986; i:1031.
- Gessain A, Gout O. Chronic Myelopathy Associated with Human T-Lymphotropic Virus Type-I (HTLV-I). *Ann Intern Med* 1992; 117(11):933.
- Wiley CA, Nerenberg M, Cros D, et al. HTLV-I Polymyositis in a Patient also Infected with the Human Immunodeficiency Virus. *N Engl J Med* 1989; 320 (15): 992.
- Blattner WA. Epidemiology of HTLV-I and Associated Diseases. In: Blattner WA, ed. *Human Retrovirology: HTLV*. New York: Raven Press, 1990; 251.
- Rosenblatt JD, Golde DW, et al. A Second Isolate of HTLV-II Associated with Atypical Hairy-cell Leukemia. *N Engl J Med* 1986; 315:372.
- Jacobson S, Lehky T, et al. Isolation of HTLV-II from a Patient with Chronic, Progressive Neurological Disease Clinically Indistinguishable from HTLV-I Associated Myelopathy/tropical Spastic Paraparesis. *Ann Neurol*. 1993; 33:392.
- The T- and B-Cell Malignancy Study Group. Statistical Analysis of Clinico-pathological, Virological and Epidemiological Data on Lymphoid Malignancies with Special Reference to Adult T-Cell Leukemia/ Lymphoma: A report of the Second Nationwide Study of Japan. *Jap. J Clin Oncol* 1985; 1593:517.
- Newton M, Cruickshank K, Miller D. et al. Antibody to Human T Lymphotropic Virus Type I in West-Indian Born UK Residents with Spastic Paraparesis. *Lancet* 1987; i:415.
- Saxinger W, Blattner WA, Levine PH, et al. Human T-Cell Leukemia Virus (HTLV-I) Antibodies in Africa. *Science* 1984; 225:1473.
- Maloney EM, Biggar RJ, Neel JV, et al. Endemic Human T-cell Lymphotropic Virus Type II Infection Among Isolated Brazilian Amerindians. *J Infect Dis* 1992; 166:100.
- Yanagihara R, Jenkins CL, Alexander SS, More CA, et al. Human T-lymphotropic Virus Type I Infection in Papua New Guinea: High Prevalence Among the Hagahai Confirmed by Western Blot Analysis. *J. Infect. Dis* 1990; 162:649.
- Bastian I, Gardner J, Webb D, Gardner I. Isolation of Human T-Lymphotropic Virus Type-I Strain from Australian Aborigines. *J. Virol.* 1993; 67(2):843.
- Wiktor SZ, Cannon RO, et al. Infection with Human T-Lymphotropic Virus Type I and II in Sexually Transmitted Disease Clinics in Baltimore and New Orleans. *J Infect Dis* 1992; 165:920.
- Zella D, Mori L, et al. HTLV-II Infection in Italian Drug Abusers. *Lancet* 1990; 336:575.
- US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
- US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
- World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
- Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4, Wayne, PA: CLSI; 2014.
- U.S. Department of Health and Human Services, Food and Drug Administration, Center for Biologics Evaluation and Research. Guidance for Industry Recommendations for Obtaining a Labeling Claim for Communicable Disease Donor Screening Tests Using Cadaveric Blood Specimens from Donors of Human Cells, Tissues, and Cellular and Tissue-Based Products (HCT/PS), November 2004. <http://www.fda.gov/BiologicsBloodVaccines/GuidanceComplianceRegulatoryInformation/Guidances/Tissue/ucm073972.htm> Accessed October 24, 2017
- Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. *Clin Chem* 1988;34(1):27-33.

Key to Symbols

	Consult instructions for use
	Manufacturer
	Sufficient for
	Temperature limitation
	Use by/Expiration date
ASSAY DILUENT	Assay Diluent
CONJUGATE	Conjugate
CONTAINS: AZIDE	Contains Sodium Azide. Contact with acids liberates very toxic gas.
CONTROL NO.	Control Number
IVD	<i>In Vitro</i> Diagnostic Medical Device
LOT	Lot Number
MICROPARTICLES	Microparticles
PRE-TRIGGER SOLUTION	Pre-Trigger Solution
PRODUCT OF GERMANY	Product of Germany
REACTION VESSELS	Reaction Vessels
REAGENT LOT	Reagent Lot
REF	List Number
REPLACEMENT CAPS	Replacement Caps
SAMPLE CUPS	Sample Cups
SEPTUM	Septum
SN	Serial number
TRIGGER SOLUTION	Trigger Solution
WARNING: SENSITIZER	Warning: May cause an allergic reaction.
WASH BUFFER	Wash Buffer

ARCHITECT and Chemiflex are trademarks of Abbott Laboratories in various jurisdictions.

ProClin is property of its respective owner.

Abbott GmbH & Co. KG
Max-Planck-Ring 2
65205 Wiesbaden
Germany
+49-6122-580

0088

Customer Service: Contact your local representative
or find country-specific contact information on
www.abbottdiagnostics.com

Revised November 2017.

©2008, 2017 Abbott Laboratories

ADDENDUM

to the Instruction for Use ARCHITECT rHTLV-I/II Reagent for Russia

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REPAIRS

Not required.

DESTRUCTION AND DISPOSAL OF A MEDICAL DEVICE

Destruct and dispose of a medical device in accordance with the requirements of SanPin 2.1.7.2790-10 dated December 12, 2010 "Sanitary-epidemiological Requirements for the Management of Medical Waste".

Transmit and dispose of a product and its components and packaging in accordance with local regulations.

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Medical device is intended for one-off usage.

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- Poiesz BJ, Ruscetti FW, Gazdar AF, et al. Detection and Isolation of Type C Retrovirus Particles from Fresh and Cultured Lymphocytes of a Patient with Cutaneous T-cell Lymphoma. *Proc Natl Aca Sci USA* 1980; 77:7415.
- Poiesz BJ, Ruscetti FW, Reitz MS, et al. Isolation of a New Type C Retrovirus (HTLV) in Primary Uncultured Cells of a Patient with Sezary T-cell Leukemia. *Nature (London)* 1981; 294:268.
- Blattner WA, Takatsuki K, Gallo RC. Human T-Cell Leukemia-Lymphoma Virus and Adult T-Cell Leukemia. *JAMA* 1983; 250:1074.
- Gessain A, Barin F, Vermant JC, et al. Antibodies to Human T-Lymphotropic Virus Type-I in Patients with Tropical Spastic Paraparesis. *Lancet* 1985; ii:407.
- Rodgers-Johnson P, Gajdusek DC, Morgan OSC, et al. HTLV-I and HTLV-II Antibodies and Tropical Spastic Paraparesis. *Lancet* 1985; ii:1247.
- Osame M, Usuku K, Izumo S, et al. HTLV-I Associated Myelopathy, A New Clinical Entity. *Lancet* 1986; i:1031.
- Gessain A, Gout O. Chronic Myelopathy Associated with Human T-Lymphotropic Virus Type-I (HTLV-I). *Ann Intern Med* 1992; 117(11):933.
- Wiley CA, Nerenberg M, Cros D, et al. HTLV-I Polymyositis in a Patient also Infected with the Human Immunodeficiency Virus. *N Engl J Med* 1989; 320 (15): 992.
- Blattner WA. Epidemiology of HTLV-I and Associated Diseases. In: Blattner WA, ed. *Human Retrovirology: HTLV*. New York: Raven Press, 1990; 251.
- Rosenblatt JD, Golde DW, et al. A Second Isolate of HTLV-II Associated with Atypical Hairy-cell Leukemia. *N Engl J Med* 1986, 315:372.
- Jacobson S, Lehky T, et al. Isolation of HTLV-II from a Patient with Chronic, Progressive Neurological Disease Clinically Indistinguishable from HTLV-I Associated Myelopathy/tropical Spastic Paraparesis. *Ann Neurol*. 1993; 33:392.
- The T- and B-Cell Malignancy Study Group. Statistical Analysis of Clinico-pathological, Virological and Epidemiological Data on Lymphoid Malignancies with Special Reference to Adult T-Cell Leukemia/ Lymphoma: A report of the Second Nationwide Study of Japan. *Jap. J Clin Oncol* 1985; 1593:517.

13. Newton M, Cruickshank K, Miller D. et al. Antibody to Human T Lymphotropic Virus Type I in West-Indian Born UK Residents with Spastic Paraparesis. *Lancet* 1987; 1:415.
14. Saxinger W, Blattner WA, Levine PH, et al. Human T-Cell Leukemia Virus (HTLV-I) Antibodies in Africa. *Science* 1984; 225:1473.
15. Maloney EM, Biggar RJ, Neel JV, et al. Endemic Human T-cell Lymphotropic Virus Type II Infection Among Isolated Brazilian Amerindians. *J Infect Dis* 1992; 166:100.
16. Yanagihara R, Jenkins CL, Alexander SS, More CA, et al. Human T-lymphotropic Virus Type I Infection in Papua New Guinea: High Prevalence Among the Hagahai Confirmed by Western Blot Analysis. *J. Infect. Dis* 1990; 162:649.
17. Bastian I, Gardner J, Webb D, Gardner I. Isolation of Human T-Lymphotropic Virus Type-I Strain from Australian Aborigines. *J. Virol.* 1993; 67(2):843.
18. Wiktor SZ, Cannon RO, et al. Infection with Human T-Lymphotropic Virus Type I and II in Sexually Transmitted Disease Clinics in Baltimore and New Orleans. *J Infect Dis* 1992; 165:920.
19. Zella D, Mori L et al. HTLV-II Infection in Italian Drug Abusers. *Lancet* 1990, 336:575.
20. US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
21. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
22. World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
23. Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.
24. U.S. Department of Health and Human Services, Food and Drug Administration, Center for Biologics Evaluation and Research. Guidance for Industry Recommendations for Obtaining a Labeling Claim for Communicable Disease Donor Screening Tests Using Cadaveric Blood Specimens from Donors of Human Cells, Tissues, and Cellular and Tissue-Based Products (HCT/PS), November 2004. <http://www.fda.gov/BiologicsBloodVaccines/GuidanceComplianceRegulatoryInformation/Guidances/Tissue/ucm073972.htm> По состоянию на: 24 октября 2017 г.
25. Boscatto LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. *Clin Chem* 1988;34(1):27-33.

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INTENDED USE

The ARCHITECT rHTLV-I/II Calibrator is for the calibration of the ARCHITECT iSystem when used for the qualitative detection of antibodies to HTLV-I and HTLV-II in human serum and plasma. Refer to the ARCHITECT rHTLV-I/II reagent insert and the ARCHITECT System Operations Manual for additional information.

CONTENTS

1 Bottle (4.0 mL) of ARCHITECT rHTLV-I/II Calibrator. Calibrator 1 contains recalcified human plasma and is reactive for anti-HTLV. Preservatives: ProClin 950 and sodium azide.

Calibrator	Color	Target Value
CAL 1	Red*	N/A

* Dye: Red D&C No. 33

STANDARDIZATION

The ARCHITECT rHTLV-I/II Calibrator 1 is standardized against internal reference material.

PRECAUTIONS

- IVD
- For In Vitro Diagnostic Use

-  **CAUTION:** This product contains human-sourced and/or potentially infectious components. Refer to the **CONTENTS** section of this package insert. No known test method can offer complete assurance that products derived from human sources or inactivated microorganisms will not transmit infection. Therefore, all human-sourced materials should be considered potentially infectious. It is recommended that these reagents and human specimens be handled in accordance with the OSHA Standard on Bloodborne Pathogens. Biosafety Level 2 or other appropriate biosafety practices should be used for materials that contain or are suspected of containing infectious agents.¹⁻⁴
- The human plasma used in Calibrator 1 is reactive for anti-HTLV and is nonreactive for HBsAg, anti-HIV-1/HIV-2, anti-HCV and HIV-1 RNA or HIV-1 Ag.

The following warnings and precautions apply to: CAL 1	
	
WARNING	Contains methylisothiazolone and sodium azide.
H317	May cause an allergic skin reaction.
EUH032	Contact with acids liberates very toxic gas.
Prevention	
P261	Avoid breathing mist / vapors / spray.
P272	Contaminated work clothing should not be allowed out of the workplace.
P280	Wear protective gloves / protective clothing / eye protection.


Response	
P302+P352	IF ON SKIN: Wash with plenty of water.
P333+P313	If skin irritation or rash occurs: Get medical advice / attention.
P362+P364	Take off contaminated clothing and wash it before reuse.
Disposal	
P501	Dispose of contents / container in accordance with local regulations.

Safety Data Sheets are available at www.abbottiagnostics.com or contact your local representative.

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STORAGE

- Calibrator is stable until the expiration date when stored and handled as directed.
- Do not use past expiration date.

PREPARATION FOR ANALYSIS

Calibrator may be used immediately after removal from 2-8°C storage.

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After each use, tightly close the cap and return the calibrator to 2-8°C storage.

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Key to Symbols

	Caution
	Consult instructions for use
	Manufacturer
	Temperature limitation
	Use by/Expiration date
CAL 1	Calibrator 1
CONTAINS: AZIDE	Contains Sodium Azide. Contact with acids liberates very toxic gas.
IVD	<i>In Vitro</i> Diagnostic Medical Device
LOT	Lot Number
PRODUCT OF GERMANY	Product of Germany
REF	List Number
WARNING: SENSITIZER	Warning: May cause an allergic reaction.

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©2008, 2014 Abbott Laboratories

ADDENDUM

to the Instruction for Use ARCHITECT rHTLV-I/II Calibrators for Russia

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en

rHTLV-I/II
6L61
G5-1059 / R04
C6L610

Read Highlighted Changes: Revised November 2014.

INTENDED USE

The ARCHITECT rHTLV-I/II Controls are for the estimation of test precision and the detection of systematic analytical deviations of the ARCHITECT iSystem when used for the qualitative detection of antibodies to HTLV-I and HTLV-II in human serum and plasma. Refer to the ARCHITECT rHTLV-I/II reagent insert and the ARCHITECT System Operations Manual for additional information.

CONTENTS

2 Bottles (8.0 mL each) of ARCHITECT rHTLV-I/II Controls: Negative Control and Positive Control. The Negative Control contains recalcified human plasma. The Positive Control contains recalcified human plasma and is reactive for anti-HTLV. Preservatives: ProClin 950 and sodium azide.

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PRECAUTIONS

- **IVD**
- For In Vitro Diagnostic Use

-  **CAUTION:** This product contains human-sourced and/or potentially infectious components. Refer to the **CONTENTS** section of this package insert. No known test method can offer complete assurance that products derived from human sources or inactivated microorganisms will not transmit infection. Therefore, all human-sourced materials should be considered potentially infectious. It is recommended that these reagents and human specimens be handled in accordance with the OSHA Standard on Bloodborne Pathogens. Biosafety Level 2 or other appropriate biosafety practices should be used for materials that contain or are suspected of containing infectious agents.¹⁻⁴
- The human plasma used in the Negative Control is nonreactive for anti-HTLV-I/II, HBsAg, anti-HIV-1/HIV-2, anti-HCV and HIV-1 RNA or HIV-1 Ag.
- The human plasma used in the Positive Control is reactive for anti-HTLV and nonreactive for HBsAg, anti-HIV-1/HIV-2, anti-HCV and HIV-1 RNA or HIV-1 Ag.

The following warnings and precautions apply to: **CONTROL -** / **CONTROL +**

WARNING	Contains methylisothiazolone and sodium azide.
H317	May cause an allergic skin reaction.
EUH032	Contact with acids liberates very toxic gas.

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PREPARATION FOR ANALYSIS

Controls may be used immediately after removal from 2-8°C storage.

Prior to each use, mix by gentle inversion.

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BIBLIOGRAPHY

1. US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
2. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
3. World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
4. Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline—Third Edition*. CLSI Document M29-A3. Wayne, PA: CLSI; 2005.

Key to Symbols

	Caution
	Consult instructions for use
	Manufacturer
	Temperature limitation
	Use by/Expiration date
CONTAINS: AZIDE	Contains Sodium Azide. Contact with acids liberates very toxic gas.
CONTROL -	Negative Control
CONTROL +	Positive Control
IVD	<i>In Vitro</i> Diagnostic Medical Device
LOT	Lot Number
PRODUCT OF GERMANY	Product of Germany
RANGE	Range
REF	List Number
WARNING: SENSITIZER	Warning: May cause an allergic reaction.

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Revised November 2014.
©2008, 2014 Abbott Laboratories

ADDENDUM

to the Instruction for Use ARCHITECT rHTLV-I/II Controls for Russia

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REF 6C20-25
REF 6C20-35

Toxo IgM
6C20
G4-7818/ R03
B6C200

REAGENTS
Kit Contents
ARCHITECT Toxo IgM 6C
NOTE: Some kit sizes are

Read Highlighted Changes: Revised October 2014.

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ARCHITECT Toxo IgM

INTENDED USE

The ARCHITECT Toxo IgM assay is a chemiluminescent microparticle immunoassay (CMIA) for the qualitative detection of IgM antibodies to *Toxoplasma gondii* in human serum and plasma.

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Toxoplasma gondii is an obligate intracellular protozoan parasite that infects most species of warm-blooded animals, including humans.¹ Toxoplasmosis is primarily acquired by ingestion of undercooked, infected meat; via oocysts from fecally contaminated hands, food and water; and maternally through transplacental transmission.² In addition, transmission associated with organ transplantation and during blood transfusion has been reported, although the risk of transmission through blood transfusion is extremely low.³

Acquired infection with *Toxoplasma gondii* in healthy individuals is commonly asymptomatic, however 10-20% of patients with acute infection may develop lymphadenopathy.⁴

Severe infections can occur in AIDS patients and adults immunocompromised by cancer chemotherapy or transplant recipients receiving immunosuppressive treatment. These infections can be fatal. Toxoplasmic encephalitis is the most common presentation and is the most frequent cause of focal central nervous system lesions in AIDS patients.⁵

Primary infection during pregnancy can result in transplacental transmission of the parasite resulting in congenital infection. The risk of congenital infection is lowest (10-25%) if acute maternal infection occurs during the first trimester and highest (60-90%) if it occurs during the third trimester.²

Severity of congenital infection is greatest when maternal infection is acquired early during pregnancy. Common outcomes of congenital toxoplasmosis include chorioretinitis, intracranial calcifications, and hydrocephalus. The majority of infants infected later in pregnancy are asymptomatic at birth, with sequelae occurring later in life.

Early treatment after prenatal diagnosis of *Toxoplasma gondii* infection has been shown to reduce the frequency and severity of congenital toxoplasmosis.⁶ Serological tests can be used to identify seronegative women, who then should be monitored during pregnancy.

The presence of IgG antibodies to *Toxoplasma gondii* indicates that infection has occurred but does not distinguish between recent and past infection. IgM antibodies are detected in individuals with a recently acquired infection, but antibodies may persist for up to 18 months post-infection.²

To differentiate between a recently acquired and a past infection, IgM and IgG positive specimens should be tested for IgG avidity. A high avidity index for IgG antibodies is a strong indication that an infection took place more than 4 months ago. Low avidity results cannot be used to diagnose an acute toxoplasmosis.

Toxo IgG	Toxo IgM	Toxo IgG Avidity	May indicate.../ Testing recommendation
nonreactive	nonreactive	N/A	no infection
nonreactive	reactive	N/A	obtain new sample 2-3 weeks after initial sample and test for Toxo IgG and Toxo IgM
reactive	nonreactive	high avidity	past infection. Strong indication that an infection took place more than 4 months ago
reactive	reactive	low avidity	obtain new sample 3 weeks after initial sample and test for Toxo IgG and Toxo IgM
reactive	reactive	high avidity	past infection. Strong indication that an infection took place more than 4 months ago

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The ARCHITECT Toxo IgM assay is a two-step immunoassay for the qualitative detection of IgM antibodies to *Toxoplasma gondii* in human serum and plasma using CMIA technology with flexible assay protocols, referred to as Chemiflex.

- In the first step, pre-diluted sample and anti-human IgM mouse monoclonal antibody coated paramagnetic microparticles are combined. Together with IgM antibodies of other specificities, anti-Toxo specific IgM present in the sample is bound by the anti-human IgM mouse monoclonal antibody coated microparticles, forming an antibody-antibody complex.
- After washing, a conjugate complex consisting of an acridinium-labeled anti-Toxo p30 antigen mouse monoclonal F(ab')₂ fragment and native *Toxoplasma gondii* lysate, containing the p30 antigen, is added to create a reaction mixture in the second step. This conjugate complex is bound by anti-Toxo specific IgM that has been captured by the anti-human IgM mouse monoclonal antibody coated microparticles in the first step, forming an antibody-antibody-conjugate complex.
- Following another wash cycle, Pre-Trigger and Trigger Solutions are added to the reaction mixture.
- The resulting chemiluminescent reaction is measured as relative light units (RLUs). There is a direct relationship between the amount of anti-Toxo IgM in the sample and the RLUs detected by the ARCHITECT iSystem optics.

The presence or absence of anti-Toxo IgM in the specimen is determined by comparing the chemiluminescent signal in the reaction to the cutoff signal determined from an active calibration. If the chemiluminescent signal in the reaction is greater than or equal to the cutoff signal, the specimen is considered reactive for anti-Toxo IgM.

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REF	6C20-25	6C20-35
	100	500
MICROPARTICLES	1 x 6.6 mL	1 x 27.0 mL
CONJUGATE	1 x 5.9 mL	1 x 26.3 mL
MICROPARTICLES Anti-human IgM (mouse, monoclonal) antibody coated microparticles in TRIS buffer with protein stabilizers and detergent. Minimum concentration: 0.08% solids. Preservatives: antimicrobial agents.		
CONJUGATE Conjugate complex consisting of acridinium-labeled anti-Toxoplasma p30 antigen (mouse, monoclonal) antibody and native <i>Toxoplasma gondii</i> lysate in phosphate buffer with protein stabilizers and detergent. Minimum concentration: 25 µg/mL. Preservative: sodium azide.		

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-  For *In Vitro* Diagnostic Use

Safety Precautions

CAUTION: This product contains human-sourced and/or potentially infectious components. Refer to the **REAGENTS** section of this package insert. No known test method can offer complete assurance that products derived from human sources or inactivated microorganisms will not transmit infection. Therefore, all human-sourced materials should be considered potentially infectious. It is recommended that these reagents and human specimens be handled in accordance with the OSHA Standard on Bloodborne Pathogens. Biosafety Level 2 or other appropriate biosafety practices should be used for materials that contain or are suspected of containing infectious agents.⁷⁻¹⁰

The following warnings and precautions apply to: CONJUGATE	
Contains sodium azide.	
EUH032	Contact with acids liberates very toxic gas.
P501	Dispose of contents / container in accordance with local regulations.

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Reagent Handling

- Do not use reagent kits beyond the expiration date.
- Do not pool reagents within a kit or between kits.**
- Before loading the reagent kit on the system for the first time, the microparticle bottle requires mixing to resuspend microparticles that may have settled during shipment. For microparticle mixing instructions, refer to the **PROCEDURE, Assay Procedure** section of this package insert.
- Septums MUST be used to prevent reagent evaporation and contamination and to ensure reagent integrity. Reliability of assay results cannot be guaranteed if septums are not used according to the instructions in this package insert.**
 - To avoid contamination, wear clean gloves when placing a septum on an uncapped reagent bottle.
 - Once a septum has been placed on an open reagent bottle, **do not invert the bottle** as this will result in reagent leakage and may compromise assay results.
 - Over time, residual liquids may dry on the septum surface. These are typically dried salts and have no effect on assay efficacy.

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INSTRUMENT PROCEDURE

The ARCHITECT Toxo IgM assay file must be installed on the ARCHITECT iSystem from an ARCHITECT iSystem Assay CD-ROM prior to performing the assay.

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- Performance has not been established for the use of cadaveric specimens or the use of body fluids other than human serum and plasma.
- Liquid anticoagulants may have a dilution effect resulting in lower concentrations for individual patient specimens.
- Sodium citrate tubes cannot be used with the ARCHITECT Toxo IgM assay.
- The instrument does not provide the capability to verify specimen type. It is the responsibility of the operator to verify that the correct specimen types are used in the assay.

Specimen Conditions

- Do not use specimens with the following conditions:
 - heat-inactivated
 - pooled
 - grossly hemolyzed (> 500 mg/dL hemoglobin)
 - obvious microbial contamination
- For accurate results, serum and plasma specimens should be free of fibrin, red blood cells, and other particulate matter. Serum specimens from patients receiving anticoagulant or thrombolytic therapy may contain fibrin due to incomplete clot formation.
- To prevent cross contamination, use of disposable pipettes or pipette tips is recommended.
- All samples (calibrator, controls, and patient specimens) should be tested within 3 hours of being placed on board the ARCHITECT iSystem.

Preparation for Analysis

- Follow the tube manufacturer's processing instructions for collection tubes. Gravity separation is not sufficient for specimen preparation.

- Mix thawed specimens thoroughly by low speed vortexing or by inverting 10 times. Visually inspect the specimens. If layering or stratification is observed, continue mixing until specimens are visibly homogeneous.
- To ensure consistency in results, specimens must be transferred to a centrifuge tube and centrifuged at $\geq 10,000$ RCF (Relative Centrifugal Force) for 10 minutes before testing if
 - they contain fibrin, red blood cells, or other particulate matter,
 - they require repeat testing, or
 - they were frozen and thawed.
- Transfer clarified specimen to a sample cup or secondary tube for testing. For centrifuged specimens with a lipid layer, transfer only the clarified specimen and not the lipemic material.
- Inspect all specimens for bubbles. Remove bubbles with an applicator stick before analysis. Use a new applicator stick for each specimen to prevent cross contamination.

Specimen Storage

Specimen Type	Storage Temperature	Maximum Storage Time
Serum/plasma	15-30°C	≤ 3 days
	2-8°C	≤ 14 days
	-20°C or colder	--

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Remove serum or plasma from the clot, red blood cells, or separator gel if stored longer than the maximum 15-30°C or 2-8°C storage time and store frozen at -20°C or colder.

No qualitative performance differences were observed between experimental controls and nonreactive or reactive specimens subjected to 6 freeze/thaw cycles; however, multiple freeze/thaw cycles should be avoided.

Specimen Shipping

- Package and label specimens in compliance with applicable state, federal, and international regulations covering the transport of clinical specimens and infectious substances.
- It is recommended that specimens be removed from the clot, red blood cells, or separator gel.
- Ship on wet ice or on dry ice.
- Do not exceed the storage time limitations listed above.

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- ARCHITECT Toxo IgM Assay file obtained from the ARCHITECT iSystem e-Assay CD-ROM found on www.abbott diagnostics.com.
- 6C20-01 ARCHITECT Toxo IgM Calibrator
- 6C20-10 ARCHITECT Toxo IgM Controls
- ARCHITECT Pre-Trigger Solution
- ARCHITECT Trigger Solution
- ARCHITECT Wash Buffer
- ARCHITECT Reaction Vessels
- ARCHITECT Sample Cups
- ARCHITECT Septum
- ARCHITECT Replacement Caps
- Pipettes or pipette tips (optional) to deliver the volumes specified on the patient or control order screen.

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Assay Procedure
Before loading the reagent
microparticle bottle
that may have
microorganisms

How speed testing or by
the specimens. If layering of
specimens until specimens are
or (Relative

Assay Procedure

- Before loading the reagent kit on the system for the first time, the microparticle bottle requires mixing to resuspend microparticles that may have settled during shipment. After the first time the microparticles have been loaded, no further mixing is required.
 - Invert the microparticle bottle 30 times.**
- Visually inspect the bottle to ensure microparticles are resuspended. If microparticles are still adhered to the bottle, continue to invert the bottle until the microparticles have been completely resuspended.
- If the microparticles do not resuspend, DO NOT USE. Contact your local Abbott representative.**
- Once the microparticles have been resuspended, discard the cap and place a septum on the bottle. For instructions on placing septums on bottles, refer to the **Reagent Handling** section of this package insert.
- Load the reagent kit on the ARCHITECT iSystem.
 - Verify that all necessary reagents are present.
 - Ensure that septums are present on all reagent bottles.
- Order calibration, if necessary.
 - For information on ordering calibrations, refer to the ARCHITECT System Operations Manual, Section 6.
- Order tests.
 - For information on ordering patient specimens and controls and for general operating procedures, refer to the ARCHITECT System Operations Manual, Section 5.
- Minimum sample cup volume is calculated by the system and printed on the Orderlist report. To minimize the effects of evaporation, verify adequate sample cup volume is present prior to running the test.

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 - Priority:

Sample volume for first test: 70 µL

Sample volume for each additional test from same sample cup: 20 µL
 - ≤ 3 hours on board:

Sample volume for first test: 150 µL

Sample volume for each additional test from same sample cup: 20 µL
 - > 3 hours on board: Additional sample volume required. For information on sample evaporation and volumes, refer to the ARCHITECT System Operations Manual, Section 5.
 - If using primary or aliquot tubes, use the sample gauge to ensure sufficient patient specimen is present.
- Prepare ARCHITECT Toxo IgM Calibrator and Controls.
 - Mix calibrator(s) and controls by gentle inversion before use.
 - Hold bottles **vertically** and dispense recommended volumes into each respective sample cup.
 - Recommended volumes:

for Calibrator 1: 6 drops

for each control: 4 drops
- Load samples.
 - For information on loading samples, refer to the ARCHITECT System Operations Manual, Section 5.
- Press RUN.
- For additional information on principles of operation, refer to the ARCHITECT System Operations Manual, Section 3.
- For optimal performance, it is important to perform routine maintenance as described in the ARCHITECT System Operations Manual, Section 9. Perform maintenance more frequently when required by laboratory procedures.

Specimen Dilution Procedures

Specimens cannot be diluted for the ARCHITECT Toxo IgM assay.

Calibration

- Test Calibrator 1 in replicates of three. Calibrator 1 should be priority loaded.

A single sample of each control level must be tested to evaluate the assay calibration. Ensure that assay control values are within the ranges specified in the respective control package insert.
- Once an ARCHITECT Toxo IgM calibration is accepted and stored, all subsequent samples may be tested without further calibration unless:
 - A reagent kit with a new lot number is used or
 - Controls are out of range.
- It is recommended that the assay be calibrated every 30 days.
- For detailed information on how to perform an assay calibration, refer to the ARCHITECT System Operations Manual, Section 6.

Quality Control Procedures

The recommended control requirement for the ARCHITECT Toxo IgM assay is that a single sample of each control level be tested once every 24 hours each day of use. If the quality control procedures in your laboratory require more frequent use of controls to verify test results, follow your laboratory-specific procedures.

The ARCHITECT Toxo IgM Control values must be within the acceptable ranges specified in the control package insert. If a control is out of its specified range, the associated test results are invalid and samples must be retested. Recalibration may be indicated.

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The ARCHITECT Toxo IgM assay belongs to method group 5 (except functional sensitivity).

RESULTS

Calculation

The ARCHITECT iSystem calculates the Calibrator 1 mean chemiluminescent signal from three Calibrator 1 replicates and stores the result. Results are reported by dividing sample result by the stored Calibrator 1 result. The default result unit for the ARCHITECT Toxo IgM assay is Index. Sample results may also be reported as sample to cutoff (S/CO). Index value divided by 0.60 equals S/CO value.

Interpretation of Results

Initial Results

Specimens with results	Interpretation	Retest Procedure
< 0.50 Index (< 0.83 S/CO)	Nonreactive for IgM antibodies to <i>Toxoplasma gondii</i>	--
0.50 ≤ x < 0.60 Index (0.83 ≤ x < 1.00 S/CO)	Grayzone	It is recommended to take a second sample within a reasonable period of time (e.g. two weeks) and repeat ARCHITECT Toxo IgM testing.
≥ 0.60 Index (≥ 1.00 S/CO)	Reactive for IgM antibodies to <i>Toxoplasma gondii</i>	--

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LIMITATIONS OF THE PROCEDURE

- If the ARCHITECT Toxo IgM results are inconsistent with clinical evidence, additional testing is suggested to confirm the result.
- For diagnostic purposes, results should be used in conjunction with other data; e.g., results of other tests (Toxo IgG, Toxo IgG Avidity), clinical impressions, etc.
- Heterophilic antibodies in human serum can react with reagent immunoglobulins, interfering with *in vitro* immunoassays. Patients routinely exposed to animals or to animal serum products can be prone to this interference, and anomalous values may be observed. Additional information may be required for diagnosis.¹¹
- Specimens from patients who have received preparations of mouse monoclonal antibodies for diagnosis or therapy may contain human anti-mouse antibodies (HAMA). Such specimens may show either falsely elevated or depressed values when tested with assay kits such as ARCHITECT Toxo IgM that employ mouse monoclonal antibodies. Additional information may be required for diagnosis.^{12, 13}
- Specimens from patients with high levels of IgM, e.g. from patients with multiple myeloma, may show depressed values when tested with μ -capture format assays.
- Human plasma collected in sodium citrate tubes cannot be used with this assay, since the ARCHITECT Toxo IgM results may be affected by this tube type.

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Precision

The ARCHITECT Toxo IgM assay is designed to have a total** precision of $\leq 10\%$ CV for a positive specimen in the range of 0.60 – 2.40 Index. The study was performed at one internal and one external evaluation site each using one instrument. Precision was assessed on a panel consisting of three different control lots and one human plasma specimen.

Panel members were tested in replicates of four across three reagent lots and one calibrator lot at the external site and three calibrator lots at the internal site, on one instrument at each site. Each combination of panel members and reagent lots was tested in four runs across several days. Data from this study are summarized in the following table*.

Member	N	Mean Index	Within Run		Total**	
			SD	%CV	SD	%CV
Negative Control	288	0.03	0.01	16.79	0.01	17.13
Positive Control	288	1.52	0.04	2.80	0.05	3.02
Human Plasma Specimen	96	1.42	0.05	3.26	0.05	3.72

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Seroconversion Sensitivity

The assay is designed to show a seroconversion sensitivity comparable to a commercially available diagnostic kit. A total of 122 bleeds from 39 different individuals seroconverting during acute toxoplasmosis infection were tested. Data from seven representative seroconverting individuals are shown in the following table*.

Sample ID	Months after last negative bleed	ARCHITECT Toxo IgM [Index]	Commercially available diagnostic kit [Index]	Isaga (Toxo-M) [Index]	Sabin-Feldman Dye Test [IU/mL]	HS Agglutination Test [IU/mL]
Grayzone		0.50 - 0.59	0.500 - 0.599	6-8	N/A	1
Reactive cutoff		0.60	0.600	9	2	2
30944001	0.0	0.05	0.061	0	< 2	< 1
30944002	1.2	0.56	0.351	4	5	1
30944003	2.1	0.89	0.588	10	800	64
30944004	2.2	0.77	0.542	10	800	64
30944005	4.3	0.38	0.193	1	1600	200
30944016	0.0	0.04	0.099	0	< 2	< 1
30944017	1.4	1.16	1.038	11	20	2
30944018	1.6	1.04	0.897	10	20	4
30944019	4.1	0.52	0.380	9	10	8
30944033	0.0	0.04	0.218	0	< 2	< 1
30944034	2.6	1.99	2.026	12	400	64
30944035	7.5	0.05	0.057	0	800	100
30944073	0.0	0.04	0.078	0	< 2	< 1
30944074	0.9	1.32	1.125	12	5	1
30944075	1.4	2.35	1.733	12	200	16
30944076	3.8	1.12	0.877	12	100	8
30944086	0.0	0.30	0.437	7	2	< 1
30944087	0.5	10.39	7.974	12	40	8
30944088	1.3	9.23	6.464	12	400	128
30944089	2.3	8.53	5.398	12	400	50
30944090	0.0	0.05	0.081	0	< 2	< 1
30944091	1.2	5.95	4.195	12	20	2
30944092	1.5	5.72	3.679	12	200	16
30944093	4.7	2.66	1.700	12	400	50
30944118	0.0	0.05	0.113	0	< 2	< 1
30944119	1.0	5.76	3.784	12	2	1
30944120	1.8	6.56	3.536	12	200	64
30944121	2.5	3.88	1.910	12	1600	100

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Resolved Relative Specificity

The ARCHITECT Toxo IgM assay is designed to have a resolved relative specificity comparable to a commercially available diagnostic kit. A study was performed at one internal and one external evaluation site. From the 2772*** specimens evaluated to assess resolved relative specificity 36 specimens were concordant reactive and additional three were confirmed positive after discordant resolution and therefore excluded from the specificity calculation.

*** NOTE: Specimens that could not be resolved or showed grayzone result interpretation on any assay being compared or used for discordant resolution were not included in the evaluation of resolved relative specificity.

Data from this study are summarized in the following table*.

Resolved Relative Specificity

Sample Type	ARCHITECT Toxo IgM		Commercially available diagnostic kit	
	Observed	Lower 95% Confidence Limit	Observed	Lower 95% Confidence Limit
Pregnant Women	99.95% (1987/1988)	99.72%	99.95% (1987/1988)	99.72%
Diagnostic / Hospital Patients	100% (451/451)	99.19%	100% (451/451)	99.19%
Blood Donors (Serum)	100% (154/154)	97.63%	100% (154/154)	97.63%
Blood Donors (Plasma)	98.57% (138/140)	94.93%	100% (140/140)	97.40%
Total	99.89% (2730/2733)	99.68%	99.96% (2732/2733)	99.80%

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Interference

No interference was observed between experimental controls and nonreactive or reactive specimens tested with elevated levels of bilirubin (20 mg/dL), triglycerides (3000 mg/dL), protein (12 g/dL), red blood cells (0.4% v/v), or hemoglobin (500 mg/dL).

Other Potential Interferants

Additional studies were performed to evaluate other potential interfering disease states on the ARCHITECT Toxo IgM assay. Sample categories were tested both unspiked and spiked with anti-Toxo IgM positive plasma.

A total of 167 unspiked specimens and a total of 165 spiked anti-Toxo IgM reactive specimens were tested from the following categories: Anti-nuclear antibodies (ANA), anti-dsDNA antibodies, Rheumatoid Factor, Herpes Simplex Virus 1 (anti-HSV-1 positive), Herpes Simplex Virus 2 (anti-HSV-2 positive), Epstein-Barr Virus (anti-EBV positive), Measles, Parvovirus B19 (anti-B19 virus positive), Varicella Zoster Virus (anti-VZV positive), Rubella Virus (anti-Rubella positive), Cytomegalovirus (anti-CMV positive), Hyperpolyclonal IgG, Hyperpolyclonal IgM, Monoclonal IgG, Monoclonal IgM, Human anti-mouse antibodies (HAMA), Influenza vaccine recipients, and Syphilis.

ARCHITECT Toxo IgM showed expected qualitative results in all categories with the exception of the categories summarized in the following table*.

Category	ARCHITECT Toxo IgM			
	N tested	Nonreactive	Reactive	Grayzone
Anti-nuclear Antibodies (ANA) Unspiked	10	8	1	1
Hyperpolyclonal IgM Spiked Reactive	10	1	8	1
Monoclonal IgM Spiked Reactive	5	4	0	1

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After discordant resolution the unspiked ANA specimen tested reactive on ARCHITECT Toxo IgM could not be resolved. Specimens from patients with high levels of IgM, e.g. from patients with multiple myeloma, may show depressed values when tested with μ -capture format assays. Refer to the **LIMITATIONS OF THE PROCEDURE** section of this package insert.

BIBLIOGRAPHY

- Willis MS, Southern P, Latimer MJ. Toxoplasma infection. Making the best use of laboratory tests. *Infect Med* 2002; 19: 522-532.
- Jones JL et al. Congenital Toxoplasmosis: A Review. *CME Review Article* Vol. 56, Number 5 2001; 296-305.
- Shulman IA. Parasitic infections and their impact on blood donor selection and testing. *Arch Pathol Lab Med.* 1994; 118: 366-370.
- Montoya J.G. Laboratory Diagnosis of Toxoplasma gondii Infection and Toxoplasmosis. *Journal of Infectious Diseases* 2002; 185 (Suppl 1) 73-82.
- Israelski DM, Remington JS. Toxoplasmosis in the non-AIDS immunocompromised host. *Curr Clin Top Infect Dis* 1993; 13: 322-356.
- Wong SY and Remington. Toxoplasmosis in Pregnancy. *Clinical Infectious Diseases* 1994; 18:853-862.
- US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
- US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories.* 5th ed. Washington, DC: US Government Printing Office; December 2009.
- World Health Organization. *Laboratory Biosafety Manual.* 3rd ed. Geneva: World Health Organization; 2004.
- Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline—Third Edition.* CLSI Document M29-A3. Wayne, PA: CLSI; 2005.
- Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. *Clin Chem* 1988;34(1):27-33.
- Primus FJ, Kelley EA, Hansen HJ, et al. "Sandwich"-type immunoassay of carcinoembryonic antigen in patients receiving murine monoclonal antibodies for diagnosis and therapy. *Clin Chem* 1988;34(2):261-264.
- Schroff RW, Foon KA, Beatty SM, et al. Human anti-murine immunoglobulin responses in patients receiving monoclonal antibody therapy. *Cancer Res* 1985;45(2):879-885.

ADDENDUM
to the Instruction for Use ARCHITECT
Toxo IgM Reagent for Russia
SHIPPING
Recommendations

■ Key to Symbols

	Caution
	Consult instructions for use
	Manufacturer
	Sufficient for
	Temperature limitation
	Use by/Expiration date
CONJUGATE	Conjugate
CONTAINS: AZIDE	Contains Sodium Azide. Contact with acids liberates very toxic gas.
CONTROL NO.	Control Number
IVD	<i>In Vitro</i> Diagnostic Medical Device
LOT	Lot Number
MICROPARTICLES	Microparticles
PRE-TRIGGER SOLUTION	Pre-Trigger Solution
PRODUCT OF GERMANY	Product of Germany
REACTION VESSELS	Reaction Vessels
REAGENT LOT	Reagent Lot
REF	List Number
REPLACEMENT CAPS	Replacement Caps
SAMPLE CUPS	Sample Cups
SEPTUM	Septum
SN	Serial number
TRIGGER SOLUTION	Trigger Solution
WASH BUFFER	Wash Buffer

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©2008, 2014 Abbott Laboratories

ADDENDUM

to the Instruction for Use ARCHITECT

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- Willis MS, Southern P, Latimer MJ. Toxoplasma infection. Making the best use of laboratory tests. *Infect Med* 2002; 19: 522-532.
- Jones JL et al. Congenital Toxoplasmosis: A Review. *CME Review Article* Vol. 56, Number 5 2001; 296-305.
- Shulman IA. Parasitic infections and their impact on blood donor selection and testing. *Arch Pathol Lab Med*. 1994; 118: 366-370.
- Montoya J.G. Laboratory Diagnosis of Toxoplasma gondii Infection and Toxoplasmosis. *Journal of Infectious Diseases* 2002; 185 (Suppl 1) 73-82.
- Israelis DM, Remington JS. Toxoplasmosis in the non-AIDS immunocompromised host. *Curr Clin Top Infect Dis* 1993; 13: 322-356.
- Wong SY and Remington. Toxoplasmosis in Pregnancy. *Clinical Infectious Diseases* 1994; 18:853-862.
- US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
- US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
- World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
- Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline—Third Edition*. CLSI Document M29-A3. Wayne, PA: CLSI; 2005.
- Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. *Clin Chem* 1988;34(1):27-33.

- Primus FJ, Kelley EA, Hansen HJ, et al. "Sandwich"-type immunoassay of carcinoembryonic antigen in patients receiving murine monoclonal antibodies for diagnosis and therapy. *Clin Chem* 1988;34(2):261-264.
- Schroff RW, Foon KA, Beatty SM, et al. Human anti-murine immunoglobulin responses in patients receiving monoclonal antibody therapy. *Cancer Res* 1985;45(2):879-885.

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INTENDED USE

The ARCHITECT Toxo IgM Calibrator is for the calibration of the ARCHITECT iSystem when used for the qualitative detection of IgM antibodies to *Toxoplasma gondii* in human serum and plasma. Refer to the reagent insert and the ARCHITECT System Operations Manual for additional information.

CONTENTS

1 Bottle (4.0 mL) of ARCHITECT Toxo IgM Calibrator contains anti-Toxoplasma p30 antigen (human, monoclonal) IgM antibody prepared in recalcified human plasma. The calibrator is reactive for anti-Toxo IgM. Preservatives: ProClin 950 and sodium azide.

STANDARDIZATION

The value of the ARCHITECT Toxo IgM Calibrator is determined against internal reference calibrators.

PRECAUTIONS

- **IVD**
- For In Vitro Diagnostic Use

- **CAUTION:** This product contains human-sourced and/or potentially infectious components. Refer to the **CONTENTS** section of this package insert. No known test method can offer complete assurance that products derived from human sources or inactivated microorganisms will not transmit infection. Therefore, all human-sourced materials should be considered potentially infectious. It is recommended that these reagents and human specimens be handled in accordance with the OSHA Standard on Bloodborne Pathogens. Biosafety Level 2 or other appropriate biosafety practices should be used for materials that contain or are suspected of containing infectious agents.¹⁻⁴
- The human plasma used in the Calibrator is nonreactive for HBsAg, HIV-1 RNA or HIV-1 Ag, anti-HIV-1/HIV-2 and anti-HCV.

The following warnings and precautions apply to: CAL 1	
WARNING	Contains methylisothiazolone and sodium azide.
H317	May cause an allergic skin reaction.
EUH032	Contact with acids liberates very toxic gas.
Prevention	
P261	Avoid breathing mist / vapors / spray.
P272	Contaminated work clothing should not be allowed out of the workplace.
P280	Wear protective gloves / protective clothing / eye protection.

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INTENDED USE

The ARCHITECT Toxo IgM Controls are for the estimation of test precision and the detection of systematic analytical deviations of the ARCHITECT iSystem when used for the qualitative detection of IgM antibodies to *Toxoplasma gondii* in human serum and plasma. Refer to the reagent insert and the ARCHITECT System Operations Manual for additional information.

CONTENTS

2 Bottles (4.0 mL each) of ARCHITECT Toxo IgM Controls: Negative Control and Positive Control prepared in recalcified human plasma. The Positive Control contains anti-Toxoplasma p30 antigen (human, monoclonal) IgM antibody and is reactive for anti-Toxo IgM. Preservatives: ProClin 950 and sodium azide.

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Control	Color	Target Anti-Toxo IgM Concentration (Index)	Control Range (Index)
CONTROL -	Natural	N/A	< 0.40
CONTROL +	Natural	1.50	0.75 - 2.25

Control	Color	Target Anti-Toxo IgM Concentration (S/CO)	Control Range (S/CO)
CONTROL -	Natural	N/A	< 0.67
CONTROL +	Natural	2.50	1.25 - 3.75

PRECAUTIONS

- **IVD**
- For *In Vitro* Diagnostic Use

-  **CAUTION:** This product contains human-sourced and/or potentially infectious components. Refer to the **CONTENTS** section of this package insert. No known test method can offer complete assurance that products derived from human sources or inactivated microorganisms will not transmit infection. Therefore, all human-sourced materials should be considered potentially infectious. It is recommended that these reagents and human specimens be handled in accordance with the OSHA Standard on Bloodborne Pathogens. Biosafety Level 2 or other appropriate biosafety practices should be used for materials that contain or are suspected of containing infectious agents.¹⁻⁴
- The human plasma used in the Negative and Positive Control is nonreactive for HBsAg, HIV-1 RNA or HIV-1 Ag, anti-HIV-1/HIV-2 and anti-HCV.

The following warnings and precautions apply to: CONTROL - / CONTROL +	
	
WARNING	Contain methylisothiazolone and sodium azide.
H317	May cause an allergic skin reaction.
EU032	Contact with acids liberates very toxic gas.
Prevention	
P261	Avoid breathing mist / vapors / spray.
P272	Contaminated work clothing should not be allowed out of the workplace.
P280	Wear protective gloves / protective clothing / eye protection.
Response	
P302+P352	IF ON SKIN: Wash with plenty of water.
P333+P313	If skin irritation or rash occurs: Get medical advice / attention.
P362+P364	Take off contaminated clothing and wash it before reuse.
Disposal	
P501	Dispose of contents / container in accordance with local regulations.

Safety Data Sheets are available at www.abbottddiagnostics.com or contact your local representative.

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PREPARATION FOR USE

Controls may be used immediately after removal from 2-8°C storage. Prior to each use, mix by gentle inversion.

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7816 / R03
6C20
C200
Toxo IgM
en

Key to Symbols

	Caution
	Consult instructions for use
	Manufacturer
	Temperature limitation
	Use by/Expiration date
CONC	Concentration
CONTAINS: AZIDE	Contains Sodium Azide. Contact with acids liberates very toxic gas.
CONTROL -	Negative Control
CONTROL +	Positive Control
INDEX	Index
IVD	<i>In Vitro</i> Diagnostic Medical Device
LOT	Lot Number
PRODUCT OF GERMANY	Product of Germany
RANGE	Range
REF	List Number
WARNING: SENSITIZER	Warning: May cause an allergic reaction.

ARCHITECT is a trademark of Abbott Laboratories in various jurisdictions.
ProClin is property of its respective owner.

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Revised October 2014.
©2008, 2014 Abbott Laboratories

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