Abbott

Abbott GmbH & Co. KG
Max-Planck-Ring 2
65205 Wiesbaden
Postfach 1303, 65011 Wiesbaden

Tel. +49 61 22 58 0
Fax +49 61 22 58 1244

I, Zaman Khan, Associate Director Regulatory Affairs EMEA and Pakistan, Abbott GmbH & Co. KG., hereby certify that the attached documents that keep the number of pages indicated below are the operational documentation. These documents are in a full accordance with actual manuals of manufacturer.

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Alinity i
HIV Ag/Ab Combo
dated August 2016.

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Signed: Zo
Zaman Khan,
Associate Director Regulatory
Affairs EMEA and Pakistan,
Abbott GmbH & Co. KG.

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
 **en**
HIV Ag/Ab
REF 08P0710
G71154R02
C8P070

NAME

Alinity i HIV Ag/Ab Combo Controls (also referred to as HIV Ag/Ab Ctrls)

INTENDED USE

The Alinity i HIV Ag/Ab Combo Controls are for the estimation of test precision and the detection of systematic analytical deviations of the Alinity i analyzer when used for the simultaneous qualitative detection of HIV p24 antigen and antibodies to human immunodeficiency virus type 1 and/or type 2 (HIV-1/HIV-2) in human serum or plasma.

For additional information, refer to the Alinity i HIV Ag/Ab Combo reagent package insert and the Alinity ci-series Operations Manual.

CONTENTS

CONTROL -, **CONTROL +1**, and **CONTROL +2** are prepared in recalcified human plasma. The Positive Control 1 (inactivated) is reactive for anti-HIV-1. The Positive Control 2 (inactivated) is reactive for anti-HIV-2. The **CONTROL +3** is purified HIV viral lysate prepared in TRIS buffered saline with protein (bovine) stabilizer.

Preservatives for Negative Control, Positive Control 1, and Positive Control 2: sodium azide and antimicrobial agent. Preservative for Positive Control 3: sodium azide.

The controls are at the following S/CO ranges:

Control	Quantity	Color	RANGE (S/CO)
CONTROL -	1 x 8.0 mL	Natural	0.00 - 0.50
CONTROL +1	1 x 8.0 mL	Blue ^a	1.20 - 11.50
CONTROL +2	1 x 8.0 mL	Yellow ^b	1.52 - 8.30
CONTROL +3	1 x 8.0 mL	Purple ^c	1.87 - 4.59

^a Dye: Acid Blue No. 9

^b Dye: Acid Yellow No. 23

^c Dye: Acid Blue No. 9 and Red D&C No. 33

NOTE: The insert ranges for the controls are not lot specific and represent the total range of values which may be generated throughout the life of the product. It is recommended that each laboratory establish its own means and acceptable ranges which should fall within the package insert ranges. Sources of variation that can be expected include:

- Calibration
- Control lot
- Reagent lot
- Calibrator lot
- Instrument

PRECAUTIONS

- **IVD**
- For In Vitro Diagnostic Use

Safety Precautions

-  **CAUTION:** This product contains human-sourced and/or potentially infectious components. Refer to the CONTENTS section of this package insert. No known test method can offer complete assurance that products derived from human sources or inactivated microorganisms will not transmit infection. Therefore, all human-sourced materials should be considered potentially infectious. It is recommended that these reagents and human specimens be handled in accordance with the OSHA Standard on Bloodborne Pathogens. Biosafety Level 2 or other appropriate biosafety practices should be used for materials that contain or are suspected of containing infectious agents.¹⁻⁴
- The human-sourced material used in the Negative Control is nonreactive for HBsAg, anti-HIV-1/HIV-2, anti-HCV and HIV-1 RNA or HIV-1 Ag.
- The human-sourced material used in the Positive Control 1 is reactive for anti-HIV-1. The human-sourced material is also tested for HIV-1 by either HIV-1 Ag and is nonreactive, or by HIV-1 RNA and may be reactive. The human-sourced material is nonreactive for anti-HIV-2, HBsAg and anti-HCV.
- The human-sourced material used in the Positive Control 2 is reactive for anti-HIV-2. The human-sourced material is also tested for HIV-1 by either HIV-1 Ag or by HIV-1 RNA and is nonreactive. The human-sourced material is nonreactive for anti-HIV-1, HBsAg, and anti-HCV.
- The Positive Control 3 is purified HIV-1 viral lysate.

The following warnings and precautions apply to:

CONTROL - , CONTROL +1 , CONTROL +2 , and CONTROL +3	
Contains sodium azide.	
EUH032	Contact with acids liberates very toxic gas.
P501	Dispose of contents / container in accordance with local regulations.

Safety Data Sheets are available at www.abbottdiagnostics.com or contact your local representative.

For a detailed discussion of safety precautions during system operation, refer to the Alinity ci-series Operations Manual, Section 8.

PREPARATION FOR USE

- This product is liquid ready-to-use.
- This product may be used immediately after removal from 2 to 8°C storage.
- Prior to each use, mix by gentle inversion (5 to 10 times).

STORAGE

- Do not use past expiration date.

	Storage Temperature	Maximum Storage Time	Additional Storage Instructions
Unopened	2 to 8°C	Until expiration date	
Opened	2 to 8°C	Until expiration date	Store tightly capped. Return to refrigerated storage after use.

INSTRUMENT PROCEDURE

- To obtain the recommended volume requirements for the controls, hold the bottle vertically, and dispense 10 drops of the Negative Control, 10 drops of the Positive Control 1, 10 drops of the Positive Control 2, and 10 drops of the Positive Control 3 into each sample cup in the assigned position.
- For instructions on ordering and loading controls on the instrument, refer to the Alinity ci-series Operations Manual, Section 5.

INDICATIONS OF INSTABILITY OR DETERIORATION

Instability or deterioration should be suspected if there are precipitates, visible signs of leakage, turbidity, or if controls do not meet the appropriate package insert and/or Alinity ci-series Operations Manual criteria.

BIBLIOGRAPHY

- US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
- US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
- World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
- Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.

Note for number formatting:

- A space is used as thousands separator (example: 10 000 specimens).
- A period is used to separate the integer part from the fractional part of a number written in decimal form (example: 3.12%).

Key to Symbols

	Caution
	Consult instructions for use
	Manufacturer
	Temperature limitation
	Use by/Expiration date
CN	Control Number
CONTAINS: AZIDE	Contains Sodium Azide. Contact with acids liberates very toxic gas.
CONTROL -	Negative Control
CONTROL + 1	Positive Control 1
CONTROL + 2	Positive Control 2
CONTROL + 3	Positive Control 3
IVD	In Vitro Diagnostic Medical Device
LOT	Lot Number
PRODUCT OF GERMANY	Product of Germany
RANGE	Range
REF	List Number

Alinity is a trademark of Abbott Laboratories in various jurisdictions. All other trademarks are property of their respective owners.

Abbott GmbH & Co. KG
Max-Planck-Ring 2
65205 Wiesbaden
Germany
+49-6122-580

Customer Service: Contact your local representative or find country-specific contact information on www.abbottdiagnostics.com

Created August 2016.

©2016 Abbott Laboratories

Instruction for Use Alinity
STORAGE AND STORAGE

recommended shipping
refrigerators/refrigerator

medical devices must
thermo containers with

Device is stable until
temperature limits
rules, applicable to

WARRANTY

Manufacturer's
Manufacturer's

Alinity i HIV

Device sh

MAINTENANCE

Not re

REP

No

D

ENDUM

Instruction for Use Alinity i HIV Ag/Ab Combo Controls for Russia

SHIPPING AND STORAGE CONDITIONS

Recommended shipping conditions for devices are 2-8°C. Medical devices must be stored in refrigerators/refrigerator chambers to ensure regulated temperature conditions with a daily temperature record.

Medical devices must be carried out in single-use thermal containers containing coolant elements, or in reusable thermo containers with automatically maintained temperature or in refrigerated trucks using thermal indicators.

Device is stable until the expiration date indicated on the package. Device is no longer suitable for use if shipping temperature limits have been violated. Ship in any kind of enclosed carrier vehicle in accordance with shipping rules, applicable to the respective type of transport.

WARRANTY

Manufacturer guarantees that device remains stable until the expiration date when stored and handled as intended. Manufacturer set the maximum shelf life of

Alinity i HIV Ag/Ab Combo Controls – 10 months from date of manufacturing.

Device should not be used if stored with deviations from manufacturer's requirements (GOST R 51088-2013).

MAINTENANCE

Not required.

REPAIRS

Not required.

DESTRUCTION AND DISPOSAL OF A MEDICAL DEVICE

Destruct and dispose of a medical device in accordance with the requirements of SanPin 2.1.7.2790-10 dated December 12, 2010 "Sanitary-epidemiological Requirements for the Management of Medical Waste". Transmit and dispose of a product and its components and packaging in accordance with local regulations.

LIST OF APPLICABLE STANDARDS

98/79/EC, 1272/2008/EC (CLP), 1907/2006/EC (REACH), 2002/364/EC, 2009/886/EC, Regulation (EC) 1069/2009, Regulation (EC) 142/2011, CLASS 2013, GHS, WHSR Part 7.1, EN ISO 9001, ISO 13485, EN ISO 14971, EN 13641, EN 13612, EN ISO 15223-1, EN 980 EN ISO 18113-1, EN ISO 18113-2, EN ISO 18113-3, EN ISO 23640*, EN ISO 20776-1, EN 13975, EN ISO 17511, EN ISO 15193, EN ISO 15194

POTENTIAL USER

The Alinity i HIV Ag/Ab Combo Controls is intended to be used in a variety of clinical and laboratory conditions, including, without limitation, hospitals, diagnostic centers, reference laboratories and blood transfusion facilities around the world, whose personnel have been trained in the use of the Alinity i analyzer.

INFORMATION ABOUT PRESENCE OF DRUGS IN COMPOSITION OF MEDICAL DEVICE

Does not contain any drugs

STERILIZATION

Medical device is not required to be sterile.

INFORMATION ABOUT POSSIBILITY OF MULTIPLE APPLICATION

Medical device is intended for one-off usage

MATERIALS OF ANIMAL ORIGIN

Materials of animal origin are tested and risks of using MD are indicated in SDS. Safety Data Sheets are available at www.abbottdiagnostics.com or contact your local representative.

INDICATIONS AND CONTRAINDICATIONS FOR USE

Not applicable.

TECHNICAL CHARACTERISTICS

- Negative Control 1: pH-value at 20 °C: 7.3
Density at 20 °C: 1.025 g/cm³
- Positive Control 1: pH-value at 20 °C: 7.3
Density at 20 °C: 1.025 g/cm³
- Positive Control 2: pH-value at 20 °C: 7.3
Density at 20 °C: 1.025 g/cm³
- Positive Control 3: pH-value at 20 °C: 7.2
Density at 20 °C: 1.007 g/cm³

COMPATIBILITY RESTRICTIONS OF MEDICAL DEVICE

No available data.

COLOR AND ODOR

Negative Control 1: Yellow and odorless
Positive Control 1: Blue and odorless
Positive Control 2: Yellow and odorless
Positive Control 2: Purple and odorless

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1. US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
2. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2008.
3. World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
4. Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.

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