
VALIDATED

Magnus Nilsson
CEO and Dr Med Sci

« _____ » 2014

Instruction For Use - Perfadex®

1

1

NOTARIUS PUBLICUS I KUNGSBACKA

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NOTARIUS PUBLICUS I KUNGSBACKA
NOTARY PUBLIC KUNGSBACKA, SWEDEN**CERTIFICATE**

I, the undersigned, Joakim Gillersten, Notary Public in the City of Kungsbacka Sweden, do hereby certify that Dr Magnus Nilsson has personally signed this enclosed document and that he is duly authorized to sign for the company, XVIVO Perfusion AB, in his normal business activities as CEO.

Furthermore, I certify that attached document is identical with original, which has been presented before me.

Kungsbacka, this the 4th day of September, 2014

Ex officio:

Joakim Gillersten

Notary Public in Kungsbacka, Sweden

Exp. No: 1954-2014

NOTARIUS PUBLICUS KUNGSBACKA


ADRESS
ENERGIGATAN 12
SE 434 37 KUNGSBACKA
SWEDEN

TELEFON
+46-300-710 18

E-MAIL
JOAKIM@GILLERSTEN.COM

VALIDATED

Magnus Nilsson
CEO and Dr Med Sci

« _____ » _____ 2014

Instruction For Use – STEEN Solution™

STEEN Solution™

xvivo
PERFUSION

EN: Description

STEEN Solution™ is a physiological salt solution containing human serum albumin and dextran.

The solution is sterile (aseptic processing) and intended for single use only.

The bottle is made of PETg and equipped with a PE screw cap lined with a silicone septum closure which facilitates aseptic transfer of the solution. The screw cap is sealed by a tamper evident plastic sleeve.

Intended purpose

STEEN Solution™ is intended for assessment of isolated lungs after removal from the donor in preparation for eventual transplantation into a recipient.

STEEN Solution™ is not intended for any other use than that stated in the section Intended purpose.

STEEN Solution™ is not intended for infusion.

Precautionary Measures

The responsibility for correct clinical use and technique rests with the user. The instructions are only provided as a suggestion for procedure, the user must on the basis of her or his medical training and experience evaluate the suitability of this procedure. Do not use STEEN Solution™ if the solution is unclear, the bottle is damaged, the flip-tear seal has been tampered with or if the 'use by' date has expired.

STEEN Solution™ is intended for single use only and MAY NOT BE REUSED. Any left over solution must be disposed of after the procedure.

Caution: All blood products should be treated as potentially infectious. Source material from which this product was derived was found negative when tested for antibodies to HIV, HBV, HCV, and HTLV I/II and non-reactive for HbsAg, HCV RNA and HIV-1 RNA and syphilis. No known test methods can offer assurance that products derived from human blood will not transmit infectious agents.

Adverse Effects

When administered systemically, human serum albumin has been associated with rare allergic reactions (<1 in 1000) such as urticaria, fever, chills, pruritus and anaphylaxis. Dextran has also been associated with rare anaphylactic reactions (<1 in 1000). However, no such reactions have been reported with either of these substances when used for ex vivo organ or tissue perfusion or preservation. The risk of rare reactions to dextran following systemic use can be virtually eliminated by prior injection of Promit; similar protection is afforded in patients who have received i.v. dextran during the preceding 48 hours.

Storage Conditions

STEEN Solution™ shall be stored dark at +2 °C to +8 °C with the bottle in an up-right position.

Lung evaluation set-up

The lungs or the heart-lung block are enclosed in a transparent plastic container to maintain suitable humidity. The 'venous' afferent side of the closed circuit is connected to a pump and a heat and gas exchanger so that the perfusate assumes physiological temperature and partial gas pressures. A leucocyte filter is connected before the inflow to prevent leucocyte-induced tissue injury.

Perfusion pressure is closely monitored and should never exceed 20 mm Hg. Since the lungs are generally cold before perfusion, initial flow rates must be very low, gradually increasing as the temperature increases. The flow should never exceed 4 l/min. Careful ventilation is begun as the temperature of the perfusate outflow reaches 32 °C.

Functional assessment can begin when the circuit has reached

the steady state at 37 °C, monitoring arterial and venous blood gases, endtidal carbon dioxide and a range of hemodynamic parameters.

Instructions

1. Read and understand all of the instructions below before proceeding.

2. Assemble a sterile extra corporal circulation system including pump, deoxygenator, fluid reservoir and facilities for measurements of gas pressure, pH, electrolyte balance and any other relevant parameter.

3. Remove the bottle from the carton and check that:

- the bottle is intact,
- the plastic sleeve around the screw cap has not been tampered with,
- and the 'use by' date has not expired.

Hang the bottle upside down in an IV stand and remove the plastic sleeve.

4. Disinfect the septum closure of the screw cap and transfer the solution aseptically via a sterile transfer set equipped with a suitable spike to the fluid reservoir in the extra corporal circulation system. Circulate the solution and de-air the system.

5. If the system is not primed properly, repeat with a new bottle of STEEN Solution™ from procedure 3.

6. Cellular or acellular STEEN Solution™ can be used in lung transplantation. If the cellular alternative is preferred aseptically add sufficient Packed Red Blood Cells (PRBC) to achieve a final hematocrit of about 15%.

7. Circulate the final solution until appropriate baseline values of pH, pCO₂, pO₂, electrolyte balance and any other relevant parameters have been established.

8. Once appropriate values have been established proceed with connecting the lungs to the extra corporal circulation system and progressively increase perfusate temperature to normothermia before commencing evaluation of lung function according to established clinical routines (see ref 1 and 2).

9. Dispose of the bottles and any left over solution after the procedure according to the hospital's standard routines. Note that the solution contains human serum albumin (and may contain human red cells).

STEEN Solution™ contains

- Human Serum Albumin – provides normal oncotic pressure preventing oedema formation.
- Dextran – a mild scavenger which coats and protects endothelium from subsequent excessive leucocyte interaction and thrombogenesis.
- Extra-cellular electrolyte composition (lowK⁺) – reduces free radical generation and avoids vascular spasm under normothermic conditions.

Specifications

Sterility	SAL 10 ⁻⁶
Bacterial endotoxin [EU/mL]	<0.25
Osmolality [mOsm/kg]	275 to 315
pH	7.1 to 7.7

Symbols

	Sterilized using aseptic processing techniques
	Containing dextran 40

	Containing human serum albumin
	Do not re-use, discard after procedure. Reuse of STEEN Solution is not allowed due to the risk for infection.
	Use by
	Temperature limitation
	Caution: Consult accompanying document
	Catalog number
	Batch code

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