Abbott GmbH & Co. KG
Max-Planck-Ring 2
65205 Wiesbaden
Postfach 1303, 65011 Wiesbaden

Tel. +49 61 22 58 0
Fax +49 61 22 58 1244

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
Signed: 
Zaman Khan,
Associate Director Regulatory
Affairs EMEA and Pakistan,
Abbott GmbH & Co. KG.

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G71284R02

C7P460

Read Highlighted Changes: Revised July 2019.

NAME

Alinity i Toxo IgG Avidity Controls (also referred to as Tox Avidity Ctrl)

INTENDED USE

The Alinity i Toxo IgG Avidity Controls are for the estimation of test precision and the detection of systematic analytical deviations of the Alinity i analyzer.

For additional information, refer to the Alinity i Toxo IgG Avidity reagent package insert and the Alinity ci-series Operations Manual.

CONTENTS

CONTROL LOW AVIDITY contains anti-*Toxoplasma gondii* monoclonal antibody in phosphate buffered saline with protein (bovine) stabilizers. Preservative: ProClin 300.

CONTROL HIGH AVIDITY contains recalcified human plasma.

Preservatives: ProClin 950 and sodium azide.

The controls are reactive for anti-Toxo IgG.

The controls are at the following ranges:

Control	Quantity	Color	Toxo IgG Avidity
			RANGE (%Avi)
CONTROL LOW AVIDITY	1 x 6.0 mL	Natural	< 50.0
CONTROL HIGH AVIDITY	1 x 6.0 mL	Natural	≥ 60.0

NOTE: The controls have been designed exclusively for use with this assay and assay format and therefore are not suitable for use with other assays. The controls cannot be run as patient specimens. If multiple control replicates are required, note that only one avidity determination can be ordered per sample cup.

The controls are manufactured and released within specific ranges and not to a specified target value. Therefore lot-to-lot differences are expected. Levey-Jennings may be utilized for data tracking, however application of Westgard Rules is not recommended due to the stated lot-to-lot variability.

STANDARDIZATION

The Alinity i Toxo IgG Avidity assay is standardized against internal reference material.

PRECAUTIONS

- IVD**
- For In Vitro Diagnostic Use

Safety Precautions

-  **CAUTION:** This product contains human-sourced and/or potentially infectious components. Refer to the CONTENTS section of this package insert. No known test method can offer complete assurance that products derived from human sources or inactivated microorganisms will not transmit infection. Therefore, all human-sourced materials should be considered potentially infectious. It is recommended that these reagents and human specimens be handled in accordance with the OSHA Standard on Bloodborne Pathogens, Biosafety Level 2 or other appropriate biosafety practices should be used for materials that contain or are suspected of containing infectious agents.¹⁻⁴

- The human-sourced material used in the Control High Avidity is reactive for anti-Toxo IgG, and nonreactive for HBsAg, HIV-1 RNA or HIV-1 Ag, anti-HIV-1/HIV-2, and anti-HCV.

The following warnings and precautions apply to:	
CONTROL LOW AVIDITY	
	
WARNING	Contains methylisothiazolones.
H317	May cause an allergic skin reaction.
Prevention	
P261	Avoid breathing mist / vapors / spray.
P272	Contaminated work clothing should not be allowed out of the workplace.
P280	Wear protective gloves / protective clothing / eye protection.
Response	
P302+P352	IF ON SKIN: Wash with plenty of water.
P333+P313	If skin irritation or rash occurs: Get medical advice / attention.
P362+P364	Take off contaminated clothing and wash it before reuse.
Disposal	
P501	Dispose of contents / container in accordance with local regulations.

The following warnings and precautions apply to:	
CONTROL HIGH AVIDITY	
	
WARNING	Contains methylisothiazolones and sodium azide.
H317	May cause an allergic skin reaction.
EUH032	Contact with acids liberates very toxic gas.
Prevention	
P261	Avoid breathing mist / vapors / spray.
P272	Contaminated work clothing should not be allowed out of the workplace.
P280	Wear protective gloves / protective clothing / eye protection.
Response	
P302+P352	IF ON SKIN: Wash with plenty of water.
P333+P313	If skin irritation or rash occurs: Get medical advice / attention.
P362+P364	Take off contaminated clothing and wash it before reuse.
Disposal	
P501	Dispose of contents / container in accordance with local regulations.

Safety Data Sheets are available at www.abbottdiagnostics.com or contact your local representative.

For a detailed discussion of safety precautions during system operation, refer to the Alinity ci-series Operations Manual, Section 8.

PREPARATION FOR USE

- This product is liquid ready-to-use.
- This product may be used immediately after removal from 2 to 8°C storage.
- Prior to each use, mix by gentle inversion.

STORAGE

- Do not use past expiration date.

	Storage Temperature	Maximum Storage Time	Additional Storage Instructions
Unopened	2 to 8°C	Until expiration date	
Opened	2 to 8°C	Until expiration date	Store tightly capped. Return to refrigerated storage after use.

INSTRUMENT PROCEDURE

- To obtain the recommended volume requirements for the controls, hold the bottle vertically, and dispense 8 drops of the control low avidity and 8 drops of the control high avidity into each sample cup in the assigned position.
- For instructions on ordering and loading controls on the instrument, refer to the Alinity ci-series Operations Manual, Section 5.

INDICATIONS OF INSTABILITY OR DETERIORATION

Instability or deterioration should be suspected if there are precipitates, visible signs of leakage, turbidity, or if controls do not meet the appropriate package insert and/or Alinity ci-series Operations Manual criteria.

BIBLIOGRAPHY

1. US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
2. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
3. World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
4. Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.

Note for number formatting:

- A space is used as thousands separator (example: 10 000 specimens).
- A period is used to separate the integer part from the fractional part of a number written in decimal form (example: 3.12%).

Key to Symbols

	Caution
	Consult instructions for use
	Manufacturer
	Temperature limitation
	Use by/Expiration date
	Control Number
	Contains Sodium Azide. Contact with acids liberates very toxic gas.
	Control High Avidity
	Control Low Avidity
	In Vitro Diagnostic Medical Device
	Lot Number
	Product of Germany
	Range
	List Number

Alinity is a trademark of Abbott Laboratories in various jurisdictions. All other trademarks are property of their respective owners.

Abbott GmbH & Co. KG
Max-Planck-Ring 2
65205 Wiesbaden
Germany
+49-6122-580

0123

Customer Service: Contact your local representative or find country-specific contact information on www.abbottdiagnostics.com

Revised July 2019.

©2016, 2019 Abbott Laboratories

ADDENDUM

to the Instruction for Use Alinity i Toxo IgG Avidity Controls for Russia

SHIPPING AND STORAGE CONDITIONS

Recommended shipping conditions for devices are 2-8°C. Medical devices must be stored in refrigerators/refrigerator chambers to ensure regulated temperature conditions with a daily temperature record.

Medical devices must be carried out in single-use thermal containers containing coolant elements, or in reusable thermo containers with automatically maintained temperature or in refrigerated trucks using thermal indicators.

Device is stable until the expiration date indicated on the package. Device is no longer suitable for use if shipping temperature limits have been violated. Ship in an any kind of enclosed carrier vehicle in accordance with shipping rules, applicable to the respective type of transport.

WARRANTY

Manufacturer guarantees that device remains stable until the expiration date when stored and handled as intended. Manufacturer set the maximum shelf life of

Alinity i Toxo IgG Avidity Controls – 12 months from date of manufacturing.

MAINTENANCE

Not required.

REPAIRS

Not required.

DESTRUCTION AND DISPOSAL OF A MEDICAL DEVICE

Destruct and dispose of a medical device in accordance with the requirements of SanPin 2.1.7.2790-10 dated December 12, 2010 "Sanitary-epidemiological Requirements for the Management of Medical Waste". Transmit and dispose of a product and its components and packaging in accordance with local regulations.

LIST OF NATIONAL STANDARDS AND REGULATORY DOCUMENTS APPLICABLE TO MEDICAL PRODUCT

GOST ISO 17511-2011, GOST R 51352-2013, GOST R ISO 23640-2015, GOST R ISO 18113-1-2015, GOST R ISO 18113-2-2015, GOST R ISO 15223-1-2014, GOST R ISO 15223 -1-2014, GOST R EN 13612-2010, GOST R 53022.2-2008, GOST R 53022.3-2008

Decree of the Government of the Russian Federation of December 27, 2012 No. 1416 "On approval of the rules of state registration of medical devices" (as amended by Decree of the Government of the Russian Federation of July 17, 2014 No. 670), Order of the Ministry of Health of Russia of 06.06.2012. No. 4n (as amended according to the Order of the Ministry of Health of the Russian Federation of September 25, 2014 N 557n), Order of the Ministry of Health of Russia of 01/09/2014. No. 2n, Order of the Ministry of Health of the Russian Federation No. 11n of January 19, 2017

POTENTIAL USER

The Alinity i Toxo IgG Avidity Controls is intended to be used in a variety of clinical and laboratory conditions, including, without limitation, hospitals, diagnostic centers, reference laboratories and blood transfusion facilities around the world, whose personnel have been trained in the use of the Alinity i analyzer.

INFORMATION ABOUT PRESENCE OF DRUGS IN COMPOSITION OF MEDICAL DEVICE

Does not contain any drugs

STERILIZATION

Medical device is not required to be sterile.

INFORMATION ABOUT POSSIBILITY OF MULTIPLE APPLICATION

Medical device is intended for one-off usage

MATERIALS OF ANIMAL ORIGIN

Materials of animal origin are tested and risks of using MD are indicated in SDS. Safety Data Sheets are available at www.abbottdiagnostics.com or contact your local representative.

INDICATIONS FOR USE

The Alinity i Toxo IgG Avidity Controls are for the estimation of test precision and the detection of systematic analytical deviations of the Alinity i analyzer.

CONTRAINDICATIONS FOR USE

Not applicable.

TECHNICAL CHARACTERISTICS

- Low Control: pH-value at 20 °C: 7.4
Density at 20 °C: 1.009 g/cm³
- High Control: pH-value at 20 °C: 7.5
Density at 20 °C: 1.025 g/cm³

COMPATIBILITY RESTRICTIONS OF MEDICAL DEVICE

No available data.

COLOR AND ODOR

Low Control: Colourless and odorless
High Control: Brown and odorless

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