Abbott

Abbott GmbH & Co. KG
Max-Planck-Ring 2
65205 Wiesbaden
Postfach 1303, 65011 Wiesbaden

Tel. +49 61 22 58 0
Fax +49 61 22 58 1244

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Signed: ZO
Zaman Khan,
Associate Director Regulatory
Affairs EMEA and Pakistan,
Abbott GmbH & Co. KG.

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en

Toxo IgG
07P45
G71279R03
B7P450

Read Highlighted Changes: Revised July 2019.

REF 07P4522

REF 07P4532

Instructions must be carefully followed. Reliability of assay results cannot be guaranteed if there are any deviations from these instructions.

NAME

Alinity i Toxo IgG Reagent Kit

INTENDED USE

The Alinity i Toxo IgG assay is a chemiluminescent microparticle immunoassay (CMIA) used for the quantitative determination of IgG antibodies to *Toxoplasma gondii* in human serum and plasma on the Alinity i analyzer.

SUMMARY AND EXPLANATION OF THE TEST

Toxoplasma gondii is an obligate intracellular protozoan parasite that infects most species of warm-blooded animals, including humans.¹

Toxoplasmosis is primarily acquired by ingestion of undercooked, infected meat; via oocysts from fecally contaminated hands, food and water; and maternally through transplacental transmission.² In addition, transmission associated with organ transplantation and during blood transfusion has been reported, although the risk of transmission through blood transfusion is extremely low.³

Acquired infection with *Toxoplasma gondii* in healthy individuals is commonly asymptomatic, however 10-20% of patients with acute infection may develop lymphadenopathy.⁴

Severe infections can occur in AIDS patients and adults immunocompromised by cancer chemotherapy or transplant recipients receiving immunosuppressive treatment. These infections can be fatal. Toxoplasmic encephalitis is the most common presentation and is the most frequent cause of focal central nervous system lesions in AIDS patients.⁵

Primary infection during pregnancy can result in transplacental transmission of the parasite resulting in congenital infection. The risk of congenital infection is lowest (10-25%) if acute maternal infection occurs during the first trimester and highest (60-90%) if it occurs during the third trimester.² Severity of congenital infection is greatest when maternal infection is acquired early during pregnancy. Common outcomes of congenital toxoplasmosis include chorioretinitis, intracranial calcifications, and hydrocephalus. The majority of infants infected later in pregnancy are asymptomatic at birth, with sequelae occurring later in life.

Early treatment after prenatal diagnosis of *Toxoplasma gondii* infection has been shown to reduce the frequency and severity of congenital toxoplasmosis.⁶ Serological tests can be used to identify seronegative women who then should be monitored during pregnancy.

The presence of IgG antibodies to *Toxoplasma gondii* indicates that infection has occurred but does not distinguish between recent and past infection. IgM antibodies are detected in individuals with a recently acquired infection, but antibodies may persist for up to 18 months post-infection.² To differentiate between a recently acquired and a past infection, IgM and IgG positive specimens should be tested for IgG avidity. A high avidity index for IgG antibodies is a strong indication that an infection took place more than 4 months ago.

Toxo IgG	Toxo IgM	Toxo IgG Avidity	May indicate... / Testing recommendation
nonreactive	nonreactive	N/A	no infection
nonreactive	reactive	N/A	obtain new sample 2-3 weeks after initial sample and test for Toxo IgG and Toxo IgM
reactive	nonreactive	high avidity	past infection. Strong indication that an infection took place more than 4 months ago
reactive	reactive	low avidity	obtain new sample 3 weeks after initial sample and test for Toxo IgG and Toxo IgM
reactive	reactive	high avidity	past infection. Strong indication that an infection took place more than 4 months ago

BIOLOGICAL PRINCIPLES OF THE PROCEDURE

This assay is a two-step immunoassay for the quantitative determination of IgG antibodies to *Toxoplasma gondii* in human serum and plasma using chemiluminescent microparticle immunoassay (CMIA) technology.

Pre-diluted sample, recombinant *Toxoplasma gondii* antigen (containing recombinant antigens P30(SAG1) and P35(GRA8)) coated paramagnetic microparticles, and assay diluent are combined and incubated. The *Toxoplasma gondii* specific antibodies present in the sample bind to the recombinant *Toxoplasma gondii* antigen coated microparticles. The mixture is washed. Murine anti-human IgG acridinium-labeled conjugate is added to create a reaction mixture and incubated. Following a wash cycle, Pre-Trigger and Trigger Solutions are added.

The resulting chemiluminescent reaction is measured as relative light units (RLUs). There is a direct relationship between the amount of anti-Toxo IgG in the sample and the RLUs detected by the system optics.

For additional information on system and assay technology, refer to the Alinity ci-series Operations Manual, Section 3.

REAGENTS

Kit Contents

Alinity i Toxo IgG Reagent Kit 07P45

NOTE: Some kit sizes are not available in all countries. Please contact your local distributor.

Volumes (mL) listed in the table below indicate the volume per cartridge.

REF	07P4522	07P4532
Tests per cartridge	100	500
Number of cartridges per kit	2	2
Tests per kit	200	1000
MICROPARTICLES	6.6 mL	27.0 mL
CONJUGATE	6.1 mL	26.5 mL
ASSAY DILUENT	10.4 mL	47.1 mL

MICROPARTICLES Recombinant *Toxoplasma gondii* antigen coated microparticles in MES buffer with protein (bovine) stabilizers. Minimum concentration: 0.03% solids. Preservative: ProClin 300.

CONJUGATE Murine anti-human IgG acridinium-labeled conjugate in MES buffer with protein (bovine) stabilizers. Minimum concentration: 0.05 µg/mL. Preservatives: antimicrobial agents.

ASSAY DILUENT TRIS buffer with protein (bovine) stabilizers. Preservative: ProClin 300.

Warnings and Precautions

- **IVD**
- For *In Vitro* Diagnostic Use

Safety Precautions

CAUTION: This product requires the handling of human specimens. It is recommended that all human-sourced materials be considered potentially infectious and handled in accordance with the OSHA Standard on Bloodborne Pathogens, Biosafety Level 2 or other appropriate biosafety practices should be used for materials that contain or are suspected of containing infectious agents.⁷⁻¹⁰

The following warnings and precautions apply to: **MICROPARTICLES** and **ASSAY DILUENT**

WARNING	Contains methylisothiazolones.
H317	May cause an allergic skin reaction.
Prevention	
P261	Avoid breathing mist / vapors / spray.
P272	Contaminated work clothing should not be allowed out of the workplace.
P280	Wear protective gloves / protective clothing / eye protection.
Response	
P302+P352	IF ON SKIN: Wash with plenty of water.
P333+P313	If skin irritation or rash occurs: Get medical advice / attention.
P362+P364	Take off contaminated clothing and wash it before reuse.
Disposal	
P501	Dispose of contents / container in accordance with local regulations.

Safety Data Sheets are available at www.abbottdiagnostics.com or contact your local representative.

For a detailed discussion of safety precautions during system operation, refer to the Alinity ci-series Operations Manual, Section 8.

Reagent Handling

- Upon receipt, gently invert the unopened reagent kit by rotating it over and back for a full 180 degrees, 5 times with green label stripe facing up and then 5 times with green label stripe facing down. This ensures that liquid covers all sides of the bottles within the cartridges. During reagent shipment, microparticles can settle on the reagent septum.
 - Place a check in the square on the reagent kit to indicate to others that the inversions have been completed.

- After mixing, place reagent cartridges in an upright position for 1 hour before use to allow bubbles that may have formed to dissipate.
- If a reagent cartridge is dropped, place in an upright position for 1 hour before use to allow bubbles that may have formed to dissipate.
- Reagents are susceptible to the formation of foam and bubbles. Bubbles may interfere with the detection of the reagent level in the cartridge and cause insufficient reagent aspiration that may adversely affect results.

For a detailed discussion of reagent handling precautions during system operation, refer to the Alinity ci-series Operations Manual, Section 7.

Reagent Storage

	Storage Temperature	Maximum Storage Time	Additional Storage Instructions
Unopened	2 to 8°C	Until expiration date	Store in upright position. If cartridge does not remain upright, gently invert the cartridge 10 times and place in an upright position for 1 hour before use.
Onboard			
100-Test Cartridge	System Temperature	30 days	Discard after 30 days.
500-Test Cartridge	System Temperature	30 days	Discard after 30 days.
Opened	2 to 8°C	Until expiration date	Store in upright position. If cartridge does not remain upright during storage, discard the cartridge. Do not reuse original reagent caps or replacement caps due to the risk of contamination and potential to compromise reagent performance.

Reagents may be stored on or off the system. If removed from the system, store reagents with new replacement caps in an upright position at 2 to 8°C. For reagents stored off the system, it is recommended that they be stored in their original trays or boxes to ensure they remain upright.

For information on unloading reagents, refer to the Alinity ci-series Operations Manual, Section 5.

Indications of Reagent Deterioration

Deterioration of the reagents may be indicated when a calibration error occurs or a control value is out of the specified range. Associated test results are invalid, and samples must be retested. Assay recalibration may be necessary.

For troubleshooting information, refer to the Alinity ci-series Operations Manual, Section 10.

INSTRUMENT PROCEDURE

The Alinity i Toxo IgG assay file must be installed on the Alinity i analyzer prior to performing the assay.

The following assay files are available for testing:

- "Toxo IgG" – which does not automatically dilute and retest specimens with an anti-Toxo IgG concentration of > 200 IU/mL.
- "Toxo IgG-R" – which automatically dilutes and retests specimens with an anti-Toxo IgG concentration of > 200 IU/mL.

For detailed information on assay file installation and viewing and editing assay parameters, refer to the Alinity ci-series Operations Manual, Section 2.

For information on printing assay parameters, refer to the Alinity ci-series Operations Manual, Section 5.

For a detailed description of system procedures, refer to the Alinity ci-series Operations Manual.

SPECIMEN COLLECTION AND PREPARATION FOR ANALYSIS

Specimen Types

The specimen types listed below were verified for use with this assay on the ARCHITECT i System.

Other specimen types and collection tube types have not been verified with this assay.

Specimen Types	Collection Tubes
Serum	Serum Serum separator
Plasma	Lithium heparin plasma separator Potassium EDTA Sodium citrate Lithium heparin Sodium heparin

- Performance has not been established for the use of cadaveric specimens or the use of bodily fluids other than human serum and plasma.
- Liquid anticoagulants may have a dilution effect resulting in lower concentration values for individual specimens.
- The instrument does not provide the capability to verify specimen types. It is the responsibility of the operator to verify that the correct specimen types are used in the assay.

Specimen Conditions

- Do not use:
 - heat-inactivated specimens
 - pooled specimens
 - grossly hemolyzed specimens
 - specimens with obvious microbial contamination
- For accurate results, serum and plasma specimens should be free of fibrin, red blood cells, and other particulate matter. Serum specimens from patients receiving anticoagulant or thrombolytic therapy may contain fibrin due to incomplete clot formation.
- To prevent cross contamination, use of disposable pipettes or pipette tips is recommended.

Preparation for Analysis

- Follow the tube manufacturer's processing instructions for collection tubes. Gravity separation is not sufficient for specimen preparation.
- Specimens should be free of bubbles. Remove bubbles with an applicator stick before analysis. Use a new applicator stick for each specimen to prevent cross contamination.

To ensure consistency in results, recentrifuge specimens prior to testing if

- they contain fibrin, red blood cells, or other particulate matter
- they require repeat testing.

NOTE: If fibrin, red blood cells, or other particulate matter are observed, mix by low speed vortex or by inverting 10 times prior to recentrifugation.

Prepare frozen specimens as follows:

- Frozen specimens must be completely thawed before mixing.
- Mix thawed specimens thoroughly by low speed vortex or by inverting 10 times.
- Visually inspect the specimens. If layering or stratification is observed, mix until specimens are visibly homogeneous.

- If specimens are not mixed thoroughly, inconsistent results may be obtained.
- Recentrifuge specimens.

Recentrifugation of Specimens

- Transfer specimens to a centrifuge tube and centrifuge at a minimum of 100 000 g-minutes.
- Examples of acceptable time and force ranges that meet this criterion are listed in the table below.

Centrifugation time using alternate RCF values can be calculated using the following formula:

$$\text{Minimum Centrifugation time (minutes)} = \frac{100\,000 \text{ g-minutes}}{\text{RCF}}$$

Recentrifugation Time (Minutes)	RCF (x g)	g-Minutes
10	10 000	100 000
20	5000	100 000
40	2500	100 000

$$\text{RCF} = 1.12 \times r_{\text{max}} (\text{rpm}/1000)^2$$

RCF - The relative centrifugal force generated during centrifugation.

rpm - The revolutions per minute of the rotor on which the specimens are being spun (usually the digital readout on the centrifuge will indicate the rpm).

Centrifugation Time - The time should be measured from the time the rotor reaches the required RCF or rpm to the time it begins decelerating.

r_{max} - Radius of the rotor in millimeters. NOTE: If custom tube adapters (i.e., adapters not defined by the centrifuge manufacturer) are used, then the radius (r_{max}) should be manually measured in millimeters and the RCF calculated.

g-minutes - The unit of measure for the product of RCF (x g) and centrifugation time (minutes).

- Transfer clarified specimen to a sample cup or secondary tube for testing. For centrifuged specimens with a lipid layer, transfer only the clarified specimen and not the lipemic material.

Specimen Storage

Specimen storage conditions were verified on the ARCHITECT i System.

Specimen Type	Temperature	Maximum Storage Time	Special Instructions
Serum/ Plasma	15 to 30°C	3 days	Specimens may be stored on or off the clot, red blood cells, or separator gel.
	2 to 8°C	14 days	Specimens may be stored on or off the clot, red blood cells, or separator gel.

Remove serum or plasma from the clot, red blood cells, or separator gel if stored longer than the maximum 15-30°C or maximum 2-8°C storage time and store frozen at -10°C or colder.

Avoid more than 6 freeze/thaw cycles.

Specimen Shipping

Package and label specimens in compliance with applicable state, federal, and international regulations covering the transport of clinical specimens and infectious substances.

PROCEDURE

Materials Provided

07P45 Alinity i Toxo IgG Reagent Kit

Materials Required but not Provided

- Alinity i Toxo IgG assay file
- 07P4501 Alinity i Toxo IgG Calibrators
- 07P4510 Alinity i Toxo IgG Controls or other control material
- Alinity Trigger Solution
- Alinity Pre-Trigger Solution
- Alinity i-series Concentrated Wash Buffer

For information on materials required for operation of the instrument, refer to the Alinity ci-series Operations Manual, Section 1.

For information on materials required for maintenance procedures, refer to the Alinity ci-series Operations Manual, Section 9.

Assay Procedure

For a detailed description of how to run an assay, refer to the Alinity ci-series Operations Manual, Section 5.

- If using primary or aliquot tubes, refer to the Alinity ci-series Operations Manual, Section 4 to ensure sufficient specimen is present.
- To minimize the effects of evaporation, verify adequate sample cup volume is present prior to running the test.
- Maximum number of replicates sampled from the same sample cup: 10
 - Priority:
 - Sample volume for first test: 75 µL
 - Sample volume for each additional test from same sample cup: 25 µL
 - ≤ 3 hours on the reagent and sample manager:
 - Sample volume for first test: 150 µL
 - Sample volume for each additional test from same sample cup: 25 µL
 - > 3 hours on the reagent and sample manager:
 - Replace with a fresh aliquot of sample.
- Refer to the Alinity i Toxo IgG calibrator package insert and/or Alinity i Toxo IgG control package insert for preparation and usage.
- For general operating procedures, refer to the Alinity ci-series Operations Manual, Section 5.
- For optimal performance, it is important to perform routine maintenance as described in the Alinity ci-series Operations Manual, Section 9. Perform maintenance more frequently when required by laboratory procedures.

Sample Dilution Procedures

Samples with an anti-Toxo IgG concentration value exceeding 200.0 IU/mL are flagged with the code "> 200.0 IU/mL" and may be diluted with the Automated Dilution Protocol.

Automated Dilution Protocol

The system performs a 1:10 dilution of the sample and automatically calculates the concentration by multiplying the result by the dilution factor.

When testing is conducted using the "Toxo IgG-R" assay file, specimens flagged as "> 200.0 IU/mL" will be automatically retested in 1:10 dilution.

For detailed information on ordering dilutions, refer to the Alinity ci-series Operations Manual, Section 5.

Calibration

For instructions on performing a calibration, refer to the Alinity ci-series Operations Manual, Section 5.

Each assay control must be tested to evaluate the assay calibration. Once a calibration is accepted and stored, all subsequent samples may be tested without further calibration unless:

- A reagent kit with a new lot number is used.
- Daily quality control results are outside of statistically-based quality control limits used to monitor and control system performance, as described in the Quality Control Procedures section of this package insert.
 - If statistically-based quality control limits are not available, then the calibration should not exceed a 30-day limit for recalibration frequency.

This assay may require recalibration after maintenance to critical parts or subsystems or after service procedures have been performed.

Quality Control Procedures

The recommended control requirement for the Alinity i Toxo IgG assay is that a single sample of each control level be tested once every 24 hours each day of use.

Additional controls may be tested in accordance with local, state, and/or federal regulations or accreditation requirements and your laboratory's quality control policy.

To establish statistically-based control limits, each laboratory should establish its own concentration target and ranges for new control lots at each clinically relevant control level. This can be accomplished by assaying a minimum of 20 replicates over several (3-5) days and using the reported results to establish the expected average (target) and variability about this average (range) for the laboratory. Sources of variation that should be included in this study in order to be representative of future system performance include:

- Multiple stored calibrations
 - Multiple reagent lots
 - Multiple calibrator lots
 - Multiple processing modules (if applicable)
 - Data points collected at different times of the day
- Refer to published guidelines for information or general control recommendation, for example Clinical and Laboratory Standards Institute (CLSI) Document C24-A3 or other published guidelines, for general quality control recommendations.¹¹
- If more frequent control monitoring is required, follow the established quality control procedures for your laboratory.
 - If quality control results do not meet the acceptance criteria defined by your laboratory, sample results may be suspect. Follow the established quality control procedures for your laboratory. Recalibration may be necessary. For troubleshooting information, refer to the Alinity ci-series Operations Manual, Section 10.
 - Review quality control results and acceptance criteria following a change of reagent or calibrator lot.

Quality Control Guidance

Refer to "Basic QC Practices" by James O Westgard, Ph.D. for guidance on laboratory quality control practices.¹²

Verification of Assay Claims

For protocols to verify package insert claims, refer to Verification of Assay Claims in the Alinity ci-series Operations Manual.

RESULTS

Calculation

The Alinity i Toxo IgG assay utilizes a 4 Parameter Logistic Curve fit data reduction method (4PLC, Y-weighted) to generate a calibration and results.

Interpretation of Results

Concentration values	Instrument Interpretation	Interpretation of Results and Retest Procedure
< 1.6 IU/mL	Nonreactive	Individuals with such results are presumed to be not infected with <i>Toxoplasma gondii</i> and susceptible to acute infection. A negative result does not always exclude the possibility of <i>Toxoplasma gondii</i> infection. Patients with negative results in suspected early disease cases should be retested in 3 weeks.
1.6 to < 3.0 IU/mL	Grayzone	Specimens that are considered grayzone may contain low levels of IgG. It is recommended to test those specimens using a Toxo IgM test, and/or a second sample should be taken within a reasonable period of time (e.g., two weeks) and used to repeat Alinity i Toxo IgG testing.
≥ 3.0 IU/mL	Reactive	Specimens that are considered reactive for IgG antibodies to <i>Toxoplasma gondii</i> indicate past or acute infection.

For details on configuring the Alinity i analyzer to use grayzone interpretations, refer to the Alinity ci-series Operations Manual, Section 2.

Flags

Some results may contain information in the Flags field. For a description of the flags that may appear in this field, refer to the Alinity ci-series Operations Manual, Section 5.

Measuring Interval

Measuring interval is defined as the range of values in IU/mL which meets the limits of acceptable performance for imprecision and Limit of Quantitation (LoQ).

The measuring interval of the Alinity i Toxo IgG assay is 0.2 to 200.0 IU/mL.

LIMITATIONS OF THE PROCEDURE

- If the Alinity i Toxo IgG results are inconsistent with clinical evidence, additional testing is suggested to confirm the result.
- For diagnostic purposes, results should be used in conjunction with other data; e.g., results of other tests (Toxo IgM, Toxo IgG Avidity), clinical impressions, etc.
- Heterophilic antibodies in human serum can react with reagent immunoglobulins, interfering with *in vitro* immunoassays. Patients routinely exposed to animals or to animal serum products can be prone to this interference, and anomalous values may be observed. Additional information may be required for diagnosis.¹³
- Specimens from patients who have received preparations of mouse monoclonal antibodies for diagnosis or therapy may contain human anti-mouse antibodies (HAMA). Such specimens may show either falsely elevated or depressed values when tested with assay kits such as Alinity i Toxo IgG that employ mouse monoclonal antibodies. Additional information may be required for diagnosis.^{14, 15}
Alinity i Toxo IgG reagents contain a component that reduces the effect of HAMA reactive specimens.
- Assay to assay variation in results: **The concentration values for Toxo IgG in a given specimen can vary based on the assay method and standardization and should not be used interchangeably.** Note: The Alinity i Toxo IgG Calibrators are referenced to the World Health Organization (WHO) First International Standard (01/600) for Anti-Toxoplasma IgG.¹⁶
- IgG antibodies to *Toxoplasma gondii* might not dilute linearly in all samples. Do not use different dilution protocols when analyzing sequential samples taken from an individual.

EXPECTED VALUES

This study was performed on the ARCHITECT i System.

Representative performance data are provided in this section. Results obtained in individual laboratories may vary.

It is recommended that each laboratory determine its own reference range based upon its particular locale and population characteristics. The prevalence of Toxo IgG antibody to *Toxoplasma gondii* will vary with age and geographic location. In this study 1270 specimens from pregnant women and 1297 specimens from random individuals were tested. Of these specimens 1115 (43.4%) were positive, 83 (3.2%) were grayzone and 1369 (53.3%) were nonreactive by the ARCHITECT Toxo IgG assay.

IU/mL	n (%) Overall Specimen Categories	n (%) Blood Donor	n (%) Diagnostic/Hospitalized	n (%) Pregnant Women
0.0 to < 1.6	1369 (53.3)	278 (45.8)	270 (39.1)	821 (64.6)
1.6 to < 3.0	83 (3.2)	16 (2.6)	38 (5.5)	29 (2.3)
3.0 to < 10.0	450 (17.5)	134 (22.1)	161 (23.3)	155 (12.2)
10.0 to < 50.0	554 (21.6)	160 (26.4)	182 (26.4)	212 (16.7)
50.0 to < 100.0	59 (2.3)	9 (1.5)	26 (3.8)	24 (1.9)
100.0 to < 150.0	22 (0.9)	5 (0.8)	6 (0.9)	11 (0.9)
150.0 to < 200.0	7 (0.3)	2 (0.3)	0 (0.0)	5 (0.4)
> 200.0	23 (0.9)	3 (0.5)	7 (1.0)	13 (1.0)
Total	2567 (100.0)	607 (100.0)	690 (100.0)	1270 (100.0)

SPECIFIC PERFORMANCE CHARACTERISTICS

Representative performance data are provided in this section. Results obtained in individual laboratories may vary.

The Alinity i analyzer and the ARCHITECT i System utilize the same reagents and sample/reagent ratios.

Unless otherwise specified, all studies were performed on the Alinity i analyzer.

Precision

Within-Laboratory Precision

A study was performed based on guidance from CLSI EP05-A2.¹⁷ Testing was conducted using 3 lots of the Alinity i Toxo IgG Reagent Kit, 3 lots of the Alinity i Toxo IgG Calibrators, and 3 lots of the Alinity i Toxo IgG Controls and 1 instrument. Three controls and 4 human plasma panels were assayed in a minimum of 2 replicates at 2 separate times per day on 20 different days.

Sample	n	Within-Run (Repeatability)			Within-Laboratory (Total) ^a	
		Mean (IU/mL)	SD	%CV	SD (Range ^b)	%CV (Range ^b)
Negative Control	357	0.0	0.03	95.8	0.05 (0.02 - 0.05)	166.5 (116.2 - 187.8)
Positive Control 1	357	5.9	0.15	2.5	0.25 (0.24 - 0.27)	4.2 (4.0 - 4.6)
Positive Control 2	360	111.2	3.74	3.4	7.62 (5.96 - 9.06)	6.8 (5.5 - 8.0)
Panel 1	356	1.2	0.04	3.3	0.07 (0.06 - 0.07)	5.3 (4.6 - 5.8)
Panel 2	355	3.0	0.08	2.6	0.13 (0.12 - 0.15)	4.4 (3.7 - 5.0)
Panel 3	357	44.6	1.07	2.4	2.09 (1.87 - 2.39)	4.7 (4.1 - 5.3)
Panel 4	359	145.6	4.97	3.4	9.99 (8.84 - 11.29)	6.9 (6.3 - 7.6)

^a Includes within-run, between-run, and between-day variability.

^b Maximum and minimum SD or %CV for each reagent lot and instrument combination.

Lower Limits of Measurement

A study was performed based on guidance from CLSI EP17-A2. Testing was conducted using 3 lots of the Alinity i Toxo IgG Reagent Kit on each of 2 instruments over a minimum of 3 days. The maximum observed Limit of Blank (LoB), Limit of Detection (LoD), and Limit of Quantitation (LoQ) values are summarized below.¹⁸

	IU/mL
LoB ^a	0.1
LoD ^b	0.2
LoQ ^c	0.2

^a The LoB represents the 95th percentile from $n \geq 60$ replicates of zero-analyte samples.

^b The LoD represents the lowest concentration at which the analyte can be detected with 95% probability based on $n \geq 60$ replicates of low-analyte level samples.

^c The LoQ was determined from $n \geq 60$ replicates of low-analyte level samples and is defined as the lowest concentration at which a maximum allowable precision of 20 %CV was met.

Resolved Relative Sensitivity, Specificity and Relative Agreement

This study was performed on the ARCHITECT i System.

Resolved relative sensitivity/specificity and relative agreement was assessed on 2464 specimens from pregnant females, diagnostic/hospitalized specimens and randomly selected volunteer blood donors.

103 specimens giving grayzone results using ARCHITECT Toxo IgG assay and/or the comparator Toxo IgG assay and/or any other Toxo assay were not included in the calculation of resolved relative sensitivity, specificity and relative agreement.

Resolved Relative Sensitivity

From 2464 specimens evaluated, 1099 were classified as positive. 1096 were reactive by ARCHITECT Toxo IgG and 3 specimens were nonreactive by ARCHITECT Toxo IgG. The resolved relative sensitivity was 99.7% (1096/1099) with a 95% confidence interval of 99.2% to 99.9%.

Resolved Relative Specificity

From 2464 specimens evaluated, 1365 were classified as negative. 1359 were nonreactive by ARCHITECT Toxo IgG and 6 specimens were reactive by ARCHITECT Toxo IgG. The resolved relative specificity was 99.6% (1359/1365) with a 95% confidence interval of 99.0% to 99.8%.

Relative Agreement

From the 2464 specimens evaluated, 12 specimens were tested discordant between ARCHITECT Toxo IgG and the comparator Toxo IgG assay resulting in a relative agreement of 99.5% (2452/2464) with a 95% confidence interval of 99.2% to 99.7%.

Comparator Toxo IgG Assay	ARCHITECT Toxo IgG			Total
	POS	REA	NR	
POS		1094	4**	1098
NEG		8*	1358	1366
Total		1102	1362	2464

REA = reactive, NR = nonreactive, POS = positive, NEG = negative

* Two specimens out of 8 specimens reactive on ARCHITECT Toxo IgG and nonreactive on the comparator Toxo IgG assay were confirmed reactive by testing with a commercially-available assay, the Sabin-Feldman Dye Test and the HS Agglutination test. Six specimens out of 8 specimens reactive on ARCHITECT Toxo IgG and nonreactive on the comparator Toxo IgG assay were confirmed nonreactive by testing with a commercially-available assay, the Sabin-Feldman Dye Test and the HS Agglutination test. Four specimens out of these 6 unconfirmed specimens showed reactivity to GRA8 (p35) on a commercially-available blot.

** One specimen out of 4 specimens reactive on the comparator Toxo IgG assay and nonreactive on ARCHITECT Toxo IgG could not be confirmed by additional testing as outlined above whereas 3 specimens were confirmed reactive.

Negative and Positive Percent Agreement

A study was performed to determine the negative and positive percent agreement between the Alinity i Toxo IgG assay and a commercially-available Toxo IgG assay by testing 899 specimens from a population consisting of pregnant women, hospitalized individuals, and blood donors. The specimens were tested in singlicate using 3 Alinity i Toxo IgG Reagent Kit lots and up to 3 lots of Alinity i Toxo IgG Calibrators. The specimens were also tested in singlicate using a commercially-available Toxo IgG assay. Specimens that showed a grayzone result interpretation on either the Alinity i Toxo IgG assay or the commercially-available Toxo IgG assay were not included in the evaluation.

The negative percent agreement ranged from 99.77% (425/426) to 100.00% (424/424; 426/426). The 95% confidence interval of the commercially-available Toxo IgG assay assumed specificity of 100% ranged from 99.14% to 100.00%. The positive percent agreement was at 100.00% (460/460) for all tested lots. The 95% confidence interval of the commercially-available Toxo IgG assay assumed sensitivity of 100% ranged from 99.20% to 100.00%.

Alinity i Toxo IgG	Commercially-Available Toxo IgG Assay	
	Reactive	Nonreactive
Reagent Lot 1		
Reactive	460	0
Nonreactive	0	424
Reagent Lot 2		
Reactive	460	1
Nonreactive	0	425
Reagent Lot 3		
Reactive	460	0
Nonreactive	0	426

Reagent Lot 1: Positive % agreement: 100.00% (460/460), 95% CI (99.20%, 100.00%); Negative % agreement: 100.00% (424/424), 95% CI (99.13%, 100.00%); 15 grayzone results were excluded

Reagent Lot 2: Positive % agreement: 100.00% (460/460), 95% CI (99.20%, 100.00%); Negative % agreement: 99.77% (425/426), 95% CI (98.70%, 99.99%); 13 grayzone results were excluded

Reagent Lot 3: Positive % agreement: 100.00% (460/460), 95% CI (99.20%, 100.00%); Negative % agreement: 100.00% (426/426), 95% CI (99.14%, 100.00%); 13 grayzone results were excluded

Seroconversion Sensitivity

This study was performed on the ARCHITECT i System.

A total of 50 bleeds from 12 different seroconversion panels were tested. ARCHITECT Toxo IgG showed comparable seroconversion sensitivity to a comparator Toxo IgG assay on the panels tested in this study.

During an acute infection, the ARCHITECT Toxo IgG assay typically shows a substantial rise in anti-Toxo IgG concentrations in consecutive draws. Data from selected seroconversion panels are shown in the following table.

Sample ID	Months after last negative bleed	ARCHITECT		Sabin-Feldman		HS Agglutination Test IU/mL	ISAGA (Toxo-M) Index
		Toxo IgG IU/mL	Comparator Toxo IgG IU/mL	Dye Test IU/mL	Dye Test IU/mL		
GZ range		1.6 - 2.9	2.0 - 2.9			1	6 - 8
Reactive cutoff		3.0	3.0	2	2	2	9
29944022	0.0	0.8	0.3	< 2	< 1		0
29944023	1.1	5.0	19.4	40	8		12
29944024	1.4	18.7	54.8	200	16		12
29944025	1.7	21.9	60.8	400	32		12
29944029	0.0	0.6	0.4	< 2	< 1		0
29944030	1.2	6.8	7.6	10	2		12
29944031	0.0	0.4	0.4	< 2	< 1		0
29944032	2.1	172.6	482.2	1600	200		12
29944033	4.1	141.2	466.5	1600	400		12
29944034	8.4	39.6	84.8	400	50		12
29944035	0.0	0.4	0.1	< 2	< 1		0
29944036	2.3	1.6	2.3	5	1		12
29944037	2.4	5.8	6.7	10	1		12

Sample ID	Months after last negative bleed	ARCHITECT Toxo IgG IU/mL	Comparator Toxo IgG IU/mL	Sabin-Feldman Dye Test IU/mL	HS Agglutination Test IU/mL	ISAGA (Toxo-M) Index
29944038	0.0	0.4	0.1	< 2	< 1	0
29944039	1.0	9.1	9.9	10	2	12
29944040	0.0	0.3	0.3	< 2	< 1	0
29944041	1.6	15.1	21.8	40	16	12
29944042	11.8	21.9	28.1	80	32	12
29944043	0.0	0.5	0.4	< 2	< 1	0
29944044	1.5	3.2	2.9	10	1	11
29944045	1.8	9.4	10.1	20	2	11
00006214	0.0	0.8	0.4	< 2	< 1	0
00006215	3.5	46.3	140.9	400	200	12
00006216	4.2	125.2	342.9	1600	400	12
00006217	0.0	0.5	0.2	< 2	< 1	0
00006218	3.2	28.4	38.1	100	16	12
00006219	7.0	20.3	18.6	100	100	12

This study was performed on the Alinity i analyzer.

One seroconversion panel of 20 bleeds was tested using 3 lots of the Alinity i Toxo IgG Reagent Kit. Testing with each reagent lot was performed on 1 of 3 Alinity i analyzers. The seroconversion panel bleeds were also tested using a commercially-available Toxo IgG assay. The Alinity i Toxo IgG and commercially-available assays exhibited equivalent Toxo IgG detection.

Panel ID	Days Since 1st Bleed	Alinity i Toxo IgG (IU/mL)			Commercially-Available Toxo IgG Assay (IU/mL)
		Reagent Lot 1	Reagent Lot 2	Reagent Lot 3	Reactive ≥ 3.0
461957152-01	0	0.1	0.1	0.2	0.1
461957152-02	7	0.1	0.1	0.1	0.1
461957152-03	9	0.1	0.1	0.2	0.1
461957152-04	14	0.2	0.2	0.3	0.1
461957152-05	16	0.3	0.2	0.3	0.2
461957152-06	23	6.7	5.8	6.1	6.2
461957152-07	27	17.6	15.4	16.8	17.5
461957152-08	30	28.5	24.9	25.8	26.5
461957152-09	49	114.5	102.4	114.8	117.0
461957152-10	53	129.2	117.8	125.1	137.5
461957152-11	86	194.5	167.6	200.6	189.2
461957152-12	88	176.4	163.2	171.0	189.8
461957152-13	93	170.2	151.1	172.2	178.3
461957152-14	95	154.0	156.9	175.3	166.7
461957152-15	100	178.1	163.3	181.4	188.7
461957152-16	104	172.5	173.5	176.9	180.6
461957152-17	107	170.4	156.6	173.3	189.4
461957152-18	112	171.0	157.0	176.9	175.3
461957152-19	114	146.1	159.7	160.0	162.3
461957152-20	119	151.4	147.1	160.9	162.5

Other Potential Interferents

This study was performed on the ARCHITECT i System.

Studies were performed to evaluate other potential interfering disease states on the ARCHITECT Toxo IgG assay. Eight specimens grayzone on either ARCHITECT Toxo IgG or a comparator Toxo IgG assay were not included in the calculation of the relative agreement.

Interfering Substance	n	Relative Agreement between ARCHITECT Toxo IgG and Comparator Toxo IgG Assay
Anti-nuclear antibody	8	100 %
Cytomegalovirus	10	100 %
Epstein-Barr Virus	10	100 %
Influenza Vaccines	10	90 %
HAMA	10	100 %
Herpes Simplex Virus 1	10	80 %
Herpes Simplex Virus 2	9	100 %
Hyperpolyclonal IgG	9	100 %
Hyperpolyclonal IgM	10	100 %
Measles	10	100 %
Parvovirus B19	10	100 %
Rheumatoid Factor	10	100 %
Rubella	8	100 %
Syphilis	10	100 %
Varicella Zoster Virus	8	100 %
anti ds DNA antibodies	7	100 %
Monoclonal IgG	10	90 %
Monoclonal IgM	9	100 %

After discordant resolution of ARCHITECT Toxo IgG nonreactive or reactive results, 2 discordant specimens from patients infected with anti-HSV-1 and 1 discordant specimen containing monoclonal IgG were confirmed by an additional commercially-available assay, the Sabin-Feldman Dye Test, HS Agglutination test and/or a commercially-available blot.

Interference

This study was performed on the ARCHITECT i System.

Potentially Interfering Endogenous Substances

No interference was observed between experimental controls and nonreactive or reactive specimens tested with elevated levels of the following potentially interfering substances.

Potentially Interfering Substance	Interferent Level
Bilirubin	≤ 20 mg/dL
Triglycerides	≤ 3000 mg/dL
Protein	≤ 12 g/dL
Hemoglobin	≤ 500 mg/dL

BIBLIOGRAPHY

- Willis MS, Southern P, Latimer MJ. Toxoplasma infection. Making the best use of laboratory tests. *Infect Med* 2002; 19: 522-532.
- Jones JL et al. Congenital Toxoplasmosis: A Review. *CME Review Article* Vol. 56, Number 5 2001; 296-305.
- Shulman IA. Parasitic infections and their impact on blood donor selection and testing. *Arch Pathol Lab Med*. 1994; 118: 366-370.
- Montoya J.G. Laboratory Diagnosis of *Toxoplasma gondii* Infection and Toxoplasmosis. *Journal of Infectious Diseases* 2002; 185 (Suppl 1) 73-82.
- Israelski DM, Remington JS. Toxoplasmosis in the non-AIDS immunocompromised host. *Curr Clin Top Infect Dis* 1993; 13: 322-356.
- Wong SY and Remington. Toxoplasmosis in Pregnancy. *Clinical Infectious Diseases* 1994; 18:853-862.
- US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
- US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
- World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
- Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.
- Clinical and Laboratory Standards Institute (CLSI). *Statistical Quality Control for Quantitative Measurement Procedures: Principles and Definitions; Approved Guideline—Third Edition*. CLSI Document C24-A3. Wayne, PA: CLSI; 2006.
- Westgard JO. *Basic QC Practices*. 3rd ed. Madison, WI: Westgard Quality Corporation; 2010.
- Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. *Clin Chem* 1988;34(1):27-33.
- Primus FJ, Kelley EA, Hansen HJ, et al. "Sandwich"-type immunoassay of carcinoembryonic antigen in patients receiving murine monoclonal antibodies for diagnosis and therapy. *Clin Chem* 1988;34(2):261-264.
- Schroff RW, Foon KA, Beatty SM, et al. Human anti-murine immunoglobulin responses in patients receiving monoclonal antibody therapy. *Cancer Res* 1985;45(2):879-885.
- Rigsby P, Rijpkema S, Guy EC, et al. Evaluation of a candidate international standard preparation for human anti-Toxoplasma immunoglobulin G. *J Clin Microbiol* 2004;42(11):5133-5138.
- Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline—Second Edition*. CLSI Document EP05-A2. Wayne, PA: CLSI; 2004.
- Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition*. CLSI Document EP17-A2. Wayne, PA: CLSI; 2012.

Note for number formatting:

- A space is used as thousands separator (example: 10 000 specimens).
- A period is used to separate the integer part from the fractional part of a number written in decimal form (example: 3.12%).

Key to Symbols

	Consult instructions for use
	Manufacturer
	Sufficient for
	Temperature limitation
	Use by/Expiration date
ASSAY DILUENT	Assay Diluent
CONJUGATE	Conjugate
INVERSIONS PERFORMED	Inversions Performed
IVD	In Vitro Diagnostic Medical Device
LOT	Lot Number
MICROPARTICLES	Microparticles
PRODUCT OF GERMANY	Product of Germany
REF	List Number
SN	Serial number

Alinity and ARCHITECT are trademarks of Abbott Laboratories in various jurisdictions. All other trademarks are property of their respective owners.

Abbott GmbH & Co. KG
Max-Planck-Ring 2
65205 Wiesbaden
Germany
+49-6122-580

CE
0123

Customer Service: Contact your local representative or find country-specific contact information on www.abbottdiagnostics.com

Revised July 2019.

©2017, 2019 Abbott Laboratories

ADDENDUM

to the Instruction for Use Alinity i Toxo IgG Reagent Kit for Russia

SHIPPING AND STORAGE CONDITIONS

Recommended shipping conditions for devices are 2-8°C. Medical devices must be stored in refrigerators/refrigerator chambers to ensure regulated temperature conditions with a daily temperature record.

Medical devices must be carried out in single-use thermal containers containing coolant elements, or in reusable thermo containers with automatically maintained temperature or in refrigerated trucks using thermal indicators.

Device is stable until the expiration date indicated on the package. Device is no longer suitable for use if shipping temperature limits have been violated. Ship in an any kind of enclosed carrier vehicle in accordance with shipping rules, applicable to the respective type of transport.

WARRANTY

Manufacturer guarantees that device remains stable until the expiration date when stored and handled as intended. Manufacturer set the maximum shelf life of

Alinity i Toxo IgG Reagent Kit – 12 months from date of manufacturing.

MAINTENANCE

Not required.

REPAIRS

Not required.

DESTRUCTION AND DISPOSAL OF A MEDICAL DEVICE

Destruct and dispose of a medical device in accordance with the requirements of SanPin 2.1.7.2790-10 dated December 12, 2010 "Sanitary-epidemiological Requirements for the Management of Medical Waste". Transmit and dispose of a product and its components and packaging in accordance with local regulations.

LIST OF APPLICABLE STANDARDS

98/79/EC, 1272/2008/EC (CLP), 1907/2006/EC (REACH), 80/181/EEC, 2009/886/EC, 2009/3/EC, 1069/2009, Regulation (EC) 142/2011, CLASS 2013, GHS, WHSR Part 7.1, EN ISO 9001, ISO 13485, EN ISO 14971, EN 13641, EN 13612, EN ISO 15223-1, EN 980 EN ISO 18113-1, EN ISO 18113-2, EN ISO 18113-3, EN ISO 23640*, EN ISO 20776-1, EN 13975, EN ISO 17511, EN ISO 15193, EN ISO 15194, ISO 80000-1

POTENTIAL USER

The Alinity i Toxo IgG Reagent Kit is intended to be used in a variety of clinical and laboratory conditions, including, without limitation, hospitals, diagnostic centers, reference laboratories and blood transfusion facilities around the world, whose personnel have been trained in the use of the Alinity i analyzer.

INFORMATION ABOUT PRESENCE OF DRUGS IN COMPOSITION OF MEDICAL DEVICE

Does not contain any drugs

STERILIZATION

Medical device is not required to be sterile.

INFORMATION ABOUT POSSIBILITY OF MULTIPLE APPLICATION

Medical device is intended for one-off usage

MATERIALS OF ANIMAL ORIGIN

Materials of animal origin are tested and risks of using MD are indicated in SDS. Safety Data Sheets are available at www.abbottdiagnostics.com or contact your local representative.

INDICATIONS AND CONTRAINDICATIONS FOR USE

Not applicable.

TECHNICAL CHARACTERISTICS

- Microparticles: pH-value at 20 °C: 6,5
Density at 20 °C: 1.018 g/cm³
- Conjugate: pH-value at 20 °C: 5.8
Density at 20 °C: 1.03 g/cm³
- Assay Diluent: pH-value at 20 °C: 7.5
Density at 20 °C: 1.026 g/cm³

COMPATIBILITY RESTRICTIONS OF MEDICAL DEVICE

No available data.

COLOUR AND ODOUR

Microparticles: Brown and odorless
Conjugate: Colourless and odourless
Assay Diluent: Colourless and odourless

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Signed: Zaman Khan
Zaman Khan,
Associate Director Regulatory
Affairs EMEA and Pakistan,
Abbott GmbH & Co. KG.

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- **IVD**
- For *In Vitro* Diagnostic Use

Safety Precautions

CAUTION: This product requires the handling of human specimens. It is recommended that all human-sourced materials be considered potentially infectious and handled in accordance with the OSHA Standard on Bloodborne Pathogens. Biosafety Level 2 or other appropriate biosafety practices should be used for materials that contain or are suspected of containing infectious agents.⁷⁻¹⁰

The following warnings and precautions apply to: MICROPARTICLES and ASSAY DILUENT	
	
WARNING	Contains methylisothiazolones.
H317	May cause an allergic skin reaction.
Prevention	
P261	Avoid breathing mist / vapors / spray.
P272	Contaminated work clothing should not be allowed out of the workplace.
P280	Wear protective gloves / protective clothing / eye protection.
Response	
P302+P352	IF ON SKIN: Wash with plenty of water.
P333+P313	If skin irritation or rash occurs: Get medical advice / attention.
P362+P364	Take off contaminated clothing and wash it before reuse.
Disposal	
P501	Dispose of contents / container in accordance with local regulations.

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- If specimens are not mixed thoroughly, inconsistent results may be obtained.

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- Multiple stored calibrations
- Multiple reagent lots
- Multiple calibrator lots
- Multiple processing modules (if applicable)
- Data points collected at different times of the day

Refer to published guidelines for information or general control recommendation, for example Clinical and Laboratory Standards Institute (CLSI) Document C24-A3 or other published guidelines, for general quality control recommendations.¹¹

- If more frequent control monitoring is required, follow the established quality control procedures for your laboratory.
- If quality control results do not meet the acceptance criteria defined by your laboratory, sample results may be suspect. Follow the established quality control procedures for your laboratory. Recalibration may be necessary. For troubleshooting information, refer to the Alinity ci-series Operations Manual, Section 10.
- Review quality control results and acceptance criteria following a change of reagent or calibrator lot.

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BIBLIOGRAPHY

- Willis MS, Southern P, Latimer MJ. Toxoplasma infection. Making the best use of laboratory tests. *Infect Med* 2002; 19: 522-532.
- Jones JL et al. Congenital Toxoplasmosis: A Review. *CME Review Article* Vol. 56, Number 5 2001; 296-305.
- Shulman IA. Parasitic infections and their impact on blood donor selection and testing. *Arch Pathol Lab Med.* 1994; 118: 366-370.
- Montoya J.G. Laboratory Diagnosis of Toxoplasma gondii Infection and Toxoplasmosis. *Journal of Infectious Diseases* 2002; 185 (Suppl 1) 73-82.
- Israelski DM, Remington JS. Toxoplasmosis in the non-AIDS immunocompromised host. *Curr Clin Top Infect Dis* 1993; 13: 322-356.
- Wong SY and Remington. Toxoplasmosis in Pregnancy. *Clinical Infectious Diseases* 1994; 18:853-862.
- US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
- US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
- World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
- Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.
- Clinical and Laboratory Standards Institute (CLSI). *Statistical Quality Control for Quantitative Measurement Procedures: Principles and Definitions; Approved Guideline—Third Edition*. CLSI Document C24-A3. Wayne, PA: CLSI; 2006.
- Westgard JO. *Basic QC Practices*. 3rd ed. Madison, WI: Westgard Quality Corporation; 2010.
- Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. *Clin Chem* 1988;34(1):27-33.
- Primus FJ, Kelley EA, Hansen HJ, et al. "Sandwich"-type immunoassay of carcinoembryonic antigen in patients receiving murine monoclonal antibodies for diagnosis and therapy. *Clin Chem* 1988;34(2):261-264.
- Schroff RW, Foon KA, Beatty SM, et al. Human anti-murine immunoglobulin responses in patients receiving monoclonal antibody therapy. *Cancer Res* 1985;45(2):879-885.
- Rigsby P, Rijpkema S, Guy EC, et al. Evaluation of a candidate international standard preparation for human anti-Toxoplasma immunoglobulin G. *J Clin Microbiol* 2004;42(11):5133-5138.
- Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline—Second Edition*. CLSI Document EP05-A2. Wayne, PA: CLSI; 2004.
- Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition*. CLSI Document EP17-A2. Wayne, PA: CLSI; 2012.

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Device is stable until the expiration date indicated on the package. Device is no longer suitable for use if shipping temperature limits have been violated. Ship in any kind of enclosed carrier vehicle in accordance with shipping rules, applicable to the respective type of transport.

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MATERIALS OF ANIMAL ORIGIN

Materials of animal origin are tested and risks of using MD are indicated in SDS. Safety Data Sheets are available at www.abbottiagnostics.com or contact your local representative.

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- Willis MS, Southern P, Latimer MJ. Toxoplasma infection. Making the best use of laboratory tests. *Infect Med* 2002; 19: 522-532.
- Jones JL et al. Congenital Toxoplasmosis: A Review. *CME Review Article* Vol. 56, Number 5 2001; 296-305.
- Shulman IA. Parasitic infections and their impact on blood donor selection and testing. *Arch Pathol Lab Med.* 1994; 118: 366-370.
- Montoya J.G. Laboratory Diagnosis of Toxoplasma gondii Infection and Toxoplasmosis. *Journal of Infectious Diseases* 2002; 185 (Suppl 1) 73-82.
- Israeli DM, Remington JS. Toxoplasmosis in the non-AIDS immunocompromised host. *Curr Clin Top Infect Dis* 1993; 13: 322-356.
- Wong SY and Remington. Toxoplasmosis in Pregnancy. *Clinical Infectious Diseases* 1994; 18:853-862.
- US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
- US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
- World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
- Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.
- Clinical and Laboratory Standards Institute (CLSI). *Statistical Quality Control for Quantitative Measurement Procedures: Principles and Definitions; Approved Guideline—Third Edition*. CLSI Document C24-A3. Wayne, PA: CLSI; 2006.
- Westgard JO. *Basic QC Practices*. 3rd ed. Madison, WI: Westgard Quality Corporation; 2010.
- Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. *Clin Chem* 1988;34(1):27-33.
- Primus FJ, Kelley EA, Hansen HJ, et al. "Sandwich"-type immunoassay of carcinoembryonic antigen in patients receiving murine monoclonal antibodies for diagnosis and therapy. *Clin Chem* 1988;34(2):261-264.
- Schroff RW, Foon KA, Beatty SM, et al. Human anti-murine immunoglobulin responses in patients receiving monoclonal antibody therapy. *Cancer Res* 1985;45(2):879-885.
- Rigsby P, Rijpkema S, Guy EC, et al. Evaluation of a candidate international standard preparation for human anti-Toxoplasma immunoglobulin G. *J Clin Microbiol* 2004;42(11):5133-5138.
- Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline—Second Edition*. CLSI Document EP05-A2. Wayne, PA: CLSI; 2004.
- Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition*. CLSI Document EP17-A2. Wayne, PA: CLSI; 2012.

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