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Instructions for Use

DESCRIPTION

RAPID RHINO™ SINU-KNIT™ Dressings are weft knitted carboxymethyl cellulose (CMC) hydro fibres. The dressing is an inverted multi layer design that allows the fabric to act as a stent. The hydrocolloid fabric creates a moist gel barrier between the device and the mucosal tissue due to its hygroscopic capability.

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CONTRAINDICATIONS

There are no known contraindications. Medical judgment is necessary when using RAPID RHINO SINU-KNIT Dressings.

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Adverse Reactions

- These include the possibility of infection, recurrent or persistent bleeding, and dislodgement of the packing necessitating medical treatment.
- Toxic shock syndrome has been reported to be associated with the use of certain nasal packing materials with fluid absorbing properties. This condition can cause death; hence all patients should be monitored for any signs, which indicate the development of this condition.
- Breathing distress is a possibility in physically small patients or patients with a shallow nasopharynx due to improper positioning of the device.

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DIRECTIONS FOR USE

1. Remove RAPID RHINO SINU-KNIT Dressing from the sterile packaging and place on a sterile surgical field.
2. Grasp the fabric with dry non-serrated forceps and gently advance the dressing into the desired location in a dry (un-hydrated state). To ensure the maximum efficacy, be sure to cover the denuded tissue and mucosa.
3. Wet the fabric with Sterile Water (not Saline). Ensure that the entire fabric from the deepest portion of the posterior cavity, to the anterior has been totally hydrated allowing the fabric conversion into its gelled state.
4. Have the patient irrigate with light saline solution until first office visit. (Usually 7-10 days post op) Do not irrigate with hypertonic solution until after first office visit.
5. Any residual gel that has not dissolved or left the surgical site through normal outflow passages may be easily removed through gentle suction.

HOW SUPPLIED

The device is supplied sterile, and is for **SINGLE USE ONLY. Do not clean, resterilize, or reuse the device, as this may damage or compromise the performance resulting in product malfunction, failure, or patient injury.** Cleaning, resterilization, or reuse of the device may also expose the patient to the risk of transmitting infectious diseases.

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The product is not made with natural rubber latex and is DEHP Free.

DISPOSAL

The device must be disposed of in accordance with applicable government regulations, and the Medical Waste Management Plan for the facility.


CUSTOMER SERVICE

Warranty Information

This product is guaranteed for materials, function, and workmanship for single patient use only. **DO NOT REUSE. THIS WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS, IMPLIED, AND/OR STATUTORY, INCLUDING BUT NOT LIMITED TO, WARRANTIES OF MERCHANTABILITY, FITNESS, AND/OR SUITABILITY FOR A PARTICULAR PURPOSE, AND OF ALL OTHER OBLIGATIONS OR LIABILITIES ON THE PART OF ARTHROCARE.**

Производитель
ArthroCare Corporation
7000 West William Cannon Drive
Austin, TX 78735 USA
Phone: (800) 343-5717
CUSTOMER SERVICE
www.arthrocare.com

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
Manufacturer
ArthroCare Corporation
7000 West William Cannon Drive
Austin, TX 78735 USA
Phone: (800) 343-5717
CUSTOMER SERVICE
www.arthrocare.com

Authorized European Representative
Smith & Nephew Orthopaedics
GmbH
Alemannenstraße 14
78532 Tuttlingen
Germany
ec.rep@smith-nephew.com

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Sterilized using ethylene oxide

Manufacturer

Authorized Representative in the
European Community

Rx only

CAUTION: Federal (USA) law
restricts this device to sale by or
on the order of a physician.

CE mark and identification
number of
Notified Body.
The product meets the essential
requirements of Medical Device
Directive (93/42/EEC).

TMTrademark of Smith & Nephew.

This product may be covered by one or more US
Patents.

See smith-nephew.com/patents for details.

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NASASTENT™

Dissolvable Nasal Dressing Instructions for Use

INDICATIONS FOR USE

NASASTENT™ dressing is intended to minimize bleeding and edema and to prevent adhesions between the septum and the nasal cavity. It is placed in the nasal cavity after surgery or trauma.

DEVICE DESCRIPTION

NASASTENT dressing is a dissolvable polysaccharide intranasal splint made from plant-based CarboxyMethyl Cellulose (CMC). As it absorbs nasal fluids, it turns into a hydrocolloidal gel that naturally drains from the nasal cavity over time.

CONTRAINDICATIONS

There are no known contraindications. Medical judgement is necessary when using NASASTENT dressing.

WARNINGS AND PRECAUTIONS

Warnings

- DO NOT attempt to use the device for anything other than the indicated use.

Precautions

- Prior to using NASASTENT dressing, ensure that all package inserts, particularly all warnings, precautions and instructions for use, are read and understood.
- NASASTENT dressing is supplied sterile and intended for SINGLE USE ONLY.
- If at any time the device becomes inflexible, do not use the device.

POSSIBLE ADVERSE EFFECTS

- Adverse effects include the possibility of infection, recurrent or persistent bleeding, and dislodgement of NASASTENT dressing necessitating medical treatment.
- Respiratory distress or decreased oxygenation can occur in patients at risk for impaired respiratory condition.
- In rare instances, the physiochemical condition associated with nasal surgery, both with and without nasal packing, may present a risk of toxic shock syndrome (TSS).
- NASASTENT dressing exhibits no antimicrobial properties. It is not bacteriostatic toward pre-existing infections, nor does it prevent the occurrence of new

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infections. In the case of pre-existing infections it is recommended that appropriate treatment be instituted.

PHYSICIAN TRAINING REQUIREMENTS

- Use of this device requires an appropriate level of medical skill and experience. Review the Instructions for Use prior to clinical use.
- Completion of a skills laboratory and/or observation of/assistance with similar procedures is recommended.
- Preoperative and operative procedures, including proper patient selection, knowledge of techniques, and proper device selection are important considerations when using the device.
- The physician should advise the patient of the known risks and complications associated with the procedure and use of this device.

INSTRUCTIONS FOR USE

1. Immediately prior to insertion, peel open the NASASTENT pouch and remove device using sterile technique to maintain product integrity.
2. NASASTENT dressing, as provided, is compressible and pliable for ease of insertion into the nasal cavity.

NOTE: NASASTENT dressing does not require hydration prior to placement.

3. Evaluate the patient's anatomy to determine if the NASASTENT dressing requires trimming prior to placement.

NOTE: NASASTENT dressing may be compressed for ease of insertion into the nasal anatomy.

4. Grasp the device with forceps and advance into the nasal cavity leading with tapered end.
5. After placement, NASASTENT dressing may be hydrated if desired.
6. Do not excessively irrigate after placement of NASASTENT dressing or use excessive force while irrigating.

NOTE: Saline nose spray may be used after placement of NASASTENT dressing.

7. Discard any unused product after opening.

Post-Operative Care

1. NASASTENT dressing is dissolvable and any residual material that has not dissolved or left the surgical site through normal outflow passages can be removed with irrigation and aspiration at the discretion of the surgeon.

2. Please advise the patient not to excessively irrigate after placement of NASASTENT dressing or to use excessive force while irrigating. (Do not irrigate with bulb style syringe until after first office visit.)

NOTE: Saline nose spray may be used as needed.

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HOW SUPPLIED

NASASTENT dressing is provided STERILE and is intended for SINGLE PATIENT USE ONLY. Do not resterilize. This product is sterilized using Ethylene Oxide.

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CUSTOMER SERVICE

Warranty Information

This product is guaranteed for materials, function and workmanship for single patient use only. DO NOT REUSE. **THIS WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS, IMPLIED, AND/OR STATUTORY, INCLUDING BUT NOT LIMITED TO, WARRANTIES OF MERCHANTABILITY, FITNESS, AND/OR SUITABILITY FOR A PARTICULAR PURPOSE, AND OF ALL OTHER OBLIGATIONS OR LIABILITIES ON THE PART OF ARTHROCARE.**

Product Complaints

All questions or concerns related for quality, reliability and/or durability of this product should be directed to Customer Service or an authorized representative.

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
Manufacturer

Rx only

Authorized Representative in the European Community
CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician. CE mark and identification number of Notified Body. The product meets the essential requirements of Medical Device Directive (93/42/EEC).

™ Trademark of Smith & Nephew.

This product may be covered by one or more US Patents. See smith-nephew.com/patents for details.

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RAPID RHINO™
SINU-FOAM™ Dressing
Instructions for Use**DESCRIPTION**

SINU-FOAM™ Dressing functions as a hydrocolloidal physical barrier that starts as a dry carboxymethylcellulose (CMC) knit within a syringe. When properly mixed, the knit gels to form a moist, viscous, and dissolvable dressing. The kit has an additional syringe for mixing and a dispensing cannula.

INDICATIONS

An intranasal splint is intended to minimize bleeding and edema and to prevent adhesions between the septum and the nasal cavity after surgery or trauma.

CONTRAINDICATIONS

There are no known contraindications. Medical judgment is necessary when using SINU-FOAM Dressing.

**WARNINGS and PRECAUTIONS****Warnings**

- Do not attempt to use the device for anything other than the indicated use.

Precautions

- Prior to mixing the SINU-FOAM Dressing, ensure that all package inserts; including warnings, precautions and instructions for use are read and understood.

- SINU-FOAM Dressing is supplied sterile and intended for **SINGLE USE ONLY**.

- Use **STERILE** Lactated Ringer Solution or Sterile Water (H₂O) **ONLY**.

- Adding excess Lactated Ringer Solution or Sterile Water (H₂O) will result in a less viscous dressing, whereas less Lactated Ringer Solution or Sterile Water (H₂O) will result in a more viscous dressing that may not perform as desired.

Adverse Reactions

- Adverse effects include the possibility of infection, recurrent or persistent bleeding and dislodgement of the dressing necessitating medical treatment.

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- Respiratory distress or decreased oxygenation can occur in patients at risk for impaired respiratory condition.

- Toxic shock syndrome has been reported to be associated with the use of certain nasal packing materials with fluid absorbing properties. This condition can cause death; hence all patients should be monitored for any signs indicating the development of this condition.

- SINU-FOAM Dressing exhibits no antimicrobial properties. It is not bacteriostatic toward pre-existing infections, nor does it prevent the occurrence of new infections. In the case of pre-existing infections it is recommended that appropriate treatment be instituted.

DIRECTIONS FOR USE

1. Using sterile techniques, remove the three pieces from the sterile packaging and place on a sterile surgical field. Do not disturb the position of the piston within the syringe containing the CMC fiber.

2. Fill a clean, empty container with STERILE Lactated Ringer Solution or Sterile Water (H_2O). Draw the sterile Lactated Ringer Solution or Sterile Water (H_2O) into the **empty small syringe** up to the fill line. Less Lactated Ringer Solution or Sterile Water (H_2O) provides a more viscous dressing, whereas more sterile Lactated Ringer Solution or Sterile Water (H_2O) will result in a less viscous dressing that may not perform as desired.

3. Remove the air from the small syringe containing the Lactated Ringer Solution or Sterile Water (H_2O).

4. Position the two syringes vertically such that the small syringe containing Lactated Ringer Solution or Sterile Water (H_2O) is above the 20 ml syringe containing the CMC knit.

5. Ensuring that the two syringes remain securely connected, transfer the Lactated Ringer Solution or Sterile Water (H_2O) from the small syringe into the syringe containing the CMC knit, then mix the resulting solution between the two syringes back and forth until an even mix is observed (about 10 – 15 transfers). Finish with the mixture in the small syringe.

6. Separate the two syringes, discard the large syringe, and connect the semi flexible dispensing cannula to the small syringe containing the prepared SINU-FOAM Dressing.

7. Instill the SINU-FOAM Dressing into the desired location.

8. To ensure maximum efficacy, cover the denuded tissue and mucosa.

9. Have the patient irrigate with light isotonic saline solution (0.65%) until first office visit (Do not irrigate with hypertonic saline solution or bulb style syringe until after first office visit).

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10. Any residual dressing that has not dissolved or left the surgical site through normal outflow passages may be easily irrigated and aspirated at the discretion of the surgeon.

HOW SUPPLIED

The device is supplied sterile, and is for **SINGLE USE ONLY. Do not clean, resterilize, or reuse the device, as this may damage or compromise the performance resulting in product malfunction, failure, or patient injury.**

Cleaning, resterilization, or reuse of the device may also expose the patient to the risk of transmitting infectious diseases.

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Product Complaints & Return Authorization

All questions, complaints, or return authorization requests should be directed to Customer Service or an authorized representative.

Manufacturer

ArthroCare Corporation
7000 West William Cannon Drive
Austin, TX 78735 USA
Phone: (800) 343-5717
CUSTOMER SERVICE
www.arthrocare.com

Authorized European Representative

Smith & Nephew Orthopaedics GmbH
Alemannenstraße 14
78532 Tuttlingen
Germany
ec.rep@smith-nephew.com

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RAPID RHINO™**RAPID PAC****Instructions for Use****DESCRIPTION**

The RAPID RHINO™ RAPID PAC™ removable nasal dressing is constructed of two layers of medical grade polyurethane foam with an open cell malleable structure that can be compressed and yet returns to its maximum memory limits. Within the dual layers of foam, is a malleable obturator which can assist with the placement. The foam is encapsulated into a welded PVC sleeve. A shroud of reinforced weft knitted hydro fiber is donned over the PVC. The fabric, which by its design wicks to form a gel, not only controls minor bleeding, but also provides a moist environment which facilitates ease of removal. The hydrocolloid fabric creates a moist gel barrier between the device and the mucosal tissue due to its hygroscopic capability.

INDICATIONS FOR USE

An intranasal splint is intended to minimize bleeding and edema and to prevent adhesions between the septum and the nasal cavity. It is placed in the nasal cavity after surgery or trauma.

CONTRAINDICATIONS

There are no known contraindications. Medical judgment is necessary when using Rapid Rhino removable nasal dressing Rapid Pac 5.5 cm in patients with severe head trauma and in patients whose facial anatomy may not accommodate the full length of the device.

**WARNINGS AND PRECAUTIONS****Warnings**

Do not attempt to use the RAPID RHINO RAPID PAC device for anything other than the indicated use.

Adverse Reactions

- These include the possibility of infection, recurrent or persistent bleeding, and dislodgement of the packing necessitating medical treatment.
- Toxic shock syndrome has been reported to be associated with the use of certain nasal packing materials with fluid absorbing properties. This condition can cause death; hence all patients should be monitored for any signs indicating the development of this condition.

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- Breathing distress is a possibility in physically small patients or patients with a shallow nasopharynx due to improper positioning of the device.

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- To ensure that migration does not occur, the Anti-Aspiration Guard (plastic butterfly) must be taped to the patient's cheek.

- To prevent secondary necrosis of the cartilaginous nasal pyramid or columella, threads should be taped to the cheek without tension to avoid necrosis.

DIRECTIONS FOR USE

1. Remove the RAPID RHINO RAPID PAC device from the packaging and place in bowl of sterile water.

2. Allow the device to soak from 30 to 45 seconds in sterile water to activate the hydrocolloid material. **DO NOT USE SALINE.**

3. Grasp and gently advance the RAPID RHINO RAPID PAC device into the anterior chamber. To ensure the maximum efficacy, be sure to completely insert the entire RAPID RHINO RAPID PAC device inside the naris.

4. **It is imperative to tape the Anti-Aspiration Guard (plastic butterfly) to the patient's cheek to avoid migration.**

5. To remove, simply un-tape the Anti-Aspiration Guard from patient's cheek, grasp the suture and gently pull out. The RAPID RHINO RAPID PAC device can usually be removed within 24 hours.

HOW SUPPLIED

The device is supplied sterile, and is for SINGLE USE ONLY. Do not clean, resterilize, or reuse the device, as this may damage or compromise the performance resulting in product malfunction, failure, or patient injury. Cleaning, resterilization, or reuse of the device may also expose the patient to the risk of transmitting infectious diseases.

The product is sterilized using ethylene oxide.

The product is not made with natural rubber latex.

DISPOSAL

The device must be disposed of in accordance with applicable government regulations, and the Medical Waste Management plan for the facility.

CUSTOMER SERVICE

Warranty information

This product is guaranteed for materials, function, and workmanship for single patient use only. DO NOT REUSE. THIS WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS, IMPLIED, AND/OR

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STATUTORY, INCLUDING BUT NOT LIMITED TO, WARRANTIES OF MERCHANTABILITY, FITNESS, AND/OR SUITABILITY FOR A PARTICULAR PURPOSE, AND OF ALL OTHER OBLIGATIONS OR LIABILITIES ON THE PART OF ARTHROCARE.

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Manufacturer

EC REP

Authorized Representative in the European Community

Rx only

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician. CE mark and identification number of Notified Body. The product meets the essential requirements of Medical Device Directive (93/42/EEC).

CE
0123

™Trademark of Smith & Nephew.
This product may be covered by one or more US Patents.
See smith-nephew.com/patents for details.
© Smith & Nephew 2015.

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RAPID RHINO™**Nasal Sinus Dressings****Instructions for Use****DESCRIPTION**

RAPID RHINO™ Nasal Sinus Dressings are constructed of medical grade polyurethane foam with an open cell malleable structure that can be compressed and yet returns to its maximum memory limits. The foam substrate acts to support and fill the nasal or sinus chambers. The foam is covered with a PVC film. A shroud of reinforced knitted carboxymethylcellulose (CMC) fabric is donned over the film. The fabric, which by its design wicks to form a gel, will not only control minor bleeding, but also provide a moist environment which facilitates ease of removal. The hydrocolloid fabric creates a moist gel barrier between the device and the mucosal tissue.

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CONTRAINDICATIONS

There are no known contraindications. Medical judgment is necessary when using RAPID RHINO Nasal Sinus Dressings.

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Do not attempt to use the device for anything other than the indicated use.

Adverse Reactions

- Infection, recurrent or persistent bleeding, and dislodgement of the packing necessitating medical treatment.
- Toxic shock syndrome has been reported to be associated with the use of certain nasal packing materials with fluid absorbing properties. This condition can cause death; hence all patients should be monitored for any signs, which indicate the development of this condition.
- Breathing distress is a possibility in physically small patients or patients with a shallow nasopharynx due to improper positioning of the device.

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DIRECTIONS FOR USE

1. Remove RAPID RHINO Nasal Sinus Dressings from the packaging and place in bowl of sterile water.
2. Allow the device to soak 30 seconds in sterile water to activate the hydrocolloid material, **DO NOT USE SALINE.**
3. Grasp with forceps and gently advance the dressing into the desired location. To ensure the maximum efficacy, be sure to cover the denuded tissue.
4. Tape the removal thread to the patient's cheek to ensure ease of removal and minimize migration.
5. To remove, simply grasp the retrieval thread and gently pull.
6. To prevent secondary necrosis of the cartilaginous nasal pyramid or columella, the removal thread should be attached to the cheek or bony nasal pyramid without tension.
7. Cat. # RR 300 and RR 400 remove within 72 hrs.
8. Cat. # RR 500 and RR 800 remove within 24 hrs.

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**Authorized European Representative**

Smith & Nephew Orthopaedics GmbH
Alemannenstraße 14
78532 Tuttlingen
Germany
ec.rep@smith-nephew.com

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Manufacturer

Rx only

Authorized Representative in
the European Community
CAUTION: Federal (USA) law
restricts this device to sale by
or on the order of a
physician.
CE mark and identification
number of
Notified Body.
The product meets the
essential
requirements of Medical
Device
Directive (93/42/EEC).

™Trademark of Smith & Nephew.

This product may be covered by one or more US
Patents. See smith-nephew.com/patents for details.

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**ADDENDUM TO INSTRUCTION FOR USE
OF MEDICAL DEVICE**

"RAPID RHINO NASAL DRESSING, VARIANTS"

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USE & OPERATION CONDITIONS OF MD

Temperature (°C)	+5.....+42
Humidity (%)	≤100
Atmospheric Pressure	Special pressure conditions are not required during operation of the product.

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MANUFACTURING FACILITIES:

1. ArthroCare Corporation, 7000 West William Cannon Drive, Austin TX 78735-8531, USA.

2. ArthroCare Corporation, B32. 1, St2, Zona Franca Coyol, Coyol, Costa Rica, 20101, COSTA RICA

AUTHORIZED REPRESENTATIVE IN RUSSIAN FEDERATION

LLC "SMITH AND NEPHEW"

105120, Moscow, Syromyatnichesky 2nd lane, 1, floor 9, prem. 1, room 1

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