State of Minnesota - Secretary of State

This Certificate is not valid for use anywhere within the United States of America, its territories or possessions.

Apostille

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1. **Country (Pays): United States of America**
This public document (Le présent acte public) **Instructions for Use**
2. **has been signed by (a été signé par) Christopher James Devet**
3. **acting in the capacity of (agissant en qualité de) Notary Public, State of Minnesota**
4. **bears the seal / stamp of (est revêtu du sceau / timbre de)**
Christopher James Devet, Notary Public, State of Minnesota

Certified Attesté

5. **at (à) St. Paul, Minnesota**
6. **Date (Date) 06/23/2021**
7. **by (par) Secretary of State, State of Minnesota**
8. **File No (Sous n°) 1240987800057**
9. **Seal / Stamp (Sceau / Timbre):**
10. **Signature (Signature):**

Steve Simon

Steve Simon
Secretary of State
State of Minnesota

This Apostille only certifies the authenticity of the signature and the capacity of the person who has signed the public document, and, where appropriate, the identity of the seal or stamp which the public document bears.

This Apostille does not certify the content of the document for which it was issued.

To verify the issuance of this Apostille, see <https://apostille.sos.state.mn.us>.

This certificate does not constitute an Apostille under the Hague Convention of 5 October 1961, when it is presented in a country which is not a party to the convention. In such cases, the certificate should be presented to the consular section of the mission representing that country.

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2510 Conway Avenue, St. Paul,
Minnesota, 55144-1000, USA

Dear Sir/Madame,

In my capacity as Regulatory Affairs Specialist, Oral Care Solutions Division, 3M ESPE
St. Paul, MN, USA, I certify that attached document (IFU) for 3M™ Filtek™ Universal
Restorative in individual packs in English is truthful and accurate

Regina Feferman
Regulatory Affairs Specialist
3M ESPE
Oral Care Solutions Division


6/16/2021

Date

State of Minnesota)
County of Ramsey) SS

Subscribed and sworn to before me this 16 day of June 2021.

Christopher James Devet, Notary Public
My Commission Expires: January 31, 2025

Instructions for use

Filtek™ Universal Restorative material in individual packs

General Information

Filtek™ Universal Restorative is a visible-light activated restorative composite optimized to create esthetic anterior and posterior restorations. The shades are body like opacity enabling up to a 2 mm depth-of-cure. The pink opaquer can be placed in 1 mm thick increments. All shades are radiopaque. Filtek Universal Restorative is offered in the following shades: A1, A2, A3, A3.5, A4, B1, B2, D3, XW, and PO. The pink opaque shade option can be used to mask discolored or stained tooth structure, metal discolorations, and amalgam stains.

The fillers are a combination of a non-agglomerated/non-aggregated 20 nm silica filler, a non-agglomerated/non-aggregated 4 to 11 nm zirconia filler, an aggregated zirconia/silica cluster filler (comprised of 20 nm silica and 4 to 11 nm zirconia particles), and a ytterbium trifluoride filler consisting of agglomerated 100 nm particles. The inorganic filler loading is about 76.5% by weight (58.4% by volume). Filtek Universal Restorative contains AUDMA, AFM, diurethane-DMA, and 1,12-dodecane-DMA. Filtek Universal Restorative is applied to the tooth following use of a methacrylate-based dental adhesive, such as manufactured by 3M ESPE, which permanently bonds the restoration to the tooth structure. Filtek Universal Restorative is packaged in traditional syringes and single-dose capsules.

Dental restorative material Filtek Universal can also be supplied in individual packs:

-Filtek Universal Syringe A1 Shade 1 - 4g, REF 6570A1, Instructions for Use, Visual technical guide; -Filtek Universal Syringe A2 Shade 1 - 4g, REF 6570A2, Instructions for Use, Visual technical guide; -Filtek Universal Syringe A3 Shade 1 - 4g, REF 6570A3, Instructions for Use, Visual technical guide; -Filtek Universal Syringe 3.5 Shade 1 - 4g, REF 6570A3.5, Instructions for Use, Visual technical guide; -Filtek Universal Syringe A4 Shade 1 - 4g, REF 6570A4, Instructions for Use, Visual technical guide; -Filtek Universal Syringe B1 Shade 1 - 4g, REF 6570B1, Instructions for Use, Visual technical guide; - Filtek Universal Syringe B2 Shade 1 - 4g, REF 6570B2, Instructions for Use, Visual technical guide; -Filtek Universal Syringe D3 Shade 1 - 4g, REF 6570D3, Instructions for Use, Visual technical guide; -Filtek Universal Syringe XW Shade 1 - 4g, REF 6570XW, Instructions for Use, Visual technical guide; -Filtek Universal Syringe PO Shade 1 - 4g, REF 6570PO, Instructions for Use, Visual technical guide; -Filtek Universal Syringe A2 Shade 1 - 1g, REF 6555A2-S, Instructions for Use, Visual technical guide; -Filtek Universal Syringe A3 Shade 1 - 1g, REF 6555A3-S, Instructions for Use, Visual technical guide; -Filtek Universal Capsules A1 Shade 20 - 0.2g, REF 6575A1, Instructions for Use, Visual technical guide; -Filtek Universal Capsules A2 Shade 20 - 0.2g, REF 6575A2, Instructions for Use, Visual technical guide; -Filtek Universal Capsules A3 Shade 20 - 0.2g, REF 6575A3, Instructions for Use, Visual technical guide; -Filtek Universal Capsules A3.5 Shade 20 - 0.2g, REF 6575A3.5, Instructions for Use, Visual technical guide; -Filtek Universal Capsules A4 Shade 20 - 0.2g, REF 6575A4, Instructions for Use, Visual technical guide; -Filtek Universal Capsules B1 Shade 20 - 0.2g, REF 6575B1, Instructions for Use, Visual technical guide; -Filtek Universal Capsules B2 Shade 20 - 0.2g, REF 6575B2, Instructions for Use, Visual technical guide; -Filtek Universal Capsules D3 Shade 20 - 0.2g, REF 6575D3, Instructions for Use, Visual technical guide; -Filtek Universal Capsules XW Shade 20 - 0.2g, REF 6575XW, Instructions for Use, Visual technical guide; -Filtek Universal Capsules PO Shade 10 - 0.2g, REF 6575PO, Instructions for Use, Visual technical guide; -Filtek Universal Capsules A2 Shade 5 - 0.2g, REF 6550A2-S, Instructions for Use, Visual technical guide; -Filtek Universal Capsules A3 Shade 5 - 0.2g, REF 6550A3-S, Instructions for Use, Visual technical guide.

Indications

- Direct anterior and posterior restorations (including occlusal surfaces)
- Core build-ups
- Splinting
- Indirect restorations including inlays, onlays, and veneers

Contraindications

Avoid use of this product in patients with known acrylate allergies

Potential Adverse Events (Potential Complications)

Possible acrylate allergies

Scope of use

Dentistry. For professional use only.

Intended Use

Used for restoration of anterior and posterior teeth.

Precautionary Information for Patients

This product contains substances that may cause an allergic reaction by skin contact in certain individuals. Avoid use of this product in patients with known acrylate allergies. If prolonged contact with oral soft tissue occurs, flush with large amounts of water. If allergic reaction occurs, seek medical attention as needed, remove the product if necessary and discontinue future use of the product.

Precautionary Information for Dental Personnel

Capsules may be warmed (Do not warm syringes).

This product contains substances that may cause an allergic reaction by skin contact in certain individuals. To reduce the risk of allergic response, minimize exposure to these materials. In particular, avoid exposure to uncured product. If skin contact occurs, wash skin with soap and water. Use of protective gloves and a no-touch technique is recommended. Acrylates may penetrate commonly used gloves. If product contacts glove, remove and discard glove, wash hands immediately with soap and water and then re-glove. If allergic reaction occurs, seek medical attention as needed.

3M SDS information can be obtained from www.3M.com or contact your local subsidiary.

Instructions for Use Preparation

1. **Prophy:** Teeth should be cleaned with pumice and water to remove surface stains and debris.
2. **Shade Selection:** Before isolating the tooth, select the appropriate shade(s) of restorative material.
3. **Isolation:** A rubber dam is the preferred method of isolation. However, other isolation systems can be used by following their appropriate instructions for use.

Direct Restorations

1. Cavity Preparation:

1.1. **Anterior restorations:** Use conventional cavity preparations for all Class III, IV, and V restorations.

1.2. **Posterior restorations:** Prepare the cavity. Line and point angles should be rounded. No residual amalgam or other base material should be left in the internal form of the preparation that would interfere with light transmission and therefore, the hardening of the restorative material.

2. **Pulp Protection:** If a pulp exposure has occurred and if the situation warrants a direct pulp capping procedure, use a minimum amount of calcium hydroxide on the exposure followed by an application of Vitrebond™ Plus Light Cure Glass Ionomer Liner/Base, manufactured by 3M ESPE. Vitrebond liner/bases may also be used to line areas of deep cavity excavation. See the Vitrebond liner/base instructions for details.

3. Placement of Matrix:

3.1. **For Anterior & Posterior restorations:** Place the matrix system of choice by following the manufacturer's instructions for use. 4. **Adhesive System:** To bond Filtek Universal Restorative to tooth structure, use of a 3M ESPE dental adhesive (for example 3M™ ESPE™ Single Bond Universal Adhesive) is recommended. Refer to adhesive system product instructions for full instructions and

precautions for the products. After curing the adhesive, continue to maintain isolation from blood, saliva, and other fluids and proceed immediately to placement of Filtek Universal Restorative.

5. Opaquer Placement (Optional): If staining or discoloration of the tooth structure has occurred, we recommend placing 3M PO. Because of the higher opacity, it is important to apply the paste in layers with a maximum thickness of 1 mm to ensure complete curing of the material. Follow the dispensing, placement, and light curing instructions below for proper use.

6. Dispensing the Composite: Follow the directions corresponding to the dispensing system chosen.

6.1. Syringe: Dispense the necessary amount of restorative material from the syringe onto the mix pad by turning the handle slowly in a clockwise manner. To prevent oozing of the restorative when dispensing is completed, turn the handle counterclockwise a half turn to stop paste flow. Immediately replace the syringe cap. If not used immediately, the dispensed material should be protected from light with an appropriate cover.

6.2. Single-Dose Capsule: Insert capsule into 3M™ ESPE™ Restorative Dispenser, manufactured for 3M ESPE. Refer to separate restorative dispenser instructions for full instructions and precautions. Extrude restorative directly into cavity.

6.3. Placement & Light Curing: Place in increments and light cure the composite as indicated in Step 8.

7. Placement:

7.1. Place and light cure restorative in increments as indicated in Step 8.

7.2. Contour and shape with appropriate composite instruments.

7.3. Avoid intense light in the working field.

7.4. Posterior placement hints:

7.4.1. To aid in adaptation, the first 1 mm layer may be placed and adapted to the proximal box.

7.4.2. A condensing instrument (or similar device) can be used to adapt the material to all of the internal cavity aspects.

8. Curing: This product is intended to be cured by exposure to a halogen or LED light with a minimum intensity of 550mW/cm² in the 400-500 nm range. Cure each increment by exposing its entire surface to a high intensity visible light source, such as a 3M ESPE curing light. Hold the light guide tip as close to the restorative as possible during light exposure.

		Cure Time	
Shades	Increment Depth	All halogen lights (with output 550-1000 mW/cm ²)	LED lights (with output 1000-2000 mW/cm ²)
Body	2.0 mm	20 sec.	10 sec.
Pink Opaquer	1.0 mm	40 sec.	20 sec.

9. Contouring: Contour restoration surfaces with fine finishing diamonds, burs or stones. Contour proximal surfaces with Sof-Lex™ Finishing Strips, manufactured for 3M ESPE.

10. Adjust Occlusion: Check occlusion with a thin articulating paper. Examine centric and lateral excursion contacts. Carefully adjust occlusion by removing material with a fine polishing diamond or stone.

11. Finishing and Polishing: Polishing with the Sof-Lex™ Diamond Polishing System, or Sof-Lex™ Finishing and Polishing System is recommended.

Indirect Procedure for Inlays, Onlays, or Veneers 1. Dental Operatory Procedure

1.1 Shade selection: Choose the appropriate shade(s) of Filtek Universal Restorative prior to isolation.

1.2 Preparation: Prepare the tooth.

1.3 Impressioning: After preparation is complete, make an impression of the prepared tooth by following the manufacturer's instructions of the impressioning material chosen. An impressioning material, such as manufactured by 3M ESPE, may be used.

1.4 Digital Scanning Systems: Alternatively, a digital scan may be taken to replace the impressioning step above. Follow the manufacturer's instructions of the scanning system chosen. A scanning system, such as one manufactured by 3M ESPE[®], may be used.

2. Laboratory Procedure

2.1 Pour the impression of the preparation with die stone. Place pins at the preparation site at this time if a "triple tray" type of impression was used.

2.2 Separate the cast from the impression after 45 to 60 minutes.

Place pins in die and base the cast as for a typical crown and bridge procedure. Mount or articulate the cast to its counter model on an adequate articulator.

2.3 If a second impression was not sent, pour a second cast using the same impression registration. This is to be used as a working cast.

2.4 Section out the preparation with a laboratory saw and trim away excess or, expose the margins so they can be easily worked. Mark the margins with a red pencil if needed. Add a spacer at this time if one is required.

2.5 Soak the die in water, then with a brush, apply a very thin coat of a separating medium to the preparation, let it dry somewhat, and then add another thin layer.

2.6 Add the first increment of composite to the floor of the preparation, stay short of the margins, and follow the cure recommendations described in the **Direct Restorations (Step 8)**.

2.7 Place and cure additional increments of composite. Allow for the last increment (incisal) to include the contact areas.

2.8 Place the die back into the articulated arch. Add the last increment of composite to the occlusal surface. Overfill very slightly mesially, distally, and occlusally. This will allow for the mesiodistal contacts and the proper occlusal contact when the opposing arch is brought into occlusion with the uncured increment. Light cure for only ten seconds, then remove the die to prevent adhering to adjacent surfaces. Finish the curing process following the cure times in the **Direct Restorations (Step 8)**.

2.9 With the occlusal contacts already established, begin removing the excess composite from around the points of contact. Develop the inclines and ridges as per remaining occlusal anatomy.

2.10 Care must be taken when removing the prosthesis from the die. Break off small amounts of the die from around the restoration, the die stone should break away cleanly from the cured restoration, until all of the restoration is recovered.

2.11 Using the master die, check the restoration for flash, undercuts, and fit. Adjust as necessary, and then polish as noted above in **Direct Restorations (Steps 9 – 11)**.

3. Dental Operatory Procedure

3.1 Roughen the interior surfaces of the indirect restoration.

3.2 Clean the prosthesis in a soap solution in an ultrasonic bath and rinse thoroughly.

3.3 Cementation: Cement the prosthesis using a 3M ESPE resin cement system, manufactured by 3M ESPE following manufacturer's instructions.

Cleaning & Disinfection

The multiple-use syringe dispenser is not intended for direct patient contact. Use new, uncontaminated gloves when handling the syringe dispenser. Directions for cleaning and low level disinfection of the syringe dispenser are provided below:

Step 1 (Cleaning):

Use a CaviWipes[™], or equivalent cleaning wipe, and wipe the entire surface of the device thoroughly for at least 30 seconds and until no visible soil remains on the device.

Step 2 (Disinfecting):

Use a new CaviWipes, or equivalent alcohol-quaternary ammonium disinfectant wipe, to disinfect the entire surface of the device by keeping wet for the contact time listed on the disinfectant label.

Storage and Use

1. This product is designed to be used at room temperature. If desired, the product may be warmed in a commercial warmer prior to use (no higher than 70°C/158°F, no longer than 1 hour); for capsules only.

2. The product is best stored at room temperature. If stored in cooler allow product to reach room temperature prior to use. Shelf life at room temperature is 36 months. Ambient temperatures routinely higher than 27°C/80°F may reduce shelf life. See outer package for expiration date.

3. Do not expose restorative materials to intense light.

4. Do not store materials in proximity to eugenol containing products.

Temperature °C	+2°C - +27°C
Humidity, (%)	до 80%
Pressure (ГПа)	840 - 1065

Specifications: syringe Filtek Universal: length (mm): $120,5 \pm 5$ %, width (mm): $26,5 \pm 5$ %; diameter (mmmm): cylinder diameter: $11,2 \pm 5$ %, weight (g): 12 ± 5 %. dead space 0,02-0,15 g
Filtek Universal cap: length (mm): $11,3 \pm 5$ %, width (mm): $14,7 \pm 5$ % (cap only).

Applicable standards: ISO 10993-1:2009, ISO 10993-3:2014, ISO 10993-5:2009, ISO 10993-6:2016, ISO 10993-10:2010, ISO 10993-11:2017, ISO 10993-12:2012, ISO 10993-18:2005, ISO/TR 10993-22, ISO 7405:2008, EN ISO 14971:2012, ISO 13485:2003 + Kop. 1:2009, EN 1041:2008+A1:2013, ISO 15223-1:2016, IEC 62366:2015/Kop. 1:2016, ASTM D4169:2016, ISO 4049:2009.

ISO 4049 classification:

Type 1, class 2; group 1 +2

Transportation conditions

Transportation conditions: Shelf life of the products is guaranteed as labeled, if stored according to the recommended storage conditions. During transportation the products may experience short term temperature variations without affecting product safety and performance.

Transportation of medical device is carried out by all type of transport.

Medicinal substances

Medical devices does not contain or deliver medicinal substances.

Human Blood Derivatives

Medical devices does not contain or deliver human blood derivatives.

Animal-Derived Substances

No materials of animal origin are intentionally being added to the products.

Disposal

See the Safety Data Sheet (available through your local subsidiary for disposal information).

Disposal methods: Dispose of content/container in accordance with the local/regional/national/international/regulations.

Medical devices in Russian Federation must be disposed of in the local public entities in accordance with SanPiN 2.1.7.2790-10 "Sanitary-epidemiological

Customer Information

No person is authorized to provide any information which deviates from the information provided in this instruction sheet.

Caution: U.S. Federal Law restricts this device to sale or use on the order of a dental professional.

Information about Medical Device Manufacturer

- *Developer*

3M ESPE Dental Products

Address: 2510 Conway Avenue, St. Paul, Minnesota, 55144-1000, USA

Phone: 651-7334365

- *Manufacturer*

3M ESPE Dental Products

Address: 2510 Conway Avenue, St. Paul, Minnesota, 55144-1000, USA

Phone: 651-7334365

– *Addresses of Manufacturing Sites*

3M ESPE Dental Products, 2111 McGaw Avenue, Irvine, California 92614, USA

Warranty

3M warrants this product will be free from defects in material and manufacture. 3M MAKES NO OTHER WARRANTIES INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and 3M's sole obligation shall be repair or replacement of the 3M product.

Limitation of Liability

Except where prohibited by law, 3M will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

Characteristics

Properties	3M Filtek Universal Restorative
Filler particle size distribution	The fillers are a combination of a non-agglomerated/non- aggregated 20 nm silica filler, a non-agglomerated/non- aggregated 4 to 11 nm zirconia filler, an aggregated zirconia/silica cluster filler (comprised of 20 nm silica and 4 to 11 nm zirconia particles), and a ytterbium trifluoride filler consisting of agglomerated 100 nm particles. The inorganic filler loading is about 76.5% by weight (58.4% by volume).
Wavelength (nm) for curing	400nm to 500nm
Intensity (mW/cm ²) for curing	Instructions provided for: • 550 to 1000 mW/cm² lights • 1000 to 2000 mW/cm² lights
Curing time recommendations (sec)	Light-cure each increment of 2 mm depth of cure (body shades) • 550 to 1000 mW/cm²: 20 sec • 1000 to 2000 mW/cm²: 10 sec Light-cure each increment of 1 mm depth of cure (Pink Opaque shade) • 550 to 1000 mW/cm²: 40 sec • 1000 to 2000 mW/cm²: 20 sec
Depth of cure recommendation (mm) .	2mm increment for body shades 1mm increment for opaque shade (see curing time recommendations above)

Labels and packaging symbols

	Use by date		Do not reuse
	Catalogue number		Lot number
	Caution	CONTENTS	Contained items

	Recycling code (Low-density polyethylene)		Product meets the requirements of the applicable EC directives
	Manufacturer		Temperature limit (storage conditions)
	Recycling code (Non-corrugated fiberboard)		Green dot

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