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IMAGE1 S 4U RUBINA OPAL1 NIRICG camera head for endoscopy

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IMAGE1 S™ 4U Rubina, OPAL1® NIR/ICG Kamerakopf TH121


INSTRUCTION MANUAL

IMAGE1 S™ 4U Rubina, OPAL1® NIR/ICG camera head TH121

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Thank you for placing your trust in the name KARL STORZ. Like all of our other products, this product is the result of years of experience and great care in manufacture. You and your organization have decided in favor of a modern, high quality piece of equipment from KARL STORZ.

This instruction manual, in conjunction with the instruction manual for the IMAGE1 S™ (96206286DER), is intended to serve as an aid in the proper connection and operation of the camera head. All required details and all actions on your part are clearly presented and explained. Please read this instruction manual carefully before proceeding with the equipment. Keep this instruction manual available for ready reference in a convenient location near the equipment.

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It is recommended to check the suitability of the product for the intended procedure prior to use.

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1 Safety instructions	3
1.1 Explanation of warnings and cautions	3
1.2 General safety information	4
2 Intended Use Camera Heads	6
2.1 Contraindications	6
2.2 User qualification	6
3 Images of the equipment	7
4 Symbols employed	8
5 Operating instructions for the IMAGE1 S™ 4U Rubina, OPAL1® NIR/ICG camera head	9
5.1 Frequency of the IMAGE1 S™ 4U Rubina, OPAL1® NIR/ICG camera head (50/60 Hz)	9
5.2 Operating the IMAGE1 S™ 4U Rubina, OPAL1® NIR/ICG camera head	9
5.2.1 Connecting an endoscope and light cable	9
5.2.2 Adjusting the focus	10
6 Technical description	11
6.1 Specifications	11
6.2 Storage and Operating conditions	11
6.3 Standard compliance	11
6.4 Directive compliance	11
7 Reprocessing (not validated for US customers)	12
7.1 Care and handling	12
7.2 Reprocessing	13
7.2.1 Reference chart	13
7.2.2 General warnings	14
7.2.3 Accessories	15
7.2.4 Manual wipe-down disinfection	15
7.2.5 Preparation for cleaning and disinfection	15
7.2.6 Manual cleaning	15
7.2.7 Manual disinfection	16
7.2.8 Assembly, inspection and care	16
7.2.9 Packaging systems	16
7.2.10 Sterilization	17
7.2.11 Reprocessing limits	18

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8 Reprocessing (validated for US customers)	19
8.1 Care and handling	19
8.2 Reprocessing	20
8.2.1 Reference chart	20
8.2.2 General warnings	21
8.2.3 Accessories	21
8.2.4 Bedside Pre-cleaning	22
8.2.5 Manual cleaning	22
8.2.6 Assembly, inspection and care	24
8.2.7 Packaging systems	24
8.2.8 Sterilization	24
8.2.9 High-level disinfection	27
8.2.10 Reprocessing limits	27
9 Maintenance	28
9.1 Maintenance	28
9.2 Servicing and repair	28
9.3 Disposal	28
9.4 Repairing devices	29
9.5 Infection prevention	29
9.6 Limitation of liability	30
9.7 Warranty	30
10 Compatible camera control unit (CCU)	31
11 Spare parts, recommended accessories	32
12 Information on electromagnetic compatibility (EMC)	33
13 Subsidiaries	49

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1 Safety instructions

Safety instructions are measures intended to protect the user and patients from the risks which could arise through the use of the system.

1.1 Explanation of warnings and cautions

Please read this manual and follow the instructions carefully. The words **Warning**, **Caution** and **Note** convey special meanings. Wherever they are used in this manual, the accompanying text should be carefully reviewed to ensure the safe and effective operation of the system. To make these words stand out more clearly, they are accompanied by a pictogram.

WARNING: A Warning indicates that the personal safety of the patient or physician is at risk.

CAUTION: A Caution indicates that particular service procedures or precautions must be followed to avoid possible damage to the device.

NOTE: A Note indicates special information about operating the device or clarifies important information.

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1.2 General safety information

-  **WARNING:** Read this instruction manual thoroughly before using the equipment.
-  **WARNING:** Check that the equipment is functioning properly prior to each surgical intervention. Ensure that the correct video image appears on all video monitors before beginning each intervention. It is recommended that a replacement system is kept available.
-  **WARNING:** Perceived electric shocks may be experienced when using the camera system in conjunction with electrosurgical units.
-  **WARNING:** When using high frequency surgical accessories, care must be taken to avert contact between the HF surgical device and the endoscope or endoscopic camera. Inadvertent contact may represent a thermal hazard to the surgeon.
-  **WARNING:** When energized endoscopes are used with energized endotherapy devices, patient leakage currents may be additive. It should be noted that this is particularly important if a Type CF applied part endoscope is used, in which case a Type CF applied part energized endotherapy device should be used in order to minimize total patient leakage current.
-  **WARNING:** Always ensure the connector is completely clean and free of debris. If moisture is present, dry thoroughly with sterile towel prior to insertion.

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-  **WARNING:** HF-generating equipment is dangerous and cannot be allowed to make contact with the camera head or system via the accessories or endoscopes connected to the camera head or system.
-  **WARNING:** Risk of blindness. Never look into the open end of a light cable.
-  **WARNING:** Risk of burns. Contact with the open end of a light cable connected to a light source may cause burns.
-  **WARNING:** Test all equipment prior to each surgical use. In the event that the image becomes unusable during surgery, the camera may be disengaged from the endoscope and the procedure continued optically. If this is not possible, it is up to the discretion of the surgeon how best to proceed. Availability of a spare system is recommended.
-  **WARNING:** High energy radiated light through endoscopes may give rise to high temperatures in front of the light outlet and to the tip of the endoscope. To minimize the risk of burns, always adjust the light source to the minimum illumination intensity necessary to achieve optimum illumination of the endoscopic scene when coupled to the video camera.

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2 Intended Use Camera Heads

KARL STORZ camera heads, in combination with the corresponding camera control unit, light source, monitor and telescopes or microscopes, are used for real-time display of the endoscopic or open operating field in diagnostic or surgical procedures.

Additionally, they are used for detection during fluorescence imaging in the visible and near-infrared spectral range.

KARL STORZ camera heads have no direct contact with the patient.

2.1 Contraindications

No contraindications relating directly to the medical device are currently known. The responsible physician must decide whether the foreseen application is admissible based on the general condition of the patient.

2.2 User qualification

The IMAGE1 S™ 4U Rubina, OPAL1® NIR/ICG camera head may only be used by physicians and medical assistants who have a corresponding specialized qualification and who have been instructed in the use of the IMAGE1 S™ 4U Rubina, OPAL1® NIR/ICG camera head.

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



- ① Instrument coupler
- ② Focusing ring/knob
- ③ User control button: for scrolling (downwards) through the menu or activating a camera function
- ④ User control button: for scrolling (upwards) through the menu or activating a camera function
- ⑤ User control button: for accessing the menu/ making a selection

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

4 Symbols employed

	Fragile, handle with care
	Keep dry
	Temperature limit
	Serial number
	Number of products in the product packaging
	CE marking
	top
	Consult instructions for use
	Manufacturer
	Federal (USA) law restricts this device to sale by or on the order of a physician.

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**5 Operating instructions
for the IMAGE1 S™ 4U
Rubina, OPAL1® NIR/ICG
camera head**

i NOTE: While operating with surgical lasers intensive laser light may occur and affect image visualization.

**5.1 Frequency of the IMAGE1 S™
4U Rubina, OPAL1® NIR/ICG
camera head (50/60 Hz)**

The IMAGE1 S™ 4U Rubina, OPAL1® NIR/ICG camera head is suitable for operation at either 50 Hz or 60 Hz. The camera head must be selected to suit the local conditions.

i NOTE: The method for switching frequency is described in the instruction manual IMAGE1 S™ – TC200, TC201, TC300, TC301, TC304 (96206286DER).

**5.2 Operating the IMAGE1 S™
4U Rubina, OPAL1® NIR/ICG
camera head**

**5.2.1 Connecting an endoscope
and light cable**

The IMAGE1 S™ 4U Rubina, OPAL1® NIR/ICG camera head has an integrated instrument coupler and integrated telescope. To connect an endoscope to the instrument coupler, turn the outer ring of the coupler clockwise and insert the endoscope eyepiece. To prevent rotation of the endoscope, the outer ring of the instrument coupler must be turned counterclockwise until tight. Connect the light cable to the endoscope eyepiece by tightening the knurled screw on the light cable through a quarter turn.

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CAUTION: Heat from the light source can cause disinfectant residue on the light cable or on the light inlet of the endoscope to burn. This can have a detrimental effect on the light conduction and thus the image quality.


NOTE: For further information regarding light cables, please refer to the manual 'Fiberoptic light cables' (97000213), which can be downloaded or requested at www.karlstorz.com.

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5.2.2 Adjusting the focus

The focus on the IMAGE1 S™ 4U Rubina, OPAL1® NIR/ICG camera head can be adjusted by turning the focusing ring/knob ② on the camera lens.

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IMAGE1 S™ 4U Rubina, OPAL1® NIR/ICG camera head	
Image sensor	1/2.5"
Image format	16:9
Image refresh rate	50 Hz/60 Hz
Focal length	19 mm
Dimensions (L x H x W)	150 mm x 55 mm x 41 mm
Weight	260 g

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6.3 Standard compliance

According to IEC 60601-1; IEC 60601-2-18.

Provides Type CF protection when used with KARL STORZ Camera Control Units that have CF symbol on connection port.

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6.4 Directive compliance

According to the Medical Device Directive (MDD):

This medical device belongs to Class I

This medical device bears a CE mark in accordance with MDD 93/42/EEC.

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**7 Reprocessing
(not validated for
US customers)**

In the following chapter, you will find information on care and cleaning of the IMAGE1 S™ 4U Rubina, OPAL1® NIR/ICG camera head.

7.1 Care and handling

To maximize the life of your camera, please follow the guidelines given below for the correct handling and care of the camera head:

DO NOT PLUG THE CAMERA HEAD CABLE INTO THE PROCESSOR RECEPTACLE IF IT IS MOIST OR NOT CLEAN. Always ensure that the connector is completely clean. Remove any moisture present with a sterile cloth before connecting.

DO NOT PICK UP THE CAMERA HEAD BY ITS CABLE. Always pick up the camera head first, and then pick up and support the cable and connector. The cable contains many tiny conductors and must therefore be handled with care.

DO NOT STORE THE CAMERA SYSTEM in direct sunlight or at high temperatures.

ENSURE THAT THE CABLE IS NOT COILED TOO TIGHTLY OR FOLDED UP WHILE THE CAMERA IS IN STORAGE. Always arrange the cable in loose coils with a diameter of 6 inches (15 cm) or more.

DO NOT PULL FORCEFULLY ON THE CABLE WHEN UNCOILING. Undesirable 'kinks' may result. Always allow the cable to uncoil gently and loosely.

HANDLE THE SYSTEM WITH CARE, avoiding any strong impacts. When transporting, use the original packaging to prevent damage to the device.

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**Reprocessing
(not validated for
US customers)**

7.2.2 General warnings

WARNING: Risk of infection: Incorrectly reprocessed medical devices expose patients, users and third parties to a risk of infection as well as the risk that the medical device may malfunction. Observe the 'Cleaning, Disinfection, Care, and Sterilization of KARL STORZ Instruments' instructions and the accompanying documentation.

WARNING: Risk of infection: These medical devices are not sterile when delivered. The use of non-sterile medical devices poses a risk of infection for patients, users and third parties. Inspect medical devices for visible contamination. Visible contamination is an indication that reprocessing has not been carried out, or has been carried out incorrectly. Reprocess the medical devices before initial use and before and after every subsequent use using validated procedures.

WARNING: When carrying out any work on contaminated medical devices, the guidelines of the Employers' Liability Insurance Association and equivalent organizations striving to ensure personal safety must be observed.

CAUTION: When preparing and using the solutions, follow the chemical manufacturer's specifications, paying close attention to proper concentration, exposure time and service life. Prolonged immersion and incorrect concentration may result in damage. Bear in mind the microbiological range of action of the chemicals used.

CAUTION: Damage to product. Only perform reprocessing using chemicals approved by KARL STORZ. A list of approved chemicals can be requested from hygiene@karlstorz.com.

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Reprocessing (not validated for US customers)

! **CAUTION:** National laws and regulations must be observed.

i **NOTE:** The manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' (art. no. 96216003) can be downloaded or requested by visiting www.karlstorz.com.

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- Brush: 27652

7.2.4 Manual wipe-down disinfection

The exterior surfaces of the medical device must be wiped clean with a disposable cloth moistened with disinfectant or a ready to use soaked disinfectant cloth. Due to their protein-fixating effect and material incompatibility, alcohol-based agents must not be used. The chemical manufacturer's specifications regarding material compatibility must be observed. At the end of the necessary exposure time, wipe the surface with a dry low lint cloth.

7.2.5 Preparation for cleaning and disinfection

Heavy soiling, corrosive solutions and pharmaceuticals must be removed from the medical device immediately after use. To this end, preclean the medical device by wiping down and rinsing, for example. As a general rule, KARL STORZ recommends manual precleaning under cold running water.

7.2.6 Manual cleaning

The medical device must be completely immersed in a cleaning solution. To ensure bubble-free wetting, the lumina in particular must be filled. After the necessary exposure time, clean the instrument mechanically with the aid of brushes or sponges. Finally, it must be rinsed with cold water to ensure neutralization.

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7.2.7 Manual disinfection

The medical device must be completely immersed in a disinfectant solution. To ensure bubble-free wetting, the lumina in particular must be filled. At the end of the necessary exposure time, the medical device must be rinsed several times in order to remove all chemical residues. To this end, water of the best possible quality must be used, taking into account the national regulations. Finally, all of the surfaces, joints, openings, channels and lumina are dried completely with (preferably sterile) compressed air in accordance with the national regulations. The cleaning gun with accessories (Art. no. 27660) is ideal for this purpose.

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i NOTE: During care procedures, use items from the catalog 'Care, Sterilization and Storage Techniques'.

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7.2.10 Sterilization

The procedures as well as the process-relevant parameters for the individually validated methods are described in detail in the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' (Art. no. 96216003). The method must be selected taking into account the respective applicable national requirements and in consultation with the device and product manufacturers.

The following sterilization methods have been validated and approved by KARL STORZ for this medical device:

**Hydrogen peroxide (H₂O₂)-
sterilization – ASP STERRAD®**

WARNING: Medical devices must be sterilized in dismantled condition.

WARNING: Remember that there are restrictions concerning the instruments which can be sterilized in the different STERRAD® sterilization systems in respect of lumen dimensions and material.

CAUTION: For detailed information, please consult the user handbook of the respective device.

NOTE: Refer to the 'STERRAD® Sterility Guide' to ensure that the medical device concerned can be sterilized using the different STERRAD® devices.

NOTE: Sterilization is not possible on surfaces which have been lubricated and oiled.

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8 Reprocessing (validated for US customers)

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 DO NOT PLUG THE CAMERA HEAD CABLE INTO THE PROCESSOR RECEPTACLE IF IT IS MOIST OR NOT CLEAN. Always ensure that the connector is completely clean. Remove any moisture present with a sterile cloth before connecting.

 DO NOT PICK UP THE CAMERA HEAD BY ITS CABLE. Always pick up the camera head first, and then pick up and support the cable and connector. The cable contains many tiny conductors and must therefore be handled with care.

 DO NOT STORE THE CAMERA SYSTEM in direct sunlight or at high temperatures.

 ENSURE THAT THE CABLE IS NOT COILED TOO TIGHTLY OR FOLDED UP WHILE THE CAMERA IS IN STORAGE. Always arrange the cable in loose coils with a diameter of 6 inches (15 cm) or more.

 DO NOT PULL FORCEFULLY ON THE CABLE WHEN UNCOILING. Undesirable 'kinks' may result. Always allow the cable to uncoil gently and loosely.

 HANDLE THE SYSTEM WITH CARE, avoiding any strong impacts. When transporting, use the original packaging to prevent damage to the device.

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- Bürsten: 27652

8.2.2 General warnings

WARNING: Risk of infection: Incorrectly reprocessed medical devices expose patients, users and third parties to a risk of infection as well as the risk that the medical device may malfunction. Observe the "Cleaning, Disinfection, Care, and Sterilization of KARL STORZ Instruments" instructions and the accompanying documentation.

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CAUTION: National laws and regulations must be observed.

8.2.3 Accessories

Accessories required for carrying out reprocessing:

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8.2.4 Bedside Pre-cleaning

1. Immediately after a procedure, at the point of use, wipe the camera head with a soft, lint-free, disposable cloth moistened with water or enzymatic cleaning solution to remove all heavy soiling.
2. Contain and transport the camera head to the decontamination area for cleaning.
3. If immediate cleaning is not feasible, place the camera head in a container and soak in neutral/mild pH enzymatic cleaning solution (e.g. Enzol), diluted with distilled/demineralized water to proper concentration per manufacturer's instructions for no longer than an hour.

i NOTE: Complete cleaning of the camera head should be started within two hours of the bedside pre-cleaning.

8.2.5 Manual cleaning

1. Prepare neutral pH enzymatic solution (e.g. Enzol) in a container large enough to fully submerge the camera head. Dilute the solution to proper concentration per manufacturer's instructions.
2. Under running cold tap water:
 - Rinse the camera head to remove all gross debris, paying close attention to the seam beneath the focus ring, for 2 minutes.
 - Brush all camera head surfaces using the soft bristle brush (P/N 27652), paying close attention to the seam beneath the focus ring. Brush for 2 minutes or until there is no visible debris on the brush.
- i NOTE:** If heavy debris collects on the brush while brushing, rinse the brush until there is no visible debris on the brush.
3. Completely immerse the camera head in a neutral pH enzymatic solution. Keep immersed for 5 minutes.

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CAUTION: Do not brush the glass window of the camera head. Doing so may cause damage to the glass. Use a soft, lint-free cloth.

4. While immersed, use the soft bristle brush (P/N 27652) to:

- Brush all camera head surfaces, paying close attention to the seam beneath the focus ring. Brush for 1 minute or until there is no visible debris on the brush.

NOTE: If heavy debris collects on the brush while brushing, rinse the brush until there is no visible debris on the brush.

5. While immersed, use the soft bristle brush (P/N 27652) to brush the following locations for 30 seconds each:

- Entire silicone cable
- Card edge connector
- Camera head buttons
- Cam-Lock mechanism, including its crevices within the coupling mechanism

NOTE: If heavy debris collects on the brush while brushing, rinse the brush until there is no visible debris on the brush.

6. While immersed, use a soft, lint-free cloth to wipe the entire camera head cable.

7. Remove the camera head from the cleaning solution.

8. Under running cold tap water, rinse the entire camera head and cable for 3 minutes including:

- Focus ring
- Camera head buttons
- Cam-Lock coupling mechanism
- Crevices and seams
- Entire cable

9. Repeat rinsing (step 8) for a total of 3 times.

10. Dry the camera head with a soft, lint-free, disposable cloth or compressed air.

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**Reprocessing
(validated for
US customers)**

8.2.6 Assembly, inspection and care

The cleaned and disinfected medical device must be visually inspected for cleanliness, completeness, damage and dryness:

- If residues or contamination are still present, the medical device must be manually cleaned and subjected to a full cleaning and disinfection procedure once more.
- Damaged or corroded medical devices must be withdrawn from use.
- Afterwards, a functional check must be carried out.

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Hydrogen peroxide (H₂O₂)-
sterilization – ASP STERRAD®

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Reprocessing (validated for US customers)

The following STERRAD® sterilization methods have been validated and approved by KARL STORZ for this medical device:

- STERRAD® 100S
- STERRAD® NX® standard cycle
- STERRAD® 100NX® standard cycle

1. Place the camera head in an FDA-cleared STERRAD® compatible instrument tray.
2. Wrap the sterilization tray with two single layers of FDA-cleared polypropylene sterilization wrap, using a sequential double wrapping technique. If using a sterilization sealed container, refer to the manufacturer's instructions for use on using the correct disposable filters for the indicated sterilization method.
3. Load the STERRAD® sterilizer, arranging the items such that the hydrogen peroxide plasma can surround them. Do not allow any items to touch the wall of the sterilizer.
4. Start the STERRAD® sterilizer.

Hydrogen peroxide (H₂O₂) sterilization – STERIS® V-PRO®

The following STERIS® V-PRO® sterilization methods have been validated and approved by KARL STORZ for this medical device:

- V-PRO® 1 „Lumen“ Cycle
- V-PRO® 1 Plus „Lumen“ and „Non-lumen“ Cycle
- V-PRO® maX „Lumen“ and „Non-lumen“ Cycle

1. Place the camera head in an FDA-cleared V-PRO® compatible instrument tray.
2. Wrap the sterilization tray with two single layers of FDA-cleared polypropylene sterilization wrap, using a sequential double wrapping technique. If using a sterilization sealed container, refer to the manufacturer's instructions for use on using the correct disposable filters for the indicated sterilization method.
3. Load the V-PRO® sterilizer, arranging the items such that the vaporized hydrogen peroxide can surround them. Do not allow any items to touch the wall of the sterilizer.
4. Start the V-PRO® sterilizer.

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Ethylene oxide (EO) sterilization

1. Place the camera head in an FDA-cleared sterilization tray.
2. Double wrap the sterilization tray with polypropylene wrap. If using a sterilization sealed container refer to the manufacturer's instructions for use on using the correct disposable filters for the indicated sterilization method.
3. Start sterilization under the following parameters:

Conditioning parameters

Temperature: 55 °C (131 °F)
Humidity: 70 % RH
Conditioning dwell time: 30 minutes

Sterilization parameters

Sterilant: 100 % ethylene oxide
Temperature: 55 °C (131 °F)
Gas concentration: 735±30 mg/l
Exposure time: 180 minutes

Aeration parameters

Time: 12 hours, minimum
Temperature: 55 °C (131 °F)

Chemical low-temperature sterilization with peracetic acid – STERIS® System 1E®

1. Place the camera head in one of the following trays:
 - C1200E SS1E General Processing Tray/Container
 - C1220E SS1E Directed Flow Processing Tray/Container
 - C1140E SS1E Flexible Endoscope Processing Tray/Container
 - C1160E SS1E Universal Flexible Processing (open) Tray
2. Load the tray into the SS1E and start the sterilizer.

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8.2.9 High-level disinfection

The following high-level disinfection methods have been validated and approved by KARL STORZ:

Revital-Ox™ RESERT® (e.g. 2.0% Accelerated Hydrogen Peroxide solution)

In order to prevent potential damage, do not place camera head in a container with other instruments.

Completely immerse the camera head in the undiluted 2.0% hydrogen peroxide solution for a minimum of 8 minutes at 20 °C (68 °F). It must be ensured that the surfaces are covered and that no air bubbles are present.

To finish, immerse the camera head three times in water. Each immersion rinse should be for a minimum of 1 minute in duration. Discard the water after each rinse, as it will be contaminated with the disinfectant solution. Use fresh water for each rinse. Thorough rinsing of the camera head with water is very important for removing any residual disinfectant.

Dry the camera head and the card-edge connector with a lint-free, sterile cloth. Ensure that the card-edge connector is completely clean and free of moisture before inserting it into the camera processor. Wipe the camera head connector with isopropyl alcohol to remove any residual disinfectant.

8.2.10 Reprocessing limits

The end of the product's service life is largely determined by wear, reprocessing methods, the chemicals used and any damage resulting from use.

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9 Maintenance

9.1 Maintenance

Preventive maintenance is not essential. Regular maintenance can, however, contribute to identifying potential problems before they become serious, thus enhancing the instrument's reliability and extending its useful service life. Maintenance services can be requested from your local representative or from the manufacturer.

9.2 Servicing and repair

Defective items of equipment must be serviced and repaired exclusively by KARL STORZ or by persons authorized by KARL STORZ; all repair work must employ original KARL STORZ parts only.

CAUTION: Do not attempt to open the camera head or lens. Opening these sealed assemblies affects the camera's soakability and voids all warranties. Do not remove the covers on the camera processor. Warranty is voided unless removed by authorized personnel.

9.3 Disposal

This unit has been marked in accordance with the European Directive on waste electrical and electronic equipment (WEEE).

At the end of its useful operating life, dispose of the unit as electronic scrap.

Please ask either KARL STORZ SE & Co. KG, a KARL STORZ subsidiary or your specialist dealer for information on your local collection point. Within the scope of application of this Directive, KARL STORZ SE & Co. KG is responsible for the proper disposal of this unit.

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9.4 Repairing devices

Only personnel authorized by the manufacturer may repair devices. Interventions described in these instructions for use are exempt from this rule.

In Germany, you can contact the following address for repairs:

KARL STORZ SE & Co. KG
Abt. Reparaturservice
Take-Off Gewerbepark 83
78579 Neuhausen

Service hotline: +49 7461 708-980
E-mail: technicalsupport@karlstorz.com

In other countries, the respective KARL STORZ branches or specialist dealers are responsible.

9.5 Infection prevention

Contaminated devices may not be shipped. To prevent contact infections and airborne infections, devices must first be decontaminated. KARL STORZ reserves the right to send back contaminated devices.

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9.6 Limitation of liability

KARL STORZ shall be liable for the safe operation, operational reliability, and performance of this equipment only subject to the following conditions:

- all assembly operations, system expansions, readjustments, modifications, or repairs have been performed by a person or persons duly authorized by KARL STORZ;
- all electrical installations at the location of use meet applicable national and local electrical codes, and
- the unit has been used in accordance with its operating instructions at all times.

9.7 Warranty

The guarantees provided can be found in the Standard Conditions of Business of KARL STORZ.

The medical device must always be sent to your local subsidiary (see "Subsidiaries" section), even during the warranty period.

Opening the equipment or performance of any repairs or modifications to the equipment by unauthorized persons shall relieve us of any liability for its performance. Any such opening, repair, or modification performed during the warranty period shall void all warranty.

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**10 Compatible
camera control unit (CCU)**

Item	Cat. No.
IMAGE1 S™ 4U-LINK in combination with IMAGE1 S CONNECT® II	TC304 TC201*

* IMAGE1 camera control units are available with various configurations and image outputs. These camera control units are compatible with IMAGE1 S™ camera heads. Please check your local region for device configuration availability.

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**11 Spare parts,
recommended accessories**

Please check your local region for recommended accessory availability.

Article	Order no.
Adaptor, autoclavable, allows telescopes to be changed under sterile conditions	533TVA
Integral rotating beamsplitter, for simultaneous viewing through the telescope and on the monitor	9530BD
Dust cap for camera heads	6349190
Sterile drapes* for KARL STORZ cameras	
Instruction Manual	96286011DER
Instruction manual IMAGE1 S™	96206286DER

* For safe and comfortable use of KARL STORZ cameras, we recommend the use of sterile drapes.

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Information on electromagnetic compatibility (EMC)

12 Information on electromagnetic compatibility (EMC)

The IMAGE1 S™ 4U Rubina, OPAL1® NIR/ICG camera head TH121 used with the IMAGE1 S™ is suitable for use within close proximity to an active HF electro-surgical device in professional healthcare facility environment. Professional healthcare facility includes physician offices, dental offices, limited care facilities, freestanding surgical centers, freestanding birth centers, multiple treatment facilities, hospitals (emergency rooms, patient rooms, intensive care, surgical rooms, outside the RF shielded room of an ME system for magnetic resonance imaging).

WARNING: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the IMAGE1 S™ 4U Rubina, OPAL1® NIR/ICG camera head, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

NOTE: The emission characteristics of this equipment make it suitable for use in industrial areas as well as in hospitals (CISPR 11 Class A) and other professional healthcare environment. If it is used in a residential environment (for which CISPR 11 Class B is normally required) this equipment might not offer adequate protection to radio communication service. The user might need to take mitigation measures, such as relocating or re-orienting equipment.

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Table 1 – Compliance level for Immunity tests

Guidance and manufacturer's declaration – electromagnetic immunity

The IMAGE1 S™ 4U Rubina, OPAL1® NIR/ICG camera head is intended for use in the electromagnetic environment specified below. The user of the IMAGE1 S™ 4U Rubina, OPAL1® NIR/ICG camera head should make sure that it is used in such an environment.

Immunity tests	EN/IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) acc. to IEC 61000-4-2	± 8 kV Contact discharge ± 15 kV Air discharge	± 8 kV Contact discharge ± 15 kV Air discharge	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Rapid transient electrical Interferences/bursts acc. to IEC 61000-4-4	± 2 kV for power lines ± 1 kV for input and output lines 100 kHz repetition	± 2 kV/1 kV ± 1 kV for input and output lines 100 kHz repetition	Mains power quality should be that of a typical commercial or hospital environment.
Surges acc. to IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and fluctuations of the supply voltage acc. to IEC 61000-4-11	Voltage Dip: Dip to 0 % for 1 cycle @ 0° phase angle Dip to 70 % for 25/30 cycles @ 0° phase angle Dropout to 0 % for 0.5 cycle @ 0°, 45°, 90°, 135°, 180°, 225°, 270° & 315° phase angles Voltage interruption: 100 % for 250/300 cycles	Voltage Dip: Dip to 0 % for 1 cycle @ 0° phase angle Dip to 70 % for 25/30 cycles @ 0° phase angle Dropout to 0 % for 0.5 cycle @ 0°, 45°, 90°, 135°, 180°, 225°, 270° & 315° phase angles Voltage interruption: 100 % for 250/300 cycles	Mains power quality should be that of a typical commercial or hospital environment. If the user of the IMAGE1 S™ 4U Rubina, OPAL1® NIR/ICG camera head require continued operation during power mains interruptions, it is recommended that the IMAGE1 S™ 4U Rubina, OPAL1® NIR/ICG camera head be powered from an uninterruptible power supply or a battery.
Magnetic field at the supply frequency • (50/60 Hz) acc. to IEC 61000-4-8	30 A/m at 50Hz/60 Hz	30 A/m at 50 Hz/60 Hz	If image distortion occurs, it may be necessary to position the IMAGE1 S™ 4U Rubina, OPAL1® NIR/ICG camera head further from sources of power frequency magnetic fields or to install magnetic shielding. The power frequency magnetic field should be measured in the intended installation location to assure that it is sufficiently low.
Radiated, radio-frequency, electromagnetic field immunity test acc. to IEC 61000-4-3	3 V/m 80 MHz to 2.7 GHz *Refer table 2 for wireless Proximity RF field test levels	3 V/m 80 MHz to 2.7 GHz	
Immunity to conducted disturbances, induced by radio-frequency fields acc. to IEC 61000-4-6	3 V _{rms} on 150 kHz to 80 MHz 1 kHz 80% AM modulation 6 V _{rms} in ISM bands	3 V _{rms} on 150 kHz to 80 MHz 1 kHz 80 % AM modulation 6 V _{rms} in ISM bands	

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Table 3 – Test levels for Radiated and conducted Immunity Tests
Guidance and manufacturer's declaration – electromagnetic immunity

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Immunity tests	EN/IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF disturbances acc. to IEC 61000-4-6	3 V _{rms} 150 kHz to 80 MHz	3 V _{rms}	Portable and mobile RF communications equipment should be used no closer to any part of the IMAGE1 S™ 4U Rubina, OPAL1® NIR/ICG camera head, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = 1.2 \sqrt{P}$ where P is the nominal power of the transmitter in watts [W] according to the information provided by the manufacturer of the transmitter and d is the recommended separation distance in meters [m]. Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a should be less than the compliance level in each frequency range ^b .
Radiated RF disturbances acc. to IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	$d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2.3 \sqrt{P}$ 800 MHz to 2.5 GHz Interference may occur in the vicinity of equipment marked with the following symbol: 

Note: At 80 MHz and 800 MHz, the higher frequency range applies.

Note: These guidelines may not apply in all situations. The propagation of electromagnetic quantities is affected by absorptions and reflections of buildings, objects, and persons.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the 4U Rubina, OPAL1® NIR/ICG camera head is used exceeds the applicable RF compliance level above, the 4U Rubina, OPAL1® NIR/ICG camera head should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the 4U Rubina, OPAL1® NIR/ICG camera head.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

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Table 4 – Emission class and group		
Guidelines and manufacturer's declaration – electromagnetic emissions		
The IMAGE1 S™ 4U Rubina, OPAL1® NIR/ICG camera head is intended for use in the electromagnetic environment specified below. The customer or user of the IMAGE1 S™ 4U Rubina, OPAL1® NIR/ICG camera head should ensure that it is used in such an environment.		
RF emissions as per CISPR 11	Group 1	The IMAGE1 S™ 4U Rubina, OPAL1® NIR/ICG camera head uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause interference in nearby electronic equipment.
RF Emissions as per CISPR 11	Class A	The IMAGE1 S™ 4U Rubina, OPAL1® NIR/ICG camera head is suitable for use in all establishments other than domestic and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions as per IEC 61000-3-2	Class A	
Voltage fluctuations/flicker Emissions as per IEC 61000-3-3	Complies	

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Table 5

Recommended separation distances between portable and mobile RF communications equipment and the IMAGE1 S™ 4U Rubina, OPAL1® NIR/ICG camera head

The IMAGE1 S™ 4U Rubina, OPAL1® NIR/ICG camera head is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the IMAGE1 S™ 4U Rubina, OPAL1® NIR/ICG camera head can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the IMAGE1 S™ 4U Rubina, OPAL1® NIR/ICG camera head as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter [W]	Separation distance d [m] as a function of the transmitter frequency		
	150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance D in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note: At 80 MHz and 800 MHz, the separation distance for higher frequency range applies.

Note: These guidelines may not apply in all situations. The propagation of electromagnetic quantities is affected by absorptions and reflections of structures, objects and people.

HF surgical instrument compatibility

This equipment has been evaluated for compatibility with high-frequency surgical equipment. It has been tested according to IEC 60601-2-2 Annex BB.

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KARL STORZ Endoscopy Canada Ltd.
7171 Millcreek Drive, Mississauga, Ontario
L5N 3R3 Canada
Phone: +1 905 816-4500, Fax: +1 905 816-4599
Toll free phone: 1-800-268-4880 (Canada only)
Toll free fax: 1-800-482-4198 (Canada only)
E-Mail: info-canada@karlstorz.com

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100, Fax: +1 424 218-8525
Toll free phone: 800 421-0837 (USA only)
Toll free fax: 800 321-1304 (USA only)
E-Mail: communications@kse.com

KARL STORZ Veterinary Endoscopy-America, Inc.
1 South Los Cameros Road
Goleta, CA 93117, USA
Phone: +1 805 968-7776, Fax: +1 805 685-2588
E-Mail: info@karlstorzvet.com

KARL STORZ Endoscopia Latino-America, Inc.
815 N. W. 57th Avenue, Suite 480
Miami, FL 33126-2042, USA
Phone: +1 305 262-8980, Fax: +1 305 262-8986
E-Mail: info@kse.com

KARL STORZ Endoscopia México S.A. de C.V.
Edificio Atlántic, Oficina 3G
Calle D e/ 1ra y 3ra
10400 Vedado, Havana, Cuba
Phone: +537 836 95 06, Fax: +537 836 97 76
E-Mail: kstorzuba@gmail.com

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TOO KARL STORZ Endoscopy Kazakhstan
Saryarka, 6, BC "Arman", off. 910
010000 Astana, Republic of Kazakhstan
Phone: +7 7172 552-549, 552-788, Fax: -444
E-Mail: info@karlstorz.kz

KARL STORZ ENDOSKOPE
East Mediterranean & Gulf (branch)
Building West Side 7A – Unit 7WA – 3008
Dubai Airport Free Zone, P.O. Box 54983
Dubai – United Arab Emirates
Phone: +971 (0)4 2958887, Fax: +971 (0)4 3205282
Service Hotline: +971 (0)4 3415882
E-Mail: info-gne@karlstorz-emg.com

KARL STORZ Endoscopy India Private Limited
11th Floor, Dr. Gopal Das Bhawan
28, Barakhamba Road
New Delhi 110001, India
Phone: +91 11 4374 3000, Fax: +91 11 4374 3010
E-Mail: corporate@karlstorz.in

KARL STORZ SE & Co. KG
Interchange 21 Tower, Level 33
309 Sukhumvit Road,
North Klongtoey, Wattana,
10110 Bangkok, Thailand
Phone: +66 2 660 3669
E-Mail: info-th@karlstorz.com

KARL STORZ SE & Co. KG
Resident Representative Office
14th Floor, M Plaza Saigon
39 Le Duan, District 1
Ho Chi Minh City, Vietnam
Phone: +84 28 3823 8000, Fax: +84 28 3823 8039
E-Mail: infovietnam@karlstorz.com

KARL STORZ Endoscopy China Ltd.
Room 2503-05, 25F AXA Tower, Landmark East,
No. 100 Hong Ming Street, Kwun Tong, Kowloon,
Hong Kong, People's Republic of China
Phone: +852 28 65 2411, Fax: +852 28 65 4114
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd., Beijing Branch
Room 1805-1807, Building B, 18F Beijing IFC
No. 8, Jianguomenwai Street, Chaoyang District,
100022, Beijing, People's Republic of China
Phone: +86 10 5638188, Fax: +86 10 5638199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Shanghai Branch
Room 701A Building 5 & Room 501 Building 7,
No. 3000 Longdong Avenue, Pilot Free Trade Zone,
201203, Shanghai, People's Republic of China
Phone: +86 21 60339888, Fax: +86 21 60339808
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Chengdu Branch
Room 803-805, 8F Jin Jiang International Building
No. 1 West Linjiang Road, Wuhou District,
6100414, Chengdu, People's Republic of China
Phone: +86 28 86587977, Fax: +86 28 86587975
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Shenyang Branch
Room 2001-2005, 20F N-MEDIA International Center,
No. 167 Youth Avenue, Shenhe District,
110014, Shenyang, People's Republic of China
Phone: +86 24 23181118, Fax: +86 24 23181119
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Guangzhou Branch
Room 02B & 03 & 04A, 35F Team Tower,
No. 208 Tianhe Road, Tianhe District,
510620, Guangzhou, People's Republic of China
Phone: +86 20 87321281, Fax: +86 20 87321286
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd.
Room 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com

KARL STORZ Endoscopy Singapore Sales Pte Ltd
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com

KARL STORZ SE & Co. KG
Representative Office Indonesia
Sinarmas MSIG Tower Level 37
Jl. Jend. Sudirman No. Kav. 21
Jakarta Selatan
DKI Jakarta 12920
E-Mail: infoindonesia@karlstorz.com

KARL STORZ Endoscopy Korea Co. Ltd.
9F Hyowon-Building
97, Jungdae-ro, Songpa-gu
05719 Seoul, Korea
Phone: +82-70-4350-7474, Fax: +82-70-8277-3299
E-Mail: infokorea@karlstorz.com

KARL STORZ Endoscopy Taiwan Ltd.
12F, No. 192, Sec. 2, Chung Hsin Rd.,
Sindian District, New Taipei City, Taiwan
Phone: +886 933 014 160, Fax: +886 2 8672 6399
E-Mail: info-tw@karlstorz.com

KARL STORZ SE & Co. KG
Representative Office Philippines
1901 Picadilly Star Bldg., 4th Avenue, BGC
Taguig City 1636, Philippines
Phone: +63 2 317 45 00, Fax: +63 2 317 45 11
E-Mail: philippines@karlstorz.com

KARL STORZ Endoscopy Japan K.K.
Stage Bldg. 8F, 2-7-2 Fujimi
Chiyoda-ku, Tokyo 102-0071, Japan
Phone: +81 3 6380-8622, Fax: +81 3 6380-8633
E-Mail: info-jp@karlstorz.com

KARL STORZ Endoscopy New Zealand Ltd.
31 Morningside Drive Mt Albert
Auckland, 1025, New Zealand
PO Box 56 511, Dominion Rd
Auckland, 1446, New Zealand
Phone: +64 9 846 6044, Fax: +64 9 846 6808
Toll free: +64 508 84 84 84 (New Zealand only)
E-Mail: sales@karlstorz.co.nz

KARL STORZ Endoscopy Australia Pty. Ltd.
68 Waterloo Road, Macquarie Park NSW 2113
P.O. Box 50 Lane Cove NSW 1595, Australia
Phone: +61 (0)2 9490 6700
Toll free: 1800 996 562 (Australia only)
Fax: +61 (0)2 9420 0695
E-Mail: karlstorz@karlstorz.com.au

KARL STORZ Endoscopy (Shanghai) Ltd., Chengdu Branch
Room 803-805, 8F Jin Jiang International Building
No. 1 West Linjiang Road, Wuhou District,
6100414, Chengdu, People's Republic of China
Phone: +86 28 86587977, Fax: +86 28 86587975
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Shenyang Branch
Room 2001-2005, 20F N-MEDIA International Center,
No. 167 Youth Avenue, Shenhe District,
110014, Shenyang, People's Republic of China
Phone: +86 24 23181118, Fax: +86 24 23181119
E-Mail: info@karlstorz.com.cn

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ENDOSCOPES FOR MEDICINE AND TECHNICAL SCIENCE
INSTRUMENTS FOR OTO-RHINO-LARYNGOLOGY

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
IMAGE1 S 4U RUBINA OPAL1 NIR\ICG camera head for endoscopy

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* Formerly trading as KARL STORZ GmbH & Co. KG

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
The information contained in this Addition takes precedence in the Russian Federation over the information specified in the corresponding instruction manual for the medical device.

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Fax: +49 7461 708-105

E-mail: info@karlstorz.com

AUTHORIZED REPRESENTATIVE OF THE MANUFACTURER:

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115114, Moscow, Derbenevskaya embankment, 7, p. 4

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Compatible Camcorder Control Units:

RU RZN 2015/3274 "Karl Storz endoscopic equipment for visualization and image processing, with accessories: Endoscopic camera head unit IMAGE1 S 4U-LINK" manufactured by Karl Storz

RU RZN 2015/3274 "Karl Storz endoscopic equipment for visualization and image processing, with accessories: Endoscopic camera head unit IMAGE1 S CONNECT II" manufactured by Karl Storz

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One IMAGE1 S CONNECT II module is at the heart of the IMAGE1 S and acts as a communication interface between the image and documentation devices. For use, one IMAGE1 S CONNECT II plus at least one LINK module is required which receives and forwards the signal. The camera head is connected to the IMAGE1 S 4U-LINK module.

Compatible light source:

Cold light Fountain Power LED Rubina, OPAL NIR/ICG with accessories manufactured by Karl Storz (in the process of registration).

Compatible endoscopes:

"Rigid endoscopes HOPKINS Rubina NIR/ICG, in versions" manufactured by Karl Storz (in the process of registration).

These control units are at the same time a power supply for the camera heads and together with the camera heads, light source and endoscopes are considered the ME System. Only the camera control units specified in this Addendum should be used with the system.

The rated voltage of the supply mains of the camera control units compatible with the product, which are the power source of the camera heads, does not exceed 250 V.

Inout electrical characteristics of CCU

Characteristics	IMAGE1 S 4U-LINK	IMAGE1 S CONNECT II
Power supply	100-240 VAC	100-240 VAC
Operating frequency	50/60 Hz	50/60 Hz
Power input	86 VA	100-135 VA

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- Legal and regulatory requirements for reprocessing medical products
- The safety of patients, users and third parties
- Lowering the status of contamination to unavoidable residual risks

The KARL STORZ material compatibility list is your guide for the reprocessing of rigid and flexible endoscopes as well as endoscopic instruments and provides the best possible value preservation of your investments. The chemicals listed here have been tested and approved by KARL STORZ with regard to material compatibility.

A differentiation is made between manual cleaning and disinfection/ultrasonic cleaning and machine cleaning and disinfection for which various chemical agents are used. Material compatibility approval is based on tests on various processing chemicals. Therefore, any possible chemical cross-reactions due to insufficient neutralization/rinsing between process steps or incorrect dosage/exposure time can not be ruled out. This particularly applies to subsequent low-temperature sterilization as this process also uses chemical agents.

To ensure efficacious reprocessing conforming to legal and statutory requirements, please observe user manuals for each specific product and the detailed user manual "Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments".

① Note:

Current situation

KARL STORZ SE & Co. KG does currently not test any reprocessing chemicals for material compatibility. This program is stopped and will be replaced by a new concept which is actually under development. At this point, no further chemicals will be added to the list of approved detergents any more until the new program is launched.

In the meantime, the options for detergent manufacturers are to either conduct individual tests on flexible endoscopes and standard instruments to confirm material compatibility, or to compare their chemicals with one that is already listed as approved for material compatibility and base a confirmation on this assessment where applicable.

1. Manual cleaning and disinfection/ultrasonic cleaning

The substances listed here have been approved and tested by KARL STORZ with regard to material compatibility. This list is not a statement concerning the microbiological effectiveness of the

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products. Please contact the chemical manufacturer regarding microbiological effectiveness.

The following KARL STORZ instruments are not suitable for ultrasonic cleaning:

- Rigid telescopes
- Flexible and semiflexible fiberscopes
- Other instruments with optical components made of glass, with chips or with electronic components

1.1 General warnings

Warning:

During all cleaning and disinfection procedures on contaminated instruments, please observe the occupational safety guidelines of the employers' liability insurance association (Berufsgenossenschaft) or equivalent organization.

Warning:

Gigasept FF from the chemical manufacturer Schülke & Mayr GmbH at a concentration of 8% and a contact time of 60 minutes is not approved for fiberscopes.

Warning:

In the case of flexible endoscopes, first vent and then remove the leakage tester so that the pressure compensation valve closes again.

Warning:

In the case of flexible endoscopes, the pressure compensation cap must be removed from the vent connection during liquid immersion/chemical disinfection.

Caution:

Please follow the instructions "Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments." These provide a detailed description of cleaning, disinfection, and sterilization procedures.

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Caution:

Damage to materials: Cleaning the outer sheath of flexible fiberscopes with ammonium chlorides will damage the material, particularly at the hose nozzle.

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Cleaning all of the materials with alcohol-based cleaning agents can lead to damage to the products as well as contamination becoming fixed to the surface of the products.

! Caution:

Please note that when preparing and using the solutions, the chemical manufacturers' instructions on concentration and contact time must be strictly followed. Excessively long immersion or incorrect concentrations may cause damage. Please contact the chemical manufacturer regarding microbiological effectiveness.

! Caution:

Using SODIUM HYDROXIDE SOLUTION may lead to surface alterations in aluminum parts (even if coated), plastics, and soldered joints, and it may reduce the service life of the instruments.

! Caution:

Before placing the flexible endoscope in the solution, check leak tightness since any liquid penetration may cause damage.

! Caution:

Regardless of the selected cleaning and disinfection procedure, flexible endoscopes should not be heated to temperatures above 60 °C.

! Caution:

Endoscopes should not be cleaned in an ultrasound bath.

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① Note:

National laws and regulations must be observed.

In the case of semi-critical medical devices that are used in a non-sterile condition, disinfectants with full virucidal efficacy should be used.

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① Note:

The use of cleaning and disinfectant agents with pH values outside the neutral range can cause color-changes in color-anodized components. These changes are cosmetic and do not affect the functionality of the components.

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Bandelin electronic GmbH	STAMMOPUR GR	✓	-	✓	✓	✓	✓	-	-
Bandelin electronic GmbH	STAMMOPUR R	✓	-	✓	✓	✓	✓	-	-
Bochemie a.s.	Chirosan	-	✓	-	✓	✓	✓	-	-
Bochemie a.s.	Chirosan plus	✓	✓	-	✓	✓	✓	-	-
Bochemie a.s.	Discleen Endo AF	-	✓	-	✓	✓	✓	-	-
Bochemie a.s.	Discleen Endo PAA	✓	✓	-	✓	✓	✓	✓	✓
Bochemie a.s.	Discleen enzyme	✓	-	✓	✓	✓	✓	✓	✓
Bochemie a.s.	Discleen Extra	✓	✓	-	✓	✓	✓	✓	✓
Bode Chemie GmbH	Bodedex® forte	✓	-	✓	✓	✓	✓	✓	✓
Bode Chemie GmbH	Bomix® plus	✓	✓	✓	✓	✓	✓	✓	✓
Bode Chemie GmbH	Dismoclean® twin basic/ Dismoclean® twin zyme	✓	-	✓	✓	✓	-	-	-
Bode Chemie GmbH	Korsolex® basic	-	✓	-	✓	✓	✓	✓	✓
Bode Chemie GmbH	Korsolex® extra	-	✓	-	✓	✓	✓	✓	✓
Bode Chemie GmbH	Korsolex® FF	-	✓	-	-	-	-	✓	✓
Bode Chemie GmbH	Korsolex® med AF	-	✓	-	✓	✓	✓	-	-
Bode Chemie GmbH	Korsolex® PAA	-	✓	-	✓	✓	✓	✓	✓
Borer Chemie AG	deconex® 50 FF	-	✓	-	✓	✓	-	-	-

Borer Chemie AG	deconex® 53 INSTRUMENT	✓	✓	✓	✓	✓	✓	-	-
Borer Chemie AG	deconex® 53 PLUS	✓	✓	✓	✓	✓	✓	✓	✓
Borer Chemie AG	deconex® 54 SPORCIDE	-	✓	-	✓	✓	✓	✓	✓
Borer Chemie AG	deconex® HLD PA/PA20	-	✓	-	✓	✓	✓	✓	✓
Borer Chemie AG	deconex® POWER ZYME	✓	-	✓	✓	✓	✓	✓	✓
Borer Chemie AG	deconex® TWIN PH10 deconex® TWIN ZYME	✓	-	✓	✓	✓	-	-	-
Calgon Vestal Laboratories	Klenzyme	✓	-	✓	-	-	-	✓	✓
DuPont	PeraSafe	-	✓	-	✓	✓	✓	-	-
DuPont	Rely+On™ PeraSafe®	✓	✓	-	✓	✓	✓	✓	✓
Deppe, Laboratorium Dr. rer. nat.	EndoStar	-	✓	-	✓	✓	✓	✓	✓
Deppe, Laboratorium Dr. rer. nat	InstruPlus®	✓	✓	✓	✓	✓	✓	✓	✓
Deppe, Laboratorium Dr. rer. nat	InstruPlus N	-	✓	✓	✓	✓	✓	-	-
Deppe, Laboratorium Dr. rer. nat	InstruStar®	✓	✓	✓	✓	✓	-	-	-
Deppe, Laboratorium Dr. rer. nat	InstruZym®	✓	-	✓	✓	✓	✓	✓	✓
DERCLEAN	DERDEVICE PLUS A2	-	✓	-	✓	✓	✓	-	-
DERCLEAN	DERDEVICE PLUS OPA	-	✓	-	✓	✓	✓	-	-
Detrox	Detro Forte®	-	✓	-	✓	✓	-	-	-

Detrox	Detro OPA®	-	✓	-	✓	✓	-	-	-
Detrox	Detro Plus®	-	✓	-	✓	✓	✓	-	-
Detrox	Detro Plus OPA®	-	✓	-	✓	✓	✓	-	-
Ecolab Deutschland GmbH	Sekusept™ aktiv	✓	✓	-	✓	✓	✓	✓	✓
Ecolab Deutschland GmbH	Sekusept™ Cleaner	✓	-	✓	✓	✓	-	✓	✓
Ecolab Deutschland GmbH	Sekusept™ easy	-	✓	-	✓	✓	✓	✓	✓
Ecolab Deutschland GmbH	Sekusept™ Extra N	✓	✓	✓	✓	✓	✓	-	-
Ecolab Deutschland GmbH	Sekusept™ forte	-	✓	✓	✓	✓	✓	✓	✓
Ecolab Deutschland GmbH	Sekusept™ MultiEnzyme	✓	-	✓	✓	✓	✓	✓	✓
Ecolab Deutschland GmbH	Sekusept™ PLUS	✓	✓	✓	✓	✓	✓	✓	✓
Ecolab Deutschland GmbH	Sekusept™ Pulver CLASSIC	✓	✓	-	✓	✓	✓	✓	✓
Esteer Pharma GmbH	ULTRA- DESMIT® AF	✓	✓	-	✓	✓	✓	-	-
Esteer Pharma GmbH	ULTRASEPTIN® aktiv	✓	✓	-	✓	✓	✓	-	-
Esteer Pharma GmbH	ULTRASEPTIN® Classic	✓	✓	✓	✓	✓	✓	-	-
Esteer Pharma GmbH	ULTRASEPTIN® PAA	✓	✓	-	✓	✓	✓	✓	✓
Fresenius AG	Afid	✓	✓	-	✓	✓	✓	-	-
Fresenius AG	Afid plus	-	✓	-	✓	✓	✓	-	-

Fresenius AG	Sporizid	-	✓	-	-	-	-	✓	-
Fresenius AG	Sporizid FF	✓	✓	-	✓	✓	✓	✓	✓
Heck Dental GmbH	Ventisept M Plus Neu	✓	✓	-	✓	✓	✓	✓	✓
Helvemed	Instrument Forte	✓	✓	✓	✓	✓	✓	✓	✓
HIMED GmbH	EndoSeptol FF	-	✓	-	✓	✓	✓	-	-
HIMED GmbH	InstruCid AF	✓	✓	✓	✓	✓	-	-	-
Holifa Poska Sp. z.o.o	Poisept Holifa	✓	✓	-	✓	✓	✓	-	-
Holifa Poska Sp. z.o.o	Poisept Plus	✓	✓	-	✓	✓	✓	✓	✓
Indeba Industria E Comercio LTDA	Letahdeido	-	✓	-	✓	✓	✓	-	-
INIBSA Laboratorios	INSTRUNET F.A. CONCEN.	✓	✓	✓	✓	✓	✓	-	-
Jose Collado S.A.	Darodor® 4000 Liquido	✓	-	-	✓	✓	✓	-	-
Jose Collado S.A.	Darodor® 9000	-	✓	-	✓	✓	✓	-	-
Jose Collado S.A.	Darodor® Sinaldehyd 2000	-	✓	-	✓	✓	✓	-	-
Lonza AG	Formulierung ID 50	-	✓	-	✓	✓	✓	-	-
Lysoform Dr. Rosemann GmbH	Almyrol®	✓	✓	✓	✓	✓	✓	-	-
Lysoform Dr. Rosemann GmbH	Desoform®	✓	✓	✓	✓	✓	✓	-	-

Lysoform Dr. Rosemann GmbH	Lysoformin® 2000	✓	✓	✓	✓	✓	✓	-	-
Lysoform Dr. Rosemann GmbH	Lysoformin® 3000	✓	✓	-	✓	✓	✓	-	-
Merz Hygiene GmbH	Edisonite Super	✓	-	✓	✓	✓	✓	-	-
Merz Hygiene GmbH	Mucadont-IS	-	✓	✓	✓	✓	✓	-	-
Merz Hygiene GmbH	Mucocit-T neu	✓	✓	✓	✓	✓	✓	✓	✓
Merz Hygiene GmbH	Mucadont- Zymaktiv	✓	-	✓	✓	✓	✓	✓	✓
METASYS Medizintechnik GmbH	ID 50	-	✓	-	✓	✓	✓	-	-
METASYS Medizintechnik GmbH	Green & Clean ID	✓	✓	✓	✓	✓	-	-	-
Metrex Research Corporation	Coldspor	-	✓	-	-	-	-	✓	✓
Metrex Research Corporation	EmPower™	✓	-	-	-	-	-	✓	✓
Metrex Research Corporation	MetriCide™	-	✓	-	-	-	-	✓	✓
Metrex Research Corporation	MetriCide™ OPA Plus Solution	-	✓	-	✓	✓	✓	✓	✓
Metrex Research Corporation	MetriZyme™	✓	-	-	-	-	-	✓	✓
Miltex, Inc.	EZ – Zyme	✓	-	-	-	-	-	✓	✓
Minntech International	Adaspor®	-	✓	-	✓	✓	✓	-	-
Mir Dezinfektsii	MIROXIDE® – 2000	-	✓	-	✓	✓	✓	-	-
Necmkimya	BIORAD GDA 2%	-	✓	-	✓	✓	-	-	-

Necmkimya	BIORAD Rayd Guanide	✓	✓	✓	✓	✓	-	-	-
Nova-Praxis-Hygiene	NOVA Biguades AF	✓	✓	✓	✓	✓	✓	-	-
PRISMAN GmbH	POW(D)ER Konz. IDP 700	✓	✓	-	✓	✓	✓	✓	✓
PRISMAN GmbH	IDP-ic 70 POWDER	✓	✓	-	✓	✓	✓	✓	✓
Profarma	Ajatin® Extra Strong	-	✓	-	✓	✓	✓	✓	✓
Promagent AB	Wavacide	-	✓	-	✓	✓	✓	✓	✓
Rheosol Wachendorff – Chemie GmbH	RHEOSEPT-ID aktiv Granulat	✓	✓	✓	✓	✓	✓	✓	✓
Rheosol Wachendorff – Chemie GmbH	RHEOSEPT-ID plus+	✓	✓	✓	✓	✓	✓	✓	✓
Ruhof Corporation	Endozime®	✓	-	-	-	-	-	✓	✓
Schumacher, Dr. GmbH	DESCOTON extra	✓	✓	-	✓	✓	✓	✓	✓
Schumacher, Dr. GmbH	DESCOTON forte	-	✓	-	✓	✓	✓	✓	✓
Schumacher, Dr. GmbH	DESCOTON 2% GDA	-	✓	-	✓	✓	✓	✓	✓
Schumacher, Dr. GmbH	INSTRU PLUS flussig	✓	-	✓	-	-	-	✓	✓
Schumacher, Dr. GmbH	MULTIZYME FORTE	✓	-	✓	✓	✓	✓	✓	✓
Schumacher, Dr. GmbH	PERFEKTAN® ENDO	✓	✓	✓	✓	✓	✓	✓	✓
Schumacher, Dr. GmbH	PERFEKTAN® NEU	✓	✓	✓	✓	✓	✓	✓	✓
Schumacher, Dr. GmbH	PERFEKTAN® TB	✓	✓	✓	✓	✓	-	-	-

Schumacher, Dr. GmbH	PLURAZYME	✓	-	✓	✓	✓	✓	✓	✓
Schulke & Mayr GmbH	gigasept AF® forte	✓	✓	-	✓	✓	✓	✓	✓
Schulke & Mayr GmbH	gigasept FF® (neu)	-	✓	-	✓	✓	✓	-	-
Schulke & Mayr GmbH	gigasept® instru AF	✓	✓	✓	✓	✓	-	-	-
Schulke & Mayr GmbH	gigasept® med forte	✓	✓	✓	✓	✓	✓	✓	✓
Schulke & Mayr GmbH	gigasept® PAA concentrate	-	✓	-	✓	✓	✓	✓	✓
Schulke & Mayr GmbH	gigazyme®	✓	-	✓	✓	✓	✓	✓	✓
Schulke & Mayr GmbH	Lysetol® V	-	✓	-	✓	✓	✓	-	-
Schulke & Mayr GmbH	mikrozyd® PAA	-	✓	-	✓	✓	✓	✓	✓
Sirmaxo Chemicals Pvt. Ltd.	Sanidex OPA	-	✓	-	✓	✓	✓	-	-
Sirmaxo Chemicals Pvt. Ltd.	Saniscope	-	✓	-	✓	✓	✓	✓	✓
Sirmaxo Chemicals Pvt. Ltd.	Sanizyme	✓	-	-	✓	✓	✓	✓	✓
Steril-4 s.r.l.	STERIL-C	-	✓	-	✓	✓	✓	-	-
STERIS Corporation	EnzyCare 2	✓	-	✓	✓	✓	✓	-	-
STERIS Corporation	Revital-Ox® Resert™ HLD	-	✓	-	✓	✓	✓	✓	✓
STERIS Corporation	Revital-Ox™ Resert™ HLD Chemosterilant	-	✓	-	✓	✓	✓	✓	✓
STERIS Corporation	Resert™ XL HLD	-	✓	-	✓	✓	✓	✓	✓

Tristel Solutions Ltd.	Medistel	-	✓	✓	✓	✓	✓	-	-
Tristel Solutions Ltd.	Enzystel	✓	-	✓	✓	✓	✓	-	-
Tristel Solutions Ltd.	Stella™ System mit: Tristel™ Clean for Stella Tristel Fuse™ for Stella	✓	✓	-	-	-	✓	✓	-
Tristel Solutions Ltd.	Perastel 3 Day	-	✓	-	✓	✓	✓	✓	✓
Tristel Solutions Ltd.	Tristel Wipes System	✓	✓	-	✓	-	✓	✓	✓
Wavi Energie System INC.	Wavicide-01	-	✓	-	-	-	-	✓	✓
Weigert, Dr. GmbH & Co.	neodisher endo® CLEAN	✓	-	✓	✓	✓	✓	✓	✓
Weigert, Dr. GmbH & Co.	neodisher endo® DIS active	✓	✓	✓	✓	✓	✓	✓	✓
Weigert, Dr. GmbH & Co.	neodisher® LM 2	✓	-	✓	✓	✓	✓	✓	✓
Weigert, Dr. GmbH & Co.	neodisher® MediClean	✓	-	✓	✓	✓	✓	✓	✓
Weigert, Dr. GmbH & Co.	neodisher® MediClean forte	✓	-	✓	✓	✓	✓	✓	✓
Weigert, Dr. GmbH & Co.	neodisher® MediZym	✓	-	✓	✓	✓	✓	✓	✓
Weigert, Dr. GmbH & Co.	neodisher® PreStop	✓	-	-	✓	✓	-	-	-
Weigert, Dr. GmbH & Co.	neodisher® Septo 3000	-	✓	✓	-	-	-	✓	✓
Weigert, Dr. GmbH & Co.	neodisher® Septo Active	-	✓	✓	✓	✓	-	-	-
Weigert, Dr. GmbH & Co.	neodisher® SeptoClean*	✓	✓	✓	✓	✓	-	-	-

Weigert, Dr. GmbH & Co.	neodisher® Septo MED	✓	✓	✓	✓	✓	✓	✓	✓
Weigert, Dr. GmbH & Co.	neodisher® Septo PreClean	✓	✓	✓	✓	✓	✓	✓	✓
Whiteley Industries Pty. Ltd.	Aidal® Plus	-	✓	-	✓	✓	✓	-	-
Whiteley Industries Pty. Ltd.	Matrix®	✓	-	-	✓	✓	✓	✓	✓

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- ✓ Suitable chemical agents
- Chemicals not tested

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2. Machine cleaning and disinfection

The substances listed here have been approved and tested by KARL STORZ with regard to material compatibility. This list is not a statement concerning the microbiological effectiveness of the products. Please contact the chemical manufacturer regarding microbiological effectiveness.

Select telescope cleaning and disinfection procedures in consultation with the manufacturers of the machine and the chemicals. Only use special processes that have been validated for this purpose.

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! Caution:

To minimize stress on materials, avoid repeated switching between processing procedures.

! Caution:

Flexible endoscopes should never be steam sterilized (autoclaved).

! Caution:

The pressure compensation cap must be removed from the vent connection during liquid immersion/chemical disinfection.

! Caution:

Please note that when preparing and using the solutions, the **chemical manufacturers' instructions on concentration and contact time must be strictly followed**. Excessively long immersion or incorrect concentrations may cause damage. Please contact the chemical manufacturer regarding microbiological effectiveness.

! Caution:

Regardless of the selected cleaning and disinfection procedure, **flexible endoscopes should not be heated to temperatures above 60 °C**.

! Caution:

The use of rinse aids is not approved for plastics. Rinse aids have a negative effect on plasticizers, which may lead to breakage.

! Caution:

The product thermosept RKF forte from Schülke & Mayr GmbH is not suitable for aluminum materials.

i Note:

The C-MAC® video laryngoscope series (840XXX) is considered to be a video endoscope.

i Note:

For detailed information on processing, please refer to the respective instructions for use.

i Note:

Please contact the chemical manufacturer regarding the range of action of the chemical agents used.

i Note:

The use of cleaning and disinfectant agents with pH

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values outside the neutral range can cause color-changes in color-anodized components. These changes are cosmetic and do not affect the functionality of the components.

① Note:

When using cleaning and disinfecting agents, chemicals should be carefully coordinated in order to avoid any chemical interactions leading to material damage. If you have any questions or doubts, please contact the chemical manufacturer.

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Bode Chemie GmbH	Dismoclean® 28 Alka one	✓	-	✓	✓	✓	✓	-	-
Bode Chemie GmbH	Dismoclean® twin basic/ Dismoclean® twin zyme	✓	-	✓	✓	✓	-	-	-
Bode Chemie GmbH	Korsolex® Endo- Cleaner	✓	-	✓	✓	✓	✓	✓	✓
Bode Chemie GmbH	Korsolex® Endo- Disinfectant	-	✓	✓	✓	✓	✓	✓	✓
Borer Chemie AG	deconex® ALKAMILD	✓	-	-	✓	✓	✓	-	-
Borer Chemie AG	deconex® ENDOMATIC	-	✓	-	-	-	✓	✓	-
Borer Chemie AG	deconex® POWER ZYME	✓	-	✓	✓	✓	✓	✓	✓
Borer Chemie AG	deconex® TWIN PH10 deconex® TWIN ZYME	✓	-	✓	✓	✓	-	-	-
B. Braun Medical AG	Helimatic® Cleaner alkaline	✓	-	-	✓	✓	✓	-	-
B. Braun Medical AG	Helimatic® Cleaner enzymatic	✓	-	-	✓	✓	✓	✓	✓
B. Braun Medical AG	Helimatic® Cleaner neutral	✓	-	-	✓	✓	-	-	-
B. Braun Medical AG	Helimatic® Disinfectant	-	✓	-	✓	✓	✓	✓	✓
Deppe, Laboratorium Dr. rer. nat.	Endomat Plus	-	✓	-	✓	✓	✓	✓	✓
Deppe, Laboratorium Dr. rer. nat.	EndomatPlus viruguard	-	✓	-	✓	✓	✓	✓	✓
Deppe, Laboratorium Dr. rer. nat.	Labomat E	✓	-	-	✓	✓	✓	✓	✓
DiverseyLever	Sumatox E	✓	-	✓	✓	✓	✓	-	-

Ecolab Deutschland GmbH	MetalClean Plus	✓	✓	-	✓	✓	✓	-	-
Ecolab Deutschland GmbH	Sekumatic™ FD	-	✓	-	✓	✓	✓	-	-
Ecolab Deutschland GmbH	Sekumatic™ FR	✓	-	-	✓	✓	✓	-	-
Ecolab Deutschland GmbH	Sekumatic™ FRE	✓	-	-	✓	✓	✓	-	-
Ecolab Deutschland GmbH	Sekumatic™ MultiClean	✓	-	-	✓	✓	✓	-	-
IMS	Adapspor® Ready to Use	-	✓	✓	✓	✓	✓	-	-
Medisafe UK, Ltd.	3E-Zyme/HS-Zyme	✓	-	✓	✓	✓	✓	-	-
MEDIVATORS Inc.	Rapicide®	-	✓	-	✓	✓	✓	✓	✓
MEDIVATORS Inc.	Rapicide® PA	-	✓	-	✓	✓	✓	✓	✓
Merz Consumer Care GmbH	mucapur® AF	✓	-	✓	✓	✓	✓	-	-
Merz Consumer Care GmbH	mucapur® ED	-	✓	✓	✓	✓	✓	✓	✓
Merz Consumer Care GmbH	mucapur® ED plus	-	✓	✓	✓	✓	✓	✓	✓
Merz Consumer Care GmbH	mucapur® ER	✓	-	✓	✓	✓	✓	✓	✓
Merz Consumer Care GmbH	mucapur® ER plus	✓	-	✓	✓	✓	✓	✓	✓
Ruhof Corporation	Endozime® AW	✓	-	✓	✓	✓	✓	-	-
Schumacher, Dr. GmbH	THERMOTON® CLEANER	✓	-	-	✓	✓	✓	-	-
Schumacher, Dr. GmbH	THERMOTON® ENDO	-	✓	✓	✓	✓	✓	✓	✓

Schumacher, Dr. GmbH	THERMOTON® DESIN- FEKTANT	-	✓	-	✓	✓	✓	✓	✓
Schumacher, Dr. GmbH	THERMOTON® NR	✓	-	-	✓	✓	✓	✓	✓
Schulke & Mayr GmbH	thermosept® alka clean	✓	-	-	✓	✓	✓	-	-
Schulke & Mayr GmbH	thermosept® alka clean forte	✓	-	-	✓	✓	✓	-	-
Schulke & Mayr GmbH	thermosept® DK	-	✓	-	✓	✓	✓	-	-
Schulke & Mayr GmbH	thermosept® ED	-	✓	-	✓	✓	✓	✓	✓
Schulke & Mayr GmbH	thermosept® ER	✓	-	-	✓	✓	✓	✓	✓
Schulke & Mayr GmbH	thermosept® PAA base thermosept® PAA additive	-	✓	-	✓	✓	✓	✓	✓
Schulke & Mayr GmbH	thermosept® RKF forte	✓	-	-	✓	✓	✓	-	-
Schulke & Mayr GmbH	thermosept® RKN- zym	✓	-	-	✓	✓	✓	-	-
Schulke & Mayr GmbH	thermosept® X-tra	✓	-	-	✓	✓	-	-	-
STERIS Corporation	Instru-Klenz	✓	-	✓	✓	✓	✓	-	-
Weigert, Dr. GmbH & Co.	neodisher endo® CLEAN	✓	-	✓	✓	✓	✓	✓	✓
Weigert, Dr. GmbH & Co.	neodisher endo® SEPT GA	-	✓	-	✓	✓	✓	✓	✓
Weigert, Dr. GmbH & Co.	neodisher endo® SEPT PAC	-	✓	-	-	-	✓	✓	✓
Weigert, Dr. GmbH & Co.	neodisher® DuoClean	✓	-	-	✓	✓	-	-	-

Weigert, Dr. GmbH & Co.	neodisher® FA	✓	-	-	✓	✓	-	-	-
Weigert, Dr. GmbH & Co.	neodisher® MediClean	✓	-	✓	✓	✓	✓	✓	✓
Weigert, Dr. GmbH & Co.	neodisher® MediClean forte	✓	-	✓	✓	✓	✓	✓	✓
Weigert, Dr. GmbH & Co.	neodisher® MediZym	✓	-	✓	✓	✓	✓	✓	✓
Weigert, Dr. GmbH & Co.	neodisher® OXIVARIO	✓	-	-	✓	✓	-	-	-
Weigert, Dr. GmbH & Co.	neodisher® SeptoClean*	✓	✓	✓	✓	✓	-	-	-
Weigert, Dr. GmbH & Co.	neodisher® Septo DN	-	✓	-	✓	✓	✓	✓	✓

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✓ Suitable chemical agents

- Chemicals not tested

Cleaning agents that are marked with*:

These cleaning agents are proven to be effective against prions

① **Note:**

If agents other than those listed above have been used, any arising damage is not covered by KARL STORZ.

Here with KARL STORZ SE & Co. KG confirms the material compatibility of the following chemicals for processing endoscopes and endoscopic instruments.

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