DECLARATION

I, Futoshi Inoue, President and Representative Director of Kitazato Corporation, do hereby solemnly and sincerely declare that the attached document herewith is the original of "Instruction for Use of Medical Device".

Futoshi INOUE
President and Representative Director
Kitazato Corporation

Date: 18th day of December, 2019


KITAZATO®

Instruction for Use
for Medical Device

Vitrification solutions
for oocytes and embryos, kit

Manufactured by Kitazato Corporation

1. Name of Medical Device

‘Vitrification Solutions for oocytes and embryos, kit’:

1. Basic Solution (BS): 1 vial 1,5 mL;
2. Equilibration Solution (ES): 1 vial 1,5 mL;
3. Vitrification Solution (VS): 2 vials 1,5 mL.

Manufacturer: Kitazato Corporation, 81 Nakajima, Fuji-city, Shizuoka 416-0907, Japan.

Place of manufacturing: Kitazato Corporation Shizuoka Office, 81 Nakajima, Fuji-city, Shizuoka 416-0907, Japan.

2. Use of Medical Device

‘Vitrification Solutions for oocytes and embryos, kit’ are used in Assisted Reproductive Technologies (ART) for vitrification of oocytes and embryos.

3. Potential customers

‘Vitrification Solutions for oocytes and embryos, kit’ are intended to be used by medical specialists trained in fertility treatment.

4. Functional Characteristics. Description of Device Operation

‘Vitrification Solutions for oocytes and embryos, kit’ ensure safe process of vitrification of good quality reproductive cells for indefinite period of time.

The Solutions are universal and suitable for vitrification of oocytes and embryos at all stages of development.

Due to the concentration of cryoprotectant agents in the content of Solutions it is possible to avoid ice formation inside the cells during the process of vitrification in liquid nitrogen at -196°C and thus allow cytoplasm to be in gel state for further recovery of cell functions.

At present, the most widely used CPAs are composed of agents that penetrate the cell (such as ethylene glycol (EG), ProH or glycerol) as well as low-molecular-weight non-penetrating cryoprotectants (such as sucrose and trehalose). The concentration of cryoprotectants in the commercially marketed vitrification solutions, including ‘Vitrification Solutions for oocytes and embryos, kit’ is very low. In order to ensure safe vitrification process for cells it is very important to strictly follow the protocol and time tracking, indicated by manufacture in present IFU.

The low temperature during vitrification process makes it possible to stop the biochemical processes of the cells and leave them for storage for a long period of time with a further possibility of their use for infertility treatment. At the same time, the shelf life of the solutions,

set by the manufacturer, means the period during which the solutions are stable under normal atmospheric conditions, and there is no growth of microorganisms. Oocytes and embryos that have been subjected to freezing using vitrification solutions can be stored indefinitely in liquid nitrogen, where there is no likelihood of the growth of microorganisms. After thawing, the cells are at the same stage of development at which they were frozen, because during storage of oocytes and embryos with stopped biochemical life processes in liquid nitrogen their biological growth and aging are absent, and, accordingly, the solutions remain sterile, as there is no microbial growth.

During vitrification procedure the Solutions are used in sequence. Each of them has specially designed content required for the reproductive cells at every stage of procedure in order to assure their safety and maximize survival rate.

Name of Solution	Function
Basic Solution (BS)	For easing equilibration of Oocytes by the ES. First of all, the Oocytes are equilibrated by the BS before the ES is added. The BS is not used for embryos.
Equilibration Solution (ES)	ES contains Ethylene Glycol that is cryoprotectant. Concentration of Ethylene Glycol contained in ES is half of concentration of those reagents contained in VS. The ES is used to equilibrate the Oocyte or Embryo with concentration of cryoprotectants of the ES before soaking the Oocyte or Embryo with cryoprotectants of the VS
Vitrification Solution (VS)	VS contains Trehalose and Ethylene Glycol that are cryoprotectants. The VS is intended to soak the Oocyte or the Embryo with concentration of cryoprotectants that can vitrify the Oocyte or the Embryo after the Oocyte or Embryo was soaked in ES.


Kitazato is recognized as one of the pioneering brands in driving and improving vitrification. Its greatest contribution in this field has been the development of the vitrification and thawing system which has become the global leader in vitrification of oocytes and embryos, in all stages of development.

Due to the minimum volume cooling, it is possible to realize highest cooling & warming rates resulting in over 90% post-thaw survival.

After over a decade on the market, Kitazato Vitrification products has been applied in over 1,500,000 clinical cases in over 90 countries and 2,200 assisted reproduction centers. Hundreds of scientific publications certify their excellent results.

5. Indications for Use

‘Vitrification Solutions for oocytes and embryos, kit’ are used for vitrification of oocytes and

embryos and provide the following advantages for the patients:

- Provide an opportunity to limit the number of embryos transferred while maximizing the usable embryo per oocyte retrieval ratio at each stimulation attempt, which is invading treatment with physically stressful and also costly for patients.
- For fertility preservation which improve their potential to achieve pregnancy at later maternal age, by cryopreservation of oocytes retrieved at younger age or before patients start oncology treatment.
- Maximize opportunity for pregnancy by oocyte donation cycle, pre-implantation genetic diagnosis or by advancing embryo transfer until the endometrium recovers to increase embryo respectability.

6. Contraindications


- Do not re-sterilize.
- To avoid contamination, do not re-use.
- Do not use solution that shows cloudiness or has turned yellow.
- Do not use if any problems are observed in the specification of the container or labeling (number, color, name, volume).
- Sterile if the package is unopened or undamaged. Do not use if the package is opened or damaged.
- The media must be received in the original unopened container and refrigerated at 2 to 8 °C.
- Do not use if the media were not stored under refrigeration (2 to 8 °C).
- Do not use the media with expired shelf life.
- Do not use for injection.

7. Potential Side Effects

There are no side effects.

8. Cautions for Use

- Read the instructions for use prior to use.
- Use the vitrification solutions to vitrify the oocytes and embryos.
- This product is intended to be used by medical specialists trained in fertility treatment.
- Aseptic technique should be used.
- Use sterilized equipment and materials only.
- Decontaminate the workroom for aseptic procedures.
- Follow procedures in an environmentally controlled room.
- Avoid eye contact or contact with skin. In case of eye or skin contact immediately wash

thoroughly with water.

- Devices and disposables used together with solutions should be intended for such use. Before use familiarize yourself with the safety instructions for the devices used together.
- Observe all federal, state and local environmental regulations when discarding the product.
- Request a qualified waste disposal contractor to dispose the products as residual and non-reusable solutions.
- Dispose the contaminated containers and packages in the same way as the containers with the product inside.
- In case of infection, dispose the product appropriately in a prescribed manner.
- The user shall be responsible for any problems caused by non-conformity to the present IFU.

All warnings and precautions necessary for the use of this device are included in the documentation accompanying the device.

9. Transportation and Storage Conditions

Transportation of 'Vitrification Solutions for oocytes and embryos, kit' should be made at +2 – +8°C.

Transportation from the manufacturing site to a customer is performed in the cold chain. For delivery to the customer, the carton boxes with the solutions are taken from the cooling chamber and are packed in insulated shipping containers from styrene with cold packs.

All the deliveries for long distances are made by air cargo and are tracked in order to avoid delays. Also the delivery by road is allowed when is made by auto truck equipped with the refrigerators adjusted for the required temperature regulation range.

During transportation it is necessary to avoid high temperature, high humidity, direct sunlight and water wet.


Upon receipt the solutions should be stored in cooling chamber or refrigerator at +2 – +8°C until the expiry date indicated on the package avoiding direct sunlight. Do not freeze.

10. Terms of Use

The use of 'Vitrification Solutions for oocytes and embryos, kit' is allowed only in the sterile laboratories with observance of the aseptic techniques. The device is not allowed for at-home use.

The procedures should be performed in the area with appropriate ventilation and air conditioning.

Temperature and % of CO₂ in the incubators should be controlled on the permanent basis with

the independent measurement devices and should be compared with the indications on their displays.

For all the manipulations with oocytes and embryos only sterile non-toxic disposables should be used in order to avoid a risk of bacterial inoculation and contamination.

During vitrification procedure the solutions are used together with the pipettes (for dosing the solutions, together with the pipettors), tips (for manipulation with cells, together with the holders), repro plates (it is recommended to use Repro Plates with six wells) as well as containers (straws) for vitrification and storage of oocytes and embryos (it is recommended to Cryotop and put up to 4 oocytes or embryos on one straw). The vitrification process is performed in cryo box filled with liquid nitrogen.


11. Safety Requirements

'Vitrification Solutions for oocytes and embryos, kit' should be used only by medical specialist of fertility treatment with observance of the aseptic techniques.

'Vitrification Solutions for oocytes and embryos, kit' should not be deliberately used for purposes other than intended. Misuse or incorrect operation may cause cell tissue injury.

During use of 'Vitrification Solutions for oocytes and embryo, kit' the strict compliance with the terms of Instruction for Use is essential.


12. Content and Characteristics of Solutions

Contains:

Component Name	Volume	Quantity	Cap Color of Tube
Basic Solution (BS)*	1.5 mL	1	White
Equilibration Solution (ES)	1.5 mL	1	Blue
Vitrification Solution (VS)	1.5 mL	2	Green

* Used only for oocyte vitrification.

Composition

Basic Solution (BS)

Ingredient of Solution	Purpose
Hydroxypropyl Cellulose	Substitution of Blood Albumin
HEPES sodium salt	Buffer
HEPES	Buffer
Sodium Bicarbonate	pH adjuster
Polyvinylpyrrolidone	Thickener
Medium	Basal Medium
Deionized, Distilled Water	Solvent

Equilibration Solution (ES)

Ingredient of Solution	Purpose
Ethylen Glycol	Cryoprotectant
Hydroxypropyl Cellulose	Substitution of Blood Albumin
HEPES sodium salt	Buffer
HEPES	Buffer
Sodium Bicarbonate	pH adjuster
Polyvinylpyrrolidone	Thickener
Medium	Basal Medium
Deionized, Distilled Water	Solvent

Vitrification Solution (VS)

Ingredient of Solution	Purpose
Trehalose	Osmotic adjuster
Ethylen Glycol	Cryoprotectant
Hydroxypropyl Cellulose	Substitution of Blood Albumin
HEPES sodium salt	Buffer
HEPES	Buffer
Sodium Bicarbonate	pH adjuster
Polyvinylpyrrolidone	Thickener
Medium	Basal Medium
Deionized, Distilled Water	Solvent

Specification

Item	Specification
Appearance (Color)	Clear, Particle Free /Rose Pink
pH	7.20 – 7.60
Osmolality	Basic Solution (BS): 240-300 mmol/kg Equilibration Solution (ES): 2300-2800 mmol/kg


	Vitrification Solution (VS): 4900-6000 mmol/kg
Sterility	No microbial growth
Endotoxin	≤ 0.25 EU/mL
Mouse Embryo Assay	One cell mouse embryo testing is performed using fresh mouse embryos. Test results indicate the percentage of mouse embryos developing fully expanded blastocyst after 96 hours exposing in the sample. Greater than or equal to 80 %.

Specification for each batch of the Solutions can be found in the Certificate of Analysis for the exact LOT.

All the solutions are intended for single use only. The shelf life of the Solutions is 3 month from date of issue of Certificate of Analysis.

One kit of solutions enables to vitrify up to 16 oocytes or embryos.

13. Instruction for Use


Preparation

1. Fill 90% of Cooling Rack with fresh liquid nitrogen.
2. Compare the width of perivitelline space with the thickness of zona pellucida and record it.
NOTE: Use a Pasteur pipette that has a suitable internal diameter for Oocyte (External diameter: 120 µm) or Embryo.
CAUTION: Equilibration procedures for Oocytes and Embryo are different.

Oocyte Equilibration

1. Drop 20 µl for BS and 300 µl each for VS1 и VS2 on the Repro Plate with micropipette.
2. Transfer the Oocyte from the culture dish to the **BOTTOM** of BS.
3. Step 1. Add ES 20 µl gently to the **TOP** of BS and leave it for 3 minutes.
Step 2. Add another ES 20 µl gently to the **TOP** of BS and leave it for 3 minutes.
Step 3. Add another ES 240 µl gently to the **TOP** of BS and leave it for 9 минут.

Embryo Equilibration

1. Drop 300 µl for **ES**, **VS1** and **VS2** on the Repro Plate using micropipette.
2. Transfer the Embryo to the **TOP** center of **ES**.
3. It spontaneously begins to shrink and then gradually returns to its original size by **ES** infiltration (within 15 minutes).

Vitrification

CAUTION: The following steps from 4 to 9 should be completed between 60 and 90 seconds.

4. Aspirate the Oocyte (Embryo) in **ES** at the tip of pasteur pipette. Transfer the Oocyte (Embryo) to the **TOP** center of **VS1**.
5. Aspirate the Oocyte (Embryo) with Pasteur pipette and blow it out. Repeat this work 3times changing the position in **VS1**.
6. Transfer the Oocyte (Embryo) to **VS2**. Move the Oocyte (Embryo) changing positions twice in **VS2** with Pasteur pipette.
7. Place the Oocyte (Embryo) by the black part of Cryotop sheet. Make a planar droplet.
8. Check if the Oocyte (Embryo) is on the sheet with minimal volume of **VS2** (less than 0.1 µl) under the microscope.
9. Plunge the Cryotop into liquid nitrogen quickly.
10. Put the Cryotop in a cane and store it in a storage tank.

For complete instruction for the vitrification procedure of oocytes and embryos, please refer to Vitrification Protocol for Cryotop Method.

13. Compliance with International Standards

For manufacture of 'Vitrification Solutions for oocytes and embryos, kit' a Quality Management System which complies with the requirements of ISO 13485:2016 is operated.

'Vitrification Solutions for oocytes and embryos, kit' meet the provisions of the MDD 93/42/ECC with the amendments as of Directive 2007/47/EC.

14. Sterilization

Hermetically closed and unbroken packaging assures sterility. Do not use solutions if sterile packaging is broken.

'Vitrification Solutions for oocytes and embryos, kit' are manufactured and filled into sterile vials in aseptic conditions according to EN ISO 13408-1. For the production of 'Vitrification Solutions for oocytes and embryos, kit' the sterile raw materials and components are used, and the filtration of components into final product and the filling into vials are performed in a LAF-unit Class A, located in Class B according to Good Manufacturing Practice (GMP). The solutions are sterilized by sterile filtration. The sterile fill of the vials is assured by media fill test.

Vials and caps are sterilized by gamma irradiation. The sterilization process meets the terms of standard ISO 11137 according to which Sterility Assurance Level (SAL) 10^{-6} is adequately assured

at the sterilization dose 15.0 kGy.

The tightly closed and sealed package assures sterility of the product. Do not use the solutions if sterile packaging is broken.

15. Stability and Shelf Life

In case of transportation and storage in accordance with the instructions of the manufacturer 'Vitrification Solutions for oocytes and embryos, kit' are stable up to their expiry date, indicated on the labels attached to the vials and carton boxes.

After warming the rest of the solutions should be discarded.

After opening it is allowed to use and store the solutions in cooling chamber or refrigerator during 2 days, after that they should be disposed.

16. Environmental Protection Requirements

'Vitrification Solutions for oocytes and embryos, kit' during use, transportation and storage do not have an adverse impact on human or environment.

17. Disposal

Used solutions should be disposed in compliance with the requirements of national legislation in regard to disposal and destruction of medical waste.

The rest of solutions and empty vials are recommended to be disposed in single use water resistant containers for dangerous medical waste with adherent cap, which avoids inadvertent opening of the container.

It is necessary to contact a qualified waste disposal contractor to dispose the solutions as residual and unsuitable for reuse.

18. Packaging and Labelling

'Vitrification Solutions for oocytes and embryos, kit' are supplied in sterile vials with screw caps made of polypropylene (primary packaging). Cap of each vial has silicone rubber pad for leak resistance.

Each solution has individual color-coding of cap. For convenience of user each cap also has letter identification of solution and number of step according to Vitrification Protocol (See Vitrification Protocol for Cryotop Method).

For filling and sealing of vials with solutions the aseptic process is used according to EN ISO 13408-1. Every time prior to the production process the sealing machine is performed is inspected under startup and also after change of the settings and sizes of the vials. Adjusted and inspected on the permanent basis sealing process allows to provide firm adherence of the cap, avoid leakage and

assure antimicrobial barrier.

For supply to the end-users the vials with the solutions are packed into the carton boxes (secondary packaging) which protect the solutions during storage and transportation.

The vials and carton boxes with the solutions are supplied with the labels. The labels contain the following information: product name, intended use, content and volume of vials, reference, batch code (LOT), product code, expiry date, single use product indication, sterility indication and sterilization method, prohibition of use in case of damaged package and re-sterilization, storage conditions as well as name and address of the manufacture and EU representative, CE mark.

The labels on the vials due to their small size contain only partial information on the product: letter identification of the solution, batch code (LOT), expiry date, single use product indication, storage conditions and number of step according to Vitrification Protocol (See Vitrification Protocol for Cryotop Method).

Symbols on Labels

	Sterilized using Aseptic Processing Techniques
	Do not reuse
	Catalogue number
	Batch Code (LOT)
	Storage Temperature
	Use by
	Consult instructions for use
	Do not use if the package is damaged
	Do not re-sterilize
	Manufacturer
	EC Representative
	CE mark

According to the requirements of national legislation in terms of labelling of medical devices the labels in Russian language are placed on the packaging of the product.

19. Contents of Delivery

'Vitrification Solutions for oocytes and embryos, kit' are supplied in kits of 4 vials of 1.5 mL packed into a carton box and accompanied by Instruction for Use.

20. Warranty

The manufacture guarantees the quality of device during its shelf life when stored and used properly.

The manufacture is not liable for any damage resulting from misuse or incorrect operation or failure to follow the cautions or instructions contained in the instruction for use or the instruction manual for any compatible medical equipment.

21. Reclamation

In case of claims concerning the quality of medical device, please contact Manufacture and / or its Authorized Representative.

Manufacturer:


Kitazato Corporation
81 Nakajima, Fuji-city, Shizuoka 416-0907, Japan
Tel.: +81-545-66-2202 / Fax: +81-545-60-5772
E-mail: info@kitazato.co.jp

Place of manufacturing:

Kitazato Corporation Shizuoka Office
81 Nakajima, Fuji-city, Shizuoka 416-0907, Japan
Tel.: +81-545-66-2202 / Fax: +81-545-60-5772
E-mail: info@kitazato.co.jp

Authorized Representative in Russian Federation:

EkoZdrav Ltd.
6, Yartsevskaya St., 1st floor, premises 2, room 11 Moscow 121351, Russia
Tel.: +7 (916) 640-24-44
E-mail: info@ekozdrav.ru

Futoshi Inoue
President and Representative Director
Kitazato Corporation
Date: 18th day of December, 2019


Registered No.'20-0006

NOTARIAL CERTIFICATE

THIS IS TO CERTIFY:

That **Mr. Hisashi HAYASHI**, a duly authorized agent of **Mr. Futoshi INOUE**, President and Representative Director of **Kitazato Corporation**, being legally constituted and existing in accordance with the laws of Japan, has stated in my presence that the said **Mr. Futoshi INOUE** has acknowledged himself that the signature on the attached document is his own.

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This public document

2. has been signed by Kikuo ASAKA

3. acting in the capacity of Notary of the Tokyo Legal Affairs Bureau

4. bears the seal/stamp of Kikuo ASAKA, Notary

Certified

5. at Tokyo

6. JAN. - 6. 2020

7. by the Ministry of Foreign Affairs

8. 20- № 003004

9. Seal/stamp:

10. Signature

T. TANAKA

Toshie TANAKA

For the Minister for Foreign Affairs

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