
DECLARATION

I, Futoshi Inoue, President and Representative Director of Kitazato Corporation, do hereby solemnly and sincerely declare that the attached document herewith is the original of "Instruction for Use of Medical Device".

Futoshi INOUE

President and Representative Director
Kitazato Corporation

Date: 18th day of December, 2019


KITAZATO®

Instruction for Use
for Medical Device

**Thawing solutions
for oocytes and embryos, kits»**

Manufactured by Kitazato Corporation

1. Name of Medical Device

'Thawing Solutions for oocytes and embryos, kit'

1. Thawing Solution (TS): 2 vials 4,0 mL;
2. Diluent Solution (DS): 1 vial 4,0 mL;
3. Washing Solution (WS): 1 vial 4,0 mL.

Manufacturer: Kitazato Corporation, 81 Nakajima, Fuji-city, Shizuoka 416-0907, Japan.

Place of manufacturing: Kitazato Corporation Shizuoka Office, 81 Nakajima, Fuji-city, Shizuoka 416-0907, Japan.

2. Use of Medical Device

'Thawing Solutions for oocytes and embryos, kit' are used in Assisted Reproductive Technologies (ART) for thawing of oocytes and embryos.

3. Potential customers

'Thawing Solutions for oocytes and embryos, kit' are intended to be used by medical specialists trained in fertility treatment.


4. Functional Characteristics. Description of Device Operation

'Thawing Solutions for oocytes and embryos, kit' ensure safe process of thawing of vitrified oocytes and embryos.

The Solutions are universal and suitable for thawing of oocytes and embryos at all stages of development.

During thawing procedures the Solutions are used in sequence. Each of them has specially designed content required for the reproductive cells at every stage of procedure in order to assure their safety and maximize survival rate after thawing.

'Solutions for Thawing of oocytes and embryos, kit' allow to remove cryoprotectants from cells. 'Solutions for Thawing of oocytes and embryos, kit' contain a lower concentration of trehalose compared to vitrification solutions. Trehalose is one of the most widely used and proven non-penetrating cryoprotectants. By gradually reducing the concentration of trehalose, thawing solutions prevent excessive cell expansion during the thawing process and help to restore normal osmotic pressure in them.

In order to ensure safe thawing process for cells it is very important to follow strictly the protocol and time tracking, indicated by manufacture in the present IFU.

Name of Solution	Function
Thawing Solution (TS)	TS contains Trehalose that regulates osmotic pressure. The TS is intended to thaw the cryopreserved Oocyte or Embryo.
Diluent Solution (DS)	DS contains Trehalose that regulates osmotic pressure. The DS is intended to dilute the thawed Oocyte or Embryo.
Washing Solution (WS)	WS is intended to wash the diluted Oocyte or Embryo. The composition of the WS is same as the composition of the BS.

Kitazato is recognized as one of the pioneering brands in driving and improving vitrification. Its greatest contribution in this field has been the development of the vitrification and thawing system which has become the global leader in vitrification of oocytes and embryos, in all stages of development.

Due to the minimum volume cooling it is possible to realize highest cooling & warming rates resulting in over 90% post-thaw survival.

After over a decade on the market, Kitazato Vitrification and Thawing products has been applied in over 1,500,000 clinical cases in over 90 countries and 2,200 assisted reproduction centers. Hundreds of scientific publications certify their excellent results.


5. Indications for Use

'Thawing Solutions for oocytes and embryos, kit' are used for warming of oocytes and embryos provide the following advantages for the patients;

- provides an opportunity to limit the number of embryos transferred while maximizing the usable embryo per oocyte retrieval ratio at each stimulation attempt, which is invading treatment with physically stressful and also costly for patients.
- For fertility preservation which improve their potential to achieve pregnancy at later maternal age, by cryopreservation of oocytes retrieved at younger age or before patients start oncology treatment.
- Maximize opportunity for pregnancy by oocyte donation cycle, pre-implantation genetic diagnosis or by advancing embryo transfer until the endometrium recovers to increase embryo respectability.

6. Contraindications

- Do not re-sterilize.
- To avoid contamination, do not re-use.

- Do not use solution that shows cloudiness or has turned yellow.
- Do not use if any problems are observed in the specification of the container or labeling (number, color, name, volume).
- Sterile if the package is unopened or undamaged. Do not use if the package is opened or damaged.
- The media must be received in the original unopened container and refrigerated at 2 to 8 °C.
- Do not use if the media were not stored under refrigeration (2 to 8 °C).
- Do not use the media with expired shelf life.
- Do not use for injection.

7. Potential Side Effects


There are no side effects.

8. Cautions for Use

The following warnings and precautions are included on the label. Warning and precautions include symbols to communicate to the user:

WARNINGS

- Read the instructions for use prior to use.
- Use the thawing solutions to thaw the oocytes or embryos vitrified with vitrification solutions.
- This product is intended to be used by medical specialists trained in fertility treatment.
- Aseptic technique should be used.
- Use sterilized equipment and materials only.
- Decontaminate the workroom for aseptic procedures.
- Follow procedures in an environmentally controlled room.
- Avoid eye contact or contact with skin. In case of eye or skin contact immediately wash thoroughly with water.
- Devices and disposables used together with solutions should be intended for such use. Before use familiarize yourself with the safety instructions for the devices used together.
-
- Observe all federal, state and local environmental regulations when discarding the product.
- Request a qualified waste disposal contractor to dispose the products as residual and non-reusable solutions.

- Dispose the contaminated containers and packages in the same way as the containers with the product inside.
- In case of infection, dispose the product appropriately in a prescribed manner.
- The user shall be responsible for any problems caused by non-conformity to the present IFU.

All warnings and precautions necessary for the use of this device are included in the documentation accompanying the device.

9. Transportation and Storage Conditions

Transportation of 'Thawing Solutions for oocytes and embryos, kit' should be made at +2 – +8°C.

Transportation from the manufacturing site to a customer is performed in the cold chain. For delivery to the customer, the carton boxes with the solutions are taken from the cooling chamber and are packed in insulated shipping containers from styrene with cold packs.

All the deliveries for long distances are made by air cargo and are tracked in order to avoid delays. Also the delivery by road is allowed when is made by auto truck equipped with the refrigerators adjusted for the required temperature regulation range.

During transportation it is necessary to avoid high temperature, high humidity, direct sunlight and water wet.

Upon receipt the solutions should be stored in cooling chamber or refrigerator at +2 – +8°C until the expiry date indicated on the package avoiding direct sunlight. Do not freeze.

10. Terms of Use


The use of 'Thawing Solutions for oocytes and embryos, kit' is allowed only in the sterile laboratories with observance of the aseptic techniques. The device is not allowed for at-home use.

The procedures should be performed in the area with appropriate ventilation and air conditioning.

Temperature and % of CO₂ in the incubators should be controlled on the permanent basis with the independent measurement devices and should be compared with the indications on their displays.

For all the manipulations with oocytes and embryos only sterile non-toxic disposables should be used in order to avoid a risk of bacterial inoculation and contamination.

During thawing procedure the solutions are used together with the pipettes (for dosing the solutions, together with the pipettors), tips (for manipulation with cells, together with the holders), Petri dishes 35 mm and repro plates (it is recommended to use Repro Plates with six wells).

11. Safety Requirements

'Thawing Solutions for oocytes and embryos, kit' should be used only by medical specialist of fertility treatment with observance of the aseptic techniques.

'Thawing Solutions for oocytes and embryos, kit' should not be deliberately used for purposes other than intended. Misuse or incorrect operation may cause cell tissue injury.

During use of Vitrification and Thawing Solutions for oocytes and embryo, kits the strict compliance with the terms of Instruction for Use is essential.

12. Content and Characteristics of Solutions


Contains:

Component Name	Volume	Quantity	Cap Color of Tube
Thawing Solution (TS)	4.0 mL	2	Red
Diluent Solution (DS)	4.0 mL	1	Yellow
Washing Solution (WS)	4.0 mL	1	White

Composition

Thawing Solution (TS)

Ingredient of Solution	Purpose
Trehalose	Osmotic adjuster
Hydroxypropyl Cellulose	Substitution of Blood Albumin
HEPES sodium salt	Buffer
HEPES	Buffer
Sodium Bicarbonate	pH adjuster
Polyvinylpyrrolidone	Thickener
Medium*	Basal Medium
Deionized, Distilled Water	Solvent

**Diluent Solution (DS)**

Ingredient of Solution	Purpose
Trehalose	Osmotic adjuster
Hydroxypropyl Cellulose	Substitution of Blood Albumin
HEPES sodium salt	Buffer
HEPES	Buffer
Sodium Bicarbonate	pH adjuster
Polyvinylpyrrolidone	Thickener
Medium*	Basal Medium
Deionized, Distilled Water	Solvent

Washing Solution (WS)

Ingredient of Solution	Purpose
Hydroxypropyl Cellulose	Substitution of Blood Albumin
HEPES sodium salt	Buffer
HEPES	Buffer
Sodium Bicarbonate	pH adjuster
Polyvinylpyrrolidone	Thickener
Medium*	Basal Medium
Deionized, Distilled Water	Solvent

Specification


Item	Specification
Appearance (Color)	Clear, Particle Free /Rose Pink
pH	7.20 – 7.60
Osmolality	Thawing Solution (TS): 1600-2000 mM/kg Diluent Solution (DS): 830-1020 mM/kg Washing Solution (WS): 240-300 mM/kg
Sterility	No microbial growth
Endotoxin	≤ 0.25 EU/mL
Mouse Embryo Assay	One cell mouse embryo testing is performed using fresh mouse embryos. Test results indicate the percentage of mouse embryos developing fully expanded blastocyst after 96 hours exposing in the sample. Greater than or equal to 80 %.

Specification for each batch of the Solutions can be found in the Certificate of Analysis for the exact LOT.

All the solutions are intended for single use only. The shelf life of the Solutions is 3 month from the date of issue of certificate of analysis.

One kit of solutions enables to thaw up to 16 oocytes or embryos.

13. Instruction for Use

Preparation

1. Warm **TS** vial with a cap and a Petri Dish in an incubator to 37 °C. Pour the full content of **TS** into the Petri Dish.
2. Drop 300 µl each of **DS**, **WS1** and **WS2** on the Repro Plate with micropipette.
NOTE: Use pasteur pipette that has a suitable internal diameter for Oocyte (External diameter: 120 µm) or Embryo.

Thawing

1. Quickly immerse the Cryotop strip completely into the **TS**. Leave it for 1 minute.
2. Aspirate the Oocyte (Embryo) with the pasteur pipette and gently place it on the **BOTTOM** of the **DS**. Leave it for 3 minutes.
3. Aspirate the Oocyte (Embryo) with the pasteur pipette and gently place it on the **BOTTOM** of the **WS1**. Leave it for 5 minutes.
- 4.
5. Aspirate the Oocyte (Embryo) with the pasteur pipette and gently place it on the **TOP** of the **WS2**. After the Oocyte (Embryo) drops to the bottom of the, repeat this process twice.
6. Transfer the Oocyte (Embryo) to a culture dish containing the appropriate culture medium. Incubate the Oocyte (Embryo) in a 37 °C incubator to complete recovery.
NOTE: Oocytes for 2 hours and Embryos for 3 hours as an idea.

For complete instruction for the thawing procedure of oocytes and embryos, please refer to Vitrification Protocol for Cryotop Method attached as a separate file.

13. Compliance with International Standards

'Thawing Solutions for oocytes and embryos, kit' meet the provisions of the MDD 93/42/ECC with the amendments as of Directive 2007/47/EC.

Hermetically closed and unbroken packaging assures sterility. Do not use solutions if sterile packaging is broken.

'Thawing Solutions for oocytes and embryos, kit' are manufactured and filled into sterile vials in aseptic conditions according to EN ISO 13408-1. For the production of 'Thawing Solutions for oocytes and embryos, kit' the sterile raw materials and components are used, and the filtration of components into final product and the filling into vials are performed in a LAF-unit Class A, located in Class B according to Good Manufacturing Practice (GMP). The solutions are sterilized by sterile filtration. The sterile fill of the vials is assured by media fill test.

Vials and caps are sterilized by gamma irradiation. The sterilization process meets the terms of standard ISO 11137 according to which Sterility Assurance Level (SAL) 10^{-6} is adequately assured at the sterilization dose 15.0 kGy.

The tightly closed and sealed package assures sterility of the product. Do not use the solutions if sterile packaging is broken.

In case of transportation and storage in accordance with the instructions of the manufacturer 'Thawing Solutions for oocytes and embryos, kit' are stable up to their expiry date, indicated on the labels attached to the vials and carton boxes.

After warming the rest of the solutions should be discarded.

After opening it is allowed to use and store the solutions in cooling chamber or refrigerator during 2 days, after that they should be disposed.

'Thawing Solutions for oocytes and embryos, kit' during use, transportation and storage do not have an adverse impact on human or environment.

Used solutions should be disposed in compliance with the requirements of national legislation in regard to disposal and destruction of medical waste.

The rest of solutions and empty vials are recommended to be disposed in single use water resistant containers for dangerous medical waste with adherent cap, which avoids inadvertent opening of the container.

It is necessary to contact a qualified waste disposal contractor to dispose the solutions as residual and unsuitable for reuse.

18. Packaging and Labelling

'Thawing Solutions for oocytes and embryos, kit' are supplied in sterile vials with screw caps made of polypropylene (primary packaging). Cap of each vial has silicone rubber pad for leak resistance.

Each solution has individual color-coding of cap. For convenience of user each cap also has letter identification of solution and number of step according to Thawing Protocol (See Vitrification Protocol for Cryotop Method).

For filling and sealing of vials with solutions the aseptic process is used according to EN ISO 13408-1. Every time prior to the production process the sealing machine is performed is inspected under startup and also after change of the settings and sizes of the vials. Adjusted and inspected on the permanent basis sealing process allows to provide firm adherence of the cap, avoid leakage and assure antimicrobial barrier.


For supply to the end-users the vials with the solutions are packed into the carton boxes (secondary packaging) which protect the solutions during storage and transportation.

The vials and carton boxes with the solutions are supplied with the labels. . The labels contain the following information: product name, intended use, content and volume of vials, reference, batch code (LOT), product code, expiry date, single use product indication, sterility indication and sterilization method, prohibition of use in case of damaged package and re-sterilization, storage conditions as well as name and address of the manufacture and EU representative, CE mark.

The labels on the vials due to their small size contain only partial information on the product: letter identification of the solution, batch code (LOT), expiry date, single use product indication, storage conditions and number of step according to Thawing Protocol (See Vitrification Protocol for Cryotop Method).

Symbols on Labels

STERILE A	Sterilized using Aseptic Processing Techniques
	Do not reuse
REF	Catalogue number

	Batch Code (LOT)
	Storage Temperature
	Use by
	Consult instructions for use
	Do not use if the package is damaged
	Do not re-sterilize
	Manufacturer
	EC Representative
	CE mark

According to the requirements of national legislation in terms of labelling of medical devices the labels in Russian language are placed on the packaging of the product.

19. Contents of Delivery

'Thawing Solutions for oocytes and embryos' are supplied in kits of 4 vials of 1.5 mL / 4.0 mL packed into a carton box and accompanied by Instruction for Use.

20. Warranty

The manufacture guarantees the quality of device during its shelf life when stored and used properly.

The manufacture is not liable for any damage resulting from misuse or incorrect operation or failure to follow the cautions or instructions contained in the instruction for use or the instruction manual for any compatible medical equipment.

21. Reclamation

In case of claims concerning the quality of medical device, please contact Manufacture and / or its Authorized Representative.

Manufacturer:

Kitazato Corporation
81 Nakajima, Fuji-city, Shizuoka 416-0907, Japan
Tel.: +81-545-66-2202 / Fax: +81-545-60-5772
E-mail: info@kitazato.co.jp

Place of manufacturing:

Kitazato Corporation Shizuoka Office
81 Nakajima, Fuji-city, Shizuoka 416-0907, Japan
Tel.: +81-545-66-2202 / Fax: +81-545-60-5772
E-mail: info@kitazato.co.jp

Authorized Representative in Russian Federation:

EkoZdrav Ltd.
6, Yartsevskaya St., 1st floor, premises 2, room 11 Moscow 121351, Russia
Tel.: +7 (916) 640-24-44
E-mail: info@ekozdrav.ru

Futoshi Inoue
President and Representative Director
Kitazato Corporation
Date: 18th day of December, 2019

Registered No. '20-0007

NOTARIAL CERTIFICATE

THIS IS TO CERTIFY:

That **Mr. Hisashi HAYASHI**, a duly authorized agent of **Mr. Futoshi INOUE**, President and Representative Director of **Kitazato Corporation**, being legally constituted and existing in accordance with the laws of Japan, has stated in my presence that the said **Mr. Futoshi INOUE** has acknowledged himself that the signature on the attached document is his own.

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


This public document

2. has been signed by Kikuo ASAKA

3. acting in the capacity of Notary of the Tokyo Legal Affairs Bureau

4. bears the seal/stamp of Kikuo ASAKA, Notary

Certified

5. at Tokyo

6. JAN. - 6. 2020

7. by the Ministry of Foreign Affairs

8. 20- № 003005

9. Seal/stamp:

10. Signature

T. TANAKA

Toshie TANAKA

For the Minister for Foreign Affairs

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