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APPROVED BY
REPROLIFE, Inc., Japan

Signature, seal/stamp, date
2021.3.19.

Instructions for use for the medical device “Manipulation plate for assisted reproductive technologies, in variants”

Prepared in accordance with the requirements of the Recommendations of International Medical Devices Regulations Forum (IMDRF) and the Order of the Ministry of Health of the Russian Federation No. 11n dated January 19th, 2017.

CONTENT:

Instructions for use for medical professionals for variant 1: Vitri Plate

Instructions for use for medical professionals for variant 2: Warm Plate

	Instructions for use for medical professionals Manipulation plate for assisted reproductive technologies, in variants Vitri Plate Draft the Russian Federation market
--	--

INSTRUCTIONS FOR USE FOR MEDICAL PROFESSIONALS

NAME OF THE MEDICAL DEVICE

Manipulation plate for assisted reproductive technologies, in variants

Variant 1:

- Vitri Plate;
- Instructions for use.

Hereinafter it may also be referred to as follows: 'Vitri plate', 'Vitrification plate', 'Plate', 'Medical device', spelling in a lowercase or capital letter, singular or plural, is considered equivalent.

RESPONSIBLE MANUFACTURER OF THE MEDICAL DEVICE

REPROLIFE, Inc., Japan
2-5-3-9F Shinjuku, Shinjuku-ku Tokyo,
160-0022,
Japan
Tel.: +813 59258931
info@reprolife.jp

MEDICAL DEVICE MANUFACTURING ADDRESS

REPROLIFE, Inc., Japan
2-5-3-9F Shinjuku, Shinjuku-ku Tokyo,
160-0022,
Japan
Tel.: +813 59258931
info@reprolife.jp

INTENDED USE OF THE MEDICAL DEVICE

Manipulation plate for assisted reproductive technologies, in variants, is intended for setting up the procedures of vitrification in the liquid nitrogen and consequent warming of patient's oocytes and embryos as part of assisted reproductive technologies (ART) implementation.

FIELD OF APPLICATION OF THE MEDICAL DEVICE

Vitri Plate is used in the field of reproductive technologies (IVF) for curing infertility associated with any disorders of the reproductive function, including uterine tubes patency disorder and commissures of female pelvic organs, as well as disorders related to male infertility factor. Medical device is also used for vitrification of the oocytes of oncology patients, vitrification taking place before the beginning of treatment (chemotherapy). Medical device can be also used in scientific studies that require vitrification or warming of biological objects.

OPERATING PRINCIPLE OF THE MEDICAL DEVICE

Vitrification is an ultrafast cell freezing method using liquid nitrogen (the cooling speed reaches 180,000 degree Celsius/min), during which intracellular water hardens, reaching a vitreous state bypassing the stage of ice crystals formation. Formation of ice crystals is the reason of cells damaging during freezing. Vitrification allows to avoid the formation of ice crystals and ensures better cell survival rates after warming. During vitrification, the water inside the cell is replaced by special cryoprotectant substances (ethylene glycol, hydroxypropyl cellulose) by the action of osmosis by means of dimethyl sulfoxide (DMSO), the universal solvent that opens cellular barriers. Vitrification is used for oocytes, blastocysts, and embryos at the cleavage stage.

	Instructions for use for medical professionals Manipulation plate for assisted reproductive technologies, in variants Vitri Plate Draft the Russian Federation market
--	--

Vitri Plate is a medical device which serves as a container for solutions: Plate wells are filled with solutions in which manipulations with objects are performed.

POTENTIAL CONSUMERS OF THE MEDICAL DEVICE

Vitri plate is intended for use in medical and healthcare facilities that are certified to provide assisted reproductive technology (ART) care, such as clinics, laboratories, and IVF centers.

Medical device must be handled by qualified embryologists with the necessary specialization and with relevant professional knowledge and skills, as well as experience in the field of assisted reproductive technologies (ART).

Specialized training in vitrification and warming procedures is recommended since the method has many technical nuances which play an important role in minimizing cell damage.

INDICATIONS FOR USE OF THE MEDICAL PRODUCT

- After hormonal stimulation of the patient, cryopreservation of embryo is recommended for creating a possibility of transfer during 'natural cycle'.
- Cryopreservation of oocytes is recommended in case obtaining the sperm of the husband is impossible, and when the sperm donorship issue is yet to be resolved.
- In case of a large number of embryos, it is recommended to perform cryopreservation on all of them except for 1–3 (by indication) if the decision to immediately transfer embryos was made.
- It is recommended for creation of donor oocyte bank.
- It is recommended for creating the bank of donor oocytes from patients in childbearing age who were diagnosed with cancer and are planning to have children, before the start of chemotherapy.

CONTRAINDICATIONS FOR USE OF THE MEDICAL PRODUCT

There are no special contraindications.

USAGE CONDITIONS FOR THE MEDICAL DEVICE

Vitri Plate can only be used in the specialized IVF laboratory with an incubator and cryostorage facility that is equipped for working with liquid nitrogen. The medical device must not be used in any other conditions.

Operating temperature 25–27 °C, RH — 45–80 %, atmospheric pressure 84.0–106.7 kPa (630–800 mm Hg).

Since the shelf-life of the medical device is limited (2 years), it is necessary to examine the label and expiry date before use. An expired medical device must not be used.

Vitri Plate must be stored and used in accordance with the instructions provided on the label and in the Instructions For Use for each version. Medical device that has been stored incorrectly must not be used.

Inspect the individual (primary) package for damage before use. Do not use medical device in case there are signs of damage to the individual (primary) package. Inspect medical device for signs of damage after opening the individual (primary) package. Do not use medical device in case there are signs of damage.

PRECAUTIONS

In vitro fertilization (IVF) and other assisted reproductive technologies departments must meet the requirements of the "Sanitary and epidemiological requirements for health-care institutions".

	Instructions for use for medical professionals Manipulation plate for assisted reproductive technologies, in variants Vitri Plate Draft the Russian Federation market
--	--

The staff of the institution should adhere to the security and safety measures generally established in laboratories, and use personal protection equipment.

The staff should be appropriately instructed and take precautions when working with liquid nitrogen and Dewar vessels.

Access of unauthorized persons and unqualified staff to the laboratory should be restricted.

Access of unauthorized persons and unqualified staff to the medical device should be restricted.

DESCRIPTION OF MEDICAL DEVICE VARIANT

Vitri Plate is a plastic medical device in the shape of a rectangular plate with three hemispheric wells and a slot for the vitrification carrier. Vitri plate is equipped with a plastic transparent lid.

Functionally, Vitri Plate is a medical device which serves as a container for solutions: solutions, in which manipulations with objects are performed, are transferred into the wells of the Plate. Each Vitri Plate is for single use.

The slot on the Vitri Plate is designed in a way, that vitrification carrier located in it will be in the same focal plane as the object, located in the wells of the Vitri Plate, which will be frozen. This allows to significantly simplify the vitrification procedure at the stage of transferring the object to be frozen on the carrier.

Vitrification carrier is a medical device recommended to use with the Plate (see Section "Description of medical products recommended for combined application").

Vitri Plate dimensions, L x W x H: (65.92 ± 0.50) mm x (26.92 ± 0.30) mm x (7.50 ± 0.20) mm.

Number of round wells: 3.

Round well diameter: (15.0 ± 0.2) mm.

Vitrification carrier device slot dimensions, L x W x H: (62.93 ± 0.50) mm x (3.25 ± 0.10) mm x (3.50 ± 0.10) mm.

Vitri Plate lid dimensions, L x W x H: (65.10 ± 0.50) mm x (26.10 ± 0.30) mm x (5.80 ± 0.20) mm.

Vitri Plate weight, with lid: (6.7 ± 0.3) g

Vitri Plate is for single use only. Do not reuse. Do not re-sterilize.

Shelf-life: two (2) years. Medical device must not be used after the expiry date. After the expiry date, an unused Vitri Plate must be discarded as a non-hazardous product in accordance with the disposal regulations adopted by the health-care institution.

MEDICAL DEVICE COMPOSITION

Polystyrene, liquid paraffin, other additives.

Medical device does not contain pharmaceuticals as well as materials of animal and/or human origin.

Medical device does not contain any substances classified as 'hazardous' or 'requiring special precautions' in handling and storage.

STERILIZATION

Medical device is sterile.

	Instructions for use for medical professionals Manipulation plate for assisted reproductive technologies, in variants Vitri Plate Draft the Russian Federation market
--	--

Sterilization type: irradiation sterilization.

Do not resterilize.

Medical device must not be used in case the integrity of the individual (primary) transparent polypropylene packaging has been compromised.

Single use only.

DESCRIPTION OF THE MEDICAL DEVICES RECOMMENDED TO USE WITH THE PLATE

Vitri Plate is recommended to use as a specialized medical device along with the following products produced by REPROLIFE, Inc., Japan, or equivalents:

- "Vitrification Kit for oocytes and embryos, variants" (hereinafter can be referred to as 'Vitrification Kit'). MA No. P3H 2019/9054 dated October 17, 2019;

- "Cryotec" container, different color variants" (hereinafter can be referred to as 'Cryotec'). Cryotec is available for combined usage with the Plate only upon receipt of the MA.

Description of the Vitrification Kit for oocytes and embryos variant "Vitrification Kit 101 for oocytes and embryos, composed"	Description of the Vitrification Kit for oocytes and embryos variant "Vitrification Kit 110 for oocytes and embryos, composed"
Vitrification Kit 101 for oocytes and embryos, composed, is intended for vitrification (freezing) of patient's oocytes and embryos and storing them in liquid nitrogen as part of assisted reproductive technologies (ART) implementation. The carton package contains: 1. Equilibration solution (ES), 1.0 ml (1 vial). 2. Vitrification solution (VS), 1.0 ml (2 vials). 3. Cryotec vitrification carrier device, 4 pcs 4. Vitrification plate, 3 pieces. 5. Instructions for use. Vitrification Kit 101 for oocytes and embryos is for single use only. Do not reuse. Do not re-sterilize. Shelf-life: one (1) year. Vitrification Kit must not be used after the expiry date. A utilized Vitrification Kit is to be disposed of as potentially contaminated waste in accordance with the State rules and regulations adopted by the health-care institution. After the expiry date, an unused Vitrification Kit must be discarded as a non-hazardous product in accordance with the disposal regulations adopted by the health-care institution. MA No. P3H 2019/9054 dated October 17, 2019	Vitrification Kit 110 is intended for vitrification (freezing) of patient's oocytes and embryos and storing them in liquid nitrogen as part of assisted reproductive technologies (ART) implementation. The carton package contains: 1. Equilibration solution (ES), 1.8 ml (2 vials); 2. Vitrification solution (VS), 1.8 ml (4 vials); 3. Instructions for use. Vitrification Kit 110 for oocytes and embryos is intended for five (5) applications Do not re-sterilize. Shelf life: 1 year. The Vitrification Kit is not to be used after the expiry date. A utilized Vitrification Kit is to be disposed of as potentially infectious waste in accordance with the public standards and regulations adopted by the health-care institution. After the expiry date, an unused Vitrification Kit is to be discarded as a non-hazardous product in accordance with the disposal regulations adopted by the health-care institution. MA No. P3H 2019/9054 dated October 17, 2019.

	Instructions for use for medical professionals Manipulation plate for assisted reproductive technologies, in variants Vitri Plate Draft the Russian Federation market
--	--

Please note that the plates included in the Vitrification Kit for oocytes and embryos, variant 'Vitrification Kit 101 for oocytes and embryos, composed', despite being named 'Vitrification plates', are equivalent to "Manipulation plate for assisted reproductive technologies, variant 1 — 'Vitri Plate'".

Description of "Cryotec" container, different color variants

"Cryotec" container, different color variants, (vitrification carrier) is intended for freezing (vitrification) in liquid nitrogen, storage and following warming of patient's oocytes and embryos during the use of assisted reproductive technology (ART), thus ensuring survival of the bioobjects.

"Cryotec" container is a special "carrier" for vitrification, consisting of thin transparent film, binded with a plastic handle (operating part). Operating part is closed by cap. All parts of "carrier" are stable to liquid nitrogen. The vitrification objects are placed on the film surface.

"Cryotec" is produced in different colors.

The design of the carriers is the same in all cases, the devices differ in handle color only: white, red, blue, green, yellow.

Shelf life: two (2) years. "Cryotec" must not be used after the expiry date. A utilized "Cryotec" must be disposed of as potentially contaminated waste in accordance with the State rules and regulations adopted by the health-care institution. After the expiry date, an unused "Cryotec" must be discarded as a non-hazardous product in accordance with the disposal regulations adopted by the health-care institution.

"Cryotec" is available for combined usage with the Plate only upon receipt of the MA.

MEDICAL PRODUCT APPLICATION PROCEDURE

Perform these actions in conditions described in Section 'Medical product application conditions'.

Open individual (primary) package of the Vitri Plate. Set the Vitri Plate on the microscope stand. Three hemispheric wells should be at the front, and the carrier slot should be in the back. Remove the lid.

It is recommended to use the Vitri Plate in combination with the Vitrification Kits. Perform the vitrification procedure in accordance with the guidelines from the Instructions For Use of the applied Vitrification Kits (see Section 'DESCRIPTION OF THE MEDICAL DEVICES RECOMMENDED TO USE WITH THE PLATE').

Perform these actions in conditions described in Section 'USAGE CONDITIONS FOR THE MEDICAL DEVICE'.

MEDICAL DEVICE STORAGE CONDITIONS

Store at 10–25 °C, keep away from direct sunlight.

Keep dry.

Do not use after expiry date.

	Instructions for use for medical professionals Manipulation plate for assisted reproductive technologies, in variants Vitri Plate Draft the Russian Federation market
--	--

Do not re-sterilize.

Do not reused.

TRANSPORT OF THE MEDICAL PRODUCT

Transport conditions: temperature 10–25 °C.

Keep dry.

Keep away from sunlight.

Proceed with caution during transfer.

Transportation is possible by any means of transport in covered vehicles, in compliance with the transportation conditions and in accordance with the shipping requirements and regulations for particular types of vehicles.

MEDICAL DEVICE DISPOSAL CONDITIONS

An utilized medical device must be disposed of as potentially contaminated waste in accordance with the State rules (Sanitary Regulations and Norms 2.1.7.2790-10 Sanitary and epidemiological requirements for handling medical waste) and regulations adopted by the health-care institution.

Do not reuse medical device.

Unused medical device must be disposed of as a non-hazardous product in accordance with the disposal regulations adopted by the health-care institution.

Medical device is not to be utilized along with household waste.

REQUIREMENTS FOR MAINTENANCE AND REPAIR

Medical device does not contain any components that could be subject to maintenance and repair.

Medical device is single use only.

AUTHORIZED REPRESENTATIVE OF THE MANUFACTURER IN THE RUSSIAN FEDERATION

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info@intergen.ru

Marketing Authorization: **** dated ****

SYMBOLS USED FOR MEDICAL PRODUCT LABELING

Symbol	Definition
REF	Catalogue number

	Instructions for use for medical professionals Manipulation plate for assisted reproductive technologies, in variants Vitri Plate Draft the Russian Federation market
--	--

	Batch code
	Manufacturer
	Date of manufacture
	Use by Date
	Sterilized using irradiation
	Caution, consult accompanying documents
	Caution, consult accompanying documents
	Do not use if the package is damaged
	Do not reuse
	Do not resterilize
	Keep dry
	Keep away from sunlight
	Storage conditions: temperature range from 10 °C to 25 °C
	Separate collection for waste

	Instructions for use for medical professionals Manipulation plate for assisted reproductive technologies, in variants Vitri Plate Draft the Russian Federation market
--	--

--	--

	Instructions for use for medical professionals Manipulation plate for assisted reproductive technologies, in variants Warm Plate Draft the Russian Federation market
--	---

INSTRUCTIONS FOR USE FOR MEDICAL PROFESSIONALS

NAME OF THE MEDICAL DEVICE

Manipulation plate for assisted reproductive technologies, in variants

Variant 2:

- Warm Plate;
- Instructions for use.

Hereinafter it may also be referred to as follows: 'Warm plate', 'Warming plate', 'Plate', 'Medical device', spelling in a lowercase or capital letter, singular or plural, is considered equivalent.

RESPONSIBLE MANUFACTURER OF THE MEDICAL DEVICE

REPROLIFE, Inc., Japan
2-5-3-9F Shinjuku, Shinjuku-ku Tokyo,
160-0022,
Japan
Tel.: +813 59258931
info@reprolife.jp

MEDICAL DEVICE MANUFACTURING ADDRESS

REPROLIFE, Inc., Japan
2-5-3-9F Shinjuku, Shinjuku-ku Tokyo,
160-0022,
Japan
Tel.: +813 59258931
info@reprolife.jp

INTENDED USE OF THE MEDICAL DEVICE

Manipulation plate for assisted reproductive technologies, in variants, is intended for setting up the procedures of vitrification in the liquid nitrogen and consequent warming of patient's oocytes and embryos as part of assisted reproductive technologies (ART) implementation.

FIELD OF APPLICATION OF THE MEDICAL DEVICE

Warm Plate is used in the field of reproductive technologies (IVF) for curing infertility associated with any disorders of the reproductive function, including uterine tubes patency disorder and commissures of female pelvic organs, as well as disorders related to male infertility factor. Medical device is also used for warming of the oocytes of oncology patients after vitrification, vitrification taking place before the beginning of treatment (chemotherapy). Medical device can be also used in scientific studies that require vitrification or warming of biological objects.

OPERATING PRINCIPLE OF THE MEDICAL DEVICE

Warming is a reverse process to vitrification. During warming the cryoprotector is removed from biological objects and is gradually being substituted by water, with the equalization of osmotic pressure in a controlled environment.

Warm Plate is a medical device which serves as a container for solutions: Plate wells are filled with solutions in which manipulations with objects are performed.

POTENTIAL CONSUMERS OF THE MEDICAL DEVICE

Warm plate is intended for use in medical and healthcare facilities that are certified to provide assisted reproductive technology (ART) care, such as clinics, laboratories, and IVF centers.

	Instructions for use for medical professionals Manipulation plate for assisted reproductive technologies, in variants Warm Plate Draft the Russian Federation market
--	---

Medical device must be handled by qualified embryologists with the necessary specialization and with relevant professional knowledge and skills, as well as experience in the field of assisted reproductive technologies (ART).

Specialized training in vitrification and warming procedures is recommended since the method has many technical nuances which play an important role in minimizing cell damage.

INDICATIONS FOR USE OF THE MEDICAL PRODUCT

- After hormonal stimulation of the patient, cryopreservation of embryo is recommended for creating a possibility of transfer during 'natural cycle'.
- Cryopreservation of oocytes is recommended in case obtaining the sperm of the husband is impossible, and when the sperm donorship issue is yet to be resolved.
- In case of a large number of embryos, it is recommended to perform cryopreservation on all of them except for 1–3 (by indication) if the decision to immediately transfer embryos was made.
- It is recommended for creation of donor oocyte bank.
- It is recommended for creating the bank of donor oocytes from patients in childbearing age who were diagnosed with cancer and are planning to have children, before the start of chemotherapy.

CONTRAINDICATIONS FOR USE OF THE MEDICAL PRODUCT

There are no special contraindications.

USAGE CONDITIONS FOR THE MEDICAL DEVICE

Warm Plate can only be used in the specialized IVF laboratory with an incubator and cryostorage facility that is equipped for working with liquid nitrogen. The medical device must not be used in any other conditions.

Operating temperature 25–27 °C, RH — 45–80 %, atmospheric pressure 84.0–106.7 kPa (630–800 mm Hg).

Since the shelf-life of the medical device is limited (2 years), it is necessary to examine the label and expiry date before use. An expired medical device must not be used.

Warm Plate must be stored and used in accordance with the instructions provided on the label and in the Instructions For Use for each version. Medical device that has been stored incorrectly must not be used.

Inspect the individual (primary) package for damage before use. Do not use medical device in case there are signs of damage to the individual (primary) package. Inspect medical device for signs of damage after opening the individual (primary) package. Do not use medical device in case there are signs of damage.

PRECAUTIONS

In vitro fertilization (IVF) and other assisted reproductive technologies departments must meet the requirements of the "Sanitary and epidemiological requirements for health-care institutions".

The staff of the institution should adhere to the security and safety measures generally established in laboratories, and use personal protection equipment.

The staff should be appropriately instructed and take precautions when working with liquid nitrogen and Dewar vessels.

Access of unauthorized persons and unqualified staff to the laboratory should be restricted.

	Instructions for use for medical professionals Manipulation plate for assisted reproductive technologies, in variants Warm Plate Draft the Russian Federation market
--	---

Access of unauthorized persons and unqualified staff to the medical device should be restricted.

DESCRIPTION OF MEDICAL DEVICE VARIANT

Warm Plate is a plastic medical device in the shape of a rectangular plate with 3 hemispheric wells and 1 rectangular well with a slightly sloping bottom. Warm Plate is equipped with a plastic transparent lid.

Functionally, Warm Plate is a medical device which serves as a container for solutions: solutions, in which manipulations with objects are performed, are transferred into the wells of the Plate. Each Warm Plate is for single use.

The '3 mm' marking on the lid of the Warm Plate is intended only for the approximate estimation of the solution volume in the pipette. It is not intended to measure anything and is not a measurement device for the medical use.

Warm Plate dimensions, L x W x H: (95.5 ± 0.5) mm x (31.0 ± 0.3) mm x (7.80 ± 0.2) mm.

Number of round wells: 3.

Round well diameter: (15.00 ± 0.20) mm.

Number of rectangular wells: 1.

Dimensions of a rectangular well: (25.23×25.23) mm ± 0.30 mm.

Warm Plate lid dimensions, L x W x H: (94.6 ± 0.5) mm x (30.1 ± 0.3) mm x (5.3 ± 0.2) mm.

Warm Plate weight, with lid: (9.5 ± 0.3) g

Warm plate is for single use only. Do not reuse. Do not re-sterilize.

Shelf-life: two (2) years. Medical device must not be used after the expiry date. After the expiry date, an unused Warm Plate must be discarded as a non-hazardous product in accordance with the disposal regulations adopted by the health-care institution.

MEDICAL DEVICE COMPOSITION

Polystyrene, liquid paraffin, other additives.

Medical device does not contain pharmaceuticals as well as materials of animal and/or human origin.

Medical device does not contain any substances classified as 'hazardous' or 'requiring special precautions' in handling and storage.

STERILIZATION

Medical device is sterile.

Sterilization type: irradiation sterilization.

Do not resterilize.

Medical device must not be used in case the integrity of the individual (primary) transparent polypropylene packaging has been compromised.

Single use only.

DESCRIPTION OF THE MEDICAL DEVICES RECOMMENDED TO USE WITH THE PLATE

Warm Plate is recommended to use as a specialized medical device along with the following products produced by REPROLIFE, Inc., Japan, or equivalents:

	Instructions for use for medical professionals Manipulation plate for assisted reproductive technologies, in variants Warm Plate Draft the Russian Federation market
--	---

- “Warming Kit for vitrified oocytes and embryos, variants” (hereinafter can be referred to as ‘Warming Kit’). Warming Kit is available for combined usage with the Plate only upon receipt of the MA;
- “Cryotec” container, different color variants” (hereinafter can be referred to as ‘Cryotec’). Cryotec is available for combined usage with the Plate only upon receipt of the MA.

Description of the Warming Kit for vitrified oocytes and embryos, variant ‘Warming Kit 102 for vitrified oocytes and embryos’	Description of the Warming Kit for vitrified oocytes and embryos, variant ‘Warming Kit 205 for vitrified oocytes and embryos’
Warming Kit 102 for vitrified oocytes and embryos is intended for warming patient’s oocytes and embryos in order to restore the activity of human cells after vitrification (freezing) as part of assisted reproductive technologies (ART) implementation. The carton package contains: 1. Warming Solution (TS), 1.8 ml (1 vial). 2. Diluent Solution (DS), 0.5ml (1 vial). 3. Washing Solution (WS), 1.0 ml (1 vial). 4. Warm plate, 1 piece. 5. Instructions for use. Warming Kit 102 for vitrified oocytes and embryos is for single use only. Do not reuse. Do not re-sterilize. Shelf-life: one (1) year. Warming Kit must not be used after the expiry date. A utilized Warming Kit is to be disposed of as potentially contaminated waste in accordance with the State rules and regulations adopted by the health-care institution. After the expiry date, an unused Warming Kit must be discarded as a non-hazardous product in accordance with the disposal regulations adopted by the health-care institution. Warming Kit is available for combined usage with the Plate only upon receipt of the MA.	Warming Kit 205 for vitrified oocytes and embryos is intended for warming patient’s oocytes and embryos in order to restore the activity of human cells after vitrification (freezing) as part of assisted reproductive technologies (ART) implementation. The carton package contains: 1. Warming Solution (TS), 1.8 ml (5 vials). 2. Diluent Solution (DS), 1.8 ml (1 vial). 3. Washing Solution (WS), 1.8 ml (2 vials). 4. Instructions for use. Warming Kit 205 for vitrified oocytes and embryos is intended for five (5) applications. Do not re-sterilize. Shelf-life: one (1) year. Warming Kit must not be used after the expiry date. A utilized Warming Kit is to be disposed of as potentially contaminated waste in accordance with the State rules and regulations adopted by the health-care institution. After the expiry date, an unused Warming Kit must be discarded as a non-hazardous product in accordance with the disposal regulations adopted by the health-care institution. Warming Kit is available for combined usage with the Plate only upon receipt of the MA.

Description of “Cryotec” container, different color variants
“Cryotec” container, different color variants, (vitrification carrier) is intended for freezing (vitrification) in liquid nitrogen, storage and following warming of patient’s oocytes and embryos during the use of assisted reproductive technology (ART), thus ensuring survival of the bioobjects.

	Instructions for use for medical professionals Manipulation plate for assisted reproductive technologies, in variants Warm Plate Draft the Russian Federation market
--	---

"Cryotec" container is a special "carrier" for vitrification, consisting of thin transparent film, binded with a plastic handle (operating part). Operating part is closed by cap. All parts of "carrier" are stable to liquid nitrogen. The vitrification objects are placed on the film surface.

"Cryotec" is produced in different colors.

The design of the carriers is the same in all cases, the devices differ in handle color only: white, red, blue, green, yellow.

Shelf life: two (2) years. "Cryotec" must not be used after the expiry date. A utilized "Cryotec" must be disposed of as potentially contaminated waste in accordance with the State rules and regulations adopted by the health-care institution. After the expiry date, an unused "Cryotec" must be discarded as a non-hazardous product in accordance with the disposal regulations adopted by the health-care institution.

"Cryotec" is available for combined usage with the Plate only upon receipt of the MA.

MEDICAL PRODUCT APPLICATION PROCEDURE

Open individual (primary) package of the Warm Plate. Keeping lid closed, place the Warm Plate in the incubator at 37 °C for not less than 3 hours before using. Take the Plate out from the incubator and set it on the microscope stand. The rectangular well should be on your left. Remove the lid.

It is recommended to use the Warm Plate in combination with the Warming Kits. Perform the warming procedure in accordance with the guidelines from the Instructions for use of the applied Warming Kits (see Section 'DESCRIPTION OF THE MEDICAL DEVICES RECOMMENDED TO USE WITH THE PLATE').

Perform these actions in conditions described in Section 'USAGE CONDITIONS FOR THE MEDICAL DEVICE'.

MEDICAL DEVICE STORAGE CONDITIONS

Store at 10–25 °C, keep away from direct sunlight.

Keep dry.

Do not use after expiry date.

Do not re-sterilize.

Do not reused.

TRANSPORT OF THE MEDICAL PRODUCT

Transport conditions: temperature 10–25 °C.

Keep dry.

Keep away from sunlight.

Proceed with caution during transfer.

Transportation is possible by any means of transport in covered vehicles, in compliance with the transportation conditions and in accordance with the shipping requirements and regulations for particular types of vehicles.

MEDICAL DEVICE DISPOSAL CONDITIONS

	Instructions for use for medical professionals Manipulation plate for assisted reproductive technologies, in variants Warm Plate Draft the Russian Federation market
--	---

An utilized medical device must be disposed of as potentially contaminated waste in accordance with the State rules (Sanitary Regulations and Norms 2.1.7.2790-10 Sanitary and epidemiological requirements for handling medical waste) and regulations adopted by the health-care institution.

Do not reuse medical device.

Unused medical device must be disposed of as a non-hazardous product in accordance with the disposal regulations adopted by the health-care institution.

Medical device is not to be utilized along with household waste.

REQUIREMENTS FOR MAINTENANCE AND REPAIR

Medical device does not contain any components that could be subject to maintenance and repair.

Medical device is single use only.

AUTHORIZED REPRESENTATIVE OF THE MANUFACTURER IN THE RUSSIAN FEDERATION

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
info@intergen.ru


Marketing Authorization: **** dated ****

SYMBOLS USED FOR MEDICAL PRODUCT LABELING

Symbol	Definition
	Catalogue number
	Batch code
	Manufacturer
	Date of manufacture
	Use by Date
	Sterilized using irradiation

2023

	Instructions for use for medical professionals Manipulation plate for assisted reproductive technologies, in variants Warm Plate Draft the Russian Federation market
--	---

	Caution, consult accompanying documents
	Caution, consult accompanying documents
	Do not use if the package is damaged
	Do not reuse
	Do not re sterilize
	Keep dry
	Keep away from sunlight
	Storage conditions: temperature range from 10 °C to 25 °C
	Separate collection for waste

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
APOSTILLE

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1. Country: JAPAN
This public document
2. has been signed by AONO Hiroshi
3. acting in the capacity of Notary of the Tokyo Legal Affairs Bureau
4. bears the seal/stamp of AONO Hiroshi . notary

Certified

5. at Tokyo
6. MAR. 29. 2021
7. by the Ministry of Foreign Affairs
8. 21- №041046
9. Seal/stamp:
10. Signature

Tanaka Toshie

TANAKA Toshie


For the Minister for Foreign Affairs

Registered No.203

NOTARIAL CERTIFICATE

This is to certify that Ms. TSUCHIYA Kaori, an agent of Mr. MUKAI Masahiro, CEO of REPROLIFE Inc., has stated in my very presence that said Mr. MUKAI Masahiro acknowledged himself to have signed and sealed the attached documents.

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