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CERTIFICATE

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China Council for the Promotion of International Trade is China Chamber of International Commerce

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China Council for the Promotion of International Trade
China Chamber of International Commerce

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THIS IS TO CERTIFY THAT: the signature of XIAOZHONG YUAN of JIANGYIN EVERISE MEDICAL EQUIPMENT CO.,LTD and the seal of the said company on the annexed Operation Manual are genuine.

China Council for the Promotion
of International Trade

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Authorized
Signature:

Shen Xiaohua

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(Date: Aug. 30, 2021)

Approved by

General Director

JIANGYIN EVERISE MEDICAL
EQUIPMENT CO.,LTD

Mr. Xiaozhong Yuan

Operation Manual

Patient Transfer Board with Carbon Fiber Frame WSX-S

Table of Contents

1.	Description, purpose and scope.....	3
2.	Principle of action.....	4
3.	Technical specifications.....	4
4.	Preventive measures.....	5
5.	Possible side effects.....	6
6.	Contraindications to use.....	6
7.	Scope of supply.....	6
8.	Working procedures.....	6
9.	Maintenance and storage.....	7
10.	Operating conditions.....	7
11.	Transportation terms and conditions.....	8
12.	Service life.....	8
13.	Recycling.....	8
14.	Marking.....	8
15.	Manufacturer.....	9
16.	Authorized Representative.....	9
17.	Warranty.....	9
18.	Regulatory Compliance	10

1. Description, purpose and scope

WSX-S patient transfer board with carbon fiber frame (hereinafter referred to as spine board) is used to transport patients with upper and lower limb injuries, back and cervical spine injuries.

The board can be used in traumatology, emergency medical services, military and extreme medicine, health care institutions, for transporting patients in the mountains, on the water, in the field.

WSX-S spine board with carbon fiber frame inside is a plastic board with holes around the perimeter, reinforced with carbon fiber rods to make it lighter and withstand heavy loads.

The holes have a dual purpose – they are used as handles for carrying the patient, as well as for attaching belts and fixing the patient. The design of the handles, raised relative to the base, allows easy lifting the board with the patient from the surface. The handles are easy of access and provide safe support when lifting or lowering and wearing the board.

The patient is fixed to the board with the help of fixing belts or belt system. The belt system consists of a central belt and 5 pairs of belts connected together. The end belts are statically sewn to the central belt, and the three belts located in the middle can move along the central belt to allow patients of any build and height to be fixed. The belts are fixed with a Velcro hook-and-loop fastener.

The set of fixing belts consists of 3 separate belts, adjustable in length and fixed with a Fastex fastener.

WSX-H1 head immobilizer a base with two side supports attached with holes for examining the patient's auricle. The head immobilizer is fixed with the help of fixing belts. The shape of the head immobilizer side supports provides an effective fixation around the ear and other head areas, limiting unintentional destabilization. Side supports can be fixed in different positions, making head immobilizer suitable for different categories of patients.

The materials used for the boards are easy to clean, resistant to mineral oil, fungal formation, and liquid penetration.

Figure 1 is a WSX-S spine board with carbon fiber frame with WSX-H1 head immobilizer and a set of fixing belts.

Figure 2 is a WSX-S spine board with carbon fiber frame with WSX-H1 head immobilizer and belt system.

2. Principle of action

In case of fractures of the spine and limbs, the patient requires immobilization. The purpose of spinal injuries immobilization is to prevent the displacement of the fractured vertebral body, to prevent compression of the spinal cord or its repeated traumatization during transportation, as well as damage to the vessels of the spinal canal and the formation of hematomas. For this purpose, the patient is placed on the spinal board. When a cervical fracture is suspected, the board is used in conjunction with the head immobilizer.

3. Technical specifications

Parameter Name	Value
Spine board with carbon fiber frame WSX-S	
Dimensions of spine board without accessories, mm	1845 x400x52
Wight of spine board without accessories, kg	4,3
Rated load, kg	250
Maximum deflection at lifting with 4 end handles, not more than, mm:	
- at nominal load 165 kg	45
- at nominal load 250 kg	90
Radiotransparency	radiotransparent, less than 1.5 mm Al
Handles	The design of the handles, raised relative to the base, allows easy lifting the board with the patient from the surface. The handles are easy of access.
Quantity of handles for carrying, pcs	14 (2 on each transverse side, 5 on each longitudinal side)
Rated load on the handles, kg	60
Head immobilizer WSX-H1	
Head immobilizer WSX-H1,dimensions,mm	395 x250x160
Weight, kg	0,6
Type of Fixation	Belt with Velcro fastener
Type of Fixation of side supports on the base	“Velcro” hook-and-loop fastener

The strength of the fastening of the belt fastener, not less than, N	200
The force required to open/close the belt fastener and remove the supports from the base, N	50
Belt system	
Total belt length of the belt system (mm)	7900
Belt width of the belt system, mm	60
Weight of the belt system, kg	0,35
Rated load, kg	55
Belt fixing type	"Velcro" hook-and-loop fasteners
Fastening strength, not less than, N	550
Force required to turn on/off fasteners, not greater than, N	50
Fixing belts set	
Maximum size of fixing belt, mm	1300 x50
Weight of full set, kg	0,27
Amount of fixing belts in a set, pcs	3
Rated load on a fixing belt, kg	85
Fastening strength, not less than, N	850
N force required to turn on/off fasteners is not greater than	50
Fixing type	Fastex
Packing specifications	
Dimensions of spine board packaging, mm	1840x430x60
Dimensions of head immobilizer packaging, mm	430x270x140
Weight of the spine board with accessories in package (including a fixing belts and belt system), kg	5,9
Weight of the Head immobilizer in package, kg	1,1

*-Allowable deviation:10%

eventive measures

Do not use boards for permanent fixation.

The board is intended for use by qualified specialists or specially trained personnel. The personnel should be physically prepared to transport the board with the patient.

- Please read the operation manual carefully before you start using the board.
- When carrying a patient on the spine board, make sure that the patient is fixed with the fixing belts.
- Control the patient and the spine board position during carrying.
- Do not load WSX-S spine board by more than 250 kg.
- At transporting a patient weighing more than 165 kg with a spinal injury use the board together with vacuum mattress.
- Do not operate the spine board with visible defects, as well as with violations of the structure of the materials.
- Do not use the spine board for a purpose other than its intended use.
- Do not leave the patient unattended.
- The product is not intended to be carried by one person.
- Check the product regularly. Perform the necessary maintenance to maintain the product in good condition and maintain its operational properties for a long time.
- Do not change or modify the product, do not repair the product.
- Avoid contact with heat sources, flame or flammable substances.

5. Possible side effects

If the product is used correctly, there will be no side effects.

6. Contraindications to use

There are no known contraindications.

7. Scope of supply *

No.	Name	Quantity, pcs
Patient Transfer Board with Carbon Fiber Frame WSX-S		
1.	Patient Transfer Board with Carbon Fiber Frame WSX-S	1
2.	A Set of Fixation belts and/or belt system	Not less than 1
3.	Head immobilizer WSX-H1	if necessary
4.	Operation Manual	1

* the specified scope of delivery is indicated on the product packaging.

8. Working procedures

8.1. Remove the shipping packaging from the board.

8.2. Check the appearance and configuration of the board. The surface must be free of mechanical defects, contamination, as well as violations of the structure of the materials.

Attention! The use of defective products is not allowed.

8.3. Lift the patient and carefully transfer him to the board, this should be done in one step with the help of 3 to 5 people. If the neck area has not been previously fixed, it is necessary to fix the patient's head with your hands with an emphasis on the shoulders. Make sure that the patient's legs and pelvic bones remain in the same position relative to the torso.

At suspicion on cervical spine injury trauma, the patient's cervical spine must be pre-fixed with a collar splint (not included).

At transporting a patient weighing more than 165 kg with a spinal injury use the board together with vacuum mattress.

Attention! Careless transfer of the patient may cause additional injuries.

8.4. Place the patient on the board exactly in the center. If the patient is more than 185 cm tall, the patient's feet may hang over the edge of the board, and the head must not hang down.

8.5. If you need to use WSX-H1 head immobilizer (suspected cervical spine injury), fix its base in the board with the help of fixation straps before placing the patient on the board. Place the patient on the board so that the head is centered on the base of the head immobilizer. Then install the side supports on both sides of the patient's head so that they fit closely to the patient's head and fix them to the base with "Velcro" hook-and-loop fasteners. Make sure that the patient's ears are located in the holes of the side supports. Secure the patient's forehead and chin with the head immobilizer straps.

Attention! Using the head immobilizer in prone position is prohibited!

Attention! Improper position of the patient's ears relative to the side supports can cause the damage of the auricle.

8.6. Fix the patient to the board with a belt system or a set of fixation belts, pass the belt through the holes in the board, and then fix it.

Attention! Do not tighten the belts too tightly – this can lead to compression and the formation of pressure sores.

Attention! Too weak fixation will not give the desired effect of immobilization.

8.7. Carefully move the patient.

Attention! Do not transfer the patient on the board unless the patient is strapped in with a belt system.

8.8. The spine board and accessories should be cleaned after use (see Paragraph 9.2.)

Additional information:

Transportation of a patient with cervical spine trauma should be carried out in the supine position. At thoracic section or lumbar spine trauma the normal position of the patient is on his back, if medically required, transportation can be carried in a prone position. To avoid dislocation of the vertebral fragments and additional damage to the spinal cord, the patient should not be transported in a side position.

9. Maintenance and storage

9.1. Before reuse, check whether the board is in proper condition: check the product for visible damage. Damaged products and accessories should be replaced.

9.2. After use, the board and accessories should be cleaned.

Spin Board, Head Immobilizer base and side parts should be cleaned with water, sponge and non-corrosive detergent cleaning.

Belt systems, fixation belts and head immobilizer belts should be washed in washing machine at 60°C with non-corrosive cleaners.

Disinfect of the board and accessories should be performed with chlorhexidine solution or 3% hydrogen peroxide solution and add non-corrosive detergents.

Attention! Alcohol solution and thermal treatment (boiling, autoclaved, etc.) are not allowed.

Avoid rigid mechanical cleaning that may damage the board surface.

After cleaning, thoroughly dry the product.

Attention! Re-use is possible only if the product is in normal technical condition, has no contamination and is disinfected.

9.3. Store the product at a temperature of -10° to +25° in well-ventilated enclosed areas without exposure to sunlight at a relative humidity of 25-80%. Do not store items in unheated ambulances.

9.4. The shelf life of the product is 3 years.

10. Operating conditions

Use the product at temperatures from *-40°C~+70°C and air relative humidity of 5-95%.

*-when using the product at too high or low temperatures, monitor the patient's condition and, if necessary, take measures to maintain optimal body temperature.

11. Transportation terms and conditions

The parts are transported by covered vehicles according to the applicable cargo transport rules for each mode of transport at temperatures of $-40^{\circ}\text{C}\sim+40^{\circ}\text{C}$ and relative humidity of 25~80. The transportation process should prevent shocks, severe fluctuations and humidity.

12. Service life

The average service life under the operating conditions of the board is 3 years, the accessories-1 year.

13. Recycling

Products and accessories at the end of their service life are subject to disposal as medical waste of class "A", are disposed of as solid household waste in a standard way according to SanPiN 2.1.7.3684-21 with preliminary disinfection.

Dispose of or recycle packaging materials in accordance with applicable local, state and regional regulations.

14. Marking

14.1. The marking is provided:

- Directly on board and head immobilizer;
- On the package.

14.2. The markings on the board and head immobilizer shall contain the following information:

- Name of the product;
- Manufacturer's logo;
- Symbol CE;
- Year of manufacturing.

14.3. The markings on the package shall include the following information*:

- Name of the product;
- Scope of supply;

- Name and address of the manufacturer;

- Symbols; 

- Symbols; 

- Lot number; 


- Date of production; 

- The name and address of the authorized representative;
- Registration Certificate number.

*- it is allowed to supplement the marking with other information that does not contradict the information specified in the technical documentation and other documents

14.4 Symbols description

Symbol	Description
--------	-------------

	Producers
	Marking in compliance with the European Quality Directive
	Note, refer to instructions
	LOT Number
	Date of production

15. Manufacturer

JIANGYIN EVERISE MEDICAL EQUIPMENT CO.,LTD(1001,Chengyang Road,
Jiangyin,Jiangsu,China) Tel :86-0510-86900998,medical@everisegroup.com www.everisegroup.com.

16. Authorized Representative

On the territory of the Russian Federation the Authorized Representative who carries out the sales, warranty obligations and repair of products is

RDM LLC

Address: Room 504, Floor 5, Grayvoronovskaya St. 23, Moscow, Russia 109518

Telephone:+7(495)109-59-55, +7(929)999-89-87 info@rdmed.ru

17. Warranty

The Manufacturer guarantees good condition of the spine board within 12 months from the date of sale while meeting the requirements of Operation Manual.

The warranty period of accessories (head immobilizer, belt system, fixation belts set) is 6 months.

Warranty is not applicable to damage caused by negligent custody, transportation and operation. If maintenance is required, please contact the manufacturer's authorized representative in the Russian Federation at the address and telephone number specified in this operation manual.

18. Regulatory compliance

Document	Description
Compliance with international standards	
EN ISO 13485:2016	Medical products. Quality management system.Requirements Regulations
EN 980:2008	Medical devices. Graphical symbols used in marking Medical devices
EN 1041:2008	Information provided by the manufacturer Medical devices

ISO 14971:2009	Medical devices. Application of Risk Management in Health Care Equipment
EN 1865-1:2010	Equipment used in ambulances to transport patients. Part 1 Specifications General stretcher system and patient transport equipment
ISO 10993	Evaluation of biological effects of medical devices
Compliance with Russian Federation standards	
GOST P 56330-2016	Medical products. Technical means of settling and transferring patients and wounded at pre-hospital stage. General technical requirements and test methods
GOST P 50444-2020	Medical devices, instruments and equipment. General technical conditions
ISO 10993 Series Standards	Evaluation of biological effects of medical devices
GOST P 52770-2016	Medical products. Safety requirements. Sanitary chemical and toxicological test methods
GOST 31214-2016	Medical products. Requirements for samples and documents submitted to toxicology, sanitary chemistry studies, aseptic and pyrogen tests
GOST P ISO 15223-1-2020	Medical products. Symbols used for medical devices, labels and marks on accompanying documents. Part 1 General requirements

Change Registration Form

Date	Version	Change
25.07.2020	001	Original version
02.08.2021	002	Changing of SanPiN number, introduction of a new technical characteristic "maximum deflection"

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