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LOT 47253901

 2021-01-31

COBAS, COBAS E and ELECSYS
are trademarks of Roche.

 Roche Diagnostics GmbH
Sandhofer Str. 116
D-68305 Mannheim

Distribution in USA by:
Roche Diagnostics,
Indianapolis, IN
Made in Germany

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REF 08814872190

For USA: **CONTENT**

- M** Streptavidin-coated microparticles 0.72 mg/mL
- R1** Two biotinylated monoclonal anti-HBsAg antibodies (mouse) > 0.5 mg/L
- R2** Monoclonal anti-HBsAg antibody (mouse), polyclonal anti-HBsAg antibodies (sheep) labeled with ruthenium complex > 1.5 mg/L
- Cal1** Human serum
- Cal2** HBsAg approximately 50 IU/mL in human serum
- Dil HepB** Human serum negative for HBsAg and anti-HBs

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



English

Precautions and warnings

For in vitro diagnostic use.
Exercise the normal precautions required for handling all laboratory reagents.
Disposal of all waste material should be in accordance with local guidelines.
Safety data sheet available for professional user on request.
This kit contains components classified as follows in accordance with the Regulation (EC) No. 1272/2008:
2-methyl-2H-isothiazol-3-one hydrochloride

EUH 208 May produce an allergic reaction.
Product safety labeling follows EU GHS guidance.
All human material should be considered potentially infectious.
The calibrators and HBSAGQ2 Dil HepB have been prepared exclusively from the blood of donors tested individually and shown to be free from HBsAg (HBSAGQ2 Cal1 and HBSAGQ2 Dil HepB only) and antibodies to HCV and HIV.

The testing methods used assays approved by the FDA or cleared in compliance with the European Directive 98/79/EC, Annex II, List A.
The serum containing HBsAg (HBSAGQ2 Cal2) was inactivated using β -propiolactone and UV-radiation.

However, as no inactivation or testing method can rule out the potential risk of infection with absolute certainty, the material should be handled with the same level of care as a patient specimen. In the event of exposure, the directives of the responsible health authorities should be followed.

Avoid foam formation in all reagents and sample types (specimens, calibrators and controls).

Method Sheet download

The box at the top of this document shows the product code together with the respective document version.
Please go to dialog.roche.com to download this document.
The website displays either the version listed here or, if available, any later version that is also valid for this product code. If you have no internet access, please contact your local Roche affiliate to obtain the document free of charge.

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