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CERTIFICATE

中国国际贸易促进委员会暨中国国际商会
China Council for the Promotion of International Trade is China Chamber of International Commerce

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China Council for the Promotion of International Trade
China Chamber of International Commerce

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THIS IS TO CERTIFY THAT: the seal of Jiangsu Suyun Medical Materials Co., Ltd. on the annexed DOCUMENT is genuine.

China Council for the Promotion
of International Trade

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证明书查询网址 Website for verifying the certificate: <http://www.rzccpit.com/validate.html>

«APPROVE»

Jiangsu Suyun Medical Materials Co., Ltd

Director, Chairman of the Board

(position)

_____/ Qin Hongping

(name)

(signature)

31. 03. 2021

«day» month, year (numerals)

Stamp

2019

User Manual

For Medical Devices

«DISPOSABLE STERILE UROGENITAL PROBES»

Version 02 from 10/10/2019

1 Medical device name

“Disposable sterile urogenital probes” in the following types:

1. Type A – “Universal”;
2. Type B – “Small curet”;
3. Type C – “Papel”: mod.1, mod.2;
4. Type D – “Cytobrush”: mod.1, mod.2;
5. Type E – “Ayre spatula”: mod.1, mod.2;
6. Type F – “Combined”: mod.1, mod. 2, mod. 3, mod. 4, mod. 5.

2 INFORMATION ABOUT THE MANUFACTURER OF THE MEDICAL DEVICE

Headquartes: Jiangsu Suyun Medical Materials Co., Ltd, No. 1 Medicine Lane, Renmin Rd., 222002 Lianyungang, People's Republic of China, Tel.: +8651885476682, E-mail: chairman@suyunmedical.com, web site: suyunmedical.com.

Manufacturer: Jiangsu Suyun Medical Materials Co., Ltd, No. 1 Medicine Lane, Renmin Rd., 222002 Lianyungang, People's Republic of China, Tel.: +8651885476682, E-mail: airman@suyunmedical.com, web site: suyunmedical.com.

Authorized representative in the Russian Federation: Limited Liability Company Unicornmed (Uni-cornmed LLC) 192019, St. Petersburg, Melnichnaya street, 18 lit. A. Tel.: +7 (812) 702-33-04, E-mail: info@unicmed.ru, website: unicmed.ru

3 MEDICAL DEVICE PURPOSE

Disposable sterile urogenital probes (hereinafter referred to as products, probe or probes) are intended for the collection of biological material from the surface of the mucous membrane of the vagina opening, cervix, cervical canal or urethra for the purpose of carrying out cytological diagnostic studies in gynecology, urology and venereology.

The product is sterile, cannot be sterilized, for single use.

1.1 Indications for use: collection of biological material from the surface of the vagina opening mucosa, cervix, cervical canal or urethra for the purpose of carrying out cytological diagnostic studies in gynecology, urology and venereology.

1.2 Side effects: none.

1.3 Contraindications: individual intolerance to materials.

1.4 Potential users: qualified medical personnel in the conditions of medical treatment-and-prophylactic institutions.

1.5 Method of application: use in accordance with the instructions.

2.6 Conditions of application: in medical institutions by qualified specialists.

2.7 Warnings:

ATTENTION! The product is intended for single use only!

ATTENTION! Reuse can lead to cross-infection.

ATTENTION! Medical personnel must wear gloves.

4 BASIC PARAMETERS AND CHARACTERISTICS OF THE MEDICAL DEVICE


Disposable sterile urogenital probe is a product intended for single use. It is packed in individual blister or polyethylene sterile packaging.

"Disposable sterile urogenital probes" in the following versions:

Name

Disposable sterile urogenital probe

№	Type	Model (mod.)	Probe length, mm (±10%)	Working part size, mm (length×diameter) (±10%)	Material	Application area
1	Type A - "Universal"	—	178.0	22.0×3.0	PP Working part of the probe with sputtering of cellulose esters	For collecting biological material from the cervical canal and urethra

Polypropylene grade PP1013H1 (PP).

The probe handle has 2 notches. The first notch at a distance of 87 mm allows the end with the working part to be broken off for placement in the container located in the hollow handle of the probe. The container is used to store and transport the working part of the probe with biomaterial. The second notch at a distance of 51 mm allows you to break off the end with the working part for placing directly into the test tube.


№	Type	Model (mod.)	Probe length, mm (±10%)	Working part size, mm (length×diameter) (±10%)	Material	Application area
2	Type B – «Small curet»	—	212.0	7.0×3.0 7.0×4.0	ABS	For collecting material from the surface of the vaginal mucosa, cervix and urethra for bacteriological study

Acrylonitrile butadiene styrene ABS-1000 (ABS).

A double-ended spoon-shaped probe with smooth or sharp edges allows efficient collection of material from the entire surface of the cervix.


Type	Model (mod.)	Probe length, mm ($\pm 10\%$)	Piston length, mm ($\pm 10\%$)	Material	Application area
Type C - "Papel"	1	277.0	220.0	PP	For intrauterine biopsy of the endometrium for histological study

Polypropylene grade PP1013H1 (PP).

The probe has a closed end with 4 side holes. It consists of two parts: a piston with a limiter and a probe. Insertion depth marking from 20 to 120 mm, graduated by 10 mm. The piston is cylindrical with a handle for holding and a depth stop for insertion into the probe.


№	Type	Model (mod.)	Probe length, mm ($\pm 10\%$)	Working part size, mm ($\pm 10\%$)	Material	Application area
4	Type C - "Papel"	2	255.0	30.0	PP	For collecting endometrial smears in the absence of the need for traction or expansion of the cervical canal for cytological study

Polypropylene grade PP1013H1 (PP)

Outside part: a transparent outline with marks of 40, 70 and 100 mm from the distal end, corresponding to the depth of insertion. Inner part: mandrel connected to a nylon brush, the tip of which is equipped with a cap similar in diameter to the outer part.


Type	Model (mod.)	Probe length, mm ($\pm 10\%$)	Working part size, mm ($\pm 10\%$)	Material	Application area
Type D - "Cytobrush"	1	200.0	23.0	ABS The working part is made of polyethylene bristles fixed on a rigid twisted stainless steel wire	For collecting material from the surface of the cervix and from the depth of the cervical canal for cytological and bacteriological studies

Acrylonitrile butadiene styrene ABS-1000 (ABS).

It has an atraumatic polymer ball at the end of the working part. On the handle of the probe at a distance of 38 mm there is a notch that allows you to break off the end with the working part for placing in a test tube.

№	Type	Model (mod.)	Probe length, mm ($\pm 10\%$)	Working part size, mm ($\pm 10\%$)	Material	Application area
6	Type D - "Cytobrush"	2	200.0	20.0	ABS The working part is made of polyethylene bristles fixed on a rigid twisted stainless steel wire	For collecting material from the surface of the cervix and from the depth of the cervical canal for cytological and bacteriological studies


Acrylonitrile butadiene styrene ABS-1000 (ABS).

On the handle of the probe at a distance of 38 mm there is a notch that allows you to break off the end with the working part for placing in a test tube.

№	Type	Model (mod.)	Probe length, mm ($\pm 10\%$)	Working part size, mm ($\pm 10\%$)	Material	Application area
7	Type E - "Ayre spatula"	1	180.0	127.0	PP	To obtain ectocervical smears and a smear of the posterior fornix, simultaneously from the surface of the cervix, from


the pharynx and the lower third of the cervical canal for cytological study

Polypropylene grade PP1013H1 (PP)

It is a bilateral scapula probe. The ends of the probe are of various shapes and sizes. Has perforation along the edge of the shoulder blades.


№	Type	Model (mod.)	Probe length, mm ($\pm 10\%$)	Working part size, mm ($\pm 10\%$)	Material	Application area
8	Type E - "Ayre spatula"	2	203.0	138.0	PP	To obtain ectocervical smears and a smear of the posterior fornix, simultaneously from the surface of the cervix, from the pharynx and the lower third of the cervical canal for cytological study

Polypropylene grade PP1013H1 (PP)

It is a bilateral scapula probe. The ends of the probe are of various shapes and sizes.


№	Type	Model (mod.)	Probe length, mm ($\pm 10\%$)	Working part size, mm (length $\times$ diameter) ($\pm 10\%$)	Material	Application area
9	Type F - "Combined"	1	200.0	30.0 $\times$ 18.0	PP	For the simultaneous collection of scrapings from the surface of the cervix and from the cervical canal (ecto- and endocervix) for the purpose of conducting cytological studies

Polypropylene grade PP1013H1 (PP)

It consists of a handle and a removable working part. The working part is a flat cervical brush in the form of a trapezoid, consisting of flexible polymer bristles of various lengths.


№	Type	Model (mod.)	Probe length, mm ($\pm 10\%$)	Working part size, mm ($\pm 10\%$)	Material	Application area
10	Type F - "Combined"	2	175.0	37.0	PP	For collecting scrapings from the surface of the cervix and cervical canal (ecto- and endocervix) for the purpose of conducting cytological studies

Polypropylene grade PP1013H1 (PP)

The working part of the probe is sprayed with cellulose esters. It consists of two parts: a hollow plastic handle and a diamond-shaped working part with variable geometry.


№	Type	Model (mod.)	Probe length, mm ($\pm 10\%$)	Working part size, mm ($\pm 10\%$)	Material	Application area
11	Type F - "Combined"	3	198.0	33.0	PP	For collecting scrapings from the surface of the cervix and cervical canal (ecto- and endocervix) for the purpose of conducting cytological studies, including liquid cytology

Polypropylene grade PP1013H1 (PP)

It consists of a handle and a removable working part. The working part is in the form of a combined cervical brush, consisting of flexible polymer bristles 13 mm long, located parallel to the probe handle and the central part of the brush 22 mm long with bristles 1 mm long perpendicularly arranged on it in two rows.


№	Type	Model (mod.)	Probe length, mm (±10%)	Working part size, mm (± 10 %)	Material	Application area
12	Type F - "Combined"	4	205.0	38.0	PP	For collecting scrapings from the surface of the cervix and cervical canal (ecto- and endocervix) for the purpose of conducting cytological studies

Polypropylene grade PP1013H1 (PP)

The working part in the form of a combined cervical brush, consisting of flexible polymer bristles 13 mm long, located parallel to the probe handle and a central part made of polyethylene bristles, fixed on a rigid twisted stainless steel wire with an atraumatic polymer ball at the end.

№	Type	Model (mod.)	Probe length, mm (±10%)	Working part size, mm (± 10 %)	Material	Application area
13	Type F - "Combined"	5	212	20.0 40.0	PP	For collecting scrapings from the surface of the cervix and cervical canal (ecto- and endocervix) for the purpose of conducting cytological studies

Polypropylene grade PP1013H1 (PP)

It consists of a comfortable handle and two working parts located at different ends of the handle.

Working part 1 - cytobrush 22 mm long. Made of polyethylene bristles attached to a tough twisted stainless steel wire.

Working part 2 - a spatula 26 mm long. It is a flat plastic spatula with a curved edge.

N o t e – the following are not used in the manufacture of the product: materials of animal and (or) human origin, medicines.

5 METHOD OF APPLICATION

- 4.1 Check the integrity of packaging and the expiration date of the product. Specify the version and model of the product, verify the method of use of the product using the flowchart shown on the product packaging.
- 4.2 Put on gloves after pretreating hands with antiseptic solutions. Carefully open the package, remove the probe without affecting its sterility.
- 4.3 Observing the rules of asepsis, carefully insert the probe without applying excessive force to avoid injury.
- 4.4 Carry out the necessary manipulations.
- 4.5 At the end of the manipulation, remove the probe.
- 4.6 Handle and dispose of the probe in accordance with the established procedure.

6 SAFETY REQUIREMENTS

The product is intended for single use only!

Reuse can lead to cross-infection.

Do not use the products if the package has been opened or damaged.

7 CLEANING AND DISINFECTION

This medical device is supplied sterile for single use. No additional cleaning and disinfection required before use.

8 STERILIZATION

This medical device is supplied sterile (sterilized using ethylene oxide) in accordance with the requirements of ISO 11135-1.

The method of validation and monitoring of the sterilization process is also carried out in accordance with ISO 11135-1.

Resterilization is not allowed.

9 LABELING, PACKAGING

Symbols used on the product and its packaging

Symbol	Description
	Date of manufacture Date when the product was manufactured
	Expiration date Date after which the product cannot be used
	Batch code Product batch code (or batch of goods)
	Prohibition on the reapplication The medical device is intended for a single use or for use on one patient during one procedure
	Sterilization with ethylene oxide It indicates that the product has been sterilized with ethylene oxide
	Do not re-sterilize It indicates that the product cannot be resterilized
	Do not use the product if the packaging is damaged It indicates that if the packaging is damaged, the medical device must not be used
	CE mark The product meets the requirements of European directives
	Keep away from moisture The product needs to be protected from moisture
	Temperature limits Temperature range at which the product should be stored or handled
	Relative humidity range Range of relative humidity at which the product should be stored or handled
	Protect from sunlight The product needs to be protected from sunlight
	Do not dispose with household waste Used products are decontaminated and disposed in accordance with the requirements of local organizations and disposal regulations in the area.

Medical device labeling includes the following information:

- trademark and (or) name of the manufacturer;
- trademark and (or) name of the manufacturer's authorized representative in the Russian Federation;
- product name;
- quantity;
- product information (technical specifications);
- lot (batch number) of the product;
- product shelf life;
- symbol "Sterilization using ethylene oxide" according to ISO 15223-1;
- symbol "Do not re-sterilize" according to ISO 15223-1;

- symbol "Do not use with damaged packaging" according to ISO 15223-1;
- flow chart for using the product.

Note – if necessary, additional signs may be placed on the product packaging.

Each product is packed in individual blister or polyethylene packaging. Individual packages with products according to the types of executions, packed in cardboard or polyethylene intermediate packages, and then in cardboard boxes.

10 USE, TRANSPORTATION AND STORAGE CONDITIONS

This medical device can be transported by all types of transport in covered vehicles in accordance with the rules for the carriage of goods in force for each type of transport.

When used, the products are resistant to the influence of biological fluid.

During use, transportation and storage, this medical device is resistant to temperatures from -5°C to $+35^{\circ}\text{C}$ and relative air humidity from 30 to 75%, without condensation, at atmospheric pressure from 200 to 1500 hPa. The temperature should not exceed $+45^{\circ}\text{C}$ in order to avoid deformation of products.

11 SERVICE LIFE (EXPIRATION DATE) AND INDICATION ON THE PROHIBITION OF THE USE OF THE MEDICAL DEVICE ON THE EXPIRATION OF THE SERVICE LIFE (EXPIRATION DATE)

The shelf life of a medical device is 5 years from the date of manufacture.

Do not use the products after the expiration date.

12 PROCEDURE FOR MEDICAL DEVICE DISPOSAL OR DESTRUCTION

11.1 The packaging is disposed at a household waste collection point.

Note. For the Russian Federation, the packaging of a new product is disposed in accordance with SanPiN (Sanitary Rules and regulations) 2.1.7.1322 "Hygienic requirements for the placement and disposal of production and consumption waste".

11.2 Used products are decontaminated and disposed in accordance with the requirements of local organizations and disposal regulations in the area. For more information on dates and locations where disposal takes place centrally, contact your local authority or waste disposal company.

Note. For the Russian Federation, the disposal of products used in the infectious ward of a medical facility is carried out in the manner prescribed by SanPiN 2.1.7.2790 "Sanitary and Epidemiological Requirements for the Management of Medical Waste" for class B waste (epidemiologically hazardous waste).

11.3 Recycling and disposal of unused products after the expiration date is carried out in accordance with the requirements for the disposal of epidemiologically safe waste.

Note. For the Russian Federation, in order to prevent environmental pollution, all waste generated during the disposal of products is subject to mandatory collection with subsequent disposal, in accordance with the established procedure and in accordance with the current requirements and norms of local regulatory documents, including SanPiN 2.1.7.2790 "Sanitary and Epidemiological Requirements for Medical Waste Management" for Class A waste.

11.4 Since the regulations for disposal may vary in different countries, if necessary, contact the manufacturer's authorized representative in the Russian Federation Unicornmed Limited Liability Company (Unicornmed LLC), 192019, St. Petersburg, Melnichnaya street, 18 lit. A, E – mail: info@unicmed.ru (tel. 8 (812) 702-33-04).

13 MAINTENANCE AND REPAIR

Not provided.

14 WARRANTY OBLIGATIONS

Jiangsu Suyun Medical Materials Co., Ltd assumes a warranty for the medical device to be free from defects and conform to the manufacturer's declared characteristics during the shelf life, subject to the conditions specified in the instructions for use.

The manufacturer is not responsible for damage or accidents resulting from misuse, user error or external emergencies. For product quality questions, please contact an authorized representative of the manufacturer.

15 RECLAMATION

In the event of a defect, the developer and manufacturer Jiangsu Suyun Medical Materials Co., Ltd at No. 1 Medicine Lane, Renmin Rd., 222002 Lianyungang, People's Republic of China guarantees the quality of the product. If you find a defect, please contact your authorized representative in the Russian Federation.

The authorized representative of the developer and manufacturer in the Russian Federation who accepts complaints is Unicornmed Limited Liability Company (Unicornmed LLC), 192019, St.

Petersburg, Melnichnaya street, 18 lit. A, E – mail: info@unicmed.ru (tel. 8 (812) 702-33-04) accepts complaints.

Defects are not considered violations of the product received during use and careless handling.

16 LIST OF INTERNATIONAL REGULATORY DOCUMENTS

ISO 15223–1:2016 "Medical devices - Symbols to be used with medical device labels, labeling, and information to be supplied - Part 1: General requirements";

ISO 11135:2014 "Sterilization of health-care products - Ethylene oxide - Requirements for the development, validation and routine control of a sterilization process for medical devices".

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