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
Flexible video cystourethroscope HD-VIEW in variants

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
Content

1	Important information	1
2	Symbols employed	3
3	Intended use	6
3.1	Indications	6
3.2	Contraindications	7
4	User qualification	8
5	Warnings and cautions	8
6	Safety instructions	9
6.1	Unpacking	18
6.2	Basic equipment	19
6.3	Operating conditions	20
7	External system components	20
8	Preparation before use	21
9	Handling	21
9.1	Operation	21
9.2	Holding the device correctly	21
9.3	Operating the suction valve	22
9.4	Handling the sheath	22
9.5	Handling the deflection mechanism	23
9.6	Handling the working channel	24
9.7	Handling the accessories (forceps, cleaning brushes, HF electrodes etc.)	26
9.8	Inspecting the surfaces	28
9.9	Inspecting the deflection	29
9.10	Connecting to IMAGE1 STM	29
9.11	Focus	29
9.12	Function buttons of the endoscope	30
9.13	Switching the flexible cystoscope HD VIEW ON/OFF	31
10	Reprocessing	32
10.1	General warnings	32
10.2	Accessories	33
10.3	Precleaning directly after use	34
10.4	Leakage test	35
10.5	Disassembly	37

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Content

10.6	Manual disinfection	40
10.6.1	Manual High Level Disinfection (HLD)	40
10.6.2	Drying	41
10.7	Machine cleaning and disinfection	41
10.7.1	Manual preparation	41
10.7.2	Connecting	42
10.7.3	Machine cleaning/chemical-thermal disinfection	42
10.8	Assembly, inspection, care and packaging	43
10.8.1	Packaging	44
10.9	Sterilization	44
10.10	Limits of reprocessing	46
10.11	Storage	46
11	Reprocessing of LUER port with stopcocks	48
11.1	Accessories	49
11.2	Preparation for cleaning and disinfection	49
11.3	Manual pre-cleaning	49
11.4	Machine cleaning and disinfection	50
11.5	Connecting	51
11.6	Assembly, inspection and care	51
11.7	Packaging systems	51
11.8	Sterilization	52
11.9	Limits of reprocessing	52
12	Technical data	53
13	Servicing and repair	55
14	Maintenance	56
15	Limitation of liability	56
16	Warranty	56
17	Directive compliance	56
18	Important information	56
19	Disposal	56
20	Electromagnetic Compatibility user information	57

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Important Information

1 Important information

Thank you for your expression of confidence in the KARL STORZ brand name. Like all of our other products, this product is the result of years of experience and great care in manufacture. You and your organization have decided in favor of a modern, high-quality item of equipment from KARL STORZ.

This instruction manual is intended to serve as an aid in the proper handling, cleaning, disinfection and, if appropriate, sterilization of the flexible cystoscope HD VIEW. All the necessary details and actions are clearly explained. We therefore ask that you read this manual carefully before proceeding to work with the instrument. Keep this manual available for easy reference.

It is recommended that the suitability of the products for the intended procedure be checked prior to use.

KARL STORZ is continuously working on the further development of all products. Please appreciate that changes to the scope of supply in form, equipment and technology are possible for this reason. Therefore, no claims may be derived from the information, figures and descriptions in this manual.

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
2 Symbols employed

	Manufacturer
	Consult instructions for use
	Catalogue number
	Number of products in the product packaging
	Serial number
	Non-sterile
	Fragile, handle with care
	Keep dry
	Temperature limit
	Federal (USA) law restricts this device to sale by or on the order of a physician.

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
Controls

- ① Connector
- ② Suction valve
- ③ Connection for leakage tester or pressure compensation cap (vent port)
- ④ Handle
- ⑤ Irrigation and working channel
- ⑥ Deflection lever
- ⑦ Distal tip
- ⑧ Working channel connector seal 11014 DS
- ⑨ Leakage tester 13242 XL
- ⑩ Suction valve, sterile, for single use 091011-20
- ⑪ Suction valve, reusable 11301 CB1
- ⑫ Pressure compensation cap 11025 E
- ⑬ LUER port with stopcocks 11014 L
- ⑭ Double LUER port 11014 LU
- ⑮ Cleaning adaptor for working channel 11014 RA
- ⑯ Cleaning adaptor 11301 CD1
- ⑰ Right navigation button
- ⑱ Left navigation button
- ⑲ Menu button

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3 Intended use

Flexible cysto-urethrosopes (fiber endoscopes and videoscopes) are used for endoscopic imaging of the urethra and the bladder for the purposes of diagnosis and treatment. Flexible cysto-urethrosopes with an integrated working channel allow for the use of flexible auxiliary instruments such as forceps, stone baskets etc., coagulation electrodes and laser fibers for lithotripsy or for the ablation of benign/malignant tissue. With flexible cysto-urethrosopes, the urethra and the entirety of the bladder can be visualized/treated. The special feature of the flexible cysto-urethroscope lies in the fact that it can be used percutaneously in the kidneys or suprapubically in the bladder via an operating sheath.

Flexible cysto-urethrosopes (fiber endoscopes and videoscopes) are designed for temporary or short-term use in endoscopic interventions performed through natural orifices or surgical invasive interventions.

3.1 Indications

Description of all known indications:

The use of cysto-urethrosopes is indicated if urethro-cystoscopy is indicated in the opinion of the attending physician.

This includes, amongst other things:

- Clarification of micro and macro hematuria
- Clarification/biopsy in the event of suspected bladder tumors
- Control urethro-cystoscopy as per TUR-B or P
- Clarification in particular in the event of suspected tumor infiltration due to adjacent organ cancer (e.g., colon carcinoma, cervical carcinoma)
- Clarification or treatment in particular in the event of suspected vesicenteric or obstetric fistulas
- Urethral or vesical foreign body removal
- Clarification of recurring urinary tract infections
- Clarification of LUTS/BES, subvesical obstructions, BOO, neurogenic micturition disorder

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Intended use

- Visual assessment and possibly treatment (e.g., injection) of the orifices of the ureters (localization, condition, etc.)
 - Visual clarification of the bladder sphincter in the event of incontinence
 - Visual assessment and possibly treatment of the urethra (tumor, stricture/stenosis, fistula)
 - Visual assistance for indwelling catheter system, e.g., with urethral trauma in the event of a pelvic fracture
 - Retrograde pyelography, ureteral splinting
- The use of cysto-urethrosopes is indicated if urethro-cystoscopy is indicated in the opinion of the attending physician.

3.2 Contraindications

The use of the cysto-urethrosopes is contraindicated if, in the opinion of the attending physician, the surgical or examination method as such is contraindicated or if the patient is not able to undergo surgery or anesthesia due to his or her general condition.

Cysto-urethrosopes must not be used for interventions in direct contact with the central nervous system (CNS) and central cardiovascular system.

The following contraindications also apply:

- Acute urethritis
- Acute epididymitis
- Acute prostatitis
- Acute/chronic cystitis/urinary tract infection
- Possible coagulation disorders
- Restless, agitated patients

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4 User qualification

The flexible cystoscope HD VIEW from KARL STORZ may only be used by persons with an appropriate medical qualification and who are acquainted with the flexible cystoscopy technique. The information given in these instructions only serves to instruct in the correct handling and reprocessing of the flexible cystoscope HD VIEW. It is not intended as an introduction to the techniques of flexible cystoscopy.

5 Warnings and cautions

 **HAZARD SYMBOL:** This is a hazard symbol. It warns of the risk of injury. To prevent injuries, all the measures marked with the hazard symbol must be followed.

Signal words:

WARNING! designates a possible imminent risk. If this is not avoided, it could lead to death or serious injuries.

ATTENTION! designates a possibly harmful situation. If it is not avoided, the medical devices could be damaged.

NOTE! indicates special information on the operation of the device or important information.

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6 Safety instructions

The flexible cystoscope HD VIEW must only be used according to the medical rules and procedures of endoscopy recognized for this method.

WARNING:

Risk of injury and damage to the product:

- It must be ensured that the electrically non-insulated parts of the endoscope do not come into contact with conductive surfaces or voltage-carrying elements of other devices at any time.

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- Never look directly into the light emitted by the flexible videoscope.

WARNING:

Risk of explosion: If flammable anesthetics are used in the immediate proximity of the flexible videoscope.

- The flexible videoscope must not be used in the presence of flammable or combustible media.

WARNING:

Risk of injury: When in operation, the integrated LED light source heats up the flexible videoscope.

- Switch off the flexible videoscope before putting it down, and do not place it in the vicinity of the patient (as per chapter 16 of the 3rd edition of IEC 60601-1).

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Safety Instructions

⚠️ WARNING:

Risk of injury: High energy light radiated through endoscopes may give rise to high temperatures in front of the light outlet, at the end faces of the fiberoptic light cable and at the tip of the endoscope.

- Prolonged direct contact of light-emitting surfaces on tissue can lead to burns; this must therefore be avoided at all costs.

⚠️ WARNING:

Risk of infection: The flexible videoscopes are not sterile when delivered, and must be cleaned and disinfected/sterilized prior to initial use and each subsequent reuse.

⚠️ WARNING:

Risk of infection: To avoid infection and maintain functionality, the chapter on reprocessing in the instruction manual must be observed. The cleaning, disinfection and sterilization procedures are explained therein in detail.

⚠️ WARNING:

Risk of infection: Before the device is used for the first time, it must be subjected to a complete reprocessing cycle, excluding precleaning.

⚠️ WARNING:

Risk from electromagnetic incompatibility: Medical electrical devices are subject to special precautions regarding Electromagnetic Compatibility (EMC). Observe the stated EMC instructions during installation and operation.

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⚠️ WARNING:

The use of accessories other than those specified in the KARL STORZ instruction manual may result in increased emissions or decreased interference immunity of the flexible videoscope. The accessories listed below have been shown to comply with the requirements of EN/IEC 60601-1-2. When using accessories other than those specified here, it is the responsibility of the operator to ensure that they comply with EN/IEC 60601-1-2.

⚠️ WARNING:

Risk of injury: The videoscope does not offer any electric insulation against high frequency (HF) voltages. Protection against HF voltage must be achieved using the appropriate HF instrument set.

- Please read the instructions for use in the instruction leaflet supplied with the HF application part for more information.

⚠️ WARNING:

Combinations of medical devices are only assured to be safe if they are listed as such in the respective instruction manuals or the intended use and interface specifications of the products used in the combination permit this.

- The instruction manuals and interface specifications of the medical devices used in combination must be observed precisely.

⚠️ WARNING:

The pressure compensation cap must be removed from the vent port for:

- endoscopic interventions,
- immersing in liquid, e.g., disinfectant solution and for machine reprocessing

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⚠ WARNING:

The pressure compensation cap must be in place on the vent port for:

- gas sterilization (EO, FO, hydrogen peroxide)
- transport and shipping

⚠ WARNING:

Only use accessories that are specifically suited for use with videoscopes in accordance with the intended use and which can be employed in line with generally accepted technical standards and safety guidelines.

⚠ WARNING:

Risk of injury: The flexible videoscope must not be used during the discharge of a cardiac defibrillator.

- Remove the flexible videoscope from the operation area before discharge takes place.

⚠ WARNING:

Risk of injury: Only KARL STORZ devices may be used in combination with the flexible videoscope. Disregarding this warning can endanger the patient and user.

ATTENTION: Risk of damage: Only use accessories (guide probes, cleaning brushes, etc.) from KARL STORZ. Accessories from other manufacturers may damage the videoscope.

ATTENTION: Do not allow any live components or instruments to come into contact with the flexible videoscope.

NOTE: Do not hit the distal tip of the flexible videoscope against hard objects, and never use a sharp object to remove dirt.

NOTE: Never place heavy objects on the flexible videoscope.

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Safety instructions

NOTE: Never clamp or fold the sheath of the flexible videoscope. This can result in damage to the videoscope.

NOTE: Constant changing between various reprocessing methods is especially hard on the materials and must therefore be avoided.

NOTE: The flexible videoscope must not be carried by the sheath. Do not pull, clamp or twist the sheath of the videoscope.

NOTE: The sheath of the flexible videoscope must never be clamped or coiled up too tightly. It is best to store the video videoscope in a suspended position or to keep it sterilized in a tray specifically designed for starting and storage (39406 AS).

WARNING:

Risk of injury from clouding of the lens.

- Do not apply lubricant to the lens, to avoid clouding the image and thereby prevent risks for the patient!

WARNING:

When inserting the flexible videoscope into the patient, the distal tip of the flexible videoscope must be in the straight, non-deflected position.

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⚠️ WARNING:

Risk of injury: Before inserting or removing the flexible videoscope, the distal tip must be put into the neutral position, as otherwise injuries to the patient or damage to the instrument may occur.

⚠️ WARNING:

Risk of injury: If the distal tip can no longer be properly controlled (sticking etc.) via the deflection lever, or if the deflection mechanism develops any other fault, the flexible videoscope must not be used.

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- Place the distal tip of the flexible videoscope in the neutral, non-deflected position and remove the flexible videoscope slowly and carefully from the patient.

NOTE: The resistance that can be felt when deflecting the distal tip using the deflection lever increases as the angle of deflection increases. Actuate the deflection gently and carefully, without any sudden or jerky movements, as far as the limit stop.

⚠️ WARNING:

Risk of injury: The flexible videoscopes from KARL STORZ are not intended for use with light sources. The use of non-compatible devices can lead to patient injury or device damage.

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⚠️ WARNING:

Risk of injury: Whenever the flexible videoscope is used, you must first check that it, and any accessories used in combination with it, are in perfect condition. Damaged flexible videoscopes or damaged accessories must not be used.

NOTE: There should be no instruments in the suction channel whilst suctioning is being performed. The instrument channel inlet should be sealed to prevent additional air from being drawn in.

NOTE: Do not squeeze the distal tip of the flexible videoscope too tightly.

⚠️ WARNING:

Risk of injury: Actuate the deflection lever gently and carefully as far as the limit stop. The deflection lever must never be operated abruptly or with force, to avoid injuring the patient or damaging the device.

➤ Move the deflection lever smoothly and slowly!

⚠️ WARNING:

Risk of injury: Before each use, the functionality of the deflection lever and the connected distal tip must be checked.

➤ Defective flexible videoscopes must not be used.

⚠️ WARNING:

Risk of injury: If the image fails or there is interference during use in the patient, put the distal tip into the neutral position and remove the device from the patient's body gently and carefully. If the flexible videoscope continues to be used in the patient after the image has failed or when there is interference, this may result in serious injury to the patient.

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Safety Instructions

⚠️ WARNING:

Risk of injury: Only insert suitable accessories, e.g., biopsy forceps, into the instrument channel of the videoscope when the distal tip is in a neutral, non-deflected position.

- Once the instrument has been introduced, carefully actuate the lever gently as far as the limit stop.
- The deflection lever must never be operated abruptly or with force, to avoid injuring the patient or damaging the device.

⚠️ WARNING:

Risk of injury: When HF electrodes or laser probes are used, it must be ensured that power is not applied until the active parts have left the instrument channel and are in the user's field of vision, and the site of the operation is clearly recognizable. Activation in the instrument channel may lead to serious injury to the patient and damage to the videoscope.

NOTE: Do not apply force to insert accessories, e.g., biopsy forceps, through the working channel, to prevent damaging the device. It should be possible to move the accessory easily, and without any hindrance, through the entire channel.

NOTE: When the distal tip is deflected, do not insert instruments into the videoscope. Only insert the instrument when the distal tip is in the neutral position and deflect the deflection lever gently and carefully up to the limit stop.

NOTE: Inserting sharp-tipped accessories into the instrument channel may damage the videoscope.

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Safety Instructions

⚠️ WARNING:

Risk of injury: Foreign bodies and items to be extracted must not be removed via the instrument channel with the aid of instruments. Instead, the entire videoscope must be taken out of the body and the foreign body must then be removed.

NOTE: For the correct settings for the camera control unit, please see the operating instructions of the devices used.

NOTE: Accessories such as biopsy forceps that are kinked or bent will not function properly, are difficult to guide through the channel and may cause severe damage to it. They must therefore be replaced.

⚠️ WARNING:

Ensure that the working channel is completely dry prior to sterilization to avoid damage to the working channel.

⚠️ WARNING:

Always unplug the flexible videoscope before carrying out any maintenance and cleaning procedures on it.

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Handling

6.2 Basic equipment

Flexible cystoscope HD VIEW 11272 VHK/VHU with following accessories:

Set 11272VHK

11272 VH	Flexible cystoscope HD VIEW
11301 P	Seal for suction valve seat
11014 L	LUER port with stopcocks
11014 DS	Working channel connector seal
11014 RA	Cleaning adaptor for working channel
11301 CD1	Irrigation adaptor

Set 11272 VHUK

11272 VHU	Flexible cystoscope HD VIEW
11301 P	Seal for suction valve seat
11014 LU	Double LUER port
11014 DS	Working channel connector seal
11014 RA	Cleaning adaptor for working channel
11301 CD1	Irrigation adaptor

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6.3 Operating conditions Flexible cystoscope HD VIEW:

- Ambient temperature: 15 °C ... 30 °C
- Humidity: 15% – 80%

7 External system components

Further information on external system components can be found in the IMAGE1 S™ (96206286) instruction manual.

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8 Preparation before use

1. Assemble LUER attachment. To do this, insert the lug of the LUER attachment into the recess on the working channel. Push and turn clockwise until it audibly clicks into place.

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2. Assemble the suction valve/valve seal.

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9 Handling

NOTE: Any damage to the flexible cystoscope HD VIEW resulting from incorrect operation is not covered by the manufacturer's warranty.

9.1 Operation

The following instructions outline the correct operation of the flexible cystoscope HD VIEW and the corresponding accessories. This is not an introduction to the techniques of cystoscopy. Instructions on endoscopic technique can be found in the relevant medical literature.

9.2 Holding the device correctly

Grip the handle of the flexible cystoscope HD VIEW in your left or right hand, as shown in the picture. With this kind of control grip, the flexible cystoscope HD VIEW can be operated with one hand.

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9.3 Operating the suction valve

When the suction valve is assembled, a suction pump must be connected to the suction valve via a tube. Suction is activated upon pressing the valve button.

9.4 Handling the sheath

ATTENTION: Risk of damage to the sheath:

- Do not carry the flexible cystoscope HD VIEW by the sheath.
- Do not pull, clamp or twist the sheath of the flexible cystoscope HD VIEW.
- The sheath of the flexible cystoscope HD VIEW must not be linked or coiled up too tightly.
- It is best to store the flexible cystoscope HD VIEW in a suspended position or to keep it sterilized in a tray specifically designed for sterilizing and storage (39403 AS).

WARNING:

Risk of injury: Do not apply lubricant to the lens, to avoid clouding the image and thereby prevent risks for the patient!

If, during use, the image falls, is disturbed, or becomes blurred, or if vision in general is impeded, put the deflection mechanism into the neutral position and slowly pull the sheath of the flexible cystoscope HD VIEW out of the patient's body.

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9.5 Handling the deflection mechanism

NOTE: The neutral position of the distal tip is straight, without any deflection whatsoever.

Deflection functions as follows with type 11272 VH endoscopes:

If the lever is pushed upwards, the tip is lowered.

If the lever is pushed downwards, the tip is deflected upwards.

This mechanism is inverted with model 11272 VHU.

ATTENTION: When the flexible videoscope is inserted into the patient, the distal tip must be in a straight, non-deflected position.

⚠️ WARNING:

The deflection lever must never be operated with force or abruptly to avoid injuring the patient or damaging the device!

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⚠️ WARNING:

Risk of injury and damage from faulty deflection mechanisms: If the distal tip can no longer be properly controlled (sticking etc.) via the deflection lever, or if the deflection mechanism develops any other fault, the flexible videoscope must not be used.

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⚠️ WARNUNG:

Risk of injury: If the flexible videoscope malfunctions during use on the patient, discontinue use immediately.

- Place the distal tip of the flexible videoscope in the neutral, non-deflected position and remove slowly and carefully from the patient.

⚠️ WARNUNG:

Risk of damage: The lens at the distal end of the flexible videoscope may become detached, e.g., due to mechanical forces.

- For this reason, it is essential to check that the lens is undamaged and present before and after every decontamination.

9.6 Handling the working channel

Do not force instruments or cleaning brushes through the channels. This can lead to damage to the working channel, especially if the tip of the instrument is deflected.

When connecting and operating suction/irrigation devices, pay attention to the instruction manuals for the particular devices. The irrigation devices recommended by KARL STORZ are set to a maximum pressure of 54 kPa (400 mmHg/0,54 bar/7,8 psi).

- The LUER adaptor 27014 Y or the seal 27550 C can be used to seal the instrument channel when inserting instruments.
- Instruments can be introduced through the instrument channel of the flexible cystoscope HD VIEW. Make sure that the accessories are of the appropriate length and diameter. Unsuitable accessories can damage the flexible cystoscopes HD VIEW. If you have any questions concerning accessories, please contact KARL STORZ directly or an authorized specialist dealer.

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Handling

- The flexible cystoscope HD VIEW must be in the straight, non-deflected position when inserting an instrument into the working channel.

Forceps, scissors and other accessories whose jaw can be opened or closed must always be inserted into the instrument channel with closed jaws. The jaws may only be opened when the distal tip of the instrument is completely clear of the channel. Opening the jaws in the instrument channel can cause damage to the videoscope.

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ATTENTION: Risk of damage: Only use
accessories (guide probes, cleaning brushes,
etc.) from KARL STORZ. Accessories from other
manufacturers may damage the videoscope.

9.7 Handling the accessories (forceps, cleaning brushes, HF electrodes etc.)

 WARNING:

Risk of injury: When HF electrodes or laser
probes are used, it must be ensured that power
is not applied until the active parts have left the
instrument channel and are in the user's field
of vision, and the site of the operation is clearly
recognizable. Activation in the instrument channel
may lead to serious injury to the patient and
damage to the videoscope.

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Handling

ATTENTION: Risk of damage: Only use
accessories (guide probes, cleaning brushes,
etc.) from KARL STORZ. Accessories from other
manufacturers may damage the videoscope.

NOTE: Do not apply force to insert accessories,
e.g., biopsy forceps, through the working channel,
to prevent damaging the device. It should be
possible to move the accessory easily, and without
any hindrance, through the entire channel.

NOTE: When the distal tip is deflected, do not
insert instruments into the videoscope. Only insert
the instrument when the distal tip is in the neutral
position and deflect the deflection lever gently and
carefully up to the limit stop.

NOTE: Inserting sharp-tipped accessories into the
instrument channel may damage the videoscope.

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9.8 Inspecting the surfaces

⚠ WARNING:

Risk of injury: The surface of the flexible sheath must be inspected visually and by carefully feeling it, to detect any cracks, cuts, dents, foreign bodies or other forms of damage.

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Before each use, verify that the device is in good condition by checking the surface of the flexible cystoscope HD VIEW, especially the sheath area and the distal tip, for sharp edges. Cystoscopes HD VIEW with sharp edges must not be used.

Otherwise they may cause serious injury to the patient or user.

Check the connecting cable for damage (such as fractures, cracks, twisted or crushed areas, etc.) and the connecting plug for corrosion.

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- Move the deflection lever smoothly and slowly!

⚠ WARNING:

Risk from damage: Faulty flexible videoscopes must not be used.

- Before each use, check the functionality of the deflection lever and the distal tip.

The distal tip is deflected by actuating the deflection lever slowly.

The distal end will deflect in accordance with the direction of actuation (forward or backward).

9.10 Connecting to IMAGE1 S™

The flexible cystoscope HD VIEW is connected to the camera control unit for operation.

For operation, an IMAGE1 S™ system with an IMAGE1 S Connect® and an IMAGE1 X-Link are required. Both camera modules must be connected to each other.

The flexible cystoscope HD VIEW is then connected to the IMAGE1 X-Link.

9.11 Focus

The image focus is adjusted automatically. Manual focusing is not possible/necessary.

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9.13 Switching the flexible cystoscope HD VIEW ON/OFF

The endoscope is switched off when the plug is pulled out, the camera control unit is switched off or when switching to another Link module.

WARNING:

Risk of burns and glare: The LED light emitted through the flexible videoscope may cause glare or high temperatures, and thereby burns, at the distal light outlet.

- > Never look directly into the light emitted by the flexible videoscope.

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WARNING:

Risk of infection: Only reprocess the medical device as per the documents provided with the device or the reprocessing instructions issued by KARL STORZ. These can be downloaded from www.karlstorz.com.

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WARNING:

Risk of infection from non-sterile delivery. Reprocess the device before each use. Check the medical device for visible soiling before use.

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WARNING:

Before reusing the device, replace any defective parts or send the flexible endoscope for repair.

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ATTENTION: Risk from damage to the device: The pressure compensation cap must be in place on the vent port for sterilization and transport to ensure pressure compensation.

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ATTENTION: Risk from damage to the device: Never clean flexible endoscopes with metal brushes.

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Reprocessing

ATTENTION: Risk from damage to the device: Flexible endoscopes must never be steam sterilized (autoclaved).

ATTENTION: Risk from damage to the device: To prevent damage to the flexible endoscope and accessories, it must be ensured that all connections are completely dry prior to use.

NOTE: When carrying out any work on contaminated medical devices, the guidelines of the Employers' Liability Insurance Association and equivalent organizations striving to ensure personal safety must be observed.

NOTE: Observe the chemical manufacturer's specifications concerning concentration and exposure times. Bear in mind the microbiological range of action of the chemicals.

NOTE: Damage to the device. Only perform reprocessing using chemicals approved by KARL STORZ. A list of approved chemicals can be found at www.karlstorz.com

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10.2 Accessories

Accessories required for carrying out reprocessing: pressure compensation cap, manual leakage tester

Place of application	Art. number brush/ adaptor
Suction channel	27651 H, 11301 CD1
Working channel	27651 H, 11014 RA
Clean surface	27652
Vent port	13242 XA, 11025 E

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10.3 Precleaning directly after use

Heavy soiling, corrosive solutions and pharmaceuticals must be removed from the endoscope immediately after use.

Wipe clean the exterior surfaces of the endoscope with a damp, single-use cloth.

Rinse out the working channel in order to prevent drying.

NOTE: To avoid damaging the sheath, do not use pushing and pulling movements.

NOTE: Do not use disinfectants for precleaning; these cause fixation of proteins and organic residues.

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


10.4 Leakage test

The endoscope should be tested to ensure it is leakproof after use and before reprocessing.

Performing the leakage test – dry

The connection tube of the leakage tester and the vent port on the endoscope must be dry. Connect the leakage tester connector securely to the supply connector vent port and lock by turning the bayonet catch clockwise 25°.

NOTE: The rubber seal in the connection adaptor of the leakage tester should be treated regularly with instrument oil.

Inflate the leakage faster until the maximum test pressure of 160–180 mmHg (blue area) is reached. The deflection rubber on the distal end may expand slightly. Using the deflection mechanism on the control body, move the deflection part in all possible directions as far as possible and observe the pointer on the manometer. The test pressure can drop by up to 10 mmHg until the endoscope is completely filled with air.

After stabilizing the pressure on the manometer, wait at least 2 minutes to see whether the pressure remains constant or drops.

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Reprocessing

Performing the leakage test – in liquid

NOTE: To avoid damage from the penetration of liquid, a test in the dry state must be performed in advance.

NOTE: Any lumina must be covered with no air bubbles by filling in a targeted manner.

Immerse the device in solution so that it is completely covered. It must be ensured that the leakage tester does not become wet. Perform the leakage test as described in "Leakage test – dry" and note whether bubbles emerge along the device. If, after an observation period of at least 3 minutes, no bubbles can be seen emerging (continuously) from the device and, in addition, the manometer does not show a fall in pressure, the device is leakproof.

Remove the endoscope from the liquid and press the vent button on the manometer. The test pressure must drop to 0 before the leakage tester can be removed.

NOTE: Never release the pressure in the liquid. Cleaning can then begin.

If a pressure loss occurs, this means that the endoscope leaks and may no longer be used. Immersed in liquid or reprocessed by machine.

ATTENTION: In this case, fill out the accompanying repair form (see Appendix or Art. no. 96096004DER), mark the endoscope, which has been cleaned by wiping to remove organic soiling/residue, as not disinfected and return it to the Service Department for repair or, alternatively, contact KARL STORZ. Line the storage/carrying case with the transport protection sheet 13990 SFN.

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10.5 Disassembly

The flexible cystoscope HD VIEW must be completely dismantled for cleaning/disinfection.

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
1. Dismantle the suction valve by pressing and turning.

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2. Dismantle the LUER attachment by pressing and turning.

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3. Dismantle the stopcock (if present) on the LUER attachment.

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
5. Clean the channels and lumen of the LUER attachment with a suitable brush.

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6. Dismantle the reusable suction valve 11301 CB1 by unscrewing the metal cap.

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
8. Brush the working channel connector.

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
9. Brush the working channel and suction channel.

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10. Connect the cleaning adaptor 11301 CD1 to
the valve seat.

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10.6 Manual disinfection

The endoscope must be completely immersed in a cleaning solution. The suction and working channel must be covered with no air bubbles by filling in a targeted manner. At the end of the necessary exposure time as specified by the manufacturer, the endoscope must be rinsed several times in order to remove all chemical residues. Sterile filtered water must be used for the final rinsing.

10.6.1 Manual High Level Disinfection (HLD)

The endoscope must be completely immersed in a high level disinfectant solution. The working and suction channel must be covered with no air bubbles by filling in a targeted manner. At the end of the necessary exposure time as specified by the manufacturer, the endoscope must be rinsed several times in order to remove all chemical residues. Sterile filtered water must be used for the final rinsing.

The endoscope can be disinfected with the following high level disinfection solutions which have been tested by way of example:

- 2.4% glutaraldehyde solution (e.g., CIDEX® Activated Dialdehyde Solution)
- 0.55% ortho-phthalaldehyde solution (e.g., CIDEX® OPA Solution)
- 2.0% accelerated hydrogen peroxide solution (e.g., Revital-Ox® Resert™ or Revital-Ox Resert XL HLD)

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Reprocessing

10.6.2 Drying

Finally, the surfaces and channels/lumina are dried completely with (preferably sterile) compressed air in accordance with the national regulations. The cleaning gun with accessories (Art. no. 27660) is ideal for this purpose.

ATTENTION: The pressure must not exceed 100 kPa (1 bar/14.5 psi). It may therefore be necessary to operate the cleaning gun with a pressure reducer.

10.7 Machine cleaning and disinfection

The following methods for machine decontamination have been validated and approved subject to compliance with the process parameters described in the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' (Art. no. 96216003DER):

10.7.1 Manual preparation

The endoscope must be completely immersed in a cleaning solution. The working and suction channel must be covered with no air bubbles by filling in a targeted manner. After the necessary exposure time as specified by the manufacturer, clean the instrument mechanically with the aid of a suitable brush. The channels are brushed until the brush head emerges and is visibly clean. Before retracting the brush head that emerges at the distal end, clean it. Never brush back and forth within the channel.

Finally, it must be rinsed with cold water to ensure neutralization.

NOTE: KARL STORZ recommends using the adaptors 11301CD1 and 11014RA for simpler filling of the lumina.

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10.7.2 Connecting

In order to ensure effective machine cleaning and disinfection, the endoscope's working and suction channel must be connected to the washer and disinfectant using the recommended adaptors 11301 CD1 and 11014 RA.

ATTENTION: The pressure compensation cap must not be in place on the vent port for cleaning and disinfection.

10.7.3 Machine cleaning/chemical-thermal disinfection

A maximum temperature of 65°C must not be exceeded. The relevant national requirements and the chemical manufacturer's instructions must be taken into account when using this method.

The selection of a suitable endoscope holder, which should ensure that the medical device is thoroughly rinsed out or through, must take place in consultation with the manufacturer of the device.

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10.8 Assembly, inspection, care and packaging

The cleaned and disinfected medical device must be visually inspected for cleanliness, completeness, damage and dryness:

- If residues or contamination are still present, the medical device must be manually cleaned and subjected to a full cleaning and disinfection procedure once more.

NOTE: Clean the distal lens with a cotton tip applicator soaked in lens cleaner or 70% isopropyl alcohol. When cleaning the distal lens, always wipe with a circular movement.

- Damaged or corroded medical devices must be withdrawn from use.

NOTE: Inspect plastic parts for discoloration, porosity and flexibility. If the parts have changed substantially from their original condition, they must be replaced. This especially applies to seals and rubber caps.

- If there is still residual moisture on or in the endoscope, it must be dried again manually. The surfaces (e.g., connection) and channels/lumina must be completely dried using (preferably sterile) compressed air in accordance with the national regulations. The cleaning gun with accessories (Art. no. 27660) is ideal for this purpose.

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- Afterwards, a functional check must be carried out.

NOTE: Check for a clear view and steerability.

NOTE: During care procedures, use items from the catalog „Care, Sterilization and Storage Techniques“.

10.8.1 Packaging

Only standardized and approved packaging materials and systems may be used (EN 868 Parts 2-10, EN ISO 11607 Parts 1 + 2, DIN 58953).

10.9 Sterilization

The procedures as well as the process-relevant parameters for the individually validated methods are described in detail in the manual „Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments“ (Art. no. 96218003). The method must be selected taking into account the respective applicable national requirements and in consultation with the device and product manufacturers.

⚠️ WARNING:

WARNING: Before sterilization, all parts of the endoscope must be thoroughly cleaned, disinfected, rinsed and then dried as otherwise sterility cannot be achieved.

⚠️ WARNING:

Risk of infection: Sterilize the medical device in a disassembled state.

ATTENTION: Flexible endoscopes must never be steam sterilized (autoclaved).

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Reprocessing

ATTENTION: The pressure compensation cap must be in place on the vent port for sterilization and transport to ensure pressure compensation.

NOTE: When using reprocessing equipment, please observe the equipment manufacturers' instructions.

NOTE: All deviations from the recommended sterilization parameters must be validated by the user.

NOTE: Sterilization is not possible on surfaces which have been lubricated and oiled.

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Hydrogen peroxide (H₂O₂) sterilization - ASP STERRAD®

ATTENTION: Please note that there are restrictions concerning the instruments which can be sterilized in the different STERRAD® sterilization systems in respect of lumen dimensions and material.

NOTE: Refer to the "STERRAD® Sterility Guide" to ensure that the medical device concerned can be sterilized using the different STERRAD® devices.

The following STERRAD® sterilization methods have been validated and approved for this medical device by KARL STORZ:

- STERRAD® 100S Long Cycle with Boosters
- STERRAD® NX® Advanced Cycle
- STERRAD® 100NX® Flex Cycle
- STERRAD® 100NX® DUO Cycle

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Reprocessing

Hydrogen peroxide (H₂O₂) sterilization - STERIS® AMSCO® V-PRO™1

NOTE: The "STERIS Low Temperature Reprocessing Device Compatibility Matrix" can be used to verify whether the medical device can be sterilized using the various STERIS V-PRO® devices.

The following STERIS V-PRO® sterilization methods have been validated and approved for this medical device by KARL STORZ:

- V-PRO® maX Flexible Cycle
- V-PRO® 60 Flexible Cycle

10.10 Limits of reprocessing

The product's service life and correct functioning are largely determined by mechanical stress and chemical influences within the scope of reprocessing and application.

10.11 Storage

The storage method must be selected taking into account the respective applicable national requirements.

KARL STORZ only recommends dust-protected, suspended storage for decontaminated endoscopes.

Horizontal storage is possible for sterilized endoscopes. KARL STORZ recommends the use of a specially designed tray for sterilization and storage.

ATTENTION: Do not expose flexible endoscopes to direct sunlight or excessive heat.

ATTENTION: The distal tip of the flexible endoscope must be in the neutral position during storage.

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Reprocessing


ATTENTION: Do not store the flexible endoscope in the case.

ATTENTION: The sheath of the flexible endoscope must never be kinked or coiled up too tightly.

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11 Reprocessing of LUER port

⚠ **WARNING:**

Risk of infection: Incorrectly reprocessed medical devices expose patients, users and third parties to a risk of infection as well as the risk that the medical device may malfunction. Observe the manual "Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments (art. no. 96216003DER)" and the accompanying documentation.

⚠ **WARNING:**

Risk of infection: These medical devices are not sterile when delivered. The use of non-sterile medical devices poses a risk of infection for patients, users and third parties. Inspect medical devices for visible contamination. Visible contamination is an indication that reprocessing has not been carried out or has been carried out incorrectly. Reprocess the medical devices before initial use and before and after every subsequent use using validated procedures.

⚠ **WARNING:**

When carrying out any work on contaminated medical devices, the guidelines of the Employers' Liability Insurance Association and equivalent organizations striving to ensure personal safety must be observed.

ATTENTION: When preparing and using solutions, follow the chemical manufacturer's instructions, paying close attention to proper concentration, exposure time and service life. Prolonged immersion and incorrect concentration may result in damage. Bear in mind the microbiological range of action of the chemicals used.

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Reprocessing

ATTENTION: Danger of damage to medical devices: The use of chemicals that have not been approved by KARL STORZ or chemical manufacturers may cause damage to the medical devices. Only use chemicals approved for the medical device for reprocessing. A helpful guide is available on the internet at www.karlstorz.com.

NOTE: National laws and regulations must be observed.

NOTE: The manual 'Cleaning, Disinfection, Care, and Sterilization of KARL STORZ Instruments' (Art. no. 96216003) can be downloaded or requested by visiting www.karlstorz.com.

11.1 Accessories

Accessories required for carrying out reprocessing:

Brushes: 27652, 27650B

Adaptor: 39107ALK (only for 11014L)

11.2 Preparation for cleaning and disinfection

Heavy soiling, corrosive solutions and pharmaceuticals must be removed from the medical device immediately after use. To this end, pre-clean the medical device by wiping down and rinsing, for example. As a general rule, KARL STORZ recommends manual pre-cleaning under cold running water.

11.3 Manual pre-cleaning

Disassembly (only valid for 11014L)

Prior to cleaning and disinfection, the medical device must be separated into its individual components as far as possible, and/or must be opened.

Brushing the surfaces

Any visible contamination or heavy soiling must be removed by cleaning the surfaces under cold running water with the aid of a brush or sponge.

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Reprocessing

Brushing the lumina

The working channels, lumina and hollow spaces must be pre-cleaned under cold running water with the aid of a suitable brush.

Rinsing out with a pressure water gun

The working channel, lumina and hollow spaces must be rinsed using a pressure water gun and a suitable rinsing attachment. For rinsing, pulse at least 3 times for 5 seconds from a jet of water with a pressure of at least 2.5 bar. The cleaning gun with accessories (Art. no. 27660) is ideal for this purpose.

11.4 Machine cleaning and disinfection

The following methods for machine decontamination have been validated and approved subject to compliance with the process parameters described in the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' (Art. no. 96216003):

Machine cleaning/chemical-thermal disinfection

A maximum temperature of 65 °C must not be exceeded. The relevant national regulations and the chemical manufacturer's instructions must be taken into account when using this method.

Machine cleaning/thermal disinfection

Thermal disinfection is preferred. The relevant national requirements and the A₀ value must be taken into account when using this method.

The selection of a suitable slide-in tray or instrument holder, which should ensure that the medical device is thoroughly rinsed out or through, must take place in consultation with the manufacturer of the device.

NOTE: If necessary, the instrument must be dried off manually afterwards.

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Reprocessing

11.5 Connecting

In order to ensure effective machine cleaning, disinfection and rinsing, the instrument must be connected in the washer and disinfectant.

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- If residues or contamination are still present, the medical device must be manually cleaned and subjected to a full cleaning and disinfection procedure once more.
- Damaged or corroded medical devices must be withdrawn from use.
- Joints, threads and glide surfaces must be treated locally with instrument oil and/or special lubricant. (only valid for 11014L)
- Dismantled medical devices must be assembled.
- Afterwards, a functional check must be carried out.

NOTE: The oil used for this purpose must be suitable for the subsequent sterilization procedure (silicone-free and paraffin- or white oil-based).

NOTE: During care procedures, use items from the catalog 'Care, Sterilization and Storage Techniques'.

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11.8 Sterilization

The procedures as well as the process-relevant parameters for the individually validated methods are described in detail in the manual "Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments" (Art. no. 96216003). The method must be selected taking into account the respective applicable national requirements and in consultation with the device and product manufacturers.

The following sterilization methods have been validated and approved by KARL STORZ for this medical device:

Steam sterilization using the fractionated prevacuum procedure

The medical device must be sterilized in its fully assembled state using the fractionated prevacuum procedure (DIN EN ISO 17665-1) at 132°C–137°C with a minimum exposure time of 3 minutes and a maximum of 18 minutes. This procedure is only suitable for thermostable instruments.

NOTE: Greased components must be disassembled for sterilization, so that steam penetration is ensured.

11.9 Limits of reprocessing

The end of the product's service life is largely determined by wear, reprocessing methods, the chemicals used and any damage resulting from use.

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Servicing and repair

13 Servicing and repair

Defective instruments must be serviced and repaired exclusively by persons authorized by KARL STORZ; all repair work must employ original parts only.

14 Maintenance

Preventive maintenance is not essential. However, regular maintenance can contribute to identifying potential problems before they become serious, enhancing the instrument's reliability and extending its useful operating life.

Maintenance services can be requested from your local representative or from the manufacturer. Regardless of the accident prevention regulations or testing intervals for medical instruments prescribed in different countries, we recommend a functional or safety test of the unit at least once a year.

15 Limitation of liability

KARL STORZ shall be liable for failure or deterioration in the safe operation, operational reliability, and performance of this instrument only subject to the conditions that all assembly, system expansions, readjustments, modifications, or repairs to same have been performed by persons authorized by KARL STORZ, and that the instrument has been used in accordance with its operating instructions at all times.

16 Warranty

The warranties provided can be found in the Standard Conditions of Business of KARL STORZ. The medical device must always be sent to your local subsidiary (see the section "Subsidiaries"), even during the warranty period.

Opening of the equipment or performance of any repairs or modifications to the equipment by unauthorized persons shall relieve us of any liability for its performance. Any such opening, repair, or modification performed during the warranty period shall void all warranty.

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
17 Directive compliance

In accordance with the Medical Device Directive (MDD):

11272 VHA/VHU:

This medical device belongs to Class IIa

This medical device bears a CE mark in accordance with the Medical Device Directive (MDD) 93/42/EEC. A code number after the CE mark indicates the responsible notified body.

18 Important information

To avoid the spread of infectious diseases among hospital staff or among KARL STORZ employees, telescopes, instruments and equipment must be cleaned and sterilized/disinfected before they are sent for repair. We reserve the right to return contaminated instruments/devices to the sender.

Repairs, modifications or expansions which are not performed by KARL STORZ or by experts authorized by KARL STORZ will invalidate all warranty rights.

KARL STORZ gives no warranty on the correct functioning of equipment or instruments which have been repaired by unauthorized third parties.

19 Disposal

This device has been marked in accordance with the European Directive 2012/19/EC on waste electrical and electronic equipment (WEEE). At the end of its useful operating life, dispose of the device as electronic scrap. Please ask KARL STORZ SE & Co. KG, a KARL STORZ subsidiary, or your authorized dealer about the collection point in your area. Within the scope of application of this Directive, KARL STORZ SE & Co. KG is responsible for the proper disposal of this device.

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20 Electromagnetic Compatibility user information

The flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™ is suitable for use within close proximity to an active HF electro-surgical device in professional healthcare facility environment.

Professional healthcare facility includes physician offices, dental offices, limited care facilities, freestanding surgical centers, freestanding birth centers, multiple treatment facilities, hospitals (emergency rooms, patient rooms, intensive care, surgical rooms, outside the HF shielded room of an ME system for magnetic resonance imaging).

WARNING:

Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

WARNING:

The use of an accessory or cable with the flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™ other than those specified in this manual may result in increased emissions or decreased immunity of the flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™. When using an accessory or cable with the flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™ other than those specified in this manual, it becomes the responsibility of the user of the flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™ to determine compliance with IEC 60601-1-2 when using this item.

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⚠️ WARNING:

Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

NOTE: The emission characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 Class A). If it is used in a residential environment (for which CISPR 11 Class B is normally required) this equipment might not offer adequate protection to radio communication service. The user might need to take mitigation measures, such as relocating or re-orienting equipment.

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Table 1 Guidelines and manufacturer's declaration – electromagnetic immunity			
The flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™ is intended for use in an electromagnetic environment specified below. The user of the flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™ should make sure that it is used in such an environment.			
Immunity tests	EN/IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) acc. to IEC 61000-4-2	± 8 kV contact discharge ± 15 kV air discharge	± 8 kV contact discharge ± 15 kV air discharge	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Rapid transient electrical Interferences/bursts acc. to IEC 61000-4-4	± 2 kV for power lines lines ± 1 kV for input and output lines 100 kHz repetition	± 2 kV/1kV ± 1 kV for input and output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surges acc. to IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and fluctuations of the supply voltage acc. to IEC 61000-4-11	Voltage Dip: Dip to 0% for 1 cycles @ 0° phase angle Dip to 70% for 25/30 cycles @ 0° phase angle Dropout to 0% for 0.5 cycles @ 0°, 45°, 90°, 135°, 180°, 225°, 270° & 315° phase angles Voltage interruption: 100% for 250/300 cycles	Voltage Dip: Dip to 0% for 1 cycles @ 0° phase angle Dip to 70% for 25/30 cycles @ 0° phase angle Dropout to 0% for 0.5 cycles @ 0°, 45°, 90°, 135°, 180°, 225°, 270° & 315° phase angles Voltage interruption: 100% for 250/300 cycles	Mains power quality should be that of a typical commercial or hospital environment. If the user of the flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™ requires continued operation during power supply interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.
Magnetic field at the supply frequency (50/60 Hz) acc. to IEC 61000-4-8	30 A/m at 50Hz/60Hz	30 A/m at 50Hz/60Hz	If image distortion occurs, it may be necessary to position the flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™ further from sources of power frequency magnetic fields or to install magnetic shielding. The power frequency magnetic field should be measured in the intended installation location to assure that it is sufficiently low.
Radiated, radio-frequency, electromagnetic field immunity test acc. to IEC 61000-4-3	3V/m 80 MHz to 2.7 GHz Refer table 2 for wireless Proximity RF field test levels	3V/m 80 MHz to 2.7 GHz	
Immunity to conducted disturbances, induced by radio-frequency fields acc. to IEC 61000-4-6	3Vrms on 150 kHz to 80MHz 1kHz 80% AM modulation 6Vrms in ISM bands	3Vrms on 150 kHz to 80MHz 1kHz 80% AM modulation 6Vrms in ISM bands	

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Table 2

Test levels for Proximity fields from RF wireless communications equipment

Test Frequency MHz	Band MHz	Service	Modulation	Immunity Test Level V/m	Compliance level V/m
385	380 – 390	TETRA 400	Pulse modulation 18 Hz	27	27
450	430 – 470	GMPRS 460 FRS 460	FM $\pm$ 5 kHz deviation 1 kHz Sine Wave	28	28
710	704 – 787	LTE Band 13 & 17	Pulse modulation 217 Hz	9	9
745					
780					
810	800 – 960	GSM 800/900 TETRA 800 IDEN 820 CDMA 850 LTE Band 5	Pulse modulation 18 Hz	28	28
870					
930					
1720	1700 – 1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation 217 Hz	28	28
1845					
1970					
2450	2400 – 2570	Bluetooth WLAN 802.11 b/g/n RFID 2450 LTE Band 7	Pulse modulation 217 Hz	28	28
5240	5100 – 5800	WLAN 802.11 a/n	Pulse modulation 217 Hz	9	9
5500					
5785					

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Table 3

**Guidance and manufacturer's declaration – electromagnetic immunity –
for equipment and systems that are not life-supporting**

The flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™ is designed for operation in an electromagnetic environment as defined below. The user of the device should make sure that it is used in such an environment.

Interference immunity tests	EN/IEC 60601 test level	Conformity level	Electromagnetic environment – guidelines
Conducted RF disturbances acc. to IEC 61000-4-6	3 V _{rms} 150 kHz to 80 MHz	3 V _{rms}	Portable and mobile RF communications equipment should be used no closer to any part of the flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = 1.2 \sqrt{P}$
Radiated RF disturbances acc. to IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	where P is the nominal power of the transmitter in watts [W] according to the information provided by the manufacturer of the transmitter and d is the recommended separation distance in meters [m]. Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a should be less than the $d = 1.2 \sqrt{P} \quad 80 \text{ MHz to } < 600 \text{ MHz}$ $d = 2.3 \sqrt{P} \quad 800 \text{ MHz to } 2.5 \text{ GHz}$ compliance level in each frequency range ^b . Interference may occur in the vicinity of equipment marked with the following symbol: 

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. The propagation of electromagnetic quantities is affected by absorptions and reflections of buildings, objects and persons.

^a Field strengths from fixed transmitters, such as base stations for radio [cellular/cordless] telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to stationary transmitters, an electromagnetic site survey should be considered. If the measured field strength at the location where the flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™ is used exceeds the above compliance levels, the flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™ should be monitored to ensure proper function. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

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Table 4 Emission class and group			
Guidelines and manufacturer's declaration – electromagnetic emissions			
The flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™ is intended for use in the electromagnetic environment specified below. The customer or user of the flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™ should ensure that it is used in such an environment.			
RF emissions as per CISPR 11		Group 1	The flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™ uses RF energy only for its internal function. The customer or user of the flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™ should make sure that it is used in such an environment. Therefore, its RF emissions are very low and are not likely to cause interference in nearby electronic equipment.
RF Emissions as per CISPR 11		Class A	The flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™ is suitable for use in all establishments other than domestic and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions as per IEC 61000-3-2		Class A	
Voltage fluctuations/flicker Emissions as per IEC 61000-3-3		Complies	

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Table 5

Recommended separation distances between portable and mobile RF telecommunication devices and the flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™

The flexible cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™ is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™ can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the cystoscope HD VIEW 11272VH/VHU together with the IMAGE1 S™ as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter [W]	Separation distance d [m] as a function of the transmitter frequency		
	150 kHz to 80 MHz $d = 1,2 \sqrt{P}$	80 MHz to 800 MHz $d = 1,2 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2,3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance D in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE: At 80 MHz and 800 MHz, the separation distance for higher frequency range applies.

NOTE: These guidelines may not apply in all situations. The propagation of electromagnetic quantities is affected by absorptions and reflections of structures, objects and people.

HF surgical instrument compatibility

This equipment has been evaluated for compatibility with high-frequency surgical equipment.

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KARL STORZ Endoscopy Canada Ltd.
7171 Millcreek Drive, Mississauga, Ontario
L5N 3R3 Canada
Phone: +1 905 816-4500, Fax: +1 905 816-4599
Toll free phone: 1-800-268-4880 (Canada only)
Toll free fax: 1-800-492-4199 (Canada only)
E-Mail: info-canada@karlstorz.com

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100, Fax: +1 424 218-8525
Toll free phone: 800 421-0837 (USA only)
Toll free fax: 800 321-1304 (USA only)
E-Mail: communications@kseaa.com

KARL STORZ Veterinary Endoscopy-America, Inc.
1 South Los Cameros Road
Goleta, CA 90317, USA
Phone: +1 805 968-7776, Fax: +1 805 685-2583
E-Mail: info@karlstorzvet.com

KARL STORZ Endoscopy Latino-America, Inc.
815 N. W. 57th Avenue, Suite 480
Miami, FL 33126-2042, USA
Phone: +1 305 262-8980, Fax: +1 305 262-8986
E-Mail: info@kseaa.com

KARL STORZ Endoscopia México S.A. de C.V.
Edificio Atlantic, Oficina 36
Calle Del Tra y 3ra
10400 Vedado, Manana, Cuba
Phone: +537 836 95 06, Fax: +537 836 97 76
E-Mail: kstorz@kscm.com

KARL STORZ Endoscopia México S.A. de C.V.
Lago Constanza No 325,
Col. Granada Del. Miguel Hidalgo
C.P. 11520 México D.F.
Phone: +52 (55) 1101 1520
E-Mail: mx-info@karlstorz.com

KARL STORZ Marketing América Do Sul Ltda.
Rua Joaquim Floriano, nº. 413, 2º andar - Itaim Bibi,
CEP-04534-011 São Paulo, Brasil
Phone: +55 11 3526-4600, Fax: +55 11 3526-4680
E-Mail: br-info@karlstorz.com

KARL STORZ Endoscopy Argentina S.A.
Zufalegui 627 6º Piso, B1638 CAA - Vicente Lopez
Provincia de Buenos Aires, Argentina
Phone: +54 11 4718 0918, Fax: +54 11 4718 2773
E-Mail: info@karlstorz.com.ar

KARL STORZ Endoskop Norge AS
Østergate 15B
0661 Oslo, Norway
Phone: +47 6380 5600, Fax: +47 6380 5601
E-Mail: post@karlstorz.no

KARL STORZ Endoskop Sverige AB
Storsågarvägen 14
127 38 Skärholmen, Sweden
Phone: +46 8 505 648 00
E-Mail: kundservice@karlstorz.se

KARL STORZ Endoscopy Suomi Oy
Jämsäentie 5
01610 Vantaa, Finland
Phone: +358 (0)96824774, Fax: +358 (0)96824775
E-Mail: asiakaspalvelu@karlstorz.fi

KARL STORZ SE & Co. KG
Representation Office
Kęstutis st. 59 / Lenkijai st. 27
06124 Vilnius, Lithuania
Phone: +370 5 272 0448
Mobile: +370 685 87 000
E-Mail: info-it-ly@karlstorz.com

KARL STORZ Endoskop Danmark A/S
Skovlystvej 33
2840 Holte, Denmark
Phone: +45 45162600, Fax: +45 45162609
E-Mail: marketing@karlstorz.dk

KARL STORZ Endoscopy (UK) Ltd.
415 Perth Avenue, Slough
Berkshire, SL1 4TG, United Kingdom
Phone: +44 1753 503500, Fax: +44 1753 578124
E-Mail: info-uk@karlstorz.com

KARL STORZ Endoscopy Nederland B.V.
Displayweg 2
3821 BT Amersfoort, Netherlands
Phone: +31 (0)33 4545890
E-Mail: info-nl@karlstorz.com

KARL STORZ Endoscopy Belgium N.V.
Phone: +31 (0)33 4545890
E-Mail: info-be@karlstorz.com

KARL STORZ Endoscopia France S.A.S.
12, rue Georges Guynemer, Quartier de l'Europe
78280 Guyancourt, France
Phone: +33 1 30484200, Fax: +33 1 30484201
E-Mail: marketing-fr@karlstorz.com

KARL STORZ Endoskop Austria GmbH
Landstraßer Hauptstr. 148/1/G1
1030 Wien, Austria
Phone: +43 1 71 56 0470, Fax: +43 1 71 56 0479
E-Mail: storz-austria@karlstorz.at

KARL STORZ Endoscopia Iberica S. A.
Parque Empresarial San Fernando
Edificio Munich - Planta Baja
28930 Madrid, Spain
Phone: +34 91 6771051, Fax: +34 91 6772981
E-Mail: info-es@karlstorz.com

KARL STORZ Endoscopia Italia S. r.l.
Via dell'Artigianato, 3
37135 Verona, Italy
Phone: +39 045 8222000, Fax: +39 045 8222001
E-Mail: info-it@karlstorz.com

KARL STORZ Adria Eos d.o.o.
Capriska 6
10000 Zagreb, Croatia
Phone: +385 1 6406 070, Fax: +385 1 6406 077
E-Mail: info@karlstorz.hr

KARL STORZ Endoscopy d.o.o.
Cesta v Gledce 34b
1000 Ljubljana, Slovenia
Phone: +386 1 620 5880, Fax: +386 1 620 5882
E-Mail: ptuama@karlstorz.si

KARL STORZ Endoskop Polska Marketing Sp. z o.o.
ul. Bojkowska 47
44-100 Gliwice, Poland
Phone: +48 32 706 13 00, Fax: +48 32 706 13 07
E-Mail: info-pl@karlstorz.com

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TOO KARL STORZ Endoscopy Kazakhstan
Saryarka, B, BC "Arman", off. 910
010000 Astana, Republic of Kazakhstan
Phone: +7 7172 552-549, 552-788, Fax: -444
E-Mail: info@karlstorz.kz

KARL STORZ ENDOSCOPE
East Mediterranean & Gulf (branch)
Building West Side 7A - Unit 7WA - 3008
Dubai Airport Free Zone, P.O. Box 54983
Dubai - United Arab Emirates
Phone: +971 (0)4 2958887, Fax: +971 (0)4 3205282
Service Hotline: +971 (0)4 3415882
E-Mail: info-gne@karlstorz-emg.com

KARL STORZ Endoscopy India Private Limited
11th Floor, Dr. Gopal Das Bhawan
28, Barakhamba Road
New Delhi 110001, India
Phone: +91 11 4374 3000, Fax: +91 11 4374 3010
E-Mail: corporate@karlstorz.in

KARL STORZ SE & Co. KG
Resident Representative Office
14th Floor, M Plaza Saigon
38 Le Quan, District 1
Ho Chi Minh City, Vietnam
Phone: +84 8 3823 8000, Fax: +84 8 3823 8039
E-Mail: infovietnam@karlstorz.com

KARL STORZ Endoscopy Australia Pty. Ltd.
15 Orion Road Lane Cove NSW 2066
P.O. Box 50 Lane Cove NSW 1595, Australia
Phone: +61 (0)2 9480 8700
Toll free: 1800 996 562 (Australia only)
Fax: +61 (0)2 9420 0685
E-Mail: info@karlstorz.au

KARL STORZ Endoscopy China Ltd., Hong Kong
Room 2503-05, 25F AXA Tower, Landmark East,
No. 100 How Ming Street, Kowloon, Kowloon,
Hong Kong, People's Republic of China
Phone: +852 28 65 2411, Fax: +852 28 65 4114
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd.
Beijing Branch,
Room 1805-1807, Building B, 18F Beijing JFC
No. 8, Jiangnanmenwai Street, Chaoyang District,
100022, Beijing, People's Republic of China
Phone: +86 10 5638188, Fax: +86 10 5638199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shanghai Branch
36F, Tower 1 Grand Gateway
No. 1 Hongqiao Road, Xuhui District,
200030, Shanghai, People's Republic of China
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Chengdu Branch
Room 805, 8F Jin Jiang International Building
No. 1 West Unjiang Road, Wuhou District,
610041, Chengdu, People's Republic of China
Phone: +86 28 86587977, Fax: +86 28 86587975
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shenyang Branch,
Room 2001-2005, 20F N-MEDIA International Center,
No. 187 Youth Avenue, Shenhe District,
110014, Shenyang, People's Republic of China
Phone: +86 24 23181118, Fax: +86 24 23181119
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Guangzhou Branch
Room 02B, 35F Teem Tower,
No. 208 Tianhe Road, Tianhe District,
510620, Guangzhou, People's Republic of China
Phone: +86 20 8732-1251, Fax: +86 20 8732-1286
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd.
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com

KARL STORZ Endoscopy Singapore Sales Pte Ltd
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com

KARL STORZ Endoscopy Korea Co. Ltd.
9F Hyowon-Building
97, Jungdae-ro, Songju-gu
05719 Seoul, Korea
Phone: +82-70-4350-7474, Fax: +82-70-8277-3299
E-Mail: infokorea@karlstorz.com


KARL STORZ Endoscopy Taiwan Ltd.
12F, No. 192, Sec. 2, Chung Hsin Rd.,
Sindian District, New Taipei City, Taiwan
Phone: +886 933 014 160, Fax: +886 2 8672 6399
E-Mail: info-tw@karlstorz.com

KARL STORZ SE & Co. KG
Representative Office Philippines
1901 Picadilly Star Bldg., 4th Avenue, BGC
Taguig City 1636, Philippines
Phone: +63 2 317 45 00, Fax: +63 2 317 45 11
E-Mail: philippines@karlstorz.com

KARL STORZ Endoscopy Japan K.K.
Stage Bldg. 8F, 2-7-2 Fujimi
Chiyoda-ku, Tokyo 102-0071, Japan
Phone: +81 3 6360-8622, Fax: +81 3 6360-8633
E-Mail: info@karlstorz.co.jp

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