
Operational documentation of the manufacturer: Viora Ltd.

**Operational documentation of the
manufacturer for medical device:**

Vtouch

CE 0344

MK-096:EN

Operational documentation of the manufacturer

Copyright © 2012 Viora. All rights reserved.

Viora reserves the right to make changes to its products or specifications to improve performance, reliability or manufacturability. Information furnished by Viora is believed to be accurate and reliable. However, Viora assumes no responsibility for its use. No license is granted by its implication or otherwise under any patent or patent rights of Viora.

No part of this document may be produced or transmitted in any form or by any means, electronic or mechanical, for any purpose, without the express written permission of Viora.

Viora has patents and pending patent applications, trademarks, copyrights, or other intellectual property rights covering subject matter in this document. The furnishing of this document does not give you any license to these patents, trademarks, copyrights, or other intellectual property rights except as expressly provided in any written agreement from Viora.

Specifications are subject to change without notice.

P/N: MK-096 EN

Release Date: December 2012

- Before using the system, please check local regulations. If any local legislation is violated, use cannot be authorized.
- In the United States, federal law restricts prescription medical devices, like V-touch, to sale by or on the order of a physician, or properly licensed practitioner. Viora makes no representations regarding federal, state, or local laws or regulations that might apply to the use and operation of this system.

Manufacturer's Address

Viora Ltd.
3 Maskit Street
Herzliya, Israel 46733
Tel: +972 9955 1344
Fax: +972 9955 1345
Website: www.vioramed.com
E-mail: info@vioramed.com

Manufacturer's Authorized Representative

CEPartner4U BV
Esdoornlaan 13
3951 DB Maarn
The Netherlands
Tel: +31 (0) 343 442 52
Fax: +31 (0) 343 442 162

This device complies with the European Medical Device Directive 93/42/EEC

TABLE OF CONTENTS

1	Before You Start	5
1.1	Conventions Used in the Manual	5
	Explanation of the Symbols Used in the System	6
2	Introduction	6
2.1	Regulating data for the construction of a medical device	7
2.1.1	System Console	8
2.1.2	Handpieces	8
2.1.3	Disposable Sterilized Tip – FR handpiece	8
2.1.4	Display Unit	9
3	Safety	9
3.1	Introduction	9
3.2	Technical requirements for medical device	10
3.2.1	System Safety Features	17
3.3	Potential Fire Hazards	18
3.4	System Labels	18
3.5	System Classifications	19
3.6	Data for development of medical device	19
3.7	Data for the manufacturing of medical device	20
3.8	Patient and Personnel Safety	21
3.9	Safety when Operating the System	21
3.10	Electrical Requirements	22
3.11	Environmental Requirements	22
3.12	Safety in Maintaining the System	22
3.13	Transport and Storage Conditions	23
4	System Installation	23
4.1	Equipment List	23
4.2	The using data for the medical device	23
4.3	Moving the System within the Facility	25
4.4	Moving the System to Another Facility	25
5	The data for the operating medical device	26
5.1	ST Treatment Menu	26
5.2	ST Basic Treatment Screen	27
5.2.1	ST Advanced Treatment Screen	29
5.3	FR Treatment Menu	29
5.4	Set Up Menu	31

Operational documentation of the manufacturer

5.5	System Messages	32
6	Clinical Guide	34
6.1	Indications for use	34
6.2	Contraindications	34
6.3	Possible Adverse Effects	35
6.4	Realistic Expectations	36
6.5	Pre-Treatment	36
6.5.1	Consult Meeting	36
6.5.2	Treatment Parameters	37
6.6	FR Clinical Guide	38
6.6.1	General	38
6.6.2	Pre-Treatment Preparation	38
6.6.3	Test Procedure	39
6.6.4	General treatment Guidelines	40
6.6.5	Treatment Schedule	41
6.6.6	Treatment Follow-Up	41
6.6.7	Post-Treatment Care	42
6.7	ST Clinical Guide	42
6.7.1	General	42
6.7.2	Test Procedure	43
6.7.3	General treatment Guidelines	43
6.7.4	Treatment Schedule	43
6.7.5	Post-Treatment Care	43
7	Data for the technical maintenance of medical device	43
7.1	Cleaning the System	44
7.2	Cleaning the ST Electrodes	44
8	Data for the service of a medical device	44
9	System Specifications	46
10	System composition	52
11	Data for recycling of medical device	52

1 BEFORE YOU START

Warning:

- Read this manual to become familiar with all safety requirements and operating procedures before attempting to operate the system.
 - Any intense radio frequency (RF) device can cause injury if used improperly.
 - High voltage is present inside the system. Always be aware of the possible dangers and take proper safeguards as described in this manual.
 - The system and its components must be serviced only by qualified personnel.
-

1.1 CONVENTIONS USED IN THE MANUAL

The following conventions in the form of notes, cautions and warning are used in this manual:

Note:

The content of this Note offers general information that is important to keep in mind.

Caution:

A Caution alerts the user to the possibility of a problem with the device associated with its use or misuse. Such problems include system malfunction, failure, and damage to the system or other property. The caution statement includes a precautionary action that should be taken to avoid the hazard.

Warning:


A Warning alerts the user to the possibility of injury, death, or serious adverse reactions associated with the use or misuse of the system.

Note:

In many places throughout this manual you will see words or phrases in the text which are bolded in a sans-serif typeface, such as this. These are words and phrases that you will always see on the LCD display of the control panel.

EXPLANATION OF THE SYMBOLS USED IN THE SYSTEM

Symbol	Description
	Attention: Consult Accompanying Document
	CE Compliance Symbol
	CSA Compliance Symbol (244308 CSA master contract number, project 2137900)
	Dispose of according to WEEE regulations
	Fuse
	Type BF Equipment
	Manufacturer
	Non-ionizing electromagnetic radiation
	Follow instructions for use
	CAUTION: Federal law restricts this device to sale by or on the order of a physician

2 INTRODUCTION

The V-touch system developed by Viora is a radiofrequency (RF) device designed for aesthetic and dermatological procedures, with applications for the face, neck, chest and hands. The system is based on well-established technology of bipolar radiofrequency, which is extensively used for the coagulation and ablation of soft tissues. Over the past decade, RF technology has become widely used in the aesthetics and dermatology field for skin tightening and body shaping.

The V-touch is the new Fractional RF device which is designed to direct RF energy exclusively to damaged areas while leaving healthy surrounding tissue intact. It takes the production of new collagen in the skin to a new dimension. Through its grid of electrodes, the V-touch treats only a fraction of skin during each session, causing the surrounding healthy tissue repair the targeted lesions using the body's natural healing processes and replacing new cells. This leads to remarkably visible results in skin texture, acne scars, fine lines and skin tone. Treatment tips are disposable, sterilized and allocated for one patient use only.

The V-touch is in Class IIb according to Directive 93/42/EEC of 14 June 1993, amended by Council Directive 2007/47/EC of 5 September 2007, concerning medical devices, Annex IX, Rule 9, "all active therapeutic devices intended to administer or exchange energy are in Class IIb".

Note

Only personnel who have undergone and completed a thorough training course are permitted to operate the system.

The intended users are aesthetic physicians or dermatologists. Eligibility of users shall be in accordance with local market regulations.

2.1 REGULATING DATA FOR THE CONSTRUCTION OF A MEDICAL DEVICE


Figure 1: System Console


Figure 2: FR Hand piece

Figure 3: FR tip

Figure 4: ST Handpiece

The System consists of the following components:

- System console
- FR (Fractional RF) handpiece
- ST handpiece
- Disposable, sterilized tips
- User Manual
- Marketing kit

2.1.1 SYSTEM CONSOLE

The system console houses the following components:

- Power Supply
- RF Generator
- Main CPU
- Display Unit – touch screen

2.1.2 HANDPIECES

The handpiece is attached to the main console through an umbilical cord.


Figure 5: System's rear panel

The FR handpiece is powered by a 50 W, 1 MHz RF generator.

2.1.3 DISPOSABLE STERILIZED TIP – FR HANDPIECE

A disposable tip, containing a matrix of bi-polar electrode pins attaches to the end of the handpiece and placed on the skin during treatment. The distance between the pins is 1 mm and different size tips have up to 8 rows and 7 columns of pins.


Figure 6: Tip Matrix

A thermoelectrically cooled bracket surrounding the grid area is designed to cool surrounding skin and reduce pain during treatment.

2.1.4 DISPLAY UNIT

The display unit's touch screen mainly presents the following information:

1. Treatment type
2. Energy
3. Pulse counter
4. Treatment time

3 SAFETY

This chapter describes the safety issues regarding the use and maintenance of the V-touch system.

3.1 INTRODUCTION

To ensure dependable, consistent and safe treatments, the following safety measures must be taken into consideration when operating and maintaining the system. The person operating the system must be familiar with the safety regulations provided in this chapter. The user manual must be read completely before operating the V-touch. It is highly important that the safety chapter is carefully read before attempting to install or operate the system. If any sections of the safety chapter are unclear, or if any questions arise prior to, or during installation and/or operation of the handpiece, contact your sales representative or distributor.

Viora assumes no liability whatsoever for any damage or injury caused as a result of actions taken which are not in strict accordance with the instructions provided in this manual. Improper use may result in injury and the primary consideration when operating the system must be to maximize safety for both patient and operator.

Warning

- It is important that the operators of this system familiarize themselves with all safety requirements and operating procedures prior to operating the

system

- Any intense RF device may cause injury if used improperly.
- High voltage is present inside the system.
- Always be aware of the possible dangers and take proper safeguards as described in Section 3.7 of this user manual.
- For complete contact information please refer to Page 2 of this user manual.

3.2 TECHNICAL REQUIREMENTS FOR MEDICAL DEVICE

Power Source

The device is powered by an AC power source, 230 V 50 Hz or 115 V 60 Hz $\pm 10\%$.

Console specification

Power entry

Power inlet

Type: IEC Medical Grade, Illuminated

Rating: 5A

Fuse: 2 X (250 VAC, 4A, 5X20 mm)

Line Filter

Grade: Medical

Rating: 10A

Power Cable

Compliance with standard: UL817

Power supply

Grade: Medical

Input Voltage: 90-254 VAC, 50-60 Hz

Output Voltage: 12V, 24V

Output Power: up to 750W

User Interface

On/Off switch on power inlet

RF release Trigger on each H/P, and on Foot Pedal

GUI: Graphical User Interface

LCD:

Size: 8"
Resolution: 800X600
Communication: LVDS
Back light: LED
Touch method: Resistive
Controller: Freescale IMX53

Main CPU:

MC Working Voltage: ~3V
MC Consumption: 71 mA @ 3V
Clock: external 25 MHz

External software plug

Type: USB Dongle
Use: Technician, Production, SW update (for GUI and MAIN)
Memory: 1 GB

RF Board

Inductor type: Flat
Input Voltage: 48 V
Output: up to 500 [W]
Pulse width range: 10-200 [mSec]
Pre-Pulse width: 80 [mSec]
RF Frequency: 0.8, 1, 1.7, 2.54 [MHz]
Calibration: via SW

Fans

Back Fan:
Dimensions: 80X80X25 mm.
CFM: 37
Operating Voltage: 12 V
Consumption: 1.4 W
Noise: 30 dBA

RF Fan:

Operational documentation of the manufacturer

Type: Heat Sink mounted

Fan Dimensions: 80X80X25 mm

CEM: 27-16

Operating Voltage: 12V

Consumption: 1.32W

Noise: <27.4 dbA

H/P Connectors:

No. of: 2

Type: Female (socket)

Type: Circular

Materials:

Side panels: Aluminium

Chassis: Aluminium

Front cover: ABS

H/P cradle: ABS

Door Cover: ABS

Feet: Rubber

Console Feet:

No. of: 4

Type: Rubber

Cart Specifications:

Wheels:

No. of: 4

Type: 2 rear, 2 front with brakes

Diameter: 7.5 cm

Door hinges:

Type: 90 Deg. Concealed

Springs: No Springs

ST H/P Specification

Materials

Shell material: ABS

Electrodes material: Electrolytic tough pitch Copper, UNS C11000, Hard Chrom plated

Cable

Type: Servo, Shielded

RF Wires: 2X20 AWG, twisted pair, shielded

Signal Wires: 4X18 AWG

Rated Voltage: 1000V

Tested: 4000V

Rating: UL

TEC

Function: Electrodes cooling

Size: 15x30X3.5 mm

Power: 8W (1 units)

Voltage: 4 V

Current: 2 A

Delta Temp: 67 Celsius max

Temperature Sensor

Range: -55° to 150° Celcius

Accuracy: $\pm 0.25^\circ$ Celcius

CPU unit

Working Voltage: 3V

Consumption: 250 mA @ 3V

Clock: internal 8 MHz

Cooling fan

Dimensions: 25X25X10 mm

CFM: 3

Operating Voltage: 12 V

Consumption: 0.45 W

Noise: 16 dBA

Operational documentation of the manufacturer

Memory chip

Type: EEPROM

Size: 1Kbit

Treatment Area: 8mmX8mm

FR H/P Specification

Operation parameters

See section 3.9

Materials

H/P Shell material: ABS

Cable

Type: Servo, Shielded

RF Wires: 2X20 AWG, twisted pair, shielded

Signal Wires: 4X18 AWG

Rated Voltage: 1000V

Tested: 4000V

Rating: UL

Disposable Tip: Regular

Electrode Matrix

No. of electrodes: 56 (7 X 8)

Electrode's pitch: 1.3 mm

Single electrode diameter: 0.3 mm

Area of treatment: 14X15.5 mm

PCB material: FR4 2.4mm thick

Connection to H/P

Connector style: contact

Identification: EEPROM

Shell Material: Santoprene TPV 281-87 med

Sterilization: EtO Gas

Vacuum Pump

Type: Diaphragm
Motor: Brushless
Voltage: 12V
Flow Rate: 1.5 L/Min @ 0 mbar, 0.6 @ -300 mbar
Max Vacuum: -500 mbar
Consumption: ~220 mA @ 12V

Solenoid

Function: build and release vacuum at Tip
Function: 3
Operation Voltage: 3V

TEC

Function: Cooling the skin via the electrodes' array
Size: 15X15 mm
Voltage: 3.8V
Current: 4A
DT max: 66 Celsius

CPU unit

Working Voltage: 3V
Consumption: 250 mA @ 3V
Clock: internal 8 MHz

Fan

Dimensions: 30X30X10 mm
CFM: 4.7
Operating Voltage: 12 V
Consumption: 1.0 W
Noise: <23 dBA

Temperature Sensor

Range: -55° to 150° Celsius
Accuracy: $\pm 0.25^\circ$ Celsius
Memory chip:
Type: EEPROM
Communication: SPI
Size: 1K

Operational documentation of the manufacturer

Foot Pedal

Type: Momentary close
Cable Termination: DB9 Male
Sealing Rating: IPX2

Shelf Life

Console shelf life: 5 years

Operation Profile

Operation limitations:
ST H/P work minutes limit: 65,000
Disposable (regular) Tip Pulse limit: 600 pulses

RF Frequencies of operation

FR H/P: 1 MHz
ST H/P: 0.8 MHz, 1.7 MHz, 2.45 MHz

RF Power output

For ST H/P: 0-130 J/cm³
For FR H/P: 0-10 Joule
Pre-Pulse ST: 8 Watt
Pre-Pulse FR: 4 Watt
FR pre-heating (Smart Heating) pulse: 9 Watt or 5 J/cm³

RF pulse width:

For ST H/P: 50-200 [mSec]

For FR H/P: 10-100 [mSec]

Pre-Pulse width: 80 [mSec]

Pre-heating pulse: 100 [mSec]

Switches

The V-touch console has an ON/OFF switch (on the inlet).

Each H/P has a trigger for releasing RF pulses.

The V-touch has an optional footswitch, with the same functionality as the trigger.

Indicators

The V-touch is provided with an LCD Graphic display, showing the different parameters and operation modes of the system.

The V-touch system has a speaker for audio indications.

The FR applicator is provided with a LED indicator to indicate RF Pulse release.

Physical Dimensions

System dimensions: 100HX54WX46D

3.2.1 SYSTEM SAFETY FEATURES

The V-touch system incorporates the following safety features:

- Default treatment parameters are automatically recommended when operator makes a treatment selection.
- The hand piece is equipped with a trigger; if the trigger is not pressed, no energy will be emitted.
- The system automatically tests for proper coupling prior to any pulses emitted by releasing a short sub pulse at low voltage. If coupling between electrodes and skin is insufficient the RF energy will be disabled.
- The hand piece and treatment tip must be properly connected to the system in order for the treatment mode to be activated.
- The system can be switched to a 'Standby' mode disabling the release of energy.

Warning

Objects such as superficial metal or conductive implants, rods, plates or pins in the treatment areas are contraindicated to the treatment.

In case of improper coupling, indicated on screen and by audio signal the trigger

should be released immediately and the hand piece should be re-oriented on the skin. In rare cases superficial burns may occur from the heat generated by the RF electrodes.

3.3 POTENTIAL FIRE HAZARDS

- Do not use the system in the presence of explosive or flammable materials.
- Do not use flammable substances when preparing the skin for treatment.
- If alcohol is used for disinfecting the hand piece, it must be completely dry before the system is used.

3.4 SYSTEM LABELS


Figure 7: Labels Affixed to the System

Figure 8: System Labels.

3.5 SYSTEM CLASSIFICATIONS:

- Electric shock protection: Class I; Type BF
- Not suitable for use in presence of flammable substances.
- Protection against ingress of liquids: Ordinary equipment
- Mode of operation.—continues

3.6 DATA FOR DEVELOPMENT OF MEDICAL DEVICE

LABELING

ISO 15223-1:2012-symbols to be used with medical device labels.

BIOCOMPATIBILITY

- a. AAMI/ ANSI / ISO 10993-1:2009 - Biological evaluation of medical devices -- Part 1: Evaluation and testing
- b. ISO:10993-5:2009 Test for in vitro cytotoxicity
- c. ISO 10993-10:2009 Tests for sensitization and irritation

CLINICAL

- a. EN ISO 14155:2011- Clinical Investigation of Medical Devices for Human Subjects

STERILIZATION (APPLICABLE FOR THE DISPOSABLE TIP ONLY).

- a. EN ISO 11135-1:2007 - Sterilization of health care products – Ethylene oxide – Requirements for development, validation and routine control of sterilization process for medical devices
- b. ISO-TS 11135-2:2008 - Sterilization of health care products – Ethylene oxide – Guidance on the application of ISO 11135-1
- c. EN 556-1:2001 - Sterilization of medical devices; Requirements for medical devices to be designated "STERILE".
- d. ANSI/AAMI ST72: 2011, Bacterial Endotoxin – Test methodologies, routine monitoring, and alternatives to batch testing.

RISK ANALYSIS

- a. EN ISO 13971:2012 Medical devices Application of risk management to medical devices
- b. EN IEC 60601-6:2010 Medical electrical equipment - part 1-6 general requirements for basic safety and essential performance - collateral standard usability

QA

- a. EN ISO 13485:2003 Medical devices Quality management systems Requirements regulatory purposes

ELECTRICAL SAFETY AND EMC

- a. IEC60601-1 - Ed 3.0: 2006 Medical electrical equipment - Part 1: General requirements for basic safety and essential performance
- b. IEC60601-1-2:2007 Ed 3.0

Medical electrical equipment - Part 1-2: General requirements for safety Collateral Standard: Electromagnetic compatibility Requirements and tests

- c. FCC part 18 - Radiation

3.7. DATA FOR THE MANUFACTURING OF MEDICAL DEVICE

The system's PCBs, as assembled by the Sub-contractor, include use of the following technologies:

- a. Surface Mount Technology
- b. Reflow soldering
- c. Through-hole

The V-touch plastic will be manufactured in the following technologies:

- a. High pressure Molds
- b. CNC

The V-touch metal sheet will be manufactured in the following technologies:

- a. Punch
- b. Laser Cut

ASSEMBLY PROCESS

The V-touch system will be assembled according to procedure by Viora.

TESTING PROCESS:

RF board, Main board, UI Board, Applicators Board shall be tested in a designated tester/Jig for each board. Assembly using customary techniques.

Each system will be tested under the Acceptance Test Procedure.

COLORS:

The V-touch visible components color are:

Console – White Shell

V-touch H/P – RAL 9003 + RAL 7024

ST H/P – RAL 9003 + RAL 7024

Trolley – Pearl

Wheel/base – White Shell

Front panel – Black

Side Panels – Golden Silver

3.8. PATIENT AND PERSONNEL SAFETY

When handling the V-touch system, the primary concern should be the safety of both operator and patient. Carefully read the following safety guidelines:

- Before beginning any treatment course, the operator should complete a detailed patient history to verify that the patient is eligible for the treatment and is in a healthy condition.
- The operator should provide the patient with information about the treatment, the treatment protocol, expected results and any risks or side effects associated with the treatment.
- Only trained operators may perform treatments.
- A test procedure must be performed before any treatment course and at any new anatomical area, in order to evaluate the treatment parameters.

Warning

- Improper treatment may cause local burns.
 - Untrained operators should not operate the system.
-

3.9 SAFETY WHEN OPERATING THE SYSTEM

- Keep all system covers closed at all times and do not touch any of the system's inner components. Only authorized service technicians are permitted to remove the panels.
- Always make sure that the disposable tip is replaced before a new treatment session begins. In no case should a disposable tip be used for more than one treatment.
- Avoid operating the system in the presence of potentially explosive or flammable materials.
- Always make sure that all operating personnel have been properly trained and are familiar with the system's controls and know the proper steps on how to shut down the system instantly when required.

Operational documentation of the manufacturer

- Only aesthetic physicians and dermatologists are allowed to operate the system.
- No part of the system should be immersed in liquid.
- Medical electrical equipment can be affected and experience interference by portable and mobile RF communications equipment.
- Only original and specified accessories and cables should be used. Using any other equipment may result in increased emissions or decreased durability of the system.

Warning

- It is forbidden to connect any third party equipment to the system without written approval from Viora.

3.10 ELECTRICAL REQUIREMENTS

The V-touch system is preset to accommodate the local line voltage according to the customer's order. Accordingly, the system will require a separate line supply of one of the following (appropriate to customer's country):

Single phase 230 Vac, 4 A, 50 Hz

or

Single phase 115 Vac, 6.3 A, 60 Hz

Warning

- To protect against risk of fire, replace the fuse with the same type and rating only.
- Proper grounding is essential for safe operation. The system is grounded through the power cable.

3.11 ENVIRONMENTAL REQUIREMENTS

- Electronic parts can be damaged by corrosive materials. System must operate in a non-corrosive environment.
- For optimal operation of the system, maintain room temperature between 15°-26°C (59°-79°F) and relative humidity of less than 80%.
- If the system has been stored in an environment cooler than the above mentioned temperatures, allow it time to warm up to room temperature before the first activation.
- Metallic dust is potentially harmful to electrical equipment.

3.12 SAFETY IN MAINTAINING THE SYSTEM

- Only authorized service technicians are qualified to open the system's panels to service the system.

- Shut down the system and disconnect the power cable before performing any of the noted maintenance procedures in Section 7 of this user manual.

3.13 TRANSPORT AND STORAGE CONDITIONS

The optimal transport and storage conditions of the system are:

- Ambient temperature: 5°C to 55°C.
- Relative humidity: from 35% to 80%, non condensing.
- The packaging protects the V-touch system from breakage or damage according to the company's packaging instructions.
- Atmospheric Pressure Range: 500hPa – 1060hPa

4. SYSTEM INSTALLATION

4.1 EQUIPMENT LIST

The V-touch system includes the following:

- System console
- ST Handpiece
- FR handpiece
- Disposable tips
- Foot pedal (optional)
- Left & right cradles with attachment screws
- Power cable
- User manual
- Marketing Kit

4.2 THE USING DATA FOR MEDICAL DEVICE

After unpacking the V-touch system, perform the following installation procedures:

- Installing the cradle: the cradle has three captive screws. To attach the cradle to the system's main unit, place the cradle in front of the holes on either side of the panel, with the opening side facing the front (left cradle on the left side; right cradle on the right side), and fasten the screws. Make sure the cradles are firmly attached to the main unit before proceeding to the next step.

Figure 9: Installing the cradle

- Connecting the handpiece: unpack the handpiece from its box. Place the handpiece in the cradle. Align the pins on the HP connector to the connector opening in the back of the system. Rotate the connector until it fits into place and lock it by rotating the connector's

Operational documentation of the manufacturer

secure ring. Place the cable inside the cable clamp on the upper side of the rear panel. Close the clamp to secure the cable.


Figure 10: Connecting the hand piece.

Figure 11: Rotating the connector.

Note

There are 2 HP connectors on the back of the system. Both connectors are identical and can serve either the FR or ST hand piece.

- c) Connecting the foot pedal: connect the foot pedal cable to the connector on the rear panel and fasten the screws to secure it into place.


Figure 12: Connecting the foot pedal.

Note

- In order to avoid a potential accident, please store the foot pedal in safe place when not in use.
- Foot pedal must not get wet. Do not keep any liquids in the vicinity of the foot pedal.

- d) Connect the power cable to the power connection port on the system's rear panel and plug the electrical cable into an appropriate electrical outlet.

Figure 13: Connecting the power cable.

- e) To connect the FR tip to the handpiece, open the supplied box and remove one sterile FR tip. Peel off the seal and place the tip on the handpiece, without removing the tip from the blister. Tip can be connected directly to the handpiece.

Figure 14: FR Tip connection.

- f) The system is equipped with wheel brakes. To complete installation – please make sure that the wheel brakes are secured and locked.

4.3 MOVING THE SYSTEM WITHIN THE FACILITY

- Make sure the system is turned off and the power cable is disconnected.
- Make sure that the wheel brakes are released to allow movement.
- Use the handle located at the back of the system to carefully move it around.

4.4 MOVING THE SYSTEM TO ANOTHER FACILITY

- Make sure the system is turned off and the power cable is disconnected.
- Make sure that the wheel brakes are released to allow movement.
- Use the handle located at the back of the system to carefully move it around.
- Slowly load the system onto the vehicle. Make sure the system is secured in place and proper padding is used to protect the system while traveling.
- To unload the system in the new facility, carefully remove it from the vehicle and use the handle on the back of the system to place it in its new location.
- Reinstall the system.

Note

It is recommended that the system only be moved when it is packaged in its

5 THE DATA FOR THE OPERATING OF MEDICAL DEVICE

To turn the system on, press the ON/OFF button located at the back of the system.

Following start-up, the system will display the main menu screen.


Figure 15: Main Menu Screen.

Note

If both handpieces are attached, choose the desired treatment by pressing the treatment selection key: ST or FR treatment on the main menu screen. When only one handpiece is attached, other handpiece selection will not be available.

Note

Please make sure that the handpiece in use corresponds with the displayed treatment screen.

5.1 ST TREATMENT MENU

The ST Treatment menu offers four different treatment programs; one for body treatment and three for facial treatments (forehead, cheeks, neck and décolleté). Choose the desired treatment area by pressing the touch screen.

From the treatment menu, the operator may access treatment tips by pressing the  Tips key. Treatment tips are divided into three categories: General, ST and FR. To return to the treatment menu, operator should press 'close'.

5.2 ST BASIC TREATMENT SCREEN


Figure 15: ST Basic Treatment Screen

A. **Mode** - The system has four different RF treatment modes that can be controlled by the operator. The modes are based on the following integrated technologies:

- Mode 1 – Low RF frequency of 0.8MHz (suitable for deep treatments).
- Mode 2 – Medium RF frequency of 1.7 MHz (suitable for mid-layer treatments).
- Mode 3 – High RF frequency of 2.45 MHz (suitable for superficial treatments).
- Mode 4 – Combines all three frequencies for multi layer treatments

To choose the desired treatment mode, press the mode number

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Note

Changing modes will move the system to standby. Press the ready/standby key to resume work.

- B. Energy:** To control the energy level, press the energy (+) or (-) keys. The smaller (+) / (-) keys will change energy by 2J. The larger (+) / (-) keys will change energy by 10J.
- C. Hour counter:** Counts total treatments time usage on the handpiece.
- D. Minute counter:** Counts time usage for current treatment. Press the  key to zero the counter.
- E. Cooling indicator:** The RF electrodes are cooled during treatment. The cooling icon  indicates that cooling is activated. If a pulse is not released within 60 seconds, cooling will stop and the system will move to Standby mode. In order to resume treatment, press the  key.
- F. Home:** By pressing the  key, operator will be prompted to main menu.

- G. Menu:** By pressing the  key, operator will be prompted to treatment menu.

- H. Basic / Advanced:** The colored key, with bold arrow facing down, shows operator's current mode of operation (basic or advanced). 

- I. Standby / Ready:**  indicates that the system is ready for work and pulses can be emitted by pressing the trigger on the handpiece or the foot pedal. System will move to  mode if:

- Operator changes mode of treatment.
- Operator changes energy by 10J.

To resume work, press the standby key. Wait for the cooling to restart (indicated by the progress bar). When pressing the key, make sure the trigger is not pressed.

To release an RF pulse, press the trigger located on the handpiece OR the foot pedal.

Note

- Pressing the trigger or foot pedal continuously will emit 5 consequent RF pulses followed by a warning screen. Operator must press OK in order to resume work.
- Pressing the trigger or foot pedal accidentally with no contact on the skin will cause bad coupling alert on the screen.

Note

Follow the treatment movements as per the arrow indications on the image.

5.2.1 ST ADVANCED TREATMENT SCREEN


Figure 16: ST Advanced Treatment Screen

The advanced treatment screen is for experienced users and allows control of additional treatment parameters.

- J. **Pulse Duration** - Determines the length of the pulse during which energy is released (On time).
- K. **Pulse Interval** - Determines the pause between pulses (Off time).

Warning

Changing pulse duration is directly related to amount of energy released onto the treatment area.

5.3 FR TREATMENT MENU

There are three pre-set treatment programs for FR treatments:

- Shallow
- Medium
- Deep

Once a treatment program is selected, the operator is prompted to the treatment screen.


Figure 17: FR Basic Treatment Screen

- A. **Vacuum** - Indicates whether vacuum feature is enabled. Operator may switch the vacuum feature on and off as desired by pressing the vacuum key.

**Note:**

Vacuum ensures optimal coupling between skin and electrodes and is recommended to use at all times.

- B. **Energy** - Control the desired RF energy using the energy (+) / (-) keys.
- C. **Tip pulse counter** - Counts total pulses with inserted tip.
- D. **Pulse counter** - Counts emitted pulses for current treatment. Press the  key to zero the counter.
- E. **Cooling** - Indicates whether cooling feature is enabled. Operator may switch the cooling feature on and off as desired by pressing the cooling key.


Note

Cooling reduces pain and makes the treatment more comfortable to the patient, hence recommended to be used at all times.

F. **Home** - By pressing the  key, operator will be prompted to main menu:

G. **Menu** - By pressing the  key, operator will be prompted to treatment menu.

H. **Basic / Advanced** - The colored key with bold arrow facing down, shows operator's current mode of operation (basic or advanced) .

I. **Standby / Ready** -  indicates that the system is ready for work and pulses can be emitted by pressing the trigger on the handpiece or foot pedal. System will move to  mode if:

- Operator changes energy of treatment.
- Pulse is not emitted for 2 minutes or longer.

To resume work, press the 'standby' key. Wait for the cooling to restart (indicated by the progress bar).

To release an RF pulse, press the trigger located on the handpiece OR the foot pedal.

5.4 SET UP MENU

- **Language** - Choose the desired language setting by pressing the language name. To confirm selection, press 'OK'.
- **Treatment preferences** - Allows operator to set a default mode for treatment screen. Operator may always switch between basic and advanced mode through the treatment screen or by resetting the default preferences.
- **Treatments History** - Shows details about past treatments:
 - Treatment date
 - Duration
 - Handpiece used
 - Handpiece serial number
 - Total pulses used
- **Events History** - Shows history of system events such as self test failures.
- **System Information** - Shows indications for:
 - Serial numbers of the system and connected hand pieces
 - Software versions
 - Usage data for connected hand pieces and tips

Setup		System Information	
Serial Number		Software Version	
System	1.0	GUI	1.2
ST Handpiece	1.0	Main	1.0
FR Handpiece	1.0	ST Handpiece	2.0
FR Tip	3.2	FR Handpiece	2.0
Usage Data		Boot loader	2.0
Remaining ST pulses	1,200	Operating system	1.0
Remaining pulses on FR Tip	300		
No. of FR Tips used	3,200		

Figure 18: System Information Screen

5.5 SYSTEM MESSAGES

The system continuously performs self-test procedures relating to the hand pieces attached and the internal operation. Status or error messages will appear on the screen and display the relevant information when required. To continue, operator must press the 'OK' key.

The following screen categories may appear:

- **Consumable time:** The system is set to remind the operator of the hand piece running time. A reminder message will appear on the screen alerting the user of the remaining time on system and will refer the user to the system information screen for more details.
- **Cooling:** system will alert once electrodes should be cooled down during the treatment. A progress bar will indicate when cooling process is completed.
- **FR tip warnings:**
 - Treatment tips should only be used for a single treatment application on a single patient. Tips should be disposed of immediately following each treatment to avoid cross contamination.
 - System will alert the user if an empty or a used tip is inserted.
 - Every 20 pulses, the operator should stop treatment and clean the FR tip of any debris that may affect the treatment. System will go into 'Standby' mode and trigger/foot pedal will be disabled until 'OK' key is pressed.

Warning

- FR tip replacements should only be purchased through Viora's authorized local distributor.

- **Low impedance warning:** This warning will appear when the lower threshold of the impedance is detected. Operator should check the integrity of the tip and its electrodes and change tip if required.
- **High impedance warning (bad coupling):** operator should make sure there is a good contact between the skin and the electrodes.

Figure 19: Warning Screen

Figure 20: Note Screen

6 CLINICAL GUIDE

6.1 INDICATIONS FOR USE

See sections 6.6.1 and 6.7.1

6.2 CONTRAINDICATIONS

The following conditions are contraindicated to the treatment:

- Pregnancy and breast feeding.
- Treatment over tattoo or permanent makeup.
- Any Permanent Implant in the designated treatment area (including but not limited to: silicon or injected chemical material, metal plates or other metal material).
- Any implanted electronic device anywhere in the body (such as cardiac pacemaker, defibrillator and cochlear implant).
- Severe concurrent conditions (such as cardiac disorders).
- Blood coagulopathy or excessive bleeding or bruising or use of blood thinning medications (anticoagulants) whether prescribed or over the counter.
- Active skin disease in the treatment area (such as psoriasis, sores, eczema and any type of rash).
- Tendency to skin disorders (such as keloids and impaired wound healing process) and extremely dry and fragile skin.
- Skin cancer (active or in the past), pre-malignant moles, malignancy (active or recent) or history of any kind of cancer.
- Immunosuppressive diseases (such as HIV-Positive) or the use of immunosuppressive medications.
- Endocrine disorder (such as diabetes and thyroid disease).
- History of disease which may be stimulated by heat, such as recurrent Herpes in the treatment area.
- Plastic surgery such as face lifting or eyelid surgery within the past 12 months.
- Use of Accutane within the past 6 months.
- Augmentation techniques with injected bio-materials (such as injections of Botulinum toxin, collagen protein or fat) within the past 6 months.
- Natural filler (such as hyaluronic acid) in the treatment area within the past 6 months.
- Surgical procedure in the treatment area within the past 6 months (or during the healing process).
- Facial controlled abrasion (such as dermabrasion, facial resurfacing or deep chemical peeling) within the past 3 months.
- Any esthetic treatment with light (such as IPL or laser device), radiofrequency or other devices in the treatment area within the past 3 months.
- Extremely tanned skin (including tanning beds or tanning creams) within the past 2 weeks.

- Any use of NonSteroidal Anti-Inflammatory Drugs (NSAIDs, such as Nurofen, Advil, and Motrin), 7 days before and 7 days after every treatment:

Caution

Though the treatment is not advised for any of the above conditions, patients may be treated in some cases after consulting with their primary physician and providing a written consent for the treatment at their physician's discretion.

6.3 POSSIBLE ADVERSE EFFECTS

Improper use of the system may result in side effects. Although these effects are highly rare and expected to pass, any adverse reaction should be immediately reported to the patient's physician and to Viora's clinical team:

Appearance of side effects may occur during the treatment or shortly after. These side effects may include the following:

- Discomfort
- Appearance of any pigmentation (hyper-pigmentation or hypo-pigmentation)
- Changes in skin texture, such as burn, blister and crust
- Excessive skin edema and /or erythema.
- Urticaria (hives).
- Scarring.

The following adverse effects relate only to FR treatment:

- Post treatment mild or strong pain
- In-grown hair
- Skin infection

Caution

In the event side effects persist or are severe than the above, the patient should consult a physician and notify Viora immediately.

Note

- Some dark-skinned patients (including Asian skin) may have a delayed response of one-to-two days after treatment and should be evaluated post-test procedure accordingly.
- Proper use of the system by following the treatment protocol and operation instructions will prevent or minimize adverse effects.

6.4 REALISTIC EXPECTATIONS

It is extremely important for the treatment's success and patient's satisfaction to set realistic expectations and advise the client of the following:

- Shortly after the treatment erythema and edema are likely to appear in the treated zone. The edema is likely to clear up within a few hours while erythema may persist up to 12-24 hours from treatment time.
- The degree of response to the treatment, and the number of treatment sessions required, will vary among patients and will depend on the clinical and physiological condition at the start of the treatment regimen.
- The outcomes of the treatment will appear progressively, while temporary improvement in the skin will be observed immediate after each treatment; the maximum long-term result can be observed 2-4 months after the last treatment.
- Maintenance sessions may be required.

Additionally, the following should be explained to the patient prior to FR treatment:

- Light ablation may be detected with gradual improvement in skin texture/flexity.
- Following treatment, and during healing, small scab-like dots may form in the treatment area, especially at the point of contact between the micro-electrodes and skin.
- After the treatment and throughout the healing stage, small pin-point crusts may appear in the area of each ablated spot, particularly at the contact area of the electrode pins with the skin. These crusts will undergo natural exfoliation.

6.5 PRE-TREATMENT

6.5.1 CONSULT MEETING

The first visit of any patient should be dedicated to reviewing treatment possibilities, contraindications and performing a test spot in order to establish the patient is an appropriate candidate for the treatment, and test for the optimal treatment parameters.

Upon the patient's initial visit, practitioners should follow the below steps during a consult meeting:

- Establish patient's needs and review symptoms.
- Set patient's expectations and desired end point for the treatment.
- Complete patient's medical history and review contraindications.
- Review with the patient the application method and operation, expected treatment results and possible side/adverse effects during or post treatment.
- Prior to commencing treatment sessions advise the patient to:
 - α Avoid skin irritation or exposure to sun and using sunscreen during daylight hours.
- Stop using any irritant topical agents for at least 2-3 days before to treatment session.
- Sign the patient on an informed consent form including the above.
- Establish a schedule of treatments.
- Take photos, etc.

- Perform a test procedure.

6.5.2 TREATMENT PARAMETERS

- Selecting the appropriate RF energy level is crucial for safe and effective operation.
- When changing parameters always take under consideration:
 - The diffusion and depth of the lesion and its anatomical location (such as proximity to bony area).
 - Patient's skin type, tolerance and tendency to scar.
- Examine the skin response to the chosen parameters before increasing the RF energy level.
- According to the practitioner's discretion and specific patient's needs, the treatment parameters can be modified during subsequent sessions.

Note

- Different temperatures cause different effects in the tissue, from activation of cellular components of the skin and cellular destruction (such as necrosis and coagulation) up to the effect of ablation.

Warning

- Patients with skin types V-VI may be more sensitive to treatment. In the case of dark skin patients, lower the RF energy by 20%.
- Additional tests should be conducted before raising the RF energy for patients with such skin type.

6.6 FR CLINICAL GUIDE

6.6.1 GENERAL

The V-touch system FR Treatment is indicated for dermatological procedures requiring ablation and resurfacing of the skin. This results in heating of isolated spots of skin that leave a healthy surrounding tissue allowing fast healing process and extra-cellular matrix remodeling mainly driven by fibroblast activity.

The V-touch uses bi-polar radio frequency technology to deliver energy in a pixelated fashion on to the skin; in order to treat dermal layers with minimal injury to the epidermis.

The procedure leads to pin point localized thermal zones inducing ablation and coagulation effect in the dermis. This results in heating of isolated spots of skin that leave a healthy surrounding tissue allowing fast healing process and extra-cellular matrix remodeling mainly driven by fibroblast activity.

Warning

Dermatologists and aesthetic physicians may treat patients only after passing appropriate training and adhere to all local and state regulations regarding the use of this system.

6.6.2 PRE-TREATMENT PREPARATION

The following steps should be taken before any treatment or test is performed:

- Hair in the treatment area should be shaved. Hair may disrupt electrodes' full coupling with the skin and as result may cause burns.
- Use a mild cleansing solution to clean the treatment area. Make sure the treatment area is dry and free of any moisture, liquids or lint.
- Apply a topical anesthetic material (water-based only, containing at least 15% active ingredients) on the entire treatment area. Follow manufacturer guidelines when applying material.
- At the end of anesthetic activity (according to manufacturer guidelines) wipe off the excess material from the skin with a tissue.
- Wash the skin with soap and lukewarm water and make sure the skin is completely clean of any residues.
- Once the skin is clean, sterilize the treatment area using an alcohol pad (at least 70% ethanol).
- Make sure the skin is completely dry before initiating the treatment.
- There is no need for the use of any medium solutions during the treatment (such as ultrasound gel or coupling gel).
- Before beginning the treatment make sure the hand piece is clean, dry and intact.

- Open the sealed pack of the treatment tip and check the integrity of the pins making sure they are all intact and that the tip's frame is undamaged.
- Connect the treatment tip to the hand piece and ensure the tip is secured in place.

Note

Areas of metal dental implants (such as dental braces, caps and crowns) may be more sensitive to RF energy. In this case, the practitioner may isolate the area using gauze, dental rolls or a tongue depressor.

6.6.3 TEST PROCEDURE.

Before receiving a patient for a V-touch treatment session, a test procedure must be performed during the initial consult meeting. The test procedure should be performed at non-exposed site of the treatment area. If more than one area is to be treated, a test should be performed on all treatment sites.

The test procedure serves several purposes:

- Establishing the client is an appropriate candidate for the treatment:
 - For skin types I-III wait 1-2 days before assessing the skin response.
 - For skin types IV-VI wait 5-7 days before assessing the skin response.
- Providing the client with an opportunity to experience the V-touch treatment.
- Finding the correct parameters according to patient's tolerance to the test procedure
- Minimizing possibilities for adverse effects.

Warning

- Patients with skin types V-VI may be more sensitive to treatment. Lower energy levels should be applied:
- When treating Asian skin classify it as one grade higher on the Fitzpatrick scale.

The test procedure should be conducted as follows:

- a. Emit a single pulse and wait 10 minutes to allow pattern to appear.
- b. In the event a pattern does not appear, increase the RF energy and release a second pulse in an adjacent location to the previous pulse.
- c. Desired end point for the test and the treatment procedures is the appearance of a grid like pattern matching the electrode pins accompanied by erythema and slight edema.

Note

- Test procedure is to be executed and evaluated prior to beginning treatment schedule.

- For complete treatment protocol please refer to the treatment protocol attached to this user manual.

At the start of the treatment course, and at every subsequent session, it is important to inquire and update the patient history portfolio, being mindful of new prescriptions, aesthetic or medical procedures, etc.

6.6.4 GENERAL TREATMENT GUIDELINES

- After completing the pre-treatment guidelines, set the treatment parameters according to the test results.
 - Be careful when treating close to bony areas; use lower RF energy levels.
- Place the tip in full contact with the skin at a 90 degree angle to the surface of the target region.
- Apply one pulse by pressing down the trigger on the hand piece or the foot pedal at the same time.

Note

Pressing the trigger/ foot pedal continuously will change the system into 'Standby' mode.

Only one pulse at a time is possible when treating with the RF hand piece.

- The LED located at the back of the hand piece will light up as an indication of the released pulse. The light has no bearing on the actual treatment and is only a visual indication that a pulse has been released.
- On Standby mode, the LED light will be OFF.
- On Ready mode, the LED light will blink.
- Once the pulse is completed a beeping sound is heard to indicate the pulse is completed. The hand piece can then be safely moved to an adjacent spot in the treatment area.
- For every pulse released, the number of pulses indicated on the display screen will increase by one.
- Do not leave untreated spaces between the pulses and do not overlap between the pulses.
- Pass over the treatment area once.
- In case the system recognizes insufficient coupling, the RF energy will be disabled and the system will not emit a pulse. This will be accompanied by a message on screen and an audio indication. The system will not count this as a pulse.
- Throughout the treatment, make sure that the electrode pins are:
 - Clean and have no cellular remains (appearing as black material).
 - Make sure that the tips are cleaned every 20 pulses while the system is on Standby mode. A note will appear on the screen every 20 pulses.

- Unharmful – ensure that no pins are missing or damaged.
- Examine the treated area for the desired end point; erythema and edema on the treatment area.
- Edema will normally appear within several minutes after the treatment.
- At the end of the treatment, examine the area to make sure there are no untreated patches of skin. If such patches exist follow the below instructions:
 - Set the system to *Standby Mode*.
 - Clean the untreated patches and the tip with a dry lint-free gauze pad.
 - Pass over the untreated patches and avoid treated areas.
- At the end of every treatment session discard the treatment tip and clean the hand piece according to the cleaning guidelines.

Warning

- Do not treat over the eyelids or the lips.
 - Do not treat different patients with the same tip, in order to prevent cross-contamination.
 - Note that in darker skin types erythema is not easily noticeable; therefore, the end point of edema will be more outstanding in these patients.
 - Before cleaning the tip make sure to set the system to *Standby mode*
-

6.6.5 TREATMENT SCHEDULE

Note that treatment sessions may be considerably longer and depend on the treatment area, physiological attributes and patient's response.

- Number of Treatments: 3-4
- Treatment Intervals: 4-6 weeks
- Maintenance: 1 treatment every 6-12 months, or as required.

Warning

It is not advisable to perform more than five consecutive sessions.

6.6.6 TREATMENT FOLLOW-UP

Post-treatment follow-ups are very important to determine the most suitable regimen for each patient.

- The patient must be evaluated 2-3 days post treatment to make sure no side/adverse effects appear.
- Always inquire and update the patient history portfolio (including medical history and contraindications).

Operational documentation of the manufacturer:

- Following completion of treatment course, patient should return after 4-6 weeks for evaluation.
- Prior to any maintenance treatment, perform a new test procedure in order to choose and verify the treatment parameters.

6.6.7

- Emollient cream should be applied to the treatment area.
- If patient notes excessive heat sensation or strong erythema or any discomfort, the area may be cooled with the following: chilled Aloe vera gel, cold (not frozen) packs or cold compresses.
- If rare blistering or ulcerated skin occurs, mild hydrocortisone ointment may be applied or a prescribed antibiotic ointment.
- Usually, 1-3 days after the treatment small scabs (less than 1 mm in diameter) will appear and will possibly remain for several days. Do not scratch or touch the scabs and allow them to shed off naturally from the skin.
- Advise the patient of the following precautions during the first 2 days after the V-touch treatment:
 - Avoid using hot water massages.
 - Avoid using of abrasive or harsh products.
 - Avoid any mechanical or thermal damage to the treated area.
 - Avoid any type of skin tanning (by sun exposure, tanning beds or artificial sunless tanning lotions) during the entire course of the treatment.
 - Keep the skin clean in order to avoid contamination or infection.
 - Use a high-factor sunscreen (at least 30 SPF) to protect the treated area for at least 1 month after every treatment.
- Moisturizing cream may be applied several hours (~½ a day) after treatment.
- Make-up may be applied 24 hours after treatment.
- Avoid using any kind of soap on the treatment area 24 hours after the treatment.

6.7 ST CLINICAL GUIDE

6.7.1 GENERAL

The V-touch system ST treatment is indicated for dermal heating for the purpose of raising skin temperature for the treatment of selected medical conditions such as relief of pain, muscle aches and spasm and for increase in local blood circulation.

6.7.2. TEST PROCEDURE

The first clinical visit of any patient should be dedicated to reviewing treatment possibilities, contraindications and performing a test procedure in order to establish the patient is an appropriate candidate for the treatment, as well as testing for the optimal treatment parameters. In addition, performing a test procedure minimizes the possibilities for side effects and adverse response to the treatment. At the start of the treatment course and at every subsequent session, it is important to inquire and update the patient history portfolio, being mindful of new prescriptions, aesthetic or medical procedure, etc.

6.7.3. GENERAL TREATMENT GUIDELINES

- Set the RF mode according to the treatment area and as indicated in the treatment protocol.
- Start with the minimal RF intensity indicated on the energy bar.
- Place the handpiece on the treatment area in a perpendicular orientation and apply slight pressure, making sure both electrodes are in good contact with the skin surface.
- Apply a layer of ultrasound gel in an amount covering half the height of the electrodes.

6.7.4. TREATMENT SCHEDULE

Note that treatment sessions may be considerably longer and depend upon the treatment area, physiological attributes and patient's response.

- Number of Treatments: 3-6
- Treatment Intervals: 2-3 weeks
- Maintenance: 1 treatment every 3-6 months, or as required

6.7.5. POST TREATMENT CARE

- If patient notes excessive heat sensation or strong erythema, the area may be cooled with the following: chilled aloe vera gel, cold (not frozen) packs or cold compresses.
- In any other case of rare blistering, refer the patient to a doctor.
- If adverse effects appear, promptly stop the treatment and apply cooling gel. Immediately refer to a doctor.
- Treatment should be resumed only after the patient has been thoroughly evaluated by a doctor and the side effects have subsided. Reduce the treatment parameters in subsequent treatments.
- Following the skin tightening treatment, patients should take special care of the treated area, making sure to avoid the use of hot water and discontinue use of abrasive or harsh products for two days.

7 DATA FOR THE TECHNICAL MAINTENANCE OF MEDICAL DEVICE

The V-touch system requires standard maintenance and cleaning procedures to ensure optimal system performance. However, please note that the system should be technically serviced only by authorized service technicians.

This section details the general routine maintenance procedures that may be performed by trained and approved operators.

Warning

DISCONNECT THE POWER SUPPLY CABLE BEFORE PERFORMING ANY MAINTENANCE OR SERVICE PROCEDURES TO AVOID ELECTRIC SHOCK.

- For continued protection against fire, replace the fuse only with one of the same type and rating.

7.1 CLEANING THE SYSTEM

- It is recommended to clean the system at least once a week.
- Turn the system OFF.
- Wipe all external surfaces with a soft, damp, non-abrasive cloth.
- To clean the LCD screen, use a suitable solution for computer screens.

7.2 CLEANING THE ST ELECTRODES

ST electrodes must be thoroughly cleaned after every treatment.

Gel accumulation on ST electrodes may cause a spark and potentially a superficial burn in the next treatment performed.

- Turn the system OFF.
- Wipe the applicator tip with a dry cloth to remove excess ultrasound gel.
- Using a cotton pad and an isopropyl-based alcohol solvent (70%), carefully wipe the applicator electrodes, making sure all corners are thoroughly cleaned.

7.3 CLEANING THE FR TIP

Throughout the treatment, the FR tip should be cleaned every 20 pulses.

- Make sure system is on 'standby' mode before cleaning the tip.
- Use a clean gauze pad soaked with alcohol (70%) to gently clean the tip.
- After cleaning, use a dry gauze pad to wipe off the electrodes. Make sure the tip is completely dry before resuming treatment.

8 DATA FOR THE SERVICE OF A MEDICAL

Manufacturer's Authorized Representative in Russia

- AlfaSpa Med. EtD, 109240 Moscovretskay nab, 2a, building 1, Moscow, Russia.
e-mail: alfaspa@alfaspa.ru, call: www.alfaspa.ru

The V-touch system is equipped with a self-testing software routine that continuously monitors system operation. The following troubleshooting guide does not attempt to list all possible system failures:

Table 1 offers a list of the most probable system errors and instructions on how to solve them. If the corrective actions listed in the table do not solve the problem, please contact your local distributor for further assistance.

Message on screen	Probable Cause	Corrective Action
Bad coupling	Bad coupling between the electrodes and the skin.	Make sure that all electrodes are in full contact with the skin.
Overheating	System is over heated	Wait for few minutes until system is cooled down and the message disappears from the screen.
No HP Connected	Improper connection of the HP	1. Disconnect and reconnect the HP 2. Contact the local distributor
No Tip Connected	When FR Treatment is selected from the main menu and there is no treatment tip connected	1. Connect the FR Tip 2. If there is already a tip connected, please replace it with a new one 3. Contact the local distributor
Foot pedal failure. Please make sure foot pedal is not pressed	Foot pedal is pressed	Release foot pedal
HP failure Please make sure trigger is released	Trigger is pressed	Release Trigger

Table 1.

- **Technical Warning screens** - Will appear under various, faulty conditions. Upon receiving one of the following warning screens, please contact your local distributor for more information and repair:
 - System failure
 - ST failure
 - FR failure
 - ST HP cooling failure
 - FR HP cooling failure
 - Communication failure

9 SYSTEM SPECIFICATIONS

	ST	FR
RF Energy	Up to 130J/cm ²	Up to 10J
RF Frequency	Frequency Mode 1 - 0.8MHz Frequency Mode 2 - 1.7MHz Frequency Mode 3 - 2.45MHz Frequency Mode 4 - 0.8MHz, 1.7MHz, 2.45MHz	1 MHz
Treatment area	8mmX8mm	8mmx7mm
Pulse duration	50-200 ms	10-100 ms
Electrical requirements	Maximum 250 watts	
Weight	22 Kg	
Dimension	100H x 54W x 46D Cm	
Electrical Requirements	100-240V, 50/60 Hz	

About Viora

Viora is a market leader in the development of technologically advanced medical aesthetic systems for physicians and aesthetic clinics worldwide. Dedicated to the research and development of safe and effective aesthetic solutions, Viora provides the highest standard of professional grade solutions while keeping your business in mind. Viora offers a wide range of applications including body contouring, cellulite reduction, skin tightening procedures as well as skin resurfacing, phototherapy and anti-aging treatments.

All of the products are designed by a team of veteran medical device developers in close cooperation with the company's advisory board including medical aesthetics practitioners, plastic surgeons, and dermatologists. Viora's aesthetic devices are clinically tested and endorsed by physicians and aesthetic specialists worldwide.

APPENDIX A – EMC Tables

Guidance and manufacturer's declaration – electromagnetic emissions		
The V-touch is intended for use in the electromagnetic environment specified below. The customer or the user of the V-touch should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The V-touch uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The V-touch is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	
RF emissions CISPR 15	Complies	The V-touch is not suitable for interconnection with other equipment.

Operational documentation of the manufacturer

Guidance and manufacturer's declaration – electromagnetic immunity			
The V-touch is intended for use in the electromagnetic environment specified below. The customer or the user of the V-touch should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines N/A	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV differential mode ± 2 kV common mode	± 1 kV differential mode ± 2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	< 5 % U_n (> 95 % dip in U_n) for 0,5 cycle 40 % U_n (60 % dip in U_n) for 5 cycles 70 % U_n (30 % dip in U_n) for 25 cycles < 5 % U_n (> 95 % dip in U_n) for 5 sec	> 95 % dip for 10 ms 60 % dip for 100 ms 30 % dip for 500 ms 95 % dip for 5000 ms	Mains power quality should be that of a typical commercial or hospital environment. If the user of the V-touch requires continued operation during power mains interruptions, it is recommended that the V-touch be powered from an uninterruptible power

			supply.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE - U_T is the AC mains voltage prior to application of the test level.			

Guidance and manufacturer's declaration – electromagnetic immunity			
The V-touch is intended for use in the electromagnetic environment specified below. The customer or the user of the V-touch should assure that it is used in such an environment.			
IMMUNITY test	IEC 60601 TEST LEVEL	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	[3] V	Portable and mobile RF communications equipment should be used no closer to any part of the V-touch, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = \left[\frac{3,5}{V_1} \right] \sqrt{P} = \left[\frac{3,5}{3} \right] \sqrt{200}$ $= 16.5 [m]$ $d = \left[\frac{3,5}{E_1} \right] \sqrt{P} = \left[\frac{3,5}{3} \right] \sqrt{200}$ $= 16.5 [m]$ 80 MHz to 800 MHz
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,5 GHz	[3] V/m	

			$d = \left[\frac{2}{E_1} \right] \sqrt{P} = \left[\frac{2}{3} \right] \sqrt{260}$ $= 33 \text{ [m]}$ 80 MHz to 2,5 GHz 800 MHz to 2,5 GHz, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b Interference may occur in the vicinity of equipment marked with the following symbol: 
NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the V-touch is used exceeds the applicable RF compliance level above, the [ME EQUIPMENT or ME SYSTEM] should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the V-touch.			
^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than [V1] V/m.			

Recommended separation distances between portable and mobile RF communications equipment and the [ME EQUIPMENT or ME SYSTEM]			
The V-touch is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the V-touch can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the V-touch as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter [W]	Separation distance according to frequency of transmitter [m]		
	150 kHz to 80 MHz $d = \left[\frac{3,5}{V_1} \right] \sqrt{P}$	80 MHz to 800 MHz $d = \left[\frac{3,5}{E_1} \right] \sqrt{P}$	800 MHz to 2,5 GHz $d = \left[\frac{7}{E_1} \right] \sqrt{P}$
0,01	0.12	0.12	0.23
0,1	0.37	0.37	0.74
1	1.17	1.17	2.33
10	3.69	3.69	7.38
100	11.67	11.67	23.33
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.			
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

10 SYSTEM COMPOSITION:

I. COMPOSITION

1. Device V-touch - 1pcs.
2. Nozzle FR for fractional procedure 1-10 pcs.
3. Nozzle FR for fractional procedure - 1pcs.
4. Network cable - 1pcs.
5. Bracket - 1-10pcs.
6. Operation manual - 1pcs.

II. ACCESSORIES

1. Nozzle ST for lifting-1-10 pcs.
2. Nozzle FR for fractional procedure 1-10 pcs.
3. Disposable sterile pipette nozzle FR - 1-10 pcs.
4. Service Kit (Line Filter-1-10pcs, Power Entry 2 pol-1-10pcs, FUSE -FST 5x20 4A 1-10pcs, POWER SUPPLY-1-10pcs, Fan 80x25-1-10pcs, RF Board Assembly-1-10pcs, UI Board ASSY-1-10pcs, PCB-CPU -1-10pcs, Touch panel filter ASSY -1-10pcs, Screen Touch LCD -1-10pcs, Wire connection- 1-10pcs,).
5. Bracket - 1-10pcs.

11 DATA FOR RECYCLING

	Dispose of according to WEEE regulations
---	--

Numbered, bound and sealed
Viora Ltd

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Copyright © 2012 Viora. Все права защищены.