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CERTIFICATE

中国国际贸易促进委员会暨中国国际商会
China Council for the Promotion of International Trade is China Chamber of International Commerce

00498434

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China Council for the Promotion of International Trade
China Chamber of International Commerce

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THIS IS TO CERTIFY THAT: the seal of SYSMED (CHINA) CO., LTD. on the annexed DOCUMENT is genuine.

China Council for the Promotion
of International Trade

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
证明书查询地址 Website for verifying the certificate: <http://www.rzccpit.com/validate.html>

SCALEO

M E D I C A L

Stationary oxygen concentrator model HO₂RIZON S5.

USER'S MANUAL

CE
0123

USER MANUAL Initial release.

Name of medical device

Stationary oxygen concentrator model HO₂RIZON S5, composed of:

1. Main unit assembly - 1 pc.
2. Humidifier - from 1 to 5 pcs. (if necessary)
- 2.1 Humidifier connecting tube - from 1 to 10 pcs. (if necessary)
3. Compressor filter - from 1 to 10 pcs. (if necessary)
4. Porous sponge filter - from 1 to 10 pcs. (if necessary)
5. Nasal cannula - from 1 to 10 pcs. (if necessary) design options:
 - 5.1 Nasal oxygen cannula MEDEREN, manufactured by Mederen Neotech Limited, Israel (RU No. RZN 2018/7623 dated 17.09.2018).
 - 5.2 Nasal oxygen cannula, manufactured by "Plasti-Med", Turkey, (RU No. FSZ 2012/12924 dated 01.02.2021)

6. Fire damper - from 1 to 10 pcs. (if necessary)
- 6.1 Fire damper tube - from 1 to 10 pcs. (if necessary)
7. User manual - 1 pc.

Intended use: Stationary oxygen concentrator, model HO₂RIZON S5 is designed to administrate the direct oxygen therapy (oxygen-air mixture of $93 \pm 3\%$) or joint work with anesthesia and respiratory equipment to provide patients with oxygen, with the exception of anesthesia and respiratory equipment operating in a "closed" (with CO₂ absorber) circuit. Electric mains (alternating current) device for enrichment with oxygen-air mixture $93 \pm 3\%$, supplied to the patient through a built-in flow meter and a nasal cannula. This device is not intended to support vital functions and does not allow monitoring the condition of patients. The equipment is intended for use at home, in various medical centers and institutions.

Stationary oxygen concentrator, model HO₂RIZON S5 is suitable for adults and children.

Scope of application: The device is intended for installation in a medical facility or at home.

Information about potential consumers of a medical device.

The products can be used for professional use in medical and treatment-and-prophylactic medical organizations and institutions, in a hospital, in an outpatient clinic, for independent use at home as prescribed by a doctor.

Medical device code in accordance with the OKPD2 Classifier.

OKPD2 32.50.21.129

Operating conditions:

- Room temperature during operation - from 10 to 37 °C
- Indoor relative humidity during operation - 20 % - 65 %
- Atmospheric pressure during operation from 80 kPa to 101 kPa
- Operation room - Dry and ventilated, no dust or any corrosive gases, no strong electromagnetic interference
- Recommended altitude for use - below 3000 meters

Type of medical device in accordance with the nomenclature classification of medical devices approved by the Ministry of Health of the Russian Federation - 191160.

Class depending on the potential risk of use – 2a according to GOST 31508, the order of the Ministry of Health of the Russian Federation dated 06.06.2012 № 4n.

Indications for use:

Stationary oxygen concentrator, model HO₂RIZON S5 is designed to administrate the direct oxygen therapy (oxygen-air mixture of $93 \pm 3\%$) or joint work with anesthesia and respiratory equipment to provide patients with oxygen, with the exception of anesthesia and respiratory equipment operating in a "closed" (with CO₂ absorber) circuit. This device is not intended to support vital functions and does not allow monitoring the condition of patients. The equipment is intended for use at home, in various medical centers and institutions.

Stationary oxygen concentrator, model HO₂RIZON S5 is suitable for adults and children.

Contraindications

Stationary oxygen concentrator, model HO₂RIZON S5, is contraindicated in patients with severe carbon monoxide poisoning.

Possible side effects, warnings

Excessively high concentrations of the oxygen-air mixture can cause poisoning, which leads to lung damage or the development of respiratory failure in patients with a predisposition to these disorders. Higher air/oxygen concentrations also increase the risk of fire, especially when smoking, and without moisturizing can also cause nasal dryness. The recommended oxygen-air saturation range depends on the condition of the patient being treated.

An oxygen-air mixture of $93 \pm 3\%$ should not be administered to patients with paraquat poisoning, as the administration of an oxygen-air mixture of $93 \pm 3\%$ may increase the degree of intoxication.

Software:

Software version: A / 1 and above

A - major version number

1 - minor version number, increases as changes are made Software class: B (Possible minor injuries).

Environmental protection requirements, including requirements for safe destruction and disposal. Environment protection

The devices must be non-toxic and meet the requirements of ISO 10993-1. During the production, storage, transportation and operation of devices, toxic substances are not released into the environment and do not have (in direct contact) harmful effects on the human body. Working with the devices requires the special precautions described above. Medical waste hazard criterion - class B.

Disposal of devices must be in accordance with 2002/96/CEE Disposal of Electrical and Electronic Equipment. In order to preserve the environment, disposal of electronic equipment may be regulated.

Dispose devices and their components at the end of their useful life in accordance with local regulations. For more information on disposal or recycling, contact your local authority.

When the device has reached its useful life, it must be disposed of. The disposal procedure includes the following steps: first, you must disinfect the device completely, and then send it to a special return department for allocation according to the type of material.

Manufacturer: SysMed (China) Co., LTD.

Address: No.17, Wensu Street, Hunnan New District, Shenyang, China, Китай

Tel. 024-24697305 info@sysmed.cn

Authorized representative of the manufacturer in the Russian Federation: Limited Liability Company "Scaleo Medical" (Scaleo Medical LLC)

Legal address: 117513, Moscow, Leninsky prospect, house 139, underground floor No. 1, premises VI, room 3. info@scaleomedical.com contact@scaleomedical.com

Mailing address: 115172, Moscow, Narodnaya str., 9, PO Box 1.

Contents

Safety instructions	5
Safety instructions	5
Warnings	6
Symbols and legend	9
General information about Stationary oxygen concentrator, model HO₂RIZON S5....	11
Stationary oxygen concentrator, model HO ₂ RIZON S5	11
Accessories to stationary oxygen concentrator, model HO ₂ RIZON S5	12
Accessories and spare parts	13
Unpacking the stationary oxygen concentrator, model HO ₂ RIZON S5	13
Instructions for use	13
Warnings and troubleshooting.....	16
Warning signals	16
Fault handling	17
Cleaning and repair	18
Specifications	20
Classification	22
Operating conditions	23
Storage and transportation conditions	23
Basic electrical diagram	24
Electromagnetic compatibility	24
Guidance and Manufacturer's Declaration - Electromagnetic Compatibility	24
Guidance and Manufacturer's Declaration - Electromagnetic Immunity	25
Warranty and liability of the parties	28
Medical device labeling and packaging	28
Complete set of a medical device.....	29
Service life and disposal of the device.....	29
Service life.....	29
Device disposal	29
Current regulations	30
Feedback	31

Safety instructions

Safety instructions

Before using the device, carefully read the user manual, save it for future use. Do not use this device if you have not read or understood the instructions completely. This can lead to injury and damage.

If the information in this user manual is difficult to understand, please contact your dealer or an authorized representative of the manufacturer Scaleo Medical LLC.

If oxygen therapy is carried out at home, the doctor's instructions must be strictly followed. Oxygen therapy should be carried out in accordance with the observance of the daily dose of the oxygen-air mixture and the number of inhalations during sleep, wakefulness and rest. In order to avoid an insufficient supply of the oxygen-air mixture or oxygen poisoning, it is strictly forbidden to change the flow of the oxygen-air mixture and change the time of the procedure.

Patients who urgently need an oxygen-air mixture or when it comes to seriously ill patients and patients who need constant oxygen inhalation, it is obligatory to have other oxygen-air mixture delivery devices when using a stationary oxygen concentrator, model HO₂RIZON S5 - cylinders with compressed oxygen or liquid oxygen systems.

If during oxygen inhalation the patient has an adverse reaction and discomfort, it is worth stopping the use of the device immediately and contacting an authorized representative of the manufacturer Scaleo Medical LLC or a doctor.

In the event of a device malfunction, immediately contact the dealer or an authorized representative of the manufacturer Scaleo Medical LLC or a certified after-sales service center, under no circumstances disassemble or repair the device by yourself.

Additional supervision and attention when using this device may be required by children, elderly patients, or patients with communication, hearing or vision impairments who are unable to respond to alarms or report discomfort. This device is not intended for infants.

Patients with severe carbon monoxide poisoning cannot use this machine.

The nasal cannula is intended for individual use and must be changed periodically. Use and clean the nasal cannula as recommended by the cannula manufacturer.

The flowmeter is a plastic flask with a ball installed inside. On the surface of the flask, a graduation is applied to determine the O₂ flow rate (from 0.5 to 5 l/min with an interval of 0.5). The position of the ball opposite a certain mark shows what kind of O₂ flow rate is set at the outlet of the oxygen concentrator. The required flow rate is set by turning the flow regulator regulator. Be careful

turning the regulator as much as possible, as the flow rate can go beyond the limit red line and up to 6.5 l/min.

Prolonged use can reduce the efficiency of the stationary oxygen concentrator, model HO₂RIZON S5, due to the fact that the molecular sieve loses its effectiveness over time if the operating conditions specified in this manual are not followed (at low and high temperatures, high humidity).

Molecular sieve is warranted for thirty-six months (36) or twenty-six (26,000) business hours after which it must be contacted by a certified after-sales service provider approved by the manufacturer.

Warnings

Before using the oxygen concentrator stationary model HO₂RIZON S5, carefully read the following information.

Particular attention should be paid to reducing the risk of fire during oxygen therapy. When the concentration of the oxygen-air mixture is high, any material becomes extremely flammable and burns quickly, whether it is flammable or not. For safety reasons, flammable objects and liquids should be kept away from a stationary oxygen concentrator, model HO₂RIZON S5, preferably in another room.

Oil, grease and other oil-based products can provoke instant ignition when in contact with an oxygen-air mixture, therefore, all hazardous materials should be kept away from a stationary oxygen concentrator, model HO₂RIZON S5, pipes, couplings, and other components of a stationary oxygen concentrator, model HO₂RIZON S5. Do not use lubricants.

It is strictly forbidden to turn over or install in a horizontal position the stationary oxygen concentrator, model HO₂RIZON S5 during use and movement. Do not use the concentrator in a smoky environment. The optimum temperature in the room or place of use should be between 10 and 37°C. When the temperature is below 10°C, the compressor may have difficulty starting, and above 37°C, the compressor overheats and reduces the life of the device.

It is strictly prohibited to use a stationary oxygen concentrator, model HO₂RIZON S5 over 3000 meters. This will lead to a decrease in the quality of patient care.

Please do not move the stationary oxygen concentrator, model HO₂RIZON S5 during use, and do not turn it on when the flow meter is closed. The oxygen concentrator makes a regular sound when the oxygen-air mixture is produced at 93 ± 3%, which is absolutely normal, and it is also normal that hot air comes out through the hole located at the bottom of the stationary oxygen concentrator, model HO₂RIZON S5. Please do not cover this hole during use.

For more efficient operation of the stationary oxygen concentrator model HO₂RIZON S5, the manufacturer recommends using it for at least 30 minutes each time, since frequent and short-term use can significantly reduce the service life of the device.

If the oxygen compressor, power cord, or plug is damaged, the stationary oxygen concentrator, model HO₂RIZON S5 does not function properly, or is covered with water, contact a Certified Technical Service immediately for inspection. Keep the power cord away from heat and high temperatures. Do not use extension cords.

Do not move the stationary oxygen concentrator, model HO₂RIZON S5, by pulling the power cord. Do not drop or insert foreign objects into the openings. Never cover the openings through which air enters and exits, do not place the device on soft surfaces such as a sofa or bed, as they can block the air inlet. Keep the air inlet of the concentrator away from soft toys, hair, and other similar items.

It is forbidden to use the device in conjunction with anesthesia and respiratory equipment operating in a "closed" (with CO₂ absorber) circuit.

It is necessary to add water to the humidifier in accordance with the calibration volume provided for the humidifier. **Don't add too much water.**

Stationary oxygen concentrator, model HO₂RIZON S5 should be located at least 40 cm away from walls, carpets, furniture and similar items. It is recommended to place the oxygen concentrator stationary in such a way that there is no difficulty in disconnecting it from the power supply and access to the mains switch.

Stationary medical oxygen concentrator, model HO₂RIZON S5 is equipped with a set of oxygen nasal cannula. Patients can also choose another officially approved oxygen cannula if needed. However, the length of the nasal cannula should not exceed 4 meters. The oxygen nasal cannula must be correctly positioned to avoid the risk of pinching the head or neck.

Do not use spare parts and accessories that are not approved by the manufacturer.

The switch is designed to automatically turn off the system in case of voltage surges to protect the electrical circuit. The shutdown sensor is automatically triggered if the current in the device exceeds 5A. If overvoltage or short circuit occurs, you have to press the switch to reset and restart the concentrator. This switch is reusable. If this switch has been activated, there is no need to replace it.

Do not connect the concentrator in parallel or with other oxygen concentrators and oxygen therapy devices.

Use only the original power cord. Unsuitable and non-recommended power cords can increase electromagnetic radiation or reduce electromagnetic immunity, which can lead to malfunction and fault.

The use of humidifiers not designed for use with this fixed oxygen concentrator, Model HO₂RIZON S5, may result in poor performance.

The use of elements of the supply equipment not intended for use with this stationary oxygen concentrator, model HO₂RIZON S5, may lead to a deterioration in its performance.

The manufacturer recommends using a stationary oxygen concentrator, model HO₂RIZON S5, strictly according to the doctor's prescription. Unauthorized flow adjustments are strictly prohibited, as this can lead to oxygen-air poisoning or carbon dioxide retention. This device is designed to compensate for a lack of oxygen concentration, but not to provide first aid or to maintain vital functions.

Avoid sparks in the vicinity of medical devices producing 93 ± 3% oxygen/air mixture, as well as static sparks caused by friction.

If the stationary oxygen concentrator model HO₂RIZON S5 is used for a long time in abnormal conditions (humidity, high/low temperature, increased dust levels), this significantly reduces its efficiency. In this case, you need to contact the seller or the authorized representative of the manufacturer Scaleo Medical LLC or a certified service center.

Regularly review and evaluate your treatment, contact your doctor if you experience any pleasant symptoms while using a stationary oxygen concentrator, model HO₂RIZON S5.

Maintenance

Stationary oxygen concentrator, model HO₂RIZON S5, is designed to minimize the number of technical inspections. Maintenance should be done once a year. Only certified after-sales service providers approved by the manufacturer can disassemble, repair and assemble the device.

Preventive inspection

Before using a stationary oxygen concentrator, model HO₂RIZON S5, it is necessary to check it for serviceability, namely the integrity of the power cord, the product itself and its main parts.

RF interference

 Most devices are susceptible to radio frequency interference. **Therefore, the use of portable communication equipment near a stationary oxygen concentrator, model HO₂RIZON S5 will generate interference.**

This product must not be used in an environment where electric shock, electrosurgery, defibrillation, X-rays, gamma rays, infrared radiation are performed, where there are transient electromagnetic fields, including magnetic resonance (MRI) and radio interference.

Portable RF communications equipment, including peripherals such as antenna cables and external antennas, must be kept more than 40 cm away from the device, or the performance of the equipment will decrease.

To reduce the risk of burns, injury, electric shock, or fire, please read the following information carefully:

- Do not use a stationary oxygen concentrator, model HO₂RIZON S5, while bathing. If the patient needs its constant use, in this case, the stationary oxygen concentrator model HO₂RIZON S5 should be no closer than 3 meters from the bathroom.
- Do not touch the stationary oxygen concentrator, model HO₂RIZON S5 with wet hands or body.
- Do not use or store Stationary Oxygen Concentrator model HO₂RIZON S5 near liquids or other electrically conductive materials.
- Do not touch the stationary oxygen concentrator, model HO₂RIZON S5, if it falls into water or any other easily conductive liquid, in case of a fall, immediately disconnect the power cable from the mains.
- When the stationary oxygen concentrator, model HO₂RIZON S5 is not used, it must be disconnected from the mains.
- If the patient lacks an oxygen-air mixture of $93 \pm 3\%$, it is necessary to immediately contact the supplier of the equipment or a doctor. Do not regulate the flow of the supplied oxygen-air mixture $93 \pm 3\%$ by yourself, unless it is prescribed by your doctor.

- Do not place the stationary oxygen concentrator, model HO2RIZON S5, where it might be where other people might fall over the power cord. Small parts or particles trapped with the air in the stationary oxygen concentrator, model HO2RIZON S5, may cause suffocation if swallowed or inhaled. Store the stationary oxygen concentrator, model HO2RIZON S5, out of the reach of children.
- Do not perform maintenance while the device is in operation.

Warning:

- The oxygen-air mixture is represented by combustion gas. Oxygen compressors should be kept away from sources of ignition and flammable materials. It is strictly forbidden to smoke, as well as to use fire near persons using an oxygen concentrator.
- It is strictly forbidden to smoke during oxygen treatment procedures. If oxygen therapy is prescribed at home, the patient should quit smoking.
- Keep the power cable away from sources of heat or fire.
- The power switch must be in the off position when oxygen therapy is not being administered.
- Do not place the nasal cannula under covers or upholstery.
- Before cleaning the housing of the oxygen concentrator, the device must be disconnected from the mains in order to avoid electric shock.

Symbols and legend

Symbols used in this manual:

This symbol is used to draw attention to important safety information during work, as well as what injuries and damage can appear if ignored. ***Be sure to follow the directions provided and always be careful and attentive.***

This symbol indicates important information regarding the use of the device, which, if not followed, could result in damage or functional malfunction.

This symbol indicates important and useful information. ***If these recommendations are taken into account, it will help the user to significantly simplify the process of work and make it more efficient.***

Tools used on products and labels:

Symbol	Value
I	Turned ON
o	Turned OFF
	Follow the instructions in the user manual
	No smoking
	Attention, refer to accompanying documents
	Class II (Double insulated)
	Working part type BF
	European Directive Compliance Mark
	AC power
Rx Only	Use as directed
	Do not use oil and grease
	Do not disassemble

Symbol	Value
	Up
	Caution, fragile
	Keep away from moisture
	Temperature limitation
	Humidity limitation
	Atmospheric pressure limitation
	Do not expose to sunlight
	Keep away from fire/Risk of fire
IP21	Splash-proof equipment
	See User Manual
	Standby mode
	Danger of electric shock
	Limited number of tiers in a stack

Legend/ symbol	Value
	Manufacturer
	Date of manufacture in YYYY-MM format
	European Authorized Representative
	Product link
	Serial number of the product
	Attention: read the manual before using
	May contain waste electrical or electronic equipment


General information about the stationary oxygen concentrator, model HO₂RIZON S5

In the event that the attending physician has concluded that an additional oxygen-air mixture is necessary to improve the patient's health and has prescribed the use of an oxygen concentrator with a specific prescribed flow setting, **DO NOT CHANGE** these flow parameters by yourself, unless the doctor himself asks you to do so. Please read the information contained in this manual carefully before using this device.

About Stationary oxygen concentrator model HO₂RIZON S5.

Stationary oxygen concentrator, model HO₂RIZON S5 produces a concentrated oxygen-air mixture of $93 \pm 3\%$ from ambient air for patients in need of low-flow oxygen therapy. The oxygen in the air is concentrated using molecular sieves and a pressure change adsorption process.


A medical specialist will definitely tell you how to use a stationary oxygen concentrator, model HO₂RIZON S5 and answer all your questions.

Accessories for Stationary oxygen concentrator, model HO₂RIZON S5


Front view

1. Oxygen-air flow meter
2. Flow rate regulator
3. USB-port
4. LCD display
5. Oxygen-air mixture outlet
6. Power indicator
7. Warning light
8. Power switch

Back view

- 9. Power cable
- 10. Switch
- 11. Frame protecting the air inlet
- 12. Humidifier
- 13. Humidifier retainer
- 14. Compressor filter cover
- 15. Product label

Description of the main parts

- Flowmeter: The level of the floating ball of the flowmeter indicates the amplitude of the flow of the oxygen-air mixture at the outlet, with a scale from 0 to 5 l/min, at 0.5 intervals.
- Oxygen-air mixture outlet: Provides an oxygen-air mixture outlet of $93 \pm 3\%$.
- Alert signal: an information code is displayed on the screen.
- Power switch: «I» ON.; «O» OFF. When the machine starts working, the green indicator light turns on.
- Humidifier: Using a humidifier prevents the dry oxygen-air mixture from entering the throat and nasal mucosa, and it also prevents sputum from drying out.
- Humidifier support retainer: The retainer allows you to secure the humidifier.
- Air inlet frame: filters the air from dust.
- Compressor filter cover: contains the compressor filter.
- The switch is designed for overvoltage and short circuit protection.

Accessories and spare parts

If you have any questions about this device and its use, please contact your doctor or distributor. Use only original accessories and spare parts:

- Compressor filter
- Sponge filter
- Humidifier connection pipe

Unpacking the stationary oxygen concentrator model HO₂RIZON S5

1. Check for visible damage to the packaging or contents. If damage is evident, notify the shipping company and your local distributor.
2. Remove the foam block from the packaging and box.
3. Carefully remove the device and all accessories from the box. Check and inspect the entire exterior of the stationary oxygen concentrator model HO₂RIZON S5 for deformation, dents, scratches, or other damage.

In the event of quality or component problems, please contact your local distributor. *Save the box and packing material for storage and shipping.*

INSTRUCTION FOR USE

USER MANUAL

Do not use extension cords or adapters.

1. Select the location of the concentrator so that it can inhaust the air in the room without any restrictions. Make sure that the device is at least 40 cm away from walls, furniture, and especially away from curtains, that may obstruct the flow of air into the device. Do not place the unit near heat sources.

2. After reading this user manual, plug the power cord into a power outlet.
3. Follow pattern A or pattern B. These are shown below.

A. If you DO not use a humidifier, follow these steps:

A.1. Using a connecting tube, connect one side of the fire damper to the $93 \pm 3\%$ oxygen / air outlet as shown in diagram "1-A1".

Fire damper


Diagram: 1-A1

A.2. Next, to the other side of the fire damper, connect the nasal cannula as shown in diagram "1-A2". The nasal cannula should be connected to the side of the valve where the arrow points are marked.

Nasal headset


Diagram: 1-A2

B. When using a humidifier, follow the steps below:

B.1. Remove the humidifier from the plastic wrap.

B.2. Unscrew the brown cap counterclockwise. Pour the required amount distilled water into the container of the humidifier, then screw the brown cap clockwise.

of

Only distilled, purified, or boiled (and then chilled) water can be used for the humidifier. It is necessary to change the water once a day. Pour water so that the level is halfway between the maximum and minimum mark, or slightly lower, but so that the outlet of the humidifier bubble is closed.

To detach the humidifier container, turn it counterclockwise. When reassembling, align the humidifier container and lid, turn the container clockwise and screw tight to avoid air leakage.


B.3. Install the humidifier with water in the special compartment on the top of the HO₂RIZON S5 as shown in diagram "2-B1".

Diagram: 2-B1

B.4. Use the retainer to secure the humidifier.

B.5. Connect PVC tubing to the outlet of the stationary oxygen concentrator, model HO₂RIZON S5, and to the inlet of the humidifier as shown in diagram "2-B2".

PVC tube

Diagram: 2-B2

B.6. Connect the fire damper to the outlet of the humidifier as shown in diagram "2-B3".


Diagram: 2-B3

B.7. Connect the nasal cannula to the fire damper as shown in diagram "2-B4".



Diagram: 2-B4

B.8. Air duct connection using removable accessories is shown in the diagram below:

**This product must be connected to the fire
 meter according to the above instructions,
 otherwise the manufacturer will not be held liable for
 any possible consequences.**

4. Turn on the power button

The power plug must be correctly and securely plugged into a safe, approved outlet.

5. Close the outlet of the humidifier with your finger for about 20 seconds. Pay attention if the float in the flow meter falls or the humidifier emits a "kyu" sound, if yes, it means that there is no air leakage.

6. Adjust the $93 \pm 3\%$ oxygen/air outlet flow following the line of sight, ball centerline and meter scale on the same horizontal line.

The oxygen-air flow rate of $93 \pm 3\%$ increases when the flow rate control is turned counterclockwise and decreases when the control is turned clockwise. *For preventive maintenance, it is recommended to set the flow rate at the level of 1-2 l / min for a period of 45 to 60 minutes per day.*

For long-term oxygen therapy, it is very important that the patient strictly follows the advice of the attending physician. *Never independently regulate either the flow of oxygen-air mixture $93 \pm 3\%$, or the time of inhalation.*


7. Place a nasal cannula in your nostrils to begin oxygen therapy.
8. At the end of oxygen therapy, press the switch.

i When the stationary oxygen concentrator, model HO₂RIZON S5 is not used, it must be disconnected from the mains.

Warnings and troubleshooting

Warning signals

LCD indications

Zone	Indicators	Possible reasons	Condition	Inspection method	Alarm activation time	Priority
------	------------	------------------	-----------	-------------------	-----------------------	----------

②	000000	Accumulate working hours	Shows the working interface for 10 seconds.	Visual		----
	High	High purity of the oxygen-air mixture	----	Visual		----
	Low	Purity of the oxygen-air mixture is low: 82 %.	Red LED is constantly on and an intense alarm is heard ▲	Bring the flow meter regulator to the maximum level until you hear an alarm.	Less than 1 minute	High
③	H01	Output flow is too weak	Periodically, the yellow LED illuminates and an audible alarm sounds ▲	Bring the flow meter regulator to the minimum level until you hear an alarm.	Less than 1 minute	Low
	H02	Abnormal temperature rise	Red LED is constantly on, intense alarms sound. The device does not work. ▲	Disconnect the fan connection	Less than 2 minutes	High

Indicators	Possible reasons	Condition	Inspection method	Alarm activation time	Priority
① H08	Abnormal pressure observed during operation	The red LED is constantly on, and a constant intense alarm sounds. The device does not work. Δ	Use a variable frequency power supply to set the supply voltage below 85% of the rated voltage	Less than 1 minute	High
③ E01	Temperature sensor malfunction	Periodically, the yellow LED illuminates and an periodic alarm sounds Δ	Disconnect the temperature sensor connector	Less than 1 minute	Low
④ "00:00"	Current working hours		Visual		----
Poweo	Power off signal	A loud beep sounds (≥ 60 dB).	Disconnect the plug from the mains during operation.		----

- In case of numerous faults, information is displayed one after the other. Signals and indicators alternate based on priority.
- All of the above signals are technical.
- 3 minutes after connecting to the network, there is no longer a signal that the concentration of the oxygen-air mixture is below 82%.
- During troubleshooting, the operator must determine the type of failure displayed on the LCD screen and contact the distributor immediately.
- Recommended interval for monitoring the alarm system: 18 month
- The audible alarm has a sound pressure of 60 to 80 dB.

Fault handling

Problem	Possible reasons	Actions for this problem
When connected to a network, the equipment does not work	Insufficient compressor capacity to start or compressor does not work	Call your distributor or healthcare professional

Problem	Possible reasons	Actions for this problem
When the equipment is connected to the network, it does not work or works intermittently	The power cord may be plugged in incorrectly or poorly connected	* Check if the power cord is damaged * Check if the power cord is plugged in correctly * Call your distributor if you still cannot solve the problem
No oxygen-air mixture is generated or the flow of oxygen-air mixture is very low.	* Nasal cannula is kinked or pinched * Humidifier bottle may be pinched	* Detach the nasal cannula * Reinstall the humidifier cover * Call your distributor or doctor if the problem persists
The ball in the flow meter is not controlled by the regulator	* Flow regulator is not tightened * Regulator was turned abruptly or quickly.	* Turn the adjuster smoothly and slowly to tighten. * If you are unable to solve the problem, contact your distributor.
Water rises from the humidifier into the cannula.	* Temperature difference is caused by the weather * Device is too close to a wall, curtains, or furniture * Location of the equipment and cannula is different.	* Dry the inside of the humidifier cover * DO NOT use hot water. * Do not overfill the humidifier * Observe the same temperature regime for both the equipment and the cannula (store everything in the same room).

Cleaning and repair

It is very important to unplug the device before carrying out any cleaning.


Excess moisture can affect the correct operation of the stationary oxygen concentrator, model HO₂RIZON S5

1. Sponge filter:

The sponge filter should be cleaned weekly to ensure uninterrupted operation of the stationary oxygen concentrator, model HO₂RIZON S5, or more often if environmental conditions require it.

To clean the sponge filter, follow these instructions:

- a. Open the frame located on the back of the concentrator. Take out the filter and shake it to get rid of dust.
- b. Rinse the filter thoroughly with clean water, remove excess water and let it dry completely.
- c. Once the filter is dry, reinstall it and close the frame.

Frame through which air enters

Sponge filter


The presence of a second filter allows you to quickly replace it while the other filter dries. Always use original hub accessories. If the filter needs to be replaced, contact your local distributor.

2. Replacing the compressor filter:

The compressor filter is used to filter the dust entering the compressor. It must be replaced regularly. It is recommended to replace the filter every 4000 hours. The compressor filter must be purchased from the distributor or authorized representative of the manufacturer Scaleo Medical LLC.


To replace the compressor filter, please follow the instructions below:

- a) The compressor filter cover is located on the back of the device (see Figure 4-1 below). Open the compressor filter cover (see Figure 4-2 below) and unscrew the filter housing clockwise (see Figure 4-3 below) to remove the compressor filter (see Figure 4-4 below).


Compressor filter cover

Figure 4-1

Filter housing

Figure 4-2

Compressor filter


Figure 4-3

Figure 4-4

b) Replace the new compressor filter, install the filter housing counterclockwise and install the filter cover.

3. Cleaning the nasal cannula:

Use and clean the nasal cannula as recommended by the cannula manufacturer.

4. Cleaning the outer case:

Disconnect the stationary oxygen concentrator, model HO₂RIZON S5 from the mains. Use a damp towel and a small amount of detergent, then dry it off with a dry towel. This procedure must be performed once or twice a month.

5. Disinfection of the outer case:

It is carried out by wiping the outer surfaces of the product with a napkin moistened with a 3% solution of hydrogen peroxide with the addition of 0.5% detergent.

6. Disinfection of the humidifier:

The water in the humidifier must be changed daily. The water used must be distilled or pre-boiled and cooled before use.

The bottle should be cleaned and disinfected at least once every three weeks. To do this, use a mild detergent (such as Dawn TM) and rinse with warm or hot water. Then put the container in the disinfectant and leave for 20-30 minutes, then rinse well in hot or warm water and let dry.

If the humidifier is not used for a long time, it must be stored in a clean plastic bag.

The following can be used as household disinfectants:

- vinegar and distilled water in a 1: 1 ratio,
- bleach (bleach) and distilled water in a ratio of 1:10,
- or commercially available household disinfectants.

It is recommended to change the humidifier once a year.

SPECIFICATIONS


SPECIFICATIONS

Rated power (W)	320
Voltage and frequency	200-240 V (50 Hz ± 1 Hz)

Flow (l/min) *	0,5 - 5
Concentration (nominal flow)	93 % $\pm$ 3 %
Sound pressure level dB	$\leq$ 45
Sound power level dBA	$\leq$ 55
Net weight (kg)	14,96 $\pm$ 10%
Dimensions (mm)	380 x 240 x 660 $\pm$ 10%
Description	Alarm when: • abnormal temperature • low oxygen concentration • power outages • drop or over pressure • weak outlet flow • thermostat malfunction.
Filter characteristics	Particle size: 0.25 - 0.5 μm Compressor filter: 0.25 μm Bacterial filter: 0.3 μm The bacterial filter and compressor filter are classified as HEPA.
Maximum recommended flow	5 liters/min
Pressure range at which the humidifier pressure limiting mechanism is triggered	21,0 kPa $\pm$ 10%
Pressure range at which the concentrator pressure limiting mechanism is triggered	48,26 kPa $\pm$ 10%
The nominal concentration of oxygen in the produced gas, expressed as a percentage by volume of oxygen, at a flow of 2 l/min;	93 $\pm$ 3% (at least 90%; maximum 96%)
Flow values in liters per minute maximum recommended flow at 7 kPa back pressure	Not less than 4.5 l/min
Maximum outlet pressure	48.26 kPa $\pm$ 10% (7 PSI $\pm$ 10%)
Time for the concentrator to reach operating mode	No more than 5 minutes

At maximum rotation, the flow rate regulator can go beyond the limiting red line and reach 6.5 l/min.

Operating the device outside specified voltage, flow, temperature, humidity, or specified altitude range can reduce the concentration of the oxygen-air mixture.

When the nominal outlet pressure of the oxygen-air mixture is zero, the concentration of the oxygen-air mixture is $93\% \pm 3\%$, depending on the operation environment and the nominal flow. See "Oxygen/Air Outlet Concentration Diagram for Various Flow Conditions".

Classification

Electric shock protection class	Class II
Degree of protection of the hub components from electric shock	BF type
Degree of protection of the components of the concentrator from water ingress during use	IP21
Mode of operation	Continuous mode
Degree of safety for the combined use of a flammable anesthetic mixture with air, oxygen or nitrogen oxide	Not suitable for this use
Estimated product life	5 years

 If the mains voltage goes beyond the specified range of 200-240 V, the device may not work correctly, not efficiently or not turn on. High voltage can lead to damage to the device, overheating, additional wear.

If the network is unstable, please install a voltage regulator before using.

Operating conditions

Room temperature during operation	from 10 to 37 °C
Indoor relative humidity during operation	20 % ± 65 %
Atmospheric pressure during operation	from 80 kPa to 101 kPa
Operation Room	Dry and ventilated, no dust or any corrosive gases, no strong electromagnetic interference
Recommended height for use	Lower than 3000 meters

 It is not recommended to use this stationary oxygen concentrator, model HO₂RIZON S5, in an environment higher than 3000 meters above sea level, which may lead to a decrease in the concentration of the oxygen-air mixture.

Recommended temperature for using a stationary oxygen concentrator, model HO₂RIZON S5 from 10 to 37 °C.

At temperatures below 10 °C, the compressor may have difficulty starting. At temperatures above 37 °C, the compressor can overheat, which will significantly reduce its service life.

The ambient temperature of a stationary oxygen concentrator, model HO₂RIZON S5, ranges from 10 °C to 37 °C.


If the device is stored in a room where the temperature is too low or too high, it must be left in the room in which it will be used for 5 hours before use, this will allow it to adapt to the ambient temperature.

Storage and transportation conditions

During transportation and storage, the product must be strictly in vertical position. It is strictly forbidden to turn the product over, transport and store it horizontally.

Storage temperature indoors and during transportation is	from -30 °C to 70 °C
Transportation and storage relative humidity	15 – 95 %, with no condensation
Atmospheric pressure during transportation and storage	500 - 1060 hPa

Basic electrical diagram

Electromagnetic compatibility

To ensure the electromagnetic compatibility of the stationary oxygen concentrator model HO₂RIZON S5, it must be installed and operated in accordance with this user manual.

Portable and mobile radio communication equipment may affect the electromagnetic compatibility of the stationary oxygen concentrator model HO₂RIZON S5. If this happens, contact the manufacturer to find a solution to the problem.

Electromagnetic compatibility refers to the ability of a device to suppress electromagnetic interference without interfering with the operation of other similar devices. Stationary oxygen concentrator model HO₂RIZON S5 electromagnetically generates interference with other equipment.

Ways to solve electromagnetically compatibility problems:

- Use the stationary oxygen concentrator model HO₂RIZON S5 in strict accordance with the instructions provided for in this manual.
- Keep the stationary oxygen concentrator model HO₂RIZON S5 away from other devices to reduce electromagnetic interference.
- The effect of electromagnetic interference can be reduced by adjusting the relative position/angle between the stationary oxygen concentrator, model HO₂RIZON S5 and other devices.
- Electromagnetic interference can be reduced by changing the position of the power cord or signal wire of the stationary oxygen concentrator, model HO₂RIZON S5, and other equipment.
- Reduce electromagnetic interference by changing the connection position of other devices to the network.

Stationary oxygen concentrator, model HO₂RIZON S5 can only be connected using the cord specified in this manual. *The use of cords and accessories that do not comply with the requirements and regulations may increase the electromagnetic interference, as well as a decrease in the electromagnetic resistance of the oxygen concentrator.*

Stationary oxygen concentrator, model HO₂RIZON S5 should not be used in close proximity to other apparatus or equipment. If, nevertheless, you are forced to use a stationary oxygen concentrator, model HO₂RIZON S5 in close proximity to other devices or equipment, in this case, you must monitor the correct operation of all devices in these conditions.

Guidance and Manufacturer's Declaration - Electromagnetic Compatibility

This device is intended for use in the electromagnetic environment specified in Table. The user must ensure that the device is used in such an environment.

Radiation check	Compliance with the requirements	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The apparatus consumes radio frequency energy for its internal function. Therefore, the radio frequency (RF) emissions from this device are very low and may not interfere with nearby electronic equipment.
RF emissions CISPR 11	Class B	The device can be used in any institution, including residential and connected to the low-voltage public network for domestic needs.
Harmonious radiation CEI 61000-3-2	Class A	
Voltage fluctuation / emission fluctuation / CEI	Corresponds	

Guidance and Manufacturer's Declaration - Electromagnetic Immunity

This device is intended for use in an electromagnetic environment (see specifications below). The user of this device should ensure that the device is used in such an environment.

Immunity test	Verification level IEC 60601	Immunity compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2.	± 8 kV contact ± 15 kV air	± 8 kV contact ± 15 kV air	Floors should be wooded, concrete or covered with ceramic tiles. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient / pulse signal IEC 61000-4-4.	± 2 kV for power lines ± 1 kV for input / output lines	± 2 kV for power lines ± 1 kV for input/output lines.	The quality of the power supply must be appropriate for use in a medical facility or at home.
Voltage surge CEI 61000-4-5.	Differential mode ± 1 kV ± 2 kV common mode	Differential mode ± 1 kV ± 2 kV common mode	The quality of the power supply must be appropriate for use in a medical facility or at home.
Voltage dips, short circuits and voltage fluctuations on power lines IEC 61000-4-11.	$< 5\% U_T$ (Drop $U_T > 95\%$) for 0,5 cycle. $40\% U_T$ (drop U_T for 60 %) for 5 cycles. $70\% U_T$ (Drop U_T for 30%) for 25 cycles. $< 5\% U_T$ (Drop $U_T > 95\%$) for 5 sec.	$< 5\% U_T$ (Drop $U_T > 95\%$) for 0,5 cycle. $40\% U_T$ (drop U_T for 60 %) for 5 cycles. $70\% U_T$ (Drop U_T for 30%) for 25 cycles. $< 5\% U_T$ (Drop $U_T > 95\%$) for 5 sec.	The quality of the power supply must be appropriate for use in a medical facility or at home. If the user of the device requires continuous operation, including during power outages, in this case it is recommended to connect the device to an uninterruptible power supply or switch to battery mode.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-4-8	30 A/m.	30 A/m.	The frequency of the power magnetic field must be suitable for use in a medical institution or at home.

U_T - AC mains voltage before applying the test level.

Immunity test	Verification level IEC 60601	Immunity compliance level	Electromagnetic environment - guidance
Radiated interference EN 1000-4-6-6 Conducted RF Interference EN 1000-4-3	3 V (rms.) 150 kHz to 80 MHz 10 V/m 80 MHz to 2.7 GHz.	3 V (rms.) 10 V/m	Portable and mobile RF communications equipment should not be used less than 30 cm from equipment, including cables. Interference may occur near the equipment marked with the following: 

Recommended separation distances between portable and mobile RF communications equipment and stationary oxygen concentrator, model HO₂RIZON S5:

This device is intended for use in an electromagnetic environment in which radiated RF energy and interference are controlled. The user can limit electromagnetic interference by maintaining a minimum distance from portable and mobile RF communications equipment (transmitters) to the oxygen concentrator in accordance with the following guidelines, which are based on the maximum output power of the communication equipment.

Maximum calculated output power of the transmitter, W	Distance to equipment according to transmitter frequency (m)		
	150 kHz to 80 MHz $d=1.2\sqrt{P}$	80MHz to 800 MHz $d=1.2\sqrt{P}$	800 MHz to 2.5 GHz $d=2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters the maximum rated output power of which is not listed in the table, the recommended distance (d) in meters (m) is determined by the equation applicable to the frequency of the transmitter, where P is the maximum rated output power of the transmitter in watts (W) according to the manufacturer's data.

Note 1: At 80 and 800 MHz, the distance applicable to the higher frequency range should be used.

Note 2: These guidelines may not apply in all situations. The propagation of electromagnetic waves is affected by their absorption by various structures, as well as their reflection from various structures, objects and people.

Warranty and liability of the parties

The manufacturer warrants that the system is free from factory defects and that the device will function properly in accordance with its specifications for a period of 3 years from the date of sale to the retailer. The molecular sieve and compressor are warranted for thirty six months (36) or twenty six (26,000) operating hours.

The manufacturer guarantees that for stationary oxygen concentrators, model HO₂RIZON S5, which are being serviced by the manufacturer or at a certified technical service center, the warranty will be extended by 90 days, starting from the day the device is returned for service.

Accessories are warranted to be free from defects for a period of 90 days from the date of purchase. If for some reason the product does not work in accordance with its specifications, then the authorized representative of the manufacturer Scaleo Medical LLC undertakes, at its choice, either to repair or replace the product.

This warranty does not cover problems caused by accident, misuse, malfunction, or defects not related to material or manufacturing process.

The warranty is automatically terminated in the following cases:

- Any tampering with the device contrary to this user manual
- Improper and inappropriate use of the device
- Any modification of the device
- Any technical intervention by an unqualified or non-certified technician.

The manufacturer does not accept any responsibility for any injury or harm, as well as their consequences, directly or indirectly caused to users, patients or third parties in the following cases:

- Failure to follow the instructions and recommendations of this user manual.
- Use of inappropriate spare parts.
- Assembly, adjustment and repair service was carried out by unqualified specialists.
- Improper use of equipment, negligence, accident, human error, or care of equipment with inappropriate materials.

This warranty supersedes any other warranty, including the seller's warranty for product quality and adaptation to use in a specific environment. In addition, the manufacturer will not be held liable for loss of profit, loss of reputation or incidental or consequential damages, even if the possibility of such a situation is advised.

Description of the medical device marking

Marking should be done in a typographic or printed manner on a label, or applied to surfaces, indicating:

1. Name or type designation (species, model);
2. Device number according to the manufacturer's numbering system, if necessary (serial number);
3. Date of manufacture;
4. Trademark (if any);
5. The nominal concentration of oxygen in the produced gas, expressed as a percentage by volume of oxygen, at a flow of 2 l/min;
6. Rated power (W);
7. Voltage and frequency;
8. Flow (l/min);
9. Concentration (nominal flow rate);
10. Sound pressure level dB.

Note - It is allowed to apply manipulation signs on the marking in accordance with this manual.

The shipping container marking must contain the following:

- name of product;
- product designation (article - if any);
- name and location of the manufacturer (legal address, including country, and, if it does not coincide with the legal address, address(s) of production) and organization in the Russian Federation authorized by the manufacturer to accept claims from consumers on its territory (if any);
- Configuring flow control;
- Rated power (W);
- Voltage and frequency;
- Flow (l/min);
- Concentration (nominal flow rate);
- Sound pressure level dBA;
- Net weight (kg)
- Dimensions (mm)
- Pressure of the oxygen-air mixture at the outlet;
- Date of manufacture (month, year);
- batch/batch number (if any);
- barcode (if any);
- storage and transportation conditions;
- registration certificate number and date of issue;
- conformity mark (if any);
- trademark (if any).

Note - It is allowed to place on the marking manipulation signs explaining the method of transportation, loading and unloading and storage of products, and others, in accordance with this manual.

Complete set of a medical device

The assembly should include:

Stationary oxygen concentrator model HO₂RIZON S5, composed of:

1. Main unit assembly - 1 pc.
2. Humidifier - from 1 to 5 pcs. (if necessary)
- 2.1 Humidifier connecting tube - from 1 to 10 pcs. (if necessary)
3. Compressor filter - from 1 to 10 pcs. (if necessary)
4. Porous sponge filter - from 1 to 10 pcs. (if necessary)
5. Nasal cannula - from 1 to 10 pcs. (if necessary) design options:

5.1 Nasal oxygen cannula MEDEREN, manufactured by Mederen Neotech Limited, Israel (RU No. RZN 2018/7623 dated 17.09.2018).

5.2 Nasal oxygen cannula, manufactured by "Plasti-Med", Turkey, (RU No. FSZ 2012/12924 dated 01.02.2021)

6. Fire damper - from 1 to 10 pcs. (if necessary)
- 6.1 Fire damper tube - from 1 to 10 pcs. (if necessary)
7. User manual - 1 pc.

Stationary oxygen concentrator, model HO₂RIZON S5, is completed with a nasal cannula (c. 5) in the Russian Federation by the authorized representative of the manufacturer, Scaleo Medical LLC .

Nasal cannula RU No. RZN 2018/7623 of 17.09.2018 (manufactured by Mederen Neotech)
Nasal oxygen cannula RU FSZ 2012/12924 of 01.02.2021 (manufacturer "Plasti-Med").

Products Included in the delivery set with the indication "if necessary" can be supplied to the customer **both in packaging** with the product, and in separate packages.

Service life and disposal of the device

Service life

The **expected service life** of the stationary oxygen concentrator, model HO₂RIZON S5 is 5 years.

Device disposal

Stationary oxygen concentrator, model HO₂RIZON S5 and its accessories are disposed of in accordance with current legislation.

Recycling directive WEEE / RoHS

If the regulations on the e-waste declaration or the RoHS recycling directive apply, the device must be properly disposed of separately from other industrial waste. Contact your certified distributor or manufacturer for device disposal instructions.

Current regulations

This apparatus has been designed in accordance with the following standards:

GOST 31214-2016 "Medical devices. Requirements for samples and documentation submitted for toxicological, sanitary-chemical tests, tests for sterility and pyrogenicity "

GOST ISO 10993-10-2011 Medical devices. Evaluation of biological effects of medical devices. Part 10. Irritation and sensitization studies (amendment)

GOST P 52770-2016 "Medical devices. Safety requirements. Methods of sanitary-chemical and toxicological tests ".

GOST ISO 10993-1-2011 "Medical devices. Evaluation of biological effects of medical devices. Part 1. Assessment of the study".

GOST ISO 10993-11-2011 Medical devices. Evaluation of biological effects of medical devices. Part 11. General toxicity studies

GOST P 50444-2020 Medical devices, apparatus and equipment. General technical conditions".

GOST R IEC 60601-1-2-2014 "Medical electrical device. Part 1-2. General safety requirements, taking into account the basic functional characteristics. Parallel standard. Electromagnetic compatibility. Requirements and tests."

GOST R IEC 60601-1-6-2011 Medical electrical device. Part 1-6. General safety requirements. Usability

GOST 30324.0-95 (IEC 601-1-88) Medical electrical device. Part 1. General safety requirements (authentic GOST R 50267.0-92 (IEC 601-1-88))

GOST IEC 60601-1-8-2011 Medical electrical device. Part 1-8. General safety requirements. General requirements, tests and guidelines for the use of alarm systems for medical electrical devices and medical electrical systems

GOST R ISO 15223-1-2014 "Medical products. Symbols used for marking on medical devices, labels and accompanying documentation. Part 1. Primary requirements".

GOST 31056-2002 (ISO 8359:2008) Oxygen concentrators for medical use. Safety requirements.

GOST R IEC 62304-2013 Medical devices. Software. Life cycle processes.

Manufacturer: SysMed (China) Co., LTD.

Address: No.17, Wensu Street, Hunnan New District, Shenyang, China, Китай

Tel. 024-24697305 info@sysmed.cn

Authorized representative of the manufacturer in the Russian Federation: Limited Liability Company "Scaleo Medical" (Scaleo Medical LLC)

Legal address: 117513, Moscow, Leninsky prospect, house 139, underground floor No. 1, premises VI, room 3. info@scaleomedical.com contact@scaleomedical.com

Mailing address: 115172, Moscow, Narodnaya, 9, PO Box 1.

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