INSTRUCTIONS FOR USE

Sperm cryopreservation medium
GM501 SpermStore

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Approved by:

12. Nov. 2019 Dr. Fabian Sell, Managing Director

GYNEMED GmbH & Co. KG, Lübecker Strasse 9, D-23738 Lensahn, Germany

Phone: +49 0 4363/9 03 29-0;

Fax: +49 0 4363/9 03 29-19

E-mail: info@gynemed.de

Web-site: <http://www.gynemed.de>

CONTENT

1. NAME OF THE MEDICAL PRODUCT	3
2. PURPOSE, SCOPE OF THE MEDICAL PRODUCT APPLICATION	3
3. COMPOSITION AND DESCRIPTION OF THE MEDICAL PRODUCT	3
4. DATA ABOUT MATERIALS OF ANIMAL AND (OR) HUMAN ORIGIN, AS WELL AS MEDICAL PRODUCTS BEING PART OF THE MEDICAL DEVICE.....	4
5. INDICATIONS.....	5
6. CONTRAINDICATIONS	5
7. POSSIBLE ADVERSE EFFECTS AND FIRST AID MEASURES IN CASE OF INCORRECT USE	6
8. SAFETY PRECAUTIONS WHEN USING.....	6
9. DESCRIPTION OF WORK PRINCIPLES. RECOMMENDATIONS FOR USE	6
10. SPECIFICATIONS.....	10
11. PRODUCT FORM	10
12. STERILITY	10
13. METHODS AND MEANS OF DISINFECTION AND PRESTERILIZATION CLEANING	11
14. TRANSPORTATION CONDITIONS	11
15. STORAGE CONDITIONS. SHELF LIFE	11
16. OPERATING CONDITIONS	11
17. LIST OF NATIONAL STANDARDS APPLIED BY THE MANUFACTURER OF THE MEDICAL DEVICE	11
18. ENVIRONMENTAL PROTECTION REQUIREMENTS	12
19. PROCEDURE AND CONDITIONS FOR DISPOSAL OF THE MEDICAL PRODUCT	12
20. PRODUCER. GUARANTEE OBLIGATIONS	12
21. AUTHORIZED REPRESENTATIVE OF THE MANUFACTURER IN THE RUSSIAN FEDERATION	13
PACKAGING SYMBOLS	13

1. NAME OF THE MEDICAL PRODUCT

Sperm cryopreservation medium GM501 SpermStore - 1 vial (20 ml) (hereinafter - "medium", "product").

Catalogue number:

GM501 SpermStore - 1 vial (20 ml)	4 SCP-20
-----------------------------------	----------

2. PURPOSE, SCOPE OF THE MEDICAL PRODUCT APPLICATION

The medium is intended for cryopreservation of human spermatozoa (testicular and epididymal) by IVF specialists, it can be used in combination with intrauterine insemination (IUI), in vitro fertilization (IVF) and intracytoplasmic sperm injection (ICSI).

Application area: obstetrics and gynecology.

The potential consumers of this medical product are women undergoing IVF, IUI, ICSI procedure, regardless of whether male or female in the couple is infertile. Potential users are IVF specialists (laboratory technicians, embryologists and doctors).

3. COMPOSITION AND DESCRIPTION OF THE MEDICAL PRODUCT

Sperm cryopreservation medium GM501 SpermStore is a ready-to-use buffer HEPES medium that also contains physiological salts, glycine, dextrose monohydrate, calcium lactate, glycerin, sucrose, and human serum albumin to protect spermatozoa from damage caused by the freezing procedure.

Name of the component	Content, %
Water (H ₂ O)	70.0 – 90.0
Glycerin (C ₃ H ₈ O ₃)	10.0 – 20.0
HEPES (4-(2-hydroxyethyl)-1-piperazine ethanesulfonic acid) (C ₈ H ₁₈ N ₂ O ₄ S)	< 1,0
Sucrose(C ₁₂ H ₂₂ O ₁₁)	1 - 5
Sodium chloride (NaCl)	< 1.0
Potassium chloride (KCl)	< 1.0
Magnesium sulfate (MgSO ₄ * 7H ₂ O)	< 1.0
Sodium dihydrogen phosphate (NaH ₂ PO ₄ * 2H ₂ O)	< 1.0
Sodium hydrocarbonate (NaHCO ₃)	< 1.0
Glycine (C ₂ H ₅ NO ₂)	< 1.0
Dextrose monohydrate (C ₆ H ₁₂ O ₆ * H ₂ O)	< 1.0
Calcium lactate (Ca (CH ₃ CHOHCOO) ₂ * 5H ₂ O)	< 1.0
Human serum albumin (HAS) (3.95 g/l)	< 0.4

Component	Functional meaning
Physiological salts	Osmotic regulation
Sucrose, lactate	Energy source
Bicarbonate	pH buffer
Glycine	Organic osmolyte
Human serum albumin (HAS)	Several functions, the main of which is stabilization of the medium and cryoprotection of cells

The medium is a stable, colorless (or yellow), transparent, odorless liquid, completely soluble in water. It has no hazardous decomposition products.

4. DATA ABOUT MATERIALS OF ANIMAL AND (OR) HUMAN ORIGIN, AS WELL AS MEDICAL PRODUCTS BEING PART OF THE MEDICAL DEVICE

Sperm cryopreservation medium GM501 SpermStore contains human serum albumin, obtained from human blood plasma. Advantage of adding HAS:

- macromolecule with an important cryoprotective role;
- prevents the accumulation of spermatozoa and their adherence to laboratory equipment;
- counteracts the harmful effects of reactive oxygen intermediates (ROI).

When producing, Plasbumin is used-25 (produced by Grifols/Talecris) or, as an alternative, Alburnorm 25% (produced by Octapharma) (human albumin solution in accordance with Article 0255 of the EP).

Protective measures during production

1. Upon receipt of raw materials, control of the following indicators: endotoxin concentration, albumin concentration and mouse embryo analysis.

2. Certification of each batch concerning that:

- the batch was manufactured in accordance with the GMP regulations and meets the criteria for production and testing of products in accordance with the European Pharmacopoeia;
- all donated plasma was tested separately, and no reaction to the hepatitis B surface antigen, HIV-1/2 Ab and HCV Ab was observed;
- each plasma reserve was tested and recognized negative for the surface antigen of the hepatitis B virus, HIV-1/2 Ab and HCV-RNA through PCR.
- in case of product recall, we will be notified by their suppliers (Grifols, Talecris or Octapharma).
- proper attention and tracking system on our part.
- in case of changes in the albumin production process and PMF modifications, we will be notified.

Plasbumin-25

- The original plasma used for production meets the requirements set forth in the article of the European Pharmacopoeia (EP: human plasma for fractionation, and is annually approved during the procedure PMF.

- All safety-related aspects of donor exclusion criteria concerning plasma donors are consistent with the criteria specified in Directive 2001/83 / EC as amended by Directive 2004/33 / EC, a series of WHO technical reports, No. 840, 1994 and EMEA / CPMP / BWP / 2879/02 / rev. 1, as well as the recommendations specified in the article of the EP [Council Recommendation No. R (15) on the preparation, use and quality assurance of the blood components, and the Council Recommendation dated June 29, 1998 on the suitability of blood and plasma donors and screening of donated blood in the European Community (98/463 / EC)].

- Donated plasma and its reserves are tested for infection markers, as required by Directive 2001/83 / EC and the relevant US regulations.

- In addition, viral safety is provided by separate stages of removal and / or deactivation of the virus, including precipitation with the solvent (acetone), deep filtration and heat treatment. The virus removal / deactivation steps provide significant clearance from the viruses transferred with the blood and increase the safety level against the risk of virus transmission. Virus safety complies with EU CPMP guidelines. A note to the guidelines on medical products obtained from the plasma, and a

note to the guidelines on the virus test study: development, participation and interpretation of studies confirming deactivation and removal of viruses.

- Grifols / Talecris production facilities operate in accordance with Annex 2 of the Directive GMP EU "Medicinal products for use in humans and in the field of veterinary medicine: good manufacturing practice."
- Warning in the Certificate of Analysis/Release that the medium contains HAS.

Albunorm 25% (alternative to Plasbumin-25)

- All plasma used as the starting material is manufactured and tested in accordance with all applicable US FDA regulations for the production of biological products, as described in Chapter 1 of Section 21 of the Code of Federal Regulations, subsection "Biological products". Each plasma unit is tested for antibodies to HIV-1/2 of human, antibodies to HCV, surface antigens of the hepatitis B virus, and HIV-1 antigen through FDA-approved tests. Each plasma unit undergoes PCR testing, and it is verified that it is not detectable concerning HIV-1, HCV, HBV and HAV. Parvovirus B19 is <1X10⁴ international units per milliliter.
- Stages of reducing viral load during the production of albumin, which include several stages of cold fractionation of ethanol and pasteurization. Researches on the reduction of viral load were carried out by Octapharma in accordance with the current EU regulations. Researches were carried out with both viruses with a shell (HIV-1, pseudorabies, Sindbis virus), and viruses without a shell (reovirus, HAV, MEV, PPV).
- EMA approved plasma Master File (annual re-approval during plasma Master File procedure)
- The human albumin solution is prepared in accordance with the detailed good manufacturing practice guidelines for the starting materials.

5. INDICATIONS

- 1) Female infertility:
 - Female infertility associated with lack of ovulation.
 - Female infertility of tubal origin (associated with congenital anomalies of the fallopian tubes or tubular obstruction).
 - Female infertility of uterine origin (associated with congenital uterine abnormality, egg implantation defects)
 - Female infertility of cervical origin.
 - Female infertility associated with male factors.
 - Other forms of female infertility (endometriosis)
 - Female infertility not otherwise specified
- 2) Male infertility (azoospermia, oligozoospermia)
- 3) Combined infertility

6. CONTRAINDICATIONS

There are no conditions under which the product cannot be used. However, the precautions when using must be considered described in section 8.

7. POSSIBLE ADVERSE EFFECTS AND FIRST AID MEASURES IN CASE OF INCORRECT USE

The product is not considered hazardous to health in accordance with the Directives 67/548/EEC and 1999/45/EC. Does not contain any ingredients that are hazardous to health in accordance with the US Occupational Safety and Health Act regulations making more than 1% of the mixture (more than 0.1% for carcinogens).

If the product is used in violation of the requirements described in the instructions for use, irritation of skin, eyes, mucous membranes and upper respiratory tract is possible.

First aid measures:

- 1) When contacting with the skin: it is necessary to treat the area with a sufficient amount of water.
- 2) When contacting with the eyes: it is necessary to rinse with a sufficient amount of water for 15 minutes. The eyelids of the affected eye when rinsing should be drawn apart.
- 3) When swallowed: it is necessary to rinse the mouth with water and consult a doctor.

If any discomfort occurs, consult a doctor.

8. SAFETY PRECAUTIONS WHEN USING

- 1) Treat all samples as potential HIV or hepatitis carriers.
- 2) Always use protective clothing, chemical protection glasses, rubber gloves when working.
- 3) Avoid contact with skin, eyes and clothing, inhalation and mouth pipetting.
- 4) The medical device does not contain antibiotics. Work under sterile conditions (laminar flow, ISO of class 5) in order to avoid possible contamination.
- 5) Only for intended use.
- 6) Do not re-sterilize.
- 7) When using of the medium during 7 days of package opening its contamination is possible. It is required to store the medium and to work with it under sterile conditions only.
- 8) The medium does not contain substances having genotoxic and carcinogenic effect. The medium does not have toxic effect on the reproductive system.

9. DESCRIPTION OF WORK PRINCIPLES. RECOMMENDATIONS FOR USE

The medium comes into direct physical contact only with human spermatozoa, including epididymal or testicular spermatozoa.

Check before use:

- 1) Do not use the medium if at the time of delivery the vial, container seal or package is open or damaged.
- 2) Do not use the medium if it has changed color, clouded, or has any signs of microbial contamination.

Additional materials:

- 1) Cryostraws for sperm.
- 2) Sterile pipettes.
- 3) Liquid nitrogen tank.

4) LAF cabin (ISO of class 5).

There are no specific recommendations on type of culture cell plates and microtools designed for the use of medium.

Attention:

- 1) Before use, make sure that all media are well mixed.
- 2) Media required for freezing and thawing should be warmed to room temperature before use (+15 - +25 °C).

Before cryopreservation:

In case of a very low sperm concentration, it is recommended to concentrate the sperm before cryopreservation. This can improve sperm quality after thawing and reduce the number of straws required for cryopreservation.

Cryopreservation:

1. Provide sperm liquefaction for 30 minutes at room temperature. (+15 до +25 °C).
2. Add to 1.0 ml of sperm 0.7 ml media GM501 SpermStore - 1. The medium should be added dropwise, gently stirring.

Attention: in order to avoid cold shock, make sure that the medium is at room temperature.

3. Leave the mixture for 10 minutes at room temperature for balancing.
4. Take a mixture of sample and medium in cryostraws, leaving approximately 1.5 cm of air at the end of the straw.
5. Seal the cryostraw.
6. Wipe the straws with a lint free cloth.
7. Shake the straws so that the air bubbles move towards the middle.
8. Freeze the cryostraws in the vertical position for 15 minutes directly above the surface of liquid nitrogen.
9. Store in liquid nitrogen.

Thawing:

1. Collect the required amount of straws from liquid nitrogen.
2. Place the straws in running water for 5 minutes.
3. Cut the end of the straw, place the open end in the test tube and shake the straw so that the whole mixture falls into the test tube.
4. Dilute concentrated sperm sample in the sperm processing medium (at least 3.0 ml per 0.5 ml of sperm), mix thoroughly.
5. Centrifuge for 15 minutes at 300-350 x g.
6. Resuspend the residue in the required volume of sperm processing medium.



After thawing:

If necessary, to remove dead cells and debris (the mass of cell pieces) after thawing, use the sperm preparation techniques. Dilute concentrated sperm sample in the sperm processing medium.

The medium can be used in combination with intrauterine insemination (IUI), in vitro fertilization (IVF) and intracytoplasmic sperm injection (ICSI).

Recommendations for use at In Vitro Fertilization (IVF).

Four-well cell plate	Cup with microdrops
----------------------	---------------------

Day "-1": - prepare a four-well cell plate for fertilization filling each well by 400 mcl of the medium for fertilization and embrio culture with Gentamicin (GM501 Cult); - completely cover the plate by mineral oil; - place the plate in the CO₂-incubator (at 5-7% of CO₂, temperature +37°C). Day «0»: - in each well place the oocytes for fertilization in the amount of not more than 6 pcs. <i>Prepare a four-well cell plate for cultivation in a manner similar to preparing a plate for fertilization.</i> Day «1»: - partially remove cumulus and granulosa cells from the fertilized oocytes; - transfer the oocytes to the plate for cultivation in the amount of not more than 6 pieces per well.	Day "-1": - prepare a cup with microdrops with the diameter of 60 mm for fertilization by applying a few drops of the medium (80-100 mcl*); - completely cover the cup with mineral oil; - place the cup in the CO₂-incubator (at 5-7% of CO₂, temperature +37°C). Day «0»: - place oocytes in each drop for fertilization in the amount of not more than 2 pcs. <i>Prepare a cup with microdrops for cultivation in a manner similar to preparing a cup with microdrops for fertilization.</i> Day «1»: - partially remove cumulus and granulosa cells from the fertilized oocytes; - transfer the oocytes into a cup with microdrops for cultivation in the amount of not more than 3 pieces per drop. <i>*For individual culture, it is recommended to apply the medium drops of 40-50 mcl.</i>
	

Simple "one-step" cultivation

Day "-1"	Day «0»	Day «1»	Day «2» or «3»
Preparation of a plate (or a cup) for fertilization (IVF) with GM501Cult	Fertilization (IVF) and cultivation in the medium GM501 Cult	Further cultivation in the medium GM501 Cult	Embryo transfer



Extended «one-step» cultivation

Day "-1"	Day «0»	Day «1»	Day «2»
Preparation of a plate (or a cup) for fertilization (IVF)	Fertilization (IVF) and cultivation in the medium GM501 Cult	Further cultivation in the medium GM501 Cult	Further cultivation in the medium GM501 Cult. Preparation of a plate (or a

with GM501 Cult			cup) with the fresh culture medium for the day «3»
-----------------	--	--	--

Day «3»	Day «4»	Day «5»
Further cultivation in the medium GM501 Cult	Further cultivation in the medium GM501 Cult	Embryo transfer

Recommendations for use at intracytoplasmic sperm injection (ICSI).

Four-well cell plate	Cup with microdrops
Day "-1": - prepare a four-well cell plate, filling each well by 400 mcl of the medium for fertilization and embryo culture with Gentamicin (GM501 Cult); - completely cover the plate with mineral oil; - place the plate in the CO₂-incubator (at 5-7% of CO₂, temperature +37°C). Day «0»: - immediately after the microinjection, place oocytes into each well for cultivation in the amount of not more than 6 pcs.	Day "-1": - prepare a cup with microdrops with a diameter of 60 mm, applying a few drops of the medium (80-100 mcl*); - completely cover the cup with mineral oil; - place the cup in the CO₂-incubator (at 5-7% of CO₂, temperature +37°C). Day «0»: - immediately after the microinjection, place oocytes into each drop for cultivation in the amount of not more than 3 pcs. <i>*For individual cultivation, it is recommended to apply the medium drops of 40-50 mcl.</i>
	

Simple "one-step" cultivation

Day "-1"	Day «0»	Day «1»	Day «2» or «3»
Preparation of a plate (or a cup) for cultivation (ICSI) with the medium GM501 Cult	Fertilization (ICSI) and cultivation in the medium GM501 Cult	Further cultivation in the medium GM501 Cult	Embryo transfer

Extended «one-step» cultivation

Day "-1"	Day «0»	Day «1»	Day «2»
Preparation of a plate (or a cup) for cultivation (ICSI) with the medium GM501 Cult	Fertilization (ICSI) and cultivation in the medium GM501 Cult	Further cultivation in the medium GM501 Cult	Further cultivation in the medium GM501 Cult. Preparation of a plate (or a cup) with the fresh culture medium for the day «3»

Day «3»	Day «4»	Day «5»
Further cultivation in the	Further cultivation in the	Embryo transfer

Recommendations for use at intrauterine insemination (IUI).

For execution of intrauterine insemination a native previously prepared and cryopreserved sperm of a donor is used.

From recipient's part the term of execution of IUI is availability of at least one passable fallopian tube.

Before execution of IUI it is necessary for donor and recipient to undergo examination determined by clinical recommendations (protocol of treatment).

IUI procedure is carried out in accordance with clinical recommendations (protocol of treatment).

10. SPECIFICATIONS

- 1) All components of the medical device are of the highest degree of purification (European Pharmacopoeia and/or American Pharmacopoeia standard).
- 2) Guarantees and quality control of media are confirmed by the certificates of the organizations MedCert and DEKRA.
- 3) Certificate of analysis for each batch can be obtained on request on the manufacturer's website (<http://www.gynemed.de>) in section "Certificates" by a lot number.
- 4) Safety data sheet (MSDS) of the medical device may be provided on request or downloaded from the manufacturer's website (<http://www.gynemed.de>).
- 5) The production of the medium is carried out in accordance with the following characteristics:

pH	7.20-7.90 (Release criteria: 7.20-7.60)
Sterility	Sterile - SAL 10^{-3} (Level of sterility assurance)
Endotoxins (EU/ml)	< 0.25
Spermatozoon survival test (%)	≥ 80 survival after 4 hours of exposure of the native ejaculate in the test medium

11. PRODUCT FORM

A cardboard pack holds:

- medium supplied in a transparent, sealed, glass vial with a vial filling volume of 20 ml with a tightly closed aluminum cap and a teflon-coated cork Flurotec® in the amount of 1 pc.
- instructions for use.

Internal additional elements additionally fix the position of the vial and protect it from damage.

12. STERILITY

Sperm cryopreservation medium GM501 SpermStore is sterile.

Sterility is ensured by a manual or automated process of aseptic filtration and the vials filling carried out in a work table with a laminar air flow, which allows to obtain SAL < 10^{-3} . Sterility of all media is validated by testing in an independent laboratory and assessing the most unfavorable case.

Vials used as primary packaging as well as stoppers are subject to sterilization by autoclaving method before aseptic filtration of the medium and capping of vials.

The medical device cannot be re-sterilized.

13. METHODS AND MEANS OF DISINFECTION AND PRESTERILIZATION CLEANING

Not applicable.

14. TRANSPORTATION CONDITIONS

Sperm cryopreservation medium GM501 SpermStore are classified as a safe product when transporting.

The medium in the transport container of the manufacturer is transported by any closed type of transport in accordance with the regulations for this type of transport, in compliance with the temperature conditions from +2°C to +8°C. When transporting, the medical devices are packed in polystyrene boxes with cooling elements and cushioning materials.

15. STORAGE CONDITIONS. SHELF LIFE

The shelf life is 18 months from the date of production.

The product is stable after transportation at elevated temperature. ($\leq 37^{\circ}\text{C}$) no more than 5 days.

Store at the temperature from +2°C to +8°C.

Protect from light.

The product should be used within 7 days after the package opening.

Do not use after expiration date.

Do not freeze.

16. OPERATING CONDITIONS

Sperm cryopreservation medium GM501 SpermStore is used at room temperature (+15°C до +25°C).

The opening and closing of the vials should be carried out under sterile conditions (laminar flow, ISO of class 5).

17. LIST OF NATIONAL AND INTERNATIONAL STANDARDS APPLIED BY THE MANUFACTURER OF THE MEDICAL DEVICE

Directive of EU 93/42/EEC	Medical devices, devices, equipment
ISO 13485	Medical devices - Quality management systems - Regulatory requirements
ISO 10993-1	Assessment of medical devices biological action - Part 1: Assessment and testing
GOST ISO 10993-3	Medical devices. Assessment of medical devices biological action. Part 3. Research of genotoxicity, carcinogenicity and toxic action on reproductive function
GOST ISO/TR 10993-33	Medical devices. Assessment of medical devices biological action. Part 33. Guidelines on testing of genotoxicity. Addendum to ISO 10993-3.

ISO 10993-10	Assessment of medical devices biological action - Part 10: Research of irritating and sensitizing action
DN 58940-7	Medical microbiology - Determination of the sensitivity of pathogenic microorganisms to antimicrobial agents - Determination of the minimum bactericidal concentration (MBC) by the method of microdilution
ISO 13408-1	Aseptic production of medical products - -Part 1: General requirements
ISO 13408-2	Aseptic production of medical products - Part 2: Filtration
ISO 13408-6	Aseptic production of medical products - Part 6: Isolating systems
ISO 14644-1	Clean premises and associated controlled media - Part 1: Air cleanliness classification
ISO 14644-3	Clean premises and associated controlled media - Part 3: Testing methods
ISO 14971	Medical devices - Application of risk management to the medical devices
ISO 15223-1	Medical devices - Symbols used when labeling on medical devices, labels and in the accompanying documentation - Part 1: General requirements
EN 556-2	Sterilization of medical devices - Requirements for medical devices to obtain the "Sterile" definition - Part 2: Requirements for medical devices treated under aseptic conditions
MEDDEV 2.1/3	Highly sensitive products, means of medical products delivery and medical devices including additional medicinal substances or human blood derivatives
MEDDEV 2.2/4	Conformity assessment of medical devices for in vitro fertilization (IVF) and assisted reproduction
Eur. Pharm.	European Pharmacopoeia
U.S. Pharmacopeia (USP)	U.S. Pharmacopeia (USP)

18. ENVIRONMENTAL PROTECTION REQUIREMENTS

Sperm cryopreservation medium GM501 SpermStore does not have a negative impact on humans and the environment during use, transportation and storage.

19. PROCEDURE AND CONDITIONS FOR DISPOSAL OF THE MEDICAL PRODUCT

Disposal of the product is carried out in accordance with SanPiN 2.1.7.2790-10 for waste of class B (epidemiologically hazardous waste). After instrument methods of disinfection using physical techniques and changing the appearance of the waste, eliminating the possibility of their reuse, waste of class B can be accumulated, temporarily stored, transported, destroyed and disposed together with waste of class A (epidemiologically safe waste, close to MSW in composition). The packaging of disinfected medical waste of class B should be marked indicating the disinfection of the waste.

The method of cleaning the medium spill site: wash the spill site with water.

20. PRODUCER. GUARANTEE OBLIGATIONS

GYNEMED GmbH & Co. KG,
Lübecker Strasse 9, D-23738 Lensahn, Germany

Phone: +49 0 4363/9 03 29-0;
Fax: +49 0 4363/9 03 29-19
E-mail: info@gynemed.de
Web-site: <http://www.gynemed.de>

Producer GYNEMED GmbH & Co. KG guarantees the compliance of the medical device specifications with all the requirements of the applicable national and international standards provided the consumer observes the conditions of operation, transportation and storage for the entire warranty period.

Guaranteed shelf life - 18 months from the date of production.

The information contained in the safety data sheet is provided without guarantees of any kind. The users should take this information into account only as an addition to other information collected by them, and should conduct an independent determination of the suitability and completeness of the information obtained from all other sources to ensure the proper use and disposal of these materials and the safety and health of the employees and the customers.

GYNEMED GmbH & Co. KG not responsible for any damage resulting from improper treatment or contact with the medical device.

21. AUTHORIZED REPRESENTATIVE OF THE MANUFACTURER IN THE RUSSIAN FEDERATION

MEDWISE LLC,
142717, Moscow Region, Leninsky District, rural settlement Razvilkovskoe, settl. Razvilka,
territory quarter 1, estate 09, office 04
Phone/fax: +7 495 660-76-28
E-mail: office@medwize.ru

PACKAGING SYMBOLS

Symbol	Meaning
	Producer
	Production date
	Use before ...
	Batch code
	Catalog Number
	Sterilization using aseptic processing techniques
	Prohibition on reuse
	Refer to the instruction for use
	Avoid exposure to sunlight
	Temperature range of storage of the medical device

2

CE	Meets the requirements European Council Directive CEE 93/42 for medical devices
----	--

On the secondary packaging of devices designed for supplies to Russian Federation the website is specified, where the instructions for use of the medical device is placed.

[Large handwritten signature in blue ink]

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25. Nov. 2019, Dr. Fabian Sell, Managing Director

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