CALIFORNIA ALL-PURPOSE CERTIFICATE OF ACKNOWLEDGMENT

A notary public or other officer completing this certificate verifies only the identity of the individual who signed the document to which this certificate is attached, and not the truthfulness, accuracy, or validity of that document.

State of California }

County of San Diego }

On Nov. 18. 2021 before me, JOUNG HEE KIM, Notary public,
(Here insert name and title of the officer)
personally appeared JEFFREY ANDERSON

who proved to me on the basis of satisfactory evidence to be the person(s) whose name(s) is/are subscribed to the within instrument and acknowledged to me that he/she/they executed the same in his/her/their authorized capacity(ies), and that by his/her/their signature(s) on the instrument the person(s), or the entity upon behalf of which the person(s) acted, executed the instrument.

I certify under PENALTY OF PERJURY under the laws of the State of California that the foregoing paragraph is true and correct.

WITNESS my hand and official seal.

[Signature]
Notary Public Signature

(Notary Public Seal)

ADDITIONAL OPTIONAL INFORMATION

DESCRIPTION OF THE ATTACHED DOCUMENT

Instruction for use
(Title or description of attached document)

(Title or description of attached document continued)

Number of Pages 1 Document Date 11-18-21

CAPACITY CLAIMED BY THE SIGNER

- ☐ Individual (s)
☒ Corporate Officer

CEO
(Title)

- ☐ Partner(s)
☐ Attorney-in-Fact
☐ Trustee(s)
☐ Other _____

INSTRUCTIONS FOR COMPLETING THIS FORM

This form complies with current California statutes regarding notary wording and, if needed, should be completed and attached to the document. Acknowledgments from other states may be completed for documents being sent to that state so long as the wording does not require the California notary to violate California notary law.

- State and County information must be the State and County where the document signer(s) personally appeared before the notary public for acknowledgment.
- Date of notarization must be the date that the signer(s) personally appeared which must also be the same date the acknowledgment is completed.
- The notary public must print his or her name as it appears within his or her commission followed by a comma and then your title (notary public).
- Print the name(s) of document signer(s) who personally appear at the time of notarization.
- Indicate the correct singular or plural forms by crossing off incorrect forms (i.e. he/she/they, is /are) or circling the correct forms. Failure to correctly indicate this information may lead to rejection of document recording.
- The notary seal impression must be clear and photographically reproducible. Impression must not cover text or lines. If seal impression smudges, re-seal if a sufficient area permits, otherwise complete a different acknowledgment form.
- Signature of the notary public must match the signature on file with the office of the county clerk.
 - ❖ Additional information is not required but could help to ensure this acknowledgment is not misused or attached to a different document.
 - ❖ Indicate title or type of attached document, number of pages and date.
 - ❖ Indicate the capacity claimed by the signer. If the claimed capacity is a corporate officer, indicate the title (i.e. CEO, CFO, Secretary).
- Securely attach this document to the signed document with a staple.

Fantom Encore Sirolimus-Eluting Bioresorbable Coronary Scaffold Reva Medical, Inc.	Instruction for use (Russia) (for registration purposes) Version 02
--	---

Approved By:

Jeffrey Anderson
CEO

Instruction for Use
(Russia)
(for registration purposes)

Fantom Encore Sirolimus-Eluting Bioresorbable
Coronary Scaffold

SEE ATTACHMENT FOR
OFFICIAL NOTARIZATION

Nov. 18. 2021 

Fantom Encore Sirolimus-Eluting Bioresorbable Coronary Scaffold	Instruction for use (Russia) (for registration purposes)
Reva Medical, Inc.	Version 02

1. Medical device name

The Fantom Encore Sirolimus-Eluting Bioresorbable Coronary Scaffold (Fantom Encore scaffold).

2. Manufacturer

REVA Medical, Inc.
5751 Copley Drive
San Diego, CA 92111, USA

Authorized representative of the Manufacturer on the territory of the Russian Federation

SP-CARDIO Ltd.
107076, Russia, Moscow,
Matrosskaya Tishina str., h. 23., bld. 1
Tel.: +7 (903) 616-15-71
e-mail: sp_cardio@rambler.ru

3. Device description

The Fantom Encore Sirolimus-Eluting Bioresorbable Coronary Scaffold (Fantom Encore scaffold) is a medical device, which is a combination product comprised of two regulated components: a bioresorbable coronary scaffold device and a drug product contained in a polymer carrier. The Fantom Encore scaffold is a fully resorbable polymer scaffold. The Fantom Encore scaffold is pre-mounted onto the balloon of a rapid-exchange delivery catheter. After implantation, and as the vessel remodels itself, the scaffold will gradually degrade. As the scaffold degrades and is resorbed, there is a seamless integration of arterial tissue into the space previously occupied by the scaffold.

Device Component Description

The Fantom Encore scaffold is made from a polycarbonate co-polymer of tyrosine analogs. The scaffold polymer material has been formulated to allow for direct visualization under standard fluoroscopy. The radiopacity of the entire scaffold structure is equivalent to currently marketed cobalt chromium devices. The scaffold is intended to be delivered in a guide catheter with a minimum inner diameter of 0.066 in/1.68mm/6Fr. Refer to Table 1 below. The scaffold is currently available in the following sizes:

Table 1. Scaffold Information

Scaffold Size	Minimum I.D. Guide Catheter (in/mm/Fr.)
2.5 mm x 12 mm	0.066/ 1.68/ 6 Fr.
3.0 mm x 12 mm	0.066/ 1.68/ 6 Fr.
3.5 mm x 12 mm	0.066/ 1.68/ 6 Fr.
2.5 mm x 18 mm	0.066/ 1.68/ 6 Fr.
3.0 mm x 18 mm	0.066/ 1.68/ 6 Fr.
3.5 mm x 18 mm	0.066/ 1.68/ 6 Fr.
2.5 mm x 24 mm	0.066/ 1.68/ 6 Fr.
3.0 mm x 24 mm	0.066/ 1.68/ 6 Fr.
3.5 mm x 24 mm	0.066/ 1.68/ 6 Fr.

Fantom Encore scaffold degradation period – 208 weeks.

Delivery System Description

The rapid exchange delivery catheter is a typical industry standard design for delivery systems. The delivery catheter includes two radiopaque markers located on the inner shaft beneath the balloon, proximal and distal to the scaffold, to aid in the accurate placement of the scaffold under fluoroscopy. This balloon catheter is designed to be used with a 0.014" guidewire. The Fantom Encore scaffold is a balloon expandable open-cell design platform that deploys in the same fashion as metal stents, thus requiring no change in clinical practice by the physician. Delivery system usable length is 139 cm.

Drug Component Description

The active pharmaceutical ingredient in the Fantom Encore scaffold is Sirolimus (also known as Rapamycin). Sirolimus is a macrocyclic lactone produced by *Streptomyces hygroscopicus*. The chemical name of Sirolimus (also known as Rapamycin) is (3S, 6R, 7E, 9R, 10R, 12R, 14S, 15E, 17E, 19E, 21S, 23S, 26R, 27R, 34aS)-9, 10, 12, 13, 14, 21, 22, 23, 24, 25, 26, 27, 32, 33, 34, 34a-hexadecahydro-9,27-dihydroxy-3-[(1R)-2-[(1S,3R,4R)-4-hydroxy-3-methoxycyclohexyl]-1-methylethyl]-10,21-dimethoxy-6, 8, 12, 14, 20, 26-hexamethyl-23,27-epoxy-3H-pyrido [2,1-c][1,4]oxaazacyclohentriacontine-1, 5, 11, 28, 29 (4H,6H,31H)-pentone. Its molecular formula is C₅₁H₇₉NO₁₃ and its molecular weight is 914.2. The structural formula of Sirolimus is illustrated as follows:


Figure 1. Sirolimus

Sirolimus is a white to off-white powder and is insoluble in water, but freely soluble in benzylalcohol, chloroform, acetone, and acetonitrile. The target dose for each device size is provided in Table 2.

Table 2. Sirolimus Dose per Scaffold Size

Scaffold Size	Target Dose (µg)
2.5 mm x 12 mm	70
3.0 mm x 12 mm	72
3.5 mm x 12 mm	72
2.5 mm x 18 mm	105
3.0 mm x 18 mm	115
3.5 mm x 18 mm	115
2.5 mm x 24 mm	140

3.0 mm x 24 mm	145
3.5 mm x 24 mm	145

Polymer Carrier

The polymer used for the scaffold structure is the same polymer used as the drug carrier. The polymer is mixed in a solvent with the drug Sirolimus and then applied to the abluminal surface of the scaffold. There is no primer or topcoat layer.

Table 3: In vitro Compliance Information, 2.5mm Diameter

Pressure (atm)	Pressure (kPa)	2.5 x 12 mm	2.5 x 18 mm	2.5 x 24 mm
6	608	2.46 mm	2.46 mm	2.46 mm
7	709	2.50 mm	2.50 mm	2.50 mm
8	811	2.54 mm	2.54 mm	2.55 mm
9	912	2.58 mm	2.58 mm	2.59 mm
10	1013	2.63 mm	2.62 mm	2.64 mm
11	1115	2.67 mm	2.66 mm	2.68 mm
12	1216	2.71 mm	2.70 mm	2.73 mm
13	1317	2.75 mm	2.74 mm	2.77 mm
14	1419	2.79 mm	2.78 mm	2.82 mm
15	1520	2.84 mm	2.82 mm	2.86 mm
16	1621	2.88 mm	2.86 mm	2.91 mm
17	1723	2.92 mm	2.90 mm	2.95 mm
18	1824	2.95 mm	2.94 mm	3.00 mm
19	1925	2.99 mm	2.98 mm	3.04 mm
20	2027	3.03 mm	3.02 mm	3.08 mm

Nominal pressure: 7 atm (709 kPa) for all sizes.

Rated burst pressure: 18 atm (1824 kPa) for all sizes.

Table 4: In vitro Compliance Information, 3.0mm diameter

Pressure (atm)	Pressure (kPa)	3.0 x 12 mm	3.0 x 18 mm	3.0 x 24 mm
6	608	2.96 mm	2.96 mm	2.96 mm
7	709	3.00 mm	3.00 mm	3.00 mm
8	811	3.05 mm	3.05 mm	3.04 mm
9	912	3.09 mm	3.10 mm	3.08 mm
10	1013	3.14 mm	3.15 mm	3.11 mm
11	1115	3.18 mm	3.20 mm	3.15 mm
12	1216	3.23 mm	3.25 mm	3.19 mm
13	1317	3.27 mm	3.30 mm	3.23 mm
14	1419	3.32 mm	3.35 mm	3.27 mm
15	1520	3.36 mm	3.40 mm	3.30 mm
16	1621	3.41 mm	3.45 mm	3.34 mm
17	1723	3.45 mm	3.50 mm	3.38 mm
18	1824	3.50 mm	3.55 mm	3.40 mm
19	1925	3.54 mm	3.60 mm	3.44 mm
20	2027	3.58 mm	3.65 mm	3.48 mm

Fantom Encore Sirolimus-Eluting Bioresorbable Coronary Scaffold	Instruction for use (Russia) (for registration purposes)
Reva Medical, Inc.	Version 02

Nominal pressure: 7 atm (709 kPa) for all sizes.

Rated burst pressure: 18 atm (1824 kPa) for all sizes.

Table 5: In vitro Compliance Information, 3.5mm diameter

Pressure (atm)	Pressure (kPa)	3.5 x 12 mm	3.5 x 18 mm	3.5 x 24 mm
6	608	3.46 mm	3.46 mm	3.46 mm
7	709	3.50 mm	3.50 mm	3.50 mm
8	811	3.55 mm	3.55 mm	3.54 mm
9	912	3.60 mm	3.59 mm	3.68 mm
10	1013	3.64 mm	3.64 mm	3.61 mm
11	1115	3.69 mm	3.69 mm	3.65 mm
12	1216	3.74 mm	3.74 mm	3.69 mm
13	1317	3.79 mm	3.78 mm	3.73 mm
14	1419	3.83 mm	3.83 mm	3.76 mm
15	1520	3.88 mm	3.88 mm	3.80 mm
16	1621	3.93 mm	3.92 mm	3.84 mm
17	1723	3.98 mm	3.97 mm	3.88 mm
18	1824	4.04 mm	4.02 mm	3.92 mm
19	1925	4.08 mm	4.05 mm	3.95 mm
20	2027	4.12 mm	4.09 mm	3.98 mm

Nominal pressure: 7 atm (709 kPa) for all sizes.

Rated burst pressure: 18 atm (1824 kPa) for all sizes.

4. How supplied

STERILE: This device is sterilized with Electron Beam Sterilization. It is intended for single use only. It is non-pyrogenic. Do not use if the package is opened or damaged.

CONTENTS:

- One (1) Fantom Encore Sirolimus-Eluting Bioresorbable Coronary Scaffold, premounted onto the balloon of a rapid exchange delivery catheter, and covered by a removable protective sleeve.
- One (1) Instructions for Use booklet
- One (1) shipping mandrel

TRANSPORTATION: The temperature at transportation should be between 0°C to +25°C with a range of ±2°C and Relative Humidity from 30% to 75% with a range of ±5%. Devices packaged by manufacturer may be transported by all modes of transport in covered vehicles in accordance with the rules of transportation of goods for each type of transport.

STORAGE: Store at 0°C to +25°C, in an enclosed space in the original manufacturer's packaging for the labeled shelf life. Excursions above 30°C must be for less than 30 minutes.

SHELF LIFE: Shelf life of the Fantom Encore Sirolimus-Eluting Bioresorbable Coronary Scaffold is 580 days.

Fantom Encore Sirolimus-Eluting Bioresorbable Coronary Scaffold	Instruction for use (Russia) (for registration purposes)
Reva Medical, Inc.	Version 02

5. Indication for use

The Fantom Encore scaffold is indicated for improving coronary luminal diameter in patients with de novo lesions in native coronary arteries ≥ 2.5 mm to ≤ 3.75 mm in diameter.

Intended use: de novo lesions in native coronary arteries ≥ 2.5 mm to ≤ 3.75 mm in diameter.

6. Contraindications

Coronary artery stenting with the Fantom Encore scaffold is contraindicated in the following:

- Patients in whom antiplatelet and/or anticoagulant therapy is contraindicated.
- Patients with a known hypersensitivity and/or contraindication to aspirin, both heparin and bivalirudin, clopidogrel, ticlopidine, prasugrel, and ticagrelor, propane-1,3-diyl bis(3-(4-hydroxy-3,5-diiodophenyl)propanoate), sirolimus or its structurally related compounds, iodine, or with contrast sensitivity, who cannot be adequately medicated.
- Patients who are judged to have a lesion that prevents complete inflation of an angioplasty balloon.
- Patients with highly calcified lesions.

7. Warnings

- The use of this device carries an associated risk of subacute thrombosis, vascular complications, and/or bleeding events. Therefore, patients should be carefully selected.
- Scaffolding should be generally avoided in those patients at heightened risk of bleeding (e.g., those patients with recently active gastritis or peptic ulcer disease) in which anticoagulation therapy would be contraindicated.
- Persons who are unable to tolerate long-term dual antiplatelet therapy (aspirin indefinitely and clopidogrel (or other P2Y₁₂ inhibitors) for a minimum of 12 months) should not be considered for implantation of the Fantom Encore scaffold.
- Consideration should be given to the risk of allergic reaction to the contrast agent and the benefit of treating the disease when using these devices in patients with a history of severe reaction to contrast agents.
- Persons allergic to propane-1,3-diyl bis(3-(4-hydroxy-3,5-diiodophenyl) propanoate), or iodine may suffer an allergic reaction to this implant.
- Persons with hypersensitivity to Sirolimus or its structurally related compounds may suffer a reaction to this implant.
- Only physicians who have received appropriate training should perform implantation of this device.
- Scaffold placement should only be performed at hospitals where emergency coronary artery bypass graft surgery can be readily performed.
- Subsequent restenosis may require repeat dilation of the arterial segment containing the scaffold. The long-term outcome following repeat dilation of coronary scaffolds is unknown at present.
- Do not expose the delivery system or the scaffold to organic solvents, such as alcohol, or detergents.
- To maintain sterility, the package should not be opened or damaged prior to use.
- The scaffold is not intended to be used in vessels with an average reference diameter < 2.50 mm or in vessels with an individual reference segment < 2.40 mm in diameter.

8. Precautions**8.1. Scaffold Handling – Precautions**

- Contents are supplied STERILE. Do not use if the sterile barrier is damaged.
- For single patient use only. Do not reuse, reprocess or resterilize. Reusing, reprocessing or resterilizing may compromise the structural integrity of the device and/or lead to device failure, which, in turn, may result in patient injury, illness or death. Reusing, reprocessing or resterilizing may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.
- Use prior to the "Use by Date". Store at controlled room temperature.
- The Fantom Encore scaffold is designed for use as a fully intact system. The scaffold is not to be removed from its delivery balloon. The scaffold cannot be crimped onto another balloon. Removing the scaffold from its delivery balloon prior to implantation may damage the scaffold and/or lead to scaffold embolization.
- Special care must be taken not to handle the scaffold in any way to disrupt the scaffold position on the delivery device. This is most important during catheter removal from packaging, removal of the protective sleeve, placement over the guidewire, and advancement through the hemostasis valve adapter and the guide catheter hub.
- Excessive handling may cause dislodgement of the scaffold from the delivery system or scaffold damage.
- Use only the appropriate balloon inflation media (see section 12 Clinical Use Information). Do not use air or any gas medium to inflate the balloon.

8.2. Scaffold Placement – Precautions

- **Guide catheter selection:** The Fantom Encore scaffold is compatible with a 6Fr. (0.066"/1.68 mm inner diameter) guide catheter or larger only. Attempting to place this device using a guide catheter with a smaller diameter than is indicated may result in a failed procedure, damage to the scaffold, scaffold dislodgement and/or a serious adverse event to the patient, including death.
- **Scaffold expansion:** Under no circumstances should an attempt be made to expand the scaffold beyond the dilatation limit (see Table 6). Expansion beyond 3.25 mm for the 2.5 mm Fantom Encore scaffold, beyond 3.75 mm for the 3.0 mm Fantom Encore scaffold, and beyond 4.0 mm for the 3.5 mm Fantom Encore scaffold may result in scaffold damage. **In the event a secondary dilatation catheter is used, REVA recommends the use of a non-compliant balloon that is no greater than 0.5 mm of the nominal scaffold diameter specified in the applicable device compliance chart.** Any attempt to expand the scaffold beyond this specification could result in a serious adverse event for the patient, including death.

Table 6. Dilatation Limits for Fantom Encore Scaffold Sizes

Nominal Scaffold Diameter	Dilatation Limit
2.5 mm	3.25 mm Maximum post-dilatation limit
3.0 mm	3.75 mm Maximum post-dilatation limit

Fantom Encore Sirolimus-Eluting Bioresorbable Coronary Scaffold	Instruction for use (Russia) (for registration purposes)
Reva Medical, Inc.	Version 02

3.5 mm	4.0 mm Maximum post-dilation limit
--------	---------------------------------------

- Do not prepare or pre-inflate balloon prior to scaffold deployment other than as directed.
- All lesions intended for scaffold treatment **MUST** be properly prepared prior to scaffold deployment. Adequate pre-dilation is essential to an optimum clinical outcome.
- Implanting a scaffold may lead to dissection of the vessel distal and/or proximal to the treated portion, and may cause acute closure of the vessel requiring additional intervention (i.e., CABG, further dilation, placement of additional stents, scaffolds, or other).
- Do not expand the scaffold if it is not properly positioned in the vessel (See section 8.3 Scaffold System Removal – Precautions).
- Placement of the scaffold has the potential to compromise side branch patency.
- Balloon pressures should be monitored during inflation. Do not exceed rated burst pressure (RBP) as indicated by the compliance chart on the device label.
- Use of pressures higher than specified on device labeling may result in a ruptured balloon and potential intimal damage and dissection. The scaffold inner diameter should approximate 1:1 times the reference diameter of the vessel.
- If unusual resistance is felt at any time during lesion access, prior to scaffold implantation, the Scaffold System and guide catheter should be removed carefully. (See section 8.3 Scaffold System Removal – Precautions).
- The scaffold system should be introduced into the coronary arteries one time only. The system should not be subsequently moved in and out through the distal end of the guide catheter as scaffold damage, balloon damage, or scaffold dislodgement from the balloon may occur.
- Scaffold retrieval methods (use of additional wires, snares and/or forceps) may result in additional trauma to the vascular site. Complications can include bleeding, hematoma or pseudoaneurysm.

8.3. Scaffold System Removal – Precautions

NOTE: The following steps should be executed under direct visualization using fluoroscopy. When removing the Scaffold System, the entire Scaffold System must be removed carefully: Maintain guidewire placement across the lesion during the entire removal process.

- Carefully pull back the Scaffold System until the proximal balloon marker of the Scaffold System is aligned with the distal tip of the guide catheter.
- Pull back Scaffold System into guide catheter while carefully observing location of the scaffold as it relates to the balloon markers. If movement of the scaffold is observed relative to the balloon marker location, immediately stop retraction into the guide catheter. If any movement of the scaffold is observed, carefully retract the Scaffold System and the guide catheter from the patient as a single unit while leaving the guidewire across the lesion.

Failure to follow these steps and/or applying excessive force to the Scaffold System can potentially result in scaffold damage, scaffold dislodgement from the balloon, damage to the delivery system, and potential injury or death to the patient.

8.4. Post Implant – Precautions

Care must be exercised when crossing a newly deployed scaffold with an intravascular ultrasound (IVUS) catheter, an optical coherence tomography (OCT) catheter, a coronary guidewire, or a

Fantom Encore Sirolimus-Eluting Bioresorbable Coronary Scaffold Reva Medical, Inc.	Instruction for use (Russia) (for registration purposes) Version 02
---	--

balloon catheter, to avoid disrupting the scaffold placement, apposition and/or geometry.

8.4.1. Magnetic Resonance Imaging (MRI)

Fantom Encore is MR Safe. Fantom Encore is comprised of material that is electrically nonconductive, nonmetallic, and nonmagnetic.

8.4.2. Mechanical Atherectomy Devices

The safety and effectiveness of using mechanical atherectomy devices (directional atherectomy catheters, rotational atherectomy catheters) or laser angioplasty catheters in conjunction with Fantom Encore scaffold implantation have not been established.

8.5. Drug Interaction

Sirolimus is known to be a substrate for both cytochrome P-450 3A4 (CYP3A4) and p-glycoprotein (P-gp). Inducers of CYP3A4 and P-gp may decrease Sirolimus concentrations whereas inhibitors of CYP3A4 and P-gp may increase Sirolimus concentrations. Considerations should be given to the potential for drug interaction when deciding to place a Fantom Encore scaffold in a patient who is taking a drug that could interact with Sirolimus, or when deciding to initiate therapy with such a drug in a patient that has recently received a Fantom Encore scaffold. Complete information regarding known drug interactions can be found in the full prescribing information for RAPAMUNE (Sirolimus) ORAL SOLUTION AND TABLETS, Wyeth Pharmaceuticals, Initial U.S. Approval: 1999.

The effect of drug interactions on the safety or efficacy of the Fantom Encore scaffold has not been determined.

9. Potential adverse events

Adverse events (alphabetically) that may be associated with the use of a coronary scaffold or stent in native coronary arteries:

- Abrupt closure
- Access site complications which may require surgical repair
- Acute myocardial infarction
- Allergic reactions to antiplatelet agents, contrast media, or propane-1,3- diyl bis(3-(4-hydroxy-3,5-diiodophenyl)propanoate), or sirolimus and its structurally related compounds
- Aneurysm, pseudoaneurysm
- Angina, stable or unstable
- Arrhythmias
- Arteriovenous fistula
- Bleeding complications, which may require transfusion
- Cardiac arrest
- Cardiac, pulmonary or renal failure or insufficiency
- Cardiac tamponade
- Coronary artery spasm
- Damage to vessel requiring CABG
- Death
- Dissection, perforation or rupture of the coronary artery
- Drug reactions to antiplatelet/anticoagulation agents or contrast media
- Embolism (air, tissue, device or thrombus)
- Emergent or non-emergent Coronary Artery Bypass Surgery

- Failure to deliver the device to the intended area
- Hemorrhage
- Hypotension/hypertension
- Infection and/or pain at the access site
- Inflammation or infection
- Myocardial ischemia and/or infarction
- Pericardial effusion
- Peripheral ischemia
- Pulmonary edema
- Restenosis of treated segment
- Scaffold embolization, breakage or misplacement
- Shock
- Stroke/transient ischemic attack
- Thrombosis
- Total occlusion of the artery
- Vascular complications, including at the entry site, which may require vessel repair
- Ventricular arrhythmias, including ventricular fibrillation or tachycardia
- Vessel spasm

The following additional side effects/complications may be associated with, but not limited to, the use of Rapamycin-type drugs:

- Acne
- Diarrhea or constipation
- Headache
- Increased blood pressure
- Increased cholesterol or triglyceride levels
- Insomnia
- Nausea
- Rash
- Sore or weak muscles or joints
- Tremor
- Upper respiratory or urinary tract infections
- Water retention

10. Drug information

10.1. Mechanism of Action

The mechanism (or mechanisms) by which the Fantom Encore scaffold affects neointimal production as seen in preclinical studies has not been established. Sirolimus inhibits T-lymphocyte activation and proliferation that occurs in response to antigenic and cytokine (Interleukin [IL]-2, IL-4, and IL-15) stimulation by a mechanism that is distinct from that of other immunosuppressants. Sirolimus also inhibits antibody production. In cells Sirolimus binds to the immunophilin, FK Binding Protein-12 (FKBP-12), to generate an immunosuppressive complex. The Sirolimus: FKBP-12 complex has no effect on calcineurin activity. This complex binds to and inhibits the activation of the mammalian Target of Rapamycin (mTOR), a key regulatory kinase. This inhibition suppresses cytokine-driven T-cell proliferation, inhibiting the progression from the G1 to the S phase of the cell cycle.

10.2. Drug Elution and Tissue Deposition

The elution of Sirolimus from the Fantom Encore scaffold has been designed to achieve comparable arterial tissue concentrations established by the Cypher® Sirolimus-Eluting Coronary Stent (a product of Cordis/J&J).

Preclinical studies have shown the amount of Sirolimus the Fantom scaffolds eluted was 53% by 7 days. Fantom scaffolds delivered and sustained a tissue concentration of $>1\mu\text{g/g}$ out to 30 days, which has been shown to be a desirable target to achieve clinical efficacy. The Fantom scaffold with a slightly lower dose of Sirolimus than the CYPHER stents posed no greater Sirolimus toxicity risk versus the commercial CYPHER product.

The Cypher Sirolimus-Eluting Coronary Stent received CE mark approval in 2002 and FDA approval in 2003, and has been evaluated in multiple clinical studies. Following is information regarding the clinical results of the Cypher Stent as published in the Instructions for Use for the CYPHER Sirolimus-Eluting Coronary Stent on RAPTORRAIL Rapid Exchange Delivery System:

The double-blinded SIRIUS (Sirolimus-Eluting Stent in De Novo Native Coronary Lesions) trial randomized 1,058 patients with de novo coronary artery lesions to either the Cypher Stent (533 patients) or a bare metal stent (525 patients).

Five-year data from the SIRIUS (Sirolimus-Eluting Stent in De Novo Native Coronary Lesions) trial confirmed the long-term safety and efficacy of the Cypher stent in patients with simple and medium complexity native coronary lesions. At 5 years, the cumulative MACE, TLR, and TVR rates significantly favored the Cypher Stent (Table 7). There were no significant differences observed in death, MI, and non-target lesion revascularization.

Table 7. SIRIUS Trial: Cumulative Events at 5 Years

	Cypher (n = 533)	BMS (n = 525)	P Value
MACE	20.3%	33.5%	< 0.0001
TLR	9.4%	24.2%	< 0.001
TVR	22.5%	33.5%	< 0.0001

The SIRIUS trial evaluated a metallic version of a sirolimus-eluting stent in a lesion subset that is similar to that indicated for use of the Fantom Sirolimus-Eluting Bioresorbable Coronary Scaffold System. As such, emulation of Cypher drug elution and tissue deposition is a reasonable benchmark for the Fantom Scaffold.

10.3. Immune Suppression Potential

Sirolimus, the active ingredient in the Fantom Encore scaffold, is an immuno-suppressive agent that is also available in oral formulations. The peak systemic blood concentration of Sirolimus following placement of one Fantom scaffold at 1 hour is $1.43 \pm 0.55 \text{ ng/mL}$, which is substantially lower than the maximum concentration (C_{max}) values usually obtained when Sirolimus oral formulations were evaluated in either healthy volunteers ($78.2 \pm 18.3 \text{ ng/mL}$) or in renal transplant patients ($12.2 \pm 6.2 \text{ ng/mL}$). While immune suppression has not been specifically evaluated as a result of implantation of the Fantom Encore scaffold, it may be possible for systemic concentrations of Sirolimus to approach immuno-suppressive levels temporarily, especially in patients who also have hepatic insufficiency or who are taking drugs that inhibit CYP3A4 or P-glycoprotein. This possibility should be considered for such patients, particularly if they are also taking oral Sirolimus, other immunosuppressive agents, or are otherwise at risk for immune suppression. The use of Sirolimus is not recommended in liver or lung transplant patients. Complete information regarding potential immune suppression risks of Sirolimus can be found in the full prescribing information for RAPAMUNE (Sirolimus) ORAL SOLUTION AND TABLETS, Wyeth Pharmaceuticals, Initial U.S. Approval: 1999.

Fantom Encore Sirolimus-Eluting Bioresorbable Coronary Scaffold Reva Medical, Inc.	Instruction for use (Russia) (for registration purposes) Version 02
--	---

10.4. Lipid Elevation Potential

The use of oral Sirolimus in renal transplant patients was associated with increased serum cholesterol and triglycerides that in some cases required treatment. The effect was seen with both low and high dose prolonged oral therapy in a dose related manner. When used according to the indications for use, the systemic Sirolimus concentrations from the Fantom Encore scaffold are expected to be lower than the concentrations usually obtained in transplant patients, but the magnitude and duration of any effect of those concentrations on lipids is not known. Complete information regarding potential lipid elevation and other potential risks of Sirolimus can be found in the full prescribing information for RAPAMUNE (Sirolimus) ORAL SOLUTION AND TABLETS, Wyeth Pharmaceuticals, Initial U.S. Approval: 1999.

10.5. Coronary Artery Surgery – Effect On Anastomoses

There have been rare reports of bronchial anastomotic dehiscence of transplant anastomoses in lung transplant patients who were receiving oral Sirolimus therapy. In a vessel that has recently been implanted with a Fantom Encore scaffold, the Sirolimus concentrations are expected to be several fold higher than systematic Sirolimus concentrations. Therefore, consideration should be given to the possibility that the presence of a Fantom Encore scaffold may compromise the healing of coronary artery vascular anastomoses. Complete information regarding potential effects of Sirolimus on anastomoses and other potential risks can be found in the full prescribing information for RAPAMUNE (Sirolimus) ORAL SOLUTION AND TABLETS, Wyeth Pharmaceuticals, Initial U.S. Approval: 1999.

10.6. Carcinogenesis, Mutagenesis, Impairment of Fertility

The genotoxicity, carcinogenicity, and reproductive toxicity of the Fantom Encore scaffold have not been evaluated. Carcinogenicity studies were conducted by others in mice and rats. In an 86-week female mouse study at Sirolimus doses 30 to 120 times higher than the 2 mg daily clinical dose (adjusted for body surface area), there was a statistically significant increase in malignant lymphoma at all dose levels compared with controls. In a second mouse study at dosages that were approximately 3 to 16 times the clinical dose (adjusted for body surface area), hepatocellular adenoma and carcinoma in males were considered Sirolimus-related. In the 104-week rat study at dosages equal to or lower than the clinical dose of 2 mg daily (adjusted for body surface area), there were no significant findings. Sirolimus was not genotoxic in the in vitro bacterial reverse mutation assay, the Chinese hamster ovary cell chromosomal aberration assay, the mouse lymphoma cell forward mutation assay, or the in vivo mouse micronucleus assay. Fertility was diminished slightly in both male and female rats following oral administration of Sirolimus at doses approximately 10 times or 2 times, respectively, the clinical dose of 2 mg daily (adjusted for body surface area). In male rats, atrophy of testes, epididymides, prostate, seminiferous tubules and/or reduction in sperm counts were observed. In female rats, reduced size of ovaries and uteri was observed. Reduction of sperm count in male rats was reversible upon cessation of dosing in one study. Testicular tubular degeneration was also seen in a 4-week intravenous study of Sirolimus in monkeys at doses that were approximately equal to the clinical dose (adjusted for body surface area).

11. Patient selection and treatment

11.1. Individualization of Treatment

- The risks and benefits should be carefully considered for each patient before use of the Fantom Encore scaffold. Patient selection factors to be assessed should include a judgment regarding risk of prolonged anticoagulation. Scaffold placement is generally

Fantom Encore Sirolimus-Eluting Bioresorbable Coronary Scaffold Reva Medical, Inc.	Instruction for use (Russia) (for registration purposes) Version 02
---	--

avoided in those patients at heightened risk of bleeding, such as those patients with recently active gastritis or peptic ulcer disease.

- Premorbid conditions that increase the risk of poor initial results or the risks of emergency referral for bypass surgery (diabetes mellitus, renal failure, and severe obesity) should be reviewed.
- Thrombosis following scaffold implantation is affected by several baseline angiographic and procedural factors. These include vessel diameter less than 2.5 mm, vessel thrombosis, poor distal flow, and/or dissection following scaffold implantation. In patients undergoing coronary stenting, the persistence of a thrombus or dissection is considered a marker for subsequent thrombotic occlusion. These patients should be monitored very carefully during the first month after scaffold implantation as thrombosis may occur during this period.
- Antiplatelet drugs should be used in combination with the Fantom Encore scaffold. Physicians should use the information from the FANTOM II clinical trial, coupled with the current literature on drug-eluting stents and scaffolds and the specific needs of individual patients, to determine the specific antiplatelet / anticoagulation dose and duration to be used for their patients in general practice. Administration of dual antiplatelet therapy is recommended pre-procedure and for a period of 12 months post-procedure for patients who are not at high risk of bleeding.
- It is very important that the patient comply with the post procedural antiplatelet recommendations. Premature discontinuation of prescribed antiplatelet medication could result in a higher risk of thrombosis, myocardial infarction, or death. Prior to percutaneous coronary intervention (PCI), if a surgical or dental procedure is anticipated that requires early discontinuation or temporary interruption of antiplatelet therapy, the interventionalist and patient should carefully consider if a sirolimus-eluting scaffold and its associated recommended antiplatelet therapy are the appropriate PCI choice. Following PCI, should a surgical or dental procedure be recommended, the risks and benefits of the procedure should be weighed against the possible risk associated with premature or temporary discontinuation of antiplatelet therapy.
- Patients that require premature discontinuation of antiplatelet therapy secondary to significant active bleeding, should be monitored carefully for cardiac events and, once, stabilized, have their antiplatelet therapy restarted as soon as possible per the discretion of their treating physician.

11.2. Use in Specific Populations

11.2.1. Pregnancy

The Fantom Encore scaffold has not been tested in pregnant or nursing women or in men intending to father children. Effects on the developing fetus have not been studied. While there is no contraindication, the risks of reproductive effects are unknown at this time.

11.2.2. Pediatric Use

Implantation of the Fantom Encore scaffold has not been tested in patients under the age of 18 years.

11.2.3. Additional Considerations

The safety and effectiveness of the Fantom Encore scaffold have not been established in patients with:

- Unresolved thrombus at the lesion site
- A coronary artery with an average reference diameter <2.5 mm or > 3.5mm

Fantom Encore Sirolimus-Eluting Bioresorbable Coronary Scaffold	Instruction for use (Russia) (for registration purposes)
Reva Medical, Inc.	Version 02

- A lesion located in the left main coronary artery, ostial lesions or lesions located at a bifurcation
- Diffuse disease, or poor outflow distal to the identified lesions
- Overlapping scaffolds
- Lesions located in saphenous vein grafts
- Previously placed stents/scaffolds or in-stent/in-scaffold restenosis
- The need for balloon dilation of any cells of a deployed scaffold

12. Clinical use information

12.1. Inspection Prior to Use

Carefully inspect the sterile package before opening. Do not use after the "Use by Date". Verify that the temperature indicators located on the product box and inner product pouch are all colored white. If the Temperature Indicator dots on the product box or inner product pouch have turned pink or red, Do Not Use the device.

If the integrity of the sterile package has been compromised prior to the device "Use by Date" (e.g., damage of the package), or the temperature indicators are not white, contact Customer Service for return information. Do not use if any defects are noted.

NOTE: At any time during use of the scaffold system, if the metal hypotube or any other section of the catheter system has been bent or kinked, do not continue to use the system.

12.2. Materials Required (not provided in Scaffold System package)

- Appropriate guide catheter(s) 6 Fr., minimum 0.066" (1.68 mm) inner diameter
- 20 ml (cc) syringe normal heparinized saline
- ≤0.014 in./0.36 mm guidewire
- Rotating hemostatic valve
- Diluted contrast medium 1:1 with normal heparinized saline
- Inflation Device with pressure gauge
- Torque Device
- Pre-deployment dilation catheter
- Three-way stopcock
- Appropriate arterial sheath

12.3. Vessel and Lesion Selection

- Quantitative imaging is recommended for the assessment of target vessel diameter at baseline for appropriate Fantom Encore size selection.
- The target vessel diameter ranges allowed to be treated in the procedure are indicated in Table 8, along with the Fantom Encore scaffold diameter to be used.

Table 8. Target Vessel Diameter Ranges and Fantom Encore Scaffold to be Used

Target Vessel Diameter Distal and Proximal	Fantom Encore Scaffold Diameter to be Used
≥ 2.5 mm and ≤ 2.75 mm	2.5 mm
≥ 2.75 mm and ≤ 3.25 mm	3.0 mm
≥ 3.25 mm and ≤ 3.75 mm	3.5 mm

- For cases where the combination of target vessel diameter and target lesion length is appropriate to be treated by more than one scaffold size, the selection of scaffold size can be made per the judgment of the physician.

WARNING: The scaffold is not intended to be used in vessels with an average reference diameter

Fantom Encore Sirolimus-Eluting Bioresorbable Coronary Scaffold	Instruction for use (Russia) (for registration purposes)
Reva Medical, Inc.	Version 02

< 2.50 mm or in vessels with an individual reference segment < 2.40 mm in diameter.

12.4. Pre-dilation

Pre-dilate the lesion/vessel with a non-compliant or semi-compliant PTCA balloon. Pre-dilation should be at an approximate balloon to artery ratio of 1:1. The expanded balloon should have a uniform appearance, without any visible indents in the profile.


12.5. Preparation

12.5.1. Packaging Removal

Step	Action
1.	Carefully remove the delivery system from its protective tubing for preparation of the delivery system. When using a rapid-exchange system, do not bend or kink hypotube during removal.

12.5.2. Removal of Protective Sleeve


Step	Action
1.	Being careful not to pinch the scaffold, remove the protective sleeve and shipping mandrel following the instructions below:

Hold the inner and outer sleeve

Slide the blue outer sleeve up until it touches the flare of the shipping mandrel

Remove the entire sleeve and shipping mandrel as one unit

12.5.3. Guidewire Lumen Flush

- | Step | Action |
|------|---|
| 1. | Flush Scaffold System guidewire lumen with normal heparinized saline to remove any air bubbles. |
| 2. | Verify that the scaffold is positioned between the proximal and distal balloon markers. Check for bends, kinks and other damage. Do not use if any defects are noted. |

12.5.4. Delivery System Preparation

Step Action

- | Step | Action |
|------|---|
| 1. | Prepare inflation device with diluted contrast medium. |
| 2. | Attach inflation device to stopcock; attach to inflation port. With rapid-exchange systems, do not bend the hypotube when connecting to inflation device. |
| 3. | Open stopcock to scaffold system and leave on neutral. |

12.5.5. Delivery Procedure

- | Step | Action |
|------|---|
| 1. | Prepare the vascular access site according to standard PTCA practice. |
| 2. | Maintain neutral pressure on inflation device attached to scaffold system. |
| 3. | Backload Scaffold System onto proximal portion of guidewire while maintaining guidewire position across target lesion. |
| 4. | Fully open rotating hemostatic valve to allow for easy passage of the scaffold and prevent damage to the scaffold. |
| 5. | Carefully advance the Scaffold System into the hub of the guide catheter. When using a rapid-exchange system, be sure to keep the hypotube straight. Ensure guide catheter stability before advancing the Scaffold System into the coronary artery. |

NOTE: If unusual resistance is felt before the scaffold exits the guide catheter, **do not force passage**. Resistance may indicate a problem and use of excessive force may result in scaffold damage or scaffold dislodgement from the balloon. Maintain guidewire placement across the lesion, and remove the Scaffold System and guide catheter carefully (See section 8.3 Scaffold System Removal – Precautions).

Fantom Encore Sirolimus-Eluting Bioresorbable Coronary Scaffold	Instruction for use (Russia) (for registration purposes)
Reva Medical, Inc.	Version 02

6. Advance the Scaffold System over the guidewire to the target lesion under direct fluoroscopic visualization. Utilize the proximal and distal radiopaque balloon markers as a reference point. If the position of the scaffold is not optimal, it should be carefully repositioned or removed (See section Scaffold System Removal – Precautions). The inside edges of the marker bands indicate both of the scaffold edges. Expansion of the scaffold should not be undertaken if the scaffold is not properly positioned in the target lesion segment of the vessel.

NOTE: If unusual resistance is felt at any time during lesion access, and before scaffold implantation, the Scaffold System and the guide catheter should be removed carefully (See section 8.3 Scaffold System Removal – Precautions).

12.6. Deployment Procedure

Step	Action
------	--------

- | | |
|----|---|
| 1. | Open stopcock to Scaffold System. Pull negative for 15 seconds, then release to neutral for contrast fill. Repeat until all air is expelled. |
| 2. | Inflate the delivery system expanding the scaffold to the target diameter per the device compliance chart. Higher pressure may be necessary to optimize scaffold apposition to the arterial wall. Refer to the compliance chart on the product label for in vitro scaffold inner diameter and rated burst pressure (RBP). |

CAUTION: Do not exceed rated burst pressure as indicated by the compliance chart on the device label.

- | | |
|----|--|
| 3. | Maintain inflation pressure for 30 seconds for full expansion of the scaffold. |
| 4. | Deflate balloon by pulling negative on inflation device until balloon is fully deflated. |
| 5. | Confirm scaffold position and deployment using standard angiographic techniques. For optimal results, the entire stenosed arterial segment should be covered by the scaffold. Fluoroscopic visualization during scaffold expansion should be used in order to properly judge the optimum expanded scaffold diameter as compared to the proximal and distal coronary artery diameter(s). Optimal expansion requires that the scaffold be in full contact with the artery wall. Reasonable effort should be taken to assure the scaffold is not under dilated. |

6. If scaffold sizing/apposition requires optimization, advance a non-compliant balloon catheter of the appropriate size, to the treated area using standard angioplasty techniques. **CAUTION:** Under no circumstances should an attempt be made to expand the scaffold beyond the dilatation limit (see Table 8 - Dilatation Limits for Fantom Encore Scaffold Sizes). Expansion beyond 3.25 mm for the 2.5 mm Fantom Encore scaffold, beyond 3.75 mm for the 3.0 mm Fantom Encore scaffold, and beyond 4.0 mm for the 3.5 mm Fantom Encore scaffold may result in scaffold damage. In the event a secondary dilatation catheter is used, REVA recommends the use of a non-compliant balloon that is no greater than 0.5 mm of the nominal scaffold diameter specified in the applicable device compliance chart. Any attempt to expand the scaffold beyond this specification could result in a serious adverse event for the patient, including death.

7. Inflate the balloon to the desired pressure while observing under fluoroscopy. Refer to the compliance chart on the product label for in vitro scaffold inner diameter and rated burst pressure.
CAUTION: Do not exceed rated burst pressure as indicated by the compliance chart on the device label.
8. Deflate the balloon.
9. Reconfirm scaffold position and angiographic result. Repeat inflations until the desired result is achieved.
10. Post-dilation of the Fantom Encore scaffold is strongly recommended for optimal scaffold apposition. When performed, inflation with a non-compliant balloon to the intended implant diameter (scaffold to artery ratio of 1:1) is recommended.

CAUTION:

Do not dilate the scaffold beyond the dilatation limit. Expansion beyond 3.25 mm for the 2.5 mm Fantom Encore scaffold, beyond 3.75 mm for the 3.0 mm Fantom Encore scaffold, and beyond 4.0 mm for the 3.5 mm Fantom Encore scaffold may result in scaffold damage.

Table 9. Dilatation Limits for Fantom Encore Scaffold Sizes

Nominal Scaffold Diameter	Dilatation Limit
2.5 mm	3.25 mm Maximum post-dilatation limit
3.0 mm	3.75 mm Maximum post-dilatation limit
3.5 mm	4.0 mm Maximum post-dilatation limit

12.7. Removal Procedure Step Action

- | <u>Step</u> | <u>Action</u> |
|-------------|---|
| 1. | Ensure balloon is fully deflated. |
| 2. | Fully open distal rotating hemostatic valve. |
| 3. | While maintaining guidewire position and negative pressure on deflation device, withdraw Delivery System. |

13. Disposal Instructions

After use, dispose of product and packaging in accordance with hospital, administrative, and/or local government policy or regulation.

Products, which had contact with blood and/or other biological liquids of patients should be disposed and eliminated in accordance with the requirements for handling medical waste, established by SanPIN 2.1.7.2790-2010 for class B and other normative legal acts of the Russian Federation.

Unused products, including the ones which expired (which had no contact with the biological liquids of patients) should be disposed in accordance with the procedure of disposal of solid domestic waste established by the legislation of the Russian Federation for class A by SanPIN 2.1.7.2790.

14. Environment protection requirements

The Fantom Encore Sirolimus-Eluting Bioresorbable Coronary Scaffold does not adversely impact the environment.

Fantom Encore Sirolimus-Eluting Bioresorbable Coronary Scaffold	Instruction for use (Russia) (for registration purposes)
Reva Medical, Inc.	Version 02

15. Essential requirements applicability

The Fantom Encore Sirolimus-Eluting Bioresorbable Coronary Scaffold is developed, manufactured and tested with reference to the following list of standards and regulations:

Table 10: List of Harmonized Standards

Sr. No.	Standards	Title
1.	MDD/93/42/EEC as amended by 2007/47/EC	Medical Devices Directive
2.	EN ISO 13485:2016	Medical Devices – Quality Management Systems – Requirements for regulatory purposes
3.	EN ISO 14971:2019	Medical Devices: Application of Risk management to Medical Devices
4.	ISO 14630:2012	Non-active surgical implants – General requirements
5.	ISO 25539-2:2012	Cardiovascular implants – Endovascular devices Part 2: Vascular stents
6.	EN ISO 11607-1:2006	Packaging for terminally sterilized medical devices. Part – 1: Requirements for materials, sterile barrier systems and packaging systems (ISO 11607-1:2006)
7.	EN ISO 11607-2:2017	Packaging for terminally sterilized medical devices. Part – 2: Validation requirements for forming, sealing and assembly processes (ISO 11607-2:2006, including AMD 1:2014)
8.	EN ISO 11137-1:2016	Sterilization of health care products – Radiation – Part 1: Requirements for development, validation and routine control of a sterilization process for medical devices
9.	ISO 11137-2:2013	Sterilization of health care products – Radiation – Part 2: Establishing the sterilization dose
10.	EN ISO 11737-1:2008 & Cor. 1:2007	Sterilization of medical devices – Microbiological methods – Part 1: Determination of a population of microorganisms of products
11.	ISO 10993-1:2009	Biological evaluation of medical devices – Part 1: Evaluation and testing within a risk management process
12.	ISO 10993-3:2014	Biological evaluation of medical devices – Part 3: Tests for Genotoxicity, carcinogenicity and reproductive toxicity
13.	ISO 10993-4:2002	Biological evaluation of medical devices – Part 4: Selection of tests for interactions with blood
14.	ISO 10993-5:2009	Biological evaluation of medical devices – Part 5: Tests for <i>in vitro</i> cytotoxicity
15.	ISO 10993-10:2010 (E)	Biological evaluation of medical devices – Part 10: Tests for irritation and skin sensitization
16.	ISO 10993-11:2006	Biological evaluation of medical devices – Part 11: Tests for systemic toxicity
17.	ISO 10555-1:2013	Intravascular catheter – Sterile, single-use catheters – Part 1: General requirements
18.	ISO 14155:2011	Clinical Investigation of Medical devices for human subjects:

Fantom Encore Sirolimus-Ekuting Bioresorbable Coronary Scaffold	Instruction for use (Russia) (for registration purposes)
Reva Medical, Inc.	Version 02

		Good clinical practice
19.	ISO 14644-1:2015 (E)	Cleanrooms and associated controlled environments – Part 1: Classification of air cleanliness
20.	ISO 14644-2:2015 (E)	Cleanrooms and associated controlled environments – Part 2: Specifications for testing and monitoring to prove continued compliance with ISO 14644-1
21.	ISO 14644-3:2005 (E)	Cleanrooms and associated controlled environments – Part 3: Test methods
22.	ISO 14644-4:2001 (E)	Cleanrooms and associated controlled environments – Part 4: Design, construction and start-up
23.	ISO 14644-5:2004 (E)	Cleanrooms and associated controlled environments – Part 5: Operations
24.	EN ISO 15223-1:2016	Medical devices. Symbols to be used with medical device labels, labelling and information to be supplied. General requirements

16. Disclaimer of warranty

ALTHOUGH REVA MEDICAL, INC. HAS MANUFACTURED THE FANTOM ENCORE SCAFFOLD UNDER CAREFULLY CONTROLLED CONDITIONS, REVA MEDICAL, INC. AND ITS AFFILIATES (COLLECTIVELY, "REVA") HAVE NO CONTROL OVER THE CONDITIONS UNDER WHICH THE PRODUCT IS USED, THE DIAGNOSIS OF THE PATIENT, THE METHODS OF ADMINISTRATION, OR ITS HANDLING AND STORAGE AFTER SALE. REVA DISCLAIMS ALL WARRANTIES, BOTH EXPRESSED AND IMPLIED, WITH RESPECT TO THE PRODUCT, INCLUDING BUT NOT LIMITED TO ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. REVA SHALL NOT BE LIABLE TO ANY PERSON OR ENTITY FOR ANY MEDICAL EXPENSES OR ANY DIRECT, INCIDENTAL, OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, FAILURE, OR MALFUNCTION OF THE PRODUCT, WHETHER A CLAIM FOR SUCH DAMAGES IS BASED UPON WARRANTY, CONTRACT, TORT, OR OTHERWISE. NO PERSON HAS ANY AUTHORITY TO BIND REVA TO ANY REPRESENTATION OR WARRANTY WITH RESPECT TO THE PRODUCT.

The exclusion and limitations set out above are not intended to and should not be construed so as to contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held to be illegal, unenforceable, or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of the Disclaimer of Warranty shall not be affected.

17. Graphical symbols for medical device labeling

	Date of manufacture		Do not use if package is damaged
	Use by	 25°C (77°F)	Store at or below 25°C

Fantom Encore Sirolimus-Eluting Bioresorbable Coronary Scaffold	Instruction for use (Russia) (for registration purposes)
Reva Medical, Inc.	Version 02

LOT	Batch code	 30°C (86°F) EXCURSIONS	Excursions above 30°C must be <30 minutes
REF	Catalogue number		MR compatible
	One unit per carton		Inner diameter
	Caution		Minimum guide catheter requirements
	Do not reuse		Legal manufacturer
STERILE R	Sterilized using irradiation	MANUFACTURED BY	Manufactured by
	Do not resterilize		Consult Instructions for Use

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