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CERTIFICATE

中国国际贸易促进委员会暨中国国际商会
China Council for the Promotion of International Trade is China Chamber of International Commerce

00646325

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China Council for the Promotion of International Trade
China Chamber of International Commerce

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THIS IS TO CERTIFY THAT: the signature of ZuoJing of
NANCHANG KAIMEI MEDICAL APPARATUS CO.,LTD. and the seal
of the said company on the annexed DOCUMENT are genuine.

China Council for the Promotion
of International Trade

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
证明查询网址 Website for verifying the certificate: <http://www.rzccpit.com/validate.html>

Signed and approved by

Nanchang Kaimed Medical Apparatus Co., Ltd

ZUO JING

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INSTRUCTION FOR USE

INSTRUCTION FOR USE

1. NAME

Disposable sterile tracheostomy tubes

Designs:

1. Disposable sterile tracheostomy tubes with cuff 3.0, 3.5, 4.0, 4.5, 5.0, 5.5, 6.0; 6.5, 7.0, 7.5, 8.0, 8.5, 9.0, 9.5, 10.0;
2. Disposable sterile tracheostomy tubes without cuff, 3.0, 3.5, 4.0, 4.5, 5.0, 5.5, 6.0; 6.5, 7.0, 7.5, 8.0, 8.5, 9.0, 9.5, 10.0;

2. DESCRIPTION

Disposable sterile tracheostomy tubes (hereinafter referred to as tracheostomy tubes) are intended for anesthesia, artificial aeration, respiratory tract clearance and restoration of the respiratory function of the lungs for tracheostomy patients.

Designed for on-site use no more than 7 days.

A disposable sterile tracheostomy tube consists of a connector for connection to the device, a tube with an X-line, a tube for cuff inflation, a flange, and an introducer. The product is a Magill tube.

Tracheostomy tubes are available in the following designs:

1. Disposable sterile tracheostomy tubes with cuff 3.0, 3.5, 4.0, 4.5, 5.0, 5.5, 6.0; 6.5, 7.0, 7.5, 8.0, 8.5, 9.0, 9.5, 10.0;
2. Disposable sterile tracheostomy tubes without cuff 3.0, 3.5, 4.0, 4.5, 5.0, 5.5, 6.0; 6.5, 7.0, 7.5, 8.0, 8.5, 9.0, 9.5, 10.0;

3. INTENDED PURPOSE/INDICATIONS

Disposable tracheostomy tubes are intended for anesthesia, artificial aeration, respiratory tract clearance and restoration of the respiratory function of the lungs for tracheostomy patients.

Designed for on-site use for no more than 7 days.

4. WARNINGS AND CONTRAINDICATIONS

1. Do not use for tension pneumothorax (you can use a computer after closed drainage intubation), weak pulse, heart failure, especially right heart failure, pulmonary bulla, pneumothorax, mediastinal emphysema before drainage, for patients with hemoptysis and myocardial infarction, cardiac edema).
2. Never use if the package is open or damaged.
3. Service life is 5 years from the sterilization date.
4. For single use only. Destroy after use. Do not reuse.

5. CONDITIONS OF USE

Medical institutions, ambulance car

6. SIDE EFFECTS ASSOCIATED WITH THE USE OF MEDICAL DEVICE

Tension pneumothorax, Hypovolemic shock, heart failure (especially the right heart failure), pulmonary bulla, pneumothorax, pneumomediastinum, big hemoptysis patients of myocardial infarction, pulmonary edema

7. APPLICATION

1. Choose the size of the tracheostomy tube suitable for the patient.
2. Inflate the cuff using the minimum volume of air for testing. Minor cuff leaks can be detected when air is inflated into the cuff. If there is a large leak, replace the cuff immediately with a similar one.
3. The product is intended for insertion into a patient's trachea by physicians who have received medical professional training and are approved for medical practice. The X-line on the outside surface of the tube is used to determine the installation depth.
4. Inject the required amount of gas according to the tube size.
5. If the tube is used within 24 hours, deflate to avoid damage, and renew the cuff after 1 hour.
6. After use remove the air, extract, dispose of and discard the product.

8. SCOPE OF SUPPLY

Disposable sterile tracheostomy tubes are supplied in a group package of 10/100 pcs/pack with instructions for use. Each tube has an individual primary packaging.

9. SHELF LIFE

Disposable sterile tracheostomy tubes are intended for use on site for no more than 7 days. The shelf life is five years from the sterilization date.

10. MANUFACTURER GUARANTEED VALUES OF BASIC PERFORMANCE, CHARACTERISTICS (PROPERTIES) OF THE MEDICAL PRODUCT.

Tracheotomy tube size range
In millimeters

Nominal inner diameter, the size	Inner diameter and tolerance
3.0	3.0±0.15
3.5	3.5±0.15
4.0	4.0±0.15
4.5	4.5±0.15
5.0	5.0±0.15
5.5	5.5±0.15
6.0	6.0±0.15
6.5	6.5±0.20
7.0	7.0±0.20
7.5	7.5±0.20
8.0	8.0±0.20
8.5	8.5±0.20
9.0	9.0±0.20
9.5	9.5±0.20
10.0	10.0±0.20

Information about the size of the connector for connecting to the equipment: 15 mm

11. MATERIALS OF MEDICAL DEVICE

Product part		Material		
		Name	Brand	Manufacturer
Inflation tube		Polyvinyl chloride	Shuanghua	Shanghai New Shanghua Polymer Material Co.,Ltd., China
		Barium sulfate (radio-opaque strip)	Cimbar	Cimbar (foshan) Trade Co.,Ltd., China
Flange		Polyvinyl chloride	Hopefinder	Hopefinder Polymer Sci. & Tech. Co., Ltd., China
Cuff		Polyvinyl chloride		Hopefinder Polymer Sci. & Tech. Co., Ltd., China
ISO 5356-1 15-mm cone fitting		Polyethylene	Kunlun	Petrochina Daqing Petrochemical Company, China
External tube		Polyvinyl chloride	Hopefinder	Hopefinder Polymer Sci. & Tech. Co., Ltd., China
Pilot balloon	Coating	Polyvinyl chloride	Shuanghua	Shanghai New Shanghua Polymer Material Co., Ltd., China
	Lid	Polypropylene	SINOPEC	Petrochina Jiujiang Branch, China
	Kernel	ABS polymer	TERLURAN	Styrolution Korea Ltd., Korea
	Sealing plug	Silikogel	Juhe	Guangdong Juhe Silicone Material Co., Ltd, China
	Spring	304 stainless steel	KOS	Shenzhen Huijun Metals Co., Ltd., China
Tube fixing tape		Polyester	Sanding Ribbon	Sanding Ribbon Group Co., Ltd., China

12. STORAGE CONDITIONS

Disposable sterile tracheostomy tubes should be stored in a non-corrosive gas, cool, dry, well-ventilated and clean room, with a relative humidity of not more than 80%

Temperature	Relative humidity	Atmospheric pressure
+4°C to +49°C	Up to 80%.	N/A

13. TRANSPORT CONDITIONS

Disposable sterile tracheostomy tubes are transported by all types of transport in accordance with the transportation rules applicable for this type of transport.

Temperature	Relative humidity	Atmospheric pressure
+4°C to +49°C	Up to 80%.	N/A

14. OPERATION CONDITIONS

Temperature	Relative humidity	Atmospheric pressure
+4°C to +49°C	Up to 80%.	N/A

15. CONDITIONS FOR DISPOSAL OR DESTROYING OF THE MEDICAL PRODUCT

In a hospital, warehouse, etc., if necessary, the medical product is disposed of and destroyed in accordance with the rules for handling medical waste regulated by SanPiN 2.1.3684-21 for class "B" (epidemiologically hazardous waste).

16. WARRANTY

The product must be used in accordance with the instruction for use. All questions related to the operation of the product can be addressed to the Authorized Representative of the manufacturer.

The manufacturer guarantees performance throughout the shelf life, subject to storage conditions.

For questions concerning product quality, contact the Authorized Representative of the manufacturer.

17. REQUIREMENTS FOR THE MAINTENANCE AND REPAIR OF THE MEDICAL PRODUCT

N/A.

18. LIST OF NATIONAL STANDARDS

EN ISO 13485	Medical products. Quality management systems. Requirements for regulatory purposes
EN ISO 14971	Medical products. Application of risk management to medical products
EN ISO 10993-1	Medical products. Assessment of the biological effect of medical devices. Part 1. Assessment and research in the risk management process
EN ISO 10993-5	Medical products. Assessment of the biological effect of medical devices. Part 5. Cytotoxicity studies: in vitro methods
EN ISO 10993-7	Medical products. Assessment of the biological effect of medical devices. Part 7. Residual ethylene oxide content after sterilization
EN ISO 10993-10	Medical products. Assessment of the biological effect of medical devices. Part 10. Investigations into irritant and sensitizing effects
EN ISO 10993-12	Medical products. Assessment of the biological effect of medical devices. Part 15. Sample preparation and reference materials
EN 980	Medical products. Graphic symbols used in the labeling of medical devices
EN 1041	Medical products. Information provided by the manufacturer
ISO 14644-1	Cleanrooms and associated controlled environments. Part 2. Requirements for control and monitoring to demonstrate continued compliance

EN ISO 11607-1	Packaging for medical devices subject to terminal sterilization. Part 1. Requirements for materials, sterilization barrier systems and packaging systems
EN ISO 11607-2	Packaging for medical devices subject to terminal sterilization. Part 2. Requirements for validation of forming, sealing and assembly processes
EN ISO 11737-1	Sterilization of medical devices. Microbiological methods. Part 1. Assessment of the population of microorganisms on products
EN ISO 11737-2	Sterilization of medical devices. Microbiological methods. Part 2. Tests for sterility carried out in the definition, verification and maintenance of the sterilization process
EN ISO 11135-1	Sterilization of medical devices. Ethylene oxide. Part 1. Requirements for the development, validation and routine control of a sterilization process for medical devices
ISO 11138-1	Sterilization of medical products. Biological indicators. Part 1. Specifications
ISO 11138-2	Sterilization of medical products. Biological indicators. Part 2. Biological indicators for sterilization with ethylene oxide
EN 62366	Medical products. Application of usability engineering to medical devices.
EN 20594-1	Conical fittings with a 6 % (Luer) taper for syringes, needles and certain other medical equipment - Part 1: General requirements
EN 20594-2	Conical fittings with a 6 % (Luer) taper for syringes, needles and certain other medical equipment - Part 2
EN ISO 5356-1	Anaesthetic and respiratory equipment – Conical connectors – Part 1. Cones and sockets
ISO 5366-1	Tracheotomy tubes. Part 1. Connectors
ISO 5366-2	Tracheotomy tubes. Part 2. Basic requirements
EN 1618	Catheters other than intravascular catheters. Test methods for common properties
EN 1782	Tracheal tubes and connectors
MEDDEV. 2.7.1 Rev.3	Clinical Assessment Guide
EN 15986	Symbols for use in the labeling medical devices. Requirements for labeling of medical devices containing phthalates

19. MANUFACTURER AND AUTHORIZED REPRESENTATIVE OF THE MANUFACTURER IN THE TERRITORY OF THE RUSSIAN FEDERATION

**Manufacturer**

Nanchang Kaimed Medical Apparatus Co., Ltd.
Jinxian Medical Device Industrial Area 331700 Nanchang, Jiangxi Province
PEOPLE'S REPUBLIC OF CHINA.

Manufacturing site:

Nanchang Kaimed Medical Apparatus Co., Ltd.
Jinxian Medical Device Industrial Area 331700 Nanchang, Jiangxi Province
PEOPLE'S REPUBLIC OF CHINA.

AUTHORIZED REPRESENTATIVE OF THE MANUFACTURER IN THE TERRITORY OF THE RUSSIAN FEDERATION

LLC INTMED

LLC INTMED and claim acceptance address:

Room 27, suite I, bld. 12, 1 Perovskaya Street, Moscow 111524

20. SYMBOLS USED ON LABELING

SYMBOL	EXPLANATION
	Lot (batch) number
	Manufactured date (YYYY-MM-DD: year-month-day)
	Manufacturer
	Best before end date (YYYY-MM-DD: year-month-day)
	Do not reuse
	Do not resterilize
	Caution!
	Do not use if the package is damaged
	Natural rubber latex free

	See instructions for use
	EO sterilized
	CE mark
Size	Size
O.D.	Outer diameter
I.D.	Inner diameter
	Temperature limit
	Fragile, handle with care
	Keep dry
	Keep away from light

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