**КОПИЯ****“APPROVE”****"Novintethical Pharma SA"****Director***(position)***Fabio Tommassi Rosso***(name)**(signature)***12/04/ 2016***«day» month (numerals)***Stamp****NOVINTETHICAL PHARMA, S.A.****Via Pian Scairolo 11
CH-6915 Pambio-Noranco
Lugano (Switzerland)****Leaflet for medical device for gut lining protection “Adiarin[®]”**

ention!

document includes information that is a trade secret of the company. These parameters are not determined and controlled by the patient when using the product.

This document is intended for professionals with a medical education and can be provided on request.

Name

Medical Device for protection of gut lining "Adiarin®".

Registration certificate: N

Information about MD Development and Manufacture.

Medical device Developer:

Novintethical Pharma SA"

Address (Location) of the Developer:

Via Pian Scairolo, 11 -6915 Lugano Pambio-Noranco, Switzerland

Medical device Manufacturer:

"Novintethical Pharma SA"

Address (Location) of the Manufacturer:

Via Pian Scairolo, 11 -6915 Lugano Pambio-Noranco, Switzerland

Manufacturing Plant:

LABOMAR, S.R.L., Istrana (TV) Via Nazario Sauro 35/l CAP 31036, Italy

D Purpose

Protection of gut lining against bacteria and toxins.

Proposed Clinical Use

For protection of gut lining in diarrhea for adults and children.

Description

Medical device for protection of gut lining "Adiarin®" is designed for protection of gut lining against bacteria and toxins.

The functional purpose of medical device for gut lining protection "Adiarin®" is not realized via pharmacologic, immunologic, genetic or metabolic effect of human body. Medical device for gut lining protection "Adiarin®" has a property of mechanical action protecting tight intercellular junctions of epithelial cells when they demonstrate incompetence or are destroyed thus preventing massive injection of bacteria and toxins to gut lining from intestinal lumen. This

ect is achieved exclusively mechanically when using the medical device for gut lining
otection "Adiarin®".

Mode of Administration

Pour out the sachet contents to a table spoon, mix with a small amount of water and use as is;
Or dissolve the sachet content in 50-100 ml of boiled water cooled to room temperature and use
the final suspension with the formed sediment.

Recommendations for Use by Different Patient Categories			
Age	From 3 months to 3 years*	From 3 to 14 years	Over 14 years
Quantity required for single use	1 sachet	From 1 to 2 sachets**	2 sachets**
Use interval	From 4 to 6 hours	From 4 to 6 hours	From 4 to 6 hours
Time to discontinue use	gut wall physiological function recovery		
*- after consulting with healthcare professional			
** subject to gut wall physiological status			

MD Use Conditions

MD can be used by individuals w/o a degree in medicine. The shelf life is 5 years from the
date of manufacture.

Contra-Indications

Individual intolerance to Medical Device components.
No other known contra-indications exist.

MD Use Precautions

Cleaning / assembly measures:
Moisten the medical device with a sprayer and collect it mechanically. Avoid dust generation.
If required, collect same with a tested and approved vacuum cleaner. Ensure proper ventilation.

General protection and hygienic measures:
Use general security facilities when handling chemical substances. Clean dirty clothes. Do not
blow off / brush off dirt. Wash hands before break and at the end of job.

Respiratory protection devices:
Not required when the room is adequately ventilated.

Eye protection devices:
In case of dust accumulation or inadequate ventilation: Tight goggles.

Before medical device for gut lining protection "Adiarin[®]" use the following shall be taken into account:

Doctor's advice is not required before medical device for gut lining protection "Adiarin[®]" use. Doctor's advice is required in case of heavy or stable symptoms, basically with children below 3 y. o. and older adults.

In case of diarrhea in children below 3 y. o., mainly when there are other symptoms, it is recommended to seek for professional medical advice.

Plenty of drinks and standard diet for diarrhea is recommended.

Though there are no known side effects it is recommended not to use medical device for gut lining protection "Adiarin[®]" across pregnancy or in first months of breast-feeding. At any rate, always consult a doctor before using the product.

Do not use medical device for gut lining protection "Adiarin[®]" after expiration date specified on medical device package.

Do not use medical device for gut lining protection "Adiarin[®]" if its package is damaged.

Inhalation:

ensure access to fresh air, in case of complaints, consult a doctor.

Skin contact:

Medical device for gut lining protection "Adiarin[®]" usually does not irritate skin.

Eyes contact:

Rinse eyes with plenty of water.

MD Classification

- 1) All-Russian Medical device Classifier Code— 93 9800
- 2) Risk class – 2a according to Russian Regulation
- 3) Scope of medical use. — Gastroenterology.
- 4) Frequency of use — Single dose MD.
- 5) Invasiveness — invasive MD.
- 6) Sterility — non-sterile MD, sterilization-free.
- 7) Contact with human body — MD of short-time (less than 24 hours) contact with mucosal linings.

Basic Characteristics of MD:

Stable and original chemical compound – gelatin tannate consisting of animal-based gelatin and tannin. Chemical characteristics of the compound which by nature may not be deemed to be a "new chemical compound" are described below:

Gelatin

Gelatin used for gelatin tannate production is of type A or acid gelatin with pH 5.5 in solution and gelling power of 40/60 Bloom degrees (viscosity 27 mPs).

Gelatin is a pure, natural, cholesterol-free, odorless and flavorless eupeptic protein. Gelatin is 14-90 % collagen protein by dry weight and includes up to 2% of water. The protein consists of 18 amino acids. Gelatin is a purified protein produced via partial hydrolysis of collagen; this is a fibrous protein produced from animal skin, connective tissue and bones. No record of major adverse events, apart from rare instances of hypersensitivity, was found in study of nutrition test and biochemical experiment results.

GRAS Select Committee draws the following conclusion: the current information about gelatin contains no evidence demonstrating or giving sound grounds to suspect danger to people when used in quantities in which it is used at the moment or can be reasonably used in future. Porcine or bovine gelatin is used for production of the medical device for gut lining protection "Adiarin®". Medical device for protection of gut lining "Adiarin®" marketed in Russia contains porcine gelatin.


Tannin:

Molecular formula: $\text{HC}_{14}\text{H}_9\text{O}_9$

Relative molecular mass: 321.22 uma

Molecular structure:

Figure 1

Tannin (gallotannic acid) is tannin normally produced from gall nuts, i.e. excrescence galls formed on new oak branches *Quercus infectoria* Oliver and allied species *Quercus* L. (Fagaceae family), from tara achenocarps (*Caesalpinia spinosa*) or excrescence galls or sumac leaves (any of *Rhus* genus). Tannin is a natural polyphenol present in all plants and thus is an integral part of human diet. Daily human food ration may provide 1-2.5 g of edible tannins. GRAS Select Committee draws the conclusion that the current information about tannin (hydrolysable gallotannins) contains no evidence demonstrating or giving sound grounds to suspect danger to people when used in quantities and in a manner in which it is used at the moment.

Presentation: this product is manufactured in form of powder (gelatin tannate) pre-packed in sachets.

Composition: 1 sachet contains 250 mg of gelatin tannate.

Physical and Chemical Properties:

Color: yellowish-white. Earthy appearance

Odor: Little to no odor.

Contents weight: $(250 \pm 10\%)$ mg per sachet.

Weight loss: less than 12% by weight Melting point: ca. 140 °C.

Boiling point: not determined.

Flashing point: not applicable.

Ignition: MD does not ignite.

Inflammation: MD does not get inflamed. However, inflammable air/dust mixture formation is possible.

Density at 20 °C: 1.03 g/cm³ (DIN 53 243).

Microorganisms Content:

Total aerobic microorganisms below 10³ CFU/g

Total revivable fungi below 10² CFU/g

Stability and Responsiveness

Thermal decomposition: subject to observance of operation, storage and transportation conditions no thermal decomposition occurs.

Labeling information

Labeling on a sachet (primary pack) contains the following information:

- product name;
- manufacturer trade mark;
- composition;
- method of use.

Labeling on a carton (secondary pack) contains the following information:

- product name;
- composition;
- indications for use;
- method of use;
- storage conditions;
- name and address of manufacturer;
- information on marketing authorization holder;

- information on marketing authorization number;
- name and address of a party accepting the claims;
- batch number;
- manufacture date;
- expiry term “Best before..”.

Labeling on a transport box contains the following information:

- handling instructions:

	Temperature limit
	Protect from a moisture
	Recyclable material
	Protect from sunshine
	Top
	Fragile. Carefully
	Manufacturer

- consignee name
- destination point name
- consignor name
- departure point name
- mass (net and gross) / quantity of product packs
- overall dimensions

Packaging Data

The medical device is packed up in sachets made of complex material, i.e. paper-aluminum-polyethylen.

Sachets are made of complex aluminum material consisting of the following:

Stucco paper 50 g/m²;

Laminated 9-micron annealed aluminum;

Polyethylen – laminated LDPE (low density polyethylene) 25 g/m² of pharmaceutical grade, manufactured by AMCOR)

Sachets are linked in pairs and divided by vertical perforated line. A carton contains 8 sachets with prescribing information. Cartons are placed in a transport cardboard box which should provide safety of product packs when transporting and storing. The box can have next dimensions, depending on quantity of supplied products: 325 x 261 x 130 mm, 245 x 210 x 190 mm, 360 x 265 x 210 mm. A packing sheet containing names, their quantities, packing date and conditional packer number is enclosed in the transport box.

List of International Regulatory Documents/Standards to Which MD Conforms

ISO 13485:2003	Medical devices -- Quality management systems -- Requirements for regulatory purposes
Guideline 93/42/CEE with Guideline 2007/47/CE	European Medical Devices Directive
UNI EN 980 (last edition)	Symbols For Use In The Labelling Of Medical Devices - Part 1: General requirement
ISO 14971:2012	Medical devices. Application of risk management to medical devices
EN 1041:2009	Information supplied by the manufacturer of medical devices
EN ISO 10993-1:2010	Biological evaluation of medical devices. Evaluation and testing- Part 1
EN ISO 10993-10:2010	Biological evaluation of medical devices. Evaluation and testing - Part 10
EN ISO 10993-5:2010	Biological evaluation of medical devices. Evaluation and testing - Part 5
EN ISO 14155 :2011	Clinical investigation of medical devices for human subjects -- Good clinical practice
MED DEV 2.12/1 – Rev. 08	GUIDELINES ON A MEDICAL DEVICES VIGILANCE SYSTEM
MED DEV 2.11/1 Rev 2	Application of Council Directive 93/42/EEC taking into account the Commission Directive 2003/32/EC for Medical Devices utilising tissues or derivatives originating from animals for which a TSE risk is suspected

DEV 2.7/1 Rev. 3	Clinical evaluation: Guide for manufacturers and notified bodies
ED DEV 2.7/4	Guidelines on Clinical investigations: a guide for manufacturers and notified bodies
ISO 22442-1:2007	Medical devices utilizing animal tissues and their derivatives — Part 1: Application of risk management
ISO 22442-2:2007	Medical devices utilizing animal tissues and their derivatives — Part 2: Controls on sourcing, collection and handling

Transportation and Storage Conditions

Environmental temperature: from +4 °C to +40 °C.

Relative humidity max. 75 %.

Atmospheric pressure from 70 to 106 kPa.

Medical device for gut lining protection "Adiarin[®]" is transportable by any transportation means in closed vehicles according to freight rules applicable to each means of transport.

Transported medical device is resistant to air-climatic factors as well as to mechanical impact.

Environmental Requirements to MD Use

The MD does not adversely affects the environment.

MD Disposal and Destruction Procedure

Upon shelf life expiration, medical devices shall be disposed according to local statutes and statutory regulations. At the end of the shelf life of a medical device is required to dispose of as household waste.

Warranty Statement

The manufacturer warrants all specified properties and characteristics of medical device for gut lining protection "Adiarin[®]" for five years from the date of manufacture subject to compliance with transportation and storage requirements.

The holder of a registration certificate:

JADRAN - GALENSKI LABORATORIJ dd, Svilno 20, 51000 Rijeka, Croatia.

Objections to quality shall be forwarded to Limited Liability Company "JADRAN"

("Jadran" LLC), 119330, Moscow, Lomonosovsky prospekt 38, office 33 Phone: +7 (495) 970-18-82.

Patient information leaflet

Device: Medical Device for protection of gut lining "Adiarin®".

Marketing Authorization: N

Medical device composition: 1 sachet contains 250 mg of gelatin tannate.

Description: medical device for protection of gut lining "Adiarin®" is designed for protection of gut lining against bacteria and toxins.

The functional purpose of medical device for gut lining protection "Adiarin®" is not realized via pharmacologic, immunologic, genetic or metabolic effect of human body. Medical device for gut lining protection "Adiarin®" has a property of mechanical action protecting tight intercellular junctions of epithelial cells when they demonstrate incompetence or are destroyed thus preventing massive injection of bacteria and toxins to gut lining from intestinal lumen. This effect is achieved exclusively mechanically when using the medical device for gut lining protection "Adiarin®".

Purpose: protection of gut lining against bacteria and toxins.

Approved indication: it is recommended for protection of gut lining in diarrhea for adults and children.

Mode of Administration: pour out the sachet contents to a table spoon, mix with a small amount of water and use as is. Or dissolve the sachet content in 50-100 ml of boiled water cooled to room temperature and use the final suspension with the formed sediment.

Recommendations for Use by Different Patient Categories			
Age	From 3 months to 3 years*	From 3 to 14 years	Over 14 years
Quantity required for single use	1 sachet	From 1 to 2 sachets**	2 sachets**
Use interval	From 4 to 6 hours	From 4 to 6 hours	From 4 to 6 hours
Time to discontinue use	gut wall physiological function recovery		
*- after consulting with healthcare professional			
** subject to gut wall physiological status			

Contra-Indications: individual intolerance to medical device components. No other known contra-indications exist.

Precautions: before medical device for gut lining protection "Adiarin®" use the following shall be taken into account:

Doctor's advice is not required before medical device for gut lining protection "Adiarin®" use. Doctor's advice is required in case of heavy or stable symptoms, basically with children below 3 y. o. and older adults.

of diarrhea in children below 3 y. o., mainly when there are other symptoms, it is recommended to seek for professional medical advice.

of drinks and standard diet for diarrhea is recommended.

gh there are no known side effects it is recommended not to use medical device for gut protection "Adiarin®" across pregnancy or in first months of breast-feeding. At any rate, consult a doctor before using the product.

ot use medical device for gut lining protection "Adiarin®" after expiration date specified on cal device package.

ot use medical device for gut lining protection "Adiarin®" if its package is damaged.

f Life: 5 years. Do not use after expiration of shelf life specified on the package.

age conditions: Store in a dry, protected from light place out of reach of children at a verature from +4 to +40 C⁰

king: 8 sachets are placed in a carton together with package leaflet.

posal and Destruction Procedure: upon shelf life expiration, medical devices shall be posed according to local statutes and statutory regulations. At the end of the shelf life of a dical device is required to dispose of as household waste.

e holder of a registration certificate:

DRAN - GALENSKI LABORATORIJ dd, Svilno 20, 51000 Rijeka, Croatia.

ldress of manufacturer:

ovintethical Pharma SA, Via Pian Scairolo, 11 -6915 Lugano Pambio-Noranco, Switzerland

anufacturing Plant:

ABOMAR, S.R.L., Istrana (TV) Via Nazario Sauro 35/1 CAP 31036, Italy

Objections to quality shall be forwarded to Limited Liability Company "JADRAN"

"Jadran" LLC), 119330, Moscow, Lomonosovsky prospekt 38, office 33 Phone: +7 (495)

70-18-82.

Attention! This document provides only basic information, all optional parameters contained in the instructions for use for specialists with medical education and also can be provided upon request.

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This Apostille only certifies the authenticity of the signature and the capacity of the person who has signed the public document, and, where appropriate, the identity of the seal or stamp which the public document bears.

This Apostille does not certify the content of the document for which it was issued.

[This Apostille is not valid for use anywhere within Spain]

[To verify the issuance of this Apostille, see <https://eregister.justicia.es>]

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