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1. GENERAL INFORMATION

Ketac Molar Easymix Filling Glass Ionomeric X-ray Contrast Lightweight Mixing Kits

Ketac Molar Easymix Glass Ionomer Filling Material - Radiopaque is a conventional, metal-free glass ionomer cement that has been specially designed as filling material for the lateral tooth area. Its improved mechanical properties, its packable consistency and greater radiopacity make it well-suited for linings under composite fillings of cavity classes I and II, for core build up under crowns, for primary tooth fillings and for cavity class I restorations in non-occluded regions. Additional areas of indication include cavity class V fillings when aesthetic aspects do not play a major role, and temporary fillings for cavity classes I and II. Due to its wide range of indications, its quality criteria and its packable consistency Ketac Molar can also be used as a temporary amalgam substitute. Ketac Molar Easymix was developed to improve hand mixing of the powder/liquid and give a high degree of dosing reproducibility. Ketac Molar Easymix provides improved wetting of the powder by the liquid for faster and easier mixing, granulated flowable powder for more precise proportions, reduced dust formation, when opening the powder bottle, dispensing and mixing. Ketac Molar Easymix is successfully used in the A.R.T. (Atraumatic Restorative Treatment) technique which employs only hand instruments with no need for electricity and therefore can be used everywhere.

All Ketac Molar products manufactured by 3M, are glass ionomer cements that on account of their chemical adhesion to the dental enamel and dentin allow a gentle preparation and an especially tight seal at the margins. The material can be applied without lining; it releases fluoride ions, and is radiopaque.

Ketac Molar is, like all conventional glass ionomer cements, a powder/liquid system and is supplied in hand-mixed (Easymix). To achieve high strength properties and a packable consistency while maintaining good bonding characteristics, the grain distribution and the pretreatment of the glass have been specially optimized for Ketac Molar. This means that 50% of the glass particles in Ketac Molar are no greater than approx. 2.8 μm in size, with 90% of the particles having a diameter of less than 9.6 μm . No polyacid is added to

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2.2. International codes and country regulatory classifications
Table 1. International codes and country regulatory classifications:

Global Nomenclature	
GMDN	16704 Cement, Dental, Glass Ionomer
UMDNS	17-619 Restorative Material, Dental, Other
Country	Regulatory Classification
EU	Class IIa per Rule 8, Annex IX, of the Medical Device Directive 93/42/EEC
Russian Federation	2a
Sterility	Non-sterile
Invasiveness	Invasive
Duration of contact	Permanent (more than 30 days)

3. INTENDED USE AND PRODUCT DESCRIPTION

3.1. Intended use

Traditional filling:

Application as a base layer when working according to the methods of open and closed sandwich

- Stump growth.
- Filling of milk teeth.
- Single-surface fillings in non-occlusive areas.
- Filling cavities of class V with a secondary value of the aesthetic effect.
- Temporary single-surface and multisurface fillings.

Dental fillings when working according to the methods of dental treatment of minimally invasive therapy (MI) without a drill:

- filling of primary teeth;
- single-surface fillings in non-occlusal areas;
- filling cavities of the V class (according to Black) with a secondary value of the aesthetic effect;

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base reaction is triggered off during which a salt gel forms. This hardens to form an amorphous material.


Figure 1.

This process is divided into three stages: The initial phase, the primary setting phase and the secondary setting phase.

During the initial phase, the carboxylic acid groups in the polycarboxylic acid dissociate into negatively charged carboxylate ions (COO^-) and positively charged protons (H^+). These hydrogen ions attack the surface of the glass filler, which releases metallic ions in the form of fluoride chelates.


Figure 2.

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
Figure 3.

Following this, the primary setting phase commences - this is where the calcium ions released from the glasses undergo complexing via the carboxylic acid functions of the copolymer chains. A gel forms. The continual attacks by hydrogen ions also cause aluminum ions such as AlF_2^+ to be released from the silicate glass. During this secondary setting reaction, Al^{3+} ions are deposited in the preformed matrix and a three-dimensional Ca-Al carboxylate gel, which is insoluble in water, forms.

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- Subsequently clear the floor of the cavity.

- Use the largest suitable excavator to treat large cavities, and exert slight pressure in order not to inadvertently expose the pulp due to excessive pressure!

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If necessary, cover the areas in close proximity to the pulp by applying spots of a hard-setting calcium hydroxide material to avoid any pulpal irritation.

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- For optimal chemical bonding to the teeth, the smear layer produced in the preparatory steps must be carefully removed by applying Ketac Conditioner or, for use with the A.R.T. technique, Ketac Molar Easymix Liquid to the prepared surfaces and allowing the substance(s) to react for 10 sec.

- Then rinse with copious amounts of water and blow dry in only 2-3 short intervals with water-free and oil-free air or dab dry with cotton pellets.

- The cavity should have a matt shiny appearance. Excessive drying can result in post-operative sensitivity after filling.

- Avoid further contamination; repeat the procedure if surfaces become recontaminated.

Times

At an ambient temperature of 23 °C/73 °F and 50% relative humidity, the following times apply:

	min:sec
Mixing	00:30
Processing from start of mixing	03:00
Setting from start of mixing	05:00

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The mixing ratio by weight is 4.5 parts powder (1 level spoonful) to 1 part liquid (1 drop).

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- Loosen the powder by shaking the bottle.
- Dip the enclosed spoon into the powder and use the plastic scraper at the bottle opening to level it off. Avoid compressing the powder.
- Dose a sufficient amount of powder on a mixing pad.
- Hold the bottle of liquid in a vertical position and dose the appropriate number of drops for the amount of powder next to the powder.
- There must not be any dried residue of the liquid on the nozzle, and the drops must not contain any air bubbles!
- Carefully re-close both bottles after use and ensure that no powder residue remains on the plastic scraper of the powder bottle.

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- Process Ketac Molar Easymix at room temperature (20-25 °C/68-77 °F).
- Use a metal or plastic cement spatula and a mixing pad or glass slab for mixing.
- Ideally, the powder should be transferred into the liquid in not more than two portions.
- Streak out the paste repeatedly until the consistency is homogeneous.
- Moisture can have a negative effect on the properties of the glass ionomer cement. If a homogeneous paste is not obtained with the given mixing ratios, discard the material.

Application

Avoid any and all contamination from water, saliva, etc. and keep the working area dry during the entire time of the application and shaping.

Conventional filling therapy:

- Apply Ketac Molar Easymix in several portions using a suitable instrument. Undercuts may need to be treated first.

- Carefully avoid trapping air bubbles under enamel undercuts.

A.R.T. technique:

- Apply Ketac Molar Easymix in several portions using a suitable instrument. Undercuts may need to be treated first.

- Carefully avoid trapping air bubbles under enamel undercuts.

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- Slightly overfill cavities; apply the material to adjacent grooves and fissures.
- Apply a small amount of Vaseline to your forefinger and use this finger to firmly press the glass ionomer cement into cavities, fissures, and grooves.
- Repeatedly move the tip of the forefinger in a mesio-distal and buccolingual direction to distribute the material over the entire occlusal surface ("finger press technique").
- Finally, remove the finger sideways to prevent the restorative material from lifting out of the cavity.

Fissure sealing with the A.R.T. technique:

The same steps of cleaning, conditioning, application, and processing as in the A.R.T. technique are performed to seal fissures, the only difference being that there is no preparation involved.

Filling Protection

Conventional filling therapy and A.R.T. technique:

The application of a protective layer on the filling is not necessary. In the event that a protective layer is desirable, vaseline and/or cotton rolls can be used as filling protection as described below:

- Apply vaseline to all of the exposed surfaces of the glass ionomer cement immediately after shaping.

Finishing

Conventional filling therapy:

- Beginning no earlier than 5 minutes after the start of mixing, use Arkansas stone, fine-grain diamonds, polishing discs of decreasing grain size, or silicone polishers for finishing.

- If desired, apply vaseline again for surface sealing.

The patient must not expose the filling to any pressure for 1 hour!

A.R.T. technique:

- Remove any excess material with a carver or large excavator as soon as possible without pulling the material from the cavity.
- Before the filling material is completely cured, test the occlusion with some articulation paper.
- Carefully remove any increase in vertical dimension with a carver.
- Especially at the edge of the filling, the occlusion should be slightly out-of-contact.

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- Apply more Vaseline to the working areas to protect the restoration. The patient must not expose the filling to any pressure for 1 hour!

Notes

The material bonds to metal instruments. Residual quantities of the cement should be removed with water before they set on the instrument.

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The glass powder used in Ketac Molar glass ionomer is a very fine, radiopaque aluminium-calcium-lanthanum fluorosilicate glass. The particle size distribution shows that 90 percent of all particles are smaller than 9µm and 10 percent are smaller than 1µm. Fifty percent of the glass particles in Ketac Molar glass ionomer are up to approx. 2.8µm.

Powder

Function	Ingredient	Quantity, %
Reactive filler	OXIDE GLASS CHEMICALS	89 - 91
Reactive component	ACRYLIC ACID, POLYMERS	9.49 - 10.49
Coating	METHACRYLIC ACID, BUTYL ESTER, POLYMER WITH 2-(DIMETHYLAMINO)ETHYL METHACRYLATE AND METHYL METHACRYLATE	< 0.27
Pigment	CI 77492	< 0.1
Pigment	Iron Oxide (Fe ₂ O ₃)	< 0.1
Pigment	CI 77499	< 0.1

Liquid

Function	Ingredient	Quantity, %
Reactive component	ACRYLIC ACID-MALEIC ACID, COPOLYMER	26 - 36
Reaction medium	Water	55 - 64
Setting regulator	Tartaric acid	8.48 - 10.00
Stabilizer	Benzoic acid	< 0.1

Ketac Conditioner

Function	Ingredient	Quantity, %
Solvent	Water	74 - 77
Conditioning component	POLYACRYLIC ACID / ACRYLIC ACID, polymers	24 - 27
Stabilizer	BENZOIC ACID	< 0.1
Stain	C.I. ACID BLUE 9, DISODIUM SALT	≤ 0.1

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Mixing pad:

Ingredient: Paper laminated with polyethylene.

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10. PACKAGING, CATALOG NUMBERS

Trial set: 1 bottle of powder (shade A3) 12.5 g, 1 bottle of liquid 8.5 ml, 1 bottle of conditioner 10 ml, 1 spoon for dosing, 1 pad for mixing.

A trial set for the ART method: 1 bottle of powder (shade A3) 12.5 g, 1 bottle of liquid 8.5 ml, 1 dosing spoon, 1 mixing pad.

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11. TECHNICAL CHARACTERISTICS OF MD

Case with powder bottle cap: Bottle for powder case: Height, mm: 66 ± 3 ; Outer diameter, mm: 33 ± 3 ; Weight, g: $12,5 \pm 10\%$. Powder Bottle Cap: Length, mm: $20 \pm 0,2$; Thread pitch, mm: $3,5 \pm 0,2$, Nominal diameter, mm: 34 ± 1

Ketac Molar Easymix Liquid: Liquid bottle body: Height, mm: 64 ± 3 ; Outer diameter, mm: 24 ± 1 ; Volume, ml: $8,5 \pm 10\%$. Liquid bottle cap: Length, mm: $32 \pm 0,5$; Thread pitch, mm: $2 \pm 0,2$; Nominal diameter, mm: $18 \pm 0,5$.

Ketac Conditioner: Ketac Conditioner bottle case: Height, mm: $52 \pm 0,5$; Outer diameter, mm: $24 \pm 0,5$; Volume, ml: $10 \pm 10\%$. Ketac Conditioner bottle cap:

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14. SHELF LIFE AND EXPLOITATION TIME

14.1. Shelf life

Shelf life: 24 months from the manufacturing date.

14.2. Exploitation time

Not applicable. Medical device is disposable.

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15. USE/OPERATION CONDITIONS OF MD

Storage conditions

The shelf-life of the product is 24 months when stored at 15-25° C/59-77° F.

To protect the powder from moisture, do not open the powder bag until immediately prior to the first application.

Replace the cap securely after use.

Do not use the product after expiration date.

Operating conditions

Ketac Molar shall be used at an ambient temperature of 23-25°C, ambient humidity and pressure values are not regulated:

At a room temperature 23°C and relative humidity 50% mixing time, working time and setting time are the following:

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Processing from start of mixing

Setting from start of mixing

At higher temperatures, the time available for processing is shortened while it is prolonged at lower temperatures (e.g., mixing on a chilled glass slab). The processing time is also shortened if the amount of powder exceeds the recommendation. Exceeding the indicated processing time may lead to reduced bonding to enamel and dentin.

Transport conditions

During transportation the products may experience short term temperature variations without affecting product safety and performance.

Transportation of medical device is carried out by all type of transport.

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17. INFORMATION ABOUT MEDICINES AND (OR) PHARMACEUTICAL SUBSTANCES INCLUDED IN MD
This medical device does not contain drug products or does not introduce medicinal substances.

18. DESCRIPTION OF STERILIZATION METHOD

Sterilization is not applicable. The products are not labelled nor otherwise represented as sterile, nor are they intended to be sterilized by the user.

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24. WARRANTY

3M Deutschland GmbH warrants this product will be free from defects in material and manufacture. 3M Deutschland GmbH MAKES NO OTHER WARRANTIES INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and 3M Deutschland GmbH's sole obligation shall be repair or replacement of the 3M Deutschland GmbH product.

Limitation of Liability

Except where prohibited by law, 3M Deutschland GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

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25. NAME AND LEGAL ADDRESS OF MANUFACTURER AND AUTHORIZED REPRESENTATIVE IN RUSSIAN FEDERATION

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Address: Carl-Schurz-Strasse 1, 41453 Neuss, Germany

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
22.09.2021 82229 Seefeld

C. Schalk
C. Schalk
Regulatory Affairs

I hereby certify that the document
affixed hereto is a true copy
of the original.

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