02

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Registered No. 2018 - 18733

NOTARIAL CERTIFICATE

HYUNDAI LAW FIRM & NOTARY OFFICE

65, Myeongdong-gil, Jung-gu, Seoul

TEL: 02)775-7085 FAX: 02)775-7080

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Collagen membrane **T-Gen**.

2. *Description*

The **T-Gen** (Absorbable collagen membrane) is made of pure collagen, to minimize antigenicity, which is without further cross-linking or chemical additives for excellent biocompatibility. **T-Gen** is a biodegradable collagen membrane derived from porcine pericardium certified by a veterinarian and is purified to prevent immunologic reactions. **T-Gen** has a porous and fibrous microstructure that is hydrophilic and prevents the ingrowth of fibrous tissue into the bone defect. **T-Gen** is sterilized by gamma radiation.

3. *Principle of action*

The membrane has the ability to hold its shape in a wet conditions, which significantly improves the comfort of use. The membrane barrier function allows to prevent connective tissue invasion, and its surface structure features facilitate bone formation process. Unlike non absorbable products, this membrane can be used with no necessity of removal because collagen is fully resorbable.

4. *Intended use*

T-Gen is indicated for guided tissue regeneration procedures in periodontal defects: guided bone regeneration, guided tissue regeneration, augmentation in the implanted area, the reconstruction of the alveolar ridge, the filling of bone defects in maxillofacial area, immediate or delayed augmentation around implants placed in extraction sockets.

5. *Indications*

T-Gen is indicated for guided tissue regeneration procedures in periodontal defects.

6. *Contraindications*

- Should not be used in patients with sensitive reaction to porcine materials.
- Do not use when allergic reaction or infection is observed during use of products. Proper treatment from the instruction by doctors or pharmacists is required.

Note - this list of contraindications is not exhaustive; Absolute contraindications determined by the doctor conducting the treatment and the doctor sent for treatment.

7. *Side effects*

T-Gen may cause allergic reactions in the presence of sensitivity to the porcine materials. Possible complications which can occur with any dental surgery include swelling of the intraoral tissue, bleeding, flap sloughing, bone loss, infection, pain or complications associated with the use of anesthesia.

8. *Information about the acceptance*

MD is supplied in ready for use condition. Installation, assembly is not required.

9. *Procedures*

1. Peel off outer and inner pouches, choose **T-Gen** to fit the wound area and trim it to fit the size of the defect.
2. The template is a convenience item to assist in shaping **T-Gen** and not implantable, which is must be discarded following modification. In order to distinguish between the membrane and the template, the template is marked with the word "TEMPLATE".

3. The wound site of patient should rinse with antimicrobial agent or sterilized saline solution. The **T-Gen** can be used regardless of front and back side.
4. The **T-Gen** can be used in dry or hydrated conditions. If the clinician prefers the handling state of the hydrated condition, the membrane can be hydrated in sterile water or sterilized saline solution for at least 3 minute prior to use.
5. Fixation of the membrane is use with glue, clip and suture.

10. Conditions of application

Item	Range
Temperature in °C	+1 - +30
Humidity in % (non-condensing)	≤ 70
Air Pressure in kPa	60 - 106

11. Precautions

1. Device is sterile if the package is unopened and undamaged. Do not use if the package seal is broken.
2. Do not re-sterilize. Discard all open and unused portions of **T-Gen**.
3. The device must be used prior to the expiration date.
4. Discard device if mishandling has caused possible damage or contamination.
5. Do not reuse. Sterility and function no longer guaranteed in the event of reuse. Cross contamination and infection may occur if reused.
6. **T-Gen** should not be used in the presence of infected wounds.
7. **T-Gen** should only be used by trained dentists and surgeons.
8. **T-Gen** should be used with special caution in patients with acute or chronic infection at the surgical site, uncontrolled metabolic diseases, prolonged corticosteroid therapy, autoimmune diseases, radiotherapy and heavy smoking.
9. If there are any changes in the performance of the membrane (such as infection, pain, any other unusual symptoms), clinician has the responsibility to inform the patient on all proper contraindications, side effects, and precautions. When these conditions occur, clinician should be instructed to see a well trained dental specialist immediately.

12. Classification of MD.

- by the multiplicity of application - disposable;
- the nature of the invasion - invasive MD;
- sterility properties - sterile MD, does not require sterilization before use;
- method of sterilization - Gamma irradiation.
- peculiar properties of the contact with human body - implantable (in wound) the constant contact device (over 30 days).

The materials of construction and special considerations of their contact with human body are presented in Table 1.

Table 1 – Materials


#	Part	Contact	Material	Brand/type	Standard*
1	Membrane	A long-term (above 30 days) contact with the wound surface	High purified collagen derived from porcine tissue**	Type I	ASTM F2212-11
2	Template	A short-term contact with the wound surface	High-density polyethylene	Tyvek / 1059B manufacture d by DuPont	1. Basis weight : ASTM D3776, EN ISO 536 2. exfoliation : ASTM D2724 3. porous : TAPPI T460, ISO 5636-5
3	Primary package	-	a. inner package: PET + easy peel LLDPE / PET + Aluminium foil + easy peel b. outer package: PET + easy peel LLDPE / PET + Aluminium foil + easy peel LLDPE	Structure : PET12 AL9 PE80&PE T12 NY15 PE-8A50 Prospect Island Amcor	1. Seal Strength : ASTM F-88 2. Bond Strength : GB8808 3. Integrity of seal area : ASTM F1929

*The standard whose requirements the materials used comply with

**Use of materials derived from multiple animals. Origin: countries with minimal or limited risk of BSE. Hog breeds: big white breed or landrace. Age: less than 24 months. Feeding: ration does not contain ingredients derived from bovine materials or derivatives

13. Photo

T-Gen collagen membrane

The left figure is a photo image of T-Gen's general view. The right one is a scanning electron microscopy of T-Gen's cross-section area.

14. Technical requirements

Parameter (characteristic)	TG-1	TG-2	TG-3
Width, mm	14.25-15.75	19.00-21.00	28.50-31.50
Length, mm	19.00-21.00	28.50-31.50	38.00-42.00
Depth, mm	< 2	< 2	< 2
Weight, mg	15.0-30.0	34.0-47.0	70.0-120.0



Quality requirements

Parameter (specification)	Requirement
Endotoxin concentration	20EU/device
pH level	basic tested liquid $\pm 1,5$
Tensile strength	no less than 5 N/mm ² (MPa)
Water binding capacity	no less than 200 %
Lipid content	no more than 1,0 %
Color	White
Smell	Odorless

15. Supply package

Components supplied as part of T-Gen package are presented in Table 3.

Table 1 – Components

#	Name, picture and description	Quantity ¹⁾
1	T-Gen collagen membrane   The left figure is a photo image of T-Gen's general view. The right one is a scanning electron microscopy of T-Gen's cross-section area.	1
2	Template	1
3	Instruction for use	1
¹⁾ Quantity of components included in supply package (pcs)		

16. Expiration date

36 months from manufacture date.

17. Term of membrane biodegradation.

Term of membrane biodegradation within 3 ~ 5 months.

18. Storage conditions

1. Store at use only between 1 and 30 °C.
2. Keep away from heat.
3. Keep dry.

Item	Range
Temperature in °C	+1 - +30
Humidity in % (non-condensing)	≤ 70
Air Pressure in kPa	60 - 106

19. Transportation

MD transported by all types of transport in covered vehicles in accordance with the rules and the transportation of cargos, operating in each type of transport. During transportation must be observed packaging integrity.

Item	Range
Temperature in °C	+1 - +30
Humidity in % (non-condensing)	≤ 70
Air Pressure in kPa	60-106

20. Warranty

In compliance with the transportation, storage requirements and use of the product manufacturer guarantees compliance with quality of the medical device.

21. Procedure of recycling and disposal

It is required to provide the mandatory implementation of all necessary procedures in accordance with local applicable laws and regulations when during disposal of MD.

22. Symbols on the Label

#	Type code	Explanation	#	Type code	Explanation
1		Do not re-use	2		Use-by date
3		Do not re-sterilize	4		Sterilized using irradiation
5		Batch code	6		Date of manufacture
7		Caution (see instructions for use)	8		CE mark and identification number of Notified Body
9		Manufacturer	10		Do not use if package is damaged
11		Temperature limit	12		Keep away from sunlight
13		Keep dry	14	Non toxic*	Non toxic

#	Type code	Explanation	#	Type code	Explanation
15	PY № *	Information about Registration Certificate in RF.	16	Уполномоченный представитель в РФ *	Information about authorized representative in Russian Federation.

* Information will be added after registration in Russian Federation.

23. Claim

SK bioland Co., Ltd.

162, Gwahaksaneop 3-ro, Ochang-eup, Cheongwon-gu,
Cheongju-si, Chungcheongbuk-do,
28125, Republic of Korea
TEL +82-43-240-8636, 8637 FAX +82-43-240-8699
www.skbioland.com

Any health care professional (e.g. patient or user) on the territory of the Russian Federation, who have claims respect for the product or is not satisfied with the reliability, safety, effectiveness or performance characteristics of this product, should notify LLC "Implant Guru" by phone, fax or send the corresponding notice by mail to the address below:

LLC "Implant Guru" – Russia, 105062, Moscow, Pokrovka street, bid .43, str. 2
phone/fax: (495) 627-637-0

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
HYUNDAI LAW FIRM
& NOTARY OFFICE

Seoul Central District Prosecutor`s Office
65, Myeongdong-gil, Jung-gu, Seoul

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This office has been authorized by the
Minister of Justice, the Republic of
Korea, to act as Notary Public Since
7, Feb. 2015 Under Law No.90.

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