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INSTRUCTION FOR USE FOR MEDICAL DEVICE

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ELASTIC NAILING IMPLANTS SURGICAL, NON-STERILE

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Official Certification:

The undersigned Notary certifies, that the preceding signature of

Büttler Joanna Zofia, 05.05.1976, von Grenchen, in Günsberg

is authentic and that the signature currently act in the position of Senior RA Program Lead, Synthes GmbH, a company currently existing and organized under the laws of Switzerland, that the address is Elmattstrasse 3, 4436 Oberdorf.

The Notary Public is aware and attests to the above facts upon examination of the corresponding documents.

Given and signed in 4500 Solothurn, Switzerland

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STEN – Stainless Steel Elastic Nail(s)

Please read these instructions for use carefully before use. Ensure that you are familiar with the appropriate surgical technique.

The manufacturer is developing a guide on surgical techniques, which describes the use of implants. Contact the manufacturer's authorized representative for surgical techniques.

Implants are single packed, and available non-sterile.

Important note for healthcare professionals: These instructions for use do not include all of the information necessary for selection and use of a device. Please see full labeling for all necessary information (corresponding surgical technique and device-specific label).

List of MD materials

Stainless steel
Titanium alloy

Intended use

Elastic Nailing Implants are intended to be used for temporary fixation and stabilization of long bones in various anatomical regions such as proximal femur, femoral shaft, tibia and humerus.

Elastic Nailing Implants (TEN and STEN) are used as single implant or in pairs for Elastic Stable Intramedullary Fixation (ESIN).

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Principle of operation

The elastic nails are bent and inserted into the medullary canal. This elastic deformation within the medullary canal creates a bending moment within the long bone that is not rigid, but that is stable enough to reduce and fix the fracture. More details of the biomechanical and biological aspects of the device are documented in the Surgical Technical Guide.

TEN and STEN Nails are used as single implant or in pairs for Elastic Stable Intramedullary Fixation (ESIN).

The nail tip allows nail insertion and sliding along the medullary canal and facilitates manipulation for fracture reduction. The height of the tip guarantees a correct relation to the medullary cavity.

The principal surgical techniques for standard indications are:

- Femur – Ascending Technique
- Femur – Descending Technique
- Tibia – Descending Technique
- Forearm
- Humerus – Ascending Technique
- Humerus – Descending Technique
- Proximal Radius – Radial Neck
- Distal Radius and Ulna Diaphyseal – Metaphyseal, Displaced Fractures
- Clavicle

Potential users of MD

The product is intended for professional use only by healthcare professionals who have passed appropriate training.

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Indications

The Elastic Stable Intramedullary Nailing subject devices have different indications for pediatric and adult populations. Both STEN and TEN variants have indications for pediatric populations, while only the TEN variant is indicated for adults.

Indications in Pediatrics

Elastic stable intramedullary nailing (ESIN) with the Titanium Elastic Nail (TEN) or Stainless Steel Nail (STEN) is indicated for the management of diaphyseal and certain metaphyseal/epiphyseal fractures of long bones in children and young adults. As follows:

- diaphyseal and certain metaphyseal fractures of long bones
- certain metaphyseal/-epiphyseal fractures (Salter Harris I and II), including but not limited to radial neck fractures
- complex clavicular fractures (significant dislocation including shortening, “floating shoulder”)
- open fractures
- threat of skin perforation at fracture ends
- pathologic fractures.

Indications in Adults

In adult patients, TEN is used for the osteosynthesis of clavicle, forearm and humerus fractures follows:

- diaphyseal fractures of long bone fractures in upper extremity
- clavicle shaft fractures.

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Contaminated implants must not be reprocessed. Any Synthes GmbH (hereinafter Synthes) implant that has been contaminated by blood, tissue, and/or bodily fluids/matter should never be used again and should be handled according to hospital protocol. Even though they may appear undamaged, the implants may have small defects and internal stress patterns that may cause material fatigue.

Precautions

Using Elastic Nailing Implants in patients with open epiphysis may impair bone growth. Unless included within the specific indications in the corresponding surgical technique, using Elastic Nailing Implants are therefore not recommended for use in skeletally immature patients.

General Cautions

- Synthes implant that has not been used, but has become soiled with blood, tissue and/or bodily fluids/matter, should be handled according to hospital protocol. Synthes does not recommend the reprocessing of soiled implants.
- Do not use a Synthes implant if the surface has been damaged.
- Synthes implants should not be processed or transported with any type of soiled or contaminated materials.
- All devices must be thoroughly cleaned and inspected prior to sterilization. During cleaning, only use detergents that are labelled for use on medical devices and in accordance with the manufacturer's instructions. Cleaning agents with a used dilution pH of within 7–9.5 are recommended. Highly alkaline conditions (pH >11) can damage devices. Do not use saline,

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environmental disinfection (including chlorine solutions) or surgical antiseptics (such as iodine- or chlorhexidine-containing products). Do not use a cleaning aid that can damage the surface of implants such as steel wool, abrasive cleaners or wire brushes. Refer to instructions below.

- Synthes implants should not be lubricated.
- Synthes implants are critical devices and must be terminally sterilized prior to use.
- The sterilization parameters are only valid for devices that are adequately cleaned.
- The parameters listed are only valid for properly installed, maintained, calibrated and compliant reprocessing equipment in accordance with standards such as the ISO 15883 and ISO 17665 series.

Limits on reprocessing

- Repeated processing cycles that include ultrasonic, mechanical washing and sterilization have minimal effects on Synthes implants.
- Synthes implants should be inspected for corrosion, damage such as scratches and notches, debris, discoloration or residue.
- A discoloration has no adverse effect on titanium or titanium alloy implants. The protective oxide layer is fully maintained.
- Any implant with corrosion, scratches, notches, residue or debris should be discarded.

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Warnings

General warnings

Avoid creating a sharp bend which may reduce the effectiveness of the nail.

- When opening the medial side, be careful not to let the drill bit slip posteriorly into the region of the femoral artery.
- The drill must be running when angling the drill bit or drill bit breakage may result
- Avoid hitting the T-piece of the inserter directly as this may result in damage to the inserter.
- Never rotate the nail more than 180°.
- Excessively long nail ends result in pseudobursa formation and prevent free flexion of the knee. They can also perforate the skin and cause infections.
- The nail end must not be bent away from the cortex if using an elastic nail end cap. Be aware of the extensor tendons and superficial radial nerve. Take care not to penetrate the contralateral cortex. The use of a hammer is not recommended since hammering may produce further fracture fragments. Be aware of the position of the radial nerve in relation to the fracture.

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The possibility and ways of integrating with other medical devices

Synthes has not tested compatibility with devices provided by other manufacturers and assumes no liability in such instances.

Magnetic Resonance environment

Torque, Displacement and Image Artifacts according to ASTM F 2213-06, ASTM F 2052-14 and ASTM F2119-07.

Non-clinical testing of worst case scenario in a 3 T MRI system did not reveal any relevant torque or displacement of the construct for an experimentally measured local spatial gradient of the magnetic field of 3.69 T/m.

The largest image artifact extended approximately 169 mm from the construct when scanned using the Gradient Echo (GE). Testing was conducted on a 3 T MRI system.

Radio-Frequency-(RF)-induced heating according to ASTM F2182-11a.

Non-clinical electromagnetic and thermal testing of worst case scenario lead to peak temperature rise of 9.5 °C with an average temperature rise of 6.6 °C (1.5 T) and a peak temperature rise of 5.9 °C (3 T) under MRI Conditions using RF Coils (whole body averaged specific absorption rate [SAR] of 2 W/kg for 6 minutes [1.5 T] and for 15 minutes [3 T]).

Precautions: The above mentioned test relies on non-clinical testing. The actual temperature rise in the patient will depend on a variety of factors beyond

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the SAR and time of RF application. Thus, it is recommended to pay particular attention to the following points:

- It is recommended to thoroughly monitor patients undergoing MR scanning for perceived temperature and/or pain sensations.
- Patients with impaired thermoregulation or temperature sensation should be excluded from MR scanning procedures.
- Generally, it is recommended to use a MR system with low field strength in the presence of conductive implants. The employed specific absorption rate (SAR) should be reduced as far as possible.
- Using the ventilation system may further contribute to reduce temperature increase in the body.

Treatment before device is used

Synthes products supplied in a non-sterile condition must be cleaned and steam-sterilized prior to surgical use. Prior to cleaning, remove all original packaging. Prior to steam-sterilization, place the product in an approved wrap or container. Follow the cleaning and sterilization instructions below.

Methods and apparatus for disinfection and sterilization cleaning

Cleaning - Manual Method

1. Rinse device under running cold tap water for a minimum of two minutes. Use a soft-bristled brush to clean the device.

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Cleaning-Disinfection Automated Method

1. Rinse devices under running cold tap water for a minimum of one minute. Use a soft-bristled brush or soft lint-free cloth to clean the device.
2. Prepare a fresh detergent solution for the ultrasonic bath using a neutral or mild alkaline detergent. Follow the detergent manufacturer's instructions for the correct dilution, temperature, water quality and exposure time. Note: fresh solution is a newly-made, clean solution.
3. Immerse the implants fully in the detergent. Clean Synthes implants ultrasonically for a minimum of 15 minutes, using a minimum frequency of 38 kHz.
4. Rinse device using DI or PURW water for a minimum of two minutes. DI or PURW water must be used for final rinse.
5. Visually inspect device. Repeat steps 2–5 until devices are visibly clean.
6. Automated washing shall be conducted in a validated washer-disinfector in compliance to ISO 15883-1 and -2, or to an equivalent standard. Load the device components in the washer-disinfector in accordance with manufacturer's instructions. Automated washing can be included as part of a validated washing, disinfection, and/or drying cycle in accordance to manufacturer's instructions:

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* see Additional Information

Additional Information

Further information regarding the use of specific cleaning agents, ultrasonic washers, washer-disinfector, packaging materials or sterilizers during validation studies are available on request. Synthes used the following during validation of these reprocessing recommendations:

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- Manual cleaning: Manual Pre-Cleaning with Prolystica 2X Concentrate Enzymatic Cleaner 1 mL/L at 14–16 °C and Ultrasonic Cleaning Prolystica 2X Concentrate Enzymatic Cleaner 1 mL/L at 12–21 °C.
- Automated Cleaning: Manual Pre-Cleaning with Prolystica 2X Concentrate Enzymatic Cleaner 1 mL/L at 14–16 °C. Water-disinfectant cleaning with (Wash 1) Prolystica 2X Concentrate Enzymatic Cleaner 1 mL/L at 23–26 °C and (Wash 2) Prolystica 2X Neutral Detergent 1 mL/L at 44–46 °C.
- Lint-Free Cloth: Berkshire Durx 670.

The cleaning and sterilization information is provided in accordance with ISO 17664.

The recommendations provided above have been validated by the medical device manufacturer as being capable of preparing a non-sterile Synthes medical device. It remains the responsibility of the processor to ensure that the processing is performed using equipment, materials and personnel in the reprocessing facility, and achieves the desired result. This requires validation and routine monitoring of the process. Likewise, any deviation by the processor from the recommendations provided should be properly evaluated for effectiveness and potential adverse consequences.

The chemical quality of the water used during reprocessing can impact device safety. Facilities should use the recommended water quality requirements for device reprocessing in accordance with local guidance (such as AAMI TIR 34, Water for the reprocessing).

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Thermal disinfection

For automated cleaning-disinfection, thermally disinfect at a minimum of 93 °C for a minimum of 5 minutes (see above section Cleaning-Disinfection Automated Method).

Recommended temperatures and pH levels

Material	Temperature	pH
Stainless steel	up to 149 °C	7–11
Titanium alloys	up to 150 °C	6 – 9,5

Description of sterilization method

Steam (moist heat) sterilization shall be performed in a locally approved, pre-vacuum (forced air removal) cycle. The steam sterilizer should be validated to the requirements of any local standards and guidance such as EN 285 or AAMI/ANSI ST8. The steam sterilizer should be installed and maintained in compliance to manufacturer's instructions and local requirements. Ensure that a steam sterilizer cycle is chosen that is designed to remove air from porous or lumened device loads in accordance to manufacturer's instructions and does not exceed the maximal sterilizer load.

Cycle Type	Minimum Sterilization Exposure Time (minutes)	Minimum Sterilization Exposure Temperature	Minimum Dry Time

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Technical characteristics of MD

Note: The parameters specified in the MD name are part of the name and presented for the convenience of the user's choice of a particular product. The table of technical characteristics contains the accurate parameters of MD

1) STEN-Stainless Steel Elastic Nail Ø 1.5mm, L 300mm (ref. 275.915)

1. Diameter: 1.6 (+0.03/-0.07) mm
2. Usable length: 300 (+/-5) mm
3. Weight: 4.72 (+/-10%) g

2) STEN-Stainless Steel Elastic Nail Ø 2.0mm, L 440mm (ref. 275.920)

1. Diameter: 2 (+0.03/-0.07) mm
2. Usable length: 440 (+/-5) mm
3. Weight: 10.8 (+/-10%) g

3) STEN-Stainless Steel Elastic Nail Ø 2.5mm, L 440mm (ref. 275.925)

1. Diameter: 2.5 (+0.03/-0.07) mm
2. Usable length: 440 (+/-5) mm
3. Weight: 16.8 (+/-10%) g

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- Name and address of authorized representative of the manufacturer in Russia;
- Number and date of Registration certificate;
- Other clarifying text and symbols.

Conditions of use/operation, transportation and storage

Elastic Nailing Implants are intended for use in Healthcare Facilities by medical professionals and operating room staff.

Elastic Nailing Implants are recommended for use under operation conditions indicated in Tables below

Table. Operation conditions.

Environmental conditions	
Temperature	from -29 to +60 °C
Relative humidity	30 - 85 %
Athmosperic pressure	No specific conditions
Operation conditions in the human body	
Temperature	from +32 to +42 °C
Relative humidity	100 %
Atmospheric pressure	No specific conditions

Products should be stored in a dry, clean environment, protected from direct sunlight, pests, and extremes of temperature and humidity.

Table. Storage conditions

Storage conditions	
Temperature	from -29 to +60 °C
Relative humidity	30 – 85 %
Atmospheric pressure	No specific conditions

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Table. Transportation conditions

Transportation conditions	
Temperature	from -29 to +60 °C
Relative humidity	30 – 85 %
Atmospheric pressure	No specific conditions

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Information about implant removal and extraction

Bone healing and implant removal in pediatric patients:

After the fracture is healed, implant removal in children is a common. Unproven long-term complications such as growth disturbance, foreign body reaction, threat of infection, or corrosion are used as arguments for removal. Implant migration or interference with another orthopedic treatment over the later course of life is also a matter of debate. On the other hand, the difficulty in removing single implants as well as possible perioperative complications has induced discussion about the retention of implants in childhood. Clear recommendations for implant removal do not exist. The current line of action still includes routine removal, as it is preferred by some authors, whereas others argue for a selective procedure. Elastic Nailing Implants are usually removed because the ends may interfere with the surrounding tissue. Bone healing in children usually lasts no longer than 1-3 months. In terms of mechanical performance, the lifetime of Elastic Nailing Implants in children is therefore seen as up to 3 months. Implant removal – if indicated – is performed after 1 year at the earliest. In terms of biological parameters, the lifetime of the implant is therefore seen to be 1-1.5 years. Benefits of removal have to outweigh the risks and complications in the individual case and the procedure should not require a more extensive procedure than insertion.

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Bone healing and implant removal in adult patients:

The implant removal if indicated is then performed at a much later stage, or omitted at all. Typically, implants are left in place for more than a year during which the progress of the fracture healing is monitored repeatedly. The risks and benefits of implant removal must be very carefully estimated in patients with poor general health, immune deficiencies (e.g. HIV, hepatitis, tuberculosis), or local circulatory disturbances (e.g. diabetes mellitus, peripheral arterial disease). Therefore, implants are usually left in place in elderly, comorbid patients. Moreover, every surgical procedure, including implant removal, involves a small risk of infection or neurovascular damage.

In summary, bone healing in adults usually lasts no longer than 6-9 months. In terms of mechanical performance, the lifetime of the Elastic Nailing Implants is therefore seen as up to 9 months. Implant removal – if indicated – is performed after 1 or 2 years at the earliest. In terms of biological parameters, the lifetime of the implant is therefore seen to be several years.

If end caps were used, the end caps are removed using the corresponding screwdriver shaft mounted on the inserter for TEN/STEN. Open the old incisions and expose the ends of the nails. Use the extraction pliers for TEN/STEN to first bend the ends of the nails so that they are lifted clear of the formed callus. With the hammer guide for TEN/STEN firmly screwed onto the extraction pliers for TEN/STEN, the nails can easily be removed

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DL-JNJRUMedicalQA@ITS.JNJ.com,

**AUTHORIZED REPRESENTATIVE OF THE MANUFACTURER IN
RUSSIA:**

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