
To
**FEDERAL SERVICE ON SURVEILLANCE
IN HEALTHCARE (ROSZDRAVNADZOR)**
4 Slavyanskaya Square, building 1, Moscow
109074 Russia

CERTIFICATE

I, Futoshi Inoue, President and Representative Director of KITAZATO Medical Co., Ltd., do hereby solemnly and sincerely certify that the following document is correct and can be used instead of the operational documentation for registration of Kitazato ET& IUI Catheters provided previously.

Futoshi Inoue

President and Representative Director

KITAZATO Medical Co., Ltd.

Dated this 01st day of March 2016.

Operational Documentation for Medical Device

【Kitazato ET Catheter】

【Kitazato IUI Catheter】

**Sterile single-use Kitazato Catheters of different types
for Assisted Reproductive Technologies**

1. Name of Medical Device

Sterile single use Kitazato Catheters of different types for Assisted Reproductive Technologies.

Manufacture: KITAZATO Medical Co., Ltd., Shibakoen Building, 1-1-8 Shibadaimon, Minato-ku, Tokyo 105-0012, Japan.

2. Intended Use

Sterile single use Kitazato Catheters of different types for Assisted Reproductive Technologies are used for embryo transfer and artificial intrauterine insemination in assisted reproduction procedures.

3. Configurations and Technical Description

A variety of configurations and parameters of Kitazato Catheters enables to choose an instrument either to match the anatomic features of a patient or to minimize trauma during the embryo transfer or insemination procedures without damage of cells.

Sterile single-use Kitazato Catheters for Assisted Reproductive Technologies are presented in the following types:

1. Kitazato ET Catheter with guide:

- catheter, outer diameter 3 Fr/1,00 mm, length 220 mm; guide, length 170 mm;
- catheter, outer diameter 3 Fr/1,00 mm, length 250 mm; guide, length 200 mm;
- catheter, outer diameter 4 Fr/1,35 mm, length 220 mm; guide, length 170 mm;
- catheter, outer diameter 4 Fr/1,35 mm, length 250 mm; guide, length 200 mm.

2. Kitazato ET Catheter with guide and obturator:

- catheter, outer diameter 3 Fr/1,00 mm, length 220 mm; guide, length 170 mm;
- catheter, outer diameter 3 Fr/1,00 mm, length 250 mm; guide, length 200 mm;
- catheter, outer diameter 3 Fr/1,00 mm, length 400 mm; guide, length 200 mm;
- catheter, outer diameter 4 Fr/1,35 mm, length 220 mm; guide, length 170 mm;
- catheter, outer diameter 4 Fr/1,35 mm, length 250 mm; guide, length 200 mm.

3. Kitazato ET Catheter with memorable guide:

- catheter, outer diameter 4,7 Fr/1,56 mm, length 180 mm; guide, length 140 mm;
- catheter, outer diameter 4,7 Fr/1,56 mm, length 230 mm; guide, length 190 mm.

4. Kitazato IUI Catheter:

- catheter, length 100 mm, hard/intermediate/soft;
- catheter, length 180 mm, hard/intermediate/soft.

5. Kitazato ET Catheter:

- catheter stylet, outer diameter 3 Fr/1,00 mm, length 170 mm for catheter with guide 170 mm;
- catheter stylet, outer diameter 3 Fr/1,00 mm, length 200 mm for catheter with guide 200 mm;
- catheter stylet, outer diameter 4 Fr/1,35 mm, length 170 mm for catheter with guide 170 mm;
- catheter stylet, outer diameter 4 Fr/1,35 mm, length 200 mm for catheter with guide 200 mm;
- catheter stylet, outer diameter 4,7 Fr/1,56 mm, length 135 mm for catheter with guide 140 mm;
- catheter stylet, outer diameter 4,7 Fr/1,56 mm, length 185 mm for catheter with guide 190 mm.

6. Kitazato IUI Catheter:

- catheter stylet, length 118 mm for catheter 100 mm;
- catheter stylet, length 198 mm for catheter 180 mm.

7. Kitazato ET Catheter for trial transfer with guide:

- catheter, outer diameter 4,7 Fr/1,56 mm, length 180 mm; guide, length 140 mm;
- catheter, outer diameter 4,7 Fr/1,56 mm, length 230 mm; guide, length 190 mm.

A description of types of sterile single-use Kitazato Catheters for assisted reproductive technologies is shown below.

3.1. Kitazato ET Catheter with Guide

Kitazato ET Catheter contains of the inner Catheter and the outer Guide (Figure 1, 2).

During use of the ET Catheter the catheter is placed inside the guide to assist during placement into the uterine cavity. The inner catheter stands proud from the guide for 5 cm.

Figure 1: Kitazato ET Catheter with Guide


Figure 2: Kitazato ET Catheter with Guide

A distal end of the guide is precurved 30° for an easy introduction into the uterine cavity. A rounded bulb tip of the guide is designed to lessen trauma during insertion into the canal of the cervix.

The guide is supplied with the adjustable guide stopper which provides a control of the depth of introduction during procedure with the help of the resin ink depth marks placed on the guide shaft. The inner catheter has soft distal end which is designed to lessen trauma during the transfer. A stainless stiffener added to the inside core of the catheter shaft provides more rigidity to its proximal part and facilitates the manipulations with the catheter.

The catheter shaft is supplied with the silicone ink depth marks for control of the depth of introduction. The catheter has a Luer Lock connector on its end which enables to connect 1 mL sterile single-use insulenic syringe for creating positive pressure and embryo transfer into the uterine cavity.

The range of Kitazato ET Catheters with guide comprises configurations of different length and outer diameter (Table 1).

Table 1: Parameters of Kitazato ET Catheters with Guide

Catheter Length	Outer Diameter Fr / mm	Guide Length	Stiffener added to the inside core	Catheter Echo Tip	Device Weight
220 mm	3 Fr/ 1.00 mm	170 mm	×	×	0.010 kg
250 mm	3 Fr/ 1.00 mm	200 mm	×	×	0.010 kg
220 mm	4 Fr/ 1.35 mm	170 mm	×	×	0.010 kg
250 mm	4 Fr/ 1.35 mm	200 mm	×	×	0.010 kg

Depth marks are placed on the catheter and guide shaft at 1 cm interval.
Echo tip is placed at 1 cm interval from the distal end of the catheter.

3.2. Kitazato ET Catheter with Guide and Obturator

Kitazato ET Catheter contains the inner Catheter, the outer Guide and Obturator (Figure 3).

During use of the ET Catheter the guide with the obturator placed inside is inserted into the uterine cavity. Obturator (intermediate rigidity) is designed to provide more rigidity to the guide while introducing it into the canal of cervix and avoids entrance of mucus in guide during placement. After reaching the uterine cavity, the obturator is removed from the guide and then catheter is placed inside the guide. The inner catheter stands proud from the guide for 5 cm.

Figure 3: Kitazato ET Catheter with Guide and Obturator

The distal end of the outer guide (Figure 4) is precurved 30° for an easy introduction into the uterine cavity.

The guide is supplied with the adjustable guide stopper which provides a control of the depth of introduction during procedure with the help of the resin ink depth marks placed on the guide shaft.


Figure 4: Kitazato ET Catheter: Guide with Obturator

The inner catheter (Figure 5) has a soft distal end which is designed to lessen trauma during the transfer.

The catheter has echo tip for for better visualization of the catheter during embryo transfer procedure under ultrasound control.

A stainless stiffener added to the inside core of the catheter shaft provides more rigidity to its proximal part and facilitates the manipulations with the catheter.

The catheter shaft is supplied with the silicone ink depth marks for control of the depth of introduction. The catheter has a Luer Lock connector on its end which enables to connect 1 mL sterile single-use insulenic syringe for creating positive pressure and embryo transfer into the uterine cavity.

Figure 5: Kitazato ET Catheter: inner catheter

The range of Kitazato ET Catheters with guide and obturator comprises configurations of different length and outer diameter (Table 2).

Table 2: Parameters of Kitazato ET Catheters with Guide and Obturator

Catheter Length	Outer Diameter Fr / mm	Guide Length	Obturator		Stiffener added to the inside core	Echo Tip	Device Weight
			Obturator Length	Color of Knob			
220 mm	3 Fr/ 1.00 mm	170 mm	170 mm	Pink	×	×	0.018 kg
250 mm	3 Fr/ 1.00 mm	200 mm	200 mm	Light Blue	×	×	0.018 kg
400 mm	3 Fr/ 1.00 mm	200 mm	200 mm	Pink	×	×	0.018 kg
220 mm	4 Fr/ 1.35 mm	170 mm	170 mm	Light Blue	×	×	0.018 kg
250 mm	4 Fr/ 1.35 mm	200 mm	200 mm	Light Blue	×	×	0.018 kg

Depth marks are placed on the catheter and guide shaft at 1 cm interval.
Echo tip is placed at 1 cm interval from the distal end of the catheter.

3.3. Kitazato ET Catheter with Memorable Guide

Kitazato ET Catheter contains the inner Catheter and outer Guide (Figure 6).
During use of the ET Catheter the catheter is placed inside the guide to assist during placement into the uterine cavity. The inner catheter stands proud from the guide for 5 cm.

Figure 6: Kitazato ET Catheter with Memorable Guide

The outer guide (Figure 7) keeps memory when giving shape and is able to take different curve level according to anatomic requirements of a patient. The distal end of the guide is supplied with the laser depth marks for control of the depth of introduction into the uterine cavity.


Figure 7: Kitazato ET Catheter: outer guide

The inner catheter (Figure 8) has a soft distal end which is designed to lessen trauma during the transfer.

A proximal part of the catheter is supplied with resin ink depth marks for control of the depth of introduction during the procedure.

The catheter has a Luer Lock connector on its end which enables to connect 1 mL sterile single-use insulenic syringe for creating positive pressure and embryo transfer into the uterine cavity.


Figure 8: Kitazato ET Catheter: inner catheter

The range of Kitazato ET Catheters with memorable guide comprises configurations of different length (Table 3).

Table 3: Parameters of Kitazato ET Catheter with Memorable Guide

Catheter Length	Outer Diameter Fr /mm	Guide Length	Device Weight
180 mm	4.7 Fr/ 1.56 mm	140 mm	0.005 kg
230 mm	4.7 Fr/ 1.56 mm	190 mm	0.005 kg

Depth marks are placed on the catheter and guide shaft at 1 cm interval.

3.4. Kitazato ET Catheter for Trial Embryo Transfer with Memorable Guide

Kitazato ET Catheter for trial embryo transfer contains inner Catheter with closed tip and outer Guide (Figure 9).

Karerep Kitazato ET Catheter for trial transfer is used for cervix examination and determination of its difficulty before the embryo transfer procedure.

During use of the ET Catheter the catheter is placed inside the guide to assist during placement into the uterine cavity. The inner catheter stands proud from the guide for 5 cm.

Figure 9: Kitazato ET Catheter for Trial Embryo Transfer with Memorable Guide

The outer guide (Figure 10) keeps memory when giving shape and is able to take different curve level according to anatomic requirements of a patient. The distal end of the guide is supplied with the laser depth marks for control of the depth of introduction into the uterine cavity.


Figure 10: Kitazato ET Catheter for Trial Embryo Transfer with Memorable Guide

The inner catheter (Figure 11) is supplied with soft rounded closed tip which lessens trauma during examination.

A proximal part of the catheter is supplied with the resin ink depth marks for control of the depth of introduction during the procedure.

Figure 11: Kitazato ET Catheter: inner catheter with blunt tip

The range of Kitazato ET Catheters for trial transfer with memorable guide comprises configurations of different length (Table 4).

Table 4: Parameters of Kitazato ET Catheter for Trial Embryo Transfer with Memorable Guide

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Depth marks are placed on the catheter and guide shaft at 1 cm interval.

3.5. Kitazato ET Catheter Stylet

Kitazato ET Catheter Stylet contains a plastic knob and hard shaft (Figure 12, 13).

Kitazato ET Catheter Stylet is used during the transfer procedure in case of difficult access to the cervix. Kitazato ET Catheter Stylet provides extra stiffness and memory to the guide catheter during insertion.

During use Kitazato ET Catheter Stylet is placed inside the guide. After reaching the uterine cavity the stylet is removed from the guide and then the catheter is placed inside the guide.

Figure 12: Kitazato ET Catheter Stylet

Figure 13: Kitazato ET Catheter Stylet

The range of Kitazato ET Catheter Stylets comprises configurations of different outer diameter and length according to the parameters of the catheter guide which is used during the transfer procedure (Table 5).

Table 5: Parameters of Kitazato ET Catheter Stylet

Catheter Stylet Length	Compatibility with Catheter		Device Weight
	Outer Diameter Fr / mm	Guide Length	
170 mm	3 Fr/ 1.00 mm	170 mm	0.003 kg
200 mm	3 Fr/ 1.00 mm	200 mm	0.003 kg
170 mm	4 Fr/ 1.35 mm	170 mm	0.003 kg
200 mm	4 Fr/ 1.35 mm	200 mm	0.003 kg
135 mm	4.7 Fr/ 1.56 mm	140 mm	0.003 kg
185 mm	4.7 Fr/ 1.56 mm	190 mm	0.003 kg

3.6. Kitazato IUI Catheter

Kitazato IUI Catheter contains Shaft and Connector (Figure 14).

The syringe is attached to the Kitazato IUI Catheter Connector end. During the intrauterine insemination procedure the washed spermatozoa are inhaled through tip of the shaft of the catheter, then the catheter is introduced into the uterine cavity through the cervix and the spermatozoa are infused by the syringe that is coupled at the connector end.

Kitazato IUI Catheter is manufactured with the three different types of rigidity offered: Hard, Intermediate and Soft.

Hard Catheters are supplied with stainless steel core inside their shaft to provide additional rigidity (Figure 15).

Intermediate and Soft Catheters are produced without the stainless steel core (Figure 16).

Figure 14: Kitazato IUI Catheter


Figure 15: Kitazato IUI Catheter: Hard


Figure 16: Kitazato IUI Catheter: Intermediate, Soft

Kitazato IUI Catheter has two side holes located bilaterally on the distal end for better sperm dispersion in the uterine cavity. Soft rounded tip of the catheter is designed to facilitate introduction of the catheter and lessens trauma during insertion.

The catheter shaft is supplied with the depth marks for control of the depth of introduction into the uterine cavity.

The catheter has a female Luer Lock connector on its end which enables to connect 3/5 mL sterile single-use insulinic syringe for creating positive pressure and sperm infusion into the uterine cavity.

The range of Kitazato IUI Catheter comprises catheters in two lengths: 10 and 18 cm. According to the rigidity level the catheter connector is presented in different colors (Table 6).

Table 6: Parameters of Kitazato IUI Catheters

Rigidity	Connector Color	Catheter Length	Outer Diameter Fr / mm	Total Volume minus hub	Stainless Steel Core	Device Weight
Hard	Dark Gray	100 mm	6 Fr/ 2.0 mm	0.015 ml	×	0.004 kg
Intermediate	Light Gray	100 mm	6 Fr/ 2.0 mm	0.05 ml	-	0.001 kg
Soft	White	100 mm	5.2 Fr/ 1.7 mm	0.05 ml	-	0.001 kg
Hard	Dark Gray	180 mm	6 Fr/ 2.0 mm	0.03 ml	×	0.007 kg
Intermediate	Light Gray	180 mm	6 Fr/ 2.0 mm	0.10 ml	-	0.004 kg
Soft	White	180 mm	5.2 Fr/ 1.7 mm	0.10 ml	-	0.004 kg

Depth marks are placed on the catheter shaft at 1 cm interval.

3.7. Kitazato IUI Catheter Stylet

Kitazato IUI Catheter Stylet contains plastic connector and the shaft made of stainless steel (Figure 17, 18). Kitazato IUI Catheter Stylet provides memory and extra stiffness in case it is difficult to insert the catheter to the uterine cavity.

Figure 17: Kitazato IUI Catheter Stylet

Figure 18: Kitazato IUI Catheter Stylet

The catheter has a female Luer Lock connector on its end which enables to connect 3/5 mL sterile single-use insulinic syringe.

During the procedure the catheter stylet is inserted into the IUI catheter. A syringe is coupled at the connector end of the catheter stylet. Catheter Stylet has an inside lumen. Without removing catheter stylet from the IUI catheter the spermatozoa are infused into the uterine cavity in one-step.

The range of Kitazato IUI Catheter Stylets comprises configurations of different length according to the parameters of Kitazato IUI Catheter which is used during the procedure (Table 7).

Table 7: Parameters of Kitazato IUI Catheter Stylet

Length	Outer Diameter Fr/mm	Compatibility with Catheter	Device Weight
		Catheter Length	
118 mm	2 Fr / 0.65 mm	100 mm	0.002 kg
198 mm	2 Fr / 0.65 mm	180 mm	0.005 kg

4. Composition and Raw Materials

For manufacture of Kitazato ET Catheter with and compatible Kitazato ET Catheter Stylets the following raw materials are used:

Name of Component		Material
Catheter	Shaft	Silicone –180705-200
	Connector	Nylon 12
	Center Core	Stainless Steel SUS 304 Composition: Carbon (C) – up to 0.8% Chrome (Cr) – 17-19% Manganese (Mn) – up to 0.2% Silicon (Si) – up to 0.8% Phosphorus (P) – up to 0.035% Sulfur (S) – up to 0.02% Nickel (Ni) – 9-11% Titanium (Ti) – up to 0.5% Cuprum (Cu) – up to 0.3% Ferrum (Fe) – ~69%
	Echo Tip	Stainless Steel SUS 304 Composition: Carbon (C) – up to 0.8% Chrome (Cr) – 17-19% Manganese (Mn) – up to 0.2% Silicon (Si) – up to 0.8% Phosphorus (P) – up to 0.035% Sulfur (S) – up to 0.02% Nickel (Ni) – 9-11% Titanium (Ti) – up to 0.5% Cuprum (Cu) – up to 0.3% Ferrum (Fe) – ~69%
	Depth Marks	Silicone Ink - Transtech B 1 black (CPN 14994)
Guide	Shaft	Fluoric resin
	Connector	Polycarbonate
	Stopper	Styrene elastomer
	Depth Marks	Laser marking
Obturator	Knob	ABS - M205FC
	Shaft	Nylon 12
Catheter Stylet	Knob	ABS - M205FC
	Shaft	Nylon 12
Adhesive Material		Epoxy resin two-component adhesive bond Loctite 3311

For manufacture of Kitazato ET Catheter with memorable guide and Kitazato ET Catheter for trial transfer with guide as well as the compatible Kitazato ET Catheter Stylets the following raw materials are used:

Name of Component		Material
Catheter	Shaft	Silicone 180705-200
	Connector	Nylon 12
	Depth Marks	Resin Ink
Guide	Shaft	Fluoric resin
	Connector	ABS - M205FC
	Depth Marks	Laser marking
Catheter Stylet	Knob	ABS - M205FC
	Shaft	Nylon 12
Adhesive Material		Epoxy resin two-component adhesion bond Loctite 3311

For manufacture of Kitazato IUI Catheter and the compatible Kitazato IUI Catheter Stylets the following raw materials are used:

Name of Component		Material
Catheter	Shaft	Silicone 180705-200
	Connector	Nylon 12
	Center Core	Stainless Steel SUS 304 Composition: Carbon (C) – up to 0.8% Chrome (Cr) – 17-19% Manganese (Mn) – up to 0.2% Silicon (Si) – up to 0.8% Phosphorus (P) – up to 0.035% Sulfur (S) – up to 0.02% Nickel (Ni) – 9-11% Titanium (Ti) – up to 0.5% Cuprum (Cu) – up to 0.3% Ferrum (Fe) – ~69%
	Depth Marks	Ink with blue Composition: Cyclohexanone (40-50%) Aromatic mixture hydrocarbon (5-20%) Synthetic resin (10-40%) Pigment (0-50%)
Catheter Stylet	Connector	ABS - M205FC
	Shaft	Stainless Steel SUS 304 Composition: Carbon (C) – up to 0.8%

		Chrome (Cr) – 17-19% Manganese (Mn) – up to 0.2% Silicon (Si) – up to 0.8% Phosphorus (P) – up to 0.035% Sulfur (S) – up to 0.02% Nickel (Ni) – 9-11% Titanium (Ti) – up to 0.5% Cuprum (Cu) – up to 0.3% Ferrum (Fe) – ~69%
Adhesive Material		Epoxy resin two-component adhesion bond Loctite 3311

5. Instruction for Use

5.1. Indications for Use

- Infertility of the unknown etiology;
- Failure of conception due to tubal blockage (for intrauterine insemination a necessary condition is tubal patency of at least one tube);
- Cervical factor infertility;
- Male factor infertility;
- Combined infertility;
- Polycystic ovarian syndrome after failure of conception during a year of the conservative treatment;
- Endometriosis.

5.2. Contraindications

Sterile single use Kitazato Catheters for assisted reproductive technologies should not be used under acute genital infection.

Sterile single use Kitazato Catheters for assisted reproductive technologies should not be reused. It may cause infection when reuse.

Sterile single use Kitazato Catheters for assisted reproductive technologies should not be re-sterilized.

5.3. Potential Side Effects

Sterile single use Kitazato Catheters for assisted reproductive technologies are produced from soft materials such as silicone and nylon and under use of qualified operator in compliance with aseptic technique cannot injure endometrium and vaginal mucous membrane, cause bleeding or infection carry.

Sterile single use Kitazato Catheters for assisted reproductive technologies do not have potential side effects.

5.4. Preparation for Use and Use Directions

- (1) Position the patient in the lithotomy position. Insert the vaginal speculum to expose cervix. Check the condition of cervix and decide the insert length and angle of the catheter and the guide before use.
Caution: Use a catheter stylet guide in case the cervix bends and it seems to be difficult to pass the cervix.
- (2) Wash the vulva and the vagina with sterilized saline and dry them.
- (3) Connect a syringe to the connector of the catheter and aspirate embryo or load catheter the catheter with washed spermatozoa.
- (4) Insert the catheter to the cervix.
- (5) Check the catheter reached to the uterus under ultrasound. Slowly press the plunger of the syringe and expel the embryo/spermatozoa into the uterine cavity.
- (6) Carefully remove the catheter.
- (7) In case of embryo transfer check for retention of embryos in the catheter.
- (8) Take the vaginal speculum out and keep the patient at rest for further 40-50 minutes.
- (9) Discard the syringe and catheter as medical waste.

5.5. Cautions for Use

- (1) Before use please familiarize yourself with the information on the label attached to the product or package and safety instructions.
- (2) Confirm the specification (outer diameter and length) of the catheter before opening it.
- (3) Confirm the adaptability with the other medical devices and equipment -to use together beforehand.
- (4) Immediately use after opening the package.
- (5) If the product has abnormalities that its package has been broken or it has been damaged before opening it, do not use it.
- (6) Do not wash the inside of the catheter with sterilized saline before infusing the catheter.
- (7) Do not use the product that it was past an expiration date.
- (8) While introducing the catheter through the cervix there is a risk of bleeding due to cervical trauma. Do not jack up forcibly or reprocess.
- (9) Use the product by operator that has seasoned knowledge about procedure for the product.
- (10) Do not use the product except the transfer of embryo or spermatozoa.
- (11) Use the product aseptically. Carry out enough precautionary measures for bacteria.
- (12) Throw the product away promptly after having used it.
- (13) Do not refurbish.
- (14) Do not reuse.
- (15) Store away from high temperature, direct sunlight, high humidity and the water wet.
- (16) The product made from polyvinyl chloride has the availability that a plasticizer may ~~liquefy out~~ .

6. Labelling and Packaging

4.1. Primary packaging

Kitazato ET Catheter and Kitazato IUI Catheter are all packaged in an individual sterile single pouch barrier system (primary packaging) made of laminated PET, polyethylene and tyvek (Figure 19, 20).


Figure 19: Individual (primary) package of Kitazato ET Catheter: front and black sides

Figure 20: Individual (primary) package of Kitazato IUI Catheter: front and back sides

Kitazato ET Catheter with guide and obturator is packaged in double pouch primary packaging: catheter and the guide with obturator are packaged separately inside the primary packaging. Inside the sterile pouch the catheter shaft is placed into protective sleeve made of polypropylene for preserving from mechanical damage during transportation and storage. (Figure 21).

6. Labelling and Packaging

4.1. Primary packaging

Kitazato ET Catheter and Kitazato IUI Catheter are all packaged in an individual sterile single pouch barrier system (primary packaging) made of laminated PET, polyethylene and tyvek (Figure 19, 20).


Figure 19: Individual (primary) package of Kitazato ET Catheter: front and black sides

Figure 20: Individual (primary) package of Kitazato IUI Catheter: front and back sides

Kitazato ET Catheter with guide and obturator is packaged in double pouch primary packaging: catheter and the guide with obturator are packaged separately inside the primary packaging. Inside the sterile pouch the catheter shaft is placed into protective sleeve made of polypropylene for preserving from mechanical damage during transportation and storage. (Figure 21).


Figure 21: Double pouch individual (primary) package of Kitazato ET Catheter: front and back sides

Kitazato ET Catheter with memorable guide, Kitazato ET Catheter for trial transfer with memorable guide, Kitazato ET Catheter Stylets and Kitazato IUI Catheter Stylets are also supplied with the protective sleeve made of polypropylene for preserving from mechanical damage during transportation and storage.

The individual single pouch (primary packaging) of Kitazato ET Catheter and Kitazato IUI Catheter has the following dimensions (Table 8).

Table 8. Dimentions of single pouch (primary packaging)

Catheter Type	Vertically Long (mm)	Horizontally Long (mm)
Kitazato ET Catheter with Guide	100	525
Kitazato ET Catheter with Guide and Obturator	100	525
- Catheter	60	485
- Guide with Obturator	60	300
Kitazato ET Catheter with Memorable Guide	60	485
Kitazato ET Catheter for trial transfer with Memorable Guide		
Kitazato ET Catheter Stylet	60	300
Kitazato IUI Catheter	60	300
Kitazato IUI Catheter Stylet		

Gross weight of devices in single pouch (primary packaging) is shown in Table 9.

Table 9. Device Weight (Single Pouch)

Catheter Type	Gross Weight (kg)
Kitazato ET Catheter with Guide	0.019
Kitazato ET Catheter with Guide and Obturator	0.027
Kitazato ET Catheter with Memorable Guide	0.008
Kitazato ET Catheter for trial transfer with Memorable Guide	0.008
Kitazato ET Catheter Stylet	0.006
Kitazato IUI Catheter	0.007
- Hard 100 cm/180 cm	0.004/0.007
- Intermediate 100 cm/180 cm	0.001/0.004
- Soft 100 cm/180 cm	0.001/0.004
Kitazato IUI Catheter Stylet 118 cm/198 cm	0.005/0.008

Each individual package of Kitazato catheters is supplied with a product label containing product name, reference, lot, expiry date, name and address of the manufacture and EU representative as well as gamma or EOG irradiation sterilization method applied and single use product indication (Figure 22, 23).

KITAZATO® Kitazato ET Catheter

Embryo Transfer Catheter / Catetere per Transferi Embrionaria /
Catheter de Transfert Embryonnaire / Embryotransferkatheter /
Καθετήρας εμβρυομεταφοράς / Katerer do Transferu Zarodków /
Катетер для Переноса Эмбриона / Cateter de Transferencia Embrionaria

REF XXXXXX

[Technical Characteristics]

(01) 0 4580303 26609 0 (17) 000000 (10) M00000

LOT	M00000	0000.00
②	STERILE	CE 0086
EC REP	Obimed-Biomedical Supply, S.L. C/ Los Berrués, 1 - Oficina Pta 4 46010 Valencia Spain TEL +34-96-305-63-95 FAX +34-96-305-63-96 KITAZATO Medical Co., Ltd. 1-1-8 Shibadaimon Minato-ku Tokyo 105-6012 Japan TEL +81-3-3434-2731 FAX +81-3-3434-2732	

11 00001

KITAZATO® Kitazato ET Catheter

Embryo Transfer Catheter / Catetere per Transferi Embrionaria /
Catheter de Transfert Embryonnaire / Embryotransferkatheter /
Καθετήρας εμβρυομεταφοράς / Katerer do Transferu Zarodków /
Катетер для Переноса Эмбриона / Cateter de Transferencia Embrionaria

REF XXXXXX

[Technical Characteristics]

(01) 0 4580303 26609 0 (17) 000000 (10) M00000

LOT	M00000	0000.00
②	STERILE	CE 0086
EC REP	Obimed-Biomedical Supply, S.L. C/ Los Berrués, 1 - Oficina Pta 4 46010 Valencia Spain TEL +34-96-305-63-95 FAX +34-96-305-63-96 KITAZATO Medical Co., Ltd. 1-1-8 Shibadaimon Minato-ku Tokyo 105-6012 Japan TEL +81-3-3434-2731 FAX +81-3-3434-2732	

11 00001

Figure 22: Kitazato ET Catheter primary packaging label

KITAZATO® Kitazato IUI Catheter

Intrauterine Insemination Catheter / Catetere per Inseminazione Intrauterina /
Catheter d'Insemination Intra uterine / Intrauterin Inseminationskatheter /
Καθετήρας ενδομήτριας σπερματεμεύσης / Katerer Inseminacyjny /
Катетер для Внутриматочной Инсеминации / Cateter Inseminacion Intrauterina

REF XXXXXX

[Technical Characteristics]

(01) 0 4580303 26607 6 (17) 000000 (10) M00000

LOT	M00000	0000.00
②	STERILE	CE 0086
EC REP	Obimed-Biomedical Supply, S.L. C/ Los Berrués, 1 - Oficina Pta 4 46010 Valencia Spain TEL +34-96-305-63-95 FAX +34-96-305-63-96 KITAZATO Medical Co., Ltd. 1-1-8 Shibadaimon Minato-ku Tokyo 105-6012 Japan TEL +81-3-3434-2731 FAX +81-3-3434-2732	

KJ 00001

Figure 23: Kitazato IUI Catheter primary packaging label

Symbols on labels have the following meaning (Figure 24):

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Figure 24: Symbols on the label

4.2. Secondary packaging

Kitazato ET Catheter and Kitazato ET Catheter Stylet in the primary packaging are sold by a cartoon boxes of unit of 10. Kitazato IUI Catheter and Kitazato IUI Catheter Stylet in the primary packaging are sold by a cartoon boxes of unit of 50.

Secondary packaging is designed to protect the units during shipping and storage.

Dimentions of transport (secondary) packaging are shown in Table 10.

Table 10. Dimentions of secondary packaging

Catheter Type	Vertically Long (mm)	Horizontally Long (mm)	Height (mm)
Kitazato ET Catheter with Guide	105	530	85
Kitazato ET Catheter with Guide and Obturator			
Kitazato ET Catheter with Memorable Guide	50	490	70
Kitazato ET Catheter for trial transfer with Memorable guide			
Kitazato ET Catheter Stylet	50	350	65
Kitazato IUI Catheter	190	310	90
Kitazato IUI Catheter Stylet			

Secondary packaging is also supplied with a label containing product name, reference, lot, expiry date, name and address of the manufacture and EU representative as well as sterilization method applied, single use product indication and number of units in the box (Figure 25, 26). Description of symbols on the label is shown in Figure 24.

KITAZATO® Kitazato ET Catheter <i>Embryo Transfer Catheter / Catetere per Transfer Embrionale / Catheter de Transfert Embryonnaire / Embryotransferkatheter / Катетер для переноса эмбрионов / Kateter do Transferu Zarodków / Катетер для Переноса Эмбрионов / Cateter de Transferencia Embrionaria</i> REF XXXXXX 10pieces/Box <i>[Technical Characteristics]</i>  (01) 0 4580303 26609 0 (17) 000000 (10) M00000		<table border="1"> <tr> <td>LOT M00000</td> <td> 0000.00</td> </tr> <tr> <td></td> <td> STERILE</td> </tr> <tr> <td>EC REP</td> <td></td> </tr> <tr> <td> Orimed-Biomedical Supply, S.L. C/ Los Buitos, 1 Oficinas Puñet 46015 Valencia Spain TFL +34-96-305-83-95 FAX +34-96-305-83-98 </td> <td> KITAZATO Medical Co., Ltd. 1-1-8 Shibadaimon Minato-ku Tokyo 105-0072 Japan TEL +81-3-3434-2731 FAX +81-3-3434-2732 </td> </tr> </table>		LOT M00000	 0000.00		 STERILE	EC REP		Orimed-Biomedical Supply, S.L. C/ Los Buitos, 1 Oficinas Puñet 46015 Valencia Spain TFL +34-96-305-83-95 FAX +34-96-305-83-98	KITAZATO Medical Co., Ltd. 1-1-8 Shibadaimon Minato-ku Tokyo 105-0072 Japan TEL +81-3-3434-2731 FAX +81-3-3434-2732
LOT M00000	 0000.00										
	 STERILE										
EC REP											
Orimed-Biomedical Supply, S.L. C/ Los Buitos, 1 Oficinas Puñet 46015 Valencia Spain TFL +34-96-305-83-95 FAX +34-96-305-83-98	KITAZATO Medical Co., Ltd. 1-1-8 Shibadaimon Minato-ku Tokyo 105-0072 Japan TEL +81-3-3434-2731 FAX +81-3-3434-2732										

Figure 25: Kitazato ET Catheter secondary packaging Label


Figure 26: Kitazato IUI Catheter secondary packaging Label

According to the requirements of national legislation in terms of labelling of medical devices an authorized representative places the labels in Russian language on primary and secondary packaging (Figure 27-34).

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Figure 28: Kitazato ET Catheter Stylet and Kitazato ET Catheter for trial transfer primary packaging Label Draft

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Figure 29: Kitazato ET Catheter secondary packaging Label Draft

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Figure 30: Kitazato ET Catheter Stylet and Kitazato ET Catheter for trial transfer secondary packaging Label Draft

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Figure 31: Kitazato IUI Catheter primary packaging Label Draft

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Figure 32: Kitazato IUI Catheter Stylet primary packaging Label Draft

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Figure 33: Kitazato IUI Catheter secondary packaging Label Draft

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Figure 34: Kitazato IUI Catheter Stylet secondary packaging Label Draft

7. Sterilization

Kitazato ET Catheter and Kitazato IUI Catheter are sterilized by radiation method (gamma irradiation) meeting standard ISO 11137.

The Sterility Assurance Level SAL 10^{-6} is adequately assured at the sterilization dose 25.0 kGy.

Kitazato ET Catheter with memorable guide, Kitazato ET Catheter for trial transfer with memorable guide and compatible Kitazato ET Catheter Stylet are sterilized by EOG irradiation meeting standard ISO 11138.

Sterilization validation reports are attached as the separate documents.

Because high performance sterilization bag is used, 5 years shelf life is adequate. Actual shelf life of 3 years is settled considering the safety of the device.

8. Quality Control

Throughout the manufacturing process of Kitazato catheters the quality control is exercised from the input of raw material to obtaining of final products.

The final products are involved into complex check including the following tests:

A- Biological Tests - device

Title of test	Acceptance criteria	Standards	Results
1- Sterility test (SAL)	10^{-6}	USP(71)	Confirm
Ethylene oxide sterilization residuals ¹	EO, ECH	ISO10993-7	Confirm
2- Mouse Embryo Assay (MEA)	One-cell assay ≥ 80 after 96 hours	ISO 17025	Confirm
3- Endotoxin	< 20 EU/device	USP33(85)	Confirm
4- Human Sperm Survival Assay (HSSA) ²	$\geq 70\%$ motility at 24 hours	ISO 17025	Confirm
5- Cytotoxicity Test	-	ISO10993-5	Confirm
6- Intracutaneous Reactivity Test	-	ISO10993-10	Confirm
7- Sensitization Test	-	ISO10993-10	Confirm

¹ Only for Kitazato ET Catheter with memorable guide, Kitazato ET Catheter for trial transfer with memorable guide and compatible Kitazato IUI Catheter Stylets.

² Only for Kitazato IUI Catheter and compatible Kitazato IUI Catheter Stylets.

B- Physical Tests

Title of test	Acceptance criteria	Standards	Results
1.Performance Test	-	JIS T 3245:2011	Confirm

9. Contents of Delivery

Kitazato ET Catheter and Kitazato IUI Catheter depending on modification are supplied in sterile single/double pouch barrier system for the packaging. The double pouch barrier system is more for convenience to have the inner Kitazato ET Catheter in a single package for the physician performing a procedure.

Kitazato ET Catheter Stylet and Kitazato IUI Catheter Stylet are supplied in sterile single pouch barrier system for the packaging.

Kitazato ET Catheter and Kitazato IUI Catheter in individual packages are sold by a box (secondary packaging) in unit of 10 / 50 accompanied by Instruction for Use.

10. Transportation and storage conditions

Kitazato ET Catheter and Kitazato IUI Catheter are shipped by all means of transport in roofed transport facilities. During transportation a type of a transport facility and distance as well as carrier and vibration frequency are considered.

During transportation and storage of Kitazato ET Catheter and Kitazato IUI Catheter it is necessary to avoid high temperature, high humidity, direct sunlight and water wet. Transportation and storage of Kitazato ET Catheter and Kitazato IUI Catheter is acceptable at + 15 - + 30°C, 65%RH.

11. Shelf Life

The shelf life of the device is 3 years.

Do not use the device after expiration date indicated on the package.


12. Disposal

Used Kitazato ET Catheter and Kitazato IUI Catheter should be disposed in compliance with the requirements of national legislation in regard to disposal and destruction of class B medical waste. For collection of used Needles single use water resistant container is used. The container should have adherent cap, which avoids inadvertent opening of the container.

13. Manufacture. Warranty

Kitazato ET Catheter and Kitazato IUI Catheter are manufactured by KITAZATO Medical Co., Ltd., Shibakoen Building, 1-1-8 Shibadaimon, Minato-ku, Tokyo, Japan.

The quality of the medical devices and the production site of KITAZATO Medical Co., Ltd. meet the international quality standards as proved by ISO 13485:2003 & EN ISO 13485:2012 and declaration of conformity for medical device.

The manufacture guarantees the quality of device during its shelf life – 3 years from sterilization.
The manufacture is not liable for any damage resulting from misuse or incorrect operation or failure to follow the cautions or instructions contained in the instruction for use or the instruction manual for any compatible medical equipment.

14. Reclamation


In case of claims concerning the quality of medical device, please contact Manufacture and / or its Authorized Representative.

Authorized Representative in Russian Federation:

EKO-zdorovie Ltd.
12, building 6, Osenniy boulevard, Moscow, 121614 Russia
Tel.: +7 (499) 391-26-70
E-mail: info@ekozdrav.ru

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