75
Years

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Thank you for your expression of confidence in the KARL STORZ brand name. Like all of our other products, this medical device is the result of years of experience and great care in manufacture. You and your organization have decided in favor of a modern, high-quality piece of equipment from KARL STORZ.

This instruction manual is intended to aid you in the proper installation, connection, and operation of this medical device.

- Failure to observe and follow this instruction manual and the instruction manuals of medical devices used in combination can result in injury to patients, users and third parties as well as damage to the medical device.
- Please read all relevant instruction manuals and always follow the instructions given. Check the functioning of the medical devices used in combination.
- Take note of all the instructions on the device.
- Keep this instruction manual available for ready reference in a clearly visible location near the device.

KARL STORZ reserves the right to make engineering modifications, in particular in the context of ongoing development and improvement of our devices.

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It is recommended that the suitability of the products for the intended procedure be checked prior to use.

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
1	Images of the unit.....3
2	Symbols employed.....5
2.1	Safety instructions.....5
2.2	Pictograms.....6
3	Intended use.....7
3.1	Intended use.....7
3.2	Contraindications.....7
3.3	User qualification.....7
4	Security.....8
5	Description of the device.....11
6	Initial operation.....14
6.1	Unpacking.....14
6.2	Scope of supply.....14
6.3	Required accessories.....14
6.4	Safety precautions at the site of installation.....15
7	Assembly.....16
7.1	Attaching the base to the OR table.....16
7.2	Securing the ARTIP BASE base unit to the base.....16
7.3	Attaching the sterile cover and assembly of the holding arm.....17
7.4	Securing the instruments.....19
7.5	Status displays.....20
8	Operation.....21
8.1	Releasing and moving the holding arm.....21
8.2	Fixing the joints of the holding arm.....22
9	Disassembly of the ARTIP BASE holding system.....23
10	Reprocessing.....24
10.1	Reprocessing overview.....24
10.2	General warnings.....25
10.3	Reprocessing of ARTIP BASE holding arms.....26

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
10.3.1	Manual wipe-down disinfection	26
10.3.2	Accessories	26
10.3.3	Preparation for cleaning and disinfection	26
10.3.4	Manual precleaning	27
10.3.5	Machine cleaning and disinfection	27
10.3.6	Assembly, inspection and care	27
10.3.7	Packaging systems	28
10.3.8	Sterilization	28
10.4	Limits of reprocessing	28
10.5	Reprocessing of the drive unit 28272 TM, rechargeable battery 28272 TA	29
10.6	Reprocessing of the charging station	29
11	Service, repair and disposal	30
11.1	Guarantee	30
11.2	Maintenance	30
11.3	Safety inspection/repair inspection according to IEC 62353	30
11.4	Repairing device parts	31
11.5	Repair program	31
11.6	Infection prevention	31
11.7	Disposal	32
11.8	Application/storage conditions	32
12	Directive compliance	32
13	Electromagnetic Compatibility user information	33
14	Subsidiaries	51

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
- ① Base column
- ② Locking screw for base column 28272 TM
- ③ Base joint
- ④ Base unit
- ⑤ Rechargeable battery
- ⑥ Illuminated ring status display
- ⑦ Actuation ring
- ⑧ Connecting cone
- ⑨ Collar for cover
- ⑩ Holding arm
- ⑪ Central joint
- ⑫ Hand joint/ball joint
- ⑬ Handle with touch plates
- ⑭ Signal cable with cable clips
- ⑮ Screw for securing the KSLOCK clamping jaws
- ⑯ KSLOCK mount, female
- ⑰ KSLOCK clamping jaw
- ⑱ Charger

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HAZARD SYMBOL: This is a hazard symbol. It warns of the risk of injury. To prevent injuries, follow all the measures marked with the hazard symbol.

Signal words:

WARNING! designates a possible imminent risk. If this is not avoided, it **could** lead to death or serious injuries.

CAUTION! designates a possible imminent risk. If this is not avoided, it could lead to minor injuries.

ATTENTION! designates a possibly harmful situation. If this is not avoided, the medical devices could be damaged.

NOTE! indicates special information on the operation of the device or important information.

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2.2 Pictograms

Symbol	Explanation
	Follow instructions for use
	This device has been marked in accordance with the European Directive on waste electrical and electronic equipment (WEEE).
	Catalogue number
	Number of products in the product packaging
	CE mark
	IP22 (splash-water protection)
	Federal (USA) law restricts this device to sale by or on the order of a physician.

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3 Intended use

3.1 Intended use

The holding system is used for the short-term, flexible positioning of KARL STORZ endoscopes, VITOM® and instruments for open, microsurgical and minimally invasive interventions.
The holding system does not come into contact with the patient and is therefore non-invasive.

3.2 Contraindications

No contraindications relating directly to the product are currently known. The use of the holding system is contraindicated if, in the opinion of the attending physician, the surgical method as such is contraindicated or if the patient is not able to undergo surgery or anesthesia due to his or her general condition.

3.3 User qualification

Only physicians and medical support staff with a relevant specialist qualification and who have received training on the device may use the ARTIP BASE holding system. Only technically suitable persons may provide training.

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4 Security

Risks of infection

WARNING: These products are not sterile when delivered. The use of non-sterile products poses a risk of infection for patients, users and third parties.

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WARNING: If CJD (Creutzfeldt Jakob disease) is suspected or indicated, KARL STORZ recommends disposing of the instrument in accordance with national regulations and not using it again.

Combination of system components

WARNING: Risk of injury and damage to the product: Failure to observe and follow this instruction manual and the instruction manuals of products used in combination can result in injury to patients, users and third parties as well as damage to the product.

- Please read all relevant manuals and always follow the instructions given precisely. Check the functioning of the products used in combination.

WARNING: Combinations of medical devices are only assured to be safe if

- they are identified as such in the respective manuals or
- the intended use and interface specifications of the products used in combination permit this (cf. IEC 60601-1-1 or chap. 16, IEC 60601-1).

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WARNING: Additional devices that are connected to electrical medical devices must conform to the respective IEC or ISO standards (e.g., IEC 60950 for data processing equipment). Furthermore, all configurations must satisfy the requirements for medical electrical systems (see IEC 60601-1-1 and IEC 60601-1, paragraph 16, 3rd edition). Any person who connects additional devices to electrical medical devices configures a medical system, and is therefore responsible for ensuring that this system complies with the requirements for medical electrical systems. It must be taken into account that local regulations take priority over the aforementioned requirements. If in doubt, please contact your local representatives or technical support.

Improper use

WARNING: Risk of injury: Incorrect application of medical instruments poses a risk of injury for patients.

- Users of medical instruments must have an appropriate medical qualification and be acquainted with the application.

WARNING: Risk of injury: Overloading the instrument by exerting too much force may cause the medical device to break, bend or malfunction, and consequently injure the patient or user.

- Do not overload the medical devices. Do not bend bent medical devices back to their original position.

ATTENTION: Liquid which penetrates the housing can damage the device.

- Do not store liquids on or above the device.
- If liquid penetrates the device, in spite of the precautions, sufficient time should be allowed for evaporation (this also applies for the formation of condensation).

WARNING: No modifications to the ARTIP BASE are allowed.

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Risks due to damaged parts

WARNING: Incorrectly assembled and damaged medical devices can lead to injuries to the patient.

- Medical devices and all the accessories used in combination with them must be checked immediately before and after every use to ensure that they are complete, free from damage, and in full working order and have no unintentional rough surfaces, sharp corners, burred edges or projecting parts.

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5 Description of the device

The KARL STORZ holding system ARTIP BASE is a holding system which can be taken apart and which comprises a holding arm, a base unit and a rechargeable battery. The system is characterized by the fact that it can be operated with one hand thanks to the electromechanical locking.

The ARTIP BASE holding system is used for the static, free positioning of telescopes, instruments and other aids in the operating field.

The maximum holding force with a horizontally extended holding arm is 3 kg for the 28272 TS and 2 kg for the 28272 TL whereby some joints yield more easily than others. This position, however, is very unlikely in real life; the extension is considerably shorter and the holding force thus greater.

The clamping of all joints is based on friction. Thus, if overloaded, they do not give abruptly but slowly slip down. As such, the device is not damaged which is why, in the event of an error, the user can remove the holding system using force from the operating field.

Components

A ready-for-use ARTIP BASE holding system comprises a base unit, a holding arm and a battery which includes an external charger.

• Base unit

The base unit must be covered with a sterile cover.

The base unit comprises a motor, gears and the complete electric control system for the ARTIP BASE holding system and creates the clamping force within fractions of a second.

The connecting cone of the holding arm is mounted in the base unit.

The electrical signals from the handle are transmitted to the electric control system via the coupling and the clamping force of the gears is transmitted to the holding arm joints via clamping bolts. The coupling itself acts a joint around its own main axis and thus provides one of the necessary degrees of freedom.

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The base column can be attached to the operating rail (DIN EN ISO 19054) of an operating table using a standard clamp (not included with the ARTIP BASE holding system).

The base unit signalizes operating states and error states via a multi-colored illuminated ring. There is no main switch; the system is immediately ready for use once the battery is inserted.

• Holding arm

The holding arms can be completely steam sterilized and can be reprocessed before each application and used directly in the operating field without a sterile cover. The distal end of the holding arm is formed by the handle with a KSLOCK mount to which numerous other articles from the KARL STORZ portfolio can be attached using the relevant connection system.

Thanks to its ergonomic triangular shape, the handle itself is suitable for holding and positioning high working loads conveniently and safely. There are three touch plates on the three outer surfaces for releasing and realigning the position. The lock is only released if the operator holds down all three touch plates. When one of the touch plates is released again, the holding system locks.

The handle is connected to the lower arm via a ball joint which enables full rotational freedom and almost right-angled bending of the handle.

The central joint connects the upper and lower arm like a hinge and restricts their reciprocal movement to a full revolution in order to prevent damage to the outer signal cable.

The upper arm ends proximally in the coupling cone which creates the mechanical and electrical connection to the base unit.

The complete system has a total of seven degrees of freedom.

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• **Rechargeable battery**

The energy for the base unit is provided by a lithium-ion battery which is inserted in the battery compartment. The battery has four 4S1P cells with 14.4 V nominal voltage and 2.5 Ah capacity. This enables more than 1,000 releasing and locking procedures.

• **Charger**

The associated battery charger is not a medical device and is intended for operation outside of the operating theater (also see instruction manual 97000250).

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6 Initial operation

6.1 Unpacking

Carefully remove the ARTIP BASE holding system and accessories from the packaging. Check for missing items and for evidence of shipping damage.

If necessary, notify the manufacturer or the supplier.

TIP: If possible, retain the original packing materials for later use; these can be used if the device has to be transported.

6.2 Scope of supply

The ARTIP BASE holding system set 28272 TLS or 28272 TSS comprising

- Base unit 28272 TM
- Locking screw for base column
- Long holding arm 28272 TL (in set 28272 TLS) or short holding arm 28272 TS (in set 28272 TSS)
- 2 x battery 28272 TA
- Charging station 28272 TC
- Sterile cover, holding systems, 041151-30, pack of 30

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
6.4 Safety precautions at the site of installation

The holding system may only be used in medical rooms whose electrical systems have been installed in accordance with applicable national regulations.

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⚠ WARNING!

The device is not intended for use in hazardous zones. Do not operate the device within the demarcated hazard zone while explosive anesthetic gases are in use.

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⚠ WARNING!

The device must not be used in oxygenated environments.

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7 Assembly

NOTE: The following components required for operation of the ARTIP BASE holding system

- Battery (Art. no. 97000244)
- Battery charger (Art. no. 97000250)
- Sterile cover (Art. no. 97000275)

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7.1 Attaching the base to the OR table

1. Fix base, e.g., 28172 HR, to an operating rail on the OR table.
2. Insert and secure the base unit column in the base.

⚠ WARNING!

Risk of injury if the ARTIP BASE holding system is not stable. The stability of the ARTIP BASE holding system is dependent on the stability of the operating table's side rails.

- Check the operating table thoroughly. If the stability of the side rails is questionable, the ARTIP BASE holding system must not be used.
- 3. Select the position of the holding system on the operating rail according to the areas of the operating field which the range of the ARTIP BASE holding system should cover during subsequent operation.

NOTE: After covering with sterile covers and commencing the operation, the attachment point of the holding arm on the OR table should no longer be altered.

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
7.2 Securing the ARTIP BASE base unit to the base

1. Remove the screw ② at the end of the base column ①.
2. After opening the base on the OR table, the base column ① of the base unit can be pushed in the base and fixed in place there.

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3. Screw the previously removed screw ② back into the end of the base column ①. This prevents the base column from being pulled out the base unintentionally.
4. Remove a battery from the charger and check the battery level.
5. Insert the fully charged battery ⑤ in the correct position into the base unit until it audibly clicks into place.

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7.3 Attaching the sterile cover and assembly of the holding arm

⚠️ WARNING!

The base unit cannot be sterilized. It must be covered with a sterile cover.

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⚠️ WARNING!

Risk of infection: Sterility has only been validated with the ARTIP BASE holding system when using the KARL STORZ cover Art. no. 041151-30.

- Do not use any other covers.
- 1. Open the sterile packaging. Pull the sheet off the paper from the top right- and left-hand corner of the packaging using your index finger.

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
2. Grasp a large area of the middle of the sheet and paper and carefully rip open the seal using both hands.
3. Remove the product from the sterile packaging.

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
4. Insert the ARTIP BASE holding arm (10) between the two ends of the cover marked with arrows and push towards the preperforated opening. This is located in one corner of the cover.

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
5. Push the connecting cone (8) of the ARTIP BASE holding arm (10) through the preperforated opening of the cover.

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
6. Push the connecting cone through the opening until the connecting cone and the black collar (9) (see arrow) project out of the cover.

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7. Insert the connecting cone (8) of the holding arm which is not covered with the sheet into the ARTIP BASE base unit (4) until it audibly clicks into place.

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⚠ WARNING!

Risk of injury if the ARTIP BASE holding system is not stable.

- Ensure that the connection between the holding arm and the base unit is sealed correctly.

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8. Pull the cover over the base unit ④ and base.
ATTENTION: Excessive pulling of the cover could cause it to tear.

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9. Secure the cover in place with the adhesive strip provided.

⚠️ WARNING!

Risk of injury to the patient. If the sterile cover is too tight, this may restrict the movement of the holding system and thus the use of the instruments.


➤ Check to ensure freedom of movement before securing the cover.

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10. In order to lock the handle, simultaneously press all 3 touch plates on the handle ⑬ once.
The lock is only activated after actuating the touch plates on the handle.
This also acts as an implicit test for functionality.

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7.4 Securing the instruments

⚠️ WARNING!

Risk of injury to the patient: If the permissible working load (max. load) printed on the holding arm is exceeded, the safe functioning of the system can no longer be guaranteed.
➤ Never exceed the max. load.

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

1. Insert the correct KSLOCK clamping jaw ⑰ in the KSLOCK mount ⑮ on the holding system. If necessary, turn the clamping jaw gently to and fro until it audibly clicks into place. Secure with knurled screw ⑯ (turn in the direction of the lock arrow).

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
2. Insert the instrument into the corresponding holder on the clamping jaw.
3. Tighten the screw (see arrow).

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7.5 Status displays

The illuminated ring ⑥ provides details on the system's operating and error states.

- green flashing: device ready, battery inserted
- blue (permanently illuminated): request to make the device ready for operation
- green (permanently illuminated for 5 seconds): device is ready for operation
- yellow (permanently illuminated): indication that the battery voltage is too low. Check the precise battery level directly at the battery charge display.
- red (permanently illuminated): error message (e.g., battery temperature too high), system can no longer be operated

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8 Operation

⚠️ WARNUNG!

The device may be subject to sudden failure during use.

- You should plan to have a replacement device at the ready or employ alternative procedures.

⚠️ WARNUNG!

Risk of injury: Before the ARTIP BASE holding system can be used on the patient, the system must be checked with the planned instrument (load).

NOTE: Before each use, check that the working area of the arm corresponds to the surgical requirements.

NOTE: The holding arm can only be operated when the battery is charged and inserted.

NOTE: Observe the status display. If the display is red, the device is blocked.

- Remove the battery and insert it again. If the error persists, please contact our Service.

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8.1 Releasing and moving the holding arm

1. Press all 3 touch plates at the same time. All joint blocking is released and the holding arm can be moved to the desired position.

NOTE: Upon actuating all three touch plates, the blocking of all holding arm joints is released and the load of the arm and the working load are placed fully in the hands of the user.

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8.2 Fixing the joints of the holding arm

1. Release the touch plates on the handle ⑨ (one touch plate is sufficient). All joints are blocked immediately.

WARNING!

Let go of the handle slowly after releasing the touch plates and check whether the applied part is indeed held. If the system doesn't act as expected, stop the procedure and swivel the holding system out of the operating field.

WARNING!

Should the holding arm no longer release due to an error (e.g., if the display on the illuminated ring is red), remove the working load and, if possible, swivel the arm out of the operating field. Then disassemble the arm from the operating rail.

NOTE: The system only requires a current for the setting processes between release and lock. If an empty battery has to be exchanged during operation, this does not affect the holding force. The control only allows the system to be released if the remaining battery charge is sufficient for subsequent locking. This prevents the battery from having to be changed when the arm is released.

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

9 Disassembly of the ARTIP BASE holding system

1. Undo the thumb screw on the clamping jaw ⑰ and remove the instrument from the clamping jaw.
2. Undo the knurled screw ⑮ (turn in the direction of the arrow) and pull the clamping jaw out of the KSLOCK mount.

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3. To remove the holding arm ⑩ from the base unit, press the touch plates on the handle ⑬ with one hand and move the holding arm into an upright position as far as possible.
4. Hold down the touch plates and with your second hand, press down the actuation ring ⑦ on the base unit coupling (see Fig.).
5. The holding arm can now be pulled out upwards.

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7. Press the two clamps on the right and left of the battery housing and pull out the battery ⑤ from the base unit.

NOTE: The device is only currentless if the battery is removed from the device.

NOTE: Recharge the battery after each use.

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- This reprocessing step is part of the validated process.
- This reprocessing step may be carried out additionally in order to aid cleaning. The sterilization method is approved in respect of material compatibility, however must be validated by the operator on site in respect of effectiveness.
- x On the grounds of potential damage, this reprocessing step must not be performed.

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10.2 General warnings

⚠️ WARNING!

Risk of infection: Incorrectly reprocessed medical devices expose patients, users and third parties to the risk that the medical device may malfunction.

- Observe the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' and the accompanying documentation.

⚠️ WARNING!

Risk of infection: These products are not sterile when delivered. The use of non-sterile products poses a risk of infection for patients, users and third parties.

- Inspect reusable products for visible contamination. Reprocess contaminated products.
- Reprocess reusable products before and after each use. When doing so, use validated procedures.

NOTE: When carrying out any work on contaminated medical devices, the guidelines of the Employers' Liability Insurance Association and equivalent organizations striving to ensure personal safety must be observed.

⚠️ ATTENTION!

When preparing and using solutions, follow the chemical manufacturer's instructions, paying close attention to proper concentration, exposure time and service life. Prolonged immersion and incorrect concentration may result in damage.

- Bear in mind the microbiological range of action of the chemicals used.

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⚠ ATTENTION!

Danger of damage to the product: The use of chemicals that have not been approved by KARL STORZ or chemical manufacturers may cause damage to the medical devices.

- Only use chemicals approved for the medical device for reprocessing. A helpful guide is available on the Internet at www.karlstorz.com.

NOTE: National laws and regulations must be observed.

10.3 Reprocessing of ARTIP BASE holding arms

10.3.1 Manual wipe-down disinfection

Wipe-clean the exterior surfaces of the medical device with a disposable cloth moistened with disinfectant or a ready-to-use soaked disinfectant cloth. Due to their protein-fixating effect and possible material incompatibility, alcohol-based agents must not be used. The chemical manufacturer's specifications must be observed. At the end of the necessary exposure time, wipe the surface with a dry low lint cloth.

10.3.2 Accessories

Accessories required for carrying out reprocessing:
Brushes

10.3.3 Preparation for cleaning and disinfection

Heavy soiling, corrosive solutions and pharmaceuticals must be removed from the medical device immediately after use. To this end, preclean the medical device by wiping down and rinsing, for example. As a general rule, KARL STORZ recommends manual precleaning under cold running water.

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10.3.4 Manual precleaning Brushing the surfaces

Any visible contamination or heavy soiling must be removed by cleaning the surfaces under cold running water with the aid of a brush or sponge.

10.3.5 Machine cleaning and disinfection

The following methods for machine decontamination have been validated and approved subject to compliance with the process parameters described in the manual, 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' (Art. no. 96216003):

Machine cleaning/thermal disinfection

Thermal disinfection is preferred. The relevant national requirements and the A_0 value must be taken into account when using this method. The selection of a suitable slide-in tray or instrument holder, which should ensure that the medical device is thoroughly rinsed out or through, must take place in consultation with the manufacturer of the device.

NOTE: If necessary, the instrument must be dried off manually afterwards.

10.3.6 Assembly, inspection and care

The cleaned and disinfected medical device must be visually inspected for cleanliness, completeness, damage and dryness:

- If residues or contamination are still present, the medical device must be manually cleaned and subjected to a full cleaning and disinfection procedure once more.
- Damaged or corroded medical devices must be withdrawn from use.
- Disassembled medical devices must be assembled.
- A functional check must be carried out before use during surgery.

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⚠ WARNING!

The joints of the holding arm must not be oiled or greased under any circumstances, as this reduces the holding force.

10.3.7 Packaging systems

Only standardized and approved packaging materials and systems may be used (EN 868 Parts 2-10, EN ISO 11607 Parts 1 + 2, DIN 58953).

10.3.8 Sterilization

The procedures as well as the process-relevant parameters for the individually validated methods are described in detail in the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' (Art. no. 96216003). The method must be selected taking into account the respective applicable national requirements and in consultation with the device and product manufacturers.

The following sterilization methods have been validated and approved by KARL STORZ for this medical device:

Steam sterilization using the fractionated prevacuum procedure

The medical device must be sterilized in its fully assembled state using the fractionated prevacuum procedure (DIN EN ISO 17665-1) at 132°C-137°C with a minimum exposure time of 3 minutes and a maximum of 18 minutes. This procedure is only suitable for thermostable instruments.

10.4 Limits of reprocessing

The end of the product's service life is largely determined by wear, reprocessing methods, the chemicals used and any damage resulting from use.

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10.5 Reprocessing of the drive unit 28272 TM, rechargeable battery 28272 TA

⚠ ATTENTION!

Not suitable for machine cleaning and sterilization. Only suitable for wipe-down disinfection.

Manual wipe-down disinfection

The exterior surfaces must be wiped clean with a disposable cloth moistened with disinfectant or a ready-to-use soaked disinfectant cloth. Due to their protein-fixating effect and possible material incompatibility, alcohol-based agents must not be used. Observe the chemical manufacturer's specifications regarding material compatibility. At the end of the necessary exposure time for the disinfectant as specified by the manufacturer, wipe the surface with a dry low lint cloth.

10.6 Reprocessing of the charging station

⚠ ATTENTION!

Not suitable for machine cleaning and sterilization. Only suitable for wipe-down disinfection.

Manual wipe-down disinfection

The exterior surfaces must be wiped clean with a disposable cloth moistened with disinfectant or a ready-to-use soaked disinfectant cloth. Due to their protein-fixating effect and possible material incompatibility, alcohol-based agents must not be used. The chemical manufacturer's specifications regarding material compatibility must be observed. At the end of the necessary exposure time for the disinfectant as specified by the manufacturer, wipe the surface with a dry low lint cloth.

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11 Service, repair and disposal

11.1 Guarantee

The warranty conditions provided can be found in the Standard Conditions of Business of KARL STORZ.

The medical device must always be sent to your competent KARL STORZ subsidiary (see the section entitled 'Subsidiaries'), even during the warranty period.

Opening of the equipment or performance of any repairs or modifications to the device by unauthorized persons shall relieve KARL STORZ of any liability for its performance. Any such opening, repair or modification performed during the guarantee period shall void all warranty.

11.2 Maintenance

No further repair measures are required beyond reprocessing.

11.3 Safety inspection/repeat inspection according to IEC 62353

WARNING!

Regardless of the accident prevention regulations or testing intervals for medical devices prescribed in different countries, safety checks/repeat inspections according to IEC 62353 must be carried out on this device and documented at least once a year by qualified personnel.

Visual inspection

1. Check the device and accessories for any mechanical damage which may impair functionality.
2. Check that inscriptions relevant to safety are legible.

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⚠️ WARNING!

*If any defects or shortcomings which could
endanger patients, personnel or others are
ascertained during these safety checks/repeat
inspections, the device must not be operated until
such time as these defects or shortcomings have
been eliminated by qualified technical servicing.*

11.4 Repairing device parts

Only personnel authorized by the manufacturer
may repair device parts. The operations described
in this manual are exempt from this rule.

The authorized personnel may only repair device
parts which the manufacturer designates as
repairable.

Upon request, the manufacturer will provide
authorized personnel with information on repairs,
e.g., circuit diagrams, spare parts lists, setting
instructions.

The fact of having such documents at one's
disposal does not correspond to authorization.

11.5 Repair program

For this medical device, individual repair is
necessary. To bridge the repair period, you will
receive a device on loan, which must then be
returned as soon as you receive the repaired
medical device.

In the case of repairs, please contact your
local KARL STORZ subsidiary (see section
“Subsidiaries”) or authorized dealer.

11.6 Infection prevention

To prevent infections, contaminated devices and
instruments must not be returned.

Decontaminate instruments and devices before
returning. KARL STORZ returns contaminated
instruments to the sender.

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11.7 Disposal

This device has been marked in accordance with the European Directive on waste electrical and electronic equipment (WEEE).

At the end of its useful operating life, dispose of the device as electronic scrap.

Please ask KARL STORZ, a KARL STORZ subsidiary (see the section "Subsidiaries"), or your authorized dealer about the collection point in your area.

Within the scope of application of this Directive, KARL STORZ is responsible for the proper disposal of this device.

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⚠ ATTENTION:

Damage to product:

- Observe the following operation/ storage conditions for the device and the accessories.

	Operation	Transport and storage
Temperature	+5 ... +40°C	-10 ... +60°C
Rel. humidity	20 – 85%	5 – 95%
Air pressure	3,000 m	500 – 1080 hPa

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12 Directive compliance

Class I medical device

This medical device bears a CE mark in accordance with MDD 93/42/EEC. A code number after the CE mark indicates the responsible notified body.

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13 Electromagnetic Compatibility user information

The ARTIP BASE holding system 28272 TLS and 28272 TSS is suitable for use in professional healthcare facility environment.

Professional healthcare facility includes physician offices, dental offices, limited care facilities, freestanding surgical centers, freestanding birth centers, multiple treatment facilities, hospitals (emergency rooms, patient rooms, intensive care, surgical rooms, outside the RF shielded room of an ME system for magnetic resonance imaging).

WARNING!

Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

WARNING!

The use of an accessory or cable with the ARTIP BASE holding system 28272 TLS and 28272 TSS other than those specified in this manual may result in increased emissions or decreased immunity of the ARTIP BASE holding system 28272 TLS and 28272 TSS. When using an accessory or cable with the ARTIP BASE holding system 28272 TLS and 28272 TSS other than those specified in this manual, it becomes the responsibility of the user of the ARTIP BASE holding system 28272 TLS and 28272 TSS to determine compliance with IEC 60601-1-2 when using this item.

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⚠️ WARNING!

Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the ARTIP BASE holding system 28272 TLS and 28272 TSS including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

NOTE: The emission characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 Class A). If it is used in a residential environment (for which CISPR 11 Class B is normally required) this equipment might not offer adequate protection to radio communication service. The user might need to take mitigation measures, such as relocating or re-orienting equipment.

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Table 1 Guidelines and manufacturer's declaration – electromagnetic immunity			
The ARTIP BASE holding system 28272 TLS and 28272 TSS is intended for use in an electromagnetic environment specified below. The user of the C-MAC® Monitor 843 ZX should make sure that it is used in such an environment.			
Immunity tests	EN/IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) acc. to IEC 61000-4-2	± 8 kV contact discharge ± 15 kV air discharge	± 8 kV contact discharge ± 15 kV air discharge	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Rapid transient electrical Interferences/bursts acc. to IEC 61000-4-4	± 2 kV for power lines ± 1 kV for input and output lines 100 kHz repetition	n/a	
Surges acc. to IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	n/a	
Voltage dips, short interruptions and fluctuations of the supply voltage acc. to IEC 61000-4-11	<u>Voltage Dip:</u> Dip to 0% for 1 cycles @ 0° phase angle Dip to 70% for 25/30 cycles @ 0° phase angle Dropout to 0% for 0.5 cycles @ 0°, 45°, 90°, 135°, 180°, 225°, 270° & 315° phase angles <u>Voltage interruption:</u> 100% for 250/300 cycles	n/a	If image distortion occurs, it may be necessary to position the ARTIP BASE holding system 28272 TLS and 28272 TSS further from sources of power frequency magnetic fields or to install magnetic shielding. The power frequency magnetic field should be measured in the intended installation location to assure that it is sufficiently low.
Magnetic field at the supply frequency (50/60 Hz) acc. to IEC 61000-4-8	30 A/m at 50Hz/60Hz	30 A/m at 50Hz/60Hz	
Radiated, radio-frequency, electromagnetic field immunity test acc. to IEC 61000-4-3	3V/m 80 MHz to 2.7 GHz *Refer table 2 for wireless Proximity RF field test levels	10 V/m 80 MHz to 2.7 GHz	
Immunity to conducted disturbances, induced by radio-frequency fields acc. to IEC 61000-4-6	3V _{rms} on 150 kHz to 80MHz 1kHz 80% AM modulation 6V _{rms} in ISM bands	n/a	

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Table 3
Guidance and manufacturer's declaration – electromagnetic immunity –
for equipment and systems that are not life-supporting

The ARTIP BASE holding system 28272 TLS and 28272 TSS is designed for operation in an electromagnetic environment as defined below. The user of the device should make sure that it is used in such an environment.

Interference immunity tests	EN/IEC 60801 test level	Conformity level	Electromagnetic environment – guidelines
			Portable and mobile RF communications equipment should be used no closer to any part of the ARTIP BASE holding system 28272 TLS and 28272 TSS, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
Conducted RF disturbances acc. to IEC 61000-4-6	3 V _{rms} 150 kHz to 80 MHz	n/a	
Radiated RF disturbances acc. to IEC 61000-4-3	3 V/m 80 MHz to 2.7 GHz	10 V/m	Recommended separation distance: $d = 1.2 \sqrt{P}$ 80 MHz to < 800 MHz $d = 2.3 \sqrt{P}$ 800 MHz to 2.7 GHz
			Interference may occur in the vicinity of equipment marked with the following symbol: 

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. The propagation of electromagnetic quantities is affected by absorptions and reflections of buildings, objects and persons.

a) Field strengths from fixed transmitters, such as base stations for radio [cellular/cordless] telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to stationary transmitters, an electromagnetic site survey should be considered. If the measured field strength at the location where the ARTIP BASE holding system 28272 TLS and 28272 TSS is used exceeds the above compliance levels, the ARTIP BASE holding system 28272 TLS and 28272 TSS should be monitored to ensure proper function. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device.

b) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

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Table 4 Guidelines and manufacturer's declaration – electromagnetic emissions		
The ARTIP BASE holding system 28272 TLS and 28272 TSS is intended for use in the electromagnetic environment specified below. The customer or user of the ARTIP BASE holding system 28272 TLS and 28272 TSS should ensure that it is used in such an environment.		
RF emissions as per CISPR 11	Group 1	The ARTIP BASE holding system 28272 TLS and 28272 TSS uses RF energy only for its internal function. The customer or user of the ARTIP BASE holding system 28272 TLS and 28272 TSS should make sure that it is used in such an environment. Therefore, its RF emissions are very low and are not likely to cause interference in nearby electronic equipment.
RF Emissions as per CISPR 11	Class A	The ARTIP BASE holding system 28272 TLS and 28272 TSS is suitable for use in all establishments including domestic and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions as per IEC 61000-3-2	Class A	
Voltage fluctuations/flicker Emissions as per IEC 61000-3-3	n/a	

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Table 5
Recommended separation distances between portable and mobile RF telecommunication devices and the ARTIP BASE holding system 28272 TLS and 28272 TSS

The ARTIP BASE holding system 28272 TLS and 28272 TSS is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the ARTIP BASE holding system 28272 TLS and 28272 TSS can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the ARTIP BASE holding system 28272 TLS and 28272 TSS as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter [W]	Separation distance d [m] as a function of the transmitter frequency		
	150 kHz to 80 MHz $d = \left[\frac{3,5}{3} \right] \sqrt{P}$	80 MHz to 800 MHz $d = \left[\frac{3,5}{3} \right] \sqrt{P}$	800 MHz to 2.5 GHz $d = \left[\frac{7}{3} \right] \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance D in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note: At 80 MHz and 800 MHz, the separation distance for higher frequency range applies.

Note: These guidelines may not apply in all situations. The propagation of electromagnetic quantities is affected by absorptions and reflections of structures, objects and people.

Compatibility with surgical HF instruments

This device was tested for compatibility with HF surgical devices.

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KARL STORZ Endoscopy Canada Ltd.
7171 Millcreek Drive, Mississauga, Ontario
L5N 3R3 Canada
Phone: +1 905 816-4500, Fax: +1 905 816-4509
Toll free phone: 1-800-268-4880 (Canada only)
Toll free fax: 1-800-482-4188 (Canada only)
E-Mail: info-canada@karlstorz.com

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100, Fax: +1 424 218-8525
Toll free phone: 800 421-0837 (USA only)
Toll free fax: 800 321-1304 (USA only)
E-Mail: communications@ksea.com

KARL STORZ Veterinary Endoscopy-America, Inc.
1 South Los Cameros Road
Goleta, CA 93117, USA
Phone: +1 805 968-7776, Fax: +1 805 685-2588
E-Mail: info@karlstorzvet.com

KARL STORZ Endoscopy Latino-America, Inc.
815 N.W. 57th Avenue, Suite 480
Miami, FL 33126-2042, USA
Phone: +1 305 262-8980, Fax: +1 305 262-8986
E-Mail: info@ksea.com

KARL STORZ Endoscopy México S.A. de C.V.
Edificio Atlantic, Oficina 3G
Calle D e/ 1ra y 3ra
10400 Vedado, Havana, Cuba
Phone: +537 836 95 06, Fax: +537 836 97 76
E-Mail: kstorzcuaba@gmail.com

KARL STORZ Endoscopy México S.A. de C.V.
Lago Constanza No 326,
Col. Granada, Del. Miguel Hidalgo
C.P. 11520 Ciudad de México
Phone: +52 (55) 1101 1520
E-Mail: mx-info@karlstorz.com

KARL STORZ Marketing America De Sul Ltda.
Rua Joaquim Floriano, nº. 413, 20º andar – Itaim Bibi,
CEP-04534-011 São Paulo, Brasil
Phone: +55 11 3526-4600, Fax: +55 11 3526-4680
E-Mail: br-info@karlstorz.com

KARL STORZ Endoscopy Argentina S.A.
Zufriategui 627 6º Piso, B1638 CAA - Vicente Lopez
Provincia de Buenos Aires, Argentina
Phone: +54 11 4718 0919, Fax: +54 11 4718 2773
E-Mail: info@karlstorz.com.ar

KARL STORZ Endoskopi Norge AS
Stamveien 1
1483 Hagan, Norway
Phone: +47 6380 5600, Fax: +47 6380 5601
E-Mail: post@karlstorz.no

KARL STORZ Endoskopi Sverige AB
Storsåtergränd 14
127 39 Skårholmen, Sweden
Phone: +46 8 505 648 00
E-Mail: kundservice@karlstorz.se
KARL STORZ Endoscopy Suomi OY
Tavallie 5
01610 Vantaa, Finland
Phone: +358 (0)96824774, Fax: +358 (0)96824775
E-Mail: asiakaspalvelu@karlstorz.fi

KARL STORZ SE & Co. KG
Representation Office
Kestutis at. 59 / Lenktolaj st. 27
08124 Vilnius, Lithuania
Phone: +370 5 272 0448
Mobile: +370 685 67 000
E-Mail: info-lt-hu@karlstorz.com

KARL STORZ Endoscopy Denmark A/S
Skovlystvej 33
2840 Holte, Denmark
Phone: +45 45162600, Fax: +45 45162609
E-Mail: marketing@karlstorz.dk

KARL STORZ Endoscopy UK Ltd.
415 Perth Avenue, Slough
Berkshire, SL1 4TO, United Kingdom
Phone: +44 1753 503500, Fax: +44 1753 578124
E-Mail: info-uk@karlstorz.com

KARL STORZ Endoscopy Nederland B.V.
Displayweg 2
3821 BT Amersfoort, Netherlands
Phone: +31 (0)33 4545890
E-Mail: info-nl@karlstorz.com

KARL STORZ Endoscopy Belgium N.V.
Phone: +31 (0)33 4545890
E-Mail: info-be@karlstorz.com

KARL STORZ Endoscopy France S.A.S.
12, rue Georges Guynemert, Quartier de l'Europe
76280 Guyancourt, France
Phone: +33 1 30484200, Fax: +33 1 30484201
E-Mail: marketing-fr@karlstorz.com

KARL STORZ Endoscopy Austria GmbH
Landskroffer Hauptstr. 148/1G1
1030 Wien, Austria
Phone: +43 1 71 56 0470, Fax: +43 1 71 56 0479
E-Mail: storz-austria@karlstorz.at

KARL STORZ Endoscopy Ibérica S.A.
Parque Empresarial San Fernando
Edificio Munich – Planta Baja
28830 Madrid, Spain
Phone: +34 91 6771051, Fax: +34 91 6772981
E-Mail: info-es@karlstorz.com

KARL STORZ Endoscopy Italia S.r.l.
Via dell'Argentario, 3
37135 Verona, Italy
Phone: +39 045 8222000, Fax: +39 045 8222001
E-Mail: info-it@karlstorz.com

KARL STORZ Endoscopy Croatia d.o.o.
Capraška 6
10000 Zagreb, Croatia
Phone: +385 1 6406 070, Fax: +385 1 6406 077
E-Mail: info@karlstorz.hr

KARL STORZ Endoskopijska d.o.o.
Cesta v Gorici 34b
1000 Ljubljana, Slovenia
Phone: +386 1 620 5880, Fax: +386 1 620 5882
E-Mail: pisarna@karlstorz.si

KARL STORZ Endoscopy Polska Marketing Sp. z o.o.
ul. Bogkowska 47
44-100 Gliwice, Poland
Phone: +48 32 706 13 00, Fax: +48 32 706 13 07
E-Mail: info-pl@karlstorz.com

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TOO KARL STORZ Endoscopy Kazakhstan
Saryarka, 8, BC "Arman", off. 910
010000 Astana, Republic of Kazakhstan
Phone: +7 7172 552-549, 552-788, Fax: +444
E-Mail: info@karlstorz.kz

KARL STORZ ENDOSKOPE
East Mediterranean & Gulf (branch)
Building West Side 7A – Unit 7WA – 3006
Dubai Airport Free Zone, P.O. Box 54983
Dubai – United Arab Emirates
Phone: +971 (0)4 2956887, Fax: +971 (0)4 3205282
Service Hotline: +971 (0)4 3416882
E-Mail: info-gne@karlstorz-emg.com

KARL STORZ Endoscopy India Private Limited
11th Floor, Dr. Gopal Das Bhawan
28, Barakhamba Road
New Delhi 110001, India
Phone: +91 11 4374 3000, Fax: +91 11 4374 3010
E-Mail: corporate@karlstorz.in

KARL STORZ SE & Co. KG
Resident Representative Office
14th Floor, Plaza Saigon
39 Le Duan, District 1
Ho Chi Minh City, Vietnam
Phone: +84 28 3823 8000, Fax: +84 28 3823 8039
E-Mail: infovietnam@karlstorz.com

KARL STORZ SE & Co. KG
Interchange 21 Tower, Level 33
390 Sukhumvit Road,
North Klongkroey, Wattana,
10110 Bangkok, Thailand
Phone: +66 2 261 442 9500, Fax: +66 2 261 442 9030

KARL STORZ Endoscopy Australia Pty. Ltd.
15 Orion Road Lane Cove NSW 2086
P.O. Box 50 Lane Cove NSW 1585, Australia
Phone: +61 (0)2 9490 6700
Toll free: 1800 996 562 (Australia only)
Fax: +61 (0)2 9420 0695
E-Mail: info@karlstorz.au

KARL STORZ Endoscopy China Ltd., Hong Kong
Room 2503-05, 25F AXA Tower, Landmark East,
No. 100 How Ming Street, Kowloon, Kowloon,
Hong Kong, People's Republic of China
Phone: +852 28 65 2411, Fax: +852 28 65 4114
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd.
Beijing Branch,
Room 1805-1807, Building B, 18F Beijing IFC
No. 8, Jianguomenwai Street, Chaoyang District,
100022, Beijing, People's Republic of China
Phone: +86 10 56381888, Fax: +86 10 56381999
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shanghai Branch
36F, Tower 1 Grand Gateway
No. 1 Hongqiao Road, Xuhui District,
200030, Shanghai, People's Republic of China
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Japan K. K.
Stage Bldg. 8F, 2-7-2 Fujimi
Chiyoda-ku, Tokyo 102-0071, Japan
Phone: +81 3 6380-8622, Fax: +81 3 6380-8633
E-Mail: info@karlstorz.co.jp

KARL STORZ Endoscopy (Shanghai) Ltd.
Shanghai Branch
36F, Tower 1 Grand Gateway
No. 1 Hongqiao Road, Xuhui District,
200030, Shanghai, People's Republic of China
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Chengdu Branch
Room 805, 8F Jin Jiang International Building
No. 1 West Linjiang Road, Wuhou District,
610041, Chengdu, People's Republic of China
Phone: +86 28 86587977, Fax: +86 28 86587975
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shenyang Branch,
Room 2001-2005, 20F M-MEDIA International Center,
No. 167 Youth Avenue, Shenyhe District,
110014, Shenyang, People's Republic of China
Phone: +86 24 23181118, Fax: +86 24 23181119
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Guangzhou Branch
Room 02B, 35F Team Tower,
No. 208 Tianhe Road, Tianhe District,
510620, Guangzhou, People's Republic of China
Phone: +86 20 8732-1281, Fax: +86 20 8732-1286
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com

KARL STORZ Endoscopy Singapore Sales Pte Ltd
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com

KARL STORZ Endoscopy Korea Co Ltd
9F Hyowon-Building
97, Jangdae-ro, Songpa-gu
05719 Seoul, Korea
Phone: +82-70-4350-7474, Fax: +82-70-8277-3299
E-Mail: info@karlstorz.com

KARL STORZ Endoscopy Taiwan Ltd.
12F, No. 192, Sec. 2, Chung Hsin Rd.,
Sindian District, New Taipei City, Taiwan
Phone: +886 933 014 160, Fax: +886 2 8672 6399
E-Mail: info-tw@karlstorz.com

KARL STORZ SE & Co. KG
Representative Office Philippines
1901 Picadilly Star Bldg., 4th Avenue, BGC
Taguig City 1636, Philippines
Phone: +83 2 317 45 00, Fax: +83 2 317 45 11
E-Mail: philippines@karlstorz.com

KARL STORZ Endoscopy Japan K. K.
Stage Bldg. 8F, 2-7-2 Fujimi
Chiyoda-ku, Tokyo 102-0071, Japan
Phone: +81 3 6380-8622, Fax: +81 3 6380-8633
E-Mail: info@karlstorz.co.jp

KARL STORZ Endoscopy (Shanghai) Ltd.
Shanghai Branch
36F, Tower 1 Grand Gateway
No. 1 Hongqiao Road, Xuhui District,
200030, Shanghai, People's Republic of China
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shanghai Branch
36F, Tower 1 Grand Gateway
No. 1 Hongqiao Road, Xuhui District,
200030, Shanghai, People's Republic of China
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199
E-Mail: info@karlstorz.com.cn

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