
December 16, 2022

To Whom It May Concern:

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
This version of the Operator's Manual Supplement is intended for use with the CorPath GRX System software version 3.0.X with PCI/PVI Application and only on the territory of the Russian Federation.

Corindus, Inc

309 Waverley Oaks Road, Suite 105

Waltham, MA, USA 02452

Phone: 800.605.9635 or 508.653.3335

Fax: +1 (508) 653-3335

info@corindus.com | www.corindus.com

Section 1. General Information	2
The name of the Medical Device.....	2
Information about Manufacturer of the Medical Device.....	2
A. Intended use, indications, contraindications, and possible side effects.....	2
B. Conditions for use.....	3
C. Classification of the Medical Device	4
Section 2. Technical Specifications of the Medical Device	6
Section 3. The list of components of the medical device	9
Section 4. The list of International Standards / Directives / Documents applicable to the Medical Device.....	12
Section 5. Labeling	13
A. Labeling on the packaging	13
B. Component Labels	16
C. Artwork of the medical device labeling in Russian language	20
Section 5. Complaints.....	21

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Type of the Medical Device in accordance with the Russian Nomenclature Classification

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Classifications	Contents
Type of protection against electrical shock	Protection against electrical shock through reliable. Protection Earth (PE) and not in basic insulation alone. Class I according to IEC 60601-1
Degree of protection against electrical shock	Type CF applied part rating for the guide wire and catheter accessories: - Robotic Drive

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Classifications	Contents
	- Extended Reach Arm - Control Console
Degree of protection against harmful ingress of water	CorPath GRX Control Console meets the vertical drip-proof fluid ingress classification of IPX1, and the Robotic Drive meets the angled drip-proof fluid ingress classification of IPX2 in accordance with IEC 60529.
Degree of safety of application in the presence of a flammable anesthetic material with air or with oxygen or nitrous oxide	Equipment not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide.
Mode of operation	Continuous operation

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Software Version

Software Version: 3.0.X with PCI/PVI Application;

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30" Midsize Display:

Dimensions: 714 mm x 524 mm x 240 mm;
Height adjustment: 100 mm;
Weight: 17.5 kg;
Resolution: 2560 x 1600.

32" Midsize Display:

Dimensions: 761 mm x 471 mm x 90 mm;
Weight: 18.6 kg;
Resolution: 3840 x 2160.

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No	Name of component	Description
Variant 1		
1.	Robotic Drive.	Receives inputs from the Control Console and mounted onto the Extended Reach Arm via lock joint. These inputs received from the Control Console actuate mechanical operations in the Single-use Cassette as movements to a loaded guidewire, RX device, and guide catheter. The Robotic Drive moves forward and backwards to advance and retract the guide catheter. A Bedside Touchscreen on the Robotic Drive provides the user interface for bedside operator during system set-up, loading and device exchanges.
2.	Extended Reach Arm.	Mounted on the rail of the procedure table and supports the Robotic Drive. Can be positioned for any access point by sliding on the rail or adjusting using buttons on the Robotic Drive.
3.	The Control Console.	The Control Console enables physicians to manipulate the guide catheter, guidewire, and RX device through a touchscreen and three joysticks. The Control Console touchscreen advances and retracts the Rapid Exchange (RX) device, guidewire, or guide catheter in discrete 1 mm movements with the touch screen arrows. In addition to the advance/retract function, the catheter and guidewire joysticks have rotational control. To rotate catheter or guidewire in a clockwise/ counterclockwise direction, turn joystick clockwise/ counterclockwise. The control console touchscreen can also rotate the guidewire in discrete increments by touching either the clockwise or counterclockwise rotation buttons. The CorPath GRX System software offers the option for continuous motion or discrete step movement of devices. It continuously monitors the joystick positions and interprets the inputs for minimum and maximum speeds. The CorPath GRX System comes with the Rotate on Retract (RoR) feature, which is a procedural automated movement that can be enabled by pressing the RoR button on the Control Console touchscreen.
4.	Robotic Drive Power Supply.	24V Robotic Drive Power Supply with AC Mains Power.
5.	Control Console Power Supply.	12V Control Console Power Supply with AC Mains Power.
6.	Communication Cable in variants: - Communication Cable 6 Meter; - Communication Cable 20 Meter; - Communication Cable 30 Meter.	Communication Interface Cable between Control Console and Extended Reach Arm/ Robotic Drive
7.	Basic Consumable Set (no more than 10 sets): - Single-use Cassette (4 pcs.); - The Drive Gear (4 pcs.); - CorPath GRX Robotic Drive Drape (4 pcs.);	The Basic Consumable Set comes with a box of 4 per consumable item, where the main item of the set is the Single-use Cassette that comes sterile and attaches to the Robotic Drive and can be loaded with commercially available guide catheters, guidewires, and RX devices. This Single-use Cassette advances, retracts, and rotates the guidewire and guide catheter, and advances and retracts the RX device.

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No	Name of component	Description
	- CorPath GRX Control Console Drape (4 pcs.).	The set come with 3 additional accessories: The Drive Gear, CorPath GRX Robotic Drive Drape, and CorPath GRX Control Console Drape. - The Drive Gear which connects between the guide catheter and the Y-Connector is a sterile single-use gear to enable the rotation of a guide catheter when placed into the Single-use Cassette. - The consumable set comes also with the CorPath GRX Robotic Drive Drape, which is a sterile, single-use plastic layer that covers the Robotic Drive and the Extended Reach Arm. The Robotic Drive and Extended Reach Arm should be draped with a CorPath GRX Robotic Drive Drape prior to every CorPath procedure. The drape is specifically designed to fit over each system component to maintain sterility. - The CorPath GRX Control Console Drape is also included, which is a sterile, single-use plastic layer that covers the Control Console. The Control Console can be covered with sterile drapes if the operator wishes to be scrubbed during the procedure. The whole Consumable Basic Set components are provided pre-sterilized via ethylene oxide (EtO), and are intended for single patient use. The whole Consumable Basic Set components are incompatible with re-sterilization.
8.	User Documentation in Russian.	User Documentation in Russian
Accessories		
1.	Connection Artis	Identifying the connection to an ARTIS Angio system
2.	Connection Other Vendor	Identifying the connection to an OEM Angio system.
3.	CorPath Cockpit Set Stand./Midsized	The Cockpit set contains the Corindus Cockpit (lead shielded workplace for the robot operator from radiation exposure while maneuvering interventional devices, which provides space to place the Corindus console as well as a mounting place for the chosen display option).
4.	CorPath Cockpit Set PowerVision Imaging System	The Cockpit set contains the Corindus Cockpit (lead shielded workplace for the robot operator from radiation exposure while maneuvering interventional devices, which provides space to place the Corindus console as well as a mounting place for the chosen display option).
5.	Chair for Interventional Cockpit	Chair for Interventional Cockpit.
6.	Standard Displays	19" monitors for Live and REF images. The part number contains 2 monitors including power cables, video cables and monitor stand.
7.	Midsized Display	30"/32" monitor displaying a 1:1 copy of the Large/Midsized display.
8.	Power Vision Display	40" 4k color monitor which displays a 1:1 copy of the large display.
9.	Sony SRG-X120 Camera	Remote controllable camera to provide visual overview of patient and X-Ray system. Recommended for set-ups where you do not have an unobstructed line of sight between the operator and the patient table.
10.	Control Room Set Power Vision	Various Hardware parts including cabling to connect the Power Vision Display for OEM systems in the control room.
11.	Video Integration with Splitter	Recommended for X-Ray systems which have a large display. This setup uses a copy of the large display at the Boom to be displayed at the Corindus control console.

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№	Name of component	Description
		The set includes a DVI splitter, which duplicates the large display at the boom for use with a 4k monitor at the Corindus control console.

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
B. Component Labels*Extended Reach Arm*

CorPath[®] GRX Extended Reach Arm

Corindus A Siemens Healthineers Company

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

Authorized Representative
MDSS GmbH
Schiffgraben 41
30175 Hannover, Germany

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REF 308
Quantity 1
LOT 12345C1234
Batch Code 2019-11-20
Date of Manufacture 2024-11-20
Use by Date

STERILE EO
Sterilized using Ethylene Oxide Single Sterile Barrier System

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EC Regulatory Matters
MDSS GmbH Schiffgraben 41
D-30175 Hannover Germany

CE 2797

Rx Only Caution: Federal law restricts this device to sale by or on the order of a Physician

California's "Prop 65" law requires the following statement on many medical devices, even though they are properly processed per all current FDA and international standards.
Warning: this product contains a chemical known to the state of California to cause cancer, birth defects or other reproductive harm.
(trace residuals from EtO sterilization)

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REF 310
Quantity 1
LOT 12345C1234
Batch Code 2019-11-20
Date of Manufacture 2024-11-20
Use by Date

STERILE EO
Sterilized using Ethylene Oxide Single Sterile Barrier System

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EC Regulatory Matters
MDSS GmbH Schiffgraben 41
D-30175 Hannover Germany

CE 2797

Rx Only Caution: Federal law restricts this device to sale by or on the order of a Physician

California's "Prop 65" law requires the following statement on many medical devices, even though they are properly processed per all current FDA and international standards.
Warning: this product contains a chemical known to the state of California to cause cancer, birth defects or other reproductive harm.
(trace residuals from EtO sterilization)

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Section 5. Complaints

The expected service life of the system is 5 years.

In case of any questions or problems with handling the system, please use the following contacts:

Contacts for questions and reclamations in Russia:

Siemens Healthcare LLC

Dubiniskaya street 96

115093 Moscow,

Russia

Tel.: +7 495 737 1252

Customer support center

8 800 550 50 22

(24 hour free line)

service.ru@siemens-healthineers.com

Contacts for questions and reclamations in the country of Legal manufacturer:

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December 16, 2022

To Whom It May Concern:


I, Robert Lavado, certify that the following Operator's Manual Supplement for Russia has been printed by the manufacturer in the USA.

Sincerely,

Robert Lavado
Head of Regulatory Affairs

On this day sixteen day of December 2022, before me, the undersigned notary public, Robert Lavado, personally appeared, proved to me through satisfactory evidence of identification, which was a Massachusetts Driver's License, to be the person whose name is signed on the preceding or attached document in my presence.

Signature of Notary

Paul L. Speidel Jr

Name of Notary Typed, Stamped or Printed
Notary Public, Commonwealth of
Massachusetts


Robotic system for carrying out transdermal coronary and peripheral interventions CorPath® GRX with accessories

Operator's Manual Supplement for Russia

For variant 2

CorPath GRX is the equivalent name of the medical device and is used below.

This version of the Operator's Manual Supplement is intended for use with the CorPath GRX System software version 3.0.X with NVI (Neurovascular) Application and only on the territory of the Russian Federation.

Corindus, Inc

309 Waverley Oaks Road, Suite 105

Waltham, MA, USA 02452

Phone: 800.605.9635 or 508.653.3335

Fax: +1 (508) 653-3335

info@corindus.com | www.corindus.com

Section 1. General Information	2
The name of the Medical Device.....	2
Information about Manufacturer of the Medical Device.....	2
A. Intended use, indications, contraindications, and possible side effects	2
B. Conditions for use.....	3
C. Classification of the Medical Device	4
Section 2. Technical Specifications of the Medical Device	6
Section 3. The list of components of the medical device	9
Section 4. The list of International Standards / Directives / Documents applicable to the Medical Device.....	12
Section 5. Labeling	13
A. Labeling on the packaging	13
B. Component Labels	16
C. Artwork of the medical device labeling in Russian language	21
Section 5. Complaints.....	22

Section 1. General Information

The name of the Medical Device

Robotic system for carrying out transdermal coronary and peripheral interventions CorPath® GRX with accessories.

CorPath GRX is the equivalent name of the medical device and is used below.

Information about Manufacturer of the Medical Device

Legal Manufacturer

Corindus, Inc.

309 Waverley Oaks Road, Suite 105,
Waltham, Massachusetts 02452, USA
Phone: 800.605.9635 or 508.653.3335

Fax: (508) 653-3355

info@corindus.com

www.corindus.com

Addresses of the manufacturing site

Corindus, Inc.

309 Waverley Oaks Road, Suite 105,
Waltham, Massachusetts 02452, USA
Phone: 800.605.9635 or 508.653.3335

Fax: (508) 653-3355

info@corindus.com

www.corindus.com

The name and address of the local representative in Russia

Siemens Healthcare LLC

Dubininskaja street 96,
115093 Moscow, Russia

Tel.: +7 495 737 1252

www.siemens-healthineers.com/ru/

A. Intended use, indications, contraindications, and possible side effects

Intended use

The CorPath GRX System is intended for use in the remote delivery and manipulation of devices during percutaneous coronary and vascular procedures, including the following devices:

- Guidewires
- Rapid Exchange Catheters
- Guide Catheters

- Microcatheters
- Embolization Coils
- Coil Assist Stents

Indications

The CorPath GRX System is intended for use in the remote delivery and manipulation of devices during percutaneous coronary and vascular procedures, including the following devices:

- Guidewires
- Rapid Exchange Catheters
- Guide Catheters
- Microcatheters
- Embolization Coils
- Coil Assist Stents

Contraindications

There are no known contraindications for CorPath GRX.

Possible side effects

There are no known possible side effects for CorPath GRX.

B. Conditions for use

For use only by qualified medical professionals with adequate training and experience in angiographic procedures and use of the CorPath GRX System.

Environmental ConditionsSystem

- Operating Temperature: 15°C to 30°C
- Operating Relative Humidity: Keep dry, protect from moisture

Cassette / Single Use

- Operating Temperature: 15°C to 30°C
- Operating Relative Humidity: Keep dry, protect from moisture

Transportation requirementsControl Console, Robotic Drive, Extend Reach Arm

- Transportation Temperature: -20°C to +60°C
- Transportation Relative Humidity: Keep dry, protect from moisture

Cassette

- Transportation Temperature: -18°C to +60°C
- Transportation Relative Humidity: Keep dry, protect from moisture

Drive Gear, Y-Connector Adapter, Microcatheter Support Clip

Operator's Manual Supplement for Russia

- Transportation Temperature: +10°C to +30°C
- Transportation Relative Humidity: Keep dry, protect from moisture

Other components

- Transportation Temperature: -18°C to +40°C
- Transportation Relative Humidity: Keep dry, protect from moisture

StorageControl Console, Robotic Drive

- Storage Temperature: -2°C to +40°C
- Storage Relative Humidity: Keep dry, protect from moisture

Extend Reach Arm

- Storage Temperature: -2°C to +40°C (+15°C to +30°)
- Storage Relative Humidity: Keep dry, protect from moisture

Cassette

- Storage Temperature: +10°C to 30°C
- Storage Relative Humidity: Keep dry, protect from moisture

Drive Gear, Y-Connector Adapter, Microcatheter Support Clip

- Storage Temperature: +15°C to 30°C
- Storage Relative Humidity: Keep dry, protect from moisture

Other components

- Storage Temperature: 15°C to 30°C
- Storage Relative Humidity: Keep dry, protect from moisture

C. Classification of the Medical Device

Type of the Medical Device in accordance with the Russian Nomenclature Classification
240330.

Class of the potential risk of the medical device application

Class IIb.

Class of the Medical Device in accordance with IEC 60601-1

Classifications	Contents
Type of protection against electrical shock	Protection against electrical shock through reliable. Protection Earth (PE) and not in basic insulation alone. Class I according to IEC 60601-1
Degree of protection against electrical shock	Type CF applied part rating for the guide wire and catheter accessories: - Robotic Drive - Extended Reach Arm

Operator's Manual Supplement for Russia

Classifications	Contents
	- Control Console
Degree of protection against harmful ingress of water	CorPath GRX Control Console meets the vertical drip-proof fluid ingress classification of IPX1, and the Robotic Drive meets the angled drip-proof fluid ingress classification of IPX2 in accordance with IEC 60529.
Degree of safety of application in the presence of a flammable anesthetic material with air or with oxygen or nitrous oxide	Equipment not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide.
Mode of operation	Continuous operation

Class of medical waste

Class B - epidemiologically hazardous waste.

Software Risk classification

Risk classification in accordance with IEC 62304: Class C.

Software Version

Software Version: 3.0.X with NVI (Neurovascular) Application

Date of release: 04/2021

Section 2. Technical Specifications of the Medical Device

Tolerances: $\pm 5\%$ unless otherwise specified.

Physical Characteristics

Robotic Drive*:

Dimensions: 690 mm x 410 mm x 250 mm;

Weight: 14.8 kg.

** Includes portion of Extended Reach Arm permanently attached to Robotic Drive*

Extended Reach Arm*:

Dimensions (lowest storage position): 610 mm x 560 mm x 790 mm;

Weight: 24.8 kg.

** Excluding portion of arm permanently attached to Robotic Drive*

Control Console:

Dimensions: 430 mm x 530 mm x 230 mm;

Weight: 7.6 kg.

Robotic Drive Power Supply:

Dimensions: 200 mm x 80 mm x 41 mm;

Output cable length: 900 mm;

Weight: 0.62 kg.

Mains Input:

Voltage: 100-240V AC;

Frequency: 50-60-Hz;

Power: 165W.

Output:

Voltage: 24V DC;

Power: 150W.

Control Console Power Supply:

Dimensions: 192 mm x 62.2 mm x 38.7 mm;

Output cable length: 900 mm;

Weight: 0.70 kg.

Mains Input:

Voltage: 100-240V AC;

Frequency: 50-60-Hz;

Power: 140W.

Output:

Voltage: 12V DC;

Power: 165W.

Communication Cable:

Length: 6000 mm/20 000 mm/30 000 mm;
Diameter: 17 mm.

Single-use Cassette – Advanced:

Dimensions: 530 mm x 150 mm x 80 mm;
Weight: 1.3 kg.

Drive Gear:

Dimensions: 19.05 mm x 19.05 mm x 17.907 mm;
Weight: 1.75 grams.

Microcatheter Support Clip:

Dimensions: 50.8 mm x 25.4 mm x 15.24 mm;
Weight: 4.6 grams.

Y-Connector Adapter:

Dimensions: 50.8 mm x 11.05 mm x 24.4 mm;
Weight: 3.4 grams.

CorPath Cockpit Set Stand./Midsize / CorPath Cockpit Set PowerVision Imaging System:

Dimensions: 813.1 mm x 1322.3 mm x 1394.7 mm
Weight: 118.8 kg;
Cockpit radiation attenuation equivalent: 0,5 mm Al;
Maximum movement force: < 46 N (horizontal flat surface).

Chair for Interventional Cockpit:

Total Height: 1092 mm;
Seat Height: 406.4 mm to 508 mm;
Width: 711.2 mm;
Depth: 457.2 mm;
Weight: 19.5 kg.

19" Standard Displays:

Dimensions: 423 mm x 208 mm x 349 mm;
Height adjustment: 82 mm;
Weight: 5.1 kg;
Resolution: 1280 x 1024.

30" Midsize Display:

Dimensions: 714 mm x 524 mm x 240 mm;

Hight adjustment: 100 mm;

Weight: 17.5 kg;

Resolution: 2560 x 1600.

32" Midsize Display:

Dimensions: 761 mm x 471 mm x 90 mm;

Weight: 18.6 kg;

Resolution: 3840 x 2160.

Power Vision Display:

Dimensions: 933.7 mm x 570.7 mm x 84 mm;

Weight: 22.5 kg;

Resolution: 3840 x 2160.

Camera:

Dimensions: 158.4 mm x 177.5 mm x 200.2 mm;

Weight: 1.8 kg;

Image Sensor: 1/2.5 type ExmorR CMOS;

Image Sensor (Number of Effective Pixels): 8.5M pixels;

Pan Angle: Pan: $\pm 170^\circ$;

Tilt Angle: $+90^\circ/-20^\circ$;

Power Requirements: DC 12V (DC 10.8 to 13.2V), PoE+ (IEEE802.3at compliant);

Power Consumption: 25.5 W;

Camera Body Color: White/Black.

Section 3. The list of components of the medical device

No	Name of component	Description
Variant 2		
1.	Robotic Drive (Neuro)	Receives inputs from the Control Console and mounted onto the Extended Reach Arm via lock joint. These inputs received from the Control Console actuate mechanical operations in the Single-use Cassette as movements to a loaded guidewire, RX device, and guide catheter. The Robotic Drive moves forward and backwards to advance and retract the guide catheter. A Bedside Touchscreen on the Robotic Drive provides the user interface for bedside operator during system set-up, loading and device exchanges.
2.	Extended Reach Arm (Neuro).	Mounted on the rail of the procedure table and supports the Robotic Drive. Can be positioned for any access point by sliding on the rail or adjusting using buttons on the Robotic Drive.
3.	The Control Console (Neuro).	The Control Console enables physicians to manipulate the guide catheter, guidewire, and RX device through a touchscreen and three joysticks. The Control Console touchscreen advances and retracts the Rapid Exchange (RX) device, guidewire, or guide catheter in discrete 1 mm movements with the touch screen arrows. In addition to the advance/retract function, the catheter and guidewire joysticks have rotational control. To rotate catheter or guidewire in a clockwise/ counterclockwise direction, turn joystick clockwise/ counterclockwise. The control console touchscreen can also rotate the guidewire in discrete increments by touching either the clockwise or counterclockwise rotation buttons. The CorPath GRX System software offers the option for continuous motion or discrete step movement of devices. It continuously monitors the joystick positions and interprets the inputs for minimum and maximum speeds. The CorPath GRX System comes with the Rotate on Retract (RoR) feature, which is a procedural automated movement that can be enabled by pressing the RoR button on the Control Console touchscreen.
4.	Robotic Drive Power Supply.	24V Robotic Drive Power Supply with AC Mains Power.
5.	Control Console Power Supply.	12V Control Console Power Supply with AC Mains Power.
6.	Communication Cable in variants: - Communication Cable 6 Meter; - Communication Cable 20 Meter; - Communication Cable 30 Meter.	Communication Interface Cable between Control Console and Extended Reach Arm/ Robotic Drive.
7.	Advanced Consumable Set (no more than 10 sets): - Single-use Cassette Advanced (4 pcs.); - The Drive Gear (4 pcs.); - Microcatheter Support Clip (4 pcs.);	The Advanced Consumable Set comes with a box of 4 per consumable item, where the main item of the set is the Single-use Cassette Advanced that comes sterile and attaches to the Robotic Drive and can be loaded with commercially available guide catheters, guidewires, and RX devices. This Single-use Cassette Advanced advances, retracts, and rotates the guidewire and guide catheter, and advances and retracts the RX device.

Operator's Manual Supplement for Russia

No	Name of component	Description
	- Y-Connector Adapter (4 pcs.); - CorPath GRX Robotic Drive Drape (4 pcs.); - CorPath GRX Control Console Drape (4 pcs.).	The set come with 5 additional accessories: The Drive Gear, Microcatheter Support Clip, Y-Connector Adapter, CorPath GRX Robotic Drive Drape, and CorPath GRX Control Console Drape. - The Drive Gear which connects between the guide catheter and the Y-Connector is a sterile single-use gear to enable the rotation of a guide catheter when placed into the Single-use Cassette. - Microcatheter Support Clip: for microcatheter use: the plastic microcatheter support clip connects to the support track when using microcatheters. - Y-Connector Adapter: for microcatheter use: the Y-Connector Adapter attaches to the copilot connected to the base catheter (guide catheter, long sheath, or intermediate catheter), or intermediate catheter. - The consumable set comes also with the CorPath GRX Robotic Drive Drape, which is a sterile, single-use plastic layer that covers the Robotic Drive and the Extended Reach Arm. The Robotic Drive and Extended Reach Arm should be draped with a CorPath GRX Robotic Drive Drape prior to every CorPath procedure. The drape is specifically designed to fit over each system component to maintain sterility. - The CorPath GRX Control Console Drape is also included, which is a sterile, single-use plastic layer that covers the Control Console. The Control Console can be covered with sterile drapes if the operator wishes to be scrubbed during the procedure. The whole Consumable Set components are provided pre-sterilized via ethylene oxide (EtO), and are intended for single patient use. The whole Consumable Set components are incompatible with re-sterilization.
8.	User Documentation in Russian.	User Documentation in Russian
Accessories		
1.	Connection Artis	Identifying the connection to an ARTIS Angio system.
2.	Connection Other Vendor	Identifying the connection to an OEM Angio system.
3.	CorPath Cockpit Set Stand./Midsize	The Cockpit set contains the Corindus Cockpit (lead shielded workplace for the robot operator from radiation exposure while maneuvering interventional devices, which provides space to place the Corindus console as well as a mounting place for the chosen display option), as well as all the necessary parts to connect it.
4.	CorPath Cockpit Set Power Vision	The Cockpit set contains the Corindus Cockpit (lead shielded workplace for the robot operator from radiation exposure while maneuvering interventional devices, which provides space to place the Corindus console as well as a mounting place for the chosen display option), as well as all the necessary parts to connect it.
5.	Chair for Interventional Cockpit	Chair for Interventional Cockpit.
6.	Standard Displays	19" monitors for Live and REF images. The part number contains 2 monitors including power cables, video cables and monitor stand.
7.	Midsize Display	30"/32" monitor displaying a 1:1 copy of the Large/Midsize display.
8.	Power Vision Display	40" 4k color monitor which displays a 1:1 copy of the large display.

Operator's Manual Supplement for Russia

No	Name of component	Description
9.	Sony SRG-X120 Camera	Remote controllable camera to provide visual overview of patient and X-Ray system. Recommended for set-ups where you do not have an unobstructed line of sight between the operator and the patient table.
10.	Control Room Set Power Vision	Various Hardware parts including cabling to connect the Power Vision Display for OEM systems in the control room.
11.	Video Integration with Splitter	Recommended for X-Ray systems which have a large display. This setup uses a copy of the large display at the Boom to be displayed at the Corindus control console. The set includes a DVI splitter, which duplicates the large display at the boom for use with a 4k monitor at the Corindus control console.

Section 4. The list of International Standards / Directives / Documents applicable to the Medical Device

Standard	Title
EN ISO 13485:2016	Medical Devices – Quality Management Systems – Requirements for Regulatory Purposes.
EN ISO 14971:2019	Medical Devices – Application of Risk Management to Medical Devices
EN 60601-1:2005	Medical Electrical Equipment – Part 1: General Requirements for Basic Safety and Essential Performance
EN 60601-1-2:2015	Medical Electrical Equipment – Part 1-2: General Requirements for Basic Safety and Essential Performance – Collateral Standard: Electromagnetic Disturbances – Requirements and Tests
EN 62304:2006	Medical Device Software – Software Life-Cycle Processes
IEC 62366-1:2015	Medical Devices – Application of Usability Engineering to Medical Devices
BS EN ISO 20417:2021	Information Supplied by the Manufacturer of Medical Devices
EN ISO 11607-1:2020	Packaging for Terminally Sterilized Medical Devices – Part 1: Requirements for Materials, Sterile Barrier Systems and Packaging Systems
EN ISO 11135:2014	Sterilization of Health Care Products – Ethylene Oxide – Part 1: Requirements for Development, Validation and Routine Control of a Sterilization Process for Medical Devices
EN 556-1:2001	Sterilization of Medical Devices – Requirements for Medical Devices to be Designated as “STERILE” – Part 1: Requirements for Terminally Sterilized Medical Devices
EN ISO 10993-1:2018	Biological Evaluation of Medical Devices – Part 1: Evaluation and Testing Within a Risk Management Process
EN ISO 15223-1:2021	Medical Devices – Symbols to be Used with Medical Device Labels, Labeling and Information to be Supplied – Part 1: General Requirements
ISO 80369-7:2021	Small-Bore Connectors for Liquids and Gases in Healthcare Applications – Part 7: Connectors for Intravascular or Hypodermic Applications

Section 5. Labeling

Given labels are examples only.

A. Labeling on the packaging

Extended Reach Arm

Shipping box:

Robotic Drive

Shipping box:

*Control Console**Shipping box:**Advanced Cassette**Shipping box:*


Drive Gear

Shipping box:

CorPath GRX Drapes

Shipping box:

B. Component Labels*Extended Reach Arm*

CorPath® GRX Extended Reach Arm

Corindus
A Siemens Healthineers Company

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
Y-Connector Adapter

CorPath® GRX Y-Connector Adapter

Corindus
A Siemens Healthineers Company

Manufacturer
Corindus, Inc.
309 Waverley Oaks Road
Suite 105
Waltham, MA 02452, USA

EC REP Authorized Representative
MDSS GmbH
Schriffgraben 41
30175 Hannover, Germany

(01) 0 0816280 02313 0

(10) A123456 (17) 000000

REF 313 LOT A123456 YYY-MM-DD YYY-MM-DD Qty: 1

150-08338 Rev C

STERILE EO 2797

MD Rx

30°C

10°C

CorPath GRX Drapes

CorPath® GRX Drape

Intended for Use With:
Robotic Drive & Arm

REF 308

QTY 1

LOT 12345C1234

2019-11-20

2024-11-20

Corindus
A Siemens Healthineers Company

PACKAGE MUST BE ASEPTICALLY OPENED

CONTENTS STERILE UNLESS SEAL IS
BROKEN OR PACKAGE IS PUNCTURED

MD Rx

EC Regulatory Matters
MDSS GmbH Schriffgraben 41
D-30175 Hannover, Germany

CE 2797

Rx Only Caution: Federal law restricts this
device to sale by or on the order of a Physician

California's "Prop 65" law requires the following statement on many medical devices, even
though they are properly processed per all current FDA and international standards.
Warning: this product contains a chemical known to the state of California to cause cancer,
birth defects or other reproductive harm.
(trace residuals from ETO sterilization)

Made in China

32443-303-004

150-04560-01 Rev C

Robotic Drive Drape

CorPath® GRX Drape

Intended for Use With:
Control Console

REF 310

QTY 1

LOT 12345C1234

2019-11-20

2024-11-20

Corindus
A Siemens Healthineers Company

PACKAGE MUST BE ASEPTICALLY OPENED

CONTENTS STERILE UNLESS SEAL IS
BROKEN OR PACKAGE IS PUNCTURED

MD Rx

EC Regulatory Matters
MDSS GmbH Schriffgraben 41
D-30175 Hannover, Germany

CE 2797

Rx Only Caution: Federal law restricts this
device to sale by or on the order of a Physician

California's "Prop 65" law requires the following statement on many medical devices, even
though they are properly processed per all current FDA and international standards.
Warning: this product contains a chemical known to the state of California to cause cancer,
birth defects or other reproductive harm.
(trace residuals from ETO sterilization)

Made in China

32443-303-004

150-04561-01 Rev C

Control Console Drape

C. Artwork of the medical device labeling in Russian language

Additionally, for all executions (variants) of medical device that are delivered on the territory of Russian Federation, a Russian label is applied with the following content:

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