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
- ① Screw base for KARL STORZ fiberoptic light cables and Olympus/Winter & Ibe
- ② Plug-in base, 9 mm diameter, for Wolf fiber optic light cables
- ③ Light inlet piece without adaptor
- ④ Color-coded ring to indicate the telescope's direction of view
- ⑤ Eyepiece
- ⑥ Sheath
- ⑦ Type designation
- ⑧ 'AUTOCLAVE' marking ring
- ⑨ Lens and distal outlet of light cable
- ⑩ ICG marking


The telescopes named are used within the scope of fluorescence imaging with the NIR/ICG system

1 Intended use

Telescopes are used for the endoscopic display and, if necessary, to create a work and irrigation channels in laparoscopic surgery, gynecology and urology, in diagnostic and operative thoracoscopy, thyroid surgery (extracervical access) as well as transanal interventions. Telescopes are designed for short-term use in surgically invasive interventions.

2 Indications

The products are used in laparoscopy. Laparoscopic intervention in the fields of surgery, urology and gynecology in which access is created into the abdominal cavity through the abdominal wall to allow intraabdominal findings to be diagnosed and treated is fundamentally seen as an indication for these products.

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3 Contraindications

The use of the telescopes is contraindicated if, in the opinion of the attending physician, the surgical method as such is contraindicated or if the patient is not able to undergo surgery or anesthesia due to his or her general condition.

Telescopes must not be used for interventions in direct contact with the central nervous system (CNS) and central cardiovascular system.

4 Safety instructions

WARNING: Risk of injury and damage to products: Failure to observe and follow this instruction manual and the instruction manuals of products used in combination can result in injury to patients, users and third parties as well as damage to the product. Please read all relevant instruction manuals and always follow the instructions given precisely. Check the functioning of the products used in combination.

WARNING: Risk of injury: Incorrectly assembled and damaged instruments can lead to injuries to the patient. Instruments and all of the accessories used in combination with them must be checked immediately before and after every use to ensure that they are complete, free from damage, and in full working order and have no unintentional rough surfaces, sharp corners, burred edges or projecting parts. Care must be taken not to leave missing or broken-off components inside the patient.

WARNING: Risk of injury: Incorrect application of medical instruments poses a risk of injury for patients. Users of medical instruments must have an appropriate medical qualification and be acquainted with the application.

WARNING: Risk of injury and risk of damaging the products: If applied parts are used outside of the field of vision, there is a risk that tissue and accessories could be damaged unintentionally. Always hold the applied parts of the active electrode, laser and other instruments that transmit energy pointed at the target in the field of vision during use.

WARNING: Risk of injury and risk of damaging the products: The highly concentrated light energy can cause the distal end, the light port, adjacent components and tissue in front of the light emission window to heat up (can exceed 41 °C). This can cause burns to the patient, user and surgical accessories.

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
- Avoid direct tissue contact with the distal end, as well as with the light port of the telescope and the light cable.
- Do not rest the endoscope or the light cable on the patient or where they are in direct contact with surgical accessories.
- Always select the lowest possible light output of the cold light source which still allows optimal illumination of the operating field.

WARNING: Risk of injury: Optical radiation poses a risk of injury to eyes. Never look into the free end of either a connected light cable or endoscope.

WARNING: Risk of injury: Patient leakage currents may accumulate if endoscopes and endotherapy devices supplied with energy are used together. This is particularly important if an endoscope with a Type CF applied part is used. In this case, an endotherapy device with a Type CF applied part must be used in order to minimize the total patient leakage current.

WARNING: Risk of injury: Overloading the instrument by exerting too much force may cause the medical device to break, bend, or malfunction, and consequently injure the patient or user. Do not overload the instruments. Do not bend bent instruments back to their original position.

5 Precautions

In order to avoid damage to the sheath and rod lens system and thus ensure the long-term quality of the telescope, the following information must be observed:

- For cleaning, disinfection, sterilization and storage, ideally store the telescopes individually or use containers in which the telescopes can be secured (telescopes 26003ACA/BCA/FCA and 26046ACA/BCA/FCA for use with wire tray 39501B1/B2, telescope 26046FCEA for use with 39501C) (see Fig.).
- Telescopes must always be held by the eyepiece and never by the distal end alone.
- Do not bend the sheath under any circumstances.
- Place the telescopes down with care. Avoid hard impacts, especially at the distal end.

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7 Handling

It is recommended that the suitability of the products for the intended procedure be checked prior to use.

NOTE: Precise information on handling and testing can be found in the instruction manual 'Telescopes and HOPKINS® telescopes' which can be requested from KARL STORZ free of charge or downloaded from our homepage (www.karlstorz.com) (Art.No. 96216001D).

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8 Reprocessing

WARNING: Risk of infection: Incorrectly reprocessed medical devices expose patients, users and third parties to a risk of infection as well as the risk that the medical device may malfunction. Observe the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' and the accompanying documentation.

WARNING: Risk of infection: These medical devices are not sterile when delivered. The use of non-sterile medical devices poses a risk of infection for patients, users and third parties. Inspect medical devices for visible contamination. Visible contamination is an indication that reprocessing has not been carried out, or has been carried out incorrectly. Reprocess the medical devices before initial use and before and after every subsequent use using validated procedures.

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-  **WARNING:** When carrying out any work on contaminated medical devices, the guidelines of the Employers' Liability Insurance Association and equivalent organizations striving to ensure personal safety must be observed.
-  **WARNING:** If CJD (Creutzfeldt Jakob disease) is suspected or indicated, KARL STORZ recommends disposing of the product in accordance with national regulations and not using it again.
-  **CAUTION:** When preparing and using solutions, follow the chemical manufacturer's specifications, paying close attention to proper concentration, exposure time and service life. Prolonged immersion and incorrect concentration may result in damage. Bear in mind the microbiological range of action of the chemicals used.
-  **CAUTION:** Risk of damage to medical devices: The use of chemicals which have not been approved by KARL STORZ may cause damage to the medical devices. Only use chemicals approved by KARL STORZ for reprocessing. A complete list can be found on the Internet at www.karlstorz.com.
-  **CAUTION:** National laws and regulations must be observed.
-  **NOTE:** The manual 'Cleaning, Disinfection, Care, and Sterilization of KARL STORZ Instruments' (Art. no. 96216003D) can be downloaded or requested by visiting www.karlstorz.com.

8.1 Accessories

Accessories required for carrying out reprocessing:
Brushes:
Surface: 27652

8.2 Preparation for cleaning and disinfection

Heavy soiling, corrosive solutions and pharmaceuticals must be removed from the medical device immediately after use. To this end, pre-clean the medical device by wiping down and rinsing, for example. As a general rule, KARL STORZ recommends manual pre-cleaning under cold running water.

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8.3 Manual precleaning

Disassembly

Prior to cleaning and disinfection, the medical device must be separated into its individual components as far as possible, and/or must be opened.

Brushing the surfaces

Any visible contamination or heavy soiling of the surfaces must be removed by cleaning under cold running water with the aid of a brush.

No use of ultrasound

For technical reasons, the medical device is not suitable for ultrasound treatment.

8.4 Machine cleaning and disinfection

The following methods for machine decontamination have been validated and approved subject to compliance with the process parameters described in the instructions, 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' (Art. no. 96216003D):

Machine cleaning/thermal disinfection

Thermal disinfection is preferred. The relevant national requirements and the A_p value must be taken into account when using this method.

The selection of a suitable slide-in tray or instrument holder, which should ensure that the medical device is thoroughly rinsed out or through, must take place in consultation with the manufacturer of the device.

NOTE: If necessary, the instrument must be dried off afterwards by hand.

8.5 Assembly, inspection and care

The cleaned and disinfected medical device must be visually inspected for cleanliness, completeness, damage and dryness:

- If residues or contamination are still present, the medical device must be manually cleaned and subjected to a full cleaning and disinfection procedure once more.
- Damaged or corroded medical devices must be withdrawn from use.
- Dismantled medical devices must be assembled.
- Afterwards, a functional check must be carried out.

NOTE: Use items from the catalog 'Care, Sterilization and Storage Techniques' during care procedures.

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8.6 Life span

The following function tests are to be carried out in order to recognize any loss of material through wear:

1. Check the surface of the medical device for mechanical integrity and changes.
2. Check that the joined components are seated correctly and, if necessary, lock the cleaning connector.
3. Check if the image is transmitted.
4. Check if light is transmitted.

8.7 Packaging systems

Only standardized and approved packaging materials or systems may be used (EN 868 Parts 2-10, EN ISO 11607 Parts 1 + 2, DIN 58953).

8.8 Sterilization

The procedures as well as the process-relevant parameters for the individually validated methods are described in detail in the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' (Art. no. 96216003D). The method must be selected taking the respective applicable national requirements into account and in consultation with the device manufacturers.

The following sterilization methods have been validated and approved by KARL STORZ for this medical device:

Steam sterilization using the fractionated prevacuum procedure

NOTE: Sterilize greased components in a disassembled state.

The medical device must be sterilized in its fully assembled state using the fractionated prevacuum procedure (DIN EN ISO 17665-1) with the following parameters:

- 132 °C–137 °C with a minimum exposure time of 4 minutes and maximum of 18 minutes.

NOTE: This validation is applicable for reusable tube sets which are open at both ends and with a maximum length of 200 cm.

or

134 °C–137 °C with a minimum exposure time of 3 minutes and maximum of 18 minutes.

8.9 Reprocessing limits

The end of the product's service life is largely determined by wear, reprocessing methods, the chemicals used and any damage resulting from use.

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9 Overview: Recommended combinations*

If there is an identification symbol on the light cable and the endoscope (see table), assign matching combinations according to the symbols shown.

Light cable diameter	Endoscope diameter
4.8-5.0 mm 	6.6-12.0 mm 
3.0-3.5 mm 	3.0-6.5 mm 
2.0-2.5 mm 	0.8-2.9 mm 

* Special endoscopes may deviate from this

10 Operating/ storage conditions

There are no special operating and storage conditions for these medical devices.

11 Complaint shipments

In the interests of our employees' health, only sterilized or disinfected products are accepted as return shipments. These should be marked as such; otherwise no further processing may be performed.

12 Disposal

To dispose of these medical devices, no special measures are necessary. National laws and regulations must be observed.

13 Directive compliance

This medical product bears the CE mark in accordance with the Medical Device Directive (MDD) 93/42/EEC. A code number after the CE mark indicates the responsible notified body.

14 Packaging symbols

For the meaning of the symbols printed on the packaging, please refer to the 'Packaging symbols' accompanying instruction leaflet, mat. no. 96216316DF. It can be downloaded from www.karlstorz.com.

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