
 BIO SYSTEMS	Operational Documentation of the Manufacturer in relation of the medical device	Document No.: ODX_BLD_2
Leica Disposable Microtome Blades		Version: 14035859588B Blades V1.0B multi 08/2020
Manufacturer: Leica Biosystems Nussloch GmbH (Heidelberger Str. 17-19, 69226 Nussloch, Germany)	Authorized Representative in the Russian Federation: BioLine LLC • Tax identification number (TIN) applied in the Russian Federation: 7813118690 • Primary state registration number (PSRN) applied in the Russian Federation: 1037804011453	Issue date: 24-08-2021

Operational Documentation

Signature

Frank Volkwein, RA/QA Engineer

Name and Position

Leica Biosystems Nussloch GmbH
 Heidelberger Strasse 17 - 19
 D-69226 Nussloch
 Tel.: 0 62 24 - 1 43 - 0

In accordance with the Order № 11H of the Ministry of Healthcare of the Russian Federation of 19 January 2017 "On approval of requirements to the content of technical and operational documentation of the producer (manufacturer) of medical products", for medical products of 1 (low) risk class it is acceptable to supply shortened operational documentation to a user, provided that such information is enough to operate with the product safely.

The operational documentation for Leica Disposable Microtome Blades contains all necessary information for safety operation.

Hazard warnings:

1. Risk of injury. Blades are very sharp. Only trained laboratory personnel may use them.
2. Never touch the blade edge!
3. Only handle the blades by holding their sides between thumb and index finger crosswise to the blade edge.
4. Wear cut-resistant safety gloves and protective eye wear when handling the blades.
5. Do not use dispenser for disposal of the blades after cutting infectious or radioactive material.

1. Handling instructions:

1. Inspect the dispenser upon receipt. If there is any damage, do not use and contact the manufacturer.
2. If the dispenser does not eject blade, do not force the blade out of the dispenser. Do not use and contact the manufacturer.
3. If the blade is not inserted directly into the blade holder, the blade dispenser must be held in the right hand to remove a blade.
4. To push out the blade, slowly push the blade slide all the way to the left until it stops.
5. The blade must then be held by the left hand and removed from the blade dispenser.
6. If present, remove the protective oil by wiping the blade between thumb and index finger.
7. Insert the blade into the knife holder according to the specifications provided by the instrument manufacturer.
8. If, after checking the additional instrument-specific, cut-relevant parameters, the cut result is not satisfactory, the blade must be moved or replaced.
9. When cutting infectious or radioactive material, change the blade with every new specimen. This prevents unwanted contamination/ transmission.
10. Remove the blade from the knife holder according to the specifications provided by the instrument manufacturer.
11. Observe the following disposal instructions.

2. Disposal:

1. After being used to cut infectious and potentially infectious material, blades belong in the cut-resistant waste container provided for this purpose (not in the waste blade box on the bottom of the blade dispenser).

The product must be disposed of according to existing applicable, local regulations

For the Russian Federation it is recommended to be disposed in accordance with the rules implemented by 'СанПиН 2.1.7.2790-10 (Sanitary and epidemiological requirements for the medical waste handling)' for class "Б (b)" as infectious and potentially infectious waste.

2. Blades that have not come in contact with infectious material, including the non-used products, can be disposed of into the waste blade box through the hole on the bottom of the blade dispenser.

-The holding capacity of the waste blade box is sufficient if all blades are fed into the hole flat.

-The product must then be disposed of according to existing applicable, local regulations

For the Russian Federation it is recommended to be disposed in accordance with the rules implemented by "СанПиН 2.1.7.2790-10 (Sanitary and epidemiological requirements for the medical waste handling)" for class "А" waste (non-hazardous).

3. Technical data:

1. Dimensions:
 - Leica 819 (low profile), DB 80 LS/LX/LC: 80.00 x 8.00 x 0.25 mm
 - Leica 818 (high profile), DB 80 HS/HC: 80.00 x 14.00 x 0.31 mm
2. Material: High-alloyed stainless steel
3. Transport/ Storage temperature range:
 - Transport temperature: - 10 to + 50 °C
 - Storage humidity: < 80 % non-condensing
4. Operating temperature range:
 - Blade dispenser: 15 to 40 °C
 - Blade: - 40 to + 40 °C
 - Operating humidity: <80 % non-condensing

4. Warranty:

Leica Biosystems Nussloch GmbH guarantees that the contractual product delivered has been subjected to a comprehensive quality control procedure based on the Leica in-house testing standards, and that the product is faultless and complies with all technical specifications and/or agreed characteristics warranted. The scope of the warranty is based on the scope of the concluded agreement. The warranty terms of your Leica sales organization or the organization from which you have purchased the contractual product shall apply exclusively.

5. Authorized representative in Russia:

BioLine LLC

- Registered office (legal address):
197022 Saint-Petersburg, municipal district Aptekarsky ostrov, Professora Popova street 23 lit.E,
structure 3H/5H, room 208, Russia
- Post address: 197022 Saint-Petersburg, Professora Popova street 23 lit.E, Russia

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