24


NOTARY PUBLIC

I, the undersigned Gabriel Rudbeck, Notary Public of the City of Stockholm, Sweden,
hereby certify,

that this document is a true and correct copy of the document presented to me.

Stockholm, Sweden 2021-05-18

Ex officio:

Gabriel Rudbeck

Sveavägen 24-26, 3rd floor, PO Box 3098, S-103 61 Stockholm, +46 (0) 8 20 59 90

Approved by
Technical Article Manager Heavy Incontinence & Skin Care
Health and Medical Solution
Essity Hygiene and Health AB

Susanne Wedel

Susanne Wedel

Essity Hygiene and Health AB
405 03 GÖTEBORG
SWEDEN

Instruction for use

Adult absorbent diapers

TENA® Slip Plus (Medium, Large, Extra Large), TENA® Slip Super (Medium, Large)

Essity Hygiene and Health AB, SE-405 03 Goteborg, Sweden

2020

1.Name of the product.

Adult absorbent diapers TENA® Slip Plus (Medium, Large, Extra Large); TENA® Slip Super (Medium, Large).

2.Manufacturer: Essity Hygiene and Health AB, SE-405 03 Goteborg, Sweden

Production site:

- 1.Essity Operations Poland Sp. z o.o., Ul. 3 Maja 30 A, 55-200 Olawa, Poland
2. Essity Operations Gennepe B.V., Hoogveld 25,6598 BL Heijen, Gennepe, The Netherlands
3. Essity Turkey Hijyen Urunleri San.Ve Tic A. S. ;Sarimese Mah. Kanal Yolu Cad, Kerime Sok. № 27 41250, Kartepe, Turkey

3. Intended use of the product

Adult absorbent diapers TENA® Slip Plus and TENA® Slip Super offer all-in-one design for heavy protection with the outstanding dryness and comfort of FeelDry technology. Featuring a soft, cloth-like material for discretion and comfort, and multi-use, velcro-like fasteners that can be simply re-adjusted. This product fits a waist size of 72-122cm for Medium size (M); 92-144cm for Large (L); 120-160 cm for Extra Large (XL). This product is of the 1st class of potential risk of use; for a short period of time it contacts with intact human skin, non-sterile.

4.Indications for use

Incontinence

5. Contraindications

No contraindications found.

6. Side Effects

No side effects found.

7. Precautions

N/A

8. Information about potential customers

This medical device is intended for people suffering from incontinence of heavy degree. The product can be used in a home-care environment or at institutions.

9.Functional characteristics of MD and indications for use

Incontinence of heavy degree.

TENA® Slip	Absorption capacity (g) (Rothwell) ±15 %
Plus Medium	2102
Plus Large	2410
Plus Extra Large	2550
Super Medium	2483
Super Large	2670

Component & and GOST targets	TENA® Slip				
	Plus Medium	Plus Large	Plus Extra Large	Super Medium	Super Large
Sorption capacity after centrifugation, g	not less 440	not less 500	not less 560	not less 440	not less 500
Rewet not more 4,4 g	not more 4,4	not more 4,4	not more 4,4	not more 4,4	not more 4,4
Inlet time cm ³ /sec not less 2,3	not less 2,3	not less 2,3	not less 2,3	not less 2,3	not less 2,3

10. Risks of using a medical product, contraindications, expected and predictable side effects associated with the use of a medical product for the intended purpose.
None


11. Technical characteristics of MD

Measurement point, mm	TENA® Slip Plus Medium &TENA® Slip Super Medium	TENA® Slip Plus Large &TENA® Slip Super Large	TENA® Slip Plus Extra Large
	±10%	±10%	±10%
Length of diaper	828	976	980
Width of diaper	650	800	840
Absorbtion area length (big core length)	630	730	720
Absorbtion area length (small core length)	414	414	414


Weight, g	TENA® Slip Plus (Medium) ±10%	TENA® Slip Plus (Large) ±10%	TENA® Slip Plus (Extra Large) ±10%	TENA® Slip Super (Medium) ±10%	TENA® Slip Super (Large) ±10%
	93	110	121	108	121

Scheme of MD

a. Model TENA® Slip Plus Medium


d. Model TENA® Slip Super Large

e. Model TENA® Slip Plus Extra Large

12. How to apply the product

- **Non - mobile**

Remove the product from the bag and open the front and back side panel. With the user lying on their side, centre and smooth out the pad in the back and slide the pad between the legs from back to front. Make sure good fit in the groin. Centre and smooth out the pad in front. Fit the tapes, one at the time. The top tapes angled slightly downwards, and the bottom tapes angled slightly upwards. Check for comfort and fit. The pad allows for multiple refastening and adjustments of the tapes. Gently press and smooth the tapes on to the textile - like back -sheet. Give a gentle pull to anchor the tapes into the fibres of the back- sheet. This ensures the secure fixation.

- **Semi-mobile**

Remove the product from the bag and open the front and back side panel. With the user standing place the pad between the legs and ensure a good fit in the groin. Centre and smooth out the pad in front and back. Fit the tapes, one at the time. The top tapes angled slightly downwards, and the bottom tapes angled slightly upwards. Check for comfort and fit. The pad allows for multiple refastening and adjustments of the tapes. Gently press and smooth the tapes on to the textile-like back- sheet. Give a gentle pull to anchor the tapes into the fibres of the back -sheet. This ensures the secure fixation.

13. Handling conditions

To handle product under normal conditions.

Requirements for environmental conditions during handling	
Air temperature, °C	+15°C - +35°C
Air relative humidity	not more 40%

14.Materials used

Plus Medium, Super Medium, Plus Large, Super Large

Component	Material	Brand of material	Supplier	MSDS (№ or date)
Top layer	nonwoven	BSNW-SB 12 g/m ²	Fitesa Germany /Italy	15/4/2014
Standing gather nw	SMS nonwoven	B1014040AB SMS NW 14 g/m ²	Fibertex Personal Care A/S	2/7/2018
Transfer layer	PE and nonwoven	NW-AIR 50 g/m ²	Berry Italy	2/2/2016
Absorbent core (double)	Fluf pulp,	Pulp BSK E, FR SS+, 8%;	International Paper	26/5/2015
	SAP	SAP - POWD,GS-301N	LG Chem Limited	3/1/2011
Backsheet	Laminate of PE film and	B-film, 15 g/m ²	RKW;	31/1/2018
	Polypropylene NW	NW-SB, SMS, 15 g/m ²	Fibertex Personal Care A/S	2/7/2018
Elastic threads (leg)	Elastic spandex threads (synthetic)	ELA-Span 800	Hyosung Corporation	July/ 2013
Waist elastic	Thin elastic film	ELA-Film 49 g/m ²	Mondi Gronau GmbH	9/3/2017
Hook/Hook laminate	Polypropylene, adhesives, silicone, nonwoven	FAS-HLAM, TENA bl,NEW	Avery Dennison	28/10/2014
Adhesive	Blend of polymers (hotmelt)	ADH- HM K Perfectcare W 2780	Colquimica S.A.	16/2/2018 20/6/2012 12/2/2016

Plus Extra Large

Component	Material	Brand of material	Supplier	MSDS (№ or date)
Top layer	nonwoven	NW-SB 12 g/m ²	Gulsan, Turkey	20/05/2009

Transfer layer	PE and nonwoven	NW-AIR 50 g/m ²	Berry Italy	10/10/2018
Absorbent core (double)	Fluf pulp,	Pulp BSK E, FR SS+, 8%;	International Paper	26/5/2015
	SAP	SAP - POWD,GS-301N	LG Chem Limited	2013/10/22
Backsheet	Laminate of PE film and	B-film, 14 g/m ²	Koroza,Turkey	29/12/2018
	Polypropylene NW	NW-LAM, 14 g/m ²	Ritas A.S.	01/01/2018
Elastic threads (leg)	Elastic spandex threads (synthetic)	ELA-Span 800	The Lycra company	21/03/2018
Waist elastic	Thin elastic film	ELA-Film 49 g/m ²	Mondi Gronau GmbH	9/3/2017
Hook/Hook laminate	Polypropylene, adhesives, silicone, nonwoven	FAS-HOOK,	Bento AS	25/04/2011
		FAS-HLAM, FAS-DISPTAP	Koroza,Turkey	29/12/2018
Adhesive	Blend of polymers (hotmelt)	ADH- HM Full -Care 5300;	HB Fuller;	05/03/2015
		ADH- HM, HYFIX 1303	Organic Kimya San. ve. Tic.A.S	03/09/2018

15.Labeling of consumer packs and MD

LOT number is usually indicated on the product.

Information on consumer pack:

- medical device name
- medical device name
- indication for use
- size of medical device
- recommended waist size
- number of items per package
- name and address of the location of the manufacturer (seller, distributor),

- trademark
- address of production site
- Lot number
- date of production
- shelf life
- storage conditions
- disposal instructions
- text describing functions.
- pictograms how to apply product
- article
- number and date of registration certificate;
- barcode

Graphic symbols on consumer pack

Symbol	Meaning
	The possibility of processing low density polyethylene (LDPE).
	Throw in the ballot box
	Prohibition of re-use
	Packaging not intended for contact with food (TENA® Slip Plus Extra Large)
	Do not dispose into drains
	Absorption capacity
	CE Marking

	Symbol is placed on products manufactured by companies that finance the waste management program and are included in its disposal system
Proskin	The name of a marketing concept that unites a group of products for people with incontinence
LOT	Number of batch
	Manufacturer
MD logo	Medical device (TENA® Slip Plus Extra Large)

16. Information about transportation of MD

Consumer pack with the diapers enclosed in it should be complete. On the surface, tears and other damage are not allowed.

Temperature -30°C to +35°C; Air relative humidity not more than 80% .

Diapers are transported only in a dry, clean transport that does not have extraneous odors.

17. Storage and Shelf life


The product in the packaging of the enterprise should be stored in a cool place - clean, dry, well ventilated areas; It is strictly forbidden to store products near sources of sharp odors (for example: paintwork materials, household chemicals, burns, solvents, etc., etc.);

When storing and operating the product, direct sunlight and sudden changes in temperature should be avoided.

Environmental conditions during storage

Requirements for environmental conditions during storage	
Air temperature, °C	+15°C - +35°C
Air relative humidity	no more 40%

Shelf life of the diapers is 3 years from the production date indicated on the package.

18.UTILIZATION

The product must be disposed of as medical waste of hazard class "B" Do not throw into drains.

19.Manufacturer's warranty

The manufacturer guarantees the quality of diapers under the conditions of operation, storage and transportation within 3 years from the date of manufacture date indicated on the package.

20.Authorized representative in Russia

During the warranty period, claims are made to the manufacturer at the address: LLC Essity, Russia, 117218 Moscow, ul.Krzhizhanovskogo, 14, building 3, phone: +7 (495) 967-33-67.

21. List of applicable national standards

Harmonized and other International Standards that are relevant to prove conformity to Medical Devices Regulation (EU) 2017/745 by Essity Hygiene and Health AB applied in full or in part.

ISO 13485:2016	Medical devices - Quality management systems – Systems requirements for regulatory purposes This international standard specifies requirements for a quality management system that can be used by an organization for the design and development, production, installation and servicing of medical devices, and the design, development, and provision of related services.
ISO 9001:2015	Quality management systems This international standard specifies requirements for a quality management system where an organization a) needs to demonstrate its ability to consistently provide product that meets customer and applicable regulatory requirements, and b) aims to enhance customer satisfaction through the effective application of the system, including processes for continual improvement of the system and the assurance of conformity to customer and applicable regulatory requirements.
ISO 14971:2012	Risk management This international standard is a framework for effective management by the manufacturer of the risks associated with the use of medical devices. It includes a procedure to investigate, using available information, the safety of a medical device, including accessories, by identifying hazards and estimating the risks associated with the device.
EN ISO 10993:2009	Biological evaluation a) the fundamental principles governing the biological evaluation of medical devices b) the definition of categories of devices based on the nature and duration of contact with the body; c) the selection of appropriate tests.
EN ISO 17025:2005	General Requirements for the competence of testing and calibrating laboratories

	This international standard specifies the general requirements for the competence to carry out tests and/or calibrations, including sampling. It covers testing and calibration performed using standard methods, non-standard methods and laboratory-developed methods.
EN 1041:2008	Information supplied by the manufacturer of medical devices This standard specifies the information to be supplied by a manufacturer for different categories of medical devices, as required by the relevant EEC Directives. It does not specify the language to be used for such information.
ISO 15223:2016	Medical devices – Symbols to be used with medical device labels, labeling and information supplied. This international standard identifies symbols conventionally used to convey information essential for proper use to the user and others for safe and effective use of medical devices.
ISO 7000:2012	Graphical symbols for use on equipment – Index and synopsis. This international standard provides a synopsis of graphical symbols, which are placed on equipment or parts of equipment of any kind in order to instruct the persons handling the equipment as to its use and operation.
EN ISO 14155:2011	Clinical investigation of medical devices for human subjects This international standard defines procedures for the conduct and performance of clinical investigations of medical devices. It specifies general requirements intended to --protect human subjects, --ensure the scientific conduct of the clinical investigation, --assist sponsors, monitors, investigators, ethics committees, regulatory authorities and bodies involved in the conformity assessment of medical devices. Part 1 of ISO 14155 a) specifies requirements for the conduct of a clinical investigation such that it establishes the performance of the medical device during the clinical investigation intended to mimic normal clinical use, reveals adverse events under normal conditions of use, and permits assessment of the acceptable risks having regard to the intended performance of the medical device, b) specifies requirements for the organization, conduct, monitoring, data collection and documentation of the clinical investigation of a medical device, c) is applicable to all clinical investigation(s) of medical devices whose clinical performance and safety is being assessed in human subjects. Part 2 of this standard provides requirements for the preparation of a Clinical Investigation Plan (CIP) for the clinical investigation of medical devices. The compilation of a CIP in accordance with the requirements of this standard and adherence to it will help in optimising the scientific validity and reproducibility of the results of a clinical investigation.

MEDDEV 2.4/1 rev. 9	Guidelines to the Classification of Medical Devices
MEDDEV 2.7.1 rev. 4	Guidelines on Medical Devices, Clinical Evaluation
MEDDEV 2.12/1 rev.8	Guidelines on a Medical Device Vigilance System
Commission Decision 2000/532/EC	European List of Waste

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