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LOT 46605001

 2021-03-31


COBAS, COBAS E
and ELECSYS are
trademarks of Roche.

 Roche Diagnostics GmbH
Sandhofer Str. 116
D-68305 Mannheim

(D) 08086630190
Distribution in USA by:
Roche Diagnostics,
Indianapolis, IN
Made in Germany

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REF 08086630190

For **USA**: **CONTENT**

- M** Streptavidin-coated microparticles 0.72 mg/mL
- R1** HAV Ag (cell culture) 28 U/mL (Roche units)
- R2** Biotinylated monoclonal anti-HAV antibody (mouse) 0.25 µg/mL; monoclonal anti-HAV antibody (mouse) labeled with ruthenium complex 0.65 µg/mL
- Cal1** Anti-HAV negative human serum
- Cal2** Anti-HAV (human) approx. 60 IU/L in human serum

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
English

Precautions and warnings

For in vitro diagnostic use.
Exercise the normal precautions required for handling all laboratory reagents.
Disposal of all waste material should be in accordance with local guidelines.
Safety data sheet available for professional user on request.

For USA: Caution: Federal law restricts this device to sale by or on the order of a physician.

This kit contains components classified as follows in accordance with the Regulation (EC) No. 1272/2008:

2-methyl-2H-isothiazol-3-one hydrochloride

EUH 208 May produce an allergic reaction.

Product safety labeling follows EU GHS guidance.

All human material should be considered potentially infectious. All products derived from human blood are prepared exclusively from the blood of donors tested individually and shown to be free from HBsAg and antibodies to HCV and HIV. The testing methods used assays approved by the FDA or cleared in compliance with the European Directive 98/79/EC, Annex II, List A.

However, as no testing method can rule out the potential risk of infection with absolute certainty, the material should be handled with the same level of care as a patient specimen. In the event of exposure, the directives of the responsible health authorities should be followed.

Avoid foam formation in all reagents and sample types (specimens, calibrators and controls).

Method Sheet download

The product code together with the respective document version (if applicable) is stated in the download information boxes listed above. Please go to the website listed in the relevant box to download this document. The website displays either the version listed here or, if available, any later version that is also valid for this product code. If you have no internet access, please contact your local Roche affiliate to obtain the document free of charge.

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2-Methyl-2H-isothiazol-3-on-Hydrochlorid

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