# INSTRUCTIONS FOR USE

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## PRODUCT NAME

**Solutions viscoelastic for intraocular use: Eyefill® H.D..**

Type of medical device in accordance with the nomenclature classification of medical devices – 250870

## COMPLETNESS OF DELIVERY

Viscoelastic solution for cataract surgery Eyefill® H.D.:

1. Syringe with viscoelastic Eyefill® H.D. solution (2.5 ml solution volume) - 1 pc.
2. Cannula – 1 pc.
3. Instructions for the use of a medical device – 1 pc.

## INFORMATION ABOUT MANUFACTURER AND REPRESENTATIVE:

Valeant Med Sp. z o.o., POLAND

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Tel: +48 22 578 16 00/ Fax +48 22 578 16 18

### Manufacturing site of the medical device:

1. CROMA-PHARMA GmbH, Industriezeile 6, 2100 Leobendorf, AUSTRIA.
2. Sterimedix Ltd, Thornhill Road, North Moons Moat, Redditch, Worcestershire, B98 9ND, UNITED KINGDOM.

### Authorized representative of the manufacturer in the RF:

Bausch Health, Limited Liability Company (Bausch Health LLC)

31 Shabolovka St., Bld. 5, Moscow, 115162, Russia

Phone: +7 (495) 510-28-79; +7 (985) 795-60-34

**INTENDED USE OF THE MEDICAL DEVICE:** Solutions viscoelastic for intraocular use: Eyefill® H.D. are intended for use in cataract surgery as an adjunct and during funduscopy and gonioscopy, as well as used as a contact fluid for diagnostic and therapeutic contact lenses. It protects the intraocular tissues and prevents adhesions and synechiae formation during the procedure.

## APPEARANCE:

**Eyefill® H.D.** - one ready-to-use syringe contains 2.5 ml of HPMC saline solution (pH 6.8-7.6). 1 ml of the solution contains 20 mg (2%) HPMC, disodium hydrogen phosphate dodecahydrate, sodium dihydrogen phosphate dihydrate, sodium chloride and water for injection. The product is delivered in a sterile form and is intended for single use.

The osmolality of the human tear film of the eye in good condition is between 285-300 mOsm/l. Osmolality of the artificial tears is between 270 – 320 mOsmol/kg (Starazolin), 267 (Blink Contacts). In case of the Dry Eye Syndrome it could increase to 360 mOsm/l. Taking into account these values it could be stated that osmolality of all Eyefill products is close to this value. The accepted ranges of this parameter for particular products are as following:

- Eyefill H.D. – 265 – 300 mOsmol/kg.

## COMPOSITION:

It is possible to use solutions as contact fluids when applying lenses and using solutions not only for cataract removal operations, but also for funduscopy and gonioscopy.

PUSTA STRONA

## ROUTE OF ADMINISTRATION AND DOSES:

**Eyefill® H.D.** - the product is administered using a narrow cannula with a blunt tip. The osmolality of the solution is 265-300 mOsmol/kg, the viscosity is approximately 3,200 MPa at a shear rate of  $5 \text{ sec}^{-1}$  ( $T=25^{\circ}\text{C}$ ). The product's amount required depends on the type of surgical procedure.


- The difference between the dynamic viscosity data given in the final specification and the instructions is due to the difference in measurement time. In the specification, the viscosity is measured at  $5 \text{ s}^{-1}$ , and in the instructions, the value is specified at a time of  $0.1 \text{ s}^{-1}$ .

## PRINCIPLE OF OPERATION:


The product provides transparency to the mechanical field during diagnostic and therapeutic intraocular procedures. Due to its special viscoelastic properties, the product creates and maintains the surgical space (including a constant depth of the anterior chamber), protects the corneal endothelium during the surgical procedure, and allows tissue structures to be displaced with less risk of injuring them. At the end of the intervention, the product must be carefully removed by irrigation and aspiration.

## INSTRUCTIONS ON CORRECT ADMINISTRATION OF THE PRODUCT:

- Hold the Luer lock adapter as shown in **fig. 1**

- To remove the tip, turn **fig. 2** and carefully pull it **fig. 3**

**! By following the instructions above, you can prevent air bubbles.**

- Hold the syringe as shown in **fig. 4**

- Firmly insert the attached cannula as shown in **fig.5** (! do not use other cannulas).
- Secure the cannula by turning it clockwise **fig. 6**

*Qm*

PUSTA STRONA

! You hold the syringe correctly if the backstop opens back towards the hand in which it is located.

## ELASTTICITY

Elasticity of the finished product has been determined using an Anton Paar MCR101 Rheometer. Rheological measurement parameters were as follows:

- Measuring System: cone/plate 50 mm  $\square$ , 50  $\mu$ m sample thickness
- Frequency Test performed: (oscillation) for absolute complex viscosity and elasticity
- Deformation: constant at 1%
- Frequency: ramp log from = 628 – 0.1 s<sup>-1</sup>, 24 measuring points
- Temperature: 25 °C

Determination of the elasticity is described in document A\_QC\_TM\_SOP\_009.

Figure 7 shows the elastic modulus  $G'$  and the viscous modulus  $G''$  versus the oscillation frequency (both logarithmic scale) of HPMC 2%.

figure 7 Elastic modulus  $G'$  and viscous modulus  $G''$  versus oscillation of 704001


Figure 7 presents the elasticity ( $= 100 \times [G'/(G'+G'')]$ ) versus the oscillation frequency of HPMC 2%.


figure 8 Elasticity of 704001

*Am*

CC Th oth fill bei int coi AI W pr rel me dis inj vo da an an W

PUSTA STRONA

### CONTRAINDICATIONS AND POSSIBLE COUNTERACTIONS:

The product should not be administered to patients who are allergic to sodium hyaluronate, HPMC or other components of the solution. The product is not indicated for intraarticular use or as a dermal filler. Sodium hyaluronate is not compatible with quaternary ammonium compounds such as benzalkonium chloride solution. Therefore, product containing sodium hyaluronate should not come into contact with surgical instruments rinsed in such solutions, and with ophthalmic devices, which contain ammonium compounds as a preservative. No interactions with other substances were found.

### ADVERSE REACTIONS AND SIDE EFFECTS:

When using viscoelastic agents with sodium hyaluronate, there have been cases of increased intraocular pressure (IOP) after surgery. Usually, IOP rises temporarily: it peaks 4 to 7 hours after surgery and returns to normal after a few days. However, since the maximum can exceed 30 mmHg, IOP should be monitored. Measures to reduce IOP may be required, especially for patients with outflow activation disorders. An isolated case of severe anaphylaxis has been reported with the use of HPMC. Inadvertent injection of OVD into corneal stroma may result in irreversible corneal opacification. Use of excessive volume, applying excessive force during injection or mishandling of the cannula can result in iris damage, posterior capsular rupture. Other rare postoperative reactions including ocular inflammation and/or infection, toxic anterior segment syndrome (TASS), viscoelastic precipitates on the IOL surface, and corneal oedema or decompensation have been reported with the use of OVD.


### WARNINGS:

The use of the product in pregnant or lactating women, as well as in children under 18, has not been studied. Reuse of the product may result in infection of the patient or user. Do not use if the sterile barrier (blister pack) is opened or damaged. Do not use the syringe if its tip has been opened, pushed, damaged or is missing into the sterile barrier (blister pack). Refrigerated product should be allowed to reach room temperature prior to use. Use immediately after opening. Do not inject too much product into the eye. The product must be injected slowly. Avoid injecting into the corneal stroma.

**INFORMATION ON THE PRESENCE OF ANY MATERIALS OF ANIMAL AND (OR) HUMAN ORIGIN IN THE MEDICAL DEVICE FOR HUMAN USE:** No animal materials were used in the manufacture of the product or as raw materials. There are no pharmaceutical additives in the product.

*Run*

PUSTA STRONA

S  
I  
C  
t  
S  
I

Even though the presented products contain substances described in the European and Polish Pharmacopoeia, these substances are not used in accordance with their pharmaceutical mechanism of action. Sodium chloride is used as a buffer component to achieve greater stability of the product. It is not used to compensate for the concentration of chlorides and/or sodium in the human body and to maintain homeostasis. Sodium hyaluronate is used for mechanical/physical purposes, such as protecting tissues during surgery. Pharmacological advantages are not stated.

Since the drugs that are part of all Eyefill® products are used not because of their pharmacological action, but because of their physical / mechanical properties, they are classified in this case as raw materials/components of medical devices. The Technical File contains information about the role of these components, and for this reason they are not classified/considered/listed in the Technical Files.


**INFORMATION ABOUT POTENTIAL CONSUMER:** The product should only be used by ophthalmic surgeons.

**STERILIZATION:** The product is sterilized by steam with the use of a sterile barrier (blister packaging). The product is intended for single use and cannot be re-sterilized after use. Information about the method of cannula sterilization is indicated on the folding box and on the cannula label.

- Sterilization with the use of ethylene oxide.

Below is an image of the label attached to the cannula:

#### STORAGE CONDITIONS:

Eyefill® H.D. - <25 °C, protect against light and freezing

#### TRANSPORT CONDITIONS:

Our stability studies confirm the following acceptable conditions for short-term periods (for example, transportation), without having a negative impact on product quality:

| Name of products | Maximum elevated temperature and duration |
|------------------|-------------------------------------------|
| Eyefill® H.D.    | 40 °C - 7 days                            |

**SHELF LIFE:** Do not use after expiry of the shelf life.

#### LEGEND

Table 1

|  |                                             |
|--|---------------------------------------------|
|  | Caution! Refer to the instructions for use. |
|--|---------------------------------------------|

*Handwritten signature or mark.*


2  
V  
E  
O  
I  
I

PUSTA STRONA

|                                                                                     |                                                                                                                                          |
|-------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------|
| CE                                                                                  | The CE symbol on this product indicates that it has been tested and complies with the provisions of Medical Devices Directive 93/42/EEC. |
|    | Manufacturer                                                                                                                             |
| LOT                                                                                 | Batch number                                                                                                                             |
|    | Do not use if the package is damaged                                                                                                     |
|    | Temperature limit                                                                                                                        |
|    | Keep away from sunlight<br>The cargo should be protected from heat and direct sunlight.                                                  |
|    | Consult instructions for use                                                                                                             |
|    | Keep dry                                                                                                                                 |
|    | Do not re-use                                                                                                                            |
| STERILIZED                                                                          | Sterilized using steam                                                                                                                   |
| STERILE                                                                             | Sterilized using ethylene oxide                                                                                                          |
|    | Do not re-sterilize                                                                                                                      |
|    | No latex was used as a raw material                                                                                                      |
|  | Use by                                                                                                                                   |

## PACKAGING:

## DISPOSAL:

Dispose of these solutions in the same way as other medical waste, such as needles and injection tubes or metal surgical instruments in a medical facility. Do not allow the infection to spread outside the facility or to pollute the environment.

Viscoelastic solutions after use are classified as medical waste of class B and are subject to disposal according to SanPiN 2.1.3684 - 21.

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## REQUIREMENTS FOR ENVIRONMENTAL PROTECTION:

The product does not harm people and the environment when used, transported, and stored following the requirements of the operating documentation.

Quality Management Representative

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