PUBLIC

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
Regulatory Affairs responsible for products of
3Shape TRIOS A/S

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OPERATIONAL DOCUMENTATION TO MEDICAL DEVICE

Intraoral scanner TRIOS® 3, in versions,
produced by 3Shape Trios A/S, Denmark

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
3shape 

TRIOS® Safety and Setup Guide

Intraoral 3D scans for dental CAD/CAM

Models S1A, S1P

Table of Contents – TRIOS® Safety and Setup Guide

1. Dear Customer	4
2. Conformity	4
3. General Information	4
4. Description of TRIOS®	5
5. Safety and Warnings	7
6. Cleaning, Disinfection, and Sterilization	10
7. Disposal	12
8. Technical Maintenance	12
9. Environmental Considerations	12
10. Consumable Accessories	12

1. Dear Customer

Thank you for purchasing your TRIOS® from 3Shape TRIOS®.

The TRIOS® system enables you to perform and digitally save oral scans of a patient's teeth for use in dental CAD/CAM systems.

This TRIOS® safety and setup guide describes safety precautions, warnings, and considerations related to your TRIOS® system. For step by step instructions regarding configuration, and operation of TRIOS®, please refer to the online help.

Before assembling or using TRIOS®, please carefully read and follow the instructions in this safety and setup guide and in the online help system. Please always keep this safety and setup guide near the system.

Please observe all safety information and warnings to prevent personal injury, material damage or damage to your TRIOS® system.

Thank you,
3Shape TRIOS®

2. Conformity

This manual is published by 3Shape TRIOS® A/S, who reserves the right to improve and modify the contents without prior notice. Modifications will, however, be published in future editions. All rights reserved.

TRIOS® is a trademark of 3Shape TRIOS® A/S.

Products manufactured by 3Shape TRIOS® A/S, Holmens Kanal 7, 1060 Copenhagen K, Denmark.

Declaration of Conformity (EU):

3Shape TRIOS® A/S hereby declares that the TRIOS® observes the provisions of: MDD 93/42/EEC with amendments until Directive 2007/47/EEC, 2011/65/EEC RoHS Directive, 2005/20/EEC Packaging and Packaging Waste Directive and 2002/96/EEC Waste electrical and electronic equipment (WEEE), and is in compliance with the following standards: EN (IEC) 60601-1, EN (IEC) 60601-1-1, EN (IEC) 60601-1-2, EN (IEC) 62471.

3. General Information

	CAUTION Be sure to observe all warnings! Please observe all safety information and warnings to prevent personal injury, material damage or damage to your TRIOS® system. Safety information and warnings are highlighted in this guide using the words WARNING, CAUTION or NOTICE.
--	---

3.1 About this Document – Symbols

The symbols used in this document imply the following:

	WARNING Warnings regarding situations where a medium risk of injury to persons exists if the information is not observed.
--	---

	CAUTION Safety information where hazards such as: slight risk of injury to persons, risk of property damage, damage to the system exist if the information is not observed.
	IMPORTANT NOTICE Important information that is not a warning/caution but must be strictly observed.
	TIPS Hints, tips, additional information to support optimal operation of the system.

3.2 Legend of Labels and symbols

	Year of manufacture (EN/ISO 15223-1)
	Company and address (EN/ISO 15223-1)
	General warning – caution (ISO 7010-W001, black on yellow)
	Stepping prohibited (ISO 7010-P019, black and red)
	Type BF equipment (IEC 60417-5333)
	WEEE Product disposal symbol (Directive 2002/96/EC)
	Caution (ISO 7000-0434A)
	CE marking
	Lot number
	Serial Number (EN/ISO 15223-1)
	Consult Instructions for use (EN/ISO 15223-1)
	AC (IEC 60417-5032)
	Protective Earth (IEC 60417-5019)
	Warning – dangerous voltage (IEC 60878) (ISO 3864-B3.6, black on yellow)
	Expiry Date (EN/ISO 15223-1)

	Standby (IEC5009)
	Pushing prohibited (ISO 7010-P017, black and red)
	Local area network, ethernet connector (RJ45)
	Refer to instruction manual / booklet (ISO 7010-M002, blue).
	Connection point for handheld scanner; the scanner has a type BF applied part.
	Identifies the protection tip
	Identifies the calibration tip
	Warning following IEC 62471-2 (black on yellow): Possibly hazardous optical radiation emitted from this product. Do not stare at operating lamp. Under real circumstances (for exposure time, eye pupils, observation distance), it is assumed that no endangerment to the eye exists. The color scanners are risk group 2, the monochrome scanners are risk group "exempt".
	The calibration patches must not be touched and must not be exposed to liquids.
	Mass of the TRIOS® System incl. safe working load (in kilogram).

4. Description of TRIOS®

4.1 Certification and Compliance Notes

CE Mark

This product bears the CE mark in accordance with the provisions of Council Directive 93/42/EEC of June 14, 1993 concerning medical devices with amendments until 2007/47/EC, 2011/65/EEC RoHS Directive, 2005/20/EEC Packaging and Packaging Waste Directive and 2002/96/EEC Waste electrical and electronic equipment (WEEE).

	CAUTION - EU only CE mark for connected products Further products which are connected to this unit must also bear the CE mark.
--	--

Compliance

This system has been tested to conform with the EN (IEC) 60601-1 Safety Standard for medical electrical devices with a Patient applied part (the standard title is: Medical electrical equipment – Part 1: General requirements for basic safety and essential performance), EN (IEC) 60601-1-1 Safety requirements for medical electrical systems, EN (IEC) 60601-1-2 Electromagnetic compatibility – Requirements and tests, and EN (IEC) 62471 Photobiological safety of lamps.

The system has also been tested to conform with CAN/CSA-C22.2 NO 60601-1 (2008) and ANSI/AAMI ES60601-1 (2005).

	WARNING Any person who assembles or modifies equipment or system complying with 60601-1 (safety requirements for medical equipment) by combining it with other equipment for conducting appropriate inspection, continued safe use of the equipment, the requirements of this regulation for the safety of the patients, the operator. Particularly for S1P, see the instruction section Electrical Safety.
--	--

4.2 Intended Use

The 3Shape TRIOS® Intraoral scanner system is a CAD/CAM system intended to obtain 3D digital images of characteristics of teeth directly from the patient's mouth or models of the teeth – for use with CAD/CAM design and production of dental indications as well as with orthodontic help for a list of supported indications.

	CAUTION Unintended use of the TRIOS® system may cause injury to patients, operators and dental staff.
--	---

4.3 Users of TRIOS®

The TRIOS® system may be operated only by trained qualified personnel. Using the TRIOS® system for a purpose other than "intended use" described above may damage the system and operators and patients. Proper use of the TRIOS® system is described in the operation and maintenance instructions, and other caution descriptions in this user guide.

	WARNING Carefully follow the instructions If the instructions for operating the TRIOS® system described in this guide are NOT observed, the protection and safety of the user and patient are at risk.
--	--

4.4 System Description

The TRIOS® Intraoral scanner system enables trained personnel to obtain 3D images of the topographical characteristics of teeth directly from their mouths.

The system is available in these variants:

S1A	Supplied with a cart with built-in PC
S1P	Supplied with a "pod" (holder), to be connected to an external PC via an Ethernet cable

Details are described in the following. You can find more information in section 4 of this guide.

All system variants comprise (Figure 1):

- **Handheld TRIOS® Intraoral scanner (S1):** Handheld scanner for scanning the surface of the patient's teeth. The scanner can scan in color. The scanner has a handle.
- **Scanner Tip (S2):** Re-usable scanner tip attached to the scanner. When the tip is detached, the scanner's body is visible. A set of scanner tips is delivered separately. The scanner tip has an extractable mirror (S4). The scanner tip is used to scan the teeth.
- **Protection Tip (S5):** Protects the optical window of the scanner from damage and dirt when the scanner is not in use.

tion Tip (S6): Used only for calibrating the scanner.

Calibration Kit (S7): Used only for color calibrating TRIOS® Color s. The kit consists of the color calibration target (S8), the color on adaptor (S9) and the sleeve (S10).

Software: Scan application/user interface for the scanner.

(Mains) cable and plug: The system is delivered with a cable late for your location. In the USA and Canada, this is a Nema 5-15P I Grade mains power cord with plug.

3 scanner tip in close proximity to the teeth, a digital image of the ayed to the software and a 3D model is created in the computer. The is stored in the computer along with the patient's case information.

System S1A furthermore comprises:

Figure 2:

scanner mount (C1) for holding the scanner while not in use.

wheels for mobility (and a brake if mobility is not desired).

uter: Contained within the cart – for processing, je and running the TRIOS® software.

1 screen (C2): Touch-enabled screen for real-time display of scan- sults generated by the TRIOS® software, and for interacting with IOS® software.

ystem standby button (C3) is underneath the screen.

anner connection socket for the scanner (C4) is also under- the screen.

3 back, there is a power socket (C5), a network (LAN) connector i fuse drawer (C7), a brake (C8), and a USB port (C9).

IMPORTANT NOTICE

Maximum current available for the USB port is 200mA (typically covering USB2 and USB3 flash keys up to 64GB, external hard discs having their own power supply do not have this limitation).

S1P furthermore comprises:

Figure 3, P1:

mount for the scanner (S1). The type of mount depends on the er having a handle or not.

several **connection sockets**. On one side, there is a single **st for the scanner**. On the other side, there are two connection ts for the Ethernet cable to the PC and for the external medical ' supply, respectively.

1 medical-grade power supply (P2) for the pod.

extension cable to allow the pod to be placed further away from the

CAUTION

The pod must be connected directly to the PC with the Ethernet cable provided with the System.

ing up TRIOS®

ended that unpacking, assembly and connection of the TRIOS® performed by TRIOS® authorized technical service personnel.

CAUTION

If the system has just been delivered from a cold environ- ment, let it adjust to room temperature to avoid internal condensation and to allow for proper calibration.

CAUTION

Examine all parts of the system. Look for physical damage, loose parts or signs of wear that could interfere with proper use and functionality. Contact your service provider if there is visible damage or malfunctioning.

System S1A (cart) only

Please follow the steps below for assembling TRIOS® (Figure 4):

Step 1: Place the cart optimally in relation to the patient chair and your work position.

Step 2: Attach the protection tip to the TRIOS® scanner if not already in place.

Step 3: Place the scanner on its designated scanner mount on the cart.

Step 4: Connect the scanner cable to the scanner connection socket on the cart just under the touch screen.

- To connect - Push the plug into socket. There is a red orientation dot on the socket and plug that should match.
- To disconnect - Pull on the outer release sleeve.

Step 5: Connect the power cable cart end to the socket at the backside of the cart's base.

Step 6: Connect the power cable to a grounded power outlet. Ensure that the cart is placed so it does not block the switch on the mains socket-outlet. To isolate the cart and thus the scanner from supply mains, unplug the flexible power cord from mains.

Step 7: Connect the LAN cable to the cart's LAN connection at the backside of the cart's base.

System S1P (pod) only

Please follow the steps below for assembling TRIOS® (Figure 5):

Before you start: Place the pod on a flat and horizontal surface opti- mally in relation to the patient chair and your work position. Attach the protection tip to the TRIOS® scanner if not already in place. It is easiest not to place the scanner onto the pod at this time yet, because you may find it easier to lift the pod while connecting cables in the following.

Step 1: Connect the Ethernet cable to the pod. The cable goes into the left connection socket on the back side of the pod (the side that has two connection sockets). If you need to release the Ethernet cable, press on the small tab on the top of the connector.

Step 2: Connect the power supply to the pod. This cable goes into the connection socket next to the Ethernet cable described in the previous step.

- To connect - Push the plug into socket. There is a red orientation dot on the socket and plug that should match.
- To disconnect - Pull on the outer release sleeve.

The indicator light next to the scanner connection socket will light green.

Step 3: Connect the scanner cable to the scanner connection socket on the front side of the pod (the side that has only one connection socket).

- To connect - Push the plug into socket. There is a red orientation dot on the plug that should point upwards.
- To disconnect - Pull on the outer release sleeve.

Now all connections to the pod are established, and you can place the scanner onto the pod.

Step 4: Connect the other end of the Ethernet cable to the PC (if no Ethernet ports are available in the PC use an Ethernet to USB adapter).

Step 5: Connect the power (mains) cord to the power supply. The other end of the power (mains) cable goes into a grounded power outlet (not shown in the Figure). To isolate the pod and thus the scanner from supply mains, unplug the power cord from mains.

4.6 Starting and Shutting Down TRIOS®

Please follow the steps below to start TRIOS® for the first time:

Step 1: Turn the system ON

S1A: Press the standby button on the cart for 1 second until it lights green. The TRIOS® software starts automatically.

S1P: Start the PC. If the PC has no TRIOS® software installed, follow the instructions in the accompanying TRIOS® Installation Guide. Then start the TRIOS® software.

Step 2: Log on using a default TRIOS® Operator.

Step 3: Configure and use your TRIOS® system.

When using TRIOS® system the user will be guided by easy-to-under- stand messages shown in the TRIOS® software. Examples of messages include:

- When the scanner has not been calibrated for the last 8 days the user is asked: "It is 8 days since last calibration. Do you want to calibrate the scanner now?"
- When calibration is completed the user is reminded to remove calibra- tion tip: "Ensure that the calibration tip is removed before proceeding."
- When entering the scan page the user is reminded to mount the scan- ner tip: "Please make sure the scanning tip is mounted on the scanner".

Calibration

A calibration of the Handheld TRIOS® scanner is required before using TRIOS® for the first time and every 8 days during regular use or when the scanner has been moved or sustained hits or strong temperature changes. The calibra- tion is done by placing the Calibration Tip on the Scanner Tube, go to the Configure page in the TRIOS® software, select Scan → 3D calibrate scanner, and then follow the on-screen instructions.

For a more detailed description of the calibration please refer to the TRIOS® online help.

TRIOS® Color Scanner only

Color Calibration

A color calibration of the handheld TRIOS® Color scanner is required before each scan in order to get optimal colors.

- Attach a cleaned and sterilized scanner tip on the scanner tube.
- Remove the sleeve from the color calibration target.
- Mount the color calibration target on the color calibration adaptor with the color side facing up.
- Attach the color calibration adaptor on the cleaned and sterilized scanner tip.
- Go to the Configure page in the TRIOS® software, select Scan → Color calibrate scanner, and then follow the on-screen instructions which will also include validation of the performed color calibration.

For a more detailed description of the color calibration please refer to the TRIOS® online help.

Proper shutdown

Inside the TRIOS® software click on the "Close" button and select "Shutdown". System S1A (cart) only: The standby button turns off to confirm that the system is shut down.

CAUTION

When switching off the system, allow enough time for it to shut down properly before disconnecting the mains cable.

System S1A (cart) only

IMPORTANT NOTICE

The standby button can also be used to force TRIOS® shut down in case the screen stops responding. To fc shutdown, push and hold the standby button pressec about 5 seconds. Then, restart the system.

4.7 Starting TRIOS® Online Help

The TRIOS® online help can be started from TRIOS® software by clic the "Help" button and choosing "online help".

4.8 TRIOS® System Environment

Operating Conditions

Ambient temperature range: 15 – 30°C
Relative humidity: 10 – 85 % (Non-condensing)
Atmospheric pressure: 800 – 1100hPa.
The scanner should be used in an environment maintaining a consistent room temperature.

Storage Conditions

Ambient temperature range: -10 – 60°C
Relative humidity: 10 – 85 % (Non-condensing)
Atmospheric pressure: 800 – 1100hPa.

Transport Conditions

Ambient temperature range: -10 – 60°C
Relative humidity: 10 – 85 % (Non-condensing)
Atmospheric pressure: 620 – 1200hPa.

4.9 Technical specifications

Light source specifications	Do not view directly with optical instrur The product complies with EN (IEC) 62 (Photobiological safety of lamps and la systems) and UL- standard US code of regulations, 21 CF 1040.10 and 1040.1
Sterilization method	Scanner tips are autoclavable.
Power input	100-240 VAC, 50/60Hz. S1A: 250VA max (internal power suppl S1P: external class I medical power su rated 30W max, 1.0 A.

5. Safety and Warnings

5.1 Prerequisites

CAUTION

Read all instructions carefully including all warnings a caution notices. You must comply with the warnings i manual to prevent injury to persons and damage to ex ment. Proper functionality and safety can only be gua if the safety precautions in this safety guide and on th TRIOS® system are observed.

Preventive inspection before use of the system

	CAUTION Please examine the TRIOS® system for any mechanical damage on: • all enclosures • all cables Safety can only be guaranteed if NO DAMAGE on the TRIOS® system is observed.
---	--

Modifications of the system

	CAUTION Modifications of the system which may affect the safety of the operator, patients or third parties are prohibited by law!
---	--

Approved software only

	CAUTION Install only approved software – to prevent interference with the runtime reliability of the system and programs within it.
---	--

Proper training

	WARNING Before you attempt to use the TRIOS® scanner equipment with patients: • You should be trained to use the system or have read and understood all sections of this guide describing correct operation. • You should also be thoroughly familiar with the safe operation of the TRIOS® scanner equipment as described in this documentation.
---	---

In case of equipment failure

	WARNING If at any time the TRIOS® system malfunctions, or if you suspect in any way that the TRIOS® system is not working correctly: • Remove the scanner from contact with the patient. • Unplug the system and make sure it cannot be used before it is checked (marked). • Contact your reseller. • DO NOT attempt to open any covers on the TRIOS® system.
---	--

5.2 Mechanical Hazards

Fragile and movable parts

	WARNING Be careful to avoid the following potential sources of injury: • Do not lean or sit on any part of the system. • Movable parts can pinch body parts.
---	--

Dropped or damaged equipment

	WARNING If you drop a TRIOS® tip on the floor, you MUST dispose of it immediately and NOT use the same tip again for scanning. There is high risk that the mirror in the tip has become dislodged and can fall out.
---	--

	WARNING If the handheld scanner is dropped or bumped it should immediately be calibrated before further use. If calibration fails, please contact your technical service provider. See this guide's instructions on Calibrating the handheld scanner.
	WARNING If the handheld scanner is dropped examine the TRIOS® system for any mechanical damage. If a mechanical damage is observed or suspected treat this as an equipment failure (see this guide's instructions on handling equipment failures).

System S1A (cart) only

	WARNING • When moving the cart across the floor, only move the cart with both hands on the upper part of the cart. • Do not attempt to lift the cart alone. When moving the system over a door step or other unsmooth surface, get an additional person's assistance to ensure stability. • When moving the cart across the floor, maintain all its wheels in contact with the floor and ensure that they are not locked. • Always place the scanner in its intended mount on the cart when not in use and when moving the cart. • Keep your feet clear of the cart's base and wheels when moving the cart. • Never drag the cart across any of its cables. • Do not place the cart on a slope greater than 10 degrees. • Keep sharp objects away from the touch screen.
---	---

System S1P (pod) only

	WARNING • Always place the scanner in its intended mount on the pod when not in use. • Do not place the pod on a slope greater than 5 degrees. • Place all cables carefully so that patient and/or user do not get caught in them and accidentally pull any part of the system. Pulling cables may damage the system, interrupt the data connection, and/or drag the system off its supporting surface.
---	---

5.3 Explosion Hazards

Environment

	WARNING The TRIOS® system is not designed to be used in environments that are potentially explosive such as in close proximity to flammable liquids or gases or in oxygen-enriched atmospheres.
---	--

Distance to the patient

	WARNING There is a potential explosion hazard if the TRIOS® system is operated in the presence of flammable anesthetic.
---	--

5.4 Electrical Safety

Ensure a grounded/earth connection

	WARNING To avoid the risk of electric shock, this equipment must only be connected to a supply mains with grounding/protective earth. The system is equipped with a three pin (US: prongs) grounding/earth type plug. If you cannot insert the plug into your outlet, contact your electrician to replace the plug or outlet to ensure a grounded power connection. Do not try to defeat or ignore the purpose of the grounding-type plug.
---	---

Electrical shock

	WARNING There is a risk of electrical shock if you attempt to access the inside of any part of the system. Only authorized and qualified service personnel may access the inside of any part of the system.
---	--

Do not use a power strip

	WARNING Do not connect the TRIOS® to an ordinary power strip or extension cord because these may not be secure in relation to grounding. If the ground connection fails, these can cause the following hazards and potential injury: • The total leakage current for all the connected equipment can exceed the limits specified in EN/IEC 60601-1. • The impedance of the ground connection could exceed the limits specified in EN/IEC 60601-1.
---	---

Stress on cables

	CAUTION All externally connected cables such as power cable, network cable, and the hand-held scanner cable connection, must never be subjected to pulling stress.
---	---

Spilled liquids

	WARNING Do not bring liquids such as beverages near the TRIOS® system. Do not spill liquids on the TRIOS® system.
---	--

Condensation

	WARNING Variations in temperature or humidity may cause water to condense inside the system which can cause it damage. Always let the system reach room temperature before you plug it in to the power source. If the system has been subject to major temperature or humidity changes, wait at least 2 hours before plugging it in to the power source. If there is visible evidence of condensation, wait at least 8 hours before plugging it in to the power source.
---	--

Disconnecting from mains

	WARNING There is no power ON/OFF switch on the system therefore the only reliable means to disconnect the system from mains is to unplug the power cord. Do not position the system so that it is difficult to unplug the power cord. Please allow enough time for the system to shut down properly before disconnecting the power cord. System S1P (pod) only: You must also disconnect the PC from mains.
---	---

PC connection options and requirements

System S1P (pod) only

	WARNING There are two allowed ways to connect the PC: (1) The PC is standard IT equipment and is located at least 1.5 m away from the scanner. You must not touch the PC at the same time! Y-cables connected with a cable patient! Peripherals using a wireless mouse may be used. (2) The PC is a medical-grade PC. See the PC's instructions for use!
	WARNING Only use the power supply that is described in the PC's instructions. OTHER POWER SUPPLIES MAY DAMAGE THE TRIOS® SYSTEM!

USB connection options and requirements

System S1A (cart) only

	WARNING Only connect bus-powered USB devices that support USB Full speed.
---	--

5.5 Eye safety

	WARNING The scanner emits bright light from the scanner. The scanner complies with EN 18154:2007 (IEC 60825-1:2007) safety of lamps and lamp systems. Exercise caution when handling the scanner. A brief glimpse of the light into the eyes may cause temporary discomfort and do not stare into beam or view directly into eyes. In general, intense light source exposure potential due to their brightness when viewing other bright light sources. Temporary reduction in visual acuity, leading to irritation, annoyance, even accidents, depending on the exposure.
---	---

5.6 Pacemaker and ICD Interference

	WARNING Do not use on patients with pacemakers or ICDs. Due to interference risks, it is prohibited to use the system on patients that use a pacemaker or ICD.
---	---

5.7 Hygiene

	WARNING To maintain a clean work environment, patient, wear clean surgical gloves: • Handling the scanner tip – carefully flipping for upper/lower scans or • Handling the protection tip. • Handling the color calibration kit. • Applying the TRIOS® for scanning. • Touching the TRIOS® system.
---	---

WARNING

The scanner unit and its optical window must always be cleaned before attaching a Protection Tip. Do not allow patient contact with the Protection Tip.

WARNING

Always have the scanner-tip mounted when inserting the scanner into the patient's mouth or in any way in contact with the patient. Mounting of the scanner tip is always required for maximal electrical safety and due to the risk of contamination.

WARNING

Always have a cleaned and sterilized scanner tip mounted before attaching the color calibration adaptor. Do not allow patient contact with the color calibration kit.

WARNING

Before each patient, you must:

- disinfect the system
- apply a new sterilized or high-level disinfected tip to the scanner

See the detailed instruction in the chapter below.

S1P (pod) only

WARNING

Before each patient, make sure any PC connected to the system and its peripherals are disinfected. It is recommended to use medical-grade peripherals (e.g., keyboards and mice) that can be disinfected easily.

Caution when scanning with TRIOS®

CAUTION

The handheld scanner is a highly precise optical instrument and must be handled carefully. Do not drop, shake or bump the handheld scanner as these actions could cause it to break or malfunction. Follow these guidelines:

- Always return the handheld scanner to its mount when not in use.
- Do not bend or kink the cord attached to the handheld scanner.
- Do not submerge the handheld scanner in any liquids.
- Do not place the handheld scanner on heated or wet surfaces.
- Grip the handheld scanner firmly during scanning and when removing/returning to its mount.

NOTICE

The scanner will get warm during use. This is normal but not part of any patient treatment.

CAUTION

The handheld scanner tip is fragile and contains a delicate mirror. To avoid damage to the tip and mirror, take care not to bump the tip onto the patient's teeth or restorations.

5.9 Protection from overheating

System S1A (cart) only

WARNING

Ventilation openings must never be blocked or obstructed. There are openings at the back of the cart (behind the touch screen), and under the base of the cart. If overheated, the system may stop functioning.

6. Cleaning, Disinfection, and Sterilization

6.1 Cleaning and sterilizing the scanner tip

The TRIOS® scanner tip is the part of the scanner that is inserted into the patient's mouth during scanning. The tip can be reused a limited number of times, but it must be cleaned and sterilized between patients to avoid cross contamination.

The tip can be removed by pulling it off with two fingers.

Cleaning and sterilization procedure

Step 1: Clean the tip manually using soap water and a brush immediately after use. Inspect the mirror of the tip after cleaning. If the mirror appears with stains, smudged or with a milky haze, repeat cleaning using a soft dish-brush and soap water. Rinse with water. Dry the mirror carefully with a paper towel.

Step 2: Put the tip in a paper sterilization pouch

The pouch should be sealed airtight. Use either a self-adhesive pouch or heat-sealed pouch.

Step 3: Sterilize the wrapped tip in an autoclave with one of the two programs:

- At 134°C for a minimum of 4 minutes.
- At 121°C for a minimum of 45 minutes.

Use an autoclave program that dries the wrapped tip before opening the autoclave.

A number of wrapped and sterilized tips can be kept at a suitable place ready for use at your convenience.

Sterilization of the scanner tip is always desirable, but if sterilization is not practicable a high-level disinfection with a high-level disinfectant may be used (i.e. Cidex OPA Solution, 0.55% ortho-phthaldehyde, which has been validated by 3Shape and is on the list of FDA approved high-level disinfectants).

Cleaning and high-level disinfection procedure

Step 1: Clean the tip manually using soap water and a brush immediately after use. Inspect the mirror of the tip after cleaning. If the mirror appears with stains, smudged or with a milky haze, repeat cleaning using a soft dish-brush and soap water. Rinse with water. Dry the mirror carefully with a paper towel.

Step 2: Disinfect the tip using Cidex OPA Solution

Immerse the tip completely in Cidex OPA Solution (0.55% ortho-phthaldehyde) for a minimum of 20 minutes. See also the Cidex OPA Solution Instructions for Use.

Step 3: Rinse the tip

After having been immersed in the Cidex OPA Solution for min. 20 minutes remove the tip and rinse it according to the Cidex OPA Solution Instructions for Use.

Step 4: Dry the tip

Use a sterile and non-abrasive paper towel or cloth to gently remove remaining drops on the mirror and to dry the tip.

The mirror can be extracted from the scanner tip and later reinserted. By doing so the scanner tip without the mirror can be sterilized and the mirror can be high-level disinfected using the two procedures described above.

Cleaning and hybrid sterilization / high-level disinfection procedure

Step 1: Clean the tip with the mirror as described above (Step 1).

Step 2: Carefully extract the mirror from the scanner tip

Be careful not to scratch or damage the mirror while extracting it.

Step 3: Sterilize the scanner tip without the mirror using the Cleaning and sterilization procedure described above (Steps 2-3).

Step 4: High-level disinfect the mirror using the Cleaning and high-level disinfection procedure described above (for Steps 2-4 only immerse, rinse and dry the mirror).

Step 5: Insert the mirror in the scanner tip. Carefully insert the mirror in the scanner tip ensuring that it is placed correctly and magnetically fixated.

CAUTION

The mirror in the tip is a delicate optical component. Its clean and undamaged surface is critical for the scan quality. The mirror should be handled with care. Do not scratch it. Be careful not to smudge it before scanning a patient.

CAUTION

Never autoclave a tip that is not wrapped, since this will leave stains on the mirror that cannot be removed. See the autoclave manual for more information.

CAUTION

New tips must be cleaned and sterilized / autoclaved before using them for the first time.

CAUTION

A class B Vacuum autoclave according to EN13060 is required for sterilization of the scanner tip.

CAUTION

The scanner tip can be re-sterilized up to 20 times and must be disposed afterwards as described in section Disposal.

CAUTION

Always clean the tip immediately after use. If biological materials (bioburden) are allowed to dry on equipment, disinfection and sterilization processes may not be effective.

CAUTION

Be sure to reinsert the mirror in the scanner tip so that it is placed correctly and magnetically fixated.

IMPORTANT NOTICE

High-level disinfection by 0.55% Cidex OPA Solution as an alternative to sterilization has been validated by 3Shape. The Cidex OPA Solution is not allowed and/or recognized as a disinfectant of medical devices worldwide, always follow procedures that have been established for your clinic or institute, as well as any national guidelines.

IMPORTANT NOTICE

Always follow cleaning, disinfection and sterilization procedures that have been established for your clinic or institute, as well as national guidelines.

IMPORTANT NOTICE

If using Cidex OPA for high-level disinfection be sure to follow the Cidex OPA Solution Instructions for Use.

6.2 Cleaning, disinfecting and sterilizing the Color Calibration Kit

The TRIOS® color calibration adaptor must always be attached on a cleaned and sterilized scanner tip and must not have contact to patients. If the instructions are not followed the color calibration kit must be cleaned, disinfected and sterilized.

Cleaning, disinfection and sterilization procedure

Step 1: Disassemble the color calibration kit (remove the sleeve from color calibration target and disassemble the color calibration target from the color calibration adaptor).

Step 2: Clean the color calibration adaptor manually using soap water and a brush. Sterilize the cleaned color calibration adaptor in an autoclave using one of the programs listed for the scanner tip and using a paper sterilization pouch.

Step 3: Clean the sleeve manually using soap water and a brush. Disinfect the sleeve using the recommended disinfection solution in section Cleaning and disinfecting other parts of the system.

Step 4: The color calibration target must only be cleaned with air. Calibration patches must not be touched and must not be exposed to liquids. See the important notice below.

After cleaning, disinfection, and sterilization the color calibration kit is reassembled before use. Ensure that all parts are dry before reassembly. Immediately place the color calibration target in the sleeve and color calibration adaptor in the sterilization pouch.

IMPORTANT NOTICE

The color calibration target is a delicate optical component. During disassembly and assembly of the color calibration kit, only the end-part of the color calibration target (opposite of the calibration patches) must be touched, and the target must only be cleaned with air.

6.3 Cleaning and disinfecting other parts of the system

After patient treatment, clean and disinfect all other surfaces of the system: the scanner tip and color calibration kit, for which the previous sections describe the procedures. Also note that the optical window is a delicate optical component and should be handled with extra care during cleaning and disinfection (see below).

For general cleaning, a soft cloth and a mild detergent can be used.

Recommended disinfection solution:

Denatured alcohol (aka. ethyl alcohol or ethanol) – typically 60-70% / v/v

To disinfect:

Step 1: Apply the disinfection solution to a soft lint-free non-abrasive cloth, the cloth must only be damp.

Step 2: Wipe the surfaces with the damp cloth.

Step 3 (optical window on scanner only): Dry the optical window with another clean, soft, lint-free, non-abrasive cloth.

	CAUTION The scanner's optical window is a delicate optical component. Its clean and undamaged surface is critical for the scan quality. The window should be handled with care. Make sure it is free from lint, stains, and other kind of dirt.
	CAUTION Never submerge the TRIOS® scanner into the cleaning solution.
	CAUTION Do not use ammonia-based or chloride based solutions or acetone on any surface.
	CAUTION (System S1A only) Do not spray directly onto the monitor as the liquid solution can accumulate at the seals and leak into the monitor.
	NOTICE The scanner shell might get yellow-toned over a long period of use.
	IMPORTANT NOTICE All surfaces that the protection and/or calibration tip comes in contact with, particularly the scanner tube and its optical window, must always be cleaned and disinfected before attaching the protection tip / the calibration tip.

6.4 Additional cleaning if you see poor scan quality

During scanning the tip mirror may become unclean, resulting in a loss of scan quality. In this situation, you can clean the mirror following these steps:

Step 1: Remove the scanner tip.

Step 2: Apply alcohol to a clean cloth or cotton-tipped swab and wipe the tip mirror until it is clean. Use alcohol that is free of impurities, as those may stain the mirror. Ethanol or propanol (ethyl-/propyl alcohol) can be used.

Step 3: Wipe the tip mirror dry with a dry, lint-free cloth to avoid leaving dust and fiber on the mirror.

	CAUTION Do not use acetone or any oxidizing solutions to clean the tip mirror.
---	--

If you still see problems, check the scanner's optical window for cleanliness, and apply the same procedure if necessary.

7. Disposal

7.1 Disposal of the scanner tip

1. Sterilize the tip as described above.
2. The scanner tip can then be disposed normally as other clinical waste.

	CAUTION Dispose the scanner tip as other clinical waste.
	CAUTION The scanner tip must be sterilized before disposal.

7.2 Disposal of the TRIOS® System

For disposal of the TRIOS® System please contact your 3Shape reseller who can inform you about the local requirements for disposal of electronic waste.

8. Technical Maintenance

	CAUTION All technical service and repair of the TRIOS® System must be performed by 3Shape or 3Shape-authorized personnel only.
---	--

The user need not perform any regular maintenance other than calibration and cleaning, disinfection and sterilization as described here. Preventive inspections and scheduled maintenance are not required.

System S1A (cart) only

Changing the fuses (Figure 6):

Shut down the TRIOS® cart by following the procedure described in "Setting up & Starting TRIOS® for the First Time" and disconnect the mains power cable from the TRIOS® cart.

	WARNING The mains power cable must be disconnected before changing fuses.
---	---

Using a pair of tweezers or similar, push in the two tabs on each side of the fuse drawer simultaneously as shown in the picture, and the drawer will begin to emerge from its slot.

It is then possible to remove the drawer and replace the two fuses. Two spare fuses are provided with the system on delivery.

	WARNING Only fuses with the correct rating and type (250V, T3.15A H) must be used.
---	--

To replace the fuses, first mount them in the drawer, then insert the drawer back into the slot and push all the way in until the tabs click back into place and the drawer is no longer loose in the slot. Do not attempt to insert the fuses directly into the slot on the cart.

The mains cable can now be reconnected and the TRIOS® system started up.

Changing the dongle (software license protection (Figure 7):

The dongle is located behind the ventilation grid. You can remove the grid by opening the screws holding the grid, as indicated by the arrows in Figure 7.

	WARNING The scanner must not be operated while the ventilation grid on the back side of the cart is removed.
---	--

9. Environmental Considerations

3Shape recognizes our common responsibility for the Earth. Thus the TRIOS® system is designed to conform with the following directives.

RoHS, Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment (2002/95/EC).

WEEE, Waste Electrical and Electronic Equipment Directive (2002/96/EC)














10. Consumable Accessories

The only consumable accessories to the TRIOS® System are the scanner tips and the color calibration kit.

When the scanner tips are disposed new tips can be ordered by contacting your 3Shape reseller. Also, a new color calibration kit can be ordered from your 3Shape reseller.

11. Complaints


On the territory of the Russian Federation consideration of complaints and appeals of consumers takes the company "BALTEKS", St. Petersburg, Baltiyskaya 31A, tel. 8 (812) 252-56-80, e-mail: spb@baltteks.ru

Technology designed the way you work

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The screenshot shows the 'Patients' interface. At the top, there are buttons for 'New session', 'Add patient', 'Edit patient', and 'Delete sessions'. Below these is a search bar labeled 'Search by name or number'. The main area is divided into two columns. The left column lists patients, with 'Hagen, Hans' selected. The right column shows the details for a specific order, including the order ID, lab name, and a list of restorations. Below the restorations, there are three 3D scans of the teeth, labeled 'Lower', 'Upper', and 'Bite'.

Restorations	Type	Material	Shade
21	Crown zircon	Vita A1	

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

A screenshot of the zshape lab software interface. It features a top navigation bar with tabs: General (selected), Details, Address, and Images. Callout 7 points to the 'zshape lab' logo. Callout 4 points to the 'lab' text. Callout 5 points to the 'Change lab' button. Callout 2 points to the 'Tools' button. The main area displays patient information: Name: Hans Hagen, Patient ID: 12345, Date of birth: 02/02/1980. Callout 3 points to the 'Change patient' button. Delivery date: 25/01/2016 14:30 and Order ID: 61158160125142737 are shown on the right. Below this is a section for 'Enable extra scans' with a 'Pre-preparation' button. A dental arch diagram shows teeth numbered 1 to 28. Callout 6 points to tooth 21. Below the arch is an 'Open shade tool' button and a 'Clear' button. At the bottom, there's an 'Indication' section with tabs: Implant, Post & Core, Implant planning, and Orthodontics. Below these are dropdown menus for 'Type' (Crown), 'Material' (Emax), and 'Shade' (Vita A1). Callout 1 points to the 'Shade' dropdown.

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The screenshot displays the 'zshape DC Lab' software interface. At the top, there's a header with the 'zshape' logo and 'DC Lab' text. Below this, there are tabs for 'General', 'Details', 'Address', and 'Images'. The 'General' tab is active, showing patient details: Name: Hans Hagen, Patient ID: 1234-123, Date of birth: 2/2/1980. There's a 'Change patient' button. To the right, it shows 'Delivery date: 6/19/2015 12:45 AM' and 'Order ID: 97728150617124030'. Below the patient info, there's a section for 'Enable extra scans' with two checkboxes: 'Pre-preparation' and 'Emergence profile'. The main part of the interface is a dental arch diagram with teeth numbered 1-28. Tooth 11 is highlighted with a scanner icon. Below the arch, there are tabs for 'Indication' (Implant selected), 'Post & Core', and 'Abutment'. At the bottom, there are dropdown menus for 'Abutment', 'Manufacturer' (Not chosen), and 'Material' (Not chosen).

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Connection

Client sharing
Wi-Fi
Proxy Settings
Diagnose Internet
Remote screen/tablet scanning
Bluetooth Devices

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External Connection Method

Communicate

3shape	Local Lab	
3shape	Local Lab2	
	ALG_Lab	
	Alla's lab	
3shape	Daria Soldatenko	
	Evers Lab	

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Order	
Dentist address	
Tooth numbering system	FDI
Shade system	Vita Classic
Remember last indication settings	☐
Calendar	
Time interval between orders (in minutes)	60
Workday start	7:00 AM
Workday end	7:00 PM
Patient names privacy	Show toggle button

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Bite alignment operation Live ▾

Island cleanup level for monochrome scanners 2014 (default) ▾

Default scan zoom level small large

Non-preparation post processing Normal ▾

Insane Mode (fast Scan Engine) ☐

Hole close color ☒ Green ▾

3D calibrate scanner 3D calibrate scanner

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Customizing your TRIOS

TRIOS setup wizard

Welcome to the TRIOS setup wizard, where you are guided to customise your TRIOS to best suit your preferences.

The wizard will guide you through the following settings:


- **Language**
- **Lab connections** (where to send orders)
- **Dentist details** (name, email, login, password)
- **Order details** (tooth numbering and color system)

You can **Exit** this wizard any time, and then default settings will be used.

You can also be reminded **Later** to complete the wizard.

All of the settings can later be changed by going into **Configure** from the home screen.

Please click **Next** to begin customization of your TRIOS.

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
Register a Dental Clinic

Please fill in all the fields and press Register

Clinic Name:

Address:

Postal Code:

Country:

Email:

Password:

Confirm Password:

☒ I have a 3Shape Dongle

3Shape Dongle number:

☒ Receive newsletter and information?

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Register Clinic Summary

Thank you for registering 'test' at 3Shape Communicate

A confirmation email will be sent to abc_dental_clinic@abc.de. You will have to press the Activation link in the email in order to activate the account.

If you have acquired a TRIOS scanner you will need some of the following info when setting up TRIOS to connect with 3Shape Communicate:

Login: abc_dental_clinic@abc.de

Password: Look in the mail sent to you

Server URL: net.tcp://www.3shapecommunicate.com:9200/3shapecommunicateservice/_30Service.svc

Website URL: https://www.3shapecommunicate.com

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
Configure Communicate

Login (email)

Password
(If different from Operator password)

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3Shape Communicate™
shkina@3Shape.com
Log Off
English

3shape
Welcome
Cases
Settings
Connections
Order templates
3shapedental.com

Cases

Filter

OrderId

Email

Created from

Order state
All

ThreeShapeOrderNo

Patient name

Created to

Get

Alla clinic

Select All
Unselect All
Delete selected

	Patient Name	Order No.	Delivery Date		
☐	Tom	08971_140311163000	-	Preview	Open Case
☐	Tom	08971_140311144517	-	Preview	Open Case
☐	Qwerty	08995_140205100826	-	Preview	Open Case

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
Welcome Cases Settings Connections Order

Connect to other Dental Clinics and Labs

Fill in email of the Clinic or Lab you want to connect to

Result:

Alek Lab Site
Alek Dentist Site
Alek Trios-Clinic
Alek's Lab

Message to Alek Lab Site (You can customize the message)

We would like to connect to you on 3Shape Communicate.

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
6.2 Cleaning, Disinfection and Sterilization

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
Information

A new scanner firmware is available, do you wish to upgrade?
The upgrade will take about 5 minutes.

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

PRODUCT DESCRIPTION— TRIOS® 3 PRODUCT DESCRIPTION

Table of Contents

1	Introduction	3
2	Basic Information	3
	2.1 Legal manufacturer and manufacturing site.....	3
	2.2 Structure and composition	3
	2.3 Sales configurations	4
	2.5 Components and accessories of S1A and S1P	5
	2.6 Components and accessories of S1A ONLY.....	6
	2.7 Components and accessories of S1P ONLY.....	6
	2.8 Packing instruction.....	7
3	Software Description	8
	3.1 Software classification	8
	3.2 Software specification.....	8
	3.2.1 TRIOS® Application SW	8
	3.2.1.1 Product ID:	8
	3.2.1.2 Architecture	9
	3.2.2 ScanSuiteTRIOS® including firmware.....	10
	3.2.2.1 Product ID	10
	3.3 Architecture Design Chart.....	10
4	Functional principle.....	11
	3.1 Scanner light.....	11
	3.2 3D Imaging principle	11
	3.2.1 Components of the scanner	11

3.2.2 How to realize 3D imaging from 2D images.....	11
3.2.3 Transforming correlation measure extrema to 3D world coordinates.....	12
3.2.4 How imaging through different spatial relations of the scanner and the object is achieved	12
3.2.5 Method of guaranteeing the imaging precision	13

1 Introduction

The TRIOS® 3 intraoral scanner system enables trained dental professionals to obtain 3D images of the topographical characteristics of their patient's teeth directly from their mouths. Moving the scanner tip in close proximity to the teeth, a digital image of the teeth is relayed to the software and a 3D model is created in the computer. The 3D model is stored in the computer along with the patient's case information.

This digital procedure replaces the manual procedure of impression and gypsum modelling in the dental clinics and laboratories.

Intended Use; The TRIOS® 3 intraoral scanner system is intended to record 3D optical digital images of the topographical characteristics of teeth directly for restoration and implant treatment as well as orthodontics.

Indications for Use; TRIOS® 3 intraoral scanner system supports the following indications:

- Up to 5-unit bridges
- Inlays
- Onlays
- Post and core
- Single crowns
- Veneers

User Profile; The TRIOS® 3 system may be operated only by trained dental professionals and qualified personnel.

2 Basic Information

2.1 Legal manufacturer and manufacturing site

Applicant name/Legal manufacturer Address	3Shape TRIOS A/S Holmens Kanal 7, 1060 Copenhagen, Denmark
Manufacturing site Production address	3Shape Poland Sp. z o.o. UL.Zapadla 8d, Szczecin, 70-033, Poland

2.2 Structure and composition

Product name	Intraoral scanner system
Trade/Sales name	TRIOS® 3 Intraoral scanner system
Product model(s)	S1A & S1P

Software included in the system

TRIOS® Software Application (*Scan application.
User interface/application for handheld scanner*)

ScanSuiteTRIOS® (*Interface between application software
and firmware in the handheld scanner*)

TRIOS® Firmware (*Scanner firmware*)

Product component(s)

Hand-held scanner (2 different options refer to section 2.5)

Scanner tip

Protection tip

Calibration tip

Color calibration kit

Power (mains) cable and plug

License USB Dongle

S1A furthermore comprises

System supplied with a cart with build-in PC

Cart with;

-scanner mount

-wheels

-computer and pre-installed software (build-in cart)

-touchscreen

-system standby button

-scanner connection socket

Back of cart consists of;

-power socket

-network (LAN) connector

-fuse drawer

-brake

-USB port

S1P furthermore comprises

*System supplied with a "pod" (holder), to be
connected to an external PC via an Ethernet
cable*

Pod with;

-scanner mount

-single connection socket

-connection socket for Ethernet cable

-connection socket for external medical power supply

External medical-grade power supply

USB memory stick with software installer

Extension cable 8m

2.3 Sales configurations

Model (option)	Handheld scanner (handle)	Handheld scanner (pen grip)	Cart	Pod (handle)	Pod (pen grip)	PC (optional)	PSU	Extension cable (8m)
S1P (1)	x			x		(x)	x	(x)

S1P (2)		x			x	(x)	x	(x)
S1A (1)	x		x					
S1A (2)		x	x					

2.5 Components and accessories of S1A and S1P



Figure 1 Handheld TRIOS® intraoral scanner (with handle)
Scanner obtains 3D images of the topographical characteristics of patient's teeth. Patient contact with distal end of tip. Weight = 350 g

OR

Figure 2 – Handheld TRIOS® intraoral scanner (pen grip)
Scanner obtains 3D images of the topographical characteristics of patient's teeth. Patient contact with distal end of tip. Weight = 325 g

Figure 3 – Detachable scanner tip
The detachable tip is a type BF applied part according to IEC 60601-1 and a multi-use device. It is supplied non-sterile to end-customer and must be sterilized or HL disinfected in between patients to avoid cross contamination.

Figure 4 – Protection tip
When scanner is not in use the protection tip must be mounted on the scanner.


Figure 5 – Calibration tip
The calibration tip must be used when calibrating the scanner.


Figure 6 – Color calibration kit
(left to right: Color calibration adaptor, color calibration target, sleeve)
The color calibration kit is applied when calibrating the color of the scanner.

2.6 Components and accessories of S1A ONLY


Figure 7 – Cart

The cart has a build-in PC with pre-installed TRIOS® Application Software and a scanner holder. The cart has wheels and is thereby mobile.

2.7 Components and accessories of S1P ONLY


Figure 8 - Pod setup with external PC

The system setup when using pod requires a PC, external power supply unit, USB with software installer


Figure 9 – Pod without holder (for scanner with handle)
Pod device with an accompanying small USB adapter for connection to an external PC. Weight = 645 g


Figure 10 – Pod with (for pen grip scanner)
Pod device with an accompanying small USB adapter for connection to an external PC. Weight = 725 g


Figure 11 - Medical grade PSU

2.8 Packing instruction

The packing of the models are as follows:

Software included in all models	TRIOS® Software Application ScanSuiteTRIOS® (embedded in scanner) TRIOS® Firmware (embedded in scanner)
Product component(s) included in all models	1) Handheld scanner with pistol grip or 2) handheld scanner with pen grip Scanner cable (2m) Scanner tip Protection tip Calibration tip Color calibration kit Power cable and plug License USB Dongle Safety and Setup Guide TRIOS® Technical Data S1AP
S1A furthermore includes Cart supplied with build-in PC	Cart Ethernet cable

The model is packed in a wooden box with foam. The scanner is also packed in a Scanner Carry Box with bottom foam and top foam to protect the equipment during transport.

S1P furthermore includes
Pod (holder) for scanner to be connected to an external PC via an Ethernet cable

1) Pod for hand-held scanner or 2) Pod for pen-grip scanner
Pod Power supply
USB memory stick with software installer
Ethernet cable 5m
USB3 to GbE Adapter
TRIOS® Pod Installation Guide
Extension cable 8m (optional)

The model is packed in a cardboard shipping box and a Scanner-Pod Carry Box with bottom foam and top foam to protect the equipment during transport.

3 Software Description

3.1 Software classification

Product Safety level: Class A

Based on software risk analysis and control the software is considered a safety class A¹. A software failure is not considered to result in injury or damage to end-user or patient.

3.2 Software specification

S1AP registration contains the following software products:

- TRIOS® Application SW
- ScanSuiteTRIOS®
 - TRIOS® Firmware (*Firmware is a part of ScanSuiteTRIOS®*)

3.2.1 TRIOS® Application SW

The user can use the application SW to control the handheld scanner. The control and feedback to the scanner is conducted through the ScanSuiteTRIOS® interface.

3.2.1.1 Product ID:

Software name : TRIOS® Application SW

¹ CEI/IEC 62304:2006

Model No. and version No. : Version number 1.3.4.5
 Manufacturer name : 3Shape TRIOS A/S
 Manufacturer address: Holmens Kanal 7, Copenhagen 1060 Denmark.
 Production address : ul. Zapadla 8d Szczecin 70-033, Poland

3.2.1.2 Architecture

The structure and relationship of the software modules:


Figure 12 – Architecture of TRIOS® Application SW

The TRIOS® Application software is structured into 3 main layers

1. Database Layer
2. Framework
3. Modules

The database is used to store orders on the TRIOS® system PC. It is being used by the framework which facilitates communication between the modules and the database. The framework also handles the support for both 3Shape and third party modules, along with configuration and user account control.

Software module and module function :

Software Module name	Module function
Home Screen	Start page for the SW
Order tracker	Manage created orders
Order Entry	Create order
Scanning	Scan

Analysis	Analysis of scanned data prior to sending
Sending	Sending order incl. scan data


3.2.2 ScanSuiteTRIOS® including firmware

3.2.2.1 Product ID

Software name:	ScanSuiteTRIOS®
Model No. and version No.:	Version number 1.4.2.0
TRIOS®3 Firmware:	Version number 44.02.00.01
Manufacturer name :	3Shape TRIOS A/S
Manufacturer address:	Holmens Kanal 7, Copenhagen 1060 Denmark.
Production address :	ul. Zapadla 8d Szczecin 70-033, Poland

3.3 Architecture Design Chart

The relationships among software and hardware units including dataflows is depicted in the figure – Hardware/Software below:

4 Functional principle

3.1 Scanner light

The scanner light is generated by LED diode and the LED class of scanner is "moderate risk" according to IEC 62471, hence a class 1 and/or 1M according to IEC 60825-1.

3.2 3D Imaging principle

3.2.1 COMPONENTS OF THE SCANNER

The scanner head comprises a tip which can be entered into a cavity, a light source, optics to collect the light from the light source, pattern generation means, a beam splitter, an image sensor and electronics, a lens system which transmits and images the light between the pattern, the object being scanned, and the image sensor (camera).

The light from the light source travels back and forth through the optical system. During this passage the optical system images the pattern onto the object being scanned and further images the object being scanned onto the image sensor. The lens system includes a focusing element adjusted to shift the focal imaging plane of the pattern on the probed object.

3.2.2 HOW TO REALIZE 3D IMAGING FROM 2D IMAGES

A 3D image (scan) is obtained by collecting a number of 2D images at different positions of the focus plane. These different 2D images are recorded while the focus element is adjusted in such a way that the image of the pattern on the scanned object is shifted along the optical axis. The in-focus regions on the object display a spatially varying pattern, while the out-of-focus regions display smaller or no contrast in the light variation. By using knowledge of the optical system used, it is possible to transform the contrast information vs. position of the focus plane into 3D surface information, on an individual pixel group basis.

Determining in-focus information relates to calculating a correlation measure of this spatially structured light signal (which we term input signal) with the variation of the pattern itself (which we term reference signal). In general, the magnitude of the correlation measure is high if the input signal coincides with the reference signal. If the input signal displays little or no variation then the magnitude of the correlation measure is low.

The 3D surface structure of the probed object is determined by finding the plane corresponding to an extremum in the correlation measure for each group of sensor element in the camera's array of sensor elements when recording the correlation measure for a range of different focus positions.

The correlation measure is defined mathematically as a dot product computed from a signal vector, $I = (I_1, \dots, I_n)$, with $n > 1$ elements representing sensor signals and a reference vector, $f = (f_1, \dots, f_n)$, of same length as said signal vector of reference weights. The correlation measure A is then given by

$$A = f \cdot I = \sum_{i=1}^n f_i I_i$$

The indices on the elements in the signal vector represent sensor signals that are recorded at different sensors elements.

3.2.3 TRANSFORMING CORRELATION MEASURE EXTREMA TO 3D WORLD COORDINATES

Relating identified focus position for pixel groups of a stack of 2D images to 3D world coordinates is done by ray tracing through the optical system. Before such ray tracing can be performed the parameters of the optical system must be known. A calibration step provides such knowledge. The calibration is performed by recording images of an object of known geometry for a plurality of focus positions. The scanner is then calibrated by generating simulated ray traced images of the calibration object and then adjusting optical system parameters as to minimize the difference between the simulated and recorded images.

With knowledge of the parameters of the optical system backward ray tracing technique is employed to estimate the 2D -> 3D mapping. The following steps have been performed:

1. From each pixel of the 2D image (at the image sensor), a certain number of rays is traced, starting from the image sensor and through the optical system (backward ray tracing).
2. From the rays that emit, the focus point is calculated, i.e. the point where all these rays substantially intersect. This point represents the 3D coordinate of where a 2D pixel will be in focus.
3. Generating a look up table for all the pixels with their corresponding 3D coordinates.

The above steps are repeated for a number of different focus lens positions covering the scanner's operation range.

3.2.4 HOW IMAGING THROUGH DIFFERENT SPATIAL RELATIONS OF THE SCANNER AND THE OBJECT IS ACHIEVED

The 3D image (scan) is obtained from a number of 2D images as described above. Such a 3D image is called a sub-scan. A sub-scan covers a limited area. If one needs to have a 3D model of an object larger than what is covered with a single sub-scan then it is necessary to move the scanner to a new position and obtain a sub-scan with a different spatial relation of the scanner and the object. The subscans are fused (with a weighted average) to a single reconstruction model, the voxel model. For each incoming subscan a voxel-subscan correspondence is calculated and used for updating the reconstruction model if

the object is even larger than the combined 3D image of two subscans then the process can be continued and a new sub-scan partially overlapping with the existing scanned area can be fused to the reconstruction model.. In this way scans of objects that are much larger than the individual sub-scans can be obtained.

3.2.5 METHOD OF GUARANTEEING THE IMAGING PRECISION

Several processes are done to ensure the imaging precision. The first step to obtain imaging precision is to record images of an object with fiducials of known geometry. The object being used for this is a plate with a checkerboard pattern. The periodicity of the pattern is known and the fiducials are the corners of the individual checkerboard squares. A set of images with different focus positions is recorded of this checkerboard reference object. In this way, a planar array of fiducial points inside the scanning volume is obtained. For an identified in-focus fiducial point, a certain number of rays are traced, starting from the image sensor and through the optical system (backward ray tracing). From the rays that emit, the focus point is calculated, i.e. the point where all these rays substantially intersect. This point represents the 3D coordinate of where a 2D pixel will be in focus. It then becomes possible to generate a look up table for all the pixels with their corresponding 3D coordinates.

The step where images of an object with fiducials of a known geometry are recorded is repeated for a number of different positions of the object inside the scanner's operation range. This allows making a lookup table for all pixels and for all focus positions to their corresponding 3D coordinates. This lookup table represents calibration data for a given scanner. When an unknown object is scanned, the scanner identifies in-focus coordinates for certain focus positions. Using the calibration data, it becomes possible to assign these in-focus coordinates to 3D surface points.

The imaging precision is tested by scanning a sphere object with a well-known diameter of 15 mm. The scanning results in a number of points corresponding to the sphere surface. A mathematical sphere is then fitted to this recorded point cloud. The diameter of the fitted sphere should be within 20µm of the real diameter for the scanner to have an acceptable precision.

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
TRIOS® Technical Data

3shape

10014018A December 2014 Copyright 3Shape TRIOS A/S

TRIOS® Intraoral Scanner System model S1A, S1P - Electromagnetic Compatibility

The TRIOS® Intraoral Scanner System (hereinafter referred to as the UNIT) complies with the requirements for electromagnetic compatibility (EMC) for medical electrical equipment according to IEC 60601-1-2. The UNIT has been tested according Medical electrical equipment – Part 1-2: General requirements for basic safety and essential performances – Collateral standard: Electromagnetic – Requirements and tests, EN/IEC 60601-1-2 edition 3.0, 2007.

WARNING: Medical electrical equipment needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in this document.

WARNING: The UNIT should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the UNIT should be observed to ensure normal operation in the used configuration.

IMPORTANT NOTICE: Portable and mobile RF communications equipment can affect medical electrical equipment.

Declaration - electromagnetic emissions

The UNIT complies with the requirements EN/IEC 60601-1-2 (Edition 3.0) clause 6.1 emission.

The UNIT is intended for use in the electromagnetic environment specified below. The customer or user of the UNIT should assure that it is used in such an environment.

Emission measurement	Conformity	Electromagnetic environment - guidelines
RF emission CISPR11	Group 1	The UNIT uses RF energy only for its internal functions. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emission CISPR11	Class B	
Harmonic emission IEC 6100-3-2	Class A	The UNIT system is intended for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings for used for domestic purposes
Voltage fluctuation and flicker IEC 61000-3-3	Complies	

Declaration - electromagnetic immunity, part 1

The UNIT complies with the requirements EN/IEC 60601-1-2 (Edition 3.0) clause 6.2 immunity. The test levels for immunity to RF electromagnetic fields are selected to "General" medical equipment.

The UNIT is intended for use in the electromagnetic environment specified below. The customer or the user of the UNIT should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidelines
Electrostatic discharge ESD IEC 61000-4-2	Contact +/-6KV. Air: +/-8KV.	Contact +/-6KV. Air: +/-8KV.	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic materials. The relative humidity should be at least 30%. During direct contact 4 kV electrostatic discharges, the UNIT is allowed to have a temporary loss of function, which requires operator intervention or system reset.

Electrical fast transient/burst IEC 61000-4-4	+/- 2 kV for power supply lines +/-1 kV for input/output lines	+/- 2 kV for power supply lines +/-1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment. During electrical fast transient/ burst the UNIT is allowed to have a temporary loss of function, which requires operator intervention or system reset.
Surge IEC 61000-4-5	+/-1 KV Line(s) to earth: +/-2 KV Line(s) to line(s).	+/-1 KV Line(s) to earth: +/-2 KV Line(s) to line(s).	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruption and voltage variations IEC 61000-4-11	<5% UT (>95% dip in UT) for 0.5 cycle 40% UT (60% dip in UT) for 5 cycles. 70% UT (30% dip in UT) for 25 cycles. <5% UT (>95% dip in UT) for 5 s.	<5% UT (>95% dip in UT) for 0.5 cycle. 40% UT (60% dip in UT) for 5 cycles. 70% UT (30% dip in UT) for 25 cycles. <5% UT (>95% dip in UT) for 5 s.	<5% UT (>95% dip in UT) for 0.5 cycle. 40% UT (60% dip in UT) for 5 cycles. 70% UT (30% dip in UT) for 25 cycles. <5% UT (>95% dip in UT) for 5 s.
Power frequency (50/60Hz) Magnet field IEC 61000-4-8	3 A/m.	3 A/m.	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical domestic or hospital environment.

Note: UT is the ac mains voltage prior to application of the test level

Declaration - electromagnetic immunity, part 2

The UNIT is intended for use in the electromagnetic environment specified below. The customer or the user of the UNIT should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidelines
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	[V1] 3 Vrms	$d = \sqrt{\frac{35}{P_1} N F}$ $d = \sqrt{\frac{35}{E_1} N F}$ 80 MHz to 800 MHz $d = \sqrt{\frac{7}{E_1} N F}$ 800 MHz to 2.5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meter (m), V1 is the conducted RF compliance level, and E1 is the RF radiated compliance level. Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, (note a and b) should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m 80 MHz to 2.5 GHz	

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, object and people.

Note a: field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitter, an electromagnetic site survey should be considered. If the measured field strength in the location in which the UNIT should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the UNIT.

Note b: Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Separation distances between the UNIT and RF communications equipment

The UNIT is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the UNIT can help prevent electromagnetic interference by maintaining a minimum distance between portable or mobile RF communications equipment (transmitters) and the UNIT as recommended below, according to the maximum output power of the communication equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = \sqrt{\frac{35}{P_1} N F}$	80 MHz to 800 MHz $d = \sqrt{\frac{35}{E_1} N F}$	800 MHz to 2.5 GHz $d = \sqrt{\frac{7}{E_1} N F}$
0.01	0.116	0.116	0.23
0.1	0.369	0.369	0.73
1	1.16	1.16	2.33
10	3.69	3.69	7.38
100	11.6	11.6	23.3

For transmitter rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation application to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer, V1 is the conducted RF compliance level and E1 is the RF radiated compliance level.

Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, object and people.

List of cables that can affect EMC compliance (emission, immunity)

Model S1P

Connection	Type	Length restriction	Shielded?
Scanner - pod	Twisted pair	< 10m	yes
Pod - PC	Twisted pair	< 3m	yes

Model S1A

Connection	Type	Length restriction	Shielded?
Scanner - Cart	Twisted pair	< 3m	yes
Network	Twisted pair	< 5m	yes

WARNING: The use of cables other than those specified above, with the exception of cables sold by 3Shape TRIOS a/s as replacement parts, may result in increased emissions or decreased immunity of the UNIT.

WARNING: The above cables are only to be used with the UNIT. Use of these cables with other medical electrical equipment and/or medical electrical systems may result in increased emissions or decreased immunity of that equipment and/or system.

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I the undersigned, N. E. Christensen, Notary Public of
Copenhagen, Denmark hereby certify that

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has approved and signed this document in my presence.

The identity was proved to me by her passport.

I furthermore certify that pursuant to the records of the Danish
Business Authority (Erhvervsstyrelsen), from May 6, 2021, and
according to an enclosed *Power of Attorney for signing regulatory
paper for device approvals in Russia*, from March 26, 2021 (signed
by Mr. Jakob Just-Bomholt, Mr. Nikolaj Hoffmann Deichmann, and
Mr. Tais Clausen) for the company 3Shape TRIOS A/S, Ms. Ewa
Agnieszka Wybraniec is the sole authorised signatory of the public
limited company

3SHAPE TRIOS A/S.

There were no obvious amendments or additions to the document,
with the exception of those denoted by my signature (initials).


In witness whereof I have hereunto set my hand and notarial seal.

The City Court of Copenhagen, Notarial Acts Division, May 6, 2021

N. E. Christensen

N. E. Christensen

Notary Public of Copenhagen, Denmark

APOSTILLE (Convention de La Haye du 5 octobre 1961)			
1. Country: Land:		Denmark Danmark	
This public document Dette offentlige dokument			
2. has been signed by er underskrevet af		Steen Hlesner	
3. acting in the capacity of i egenskab af		Notary Public Notar	
4. bears the seal/stamp of er forsynet med segl/stempel af		Copenhagen City Court Københavns Hæret	
Certified Attesteret			
5. at i	Copenhagen København	6. the den	06 May 2021 06 maj 2021
7. by af	Ministry of Foreign Affairs of Denmark Udenrigsministeriet		
8. No nr.	9665586A		
9. Seal/stamp: Segl/stempel:		10. Signature: Underskrift:	 Michael Madsen

This Apostille only certifies the authenticity of the signature and the capacity of the person who has signed the public document, and, where appropriate, the identity of the seal or stamp which the public document bears. This Apostille does not certify the content of the document for which it was issued.

To verify the issuance of this Apostille, scan the QR code or visit the following website:

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