«_23_» June 2021 r.

To be submitted to the Federal Service for Surveillance in Healthcare (Roszdravnadzor) in Russian Federation.

We, Olympus Winter & Ibe GmbH, located at Kuehnstr. 61, 22045 Hamburg, Germany, as the legal manufacturer of the Medical device " Hysteroscopic sheath with accessories", hereby state that attached documents is true:

1. Sheath WA47777A instructions for use with addendum
2. Sealing cap A0313 instructions for use with addendum
3. Stopcock accessory WA47778A access information sheet

Best regards,

Dr. Christian Meyer
Managing Director
Olympus Winter & Ibe GmbH

OLYMPUS SURGICAL TECHNOLOGIES EUROPE

Olympus Winter & Ibe GmbH, Kuehnstraße 61, 22045 Hamburg, Germany, P.O. Box 70 17 09, 22017 Hamburg, Germany

Phone: +49 40 669 66-0, Fax: +49 40 669 66-2109, www.olympus-oste.eu

Managing Directors: Dr. André Roggan (Executive Managing Director), Kazutaka Eguchi, Dr. Christian Meyer, Tomohisa Sakurai,

Akihiro Taguchi, Carl Constantin Zangemeister, Reinhard Zentner

Trade Registration Court: Amtsgericht Hamburg HRB 16 328

I, the undersigned Civil Law Notary Dr Matthias Kleiser, Civil Law Notary of this Free and Hanseatic City of Hamburg, Federal Republic of Germany, hereby certify and attest that

Mr. Dr. Christian Martin Meyer
date of birth: 17th day of June 1977
residing: Kuenhstaße 61, 22045 Hamburg, F.R. of Germany
personally known to me, the Notary

personally signed before me.

IN WITNESS WHEREOF, I have hereonto set my hand and affixed my seal of office.

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
SHEATH, STOPCOCK ACCESSORY

Contents

1 General information	5
1.1 User instructions	5
1.2 Signal words.....	5
1.3 Conventions throughout this document.....	5
1.4 Manufacturer	6
2 Safety information.....	7
2.1 Intended use	7
2.2 Contraindications	7
2.3 User qualification	8
2.4 Environment of use	9
2.5 General dangers, warnings and cautions.....	9
3 Product description	11
3.1 Scope of delivery	11
3.2 Sheath.....	11
3.3 Stopcock accessory	13
3.4 Symbols	14
3.5 Warranty.....	14
4 Preparation	15
4.1 Safety information for preparation	15
4.2 Inspection.....	15
4.3 Assembly and connection	17
5 Use.....	20
5.1 Safety information for use	20
5.2 Inserting surgical instruments	21
6 After use	22
6.1 Removal and disassembly	22
6.2 Preparation for reprocessing at the point of use	24
7 Reprocessing	26
7.1 Safety information for reprocessing.....	26

7.2 General information for reprocessing	27
7.3 Manual cleaning	31
7.4 Manual disinfection	33
7.5 Automated cleaning and thermal disinfection	35
7.6 General information for sterilization	38
7.7 Steam sterilization	39
8 Repair, shipment and disposal	41
8.1 Repair	41
8.2 Shipment	41
8.3 Disposal	42
9 Ordering information	43
10 Compatible equipment	44
11 Technical data	45
11.1 General dimensions	45
11.2 Dimensions for reprocessing	46
11.3 Ambient conditions	46

1 General information

1.1 User instructions

- Before use, thoroughly read these instructions for use and the instructions for use of all other products that will be used during the procedure.
- If the required instructions for use are missing, immediately contact an Olympus representative.
- Keep the instructions for use in a safe, accessible location.

1.2 Signal words

The following signal words are used throughout this document.

WARNING

Indicates a potentially hazardous situation which, if not avoided, can result in death or serious injury.

CAUTION

Indicates a potentially hazardous situation which, if not avoided, can result in minor or moderate injury.

NOTICE

Indicates a property damage message.

1.3 Conventions throughout this document

	This is the safety alert symbol. It is used to alert the user to potential physical injury hazards. Observe all safety messages that follow this symbol to avoid possible injury.
	This symbol indicates additional helpful information.

General information

1. 2. 3.	A numeration indicates a sequence of actions.
•	Bullet points indicate individual actions or different options for action.
- - -	Dashes indicate the listing of data, options or objects.
1) 2) 3)	Numbers with right parenthesis name elements in illustrations.

1.4 Manufacturer

Olympus Winter & Ibe GmbH
Kuehnstr. 61
22045 Hamburg
Germany

2 Safety information

2.1 Intended use

2.1.1 General intended use

WA47777A

Hysteroscope sheath for endoscopic diagnosis and treatment in gynecologic applications.

WA47778A

Stopcock accessory to be used in combination with a compatible hysteroscope or ureteroscope. The stopcock has a sealing function, supports the insertion of semi-flexible and flexible instruments and supports the flow regulation of the irrigation fluid during gynecological or urological procedures.

2.1.2 Specific indications

WA47777A

- transcervical diagnosis and treatment of abnormal uterine bleeding, endometrial polypectomy, submucous and intramural myoma, endometrial ablation, lost intrauterine devices, septa and synechiae
- enabeling transcervical sterilization

2.2 Contraindications

Absolute contraindications

- Acute pelvic inflammatory disease

Relative contraindications

- Electrotherapeutical interventions are contraindicated, if, in the judgement of the physician, tissue resection, vaporization, ablation, incision, cutting and coagulation could have a negative effect on the state of the patient.

Safety information

- Electrosurgical treatment of tissue may be contraindicated for patients with cardiac pacemakers, weakened immune systems or blood coagulation disorders.

Hysteroscopy may be contraindicated by the following conditions, depending on the severity or extent:

- Inability to distend the uterus
- Cervical stenosis
- Cervical/vaginal infection
- Uterine bleeding or menses
- Known pregnancy
- Invasive carcinoma of the cervix
- Recent uterine perforation

The intended use and contraindications of the equipment used in combination with the product must be observed.

2.3 User qualification

Medical use

This product is only intended to be used by a trained physician. These instructions for use do not explain or discuss clinical procedures.

Reprocessing

Reprocessing of reusable products may only be performed by qualified hygiene personnel.

Repair

Repair of the product may only be performed by trained qualified servicing personnel that has been authorized by Olympus. Otherwise, Olympus cannot be held responsible for the safety and performance of the product.

2.4 Environment of use

Medical use

This product is only intended to be used in hospitals and medical rooms with appropriate endoscopic equipment.

Reprocessing

Reusable products must be reprocessed according to the applicable national and local standards and regulations.

2.5 General dangers, warnings and cautions

The following warnings and cautions apply to the general handling of the product. This information is to be supplemented by the dangers, warnings and cautions given in each chapter in this document and in the instructions for use of any product being used with this product.

WARNING

Risk of injury to the patient and the user

An insufficient understanding of the warnings, cautions, and instructions can result in serious injury.

- Make sure to have the complete set of instructions for use.
- Before use, thoroughly read these instructions for use and the instructions for use of all other products that will be used during the procedure.

CAUTION

Risk of injury to the patient and the user

Using incompatible equipment and accessories can lead to injury of the patient and the user.

- Only use compatible equipment as listed in the chapter "Compatible equipment" on page 44.

CAUTION

Risk of Injury of the patient

Failure to observe the information given in the section "Intended Use", "User Qualification" and "Environment of Use" amounts to use of the product beyond the use for which it is intended.

- Observe the information given in the section "Intended use", "User qualification" and "Environment of use".

NOTICE

Risk of damage to the product

If instruments touch each other during reprocessing, transport or storage, the instruments can be damaged.

- Avoid that instruments touch each other during reprocessing, transport or storage.

3 Product description

3.1 Scope of delivery

- Before use, check that all items listed below are available.
- Contact an Olympus representative or an authorized service center if any items are missing or damaged.

WA47777A

- Sheath
- Instructions for use

WA47778A

- Stopcock accessory
- Access information sheet

3.2 Sheath

The sheath is a rigid instrument. A telescope can be inserted into the telescope channel of the sheath from the proximal end to provide an endoscopic image during the procedure. The endoscopic image can be viewed using the ocular or by connecting a compatible camera head.


Compatible surgical instruments can be inserted into the instrument channel of the sheath.

The irrigation stopcock can be used to control the inflow of the irrigation fluid. The instrument stopcock controls the passage of the surgical instruments.

For information on compatibility, refer to the chapter "Compatible equipment" on page 44.

As there is no outflow channel for the irrigation fluid, the irrigation fluid flows out between the sheath and the cervix channel.

- 1) Distal end
- 2) Irrigation stopcock
- 3) Irrigation port
- 4) Instrument stopcock
- 5) Instrument channel
- 6) Telescope channel

Irrigation stopcock

The irrigation stopcock controls the inflow of the irrigation fluid. Depending on the required inflow, the irrigation stopcock can be opened or closed.


Instrument stopcock

The instrument stopcock controls the passage of surgical instruments through the instrument channel.


The stopcocks can be disassembled. If the components of the stopcocks are lost or damaged, replace the components with a compatible stopcock accessory.


Stopcock open

Stopcock closed


3.3 Stopcock accessory

The stopcock accessory can replace the stopcocks of the sheath if the stopcock components are lost or damaged.
For more information, refer to the chapter "Ordering Information" on page 43.

3.4 Symbols

This section gives an explanation for each symbol used on the product and on the packaging of the product.

Symbol	Explanation	Symbol	Explanation
	Catalog number		Non-sterile
	Batch code		Consult instructions for use
	Quantity of content	Rx only	Federal (USA) law restricts this device to sale by or on the order of a physician
	Manufacturer		Indicates a recovery/recyclable package or package material
	Date of manufacture		Green Dot-symbol for dual recycling system
	Storage conditions		Indicates the range of humidity to which the medical device can be safely exposed
	Transport conditions		Indicates the temperature limits to which the medical device can be safely exposed
	CE certification mark – symbol for the compliance with the Medical Device Directive 93/42/EEC		

3.5 Warranty

Any warranty claims towards Olympus are forfeited if the user or unauthorized persons attempt repair or modification of the product. No warranty is provided for any damage due to misuse of the product.

4 Preparation

4.1 Safety information for preparation

WARNING

Risk of injury to the patient and the medical personnel
Improper or incomplete reprocessing can cause infection of the patient.

- Reprocess the product before first and each subsequent use.
- For reprocessing the product, observe the instructions in the chapter "Reprocessing" on page 26.
- Do not use a product that has not been reprocessed.
- Before each use, inspect the product as described in the section "Inspection" on page 15.
- Do not use a damaged product.

CAUTION

Risk of injury to the patient and the user

The use of a damaged product or of a product with improper functioning can cause an electric shock, mechanical injury, infection and thermal injury.

- Before each use, observe the instructions in the section "Inspection" on page 15.
- Do not use damaged equipment, equipment that does not function properly or equipment with illegible or missing markings.
- Replace damaged equipment, equipment that does not function properly or equipment with illegible or missing markings.

4.2 Inspection

Inspection regarding reprocessing

- Make sure that the product has been properly reprocessed.
- Visually inspect the product thoroughly. The product must be visually clean.

General inspection


- Check that the product has:
 - no dents, cracks, kinks, or deformations
 - no deep scratches
 - no corrosion
 - no missing or loose parts
- Check all markings on the product for clear visibility.

Even products designed to be reused have a limited service life. A number of factors connected with handling and some reprocessing methods may lead to increased wear of the product. The service life can be markedly shortened. The product must be replaced if signs of wear become visible.

4.3 Assembly and connection


4.3.1 Assembling the stopcocks

- Assemble the irrigation stopcock and the instrument stopcock.
- Securely fasten the components of the stopcocks until they stop.

Preparation

4.3.2 Connecting the irrigation tubing


1. Close the irrigation stopcock.
2. Close the instrument stopcock.
3. Connect the irrigation tubing to the irrigation port.

4.3.3 Attaching the sealing cap

- Attach a sealing cap A0313 to the instrument channel.

4.3.4 Inserting the telescope

1. Insert the telescope into the telescope channel of the sheath.
2. Advance the telescope until it stops.
Make sure that the connector pin of the telescope clicks securely into position.

For information on connecting the telescope, refer to the instructions for use of compatible telescopes.

5 Use

5.1 Safety information for use

WARNING

Risk of injury to the patient

During prolonged procedures, the combination of open blood vessels and certain pressure conditions can cause the irrigation fluid to be absorbed into the circulatory system creating a water intoxication condition called fluid overload-syndrome.

- Perform the procedure with utmost care.
- Carefully monitor fluid inflow and fluid outflow to ensure that the fluid deficit is known at all times.
- If excessive fluid deficit occurs, consider stopping further inflow and ending the procedure.

CAUTION

Risk of injury to the patient

There is a risk of injury if the temperature of the irrigation fluid reaches 43 °C (109 °F).

- Do not preheat the irrigation fluid above body temperature (37 °C/99 °F).

CAUTION

Risk of injury to the patient and the user

There is a risk of injury to the patient and the user due to malfunction of the equipment.

- Always have spare equipment available.

CAUTION

Risk of injury to the patient

When using this product in combination with HF electrodes, arc-over between the endoscopic equipment and the HF electrode can occur.

- Keep a distance between the distal end of the endoscopic equipment and the tip of the HF electrode as specified by the manufacturer of the HF electrode.
For information on the minimum distance, refer to the instructions for use of the HF electrode.

The surface temperature of the product can exceed 41 °C but will not reach 43 °C.

As there is no outflow channel for the irrigation fluid, the irrigation fluid flows out between the sheath and the cervix channel.

5.2 Inserting surgical instruments

1. Close the instrument stopcock.
2. Select a compatible surgical instrument.
3. In the case of a surgical instrument with jaws, close the jaws.
4. Insert the distal end of the surgical instrument into the sealing cap on the instrument channel.
5. Advance the surgical instrument until it stops.
6. Open the instrument stopcock.
7. Advance the surgical instrument until the distal end becomes visible on the endoscopic image.

6 After use


6.1 Removal and disassembly

6.1.1 Removing the surgical instrument


1. In the case of a surgical instrument with jaws, close the jaws.
2. Remove the surgical instrument.
3. Close the instrument stopcock.

6.1.2 Disconnecting the sheath

1. Close the irrigation stopcock.


2. Remove the telescope.

3. Disconnect the irrigation tubing.
4. **Remove the sealing cap.**

After use

6.1.3 Disassembling the stopcocks

- Disassemble the irrigation stopcock and the instrument stopcock.

6.2 Preparation for reprocessing at the point of use

WARNING

Risk of injury to the patient and the medical personnel

If instruments are left in dry condition for a longer period, heavy debris may dry to the instrument, leading to incrustations that may be difficult to remove.

- Do not leave used instruments overnight before reprocessing.
- Reprocess this product immediately after use.

Immediately after use

1. Immediately after use, prepare reusable instruments for reprocessing as described in this section.

2. Remove heavy debris from the product by wiping with an appropriate, single-use lint-free cloth or sponge.
3. Bring the product to the reprocessing area.

6.2.1 Transport to the reprocessing area

WARNING

Risk of Injury to the patient and the medical personnel
Insufficient detergents or disinfectants can lead to protein fixation on the instrument.

- Olympus recommends to return the product from the point of use to the reprocessing area in dry condition to avoid protein fixation caused by some detergents or disinfectants.

Transport in dry condition

1. Make sure that gross debris does not dry.
2. Contain the product during transport to avoid possible environmental or personnel contamination.
3. Reprocess the product according to the instructions in this document.

Transport in wet condition

1. Make sure that the detergents or disinfectants do not cause protein fixation.
2. Contain the product during transport to avoid possible environmental or personnel contamination.
3. Reprocess the product according to the instructions in this document.

NOTICE

Risk of damage to the product

If instruments touch each other during reprocessing, transport or storage, the instruments can be damaged.

- Avoid that instruments touch each other during reprocessing, transport or storage.

7 Reprocessing

7.1 Safety information for reprocessing

WARNING

Risk of Injury to the patient and the medical personnel

Improper or incomplete reprocessing can cause infection of the patient and the medical personnel as well as damage to the product.

- Reprocess the product before first and each subsequent use.
- Reprocess the product according to the instructions in this chapter.
- Only use reprocessing methods which have been validated for the product.

WARNING

Risk of Injury to the medical personnel

Patient debris and reprocessing chemicals are hazardous. There is a risk of infection or skin irritation when they come into contact with exposed skin.

- Wear personal protective equipment to guard against dangerous chemicals and potentially infectious material in accordance with the applicable local guidelines.
- During cleaning and disinfection or sterilization, wear appropriate personal protective equipment, such as eye wear, face mask, moisture-resistant clothing and chemical-resistant gloves that fit properly and are long enough so that no part of the skin is exposed.

WARNING

Risk of Injury to the patient

If the stopcocks are not disassembled before cleaning and disinfection, the stopcocks cannot be reprocessed properly.

- Before cleaning and disinfection, disassemble the stopcocks as described in these instructions for use.

WARNING

Risk of Injury to the patient or the medical personnel

Residual water can affect the gas sterilization process and lead to insufficient reprocessing.

- Thoroughly clean and dry the instrument before gas sterilization.

7.2 General information for reprocessing

Brand-new instruments

- Treat brand-new instruments as if they have been used.
- Reprocess brand-new instruments using a complete reprocessing cycle.

The reprocessing cycle

- Reprocess the product according to the instructions in this chapter.
- To minimize the risk of infecting patients, users or third parties, reprocess the product before each use.
- Contact your local hygiene representative for local standards and regulations.

Manual cleaning or automated cleaning

In general, this product can be cleaned manually or automatically with sufficient results. Manual cleaning methods involve infection risks for the personnel responsible for reprocessing. Automated methods reduce these risks and provide the advantage of standardized, reproducible and validated procedures. Therefore, Olympus recommends to generally prefer automated cleaning methods.

Validated for efficacy

Validated for efficacy means that the efficacy of the method or agent has been validated for reprocessing this product as described in this document.

Reprocessing

Verified for material compatibility

Verified for material compatibility means that according to current knowledge the reprocessing method or agent does not negatively affect materials or functional performance of the product.

Verified for material compatibility does not mean that microbiological efficacy can be guaranteed.

Compatible reprocessing methods (according to ISO 17664)

The microbiological efficacy and/or material compatibility of the reprocessing methods listed in the table have been validated or verified with this product. For detailed information on validated reprocessing methods, refer to the instructions following the table.

- ++** Method validated for microbiological/cleaning efficacy and verified for material compatibility
- +** Method verified for material compatibility
- Incompatible
- o** Contact an Olympus representative for additional information.
- #** Not tested

Process				All products covered in this document
Preparation at the point of use	Deposition	Wet		+
		Dry		+
Decontamination	Preparation	Selection according to the cleaning and disinfection procedures		
	Cleaning	Manual	Ultrasonic	+
			Alkaline detergent	+
			Acidic detergent	#
			Neutral detergent	++
		Machine	Alkaline detergent	++
			Acidic detergent	0
			Neutral detergent	+
	Rinsing			+
	Disinfection	Manual		++
		Machine	Chemical	#
			Thermal	++
	Drying (dry heat)	T _{max}		95 °C (203 °F) (max. 10 min)
Maintenance				
Sterilization	Steam (preferred method)	Autoclave (prevacuum)		++
	Low temperature: steam and formaldehyde			#
	Ethylene oxide			#
	Gas Plasma	STERRAD® 100S		~
		STERRAD® NX™/100NX™		~

Increased wear

Some reprocessing methods may lead to increased wear of the products.

- Inspect the products thoroughly before use for traces of wear.

Cleaning brush

In the following reprocessing instructions, Olympus uses the terms "appropriate brush" and "appropriate surface brush". Select an appropriate brush/surface brush according to the following information:

- The brush diameter must be greater than or equal to the diameter of the inner lumen being brushed.
- The brush bristles must fully extend and contact the surface of the inner lumen.
- The brush must move easily in and out of the instrument.
- The brush must be at least 50 mm longer than the lumen of the instrument.
- Tapered or irregular lumens can require the use of different brushes with different brush diameter.
- A surface brush is for outer surfaces only.
- Do not use a lumen brush for surfaces.
- Do not use a brush with metal bristles or any other types of bristle which can scratch and damage the instrument.
- Only use brushes which the manufacturer has designated for reprocessing of medical devices.

For the dimensions of the product, refer to the section "Dimensions for reprocessing" on page 46.

Instructions for using detergents

- All components of the product must be completely immersed in the detergent.
- Select the immersion time and the concentration of the solution in accordance with the instructions given by the manufacturer of the detergent.
- Make sure that there are no air bubbles.

7.3 Manual cleaning

7.3.1 Safety information for manual cleaning

WARNING

Risk of injury to the patient and the medical personnel

There is a risk of infection when using insufficient detergents.

- Only use detergents which have been validated for the product.

WARNING

Risk of injury to the patient and the medical personnel

The cleaning efficacy can be compromised if the detergent fixes proteins to the product.

- Avoid any protein fixation prior to or during cleaning caused by the reprocessing process.

NOTICE

Risk of damage to the product

Incompatible detergents can considerably damage Olympus endoscopes and accessories.

- Use only solutions that are approved and verified (for example by FDA, DGHM or a comparable institution) and that are certified by their manufacturers as safe for endoscopic instrument cleaning and disinfection.
- Select the immersion time and the concentration of the solution in accordance with the instructions given by the manufacturer of the detergent.

NOTICE

Risk of damage to the product


Detergents can contain various aggressive compounds which can damage the product.

- To remove residues, rinse the product thoroughly with deionized water as described by the manufacturer of the detergent.
- Do not use tap water for rinsing because it might be chlorinated.

7.3.2 Manual cleaning procedure

Manual cleaning of the product has been validated with the following detergents:

- Enzyme-based detergent:
Neodisher MediZym 0.5% manufactured by Chemische Fabrik Dr. Weigert GmbH und Co. KG
- Alkaline detergent:
Neodisher MediClean forte 2% manufactured by Chemische Fabrik Dr. Weigert GmbH und Co. KG

1. Disassemble the irrigation stopcock and the instrument stopcock.
2. Immerse the product in enzyme-based detergent for at least 20 min.
3. Thoroughly rinse the product and all gaps and lumens at least 10 times with enzyme-based detergent using a syringe of at least of 10 ml.
4. Thoroughly brush all gaps and lumens at least 10 times with an appropriate surface brush.

5. Thoroughly flush all gaps and lumens at least 10 times with fresh enzyme-based detergent using a syringe of at least 10 ml.
6. Immerse the product in alkaline detergent and ultrasonicate for at least 15 min.
7. Thoroughly rinse the product with deionized water using a cleaning pistol or other rinse device.
The cleaning pistol or other rinse device must be suitable for cleaning medical devices and deliver a minimum pressure of 1 bar (14.5 psi).
8. Dry the product.

7.4 Manual disinfection

7.4.1 Safety information for manual disinfection

NOTICE

Risk of damage to the product

Incompatible disinfectants can considerably damage Olympus endoscopes and accessories.

- Use only solutions that are approved and verified (for example by FDA, DGHM or a comparable institution) and that are certified by their manufacturers as safe for endoscopic instrument cleaning and disinfection.
- Select the immersion time and the concentration of the solution in accordance with the instructions given by the manufacturer of the disinfectant.

NOTICE

Risk of damage to the product

Disinfectant solutions or tap water can contain various aggressive compounds that can damage the product.

- To remove residues, rinse the product thoroughly with sterile water.
- Do not use tap water for rinsing because it might be chlorinated.


Reprocessing

7.4.2 Precondition for manual disinfection

The product must be cleaned thoroughly before disinfection. Thorough cleaning removes both micro-organisms and organic material. Failure to remove organic material decreases the effectiveness of the disinfection process.

- To avoid damaging the product, it is recommended to use commercially available grasping forceps with rubber-padded jaws to handle the product once it is immersed.
- Grasp the product at its main body.

7.4.3 Manual disinfection procedure

1. Disassemble the irrigation stopcock and the instrument stopcock.
2. Completely immerse the product in Cidex OPA® Solution by Johnson & Johnson for at least 12 minutes at 20 °C (68 °F) or higher, filling all lumens and eliminating all air pockets.

3. Remove the product from the solution.
4. Thoroughly rinse the product by completely immersing it in a large volume (e.g. 2 gallons) of sterile water.
5. Keep the product completely immersed for at least 1 minute.
6. Flush all lumens with large volumes (at least 100 ml) of rinse water.
7. Remove the product and discard the rinse water.
8. To remove disinfectant residues, repeat the rinsing procedure 2 additional times, for a total of 3 rinses, with large volumes of fresh water.
Always use fresh volumes of water for each rinse.
9. Dry the product.

7.5 Automated cleaning and thermal disinfection

Appropriate washer-disinfector with thermal disinfection programs

- Use only a washer-disinfector with thermal disinfection programs which are specifically intended for cleaning and disinfection of reusable medical instruments and accessories by the manufacturer.
- Use only a washer-disinfector with thermal disinfection programs in accordance with international standard series ISO 15883.
- Refer to the instructions for use of the washer-disinfector.
- For automated cleaning, use alkaline detergent.

Detergents in automated cleaning

- Only use detergents which have been validated for the product.
- Select the immersion time and the concentration of the solution in accordance with the instructions given by the manufacturer of the detergent.


Automated cleaning of the product has been validated with the following detergents:

- Neodisher MediClean forte manufactured by Chemische Fabrik Dr. Weigert GmbH und Co. KG
- 0.5% solution

Reprocessing

- Neodisher Z manufactured by Chemische Fabrik Dr. Weigert GmbH und Co. KG
- 0.1% solution

7.5.1 Preparation for automated cleaning and thermal disinfection

1. Disassemble the irrigation stopcock and the instrument stopcock.
2. Place the components of the stopcocks in an instrument basket for small items.
3. Make sure that instruments are placed in the washer-disinfector so that they can be appropriately cleaned and disinfected.
4. Do not overload the washer-disinfector.

The components of the stopcocks can get lost or be damaged if reprocessed without using an instrument basket for small items.

7.5.2 Selection of programs

CAUTION

Risk of Injury to the patient

High inflow temperatures during automated cleaning and thermal disinfection can cause thermal coagulation of proteins.

- Make sure that the first water inflow into the washer-disinfector during automated cleaning and thermal disinfection starts with cold temperatures (<20 °C) to avoid thermal coagulation of proteins.

Select the program according to the following characteristics:

- The program must at least reach an A0 value of 3000.
- The temperature of the first inlet water must not exceed 20 °C.
- A temperature of 55 °C must not be exceeded during cleaning.
- Do not use programs starting with high inflow temperatures (e.g., 93 °C), because this causes protein fixation.
- Make sure that the program does not include any sudden changes in temperature.

7.5.3 After automated cleaning and thermal disinfection

1. To prevent corrosion, remove all instrument parts from the washer-disinfector after the automated procedure is completed.
2. Be careful as the instrument parts can be hot.
3. If the instruments are wet, dry the instruments.

7.6 General information for sterilization

Precondition for sterilization

The product must be cleaned thoroughly before sterilization. Thorough cleaning removes both micro-organisms and organic material. Failure to remove organic material decreases effectiveness of the sterilization process. After cleaning, make sure that the instrument is carefully dried.

Storage life of sterilized instruments


The storage life of sterilized instruments depends on the type of packaging and the storage conditions. Refer to national and local laws and guidelines.

Sterilization accessories

Use only sterilization accessories that have been cleared by the local authorities and that have been validated for the sterilization processes of the product.

7.7 Steam sterilization

Preparing the sheath for steam sterilization

1. Assemble the irrigation stopcock and the instrument stopcock.
2. Securely fasten the components of the stopcocks until they stop.
3. Seal the product in packaging appropriate for steam sterilization according to local guidelines and your hospital's protocol.

Steam sterilization conditions

Steam sterilization of the product has been validated with fractionated prevacuum at 132 °C for 3 min.

Other possible steam sterilization conditions are:

Temperature 132 to 136 °C (269.6 to 276.8 °F)

Holding time 3 to 20 min.

or

Reprocessing

Temperature 135 °C (275 °F)

Holding time 20 min.

Steam sterilization procedure

For steam sterilization, Olympus recommends a drying time of at least 15–30 minutes.

- Choose an appropriate drying time, depending on the autoclave.
- Refer to the autoclave's instructions for use.
- Be careful when unloading the autoclave, the contents may be hot.
- Make sure that the packaging of the product is not damaged.
- Resterilize the product if 1 of the following applies:
 - the sealing has been opened
 - the packaging is wet
 - the packaging is damaged

8 Repair, shipment and disposal

8.1 Repair

WARNING

Risk of injury to the patient and the user

There is a risk of injury to the patient and the user caused by unauthorized repairs and product modification.

- Do not attempt to repair or modify the product.

Repairs may only be performed by qualified servicing personnel that has been authorized by Olympus. Contact an Olympus representative or an authorized service center for repair and warranty information.

8.2 Shipment

WARNING

Risk of infection for the personnel

Sending used products to Olympus represents a risk of infection for the personnel. Olympus is entitled to refuse soiled or contaminated products for reasons of safety.

- Reprocess the product before sending in.
- If reprocessing will damage the product completely, clean the product as thoroughly as possible and mark it accordingly.

When sending products to an authorized service center for repair or warranty claims, note the following.

- Use the original cardboard packaging for the transport of the product.
- If this is not possible, wrap each component individually in sufficient paper or sheets of foamed material and place them in a cardboard box.
- If applicable, use the respective transport protection (e.g., instrument tray, protective tube) of the product.

Service centers do not accept warranty claims for damage caused by inadequate packaging.

8.3 Disposal

WARNING

Risk of injury to medical personnel and third party

Disposing of the product without prior decontamination can cause infection of the medical personnel.


- Decontaminate the product before disposal.
- When disposing of the product or any of its components, follow all applicable national and local laws and guidelines.

9 Ordering information

Product description	Catalog number
Sealing cap	A0313
Stopcock accessory	WA47778A
Sheath WA47777A with sealing cap A0313	WA47777S


10 Compatible equipment

Olympus recommends the use of equipment as listed in this chapter. If combinations are used that are not listed in this chapter, the user takes the full responsibility. Future equipment can also be compatible. For more information, contact an Olympus representative. Some of the products listed in this chapter may not be available in all sales territories.

11 Technical data

11.1 General dimensions

Maximum working length

L1: WA47777A 206.85 mm

Maximum insertion portion width


D1: WA47777A 5.25 mm

Minimum instrument channel width

WA47777A 5 Fr.

There is no guarantee that instruments selected solely using this minimum instrument channel width will be compatible in combination.


Maximum width

W1: WA47777A 3.54 mm

11.2 Dimensions for reprocessing

This section lists all relevant dimensions for reprocessing the product. For reprocessing the product, refer to the chapter "Reprocessing" on page 26.

Minimum channel diameter

D1.....	4.83 mm
D2.....	3.16 mm
D3.....	3.16 mm

Minimum channel length

L1	280.2 mm
L2	48.0 mm
L3	73.7 mm

11.3 Ambient conditions

Operating conditions

Ambient temperature.....	10 to 35 °C (50 to 95 °F)
Relative humidity.....	30 to 85%

Storage conditions

Temperature.....	10 to 40 °C (50 to 104 °F)
Relative humidity.....	30 to 85%

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
**Addendum to operating instructions
for medical device: “Hysteroscopic sheath with accessories”**

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



Specifications for the Sealing cap (catalog number A0313):


Characteristics	Value
<i>Tolerance $\pm 5\%$ if not specified</i>	
Length (L)	24 mm
Maximum diameter (D1)	10 mm
Diameter (D2)	5 mm
Diameter (D3)	3.75 mm
Weight	1 g

Specifications for the Stopcock accessory (catalog number WA47778A):

Characteristics	Value
<i>Tolerance $\pm 5\%$ if not specified</i>	
General length (L)	19.1 $\pm$ 0.5 mm
Width (D)	9.5 $\pm$ 0.5 mm

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for single use

A0282
A0313

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

Contents

1 General information	5
1.1 User instructions	5
1.2 Signal words	5
1.3 Conventions throughout this document	5
1.4 Manufacturer	6
2 Safety information	7
2.1 Intended use	7
2.2 User qualification	7
2.3 Environment of use	7
2.4 General warnings and cautions	8
3 Product description	10
3.1 Scope of delivery	10
3.2 Sealing cap	10
3.3 Symbols	11
4 Preparation	12
4.1 Safety information for preparation	12
4.2 Inspection	12
4.3 Steam sterilization	12
4.4 Assembly and connection	13
5 After use	14
5.1 Disposal	14
5.2 Shipment	14
6 Technical data	16
6.1 General dimensions	16
6.2 Ambient conditions	16

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	This is the safety alert symbol. It is used to alert the user to potential physical injury hazards. Observe all safety messages that follow this symbol to avoid possible injury.
	This symbol indicates additional helpful information.
1. 2. 3.	A numeration indicates a sequence of actions.

General information

•	Bullet points indicate individual actions or different options for action.
- - -	Dashes indicate the listing of data, options or objects.
1) 2) 3)	Numbers with right parenthesis name elements in illustrations

1.4 Manufacturer

Olympus Winter & Ibe GmbH
Kuehnstr. 61
22045 Hamburg
Germany

2 Safety information

2.1 Intended use

Sealing cap intended to be used with compatible instruments designed for endoscopic diagnosis and treatment in urologic and gynecologic applications.

- Do not use it for any other purpose.

2.2 User qualification

Medical use

This product is only intended to be used by a trained physician or trained qualified medical personnel under the supervision of a physician.

These instructions for use do not explain or discuss clinical procedures.

Reprocessing

Reprocessing may only be performed by qualified hygiene personnel.

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Reprocessing

The products must be reprocessed according to the applicable national and local standards and regulations.

2.4 General warnings and cautions

The following warnings and cautions apply to the general handling of the product. This information is to be supplemented by the dangers, warnings and cautions given in each chapter in this document or in the instructions for use of any product being used with this product.

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An insufficient understanding of the dangers, warnings, cautions, and instructions can result in serious injury or damage to the product.

- Make sure to have the complete set of instructions for use.

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The use of damaged equipment, equipment that does not function properly or equipment with illegible or missing markings may cause infection, toxic shock, electrical, mechanical and thermal injury or unintended nerve stimulation.

- Before each use, observe the instructions in the section "Inspection" on page 12.
- Do not use damaged equipment, equipment that does not function properly or equipment with illegible or missing markings.
- Replace damaged equipment, equipment that does not function properly or equipment with illegible or missing markings.

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There is a risk of injury to the patient and/or the user caused by unauthorized repairs and product modification.

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CAUTION

Risk of Injury of the patient

Failure to observe the information given in the section "Intended use", "User qualification" and "Environment of use" amounts to use of the product beyond the use for which it is intended.

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A0282, A0313:

- 10 sealing caps
- Instructions for use

If the sealing caps are delivered in combination with instruments, the scope of delivery can include less than 10 sealing caps.


3.2 Sealing cap

The sealing cap can be attached to the instrument ports of compatible equipment. Compatible surgical instruments can be inserted through the sealing cap. The sealing cap prevents the irrigation fluid from leaking out of the instrument.


The sealing caps are for single use only and may only be sterilized once.

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Symbol	Explanation	Symbol	Explanation
	Catalog number		Non-sterile
	Batch code		Do not reuse
	Quantity of content		Consult instructions for use
	Manufacturer		Green Dot-symbol for dual recycling system
	Storage conditions		Transport conditions
	Indicates the range of humidity to which the medical device can be safely exposed		Indicates the temperature limits to which the medical device can be safely exposed
	Indicates a recovery/recyclable package or package material		CE certification mark – symbol for the compliance with the Medical Device Directive 93/42/EEC

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WARNING

Risk of injury to the patient

The product is for single use only and delivered in non-sterile condition. Reuse or resterilization can lead to malfunction of the product and injury of the patient.

- Sterilize the product before first use.
Observe the instructions in the section "Steam sterilization" on page 12.
- Dispose of the product after use.
- Do not reuse or attempt to reprocess after use.

4.2 Inspection

General Inspection

- Check that the product has:
 - no dents, cracks
 - no deep scratches
- Check all markings on the product for clear visibility.

4.3 Steam sterilization

WARNING

Risk of injury to the patient and/or the medical personnel
Improper or incomplete reprocessing can cause infection of the patient and the medical personnel as well as damage to the product.

- Reprocess the sealing cap before first use using steam sterilization.

Preparing the sealing cap for steam sterilization

- If the sealing caps are attached to an instrument at delivery, remove the sealing caps from the instruments.

Steam sterilization conditions

Steam sterilization of the sealing cap has been validated with fractionated prevacuum at 132 °C for 3 min.

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- Choose an appropriate drying time, depending on the autoclave.
- Refer to the autoclave's instructions for use.
- Be careful when unloading the autoclave, the contents may be hot.
- Make sure that the sterile packaging of the product is not damaged.

4.4 Assembly and connection

- Attach the distal end of the sealing cap to the instrument channel of a compatible instrument.

5 After use

5.1 Disposal

WARNING

Risk of injury to the medical personnel

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-
- Dispose of the product after use.
 - Do not reuse or attempt to reprocess.

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- Reprocess the product as thoroughly as possible before sending in.
- Clean the product as thoroughly as possible and mark it accordingly.

When sending products to an authorized service center for warranty claims, note the following.

- Use the original blister and cardboard packaging to transport the product.
- If this is not possible, wrap each component individually in sufficient paper or sheets of foamed material and place them in a cardboard box.

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6.1 General dimensions

A0282

Color.....transparent
Passage for instruments 6-8 Fr. instruments

A0313

Color.....black
Passage for instruments 4-5 Fr. instruments

6.2 Ambient conditions

Operating conditions

Ambient temperature.....10 to 35 °C (50 to 95 °F)
Relative humidity30 to 85%

Storage conditions

Temperature10 to 40 °C (50 to 104 °F)
Relative humidity30 to 85%

Transport conditions

Temperature..... -40 to 70 °C (-40 to 158 °F)
Relative humidity 10 to 95%

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
**Addendum to operating instruction for
device: “Sealing cap A0313”**

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
Reclamation

All complaints related to the circulation of a medical device in the Russian Federation should be sent:
 Olympus Moscow LLC
 Elektrovzavodskaya str. 27 bld. 8, 107023 moscow, Russia
 Tel.: +7 (495) 926-7077
 info@olympus.com

en

Access information sheet

When using this product, follow the instructions for use for WA47777A or WA2UR11A, WA2UR12A, WA2UR13A, WA2UR14A, WA2UR21A, WA2UR22A, WA2UR23A, WA2UR31A, WA2UR32A, WA2URADP.

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© Copyright 2017 Olympus Winter & Ibe GmbH
© Copyright 2019 Olympus Winter & Ibe GmbH