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NOTARIAL CERTIFICATE

LAW FIRM IDEA & NOTARY OFFICE INC.


(Cheongjin-dong, Samsong Bldg.) 701, 702,
34, Jong-ro 3-gil, Jongno-gu, Seoul, Republic of Korea

Instruction for Use
DS MED Intradermal cannula in the set with needles for puncture, in kits

«I certify»

"26" January 2022

Mr. Kim Keun Sik

CEO of FEEL TECH CO., LTD.



3,4 Floor, Standard Factory 2-dong, 15,

Jayumuyeok 2-gil, Gunsan-si, Jeollabuk-do, Korea

Instruction for Use


for medical device

DS MED Intradermal cannula in the set with needles for puncture, in kits

 KIM KEUN SIK
FEEL TECH CO., LTD.

PRESIDENT KIM KEUN SIK

- 2022 -

Instruction for Use
DS MED Intradermal cannula in the set with needles for puncture, in kits

Instruction for Use

DS MED Intradermal cannula in the set with needles for puncture, in kits

Content

1. Identification of medical device.....	4
1.1. Name of the medical device.....	4
1.2. Product classification and risk class.....	4
1.3. Information about Developer.....	5
1.4. Information about Manufacturer.....	5
1.5. Embodiments of medical device.....	6
2. Appointment of medical device.....	6
2.1. Principles of Operation.....	7
2.2. Application area.....	7
2.3. Potential users.....	7
3. Use of medical device and precautions.....	7
3.1. Mode of application.....	7
3.2. Application conditions.....	10
3.3. Indications.....	10
3.4. Contraindications.....	10
3.4.1. General contraindications.....	10
3.4.2. Local contraindications.....	11
3.5. Side effects.....	11
3.6. Safety requirements.....	11
4. Description of medical device.....	12
4.1. Overview.....	12
4.2. Parts and components.....	13
4.3. Materials in contact with the human body.....	14
4.4. Specifications.....	18
5. Packaging and labeling.....	23
5.1. Packaging of medical device.....	23
5.2. Labelling of medical device.....	24
6. Aftermarket Consumer Information.....	30
6.1. Shelf life.....	30
6.2. Terms of transportation.....	30
6.3. Storage conditions.....	31
6.4. Maintenance and repairs.....	31
6.5. Procedure and conditions for disposal.....	31
6.6. Environmental requirements.....	31
6.7. Warranty.....	31
7. Authorized representative of the manufacturer in the Russian Federation.....	32

1. Identification of medical device

1.1. Name of the medical device

DS MED Intradermal Cannula in the set with needles for puncture, in kits (also hereinafter referred to for convenience "cannula (s) DSMED" and / or " DSMED cannula (s) " and / or "cannula (s)" and / or "the device(s)").

1.2. Product classification and risk class

Risk Class IIa in accordance with Council Directive 93/42/EEC+2007/47/EC

Annex IX - Classification Criteria

Chapter III - Classification

Section 2 - Invasive Devices

Sub-section 2.2. - Rule 6.

All surgically invasive devices intended for transient use are in Class IIa unless they are:

- intended specifically to control, diagnose, monitor or correct a defect of the heart or of the central circulatory system through direct contact with these parts of the body, in which case they are in Class III,
- reusable surgical instruments, in which case they are in Class I,
- intended specifically for use in direct contact with the central nervous system, in which case they are in Class III,
- intended to supply energy in the form of ionising radiation in which case they are in Class IIb,
- intended to have a biological effect or to be wholly or mainly absorbed in which case they are in Class IIb,

Instruction for Use

DS MED Intradermal cannula in the set with needles for puncture, in kits

— intended to administer medicines by means of a delivery system, if this is done in a manner that is potentially hazardous taking account of the mode of application, in which case they are in Class IIb.

1.3. Information about Developer

Developer:


FEEL TECH CO., LTD.

3,4 Floor, Standard Factory 2-dong, 15, Jayumuyeok 2-gil, Gunsan-si, Jeollabuk-do, Korea / 3,4 floor Standard Factory 2-dong, 15, Dzhayumuyok 2-gil, Gunsan- si Dzeollabuk-do, Korea

1.4. Information about Manufacturer

Feel Tech Co., Ltd. was established on June 15 in 1991. On January 29 in 1999, we were acquired certification of ISO9001:2008, ISO13485:2003

1.4.1. Homepage : www.feeltech.co.kr

1.4.2. Manufacturer : Feel Tech Co., Ltd

1.4.13. Manufacturer address: 3,4 Floor, Standard Factory 2-dong, 15, Jayumuyeok2-gil, Gunsan-si, Jeollabuk-do, Korea

1.4.4. Production place: 3,4 Floor, Standard Factory 2-dong, 15, Jayumuyeok2-gil, Gunsan-si, Jeollabuk-do, Korea

1.4.5. Storage place: 3,4 Floor, Standard Factory 2-dong, 15, Jayumuyeok2-gil, Gunsan-si, Jeollabuk-do, Korea

1.4.6. Manufacturer representative : Geun Sik, Kim

1.4.7. Phone TEL : +82-63-468-6626~8 FAX : +82-63-468-6623

1.4.8. Contact person : Cheol Ho, Choi

1.4.9. E-mail : info@feeltech.co.kr

Instruction for Use
DS MED Intradermal cannula in the set with needles for puncture, in kits

1.5. Embodiments of medical device

DS MED Intradermal Cannula in in the set with needles for puncture, in kits:

№№	Kit	Content
1	FTMN 22-70-21	1.1. Cannula 22G, length 70 mm (0,70 x 70 mm) 1.2. Needle 21G, length of 1" (0,80 x 25,4 mm)
2	FTMN 22-50-21	2.1. Cannula 22G, 50 mm length (0,70 x 50 mm) 2.2. Needle 21G, length of 1" (0,80 x 25,4mm)
3	FTMN 25-70-23	3.1. Cannula 25 G, 70 mm length (0,50 x 70 mm) 3.2. Needle 23G, length 1" (0,60 x 25,4 mm)
4	FTMN 25-50-23	4.1. Cannula 25 G, 50 mm length (0,50 x 50 mm) 4.2. Needle 23G, length of 1" (0,60 x 25,4 mm)
5	FTMN 25-40-23	5.1. Cannula 25 G, 40 mm length (0,50 x 40 mm) 5.2. Needle 23G, length 1" (0,60 x 25,4 mm)
6	FTMN 27-50-26	6.1. Cannula 27G, 50 mm length (0,40 x 50 mm) 6.2. Needle 26G, length 1/2" (0,45 x 12,7 mm)
7	FTMN 27-40-26	7.1. Cannula 27G, 40 mm length (0,40 x 40 mm) 7.2. Needle 26G, length 1/2" (0,45 x 12,7 mm)
8	FTMN 27-25-26	8.1. Cannula 27G, length 25,4 mm (0,40 x 25,4 mm) 8.2. Needle 26G, length 1/2" (0,45x12,7 mm)
9	FTMN 30-25-27	9.1. Cannula 30G, length 25,4 mm (0,30 x 25,4 mm) 9.2. Needle 27G, length 1/2" (0,40 x 12,7 mm)

2. Appointment of medical device

The product is intended for administration to a dermal filler or subcutaneous layers to impart additional volume of tissue in cosmetic or reconstructive purposes.

2.1. Principles of Operation

Injecting the medicine with and/or without syringes. This device can be inserted into the soft tissue with a little injury of tissue. The device is inserted into vascular or coelom in general-purpose. And this is composed of cannula, hub, protect cap and is connected to connector or tube. The end of cannula is closed and has a hole on the side of cannula which has the function to inject or suck the body fluid or liquid medicine

2.2. Application area


Aesthetic Medicine, cosmetology.

2.3. Potential users

The product is intended for professional use by healthcare professionals.

3. Use of medical device and precautions

3.1. Mode of application

Selection of administration techniques (linear, fan-shaped, and the like), filler, as well as instruments for its administration in accordance with the specific indication.

Selecting the length and outer diameter (G) of the cannula is dependent on the anatomical parameters of the patient and the filler used.

Before using the product, check the packaging for damage and make sure the expiration date has not expired.

The injection site to be treated by the usual method antiseptic (recommended to use 70% alcohol solution) and, if necessary, a means for analgesia.

Analyze the face, skin, muscle, if necessary, to make the marking areas with a pencil (to determine the location and amount of filler administration). Determine the injection point depending on the required correction zones.



Before the introduction of the filler via a cannula pre-made pointed puncture needle supplied in the package with the cannula.

Instruction for Use

DS MED Intradermal cannula in the set with needles for puncture, in kits

The rounded tip of the cannula, previously fixed on the syringe cone with fillers, accurately placed in a puncture made by the needle with simultaneous stretching of the skin. Then smooth translational motion correction is carried out to place the cannula and the filler is introduced in a retrograde fashion.

After use the needle and cannula must be discarded into the container for medical waste in compliance with safety measures: to not touch the needle and is not suitable to put on the cap a used needle.


Selection Table DSMED cannula in accordance with the field of application:

dimensions	Field of application (recommended for delicate zones correction) *	surfactants subcutaneously
30G	- lips and paraoral zone; - purse string wrinkles; - paroorbital zone; - glabella region; - bridge; - forehead	shallow subcutaneous

Instruction for Use

DS MED Intradermal cannula in the set with needles for puncture, in kits

27G	- wrinkles puppets; - lacrimal groove; - forehead, - back of hands - glabella region, - wrinkles, "crow's feet" - the area around the eyes, - purse string wrinkles - nasolabial folds, - lips	The surface layer VLS (subcutaneous fat) in a supra- periosteal way.
25G	- frown wrinkles; - nasolabial folds; - the neck and chest; - corners of the mouth. - lacrimal groove, - nose, - the middle third of the face, cheeks, cheekbones, nososkulovaya furrow - facial contour, - lips - the chin	 The surface layer VLS (subcutaneous fat) in a supra- periosteal way.
22G	- cheeks / cheek / chin; - nasolabial and glabellar wrinkles; - nose	The VLS (subcutaneous fat), nanadkostnichno.

*Note:

The diameter and length of the cannula and needle are selected in accordance with the recommendations set forth in the documentation on the gels and the solutions to be injected or physician preferences.

3.2. Application conditions


Temperature $+ 5^{\circ} \text{C} - + 25^{\circ} \text{C}$;

Relative humidity (less than 70) $\pm 2\%$;

Atmospheric pressure is $755 \pm 20 \text{ mm Hg. Art.}$

3.3. Indications

- the introduction of fillers into the skin and subcutaneous layers.

3.4. Contraindications

Contraindications to the use of products include fillers contraindications for administration as a whole, namely:


3.4.1. General contraindications

- allergic reactions or intolerance of individual components of the product or products;
- autoimmune disease or systemic connective tissue disorders;
- mental disorders;
- infectious diseases;
- severe somatic diseases;
- any exacerbation of chronic diseases;
- Pregnancy and lactation;
- anaphylactic shock cases or "rejection" foreign body in the history;
- Hemophilia, anticoagulation and antiplatelet agents;

- bleeding disorder;
- recurrent sore throats and fevers.

3.4.2. Local contraindications

- inflammation, violation of the integrity of the skin and the presence of nerezorbiruemogo implant in the area where you plan to do a puncture;
- presence in areas of prospective injections of other spies, biodegradable materials based on polylactic acid, collagen, hydroxyapatite or calcium;
- malignancies in the proposed injection zone;
- recently conducted procedure chemical peeling, laser or mechanical grinding;
- tendency to keloid and hypertrophic scars;
- viral, fungal diseases of the skin.

3.5. Side effects

- Pain.
- Hematoma.
- Local edema.
- Erythema.
- Seal.
- Telangiectasia.
- Infectious complications.

3.6. Safety requirements

- Prior to use, make sure that the product is sterile (individual sterile packaging is not damaged).
- Before operating the read instructions for use.

- Do not use products that have expired.
- The product is for single use, do not use the product again, not to re-sterilization.
- Any use of the product other than its intended use is prohibited.
- It can only be used by trained and qualified medical personnel.

4. Description of medical device

4.1. Overview

DS MED Intradermal Cannula in in the set with needles for puncture, in kits produced by "Feel Tech Co., Ltd.," Korea are thin flexible hollow tubes made of stainless steel with a rounded end at the distal end which is the side opening and at the proximal end - a connector for connection to device with fillers.

A feature of the cannula DSMED is the minimum distance from the tip opening, which reduces the flow of the filler by reducing its volume in the distal portion of the cannula (with the use of conventional cannulae drug consumption increases by 10-15%).

Furthermore, DSMED cannula possess sufficient flexibility and elasticity for their passage through the tissue and return to their original shape as retrograde excretion outwardly.


Labels applied via anodic oxidation at regular intervals on the tube of the cannula, can control depth of its administration, and on the connector cannula point indicates the location of the outlet cannula. Silicone coating cannula tube greatly facilitates the introduction of the process and increases the flexibility of the cannula.

INSTRUCTION FOR USE
DS MED Intradermal cannula in the set with needles for puncture, in kits


Due to the fact that the rounded end of the cannula prevents disrupt the integrity of the skin, for skin puncture used disposable sterile needle (the needle for puncture the primary), supplied with the cannula.

Cannulae and needles for puncture with safety caps to protect against external influences.

4.2. Parts and components


Cannula:

- 1 - connector
- 2 - the cannula tube
- 3 - protective cap

Needle for primary puncture:

INSTRUCTION FOR USE
DS MED Intradermal cannula in the set with needles for puncture, in kits

1 - needle tube 2 - connector 3 - protective cap


4.3. Materials in contact with the human body

	Name	Material
1.	Cannula	Stainless steel: (STS 304L) Manufacturer: GOKA CO, LTD TESHIMA KOREA CO, LT).. Polypropylene: HI-PRENE H551 (CAS №: 9003-07-0), Manufacturer: LG-Caltex Oil Corporation, KOREA Polyethylene GRADE: XJ700 (CAS №: 9002-88-4), manufacturer LG-Caltex Oil Corporation, KOREA Epoxy resin: (CAS №: 25068-38-60), Manufacturer: SEGO ENGINEERING COMPANY, Korea Silicone: polydimethylsiloxane (PDMS) (SA-1A),

	Manufacturer: SEGO ENGINEERING COMPANY, Korea Dyes: 1. Black {Raw materials(quantity): PANAX YELLOW 3069 0.75g, CAS No:71832-85-4 / PANAX BLUE B 0.2g, CAS No: 147-14-8 / Carbon Black 1g, CAS No: 1333-86-4} Mixing rate (Color+PP) in Total quantity 100%: 0.075% + 0.002% + 0.1% + 99.823% Manufacturer: LG-Caltex Oil Corporation, Korea 2. Orange {Raw materials(quantity): PIGMENT YELLOW 180 2g, CAS No: 77804-81-0 / PANAX RED 2R 2g, CAS No: 7023-61-2} Mixing rate (Color+PP) in Total quantity 100%: 0.2%+0.2%+99.6% Manufacturer: LG-Caltex Oil Corporation, Korea 3. Medium-Grey {Raw materials(quantity): PANAX BLUE B 3.8g, CAS No: 147-14-8 / PANAX YELLOW 3069 0.5g, CAS No: 71832-85-4} Mixing rate (Color+PP) in Total quantity 100%: 0.38%+0.05%+99.57% Manufacturer: LG-Caltex Oil Corporation, Korea 4. Yellow {Raw materials(quantity): PIGMENT YELLOW 180 1.25g, CAS No: 77804-81-0} Mixing rate (Color+PP) in Total quantity 100%: 0.125%+99.875% Manufacturer: LG-Caltex Oil Corporation, Korea
--	---

2.	Needle for primary puncture	Stainless steel: (STS 304L) Manufacturer: GOKA CO, LTD TESHIMA KOREA CO, LT).. Polypropylene: HI-PRENE H551 (CAS №: 9003-07-0), Manufacturer: LG-Caltex Oil Corporation, KOREA Polyethylene GRADE: XJ700 (CAS №: 9002-88-4), manufacturer LG-Caltex Oil Corporation, KOREA Epoxy resin: (CAS №: 25068-38-60), Manufacturer: SEGO ENGINEERING COMPANY, Korea Silicone: polydimethylsiloxane (PDMS) (SA-1A), Manufacturer: SEGO ENGINEERING COMPANY, Korea Dyes: 1. Deep-Green {Raw materials(quantity): PANAX GREEN 4059 3.5g, CAS No: 1328-53-6 / PIGMENT YELLOW 180 0.3g, CAS No:77804-81-0} Mixing rate (Color+PP) in Total quantity 100%: 0.35%+0.08%+99.57% Manufacturer: LG-Caltex Oil Corporation, Korea 2. Deep-Blue {Raw materials(quantity): PANAX BLUE B 4g, CAS No: 147-14-8} Mixing rate (Color+PP) in Total quantity 100%: 0.4%+99.6% Manufacturer: LG-Caltex Oil Corporation, Korea
----	-----------------------------	---

		3. Brown {Raw materials(quantity): PIGMENT YELLOW 180 2g, CAS No: 77804-81-0 / PANAX RED 112 1.5g, CAS No: 6535-46-2} Mixing rate (Color+PP) in Total quantity 100%: 0.2%+0.15%+99.65% Manufacturer: LG-Caltex Oil Corporation, Korea 4. Medium-Grey {Raw materials(quantity): PANAX BLUE B 3.8g, CAS No: 147-14-8 / PANAX YELLOW 3069 0.5g, CAS No: 71832-85-4} Mixing rate (Color+PP) in Total quantity 100%: 0.38%+0.05%+99.57% Manufacturer: LG-Caltex Oil Corporation, Korea
--	--	---

4.4. Specifications

DS MED Intradermal Cannula:

	Option			
Gauge	22G	25G	27G	30G
Wall thickness (RW – regular wall, TW – thin wall, ETW – extra thin wall, UTW – ultra thin wall)	TW	TW	TW	TW
Hub color	Black	Orange	Medium-Grey	Yellow
Width of the label mm	0.80 ± 0.05	0.80 ± 0.05	0.80 ± 0.05	0.80 ± 0.05
Distance between the marks, mm	8.5 ± 0.5	8.5 ± 0.5	8.5 ± 0.5	8.5 ± 0.5
Distance from the distal end to the far double mark, mm	$50 \pm 0,02$	$50 \pm 0,02$	$50 \pm 0,02$	$50 \pm 0,02$
Outer tube diameter (OD), mm	0,698-0,730	0,500-0,530	0,400-0,420	0,298-0,320
Inner diameter, ($\pm 0,01$) mm	Min. 0,440	Min. 0,292	Min. 0,241	Min. 0,165
Patency of lumen: the diameter of stainless steel stylet to pass through the needle, mm	0.35	0.23	0.19	0.13
Nominal tube length, mm	70 (+ 1.5 / -2.5) 50 (+ 1.5 / -2.5)	70 (+ 1.5 / -2.5) 50 (+ 1.5 / -2.5) 40 (0 / -4)	50 (+ 1.5 / -2.5) 40 (0 / -4) 25,4(+ 1.5 / -2.5)	25,4 (+ 1.5 / -2.5)


DS MED Intradermal cannula in the set with needles for puncture, in kits

Dimensions of exit hole (length x width), mm	$1,70 \pm 0,2$ x $0,9 * OD \sim 1 * OD$	$1,50 \pm 0,2$ x $0,9 * OD \sim 1 * OD$	$1,30 \pm 0,2$ x $0,9 * OD \sim 1 * OD$	$0,90 \pm 0,2$ x $0,9 * OD \sim 1 * OD$
Distance from the opening to the distal end, mm $\pm 0,05$	0.80	0.70	0.50	0.40
Force to test bond between hub and needle tube, N	40	22	22	22
Radius of cannula blunting, mm $\pm 10\%$	0.35	0.25	0.20	0.15
Cleanliness of outer surface of cannula tube	When inspected by normal or corrected-to-normal vision without magnification under an illuminance of 300 lx to 700 lx, the surface of the cannula tube shall appear free from particles and extraneous matter.			
Cleanliness of hub socket	When examined under $2,5\times$ magnification, the hub socket (fluid path surface) shall appear free from particles and extraneous matter			
Limits for acidity or alkalinity	When determined with a laboratory pH meter and using a general purpose electrode, the pH value of a prepared extract shall be within one unit of pH of that of the control fluid.			
Limits for extractable metals	Cannula extract shall, when corrected for the metals content of the control fluid, contain not greater than a combined total of 5 mg/l of lead, tin, zinc and iron. The cadmium content of the extract shall, when corrected for the cadmium content of the control fluid, be lower than 0,1 mg/l			

DS MED Intradermal cannula in the set with needles for puncture, in kits

The alignment of the connection	Axis of the tube and the hub must be the same. Maximum allowable deviation from the axis of the hub is $\pm 3^\circ$
The roughness of the outer surface of the tube	Ra less than 0.20 μm
Conical fitting type of needle hub	"Luer" / "luer-lock" taper 6%
Corrosion resistance	Must be corrosion resistant
Lubricant	Amount of lubricant should be minimal and not form liquid droplets on the outer surface and the inner needle when viewed normal or corrected to normal vision
Balance of ethylene oxide	Ethylene oxide (EO) and ethylene chlorohydrin (ECH) <25 p. / million. Ethylene glycol <250 p. / million


Needles for puncture:

	option			
Gauge	21G	23G	26G	27G
Wall thickness (RW – regular wall, TW – thin wall, ETW – extra thin wall, UTW – ultra thin wall)	TW	TW	TW	TW
Hub Color *	Deep-Green	Deep-Blue	Brown	Medium- Grey
The nominal outer diameter of the needle tube, ($\pm 0,01$) mm	0,800-0,830	0,600-0,673	0.440-0.470	0,400-0,420
Inner diameter, ($\pm 0,01$) mm	Min. 0,547	Min. 0,370	Min. 0.292	Min. 0.241

DS MED Intradermal cannula in the set with needles for puncture, in kits

Patency of lumen: the diameter of stainless steel stylet to pass through the needle, mm	0.42	0.29	0.23	0.19
Nominal length of the needle tube, mm	25,4 (+1,5 / -2,5)	25,4 (+1,5 / -2,5)	12,7 (+ 1 / - 2)	12,7 (+ 1 / - 2)
Angle of the combined secondary slice	110 ± 10 °			
Nominal length of secondary bevel, mm	1,89-2,12	1.40-1.58	0.95-1,08	0,77-0,99
Point length, mm	4,20-4,70	3.10-3.50	2.10-2.40	1.70-2.20
Force to test bond between hub and needle tube, N	44	34	22	22
Stiffness test, maximum deflection, mm	No more than 0.45	Not more than 0.43	No more than 0.51	No more than 0.34
Cleanliness of outer surface of needle tube	When inspected by normal or corrected-to-normal vision without magnification under an illuminance of 300 lx to 700 lx, the surface of the hypodermic needle tube shall appear free from particles and extraneous matter.			
Cleanliness of hub socket	When examined under 2,5× magnification, the hub socket (fluid path surface) shall appear free from particles and extraneous matter			
Limits for acidity or alkalinity	When determined with a laboratory pH meter and using a general purpose electrode, the pH value of a prepared extract shall be within one unit of pH of that of the control fluid.			

Limits for extractable metals	Needle extract shall, when corrected for the metals content of the control fluid, contain not greater than a combined total of 5 mg/l of lead, tin, zinc and iron. The cadmium content of the extract shall, when corrected for the cadmium content of the control fluid, be lower than 0,1 mg/l
Test on resistance to breakage	It should withstand 20 cycles bends at an angle of 25 ° in opposite directions (An angle of $(25 \pm 1)^\circ$ for regular-walled tubing, $(20 \pm 1)^\circ$ for thin-walled tubing, or $(15 \pm 1)^\circ$ for extra-thin-walled and ultra-thin walled tubing.)
The alignment of the connection, mm	Axis of the tube and the hub must be the same. Maximum allowable deviation from the axis of the hub is $\pm 3^\circ$
The roughness of the outer surface of the tube	Ra less than 0.20 μm 
Needle sharpening surface roughness	Ra less than 0.63 μm
Needle point	When examined under 2,5× magnification the needle point shall appear sharp and free from feather edges, burrs and hooks.
Conical fitting type of needle hub	6 % Luer conical fitting
Lubricant	The lubricant shall not be visible, under normal or corrected vision, as droplets of fluid on the outer or inner surfaces of the needle tube. The quantity of lubricant used should not exceed 0,25 mg/cm ² of the lubricated surface area of the needle tube.
Balance of ethylene oxide	Ethylene oxide (EO) and ethylene chlorohydrin (ECH) <25 p. / million. Ethylene glycol <250 p. / million

* - the color of the connector corresponds to the manufacturer's specifications.

5. Packaging and labeling

5.1. Packaging of medical device

5.1.1. Packaging design

DS MED Intradermal cannula in the set with needles for puncture, in kits is sealed in individual sterile packaging (blister). Material and package design enable visual determination of the color of the needle head, not adversely affecting the contents.


Design and packaging material provides:

- a) preserving sterility during storage in the dry, clean, well-ventilated area;
- b) minimal risk of contamination of the contents on opening the package;
- c) protection of transportation and storage of the contents under normal conditions;
- d) inability to reseal the package needle if its integrity has been violated, opening the package will be obvious.

1st Packaging : Blister Packaging

2nd Packaging : Inner Box

3rd Packaging : Out Box

DS MED Intradermal cannula in the set with needles for puncture, in kits are multipacked in an inbox in an amount of 24 pieces together with instruction for use (1 pc.). Freight Packaging – cartonbox (outbox). Number of inboxes in the cartonbox - 50 pieces.

5.1.2. Packaging material

Packaging materials	Material	Manufacturer
Individual / unit packaging (blister)		
Medical Grade Film K-CPP90	Polypropylene 65%, Polyethylene 35%	Company Name: Kyung Il Industrial Co. Address : #79-1, Ssangsong- Ri, Mado-Mteon, Hwasung- Gun, Kyungki-Do, Korea Tel: +82 031-356-0992
Medical Grade Paper K-CP66	Virgin papers, Cellulose derirerties, Agheous resin polymer	Company Name: Kyung Il Industrial Co.

		Address : #79-1, Ssangsong-Ri, Mado-Miteon, Hwasung-Gun, Kyungki-Do, Korea Tel: +82 031-356-0992
User packaging (inbox)		
Brand paperboard Pottery white board (HI-Q AB Plus)	Paper, pulp, limestone, water	Company Name: Hanlimjigi Address : 348, Cheongnyong-ro, Seosu-myeon, Gunsan-si, Jeollabuk-do, Korea TEL: +82 063-451-7034
Transport wrapping (cartonbox)		
Mark cardboard K2	Paper, cellulose, water	Company Name: Samyang Packing Company Address : 43, Gwangworan-gil, Okgu-eup, Gunsan-si, Jeollabuk-do, Korea TEL: +82 063-464-5522

5.1.3. Features of packaging

Individual / unit packaging (blister)	
The width of the adhesive joint, mm	3-7mm
The force required for opening bondline, H	2,5-5N
Size mm $\pm$ 5mm	140x30mm
Weight, g	Not more than 1
User packaging (inbox)	
Size mm (LxWxH)	195 $\pm$ 5 x 150 $\pm$ 5 x 55 $\pm$ 1
Weight, g	Not more than 150
Transport wrapping (cartonbox)	
Size (LxWxH) mm	760 $\pm$ 5 x 400 $\pm$ 5 x 300 $\pm$ 5

5.2. Labelling of medical device


5.2.1. Individual / unit packaging (blister)

- The name and logo of the device
- Name of the manufacturer and its address
- Country of Origin

- Model
- The size of the device
- The inscription "peel open"
- symbols containing information important for the proper use of medical device.

Deciphering the symbols used in the unit packaging labeling:

	date of manufacture
	manufacturer
	Logo of distributor
	Trademark 
	batch code
	The prohibition on re-use
	Do not re-sterilize
	Sterilization by ethylene oxide
	Do not use if packaging is damaged

	Avoid exposure to sunlight
	Refer to the instructions for use
	Protect from moisture
	use up to
	Temperature Range
	Caution
FTMN	Starting symbols of model letter code


5.2.2. User packaging (inbox)

The user packaging contains the following information:

- the name and logo of the device
- device model,
- device size
- needle gauge and length (as a sticker)
- quantity in package
- indication of sterility and method of sterilization;
- indication of a single application;
- country of origin

- manufacturer and its address;
- authorized representative in Russia and its address
- Information about the conditions of use, storage, disposal and transport of medical device
- date of expiry
- information sticker in Russian
- symbols containing information important for the proper use of medical devices.

Deciphering the symbols used in the user packaging labeling:

	batch code
	Trademark
	Logo of distributor 
	date of manufacture
	manufacturer
	use up to
	Refer to the instructions for use
	The prohibition on re-use

	Do not re-sterilize
STERILE EO	Sterilization by ethylene oxide
	Do not use if packaging is damaged
	Avoid exposure to sunlight
	Protect from moisture
	Temperature Range
	Caution
	Quantity in package



5.2.3. Transport wrapping (cartonbox)


The transport wrapping contains the following information:


- the name and logo of the device
- device size
- quantity in package
- indication of sterility and method of sterilization;
- indication of a single application;

- country of origin
- manufacturer and its address;
- Information about the conditions of use, storage, disposal and transport of medical products
- date of expiry
- number and date of registration certificate
- symbols containing information important for the proper use of medical devices.

Deciphering the symbols used in the transport wrapping labeling:

	manufacturer
	Trademark
	Logo of distributor
	Refer to the instructions for use 
	The prohibition on re-use
	Do not re-sterilize
	Sterilization by ethylene oxide
	Do not use if packaging is damaged

	Avoid exposure to sunlight
	Protect from moisture
	Temperature Range
	Caution
	handle with care
	side to face upward
	protect from sunlight

6. Aftermarket Consumer Information

6.1. Shelf life

Shelf life when stored - 5 years from date of manufacture, provided that all the storage requirements specified by the manufacturer.

6.2. Terms of transportation

Temperature: $+5^{\circ}\text{C} - +25^{\circ}\text{C}$

Relative Humidity: $70 \pm 2\%$;

Atmospheric pressure: $755 \pm 20 \text{ mm Hg}$

6.3. Storage conditions


Store in original container, in closed, dry, cool place at temperature of: $+5^{\circ}\text{C} - +25^{\circ}\text{C}$; relative humidity of: $70 \pm 2\%$ and atmospheric pressure of: $755 \pm 20 \text{ mm Hg}$. Avoid exposure to direct sunlight, high temperatures and high humidity. Packages must be protected from damage.

6.4. Maintenance and repairs

DS MED Intradermal cannula in the set with needles for puncture, in kits is a disposable medical device, in case of damage it can not be repaired.

6.5. Procedure and conditions for disposal

The class of hazardous medical waste in the disposal medical device DSMED intradermal cannulas in the set with needles for puncture, in kits can be attributed to the epidemiologically hazardous waste (infectious and potentially infectious waste contaminated with blood and / or body fluids).

Disposal should be in accordance with the laws and local regulations, as well as the protocols of the medical institution.

6.6. Environmental requirements

Medical device during use, transportation and storage, has no negative effects on humans and the environment.

6.7. Warranty


Warranty valid all the declared shelf life products and extends to articles subject to the conditions of storage, transportation and use.

The company FEEL TECH CO., LTD., Korea, guarantees that the product has been produced and packed with due diligence required in the EU Directive 93/42 / EEC for medical devices. The company FEEL TECH CO, LTD, Korea has the CE mark,

indicating that the company's quality management system complies with the requirements of Annex II (excluding section 4) of Directive 93/42 / EEC, ISO 13485: 2003.

Products have been individually tested and checked carefully. The company FEEL TECH CO., LTD, Korea will replace any product in the event of a defect, admitted during the production and packaging (evidence can be provided), and reports of such defects in the company FEEL TECH CO., LTD, Korea or its authorized representative or distributors. This warranty is in lieu of any other warranty, oral or written, express or implied, and excludes it.

The user is responsible for determining the suitability of the product for use on patients. The company FEEL TECH CO., LTD, Korea and its authorized representative or distributors are not liable for any loss or damage of any nature arising directly or indirectly from the use of the product.

If this product was defective within the warranty period, your exclusive options and sole responsibility of the company FEEL TECH CO., LTD, Korea should be the correction or replacement of the product DS MED Intradermal cannula in the set with needles for puncture, in kits. Liability in respect of the product is stopped, if there was any damage caused by transportation, storage or improper handling user, as well as other factors beyond the control of the manufacturer.

7. Authorized representative of the manufacturer in the Russian Federation

Limited Liability Company "DC MEDICAL"

(OOO "DC MEDICAL"); INN / KPP: 9715354188/771501001

127273, Moscow, Signalny passage, Building 18, room. 109, room 21

Phone: +7 (495) 902-77-78.

DS MED Intradermal cannula in the set with needles for puncture, in kits

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NOTARIAL CERTIFICATE

AHN, JIN WOO -----

attorney-in-fact of
FEEL TECH CO., LTD.
PRESIDENT KIM KEUN SIK -----

appeared before me and admitted
said principal`s subscription to
the attached CERTIFICATE. -----

This is hereby attested on this
26th day of Jan. 2022 at this office.

LAW FIRM IDEA & NOTARY OFFICE INC.

Belong to Seoul Central
District Prosecutor's Office
(Cheongjin-dong, Samsong Bldg.) 701, 702,
34, Jong-ro 3-gil, Jongno-gu, Seoul,
Republic of Korea

Signature of the Notary Public
Park Jong Cheil

This office has been authorized by the
Minister of Justice, the Republic of
Korea, to act as Notary Public Since
7, Feb. 2020 Under Law No.243.

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